

---

# **pylablib Documentation**

***Release 1.4.1***

**Alexey Shkarin**

**Oct 06, 2022**



---

## Contents:

---

<b>1</b>	<b>Related projects</b>	<b>3</b>
1.1	Installation . . . . .	3
1.2	Devices overview . . . . .	6
1.3	Data processing . . . . .	85
1.4	Data storage . . . . .	89
1.5	Various utilities . . . . .	93
1.6	Change log . . . . .	96
1.7	pylablib . . . . .	103
<b>2</b>	<b>Indices and tables</b>	<b>877</b>
	<b>Python Module Index</b>	<b>879</b>
	<b>Index</b>	<b>883</b>



PyLabLib aims to provide support for device control and experiment automation. It interfaces with lots of different *devices*, including several different *camera interfaces*, *translational stages*, *oscilloscopes*, *AWGs*, *sensors*, and more. The interface is implemented in a natural way through Python objects, and is easy to understand. For example, here is a complete script which steps *Thorlabs KDC101* stage by 10000 steps ten times, and each time grabs a frame with *Andor iXon camera*:

```
from pylablib.devices import Thorlabs, Andor # import the device libraries
import numpy as np # import numpy for saving

# connect to the devices
with Thorlabs.KinesisMotor("27000000") as stage, Andor.AndorSDK2Camera() as cam:
    # change some camera parameters
    cam.set_exposure(50E-3)
    cam.set_roi(0, 128, 0, 128, hbin=2, vbin=2)
    # start the stepping loop
    images = []
    for _ in range(10):
        stage.move_by(10000) # initiate a move
        stage.wait_move() # wait until it's done
        img = cam.snap() # grab a single frame
        images.append(img)

np.array(images).astype("<u2").tofile("frames.bin") # save frames as raw binary
```

The list of the devices is constantly expanding.

Additional utilities are added to simplify data acquisition, storage, and processing:

- Simplified data processing utilities: convenient *fitting*, *filtering*, *feature detection*, *FFT* (mostly wrappers around NumPy and SciPy).
- Universal multi-level *dictionaries* which are convenient for *storing* heterogeneous data and settings in human-readable format.
- Assorted functions for dealing with *file system* (creating, moving and removing folders, zipping/unzipping, path normalization), *network* (simplified interface for client and server sockets), *strings* (conversion of various Python objects to and from string), and more.
- Tools for GUI generation and advanced multi-threading built on top of Qt5 (*still in development stage: the documentation is incomplete, and the interfaces can change in later versions*)

The library only works on Python 3, and has been most extensively tested on Windows 10 with 64-bit Python. Linux is, in principle, supported, but devices which require manufacturer-provided DLLs (mostly cameras) might, potentially, have problems.

---

**Note:** This is documentation for the newer **1.x** version of the library. The older **0.x** documentation can be found at <https://pylablib-v0.readthedocs.io/en/latest/>.

---



[Pylablib cam-control](#) - software for universal camera control and camera data acquisition.

## 1.1 Installation

### 1.1.1 Standard install

You can install the library from PyPi:

```
pip install pylablib
```

If you already have it installed, you can upgrade it to get the newest version:

```
pip install -U pylablib
```

This will install the full set of dependencies: basic dependencies and computing packages (`numpy`, `scipy`, `pandas`, `numba`, `rpyc`), basic device communication packages (`pyft232`, `pyvisa`, `pyserial`, `pyusb`), and PyQt5-based GUI (`pyqt5` and `pyqtgraph`). You can also install additional device library dependencies (`nidaqmx` and `websocket-client`) using the extra requirements feature of pip:

```
pip install -U pylablib[devio-full]
```

### 1.1.2 Minimal install

In case you do not want some of these packages installed, or they are unavailable on your platform, you can install a lightweight version of pylablib called `pylablib-lightweight`. It contains exactly the same code, but has only the most basic dependencies (`numpy`, `scipy`, and `pandas`):

```
pip install -U pylablib-lightweight
```

With this, the basic functionality (such as data processing or file IO) will work, but more advanced features such as device communication and GUI, will require additional packages. In most cases, the raised errors will notify which packages are missing. These can be installed either manually, or using the extra requirements:

- [extra] extra packages used in some situations: numba (speeds up some data processing) and rpyc (communication between different PCs)
- [devio] basic devio packages: pyft232, pyvisa, pyserial, and pyusb
- [devio-extra] additional devio packages: nidaqmx and websocket-client
- [gui-pyqt5] **PyQt5**-based GUI: pyqt5 and pyqtgraph. Should not be used together with [gui-pyside2]
- [gui-pyside2] **PySide2**-based GUI: pyside2 and pyqtgraph. Should not be used together with [gui-pyqt5]

The options can be combined. For example,

```
pip install pylablib-lightweight[extra,devio,gui-pyside2]
```

installs the dependencies as the usual pylablib distribution, but with PySide2 Qt5 backend instead of PyQt5.

### 1.1.3 Anaconda install

The package is also available on Anaconda via conda-forge channel. To install it, run

```
conda install -c conda-forge pylablib
```

in the Anaconda prompt.

The Anaconda version of pylablib comes with all the standard dependencies except for pyft232, nidaqmx and websocket-client, which are not available on conda-forge channel. This means, that *Thorlabs APT/Kinesis*, *NI DAQs*, and some functionality of *M2 Solstis laser* are not accessible. To use those, it is recommended to either install those packages explicitly via pip (keep in mind that it can break Anaconda environment), or use a standalone Python distribution.

### 1.1.4 Usage

To access to the most common features simply import the library:

```
import pylablib as pll
# Create a parameter dictionary (e.g., for some processing script)
parameters = pll.Dictionary({"par/x":1, "par/y":2, "par/z":[3,4,5], "out":"result"})
pll.save_dict(parameters, "parameters.dat") # save parameters to a text file
```

More advanced features (e.g., *device communication*) should be imported directly:

```
from pylablib.devices import Andor # import Andor devices module
cam = Andor.AndorSDK2Camera() # connect to Andor SDK2 camera (e.g., iXon)
cam.set_exposure(0.1) # set exposure to 100ms
frame = cam.snap() # grab a single frame
cam.close() # close the connection
```



### 1.1.5 Dependencies and requirements

The basic package dependencies are `NumPy` for basic computations and overall array interface, `SciPy` for advanced computations (interpolation, optimization, special functions), and `pandas` for heterogeneous tables (`DataFrame`). In addition, it is recommended to have `Numba` package to speed up some computations. Finally, if you use options for remote computing and communication between different PCs, you need to install `RPyC`. Note that when installed directly from pip, `numpy` comes with the OpenBLAS version of the linear algebra library; if other version (e.g., Intel MKL) is preferred, it is a good idea to have `numpy` already installed before installing `pylablib`.

The main device communication packages are `PyVISA` and `pySerial`, which cover the majority of devices. Several devices (e.g., *Thorlabs Kinesis* and *Attocube ANC 350*) require additional communication packages: `pyft232` and `PyUSB`. Finally, some particular devices completely or partially rely on specific packages: `NI-DAQmx` for *NIDAQ* and `websocket-client` for additional *M2 Solstis* functionality.

Finally, GUI and advanced multi-threading relies on Qt5, which has two possible options. The first (default) option is `PyQt5` with `sip` for some memory management functionality. Note that while newer `PyQt5` versions  $\geq 5.11$  already come with `PyQt5-sip`, older versions require a separate `sip` installation. Hence, if you use an older `PyQt5` version, you need to install `sip` separately. The second possible Qt5 option is `PySide2` with `shiboken2`. Both `PyQt5` and `PySide2` should work equally well, and the choice mostly depends on what is already installed, because having both `PyQt5` and `PySide2` might lead to conflicts. Finally, plotting relies on `pyqtgraph`, which, starting with version 0.11m is compatible with both `PySide2` and `PyQt5`.

The package has been tested with Python 3.6 through 3.9, and is incompatible with Python 2. The last version officially supporting Python 2.7 is 0.4.0. Furthermore, testing has been mostly performed on 64-bit Python. This is the recommended option, as 32-bit version limitations (most notably, limited amount of accessible RAM) mean that it should only be used when absolutely necessary, e.g., when some required packages or libraries are only available in 32-bit version.

### 1.1.6 Installing from GitHub

The most recent and extensive, but less tested and documented, version of this library is available on GitHub at <https://github.com/AlexShkarin/pyLabLib/>. There are several versions of installing it:

- Install using pip using GitHub as a library source:

```
pip install -U git+https://github.com/AlexShkarin/pyLabLib.git
```

- Download it as a zip-file and unpack it into any appropriate place (can be folder of the project you're working on, Python `site-packages` folder, or any folder added to `PATH` or `PYTHONPATH` variable).

To download the code of a specific version, you can choose it in the dropdown *Branch* menu under *Tags* tab. This is the same code as available on PyPi.

Keep in mind that, unlike the first method, the required packages will not be automatically installed, so this has to be done manually:

```
pip install numpy scipy pandas numba rpyc
pip install pyft232 pyvisa pyserial pyusb nidaqmx websocket-client
pip install pyqt5 pyqtgraph
```

- Clone the repository to your computer In order to easily get updates in order to easily get updates. For that, you need to install Git (<https://git-scm.com/>), and use the following commands in the command line (in the folder where you want to store the library):

```
git clone https://github.com/AlexShkarin/pyLabLib
cd ./pyLabLib
```

Whenever you want to update to the most recent version, simply type

```
git pull
```

in the library folder. Keep in mind that any changes that you make to the library code might conflict with the new version that you pull from GitHub, so you should not modify anything in this folder if possible.

### 1.1.7 Support and feedback

If you have any issues, suggestions, or feedback, you can either raise an issue on GitHub at <https://github.com/AlexShkarin/pyLabLib/issues>, or send an e-mail to [pylablib@gmail.com](mailto:pylablib@gmail.com).

## 1.2 Devices overview

Basic concepts are described at the *general device communication page*.

Currently supported devices:

- *Cameras*
  - *Andor SDK2* and *Andor SDK3*: variety of Andor (currently part of Oxford Instruments) cameras. Tested with Andor iXon, Luca, Zyla, Neo and Marana.
  - *Allied Vision Bonito* cameras: CameraLink-interfaced cameras. Tested with Bonito CL-400B/C and NI IMAQ frame grabber.
  - *Basler*: Basler pylon-compatible cameras. Tested with an emulated Basler camera.
  - *BitFlow*: BitFlow Axion family frame grabbers. Tested with BitFlow Axion 1xB frame grabber together with PhotonFocus MV-D1024E camera.
  - *DCAM*: Hamamatsu cameras. Tested with Hamamatsu Orca Flash 4.0 and ImagEM.
  - *NI IMAQ*: National Instruments frame grabbers. Tested with NI PCI-1430 and PCI-1433 frame grabbers together with PhotonFocus MV-D1024E camera.
  - *NI IMAQdx*: National Instruments universal camera interface. Tested with Ethernet-connected PhotonFocus HD1-D1312 camera.
  - *Photon Focus*: Photon Focus pfcam interface. Tested with PhotonFocus MV-D1024E camera connected through either NI frame grabbers (PCI-1430 and PCI-1433) or Silicon Software frame grabbers (microEnable IV AD4-CL).
  - *PCO SC2*: PCO cameras. Tested with pco.edge cameras with CLHS and regular CameraLink interfaces.
  - *Picam*: Princeton Instruments cameras. Tested with a PIXIS 400 camera.
  - *PVCAM*: Photometrics cameras. Tested with a Prime 95B camera.
  - *Silicon Software*: Silicon Software frame grabbers. Tested with microEnable IV AD4-CL frame grabbers together with PhotonFocus MV-D1024E camera.
  - *Thorlabs Scientific Cameras*: Thorlabs sCMOS cameras. Tested with Thorlabs Kiralux camera.
  - *Uc480/uEye*: multiple cameras, including simple Thorlabs and IDS cameras. Tested with IDS SC2592R12M and Thorlabs DCC1545M.
- *Stages*
  - *Attocube ANC300* and *Attocube ANC350*: most common Attocube positioner controllers. Tested with Ethernet and USB connection for ANC300, and USB connection for ANC350.

- *Thorlabs APT/Kinesis*: basic Thorlabs motorized stages and optomechanics devices. Tested with KDC101, K10CR1, and BSC201 motor controllers, KIM101 piezo motor controller, as well as MFF101 and FW102 (described at a [different page](#))
- *Thorlabs Elliptec*: resonant piezoelectric Thorlabs stages. Tested with ELL18 and ELL14 rotational mounts.
- *Newport Picomotor*: precision piezo-actuated screws based on slip-stick principle. Tested with Newport 8742 Picomotor driver using Ethernet or USB connection.
- *Arcus Performax*: fairly common single- and multi-axis motor controllers sold under different brands: Arcus, Nippon Pulse America, or Newmark Systems. Tested with PMX-4EX device with USB connection.
- *Trinamic*: universal motor controllers and drivers. Tested with a single-axis TMC2130 controller with USB connection.
- *SmarAct*: high-performance piezo sliders. So far only simple open-loop *SCU controllers* are supported. Tested with a standard HCU controller unit.
- *Physik Instrumente*: piezo controllers. So far only PI E-516 is supported and tested.
- *Basic sensors*
  - *HighFinesse*: laser wavelength meters. Tested with WS6 and WS7 USB-controlled devices.
  - *Ophir*: optical power and energy meters. Tested with Ophir Vega.
  - *Lakeshore*: temperature sensors. Tested with Lakeshore 218.
  - *Cryocon*: temperature sensors. Tested with CryoCon 14C.
  - *Cryomagnetics*: liquid nitrogen or helium level sensor. Tested with LM-500 and LM-510 sensors.
  - *Pfeiffer*: pressure gauges. Tested with TPG261 and DPG202 controllers.
  - *Leybold*: pressure gauges. Tested with ITR90 gauge.
  - *Kurt J. Lesker*: pressure gauges. Tested with KJL300 gauge.
  - *Thorlabs quadrature detector controller*. Tested with TPA101.
- Lasers
  - *Basic lasers*
    - \* *Lighthouse Photonics SproutG*
    - \* *Laser Quantum Finesse*
  - *M2 Solstis laser*
  - *Toptica iBeam Smart laser*
  - *Sirah Matisse laser*
- *Tektronix oscilloscopes*. Tested with TDS2002B, TDS2004B, and DPO2004B.
- *NI DAQs*. Tested with NI USB-6008, NI USB-6343, and NI PCIe-6323.
- *Generic AWGs*. Tested with Agilent 33500 and 33220A, Rigol DG1022, Tektronix AFG1022, GW Instek AFG2225 and AFG2115, and RS Comp AFG21005.
- *Miscellaneous Thorlabs devices*: *MFF101/102* motorized flip mirror mount, *FW102/212* motorized filter wheel, and *MDT693/694* high-voltage source.
- *Miscellaneous OZOptics devices*: *EPC04* fiber polarization controller, *DD100* motorized fiber attenuator, and *TF100* motorized fiber filter.

- *Miscellaneous devices*
  - *Conrad relay board*
  - *Basic Arduino communication*

## 1.2.1 Basics of device communication

The devices are represented as Python objects. In most cases, one object controls one device, although sometimes one object can be responsible for multiple interconnected devices (e.g., when daisy-chaining of several devices is used, as in *Picomotor stage*). All the device control functions are contained within the class. Occasionally, there are auxiliary function present for listing available devices, dealing with data generated by the device, or adjusting global parameters.

---

**Note:** Some specific devices functionality might not be completely covered in the current release. If this is the case for your device, you can let the developers know by raising an [issue on GitHub](#), or sending an e-mail to [pylablib@gmail.com](mailto:pylablib@gmail.com).

---

### Connection

The device identifier or address needs to be provided upon the device object creation, after which it is automatically connected. Getting the address usually depends on the kind of device:

- Simple message-style devices, such as AWG, oscilloscopes, sensors and gauges, require an address which depends on the exact connection protocol. For example, serial devices addresses look like "COM1", Visa addresses as "USB0::0x1313::0x8070::000000::INSTR", and network addresses take IP and, possibly, port "192.168.1.3:7230". To get the list of all connected devices, you can run `comm_backend.list_backend_resources()`:

```
>> import pylablib as pll
>> pll.list_backend_resources("serial") # list serial port resources
['COM38', 'COM1', 'COM36', 'COM3']
>> pll.list_backend_resources("visa") # note that, by default, visa also
↳ includes all the COM ports
('USB0::0x1313::0x8070::000000::INSTR',
'ASRL1::INSTR',
'ASRL3::INSTR',
'ASRL10::INSTR',
'ASRL36::INSTR',
'ASRL38::INSTR')
```

Network devices do not easily provide such functionality (and there are, in principle, many unrelated devices connected to the network), so you might need to learn the device IP elsewhere. Usually, it is set on the device front panel or using some kind of configuration tool and a different connection, such as serial or USB.

In most cases, the connection address is all you need. However, sometimes the connection might require some additional information. The most common situations are ports for the network connection and baud rates for the serial connections. Ports can be supplied either as a part of the string "192.168.1.3:7230", or as a tuple ("192.168.1.3", 7230). The baud rates are, similarly, provided as a tuple: ("COM1", 19200). By default, the devices would use the baud rate which is most common for them, but in some cases (e.g., if the device baud rate can be changed), you might need to provide it explicitly. If it is provided incorrectly, then no communication can be done, and requests will typically return a timeout error:

```
>> from pylablib.devices import Ophir
>> meter = Ophir.VegaPowerMeter("COM3") # for this power meter 9600 baud are
↳used by default
>> meter.get_power() # let us assume that the devices is currently set up with
↳38400 baud
...
OphirBackendError: backend exception: 'timeout during read'
>> meter.close() # need to close the connection before reopening
>> meter = Ophir.VegaPowerMeter(("COM3",38400)) # explicitly specifying the
↳correct baud rate
>> meter.get_power()
1E-6
```

- More complicated devices using custom DLLs (usually cameras or some translation stages) will have more unique methods of addressing individual devices: serial number, device index, device ID, etc. In most cases such devices come with `list_devices` or `get_devices_number` functions, which give the necessary information.

After communication is done, the connection needs to be closed, since in most cases it can only be opened in one program or part of the script at a time. It also implies that usually it's impossible to connect to the device while its manufacturer software is still running.

The devices have `open` and `close` methods, but they can also work in together with Python with statements:

```
# import Thorlabs device classes
from pylablib.devices import Thorlabs

# connect to FW102 motorized filter wheel
wheel = Thorlabs.FW("COM1")
# set the position
wheel.set_position(1)
# close the connection (until that it's impossible to establish a different
↳connection to this device)
wheel.close()

# a better approach
with Thorlabs.FW("COM1") as wheel: # connection is closed automatically when leaving
↳the with-block
    wheel.set_position(1)
```

Because the devices are automatically connected on creation, `open` method is almost never called explicitly. It is generally only used to reconnect to the device after the connection has been previously closed, although in this case creating a new device object would work just as well.

## Operation

The devices are controlled by calling their methods; attributes and properties are very rarely used. Effort is made to maintain consistent naming conventions, e.g., most getter-methods will start with `get_` and setter methods with `set_` or `setup_` (depending on the complexity of the method). It is also common for setter methods to return the new value as a result, which is useful in CLI operation and debugging. Devices of the same kind have the same names for similar or identical functions: most stages have `move_by`, `jog` and `stop` methods, and cameras have `wait_for_frame` and `read_multiple_images` methods. Whenever it makes sense, these methods will also have the same signatures.

## Asynchronous operation and multi-threading

For simplicity of usage and construction, devices interfaces are designed to be synchronous and single-threaded. Asynchronous operation can be achieved by explicit usage of Python multi-threading. Furthermore, the device classes are not designed to be thread safe, i.e., it is not recommended to use the same device simultaneously from two separate threads. However, non-simultaneous calling of device methods from different threads (synchronized, e.g., using locks) or simultaneous usage of several separate devices of the same class is supported.

## Error handling

Errors raised by the devices are usually specific to the device and manufacturer, e.g., `AttocubeError` or `TrinamicError`. These can be obtained from the module containing the device class, or from the class itself as `Error` attribute:

```
>> from pylablib.devices import Attocube
>> atc = Attocube.ANC300("192.168.1.1")
>> atc.disable_axis(1)
>> atc.move_by(1,10) # move on a disabled axis raises an error for ANC300
...
AttocubeError: Axis in wrong mode
>> try:
..     atc.move_by(1,10)
.. except atc.Error: # could also write "except Attocube.AttocubeError"
..     print("Can not move")
Can not move
```

All of the device errors inherit from `DeviceError`, which in turn is a subclass of `RuntimeError`. Therefore, one can also use those exception classes instead:

```
>> import pylablib as pll
>> try:
..     atc.move_by(1,10)
.. except pll.DeviceError:
..     print("Can not move")
Can not move
```

## Getting more information

A lot of information about the devices can be gained just from their method names and descriptions (docstrings). There are several ways of getting these:

- In many cases your IDE (PyCharm, Spyder, VS Code with installed Python extension) supports code inspection. In this case, the list of methods will usually pop up after you type the device object name and a dot (such as `cam.`), and the method docstring will show up after you type the method name and parenthesis (such as `cam.get_roi()`). However, sometimes it might take a while for these pop-ups to show up.
- You can use console, such as Jupyter QtConsole, Jupyter Notebook, or a similar console built into the IDE. Here the list of methods can be obtained using the autocomplete feature: type name of the class or object with a dot (such as `cam.`) and then press `Tab`. The list of all methods should appear. To get the description of a particular class or method, type it with a question mark (such as `cam?` or `cam.get_roi?`) and execute the result (`Enter` or `Shift-Enter`, depending on the console). A description should appear with the argument names and the description.

- You can also use the auto-generated documentation within this manual through the search bar: simply type the name of the class or the method (such as `AndorSDK3Camera` or `AndorSDK3Camera.get_roi`) and look through the results. However, the formatting of the auto-generated documentation might be a bit overwhelming.

## Universal settings access

All devices have `get_settings` and `apply_settings` methods which, correspondingly, return Python dictionaries with the most common settings or take these dictionaries and apply the contained settings. These can be used to easily store and re-apply device configuration within a script.

Additionally, there is `get_full_info` method, which returns as complete information as possible. It is particularly useful to check the device status and see if it is connected and working properly, and to save the devices configuration when acquiring the data. Finally, the settings can also be accessed through `.dv` attribute, which provides dictionary-like interface:

```
>>> wheel = Thorlabs.FW("COM1") # connect to FW102 motorized filter wheel
>>> wheel.get_position()
1
>>> wheel.get_settings()
{'pcount': 6,
 'pos': 1,
 'sensors_mode': 'off',
 'speed_mode': 'high',
 'trigger_mode': 'in'}
>>> wheel.dv["pos"]
1
>>> wheel.apply_settings({"pos":2})
>>> wheel.get_position()
2
>>> wheel.dv["pos"] = 3
>>> wheel.get_position()
3
>>> wheel.close()
```

By default not all information is shown, as it can take long time (up to several seconds) to obtain it, and it takes a lot of space on the screen. To get a full set of parameters, you can call `get_full_info("all")`:

```
>> cam = IMAQdx.IMAQdxCamera()
>> cam.get_full_info()
{ 'roi': (0, 1312, 0, 1082),
  'acquisition_in_progress': False,
  'frames_status': TFramesStatus(acquired=0, unread=0, skipped=0, buffer_size=0),
  'cls': 'IMAQdxCamera',
  'conn': 'cam0',
  'detector_size': (1312, 1082),
  'device_info': TDeviceInfo(vendor='Photonfocus AG', model='HD1-D1312-80-G2-12',
↳serial_number='0000000000000000', bus_type='Ethernet') }
>>
>> cam.get_full_info("all")
{ 'roi': (0, 1312, 0, 1082),
  'acquisition_in_progress': False,
  'frames_status': TFramesStatus(acquired=0, unread=0, skipped=0, buffer_size=0),
  'camera_attributes': Dictionary('AcquisitionAttributes/AdvancedEthernet/
↳BandwidthControl/ActualPeakBandwidth': 1000.0
  ... lots and lots of attributes
  'OffsetX': 0
```

(continues on next page)



(continued from previous page)

```

'OffsetY': 0
'PayloadSize': 1419584
'PixelFormat': Mono8
'Width': 1312),
'cls': 'IMAQdxCamera',
'conn': 'cam0',
'detector_size': (1312, 1082),
'device_info': TDeviceInfo(vendor='Photonfocus AG', model='HD1-D1312-80-G2-12',
↪serial_number='0000000000000000', bus_type='Ethernet') }
```

## Dependencies and external software

Many devices require external software not provided with this package.

The simpler devices using serial connection (either with an external USB-to-Serial adapter, or with a similar built-in chip) only need the corresponding drivers: either standard adapter drivers or the ones supplied by the manufacturer, e.g., via Thorlabs APT software. If the device already shows up as a serial communication port in the OS, no additional software is normally needed. Similarly, devices using Ethernet connection do not need any external software, as long as they are properly connected to the network. Finally, devices using Visa connection require NI VISA Runtime, which is freely available from the [National Instruments website](#). See also [PyVISA documentation](#) for details.

Devices which require manufacturer DLLs are harder to set up. For most of them, at the very least, you need to install the manufacturer-provided software for communication. Frequently it already includes the necessary libraries, which means that nothing else is required. However, sometimes you would need to download either an additional SDK package, or DLLs directly from the website. Since these libraries take a lot of space and are often proprietary, they are not distributed with the pylablib.

Note that DLLs can have 32-bit and 64-bit version, and this version should agree with the Python version that you use. Unless you have a really good reason to do otherwise, it is strongly recommended to use 64-bit Python, which means that you would need 64-bit DLLs, which is the standard in most cases these days. To check your Python bitness, you can read the prompt when running the Python console, or run `python -c "import platform; print(platform.architecture()[0])"` in the command line.

In addition, you need to provide pylablib with the path to the DLLs. In many cases it checks the standard locations such as the default `System32` folder (used, e.g., in DCAM or IMAQ cameras), paths contained on the `PATH` environment variable, or defaults paths for manufacturer software (such as `C:/Program Files/Andor SOLIS` for Andor cameras). If the software path is different, or if you choose to obtain DLLs elsewhere, you can also explicitly provide path by setting the `library` parameter:

```

import pylablib as pll
pll.par["devices/dlls/andor_sdk3"] = "D:/Program Files/Andor SOLIS"
from pylablib.devices import Andor
cam = Andor.AndorSDK3Camera()
```

All of these requirements are described in detail for the specific devices.

Starting from Python 3.8 the DLL search path is changed to not include the files contained in `PATH` environment variable and in the script folder. By default, this behavior is still emulated when pylablib searches for the DLLs, since it is required in some cases (e.g., Photon Focus pfcam interface). If needed, it can be turned off (i.e., switched to the new default behavior of Python 3.8+) by setting `pll.par["devices/dlls/add_envIRON_paths"]=False`.

## Advanced examples

Connecting to a Cryomagnetics LM500 level meter and reading out the level at the first channel:



```
from pylablib.devices import Cryomagnetics # import the device library
with Cryomagnetics.LM500("COM1") as lm:
    level = lm.get_level(1) # read the level
```

Stepping the M Squared laser wavelength and recording an image from the Andor iXon camera at each step:

```
with M2.Solstis("192.168.1.2", 34567) as laser, Andor.AndorSDK2Camera() as cam: #_
    ↪connect to the devices
    # change some camera parameters
    cam.set_exposure(50E-3)
    cam.set_roi(0, 128, 0, 128, hbin=2, vbin=2)
    cam.setup_shutter("open")
    # start camera acquisition
    wavelength = 770E-9 # initial wavelength (in meters)
    images = []
    cam.start_acquisition()
    while wavelength < 780E-9:
        laser.coarse_tune_wavelength(wavelength) # tune the laser frequency (using_
    ↪coarse tuning)
        time.sleep(0.5) # wait until the laser stabilizes
        cam.wait_for_frame() # ensure that there's a frame in the camera queue
        img = cam.read_newest_image()
        images.append(img)
        wavelength += 0.5E-9
```

## Available devices

- *Cameras*

- *Andor SDK2* and *Andor SDK3*: variety of Andor (currently part of Oxford Instruments) cameras. Tested with Andor iXon, Luca, Zyla, Neo and Marana.
- *Allied Vision Bonito* cameras: CameraLink-interfaced cameras. Tested with Bonito CL-400B/C and NI IMAQ frame grabber.
- *Basler*: Basler pylon-compatible cameras. Tested with an emulated Basler camera.
- *BitFlow*: BitFlow Axion family frame grabbers. Tested with BitFlow Axion 1xB frame grabber together with PhotonFocus MV-D1024E camera.
- *DCAM*: Hamamatsu cameras. Tested with Hamamatsu Orca Flash 4.0 and ImagEM.
- *NI IMAQ*: National Instruments frame grabbers. Tested with NI PCI-1430 and PCI-1433 frame grabbers together with PhotonFocus MV-D1024E camera.
- *NI IMAQdx*: National Instruments universal camera interface. Tested with Ethernet-connected PhotonFocus HD1-D1312 camera.
- *Photon Focus*: Photon Focus pfcam interface. Tested with PhotonFocus MV-D1024E camera connected through either NI frame grabbers (PCI-1430 and PCI-1433) or Silicon Software frame grabbers (microEnable IV AD4-CL).
- *PCO SC2*: PCO cameras. Tested with pco.edge cameras with CLHS and regular CameraLink interfaces.
- *Picam*: Princeton Instruments cameras. Tested with a PIXIS 400 camera.
- *PVCAM*: Photometrics cameras. Tested with a Prime 95B camera.
- *Silicon Software*: Silicon Software frame grabbers. Tested with microEnable IV AD4-CL frame grabbers together with PhotonFocus MV-D1024E camera.

- *Thorlabs Scientific Cameras*: Thorlabs sCMOS cameras. Tested with Thorlabs Kiralux camera.
- *Uc480/uEye*: multiple cameras, including simple Thorlabs and IDS cameras. Tested with IDS SC2592R12M and Thorlabs DCC1545M.
- *Stages*
  - *Attocube ANC300* and *Attocube ANC350*: most common Attocube positioner controllers. Tested with Ethernet and USB connection for ANC300, and USB connection for ANC350.
  - *Thorlabs APT/Kinesis*: basic Thorlabs motorized stages and optomechanics devices. Tested with KDC101, K10CR1, and BSC201 motor controllers, KIM101 piezo motor controller, as well as MFF101 and FW102 (described at a [different page](#))
  - *Thorlabs Elliptec*: resonant piezoelectric Thorlabs stages. Tested with ELL18 and ELL14 rotational mounts.
  - *Newport Picomotor*: precision piezo-actuated screws based on slip-stick principle. Tested with Newport 8742 Picomotor driver using Ethernet or USB connection.
  - *Arcus Performax*: fairly common single- and multi-axis motor controllers sold under different brands: Arcus, Nippon Pulse America, or Newmark Systems. Tested with PMX-4EX device with USB connection.
  - *Trinamic*: universal motor controllers and drivers. Tested with a single-axis TMCM-1110 controller with USB connection.
  - *SmarAct*: high-performance piezo sliders. So far only simple open-loop *SCU controllers* are supported. Tested with a standard HCU controller unit.
  - *Physik Instrumente*: piezo controllers. So far only PI E-516 is supported and tested.
- *Basic sensors*
  - *HighFinesse*: laser wavelength meters. Tested with WS6 and WS7 USB-controlled devices.
  - *Ophir*: optical power and energy meters. Tested with Ophir Vega.
  - *Lakeshore*: temperature sensors. Tested with Lakeshore 218.
  - *Cryocon*: temperature sensors. Tested with CryoCon 14C.
  - *Cryomagnetics*: liquid nitrogen or helium level sensor. Tested with LM-500 and LM-510 sensors.
  - *Pfeiffer*: pressure gauges. Tested with TPG261 and DPG202 controllers.
  - *Leybold*: pressure gauges. Tested with ITR90 gauge.
  - *Kurt J. Lesker*: pressure gauges. Tested with KJL300 gauge.
  - *Thorlabs quadrature detector controller*. Tested with TPA101.
- *Lasers*
  - *Basic lasers*
    - \* *Lighthouse Photonics SproutG*
    - \* *Laser Quantum Finesse*
  - *M2 Solstis laser*
  - *Toptica iBeam Smart laser*
  - *Sirah Matisse laser*
- *Tektronix oscilloscopes*. Tested with TDS2002B, TDS2004B, and DPO2004B.
- *NI DAQs*. Tested with NI USB-6008, NI USB-6343, and NI PCIe-6323.

- *Generic AWGs*. Tested with Agilent 33500 and 33220A, Rigol DG1022, Tektronix AFG1022, GW Instek AFG2225 and AFG2115, and RS Comp AFG21005.
- *Miscellaneous Thorlabs devices*: *MFF101/102* motorized flip mirror mount, *FW102/212* motorized filter wheel, and *MDT693/694* high-voltage source.
- *Miscellaneous OZOptics devices*: *EPC04* fiber polarization controller, *DD100* motorized fiber attenuator, and *TF100* motorized fiber filter.
- *Miscellaneous devices*
  - *Conrad relay board*
  - *Basic Arduino communication*

## 1.2.2 Cameras

Basic concepts are described at the [cameras communication page](#).

Currently supported cameras:

- *Andor SDK2* and *Andor SDK3*: variety of Andor (currently part of Oxford Instruments) cameras. Tested with Andor iXon, Luca, Zyla, Neo and Marana.
- *Allied Vision Bonito* cameras: CameraLink-interfaced cameras. Tested with Bonito CL-400B/C and NI IMAQ frame grabber.
- *Basler*: Basler pylon-compatible cameras. Tested with an emulated Basler camera.
- *BitFlow*: BitFlow Axion family frame grabbers. Tested with BitFlow Axion 1xB frame grabber together with PhotonFocus MV-D1024E camera.
- *DCAM*: Hamamatsu cameras. Tested with Hamamatsu Orca Flash 4.0 and ImagEM.
- *NI IMAQ*: National Instruments frame grabbers. Tested with NI PCI-1430 and PCI-1433 frame grabbers together with PhotonFocus MV-D1024E camera.
- *NI IMAQdx*: National Instruments universal camera interface. Tested with Ethernet-connected PhotonFocus HD1-D1312 camera.
- *Photon Focus*: Photon Focus pfcam interface. Tested with PhotonFocus MV-D1024E camera connected through either NI frame grabbers (PCI-1430 and PCI-1433) or Silicon Software frame grabbers (microEnable IV AD4-CL).
- *PCO SC2*: PCO cameras. Tested with pco.edge cameras with CLHS and regular CameraLink interfaces.
- *Picam*: Princeton Instruments cameras. Tested with a PIXIS 400 camera.
- *PVCAM*: Photometrics cameras. Tested with a Prime 95B camera.
- *Silicon Software*: Silicon Software frame grabbers. Tested with microEnable IV AD4-CL frame grabbers together with PhotonFocus MV-D1024E camera.
- *Thorlabs Scientific Cameras*: Thorlabs sCMOS cameras. Tested with Thorlabs Kiralux camera.
- *Uc480/uEye*: multiple cameras, including simple Thorlabs and IDS cameras. Tested with IDS SC2592R12M and Thorlabs DCC1545M.

---

**Note:** General device communication concepts are described on the corresponding [page](#).

---

## Cameras control basics

### Basic examples

Basic camera usage is fairly straightforward:

```
from pylablib import Andor
cam = Andor.AndorSDK3Camera() # connect to the camera
cam.set_exposure(10E-3) # set 10ms exposure
cam.set_roi(0,128,0,128) # set 128x128px ROI in the upper left corner
images = cam.grab(10) # grab 10 frames
cam.close()
```

In case you need to grab and process frames continuously, the example is a bit more complicated:

```
with Andor.AndorSDK2Camera() as cam: # to close the camera automatically
    cam.start_acquisition() # start acquisition (automatically sets it up as well)
    while True: # acquisition loop
        cam.wait_for_frame() # wait for the next available frame
        frame=cam.read_oldest_image() # get the oldest image which hasn't been read
    yet
    # ... process frame ...
```

Some concepts are explained below in more detail.

### Basic concepts

#### Frames buffer

In most cases, the frames acquired by the camera are first temporarily stored in the local camera and / or frame grabber memory, from which they are transferred to the PC RAM by the camera drivers. Afterwards, this memory is made available to all other applications. In principle, it should be enough to store only the most recent frame in RAM, and for the user software to continuously wait for a new frame, immediately read it from RAM and process it. However, such approach is very demanding to the user code: if the new frame is acquired before the previous one is processed or copied, then the RAM data is overwritten, and the old frame is lost. Hence, it is more practical to have a *buffer* of several most recently acquired frames to account for inevitable interruptions in the user wait-read-process loop caused by OS scheduling and by other jobs. In this case, the frames get lost only when the buffer is completely filled, and the oldest frames starts getting overwritten.

When using the camera classes provided by pylablib, you do not need to worry about setting up the buffer yourself, since it is done behind the scene either by the manufacturer's code or by the device class. However, it is important to keep in mind the existence of the buffer when setting up the acquisition, interpreting the buffer and acquired frames status, or identifying the skipped frames.

The size of the buffer can almost always be selected by the user. Typically it is a good idea to have at least 100ms worth of frames there, although, depending on the other jobs performed by the software, it can be larger.

#### Acquisition setup

Setting up an acquisition process might take a lot of time (up to 10s in more extreme cases). This happens mostly because of the buffer allocation and setting up internal API structures; initiating the acquisition process itself is fairly fast. Hence, it is useful to separate setting up / cleaning up and starting / stopping.

The first two procedures correspond to `setup_acquisition` and `clear_acquisition` method, which are slow, but rarely called. Usually, they only need to be invoked right after connecting to the camera, or when the acquired image size is changed (e.g., due to a change in binning or ROI). Since these methods deal with buffer allocation, in almost all cases they take a parameter specifying buffer size (typically called `nframes`).

The other two procedures correspond to `start_acquisition` and `stop_acquisition` methods. These try to be as fast as possible, as they need to be called any time the acquisition is started or stopped, or when minor parameters (frame rate, exposure, trigger mode) are called.

## Region of interest (ROI) and binning

Most cameras allow the user to select only a part of the whole sensor for readout and transfer. Since the readout speed is usually the factor limiting the frame rate, selecting smaller ROI frequently lets you achieve higher frame rate. In addition, it also reduces the size of the frame buffer and the data transfer load. Same goes for binning: many cameras can combine values of several consecutive pixels in the same row or column (or both), which results in smaller images and, depending on the camera architecture, higher signal-to-noise ratio compared to binning in post-processing. Much less frequently you can set up subsampling instead of binning, which skips pixels instead of averaging them together.

Both operations depend very strongly on the exact hardware, so there are typically many associated restriction. The most common are minimal sizes in width and height, positions and sizes being factors of some power of 2 (up to 32 for some cameras), or equal binning for both axes. Device classes will typically round the ROI to the nearest allowed value. Furthermore, the scaling of the maximal frame rate with the ROI size is also hardware-dependent; for example, in many sCMOS chips readout speed only depends on the vertical extent, since the readout is done simultaneously for the whole row. In most cases, it takes some experiments to get a hang of the camera behavior.

## Exposure and frame rate

Almost all scientific cameras let user change the exposure, typically in a wide range (down to sub-ms). Frequently they also allow to separately change the frame period (inverse of the frame rate). Usually (but not always) the minimal frame period is set by the exposure plus some readout time, which depends on the ROI and some additional parameters such as pixel clock or simultaneous readout mode. Usually exposure takes priority over the frame period, i.e., if the frame period is set too short, it is automatically adjusted. Notable exception from this rule is [Uc480](#) interface, where this dependence is reversed.

## Triggering

Usually the cameras will have several different options for triggering, i.e., choosing when to start acquiring a new frame or a new batch of frames. The default option is the internal trigger, which means that the internal timer generates trigger event at a constant rate (frame rate). Many cameras will also take an external trigger signal to synchronize acquisition to external events or other cameras. Typically, a rising edge from 0 to 5V on the input will initiate the frame acquisition, but more exotic options (different polarities or levels, exposure control with pulse width, line-readout trigger) can be present.

## Application notes and examples

Here we talk more practically about performing tasks common to most cameras.

## Simple acquisition

Frame acquisition is, understandably, the most important part of the camera. Basic acquisition can be done without explicitly setting up the acquisition loop, simply by using `ICamera.snap()` and `ICamera.grab()` methods which, correspondingly, grab a single frame or a given number of frames:

```
from pylablib import Andor
cam = Andor.AndorSDK3Camera() # connect to the camera
img = cam.snap() # grab a single frame
images = cam.grab(10) # grab 10 frames (return a list of frames)
cam.close()
```

These allow for quick tests of whether the camera works properly, and for occasional frames acquisition. However, these methods have to start and stop acquisition every time they are called, which for some cameras can take about a second. Hence, if continuous acquisition and high frame rate are required, you would need to set up the acquisition loop.

## Acquisition loop

A typical simple acquisition loop has already been shown above:

```
# nframes=100 relates to the size of the frame buffer; the acquisition will continue_
↳ indefinitely
cam.setup_acquisition(mode="sequence", nframes=100) # could be combined with start_
↳ acquisition, or kept separate
cam.start_acquisition()
while True: # acquisition loop
    cam.wait_for_frame() # wait for the next available frame
    frame = cam.read_oldest_image() # get the oldest image which hasn't been read yet
    # ... process frame ...
    if time_to_stop:
        break
cam.stop_acquisition()
```

It relies on 3 sets of methods. First, starting and stopping acquisition using `start_acquisition` and `stop_acquisition`. As explained above, one also has an option to setup the acquisition first using `setup_acquisition`, which makes the subsequent `start_acquisition` call faster. However, one can also supply the same setup parameters to `start_acquisition` method, which automatically sets up the acquisition if it is not set up yet, or if any parameters are different from the current ones.

Second are the methods for checking on the acquisition process. The method used above is `wait_for_frame`, which by default waits until there is at least one unread frame in the buffer (i.e., it exits immediately if there is already a frame available). Its arguments modify this behavior by changing the point from which the new frame is acquired (e.g., from the current call), or the minimal required number of frames. Alternatively, there is a method `get_new_images_range`, which returns a range of the frame indices which have been acquired but not read. This method allows for a quick check of a number of unread frames without pausing the acquisition.

Finally, there are methods for reading out the frames. The simplest method is `read_oldest_image`, which return the oldest image which hasn't been read yet, and marks it as read. A more powerful is the `read_multiple_images` method, which can return a range of images (by default, all unread images). Both of these methods also take a `peek` argument, which allows one to read the frames without marking them as read.

## Returned frame format

`ICamera.read_multiple_images()` method described above has several different formats for return-

ing the frames, which can be controlled using `ICamera.set_frame_format()` and checked `ICamera.get_frame_format()`. The default format is "list", which returns a list of individual frames. The second possibility is "array", which returns a single 3D numpy array with all the frames. Finally, "chunks" returns a list of 3D arrays, each containing several consecutive frames.

While "chunks" format is the hardest to work with, it provides the best performance. First, it does not require any extra memory copies, which negatively affect performance at very high data rates, above ~1Gb/s. Second, it can combine multiple small frames together into a single array, which makes further processing faster, as it does not require explicit Python loop over every frame. This usually becomes important at frames rates above ~10kFPS, where treating each frame as an individual 2D array leads to significant overhead.

## Frame indexing

Different areas and libraries adopt different indexing convention for 2D arrays. The two most common ones are coordinate-like `xy` (the first index is the `x` coordinate, the second is `y` coordinate, and the origin is in the lower left corner) and matrix-like `ij` (the first index is row, the second index is column, the origin is in the upper right corner). Almost all cameras adopt the `ij` convention. The only exception is Andor SDK2, which uses similar row-column indexing, but counting from the bottom.

By default, the frames returned by the camera are indexed in the preferred convention, to reduce the overhead on re-indexing the frames. It is possible to check and change it using `ICamera.get_image_indexing()` and `ICamera.set_image_indexing()` methods:

```
>> cam.set_roi(0,256,0,128) # 256px horizontally, 128 vertically
>> cam.snap().shape # 128 rows, 256 columns
(128, 256)
>> cam.set_image_indexing("xyb") # standard xy indexing, starting from the bottom
>> cam.snap().shape
(256, 128)
```

## ROI, detector size and frame shape

Both ROI and binning are controlled by one pair of methods `get_roi` and `set_roi` which, depending on whether camera supports binning, take (and return) 4 or 6 arguments: start and stop positions of ROI along both axes and, optionally, binning along the axes:

```
cam.set_roi(0,128,0,256) # set 128x256px (width x height) ROI in the (typically)
↳upper left controller
cam.set_roi(0,128) # set roi with 128px width and full height (non-supplied
↳arguments take extreme values)
cam.set_roi(0,128,0,128,2,2) # set 128x128px ROI with 2x2 binning; the resulting
↳image size is 64x64
```

Regardless of the frame indexing, the first pair of arguments always controls horizontal span, the second pair controls vertical span, and the last pair controls horizontal and vertical binning (if applicable).

In addition, there is a couple of methods to acquire the detector and frame size. The first method is `get_detector_size`. It always returns the full camera detector size as a tuple (`width`, `height`) and, therefore, is not affected by ROI, binning, and indexing. The second method is `get_data_dimensions`, which returns the shape of the returned frame given the currently set up indexing. The results of this method do depend on the ROI, binning, and indexing:

```
>> cam.get_detector_size() # (width, height)
(2560, 1920)
```

(continues on next page)



(continued from previous page)

```
>> cam.get_data_dimensions() # (rows, columns), i.e., (height, width)
(1920, 2560)

>> cam.set_roi(0,256,0,128,2,2) # 256px horizontally, 128 vertically, 2x2 binning
>> cam.get_detector_size() # unaffected
(2560, 1920)
>> cam.get_data_dimensions() # depends on ROI
(64, 128)

>> cam.set_image_indexing("xyb")
>> cam.get_detector_size() # unaffected
(2560, 1920)
>> cam.get_data_dimensions() # depends on indexing
(128, 64)
```

## Exposure and frame period

In pylablib these parameters are normally controlled by `get_exposure/set_exposure` and, correspondingly `get_frame_period/set_frame_period` methods. In addition, `get_frame_timings` method provides an overview of all the relevant times. Exposure typically takes priority over frame period: if the frame period is set too small, it becomes the smallest possible for the given exposure; at the same time, if the exposure is set too big, it is still applied, and the frame period becomes the smallest possible with this exposure:

```
>> cam.get_frame_timings() # frame period is a usually bit larger due to the readout_
↳time
TAcqTimings(exposure=0.1, frame_period=0.12)

>> cam.set_exposure(0.01)
>> cam.get_frame_timings() # smaller exposure is still compatible with this frame_
↳period
TAcqTimings(exposure=0.01, frame_period=0.12)

>> cam.set_frame_period(0) # effectively means "set the highest possible frame rate"
>> cam.get_frame_timings()
TAcqTimings(exposure=0.01, frame_period=0.03)

>> cam.set_exposure(0.2)
>> cam.get_frame_timings() # frame period is increased accordingly
TAcqTimings(exposure=0.2, frame_period=0.22)
```

There are exceptions for some camera types, which are discussed separately.

## Camera attributes

Some camera interfaces, e.g., *Thorlabs Scientific Cameras*, *PCO SC2*, or *NI IMAQ* are fairly specific, and only apply to a handful of devices with very similar capabilities. In this case, pylablib usually attempts to implement as much of the functionality as possible given the available hardware, and to present it via the camera object methods.

In other cases, e.g., *NI IMAQdx*, *Andor SDK3*, or *DCAM*, the same interface deals with many fairly different cameras. This is especially true for *IMAQdx*, which covers hundreds of cameras from dozens of manufacturers, all with very different capabilities and purpose. Since managing such cameras can not usually be conformed to a small set of functions, it is implemented through camera attributes mechanism. That is, for each camera the interface defines a set of attributes (sometimes also called properties or features), which can be queried or set by their names, and whose exact meaning and possible values depend on the specific camera.



Typically, cameras dealing with attributes will implement `IAttributeCamera.get_attribute_value()` and `IAttributeCamera.set_attribute_value()` for querying and setting the attributes, as well as dictionary-like `.cav` (stands for “camera attribute value”) interface to do the same thing:

```
>> cam = Andor.AndorSDK3Camera()
>> cam.get_attribute_value("CameraAcquiring") # check if the camera is acquiring
0
>> cam.set_attribute_value("ExposureTime", 0.1) # set the exposure to 100ms
>> cam.cav["ExposureTime"] # get the exposure; could also use cam.get_attribute_
↪value("ExposureTime")
0.1
```

Additionally, there are `IAttributeCamera.get_all_attribute_values()` and `IAttributeCamera.set_all_attribute_values()` which get and set all camera attributes (possibly only within the given branch, if camera attributes form a hierarchy). Finally, methods `IAttributeCamera.get_attribute()` and `IAttributeCamera.get_all_attributes()`, together with the corresponding `.ca` interface, allow to query specific attribute objects, which provide additional information about the attributes: whether they are writable or readable, their range, description, possible values, types, etc.:

```
>> cam = DCAM.DCAMCamera()
>> attr=cam.ca["EXPOSURE TIME"] # get the exposure attribute
DCAMAttribute(name='EXPOSURE TIME', id=2031888, min=0.001, max=10.0, unit=1)
>> attr.max
10.0
>> attr.set_value(0.1) # same as cam.set_attribute_value("EXPOSURE TIME", 0.1)
```

Note that, depending on the camera, the attribute properties (especially minimal and maximal value) can depend on the other camera attributes. For example, minimal exposure can depend on the frame size:

```
>> cam = DCAM.DCAMCamera()
>> attr=cam.ca["EXPOSURE TIME"] # get the exposure attribute
DCAMAttribute(name='EXPOSURE TIME', id=2031888, min=0.001, max=10.0, unit=1)
>> attr.min
0.001
>> cam.set_roi(0, 0, 0, 0) # set the minimal possible ROI
(0, 4, 0, 4, 1, 1)
>> attr.min # minimal value hasn't been updated yet
0.001
>> attr.update_limits() # update the property limits
>> attr.min # now the minimal possible exposure is smaller
7.795e-05
```

If the documentation is not available (as is the case for, e.g., some IMAQdx cameras), the best way to learn about the attributes is to use the native software (whenever available) to modify camera settings and then check how the attributes change. Besides that, it is always useful to check attribute description (available for IMAQdx parameter), their range, and the available values for enum attributes.

## Trigger setup

The trigger is usually set up using `set_trigger_mode` method, although it might be different if more specialized modes are used. When external trigger is involved, most of the code (such as acquisition set up and start) stays the same. The only difference is the rate at which the frames are generated:

```
frame = cam.snap() # starts acquiring immediately, returns after a single frame
cam.set_trigger_mode("ext") # set up the trigger mode
```

(continues on next page)

(continued from previous page)

```
frame = cam.snap()
# after cam.snap() is called, the execution will wait
# for an external trigger pulse to acquire the frame and return
```

## Frame metainfo

Many cameras supply additional information together with the frames. Most frequently it contains the internal frames-tamp and timestamp (which are useful for tracking missing frames), but sometimes it also includes additional information such as frame size or location, status, or auxiliary input bits. To get this information, you can supply the argument `return_info=True` to the `read_multiple_images` method. In this case, instead of a single list of frames, it will return a tuple of two lists, where the second list contains this metainfo.

There are several slightly different metainfo formats, which can be set using `ICamera.set_frame_info_format()` method. The default representation is a (possibly nested) named tuple, but it is also possible to represent it as a flat list, flat dictionary, or a numpy array. The exact structure and values depend on the camera.

Keep in mind, that for some camera interfaces (e.g., *Uc480* or *Silicon Software*) obtaining the additional information might take relatively long, even longer than the proper frame readout. Hence, at higher frame rates it might become a bottleneck, and would need to be turned off.

## Related projects

Pylabelib *cam-control* is a standalone software package which builds on camera classes included in pylablib. It provides an easy way to detect and control many different cameras and acquire their data. In addition, it supports custom on-line image processing, flexible data acquisition, and control by external software using a TCP/IP server.

## Currently supported cameras

- *Andor SDK2* and *Andor SDK3*: variety of Andor (currently part of Oxford Instruments) cameras. Tested with Andor iXon, Luca, Zyla, Neo and Marana.
- *Allied Vision Bonito* cameras: CameraLink-interfaced cameras. Tested with Bonito CL-400B/C and NI IMAQ frame grabber.
- *Basler*: Basler pylon-compatible cameras. Tested with an emulated Basler camera.
- *BitFlow*: BitFlow Axion family frame grabbers. Tested with BitFlow Axion 1xB frame grabber together with PhotonFocus MV-D1024E camera.
- *DCAM*: Hamamatsu cameras. Tested with Hamamatsu Orca Flash 4.0 and Imagem.
- *NI IMAQ*: National Instruments frame grabbers. Tested with NI PCI-1430 and PCI-1433 frame grabbers together with PhotonFocus MV-D1024E camera.
- *NI IMAQdx*: National Instruments universal camera interface. Tested with Ethernet-connected PhotonFocus HD1-D1312 camera.
- *Photon Focus*: Photon Focus pfcam interface. Tested with PhotonFocus MV-D1024E camera connected through either NI frame grabbers (PCI-1430 and PCI-1433) or Silicon Software frame grabbers (microEnable IV AD4-CL).
- *PCO SC2*: PCO cameras. Tested with pco.edge cameras with CLHS and regular CameraLink interfaces.
- *Picam*: Princeton Instruments cameras. Tested with a PIXIS 400 camera.

- *PVCAM*: Photometrics cameras. Tested with a Prime 95B camera.
- *Silicon Software*: Silicon Software frame grabbers. Tested with microEnable IV AD4-CL frame grabbers together with PhotonFocus MV-D1024E camera.
- *Thorlabs Scientific Cameras*: Thorlabs sCMOS cameras. Tested with Thorlabs Kiralux camera.
- *Uc480/uEye*: multiple cameras, including simple Thorlabs and IDS cameras. Tested with IDS SC2592R12M and Thorlabs DCC1545M.

---

**Note:** General camera communication concepts are described on the corresponding [page](#)

---

## Andor cameras

Andor implements two completely separate interfaces for different cameras. The older one, called SDK2, or simply SDK, provides interface for the older cameras: iXon, iKon, iStart, iDus, iVac, Luca, Newton. The details of this SDK are available in the [manual](#).

The newer SDK, called SDK3, covers newer cameras: Zyla, Neo, Apogee, Sona, Marana, and Balor. The [manual](#) describes the cameras and capabilities in more details.

The required DLLs are distributed with [Andor Solis](#) or the corresponding [Andor SKD](#). In most cases, you have Andor Solis already installed to provide the drivers and to communicate with the cameras to begin with.

## Andor SDK 2

This is an older SDK, which mainly involves older cameras. It has been tested with Andor iXon and Andor Luca.

The code is located in `pylablib.devices.Andor`, and the main camera class is `pylablib.devices.Andor.AndorSDK2Camera`.

## Software requirements

The required DLL can have different names depending on the Solis version and SDK bitness. For 64-bit version it will be called `atmcd64d.dll` or `atmcd64d_legacy.dll`. For 32-bit version, correspondingly, `atmcd32d.dll` or `atmcd32d_legacy.dll`. By default, library searches for DLLs in Andor Solis and Andor SDK folder in Program Files folder (or Program files (x86), if 32-bit version of Python is running), as well as in the folder containing the script. If the DLLs are located elsewhere, the path can be specified using the library parameter `devices/dlls/andor_sdk2`:

```
import pylablib as pll
pll.par["devices/dlls/andor_sdk2"] = "path/to/dlls"
from pylablib.devices import Andor
cam = Andor.AndorSDK2Camera()
```

## Connection

The cameras are identified by their index, starting from zero. To get the total number of cameras, you can run `Andor.get_cameras_number_SDK2`:

```
>> from pylablib.devices import Andor
>> Andor.get_cameras_number_SDK2()
2
>> cam1 = Andor.AndorSDK2Camera(idx=0)
>> cam2 = Andor.AndorSDK2Camera(idx=1)
>> cam1.close()
>> cam2.close()
```

**Warning:** It is important to close all camera connections before finishing your script. Otherwise, DLL resources might become permanently blocked, and the only way to solve it would be to restart the PC.

## Operation

The operation of these cameras is relatively standard. They support all the standard methods for dealing with ROI and exposure, starting and stopping acquisition, and operating the frame reading loop. However, there's a couple of differences from the standard libraries worth highlighting:

- Since the manufacturer DLLs do not provide methods to get most of the camera parameters (such as exposure or ROI), it is impossible to know them when connecting the camera. To get around it, the camera is put into a “default” state any time the connection is opened.
- When applicable, it is important to properly set the cooling setpoint and the fan mode. By default, the fan is turned off, and the cooling is set to the 20'th percentile of the whole range (e.g., -80C for Andor iXon). It is possible to pass these parameters on camera creation:

```
cam = Andor.AndorSDK2Camera(temperature=-50, fan_mode="on")
```

- Often cameras have a lot of different readout parameters: channel, amplifier, vertical and horizontal scan speed, etc. These parameters greatly affect the camera sensitivity and readout speed. Upon the connection, the parameter are typically set to the slowest mode. To get the list of all possible parameter combinations, you can use `AndorSDK2Camera.get_all_amp_modes()` and `AndorSDK2Camera.get_max_vsspeed()`. Afterwards, you can set them using `AndorSDK2Camera.set_amp_mode()` and `AndorSDK2Camera.set_vsspeed()`.
- The default shutter parameter is "closed". This preserves camera from possible high illumination, but can lead to confusion, if you expect to see some image.
- This SDK does not allow for specifying number of frames in the frames buffer. However, the parameters chosen by the SDK are usually reasonable (at least a second worth of acquisition).
- Some cameras (e.g., iXon) have lots of readout (full frame, ROI, full vertical binning, etc.) and acquisition modes (single, continuous, accumulating, kinetic cycle, etc.). They are described in details in the [manual](#).

## Andor SDK 3

This is a newer SDK, which covers the newer cameras. It has been tested with Andor Zyla, Neo and Marana.

The code is located in `pylablib.devices.Andor`, and the main camera class is `pylablib.devices.Andor.AndorSDK3Camera`.

## Software requirements

This library requires several DLLs all located in the same folder: `atcore.dll`, `atblkbx.dll`, `atcl_bitflow.dll`, `atdevapogee.dll`, `atdevregcam.dll`, `atusb_libusb.dll`, `atusb_libusb10.dll`. Same as for SDK2, pylablib looks for DLLs in Andor Solis and Andor SDK3 folders in Program Files folder (or Program files (x86), if 32-bit version of Python is running), as well as in the folder containing the script. A custom DLLs path can be specified using the library parameter `devices/dlls/andor_sdk3`:

```
import pylablib as pll
pll.par["devices/dlls/andor_sdk3"] = "path/to/SDK3/dlls"
from pylablib.devices import Andor
cam = Andor.AndorSDK3Camera()
```

## Connection

The cameras are identified by their index, starting from zero. To get the total number of cameras, you can run `Andor.get_cameras_number_SDK3`:

```
>> from pylablib.devices import Andor
>> Andor.get_cameras_number_SDK3()
2
>> cam1 = Andor.AndorSDK3Camera(idx=0)
>> cam2 = Andor.AndorSDK3Camera(idx=1)
>> cam1.close()
>> cam2.close()
```

## Operation

The operation of these cameras is also relatively standard. They support all the standard methods for dealing with ROI and exposure, starting and stopping acquisition, and operating the frame reading loop. However, there's a couple of differences from the standard libraries worth highlighting:

- The SDK also provides a universal interface for getting and setting various *camera attributes* (called “features” in the documentation) using their name. You can use `AndorSDK3Camera.get_attribute_value()` and `AndorSDK3Camera.set_attribute_value()` for that, as well as `.cav` attribute which gives a dictionary-like access:

```
>> cam = Andor.AndorSDK3Camera()
>> cam.get_attribute_value("CameraAcquiring") # check if the camera is acquiring
0
>> cam.set_attribute_value("ExposureTime", 0.1) # set the exposure to 100ms
>> cam.cav["ExposureTime"] # get the exposure; could also use cam.get_attribute_
↪value("ExposureTime")
0.1
```

Some values serve as commands; these can be invoked using `AndorSDK3Camera.call_command()` method. To see all available attributes, you can call `AndorSDK3Camera.get_all_attributes()` to get a dictionary with attribute objects, and `AndorSDK3Camera.get_all_attribute_values()` to get the dictionary of attribute values. The attribute objects provide additional information: their kind, whether they are implemented, readable, or writable, what are their limits or possible values, etc:

```
>> cam = Andor.AndorSDK3Camera()
>> attr = cam.get_attribute("SensorTemperature")
```

(continues on next page)

(continued from previous page)

```
>> attr.readable
True
>> attr.writable
False
>> (attr.min, attr.max)
(-100.0, 50.0)
```

The description of the attributes is given in [manual](#).

- USB cameras can, in principle, generate data at higher rate than about 320Mb/s that the USB3 bus supports. For example, Andor Zyla with 16 bit readout has a single full frame size of 8Mb, which puts the maximal USB throughput at about 40FPS. At the same time, the camera itself is capable of reading up to 100FPS at the full frame. Hence, it is possible to overflow the camera internal buffer (size on the order of 1Gb) regardless of the PC performance. If this happens, the acquisition process halts and needs to be restarted. You can check the number of buffer overflows using `AndorSDK3Camera.get_missed_frames_status()`, and reset this counter using `AndorSDK3Camera.reset_overflows_counter()`; the counter is also automatically resets on acquisition clearing, but not stopping.

Furthermore, the class implements different strategies when encountering overflow while waiting for a new frame. The specific strategy is selected using `AndorSDK3Camera.set_overflow_behavior()`, and it can be "error" (raise `AndorFrameTransferError`, which is the default behavior), "restart" (restart the acquisition and immediately raise timeout error), or "ignore" (ignore the overflow, which will eventually lead to a timeout error, as the new frames are no longer generated).

---

**Note:** General camera communication concepts are described on the corresponding [page](#).

---

## Allied Vision Bonito cameras

Allied Vision manufactures a variety of cameras with different interfaces: USB, GigE, and CameraLink. Currently, only CameraLink Bonito cameras using NI IMAQ frame grabber are supported. It has been tested with Bonito CL-400B/C and NI IMAQ frame grabber.

The code is located in `pylablib.devices.AlliedVision`, and the main camera class is `pylablib.devices.AlliedVision.BonitoIMAQCamera`.

## Software requirements

Since the camera control is done purely through the frame grabber interface, the requirements are the same as for generic [IMAQ cameras](#). However, the correct camera file still needs to be specified to determine the correct serial communication parameters (especially the termination character)

## Connection

The cameras are identified by their name, which usually looks like "img0". To get the list of all cameras, you can use NI MAX (Measurement and Automation Explorer), or `IMAQ.list_cameras()`:

```
>> from pylablib.devices import IMAQ, AlliedVision
>> IMAQ.list_cameras()
['img0']
>> cam = AlliedVision.BonitoIMAQCamera('img0')
>> cam.close()
```

## Operation

The operation of these cameras is relatively standard. They support all the standard methods for dealing with ROI and exposure, starting and stopping acquisition, and operating the frame reading loop. However, there's a couple of differences from the standard libraries worth highlighting:

- `Bonito.BonitoIMAQCamera` supports all of `IMAQ.IMAQCamera` features, such as trigger control and fast buffer acquisition. Some methods have been modified to make them more convenient: e.g., `Bonito.BonitoIMAQCamera.set_roi()` method sets the camera ROI and automatically adjusts the frame grabber ROI to match.
- Internally the camera only supports vertical ROI (number of rows), so the horizontal ROI is set via the frame grabber. This means that regardless of the horizontal ROI settings the whole rows are always transmitted between the camera and the frame grabber, so it does not affect, e.g., the maximal frame rate.
- The camera supports a status line, which replaces the first 8 pixels in the upper row encoded frame number. You can use `AlliedVision.Bonito.get_status_lines()` function to identify and extract the data in the status lines from the supplied frames. Note that due to the full row transfer mentioned earlier, the status line is only available if the horizontal ROI span starts from zero; otherwise, it will be partially or completely cut off.
- You can use the function `AlliedVision.Bonito.check_grabber_association()` to check if the given IMAQ camera is a Bonito model by sending several standard Bonio commands and checking replies.

---

**Note:** General camera communication concepts are described on the corresponding [page](#)

---

## Basler cameras interface

Basler manufactures a wide variety of cameras, which implement GenICam-based interface through its pylon API. It has been tested with pylon-provided emulated camera camera.

The code is located in `pylablib.devices.Basler`, and the main camera class is `pylablib.devices.Basler.BaslerPylonCamera`.

## Software requirements

These cameras require `PylonC_vX_Y.dll` (where X and Y is the pylon version, e.g., `PylonC_V7_1.dll`), which is installed with the freely available upon registration [Basler pylon Camera Software Suite](#) (the current latest version is 7.1.0). After installation, the path to the DLL (for pylon 7.1.0 located by default in `Basler/pylon 7/Runtime/x64` folder in `Program Files`) is automatically added to system `PATH` variable, which is one of the places where `pylablib` looks for it by default. If the DLLs are located elsewhere, the path (either to the DLL file, or to the containing folder) can be specified using the library parameter `devices/dlls/basler_pylon`:

```
import pylablib as pll
pll.par["devices/dlls/basler_pylon"] = "path/to/dlls"
from pylablib.devices import Basler
cam = Basler.BaslerPylonCamera()
```

## Connection

The cameras are identified either by their index among the present cameras (starting from 0), or by their name. To get the list of all cameras, you can use `pylon Viewer`, or `Basler.list_cameras`:

```
>> from pylablib.devices import Basler
>> Basler.list_cameras()
[TCameraInfo(name='Emulation (0815-0000)', model='Emulation', serial='0815-0000',
↳ devclass='BaslerCamEmu', devversion='', vendor='Basler', friendly_name='Basler_
↳ Emulation (0815-0000)', user_name='', props={'DeviceFactory': 'CamEmu/BaslerCamEmu_
↳ 7.1.0.19126', 'InterfaceID': 'DefaultInterface', 'TLType': 'CamEmu'})]
>> cam = Basler.BaslerPylonCamera() # by default, connect to the first available_
↳ camera
>> cam.close()
>> cam = Basler.BaslerPylonCamera(name="Emulation (0815-0000)")
>> cam.close()
```

## Operation

The operation of these cameras is relatively standard. They support all the standard methods for dealing with ROI, starting and stopping acquisition, and operating the frame reading loop. The SDK also provides a universal interface for getting and setting various *camera attributes* using their name. You can use `BaslerPylonCamera.get_attribute_value()` and `BaslerPylonCamera.set_attribute_value()` for that, as well as `.cav` attribute which gives a dictionary-like access:

```
>> cam = Basler.BaslerPylonCamera()
>> cam.get_attribute_value("StatusInformation/AcqInProgress") # check if the camera_
↳ is acquiring
0
>> cam.set_attribute_value("Width", 512) # set the ROI width to 512px
>> cam.cav["Width"] # get the exposure; could also use cam.get_attribute_value("Width
↳ ")
512
```

To see all available attributes, you can call `BaslerPylonCamera.get_all_attributes()` to get a dictionary with attribute objects, and `BaslerPylonCamera.get_all_attribute_values()` to get the dictionary of attribute values. The attribute objects provide additional information: attribute kind (integer, enum, string, etc.), range (either numerical range, or selection of values for enum attributes), description string, etc.:

```
>> cam = Basler.BaslerPylonCamera()
>> attr = cam.get_attribute("Width")
>> attr.description
'This value sets the width of the area of interest in pixels.'
>> attr.writable
True
>> (attr.min, attr.max)
(1, 4096)
```

Since these properties vary a lot between different cameras, it is challenging to write a universal class covering a large range of cameras. Hence, currently the universal class only has the basic camera parameter control such as ROI (without binning), acquisition status, and exposure (if present). For many specific cameras you might need to explore the attributes tree (either using the Python class and, e.g., a console, or via pylon Viewer) and operate them directly in your code.

## Known issues

- Currently only the basic unpacked monochrome pixel formats are supported: Mono8, Mono10, Mono12, Mono16, and Mono32. The reason is that even nominally well-defined types (e.g., Mono12Packed) have different formats for different cameras. Currently any unsupported format will raise an error on readout by



default. It is still possible to read these out as raw frame data in the form of 1D or 2D numpy 'u1' array by enabling raw frame readout using `BaslerPylonCamera.enable_raw_readout()` method:

```
>> cam = Basler.BaslerPylonCamera()
>> cam.get_detector_size() # 1024px x 1024px frame
(1024, 1024)
>> cam.set_attribute_value("PixelFormat", "BGRA8Packed") # unsupported format
>> cam.snap().shape
...
BaslerError: pixel format BGRA8Packed is not supported
>> cam.enable_raw_readout("frame") # frame data is returned as a flat array
>> cam.snap().shape # 1024 * 1024 * 4 = 4194304 bytes
(1, 4194304)
```

---

**Note:** General camera communication concepts are described on the corresponding [page](#)

---

## BitFlow Axion frame grabbers interface

BitFlow manufactures several kinds of camera interface cards, including CameraLink. Currently, only newer CameraLink Axion family is supported. It has been tested with NI BitFlow Axion 1xB frame grabbers together with PhotonFocus MV-D1024E camera.

The code is located in `pylablib.devices.BitFlow`, and the main camera class is `pylablib.devices.BitFlow.BitFlowCamera`.

## Software requirements

This interface requires two pieces of software, both freely available on the [BitFlow website](#). First, you need [BitFlow SDK 6.5](#), which also includes all the necessary drivers. The free version does not provide any headers and documentation to the DLLs, so you use it you also need to install the manufacturer-provided Python packages, either for [Python 3.6.6](#), or for [Python 3.8.10](#). Note that only these two Python versions are officially supported.

After installation, the DLL locations are automatically added to the `PATH` environment variable. To facilitate proper package import and DLL loading on Python 3.8, it is recommended to install BitFlow SDK into its default library, or at least leave BitFlow in the folder name.

## Connection

The cameras are identified by their index, starting from 0. To get the list of all cameras, you can use `BitFlow.list_cameras()`:

```
>> from pylablib.devices import BitFlow
>> BitFlow.list_cameras()
['img0', 'img1']
>> cam = BitFlow.BitFlowCamera('img0')
>> cam.close()
```

## Operation

Unlike most camera classes, the frame grabber interface only deals with the frame transfer between the camera and the PC over the CameraLink interface. Therefore, it can not directly control camera parameters such as exposure, frame

rate, triggering, ROI, etc. Some similar-looking parameters are still present, but they have a different meaning:

- ROI is defined within the transferred image, whose size itself is determined by the camera ROI. Hence, e.g., if the camera chip is 1024x1024px and its roi is 512x512, then the frame grabber ROI can go only up to 512x512. Any attempts to set it higher result in the frozen acquisition, as the frame grabber expects a larger frame than it receives, and waits forever to get the rest.

### Fast buffer readout and frames merging

At high frame rates (above ~10kFPS) dealing with each frame individually becomes too slow for Python. Hence, there is an option to read out and process frames in larger ‘chunks’, which are 3D numpy arrays with the first axis enumerating the frame index. This approach leverages the ability to store several frame buffers in the contiguous memory locations (resulting in a single 3D array), and it essentially eliminates the overhead for dealing with multiple frames at high frame rates, as long as the total data rate is manageable (typically below 600Mb/s).

This option can be accessed by calling using `BitFlowCamera.set_frame_format()` method to set frames format to "chunks". In this case, instead of a list of individual frames (which is the standard behavior), the method returns list of chunks of varying size, which contain several consecutive frames.

On top of that, due to unavoidable Python loop required by the BitFlow Python interface, the frame rate is usually limited to about 2-4kFPS. However, there is a way to overcome this by merging  $n$  consecutive frames to a single “super-frame” with  $n$  times larger height. This merging can be specified by `frame_merge` parameter in the `BitFlowCamera.setup_acquisition()` or `BitFlowCamera.start_acquisition()` methods (by default it is 1, meaning no merging). Adjusting the frame grabber ROI and splitting the resulting files is done transparently for the user; the only difference is that frames always arrive in batches, e.g., with `frame_merge=10` and 10FPS rate the frames will arrive once a second in batches of 10. Therefore, it makes sense to adjust the merging to keep the “merged” frame rate high enough for real-time operations but lower than the 2kFPS limit (e.g., around 100FPS).

### Communication with the camera and camera files

The frame grabber needs some basic information about the camera: sensor size, bit depth, data transfer format, time-outs, aux lines mapping, etc. This information is contained in the so-called camera files, which for Axion cameras have `.bfml` extension. These files can be assigned to cameras using `SysReg` utility located in the `Bin64` folder of your BitFlow installation (by default, `C:\BitFlow SDK 6.5\Bin64`).

In addition, due to limitations of the provided Python interfaces, some operations such as changing ROI and bitness can only be done by altering the camera file. Hence, there is an option to create a temporary camera file and alter it to control these parameters. However, it needs the original camera file to serve as a template (this original file is only used as source and not modified). Since there is no possibility to get a path to this file within the Python interface, it should be provided using `camfile` parameter upon creation.

### Known issues

- As mentioned above, ROI is defined within a frame transferred by the camera. Hence, if it includes pixels with positions outside of the transferred frame, the acquisition will time out. For example, suppose the camera sensor is 1024x1024px, and the *camera* ROI is selected to be central 512x512 region. As far as the frame grabber is concerned, now the camera sensor size is 512x512px. Hence, if you try to set the same *frame grabber* ROI (i.e., 512x512 starting at 256,256), it will expect at least 768x768px frame. Since the frame is, actually, 512x512px, the acquisition will time out. The correct solution is to set frame grabber ROI from 0 to 512px on both axes. In general, it is a good idea to always follow this pattern: control ROI only on camera, and always set frame grabber ROI to cover the whole transfer frame.

---

**Note:** General camera communication concepts are described on the corresponding [page](#).

---

## DCAM cameras interface

DCAM is the interface used in Hamamatsu cameras. It has been tested with Hamamatsu Orca Flash and ImagEM.

The code is located in `pylablib.devices.DCAM`, and the main camera class is `pylablib.devices.DCAM.DCAMCamera`.

## Software requirements

These cameras require `dcamapi.dll`, which is installed with most of Hamamatsu software (such as HoKaWo or HiPic), as well as with the freely available [DCAM API](#), which also includes all the necessary drivers. Keep in mind, that you also need to install the drivers for required corresponding camera type (USB, Ethernet, IEEE 1394). These drivers are in the same installer, but need to be installed separately. You should also pay attention to the cameras supported by the given DCAM driver version, since newer version do not support older cameras (e.g., ImageEM C9100 cameras are only supported up to version 15). After installation, the DLL is automatically added to the `System32` folder, where `pylablib` looks for it by default. If the DLL is located elsewhere, the path can be specified using the library parameter `devices/dlls/dcamapi`:

```
import pylablib as pll
pll.par["devices/dlls/dcamapi"] = "path/to/dlls"
from pylablib.devices import DCAM
cam = DCAM.DCAMCamera()
```

## Connection

The cameras are identified by their index, starting from zero. To get the total number of cameras, you can run `DCAM.get_cameras_number()`:

```
>> from pylablib.devices import DCAM
>> DCAM.get_cameras_number()
2
>> cam1 = DCAM.DCAMCamera(idx=0)
>> cam2 = DCAM.DCAMCamera(idx=1)
>> cam1.close()
>> cam2.close()
```

## Operation

The operation of these cameras is relatively standard. They support all the standard methods for dealing with ROI and exposure, starting and stopping acquisition, and operating the frame reading loop. The SDK also provides a universal interface for getting and setting various *camera attributes* (called “properties” in the documentation) using their name. You can use `DCAMCamera.get_attribute_value()` and `DCAMCamera.set_attribute_value()` for that, as well as `.cav` attribute which gives a dictionary-like access:

```
>> cam = DCAM.DCAMCamera()
>> cam.get_attribute_value("BINNING") # get the camera binning (no binning, by
↪ default)
```

(continues on next page)

(continued from previous page)

```
1
>> cam.set_attribute_value("EXPOSURE TIME", 0.1) # set the exposure to 100ms
>> cam.cav["EXPOSURE TIME"] # get the exposure; could also use cam.get_attribute_
↪value("EXPOSURE TIME")
0.1
```

To see all available attributes, you can call `DCAMCamera.get_all_attributes()` to get a dictionary with attribute objects, and `DCAMCamera.get_all_attribute_values()` to get the dictionary of attribute values, with an option of representing enum attributes either as text or as integer values. The attribute objects provide additional information: attribute range, step, and units:

```
>> cam = DCAM.DCAMCamera()
>> attr = cam.get_attribute("EXPOSURE TIME")
>> (attr.min, attr.max)
(0.001, 10.0)
```

Additionally, there's a couple of differences from the standard libraries worth highlighting:

- The library supports only symmetric binning, i.e., the binning factor is the same in both directions. For compatibility `DCAMCamera.get_roi()` and `DCAMCamera.set_roi()` still return and accept both binning parameters independently, but they are always the same when returned, and `vbin` is ignored when set.
- By default, the SDK does not provide independent control of the frame period and the exposure. Hence, `set_frame_period` method is unavailable, and the frame rate is defined solely by the exposure.

---

**Note:** General camera communication concepts are described on the corresponding [page](#)

---

## NI IMAQ frame grabbers interface

NI IMAQ is the interface from National Instruments, which is used in a variety of frame grabbers. It has been tested with NI PCI-1430 and PCI-1433 frame grabbers together with PhotonFocus MV-D1024E camera.

The code is located in `pylablib.devices.IMAQ`, and the main camera class is `pylablib.devices.IMAQ.IMAQCamera`.

## Software requirements

This interfaces requires `imaq.dll`, which is installed with the freely available [Vision Acquisition Software](#), which also includes all the necessary drivers. After installation, the DLL is automatically added to the `System32` folder, where `pylablib` looks for it by default. If the DLL is located elsewhere, the path can be specified using the library parameter `devices/dlls/niimaq`:

```
import pylablib as pll
pll.par["devices/dlls/niimaq"] = "path/to/dlls"
from pylablib.devices import IMAQ
cam = IMAQ.IMAQCamera()
```

## Connection

The cameras are identified by their name, which usually looks like `"img0"`. To get the list of all cameras, you can use NI MAX (Measurement and Automation Explorer), or `IMAQ.list_cameras()`:

```
>> from pylablib.devices import IMAQ
>> IMAQ.list_cameras()
['img0', 'img1']
>> cam1 = IMAQ.IMAQCamera('img0')
>> cam2 = IMAQ.IMAQCamera('img1')
>> cam1.close()
>> cam2.close()
```

## Operation

Unlike most camera classes, the frame grabber interface only deals with the frame transfer between the camera and the PC over the CameraLink interface. Therefore, it can not directly control camera parameters such as exposure, frame rate, triggering, ROI, etc. Some similar-looking parameters are still present, but they have a different meaning:

- External trigger controls frame *transfer*, not frame *acquisition*, which is defined by the camera. By default, when the internal frame grabber trigger is used, the frame grabber transfer rate is synchronized to the camera, so every frame gets transferred. However, if the external transfer trigger is used and it is out of sync with the camera, it can result in duplicate or missing frames.
- ROI is defined within the transferred image, whose size itself is determined by the camera ROI. Hence, e.g., if the camera chip is 1024x1024px and its roi is 512x512, then the frame grabber ROI can go only up to 512x512. Any attempts to set it higher result in the frozen acquisition, as the frame grabber expects a larger frame than it receives, and waits forever to get the rest.

The SDK also provides a universal interface for getting and setting various camera attributes using their name. You can use `IMAQCamera.get_grabber_attribute_value()` and `IMAQCamera.set_grabber_attribute_value()` for that:

```
>> cam = IMAQ.IMAQCamera()
>> cam.get_grabber_attribute_value("FRAMEWAIT_MSEC") # frame read request timeout
1000
```

To get all available attributes as a dictionary, you can call `IMAQCamera.get_all_grabber_attribute_values()`. Their meaning, as well as descriptions of trigger modes and other settings, is explained in the manual supplied with the [Vision Acquisition Software](#).

## Fast buffer readout mode

At high frame rates (above ~10kFPS) dealing with each frame individually becomes too slow for Python. Hence, there is an option to read out and process frames in larger ‘chunks’, which are 3D numpy arrays with the first axis enumerating the frame index. This approach leverages the ability to store several frame buffers in the contiguous memory locations (resulting in a single 3D array), and it essentially eliminates the overhead for dealing with multiple frames at high frame rates, as long as the total data rate is manageable (typically below 600Mb/s).

This option can be accessed by calling using `IMAQCamera.set_frame_format()` method to set frames format to "chunks" (former way of supplying `fastbuff=True` in `IMAQCamera.read_multiple_images()` is now deprecated). In this case, instead of a list of individual frames (which is the standard behavior), the method returns list of chunks about 1Mb in size, which contain several consecutive frames.

## Communication with the camera and camera files

The frame grabber needs some basic information about the camera: sensor size, bit depth, data transfer format, timeouts, aux lines mapping, etc. In NI MAQ this information is contained in the so-called camera files. These files can be

assigned to cameras in the NI MAX, and are usually supplied by NI or by the camera manufacturer. In addition, NI MAX allows one to adjust some settings within these files, which are read-only within the NI IMAQ software. These include frame timeout and camera bit depth.

The communication with the camera itself greatly varies between different cameras. Some will have additional connection to control the parameters. Others use serial communication built into the CameraLink interface. This communication can be set up with `IMAQCamera.setup_serial_params()` and used via `IMAQCamera.serial_read()` and `IMAQCamera.serial_write()`. The communication protocols are camera-dependent, and are frequently described in the camera manual. However, some other cameras (e.g., Photon Focus) use proprietary communication protocol. In this case, they provide their own DLLs, which independently use NI-provided DLLs for serial communication (most notably, `clallserial.dll`) to communicate with the camera. In this case, one needs to maintain two independent connections: one directly to the NI frame grabber to obtain the frame data, and one to the manufacturer library to control the camera. This is the way it is implemented in *PhotonFocus* camera interface.

### Known issues

- Sometimes when the acquisition is stopped and restarted without being cleared, the acquired frame counter does not refresh. This might show up as the software not reporting any new frames. It has been tracked down to a very low (~1ms) frame read timeout. Hence, it is recommended to keep this timeout at least at 500ms.
- If you are unable to access full camera sensor size, check the camera file (it can be opened in the text editor). `MaxImageSize` parameter defines the maximal allowed image size, and it should be equal to the camera sensor size.
- Same goes for bitness. If the camera bitness is higher than set up in the frame grabber, a single camera pixel gets treated as several pixels by the frame grabber, typically resulting in 1px-wide vertical stripes on the image. In the opposite case, the frame grabber expects more bytes than the camera sends, it never receives the full frame, and the acquisition times out.
- Keep in mind that as long as the frame grabber is accessed in NI MAX, it is blocked from use in any other software. Hence, you need to close NI MAX before running your code.
- As mentioned above, ROI is defined within a frame transferred by the camera. Hence, if it includes pixels with positions outside of the transferred frame, the acquisition will time out. For example, suppose the camera sensor is 1024x1024px, and the *camera* ROI is selected to be central 512x512 region. As far as the frame grabber is concerned, now the camera sensor size is 512x512px. Hence, if you try to set the same *frame grabber* ROI (i.e., 512x512 starting at 256,256), it will expect at least 768x768px frame. Since the frame is, actually, 512x512px, the acquisition will time out. The correct solution is to set frame grabber ROI from 0 to 512px on both axes. In general, it is a good idea to always follow this pattern: control ROI only on camera, and always set frame grabber ROI to cover the whole transfer frame.
- Some frame grabbers have a limit on the data transfer rate (for one model observed to be about 200 Mb/s). If the camera data generation rate exceeds it (e.g., it produces 1024x1024px 16-bit frames at >100FPS), then the camera will raise `IMG_ERR_FIFO` error shortly after the acquisition start. In this case, you will need to reduce the data rate by reducing the frame rate or frame size (through ROI, binning, or bitness).

---

**Note:** General camera communication concepts are described on the corresponding [page](#)

---

### NI IMAQdx cameras interface

NI IMAQdx is the interface provided by National Instruments and which supports a wide variety of cameras. It is completely separate from IMAQ, and it supports different communication interfaces: USB, Ethernet and FireWire. It has been tested with Ethernet-connected PhotonFocus HD1-D1312 camera.

The code is located in `pylablib.devices.IMAQdx`, and the main camera class is `pylablib.devices.IMAQdx.IMAQdxCamera`.

## Software requirements

These cameras require `imaqdx.dll`, which is installed with the freely available [Vision Acquisition Software](#). However, the IMAQdx part of the software is proprietary, and requires purchase to use. If the software license is invalid, then any attempt to communicate with cameras will result in `License not activated` error (although simply listing the cameras still works). After installation, the DLL is automatically added to the `System32` folder, where `pylablib` looks for it by default. If the DLL is located elsewhere, the path can be specified using the `library` parameter `devices/dlls/niimaqdx`:

```
import pylablib as pll
pll.par["devices/dlls/niimaqdx"] = "path/to/dlls"
from pylablib.devices import IMAQdx
cam = IMAQdx.IMAQdxCamera()
```

## Connection

The cameras are identified by their name, which usually looks like `"cam0"`. To get the list of all cameras, you can use `NI MAX` (Measurement and Automation Explorer), or `IMAQdx.list_cameras()`:

```
>> from pylablib.devices import IMAQdx
>> IMAQdx.list_cameras()
['cam0', 'cam1']
>> cam1 = IMAQdx.IMAQdxCamera('cam0')
>> cam2 = IMAQdx.IMAQdxCamera('cam1')
>> cam1.close()
>> cam2.close()
```

## Operation

The operation of these cameras is relatively standard. They support all the standard methods for dealing with ROI, starting and stopping acquisition, and operating the frame reading loop. The SDK also provides a universal interface for getting and setting various *camera attributes* using their name. You can use `IMAQdxCamera.get_attribute_value()` and `IMAQdxCamera.set_attribute_value()` for that, as well as `.cav` attribute which gives a dictionary-like access:

```
>> cam = IMAQdx.IMAQdxCamera()
>> cam.get_attribute_value("StatusInformation/AcqInProgress") # check if the camera_
↳ is acquiring
0
>> cam.set_attribute_value("Width", 512) # set the ROI width to 512px
>> cam.cav["Width"] # get the exposure; could also use cam.get_attribute_value("Width
↳ ")
512
```

To see all available attributes, you can call `IMAQdxCamera.get_all_attributes()` to get a dictionary with attribute objects, and `IMAQdxCamera.get_all_attribute_values()` to get the dictionary of attribute values. The attribute objects provide additional information: attribute kind (integer, enum, string, etc.), range (either numerical range, or selection of values for enum attributes), description string, etc.:



```
>> cam = IMAQdx.IMAQdxCamera()
>> attr = cam.get_attribute("Width")
>> attr.description
'Width of the Image provided by the device (in pixels).'
>> attr.writable
True
>> (attr.min, attr.max)
(448, 1312)
```

Since these properties vary a lot between different cameras, it is challenging to write a universal class covering a large range of cameras. Hence, currently the universal class only has the basic camera parameter control such as ROI (without binning) and acquisition status. For many specific cameras you might need to explore the attributes tree (either using the Python class and, e.g., a console, or via NI MAX) and operate them directly in your code.

### Known issues

- It seems like sometimes the camera communication settings might be interfering with its operation. It can show up in an unexpected way, e.g., as an `Attribute` value is out of range error when starting acquisition. If it looks like this might be the case, it is a good idea to open the camera in NI MAX (note that Ethernet cameras are listed under `Network Devices`, not in the general device list) and try to snap a single frame. NI MAX might report some problems with the settings and suggest resolution methods. Once the camera is operational, you can close NI MAX and save the camera settings (request is shown upon closing).
- In general, Ethernet cameras work better with larger packet sizes. However, packets above 1500 bits (so-called jumbo packets) are not supported by all network adapters by default. If this is the case, any attempt to acquire images causes `IMAQdxErrorTestPacketNotReceived` error. One way to deal with that is to set the packet size to 1500, which is done automatically when `small_packet=True` is supplied upon the camera creation. The other is to enable jumbo packets in the adapter properties (in Windows this is done in Device Manager).
- Currently only the basic unpacked monochrome pixel formats are supported: `Mono8`, `Mono10`, `Mono12`, `Mono16`, and `Mono32`. The reason is that even nominally well-defined types (e.g., `Mono12Packed`) have different formats for different cameras. Currently any unsupported format will raise an error on readout by default. It is still possible to read these out as raw frame data in the form of 1D or 2D numpy 'u1' array by enabling raw frame readout using `IMAQdxCamera.enable_raw_readout()` method:

```
>> cam = IMAQdx.IMAQdxCamera()
>> cam.get_detector_size() # 1280px x 1024px frame
(1280, 1024)
>> cam.set_attribute_value("PixelFormat", "BGRA 8 Packed") # unsupported format
>> cam.snap().shape
...
IMAQdxError: pixel format BGRA 8 Packed is not supported
>> cam.enable_raw_readout("frame") # frame data is returned as a flat array
>> cam.snap().shape # 1280 * 1024 * 4 = 5242880 bytes
(5242880,)
```

---

**Note:** General camera communication concepts are described on the corresponding [page](#).

---

### Photon Focus pfcam interface

Photon Focus CameraLink cameras transfer their data to the PC using frame grabbers (e.g., via *NI IMAQ* or *Silicon Software* interfaces). Hence, the camera control is done through the serial port built into the CameraLink interface.



However, the cameras use a closed binary protocol, so all the control is done through the pfcam library provided by Photon Focus. It relies on the libraries exposed by the frame grabber manufacturers (e.g., the standard `cl*serial.dll`) to communicate with the camera directly, meaning that the pfcam user simply calls its method, and all the communication happens behind the scenes.

In principle, pfcam can work with any frame grabber. Because of that, there are two different kinds of classes for this camera. To start with, there is `.PhotonFocus.IPhotonFocusCamera`, which provides interface for addressing camera properties, but can not handle actual frame acquisition. Using this class directly leads to errors in any frame data related methods (e.g., `wait_for_frame`, or `read_multiple_images`), and it is mostly intended to serve as a base class to be combined with the actual frame grabber. Two such combined classes are already provided: `.PhotonFocus.PhotonFocusIMAQCamera` for National Instruments frame grabbers using the *NI IMAQ* interface, `.PhotonFocus.PhotonFocusSiSoCamera` for *Silicon Software* frame grabbers, and `.PhotonFocus.PhotonFocusBitFlowCamera` for *BitFlow* frame grabbers. All classes are complete and ready to use. In addition to combining camera and frame grabber control, they also implement basic consistency support, such as automatic adjustment of frame grabber ROI and data transfer format.

## Software requirements

These cameras require `pfcam.dll`, which is installed with freely available (upon registration) **PFInstaller**. In addition, this DLL requires `comdll.dll` and the DLLs referring to a particular camera, e.g., `mv_d1024e_160.dll`. After installation, the path to the DLLs (all located by default in `Photonfocus/PFRemote/bin` folder in Program Files) is automatically added to system PATH variable, which is one of the places where pylablib looks for it by default. If the DLLs are located elsewhere, the path can be specified using the library parameter `devices/dlls/pfcam`:

```
import pylablib as pll
pll.par["devices/dlls/pfcam"] = "path/to/dlls"
from pylablib.devices import PhotonFocus
cam = PhotonFocus.PhotonFocusIMAQCamera()
```

## Connection

The camera class requires two pieces of information. First is the frame grabber interface connection, e.g., NI IMAQ interface name (e.g., `"img0"`) identified as described in the *NI IMAQ* documentation, or Silicon Software board and applet described in *Silicon Software* documentation. The second piece of information is the pfcam port, which is either a number starting from zero indexing the port in the ports list, or a tuple (manufacturer, port), e.g., `("National Instruments", "port0")`. To list all of the connected pfcam-compatible cameras, you can use the PFRemote software (the interface number is given in parentheses after every connection option in the list) or run `PhotonFocus.list_cameras()`:

```
>> from pylablib.devices import PhotonFocus, IMAQ
>> IMAQ.list_cameras() # get all IMAQ frame grabber devices
['img0.iid']
>> PhotonFocus.list_cameras() # by default, get only the ports which support pfcam_
↪interface
[(1, TCameraInfo(manufacturer='National Instruments', port='port0', version=5,
↪type=0))]
>> cam = PhotonFocus.PhotonFocus.PhotonFocusIMAQCamera(imaq_name="img0.iid", pfcam_
↪port=("National Instruments", "port0"))
>> cam.close()
>> cam = PhotonFocus.PhotonFocus.PhotonFocusIMAQCamera(imaq_name="img0.iid", pfcam_
↪port=1) # same effect as above
>> cam.close()
```

## Operation

The operation of these cameras is relatively standard. They support all the standard methods for dealing with ROI and exposure, starting and stopping acquisition, and operating the frame reading loop. However, there's a couple of differences from the standard libraries worth highlighting:

- The SDK also provides a universal interface for getting and setting various *camera attributes* (called “properties” in the documentation) using their name. You can use *IPhotonFocusCamera.get\_attribute\_value()* and *IPhotonFocusCamera.set\_attribute\_value()* for that, as well as *.cav* attribute which gives a dictionary-like access:

```
>> cam = PhotonFocus.PhotonFocusIMAQCamera()
>> cam.get_attribute_value("Window/W") # get the ROI width
256
>> cam.set_attribute_value("ExposureTime", 0.1) # set the exposure to 100ms
>> cam.cav["ExposureTime"] # get the exposure; could also use cam.get_attribute_
↪value("ExposureTime")
0.1
```

Some values (e.g., *Window.Max* or *Reset*) serve as commands; these can be invoked using *PhotonFocusIMAQCamera.call\_command()* method. To see all available attributes, you can call *IPhotonFocusCamera.get\_all\_attributes()* to get a dictionary with attribute objects, and *IPhotonFocusCamera.get\_all\_attribute\_values()* to get the dictionary of attribute values. The attribute objects provide additional information: attribute range, step, and units:

```
>> cam = PhotonFocus.PhotonFocusIMAQCamera()
>> attr = cam.get_attribute("Window/W")
>> attr.writable
True
>> (attr.min, attr.max)
(16, 1024)
```

- *PhotonFocus.PhotonFocusIMAQCamera* supports all of *IMAQ. IMAQCamera* features, such as trigger control and fast buffer acquisition. Some methods have been modified to make them more convenient: e.g., *PhotonFocusIMAQCamera.set\_roi()* method sets the camera ROI and automatically adjusts the frame grabber ROI to match.
- Same is true for *PhotonFocus.PhotonFocusSiSoCamera*, which, e.g., provides access to all of the frame grabber variables.
- The camera supports a status line, which replaces the bottom one or two rows of the frame with encoded frame-related data such as frame number and timestamp. You can use *PhotonFocus.get\_status\_lines()* function to identify and extract the data in the status lines from the supplied frames. In addition, you can use *PhotonFocus.remove\_status\_line()* to remove the status lines in several possible ways: zeroing out, masking with the previous frame, cutting off entirely, etc.
- If several PhotonFocus cameras are connected, you need to correctly associate different PFCam ports with the corresponding frame grabbers. To do that, you can use the function *PhotonFocus.check\_grabber\_association()*.

---

**Note:** General camera communication concepts are described on the corresponding [page](#)

---

## PCO SC2 cameras interface

SC2 is the interface used with PCO cameras. It has been tested with pco.edge cameras with CLHS and regular CameraLink interfaces. A detailed description of the interface is given in the [manual](#).

The code is located in `pylablib.devices.PCO`, and the main camera class is `pylablib.devices.PCO.PCOS2Camera`.

## Software requirements

These cameras require SC2\_Cam.dll, which is installed with the freely available [pco.camware](#) and [pco.sdk](#) tools. By default, the library searches for DLLs in Digital Camera Toolbox/Camware4 and PCO Digital Camera Toolbox/pco.sdk/bin folder in Program Files folder (or Program files (x86), if 32-bit version of Python is running), as well as in the folder containing the script. If the DLLs are located elsewhere, the path can be specified using the library parameter `devices/dlls/pco_sc2`:

```
import pylablib as pll
pll.par["devices/dlls/pco_sc2"] = "path/to/dlls"
from pylablib.devices import PCO
cam = PCO.PCOS2Camera()
```

## Connection

The cameras are identified by their index, starting from zero, and, possibly, by their interface. To get the total number of connected cameras, you can run `PCO.get_cameras_number`:

```
>> from pylablib.devices import PCO
>> PCO.get_cameras_number()
2
>> cam1 = PCO.PCOS2Camera(idx=0)
>> cam2 = PCO.PCOS2Camera(idx=1)
>> cam1.close()
>> cam2.close()
```

## Operation

The operation of these cameras is relatively standard. They support all the standard methods for dealing with ROI and exposure, starting and stopping acquisition, and operating the frame reading loop. The class also provides read-access to all of the relevant camera data using `PCOS2Camera.get_full_camera_data()`. This method returns data in the internal manufacturer format; to interpret it, you should consult the [manual](#).

## Known issues

- Some cameras support only ROIs which are symmetric with respect to vertical flip. In other words, if the camera detector has vertical size of 2160px, the vertical ROI should always have the form `(x0, 2160-x0)`. It is still possible to set non-symmetric ROI, but it is achieved by the software clipping, while the camera still reads out the smallest symmetric ROI contained the selected one. As a result, the readout time for the same ROI size strongly depends on the ROI position. For example, while vertical ROI of `(0, 8)` has only 8 pixel rows, it is not symmetric, and requires reading the whole frame; hence, it will be as slow as the full-frame acquisition. On the other hand, ROI of `(1076, 1084)` is symmetric, so the camera does read out only 8 rows. This results in

vastly faster readout time. You can use `PCOSC2Camera.requires_symmetric_roi()` to check if the symmetric ROI is required.

---

**Note:** General camera communication concepts are described on the corresponding [page](#)

---

## Princeton Instruments Picam cameras

Picam is the interface provided by Teledyne Princeton Instruments and which supports a set of their cameras. It has been tested with PIXIS 400 camera.

The code is located in `pylablib.devices.PrincetonInstruments`, and the main camera class is `pylablib.devices.PrincetonInstruments.PicamCamera`.

## Software requirements

These cameras require `picam.dll`, which is installed with the freely available [PICam software](#). By default, the library searches for DLLs in Princeton Instruments/PICam/Runtime folder in Program Files folder (or Program files (x86), if 32-bit version of Python is running), as well as in the folder containing the script. If the DLL is located elsewhere, the path can be specified using the library parameter `devices/dlls/picam`:

```
import pylablib as pll
pll.par["devices/dlls/picam"] = "path/to/dlls"
from pylablib.devices import PrincetonInstruments
cam = PrincetonInstruments.PicamCamera()
```

## Connection

The cameras are identified by their serial number, which can look like "2800000001". To get the list of all cameras, you can use `.PrincetonInstruments.list_cameras`:

```
>> from pylablib.devices import PrincetonInstruments
>> PrincetonInstruments.list_cameras()
[TCameraInfo(name='E2V 1340 x 400 (CCD 36) (B) (R)', serial_number='2800000001', model=
↳ 'PIXIS: 400BR', interface='USB 2.0'),
  TCameraInfo(name='E2V 1340 x 400 (CCD 36) (B) (R)', serial_number='2800000002', model=
↳ 'PIXIS: 400BR', interface='USB 2.0')]
>> cam1 = PrincetonInstruments.PicamCamera('2800000001')
>> cam2 = PrincetonInstruments.PicamCamera('2800000002')
>> cam1.close()
>> cam2.close()
```

If no serial number is supplied, the first available camera is connected.

## Operation

The operation of these cameras is relatively standard. They support all the standard methods for dealing with ROI, starting and stopping acquisition, and operating the frame reading loop. The SDK also provides a universal interface for getting and setting various *camera attributes* using their name. You can use `PicamCamera.get_attribute_value()` and `PicamCamera.set_attribute_value()` for that, as well as `.cav` attribute which gives a dictionary-like access:

```
>> cam = PrincetonInstruments.PicamCamera()
>> cam.get_attribute_value("Pixel Format") # get the current pixel format
'Monochrome 16-bit'
>> cam.set_attribute_value("Exposure Time", 10) # set the exposure time to 10 ms
>> cam.cav["Exposure Time"] # get the exposure; could also use cam.get_attribute_
↪value("Exposure Time")
10.0
```

To see all available attributes, you can call `PicamCamera.get_all_attributes()` to get a dictionary with attribute objects, and `PicamCamera.get_all_attribute_values()` to get the dictionary of attribute values. The attribute objects provide additional information: attribute kind (integer, enum, float, etc.), range (either numerical range, or selection of values for enum attributes), default value, etc.:

```
>> cam = PrincetonInstruments.PicamCamera()
>> attr = cam.get_attribute("Exposure Time")
>> attr.default
100.0
>> attr.writable
True
>> (attr.min, attr.max)
(0.0, 8355840.0)
```

Since these properties vary a lot between different cameras, it is challenging to write a universal class covering a large range of cameras. Hence, currently the universal class only has the basic camera parameter control such as ROI (without binning), exposure, and acquisition status. For many specific cameras you might need to explore the attributes tree using the Python class and operate them directly in your code.

## Known issues

- Frame period obtained using `PicamCamera.get_frame_period()` can be an underestimate (i.e., it can overestimate the frame rate).
- While the cameras support multiple ROIs, only single-ROI readout is currently supported.
- Changing readout mode ("Readout Control Mode") to "Kinetics" might invalidate the current ROI, if it was originally too large. Therefore, you would need to call `set_roi` again after setting this mode.
- In principle, the cameras support a variety of different metainfos which can be enabled or disabled separately. However, for simplicity only two modes are supported in the camera class: either no metainfo, or full “standard” metainfo (frame stamp, and start and stop timestamps). Any time the metainfo is enabled, disabled, or queried, it is automatically “truncated” to one of these two modes.

---

**Note:** General camera communication concepts are described on the corresponding [page](#)

---

## Photometrics PVCAM cameras

PVCAM is the interface provided by Teledyne Photometrics and which supports a set of their cameras. It has been tested with Prime 95B camera.

The code is located in `pylablib.devices.Photometrics`, and the main camera class is `pylablib.devices.PvcamCamera`.

## Software requirements

These cameras require `pvcam32.dll` or `pvcam64.dll`, which is installed with the freely available (upon registration) [PVCAM software](#). By default, the library searches for DLL is automatically added to the `System32` folder, where `pylablib` looks for them by default. If the DLL is located elsewhere, the path can be specified using the library parameter `devices/dlls/pvcam`:

```
import pylablib as pll
pll.par["devices/dlls/pvcam"] = "path/to/dlls"
from pylablib.devices import Photometrics
cam = Photometrics.PvcamCamera()
```

## Connection

The cameras are identified by their name, which can look like `"PMUSBCam00"`. To get the list of all cameras, you can use `.Photometrics.list_cameras`:

```
>> from pylablib.devices import Photometrics
>> Photometrics.list_cameras()
['PMUSBCam00', 'PMUSBCam01']
>> cam1 = Photometrics.PvcamCamera('PMUSBCam00')
>> cam2 = Photometrics.PvcamCamera('PMUSBCam01')
>> cam1.close()
>> cam2.close()
```

If no name is supplied, the first camera in the list is connected.

## Operation

The operation of these cameras is relatively standard. They support all the standard methods for dealing with ROI, starting and stopping acquisition, and operating the frame reading loop. The SDK also provides a universal interface for getting and setting various *camera attributes* using their name. You can use `PvcamCamera.get_attribute_value()` and `PvcamCamera.set_attribute_value()` for that, as well as `.cav` attribute which gives a dictionary-like access:

```
>> cam = Photometrics.PvcamCamera()
>> cam.get_attribute_value("EXPOSURE_MODE") # get the current exposure mode
'Internal Trigger'
>> cam.set_attribute_value("METADATA_ENABLED", True) # enable frame metadata
>> cam.cav["METADATA_ENABLED"] # check if metadata is enabled; could also use cam.
    ↪ get_attribute_value("METADATA_ENABLED")
True
```

To see all available attributes, you can call `PvcamCamera.get_all_attributes()` to get a dictionary with attribute objects, and `PvcamCamera.get_all_attribute_values()` to get the dictionary of attribute values. The attribute objects provide additional information: attribute kind (integer, enum, float, etc.), range (either numerical range, or selection of values for enum attributes), default value, etc.:

```
>> cam = Photometrics.PvcamCamera()
>> attr = cam.get_attribute("EXPOSURE_TIME")
>> attr.default
0
>> attr.readable
```

(continues on next page)

(continued from previous page)

```
True
>> (attr.min, attr.max)
(0, 10000)
```

Since these properties vary a lot between different cameras, it is challenging to write a universal class covering a large range of cameras. Hence, currently the universal class only has the basic camera parameter control such as ROI (without binning), exposure, and acquisition status. For many specific cameras you might need to explore the attributes tree using the Python class and operate them directly in your code.

### Fast buffer readout mode

At high frame rates (above ~10kFPS) dealing with each frame individually becomes too slow for Python. Hence, there is an option to read out and process frames in larger ‘chunks’, which are 3D numpy arrays with the first axis enumerating the frame index. This approach leverages the ability to store several frame buffers in the contiguous memory locations (resulting in a single 3D array), and it essentially eliminates the overhead for dealing with multiple frames at high frame rates, as long as the total data rate is manageable (typically below 600Mb/s).

This option can be accessed by calling using `PvcamCamera.set_frame_format()` method to set frames format to "chunks". In this case, instead of a list of individual frames (which is the standard behavior), the method returns list of chunks, which contain several consecutive frames.

### Known issues

- Frame period obtained using `PvcamCamera.get_frame_period()` can be an underestimate (i.e., it can overestimate the frame rate), especially for USB-connected devices.
- While the cameras support multiple ROIs, only single-ROI readout is currently supported.
- Exposure time, exposure mode, and ROI are configured using special methods separately from other camera attributes. Therefore, their corresponding attributes are read-only.
- Not all horizontal and vertical binning combinations are supported. The allowed combinations can be queries using `PvcamCamera.get_supported_binning_modes()`. If the combination is not supported, it is truncated down to the smallest supported one.

---

**Note:** General camera communication concepts are described on the corresponding [page](#)

---

### Silicon Software frame grabbers interface

Silicon Software produces a range of frame grabbers, which can be used to control different cameras with a CameraLink interface. It has been tested with microEnable IV AD4-CL frame grabber together with PhotonFocus MV-D1024E camera.

The code is located in `pylablib.devices.SiliconSoftware`, and the main camera class is `pylablib.devices.SiliconSoftware.SiliconSoftwareCamera`.

### Software requirements

This interfaces requires `fglib5.dll`, which is installed with the freely available (upon registration) [Silicon Software Runtime Environment](#) (the newest version for 64-bit Windows is [5.7.0](#)), which also includes all the necessary drivers.



After installation, the path to the DLL (located by default in `SiliconSoftware/Runtime5.7.0/bin` folder in Program Files) is automatically added to system PATH variable, which is one of the places where pylablib looks for it by default. If the DLL is located elsewhere, the path can be specified using the library parameter `devices/dlls/sisofgrab`:

```
import pylablib as pll
pll.par["devices/dlls/sisofgrab"] = "path/to/dlls"
from pylablib.devices import SiliconSoftware
cam = SiliconSoftware.SiliconSoftwareCamera()
```

## Connection

Figuring out the connection parameters is a multi-stage process. First, one must identify one of several boards. The boards can be identified using `SiliconSoftware.list_boards` function. Second, one must select an applet. These provide different board readout modes and, for Advanced Applets, various post-processing capabilities. These applets can be identified using `SiliconSoftware.list_applets` method, or directly from the Silicon Software RT microDisplay software supplied with the runtime. The choice depends on the color mode (color vs. gray-scale and different bitness), readout mode (area or line), and camera connection (single, double, or quad). Finally, depending on the board and the camera connection, one of several ports must be selected. For example, if the frame grabber has two connectors, but the camera only uses a single interface, then the double camera applet (e.g., `DualAreaGray16`) must be selected, and the port should specify the board connector (0 for A, 1 for B):

```
>> from pylablib.devices import SiliconSoftware
>> SiliconSoftware.list_boards() # first list the connected boards
[BoardInfo(name='mE4AD4-CL', full_name='microEnable IV AD4-CL')]
>> SiliconSoftware.list_applets(0) # list all applets on the first board
[ ...,
  TAppletInfo(name='DualAreaGray16', file='DualAreaGray16.dll'),
  ... ]
>> cam = SiliconSoftware.SiliconSoftwareCamera(0, 'DualAreaGray16') # connect to the
↳first board (port 0 by default)
>> cam.close()
```

Note that currently the code is organized in such a way, that only one port on a single board can be in use at one time.

## Operation

Unlike most camera classes, the frame grabber interface only deals with the frame transfer between the camera and the PC over the CameraLink interface. Therefore, in can not directly control camera parameters such as exposure, frame rate, triggering, ROI, etc. Some similar-looking parameters are still present, but they have a different meaning:

- External trigger controls frame *transfer*, not frame *acquisition*, which is defined by the camera. By default, when the internal frame grabber trigger is used, the frame grabber transfer rate is synchronized to the camera, so every frame gets transferred. However, if the external transfer trigger is used and it is out of sync with the camera, it can result in duplicate or missing frames.
- ROI is defined within the transferred image, whose size itself is determined by the camera ROI. Hence, e.g., if the camera chip is 1024x1024px and its roi is 512x512, then the frame grabber ROI can go only up to 512x512. Any attempts to set it higher result in frame being misshapen or having random data outside of the image area.

The SDK also provides a universal interface for getting and setting various *attributes* using their name. You can use `SiliconSoftwareCamera.get_grabber_attribute_value()` and `SiliconSoftwareCamera.set_grabber_attribute_value()` for that, as well as `.gav` attribute which gives a dictionary-like access:



```
>> cam = SiliconSoftware.SiliconSoftwareCamera()
>> cam.get_grabber_attribute_value("CAMERA_LINK_CAMTYP") # get the camera data format
'FG_CL_SINGLETAP_8_BIT'
>> cam.set_grabber_attribute_value("WIDTH", 512) # set the readout frame width to_
↪ 512px
>> cam.gav["WIDTH"] # get the width; could also use cam.get_grabber_attribute_value(
↪ "WIDTH")
512
```

To see all available attributes, you can call `SiliconSoftwareCamera.get_all_grabber_attributes()` to get a dictionary with attribute objects, and `SiliconSoftwareCamera.get_all_grabber_attribute_values()` to get the dictionary of attribute values. The attribute objects provide additional information: attribute kind (integer, string, etc.), range (either numerical range, or selection of values for enum attributes), description string, etc.:

```
>> cam = SiliconSoftware.SiliconSoftwareCamera()
>> attr = cam.get_grabber_attribute("BITALIGNMENT")
>> attr.values
{1: 'FG_LEFT_ALIGNED', 0: 'FG_RIGHT_ALIGNED'}
```

The parameter can also be inspected in the Silicon Software RT microDisplay software.

## Fast buffer readout mode

At high frame rates (above ~10kFPS) dealing with each frame individually becomes too slow for Python. Hence, there is an option to read out and process frames in larger ‘chunks’, which are 3D numpy arrays with the first axis enumerating the frame index. This approach leverages the ability to store several frame buffers in the contiguous memory locations (resulting in a single 3D array), and it essentially eliminates the overhead for dealing with multiple frames at high frame rates, as long as the total data rate is manageable (typically below 600Mb/s).

This option can be accessed by calling using `SiliconSoftwareCamera.set_frame_format()` method to set frames format to "chunks" (former way of supplying `fastbuff=True` in `SiliconSoftwareCamera.read_multiple_images()` is now deprecated). In this case, instead of a list of individual frames (which is the standard behavior), the method returns list of chunks about 1Mb in size, which contain several consecutive frames.

## Communication with the camera

The frame grabber needs some basic information about the camera: sensor size, bit depth, data transfer format, timeouts, aux lines mapping. This information can be specified using the grabber attributes. The most important transfer parameters are the number of taps and the bitness of the transferred data, which can be set up using `SiliconSoftwareCamera.setup_camlink_pixel_format()`. The values for this parameters can usually be obtained from the camera manuals.

## Known issues

- The maximal frame rate is limited for some boards (at least for the tested microEnable IV AD4-CL board) by about 20kFPS. It seems to be relatively independent of the frame size, i.e., it is not the data transfer rate issue. One possible way to get around it is to use line readout applet, e.g., `DualLineGray16`, and set the frame height to be the integer multiple of the camera frame. This will combine several camera frames into a single frame-grabber frame, effectively lowering the frame rate at avoiding the issue. However, this sometimes leads to incorrect frame splitting: the top line of the “combined” frame does not coincide with the top line of the

original camera frame, so all frames are shifted cyclically by some number of rows. Hence, it might require some post-processing with frames merging and re-splitting.

- As mentioned above, ROI is defined within a frame transferred by the camera. Therefore, if it includes pixels with positions outside of the transferred frame, the acquisition will be faulty. For example, suppose the camera sensor is 1024x1024px, and the *camera* ROI is selected to be central 512x512 region. As far as the frame grabber is concerned, now the camera sensor size is 512x512px. Hence, if you try to set the same *frame grabber* ROI (i.e., 512x512 starting at 256,256), it will expect 768x768px frame. Since the frame is, actually, 512x512px, the returned frame will partially contain random data. The correct solution is to set frame grabber ROI from 0 to 512px on both axes. In general, it is a good idea to always follow this pattern: control ROI only on camera, and always set frame grabber ROI to cover the whole transfer frame.

---

**Note:** General camera communication concepts are described on the corresponding [page](#)

---

## Thorlabs Scientific Cameras interface

This is the interface used in Thorlabs scientific sCMOS cameras such as Kiralux or Zelux. It has been tested with Thorlabs Kiralux camera.

The code is located in `pylablib.devices.Thorlabs`, and the main camera class is `pylablib.devices.Thorlabs.ThorlabsTLCamera`.

## Software requirements

These cameras require `thorlabs_tsi_camera_sdk.dll`, as well as several additional DLLs: `thorlabs_unified_sdk_kernel.dll`, `thorlabs_unified_sdk_main.dll`, `thorlabs_tsi_usb_driver.dll`, `thorlabs_tsi_usb_hotplug_monitor.dll`, `thorlabs_tsi_cs_camera_device.dll`, `tsi_sdk.dll`, `tsi_usb.dll`. All of them is automatically installed with the freely available [ThorCam](#) tools. By default, the library searches for DLLs in Thorlabs/Scientific Imaging/ThorCam folder in Program Files folder (or Program files (x86), if 32-bit version of Python is running), as well as in the folder containing the script. If the DLLs are located elsewhere, the path can be specified using the library parameter `devices/dlls/thorlabs_tlcam`:

```
import pylablib as pll
pll.par["devices/dlls/thorlabs_tlcam"] = "path/to/dlls"
from pylablib.devices import Thorlabs
cam = Thorlabs.ThorlabsTLCamera()
```

## Connection

The cameras are identified by their serial number. To list all of the connected cameras, you can run `Thorlabs.list_cameras_tlcam`:

```
>> from pylablib.devices import Thorlabs
>> Thorlabs.list_cameras_tlcam()
['12001', '12002']
>> cam1 = Thorlabs.ThorlabsTLCamera(serial="12001")
>> cam2 = Thorlabs.ThorlabsTLCamera(serial="12002")
>> cam1.close()
>> cam2.close()
```

If no serial is provided, the software connects to the first available camera.

## Operation

The operation of these cameras is relatively standard. They support all the standard methods for dealing with ROI and exposure, starting and stopping acquisition, and operating the frame reading loop.

**Warning:** The library appears to be not entirely stable: every time acquisition start is issued, there is small (0.1-1%) chance that it will not actually start, which results in timeout errors. Furthermore, there are occasional crashes on the SDK unloading (i.e., camera closing), especially when acquisition has been started and stopped multiple times. It is unclear, what is the cause of this behavior, but it seems to originate from the manufacturer's DLL (bare-bones example and the native Python library reproduce this behavior). Hence, it might be different with different DLL versions.

**Note:** The DLL prints some debug information in the console when camera list is requested and when the camera is opened. At the moment, it is unclear how to get rid of it.

**Note:** General camera communication concepts are described on the corresponding [page](#)

## Uc480/uEye camera interface

This is the interface used in multiple cameras, including many simple Thorlabs and IDS cameras. It has been tested with IDS SC2592R12M and Thorlabs DCC1545M.

Essentially identical interface is available under two different implementations: either as Thorlabs uc480 or as IDS uEye. Both of these seem to cover exactly the same cameras, both are freely available from the manufacturers, and both implement exactly the same functionality. However, these interfaces are not interchangeable, and each camera will only interact with one of them depending on which driver it happens to use (usually based on which of the software packages was installed last). Hence, if you have both [ThorCam](#) and [IDS Software Suite](#) installed, you would need to check both interfaces. Normally, the interface should correspond to the software which can connect to the camera (either ThorCam or uEye Cockpit).

The code is located in `pylablib.devices.uc480`, and the main camera class is `pylablib.devices.uc480.UC480Camera`. Note that while the names only refer to uc480, the same functions and classes equally cover IDS uEye interface if the appropriate `backend` argument is provided.

## Software requirements

Depending on the interface, these cameras require either `uc480.dll`, or `ueye_api.dll`. These are automatically installed with, correspondingly, the freely available [ThorCam](#) software or with [IDS Software Suite](#) (upon registration; note that you need specifically IDS Software Suite, and not IDS peak). By default, the library searches for DLLs in the corresponding Program Files folder (Thorlabs/Scientific Imaging/ThorCam or IDS/uEye), in the locations placed in PATH during the installation, as well as in the folder containing the script. If the DLLs are located elsewhere, the path can be specified using the library parameter `devices/dlls/uc480` or `devices/dlls/ueye`:

```
import pylablib as pll
pll.par["devices/dlls/uc480"] = "path/to/uc480/dlls"
from pylablib.devices import uc480
cam = uc480.UC480Camera()
```

(continues on next page)

(continued from previous page)

```
pll.par["devices/dlls/ueye"] = "path/to/ueye/dlls"
cam = uc480.UC480Camera(backend="ueye")
```

## Connection

The cameras are identified by their camera ID or device ID (both starting from 1). Device ID corresponds to the connection order of the cameras: it is guaranteed to be unique, but will change if the camera is disconnected and reconnected again. On the other hand, camera ID is tied to the camera, but it is set to 1 by default for all cameras, and needs to be manually assigned using `UC480Camera.set_camera_id()`. Alternatively, one can use other characteristics (model or serial number) as a unique identifier. To list all of the connected cameras together with their basic information, you can run `uc480.list_cameras()`:

```
>> from pylablib.devices import uc480
>> uc480.list_cameras()
[TCameraInfo(cam_id=4, dev_id=1, sens_id=11, model='SC2592R12M', serial_number=
↳ '1234567890', in_use=False, status=0)]
>> cam = uc480.UC480Camera(cam_id=4) # connect to the camera using cam_id
>> img = cam.snap()
>> cam.close()
>> cam = uc480.UC480Camera(dev_id=1) # connecting to the same camera using dev_id
>> cam.close()
>> cam = uc480.UC480Camera() # connecting to the first available camera
>> cam.close()
```

If `cam_id = 0` is provided (default), the software connects to the first available camera.

By default, the code above uses Thorlabs uc480 interface. If you want to use ueye interface, you need to specify `backend="ueye"` argument to the corresponding functions and to the camera class upon creation. With that, the example above becomes:

```
>> from pylablib.devices import uc480
>> uc480.list_cameras(backend="ueye") # list all cameras for uEye backend
[TCameraInfo(cam_id=4, dev_id=1, sens_id=11, model='SC2592R12M', serial_number=
↳ '1234567890', in_use=False, status=0)]
>> cam = uc480.UC480Camera(cam_id=4, backend="ueye") # connect to the camera using_
↳ cam_id and ueye backend
>> img = cam.snap()
>> cam.close()
```

## Operation

The operation of these cameras is relatively standard. They support all the standard methods for dealing with ROI and exposure, starting and stopping acquisition, and operating the frame reading loop. However, there's a couple of differences from the standard libraries worth highlighting:

- Some cameras support both binning (adding several pixels together) and subsampling (skipping some pixels). However, only one can be enabled at a time. They can be set independently using, correspondingly, `UC480Camera.get_binning()/UC480Camera.set_binning()` and `UC480Camera.get_subsampling()/UC480Camera.set_subsampling()`. They can also be set as binning factors in `UC480Camera.get_roi()/UC480Camera.set_roi()`. Whether binning or subsampling is set there can be determined by the `roi_binning_mode` parameter supplied on creation.
- Uc480 API supports many different pixel modes, including packed ones. However, pylablib currently supports only monochrome unpacked modes.

- Occasionally (especially at high frame rates) frames get skipped during transfer, before they are placed into the frame buffer by the camera driver. This can happen in two different ways. First, the frame is simply dropped without any indication. This typically can not be detected without using the framestamp contained in the frame info, as the frames flow appear to be uninterrupted. In the second way, the acquisition appears to get “restarted” (the internal number of acquired frames is dropped to zero), which is detected by the library. In this case there are several different ways the software can react, which are controlled using `UC480Camera.set_frameskip_behavior()`.

The default way to address this “restart” event (“ignore”) is to ignore it and only adjust the internal acquired frame counter; this manifests as quietly dropped frames, exactly the same as the first kind of event. In the other method (“skip”), some number of frames are marked as skipped, so that the difference between the number of acquired frames and the internal framestamp is kept constant. This makes the gap explicit in the camera frame counters. Finally (“error”), the software can raise `uc480FrameTransferError` when such event is detected, which can be used to, e.g., restart the acquisition.

One needs to keep in mind, that while the last two methods make “restarts” more explicit, they do not address the first kind of events (quiet drops). The most direct way to deal with them is to use frame information by setting `return_info=True` in frame reading methods like `read_multiple_images`. This information contains the internal camera framestamp, which lets one detect any skipped frames.

### 1.2.3 Stages

Basic concepts are described at the [stages communication page](#).

Currently supported stages:

- Attocube ANC300* and *Attocube ANC350*: most common Attocube positioner controllers. Tested with Ethernet and USB connection for ANC300, and USB connection for ANC350.
- Thorlabs APT/Kinesis*: basic Thorlabs motorized stages and optomechanics devices. Tested with KDC101, K10CR1, and BSC201 motor controllers, KIM101 piezo motor controller, as well as MFF101 and FW102 (described at a [different page](#))
- Thorlabs Elliptec*: resonant piezoelectric Thorlabs stages. Tested with ELL18 and ELL14 rotational mounts.
- Newport Picomotor*: precision piezo-actuated screws based on slip-stick principle. Tested with Newport 8742 Picomotor driver using Ethernet or USB connection.
- Arcus Performax*: fairly common single- and multi-axis motor controllers sold under different brands: Arcus, Nippon Pulse America, or Newmark Systems. Tested with PMX-4EX device with USB connection.
- Trinamic*: universal motor controllers and drivers. Tested with a single-axis TMC2130 controller with USB connection.
- SmarAct*: high-performance piezo sliders. So far only simple open-loop *SCU controllers* are supported. Tested with a standard HCU controller unit.
- Physik Instrumente*: piezo controllers. So far only PI E-516 is supported and tested.

---

**Note:** General device communication concepts are described on the corresponding [page](#).

---

### Stages control basics

#### Basic example

Almost all stages implement the same basic functionality for moving, stopping, homing, and querying the status:

```
stage = Thorlabs.KinesisMotor("27000001") # connect to the stage
stage.home() # home the stage
stage.wait_for_home() # wait until homing is done
stage.move_by(1000) # move by 1000 steps
stage.wait_move() # wait until moving is done
stage.jog("+") # initiate jog (continuous move) in the positive direction
time.sleep(1.) # wait for 1 second
stage.stop() # stop the motion
stage.close()
```

Some stages will miss some of this functions (e.g., no homing), but if it's present, it works roughly in the same manner.

Some concepts are explained below in more detail.

## Basic concepts

### Counters, encoders, homing, and limit switches

Stages have two basic strategies for keeping track of the position. The first one is counting the steps. The problem with it is that once the device is powered up, its position is unknown. Hence, it requires some kind of homing procedure, which usually involves moving to a predefined position and zeroing out the step counter there. This position is defined by the hardware, usually in the form of a limit switch: a physical switch located at the end of the stage travel range, which changes the state when the stage reaches its position. It also usually automatically turns off the motion when tripped, to prevent the motor from overheating or the stage from breaking.

When stepper motors are used, the size of each step (or microstep, if used) is a reasonably well-defined fraction of a turn, so counting them gives fairly reproducible results. On the other hand, piezo slip-stick sliders (such as Attocube, SmarAct, or Picomotor) have inherently unreliable steps size which depends on, e.g., load, direction, position, temperature, or other environmental factors. In this case steps counting, while possible, usually leads to long-term drifts.

If the reliable counting is impossible, like in the case of sliders or regular DC (as opposed to stepper) motors, the manufacturer might add a hardware position readout. It can be digital (encoder) or analog (e.g., resistive, capacitive, or optical readout). The first kind is generally simpler, cheaper and more reliable, but the second one can provide much higher resolution, and can work in more extreme environments (high vacuum, cryogenics). In both cases, the controllers would typically have some kind of feedback loop to smoothly control the motion speed and direction to approach a given position.

### Steps and real coordinates

Almost all stages allow control or readout of position in motor steps, encoder steps, or some other internal units. It is usually not straightforward, or sometimes even impossible, to convert those to real units. In cases where it is possible, it is defined by the motor gearbox and the screw pitch (for linear stages); in most cases, this ratio is provided in the motor or translation stage manual (which can be different from the motor controller manual, and the two might even be completely independent). Sometimes, one even has to do explicit calculations, e.g., getting the number of microsteps per revolution from the controller and motor manufacturer, and the displacement per step from the stage manufacturer.

### Speed control

In many cases, the motor speed is ramped up and down linearly rather than abruptly; hence, both the “cruising” speed and the ramping acceleration can, in principle, be configured. Usually they are defined in, respectively, steps/s and steps/s<sup>2</sup>, although sometimes internal units have to be used.

## Application notes and examples

Here we talk more practically about using pylablib to perform common tasks.

### Motion

The most standard motion methods are `move_to`, which moves to a specified position, `move_by`, which moves by a specified distance or number of steps, and `jog`, which moves continuously in a given direction until stopped or run into a limit switch. If both `move_to` and `move_by` are present, they usually perform the same operation under the hood: `stage.move_by(s)` and `stage.move_to(stage.get_position()+s)` yield the same result.

In almost all cases these commands are asynchronous, in the sense that they simply initialize the motion and continue immediately:

```
>> stage.move_by(1000)
>> stage.is_moving()  # the stage is moving, but the execution continues
True
>> time.sleep(1.)
>> stage.is_moving()  # after 1s the motion is done
False
```

To stop immediately (which is usually only used with `jog` commands) you can use the `stop` method. In some cases, there are two different stop kinds: “soft” with a ramp-down, or “hard” which immediately ceases motion.

### Status and synchronization

Since the motion commands are asynchronous, the devices provide two methods to synchronize it with the script execution. The first one, `is_moving`, checks if the stage is currently in motion. The second one, `wait_move`, pauses the execution until the stage motion is finished.

In addition, many stages provide methods to obtain additional information, e.g., `get_status` (which, usually, returns state of motion, limit switches, possible errors, etc.), or `get_current_speed`.

### Position readout

If a stage has position readout (either hardware sensor, or step counting), it is implemented with the `get_position` method. In most cases, it will be accompanied with the `set_position_reference` method, which lets one change the currently stored position, effectively adding an offset to all further position readings:

```
>> stage.get_position()
10000
>> stage.set_position_reference(20000)  # change current reference to
>> stage.get_position()  # note that it reacts immediately, unlike move_to; no_
↪physical motion happened
20000
>> stage.move_to(21000)  # move by 1000 steps; equivalent to .move_by(1000), or .move_
↪to(11000) before the reference change
```

Note that it only changes the internal counter state, and does not cause any stage motion (which is performed by `move_to`).



## Axis selection

Many controllers support simultaneous control of several different motors. In this case, all of their methods take an additional `axis` (in most cases) or `channel` argument, which specify the exact motor. In cases where usually only one motor is controlled (e.g., TCM1110 or Thorlabs KDC101), this parameter is set to the default value, and is closer to the end of the parameter list. If having multiple controlled stages is the default (e.g., Attocube ANC350 or Arcus Performax), this parameter is usually the first one, and it has to be specified. In this cases, the methods frequently allow to set this parameter to `"all"`, which means that the action is performed for all axes, or the results is returned for all axes (usually in a form of a list or a dictionary).

The channels are usually specified by their index starting from 0 or 1, although some stages adopt a different labeling (e.g., Arcus Performax labels them as X, Y, Z, and U). The exact specification is given in the specific class description.

## Homing

As mentioned above, often stages require homing to get absolute position readings. It needs to be done every time the stage is power-cycled, but the homing parameters usually persist between different re-connections.

If homing is implemented, it is done using the `home` method. In addition, there can also be an `is_homed` method, which checks if the homing has already been performed. If the method is present, then by default `home` will not execute if `is_homed` returns `True`, unless forced.

Some stages do not have an explicit homing method, but can be manually homed by, e.g., running the stage to the limit switch and setting the position reference to 0.

---

**Note:** General stage communication concepts are described on the corresponding [page](#)

---

## Attocube positioners

Attocube has two main positioner controllers: ANC300 and ANC350. These cover different but somewhat overlapping positioner classes, and have fairly different programming interfaces.

### Attocube ANC300

This controller is aimed at open-loop (i.e., no position readout) positioners. It is a chassis with a single PC communication module and up to 7 individual piezo control modules: ANM150 (only stepping), ANM200 (only scanning), or ANM250 (stepping and scanning).

The device class is `pylablib.devices.Attocube.ANC300`.

## Software requirements

The controller has several communication modes: USB, RS232, and Ethernet. USB mode requires a driver supplied with the controller (or downloaded from the controller itself using its Ethernet connection and HTTP port), which makes ANC300 appear as a virtual COM port. RS232 requires a USB-to-RS232 adapter, which usually manifests in the same way. Finally, Ethernet connection works like any other networks device. The controller has been tested with USB and Ethernet communication modes (RS232 is identical to USB, so it should operate as well).

Of all of these modes only USB requires specialized drivers, and the other two are usually available purely through the built-in OS capabilities.



## Connection

The device is identified by its communication address. It can be either a serial port (e.g., "COM5"), or an IP address (e.g., "192.168.1.100"); see [connection description](#) for more information. The backend is chosen automatically based on the connection parameter. Additionally, Ethernet connection requires a password; by default, the standard Attocube password "123456" is used, but if you specified a custom password, you need to provide it upon connection:

```
>> from pylablib.devices import Attocube
>> atc1 = Attocube.ANC300("COM5") # USB or RS232 connection
>> atc2 = Attocube.ANC300("192.168.1.1", pwd="root") # Ethernet connection; no need
↳to provide a password, if it is default
>> atc1.close()
>> atc2.close()
```

Note that since Ethernet inherently supports multiple connections, it is possible to control the same devices in multiple scripts at the same time.

## Operation

This controller has several features and differences compared to most other stages and sliders:

- The controller is inherently multi-axis, hence it always take the axis as the first argument. The axes are numbered starting from 1, and are addressed according to the chassis spaces, so some can be skipped or missing. To update the list of connected axes, use `ANC300.update_available_axes()` (called automatically on connection).
- Different control modules provide different functionality. Hence, not all methods would work for all axes: offset voltage commands such as `ANC300.set_offset()` do not work with ANM150 module, while stepping commands such as `ANC300.move_by()` do not work with ANM200 module. To get the module kinds and serial numbers, use `ANC300.get_axis_serial()`.
- The most important stepping parameters are step voltage amplitude and step frequency (number of steps per second). These can be controlled with, correspondingly, `ANC300.get_voltage()/ANC300.set_voltage()` and `ANC300.get_frequency()/ANC300.set_frequency()`.
- Different axes can be enabled and disabled (i.e., connected or grounded) using `ANC300.enable_axis()` and `ANC300.disable_axis()`. Disabling an axis completely shuts off the connection to the positioner, which usually reduces the noise. In addition, there can be different operation modes for only offset, only stepping, or combination of the two.
- It is possible to measure the positioner capacitance using `ANC300.get_capacitance()`, which is useful in identifying breaks or shorts in the wiring or faults in the piezos. By default, this method simply returns the last measured value. To re-measure, call it with `measure=True`. Note that after the measurement is done, the axis is automatically disabled, and needs to be enabled explicitly:

```
>> atc = ANC300("COM5")
>> atc.get_capacitance(1, measure=True) # get the capacitance (in F) on the
↳first axis; the method waits until the measurement is done (about 1s)
200E-9
>> atc.is_enabled(1)
False
```

Note that this is also the only way to know if there is an actual positioner connected to the given control module.

## Attocube ANC350

This controller is aimed at closed-loop (i.e., with position readout) positioners. It can control up to 3 positioners.

The device class is `pylablib.devices.Attocube.ANC350`.

## Software requirements

The controller has USB and Ethernet modes. USB mode requires a driver supplied with the controller. The communication is done via [PyUSB](#), which means that it does not require any additional Attocube DLLs, although you might need to install libusb (see [PyUSB](#) for more details). Ethernet control is supplied as an additional purchaseable option and can be configured using the supplied Daisy control software.

This device has only been tested with a USB connection.

## Connection

When using a USB connection, the device is identified by its index among all the connected ANC350 devices. To get the total number of devices, you can use `Attocube.get_usb_devices_number_ANC350()`:

```
>> from pylablib.devices import Attocube
>> Attocube.get_usb_devices_number_ANC350()
2
>> atc1 = Attocube.ANC350() # use 0 index by default
>> atc2 = Attocube.ANC350(1)
>> atc1.close()
>> atc2.close()
```

Ethernet connection should work in the same manner as any other similar devices, i.e., the address and, possibly, the port should be provided.

## Operation

This controller has several features and differences compared to most other stages and sliders:

- The controller is inherently multi-axis, hence it always take the axis as the first argument. The axes are numbered 0 through 2. You can check if the slide is connected to the given axis using `ANC350.is_connected()`.
- Different axes can be enabled and disabled (i.e., connected or grounded) using `ANC300.enable_axis()` and `ANC300.disable_axis()`. Disabling an axis completely shuts off the connection to the positioner, which usually reduces the noise.
- It is also possible to control the sensor voltage using `ANC350.get_sensor_voltage()/ANC350.set_sensor_voltage()` methods. Reducing this voltage lowers the heating produced by the sensor, which becomes especially important at very low (<1K) temperatures.
- The most important stepping parameters are step voltage amplitude and step frequency (number of steps per second). These can be controlled with, correspondingly, `ANC350.get_voltage()/ANC350.set_voltage()` and `ANC350.get_frequency()/ANC350.set_frequency()`.
- It is possible to measure the positioner capacitance using `ANC350.get_capacitance()`, which is useful in identifying breaks or shorts in the wiring. By default, this method simply returns the last measured value. To re-measure, call it with `measure=True`.

- Fine positioning is performed using the position readout and the feedback loop. Then a `move_to/move_by` command is issued, this feedback loop is activated, and the positioner tries to reach and stay at the current position. You can use `ANC350.is_target_reached()` to check if the target is reached, `ANC350.get_target_position()` to get the target, and `ANC350.get_precision()/ANC350.set_precision()` to control the target precision.
- In addition, there is a method `ANC350.move_by_steps()`, which mimics `ANC300.move_by()` by moving for a given number of steps instead of a given distance. However, due to implementation limitations, this method is synchronous, i.e., it waits until all steps are performed. Nevertheless, `ANC350.jog()` is still asynchronous.

---

**Note:** General stage communication concepts are described on the corresponding [page](#)

---

## Thorlabs APT/Kinesis devices

Thorlabs has a variety of APT/Kinesis devices for various motion-related functionality (mostly motor controllers and piezo drivers), which share the same API. The library uses an older and more low-level APT protocol to communicate with these devices. So far it has been only implemented for motor controllers and some *specialized devices* and tested with KDC101, K10CR1, and BSC201 motor controllers, KIM101 piezo motor controller, and TPA101 quadrature sensor controller.

The main device classes are `pylablib.devices.Thorlabs.BasicKinesisDevice` for a generic Kinesis/APT devices `pylablib.devices.Thorlabs.KinesisMotor` aimed at motor controllers such as K10CR1 or KDC101, and `pylablib.devices.Thorlabs.KinesisPiezoMotor` for piezo drivers such as KIM and TIM.

## Software requirements

The connection is done using Thorlabs APT protocol, so it needs the corresponding APT [drivers](#). Pylablib communicates directly with the FTDI USB-to-RS232 using `pyft232` chip inside the controller, so it bypasses most of the Thorlabs software. This means that it does not need any Thorlabs-supplied DLLs, but it also means that it can not work with the simulated devices, since these are simulated on a level above the direct serial communication.

In some cases `pyft232` library can not find the required `ftd2xx.dll` library, which leads to an error. There are several ways to get around this. First, you can install the FTDI drivers from the [manufacturer's website](#). Setup executable for Windows automatically places the necessary DLL into the `System32` folder, where `pyft232` can discover them. Alternatively, you can copy the DLLs there yourself from the Thorlabs APT installation. Their default location is `Program Files\Thorlabs\APT\Drivers\APT\USB Driver\amd64` for 64-bit version or `Program Files\Thorlabs\APT\Drivers\APT\USB Driver\i386` for 32-bit version. Note that in the first case the file is called `ftd2xx64.dll`, and you will need to rename it to `ftd2xx.dll` when copying to the `System32` folder.

## Connection

The devices are identified by their address, which correspond to their serial numbers. To get the list of all the connected devices, you can use `Thorlabs.list_kinesis_devices`:

```
>> from pylablib.devices import Thorlabs
>> Thorlabs.list_kinesis_devices()
[('27500001', 'Kinesis K-Cube DC Driver')]
>> stage = Thorlabs.KinesisMotor("27500001")
>> stage.close()
```

## Operation

### Standard motors

This controller has several features and differences compared to most other stages and sliders:

- There are two different classes of devices which require slightly different communication approach: generic USB devices and rack-bay devices. These are hard to detect a priori, so by default generic USB device (which covers the majority of equipment) is assumed. If this assumption is incorrect, the communication becomes impossible, and an attempt to connect to the device raises a communication error `ThorlabsBackendError: backend exception: 'read returned less data than expected' ('read returned less data than expected')`. If you experience this error, you should first power-cycle the device, as it often gets stuck in a non-communicable state, and then double-check that the standard Thorlabs software (Kinesis or APT) can detect and control it. If this is the case, you should supply `is_rack_system=True` to the controller:

```
stage = Thorlabs.KinesisMotor("70000001", is_rack_system=True)
```

- There are several different ways to specify the stage calibration, which are controlled by the `scale` parameter supplied upon the connection. By default (`scale = "step"`), it accepts and returns position in motor steps, velocity in steps/s and acceleration in steps/s<sup>2</sup> (scaling coefficients for the latter two are determined from the controller model). If `scale = "stage"`, the class attempts to autodetect the stage and use meters or degrees instead of steps; in addition you can supply the stage name (e.g., "MTS25-Z8") as a scale instead of relying on the autodetection. If there is no calibration for the stage that you have, you can instead supply a single scaling factor, which specifies the number of steps per physical unit (e.g., for "MTS25-Z8" stage and mm units, one would supply `scale = 34304`). The stage scaling can be obtained from the [APT manual](#). Finally, one can supply a 3-tuple of scales for position, velocity and acceleration (all relative to the internal units). The details are given in the APT manual. To ensure that the units have been applied and/or autodetected correctly, you can use `KinesisMotor.get_scale()`, `KinesisMotor.get_scale_units()` and `KinesisMotor.get_stage()` methods.
- By default, the controllers are treated as single-axis. If several axes are supported, they can be specified using `channel` argument in the corresponding methods such as `move_to` or `get_status`. In addition, you can specify the number of channels using `KinesisMotor.set_supported_channels()` method, in which case settings `channel="all"` in the method would act on all the channels.
- The motor power-up parameters for homing, jogging, limit switches, etc., can be different from the parameters showing up in the APT/Kinesis controller. This can lead to problems if, e.g., homing speed is too low, so the motor appears stationary while homing. You should make sure to check those parameters using `KinesisMotor.get_velocity_parameters()`, `KinesisMotor.get_jog_parameters()`, `KinesisMotor.get_homing_parameters()`, `KinesisMotor.get_gen_move_parameters()`, and `KinesisMotor.get_limit_switch_parameters()`.

### Piezo motors

This controller has several features and differences compared to most other stages and sliders:

- The controllers are treated as multi-axis. However, to be compatible with other Kinesis motor, the `channel` argument is not required, and it defaults to the currently selected “default” channel (1 in the beginning). To control different channels, you can either supply `channel` argument explicitly, or specify a different default channel using `KinesisPiezoMotor.set_default_channel()` or `KinesisPiezoMotor.using_channel()`.
- The motor power-up parameters for jogging and drive can be different from the parameters showing up in the APT/Kinesis controller. This can lead to problems if, e.g., speed is too low. You should make sure to check those

parameters using `KinesisPiezoMotor.get_drive_parameters()` and `KinesisPiezoMotor.get_jog_parameters()`.

- Even open-loop controllers support absolute positioning, which is achieved simply by counting steps in both directions. However, unlike stepper motors or encoders, these steps can be different depending on the direction, position, instantaneous load, speed, etc. Hence, the absolute positions quickly become unreliable. It is, therefore, recommended to generally use relative positioning using `KinesisPiezoMotor.move_by()` method.

## Quadrature detector

These are fairly different from the other discussed devices, since they are more related to sensors than to motors. This controller takes signal from a quadrature photodetector and implements a PI control loop to feed back to some control device (e.g., a piezo driver or a galvo mirror). Hence, all of its methods are fairly distinct from the usual motors. Nevertheless, it is described here, since it still belongs to the APT/Kinesis family of devices and shares their detection and connection approach. The device is implemented in the `pylablib.devices.Thorlabs.KinesisQuadDetector` class.

The operation is fairly straightforward: it implements control of PID parameters, output parameters (such as limits), operation mode (open/close loop), allows for reading current state and setting outputs in the open-loop mode.

---

**Note:** General stage communication concepts are described on the corresponding [page](#)

---

## Thorlabs Elliptec devices

Thorlabs has a line of basic resonant piezoelectric motor stages from Elliptec, which include several rotational and linear stages and feature step-motion and position readout. The library has been tested with ELL18 and ELL14 rotational mounts.

The main device class is `pylablib.devices.Thorlabs.ElliptecMotor`.

## Software requirements

The connection is done using a USB connection together with a built-in USB-to-RS232 chip. It is automatically recognized as a serial port, and no additional software is required. In case the device is not recognized as a serial port, you can fix it by installing freely available [Thorlabs Elliptec software](#).

## Connection

Since the devices are identified as COM ports, they use the standard [connection method](#), and all you need to know is their COM-port address (e.g., COM5):

```
>> from pylablib.devices import Thorlabs
>> stage = Thorlabs.ElliptecMotor("COM5")
>> stage.close()
```

## Operation

These devices have several features and differences compared to most other stages and sliders:

- There is a possibility to have several (up to 16) devices connected to the same controller board (i.e., the same serial port address) using bus distributor. However, since they all use the same serial port, they are all controlled from a single `ElliptecMotor` instance. Hence, in order to refer to specific devices, each communication requires an address (integer from 0 to 15), which is specified by `addr` argument available in almost all methods. When this argument is `None` (which is the default value), the so-called default address is used, which can be accessed via `ElliptecMotor.get_default_addr()` and `ElliptecMotor.set_default_addr()` methods. By default, all connected devices are discovered up the connection, and the first available devices is used as default; therefore, if only a single device is connected, `addr` argument does not have to be used.
- Compared to most motor controllers, Elliptec devices have some limitation related to their inability to communicate while the motor is moving. Therefore, there are no methods to query whether the motor is moving, or stop the motion once initiated. To address that and to simplify the library and the user code, all motion-related methods (`ElliptecMotor.move_to()`, `ElliptecMotor.move_by()`, and `ElliptecMotor.home()`) are made synchronous, i.e., the execution is paused until the motion is complete. Note that this is true even when several devices are connected to the same port.
- There are several different ways to specify the stage calibration, which are controlled by the `scale` parameter supplied upon the connection. By default (`scale = "stage"`), the internal device calibration is used, so all of the positions are expressed in device-specific units (deg or mm). If `scale = "step"`, all of the position are specified in internal device steps instead. Finally, if `scale` is a number, it is the proportionality coefficient between the position units and the internal steps, i.e., the position in user-defined units is multiplied by it to specify the position in steps. The scale for individually addressed devices can be set using `ElliptecMotor.get_scale()` and `ElliptecMotor.set_scale()` methods.

---

**Note:** General stage communication concepts are described on the corresponding [page](#)

---

## Newport Picomotor controller

Newport Picomotor is a series of actuators, usually in a screw format, based on the slip-stick piezo actuation mechanism (similar to, e.g., Attocubes). Operating them requires a driver/controller to output specific voltage pulses. The basic modern open-loop controller is Newport 8742, which can drive up to 4 actuators (but only one at a time), supports connection via USB or Ethernet, and can be daisy-chained to communicate with several controllers through one connection. The class has been tested with this controller and a single standard actuator.

The device class is `pylablib.devices.Newport.Picomotor8742`.

## Software requirements

The controller has two communication modes: USB, and Ethernet. USB mode requires a driver supplied with the freely available [PicomotorApp software](#), while Ethernet connection works like any other networks device and does not require any additional software. The controller has been tested both with USB and Ethernet communication modes.

## Connection

When using the USB connection, the device is identified by its index, starting from 0. To get the number of connected devices, you can use `Newport.get_usb_devices_number_picomotor`:

```
>> from pylablib.devices import Newport
>> Newport.get_usb_devices_number_picomotor()
2
>> stage1 = Newport.Picomotor8742()
```

(continues on next page)

(continued from previous page)

```
>> stage2 = Newport.Picomotor8742(1)
>> stage1.close()
>> stage2.close()
```

Ethernet connection requires a host name or an IP address. Both can be set up by first connecting the device via USB or by using the PicomotorApp software (in the Setup -> Ethernet menu). After that, they can be supplied to the class instead of index:

```
>> from pylablib.devices import Newport
>> stage1 = Newport.Picomotor8742("8742-12345") # by default, all host names start_
↳with 8742
>> stage1.close()
```

## Operation

This controller has several features and differences compared to most other stages and sliders:

- The controller is inherently multi-axis, hence it always take the axis as the first argument. The axes are labeled numerically starting from 1 (i.e., 1, 2, 3, and 4). The list of all axes is related to the exact controller, and can be obtained using `Picomotor8742.get_all_axes()`.
- There is an option to auto-detect motors and their kind using `Picomotor8742.autodetect_motors()` method. However, since it involves stepping the motor, it usually makes more sense to detect them once and then store them into the non-volatile (i.e., power-independent) memory using `Picomotor8742.save_parameters()`.
- Even open-loop controllers support absolute positioning, which is achieved simply by counting steps in both directions. However, unlike stepper motors or encoders, these steps can be different depending on the direction, position, instantaneous load, speed, etc. Hence, the absolute positions quickly become unreliable. It is, therefore, recommended to generally use relative positioning using `Picomotor8742.move_by()` method.
- As mentioned above, the controller support daisy-chaining using RS-485 connections. It allows to connect several controllers together while still only using a single PC connection. In this case, it is recommended to supply `multiaddr=True` upon connecting to the device. If, in addition `scan=True` is set (default), then upon connection the controller scans for all other connected devices, resolves their address conflicts, and builds the list of the available addresses (address is a number between 1 and 31). The list can later be read using `Picomotor8742.get_addr_map()`, and the network rescanned using `Picomotor8742.scan_devices()`. To refer to a specific device, its address should be specified using `addr` parameter of a method; by default it is set to `None`, which selects the device connected to the PC.

---

**Note:** General stage communication concepts are described on the corresponding [page](#)

---

## Arcus Performax positioners

Arcus has several motor controllers and drivers, which are mainly different in their number of axes, communication possibilities, and driving function. They are also distributed under different names, e.g., Nippon Pulse America (NPA) or Newmark Systems. However, the models nomenclature is the same: there is 4EX for 4-axis controllers with USB and RS485 connection, 2EX/2ED for 2-axis controllers with USB and RS485 connections, and 4ET for 4-axis controllers with Ethernet connection. The class has been tested with 4EX and (partially) 2ED controllers with USB and RS-485 connectivity mode, but other controllers mentioned above should also work.



The main device classes are `pylablib.devices.Arcus.Performax4EXStage` or 4-axis controllers, `pylablib.devices.Arcus.Performax2EXStage` for 2-axis controllers, and `pylablib.devices.Arcus.PerformaxDMXJSASStage` for simple single-axis controller (DMX-J-SA). In addition to a different number of axes, they have several syntax differences, so one can not substitute for the other.

In addition, there is also a generic Performax stage class `pylablib.devices.Arcus.GenericPerformaxStage`, which implements only the most basic functions: ASCII communication with the device and basic methods such as device name request. It can be used with new or not currently supported Arcus stages to directly control them using the ASCII control language (usually described in the stage manual).

## Software requirements

The controller has several communication modes: USB, RS485, and Ethernet. USB mode requires a driver supplied with the operation software: [Arcus Drivers and Tools](#), [Performax Series Installer](#), and [Performax USB Setup](#) (all obtained at [Arcus website](#)). Installing all three seem to be sufficient. Once the appropriate USB drivers are installed, one can connect the device directly via its USB port and use the manufacturer DLLs `PerformaxCom.dll` and `SiUSBXp.dll` to communicate with the device. They can be obtained on the manufacturer's [website](#) and placed in the folder with the script, or in the System32 Windows folder. If the DLL is located elsewhere, the path can be specified using the library parameter `devices/dlls/arcus_performax`:

```
import pylablib as pll
pll.par["devices/dlls/arcus_performax"] = "path/to/dll"
from pylablib.devices import Arcus
stage = Arcus.Performax4EXStage()
```

**Warning:** There appear to be some issues for USB-controlled devices with Python 3.6 which result in out-of-bounds write, memory corruption, and undefined behavior. Hence, Python 3.7+ is required to work with this device.

RS-485 connection does not require any device-specific drivers or DLLs, but it does need RS-485 controller connected to the PC. Such controllers usually show up as virtual COM ports, and they typically do not need any additional drivers.

## Connection

When using the USB connection, the device is identified by its index, starting from 0. To get the list of all the connected devices, you can use `Arcus.list_usb_performax_devices`:

```
>> from pylablib.devices import Arcus
>> Arcus.list_usb_performax_devices()
[(0, '4ex01', 'Performax USB',
  '\\\\?\\usb#vid_1589&pid_a101#4ex01#{xxxxxxxx-xxxx-xxxx-xxxx-xxxxxxxxxxxx}', '1589',
  ↪ 'a101'),
 (1, '4ex21', 'Performax USB',
  '\\\\?\\usb#vid_1589&pid_a101#4ex21#{xxxxxxxx-xxxx-xxxx-xxxx-xxxxxxxxxxxx}', '1589',
  ↪ 'a101')]
>> stage1 = Arcus.Performax4EXStage()
>> stage2 = Arcus.Performax2EXStage(idx=1)
>> stage1.close()
>> stage2.close()
```

When using the RS-485 connection, you need to specify the serial port corresponding to your RS-485 connection and, possibly, its baud rate:



```
stage = Arcus.Performax4EXStage(conn = "COM5")
stage2 = Arcus.Performax4EXStage(conn = ("COM5", 38400)) # specify a baud rate
```

The baud rate is 9600 by default, which is the standard value for the controllers. However, it can be changed using `Performax4EXStage.set_baudrate()` method, in which case you would need to explicitly specify it during the next connection.

In RS-485 mode `idx` parameter is still used, and it specifies the device number connected to this controller. By default this number is 0, and it can be queried (using USB connection) via `Performax4EXStage.get_device_number()`. It can also be set using `Performax4EXStage.set_device_number()`, although the changes takes effect only after the device is power cycled. Although in principle `idx` can be used to distinguish several Arcus controllers connected to the same bus (i.e., sharing the same RS-485 COM port), currently only single device connection is supported.

To switch between USB and RS-485 control modes, you need to plug or unplug USB connection. It is strongly recommended to power cycle the device after that, since otherwise it might stop responding to RS-485 commands.

## Operation

This controller has several features and differences compared to most other stages and sliders:

- The 4-axis and 2-axis controllers are inherently multi-axis, hence they always take the axis as the first argument. The axes are labeled with letters "x", "y" for a 2-axis version, or "x", "y", "z", "u" for a 4-axis one. The list of all axes is related to the exact controller, and can be obtained using `Performax4EXStage.get_all_axes()`. A single-axis controller does not take an axis argument.
- Different axes can be enabled and disabled using `Performax4EXStage.enable_axis()`. Note that disabled axes still behave the same as the enabled ones; e.g., their position will increment as usual, when `move_to` is called. This can lead to some confusion, as the axis appears mostly operational, but the motor does not move.
- In the default controller configuration the limit errors are enabled. In this case, once a single axes reaches the limit switch during motion, it is put into an error state, which immediately stops this and all other axes. Any further motion command on this axis will raise an error, although it is still possible to restart motion on other axes. The axis motion can only be resumed by calling `Performax4EXStage.clear_limit_error()`. If, however, limit errors are disabled, then only the axis which reached the limit is stopped, and all other axes are unaffected. Furthermore, the motion on the offending axis can be resumed without clearing its error status. In many cases the default limit error behavior is undesirable, so the class turns it off upon connection. It can be subsequently turned on and off using `Performax4EXStage.enable_limit_errors()`, and checked using `Performax4EXStage.limit_errors_enabled()`.
- Since simplified single-axis controller (DMX-J-SA) always has limit errors disabled, its behavior is specified a bit differently. Upon connection you can specify `autoclear` argument (`True` by default), which indicates that before every movement command the limit error should be automatically cleared.
- The controllers also have analog and digital inputs and digital outputs, which can be queried and set with the corresponding commands.
- The controller has an option to connect an encoder for a separate position readout. By default, all of the commands (e.g., for moving, getting position, getting current speed, etc.) still work in the step-counting mode, and the encoder values are only accessed via `Performax4EXStage.get_encoder()/Performax4EXStage.set_encoder_reference()`. In principle, there is a closed-loop mode call `StepNLoop`, but it is not currently supported in the code.
- The built-in motion command has 2 modes: relative and absolute. The code sets the absolute mode on connection and assumes it in all commands. However, if the mode changes for any reason, the move commands will stop working properly.

---

**Note:** General stage communication concepts are described on the corresponding [page](#)

---

## Trinamic TMC1110 controller

TMC1110 is a universal single-axis stepper motor controller from Trinamic. It provides multiple connection options, but so far has only been tested with USB connection.

The main device class is `pylablib.devices.Trinamic.TMC1110`.

## Software requirements

USB connection needs drivers, which are supplied with the freely-available [TMCL-IDE](#), or [TMCL-LITE](#). With those drivers installed, the controllers show up as virtual COM ports. Note that when several devices are connected, they sometimes get assigned conflicting (i.e., overlapping) COM ports. In this case, you might need to manually reassign these in the Device Manager.

## Connection

Since the devices are identified as virtual COM ports, they use the standard [connection method](#), and all you need to know is their COM-port address (e.g., COM5):

```
>> from pylablib.devices import Trinamic
>> stage1 = Trinamic.TMC1110("COM5")
>> stage2 = Trinamic.TMC1110("COM8")
>> stage1.close()
>> stage2.close()
```

## Operation

This controller has several features and differences compared to most other stages and sliders:

- The controller allows one to control the number of microsteps per step using `TMC1110.get_microstep_resolution()` and `TMC1110.set_microstep_resolution()`. Hence, the calibration of the real position to the controller readout position depends on this resolution. Furthermore, changing this resolution does not affect the step counter, meaning that changing it, performing a move, and changing it back will result in a different position. Hence, it is not recommended to change it after homing or referencing the position.
- Similarly, the controller has variable frequency divisors, which control the ratio between internal and real units for the velocity and the acceleration. They are set up together with the maximal velocity and acceleration using `TMC1110.setup_velocity()` and `TMC1110.get_velocity_parameters()`, and the conversion factors can be obtained using `TMC1110.get_acceleration_factor()` and `TMC1110.get_velocity_factor()`.
- The device has an option of controlling maximal output current using `TMC1110.setup_current()` and `TMC1110.get_current_parameters()`. Change them carefully, since the values which are too large can damage the motor. Also take into account, that the currents are defined relative to the maximal output current, which is controlled using the physical jumper on the board.

---

**Note:** General stage communication concepts are described on the corresponding [page](#)

---

## SmarAct positioners

SmarAct has multiple different controller covering different slider kinds. So far only simple controllers (CU/HCU/SCU) are implemented.

### SmarAct CU/HCU/SCU

This is a simple controller, which is mostly aimed at open-loop (i.e., no position readout) positioners. It can control up to 3 axes, and connects to the PC via the USB port.

The device class is `pylablib.devices.SmarAct.SC3D`. Currently only open-loop controllers are supported.

## Software requirements

The controller shows up as a virtual COM port, and it has a standard FTDI chip, so it does not need any special drivers. However, to communicate with the device, it still needs `SC3DControl.dll` library. It is supplied on a CD together with the device, although it might also be possible to request it from SmarAct.

## Connection

The devices are identified by their index starting from 0. To get the list of all the connected devices, you can use `SmarAct.list_scu_devices`:

```
>> from pylablib.devices import SmarAct
>> SmarAct.list_scu_devices()
[TDDeviceInfo(device_id=0, firmware_version='1.3.0.0', dll_version='4.3.0.0')]
>> stage = SmarAct.SC3D(idx=0) # connect to the first device in the list
>> stage.close()
```

Due to the manufacturer’s API organization, it is currently only possible to “reserve” all connected stages of the same type simultaneously in one application. This means that no other application can connect to any of the stages as long as at least one stage is being controlled (though it does not make any difference if only one stage is connected).

In addition, currently there is no check on whether the stage is already controlled in the other part of the code. This is in contrast with the vast majority of the devices, which issue a unique handle making it impossible to create two different device objects even within the same application. Hence, one needs to be careful to not connect to the same device twice, which can lead to confusing behavior.

## Operation

This controller has several features and differences compared to most other stages and sliders:

- The motion is generally executed in “macrosteps”, which is a sequence of several “microsteps” with a given amplitude, frequency, and number. A single macrostep with the defined parameters can be performed with `SC3D.move_macrostep()`, while `SC3D.move_by()` executes a series of these macrosteps with one of the predefined sizes (from 0 to 20). These sizes are configured to roughly correspond to the step sizes selectable by the controller, although the agreement is not exact.

---

**Note:** General stage communication concepts are described on the corresponding [page](#)

---

## Physik Instrumente (PI) controllers

Physik Instrumente produces a variety of piezo, servo, and slider controller. So far, only PI E-516 is supported and tested via a standard serial connection.

The main device class is `pylablib.devices.PhysikInstrumente.PIE516`.

## Software requirements

The devices provide a bare RS232 interface, so any appropriate USB-to-RS232 adapter should work.

## Connection

Since the devices are identified as COM ports, they use the standard [connection method](#), and all you need to know is their COM-port address (e.g., COM5):

```
>> from pylablib.devices import PhysikInstrumente
>> stage1 = PhysikInstrumente.PIE516("COM5")
>> stage2 = PhysikInstrumente.PIE516("COM8")
>> stage1.close()
>> stage2.close()
```

## Operation

This controller has several features and differences compared to most other stages and sliders:

- The controller supports either servo (position feedback) or direct voltage output modes, controlled with `PIE516.enable_servo()` method. In the servo mode it is more similar to a stage controller, and you can use, e.g., `PIE516.move_to()` and `PIE516.stop()` methods. In the direct voltage mode you can use `PIE516.set_voltage()` to set the voltage directly.
- The controller only accepts commands from the PC when it is in the “online” (i.e., remote) mode, in which case external voltage controls are ignored. This mode is enabled automatically upon connection if `auto_online=True` is supplied upon creation (default), and can be connected via `PIE516.enable_online()` method. Note that in this case manual servo switches should be turned off, since otherwise the device is permanently in the servo mode.

### 1.2.4 Basic sensors

Basic concepts are described at the [basic sensors communication page](#).

Currently supported sensors:

- *HighFinesse*: laser wavelength meters. Tested with WS6 and WS7 USB-controlled devices.
- *Ophir*: optical power and energy meters. Tested with Ophir Vega.
- *Lakeshore*: temperature sensors. Tested with Lakeshore 218.
- *Cryocon*: temperature sensors. Tested with CryoCon 14C.

- *Cryomagnetics*: liquid nitrogen or helium level sensor. Tested with LM-500 and LM-510 sensors.
- *Pfeiffer*: pressure gauges. Tested with TPG261 and DPG202 controllers.
- *Leybold*: pressure gauges. Tested with ITR90 gauge.
- *Kurt J. Lesker*: pressure gauges. Tested with KJL300 gauge.
- *Thorlabs quadrature detector controller*. Tested with TPA101.

---

**Note:** General device communication concepts are described on the corresponding [page](#).

---

## Basics of sensors communication

### Basic example

Basic sensors usually only implement a handful of functions related to reading out the measurements (possibly on different channels) and setting up measurements modes:

```
>> gauge = Pfeiffer.TPG260("COM1") # connect to the gauge
>> gauge.enable(1) # enable the first channel (usually it's already enabled)
>> gauge.get_pressure() # read pressure at the default channel (1)
100E3
>> gauge.close()
```

### Application notes and examples

Here we talk more practically about using pylablib to perform commons sensor tasks.

### Readout

The main readout methods almost always start with `get_` prefix, e.g., `get_pressure`, `get_temperature`, or `get_level`. In some cases there would be two different measurement modes: one which just reads the latest measurement result, and one which initializes the measurement, waits until it's done, and returns the result. These two approaches may be implemented differently in different devices, and it is addressed in their description:

```
>> meter = Cryomagnetics.LM500("COM1")
>> meter.get_level(1) # immediately return the latest reading
20.0
>> meter.get_level(1) # return the same reading
20.0
>> meter.measure_level(1) # initialize a new measurement; takes some time
19.8
```

### Non-numerical values

In some cases the readout method would return a non-numerical values. This usually happens when the sensor readings are outside of its range, or if it is in a wrong state (off, warming up, error, etc.) These cases are documented in the querying method description:

```
>> meter = Ophir.VegaPowerMeter("COM1")
>> meter.get_power() # power is higher than the current range
'over'
>> meter.set_range_idx(0) # set the maximal power range
>> meter.get_power() # now the reading is numerical
10E-3
```

## Units

Unless absolutely necessary and obvious, all the readout values are specified in SI units (even, e.g., laser frequency in Hz, or pressure in Pa). In rare cases when the devices allows for selection of readout units (e.g., Pfeiffer TPG260 gauges), it only affects the displayed value, but not the results returned by the corresponding methods:

```
>> gauge = Pfeiffer.TPG260("COM1")
>> gauge.set_units("pa")
>> gauge.get_pressure()
100E3
>> gauge.set_units("mbar")
>> gauge.get_pressure() # pressure still in Pa
100E3
>> gauge.get_pressure(display_units=True) # pressure in display units
1000
```

## Channel selection

Some gauges support simultaneous readout on several channels. In this case, all of their methods take an additional `channel` (in most cases) argument, which specify the read channel.

The channels are usually specified by their index starting from 0 or 1, although some devices adopt more complicated labeling schemes (e.g., Lakeshore 218 temperature sensor can only assign a sensor type to a group of 4 sensors, which is labeled "A" or "B"). The exact specification is given in the specific class description.

## Currently supported sensors

- *HighFinesse*: laser wavelength meters. Tested with WS6 and WS7 USB-controlled devices.
- *Ophir*: optical power and energy meters. Tested with Ophir Vega.
- *Lakeshore*: temperature sensors. Tested with Lakeshore 218.
- *Cryocon*: temperature sensors. Tested with CryoCon 14C.
- *Cryomagnetics*: liquid nitrogen or helium level sensor. Tested with LM-500 and LM-510 sensors.
- *Pfeiffer*: pressure gauges. Tested with TPG261 and DPG202 controllers.
- *Leybold*: pressure gauges. Tested with ITR90 gauge.
- *Kurt J. Lesker*: pressure gauges. Tested with KJL300 gauge.
- *Thorlabs quadrature detector controller*. Tested with TPA101.

---

**Note:** General sensor communication concepts are described on the corresponding [page](#)

---

## HighFinesse wavemeters

HighFinesse produces a variety of fiber-coupled wavelength meters. Currently pylablib only deals with WS series which uses a USB connection. The code has been tested with several WS6 and WS7 wavemeters.

The main device class is `pylablib.devices.HighFinesse.WLM`.

## Software requirements

HighFinesse employs a fairly unique control system.

First, one needs to install the control software, which is uniquely tied to a particular wavemeter and is supplied with it. In theory, software from another wavemeter might still work, but the results are not guaranteed.

Second, this control software runs an application server which processes all requests from third-party software. This means, that the main application needs to be running to perform any device communication from the code. The code has an option of automatically starting it, but on some occasions it might fail, in which case it is necessary to either manually start it, or supply the location of the executable file.

---

**Note:** The control software should keep running the whole time. As soon as it is closed, the device will raise an error.

---

Finally, one needs the DLL to communicate with this software. It is usually named `wlmData.dll`, and it is located in the main controller software folder either in `Com-Test` (for 32-bit applications) or `Projects/64` (for 64-bit applications).

## Connection

The device class makes an attempt to search for the DLL and executable in the standard installation folders, as well as use the DLL in the standard location and its executable auto-detection capabilities. However, depending on the number of installed wavemeters and their installation locations, one needs to provide up to 3 arguments on connection. First, the wavemeter ID, which simply a 1 to 5-digit number (e.g. 1234). It is used to identify the correct instance of the control software, either by searching for the correct folder, or via DLL autostart capabilities. Second, one might need to provide the path to `wlmData.dll` (either including the name, or simply the containing folder). Its location is described in the above section. Finally, you might also need to give the path to the application executable, which is located in the main installation folder and is named `wlm_ws*.exe`, where `*` is the wavemeter generation (e.g., `wlm_ws7.exe` for WS7 wavemeters). Hence, the fully qualified (and, therefore, most robust) instantiation looks like this:

```
>> import os
>> from pylablib.devices import HighFinesse
>> app_folder = r"C:\Program Files\HighFinesse\Wavelength Meter WS7 1234"
>> dll_path = os.path.join(app_folder, "Projects", "64")
>> app_path = os.path.join(app_folder, "wlm_ws7.exe")
>> wm = HighFinesse.WLM(1234, dll_path=dll_path, app_path=app_path)
>> wm.close()
```

A unique property of this device is the ability to control it simultaneously from several applications. Keep this in mind, since it might cause confusion or strange results if the control attempts are not synchronized.

**Warning:** Communication with several simultaneously running wavemeters from a single application has not been tested, and might not work correctly.

## Operation

The operation of the wavemeter is fairly straightforward, but there is a couple of points to keep in mind:

- By default, the main measurement functions (`WLM.get_frequency()` and `WLM.get_wavelength()`) raise an error on over- or under-exposure. If this is undesirable (e.g., the laser has power jumps), one can instead make it return "over" or "under" on these occasions.
- The measurement result is returned immediately, but it is updated only about every 15-30ms (+ exposure time). Hence, fast consecutive calls to `WLM.get_frequency()` and `WLM.get_wavelength()` will return the same value.
- Multi-channel devices have two working modes: single-channel (when only one channel is enabled at a time) and cycling (the wavemeter constantly cycles through several channels for quasi-simultaneous measurements). Some methods only make sense in one of this modes, e.g., `WLM.set_active_channel()` only works in the single-channel mode, while `WLM.enable_switcher_channel()` only in the multi-channel mode. By default, these methods will automatically switch to the corresponding mode.
- Due to a minor control software bug, change in the exposure on some channels might not be reported until the control software is switched to the corresponding channel's exposure control tab (in the upper right corner). By default, the device class performs this switching any time the exposure value is queried, which solves the issue. However, it does take about 10ms. If it is critical, it's possible to turn off this behavior by setting `auto_channel_tab` attribute to `False`.

---

**Note:** General sensor communication concepts are described on the corresponding [page](#)

---

## Ophir power meters

Ophir produces a variety of power and energy meters with different controllers and measurement heads. The class has been tested with Ophir Vega controller with a photodiode head.

The main device classes are `pylablib.devices.Ophir.OphirDevice` for a generic device and `pylablib.devices.Ophir.VegaPowerMeter` for Vega power meter.

## Software requirements

The device provides a bare RS232 interface, so any appropriate USB-to-RS232 adapter should work.

## Connection

Since the devices are identified as COM ports, they use the standard [connection method](#), and all you need to know is their COM-port address (e.g., COM5) and the baudrate, if it is different from the standard one (9600 baud):

```
>> from pylablib.devices import Ophir
>> meter1 = Ophir.VegaPowerMeter("COM5") # default connection assumes 9600 baud
>> meter2 = Ophir.VegaPowerMeter("COM6", 19200) # if the second power meter has a
↪different baudrate
>> meter1.close()
>> meter2.close()
```



## Operation

The operation of the power meter is fairly straightforward, but there is a couple of points to keep in mind:

- On the Vega controller the results can be sent at most 15 times a second. However, they are not necessarily updated at this rate, so several consecutive request might yield the same result.
- The device provides the way to change the communication baud rate. If the rate is changed, the device is automatically disconnected, and the new object needs to be instantiated with the updated baudrate.
- The device might return "over" instead of the power reading on overexposure. To fix that, you can adjust the measurement range using `VegaPowerMeter.set_range_idx()`.

---

**Note:** Basic sensors communication concepts are described on the corresponding [page](#)

---

## Lakeshore temperature sensors

Lakeshore manufactures a range of temperature sensor controllers and resistance bridges, which are also used for temperature sensing. There is some overlap between different products, but they still use fairly distinct interfaces and interaction patterns. The code has been tested with Lakeshore 218 temperature controller.

The main device class is `pylablib.devices.Lakeshore.Lakeshore218`.

## Software requirements

The device provides a bare RS232 interface, so any appropriate USB-to-RS232 adapter should work.

## Connection

Since the devices are identified as COM ports, they use the standard [connection method](#), and all you need to know is their COM-port address (e.g., COM5):

```
>> from pylablib.devices import Lakeshore
>> sensor = Lakeshore.Lakeshore218("COM5")
>> sensor.close()
```

Note that the connection uses the standard which is fairly different from most RS232 controllers: 7 data bits, 1 parity bit, and 1 stop bit (as opposed to 8 data bits and no parity bit for most controllers). Hence, it is possible that not all RS232 controllers can communicate with it. In addition, they might need a null-modem (crossed Rx and Tx lines) RS232 cable.

## Operation

The operation of this temperature sensor is fairly straightforward, but there is a couple of points to keep in mind:

- Like most similar devices, querying temperature using `Lakeshore218.get_temperature()` immediately returns the most recently measured value. Re-measurement is periodically initiated by the devices itself.
- It is possible to specify custom response curves by using `Lakeshore218.set_curve_header()` and `Lakeshore218.set_curve()`. However, you need to be careful, as it overwrites the stored user curves.

---

**Note:** Basic sensors communication concepts are described on the corresponding [page](#)

---

## CryoCon temperature sensors

CryoCon manufactures a range of temperature sensor controllers and resistance bridges, which are also used for temperature sensing. The code has been tested with CryoCon 14C temperature controller.

The main device class is `pylablib.devices.Cryocon.Cryocon1x`.

## Software requirements

The device provides a bare RS232 interface, so any appropriate USB-to-RS232 adapter should work.

## Connection

Since the devices are identified as COM ports, they use the standard [connection method](#), and all you need to know is their COM-port address (e.g., COM5):

```
>> from pylablib.devices import Cryocon
>> sensor = Cryocon.Cryocon1x("COM5")
>> sensor.close()
```

## Operation

The operation of this temperature sensor is fairly straightforward, but there is a couple of points to keep in mind:

- Like most similar devices, querying temperature using `Cryocon1x.get_temperature()` immediately returns the most recently measured value. Re-measurement is periodically initiated by the devices itself.

---

**Note:** Basic sensors communication concepts are described on the corresponding [page](#)

---

## Cryomagnetics level monitor

Cryomagnetics manufactures cryogenic liquid level monitors, which are used for mintoring liquid nitrogen or helium levels inside cryostats. The two level meters supported in the package are LM-500 and LM-510; despide difference in appearance, their functionalities are very similar, so their interfaces are nearly identical.

The main device classes are `pylablib.devices.Cryomagnetics.LM500` and `pylablib.devices.Cryomagnetics.LM500`.

## Software requirements

LM-500 provides a bare RS232 interface, so any appropriate USB-to-RS232 adapter should work. LM-510 has a USB interface with a built-in USB-to-RS232 adapter, which is automatically recognized as a serial port, so no additional software is required.

## Connection

Since the devices are identified as COM ports, they use the standard *connection method*, and all you need to know is their COM-port address (e.g., COM5):

```
>> from pylablib.devices import Cryomagnetics
>> sensor = Cryomagnetics.LM510("COM5")
>> sensor.close()
```

## Operation

The operation of this temperature sensor is fairly straightforward, but there is a couple of points to keep in mind:

- Upon connection the devices are automatically switched into the remote mode, which disables manual controls. If this mode is manually switched off (e.g., using `Local` button in LM-510), the device will no longer obey the remote commands, even though the readout would still work.
- There are no specific commands for stopping a refill or resetting the timeout state after a timed-out refill. However, both can be achieved using `LM500.reset()` method.
- Only LM-510 supports switching the automated refill option on and off using `LM510.set_control_mode()` method.
- Like most similar devices, querying the level using `LM500.get_level()` immediately returns the most recently measured value. Re-measurement is periodically initiated by the devices itself, or can be initiated manually using `LM500.start_measurement()` or `LM500.measure_level()`.

---

**Note:** Basic sensors communication concepts are described on the corresponding *page*

---

## Pfeiffer pressure gauges

Pfeiffer manufactures a range of pressure gauges and controllers with several different standards and communication protocols. The code has been tested with Pfeiffer TPG260 series controller (specifically, TPG261) and Pfeiffer DPG202 controller.

The main device classes are `pylablib.devices.Pfeiffer.TPG260` and `pylablib.devices.Pfeiffer.DPG202`.

## Software requirements

The devices provide a bare RS232 interface, so any appropriate USB-to-RS232 adapter should work.

## Connection

Since the devices are identified as COM ports, they use the standard *connection method*, and all you need to know is their COM-port address (e.g., COM5):

```
>> from pylablib.devices import Pfeiffer
>> gauge = Pfeiffer.TPG260("COM5")
>> gauge.close()
```

## Operation

### TPG260 series

The operation of this gauge is fairly straightforward, but there is a couple of points to keep in mind:

- On measurement error `TPG260.get_pressure()` returns `None`. To get the underlying issue, you can use `TPG260.get_channel_status()`
- By default, the pressure is always returned in Pa regardless of the display units. This behavior can be overridden by setting `display_units=True` in `TPG260.get_pressure()`.
- In case an error occurs, you can use `TPG260.get_current_errors()` to get the list of currently active errors and `TPG260.reset_error()` to reset them.
- This communication protocol for 350-series gauges (361, 362 and 366) is similar, so the device class should also be able to work with them. However, it has not been tested.

### DPG202/TPG202 controller

There is a variety of different controllers which implement a similar protocol: DPG202 and TPG202, as well as a variety of RS485-controlled gauges (e.g., CPT200). It is based on requesting parameters with certain 3-digit numbers. These are fairly consistent between the devices, for example, 312 stands for the software version, 740 for pressure, and 349 for the device name. However, different devices implement different subsets of these parameters. The supplied class provides a generic interface through `DPG202.get_value()` and `DPG202.comm()` methods, which, correspondingly, request or set a value of a given parameter given its number (e.g., 740) and datatype (e.g., "string", "u\_expo\_new", or "u\_short\_int"). Both of these pieces of information are usually provided in the controller or gauge manual in the `Parameter overview` (or similar-named) section. Currently the device class provides only the most basic functionality:

```
>> from pylablib.devices import Pfeiffer
>> gauge = Pfeiffer.DPG202("COM5")
>> gauge.get_pressure() # pressure in Pa
9.78E4
>> gauge.get_value(740, "u_expo_new") # request the parameter directly, yields_
↪ pressure in mBar
9.78E2
>> gauge.close()
```

---

**Note:** Basic sensors communication concepts are described on the corresponding [page](#)

---

### Leybold pressure gauges

Leybold manufactures a range of pressure gauges and controllers with several different standards and communication protocols. The code has been tested with Leybold ITR90 pressure gauge using its built-in RS232 connection.

The main device classes are `pylablib.devices.Leybold.ITR90`.

### Software requirements

The devices provide a bare RS232 interface, so any appropriate USB-to-RS232 adapter should work.

## Connection

Since the devices are identified as COM ports, they use the standard *connection method*, and all you need to know is their COM-port address (e.g., COM5):

```
>> from pylablib.devices import Leybold
>> gauge = Leybold.ITR90("COM5")
>> gauge.close()
```

## Operation

### ITR90

The operation of this gauge is fairly straightforward, but there is a couple of points to keep in mind:

- Device operates by constantly streaming its status updates. To get the most recent and most consistent data, you can use `ITR90.get_update()`. This is also how you access the gauge status and error states.
- By default, the pressure is always returned in Pa regardless of the display units. This behavior can be overridden by setting `display_units=True` in `ITR90.get_pressure()`.

---

**Note:** Basic sensors communication concepts are described on the corresponding *page*

---

## Kurt J. Lesker pressure gauges

KJL manufactures a range of pressure gauges and controllers with several different standards and communication protocols. The code has been tested with KJL300 pressure gauge using its built-in RS232 connection.

The main device classes are `pylablib.devices.KJL.KJL300`.

## Software requirements

The devices provide a bare RS232 interface, so any appropriate USB-to-RS232 adapter should work.

## Connection

Since the devices are identified as COM ports, they use the standard *connection method*, and all you need to know is their COM-port address (e.g., COM5):

```
>> from pylablib.devices import KJL
>> gauge = KJL.KJL300("COM5")
>> gauge.close()
```

## Operation

### KJL300

The operation of this gauge is fairly straightforward, but there is a couple of points to keep in mind:

- Even standard RS232 operation requires specifying the device RS485 address. IT can be specified using `addr` parameter on creation. By default, the class assumes the factory default of 1, but if it is ever changed on the device, it needs to be specified correctly.
- By default, the pressure is always returned and set in Pa regardless of the display units.

---

**Note:** General device communication concepts are described on the corresponding [page](#).

---

## 1.2.5 Basic lasers

### Basic example

Basic lasers (such as pump lasers) usually only have very basic power-related functionality: turning it on and off, setting power, and controlling and/or requesting the shutter state:

```
>> laser = LaserQuantum.Finesse("COM1") # connect to the laser
>> laser.set_output_power(10.) # set 10W output power
>> laser.enable() # enable the laser
>> laser.get_output_power() # laser hasn't ramped up yet
0.1
>> time.sleep(30.) # wait until the ramp up is done
>> laser.get_output_power()
10.0
>> laser.enable(False)
>> laser.close()
```

### Lighthouse Photonics Sprout

Lighthouse Photonics Sprout laser implements the same basic functionality, with some small additions like reading the interlock status, output mode, temperatures, etc.

The device class is `pylablib.devices.LighthousePhotonics.SproutG`.

Since the device shows up as a COM port, it uses the standard *connection method*, and all you need to know to connect is its COM-port address:

```
from pylablib.devices import LighthousePhotonics
laser = LighthousePhotonics.SproutG("COM1")
laser.close()
```

### Laser Quantum Finesse

Laser Quantum Finesse laser implements the same basic functionality, with some small additions like controlling the shutter, reading the driving current, temperatures, etc.

The device class is `pylablib.devices.LaserQuantum.Finesse`.

Since the device shows up as a COM port, it uses the standard *connection method*, and all you need to know to connect is its COM-port address:

```
from pylablib.devices import LaserQuantum
laser = LaserQuantum.Finesse("COM1")
laser.close()
```

---

**Note:** General device communication concepts are described on the corresponding [page](#).

---

### 1.2.6 M2 Solstis laser

Solstis is a Ti:Saph laser produces by M2. It is controlled via IceBloc controller unit, which communicates with the PC via a network connection.

The main laser class is `pylablib.devices.M2.Solstis`.

#### Software requirements

The device provides a bare network interface, so no additional software is required. However, the device and the local network need to be appropriately configured, such that the PC is and the laser are in the same local network and have static IPs.

In order to access some advanced features, you will need a `websocket-client` package, which is not installed with pylablib by default. You can obtain it from PyPi either separately as

```
pip install websocket-client
```

or with the expanded pylablib version

```
pip install pylablib[devio-full]
```

#### Connection

The laser is identified by its IP address (typically starting with `192.168.1`, if it is on the local network) and the port:

```
>> from pylablib.devices import M2
>> laser = M2.Solstis("192.168.1.2", 34567)
>> laser.close()
```

The port is set up in the `Remote` interface row of the `Network Settings` menu of the laser web interface. There you also need to provide the correct IP address of the controlling PC and enable the remote interface; otherwise the connection will be rejected by the laser.

In addition, you can enable websocket interface option, which is used to send request directly though the device web interface. It is used for some options which are unavailable otherwise, such as enabling or disable the wavemeter connection, receiving some additional status information, and performing more robust control. Note that for proper operation the web interfaces should be opened in the browser and logged in.

#### Operation

The method names are pretty self-explanatory, and mostly correspond directly to the operations in the web interface. Note that, due to the remote interface organization, `terascan` requires two methods to start: first `Solstis.setup_terascan()` to specify parameters, and then `Solstis.start_terascan()` to start it.

One should note, that the device operation is not very stable, and occasionally some errors and crashes arise. These can range from failed wavelength tuning and `terascan`, to `terascan`s failing in exotic ways (e.g., the remote interface suggests that the scan is in progress while the web interface reports a crash), to complete device failure requiring Ice Bloc power cycling.

The device class attempts to somewhat mitigate it by providing relatively a robust stopping method `Solstis.stop_all_operation()`, which tries to set the devices to the default idle state. It uses web interface to get a better information about the laser crashing and send additional stopping commands. It also performs additional steps to stop scans and put the laser in an operation state after a failure, such as starting quick small fine and terascans, and tuning to a nearby frequency.

---

**Note:** General device communication concepts are described on the corresponding [page](#).

---

## 1.2.7 Toptica iBeam Smart laser

Toptica iBeam Smart is a series of CW diode lasers from Toptica. The software has been tested with the standard 633nm laser.

The main device class is `pylablib.devices.Toptica.TopticaIBeam`.

### Software requirements

The device is connected to the PC via RS232 or USB. RS232 simply requires a COM-port controller on the PC, which in most cases is a USB-to-Serial adapter. Such adapters normally come with their standard drivers. The USB version simply involves a built-in USB-to-Serial converter (e.g., a standard FTDI chip), so it also shows up as a virtual COM port. Hence, it requires relatively standard drivers, which are either included with the laser, or can be download from the [manufacturer's website](#), for example, together with the TOPAS control software.

### Connection

Since the devices are identified as virtual COM ports, they use the standard [connection method](#), and all you need to know is their COM-port address (e.g., COM5) and, possibly, baud rate, if it is different from the standard 115200 baud:

```
>> from pylablib.devices import Toptica
>> laser1 = Toptica.TopticaIBeam("COM5")
>> laser2 = Toptica.TopticaIBeam(("COM10", 38400)) # in case of 38400 baud connection
>> laser1.close()
>> laser2.close()
```

### Operation

#### Power and output control

Usually the laser has the main power control and one or several (up to 5) output channels, which can be controlled separately. To turn the whole laser on or off, you can use `TopticaIBeam.enable()`, while each channel is controlled using `TopticaIBeam.enable_channel()`. The power is set independently for each channel via `TopticaIBeam.set_channel_power()`. The actual output power can be queried using `TopticaIBeam.get_output_power()`.

### Detailed info

The most detailed information about the laser can be obtained using `TopticaIBeam.get_full_data()` method. It outputs a detailed report generated by the laser, which contains most of the adjustable parameters.



## Notes and issues

Occasionally the laser communication falls into an error state, where replies are lagging behind the requests (i.e., instead of replying to the issued command, the devices replies to the previous one). This is especially likely if several commands are issued in a rapid succession. If this happens, the laser should be rebooted using `TopicalBeam.reboot()` method.

---

**Note:** General device communication concepts are described on the corresponding [page](#).

---

### 1.2.8 Sirah Matisse laser

Matisse is a family of Ti:Saph and dye ring lasers produces by Sirah.

The main laser class is `pylablib.devices.Sirah.SirahMatisse`.

#### Software requirements

The device requires Matisse Commander software supplied by the manufacturer. When it is installed, it shows up as a VISA resource and can be accessed without further requirements.

#### Connection

The laser is identified by its VISA address, typically looking like "USB0::0x17E7::0x0102::01-01-10::INSTR":

```
>> import pylablib as pll
>> pll.list_backend_resources("visa")
('USB0::0x17E7::0x0102::01-01-10::INSTR',)
>> from pylablib.devices import Sirah
>> laser = Sirah.SirahMatisse("")
>> laser.close()
```

#### Operation

The method names are pretty self-explanatory, and mostly correspond directly to the operations in the Matisse Commander. However, only the basic tuning and scanning functions supplied by the interface are provided, and the more advanced ones like scanning BRF/etalon or interfacing with a wavemeter need to be implemented by the user based on the defined methods.

Note that depending on the specific model not all methods are available, e.g., reference cell locking is not available in TR/DR configuration.

---

**Note:** General device communication concepts are described on the corresponding [page](#).

---

### 1.2.9 Tektronix oscilloscopes

Tektronix produces a large number of very widespread oscilloscopes. They have strongly overlapping, though not entirely identical, interfaces. The library has been tested with TDS2002B, TDS2004B, and DBO2014B.

The generic oscilloscope class is `pylablib.devices.Tektronix.ITektronixScope`, and the derived classes for specific devices are `pylablib.devices.Tektronix.TDS2000` of TDS2000 series and `pylablib.devices.Tektronix.DPO2000` for DPO2000/MSO2000 series.

## Software requirements

These oscilloscopes use NI VISA communication interface. Hence, it requires NI VISA Runtime, which is freely available from the [National Instruments website](#)

## Connection

The devices are identified by their VISA connection strings, which typically start with `USB0::0x0699`, e.g., `"USB0::0x0699::0x0364::C000001::INSTR"`. To get a list of all connected VISA-enabled devices, you can run `pylablib.list_backend_resources("visa")`:

```
>> import pylablib as pll
>> pll.list_backend_resources("visa")
('USB0::0x0699::0x0364::C000001::INSTR',)
>> from pylablib.devices import Tektronix
>> osc = Tektronix.TDS2000("USB0::0x0699::0x0364::C000001::INSTR")
>> osc.close()
```

## Operation

The method names are usually pretty self-explanatory. A typical operation involves setting up channels, scales, and trigger options, acquiring a waveform, and reading the result:

```
from pylablib.devices import Tektronix
osc = Tektronix.TDS2000("USB0::0x0699::0x0364::C000001::INSTR") # connect to the_
↪oscilloscope
osc.enable_channel([1,2]) # enable channels
osc.set_horizontal_span(0.1) # set up horizontal and vertical spans
osc.set_vertical_span("CH1", 1)
osc.set_vertical_span("CH2", 0.1)
osc.setup_edge_trigger("CH1", 0., "dc", "rise") # set up edge trigger on channel 1_
↪at 0V threshold
osc.grab_single(wait_timeout=10.) # grab a single waveform and wait for up to 10s to_
↪finish acquisition
sweeps = osc.read_multiple_sweeps([1,2]) # read out the waveforms
osc.close()
```

However, there is a couple of points to keep in mind:

- The acquisition is controlled using `grab_` methods. Generally, the most convenient way is to use `ITektronixScope.grab_single()` to acquire a single waveform (analogous to pressing a Single button on the oscilloscope panel). By default, this method waits until the acquisition is complete (i.e., the oscilloscope is triggered and the waveform is completely acquired) before continuing. You can also set `wait=False` to perform other operations in the meantime. The acquisition status can be queried via `ITektronixScope.is_grabbing()`, which returns `True` while the trigger is armed or while the data is recording, and `False` after the acquisition is done.
- It appears that the software trigger does not work some time (~500 ms) after the acquisition is set up. If it is invoked in `ITektronixScope.grab_single()` method by supplying `software_trigger=True`, a 300ms delay is added automatically. However, if you invoke it manually using `ITektronixScope.force_trigger()`, you should keep it in mind.

- The waveform transfer is usually performed via `ITektronixScope.read_sweep()` or `ITektronixScope.read_multiple_sweeps()` methods. Since the waveform is transferred in raw form, it requires a preamble data (vertical and horizontal scales and offsets, data format, etc.) to translate into physical units. By default, it is acquired every time before the waveform transfer, which takes some time (up to ~200ms). Alternatively, one can acquire a preamble once and use it in subsequent reading. This method is faster, but will result in an incorrect scaling if the parameters are changed in the meantime (either remotely, or directly on the oscilloscope):

```
>> wfmpres = osc.osc.get_wfmpre([1,2])
>> %time sweeps = osc.read_multiple_sweeps([1,2])
Wall time: 2.2 s
>> %time sweeps = osc.read_multiple_sweeps([1,2], wfmpres=wfmpres)
Wall time: 450 ms
```

- The device class attempts to determine the number of channels automatically on connection, based on which requests raise device errors. However, this process takes some time, and sometimes can raise errors on not fully SCPI-compliant devices. If that is the case, it is always possible to supply the number of channels on construction:

```
>> osc = Tektronix.TDS2000("USB0::0x0699::0x0364::C000001::INSTR") # use ↪
↪ autodetection
>> osc.get_channels_number()
2
>> osc.close()
>> osc = Tektronix.TDS2000("USB0::0x0699::0x0364::C000001::INSTR", nchannels=2)
↪ # specify manually
```

---

**Note:** General device communication concepts are described on the corresponding [page](#).

---

## 1.2.10 NI DAQmx interface

National Instruments produces lots of different data acquisition devices, which support digital and analog input and output, both immediate and clocked (depending on the exact device). They are controlled via a very universal **NI DAQmx** interface. This interface is implemented in `python-nidaqmx` package, which provides a fairly close to original functionality, but with much more convenient Python wrappers. Pylablib implements a relatively thin wrapper around this package to present it in a way similar to the other device classes, and to simplify common tasks such as setting up voltage and counter input channels.

The main daq class is `pylablib.devices.NI.NIDAQ`. It has been tested with NI PCIe-6323, NI USB-6008, and NI USB-6363.

### Software requirements

This interface uses NI DAQmx library, which is freely available on the [National Instruments website](#). Additionally, it needs `python-nidaqmx` package (not to be confused with `pydaqmx`). It is not automatically installed with the base version of `pylablib`, and can be obtained from PyPi either separately as

```
pip install nidaqmx
```

or with the expanded `pylablib` version

```
pip install pylablib[devio-full]
```

## Connection

The devices are identified by their name, such as "Dev1". To list all of the connected devices together with their basic information, you can run `NI.list_nidaqmx_devices`:

```
>> from pylablib.devices import NI
>> NI.list_nidaqmx_devices()
[TDDeviceInfo(name='Dev1', model='USB-6008', serial_number='01234567')]
>> daq = NI.NIDAQ("Dev1")
>> daq.close()
```

## Operation

The typical use case involves setting up different input and output channels, starting acquisition, and acquiring some number of samples:

```
from pylablib.devices import NI
daq = NI.NIDAQ("Dev1")
daq.add_voltage_input("vin", "ai0") # add voltage input named "vin" on the terminal
    ↳ "ai0"
daq.add_voltage_input("vin2", "ai1", rng=(-1,1)) # add second channel with a smaller
    ↳ range of +/- 1V
daq.add_digital_input("din", "port0/line0")
daq.setup_clock(100) # setup 100Hz sampling clock
trace = daq.read(100) # start acquisition, read finite number of samples, and stop
    ↳ it again
# now do continuous acquisition + processing loop
nsamples = 0
daq.start() # start continuous acquisition
while nsamples<1000:
    sample = daq.read()
    ... process sample
    nsamples+=1
daq.stop()
```

The class provide basic methods to set up analog, digital, and counter inputs, and analog and digital outputs. All the analog and digital inputs are synchronized to the same clock, which is the default analog input sample clock (ai/SampleClock) by default. It is also possible to set up the external clock via `NIDAQ.setup_clock()` and export the sampling clock via `NIDAQ.export_clock()`. Not that not all devices support clocked digital inputs, which means that setting up digital inputs there would raise an error.

By default, the counter inputs are synchronized to the same clock, although it is possible to change that. The counter inputs have 3 modes for output values: bare counter (accumulates the number of counts), differential (number of new counts between the two sampling points), and rate (same as differential, but normalized by the sampling rate). In case of external clock, when the sampling rate is a priori unknown, it might be useful to setup a clock rate counter input to determine this clock rate via `NIDAQ.add_clock_period_input()`.

Acquisition is controlled with `NIDAQ.start()` and `NIDAQ.stop()` methods, and the readout is performed via `NIDAQ.read()`. The result of this is always a 2D numpy array, where the first index corresponds to samples and the second to channels. The order of channels can be obtained from `NIDAQ.get_input_channels()`.

The outputs can be either analog or digital. The digital outputs are always immediate, i.e., they immediately produce and hold the latest output value. The analog outputs can work in two modes: either immediate, or clocked. The mode is set up via `NIDAQ.setup_voltage_output_clock()`. In this case, it is possible to output a list of values, which produces a waveform clocked according to the specified clock: either a separate clock source (default), or the analog input clock, which makes voltage input and output synchronized.

---

**Note:** General device communication concepts are described on the corresponding [page](#).

---

## 1.2.11 Generic AWGs

There is a large variety of Arbitrary Waveform Generators, which have very similar characteristics and communication interface.

The generic AWG class is `pylablib.devices.AWG.GenericAWG`, and the derived classes for specific devices are `pylablib.devices.AWG.Agilent33500` and `pylablib.devices.AWG.Agilent33220A` for two different Agilent AWGs, `pylablib.devices.AWG.RigolDG1000` for Rigol DG1000 series, `pylablib.devices.AWG.TektronixAFG1000` for Tektronix AFG1000 series, `pylablib.devices.AWG.InstekAFG2000` for Instek GW 2000 series, `pylablib.devices.AWG.RSInstekAFG21000` for Iso-Tech 21000 series (a clone of Instek AFG2000, but with a couple of bugs which needs to be worked around), and `pylablib.devices.AWG.InstekAFG2225` for Instek GW 2225 (slightly advanced two-channel version of Instek AFG2000).

### Software requirements

Most of these AWGs use NI VISA communication interface. Hence, they require NI VISA Runtime, which is freely available from the [National Instruments website](#). However, Instek and Iso-Tech AWGs show up as virtual COM ports, so they require no additional software.

### Connection

The devices are identified by their VISA connection strings, (e.g., "USB0::0x0699::0x0364::C000001::INSTR") or COM-port (e.g., "COM5"). To get a list of all connected VISA-enabled devices, you can run `pylablib.list_backend_resources("visa")`:

```
>> import pylablib as pll
>> pll.list_backend_resources("visa")
('USB0::0x09C4::0x0400::DG1D150200000::INSTR',)
>> from pylablib.devices import AWG
>> dev = AWG.RigolDG1000("USB0::0x09C4::0x0400::DG1D150200000::INSTR")
>> dev.close()
```

### Operation

The method names are usually pretty self-explanatory. A typical operation involves setting up the function, its parameters, and controlling output:

```
from pylablib.devices import AWG
dev = AWG.RigolDG1000("USB0::0x09C4::0x0400::DG1D150200000::INSTR") # connect to the
↪device
dev.set_function("square", 2) # set up square waveform on the second channel
dev.set_duty_cycle(20, 2)
dev.set_output_range((-1, 1), 2) # set output span from -1V to 1V
dev.enable_output(channel=2) # enable output
dev.close()
```

However, there is a couple of points to keep in mind:

- Since the same general class architecture supports both single-channel and multichannel devices, the channel argument is usually close to the end of the argument list and is not mandatory. If it is not supplied, it is chosen to be the current default channel (1 upon creation), which can be set using *GenericAWG.select\_current\_channel()*. Hence, in the example above we can write:

```
dev.select_current_channel(2)  # now all methods assume channel 2
dev.set_function("square")
dev.set_duty_cycle(20)
dev.set_output_range((-1, 1))
dev.enable_output()
```

- Similarly, some methods can be present but not applicable to the particular AWG (e.g., burst trigger related methods, phase synchronization methods, etc.) If this is the case, they will cause an error when called.

---

**Note:** General device communication concepts are described on the corresponding [page](#).

---

## 1.2.12 Miscellaneous Thorlabs devices

Thorlabs has a variety of devices implementing different serial communication protocols, mostly related to optomechanics. Their requirements and general approach are still fairly similar, so they are all collected here.

### Software requirements

Most devices provide either a bare RS232 interface, or a USB connection with a built-in USB-to-RS232 chip. In either case, they are automatically recognized as serial ports, and no additional software is required. The only exception on this page is MFF101/102 motorized flip mount, which belongs to the [Kinesis devices](#) and requires APT software.

### Connection

Most of the devices are identified as COM ports, so they use the standard [connection method](#), and all you need to know is their COM-port address (e.g., COM5):

```
>> from pylablib.devices import Thorlabs
>> wheel = Thorlabs.FW102("COM5")
>> wheel.close()
```

The only exception is MFF101/102, which is identified by its serial number (more details are given at the [Kinesis devices page](#)).

### Operation

#### MFF101/102 flip mount

The class is provided as `pylablib.devices.Thorlabs.MFF`. It allows for control of the flip mirror position, as well as changing its motion parameters and designations of its digital input and output.

## FW102/212 filter wheel

The class is proved as `pylablib.devices.Thorlabs.FW`.

In addition to setting the position, it allows to adjust speed settings and turn the indicator LED off to minimize light contamination. By default, the wheel also “respects bound” between the first and the last position. Usually, when one orders a move from, e.g., position 2 to 6 on a 6-position wheel, it would go along the shortest route, i.e., position 1. If this is an ND filter wheel (e.g., FW102CNEB), this leads to momentary increase of the transmitted power by ND0.5 (about factor of 3) compared to start and stop positions. To avoid that, the class breaks this move into several shorter (no longer than 1/3 of the wheel) moves, which never cross the boundary between the first and the last position. This takes a bit longer (as it requires several consecutive moves), but is generally safer. This behavior can be turned off by setting `respect_bound=False` on class creation.

Note that older version (1.0) of the filter wheel do not support the full range of options and operate on a slightly different protocol. This leads to crashes on at least some of the methods, e.g., `FW.get_position()`. If this is the case, you can try `pylablib.devices.Thorlabs.FWv1` instead.

## MDT693/694 high-voltage source

The class is proved as `pylablib.devices.Thorlabs.MDT69xA`.

The class provides the ability to set and query the voltage on the three channels, as well as to query the total voltage range (it is set by a physical switch on the back panel, and can not be altered remotely).

---

**Note:** General device communication concepts are described on the corresponding [page](#).

---

## 1.2.13 OZ Optics devices

OZ Optics provides a variety of mostly fiber-optics related devices. Pyalbib covers some of its fiber optomechanics solutions: polarization controller, tunable filter and variable attenuator. Their requirements and general approach are fairly similar, so they are all collected here.

### Software requirements

All the devices provide either a bare RS232 interface, or a USB connection with built-in USB-to-RS232 chip. In either case, they are automatically recognized as serial ports, and no additional software is required.

### Connection

The devices are identified as COM ports, so they use the standard [connection method](#), and all you need to know is their COM-port address (e.g., COM5):

```
>> from pylablib.devices import OZOptics
>> ctl = OZOptics.EPC04("COM5")
>> ctl.close()
```

## Operation

### EPC04 fiber polarization controller

The class is proved as `pylablib.devices.OZOptics.EPC04`. It lets the user change the 4 control voltages, switch between DC and AC (scrambling) modes, and change the AC frequency.

### DD100 fiber attenuator

The class is proved as `pylablib.devices.OZOptics.DD100`. It simply lets the user query and change the attenuation, as well as home the device. Note that homing is required once after the device power up, and it might in general sweep over the whole range of attenuations.

### TF100 fiber filter

The class is proved as `pylablib.devices.OZOptics.TF100`. It simply lets the user query and change the central wavelength, as well as home the device. Note that homing is required once after the device power up, and it might in general sweep over the whole range of wavelengths.

---

**Note:** General device communication concepts are described on the corresponding [page](#).

---

## 1.2.14 Miscellaneous devices

There are several miscellaneous device classes, which are collected in this page. All of them implement straightforward serial communication protocol, so the software requirements and the connection approach is the same for all of them.

### Software requirements

All the devices provide either a bare RS232 interface, or a USB connection with a built-in USB-to-RS232 chip. In either way, they are automatically recognized as serial ports, and no additional software is required.

### Connection

The devices are identified as COM ports, so they use the standard [connection method](#), and all you need to know is their COM-port address (e.g., COM5):

```
>> from pylablib.devices import Conrad
>> dev = Conrad.RelayBoard("COM5")
>> dev.close()
```

## Operation

### Conrad relay board

This is a board, which has several externally-controlled relays.



The class is proved as `pylablib.devices.Conrad.RelayBoard`. It simply lets the user query and set the relay states. It also in principle supports communication with several daisy-chained boards, but it has never been tested.

### Generic Arduino class

The class is proved as `pylablib.devices.Arduino.IArduinoDevice`. It implements basic serial communication; the exact command protocol depends on the particular Arduino software written and uploaded by the user.

The main difference from directly using a serial backend is in handling of DTR line, which signal reset to the Arduino board. Unlike the standard backend, connection will not restart the board; instead, there is an explicit `IArduinoDevice.reset_board()` which pulses the DTR line to reset the board.

## 1.3 Data processing

### 1.3.1 Fitting

Class `fitting.Fitter` is a user-friendly wrapper around `scipy.optimize.least_squares()` routine. Dealing with fitting is made more convenient in a couple of ways:

- It is easy to specify the x-parameter name (in the case it is not the first parameter), or specify multiple x-parameters;
- All of the fit and fixed parameters are specified by name; it is easy to switch between any parameter being fit or fixed;
- The wrapper automatically handles complex parameters (split into real and imaginary parts), numpy arrays, lists, or tuples (including nested structures);
- The final parameters (fit and fixed) are returned in a single dictionary indexed by their names;
- The wrapper also returns the fit function with all of the parameters bound to the final fit and fixed values;
- The fit function result is flattened during fitting, so it works for functions returning multi-dimensional (for example, 2D) arrays.

### Examples

Fitting a Lorentzian:

```
def lorentzian(frequency, position=0., width=1., height=1.):
    return height/(1.+4.*(frequency-position)**2/width**2)

## creating the fitter
# fit_parameters dictionary specifies the initial guess
fit_par = {"position":0.5, "height":1.}
fitter = pll.Fitter(lorentzian, xarg_name="frequency", fit_parameters=fit_par)
# additional fit parameter is supplied during the call
fit_par, fit_func = fitter.fit(xdata, ydata, fit_parameters={"width":1.0})
plot(xdata, ydata) # plot the experimental data
plot(xdata, fit_func(xdata)) # plot fit result
```

Fitting a sum of complex Lorentzians with the same width:

```
def lorentzian_sum(frequency, positions, width, amplitudes):
    # list of complex lorentzians
    # positions and amplitudes are lists, one per peak
    lorentzians = [a/(1.+2j*(frequency-p)/width) for (a,p) in zip (amplitudes,
↪positions)]
    return np.sum(lorentzians, axis=0)

## creating the fitter
# fit_parameters dictionary specifies the initial guess
# (complex initial guess for the "amplitude" parameter hints that this parameter_
↪is complex)
fit_par = {"positions":[0.,0.5,1.], "amplitudes":[1.+0.j]*3}
fitter = pll.Fitter(lorentzian_sum, xarg_name="frequency", fit_parameters=fit_par)
# fixed parameter is supplied during the call (could have also been supplied on_
↪Fitter initialization)
fit_par, fit_func = fitter.fit(xdata, ydata, fixed_parameters = {"width":0.3})
plot(xdata, ydata.real) # plot the experimental data
plot(xdata, fit_func(xdata).real) # plot fit result
```

Fitting 2D Gaussian and getting the parameter estimation errors:

```
def gaussian(x, y, pos, width, height):
    return np.exp( -(x-pos[0])**2+(y-pos[1])**2)/(2*width**2) )*height

## creating the fitter
# fit_parameters dictionary specifies the initial guess
fit_par = {"pos":(100,100), "width":10., "height":5.}
fitter = pll.Fitter(gaussian, xarg_name=["x","y"], fit_parameters=fit_par)
xs, ys = np.meshgrid(np.arange(img.shape[0]), np.arange(img.shape[1]), indexing="ij")
↪# building x and y coordinates for the image
# fit_stderr is a dictionary containing the fit error for the corresponding parameters
fit_par, fit_func, fit_stderr = fitter.fit([xs,ys], img, return_stderr=True)
imshow(fit_func(xs, ys)) # plot fit result
```

The full module documentation is given at [pylablib.core.dataproc.fitting](http://pylablib.core.dataproc.fitting).

## 1.3.2 Filtering and decimation

There are several functions present for filtering the data to smooth it or reduce its size. Most of them are thin wrapper around standard numpy or scipy method, but they provide more universal interface which work both with numpy arrays and pandas DataFrames:

- First are the decimation functions: `filters.decimate()` (and its special case `filters.binning_average()`), `filters.decimate_full()` and `filters.decimate_datasets()`. The first one splits the supplied trace into consecutive segments of  $n$  points and compresses them into a single value using the supplied method (e.g., "mean" will average them together, which is used for `filters.binning_average()`). The second one completely decimates the dataset along the given axis (which is essentially identical to using the standard numpy methods such as `np.mean` or `np.max`). The last one decimates several datasets together, which is similar to combining them into a large  $(n+1)$ D array and fully decimating along the given axis:

```
>> trace = np.arange(10)
>> pll.binning_average(trace, 3) # average every block of 3 points to a single_
↪value
array([1., 4., 7.])
>> pll.decimate(trace, 3, dec="max")
```

(continues on next page)

(continued from previous page)

```

array([2, 5, 8])
>> pll.decimate_full(trace, "mean") # same as np.mean(trace)
4.5
>> trace2 = np.arange(10)**2
>> pll.decimate_datasets([trace, trace2], "sum") # same as np.sum([trace,
↪trace2],axis=0)
array([ 0,  2,  6, 12, 20, 30, 42, 56, 72, 90])

```

- Sliding decimation methods `filters.sliding_average()`, `filters.median_filter()` and `filters.sliding_filter()` are related, but use a sliding window of  $n$  points instead of complete decimation of  $n$  points together. They only work for 1D traces or 2D multi-column datasets. Note that `filters.sliding_filter()` is implemented through a simple Python loop, so it is fairly inefficient:

```

>> trace = np.arange(10)
>> pll.sliding_average(trace, 4) # average points in 4-point window (by default,
↪use "reflect" boundary conditions)
array([0.75, 1.5 , 2.5 , 3.5 , 4.5 , 5.5 , 6.5 , 7.5 , 8.25, 8.5 ])
>> pll.sliding_filter(trace, 4, "max") # find maximum of points in 4-point window
array([2, 3, 4, 5, 6, 7, 8, 9, 9, 9])

```

- Next are convolution filters which operate by convolving the trace with a given kernel function. These involve `filters.gaussian_filter()` (and `filters.gaussian_filter_nd()`, which is simply a wrapper around `scipy.ndimage.gaussian_filter()`), and a more generic `filters.convolution_filter()`. Related are infinite impulse response (IIR) filter `filters.low_pass_filter()` and `filters.high_pass_filter()`, which mimic standard single-pole low-pass and high-pass filters. In principle, they can be modelled as a convolution with an exponential decay, but the implementation using the recursive filters is more efficient for large widths.
- Finally, there are Fourier filters, which Fourier-transform the trace, scale the transform values, and transform it back to the real domain. These involve the main function `filters.fourier_filter()`, which takes a generic frequency response function, as well as two specific response function generators `filters.fourier_filter_bandpass()` and `filters.fourier_filter_bandstop()`, both generating hard frequency cutoff filters.
- In addition to “post-processing” filters described above, there are also “real-time” filters which serve to filter data as it is acquired, e.g., to filter out temporary noise or spikes. There are two filters of this kind: `filters.RunningDecimationFilter` and `filters.RunningDebounceFilter`. They are implemented as classes, and both have methods to add a new datapoint, to get the current filter value, and to reset the filter.

### 1.3.3 Fourier transform

There is a couple of methods to work with Fourier transform. They are built around `numpy.fft.fft()`, but allow more convenient normalization (e.g., in units of power spectral density), and work better with pandas DataFrames. They also have an option to automatically trim the trace length to the nearest “good” size, which is a product of small primes. This can have fairly strong (up to a factor of several) effect on the transform runtime, while typically trimming off less than 1% of the data.

The main methods are `fourier.fourier_transform()` for the direct transform, `fourier.inverse_fourier_transform()` for the inverse transform, and `fourier.power_spectral_density()` for the power spectral density:

```

>> x = np.random.normal(size=10**5) # normal distribution centered at 0 with a width,
↪of 1
>> PSD = pll.power_spectral_density(x, dt=1E-3) # by default, use density,
↪normalization; assume time step of 1ms

```

(continues on next page)

(continued from previous page)

```
>> df = PSD[1,0] - PSD[0,0]
>> df # total span is 1kHz with 10**5 points, resulting in 0.01Hz step
0.01
>> np.sum(PSD[:,1]) * df # integrated PSD is equal to the original trace RMS squared,
→ which is about 1 for the normal distribution
1.005262206692361
>> np.mean(x**2)
1.005262206692361
```

### 1.3.4 Feature detection

There are several methods for simple feature detection:

- The peak detection, which is usually achieved by the combination of `feature.multi_scale_peakdet()` and `feature.find_peaks_cutoff()`. The first applies difference-of-Lorentzians or difference-of-Gaussians filter, which detects peaks of a particular width. The second finds peaks using a cutoff.
- Another way to find peaks is using `feature.find_local_extrema()`, which finds local minima or maxima in a sliding window of a given width.
- Switching between two states with a noisy trace can be detected using `feature.latching_trigger()`. It implements a more robust approach to find when the trace is above/below threshold by considering two thresholds: a higher “on” thresholds and a lower “off” threshold. It makes the on/off state “latch” to its current value and is robust to small trace fluctuations around the threshold, which would lead to rapid on/off switches in a single-threshold scheme.

### 1.3.5 Miscellaneous utilities

Additionally, there is a variety of small functions to simplify some data analyses and transforms:

- Checking trace properties: `dataproc.utils.is_ascending()`, `dataproc.utils.is_descending()`, `dataproc.utils.is_ordered()`, `dataproc.utils.is_linear()`.
- Sorting by a given column: `dataproc.utils.sort_by()`; work both on pandas and numpy arrays
- Filtering: `dataproc.utils.filter_by()` and `dataproc.utils.unique_slices()` (a simple analog of pandas `pandas.DataFrame.groupby()`, which works on numpy arrays)
- Binary search (both in ordered and unordered 1D arrays): `dataproc.utils.find_closest_arg()`, `dataproc.utils.find_closest_value()`, and `dataproc.utils.get_range_indices()`.
- Traces step analysis and unwrapping: `dataproc.utils.find_discrete_step()` tries to find a single number which divides all values within a reasonable precision, and `dataproc.utils.unwrap_mod_data()` “unwraps” modulo data (e.g., phase, which is defined mod  $2\pi$ ) provided that the steps between two consecutive points are less than  $1/2$  of the module.
- Cutting the trace to the given range, or cutting out a given range: `dataproc.utils.cut_to_range()` and `dataproc.utils.cut_out_regions()`.
- Converting between 2-column “XY” and complex representations: `dataproc.utils.xy2c()` and `dataproc.utils.c2xy()`
- Scalar numerical utilities: `utils.numerical.limit_to_range()` (limit a value to lie in a given range, including option for no limits in one or both directions), `utils.numerical.gcd()` and `utils.numerical.gcd_approx()` (greatest common divisor or its approximate version for non-integer values)

## 1.4 Data storage

Complex data storage in pylablib centers around 2 main components: the multi-level dictionary for representing hierarchical data within the code, and file IO to (among other things) load and store it in a human-readable format.

### 1.4.1 Multi-level dictionary

`dictionary.Dictionary` is an expansion of the standard `dict` class which supports tree structures (nested dictionaries). The extensions include:

- handling multi-level paths and nested dictionaries, with several different indexing methods
- iteration over the immediate branches, or over the whole tree structure
- some additional methods: mapping, filtering, finding difference between two dictionaries
- combined with `pylablib.core.fileio` allows to save and load the content in a human-readable format.

Creating and indexing:

```
>>> d = pll.Dictionary()
>>> d['d/0/x'] = 5
>>> d
Dictionary('d/0/x': 5)
>>> d['d/0/x'] # string path indexing
5
>>> d['d']['0']['x'] # nested indexing
5
>>> d['d','0','x'] # multi-level path indexing
5
>>> d['d',0,'x'] # all path elements are converted into strings
5
>>> d['d/0']['x'] # indexing styles can be freely mixed
5
>>> d['d','0/x']
5
>>> b = d['d'] # indexing a branch yields another Dictionary object
>>> b
Dictionary('0/x': 5)
>>> b['0/x'] = 10 # the branch shares the data with the main dictionary
>>> d
Dictionary('d/0/x': 10)
```

A dictionary can be build from a Python `dict`, which automatically normalizes paths and nested dictionaries:

```
>>> d = pll.Dictionary({ 'a':1, 'b/i':2, 'c':{'i':3, 'ii':4}, 'd/0/x':5 })
>>> d
Dictionary('b/i': 2
'c/i': 3
'c/ii': 4
'd/0/x': 5
'a': 1)
```

**Note:** There are several limitations on the dictionary structure (mostly they involve possible paths and keys):

- As mentioned above, the keys are converted into strings to get the path; therefore, different Python object can merge together (e.g., number 0 and string literal '0'). This also discourages use of some of the objects with

“underdefined” (implementation dependent) representations, for example, floating point numbers.

- Since the `'/'` symbol is used to split different path entries, it can’t be used inside a single-level key. It is possible to re-define this symbol on dictionary creation; however, it might lead to compatibility issues.
- Empty keys are not allowed. When building a path, they are automatically dropped, so `'a/b'`, `'a/b/'`, `'a///b///'` all correspond to the same path.
- One path can either correspond to a branch node, or a leaf node. In other words, one path can’t be a prefix of other paths and also contain data: structures like `pll.Dictionary({'a':1, 'a/b':2})` are not allowed. To get around this, one can define a specific “data key” not used anywhere else, and store data in a node under that key (e.g., with the data key `'#'` the example before turns into a valid structure `pll.Dictionary({'a/#':1, 'a/b/#':2})`).

Thus, it is generally recommended to only use strings or non-negative integers as keys, and apply the same restrictions to them as to the Python variable names (with the addition of names starting with a digit).

---

## 1.4.2 File IO

`pylablib.core.fileio` contains several function for saving and loading data into different kinds of files: binary (`loadfile.load_bin()` and `savefile.save_bin()`), CSV (`loadfile.load_csv()` and `savefile.save_csv()`), or dictionary (`loadfile.load_dict()` and `savefile.save_dict()`).

### Binary files

The first (binary files) closely corresponds to `numpy.fromfile`. In addition, it also allows automatic conversion into pandas arrays, setting column names, and skipping some number of bytes from the start:

```
>> table = np.arange(6).reshape((3,2))
>> pll.save_bin(table, "table.dat", dtype="<f8")
>> pll.load_bin("table.dat", columns=["Column1", "Column2"], dtype="<f8")
  Column1  Column2
0        0.0      1.0
1        2.0      3.0
2        4.0      5.0
```

Furthermore, there is an option to save the binary data with a preamble dictionary file, which describes its structure (columns, dtype, etc.) This way, one does not have specify these parameter in the loading code:

```
>> table = pd.DataFrame({"C1":range(3), "C2":range(3)**2/3})
>> table
   C1      C2
0  0  0.000000
1  1  0.333333
2  2  1.333333
>> pll.save_bin_desc(table, "table.dat")
>> pll.load_bin_desc("table.dat")
   C1      C2
0  0.0  0.000000
1  1.0  0.333333
2  2.0  1.333333
>> np.fromfile("table_data.bin", "<f8").reshape((3, 2)) # the data is still stored_
↪ in the regular binary format
array([[0.          , 0.          ],
```

(continues on next page)

(continued from previous page)

```
[1.      , 0.33333333],
[2.      , 1.33333333]])
```

Note that only homogeneous data (i.e., all columns having the same type) is currently supported. That's why the first column got converted from integers into reals.

## CSV files

The functionality of the second one mimics pandas `read_csv`, but offers a bit more flexibility with more complicated values in columns, such as tuples or binary strings:

```
>> table = pd.DataFrame({ "C1":np.arange(3), "C2":[(i**2,i**3) for i in range(3)] })
>> table # the second columns contains tuples
   C1      C2
0  0  (0, 0)
1  1  (1, 1)
2  2  (4, 8)
>> pll.save_csv(table, "table.csv")
>> pll.load_csv("table.csv", dtype="generic") # need to specify generic values type,
↳ which handle complicated cases, but is somewhat slower
   C1      C2
0  0  (0, 0)
1  1  (1, 1)
2  2  (4, 8)
```

In addition, its default settings are a bit different: the column separator is a whitespace, the column names are contained in the comment string (which removes occasional ambiguity), and the creation date string is appended by default. Hence, the content of the file created above is

```
# C1      C2
0  (0, 0)
1  (1, 1)
2  (4, 8)

# Saved on 2021/01/01 12:00:00
```

Note that currently it operates only with simple flat tables and does not support advanced pandas features such as index or multi-index. If these are required, you can use `savefile.save_csv_desc()` and `loadfile.load_csv_desc()`. Similarly to `savefile.save_bin_desc()` and `loadfile.load_bin_desc()`, it saves a dictionary containing additional description; however, the table is inlined by default, so only one file is generated:

```
>> table = pd.DataFrame({ "C1":np.arange(3), "C2":[(i**2,i**3) for i in range(3)] },
↳ index=np.arange(3)+10)
>> table # non-trivial index column
   C1      C2
10  0  (0, 0)
11  1  (1, 1)
12  2  (4, 8)
>> pll.save_csv(table, "table.csv")
>> pll.load_csv("table.csv", dtype="generic") # index is lost
   C1      C2
0  0  (0, 0)
1  1  (1, 1)
2  2  (4, 8)
```

(continues on next page)

(continued from previous page)

```
>> pll.save_csv_desc(table, "table.dat")
>> pll.load_csv_desc("table.dat")  # index is preserved (also note that here dtype is
→ "generic" by default)
      C1      C2
10    0  (0, 0)
11    1  (1, 1)
12    2  (4, 8)
```

## Dictionary files

Finally, dictionary saving and loading operates with *dictionary* objects. It is generally useful to load or save various heterogeneous settings or parameters, such as device parameters, data processing parameters, and GUI or device state. It supports most basic Python data types as values: standard scalar types (integers, reals, complex numbers, strings, booleans, None), containers (tuples, lists, dictionaries, sets, including nested ones), binary and raw string representation (e.g., `b"\x00"` or `r"m\n\o"`), short numpy arrays (represented as, e.g., `"array([1, 2, 3])"`), and inline tables (which are interpreted as pandas table by default). The only common data type not included is named tuples; they get automatically converted to regular tuples on saving.

The dictionary files have the `key value` line formats and typically use full paths (as opposed to, say, XML hierarchy), which makes them easier to inspect and parse without pylablib. For example, the dictionary from the previous section will be saved as

```
b/i 2
c/i 3
c/ii 4
d/0/x 5
a 1
```

With more complicated data types, it might look more like

```
process/points  array([1., 2., 3.])
process/default/frequency  10+2.j
# Lines starting with # are treated as comments
plot/position  [(0,0), (1,1), (2,3)]
plot/label  r"$\nu_0$"
# Keys do not have to be in any particular order
process/default/amplitude  5.
```

which results in a dictionary

```
Dictionary('plot/label': $\nu_0$
'plot/position': [(0, 0), (1, 1), (2, 3)]
'process/default/amplitude': 5.0
'process/default/frequency': (10+2j)
'process/points': [1. 2. 3.]
```

The format also supports hierarchy using `//branch` to mark a start of sub-branch and `///` to mark its end. For example, the dictionary above can be also saved as

```
//process
# indentation is not required, but helps to see the structure
points  array([1., 2., 3.])
default/frequency  10+2.j
default/amplitude  5.
```

(continues on next page)



(continued from previous page)

```

///

//plot
    position    [(0,0), (1,1), (2,3)]
    label      r"$\nu_0$"
///

```

Finally, it is possible to specify inline tables using special comment lines. For example,

```

# The key without the value marks the path to the table within the dictionary
data/table
## Begin table
1    1.j
2    4.j
3    9.j
## End table

```

produces a dictionary containing pandas DataFrame:

```

Dictionary('data/table':
  0          1
0 1  0.000000+1.000000j
1 2  0.000000+4.000000j
2 3  0.000000+9.000000j )

```

## 1.5 Various utilities

### 1.5.1 File system

There is a number of methods which are minor expansions of the built-in file utilities:

- Accessing and changing file times: `utils.files.get_file_creation_time()`, `utils.files.get_file_modification_time()`, `utils.files.touch()` (update the modification date).
- Generating new file names (e.g., for storing a new dataset): `utils.files.generate_indexed_filename()` and `utils.files.generate_prefixed_filename()`.
- Some path analysis methods: `utils.files.fullsplit()`, `utils.files.normalize_path()`, `utils.files.paths_equal()`, `utils.files.relative_path()`; a lot of these have also been implemented in `pathlib` module, and are kept for backwards compatibility.
- Checking if a string is a valid path: `utils.files.is_path_valid()`.
- File copying and moving, which also creates containing folders if necessary: `utils.files.copy_file()`, `utils.files.move_file()`.
- Folder creation and cleaning: `utils.files.ensure_dir()`, `utils.files.remove_dir()`, `utils.files.remove_dir_if_empty()`, `utils.files.clean_dir()`.
- Analyzing folder content: `utils.files.list_dir()`, `utils.files.list_dir_recursive()`, `utils.files.dir_empty()`, `utils.files.walk_dir()`. Compared to the built-in methods, allows for more complicated (e.g., regex) filters for listed files and folders, as well as for visited folders.
- Copying, moving, and comparing folders: `utils.files.copy_dir()`, `utils.files.move_dir()`, `utils.files.cmp_dirs()`; like methods above, allows for regex filters for files and folders.

- Retrying versions of most of the above methods: e.g., `utils.files.retry_move()` or `utils.files.retry_clean_dir()`. These functions try to copy/move/remove files or folders several times if errors arise, in case the files or folders are only temporarily blocked. Useful when, e.g., using network shares or some software which makes files or folders unavailable for a short period of time.
- Wrapping methods for working with zip files: `utils.files.zip_folder()`, `utils.files.zip_file()`, `utils.files.zip_multiple_files()`, `utils.files.unzip_folder()`, `utils.files.unzip_file()`.

## 1.5.2 Network

There is a simple wrapper class `utils.net.ClientSocket`, which simplifies some operations with the built-in `socket` module. In addition, it also implements a couple of higher-level ways to send the data: either fixed length (as in the usual socket), with the length prepended (in case the total length is initially unknown at the receiving end), or using a delimiter to mark the end of the message.

In addition, there are several methods for gaining local or remote host information (`utils.net.get_local_addr()`, `utils.net.get_all_local_addr()`, `utils.net.get_local_hostname()`, `utils.net.get_all_remote_addr()`, `utils.net.get_remote_hostname()`), receiving JSON-formatted values (`utils.net.recv_JSON()`), and listening on a given port (`utils.net.listen()`).

## 1.5.3 Strings

There are several string manipulation functions present:

- Powerful to/from string conversion. The main function are `utils.string.to_string()` and `utils.string.from_string()`, which can convert a large variety of values: simple scalar values (numbers, strings, bools, None), containers (lists, tuples, sets, dictionaries), escaped and byte strings (e.g., `b"\x00"`), complex types such as numpy arrays (represented as, e.g., `"array([0, 1, 2, 3, 4])"`). The latter version requires setting `use_classes=True` in `utils.string.to_string()`, which is not enabled by default to make the results more compatible with other parsers:

```
>> pll.to_string(np.arange(5)) # by default, use the standard str method, which
↳ makes array look like a list
'[0, 1, 2, 3, 4]'
>> pll.from_string('[0, 1, 2, 3, 4]') # gets converted back into a list
[0, 1, 2, 3, 4]
>> pll.to_string(np.arange(5), use_classes=True) # use representation class
'array([0, 1, 2, 3, 4])'
>> pll.from_string('array([0, 1, 2, 3, 4])') # get converted back into an array
array([0, 1, 2, 3, 4])
```

More complex data classes can be added using `utils.string.add_conversion_class()` and `utils.string.add_namedtuple_class()`:

```
>> NamedTuple = collections.namedtuple("NamedTuple", ["field1", "field2"])
>> nt = NamedTuple(1,2)
>> nt
NamedTuple(field1=1, field2=2)
>> pll.to_string(nt, use_classes=True) # class is not registered, so use the
↳ default tuple representation
'(1, 2)'
>> pll.add_namedtuple_class(NamedTuple)
>> pll.to_string(nt, use_classes=True) # now the name marker is added
'NamedTuple(1, 2)'
```

(continues on next page)

(continued from previous page)

```
>> pll.from_string('NamedTuple(1, 2)')
NamedTuple(field1=1, field2=2)
>> DifferentNamedTuple = collections.namedtuple("DifferentNamedTuple", ["field1",
↪ "field2"])
>> pll.from_string('DifferentNamedTuple(1, 2)') # note that if the class is not_
↪ registered, it can't be parsed, so the string is returned back
'DifferentNamedTuple(1, 2)'
```

Furthermore, there is a couple of auxiliary string functions to parse more complicated situations: `utils.string.escape_string()` and `utils.string.unescape_string()` for escaping and unescaping string with potentially confusing or unprintable characters (e.g., quotation marks, spaces, new lines); `utils.string.from_string_partial()`, `utils.string.from_row_string()`, `utils.string.extract_escaped_string()` to determine and extract the first value in a string which potentially has several values.

- Comparing and searching string: `utils.string.string_equal()` (compare string using different rules such as case sensitivity), `utils.string.find_list_string()`, `utils.string.find_dict_string()` (find string in a list or a dictionary using different comparison rules).
- Filtering strings: `utils.string.get_string_filter()`, `utils.string.sfglob()`, and `utils.string.sfregeex()`. Creates filter functions which may include or exclude certain string patterns; these filter functions can be later used in, e.g., file-related methods such as `utils.files.list_dir()`.

## 1.5.4 Misc utilities

A variety of small useful methods and classes:

- Dictionary manipulation functions: `utils.general.any_item()` (get a random dict key-value pair), `utils.general.merge_dicts()` (merge several dictionaries together), `utils.general.filter_dict()` (filter dictionary according to key or value), `utils.general.map_dict_keys()`, `utils.general.map_dict_values()`, `utils.general.to_dict()` (convert a dict or a list of pairs into a dictionary, using a default value for a non-pair list elements), `utils.general.invert_dict()` (turn keys into values and vice versa).
- List manipulation functions: `utils.general.flatten_list()` (flatten a nested list structure), `utils.general.partition_list()` (split a list into two lists according to a predicate), `utils.general.split_in_groups()` (split list into several groups according to a key function), `utils.general.sort_set_by_list()` (convert set into a list, whose values are sorted according to a second supplied list), `utils.general.compare_lists()` (compare two lists and return their intersection and differences).
- `utils.general.DummyResource`: a “dummy” resource class, which can be used in a `with` block but does nothing; can be used to, e.g., replace multi-threading resources such as locks to turn them off.
- Unique ID generators: `utils.general.UIDGenerator` and `utils.general.NamedUIDGenerator`, which generate unique names (based on a counter), with a thread-safe option (useful to create, e.g., unique data markers).
- Timekeeping: `utils.general.Countdown` for single shot and `utils.general.Timer` for repeating tasks. Simplifies dealing with operation timeouts: checking how much time is left (including options for infinite timeout), checking if timeout is passed, resetting, etc.
- Script restarting via `utils.general.restart()` (thread-controller style applications can also use `thread.controller.restart_app()` for a more managed restart).
- `utils.general.StreamFileLogger`, which can be set up to log all outputs into a stream (e.g., `stdout`):

```
from pylablib import StreamFileLogger
import sys
sys.stderr = StreamLogger("logerr.txt", sys.stderr) # replace stderr stream with
↳ a logged version
# perform some tasks ...
sys.stderr = sys.stderr.stream # revert back, if necessary
```

With the code above, all output to `stderr` will be logged into `logerr.txt` to be analyzed later. It can also be set with `autoflush=True` to automatically flush the printed text, which helps with identifying crushing bugs, and it can be supplied with a lock to help separate printouts from different threads.

## 1.6 Change log

This is a list of changes between each version.

### 1.6.1 Version 1.x

Transitioning from version 0.x to version 1.x saw lots of interface changes which break backward compatibility. The previous version of the library can be either obtained on PyPi using `pip install "pylablib<1"`, or by using `legacy` module. Hence, instead of

```
import pylablib as pll
from pylablib.aux_libs.devices import Lakeshore
```

you can write

```
import pylablib.legacy as pll
from pylablib.legacy.aux_libs.devices import Lakeshore
```

#### 1.4.1

- Devices
  - Added Basler pylon-compatible cameras, BitFlow frame grabbers, AlliedVision Bonito cameras, Thorlabs Elliptec stages, PI-E516 piezo controller, and Sirah Matisse laser.
  - Minor additions to Cryocon temperature controller, Cryomagnetics LM510 level meters, and NI DAQmx DAQs. Improved performance of PCO cameras at high frame rates.
  - Multiple minor bug fixes and improved support of specific models.
- Added encoding argument to file loading.
- Improved color images support in image plotter, minor additions to trace plotter.
- Added real-time binning and debounce filters.

#### 1.4.0

- Added Photometrics cameras and Cryocon temperature controllers.
- More consistent cameras interface: attributes properties, fast chunks (former `fastbuff`) readout, frame info formats.

- Added new simple GUI elements: multiline edits, enum labels.
- Expanded image and trace plotting widgets.
- Added linear transforms to data processing.
- Minor bugfixes in threading, GUI, devices.

### 1.3.3

- Numpy loads bugfix (fixes compatibility with `numpy>=1.22`).

### 1.3.2

- Added Leybold ITR90 and KJL300 pressure gauges.
- Minor bugfixes in threading and devices.

### 1.3.1

- Added expandable edit boxes and dialog containers.
- Improved Thorlabs devices compliance.
- Additional minor bugfixes in threading, GUI, devices.

### 1.3.0

- General
  - Minor speedups through calls caching.
  - Changed `muxcall` signature to allow multiple special argument values.
- Devices
  - Added Princeton Instruments cameras, IDS uEye cameras (as an option in uc480 cameras backend), Thorlabs Kinesis piezo motor controllers (e.g., KIM101) and quadrature photo-detector controllers (e.g., KPA101).
  - Added RS485 Arcus connection and a simple single-motor stage (DMX-J-SA).
  - Improved reliability if errors are encountered upon connection.
  - Multiple minor bug fixes and improved support of specific models.
- GUI
  - Added widgets: menu dropdown button, scroll area container, area highlighter.
  - Added querying element position and layout shape in layout widgets.
  - Added more utilities methods: querying containing layout, querying top-level parent, deleting widget.
- Threading
  - Added simple profiling through `yappi`.

### 1.2.1

- General
  - Added restarting methods for regular and threaded applications.
- Threading
  - Bugfixes in cameras and camera threads.
  - Bugfixes in streaming.

### 1.2.0

- General
  - Added `timing` context manager for simple code timing checks.
  - Improved RPyC wrapper logging and reliability.
  - Added Anaconda support.
  - Added minor network and file functions.
- Devices
  - Added Newport Picomotor 8742 motor controller, Toptica iBeam Smart laser, older version of Thorlabs FW motorized filter wheel.
  - Added camera frame output format (list or array).
  - Added `use_cavity` option to M2 Solstis laser.
  - Added method for auto-detecting associations between PhotonFocus cameras and frame grabbers.
  - Updated some generic classes (DCAM cameras, Thorlabs TLCamera cameras).
  - Updated SCPI failsafe operation, improved Thorlabs FW reliability.
  - Fixed several minor bugs.
- GUI
  - Rewritten GUI values handling to pass calls in a hierarchical manner. This makes the operation more predictable and overloading the behavior a bit easier.
  - Added out-of-range value action for combo boxes.
  - Fixed `ImagePlotter` incompatibility with the newer `pyqtgraph` versions, added separate x and y axis line cuts selection.
  - Minor layout handling bugfixes.
- Threading
  - Released advanced threading functionality: table/frame streaming, device threads, basic frame processing.
  - Task thread additions: delayed batch job stopping, context manager for task loop pausing.
  - Added argument-dependent call queue limit.
  - Improved threading speed and stability.

### 1.1.0

- General
  - Reorganized the core modules import structure: now `__init__.py` modules are mostly empty, and all the necessary imports are either exposed directly in `pylablib` (e.g., `pylablib.Fitter`), or should be accessed directly by the module (e.g. `pll.core.dataproc.fitting.Fitter`). Intermediate access (e.g., `pll.core.dataproc.Fitter`) is no longer supported.
  - File IO functions (e.g., `read_csv`) can now take file-like objects in addition to paths.
- Devices
  - Added Silicon Software frame grabbers interface and rearranged PhotonFocus code to include both IMAQ and SiliconSoftware frame grabbers.
  - Fixed various compatibility bugs arising for specific versions of Python or dependency modules: Kinesis error with specific pyft232 versions, some DLL-dependent devices errors with Python 3.8+, DLL types in 32-bit Python.
  - Addressed issue with occasional uc480 acquisition restarts, fixed M2 communication report errors.
- GUI and threading
  - Added container and layout management classes in addition to parameter tables for more consistent GUI structure organization.
  - Added `pylablib.widgets` module which combines all custom widgets for the ease of using in layout managers or custom applications.
  - Fixed support for PySide2 Qt5 backed.
  - Renamed `setupUi` -> `setup` for all widgets, and changed the GUI setup organization for many of them (the functioning stayed the same).
  - Reorganized scheduling in `QTaskThread` to treat jobs, commands, and subscriptions more consistently.
  - Added basic data stream management.

### 1.0.0

There have been too many alterations to list here comprehensively. Below is the list of the largest changes.

- General
  - Removed built-in `DataTable` class (together with `core.datatable` subpackage) in favor of `pandas`.
  - Renamed file IO functions: instead of generic `load` and `save` methods there are now more specific `loadfile.load_csv()`, `loadfile.load_dict()`, etc.
  - Removed some legacy modules which are not used in the rest of the library.
  - Renamed or moved certain modules: `core.utils.rpyc` -> `core.utils.rpyc_utils`, `core.fileio.logfile` -> `core.fileio.table_stream`, `core.fileio.binio` -> `core.utils.binio`, `core.devio.backend` -> `core.devio.backencd_comm`, `core.devio.untis` -> `core.utils.units`, `core.dataproc.waveforms` -> `core.dataproc.utils`
- Devices
  - Some legacy devices have been removed, since without access to the hardware it is hard to maintain and expand them. These include most of Agilent devices (33502A amplifier, N9310A microwave generator, HP 8712B and HP 8722D network analyzers, HP 8168F laser), Rigol DSA1030A spectrum analyzer, Tektronix MDO3000 oscilloscope, Vaunix LabBrick generators, Zurich Instruments HF2 and UHF, Andor

Shamrock spectrographs (should be restored in future releases), NuPhoton NP2000 EDFA, PurePhotonics PPCL200 laser, Sirah Matisse laser (should be restored in future releases), Thorlabs PM100 power meter (should be restored in future releases), Lakeshore 370 resistance bridge (should be restored in future releases), MKS 900-series pressure gauges, and some custom devices (Arduino and Olimex AVR boards and Janis-related hardware).

- The main devices package has been moved from `pylablib.aux_libs.devices` (which now refers to the legacy code) to `pylablib.devices`. Module organization has also changed slightly. To find the required modules and device class names, see the *devices list*.
  - Lots of devices' interface has varied slightly, to make the interface more uniform and compatible between different kinds of devices. The changes are usually fairly straightforward (e.g., `move_to` instead of `move`). In many cases the interface was also expanded to include additional available methods.
  - Several devices have been added, generalized, or restructured:
    - \* Combined Thorlabs KDC101 and K10CR1 into a single class `pylablib.devices.Thorlabs.BasicKinesisDevice`, which also accommodates similar kinds of devices.
    - \* Added Arcus Performax2EXStage device for 2-axis controller with a slightly different interface (`pylablib.devices.Arcus.Performax2EXStage`)
    - \* Added *several more AWGs* with similar interfaces
  - Simplified the way external DLLs are *handled*
  - Unified the *error handling*
- GUI and threading
    - Changed module structure
      - \* threading and GUI are now separate sub-packages `core.thread` and `core.gui`
      - \* all widgets are available simply through `pylablib.widgets` (simplifies integration with Qt Designer)
      - \* moved parameter tables widgets to the core library
    - Renamed some widgets to remove the LV prefix.
    - Interfaces changes in some of the classes: thread controllers, parameter tables, value tables. The changes are mostly cosmetics and involve names and parameters order. Most important changes:
      - \* thread controller methods: `subscribe` -> `subscribe_sync`, `sync_exec` -> `sync_exec_point`,
      - \* thread controller command/query shortcut: `.c` -> `.ca`, `.q` -> `.cs`, `.qi` -> `.csi`, `.qs` -> `.css`
      - \* thread controller variable access uses `.v` shortcut, i.e., instead of `ctl[name]` it is now `ctl.v[name]`
      - \* GUI value storage `ValuesTable/IndicatorValuesTable` are now combined and named as `GUIValues`
      - \* `ParamTable` and `GUIValues` uses `.h` shortcut to access value handlers, i.e., instead of `table[name]` it is now `table.h[name]`
      - \* `ParamTable`, `ImagePlotterCtl`, `TracePlotterCtl` constructor arguments: `display_table` -> `gui_values`, `display_table_root` -> `gui_values_root`
      - \* value-changed signal names in `ParamTable` and `GUIValues`: `changed_event` -> `get_value_changed_signal`



- \* value-changed signal names in value handlers: `value_changed_signal` -> `get_value_changed_signal`
- \* ParamTable methods: `lock` -> `set_enabled`, `add_button(checkable=True)` -> `add_toggle_button`
- \* NumEdit and NumLabel methods: `set_number_format` -> `set_formatter`, `set_number_limit` -> `set_limiter` (the call signature also changed)
- \* renamed signals to multicasts to avoid confusion with built-in Qt signals. Leads to `ThreadController.send_signal` -> `send_multicast`, `ThreadController.process_signal` -> `process_multicast`, `ThreadController` constructor argument `signal_pool` -> `multicast_pool`, class `SignalPool` -> `MulticastPool`, `QSignalThreadCallScheduler` -> `QMulticastThreadCallScheduler`.

## 1.6.2 Version 0.x

### 0.4.1

#### Interface changes

- Slightly changed representations of complex number in to-string conversions depending on the conversion rules ("python" vs "text").

#### Additions

- Devices
  - Added Thorlabs K10CR1 rotational stage (`legacy.aux_libs.devices.Thorlabs.K10CR1`)
  - Added Andor Shamrock spectrographs (`legacy.aux_libs.devices.AndorShamrock`)
  - Expanded Agilent AWG class
  - Added more 32bit dlls
  - Added `list_resources` method to every backend class, which lists available connections for this backend (not available for every backend; so far only works in `legacy.core.devio.backed.VisaDeviceBackend`, `legacy.core.devio.backed.SerialDeviceBackend`, and `legacy.core.devio.backed.FT232BackendOpenError`).
- GUI and threading
  - Added `legacy.aux_libs.gui.helpers.TableAccumulatorThread.preprocess_data` method to pre-process incoming data before adding it to the table
  - Added `update_only_on_visible` argument to `legacy.aux_libs.gui.widgets.trace_plotter.TracePlotter.setupUi` method, and `legacy.aux_libs.gui.widgets.trace_plotter.TracePlotter.get_required_channels` method.

### 0.4.0

#### Interface changes

- Dictionary entries (`legacy.core.fileio.dict_entry`) system has been slightly redesigned: building entries from stored objects has been moved from `legacy.core.fileio.dict_entry.IDictionaryEntry.build_entry` class method to a dedicated function `legacy.core.fileio.dict_entry.build_entry`, and entry classes have been added.

- `legacy.aux_libs.gui.helpers.StreamFormerThread` architecture changes, so that it can accumulate several rows before adding them into the storage; this lead to replacement of `legacy.aux_libs.gui.helpers.StreamFormerThread.prepare_new_row` method by `legacy.aux_libs.gui.helpers.StreamFormerThread.prepare_new_data`.

## Additions

- General
  - Added pandas support in a bunch of places: loading/saving tables and dictionaries; data processing routines in `legacy.core.dataproc`; conversion of `legacy.core.dataproc.datatable.DataTable` and `legacy.core.utils.dictionary.Dictionary` object to/from pandas dataframes.
  - Expanded string conversion to support more explicit variable classes. For example, a numpy array `np.array([1,2,3])` can be converted into a string `'array([1, 2, 3])'` instead of a more ambiguous string `'[1, 2, 3]'` (which can also be a list). This behavior is controlled by a new argument `use_classes` in string conversion functions (such as `legacy.core.utils.string.to_string` and `legacy.core.utils.string.from_string`) and an argument `use_rep_classes` in file saving (`legacy.core.fileio.savefile.save`)
  - Added general library parameters, which can be accessed via `pylablib.par` (works as a dictionary object). So far there's only one supported parameter: the default return type of the CSV file reading (can be "pandas" for pandas dataframe, "table" for `legacy.core.dataproc.datatable.DataTable` object, or "array" for raw numpy array).
- Devices
  - Added `LaserQuantum Finesse` device class (`legacy.aux_libs.devices.devices.LaserQuantum`)
  - NI DAQ now supports output of waveforms
  - Added `legacy.aux_libs.devices.PCO_SC2.reset_api` and `legacy.aux_libs.devices.PCO_SC2.PCOSC2Camera.reboot` methods for resetting API and cameras
  - Added `legacy.aux_libs.devices.Thorlabs.list_kinesis_devices` function to list connected Kinesis devices
  - Added serial communication methods for IMAQ cameras (`legacy.aux_libs.devices.IMAQ.IMAQCamera`)
- GUI and threading
  - Added line plotter (`legacy.aux_libs.gui.widgets.line_plotter`) and trace plotter (`legacy.aux_libs.gui.widgets.trace_plotter`) widgets
  - Added virtual elements to value tables and parameter tables
  - Added `gui_thread_safe` parameter to value tables and parameter tables. Enabling it make most common methods thread-safe (i.e., transparently called from the GUI thread)
  - Added a corresponding `legacy.core.gui.qt.thread.controller.gui_thread_method` wrapper to implement the change above
  - Added functional thread variables (`legacy.core.gui.qt.thread.controller.QThreadController.set_func_variable`)
- File saving / loading
  - Added notation for dictionary files to include nested structures ('prefix blocks'). This lets one avoid common path prefix in stored dictionary files. For example, a file

```
some/long/prefix/x 1
some/long/prefix/y 2
some/long/prefix/y 3
```

can be represented as

```
//some/long/prefix
    x  1
    y  2
    z  3
///
```

The meaningful elements are `//some/long/prefix` line denoting that following elements have the given prefix, and `///` line denoting that the prefix block is done (indentation is only added for clarity).

- New dictionary entries: `dict_entry.ExternalNumpyDictionaryEntry` (external numpy array, can have arbitrary number of dimensions) and `dict_entry.ExpandedContainerDictionaryEntry` (turns lists, tuples and dicts into dictionary branches, so that their content can benefit from the automatic table inlining, dictionary entry classes, etc.).
- Data processing
  - `legacy.core.dataproc.fitting.Fitter` now takes default scale and limit as constructor arguments.
  - `legacy.core.dataproc.feature.multi_scale_peakdet` has new `norm_ratio` argument.
  - `legacy.core.dataproc.image.get_region` and `legacy.core.dataproc.image.get_region_sum` take axis argument.
- Miscellaneous
  - Functions introspection module now supports Python 3 - style functions, and added a new function `legacy.core.utils.functions.funcsig`
  - `legacy.core.utils.general.StreamFileLogger` supports multiple destination paths
  - New network function `legacy.core.utils.net.get_all_local_addr` (return list of all local addresses on all interfaces) and `legacy.core.utils.net.get_local_hostname`

## 1.7 pylablib

### 1.7.1 pylablib package

#### Subpackages

#### pylablib.core package

#### Subpackages

#### pylablib.core.dataproc package

#### Submodules

**pylablib.core.dataproc.callable module****class** pylablib.core.dataproc.callable.ICallableBases: `object`

Fit function generalization.

Has a set of mandatory argument with no default values and a set of parameters with default values (there may or may not be an explicit list of them).

All the arguments are passed explicitly by name. Passed value supersede default values. Extra arguments (not used in the calculations) are ignored.

Assumed (but not enforced) to be immutable: changes after creation can break the behavior.

Implements (possibly; depends on subclasses) call namelist binding boosting: if the function is to be called many times with the same parameter names list, one can first bind parameters list, and then call bound function with the corresponding arguments. This way, `callable(**p)` should be equivalent to `callable.bind(p.keys())(*p.values())`.

**has\_arg** (*arg\_name*)Determine if the function has an argument *arg\_name* (of all 3 categories)**filter\_args\_dict** (*args*)

Filter argument names dictionary to leave only the arguments that are used

**get\_mandatory\_args** ()

Return list of mandatory arguments (these are the ones without default values)

**is\_mandatory\_arg** (*arg\_name*)Check if the argument *arg\_name* is mandatory**get\_arg\_default** (*arg\_name*)Return default value of the argument *arg\_name*.

Raise `KeyError` if the argument is not defined or `ValueError` if it has no default value.

**bind** (*arg\_names*, *\*\*bound\_params*)Bind function to a given parameters set, leaving *arg\_names* as free parameters (in the given order)**class** NamesBoundCall (*func*, *names*, *bound\_params*)Bases: `object`**bind\_namelist** (*arg\_names*, *\*\*bound\_params*)

Bind namelist to boost subsequent calls.

Similar to `bind(arg_names)`, but bound function doesn't accept additional parameters and can be boosted.

**class** pylablib.core.dataproc.callable.MultiplexedCallable (*func*, *multiplex\_by*,  
*join\_method='stack'*)Bases: `pylablib.core.dataproc.callable.ICallable`

Multiplex a single callable based on a single parameter.

If the function is called with this parameter as an iterable, then the underlying callable will be called for each value of the parameter separately, and the results will be joined into a single array (if return the values are scalar, they're joined in 1D array; otherwise, they're joined using *join\_method*).

**Parameters**

- **func** (*callable*) – Function to be parallelized.
- **multiplex\_by** (*str*) – Name of the argument to be multiplexed by.

- **join\_method** (*str*) – Method for combining individual results together if they're non-scalars. Can be either 'list' (combine the results in a single list), 'stack' (combine using `numpy.column_stack()`, i.e., add dimension to the result), or 'concatenate' (concatenate the return values; the dimension of the result stays the same).

Multiplexing also makes use of call signatures for underlying function even if `__call__` is used.

Note that this operation is slow, and should be used only for high-dimensional multiplexing; for 1D case it's much better to just use numpy arrays as arguments and rely on numpy parallelizing.

**has\_arg** (*arg\_name*)

Determine if the function has an argument *arg\_name* (of all 3 categories)

**get\_mandatory\_args** ()

Return list of mandatory arguments (these are the ones without default values)

**get\_arg\_default** (*arg\_name*)

Return default value of the argument *arg\_name*.

Raise `KeyError` if the argument is not defined or `ValueError` if it has no default value.

**class NamesBoundCall** (*func, names, bound\_params*)

Bases: `object`

**bind** (*arg\_names, \*\*bound\_params*)

Bind function to a given parameters set, leaving *arg\_names* as free parameters (in the given order)

**bind\_namelist** (*arg\_names, \*\*bound\_params*)

Bind namelist to boost subsequent calls.

Similar to `bind(arg_names)`, but bound function doesn't accept additional parameters and can be boosted.

**filter\_args\_dict** (*args*)

Filter argument names dictionary to leave only the arguments that are used

**is\_mandatory\_arg** (*arg\_name*)

Check if the argument *arg\_name* is mandatory

**class** `pylablib.core.dataproc.callable.JoinedCallable` (*funcs, join\_method='stack'*)

Bases: `pylablib.core.dataproc.callable.ICallable`

Join several callables sharing the same arguments list.

The results will be joined into a single array (if return the values are scalar, they're joined in 1D array; otherwise, they're joined using *join\_method*).

#### Parameters

- **funcs** (*[callable]*) – List of functions to be joined together.
- **join\_method** (*str*) – Method for combining individual results together if they're non-scalars. Can be either 'list' (combine the results in a single list), 'stack' (combine using `numpy.column_stack()`, i.e., add dimension to the result), or 'concatenate' (concatenate the return values; the dimension of the result stays the same).

**has\_arg** (*arg\_name*)

Determine if the function has an argument *arg\_name* (of all 3 categories)

**get\_mandatory\_args** ()

Return list of mandatory arguments (these are the ones without default values)

**get\_arg\_default** (*arg\_name*)

Return default value of the argument *arg\_name*.

Raise `KeyError` if the argument is not defined or `ValueError` if it has no default value.

**class** `NamesBoundCall` (*func*, *names*, *bound\_params*)

Bases: `object`

**bind** (*arg\_names*, *\*\*bound\_params*)

Bind function to a given parameters set, leaving *arg\_names* as free parameters (in the given order)

**bind\_namelist** (*arg\_names*, *\*\*bound\_params*)

Bind namelist to boost subsequent calls.

Similar to `bind(arg_names)`, but bound function doesn't accept additional parameters and can be boosted.

**filter\_args\_dict** (*args*)

Filter argument names dictionary to leave only the arguments that are used

**is\_mandatory\_arg** (*arg\_name*)

Check if the argument *arg\_name* is mandatory

**class** `pylablib.core.dataproc.callable.FunctionCallable` (*func*, *function\_signature=None*, *defaults=None*, *alias=None*)

Bases: `pylablib.core.dataproc.callable.ICallable`

Callable based on a function or a method.

#### Parameters

- **func** – Function to be wrapped.
- **function\_signature** – A `functions.FunctionSignature` object supplying information about function's argument names and default values, if they're different from what's extracted from its signature.
- **defaults** (*dict*) – A dictionary {*name*: *value*} of additional default parameters values. Override the defaults from the signature. All default values must be pass-able to the function as a parameter
- **alias** (*dict*) – A dictionary {*alias*: *original*} for renaming some of the original arguments. Original argument names can't be used if aliased (though, multi-aliasing can be used explicitly, e.g., `alias={'alias':'arg', 'arg':'arg'}`). A name can be blocked (its usage causes error) if it's aliased to None (`alias={'blocked_name':None}`).

Optional non-named arguments in the form `*args` are not supported, since all the arguments are passed to the function by keywords.

Optional named arguments in the form `**kwargs` are supported only if their default values are explicitly provided in defaults (otherwise it would be unclear whether argument should be added into `**kwargs` or ignored altogether).

**has\_arg** (*arg\_name*)

Determine if the function has an argument *arg\_name* (of all 3 categories)

**get\_mandatory\_args** ()

Return list of mandatory arguments (these are the ones without default values)

**get\_arg\_default** (*arg\_name*)

Return default value of the argument *arg\_name*.

Raise `KeyError` if the argument is not defined or `ValueError` if it has no default value.

```
class NamesBoundCall (func, names, bound_params)
```

Bases: `object`

```
bind (arg_names, **bound_params)
```

Bind function to a given parameters set, leaving *arg\_names* as free parameters (in the given order)

```
bind_namelist (arg_names, **bound_params)
```

Bind namelist to boost subsequent calls.

Similar to `bind(arg_names)`, but bound function doesn't accept additional parameters and can be boosted.

```
filter_args_dict (args)
```

Filter argument names dictionary to leave only the arguments that are used

```
is_mandatory_arg (arg_name)
```

Check if the argument *arg\_name* is mandatory

```
class pylablib.core.dataproc.callable.MethodCallable (method, func-  
 tion_signature=None, de-  
 faults=None, alias=None)
```

Bases: `pylablib.core.dataproc.callable.FunctionCallable`

Similar to `FunctionCallable`, but accepts class method instead of a function.

The only addition is that now object's attributes can also parameters to the function: all the parameters which are not explicitly mentioned in the method signature are assumed to be object's attributes.

The parameters are affected by alias, but NOT affected by defaults (since it's impossible to ensure that all object's attributes are kept constant, and it's impractical to reset them all to default values at every function call).

#### Parameters

- **method** – Method to be wrapped.
- **function\_signature** – A `functions.FunctionSignature` object supplying information about function's argument names and default values, if they're different from what's extracted from its signature. If it's assumed that the first self argument is already excluded.
- **defaults** (*dict*) – A dictionary {name: value} of additional default parameters values. Override the defaults from the signature. All default values must be pass-able to the function as a parameter
- **alias** (*dict*) – A dictionary {alias: original} for renaming some of the original arguments. Original argument names can't be used if aliased (though, multi-aliasing can be used explicitly, e.g., `alias={'alias':'arg', 'arg':'arg'}`). A name can be blocked (its usage causes error) if it's aliased to None (`alias={'blocked_name':None}`).

This callable is implemented largely to be used with `TheoryCalculator` class (currently deprecated).

```
has_arg (arg_name)
```

Determine if the function has an argument *arg\_name* (of all 3 categories)

```
get_arg_default (arg_name)
```

Return default value of the argument *arg\_name*.

Raise `KeyError` if the argument is not defined or `ValueError` if it has no default value.

```
class NamesBoundCall (func, names, bound_params)
```

Bases: `object`

**bind** (*arg\_names*, *\*\*bound\_params*)

Bind function to a given parameters set, leaving *arg\_names* as free parameters (in the given order)

**bind\_namelist** (*arg\_names*, *\*\*bound\_params*)

Bind namelist to boost subsequent calls.

Similar to `bind(arg_names)`, but bound function doesn't accept additional parameters and can be boosted.

**filter\_args\_dict** (*args*)

Filter argument names dictionary to leave only the arguments that are used

**get\_mandatory\_args** ()

Return list of mandatory arguments (these are the ones without default values)

**is\_mandatory\_arg** (*arg\_name*)

Check if the argument *arg\_name* is mandatory

`pylablib.core.dataproc.callable.to_callable` (*func*)

Convert a function to an *ICallable* instance.

If it's already *ICallable*, return unchanged. Otherwise, return *FunctionCallable* or *MethodCallable* depending on whether it's a function or a bound method.

## pylablib.core.dataproc.feature module

Traces feature detection: peaks, baseline, local extrema.

**class** `pylablib.core.dataproc.feature.Baseline`

Bases: `pylablib.core.dataproc.feature.Baseline`

Baseline (background) for a trace.

*position* is the background level, and *width* is its noise width.

**count** ()

Return number of occurrences of value.

**index** ()

Return first index of value.

Raises `ValueError` if the value is not present.

**position**

**width**

`pylablib.core.dataproc.feature.get_baseline_simple` (*trace*, *find\_width=True*)

Get the baseline of the 1D trace.

If *find\_width==True*, calculate its width as well.

`pylablib.core.dataproc.feature.subtract_baseline` (*trace*)

Subtract baseline from the trace (make its background zero).

**class** `pylablib.core.dataproc.feature.Peak`

Bases: `pylablib.core.dataproc.feature.Peak`

A trace peak.

*kernel* defines its shape (for, e.g., generation purposes).

**count** ()

Return number of occurrences of value.



**height**

**index()**

Return first index of value.

Raises ValueError if the value is not present.

**kernel**

**position**

**width**

```
pylablib.core.dataproc.feature.find_peaks_cutoff(trace, cutoff, min_width=0,
                                                kind='peak', subtract_bl=True)
```

Find peaks in the data using cutoff.

#### Parameters

- **trace** – 1D data array.
- **cutoff** (*float*) – Cutoff value for the peak finding.
- **min\_width** (*int*) – Minimal uninterrupted width (in datapoints) of a peak. Any peaks this width are ignored.
- **kind** (*str*) – Peak kind. Can be 'peak' (positive direction), 'dip' (negative direction) or 'both' (both directions).
- **subtract\_bl** (*bool*) – If True, subtract baseline of the trace before checking cutoff.

**Returns** List of *Peak* objects.

```
pylablib.core.dataproc.feature.rescale_peak(peak, xoff=0.0, xscale=1.0, yoff=0,
                                           yscale=1.0)
```

Rescale peak's position, width and height.

*xscale* rescales position and width, *xoff* shifts position, *yscale* and *yoff* affect peak height.

```
pylablib.core.dataproc.feature.peaks_sum_func(peaks, peak_func='lorentzian')
```

Create a function representing sum of *peaks*.

*peak\_func* determines default peak kernel (used if `peak.kernel=="generic"`). Kernel is either a name string or a function taking 3 arguments (*x*, *width*, *height*).

```
pylablib.core.dataproc.feature.get_kernel(width, kernel_width=None, kernel='lorentzian')
```

Get a finite-sized kernel.

Return 1D array of length  $2 \times \text{kernel\_width} + 1$  containing the given kernel. By default, `kernel_width=int(width*3)`.

```
pylablib.core.dataproc.feature.get_peakdet_kernel(peak_width, background_width,
                                                  kernel_width=None, kernel='lorentzian')
```

Get a peak detection kernel.

Return 1D array of length  $2 \times \text{kernel\_width} + 1$  containing the kernel. The kernel is a sum of narrow positive peak (with the width *peak\_width*) and a broad negative peak (with the width *background\_width*); both widths are specified in datapoints (index). Each peak is normalized to have unit sum, i.e., the kernel has zero total sum. By default, `kernel_width=int(background_width*3)`.

```
pylablib.core.dataproc.feature.multi_scale_peakdet(trace, widths, background_ratio,
                                                    kind='peak', norm_ratio=None,
                                                    kernel='lorentzian')
```

Detect multiple peak widths using `get_peakdet_kernel()` kernel.

### Parameters

- **trace** – 1D data array.
- **widths** (*[float]*) – Array of possible peak widths.
- **background\_ratio** (*float*) – ratio of the *background\_width* to the *peak\_width* in `get_peakdet_kernel()`.
- **kind** (*str*) – Peak kind. Can be 'peak' (positive direction) or 'dip' (negative direction).
- **norm\_ratio** (*float*) – if not None, defines the width of the “normalization region” (in units of the kernel width, same as for the background kernel); it is then used to calculate a local trace variance to normalize the peaks magnitude.
- **kernel** – Peak matching kernel.

**Returns** Filtered trace which shows peak ‘affinity’ at each point.

```
pylablib.core.dataproc.feature.find_local_extrema(wf, region_width=3, kind='max',  
                                                  min_distance=None)
```

Find local extrema (minima or maxima) of 1D trace.

*kind* can be "min" or "max" and determines the kind of the extrema. Local minima (maxima) are defined as points which are smaller (greater) than all other points in the region of width *region\_width* around it. *region\_width* is always round up to an odd integer. *min\_distance* defines the minimal distance between the extrema (*region\_width*//2 by default). If there are several extrema within *min\_distance*, their positions are averaged together.

```
pylablib.core.dataproc.feature.latching_trigger(wf, threshold_on, thresh-  
                                                old_off, init_state='undef', re-  
                                                sult_kind='separate')
```

Determine indices of rise and fall trigger events with hysteresis (latching) thresholds.

Return either two arrays (*rise\_trig*, *fall\_trig*) containing trigger indices (if *result\_kind*=="separate"), or a single array of tuples [(*dir*,*pos*)], where *dir* is the trigger direction (+1 or -1) and *pos* is its index (if *result\_kind*=="joined"). Triggers happen when a state switch from 'high' to 'low' (rising) or vice versa (falling). The state switches from 'low' to 'high' when the trace value goes above *threshold\_on*, and from 'high' to 'low' when the trace value goes below *threshold\_off*. For a stable hysteresis effect, *threshold\_on* should be larger than *threshold\_off*, which means that the trace values between these two thresholds can not change the state. *init\_state* specifies the initial state: "low", "high", or "undef" (undefined state).

## pylablib.core.dataproc.filters module

Routines for filtering arrays (mostly 1D data).

```
pylablib.core.dataproc.filters.convolve1d(trace, kernel, mode='reflect', cval=0.0)
```

Convolution filter.

Convolves *trace* with the given *kernel* (1D array). *mode* and *cval* determine how the endpoints are handled. Simply a wrapper around the standard `scipy.ndimage.convolve1d()` that handles complex arguments.

```
pylablib.core.dataproc.filters.convolution_filter(a, width, kernel='gaussian', ker-  
                                                  nel_span='auto', mode='reflect',  
                                                  cval=0.0, kernel_height=None)
```

Convolution filter.

### Parameters

- **a** – array for filtering.

- **width** (*float*) – kernel width (second parameter to the kernel function).
- **kernel** – either a string defining the kernel function (see `specfunc.get_kernel_func()` for possible kernels), or a function taking 3 arguments (*pos*, *width*, *height*), where *height* can be `None` (assumes normalization by area).
- **kernel\_span** – the cutoff for the kernel function. Either an integer (number of points) or 'auto' (autodetect for "gaussian", "rectangle" and "exp\_decay", full trace width for all other kernels).
- **mode** (*str*) – convolution mode (see `scipy.ndimage.convolve()`).
- **cval** (*float*) – convolution fill value (see `scipy.ndimage.convolve()`).
- **kernel\_height** – height parameter to be passed to the kernel function. `None` means normalization by area.

`pylablib.core.dataproc.filters.gaussian_filter(a, width, mode='reflect', cval=0.0)`

Simple gaussian filter. Can handle complex data.

Equivalent to a convolution with a gaussian. Equivalent to `scipy.ndimage.gaussian_filter1d()`, uses `convolution_filter()`.

`pylablib.core.dataproc.filters.gaussian_filter_nd(a, width, mode='reflect', cval=0.0)`

Simple gaussian filter. Can't handle complex data.

Equivalent to a convolution with a gaussian. Wrapper around `scipy.ndimage.gaussian_filter()`.

`pylablib.core.dataproc.filters.low_pass_filter(trace, t, mode='reflect', cval=0.0)`

Simple single-pole low-pass filter.

*t* is the filter time constant, *mode* and *cval* are the trace expansion parameters (only from the left). Implemented as a recursive digital filter, so its performance doesn't depend strongly on *t*. Works only for 1D arrays.

`pylablib.core.dataproc.filters.high_pass_filter(trace, t, mode='reflect', cval=0.0)`

Simple single-pole high-pass filter (equivalent to subtracting a low-pass filter).

*t* is the filter time constant, *mode* and *cval* are the trace expansion parameters (only from the left). Implemented as a recursive digital filter, so its performance doesn't depend strongly on *t*. Works only for 1D arrays.

`pylablib.core.dataproc.filters.integrate(trace)`

Calculate the integral of the trace.

Alias for `numpy.cumsum()`.

`pylablib.core.dataproc.filters.differentiate(trace)`

Calculate the differential of the trace.

Note that since the data dimensions are changed (length is reduced by 1), the index is not preserved for pandas DataFrames.

`pylablib.core.dataproc.filters.sliding_average(a, width, mode='reflect', cval=0.0)`

Simple sliding average filter

Equivalent to convolution with a rectangle peak function.

`pylablib.core.dataproc.filters.median_filter(a, width, mode='reflect', cval=0.0)`

Median filter.

Wrapper around `scipy.ndimage.median_filter()`.

`pylablib.core.dataproc.filters.sliding_filter(trace, n, dec='bin', mode='reflect', cval=0.0)`

Perform sliding filtering on the data.

#### Parameters

- **trace** – 1D array-like object.
- **n** (*int*) – bin width.
- **dec** (*str*) – decimation method. Can be - 'bin' or 'mean' - do a binning average; - 'sum' - sum points; - 'min' - leave min point; - 'max' - leave max point; - 'median' - leave median point (works as a median filter). - a function which takes a single 1D array and compresses it into a number
- **mode** (*str*) – Expansion mode. Can be 'constant' (added values are determined by *cval*), 'nearest' (added values are end values of the trace), 'reflect' (reflect trace with respect to its endpoint) or 'wrap' (wrap the values from the other size).
- **cval** (*float*) – If mode=='constant', determines the expanded values.

`pylablib.core.dataproc.filters.decimate(a, n, dec='skip', axis=0, mode='drop')`  
Decimate the data.

Note that since the data dimensions are changed, the index is not preserved for pandas DataFrames.

#### Parameters

- **a** – data array.
- **n** (*int*) – decimation factor.
- **dec** (*str*) – decimation method. Can be - 'skip' - just leave every n'th point while completely omitting everything else; - 'bin' or 'mean' - do a binning average; - 'sum' - sum points; - 'min' - leave min point; - 'max' - leave max point; - 'median' - leave median point (works as a median filter). - a function which takes two arguments (nD numpy array and an axis) and compresses the array along the given axis
- **axis** (*int*) – axis along which to perform the decimation; can also be a tuple, in which case the same decimation is performed sequentially along several axes.
- **mode** (*str*) – determines what to do with the last bin if it's incomplete. Can be either 'drop' (omit the last bin) or 'leave' (keep it).

`pylablib.core.dataproc.filters.binning_average(a, width, axis=0, mode='drop')`  
Binning average filter.

Equivalent to `decimate()` with `dec=='bin'`.

`pylablib.core.dataproc.filters.decimate_full(a, dec='skip', axis=0)`  
Completely decimate the data along a given axis.

#### Parameters

- **a** – data array.
- **dec** (*str*) – decimation method. Can be - 'skip' - just leave every n'th point while completely omitting everything else; - 'bin' or 'mean' - do a binning average; - 'sum' - sum points; - 'min' - leave min point; - 'max' - leave max point; - 'median' - leave median point (works as a median filter). - a function which takes two arguments (nD numpy array and an axis) and compresses the array along the given axis
- **axis** (*int*) – axis along which to perform the decimation; can also be a tuple, in which case the same decimation is performed along several axes.

`pylablib.core.dataproc.filters.decimate_datasets(arrs, dec='mean')`  
Decimate datasets with the same shape element-wise (works only for 1D or 2D arrays).

Note that the index data is taken from the first array in the list.

*dec* has the same values and meaning as in `decimate()`. The format of the output (numpy or pandas, and the name of columns in pandas DataFrame) is determined by the first array in the list.

```
pylablib.core.dataproc.filters.collect_into_bins(values, distance, pre-
                                                serve_order=False,
                                                to_return='value')
```

Collect all values into bins separated at least by *distance*.

Return the extent of each bin. If `preserve_order==False`, values are sorted before splitting. If `to_return="value"`, the extent is given in values; if `to_return="index"`, it is given in indices (only useful if `preserve_order=True`, as otherwise the indices correspond to a sorted array). If *distance* is a tuple, then it denotes the minimal and the maximal separation between consecutive elements; otherwise, it is a single number denoting maximal absolute distance (i.e., it corresponds to a tuple `(-distance, distance)`).

```
pylablib.core.dataproc.filters.split_into_bins(values, max_span, max_size=None)
Split values into bins of the span at most max_span and number of elements at most max_size.
```

If *max\_size* is None, it's assumed to be infinite. Return array of indices for each bin. Values are sorted before splitting.

```
pylablib.core.dataproc.filters.fourier_filter(trace, response, dt=1, pre-
                                              serve_real=True)
```

Apply filter to a trace in the frequency domain.

*response* is a (possibly) complex function with single 1D real numpy array as a frequency argument. *dt* specifies time step between consecutive points. Note that in case of a multi-column data the filter is applied column-wise; this is in contrast with the Fourier transform methods, which would assume the first column to be times.

If `preserve_real==True`, then the *response* for negative frequencies is automatically taken to be complex conjugate of the *response* for positive frequencies (so that the real trace stays real).

```
pylablib.core.dataproc.filters.fourier_make_response_real(response)
Turn a frequency filter function into a real one (in the time domain).
```

Done by reflecting and complex conjugating positive frequency part to negative frequencies. *response* is a function with a single argument (frequency), return value is a modified function.

```
pylablib.core.dataproc.filters.fourier_filter_bandpass(pass_range_min,
                                                       pass_range_max)
```

Generate a bandpass filter function (hard cutoff).

The function is symmetric, so that it corresponds to a real response in time domain.

```
pylablib.core.dataproc.filters.fourier_filter_bandstop(stop_range_min,
                                                       stop_range_max)
```

Generate a bandstop filter function (hard cutoff).

The function is symmetric, so that it corresponds to a real response in time domain.

```
class pylablib.core.dataproc.filters.RunningDecimationFilter(n, mode='mean',
                                                            on_incomplete='none')
```

Bases: `object`

Running decimation filter.

Remembers last *n* samples and returns their averages, median, etc.

#### Parameters

- **n** – decimation length
- **mode** – decimation mode ("mean", "median", "min", or "max")

- **on\_incomplete** – determines what to return while the filter window is not yet full; can be "none" (default, return None), or "partial" (operate on the partial accumulated data)

**get** ()  
Get the filtered result

**add** (x)  
Add a new sample

**reset** ()  
Reset the filter

**class** pylablib.core.dataproc.filters.**RunningDebounceFilter** (n, precision=None, initial=None)

Bases: `object`

Running debounce filter.

“Sticks” to the current value and only switches when a new value remains constant (withing a given precision) for a given number of samples. Filters out temporary spikes and short changes, conceptually similar to a running median filter.

#### Parameters

- **n** – length of the required constant period
- **precision** – comparison precision (None means that the values should be exactly equal)
- **initial** – initial value; None means that the first sample sets this value

**get** ()  
Get the filtered result

**add** (x)  
Add a new sample

**reset** ()  
Reset the filter

## pylablib.core.dataproc.fitting module

Universal function fitting interface.

**class** pylablib.core.dataproc.fitting.**Fitter** (func, xarg\_name=None, fit\_parameters=None, fixed\_parameters=None, scale=None, limits=None, weights=None)

Bases: `object`

Fitter object.

Can handle variety of different functions, complex arguments or return values, array arguments.

#### Parameters

- **func** (*callable*) – Fit function. Can be anything callable (function, method, object with `__call__` method, etc.).
- **xarg\_name** (*str or list*) – Name (or multiple names) for x arguments. These arguments are passed to *func* (as named arguments) when calling for fitting. Can be a string (single argument) or a list (arbitrary number of arguments, including zero).

- **fit\_parameters** (*dict*) – Dictionary {name: value} of parameters to be fitted (*value* is the starting value for the fitting procedure). If *value* is None, try and get the default value from the *func*.
- **fixed\_parameters** (*dict*) – Dictionary {name: value} of parameters to be fixed during the fitting procedure. If *value* is None, try and get the default value from the *func*.
- **scale** (*dict*) – Defines typical scale of fit parameters (used to normalize fit parameters supplied of `scipy.optimize.least_squares()`). Note: for complex parameters scale must also be a complex number, with re and im parts of the scale variable corresponding to the scale of the re and im part.
- **limits** (*dict*) – Boundaries for the fit parameters (missing entries are assumed to be unbound). Each boundary parameter is a tuple (*lower*, *upper*). *lower* or *upper* can be None, `numpy.nan` or `numpy.inf` (with the appropriate sign), which implies no bounds in the given direction. Note: for compound data types (such as lists) the entries are still tuples of 2 elements, each of which is either None (no bound for any sub-element) or has the same structure as the full parameter. Note: for complex parameters limits must also be complex numbers (or None), with re and im parts of the limits variable corresponding to the limits of the re and im part.
- **weights** (*list or numpy.ndarray*) – Determines the weights of y-points. Can be either an array broadcastable to *y* (e.g., a scalar or an array with the same shape as *y*), in which case it's interpreted as list of individual point weights (which multiply residuals before they are squared). Or it can be an array with number of elements which is square of the number of elements in *y*, in which case it's interpreted as a weights matrix (which matrix-multiplies residuals before they are squared).

**set\_xarg\_name** (*xarg\_name*)

Set names of x arguments.

Can be a string (single argument) or a list (arbitrary number of arguments, including zero).

**use\_xarg** ()

Return True if the function requires x arguments

**set\_fixed\_parameters** (*fixed\_parameters*)

Change fixed parameters

**update\_fixed\_parameters** (*fixed\_parameters*)

Update the dictionary of fixed parameters

**del\_fixed\_parameters** (*fixed\_parameters*)

Remove fixed parameters

**set\_fit\_parameters** (*fit\_parameters*)

Change fit parameters

**update\_fit\_parameters** (*fit\_parameters*)

Update the dictionary of fit parameters

**del\_fit\_parameters** (*fit\_parameters*)

Remove fit parameters

**fit** (*x=None, y=0, fit\_parameters=None, fixed\_parameters=None, scale='default', limits='default', weights=1.0, parscore=None, return\_stderr=False, return\_residual=False, \*\*kwargs*)  
Fit the data.

#### Parameters

- **x** – x arguments. If the function has single x argument, *x* is an array-like object; otherwise, *x* is a list of array-like objects (can be None if there are no x parameters).



- **y** – Target function values.
- **fit\_parameters** (*dict*) – Adds to the default *fit\_parameters* of the fitter (has priority on duplicate entries).
- **fixed\_parameters** (*dict*) – Adds to the default *fixed\_parameters* of the fitter (has priority on duplicate entries).
- **scale** (*dict*) – Defines typical scale of fit parameters (used to normalize fit parameters supplied of `scipy.optimize.least_squares()`). Note: for complex parameters scale must also be a complex number, with re and im parts of the scale variable corresponding to the scale of the re and im part. If value is "default", use the value supplied on the fitter creation (by default, no specific scales).
- **limits** (*dict*) – Boundaries for the fit parameters (missing entries are assumed to be unbound). Each boundary parameter is a tuple (*lower*, *upper*). *lower* or *upper* can be `None`, `numpy.nan` or `numpy.inf` (with the appropriate sign), which implies no bounds in the given direction. Note: for compound data types (such as lists) the entries are still tuples of 2 elements, each of which is either `None` (no bound for any sub-element) or has the same structure as the full parameter. Note: for complex parameters limits must also be complex numbers (or `None`), with re and im parts of the limits variable corresponding to the limits of the re and im part. If value is "default", use the value supplied on the fitter creation (by default, no limits).
- **weights** (*list or numpy.ndarray*) – Determines the weights of y-points. Can be either an array broadcastable to *y* (e.g., a scalar or an array with the same shape as *y*), in which case it's interpreted as list of individual point weights (which multiply residuals before they are squared). Or it can be an array with number of elements which is square of the number of elements in *y*, in which case it's interpreted as a weights matrix (which matrix-multiplies residuals before they are squared). If value is "default", use the value supplied on the fitter creation (by default, no weights)
- **parscore** (*callable*) – parameter score function, whose value is added to the mean-square error (sum of all residuals squared) after applying weights. Takes the same parameters as the fit function, only without the x-arguments, and return an array-like value. Can be used for, e.g., 'soft' fit parameter constraining.
- **return\_stderr** (*bool*) – If `True`, append *stderr* to the output.
- **return\_residual** – If not `False`, append *residual* to the output.
- **\*\*kwargs** – arguments passed to `scipy.optimize.least_squares()` function.

## Returns

(*params*, *bound\_func*[, *stderr*][, *residual*]):

- *params*: a dictionary {*name*: *value*} of the parameters supplied to the function (both fit and fixed).
- *bound\_func*: the fit function with all the parameters bound (i.e., it only requires x parameters).
- *stderr*: a dictionary {*name*: *error*} of standard deviation for fit parameters to the return parameter. If the fitting routine returns no residuals (usually for a bad or an under-constrained fit), all residuals are set to NaN.
- *residual*: either a full array of residuals `func(x, **params) - y` (if `return_residual == 'full'`), a mean magnitude of the residuals `mean(abs(func(x, **params) - y) ** 2)` (if `return_residual == True` or `return_residual == 'mean'`),



or the total residuals including weights `mean(abs((func(x, **params)-y)*weights)**2)` (if `return_residual=='weighted'`).

#### Return type `tuple`

**initial\_guess** (*fit\_parameters=None*, *fixed\_parameters=None*, *return\_stderr=False*, *return\_residual=False*)

Return the initial guess for the fitting.

#### Parameters

- **fit\_parameters** (*dict*) – Overrides the default *fit\_parameters* of the fitter.
- **fixed\_parameters** (*dict*) – Overrides the default *fixed\_parameters* of the fitter.
- **return\_stderr** (*bool*) – If `True`, append *stderr* to the output.
- **return\_residual** – If not `False`, append *residual* to the output.

#### Returns

(*params*, *bound\_func*).

- *params*: a dictionary {*name*: *value*} of the parameters supplied to the function (both fit and fixed).
- *bound\_func*: the fit function with all the parameters bound (i.e., it only requires *x* parameters).
- *stderr*: a dictionary {*name*: *error*} of standard deviation for fit parameters to the return parameter. Always zero, added for better compatibility with `fit()`.
- *residual*: either a full array of residuals `func(x, **params)-y` (if `return_residual=='full'`) a mean magnitude of the residuals `mean(abs(func(x, **params)-y)**2)` (if `return_residual==True` or `return_residual=='mean'`). Always zero, added for better compatibility with `fit()`.

#### Return type `tuple`

`pylablib.core.dataproc.fitting.huge_error(x, factor=100.0)`

`pylablib.core.dataproc.fitting.get_best_fit(x, y, fits)`

Select the best (lowest residual) fit result.

*x* and *y* are the argument and the value of the bound fit function. *fits* is the list of fit results (tuples returned by `Fitter.fit()`).

## pylablib.core.dataproc.fourier module

Routines for Fourier transform.

`pylablib.core.dataproc.fourier.get_prev_len(l, maxprime=7)`

Get the largest number less or equal to *l*, which is composed of prime factors up to *maxprime*.

So far, only *maxprime* of 2, 3, 5, 7 and 11 are supported. *maxprime* of 5 gives less than 15% length reduction (and less than 6% for lengths above 400). *maxprime* of 11 gives less than 8% length reduction (and less than 4% for lengths above 300).

`pylablib.core.dataproc.fourier.truncate_trace(trace, maxprime=7)`

Truncate trace length to the nearest smaller length which is composed of prime factors up to *maxprime*.

So far, only *maxprime* of 2, 3, 5, 7 and 11 are supported. *maxprime* of 5 gives less than 15% length reduction (and less than 6% for lengths above 400). *maxprime* of 11 gives less than 8% length reduction (and less than 4% for lengths above 300).

```
pylablib.core.dataproc.fourier.normalize_fourier_transform(ft,          normaliza-  
                                                           tion='none',  
                                                           df=None,  
                                                           copy=False)
```

Normalize the Fourier transform data.

*ft* is a 1D trace or a 2D array with 2 columns: frequency and complex amplitude. *normalization* can be 'none' (standard numpy normalization), 'sum' (the power sum is preserved:  $\text{sum}(\text{abs}(ft)**2) == \text{sum}(\text{abs}(\text{trace})**2)$ ), 'rms' (the power sum is equal to the trace RMS power:  $\text{sum}(\text{abs}(ft)**2) == \text{mean}(\text{abs}(\text{trace})**2)$ ), 'density' (power spectral density normalization,  $\text{sum}(\text{abs}(ft[:,1])**2)*df == \text{mean}(\text{abs}(\text{trace}[:,1])**2)$ ), or 'dBc' (same as 'density', but normalized by the mean of the trace) If *normalization*=='density', then *df* can specify the frequency step between two consecutive bins; if *df* is None, it is extracted from the first two points of the frequency axis (or set to 1, if *ft* is a 1D trace)

```
pylablib.core.dataproc.fourier.apply_window(trace_values, window='rectangle', win-  
                                           dow_power_compensate=True)
```

Apply FT window to the trace.

If *window\_power\_compensate*==True, multiply the data is multiplied by a compensating factor to preserve power in the spectrum.

```
pylablib.core.dataproc.fourier.fourier_transform(trace,          dt=None,          trun-  
                                                cate=False,          normaliza-  
                                                tion='none',          single_sided=False,  
                                                window='rectangle',          win-  
                                                dow_power_compensate=True,  
                                                raw=False)
```

Calculate a fourier transform of the trace.

#### Parameters

- **trace** – Time trace to be transformed. It can be a 1D trace of values, a 2-column trace, or a 3-column trace. If *dt* is None, then the first column is assumed to be time (only support uniform time step), and the other columns are either the trace values (for a single data column) or real and imaginary parts of the trace (for two data columns). If *dt* is not None, then the time column is assumed to be missing, so the two columns are assumed to be the real and the imaginary parts.
- **dt** – if not None, can specify the time step between the consecutive samples, in which case it is assumed that the time column is missing from the trace; otherwise, try to get it from the time column of the trace if it exists, or set to 1 otherwise.
- **truncate** (*bool* or *int*) – Determines whether to truncate the trace to the nearest product of small primes (speeds up FFT algorithm); can be False (no truncation), an integer 2, 3, 5, 7, or 11 (truncate to a product of primes up to and including this number), or True (default prime factorization, currently set to 7)
- **normalization** (*str*) – Fourier transform normalization: - 'none': no (i.e., default numpy) normalization; - 'sum': the norm of the data is conserved ( $\text{sum}(\text{abs}(ft[:,1])**2) == \text{sum}(\text{abs}(\text{trace}[:,1])**2)$ ); - 'rms': sum of the PSD is equal to the RMS trace amplitude squared ( $\text{sum}(\text{abs}(ft[:,1])**2) == \text{mean}(\text{abs}(\text{trace}[:,1])**2)$ ); - 'density': power spectral density normalization, in x/rtHz ( $\text{sum}(\text{abs}(ft[:,1])**2)*df == \text{mean}(\text{abs}(\text{trace}[:,1])**2)$ ); - 'dBc': like 'density', but normalized to the mean trace value.
- **single\_sided** (*bool*) – If True, only leave positive frequency side of the transform.

- **window** (*str*) – FT window. Can be 'rectangle' (essentially, no window), 'hann' or 'hamming'.
- **window\_power\_compensate** (*bool*) – If True, the data is multiplied by a compensating factor to preserve power in the spectrum.
- **raw** (*bool*) – if True, return a simple 1D trace with the result.

**Returns** a two-column array of the same kind as the input, where the first column is frequency, and the second is complex FT data.

```
pylablib.core.dataproc.fourier.flip_fourier_transform(ft)
```

Flip the fourier transform (analogous to making frequencies negative and flipping the order).

```
pylablib.core.dataproc.fourier.inverse_fourier_transform(ft, df=None,
                                                         truncate=False,
                                                         zero_loc=None, symmetric_time=False,
                                                         raw=False)
```

Calculate an inverse fourier transform of the trace.

#### Parameters

- **ft** – Fourier transform data to be inverted. It can be a 1D trace of values, a 2-column trace, or a 3-column trace. If *df* is None, then the first column is assumed to be frequency (only support uniform frequency step), and the other columns are either the trace values (for a single data column) or real and imaginary parts of the trace (for two data columns). If *df* is not None, then the frequency column is assumed to be missing, so the two columns are assumed to be the real and the imaginary parts.
- **df** – if not None, can specify the frequency step between the consecutive samples; otherwise, try to get it from the frequency column of the trace if it exists, or set to 1 otherwise.
- **truncate** (*bool or int*) – Determines whether to truncate the trace to the nearest product of small primes (speeds up FFT algorithm); can be False (no truncation), an integer 2, 3, 5, 7, or 11 (truncate to a product of primes up to and including this number), or True (default prime factorization, currently set to 7)
- **zero\_loc** (*bool*) – Location of the zero frequency point. Can be None (the one with the value of f-axis closest to zero, or the first point if the frequency column is missing), 'center' (mid-point), or an integer index.
- **symmetric\_time** (*bool*) – If True, make time axis go from  $(-0.5/df, 0.5/df)$  rather than  $(0, 1./df)$ .
- **raw** (*bool*) – if True, return a simple 1D trace with the result.

**Returns** a two-column array, where the first column is frequency, and the second is the complex-valued trace data.

```
pylablib.core.dataproc.fourier.power_spectral_density(trace, dt=None, truncate=False,
                                                         normalization='density', single_sided=False,
                                                         window='rectangle', window_power_compensate=True,
                                                         raw=False)
```

Calculate a power spectral density of the trace.

#### Parameters

- **trace** – Time trace to be transformed. It can be a 1D trace of values, a 2-column trace, or a 3-column trace. If *dt* is `None`, then the first column is assumed to be time (only support uniform time step), and the other columns are either the trace values (for a single data column) or real and imaginary parts of the trace (for two data columns). If *dt* is not `None`, then the time column is assumed to be missing, so the two columns are assumed to be the real and the imaginary parts.
- **dt** – if not `None`, can specify the time step between the consecutive samples; otherwise, try to get it from the time column of the trace if it exists, or set to 1 otherwise.
- **truncate** (*bool or int*) – Determines whether to truncate the trace to the nearest product of small primes (speeds up FFT algorithm); can be `False` (no truncation), an integer 2, 3, 5, 7, or 11 (truncate to a product of primes up to and including this number), or `True` (default prime factorization, currently set to 7)
- **normalization** (*str*) – Fourier transform normalization: - `'none'`: no (i.e., default numpy) normalization; - `'sum'`: the norm of the data is conserved (`sum(PSD[:,1])==sum(abs(trace[:,1])**2)`); - `'rms'`: sum of the PSD is equal to the RMS trace amplitude squared (`sum(PSD[:,1])==mean(abs(trace[:,1])**2)`); - `'density'`: power spectral density normalization, in `x/rHz` (`sum(PSD[:,1])*df==mean(abs(trace[:,1])**2)`); - `'dBc'`: like `'density'`, but normalized to the mean trace value.
- **single\_sided** (*bool*) – If `True`, only leave positive frequency side of the PSD.
- **window** (*str*) – FT window. Can be `'rectangle'` (essentially, no window), `'hann'` or `'hamming'`.
- **window\_power\_compensate** (*bool*) – If `True`, the data is multiplied by a compensating factor to preserve power in the spectrum.
- **raw** (*bool*) – if `True`, return a simple 1D trace with the result.

**Returns** a two-column array, where the first column is frequency, and the second is positive PSD.

`pylablib.core.dataproc.fourier.get_real_part_ft(ft)`

Get the fourier transform of the real part only from the fourier transform of a complex variable.

`pylablib.core.dataproc.fourier.get_imag_part_ft(ft)`

Get the fourier transform of the imaginary part only from the fourier transform of a complex variable.

`pylablib.core.dataproc.fourier.get_correlations_ft(ft_a, ft_b, zero_mean=True, normalization='none')`

Calculate the correlation function of the two variables given their fourier transforms.

#### Parameters

- **ft\_a** – first variable fourier transform
- **ft\_b** – second variable fourier transform
- **zero\_mean** (*bool*) – If `True`, the value corresponding to the zero frequency is set to zero (only fluctuations around means of *a* and *b* are calculated).
- **normalization** (*str*) – Can be `'whole'` (correlations are normalized by product of PSDs derived from *ft\_a* and *ft\_b*) or `'individual'` (normalization is done for each frequency individually, so that the absolute value is always 1).

## pylablib.core.dataproc.iir\_transform module

Digital recursive infinite impulse response filter.

Implemented using Numba library (JIT high-performance compilation) if possible.

```
pylablib.core.dataproc.iir_transform.iir_apply_complex(trace, xcoeff, ycoeff)
```

Apply digital, (possibly) recursive filter with coefficients *xcoeff* and *ycoeff* along the first axis.

Result is filtered signal *y* with  $y[n] = \sum_j x[n-j] * xcoeff[j] + \sum_k y[n-k-1] * ycoeff[k]$ .

## pylablib.core.dataproc.image module

```
pylablib.core.dataproc.image.convert_shape_indexing(shape, src, dst, axes=(0, 1))
```

Convert image indexing style.

*shape* is the source image shape (2-tuple), *src* and *dst* are current format and desired format. Formats can be "rcb" (first index is row, second is column, rows count from the bottom), "rct" (same, but rows count from the top), "xyb" (first index is column, second is row, rows count from the bottom), or "xyt" (same but rows count from the top). "rc" is interpreted as "rct", "xy" as "xyt"

```
pylablib.core.dataproc.image.convert_image_indexing(img, src, dst, axes=(0, 1))
```

Convert image indexing style.

*img* is the source image (ND numpy array with  $N \geq 2$ ), *src* and *dst* are current format and desired format, *axes* specify correspondingly the row and the column axes (by default, the first two array axes). Formats can be "rcb" (first index is row, second is column, rows count from the bottom), "rct" (same, but rows count from the top), "xyb" (first index is column, second is row, rows count from the bottom), or "xyt" (same but rows count from the top). "rc" is interpreted as "rct", "xy" as "xyt"

```
class pylablib.core.dataproc.image.ROI (imin=0, imax=None, jmin=0, jmax=None)
```

Bases: `object`

**copy** ()

**center** (*shape=None*)

**size** (*shape=None*)

**area** (*shape=None*)

**tup** (*shape=None*)

**ispan** (*shape=None*)

**jspan** (*shape=None*)

**classmethod from\_centersize** (*center, size, shape=None*)

**classmethod intersect** (\*args)

**limit** (*shape*)

```
pylablib.core.dataproc.image.get_region(image, center, size, axis=(-2, -1))
```

Get part of the image with the given center and size (both are tuples (*i*, *j*)).

The region is automatically reduced if a part of it is outside of the image.

```
pylablib.core.dataproc.image.get_region_sum(image, center, size, axis=(-2, -1))
```

Sum part of the image with the given center and size (both are tuples (*i*, *j*)).

The region is automatically reduced if a part of it is outside of the image. Return tuple (*sum*, *area*), where *area* is the actual summer region are (in pixels).

## pylablib.core.dataproc.interpolate module

```
pylablib.core.dataproc.interpolate.interpolate1D_func(x, y, kind='linear',
axis=-1, copy=True,
bounds_error=True,
fill_values=nan, assume_sorted=False)
```

1D interpolation.

Simply a wrapper around `scipy.interpolate.interp1d`.

### Parameters

- **x** – 1D arrays of x coordinates for the points at which to find the values.
- **y** – array of values corresponding to x points (can have more than 1 dimension, in which case the output values are (N-1)-dimensional)
- **kind** – Interpolation method.
- **axis** – axis in y-data over which to interpolate.
- **copy** – if True, make internal copies of x and y.
- **bounds\_error** – if True, raise error if interpolation function arguments are outside of x bounds.
- **fill\_values** – values to fill the outside-bounds regions if `bounds_error==False`.
- **assume\_sorted** – if True, assume that *data* is sorted.

**Returns** A 1D array with interpolated data.

```
pylablib.core.dataproc.interpolate.interpolate1D(data, x, kind='linear',
bounds_error=True, fill_values=nan,
assume_sorted=False)
```

1D interpolation.

### Parameters

- **data** – 2-column array [(x,y)], where y is a function of x.
- **x** – Arrays of x coordinates for the points at which to find the values.
- **kind** – Interpolation method.
- **bounds\_error** – if True, raise error if x values are outside of *data* bounds.
- **fill\_values** – values to fill the outside-bounds regions if `bounds_error==False`
- **assume\_sorted** – if True, assume that *data* is sorted.

**Returns** A 1D array with interpolated data.

```
pylablib.core.dataproc.interpolate.interpolate2D(data, x, y, method='linear',
fill_value=nan)
```

Interpolate data in 2D.

Simply a wrapper around `scipy.interpolate.griddata()`.

### Parameters

- **data** – 3-column array [(x,y,z)], where z is a function of x and y.
- **x/y** – Arrays of x and y coordinates for the points at which to find the values.

- **method** – Interpolation method.

**Returns** A 2D array with interpolated data.

```
pylablib.core.dataproc.interpolate.interpolateND(data, xs, method='linear')
```

Interpolate data in N dimensions.

Simply a wrapper around `scipy.interpolate.griddata()`.

#### Parameters

- **data** – (N+1)-column array [ (*x*<sub>1</sub>, ..., *x*<sub>N</sub>, *y*) ], where *y* is a function of *x*<sub>1</sub>, ..., *x*<sub>N</sub>.
- **xs** – N-tuple of arrays of coordinates for the points at which to find the values.
- **method** – Interpolation method.

**Returns** An ND array with interpolated data.

```
pylablib.core.dataproc.interpolate.regular_grid_from_scatter(data, x_points,  
                                                           y_points,  
                                                           x_range=None,  
                                                           y_range=None,  
                                                           method='nearest')
```

Turn irregular scatter-points data into a regular 2D grid function.

#### Parameters

- **data** – 3-column array [ (*x*, *y*, *z*) ], where *z* is a function of *x* and *y*.
- **x\_points/y\_points** – Number of points along *x/y* axes.
- **x\_range/y\_range** – If not None, a tuple specifying the desired range of the data (all points in *data* outside the range are excluded).
- **method** – Interpolation method (see `scipy.interpolate.griddata()` for options).

**Returns** A nested tuple (*data*, (*x\_grid*, *y\_grid*)), where all entries are 2D arrays (either with data or with gridpoint locations).

```
pylablib.core.dataproc.interpolate.interpolate_trace(trace, step, rng=None,  
                                                    x_column=0, select=  
                                                    columns=None,  
                                                    kind='linear', assume=  
                                                    sorted=False)
```

Interpolate trace data over a regular grid with the given step.

*rng* specifies interpolation range (by default, whole data range). *x\_column* specifies column index for *x*-data. *select\_column* specifies which columns to interpolate and keep at the output (by default, all data). If *assume\_sorted*==True, assume that *x*-data is sorted. *kind* specifies interpolation method.

```
pylablib.core.dataproc.interpolate.average_interpolate_1D(data, step, rng=None,  
                                                         avg_kernel=1,  
                                                         min_weight=0,  
                                                         kind='linear')
```

1D interpolation combined with pre-averaging.

#### Parameters

- **data** – 2-column array [(*x*,*y*)], where *y* is a function of *x*.
- **step** – distance between the points in the interpolated data (all resulting *x*-coordinates are multiples of *step*).

- **rng** – if not `None`, specifies interpolation range (by default, whole data range).
- **avg\_kernel** – kernel used for initial averaging. Can be either a 1D array, where each point corresponds to the relative bin weight, or an integer, which specifies simple rectangular kernel of the given width.
- **min\_weight** – minimal accumulated weight in the bin to consider it ‘valid’ (if the bin is invalid, its accumulated value is ignored, and its value is obtained by the interpolation step). *min\_weight* of 0 implies any non-zero weight; otherwise, *weight*  $\geq$  *min\_weight*.
- **kind** – Interpolation method.

**Returns** A 2-column array with the interpolated data.

## pylablib.core.dataproc.specfunc module

Specific useful functions.

`pylablib.core.dataproc.specfunc.gaussian_k(x, sigma=1.0, height=None)`  
Gaussian kernel function.

Normalized by the area if *height* is `None`, otherwise *height* is the value at 0.

`pylablib.core.dataproc.specfunc.rectangle_k(x, width=1.0, height=None)`  
” Symmetric rectangle kernel function.

Normalized by the area if *height* is `None`, otherwise *height* is the value at 0.

`pylablib.core.dataproc.specfunc.lorentzian_k(x, gamma=1.0, height=None)`  
Lorentzian kernel function

Normalized by the area if *height* is `None`, otherwise *height* is the value at 0.

`pylablib.core.dataproc.specfunc.complex_lorentzian_k(x, gamma=1.0, amplitude=1j)`  
Complex Lorentzian kernel function.

`pylablib.core.dataproc.specfunc.exp_decay_k(x, width=1.0, height=None, mode='causal')`  
Exponential decay kernel function

Normalized by area if *height*=`None` (if possible), otherwise *height* is the value at 0.

**Mode determines value for  $x < 0$ :**

- 'causal' - it's 0 for  $x < 0$ ;
- 'step' - it's constant for  $x \leq 0$ ;
- 'continue' - it's a continuous decaying exponent;
- 'mirror' - function is symmetric:  $\exp(-|x|/\text{width})$ .

`pylablib.core.dataproc.specfunc.get_kernel_func(kernel)`  
Get a kernel function by its name.

Available functions are: 'gaussian', 'rectangle', 'lorentzian', 'exp\_decay', 'complex\_lorentzian'.

`pylablib.core.dataproc.specfunc.rectangle_w(x, N, ft_compensated=False)`  
Rectangle FT window function

`pylablib.core.dataproc.specfunc.gen_hamming_w(x, N, alpha, beta, ft_compensated=False)`  
Generalized Hamming FT window function.



If `ft_compensated==True`, multiply the window function by a compensating factor to preserve power in the spectrum.

```
pylablib.core.dataproc.specfunc.hann_w(x, N, ft_compensated=False)
```

Hann FT window function.

If `ft_compensated==True`, multiply the window function by a compensating factor to preserve power in the spectrum.

```
pylablib.core.dataproc.specfunc.hamming_w(x, N, ft_compensated=False)
```

Specific Hamming FT window function.

If `ft_compensated==True`, multiply the window function by a compensating factor to preserve power in the spectrum.

```
pylablib.core.dataproc.specfunc.get_window_func(window)
```

Get a window function by its name.

Available functions are: 'hamming', 'rectangle', 'hann'.

```
pylablib.core.dataproc.specfunc.gen_hamming_w_ft(f, t, alpha, beta)
```

Get Fourier Transform of a generalized Hamming FT window function.

*f* is the argument, *t* is the total window size.

```
pylablib.core.dataproc.specfunc.rectangle_w_ft(f, t)
```

Get Fourier Transform of the rectangle FT window function.

*f* is the argument, *t* is the total window size.

```
pylablib.core.dataproc.specfunc.hann_w_ft(f, t)
```

Get Fourier Transform of the Hann FT window function.

*f* is the argument, *t* is the total window size.

```
pylablib.core.dataproc.specfunc.hamming_w_ft(f, t)
```

Get Fourier Transform of the specific Hamming FT window function.

*f* is the argument, *t* is the total window size.

```
pylablib.core.dataproc.specfunc.get_window_ft_func(window)
```

Get a Fourier Transform of a window function by its name.

Available functions are: 'hamming', 'rectangle', 'hann'.

## pylablib.core.dataproc.table\_wrap module

Utilities for uniform treatment of pandas tables and numpy arrays for functions which can deal with them both.

```
class pylablib.core.dataproc.table_wrap.IGenWrapper(container)
```

Bases: `object`

The interface for a wrapper that gives a uniform access to basic methods of wrapped objects'.

```
get_type()
```

Get a string representing the wrapped object type

```
copy(wrapped=False)
```

Copy the object.

If `wrapped==True`, return a new wrapper containing the object copy; otherwise, just return the copy.

```
ndim()
```

**shape()**

**class** pylablib.core.dataproc.table\_wrap.**1DWrapper**(*container*)

Bases: [pylablib.core.dataproc.table\\_wrap.IGenWrapper](#)

A wrapper containing a 1D object (a 1D numpy array or a pandas Series object).

Provides a uniform access to basic methods of a wrapped object.

**class** **Accessor**(*wrapper*)

Bases: [object](#)

An accessor: creates a simple uniform interface to treat the wrapped object element-wise (get/set/iterate over elements).

Generated automatically for each table on creation, doesn't need to be created explicitly.

**subcolumn**(*idx*, *wrapped=False*)

Return a subcolumn at index *idx*.

If *wrapped==True*, return a new wrapper containing the column; otherwise, just return the column.

**static from\_array**(*array*, *index=None*, *force\_copy=False*, *wrapped=False*)

Build a new object of the type corresponding to the wrapper from the supplied *array* (a 1D numpy array or a list).

If *force\_copy==True*, make a copy of supplied array. If *wrapped==True*, return a new wrapper containing the column; otherwise, just return the column.

**classmethod from\_columns**(*columns*, *column\_names=None*, *index=None*, *wrapped=False*)

Build a new object of the type corresponding to the wrapper from the supplied *columns* (a list of columns; only length-1 lists is supported).

*column\_names* parameter is ignored. If *wrapped==True*, return a new wrapper containing the column; otherwise, just return the column.

**array\_replaced**(*array*, *force\_copy=False*, *preserve\_index=False*, *wrapped=False*)

Return a copy of the column with the data replaced by *array*.

All of the parameters are the same as in [from\\_array\(\)](#).

**get\_index()**

Get index of the given 1D trace, or *None* if none is available

**get\_type()**

Get a string representing the wrapped object type

**copy**(*wrapped=False*)

Copy the object.

If *wrapped==True*, return a new wrapper containing the object copy; otherwise, just return the copy.

**ndim()**

**shape()**

**class** pylablib.core.dataproc.table\_wrap.**Array1DWrapper**(*container*)

Bases: [pylablib.core.dataproc.table\\_wrap.1DWrapper](#)

A wrapper for a 1D numpy array.

Provides a uniform access to basic methods of a wrapped object.

**get\_deleted**(*idx*, *wrapped=False*)

Return a copy of the column with the data at index *idx* deleted.

If `wrapped==True`, return a new wrapper containing the column; otherwise, just return the column.

**get\_inserted** (*idx, val, wrapped=False*)

Return a copy of the column with the data *val* added at index *idx*.

If `wrapped==True`, return a new wrapper containing the column; otherwise, just return the column.

**insert** (*idx, val*)

Add data *val* to index *idx*

**get\_appended** (*val, wrapped=False*)

Return a copy of the column with the data *val* appended at the end.

If `wrapped==True`, return a new wrapper containing the column; otherwise, just return the column.

**append** (*val*)

Append data *val* to the end

**subcolumn** (*idx, wrapped=False*)

Return a subcolumn at index *idx*.

If `wrapped==True`, return a new wrapper containing the column; otherwise, just return the column.

**static from\_array** (*array, index=None, force\_copy=False, wrapped=False*)

Build a new object of the type corresponding to the wrapper from the supplied *array* (a 1D numpy array or a list).

If `force_copy==True`, make a copy of supplied array. If `wrapped==True`, return a new wrapper containing the column; otherwise, just return the column.

**get\_type** ()

Get a string representing the wrapped object type

**copy** (*wrapped=False*)

Copy the object.

If `wrapped==True`, return a new wrapper containing the object copy; otherwise, just return the copy.

**class Accessor** (*wrapper*)

Bases: `object`

An accessor: creates a simple uniform interface to treat the wrapped object element-wise (get/set/iterate over elements).

Generated automatically for each table on creation, doesn't need to be created explicitly.

**array\_replaced** (*array, force\_copy=False, preserve\_index=False, wrapped=False*)

Return a copy of the column with the data replaced by *array*.

All of the parameters are the same as in `from_array()`.

**classmethod from\_columns** (*columns, column\_names=None, index=None, wrapped=False*)

Build a new object of the type corresponding to the wrapper from the supplied *columns* (a list of columns; only length-1 lists is supported).

*column\_names* parameter is ignored. If `wrapped==True`, return a new wrapper containing the column; otherwise, just return the column.

**get\_index** ()

Get index of the given 1D trace, or `None` if none is available

**ndim** ()

**shape** ()

**class** `pylablib.core.dataproc.table_wrap.Series1DWrapper` (*container*)

Bases: `pylablib.core.dataproc.table_wrap.I1DWrapper`

A wrapper for a pandas Series object.

Provides a uniform access to basic methods of a wrapped object.

**get\_deleted** (*idx*, *wrapped=False*)

Return a copy of the column with the data at index *idx* deleted.

If *wrapped==True*, return a new wrapper containing the column; otherwise, just return the column.

**get\_inserted** (*idx*, *val*, *wrapped=False*)

Return a copy of the column with the data *val* added at index *idx*.

If *wrapped==True*, return a new wrapper containing the column; otherwise, just return the column.

**get\_appended** (*val*, *wrapped=False*)

Return a copy of the column with the data *val* appended at the end.

If *wrapped==True*, return a new wrapper containing the column; otherwise, just return the column.

**subcolumn** (*idx*, *wrapped=False*)

Return a subcolumn at index *idx*.

If *wrapped==True*, return a new wrapper containing the column; otherwise, just return the column.

**static from\_array** (*array*, *index=None*, *force\_copy=False*, *wrapped=False*)

Build a new object of the type corresponding to the wrapper from the supplied *array* (a 1D numpy array or a list).

If *force\_copy==True*, make a copy of supplied array. If *wrapped==True*, return a new wrapper containing the column; otherwise, just return the column.

**get\_index** ()

Get index of the given 1D trace, or *None* if none is available

**get\_type** ()

Get a string representing the wrapped object type

**copy** (*wrapped=False*)

Copy the object.

If *wrapped==True*, return a new wrapper containing the object copy; otherwise, just return the copy.

**class Accessor** (*wrapper*)

Bases: `object`

An accessor: creates a simple uniform interface to treat the wrapped object element-wise (get/set/iterate over elements).

Generated automatically for each table on creation, doesn't need to be created explicitly.

**array\_replaced** (*array*, *force\_copy=False*, *preserve\_index=False*, *wrapped=False*)

Return a copy of the column with the data replaced by *array*.

All of the parameters are the same as in `from_array()`.

**classmethod from\_columns** (*columns*, *column\_names=None*, *index=None*, *wrapped=False*)

Build a new object of the type corresponding to the wrapper from the supplied *columns* (a list of columns; only length-1 lists is supported).

*column\_names* parameter is ignored. If *wrapped==True*, return a new wrapper containing the column; otherwise, just return the column.

**ndim** ()

**shape()**

**class** `pylablib.core.dataproc.table_wrap.I2DWrapper` (*container*, *r=None*, *c=None*,  
*t=None*)

Bases: `pylablib.core.dataproc.table_wrap.IGenWrapper`

A wrapper containing a 2D object (a 2D numpy array or a pandas DataFrame object).

Provides a uniform access to basic methods of a wrapped object.

**classmethod** `from_columns` (*columns*, *column\_names=None*, *index=None*, *wrapped=False*)

Build a new object of the type corresponding to the wrapper from the supplied *columns* (a list of columns).

*column\_names* supplies names of the columns (only relevant for `DataFrame2DWrapper`). If *wrapped==True*, return a new wrapper containing the table; otherwise, just return the table.

**columns\_replaced** (*columns*, *preserve\_index=False*, *wrapped=False*)

Return copy of the object with the data replaced by *columns*.

If *wrapped==True*, return a new wrapper containing the table; otherwise, just return the table.

**static** `from_array` (*array*, *column\_names=None*, *index=None*, *force\_copy=False*,  
*wrapped=False*)

Build a new object of the type corresponding to the wrapper from the supplied *array* (a list of rows or a 2D numpy array).

*column\_names* supplies names of the columns (only relevant for `DataFrame2DWrapper`). If *wrapped==True*, return a new wrapper containing the table; otherwise, just return the table.

**array\_replaced** (*array*, *preserve\_index=None*, *force\_copy=False*, *wrapped=False*)

Return a copy of the column with the data replaced by *array*.

All of the parameters are the same as in `from_array()`.

**get\_index()**

Get index of the given 2D table, or `None` if none is available

**get\_type()**

Get a string representing the wrapped object type

**copy** (*wrapped=False*)

Copy the object.

If *wrapped==True*, return a new wrapper containing the table; otherwise, just return the table.

**column** (*idx*, *wrapped=False*)

Get a column at index *idx*.

Return a 1D numpy array for a 2D numpy array object, and an Series object for a pandas DataFrame. If *wrapped==True*, return a new wrapper containing the column; otherwise, just return the column.

**subtable** (*idx*, *wrapped=False*)

Return a subtable at index *idx*.

If *wrapped==True*, return a new wrapper containing the table; otherwise, just return the table.

**ndim()**

**shape()**

**class** `pylablib.core.dataproc.table_wrap.Array2DWrapper` (*container*)

Bases: `pylablib.core.dataproc.table_wrap.I2DWrapper`

A wrapper for a 2D numpy array.

Provides a uniform access to basic methods of a wrapped object.

**set\_container** (*cont*)

**class RowAccessor** (*wrapper, storage*)

Bases: `object`

A row accessor: creates a simple uniform interface to treat the wrapped object row-wise (append/insert/delete/iterate over rows).

Generated automatically for each table on creation, doesn't need to be created explicitly.

**get\_deleted** (*idx, wrapped=False*)

Return a new table with the rows at *idx* deleted.

If *wrapped*==True, return a new wrapper containing the table; otherwise, just return the table.

**get\_inserted** (*idx, val, wrapped=False*)

Return a new table with new rows given by *val* inserted at *idx*.

If *wrapped*==True, return a new wrapper containing the table; otherwise, just return the table.

**insert** (*idx, val*)

Insert new rows given by *val* at index *idx*.

**get\_appended** (*val, wrapped=False*)

Return a new table with new rows given by *val* appended to the end of the table.

If *wrapped*==True, return a new wrapper containing the table; otherwise, just return the table.

**append** (*val*)

Insert new rows given by *val* to the end of the table

**class ColumnAccessor** (*wrapper, storage*)

Bases: `object`

A column accessor: creates a simple uniform interface to treat the wrapped object column-wise (append/insert/delete/iterate over columns).

Generated automatically for each table on creation, doesn't need to be created explicitly.

**get\_deleted** (*idx, wrapped=False*)

Return a new table with the columns at *idx* deleted.

If *wrapped*==True, return a new wrapper containing the table; otherwise, just return the table.

**get\_inserted** (*idx, val, wrapped=False*)

Return a new table with new columns given by *val* inserted at *idx*.

If *wrapped*==True, return a new wrapper containing the table; otherwise, just return the table.

**insert** (*idx, val*)

Insert new columns given by *val* at index *idx*.

**get\_appended** (*val, wrapped=False*)

Return a new table with new columns given by *val* appended to the end of the table.

If *wrapped*==True, return a new wrapper containing the table; otherwise, just return the table.

**append** (*val*)

Insert new columns given by *val* to the end of the table

**set\_names** (*names*)

Set column names (does nothing)

**get\_names** ()

Get column names (all names are None)

**get\_column\_index** (*idx*)

Get number index for a given column index

**class TableAccessor** (*storage*)

Bases: `object`

A table accessor: accessing the table data through this interface returns an object of the appropriate type (numpy array for numpy wrapped object, and a DataFrame for a pandas DataFrame wrapped object).

Generated automatically for each table on creation, doesn't need to be created explicitly.

**subtable** (*idx*, *wrapped=False*)

Return a subtable at index *idx* of the appropriate type (2D numpy array).

If *wrapped==True*, return a new wrapper containing the table; otherwise, just return the table.

**column** (*idx*, *wrapped=False*)

Get a column at index *idx* as a 1D numpy array.

If *wrapped==True*, return a new wrapper containing the column; otherwise, just return the column.

**classmethod from\_columns** (*columns*, *column\_names=None*, *index=None*, *wrapped=False*)

Build a new object of the type corresponding to the wrapper from the supplied *columns* (a list of columns).

If *wrapped==True*, return a new wrapper containing the table; otherwise, just return the table. *column\_names* parameter is ignored.

**static from\_array** (*array*, *column\_names=None*, *index=None*, *force\_copy=False*, *wrapped=False*)

Build a new object of the type corresponding to the wrapper from the supplied *array* (a list of rows or a 2D numpy array).

If *wrapped==True*, return a new wrapper containing the table; otherwise, just return the table. *column\_names* parameter is ignored.

**get\_type** ()

Get a string representing the wrapped object type

**copy** (*wrapped=False*)

Copy the object.

If *wrapped==True*, return a new wrapper containing the table; otherwise, just return the table.

**array\_replaced** (*array*, *preserve\_index=None*, *force\_copy=False*, *wrapped=False*)

Return a copy of the column with the data replaced by *array*.

All of the parameters are the same as in `from_array()`.

**columns\_replaced** (*columns*, *preserve\_index=False*, *wrapped=False*)

Return copy of the object with the data replaced by *columns*.

If *wrapped==True*, return a new wrapper containing the table; otherwise, just return the table.

**get\_index** ()

Get index of the given 2D table, or None if none is available

**ndim** ()

**shape** ()

**class** `pylablib.core.dataproc.table_wrap.DataFrame2DWrapper` (*container*)

Bases: `pylablib.core.dataproc.table_wrap.I2DWrapper`

A wrapper for a pandas DataFrame object.

Provides a uniform access to basic methods of a wrapped object.

```
class RowAccessor (wrapper, storage)
```

Bases: `object`

A row accessor: creates a simple uniform interface to treat the wrapped object row-wise (append/insert/delete/iterate over rows).

Generated automatically for each table on creation, doesn't need to be created explicitly.

```
get_deleted (idx, wrapped=False)
```

Return a copy of the column with the data at index *idx* deleted.

If `wrapped==True`, return a new wrapper containing the column; otherwise, just return the column.

```
get_inserted (idx, val, wrapped=False)
```

Return a new table with new rows given by *val* inserted at *idx*.

If `wrapped==True`, return a new wrapper containing the table; otherwise, just return the table.

```
insert (idx, val)
```

Insert new rows given by *val* at index *idx*.

```
get_appended (val, wrapped=False)
```

Return a new table with new rows given by *val* appended to the end of the table.

If `wrapped==True`, return a new wrapper containing the table; otherwise, just return the table.

```
append (val)
```

Insert new rows given by *val* to the end of the table

```
class ColumnAccessor (wrapper, storage)
```

Bases: `object`

A column accessor: creates a simple uniform interface to treat the wrapped object column-wise (append/insert/delete/iterate over columns).

Generated automatically for each table on creation, doesn't need to be created explicitly.

```
get_deleted (idx, wrapped=False)
```

Return a new table with the columns at *idx* deleted.

If `wrapped==True`, return a new wrapper containing the table; otherwise, just return the table.

```
get_inserted (idx, val, column_name=None, wrapped=False)
```

Return a new table with new columns given by *val* inserted at *idx*.

If `wrapped==True`, return a new wrapper containing the table; otherwise, just return the table.

```
insert (idx, val, column_name=None)
```

Insert new columns given by *val* at index *idx*

```
get_appended (val, column_name=None, wrapped=False)
```

Return a new table with new columns given by *val* appended to the end of the table.

If `wrapped==True`, return a new wrapper containing the table; otherwise, just return the table.

```
append (val, column_name=None)
```

Insert new columns given by *val* to the end of the table

```
set_names (names)
```

Set column names

```
get_names ()
```

Get column names



**get\_column\_index** (*idx*)

Get number index for a given column index

**class TableAccessor** (*storage*)

Bases: `object`

A table accessor: accessing the table data through this interface returns an object of the appropriate type (numpy array for numpy wrapped object, and a DataFrame for a pandas DataFrame wrapped object).

Generated automatically for each table on creation, doesn't need to be created explicitly.

**subtable** (*idx*, *wrapped=False*)

Return a subtable at index *idx* of the appropriate type (pandas DataFrame).

If *wrapped==True*, return a new wrapper containing the table; otherwise, just return the table.

**column** (*idx*, *wrapped=False*)

Get a column at index *idx* as a pandas Series object.

If *wrapped==True*, return a new wrapper containing the column; otherwise, just return the column.

**classmethod from\_columns** (*columns*, *column\_names=None*, *index=None*, *wrapped=False*)

Build a new object of the type corresponding to the wrapper from the supplied *columns* (a list of columns).

*column\_names* supplies names of the columns (only relevant for `DataFrame2DWrapper`). If *wrapped==True*, return a new wrapper containing the table; otherwise, just return the table.

**static from\_array** (*array*, *column\_names=None*, *index=None*, *force\_copy=False*, *wrapped=False*)

Build a new object of the type corresponding to the wrapper from the supplied *array* (a list of rows or a 2D numpy array).

*column\_names* supplies names of the columns (only relevant for `DataFrame2DWrapper`). If *wrapped==True*, return a new wrapper containing the table; otherwise, just return the table.

**get\_index** ()

Get index of the given 2D table, or None if none is available

**get\_type** ()

Get a string representing the wrapped object type

**copy** (*wrapped=False*)

Copy the object. If *wrapped==True*, return a new wrapper containing the table; otherwise, just return the table

**array\_replaced** (*array*, *preserve\_index=None*, *force\_copy=False*, *wrapped=False*)

Return a copy of the column with the data replaced by *array*.

All of the parameters are the same as in `from_array()`.

**columns\_replaced** (*columns*, *preserve\_index=False*, *wrapped=False*)

Return copy of the object with the data replaced by *columns*.

If *wrapped==True*, return a new wrapper containing the table; otherwise, just return the table.

**ndim** ()

**shape** ()

`pylablib.core.dataproc.table_wrap.wrap1d` (*container*)

Wrap a 1D container (a 1D numpy array or a pandas Series) into an appropriate wrapper

`pylablib.core.dataproc.table_wrap.wrap2d` (*container*)

Wrap a 2D container (a 2D numpy array or a pandas DataFrame) into an appropriate wrapper

`pylablib.core.dataproc.table_wrap.wrap(container)`

Wrap container (a numpy array, a pandas Series or a pandas DataFrame) into an appropriate wrapper

## pylablib.core.dataproc.transform module

**class** `pylablib.core.dataproc.transform.LinearTransform` (*tmatr=None, shift=None, ndim=2*)

Bases: `object`

A generic linear transform which combines an affine transform with a given matrix and an additional shift.

### Parameters

- **tmatr** – translational matrix (if `None`, use a unity matrix)
- **shift** – added shift (if `None`, use a zero shift)
- **ndim** – if both *tmatr* and *shift* are `None`, specifies the dimensionality of the transform; otherwise, ignored

**i** (*coord, shift=True*)

**inverted**()

Return inverted transformation

**preceded**(*trans*)

Return a combined transformation which result from applying this transformation followed by *trans*

**followed**(*trans*)

Return a combined transformation which result from applying *trans* followed by this transformation

**shifted**(*shift, preceded=False*)

Return a transform with an added shift before or after (depending of *preceded*) the current one

**multiplied**(*mult, preceded=False*)

Return a transform with an added scaling before or after (depending of *preceded*) the current one.

*mult* can be a single number (scale), a 1D vector (scaling for each axis independently), or a matrix.

**rotated2d**(*deg, preceded=False*)

Return a transform with an added rotation before or after (depending of *preceded*) the current one.

Only applies to 2D transforms.

**class** `pylablib.core.dataproc.transform.Indexed2DTransform` (*tmatr=None, shift=None, rigid=False*)

Bases: `pylablib.core.dataproc.transform.LinearTransform`

A restriction of `LinearTransform` which only applies to 2D and only allows rotations by multiples of 90 degrees.

### Parameters

- **tmatr** – translational matrix (if `None`, use a unity matrix)
- **shift** – added shift (if `None`, use a zero shift)
- **rigid** – if `True`, only allow orthogonal transforms, i.e., no scaling

**rotated2d**(*deg, preceded=False*)

Return a transform with an added rotation before or after (depending of *preceded*) the current one.

Only applies to 2D transforms.

**followed** (*trans*)

Return a combined transformation which result from applying *trans* followed by this transformation

**i** (*coord*, *shift=True*)

**inverted** ()

Return inverted transformation

**multiplied** (*mult*, *preceded=False*)

Return a transform with an added scaling before or after (depending of *preceded*) the current one.

*mult* can be a single number (scale), a 1D vector (scaling for each axis independently), or a matrix.

**preceded** (*trans*)

Return a combined transformation which result from applying this transformation followed by *trans*

**shifted** (*shift*, *preceded=False*)

Return a transform with an added shift before or after (depending of *preceded*) the current one

## pylablib.core.dataproc.utils module

Generic utilities for dealing with numerical arrays.

`pylablib.core.dataproc.utils.is_ascending` (*trace*)

Check the if the trace is ascending.

If it has more than 1 dimension, check all lines along 0'th axis.

`pylablib.core.dataproc.utils.is_descending` (*trace*)

Check if the trace is descending.

If it has more than 1 dimension, check all lines along 0'th axis.

`pylablib.core.dataproc.utils.is_ordered` (*trace*)

Check if the trace is ordered (ascending or descending).

If it has more than 1 dimension, check all lines along 0'th axis.

`pylablib.core.dataproc.utils.is_linear` (*trace*)

Check if the trace is linear (values go with a constant step).

If it has more than 1 dimension, check all lines along 0'th axis (with the same step for all).

`pylablib.core.dataproc.utils.get_x_column` (*t*, *x\_column=None*, *idx\_default=False*)

Get x column of the table.

*x\_column* can be

- an array: return as is;
- '#': return index array;
- None: equivalent to '#' for 1D data if `idx_default==False`, or to 0 otherwise;
- integer: return the column with this index.

`pylablib.core.dataproc.utils.get_y_column` (*t*, *y\_column=None*)

Get y column of the table.

*y\_column* can be

- an array: return as is;
- '#': return index array;

- None: return *t* for 1D data, or the column 1 otherwise;
- integer: return the column with this index.

`pylablib.core.dataproc.utils.sort_by(t, x_column=None, reverse=False, stable=False)`

Sort a table using selected column as a key and preserving rows.

If `reverse==True`, sort in descending order. `x_column` values are described in `get_x_column()`. If `stable==True`, use stable sort (could be slower and uses more memory, but preserves the order of elements for the same key)

`pylablib.core.dataproc.utils.filter_by(t, columns=None, pred=None, exclude=False)`

Filter 1D or 2D array using a predicate.

If the data is 2D, `columns` contains indices of columns to be passed to the `pred` function. If `exclude==False`, drop all of the rows satisfying `pred` rather than keep them.

`pylablib.core.dataproc.utils.unique_slices(t, u_column)`

Split a table into subtables with different values in a given column.

Return a list of *t* subtables, each of which has a different (and equal among all rows in the subtable) value in `u_column`.

`pylablib.core.dataproc.utils.merge(ts, idx=None, as_array=True)`

Merge several tables column-wise.

If `idx` is not None, then it is a list of index columns (one column per table) used for merging. The rows that have the same value in the index columns are merged; if some values aren't contained in all the *ts*, the corresponding rows are omitted. If `idx` is None, just join the tables together (they must have the same number of rows).

If `as_array==True`, return a simple numpy array as a result; otherwise, return a pandas DataFrame if applicable (note that in this case all column names in all tables must be different to avoid conflicts)

**class** `pylablib.core.dataproc.utils.Range(start=None, stop=None)`

Bases: `object`

Single data range.

If `start` or `stop` are None, it's implied that they're at infinity (i.e., `Range(None, None)` is infinite). If the range object is None, it's implied that the range is empty

**start**

**stop**

**contains** (*x*)

Check if *x* is in the range

**intersect** (*\*rngs*)

Find an intersection of multiple ranges.

If the intersection is empty, return None.

**rescale** (*mult=1.0, shift=0.0*)

**tup** ()

`pylablib.core.dataproc.utils.find_closest_arg(xs, x, approach='both', ordered=False)`

Find the index of a value in *xs* that is closest to *x*.

*approach* can take values 'top', 'bottom' or 'both' and denotes from which side should array elements approach *x* (meaning that the found array element should be  $>x$ ,  $<x$  or just the closest one). If there are no elements lying on the desired side of *x* (e.g. `approach=='top'` and all elements of *xs* are less than *x*), the function returns None. if `ordered==True`, then *xs* is assumed to be in ascending or descending order, and binary search is implemented (works only for 1D arrays). if there are recurring elements, return any of them.

`pylablib.core.dataproc.utils.find_closest_value(xs, x, approach='both', ordered=False)`

`pylablib.core.dataproc.utils.get_range_indices(xs, xs_range, ordered=False)`

Find trace indices corresponding to the given range.

The range is defined as `xs_range[0]:xs_range[1]`, or infinite if `xs_range=None` (so the data is returned unchanged in that case). If `ordered==True`, then the function assumes that `xs` in ascending or descending order.

`pylablib.core.dataproc.utils.cut_to_range(t, xs_range, x_column=None, ordered=False)`

Cut the table to the given range based on `x_column`.

The range is defined as `xs_range[0]:xs_range[1]`, or infinite if `xs_range=None`. `x_column` is used to determine which column's values to use to check if the point is in range (see `get_x_column()`). If `ordered_x==True`, then the function assumes that `x_column` in ascending order.

`pylablib.core.dataproc.utils.cut_out_regions(t, regions, x_column=None, ordered=False, multi_pass=True)`

Cut the regions out of the `t` based on `x_column`.

`x_column` is used to determine which column's values to use to check if the point is in range (see `get_x_column()`). If `ordered_x==True`, then the function assumes that `x_column` in ascending order. If `multi_pass==False`, combine all indices before deleting the data in a single operation (works faster, but only for non-intersecting regions).

`pylablib.core.dataproc.utils.find_discrete_step(trace, min_fraction=1e-08, tolerance=1e-05)`

Try to find a minimal divisor of all steps in a 1D trace.

`min_fraction` is the minimal possible size of the divisor (relative to the minimal non-zero step size). `tolerance` is the tolerance of the division. Raise an `ArithmeticError` if no such value was found.

`pylablib.core.dataproc.utils.unwrap_mod_data(trace, wrap_range)`

Unwrap data given `wrap_range`.

Assume that every jump greater than `0.5*wrap_range` is not real and is due to value being restricted. Can be used to, e.g., unwrap the phase data.

`pylablib.core.dataproc.utils.pad_trace(trace, pad, mode='constant', cval=0.0)`

Expand 1D trace or a multi-column table for different convolution techniques.

Wrapper around `numpy.pad()`, but can handle pandas dataframes or multi-column arrays. Note that the index data is not preserved.

#### Parameters

- **trace** – 1D array-like object.
- **pad** (*int or tuple*) – Expansion size. Can be an integer, if pad on both sides is equal, or a 2-tuple (`left`, `right`) for pads on opposite sides.
- **mode** (*str*) – Expansion mode. Takes the same values as `numpy.pad()`. Common values are 'constant' (added values are determined by `cval`), 'edge' (added values are end values of the trace), 'reflect' (reflect trace with respect to its endpoint) or 'wrap' (wrap the values from the other size).
- **cval** (*float*) – If `mode=='constant'`, determines the expanded values.

`pylablib.core.dataproc.utils.xy2c(t)`

Convert a trace or a table from xy representation to a single complex data.

`t` is a 2D array with either 2 columns (x and y) or 3 columns (index, x and y). Return 2D array with either 1 column (c) or 2 columns (index and c).

`pylablib.core.dataproc.utils.c2xy(t)`

Convert the a trace or a table from complex representation to a split x and y data.

*t* is either 1D array (c data) or a 2D array with either 1 column (c) or 2 columns (index and c). Return 2D array with either 2 column (x and y) or 3 columns (index, x and y).

## Module contents

### pylablib.core.devio package

#### Submodules

### pylablib.core.devio.SCPi module

```
class pylablib.core.devio.SCPi.SCPiDevice(conn, term_write=None, term_read=None,
                                         wait_callback=None, backend='auto', back-
                                         end_defaults=None, failsafe=None, time-
                                         out=None, backend_params=None)
```

Bases: `pylablib.core.devio.comm_backend.ICommBackendWrapper`

A base class for a device controlled with the usual SCPI syntax.

#### Implements two functions:

- deals with composing and parsing of standard SCPI commands and simplifying repetitive property access routines
- implements automatic re-sending and reconnecting on communication failures (fail-safe mode)

#### Parameters

- **conn** – Connection parameters (depend on the backend). Can also be an opened `comm_backend.IDeviceCommBackend` class for a custom backend.
- **term\_write** (*str*) – Line terminator for writing operations.
- **wait\_callback** (*callable*) – A function to be called periodically (every 300ms by default) while waiting for operations to complete.
- **backend** (*str*) – Connection backend (e.g., 'serial' or 'visa').
- **backend\_defaults** – if not None, specifies a dictionary {backend: params} with default connection parameters (depending on the backend), which are added to *conn*
- **failsafe** (*bool*) – If True, the device is working in a fail-safe mode: if an operation times out, attempt to repeat it several times before raising error. If None, use the class value `_default_failsafe` (False by default).
- **timeout** (*float*) – Default timeout (in seconds).

#### Error

alias of `pylablib.core.devio.base.DeviceError`

**ReraiseError** = None

#### BackendError

alias of `pylablib.core.devio.comm_backend.DeviceBackendError`

**reconnect** (*new\_instrument=True, ignore\_error=True*)

Remake the connection.

If *new\_instrument==True*, create a new backend instance. If *ignore\_error==True*, ignore errors on closing.

**sleep** (*delay*)

Wait for *delay* seconds

**using\_write\_buffer** ()

Context manager for using a write buffer.

While it's active, all the consecutive *write()* operations are bundled together with ; delimiter. The actual write is performed at the *read()/ask()* operation or at the end of the block.

**get\_id** (*timeout=None*)

Get the device IDN. (query SCPI '\*IDN?' command)

**reset** ()

Reset the device (by default, "\*RST" command)

**get\_esr** (*timeout=None*)

Get the device status register (by default, "\*ESR?" command)

**wait\_sync** (*timeout=None, wait\_callback=None*)

Pause execution of the script until device overlapped commands (e.g., taking sweeps) are complete.

*timeout* and *wait\_callback* override default constructor parameters.

**wait\_dev** ()

Pause execution of the device commands until device overlapped commands (e.g., taking sweeps) are complete.

Note that the code execution is not paused.

**wait** (*wait\_type='sync', timeout=None, wait\_callback=None*)

Pause execution until device overlapped commands are complete.

*wait\_type* is either 'sync' (perform *wait\_sync()*), 'dev' (perform *wait\_dev()*) or 'none' (do nothing).

**static get\_arg\_type** (*arg*)

Autodetect argument type

**write** (*msg, arg=None, arg\_type=None, unit=None, bool\_selector=None, wait\_sync=None, read\_echo=False, read\_echo\_delay=0.0*)

Send a command.

#### Parameters

- **msg** (*str*) – Text message.
- **arg** – Optional argument to append in the end. If a list of arguments is supplied, the result is joined with ", ".
- **arg\_type** (*str*) – Argument type. Can be 'raw' (in which case data is sent raw), 'string', 'int', 'float', 'bool', a format string (such as '{:.3f}'), or a list of argument types (for an iterable argument); if format string is used and the argument is a list or a tuple, then it is expanded as a list of arguments (e.g., *arg\_type='{0};{1}'* with *arg=[1, 2]* will produce a string '1; 2'); if a list of types is used, each element of *arg* is converted using the corresponding type, and the result is joined with ", ".
- **unit** (*str*) – If not None, use it as a unit to append after the value.

- **bool\_selector** (*tuple*) – A tuple (*false\_value*, *true\_value*) of two strings to represent bool argument; by default, use `._bool_selector` attribute.
- **wait\_sync** – if `True`, append the sync command (specified as `._wait_sync_comm` attribute, `"*OPC?"` by default) after the message and pause the execution command is complete; useful in long set operations, where the device might ignore later inputs until the current command is complete; if `None`, use the class default `._default_write_sync` attribute (`False` by default).
- **read\_echo** (*bool*) – If `True`, read a single line after write.
- **read\_echo\_delay** (*float*) – The delay between write and read if `read_echo==True`.

**read** (*data\_type='string', timeout=None*)

Read data from the device.

*data\_type* determines the type of the data. Can be `'raw'` (just raw data), `'string'` (with trailing and leading spaces stripped), `'int'`, `'float'`, `'bool'` (interprets 0 or 'off' as `False`, anything else as `True`), `'value'` (returns tuple (*value*, *unit*), where *value* is float), a callable (return the result of this callable applied to the string value), a dictionary (return the stored value corresponding to the string value, or to the value converted into integer if the string value is not present), or a list of data types (the result is treated as a list of values with the given types separated by commas). *timeout* overrides the default value.

**class NoParameterCaller** (*device, kind*)

Bases: `object`

Class to simplify calling functions without a parameter

**ask** (*msg, data\_type='string', delay=0.0, timeout=None, read\_echo=False*)

Write a message and read a reply.

*msg* is the query message, *delay* is the delay between write and read. Other parameters are the same as in `read()`. If `read_echo==True`, assume that the device first echoes the input and skip it.

**close** ()

Close the backend

**get\_device\_variable** (*key*)

Get the value of a settings, status, or full info parameter

**get\_full\_info** (*include=0*)

Get dict {*name*: *value*} containing full device information (including status and settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or `"all"` (all available variables). Since the lowest priority is -10, setting `include=-10` queries all available variables, which is equivalent to `include="all"`.

**get\_full\_status** (*include=0*)

Get dict {*name*: *value*} containing the device status (including settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or `"all"` (all available variables). Since the lowest priority is -10, setting `include=-10` queries all available variables, which is equivalent to `include="all"`.

**get\_settings** (*include=0*)

Get dict {*name*: *value*} containing all the device settings.



*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting `include=-10` queries all available variables, which is equivalent to `include="all"`.

**is\_opened()**

Check if the device is connected

**lock** (*timeout=None*)

Lock the access to the device from other threads/processes (isn't necessarily implemented)

**locking** (*timeout=None*)

Context manager for lock & unlock

**open()**

Open the backend

**set\_device\_variable** (*key, value*)

Set the value of a settings parameter

**unlock()**

Unlock the access to the device from other threads/processes (isn't necessarily implemented)

**flush** (*one\_line=False*)

Flush the read buffer (read all the available data and return the number of bytes read).

If `one_line==True`, read only a single line.

**read\_binary\_array\_data** (*include\_header=False, timeout=None, flush\_term=True*)

Read a binary data in the from the device.

The data assumes the standard binary transfer header consisting of "#" symbol, then a single digit with the size of the length string, then the length string containing the length of the binary data (in bytes). If `include_header==True`, return the data with the header; otherwise, return only the content. If `flush_term==True`, flush the following line to skip terminator characters after the binary data, which are added by some devices. *timeout* overrides the default value.

**static parse\_array\_data** (*data, fmt, include\_header=False*)

Parse the data returned by the device. *fmt* is *DataFormat* description in numpy format (e.g., "<u2").

If `include_header==True`, the data is assumed to be in a (somewhat) standard SCPI format: b'#', then a single digit *s* denoting length of the size block, then *s* digits denoting length of the data (in bytes) followed by the actual data. Otherwise (`include_header==False`), assume that the header is already removed.

**apply\_settings** (*settings*)

Apply the settings.

*settings* is a dict {name: value} of the available device settings. Non-applicable settings are ignored.

## pylablib.core.devio.backend\_logger module

**class** pylablib.core.devio.backend\_logger.**BackendLogger** (*path*)

Bases: *object*

Backend logger.

Receives log requests from backends and stores them in a predefined file.

**Parameters** *path* – path to save the log

**start** (*header*)  
Start logging section

**stop** ()  
Stop logging section

**section** (*header*)  
Context manager for operations within a header

**log** (*operation, value*)  
Log the operation

`pylablib.core.devio.backend_logger.load_logfile(path)`  
Load backend log file.

Return a list of tuples [(header, section)], where header is the header name, and section is the list [(op, value)] with operations ("r", "w", or "e") and corresponding values.

### pylablib.core.devio.base module

**exception** `pylablib.core.devio.base.DeviceError`  
Bases: `RuntimeError`

Generic device communication error

**args**

**with\_traceback** ()  
Exception.with\_traceback(tb) – set self.\_\_traceback\_\_ to tb and return self.

### pylablib.core.devio.comm\_backend module

Routines for defining a unified interface across multiple backends.

**exception** `pylablib.core.devio.comm_backend.DeviceBackendError` (*exc*)  
Bases: `pylablib.core.devio.base.DeviceError`

Generic exception relaying a backend error

**args**

**with\_traceback** ()  
Exception.with\_traceback(tb) – set self.\_\_traceback\_\_ to tb and return self.

`pylablib.core.devio.comm_backend.reraise(func)`  
Wrapper for a backend method which intercepts backend exceptions and re-emits them as a subclass of `DeviceBackendError` defined in the class

`pylablib.core.devio.comm_backend.logerror(func)`  
Wrapper for a backend method which logs if any errors escaped

**class** `pylablib.core.devio.comm_backend.IDeviceCommBackend` (*conn*, *timeout=None*,  
*term\_write=None*,  
*term\_read=None*,  
*datatype='auto'*,  
*reraise\_error=None*)

Bases: `object`

An abstract class for a device communication backend.

Connection is automatically opened on creation.

## Parameters

- **conn** – Connection parameters (depend on the backend).
- **timeout** (*float*) – Default timeout (in seconds).
- **term\_write** (*str*) – Line terminator for writing operations.
- **term\_read** (*str*) – Line terminator for reading operations.
- **datatype** (*str*) – Type of the returned data; can be "bytes" (return bytes object), "str" (return str object), or "auto" (default Python result: str in Python 2 and bytes in Python 3)
- **reraise\_error** – if not None, specifies an error to be re-raised on any back-end exception (by default, use backend-specific error); should be a subclass of *DeviceBackendError*.

## BackendError = None

Base class for the errors raised by the backend operations

## Error

alias of *DeviceBackendError*

## classmethod combine\_conn (conn1, conn2)

Combined two connection parameters into a single dictionary (*conn1* overrides *conn2*)

## classmethod get\_backend\_name ()

Get string representation of the backend (e.g., "serial", "visa", or "network")

## open ()

Open the connection

## close ()

Close the connection

## is\_opened ()

Check if the device is connected

## lock (timeout=None)

Lock the access to the device from other threads/processes (isn't necessarily implemented)

## unlock ()

Unlock the access to the device from other threads/processes (isn't necessarily implemented)

## locking (timeout=None)

Context manager for lock & unlock

## setup\_cooldown (\*\*kwargs)

Setup cooldown times for various operations.

The arguments are of the form *kind=value*, where *value* is the cooldown time (in seconds), and *kind* is the operation kind (common kinds are *open*, *close*, *read*, *write*, *timeout*, and *flush*). *kind* can also be *default* (default value for all *kind*), or *all* (reset all cooldown values to this value). The cooldowns of the given kinds are usually called after the corresponding operation (it is necessary for some devices, otherwise the communication can freeze or crush). Default cooldown values are specified by *\_default\_operation\_cooldown* class attribute dictionary.

## cooldown (kind='default')

Cooldown between the operations.

*kind* specifies the operation kind (common kinds are *open*, *close*, *read*, *write*, *timeout*, and *flush*); "default" corresponds to the default cooldown (usually, specified as 0). Called automatically by various backend operations, so usually there is no need to call explicitly.

**set\_timeout** (*timeout*)

Set operations timeout (in seconds)

**get\_timeout** ()

Get operations timeout (in seconds)

**using\_timeout** (*timeout=None*)

Context manager for usage of a different timeout inside a block

**readline** (*remove\_term=True, timeout=None, skip\_empty=True*)

Read a single line from the device.

#### Parameters

- **remove\_term** (*bool*) – If `True`, remove terminal characters from the result.
- **timeout** – Operation timeout. If `None`, use the default device timeout.
- **skip\_empty** (*bool*) – If `True`, ignore empty lines (works only for `remove_term==True`).

**readlines** (*lines\_num, remove\_term=True, timeout=None, skip\_empty=True*)

Read multiple lines from the device.

Parameters are the same as in `readline()`.

**read** (*size=None*)

Read data from the device.

If *size* is not `None`, read *size* bytes (the standard timeout applies); otherwise, read all available data (return immediately).

**flush\_read** ()

Flush the device output (read all the available data; return the number of bytes read)

**write** (*data, flush=True, read\_echo=False, read\_echo\_delay=0, read\_echo\_lines=1*)

Write data to the device.

If `flush==True`, flush the write buffer. If `read_echo==True`, wait for *read\_echo\_delay* seconds and then perform `readline()` (*read\_echo\_lines* times).

**ask** (*query, delay=0.0, read\_all=False*)

Perform a write followed by a read, with *delay* in between.

If `read_all==True`, read all the available data; otherwise, read a single line.

**static list\_resources** (*desc=False*)

List all available resources for this backend.

If `desc==False`, return list of connections (usually strings or tuples), which can be used to connect to the device. Otherwise, return a list of descriptions, which have more info, but can be backend-dependent.

Might not be implemented (depending on the backend), in which case returns `None`.

`pylablib.core.devio.comm_backend.remove_longest_term(msg, terms)`

Remove the longest terminator among *terms* from the end of the message.

**exception** `pylablib.core.devio.comm_backend.DeviceVisaError` (*exc*)

Bases: `pylablib.core.devio.comm_backend.DeviceBackendError`

Visa backend operation error

**args**

**with\_traceback** ()

Exception.with\_traceback(tb) – set self.\_\_traceback\_\_ to tb and return self.

```
class pylablib.core.devio.comm_backend.VisaDeviceBackend(conn, timeout=10.0,
                                                         term_write=None,
                                                         term_read=None,
                                                         do_lock=None,
                                                         datatype='auto',
                                                         reraise_error=None)
```

Bases: `pylablib.core.devio.comm_backend.IDeviceCommBackend`

NIVisa backend (via pyVISA).

Connection is automatically opened on creation.

#### Parameters

- **conn** (*str*) – Connection string.
- **timeout** (*float*) – Default timeout (in seconds).
- **term\_write** (*str*) – Line terminator for writing operations; appended to the data
- **term\_read** (*str*) – Line terminator for reading operations (specifies when `readline()` stops).
- **do\_lock** (*bool*) – If True, employ locking operations; otherwise, locking function does nothing.
- **datatype** (*str*) – Type of the returned data; can be "bytes" (return *bytes* object), "str" (return *str* object), or "auto" (default Python result: *str* in Python 2 and *bytes* in Python 3)
- **reraise\_error** – if not None, specifies an error to be re-raised on any back-end exception (by default, use backend-specific error); should be a subclass of `DeviceBackendError`.

#### BackendError

alias of `builtins.object`

#### Error

alias of `DeviceVisaError`

#### static list\_resources(*desc=False*)

List all available resources for this backend.

If `desc==False`, return list of connections (usually strings or tuples), which can be used to connect to the device. Otherwise, return a list of descriptions, which have more info, but can be backend-dependent.

Might not be implemented (depending on the backend), in which case returns None.

#### open()

Open the connection

#### close()

Close the connection

#### is\_opened()

Check if the device is connected

#### lock(*timeout=None*)

Lock the access to the device from other threads/processes

#### unlock()

Unlock the access to the device from other threads/processes

#### locking(*timeout=None*)

Context manager for lock & unlock

**set\_timeout** (*timeout*)

Set operations timeout (in seconds)

**get\_timeout** ()

Get operations timeout (in seconds)

**readline** (*remove\_term=True, timeout=None, skip\_empty=True*)

Read a single line from the device.

#### Parameters

- **remove\_term** (*bool*) – If `True`, remove terminal characters from the result.
- **timeout** – Operation timeout. If `None`, use the default device timeout.
- **skip\_empty** (*bool*) – If `True`, ignore empty lines (works only for `remove_term==True`).

**read** (*size=None*)

Read data from the device.

If *size* is not `None`, read *size* bytes (the standard timeout applies); otherwise, read all available data (return immediately).

**write** (*data, flush=True, read\_echo=False, read\_echo\_delay=0, read\_echo\_lines=1*)

Write data to the device.

If `read_echo==True`, wait for `read_echo_delay` seconds and then perform `readline()` (`read_echo_lines` times). `flush` parameter is ignored.

**ask** (*query, delay=0.0, read\_all=False*)

Perform a write followed by a read, with *delay* in between.

If `read_all==True`, read all the available data; otherwise, read a single line.

**classmethod combine\_conn** (*conn1, conn2*)

Combined two connection parameters into a single dictionary (*conn1* overrides *conn2*)

**cooldown** (*kind='default'*)

Cooldown between the operations.

*kind* specifies the operation kind (common kinds are `open`, `close`, `read`, `write`, `timeout`, and `flush`); `"default"` corresponds to the default cooldown (usually, specified as 0). Called automatically by various backend operations, so usually there is no need to call explicitly.

**flush\_read** ()

Flush the device output (read all the available data; return the number of bytes read)

**classmethod get\_backend\_name** ()

Get string representation of the backend (e.g., `"serial"`, `"visa"`, or `"network"`)

**readlines** (*lines\_num, remove\_term=True, timeout=None, skip\_empty=True*)

Read multiple lines from the device.

Parameters are the same as in `readline()`.

**setup\_cooldown** (*\*\*kwargs*)

Setup cooldown times for various operations.

The arguments are of the form `kind=value`, where *value* is the cooldown time (in seconds), and *kind* is the operation kind (common kinds are `open`, `close`, `read`, `write`, `timeout`, and `flush`). *kind* can also be `default` (default value for all *kind*), or `all` (reset all cooldown values to this value). The cooldowns of the given kinds are usually called after the corresponding operation (it is necessary for

some devices, otherwise the communication can freeze or crush). Default cooldown values are specified by `_default_operation_cooldown` class attribute dictionary.

**using\_timeout** (*timeout=None*)

Context manager for usage of a different timeout inside a block

**exception** `pylablib.core.devio.comm_backend.DeviceSerialError` (*exc*)

Bases: `pylablib.core.devio.comm_backend.DeviceBackendError`

Serial backend operation error

**args**

**with\_traceback** ()

Exception.with\_traceback(tb) – set self.\_\_traceback\_\_ to tb and return self.

**class** `pylablib.core.devio.comm_backend.SerialDeviceBackend` (*conn, timeout=10.0, term\_write=None, term\_read=None, connect\_on\_operation=False, open\_retry\_times=3, no\_dtrrts=False, datatype='auto', reraise\_error=None*)

Bases: `pylablib.core.devio.comm_backend.IDeviceCommBackend`

Serial backend (via pySerial).

Connection is automatically opened on creation.

#### Parameters

- **conn** – Connection parameters. Can be either a string (for a port), or a list/tuple (port, baudrate, bytesize, parity, stopbits, xonxoff, rtscts, dsrdtr) supplied to the serial connection (default is ('COM1', 19200, 8, 'N', 1, 0, 0, 0)), or a dict with the same parameters.
- **timeout** (*float*) – Default timeout (in seconds).
- **term\_write** (*str*) – Line terminator for writing operations; appended to the data
- **term\_read** (*str*) – List of possible single-char terminator for reading operations (specifies when `readline()` stops).
- **connect\_on\_operation** (*bool*) – If `True`, the connection is normally closed, and is opened only on the operations (normally two processes can't be simultaneously connected to the same device).
- **open\_retry\_times** (*int*) – Number of times the connection is attempted before giving up.
- **no\_dtrrts** (*bool*) – If `True`, turn off DTR and RTS status lines before opening (e.g., turns off reset-on-connection for Arduino controllers).
- **datatype** (*str*) – Type of the returned data; can be "bytes" (return *bytes* object), "str" (return *str* object), or "auto" (default Python result: *str* in Python 2 and *bytes* in Python 3)
- **reraise\_error** – if not `None`, specifies an error to be re-raised on any backend exception (by default, use backend-specific error); should be a subclass of `DeviceBackendError`.

**BackendError**

alias of `builtins.object`

**Error**

alias of `DeviceSerialError`

**open()**

Open the connection

**close()**

Close the connection

**is\_opened()**

Check if the device is connected

**single\_op()**

Context manager for a single operation.

If `connect_on_operation==True` during creation, wrapping several command in `single_op` prevents the connection from being closed and reopened between the operations (only opened in the beginning and closed in the end).

**set\_timeout(timeout)**

Set operations timeout (in seconds)

**get\_timeout()**

Get operations timeout (in seconds)

**readline(remove\_term=True, timeout=None, skip\_empty=True, error\_on\_timeout=True)**

Read a single line from the device.

**Parameters**

- **remove\_term** (*bool*) – If `True`, remove terminal characters from the result.
- **timeout** – Operation timeout. If `None`, use the default device timeout.
- **skip\_empty** (*bool*) – If `True`, ignore empty lines (works only for `remove_term==True`).
- **error\_on\_timeout** (*bool*) – If `False`, return an incomplete line instead of raising the error on timeout.

**read(size=None)**

Read data from the device.

If `size` is not `None`, read `size` bytes (usual timeout applies); otherwise, read all available data (return immediately).

**read\_multichar\_term(term, remove\_term=True, timeout=None, error\_on\_timeout=True)**

Read a single line with multiple possible terminators.

**Parameters**

- **term** – Either a string (single multi-char terminator) or a list of strings (multiple terminators).
- **remove\_term** (*bool*) – If `True`, remove terminal characters from the result.
- **timeout** – Operation timeout. If `None`, use the default device timeout.
- **error\_on\_timeout** (*bool*) – If `False`, return an incomplete line instead of raising the error on timeout.



**write** (*data*, *flush=True*, *read\_echo=False*, *read\_echo\_delay=0*, *read\_echo\_lines=1*)

Write data to the device.

If *flush==True*, flush the write buffer. If *read\_echo==True*, wait for *read\_echo\_delay* seconds and then perform *readline()* (*read\_echo\_lines* times).

**static list\_resources** (*desc=False*)

List all available resources for this backend.

If *desc==False*, return list of connections (usually strings or tuples), which can be used to connect to the device. Otherwise, return a list of descriptions, which have more info, but can be backend-dependent.

Might not be implemented (depending on the backend), in which case returns *None*.

**ask** (*query*, *delay=0.0*, *read\_all=False*)

Perform a write followed by a read, with *delay* in between.

If *read\_all==True*, read all the available data; otherwise, read a single line.

**classmethod combine\_conn** (*conn1*, *conn2*)

Combined two connection parameters into a single dictionary (*conn1* overrides *conn2*)

**cooldown** (*kind='default'*)

Cooldown between the operations.

*kind* specifies the operation kind (common kinds are *open*, *close*, *read*, *write*, *timeout*, and *flush*); "default" corresponds to the default cooldown (usually, specified as 0). Called automatically by various backend operations, so usually there is no need to call explicitly.

**flush\_read** ()

Flush the device output (read all the available data; return the number of bytes read)

**classmethod get\_backend\_name** ()

Get string representation of the backend (e.g., "serial", "visa", or "network")

**lock** (*timeout=None*)

Lock the access to the device from other threads/processes (isn't necessarily implemented)

**locking** (*timeout=None*)

Context manager for lock & unlock

**readlines** (*lines\_num*, *remove\_term=True*, *timeout=None*, *skip\_empty=True*)

Read multiple lines from the device.

Parameters are the same as in *readline()*.

**setup\_cooldown** (*\*\*kwargs*)

Setup cooldown times for various operations.

The arguments are of the form *kind=value*, where *value* is the cooldown time (in seconds), and *kind* is the operation kind (common kinds are *open*, *close*, *read*, *write*, *timeout*, and *flush*). *kind* can also be *default* (default value for all *kind*), or *all* (reset all cooldown values to this value). The cooldowns of the given kinds are usually called after the corresponding operation (it is necessary for some devices, otherwise the communication can freeze or crush). Default cooldown values are specified by *\_default\_operation\_cooldown* class attribute dictionary.

**unlock** ()

Unlock the access to the device from other threads/processes (isn't necessarily implemented)

**using\_timeout** (*timeout=None*)

Context manager for usage of a different timeout inside a block

**exception** *pylablib.core.devio.comm\_backend.DeviceFT232Error* (*exc*)

Bases: *pylablib.core.devio.comm\_backend.DeviceBackendError*

FT232 backend operation error

**args**

**with\_traceback()**

Exception.with\_traceback(tb) – set self.\_\_traceback\_\_ to tb and return self.

```
class pylablib.core.devio.comm_backend.FT232DeviceBackend(conn, timeout=10.0,
                                                         term_write=None,
                                                         term_read=None,
                                                         open_retry_times=3,
                                                         datatype='auto',
                                                         reraise_error=None)
```

Bases: `pylablib.core.devio.comm_backend.IDeviceCommBackend`

FT232 backend (via pyft232).

Connection is automatically opened on creation.

#### Parameters

- **conn** – Connection parameters. Can be either a string (for a port), or a list/tuple (port, baudrate, bytesize, parity, stopbits, xonxoff, rtscts) supplied to the serial connection (default is ('COM1', 19200, 8, 'N', 1, 0, 0, 0)), or a dict with the same parameters.
- **timeout** (*float*) – Default timeout (in seconds).
- **term\_write** (*str*) – Line terminator for writing operations; appended to the data
- **term\_read** (*str*) – List of possible single-char terminator for reading operations (specifies when `readline()` stops).
- **open\_retry\_times** (*int*) – Number of times the connection is attempted before giving up.
- **datatype** (*str*) – Type of the returned data; can be "bytes" (return *bytes* object), "str" (return *str* object), or "auto" (default Python result: *str* in Python 2 and *bytes* in Python 3)
- **reraise\_error** – if not None, specifies an error to be re-raised on any backend exception (by default, use backend-specific error); should be a subclass of `DeviceBackendError`.

#### BackendError

alias of `builtins.object`

#### Error

alias of `DeviceFT232Error`

**open()**

Open the connection

**close()**

Close the connection

**is\_opened()**

Check if the device is connected

**single\_op()**

Context manager for a single operation.

Does nothing.

**set\_timeout** (*timeout*)  
Set operations timeout (in seconds)

**get\_timeout** ()  
Get operations timeout (in seconds)

**readline** (*remove\_term=True, timeout=None, skip\_empty=True, error\_on\_timeout=True*)  
Read a single line from the device.

#### Parameters

- **remove\_term** (*bool*) – If `True`, remove terminal characters from the result.
- **timeout** – Operation timeout. If `None`, use the default device timeout.
- **skip\_empty** (*bool*) – If `True`, ignore empty lines (works only for `remove_term==True`).
- **error\_on\_timeout** (*bool*) – If `False`, return an incomplete line instead of raising the error on timeout.

**read** (*size=None*)  
Read data from the device.

If *size* is not `None`, read *size* bytes (usual timeout applies); otherwise, read all available data (return immediately).

**read\_multichar\_term** (*term, remove\_term=True, timeout=None, error\_on\_timeout=True*)  
Read a single line with multiple possible terminators.

#### Parameters

- **term** – Either a string (single multi-char terminator) or a list of strings (multiple terminators).
- **remove\_term** (*bool*) – If `True`, remove terminal characters from the result.
- **timeout** – Operation timeout. If `None`, use the default device timeout.
- **error\_on\_timeout** (*bool*) – If `False`, return an incomplete line instead of raising the error on timeout.

**write** (*data, flush=True, read\_echo=False, read\_echo\_delay=0, read\_echo\_lines=1*)  
Write data to the device.

If `flush==True`, flush the write buffer. If `read_echo==True`, wait for *read\_echo\_delay* seconds and then perform `readline()` (*read\_echo\_lines* times).

**static list\_resources** (*desc=False*)  
List all available resources for this backend.

If `desc==False`, return list of connections (usually strings or tuples), which can be used to connect to the device. Otherwise, return a list of descriptions, which have more info, but can be backend-dependent.

Might not be implemented (depending on the backend), in which case returns `None`.

**ask** (*query, delay=0.0, read\_all=False*)  
Perform a write followed by a read, with *delay* in between.

If `read_all==True`, read all the available data; otherwise, read a single line.

**classmethod combine\_conn** (*conn1, conn2*)  
Combined two connection parameters into a single dictionary (*conn1* overrides *conn2*)

**cooldown** (*kind='default'*)

Cooldown between the operations.

*kind* specifies the operation kind (common kinds are `open`, `close`, `read`, `write`, `timeout`, and `flush`); "default" corresponds to the default cooldown (usually, specified as 0). Called automatically by various backend operations, so usually there is no need to call explicitly.

**flush\_read** ()

Flush the device output (read all the available data; return the number of bytes read)

**classmethod get\_backend\_name** ()

Get string representation of the backend (e.g., "serial", "visa", or "network")

**lock** (*timeout=None*)

Lock the access to the device from other threads/processes (isn't necessarily implemented)

**locking** (*timeout=None*)

Context manager for lock & unlock

**readlines** (*lines\_num, remove\_term=True, timeout=None, skip\_empty=True*)

Read multiple lines from the device.

Parameters are the same as in `readline` ().

**setup\_cooldown** (*\*\*kwargs*)

Setup cooldown times for various operations.

The arguments are of the form `kind=value`, where *value* is the cooldown time (in seconds), and *kind* is the operation kind (common kinds are `open`, `close`, `read`, `write`, `timeout`, and `flush`). *kind* can also be `default` (default value for all *kind*), or `all` (reset all cooldown values to this value). The cooldowns of the given kinds are usually called after the corresponding operation (it is necessary for some devices, otherwise the communication can freeze or crush). Default cooldown values are specified by `_default_operation_cooldown` class attribute dictionary.

**unlock** ()

Unlock the access to the device from other threads/processes (isn't necessarily implemented)

**using\_timeout** (*timeout=None*)

Context manager for usage of a different timeout inside a block

**exception** `pylablib.core.devio.comm_backend.DeviceNetworkError` (*exc*)

Bases: `pylablib.core.devio.comm_backend.DeviceBackendError`

Network backend operation error

**args**

**with\_traceback** ()

Exception.with\_traceback(tb) – set self.\_\_traceback\_\_ to tb and return self.

**class** `pylablib.core.devio.comm_backend.NetworkDeviceBackend` (*conn, timeout=10.0,*  
*term\_write=None,*  
*term\_read=None,*  
*datatype='auto',*  
*reraise\_error=None*)

Bases: `pylablib.core.devio.comm_backend.IDeviceCommBackend`

Serial backend (via pySerial).

Connection is automatically opened on creation.

**Parameters**

- **conn** – Connection parameters. Can be either a string "IP:port" (e.g., "127.0.0.1:80"), or a tuple (IP, port), where *IP* is a string and *port* is a number.
- **timeout** (*float*) – Default timeout (in seconds).
- **term\_write** (*str*) – Line terminator for writing operations; appended to the data
- **term\_read** (*str*) – List of possible single-char terminator for reading operations (specifies when *readline()* stops).
- **datatype** (*str*) – Type of the returned data; can be "bytes" (return *bytes* object), "str" (return *str* object), or "auto" (default Python result: *str* in Python 2 and *bytes* in Python 3)
- **reraise\_error** – if not *None*, specifies an error to be re-raised on any back-end exception (by default, use backend-specific error); should be a subclass of *DeviceBackendError*.

---

**Note:** If *term\_read* is a string, its behavior is different from the VISA backend: instead of being a multi-char terminator it is assumed to be a set of single-char terminators. If multi-char terminator is required, *term\_read* should be a single-element list instead of a string.

---

#### BackendError

alias of `builtins.OSError`

#### Error

alias of *DeviceNetworkError*

#### open()

Open the connection

#### close()

Close the connection

#### is\_opened()

Check if the device is connected

#### set\_timeout(timeout)

Set operations timeout (in seconds)

#### get\_timeout()

Get operations timeout (in seconds)

#### readline(remove\_term=True, timeout=None, skip\_empty=True)

Read a single line from the device.

#### Parameters

- **remove\_term** (*bool*) – If *True*, remove terminal characters from the result.
- **timeout** – Operation timeout. If *None*, use the default device timeout.
- **skip\_empty** (*bool*) – If *True*, ignore empty lines (works only for *remove\_term==True*).

#### read(size=None)

Read data from the device.

If *size* is not *None*, read *size* bytes (usual timeout applies); otherwise, read all available data (return immediately).

**read\_multichar\_term** (*term*, *remove\_term=True*, *timeout=None*)

Read a single line with multiple possible terminators.

**Parameters**

- **term** – Either a string (single multi-char terminator) or a list of strings (multiple terminators).
- **remove\_term** (*bool*) – If `True`, remove terminal characters from the result.
- **timeout** – Operation timeout. If `None`, use the default device timeout.

**write** (*data*, *flush=True*, *read\_echo=False*, *read\_echo\_delay=0*, *read\_echo\_lines=1*)

Write data to the device.

If `read_echo==True`, wait for `read_echo_delay` seconds and then perform `readline()` (`read_echo_lines` times). `flush` parameter is ignored.

**ask** (*query*, *delay=0.0*, *read\_all=False*)

Perform a write followed by a read, with `delay` in between.

If `read_all==True`, read all the available data; otherwise, read a single line.

**classmethod combine\_conn** (*conn1*, *conn2*)

Combined two connection parameters into a single dictionary (`conn1` overrides `conn2`)

**cooldown** (*kind='default'*)

Cooldown between the operations.

`kind` specifies the operation kind (common kinds are `open`, `close`, `read`, `write`, `timeout`, and `flush`); `"default"` corresponds to the default cooldown (usually, specified as 0). Called automatically by various backend operations, so usually there is no need to call explicitly.

**flush\_read** ()

Flush the device output (read all the available data; return the number of bytes read)

**classmethod get\_backend\_name** ()

Get string representation of the backend (e.g., `"serial"`, `"visa"`, or `"network"`)

**static list\_resources** (*desc=False*)

List all available resources for this backend.

If `desc==False`, return list of connections (usually strings or tuples), which can be used to connect to the device. Otherwise, return a list of descriptions, which have more info, but can be backend-dependent.

Might not be implemented (depending on the backend), in which case returns `None`.

**lock** (*timeout=None*)

Lock the access to the device from other threads/processes (isn't necessarily implemented)

**locking** (*timeout=None*)

Context manager for lock & unlock

**readlines** (*lines\_num*, *remove\_term=True*, *timeout=None*, *skip\_empty=True*)

Read multiple lines from the device.

Parameters are the same as in `readline()`.

**setup\_cooldown** (*\*\*kwargs*)

Setup cooldown times for various operations.

The arguments are of the form `kind=value`, where `value` is the cooldown time (in seconds), and `kind` is the operation kind (common kinds are `open`, `close`, `read`, `write`, `timeout`, and `flush`). `kind` can also be `default` (default value for all kind), or `all` (reset all cooldown values to this value). The cooldowns of the given kinds are usually called after the corresponding operation (it is necessary for

some devices, otherwise the communication can freeze or crush). Default cooldown values are specified by `_default_operation_cooldown` class attribute dictionary.

**unlock()**

Unlock the access to the device from other threads/processes (isn't necessarily implemented)

**using\_timeout** (*timeout=None*)

Context manager for usage of a different timeout inside a block

**exception** `pylablib.core.devio.comm_backend.DeviceUSBError` (*exc*)

Bases: `pylablib.core.devio.comm_backend.DeviceBackendError`

USB backend operation error

**args**

**with\_traceback()**

Exception.with\_traceback(tb) – set self.\_\_traceback\_\_ to tb and return self.

```
class pylablib.core.devio.comm_backend.PyUSBDeviceBackend(conn, timeout=10.0,
                                                         term_write=None,
                                                         term_read=None,
                                                         check_read_size=True,
                                                         datatype='auto',
                                                         reraise_error=None)
```

Bases: `pylablib.core.devio.comm_backend.IDeviceCommBackend`

USB backend (via PyUSB package).

Connection is automatically opened on creation.

#### Parameters

- **conn** – Connection parameters. Can be either a string (for a port), or a list/tuple (vendorID, productID, index, endpoint\_read, endpoint\_write, backend) supplied to the connection (default is (0x0000, 0x0000, 0, 0x00, 0x01, 'libusb1'), which is invalid for most devices), or a dict with the same parameters. vendorID and productID specify device kind, index is an integer index (starting from zero) of the device among several identical (i.e., with the same ids) ones, and endpoint\_read and endpoint\_write specify connection endpoints for the specific device.
- **timeout** (*float*) – Default timeout (in seconds).
- **term\_write** (*str*) – Line terminator for writing operations; appended to the data
- **term\_read** (*str*) – List of possible single-char terminator for reading operations (specifies when `readline()` stops).
- **datatype** (*str*) – Type of the returned data; can be "bytes" (return *bytes* object), "str" (return *str* object), or "auto" (default Python result: *str* in Python 2 and *bytes* in Python 3)
- **reraise\_error** – if not None, specifies an error to be re-raised on any backend exception (by default, use backend-specific error); should be a subclass of `DeviceBackendError`.

#### BackendError

alias of `usb.core.USBError`

#### Error

alias of `DeviceUSBError`

**open()**  
Open the connection

**close()**  
Close the connection

**is\_opened()**  
Check if the device is connected

**set\_timeout(timeout)**  
Set operations timeout (in seconds)

**get\_timeout()**  
Get operations timeout (in seconds)

**readline(remove\_term=True, timeout=None, skip\_empty=True, error\_on\_timeout=True)**  
Read a single line from the device.

#### Parameters

- **remove\_term** (*bool*) – If `True`, remove terminal characters from the result.
- **timeout** – Operation timeout. If `None`, use the default device timeout.
- **skip\_empty** (*bool*) – If `True`, ignore empty lines (works only for `remove_term==True`).
- **error\_on\_timeout** (*bool*) – If `False`, return an incomplete line instead of raising the error on timeout.

**read(size=None, max\_read\_size=65536)**  
Read data from the device.

If *size* is not `None`, read *size* bytes (usual timeout applies); otherwise, read all available data (return immediately).

**read\_multichar\_term(term, remove\_term=True, timeout=None, error\_on\_timeout=True)**  
Read a single line with multiple possible terminators.

#### Parameters

- **term** – Either a string (single multi-char terminator) or a list of strings (multiple terminators).
- **remove\_term** (*bool*) – If `True`, remove terminal characters from the result.
- **timeout** – Operation timeout. If `None`, use the default device timeout.
- **error\_on\_timeout** (*bool*) – If `False`, return an incomplete line instead of raising the error on timeout.

**write(data, flush=True, read\_echo=False, read\_echo\_delay=0, read\_echo\_lines=1)**  
Write data to the device.

If `read_echo==True`, wait for *read\_echo\_delay* seconds and then perform `readline()` (*read\_echo\_lines* times). *flush* parameter is ignored.

**static list\_resources(desc=False, \*\*kwargs)**  
List all available resources for this backend.

If `desc==False`, return list of connections (usually strings or tuples), which can be used to connect to the device. Otherwise, return a list of descriptions, which have more info, but can be backend-dependent.

Might not be implemented (depending on the backend), in which case returns `None`.



**ask** (*query*, *delay=0.0*, *read\_all=False*)  
 Perform a write followed by a read, with *delay* in between.

If *read\_all==True*, read all the available data; otherwise, read a single line.

**classmethod combine\_conn** (*conn1*, *conn2*)  
 Combined two connection parameters into a single dictionary (*conn1* overrides *conn2*)

**cooldown** (*kind='default'*)  
 Cooldown between the operations.

*kind* specifies the operation kind (common kinds are *open*, *close*, *read*, *write*, *timeout*, and *flush*); "default" corresponds to the default cooldown (usually, specified as 0). Called automatically by various backend operations, so usually there is no need to call explicitly.

**flush\_read** ()  
 Flush the device output (read all the available data; return the number of bytes read)

**classmethod get\_backend\_name** ()  
 Get string representation of the backend (e.g., "serial", "visa", or "network")

**lock** (*timeout=None*)  
 Lock the access to the device from other threads/processes (isn't necessarily implemented)

**locking** (*timeout=None*)  
 Context manager for lock & unlock

**readlines** (*lines\_num*, *remove\_term=True*, *timeout=None*, *skip\_empty=True*)  
 Read multiple lines from the device.

Parameters are the same as in *readline* ().

**setup\_cooldown** (*\*\*kwargs*)  
 Setup cooldown times for various operations.

The arguments are of the form *kind=value*, where *value* is the cooldown time (in seconds), and *kind* is the operation kind (common kinds are *open*, *close*, *read*, *write*, *timeout*, and *flush*). *kind* can also be *default* (default value for all kind), or *all* (reset all cooldown values to this value). The cooldowns of the given kinds are usually called after the corresponding operation (it is necessary for some devices, otherwise the communication can freeze or crush). Default cooldown values are specified by *\_default\_operation\_cooldown* class attribute dictionary.

**unlock** ()  
 Unlock the access to the device from other threads/processes (isn't necessarily implemented)

**using\_timeout** (*timeout=None*)  
 Context manager for usage of a different timeout inside a block

**exception** *pylablib.core.devio.comm\_backend.DeviceRecordedError* (*exc*)  
 Bases: *pylablib.core.devio.comm\_backend.DeviceBackendError*

Recorded backend operation error

**args**

**with\_traceback** ()  
 Exception.with\_traceback(tb) – set self.\_\_traceback\_\_ to tb and return self.

**class** *pylablib.core.devio.comm\_backend.RecordedDeviceBackend* (*conn*,  
*datatype='auto'*,  
*reraise\_error=None*)

Bases: *pylablib.core.devio.comm\_backend.IDeviceCommBackend*

Recorded backend.

Connection is automatically opened on creation.

#### Parameters

- **conn** – connection parameters (recorded log path)
- **datatype** (*str*) – Type of the returned data; can be "bytes" (return *bytes* object), "str" (return *str* object), or "auto" (default Python result: *str* in Python 2 and *bytes* in Python 3)
- **reraise\_error** – if not `None`, specifies an error to be re-raised on any back-end exception (by default, use backend-specific error); should be a subclass of *DeviceBackendError*.

#### BackendError

alias of `builtins.OSError`

#### Error

alias of *DeviceRecordedError*

#### `open()`

Open the connection

#### `close()`

Close the connection

#### `is_opened()`

Check if the device is connected

#### `start(header)`

Start recorded section

#### `stop()`

Stop logging section

#### `section(header)`

#### `readline(remove_term=True, timeout=None, skip_empty=True)`

Read a single line from the device.

#### Parameters

- **remove\_term** (*bool*) – If `True`, remove terminal characters from the result.
- **timeout** – Operation timeout. If `None`, use the default device timeout.
- **skip\_empty** (*bool*) – If `True`, ignore empty lines (works only for `remove_term==True`).

#### `read(size=None)`

Read data from the device.

If *size* is not `None`, read *size* bytes (usual timeout applies); otherwise, read all available data (return immediately).

#### `write(data, flush=True, read_echo=False, read_echo_delay=0, read_echo_lines=1)`

Write data to the device.

If `flush==True`, flush the write buffer. If `read_echo==True`, wait for *read\_echo\_delay* seconds and then perform *readline()* (*read\_echo\_lines* times).

#### `ask(query, delay=0.0, read_all=False)`

Perform a write followed by a read, with *delay* in between.

If `read_all==True`, read all the available data; otherwise, read a single line.

**classmethod combine\_conn** (*conn1*, *conn2*)

Combined two connection parameters into a single dictionary (*conn1* overrides *conn2*)

**cooldown** (*kind*='default')

Cooldown between the operations.

*kind* specifies the operation kind (common kinds are `open`, `close`, `read`, `write`, `timeout`, and `flush`); "default" corresponds to the default cooldown (usually, specified as 0). Called automatically by various backend operations, so usually there is no need to call explicitly.

**flush\_read** ()

Flush the device output (read all the available data; return the number of bytes read)

**classmethod get\_backend\_name** ()

Get string representation of the backend (e.g., "serial", "visa", or "network")

**get\_timeout** ()

Get operations timeout (in seconds)

**static list\_resources** (*desc*=False)

List all available resources for this backend.

If *desc*==False, return list of connections (usually strings or tuples), which can be used to connect to the device. Otherwise, return a list of descriptions, which have more info, but can be backend-dependent.

Might not be implemented (depending on the backend), in which case returns None.

**lock** (*timeout*=None)

Lock the access to the device from other threads/processes (isn't necessarily implemented)

**locking** (*timeout*=None)

Context manager for lock & unlock

**readlines** (*lines\_num*, *remove\_term*=True, *timeout*=None, *skip\_empty*=True)

Read multiple lines from the device.

Parameters are the same as in `readline()`.

**set\_timeout** (*timeout*)

Set operations timeout (in seconds)

**setup\_cooldown** (*\*\*kwargs*)

Setup cooldown times for various operations.

The arguments are of the form *kind*=*value*, where *value* is the cooldown time (in seconds), and *kind* is the operation kind (common kinds are `open`, `close`, `read`, `write`, `timeout`, and `flush`). *kind* can also be `default` (default value for all *kind*), or `all` (reset all cooldown values to this value). The cooldowns of the given kinds are usually called after the corresponding operation (it is necessary for some devices, otherwise the communication can freeze or crush). Default cooldown values are specified by `_default_operation_cooldown` class attribute dictionary.

**unlock** ()

Unlock the access to the device from other threads/processes (isn't necessarily implemented)

**using\_timeout** (*timeout*=None)

Context manager for usage of a different timeout inside a block

`pylablib.core.devio.comm_backend.autodetect_backend` (*conn*, *default*='visa')

Try to determine the backend by the connection.

*default* specifies the default backend which is returned if the backend is unclear.

```
pylablib.core.devio.comm_backend.new_backend(conn, backend='auto', defaults=None,
                                              **kwargs)
```

Build new backend with the supplied parameters.

#### Parameters

- **conn** – Connection parameters (depend on the backend). Can be simply connection parameters (tuple or dict) for the given backend (e.g., "192.168.0.1" or ("COM1", 19200)), a tuple (backend, conn) which specifies both backend and connection (in which case it overrides the supplied backend), or an already opened backend (in which case it is returned as is)
- **backend** (*str*) – Backend type. Available backends are 'auto' (try to autodetect based on the connection), 'visa', 'serial', 'ft232', 'network', and 'pyusb'. Can also be directly a backend class (more appropriate for custom backends), or a tuple ('auto', backend), which is analogous to 'auto', but it returns the specified backend if the autodetection fails; by default, the fallback backend is 'visa', so 'auto' is exactly the same as ('auto', 'visa').
- **defaults** – if not None, specifies a dictionary {backend: params} with default connection parameters (depending on the backend), which are added to the connection parameters
- **\*\*kwargs** – parameters sent to the backend.

```
pylablib.core.devio.comm_backend.backend_error(backend, conn=None)
```

Return error class corresponding to the current backend.

Like `new_backend()`, allows setting `backend="auto"`, in which case `conn` is used to try and autodetect the backend kind (not completely reliable, should be avoided).

```
pylablib.core.devio.comm_backend.list_backend_resources(backend=None,
                                                         desc=False)
```

List all resources for the given backend.

If `backend` is None, return dictionary {backend: resources} for all available backends. If `desc==False`, return list of connections (usually strings or tuples), which can be used to connect to the device. Otherwise, return a list of descriptions, which have more info, but can be backend-dependent.

```
class pylablib.core.devio.comm_backend.ICommBackendWrapper(instr)
```

Bases: `pylablib.core.devio.interface.IDevice`

A base class for an instrument using a communication backend.

**Parameters** `instr` – Backend (assumed to be already opened).

```
class NoParameterCaller(device, kind)
```

Bases: `object`

Class to simplify calling functions without a parameter

```
apply_settings(settings)
```

Apply the settings.

`settings` is the dict {name: value} of the device available settings. Non-applicable settings are ignored.

```
get_device_variable(key)
```

Get the value of a settings, status, or full info parameter

```
get_full_info(include=0)
```

Get dict {name: value} containing full device information (including status and settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting `include=-10` queries all available variables, which is equivalent to `include="all"`.

**get\_full\_status** (*include=0*)

Get dict {name: value} containing the device status (including settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting `include=-10` queries all available variables, which is equivalent to `include="all"`.

**get\_settings** (*include=0*)

Get dict {name: value} containing all the device settings.

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting `include=-10` queries all available variables, which is equivalent to `include="all"`.

**set\_device\_variable** (*key, value*)

Set the value of a settings parameter

**open** ()

Open the backend

**close** ()

Close the backend

**is\_opened** ()

Check if the device is connected

**lock** (*timeout=None*)

Lock the access to the device from other threads/processes (isn't necessarily implemented)

**unlock** ()

Unlock the access to the device from other threads/processes (isn't necessarily implemented)

**locking** (*timeout=None*)

Context manager for lock & unlock

## pylablib.core.devio.data\_format module

Library for binary data encoding/decoding for device communication and dealing with different data format representations in different contexts (numpy, SCPI, etc.).

**class** `pylablib.core.devio.data_format.DataFormat` (*size, kind, byteorder*)

Bases: `object`

Describes data encoding for device communications.

### Parameters

- **size** (*int*) – Size of a single element (in bytes).
- **kind** (*str*) – Kind of the element. Can be 'i' (integer), 'u' (unsigned integer), 'f' (floating point) or 'ascii' (text representation).
- **byteorder** (*str*) – Byte order: '>' is big-endian (MSB first), '<' is little-endian (LSB first).

**flip\_byteorder()**

Flip byteorder of the description

**is\_ascii()**

Check if the format is textual

**static from\_desc**(*desc*, *str\_type*='numpy')

Build the format from the string description.

*str\_type* is the description format. Can be 'numpy' (numpy dtype description), 'struct' (struct description) or 'SCPI' (the standard SCPI description).

**static from\_desc\_SCPI**(*desc*, *border*='norm')

Build the format from the string SCPI description.

*border* describes byte order (either 'norm' or 'swap').

**to\_desc**(*str\_type*='auto')

Build a description string of this format.

*str\_type* can be 'auto' (similar to 'numpy', but also accepts 'ascii'), 'numpy', 'struct' or 'SCPI' (return tuple (*desc*, *border*)).

**convert\_from\_str**(*data*)

Convert the string data into an array

**convert\_to\_str**(*data*, *ascii\_format*='.5f')

Convert the array into a string data.

*ascii\_format* is the `str.format()` string for textual representation.

## pylablib.core.devio.interface module

**class** pylablib.core.devio.interface.IDevice

Bases: `object`

A base class for an instrument.

Contains some useful functions for dealing with device settings.

**open()**

Open the connection

**close()**

Close the connection

**is\_opened()**

Check if the device is connected

**class** NoParameterCaller(*device*, *kind*)

Bases: `object`

Class to simplify calling functions without a parameter

**get\_settings**(*include*=0)

Get dict {*name*: *value*} containing all the device settings.

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include*=-10 queries all available variables, which is equivalent to *include*="all".

**get\_full\_status** (*include=0*)

Get dict {name: value} containing the device status (including settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_full\_info** (*include=0*)

Get dict {name: value} containing full device information (including status and settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**apply\_settings** (*settings*)

Apply the settings.

*settings* is the dict {name: value} of the device available settings. Non-applicable settings are ignored.

**get\_device\_variable** (*key*)

Get the value of a settings, status, or full info parameter

**set\_device\_variable** (*key, value*)

Set the value of a settings parameter

**class** pylablib.core.devio.interface.IParameterClass (*name*)

Bases: `object`

A generic parameter class.

Deals with converting device interface representation and the ‘internal’ representation (e.g., names used in SCPI commands or integer indices). Also responsible for validating the user-passed and device-returned parameters.

Needs to define to methods: `__call__` for converting user parameters (‘alias’) into the device parameters (‘value’) and `i()` for the opposite conversion. In addition, it provides `using_device()` context manager to temporarily change the `device` attribute, which can be used by some parameter classes for device-dependent conversions.

**Parameters** *name* – parameter class name; used to match method arguments with corresponding classes.

**using\_device** (*device*)

Context manager for temporarily changing the `device` attribute to the given device instance

**docstring** ()

Get a parameter docstring

**i** (*value, device=None*)

Convert device parameter value into a corresponding use parameter value

If not None, *device* specifies the corresponding device instance for device-dependent conversion.

**class** pylablib.core.devio.interface.ICheckingParameterClass (*name*)

Bases: `pylablib.core.devio.interface.IParameterClass`

Parameter class which separately handles checking and conversion.

Specifies six methods: `check_value()`, `to_alias()` and `_value_error_str` for handling value-to-alias conversion, and `check_alias()`, `to_value()` and `_alias_error_str` for handling alias-to-value conversion.

**check\_alias** (*alias*)

Check if the alias is valid

**check\_value** (*value*)

Check if the device value is valid

**to\_value** (*alias*)

Convert the alias into a device value

**to\_alias** (*value*)

Convert the device value into an alias

**i** (*value*, *device=None*)

Convert device parameter value into a corresponding use parameter value

If not `None`, *device* specifies the corresponding device instance for device-dependent conversion.

**docstring** ()

Get a parameter docstring

**using\_device** (*device*)

Context manager for temporarily changing the `device` attribute to the given device instance

```
class pylablib.core.devio.interface.RangeParameterClass (name,          minval=None,  
                                                         maxval=None,  
                                                         out_of_range='error')
```

Bases: `pylablib.core.devio.interface.ICheckingParameterClass`

Parameter class for numerical values constrained to a certain range.

#### Parameters

- **name** – parameter class name
- **minval** – minimal allowed value (inclusive); `None` means no lower limit
- **maxval** – maximal allowed value (inclusive); `None` means no upper limit
- **out\_of\_range** – action if an out-of-range value is supplied; can be either `"error"` (raise an error), or `"truncate"` (truncate to the nearest limit).

**check\_value** (*value*)

Check if the device value is valid

**check\_alias** (*alias*)

Check if the alias is valid

**to\_value** (*alias*)

Convert the alias into a device value

**docstring** ()

Get a parameter docstring

**i** (*value*, *device=None*)

Convert device parameter value into a corresponding use parameter value

If not `None`, *device* specifies the corresponding device instance for device-dependent conversion.

**to\_alias** (*value*)

Convert the device value into an alias

**using\_device** (*device*)

Context manager for temporarily changing the `device` attribute to the given device instance



```
class pylablib.core.devio.interface.IEnumParameterClass (name,
                                                         allowed_alias='device_values',
                                                         allowed_value='exact',
                                                         alias_case=None,
                                                         value_case=None,
                                                         match_prefix=False)
```

Bases: `pylablib.core.devio.interface.ICheckingParameterClass`

Parameter class for a generic enum (i.e., predefined values) parameter.

**Defines two methods for handling conversion:**

- `_get_value_map` which returns a dictionary for converting device values into aliases,
- `_get_alias_map` which returns a dictionary for converting aliases into device values.

These methods need to be redefined in subclasses.

#### Parameters

- **name** – parameter class name
- **allowed\_alias** – specifies a range of allowed aliases; can be "exact" (only exact map matches are allowed), "device\_value" (exact map matches and raw device values are allowed), or "all" (all values are allowed); in the latter two cases the value not in the map are passed as is.
- **allowed\_value** – specifies a range of allowed device values; can be "exact" (only exact map matches are allowed), or "all" (all values are allowed); in the latter case the value not in the map is passed as is.
- **alias\_case** – default alias parameter case for string values; can be None (no case normalization), or "lower" or "upper" (any received or returned alias will be normalized into this case)
- **value\_case** – default value parameter case for string values; can be None (no case normalization), or "lower" or "upper" (any received or returned device value will be normalized into this case)
- **match\_prefix** – if True, then the keys in the value map (returned by `_get_value_map` method) are interpreted as prefixes, so in the value-to-alias conversion the converted value matches the map value if it just starts with it; in the case of ambiguity (several map values are prefixes for the same converted value), the exact match takes priority; useful for some SCPI devices, where the shorter version of the value can sometimes be returned.

**check\_value** (value)

Check if the device value is valid

**check\_alias** (alias)

Check if the alias is valid

**to\_value** (alias)

Convert the alias into a device value

**to\_alias** (value)

Convert the device value into an alias

**docstring** ()

Get a parameter docstring

**i** (value, device=None)

Convert device parameter value into a corresponding use parameter value

If not `None`, *device* specifies the corresponding device instance for device-dependent conversion.

**using\_device** (*device*)

Context manager for temporarily changing the `device` attribute to the given device instance

```
class pylablib.core.devio.interface.EnumParameterClass (name, alias_map,
                                                         value_map=None, allowed_alias='device_values',
                                                         allowed_value='exact',
                                                         alias_case=None,
                                                         value_case=None,
                                                         match_prefix=False)
```

Bases: `pylablib.core.devio.interface.IEnumParameterClass`

Parameter class for a enum (i.e., predefined values) parameter with the specified mapping.

#### Parameters

- **name** – parameter class name
- **alias\_map** – mapping of aliases to device values; can be a dictionary, or a list of (*alias*, *value*) tuples (in the latter case non-tuple values are also allowed, indicating that value is the same as the alias); the list representation is useful in cases where the same alias maps to more than one value, so the map inversion is impossible
- **value\_map** – mapping of device values to aliases; can only be a dictionary or `None`, which means that the alias map is automatically inverted
- **allowed\_alias** – specifies a range of allowed aliases; can be "exact" (only exact map matches are allowed), "device\_value" (exact map matches and raw device values are allowed), or "all" (all values are allowed); in the latter two cases the value not in the map are passed as is.
- **allowed\_value** – specifies a range of allowed device values; can be "exact" (only exact map matches are allowed), or "all" (all values are allowed); in the latter case the value not in the map is passed as is.
- **alias\_case** – default alias parameter case for string values; can be `None` (no case normalization), or "lower" or "upper" (any received or returned alias will be normalized into this case)
- **value\_case** – default value parameter case for string values; can be `None` (no case normalization), or "lower" or "upper" (any received or returned device value will be normalized into this case)
- **match\_prefix** – if `True`, then the keys in the value map (or values in the alias map, if only it is provided) are assumed to be prefixes, so in the value-to-alias conversion the converted value matches the map value if it just starts with it; useful for some SCPI devices, where the shorter version of the value can sometimes be returned.

**check\_alias** (*alias*)

Check if the alias is valid

**check\_value** (*value*)

Check if the device value is valid

**docstring** ()

Get a parameter docstring

**i** (*value*, *device*=`None`)

Convert device parameter value into a corresponding use parameter value

If not `None`, *device* specifies the corresponding device instance for device-dependent conversion.

**to\_alias** (*value*)  
Convert the device value into an alias

**to\_value** (*alias*)  
Convert the alias into a device value

**using\_device** (*device*)  
Context manager for temporarily changing the *device* attribute to the given device instance

```
class pylablib.core.devio.interface.FunctionParameterClass (name, to_alias=None,
                                                         to_value=None,
                                                         check_value=None,
                                                         check_alias=None,
                                                         alias_err=None,
                                                         value_err=None)
```

Bases: `pylablib.core.devio.interface.ICheckingParameterClass`

Parameter class which uses supplied methods for checking, conversion, and generating error messages.

The arguments correspond to the parameter methods with the same names. When not supplied, checking methods always return `True`, conversion methods leave value intact, and error string methods generate the default error messages.

**check\_value** (*value*)  
Check if the device value is valid

**check\_alias** (*alias*)  
Check if the alias is valid

**to\_alias** (*value*)  
Convert the device value into an alias

**to\_value** (*alias*)  
Convert the alias into a device value

**docstring** ()  
Get a parameter docstring

**i** (*value*, *device*=None)  
Convert device parameter value into a corresponding use parameter value  
If not `None`, *device* specifies the corresponding device instance for device-dependent conversion.

**using\_device** (*device*)  
Context manager for temporarily changing the *device* attribute to the given device instance

```
class pylablib.core.devio.interface.CombinedParameterClass (name, parameters)
```

Bases: `pylablib.core.devio.interface.IParameterClass`

A multi-stage combined parameter class, which performs several conversion/check stages.

#### Parameters

- **name** – parameter class name
- **parameters** – list of parameters classes which are combined; the order is from the ‘most alias’ to the ‘most device parameter’, i.e., when converting an alias to a device parameter, it is first passed to the first class, then the second, etc. (the reverse is done when converting device values into aliases)

**docstring** ()  
Get a parameter docstring

**i** (*value*, *device=None*)

Convert device parameter value into a corresponding use parameter value

If not `None`, *device* specifies the corresponding device instance for device-dependent conversion.

**using\_device** (*device*)

Context manager for temporarily changing the `device` attribute to the given device instance

**class** `pylablib.core.devio.interface.TRawParameterValue` (*value*)

Bases: `tuple`

**count** ()

Return number of occurrences of value.

**index** ()

Return first index of value.

Raises `ValueError` if the value is not present.

**value**

`pylablib.core.devio.interface.pval` (*value*)

Mark that the value has already been treated by the parameter class

`pylablib.core.devio.interface.use_parameters` (\**args*, \*\**kwargs*)

Wrapper to indicate that a device class method uses device parameter classes.

The corresponding parameters classes are automatically determined if the argument name matches the parameter class name. The parameters classes can also be defined explicitly using keywords arguments `arg=parameter` supplied to the wrapper, where `arg` is the argument, and `parameter` is either a parameter class instance, or a parameter class name (the more preferable way). In addition, an argument `_returns` can be used to define the parameter class for the return value; it can also be a list or a tuple of parameter classes, indicating that the returned value is also a list or a tuple.

## Module contents

### pylablib.core.fileio package

#### Submodules

#### pylablib.core.fileio.datafile module

**class** `pylablib.core.fileio.datafile.DataFile` (*data*, *filepath=None*, *filetype=None*,  
*creation\_time=None*, *comments=None*,  
*props=None*)

Bases: `object`

Describes a single datafile.

#### Parameters

- **data** – the main content of the file (usually a numpy array, a pandas `DataFrame` or a `Dictionary`).
- **filepath** (*str*) – absolute path from which the file was read
- **filetype** (*str*) – a source type (e.g., "csv" or "bin")
- **creation\_time** (*datetime.datetime*) – File creation time
- **props** (*dict*) – all the metainfo about the file (extracted from comments, filename etc.)

- **comments** (*list*) – all the comments excluding the ones containing props

**get** (*name*, *default=None*)

Get a property from the dictionary. Use default value if it's not found

## pylablib.core.fileio.dict\_entry module

Classes for dealing with the *Dictionary* entries with special conversion rules when saved or loaded. Used to redefine how certain objects (e.g., tables) inside dictionaries are written into files and read from files.

`pylablib.core.fileio.dict_entry.is_dict_entry_branch(branch)`

Check if the dictionary branch contains a dictionary entry which needs to be specially converted.

**class** `pylablib.core.fileio.dict_entry.DictEntryBuilder` (*entry\_cls*, *pred=None*, *\*\*kwargs*)

Bases: `object`

Object for building dictionary entries from objects.

### Parameters

- **entry\_cls** – dictionary entry class
- **pred** – method used to check if an object can be turned into the corresponding entry; if None, use the default entry class checker (`entry_class.is_data_valid`)
- **kwargs** – keyword arguments passed to the entry constructor along with the data

**is\_data\_valid** (*data*)

Check if a data object can be wrapped by the current entry class

**from\_data** (*data*)

Build a dictionary entry from the data

**class** `pylablib.core.fileio.dict_entry.DictEntryParser` (*entry\_cls*, *pred=None*, *\*\*kwargs*)

Bases: `object`

Object for building dictionary entries from dictionary branches.

### Parameters

- **entry\_cls** – dictionary entry class
- **pred** – method used to check if a dictionary branch can be turned into the corresponding entry; if None, use the default entry class checker (`entry_class.is_branch_valid`)
- **kwargs** – keyword arguments passed to the entry `from_dict` class method along with the branch

**is\_branch\_valid** (*branch*)

Check if a branch can be parsed by the current entry class

**from\_dict** (*dict\_ptr*, *loc*)

Build a dictionary entry from the branch and the file location

`pylablib.core.fileio.dict_entry.add_dict_entry_builder(builder)`

Add an entry builder to the global list of builders

`pylablib.core.fileio.dict_entry.add_dict_entry_parser(parser)`

Add an entry parser to the global list of parsers

`pylablib.core.fileio.dict_entry.add_dict_entry_class(cls)`

Add an entry class.

Automatically registers builder and parser, which take no additional arguments and use default class method to determine if an object/branch can be converted into an entry.

`pylablib.core.fileio.dict_entry.from_data(data, builders=None)`

Build a dictionary entry from the data.

*builders* can contain an additional list of builder to try before using the default ones.

`pylablib.core.fileio.dict_entry.from_dict(dict_ptr, loc, parsers=None)`

Build a dictionary entry from the dictionary branch and the file location.

*parsers* can contain an additional list of parsers to try before using the default ones.

**class** `pylablib.core.fileio.dict_entry.IDictionaryEntry(data)`

Bases: `object`

A generic *Dictionary* entry.

Contains data represented by the node, as well as the way to represent this data as a dictionary branch.

**Parameters** *data* – data to be wrapped

**is\_data\_valid** (*class method*)

check if a data object can be wrapped by the current entry class

**is\_branch\_valid** (*class method*)

check if a branch can be parsed by the current entry class

**from\_dict** (*class method*)

create a dictionary entry of a given class from the dictionary branch

**to\_dict** ()

convert the entry to a dictionary branch

**classmethod is\_data\_valid** (*data*)

Check if a data object can be wrapped by the current entry class

**classmethod is\_branch\_valid** (*branch*)

Check if a branch can be parsed by the current entry class

**classmethod from\_dict** (*dict\_ptr, loc*)

Convert a dictionary branch to a specific *IDictionaryEntry* object.

**Parameters**

- **dict\_ptr** (`dictionary.DictionaryPointer`) – Pointer to the dictionary location for the entry.
- **loc** – Location for the data to be loaded.

**to\_dict** (*dict\_ptr, loc*)

Convert data to a dictionary branch on saving.

**Parameters**

- **dict\_ptr** (`dictionary.DictionaryPointer`) – Pointer to the dictionary location for the entry.
- **loc** – File location for the data to be saved.

`pylablib.core.fileio.dict_entry.parse_stored_table_data(desc=None, data=None, out_type='pandas')`

Parse table data corresponding to the given description dictionary and data.

**Parameters**

- **desc** – description dictionary; can be None, if no description is given
- **data** – separately loaded data; can be None, if no data is given (in this case assume that it is stored in the description dictionary); can be a tuple (column\_data, column\_names) (such as the one returned by `parse_csv.read_table()`), or a an *InlineTable* object containing such tuple.
- **out\_type** (*str*) – Output format of the data ('array' for numpy arrays or 'pandas' for pandas DataFrame objects).

**Returns** tuple (data, columns), where data is the data table in the specified format, and columns is the list of columns

```
class pylablib.core.fileio.dict_entry.ITableDictionaryEntry (data,
                                                             columns=None)
```

Bases: `pylablib.core.fileio.dict_entry.IDictionaryEntry`

A generic table Dictionary entry.

**Parameters**

- **data** – Table data.
- **columns** (*list*) – If not None, list of column names (if None and data is a pandas DataFrame object, get column names from that).

```
classmethod is_data_valid (data)
```

Check if a data object can be wrapped by the current entry class

```
classmethod from_dict (dict_ptr, loc, out_type='pandas')
```

Convert a dictionary branch to a specific DictionaryEntry object.

**Parameters**

- **dict\_ptr** (`dictionary.DictionaryPointer`) – Pointer to the dictionary location for the entry.
- **loc** – Location for the data to be loaded.
- **out\_type** (*str*) – Output format of the data ('array' for numpy arrays or 'pandas' for pandas DataFrame objects), used only if the dictionary doesn't provide the format.

```
classmethod is_branch_valid (branch)
```

Check if a branch can be parsed by the current entry class

```
to_dict (dict_ptr, loc)
```

Convert data to a dictionary branch on saving.

**Parameters**

- **dict\_ptr** (`dictionary.DictionaryPointer`) – Pointer to the dictionary location for the entry.
- **loc** – File location for the data to be saved.

```
class pylablib.core.fileio.dict_entry.InlineTableDictionaryEntry (data,
                                                                    columns=None)
```

Bases: `pylablib.core.fileio.dict_entry.ITableDictionaryEntry`

An inlined table Dictionary entry.

**Parameters**

- **data** – Table data.
- **columns** (*list*) – If not `None`, a list of column names (if `None` and data is a pandas DataFrame object, get column names from that).

**to\_dict** (*dict\_ptr*, *loc*)

Convert the data to a dictionary branch and write the table to the file.

**classmethod from\_dict** (*dict\_ptr*, *loc*, *out\_type*='pandas')

Build an *InlineTableDictionaryEntry* object from the dictionary and read the inlined data.

#### Parameters

- **dict\_ptr** (*dictionary.DictionaryPointer*) – Pointer to the dictionary location for the entry.
- **loc** – Location for the data to be loaded.
- **out\_type** (*str*) – Output format of the data ('array' for numpy arrays or 'pandas' for pandas DataFrame objects).

**classmethod is\_branch\_valid** (*branch*)

Check if a branch can be parsed by the current entry class

**classmethod is\_data\_valid** (*data*)

Check if a data object can be wrapped by the current entry class

```
class pylablib.core.fileio.dict_entry.IExternalTableDictionaryEntry (data,  
                                                                    file_format,  
                                                                    name,  
                                                                    columns,  
                                                                    force_name=True)
```

Bases: *pylablib.core.fileio.dict\_entry.ITableDictionaryEntry*

**classmethod from\_dict** (*dict\_ptr*, *loc*, *out\_type*='pandas')

Convert a dictionary branch to a specific DictionaryEntry object.

#### Parameters

- **dict\_ptr** (*dictionary.DictionaryPointer*) – Pointer to the dictionary location for the entry.
- **loc** – Location for the data to be loaded.
- **out\_type** (*str*) – Output format of the data ('array' for numpy arrays or 'pandas' for pandas DataFrame objects), used only if the dictionary doesn't provide the format.

**classmethod is\_branch\_valid** (*branch*)

Check if a branch can be parsed by the current entry class

**classmethod is\_data\_valid** (*data*)

Check if a data object can be wrapped by the current entry class

**to\_dict** (*dict\_ptr*, *loc*)

Convert data to a dictionary branch on saving.

#### Parameters

- **dict\_ptr** (*dictionary.DictionaryPointer*) – Pointer to the dictionary location for the entry.
- **loc** – File location for the data to be saved.



```
class pylablib.core.fileio.dict_entry.ExternalTextTableDictionaryEntry (data=None,
                                                                    file_format='csv',
                                                                    name="",
                                                                    columns=None,
                                                                    force_name=True)
```

Bases: `pylablib.core.fileio.dict_entry.IExternalTableDictionaryEntry`

An external text table Dictionary entry.

#### Parameters

- **data** – Table data.
- **file\_format** (*str*) – Output file format.
- **name** (*str*) – Name template for the external file (default is the full path connected with "\_" symbol).
- **columns** (*list*) – If not None, a list of column names (if None and data is a pandas DataFrame object, get column names from that).
- **force\_name** (*bool*) – If False and the target file already exists, generate a new unique name; otherwise, overwrite the file.

```
to_dict (dict_ptr, loc)
```

Convert the data to a dictionary branch and save the table to an external file.

```
classmethod from_dict (dict_ptr, loc, out_type='pandas')
```

Build an `ExternalTextTableDictionaryEntry` object from the dictionary and load the external data.

#### Parameters

- **dict\_ptr** (`dictionary.DictionaryPointer`) – Pointer to the dictionary location for the entry.
- **loc** – Location for the data to be loaded.
- **out\_type** (*str*) – Output format of the data ('array' for numpy arrays or 'pandas' for pandas DataFrame objects).

```
classmethod is_branch_valid (branch)
```

Check if a branch can be parsed by the current entry class

```
classmethod is_data_valid (data)
```

Check if a data object can be wrapped by the current entry class

```
class pylablib.core.fileio.dict_entry.ExternalBinTableDictionaryEntry (data=None,
                                                                    file_format='bin',
                                                                    name="",
                                                                    columns=None,
                                                                    force_name=True)
```

Bases: `pylablib.core.fileio.dict_entry.IExternalTableDictionaryEntry`

An external binary table Dictionary entry.

#### Parameters

- **data** – Table data.
- **file\_format** (*str*) – Output file format.
- **name** (*str*) – Name template for the external file (default is the full path connected with "\_" symbol).

- **columns** (*list*) – If not `None`, a list of column names (if `None` and data is a pandas DataFrame object, get column names from that).
- **force\_name** (*bool*) – If `False` and the target file already exists, generate a new unique name; otherwise, overwrite the file.

**to\_dict** (*dict\_ptr*, *loc*)

Convert the data to a dictionary branch and save the table to an external file.

**classmethod from\_dict** (*dict\_ptr*, *loc*, *out\_type*='pandas')

Build an `ExternalBinTableDictionaryEntry` object from the dictionary and load the external data.

#### Parameters

- **dict\_ptr** (`dictionary.DictionaryPointer`) – Pointer to the dictionary location for the entry.
- **loc** – Location for the data to be loaded.
- **out\_type** (*str*) – Output format of the data ('array' for numpy arrays or 'pandas' for pandas DataFrame objects).

**classmethod is\_branch\_valid** (*branch*)

Check if a branch can be parsed by the current entry class

**classmethod is\_data\_valid** (*data*)

Check if a data object can be wrapped by the current entry class

`pylablib.core.fileio.dict_entry.table_entry_builder` (*table\_format*='inline')

Make an entry builder for tables depending on the table format.

**Parameters table\_format** (*str*) – Default format for table (numpy arrays or pandas DataFrames) entries. Can be 'inline' (table is written inside the file), 'csv' (external CSV file) or 'bin' (external binary file).

**class** `pylablib.core.fileio.dict_entry.IExternalFileDictionaryEntry` (*data*,  
name="",  
force\_name=True)

Bases: `pylablib.core.fileio.dict_entry.IDictionaryEntry`

Generic dictionary entry for data in an external file.

#### Parameters

- **data** – Stored data.
- **name** (*str*) – Name template for the external file (default is the full path connected with "\_" symbol).
- **force\_name** (*bool*) – If `False` and the target file already exists, generate a new unique name; otherwise, overwrite the file.

**file\_format** = `None`

**static add\_file\_format** (*subclass*)

Register an `IExternalFileDictionaryEntry` as a possible stored file format.

Used to automatically invoke a correct loader when loading the dictionary file. Only needs to be done once after the subclass declaration.

**to\_dict** (*dict\_ptr*, *loc*)

Convert the data to a dictionary branch and save the data to an external file

**classmethod** `from_dict(dict_ptr, loc)`

Build an *IEternalFileDictionaryEntry* object from the dictionary and load the external data.

#### Parameters

- **dict\_ptr** (*dictionary.DictionaryPointer*) – Pointer to the dictionary location for the entry.
- **loc** – Location for the data to be loaded.

**get\_preamble()**

Generate preamble (dictionary with supplementary data which allows to load the data from the file)

**save\_file(location\_file)**

Save stored data into the given location.

Virtual method, should be overloaded in subclasses

**classmethod load\_file(location\_file, preamble)**

Load stored data from the given location, using the supplied preamble.

Virtual method, should be overloaded in subclasses

**classmethod is\_branch\_valid(branch)**

Check if a branch can be parsed by the current entry class

**classmethod is\_data\_valid(data)**

Check if a data object can be wrapped by the current entry class

```
class pylablib.core.fileio.dict_entry.ExternalNumpyDictionaryEntry(data,
                                                                    name="",
                                                                    force_name=True,
                                                                    dtype=None)
```

Bases: *pylablib.core.fileio.dict\_entry.IEternalFileDictionaryEntry*

A dictionary entry which stores the numpy array data into an external file in binary format.

#### Parameters

- **data** – Numpy array data.
- **name** (*str*) – Name template for the external file (default is the full path connected with "\_" symbol).
- **force\_name** (*bool*) – If *False* and the target file already exists, generate a new unique name; otherwise, overwrite the file.
- **dtype** – numpy dtype to load/save the data (by default, dtype of the supplied data).

**file\_format = 'numpy'**

**get\_preamble()**

Generate preamble (dictionary with supplementary data which allows to load the data from the file)

**save\_file(location\_file)**

Save stored data into the given location

**classmethod load\_file(location\_file, preamble)**

Load stored data from the given location, using the supplied preamble

**static add\_file\_format(subclass)**

Register an *IEternalFileDictionaryEntry* as a possible stored file format.

Used to automatically invoke a correct loader when loading the dictionary file. Only needs to be done once after the subclass declaration.

**classmethod** `from_dict(dict_ptr, loc)`

Build an *IExternalFileDictionaryEntry* object from the dictionary and load the external data.

**Parameters**

- **dict\_ptr** (*dictionary.DictionaryPointer*) – Pointer to the dictionary location for the entry.
- **loc** – Location for the data to be loaded.

**classmethod** `is_branch_valid(branch)`

Check if a branch can be parsed by the current entry class

**classmethod** `is_data_valid(data)`

Check if a data object can be wrapped by the current entry class

**to\_dict** (*dict\_ptr, loc*)

Convert the data to a dictionary branch and save the data to an external file

**class** `pylablib.core.fileio.dict_entry.ExpandedContainerDictionaryEntry(data)`

Bases: *pylablib.core.fileio.dict\_entry.IDictionaryEntry*

A dictionary entry which expands containers (lists, tuples, dictionaries) into subdictionaries.

Useful when the data in the containers is complex, so writing it into one line (as is default for lists and tuples) wouldn't work.

**Parameters** **data** – Container data.

**to\_dict** (*dict\_ptr, loc*)

Convert the stored container to a dictionary branch

**classmethod** `from_dict(dict_ptr, loc)`

Build an *ExpandedContainerDictionaryEntry* object from the dictionary

**classmethod** `is_branch_valid(branch)`

Check if a branch can be parsed by the current entry class

**classmethod** `is_data_valid(data)`

Check if a data object can be wrapped by the current entry class

## pylablib.core.fileio.loadfile module

Utilities for reading data files.

**class** `pylablib.core.fileio.loadfile.IInputFileFormat`

Bases: *object*

Generic class for an input file format.

Based on *file\_format* or autodetection, calls one of its subclasses to read the file.

Defines a single static method

**static** `detect_file_format(location_file)`

**read** (*location\_file*)

Read a file at a given location

**class** `pylablib.core.fileio.loadfile.ITextInputFileFormat`

Bases: *pylablib.core.fileio.loadfile.IInputFileFormat*

Generic class for a text input file format.

Based on *file\_format* or autodetection, calls one of its subclasses to read the file.

**static detect\_file\_format** (*location\_file*)

**read** (*location\_file*)

Read a file at a given location

```
class pylablib.core.fileio.loadfile.CSVTableInputFileFormat (out_type='default',
                                                             dtype='numeric',
                                                             columns=None,
                                                             delimiters=None,
                                                             empty_entry_substitute=None,
                                                             ig-
                                                             nore_corrupted_lines=True,
                                                             skip_lines=0)
```

Bases: *pylablib.core.fileio.loadfile.ITextInputFileFormat*

Class for CSV input file format.

#### Parameters

- **out\_type** (*str*) – type of the result: 'array' for numpy array, 'pandas' for pandas DataFrame, or 'default' (determined by the library default; 'pandas' by default)
- **dtype** – dtype of entries; can be either a single type, or a list of types (one per column). Possible dtypes are: 'int', 'float', 'complex', 'numeric' (tries to coerce to minimal possible numeric type, raises error if data can't be converted to *complex*), 'generic' (accept arbitrary types, including lists, dictionaries, escaped strings, etc.), 'raw' (keep raw string).
- **columns** – either a number if columns, or a list of columns names.
- **delimiters** (*str*) – Regex string which recognizes entries delimiters (by default `r"\s*,\s*|\s+"`, i.e., commas and whitespaces).
- **empty\_entry\_substitute** – Substitute for empty table entries. If None, all empty table entries are skipped.
- **ignore\_corrupted\_lines** (*bool*) – If True, skip corrupted (e.g., non-numeric for numeric dtype, or with too few entries) lines; otherwise, raise *ValueError*.
- **skip\_lines** (*int*) – Number of lines to skip from the beginning of the file.

**read** (*location\_file*)

Read a file at a given location

**static detect\_file\_format** (*location\_file*)

```
class pylablib.core.fileio.loadfile.DictionaryInputFileFormat (case_normalization=None,
                                                                in-
                                                                line_dtype='generic',
                                                                in-
                                                                line_out_type='default',
                                                                en-
                                                                try_format='value',
                                                                al-
                                                                low_duplicate_keys=False,
                                                                skip_lines=0)
```

Bases: *pylablib.core.fileio.loadfile.ITextInputFileFormat*

Class for Dictionary input file format.

### Parameters

- **location\_file** – Location of the data.
- **case\_normalization** (*str*) – If *None*, the dictionary paths are case-sensitive; otherwise, defines the way the entries are normalized ('lower' or 'upper').
- **inline\_dtype** (*str*) – dtype for inlined tables.
- **inline\_out\_type** (*str*) – type of the result of the inline table: 'array' for numpy array, 'pandas' for pandas DataFrame, 'raw' for raw *InlineTable* data containing tuple (column\_data, column\_names), or 'default' (determined by the library default; 'pandas' by default).
- **entry\_format** (*str*) – Determines the way for dealing with *dict\_entry.IDictionaryEntry* objects (objects transformed into dictionary branches with special recognition rules). Can be 'branch' (don't attempt to recognize those object, leave dictionary as in the file), 'dict\_entry' (recognize and leave as *dict\_entry.IDictionaryEntry* objects) or 'value' (recognize and keep the value).
- **allow\_duplicate\_keys** (*bool*) – if *False* and the same key is mentioned twice in the file, raise an error
- **skip\_lines** (*int*) – Number of lines to skip from the beginning of the file.

**read** (*location\_file*)

Read a file at a given location

**static detect\_file\_format** (*location\_file*)

```
class pylablib.core.fileio.loadfile.BinaryTableInputFileFormatter (out_type='default',  
                                                                    dtype='<f8',  
                                                                    columns=None,  
                                                                    packing='flatten',  
                                                                    preamble=None,  
                                                                    skip_bytes=0)
```

Bases: *pylablib.core.fileio.loadfile.IInputFileFormat*

Class for binary input file format.

### Parameters

- **location\_file** – Location of the data.
- **out\_type** (*str*) – type of the result: 'array' for numpy array, 'pandas' for pandas DataFrame, or 'default' (determined by the library default; 'pandas' by default)
- **dtype** – *numpy.dtype* describing the data.
- **columns** – either number of columns, or a list of columns names.
- **packing** (*str*) – The way the 2D array is packed. Can be either 'flatten' (data is stored row-wise) or 'transposed' (data is stored column-wise).
- **preamble** (*dict*) – If not *None*, defines binary file parameters that supersede the parameters supplied to the function. The defined parameters are 'dtype', 'packing', 'ncols' (number of columns) and 'nrows' (number of rows).
- **skip\_bytes** (*int*) – Number of bytes to skip from the beginning of the file.

**read** (*location\_file*)

Read a file at a given location

**static detect\_file\_format** (*location\_file*)

`pylablib.core.fileio.loadfile.build_file_format` (*location\_file*, *file\_format*='generic',  
\*\**kwargs*)

Create file format (*IInputFileFormat* instance) for given parameters and file locations.

If *file\_format* is already an instance of *IInputFileFormat*, return unchanged. If *file\_format* is generic (e.g., "generic" or "test"), attempt to autodetect it from the file. \*\**kwargs* are passed to the file format constructor.

`pylablib.core.fileio.loadfile.load_csv` (*path*=None, *out\_type*='default',  
*dtype*='numeric', *columns*=None, *delim-*  
*iters*=None, *empty\_entry\_substitute*=None,  
*ignore\_corrupted\_lines*=True, *skip\_lines*=0,  
*loc*='file', *encoding*=None, *return\_file*=False)

Load data table from a CSV/table file.

#### Parameters

- **path** (*str*) – path to the file of a file-like object
- **out\_type** (*str*) – type of the result: 'array' for numpy array, 'pandas' for pandas DataFrame, or 'default' (determined by the library default; 'pandas' by default)
- **dtype** – dtype of entries; can be either a single type, or a list of types (one per column). Possible dtypes are: 'int', 'float', 'complex', 'numeric' (tries to coerce to minimal possible numeric type, raises error if data can't be converted to *complex*), 'generic' (accept arbitrary types, including lists, dictionaries, escaped strings, etc.), 'raw' (keep raw string).
- **columns** – either a number if columns, or a list of columns names
- **delimiters** (*str*) – regex string which recognizes entries delimiters (by default `r"\s*,\s*|\s+"`, i.e., commas and whitespaces)
- **empty\_entry\_substitute** – substitute for empty table entries. If None, all empty table entries are skipped
- **ignore\_corrupted\_lines** (*bool*) – if True, skip corrupted (e.g., non-numeric for numeric dtype, or with too few entries) lines; otherwise, raise *ValueError*
- **skip\_lines** (*int*) – number of lines to skip from the beginning of the file
- **loc** (*str*) – location type ("file" means the usual file location; see `location.get_location()` for details)
- **encoding** – if a new file location is opened, this specifies the encoding
- **return\_file** (*bool*) – if True, return *DataFile* object (contains some metainfo); otherwise, return just the file data

`pylablib.core.fileio.loadfile.load_csv_desc` (*path*=None, *loc*='file', *encoding*=None, *re-*  
*turn\_file*=False)

Load data from the extended CSV table file.

Analogous to `load_dict()`, but doesn't allow any additional parameters (which don't matter in this case).

#### Parameters

- **path** (*str*) – path to the file of a file-like object

- **loc** (*str*) – location type ("file" means the usual file location; see `location.get_location()` for details)
- **encoding** – if a new file location is opened, this specifies the encoding
- **return\_file** (*bool*) – if True, return `DataFile` object (contains some metainfo); otherwise, return just the file data

```
pylablib.core.fileio.loadfile.load_bin(path=None, out_type='default', dtype='<f8',
                                       columns=None, packing='flatten', preamble=None,
                                       skip_bytes=0, loc='file', encoding=None, re-
                                       turn_file=False)
```

Load data from the binary file.

#### Parameters

- **path** (*str*) – path to the file of a file-like object
- **out\_type** (*str*) – type of the result: 'array' for numpy array, 'pandas' for pandas DataFrame, or 'default' (determined by the library default; 'pandas' by default)
- **dtype** – `numpy.dtype` describing the data.
- **columns** – either number of columns, or a list of columns names.
- **packing** (*str*) – The way the 2D array is packed. Can be either 'flatten' (data is stored row-wise) or 'transposed' (data is stored column-wise).
- **preamble** (*dict*) – If not None, defines binary file parameters that supersede the parameters supplied to the function. The defined parameters are 'dtype', 'packing', 'ncols' (number of columns) and 'nrows' (number of rows).
- **skip\_bytes** (*int*) – Number of bytes to skip from the beginning of the file.
- **loc** (*str*) – location type ("file" means the usual file location; see `location.get_location()` for details)
- **encoding** – if a new file location is opened, this specifies the encoding
- **return\_file** (*bool*) – if True, return `DataFile` object (contains some metainfo); otherwise, return just the file data

```
pylablib.core.fileio.loadfile.load_bin_desc(path=None, loc='file', encoding=None, re-
                                             turn_file=False)
```

Load data from the binary file with a description.

Analogous to `load_dict()`, but doesn't allow any additional parameters (which don't matter in this case).

#### Parameters

- **path** (*str*) – path to the file of a file-like object
- **loc** (*str*) – location type ("file" means the usual file location; see `location.get_location()` for details)
- **encoding** – if a new file location is opened, this specifies the encoding
- **return\_file** (*bool*) – if True, return `DataFile` object (contains some metainfo); otherwise, return just the file data



```
pylablib.core.fileio.loadfile.load_dict (path=None, case_normalization=None, inline_dtype='generic', entry_format='value', inline_out_type='default', skip_lines=0, allow_duplicate_keys=False, loc='file', encoding=None, return_file=False)
```

Load data from the dictionary file.

#### Parameters

- **path** (*str*) – path to the file of a file-like object
- **case\_normalization** (*str*) – If *None*, the dictionary paths are case-sensitive; otherwise, defines the way the entries are normalized ('lower' or 'upper').
- **inline\_dtype** (*str*) – dtype for inlined tables.
- **inline\_out\_type** (*str*) – type of the result of the inline table: 'array' for numpy array, 'pandas' for pandas DataFrame, 'raw' for raw *InlineTable* data containing tuple (column\_data, column\_names), or 'default' (determined by the library default; 'pandas' by default).
- **entry\_format** (*str*) – Determines the way for dealing with *dict\_entry.IDictionaryEntry* objects (objects transformed into dictionary branches with special recognition rules). Can be 'branch' (don't attempt to recognize those object, leave dictionary as in the file), 'dict\_entry' (recognize and leave as *dict\_entry.IDictionaryEntry* objects) or 'value' (recognize and keep the value).
- **allow\_duplicate\_keys** (*bool*) – if *False* and the same key is mentioned twice in the file, raise an error
- **skip\_lines** (*int*) – Number of lines to skip from the beginning of the file.
- **loc** (*str*) – location type ("file" means the usual file location; see *location.get\_location()* for details)
- **encoding** – if a new file location is opened, this specifies the encoding
- **return\_file** (*bool*) – if *True*, return *DataFile* object (contains some meta-info); otherwise, return just the file data

```
pylablib.core.fileio.loadfile.load_generic (path=None, file_format=None, loc='file', encoding=None, return_file=False, **kwargs)
```

Load data from the file.

#### Parameters

- **path** (*str*) – path to the file of a file-like object
- **file\_format** (*str*) – input file format; if *None*, attempt to auto-detect file format (same as 'generic'); can also be an *IInputFileFormat* instance for specific reading method
- **loc** (*str*) – location type ("file" means the usual file location; see *location.get\_location()* for details)
- **encoding** – if a new file location is opened, this specifies the encoding
- **return\_file** (*bool*) – if *True*, return *DataFile* object (contains some meta-info); otherwise, return just the file data

*\*\*kwargs* are passed to the file formatter used to read the data (see *CSVTableInputFileFormat*, *DictionaryInputFileFormat* and *BinaryTableInputFileFormatter* for the possible arguments). The default format names are:

- 'generic': Generic file format. Attempt to autodetect, raise `IOError` if unsuccessful;
- 'txt': Generic text file. Attempt to autodetect, raise `IOError` if unsuccessful
- 'csv': CSV file, corresponds to `CSVTableInputFileFormat`;
- 'dict': Dictionary file, corresponds to `DictionaryInputFileFormat`;
- 'bin': Binary file, corresponds to `BinaryTableInputFileFormatter`

## pylablib.core.fileio.loadfile\_utils module

Miscellaneous utilities for reading data files.

`pylablib.core.fileio.loadfile_utils.is_unprintable_character(ch)`

`pylablib.core.fileio.loadfile_utils.detect_binary_file(stream)`

Check if the opened file is binary

`pylablib.core.fileio.loadfile_utils.test_row_type(line)`

Try to determine whether the line is a comment line, a numerical data row, a dictionary row or an unrecognized row.

Doesn't distinguish with a great accuracy; useful only for trying to guess file format.

`pylablib.core.fileio.loadfile_utils.detect_textfile_type(stream)`

Try to autodetect text file type: dictionary or table

`pylablib.core.fileio.loadfile_utils.test_savetime_comment(line)`

Test if the comment resembles a savetime line

`pylablib.core.fileio.loadfile_utils.find_savetime_comment(comments)`

Try to find savetime comment

`pylablib.core.fileio.loadfile_utils.test_columns_line(line, cols_num)`

Test if the line looks like a list of columns for a given columns number

`pylablib.core.fileio.loadfile_utils.find_columns_lines(corrupted, comments, cols_num)`

Try to find a column line (for a given columns number) among the comment and corrupted lines

`class pylablib.core.fileio.loadfile_utils.InlineTable(table)`

Bases: `object`

Simple marker class that denotes that the wrapped numpy 2D array should be written inline

`pylablib.core.fileio.loadfile_utils.parse_dict_line(line)`

Parse stripped dictionary file line

`pylablib.core.fileio.loadfile_utils.read_dict_and_comments(f, case_normalization=None, in-line_dtype='generic', allow_duplicate_keys=False)`

Load dictionary entries and comments from the file stream.

### Parameters

- **f** – file stream
- **case\_normalization** – case normalization for the returned dictionary; `None` means that it's case sensitive, "upper" and "lower" determine how they are normalized

- **inline\_dtype** – dtype for inline tables; by default, use the most generic type (can include Python objects such as lists or strings)
- **allow\_duplicate\_keys** – if `False` and the same key is listed twice, raise an error

Return tuple (`data`, `comment_lines`), where `data` is a dictionary with parsed entries (tables are still represented as ‘raw’, i.e., as a tuple of columns list and column names list), and `comment_lines` is a list of comment lines

## pylablib.core.fileio.location module

Classes for describing a generic file location.

**class** `pylablib.core.fileio.location.LocationName` (`path=None`, `ext=None`)

Bases: `object`

File name inside a location.

### Parameters

- **path** – Path inside the location. Gets normalized according to the Dictionary rules (not case-sensitive; ‘/’ and ‘\’ are the delimiters).
- **ext** (`str`) – Name extension (`None` is default).

**get\_path** (`default_path=""`, `sep='/'`)

Get the string path.

If the object’s `path` is `None`, use `default_path` instead. If `sep` is not `None`, use it to join the path entries; otherwise, return the path in a list form.

**get\_ext** (`default_ext=""`)

Get the extension.

If the object’s `ext` is `None`, use `default_ext` instead.

**to\_string** (`default_path=""`, `default_ext=""`, `path_sep='/'`, `ext_sep='|'`, `add_empty_ext=True`)

Convert the path to a string representation.

### Parameters

- **default\_path** (`str`) – Use it as path if the object’s `path` is `None`.
- **default\_ext** (`str`) – Use it as path if the object’s `ext` is `None`.
- **path\_sep** (`str`) – Use it to join the path entries.
- **ext\_sep** (`str`) – Use it to join path and extension.
- **add\_empty\_ext** (`str`) – If `False` and the extension is empty, don’t add `ext_sep` in the end.

**to\_path** (`default_path=""`, `default_ext=""`, `ext_sep='|'`, `add_empty_ext=True`)

Convert the path to a list representation.

Extension is added with `ext_sep` to the last entry in the path.

### Parameters

- **default\_path** (`str`) – Use it as path if the object’s `path` is `None`.
- **default\_ext** (`str`) – Use it as path if the object’s `ext` is `None`.
- **ext\_sep** (`str`) – Use it to join path and extension.

- **add\_empty\_ext** (*str*) – If `False` and the extension is empty, don't add *ext\_sep* in the end.

**static from\_string** (*expr*, *ext\_sep*='|')

Create a *LocationName* object from a string representation.

*ext\_sep* defines extension separator; the path separators are '/' and '\'. Empty path or extension translate into `None`.

**static from\_object** (*obj*)

Create a *LocationName* object from an object.

*obj* can be a *LocationName* (return unchanged), tuple or list (use as construct arguments), string (treat as a string representation) or `None` (return empty name).

**copy** ()

**class** `pylablib.core.fileio.location.LocationFile` (*loc*, *name*=`None`)

Bases: `object`

A file at a location.

Combines information about the location and the name within this location. Can be opened for reading or writing.

#### Parameters

- **loc** – File location.
- **name** – File's name inside the location.

**loc**

File location.

**name**

File's name inside the location.

**opened**

Whether the file is currently opened.

**open** (*mode*='read', *data\_type*='text')

Open the file.

#### Parameters

- **mode** (*str*) – Opening mode. Can be 'read', 'write' or 'append', as well as standard abbreviation (e.g., "r" or "wb").
- **data\_type** (*str*) – Either 'text' or 'binary'; if *mode* is an abbreviation, this parameter is ignored (i.e., `open("r", "binary")` still opens file as text).

**close** ()

Close the file

**class** `pylablib.core.fileio.location.IDataLocation`

Bases: `object`

Generic location.

**is\_free** (*name*=`None`)

Check if the name is unoccupied

**generate\_new\_name** (*prefix\_name*, *idx*=0)

Generate a new name inside the location using the given prefix and starting index.

If *idx* is `None`, check just the *prefix\_name* first before starting to append indices.

**open** (*name=None, mode='read', data\_type='text'*)

Open a location file.

#### Parameters

- **name** – File name inside the location (None means 'default' location),
- **mode** (*str*) – Opening mode. Can be 'read', 'write' or 'append', as well as standard abbreviation (e.g., "r" or "wb").
- **data\_type** (*str*) – Either 'text' or 'binary'; if *mode* is an abbreviation, this parameter is ignored (i.e., open ("r", "binary") still opens file as text).

**close** (*name*)

Close a location file.

**list\_opened\_files** ()

Get a dictionary {string\_name: location\_file} of all files opened in this location

```
class pylablib.core.fileio.location.OpenedFileLocation (f, open_error=False,
check_mode=False,
check_data_type=True)
```

Bases: `object`

File location which corresponds to an already opened file.

**is\_free** (*name=None*)

**generate\_new\_name** (*prefix\_name, idx=0*)

**open** (*name=None, mode='read', data\_type='text'*)

**close** (*name*)

**list\_opened\_files** ()

```
class pylablib.core.fileio.location.IFileSystemDataLocation (encoding=None)
```

Bases: `pylablib.core.fileio.location.IDataLocation`

A generic filesystem data location.

A single file name describes a single file in the filesystem.

**get\_filesystem\_path** (*name=None, path\_type='absolute'*)

Get the filesystem path corresponding to a given name.

*path\_type* can be 'absolute' (return absolute path), 'relative' (return relative path; level depends on the location) or 'name' (only return path inside the location).

**is\_free** (*name=None*)

Check if the name is unoccupied

**open** (*name=None, mode='read', data\_type='text'*)

Open a location file.

#### Parameters

- **name** – File name inside the location (None means 'default' location),
- **mode** (*str*) – Opening mode. Can be 'read', 'write' or 'append', as well as standard abbreviation (e.g., "r" or "wb").
- **data\_type** (*str*) – Either 'text' or 'binary'; if *mode* is an abbreviation, this parameter is ignored (i.e., open ("r", "binary") still opens file as text).

**close** (*name*)

Close a location file

**list\_opened\_files()**

Get a dictionary {string\_name: location\_file} of all files opened in this location

**generate\_new\_name** (prefix\_name, idx=0)

Generate a new name inside the location using the given prefix and starting index.

If *idx* is None, check just the *prefix\_name* first before starting to append indices.

**class** pylablib.core.fileio.location.**SingleFileSystemDataLocation** (file\_path,  
encod-  
ing=None)

Bases: `pylablib.core.fileio.location.IFileSystemDataLocation`

A location describing a single file.

Any use of a non-default name raises `ValueError`.

**Parameters** *file\_path* (*str*) – The path to the file.

**get\_filesystem\_path** (name=None, path\_type='absolute')

Get the filesystem path corresponding to a given name.

*path\_type* can be 'absolute' (return absolute path), 'relative' (return relative path; level depends on the location) or 'name' (only return path inside the location).

**close** (name)

Close a location file

**generate\_new\_name** (prefix\_name, idx=0)

Generate a new name inside the location using the given prefix and starting index.

If *idx* is None, check just the *prefix\_name* first before starting to append indices.

**is\_free** (name=None)

Check if the name is unoccupied

**list\_opened\_files()**

Get a dictionary {string\_name: location\_file} of all files opened in this location

**open** (name=None, mode='read', data\_type='text')

Open a location file.

#### Parameters

- **name** – File name inside the location (None means 'default' location),
- **mode** (*str*) – Opening mode. Can be 'read', 'write' or 'append', as well as standard abbreviation (e.g., "r" or "wb").
- **data\_type** (*str*) – Either 'text' or 'binary'; if *mode* is an abbreviation, this parameter is ignored (i.e., open ("r", "binary") still opens file as text).

**class** pylablib.core.fileio.location.**PrefixedFileSystemDataLocation** (file\_path,  
pre-  
fix\_template='{0}\_{1}',  
encod-  
ing=None)

Bases: `pylablib.core.fileio.location.IFileSystemDataLocation`

A location describing a set of prefixed files.

#### Parameters

- **file\_path** (*str*) – A master path. Its name is used as a prefix, and its extension is used as a default.

- **prefix\_template** (*str*) – A `str.format()` string for generating prefixed files. Has two arguments: the first is the master name, the second is the sub\_location.

Multi-level paths translate into nested folders (the top level folder is combined from the *file\_path* prefix and the first path entry).

**get\_filesystem\_path** (*name=None, path\_type='absolute'*)

Get the filesystem path corresponding to a given name.

*path\_type* can be 'absolute' (return absolute path), 'relative' (return relative path; level depends on the location) or 'name' (only return path inside the location).

**close** (*name*)

Close a location file

**generate\_new\_name** (*prefix\_name, idx=0*)

Generate a new name inside the location using the given prefix and starting index.

If *idx* is None, check just the *prefix\_name* first before starting to append indices.

**is\_free** (*name=None*)

Check if the name is unoccupied

**list\_opened\_files** ()

Get a dictionary {string\_name: location\_file} of all files opened in this location

**open** (*name=None, mode='read', data\_type='text'*)

Open a location file.

#### Parameters

- **name** – File name inside the location (None means 'default' location),
- **mode** (*str*) – Opening mode. Can be 'read', 'write' or 'append', as well as standard abbreviation (e.g., "r" or "wb").
- **data\_type** (*str*) – Either 'text' or 'binary'; if *mode* is an abbreviation, this parameter is ignored (i.e., `open("r", "binary")` still opens file as text).

```
class pylablib.core.fileio.location.FolderFileSystemDataLocation (folder_path,
                                                                    de-
                                                                    fault_name='content',
                                                                    de-
                                                                    fault_ext="",
                                                                    encod-
                                                                    ing=None)
```

Bases: `pylablib.core.fileio.location.IFileSystemDataLocation`

A location describing a single folder.

#### Parameters

- **folder\_path** (*str*) – A path to the folder. Can also have one or two '|' symbols in the end (e.g., 'folder|file|dat'), which separate default name and extension (overrides *default\_name* and *default\_ext* parameters)
- **default\_name** (*str*) – The default file name.
- **default\_ext** (*str*) – The default file extension.

Multi-level paths translate into nested subfolders.

**get\_filesystem\_path** (*name=None, path\_type='absolute'*)

Get the filesystem path corresponding to a given name.

*path\_type* can be 'absolute' (return absolute path), 'relative' (return relative path; level depends on the location) or 'name' (only return path inside the location).

**close** (*name*)

Close a location file

**generate\_new\_name** (*prefix\_name*, *idx=0*)

Generate a new name inside the location using the given prefix and starting index.

If *idx* is None, check just the *prefix\_name* first before starting to append indices.

**is\_free** (*name=None*)

Check if the name is unoccupied

**list\_opened\_files** ()

Get a dictionary {*string\_name*: *location\_file*} of all files opened in this location

**open** (*name=None*, *mode='read'*, *data\_type='text'*)

Open a location file.

#### Parameters

- **name** – File name inside the location (None means 'default' location),
- **mode** (*str*) – Opening mode. Can be 'read', 'write' or 'append', as well as standard abbreviation (e.g., "r" or "wb").
- **data\_type** (*str*) – Either 'text' or 'binary'; if *mode* is an abbreviation, this parameter is ignored (i.e., `open("r", "binary")` still opens file as text).

`pylablib.core.fileio.location.get_location` (*path*, *loc*, \**args*, \*\**kwargs*)

Build a location.

If *path* or *loc* are instances of *IDataLocation*, return them unchanged. If *loc* is a string, it describes location kind:

- 'single\_file': *SingleFileSystemDataLocation* with the given *path*.
- 'file' or 'prefixed\_file': *PrefixedFileSystemDataLocation* with the given *path* as a master path.
- 'folder': *FolderFileSystemDataLocation* with the given folder *path*.

Any additional arguments are relayed to the constructors.

## pylablib.core.fileio.parse\_csv module

Utilities for parsing CSV files.

**class** `pylablib.core.fileio.parse_csv.ChunksAccumulator` (*dtype='numeric'*, *ignore\_corrupted\_lines=True*, *trim\_rows=False*)

Bases: `object`

Class for accumulating data chunks into a single array.

#### Parameters

- **dtype** – dtype of entries; can be either a single type, or a list of types (one per column). Possible dtypes are: 'int', 'float', 'complex', 'numeric' (tries to coerce to minimal possible numeric type, raises error if data can't be converted to complex), 'generic' (accept arbitrary types, including lists, dictionaries, escaped strings, etc.), 'raw' (keep raw string).



- **ignore\_corrupted\_lines** – if `True`, skip corrupted (e.g., non-numeric for numeric dtype, or with too few entries) lines; otherwise, raise `ValueError`.
- **trim\_rows** – if `True` and the row length is larger than expected, drop extra entries; otherwise, treat the row as corrupted

**corrupted\_number** ()

**convert\_columns** (*raw\_columns*)

Convert raw columns into appropriate data structure (numpy array for numeric dtypes, list for generic and raw).

**add\_columns** (*columns*)

Append columns (lists or numpy arrays) to the existing data.

**add\_chunk** (*chunk*)

Add a chunk (2D list) to the pre-existing data.

```
pylablib.core.fileio.parse_csv.read_columns(f, dtype, delimiters='\s*', \s*\|s+',
                                             empty_entry_substitute=None,
                                             ignore_corrupted_lines=True,
                                             trim_rows=False, stop_comment=None)
```

Load columns from the file stream *f*.

#### Parameters

- **dtype** – dtype of entries; can be either a single type, or a list of types (one per column). Possible dtypes are: 'int', 'float', 'complex', 'numeric' (tries to coerce to minimal possible numeric type, raises error if data can't be converted to complex), 'generic' (accept arbitrary types, including lists, dictionaries, escaped strings, etc.), 'raw' (keep raw string).
- **delimiters** (*str*) – Regex string which recognizes delimiters (by default `r"\s*, \s*\|s+"`, i.e., commas and whitespaces).
- **empty\_entry\_substitute** – Substitute for empty table entries. If `None`, all empty table entries are skipped.
- **ignore\_corrupted\_lines** – If `True`, skip corrupted (e.g., non-numeric for numeric dtype, or with too few entries) lines; otherwise, raise `ValueError`.
- **trim\_rows** – if `True` and the row length is larger than expected, drop extra entries; otherwise, treat the row as corrupted
- **stop\_comment** (*str*) – Regex string for the stopping comment. If not `None`, the function will stop if comment satisfying *stop\_comment* regex is encountered.

#### Returns

(*columns*, *comments*, *corrupted\_lines*).

*columns* is a list of columns with data.

*comments* is a list of comment strings.

*corrupted\_lines* is a dict {'size':list, 'type':list} of corrupted lines (already split into entries), based on the corruption type ('size' means too small size, 'type' means it couldn't be converted using provided dtype).

**Return type** `tuple`

```
pylablib.core.fileio.parse_csv.columns_to_table(data, columns=None,
                                                dtype='numeric',
                                                out_type='columns')
```

Convert *data* (columns list) into a table.

#### Parameters

- **columns** – either number if columns, or a list of columns names.
- **out\_type** (*str*) – type of the result: 'array' for numpy array, 'pandas' for pandas DataFrame, 'columns' for tuple (data, columns)

```
pylablib.core.fileio.parse_csv.read_table(f, dtype='numeric', columns=None,
                                          out_type='columns', delimiters='\s*',
                                          '\s*\s+', empty_entry_substitute=None, ignore_corrupted_lines=True, trim_rows=False,
                                          stop_comment=None)
```

Load table from the file stream *f*.

Arguments are the same as in `read_columns()` and `columns_to_table()`.

#### Returns

(table, comments, corrupted\_lines).

*table* is a table of the format *out\_type*.

*corrupted\_lines* is a dict {'size':list, 'type':list} of corrupted lines (already split into entries), based on the corruption type ('size' means too small size, 'type' means it couldn't be converted using provided dtype).

*comments* is a list of comment strings.

**Return type** tuple

## pylablib.core.fileio.savefile module

Utilities for writing data files.

```
class pylablib.core.fileio.savefile.IOutputFileFormat(format_name)
```

Bases: `object`

Generic class for an output file format.

**Parameters** **format\_name** (*str*) – The name of the format (to be defined in subclasses).

```
write_file(location_file, to_save)
```

```
write_data(location_file, data)
```

```
write(location_file, data)
```

```
class pylablib.core.fileio.savefile.ITextOutputFileFormat(format_name,
                                                           save_props=True,
                                                           save_comments=True,
                                                           save_time=True,
                                                           new_time=True)
```

Bases: `pylablib.core.fileio.savefile.IOutputFileFormat`

Generic class for a text output file format.

#### Parameters

- **format\_name** (*str*) – The name of the format (to be defined in subclasses).

- **save\_props** (*bool*) – If True and saving *datafile.DataFile* object, save its props metainfo.
- **save\_comments** (*bool*) – If True and saving *datafile.DataFile* object, save its comments metainfo.
- **save\_time** (*bool*) – If True, append the file creation time in the end.
- **new\_time** (*bool*) – If saving *datafile.DataFile* object, determines if the time should be updated to the current time.

**make\_comment\_line** (*comment*)

**make\_prop\_line** (*name, value*)

**make\_savetime\_line** (*time*)

**static write\_line** (*stream, line*)

**write\_comments** (*stream, comments*)

**write\_props** (*stream, props*)

**write\_savetime** (*stream, time*)

**write\_file** (*location\_file, to\_save*)

**write** (*location\_file, data*)

**write\_data** (*location\_file, data*)

```
class pylablib.core.fileio.savefile.CSVTableOutputFileFormat (delimiters='t',
                                                                value_formats=None,
                                                                use_rep_classes=False,
                                                                save_columns=True,
                                                                save_props=True,
                                                                save_comments=True,
                                                                save_time=True)
```

Bases: *pylablib.core.fileio.savefile.ITextOutputFileFormat*

Class for CSV output file format.

#### Parameters

- **delimiters** (*str*) – Used to separate entries in a row.
- **value\_formats** (*str*) – If not None, defines value formats to be passed to *utils.string.to\_string()* function.
- **use\_rep\_classes** (*bool*) – If True, use representation classes for Dictionary entries (e.g., numpy arrays will be represented as "array ([1, 2, 3])" instead of just "[1, 2, 3]"); This improves storage fidelity, but makes result harder to parse (e.g., by external string parsers).
- **save\_columns** (*bool*) – If True, save column names as a comment line in the beginning of the file.
- **save\_props** (*bool*) – If True and saving *datafile.DataFile* object, save its props metainfo.
- **save\_comments** (*bool*) – If True and saving *datafile.DataFile* object, save its comments metainfo.
- **save\_time** (*bool*) – If True, append the file creation time in the end.

**get\_table\_line** (*line*)

**get\_columns\_line** (*columns*)

**write\_data** (*location\_file*, *data*)

Write data to a CSV file.

#### Parameters

- **location\_file** – Location of the destination.
- **data** – Data to be saved. Can be a pandas DataFrame or an arbitrary 2D array (numpy array, 2D list, etc.); if the data is not DataFrame or numpy 2D array, it gets converted into a DataFrame using the standard constructor (i.e., 2D list is interpreted as a list of rows)

**make\_comment\_line** (*comment*)

**make\_prop\_line** (*name*, *value*)

**make\_savetime\_line** (*time*)

**write** (*location\_file*, *data*)

**write\_comments** (*stream*, *comments*)

**write\_file** (*location\_file*, *to\_save*)

**static write\_line** (*stream*, *line*)

**write\_props** (*stream*, *props*)

**write\_savetime** (*stream*, *time*)

```
class pylablib.core.fileio.savefile.DictionaryOutputFileFormat (param_formats=None,  
                                                                use_rep_classes=False,  
                                                                table_format='inline',  
                                                                inline_delimiters='t',  
                                                                inline_formats=None,  
                                                                save_props=True,  
                                                                save_comments=True,  
                                                                save_time=True)
```

Bases: `pylablib.core.fileio.savefile.ITextOutputFileFormat`

Class for Dictionary output file format.

#### Parameters

- **param\_formats** (*str*) – If not None, defines value formats to be passed to `utils.string.to_string()` function when writing Dictionary entries.
- **use\_rep\_classes** (*bool*) – If True, use representation classes for Dictionary entries (e.g., numpy arrays will be represented as "array([1, 2, 3])" instead of just "[1, 2, 3]"); This improves storage fidelity, but makes result harder to parse (e.g., by external string parsers).
- **table\_format** (*str*) – Default format for table (numpy arrays or pandas DataFrames) entries. Can be 'inline' (table is written inside the file), 'csv' (external CSV file) or 'bin' (external binary file).
- **inline\_delimiters** (*str*) – Used to separate entries in a row for inline tables.
- **inline\_formats** (*str*) – If not None, defines value formats to be passed to `utils.string.to_string()` function when writing inline tables.

- **save\_props** (*bool*) – If True and saving *datafile.DataFile* object, save its props metainfo.
- **save\_comments** (*bool*) – If True and saving *datafile.DataFile* object, save its comments metainfo.
- **save\_time** (*bool*) – If True, append the file creation time in the end.

**get\_dictionary\_line** (*path, value*)

**write\_data** (*location\_file, data*)

Write data to a Dictionary file.

#### Parameters

- **location\_file** – Location of the destination.
- **data** – Data to be saved. Should be object of class *Dictionary*.

**make\_comment\_line** (*comment*)

**make\_prop\_line** (*name, value*)

**make\_savetime\_line** (*time*)

**write** (*location\_file, data*)

**write\_comments** (*stream, comments*)

**write\_file** (*location\_file, to\_save*)

**static write\_line** (*stream, line*)

**write\_props** (*stream, props*)

**write\_savetime** (*stream, time*)

**class** *pylablib.core.fileio.savefile.IBinaryOutputFileFormat* (*format\_name*)

Bases: *pylablib.core.fileio.savefile.IOutputFileFormat*

**get\_preamble** (*location\_file, data*)

**write** (*location\_file, data*)

**write\_data** (*location\_file, data*)

**write\_file** (*location\_file, to\_save*)

**class** *pylablib.core.fileio.savefile.TableBinaryOutputFileFormat* (*dtype=None, transposed=False*)

Bases: *pylablib.core.fileio.savefile.IBinaryOutputFileFormat*

Class for binary output file format.

#### Parameters

- **dtype** – a string with numpy dtype (e.g., "<f8") used to save the data. By default, use little-endian ("<") variant kind of the supplied data array dtype
- **transposed** (*bool*) – If False, write the data row-wise; otherwise, write it column-wise.

**get\_dtype** (*table*)

**get\_preamble** (*location\_file, data*)

Generate a preamble (dictionary describing the file format).

The parameters are 'dtype', 'packing' ('transposed' or 'flatten', depending on the *transposed* attribute), 'ncol' (number of columns) and 'nrows' (number of rows).

**write\_data** (*location\_file*, *data*)

Write data to a binary file.

#### Parameters

- **location\_file** – Location of the destination.
- **data** – Data to be saved. Can be a pandas DataFrame or an arbitrary 2D array (numpy array, 2D list, etc.) Converted to numpy array before saving.

**write\_file** (*location\_file*, *to\_save*)

**write** (*location\_file*, *data*)

`pylablib.core.fileio.savefile.get_output_format` (*data*, *output\_format*, *\*\*kwargs*)

`pylablib.core.fileio.savefile.save_csv` (*data*, *path*, *delimiters*='\t', *value\_formats*=None, *use\_rep\_classes*=False, *save\_columns*=True, *save\_props*=True, *save\_comments*=True, *save\_time*=True, *loc*='file', *encoding*=None)

Save data to a CSV file.

#### Parameters

- **data** – Data to be saved (2D numpy array, pandas DataFrame, or a *datafile.DataFile* object containing this data).
- **path** (*str*) – Path to the file or a file-like object.
- **delimiters** (*str*) – Used to separate entries in a row.
- **value\_formats** (*str*) – If not None, defines value formats to be passed to *utils.string.to\_string()* function.
- **use\_rep\_classes** (*bool*) – If True, use representation classes for Dictionary entries (e.g., numpy arrays will be represented as "array([1, 2, 3])" instead of just "[1, 2, 3]"); This improves storage fidelity, but makes result harder to parse (e.g., by external string parsers).
- **save\_columns** (*bool*) – If True, save column names as a comment line in the beginning of the file.
- **save\_props** (*bool*) – If True and saving *datafile.DataFile* object, save its props metainfo.
- **save\_comments** (*bool*) – If True and saving *datafile.DataFile* object, save its comments metainfo.
- **save\_time** (*bool*) – If True, append the file creation time in the end.
- **loc** (*str*) – Location type.
- **encoding** – if a new file location is opened, this specifies the encoding.

`pylablib.core.fileio.savefile.save_csv_desc` (*data*, *path*, *loc*='file', *encoding*=None)

Save data table to a dictionary file with an inlined table.

Compared to *save\_csv()*, supports more pandas features (index, column multi-index), but can only be directly read by pylablib.

#### Parameters

- **data** – Data to be saved (2D numpy array, pandas DataFrame, or a *datafile.DataFile* object containing this data).
- **path** (*str*) – Path to the file or a file-like object.
- **loc** (*str*) – Location type.
- **encoding** – if a new file location is opened, this specifies the encoding.

```
pylablib.core.fileio.savefile.save_bin(data, path, dtype=None, transposed=False,
                                         loc='file', encoding=None)
```

Save data to a binary file.

#### Parameters

- **data** – Data to be saved (2D numpy array, pandas DataFrame, or a *datafile.DataFile* object containing this data).
- **path** (*str*) – Path to the file or a file-like object.
- **dtype** – *numpy.dtype* describing the data. By default, use little-endian ("*<*") variant kind of the supplied data array dtype.
- **transposed** (*bool*) – If *False*, write the data row-wise; otherwise, write it column-wise.
- **loc** (*str*) – Location type.
- **encoding** – if a new file location is opened, this specifies the encoding.

```
pylablib.core.fileio.savefile.save_bin_desc(data, path, loc='file', encoding=None)
```

Save data to a binary file with an additional description file, which contains all of the data related to loading (shape, dtype, columns, etc.)

#### Parameters

- **data** – Data to be saved (2D numpy array, pandas DataFrame, or a *datafile.DataFile* object containing this data).
- **path** (*str*) – Path to the file or a file-like object.
- **loc** (*str*) – Location type.
- **encoding** – if a new file location is opened, this specifies the encoding.

```
pylablib.core.fileio.savefile.save_dict(data, path, param_formats=None,
                                         use_rep_classes=False, table_format='inline',
                                         inline_delimiters='\\t', inline_formats=None,
                                         save_props=True, save_comments=True,
                                         save_time=True, loc='file', encoding=None)
```

Save dictionary to a text file.

#### Parameters

- **data** – Data to be saved.
- **path** (*str*) – Path to the file or a file-like object.
- **param\_formats** (*str*) – If not *None*, defines value formats to be passed to *utils.string.to\_string()* function when writing Dictionary entries.
- **use\_rep\_classes** (*bool*) – If *True*, use representation classes for Dictionary entries (e.g., numpy arrays will be represented as "array([1, 2, 3])" instead of just "[1, 2, 3]"); This improves storage fidelity, but makes result harder to parse (e.g., by external string parsers).

- **table\_format** (*str*) – Default format for table (numpy arrays or pandas DataFrames) entries. Can be 'inline' (table is written inside the file), 'csv' (external CSV file) or 'bin' (external binary file).
- **inline\_delimiters** (*str*) – Used to separate entries in a row for inline tables.
- **inline\_formats** (*str*) – If not None, defines value formats to be passed to *utils.string.to\_string()* function when writing inline tables.
- **save\_props** (*bool*) – If True and saving *datafile.DataFile* object, save its props metainfo.
- **save\_comments** (*bool*) – If True and saving *datafile.DataFile* object, save its comments metainfo.
- **save\_time** (*bool*) – If True, append the file creation time in the end.
- **loc** (*str*) – Location type.
- **encoding** – if a new file location is opened, this specifies the encoding.

```
pylablib.core.fileio.savefile.save_generic(data, path, output_format=None, loc='file',  
                                             encoding=None, **kwargs)
```

Save data to a file.

#### Parameters

- **data** – Data to be saved.
- **path** (*str*) – Path to the file or a file-like object.
- **output\_format** (*str*) – Output file format. Can be either None (defaults to 'csv' for table data and 'dict' for Dictionary data), a string with one of the default format names, or an already prepared *IOutputFileFormat* object.
- **loc** (*str*) – Location type.
- **encoding** – if a new file location is opened, this specifies the encoding.

*\*\*kwargs* are passed to the file formatter constructor (see *CSVTableOutputFileFormat*, *DictionaryOutputFileFormat* and *TableBinaryOutputFileFormat* for the possible arguments). The default format names are:

- 'csv': CSV file, corresponds to *CSVTableOutputFileFormat* and *save\_csv()*;
- 'csv': CSV file with an additional dictionary containing format description, corresponds to *DictionaryOutputFileFormat* and *save\_csv\_desc()*;
- 'bin': Binary file, corresponds to *TableBinaryOutputFileFormat* and *save\_bin()*;
- 'bin\_desc': Binary file with an additional dictionary containing format description, corresponds to *DictionaryOutputFileFormat* and *save\_bin\_desc()*;
- 'dict': Dictionary file, corresponds to *DictionaryOutputFileFormat* and *save\_dict()*

#### pylablib.core.fileio.table\_stream module

```
class pylablib.core.fileio.table_stream.TableStreamFile(path,      columns=None,  
                                                         delimiter='t', fmt=None,  
                                                         add_timestamp=False,  
                                                         header_prepend='# ')
```

Bases: *object*

Expanding table file.



Can define column names and formats for different columns, and repeatedly write data into the same file. Useful for, e.g., continuous log files.

#### Parameters

- **path** (*str*) – Path to the destination file.
- **columns** (*list*) – If not `None`, it's a list of column names to be added as a header on creation.
- **delimiter** (*str*) – Values delimiter.
- **fmt** (*str*) – If not `None`, it's a list of format strings for the line entries (e.g., "%.3f"); instead of format string one can also be `None`, which means using the standard `to_string()` conversion function
- **add\_timestamp** (*bool*) – If `True`, add the UNIX timestamp in the beginning of each line (columns and format are expanded accordingly)
- **header\_prepend** – the string to prepend to the header line; by default, a comment symbol, which is best compatibly with `loadfile.load_csv()` function

#### **write\_text\_lines** (*lines*)

Write several text lines into the file.

Create the file if it doesn't exist (in which case the header is automatically added).

**Parameters** **lines** (*[str]*) – List of lines to write.

#### **write\_row** (*row*)

Write a single data row into the file.

Create the file if it doesn't exist (in which case the header is automatically added).

**Parameters** **data** (*list or numpy.ndarray*) – Data row to be added.

#### **write\_multiple\_rows** (*rows*)

Write a multiple data lines into the file.

Create the file if it doesn't exist (in which case the header is automatically added).

**Parameters** **rows** (*[list or numpy.ndarray]*) – Data rows to be added.

## Module contents

### pylablib.core.gui package

#### Subpackages

### pylablib.core.gui.widgets package

#### Submodules

### pylablib.core.gui.widgets.button module

**class** `pylablib.core.gui.widgets.button.ToggleButton` (*parent=None*)

Bases: `object`

Expanded toggle button.

Maintains internally stored consistent value (which can be, e.g., accessed from different threads). Allows setting different captions of pressed/unpressed, and uses those to represent values.

**set\_value\_labels** (*labels*)

Set a list of values corresponding to combo box indices.

Can be either a list of values, whose length must be equal to the number of options, or `None` (don't change the button label on toggle).

**value\_changed** = `<Mock name='mock.QtCore.pyqtSignal()' id='140096333774352'>`

Signal emitted when value is changed

**get\_value** ()

Get current value

**set\_value** (*value*, *notify\_value\_change=True*)

Set current value.

If *notify\_value\_change==True*, emit the *value\_changed* signal; otherwise, change value silently.

**repr\_value** (*value*)

Return representation of *value* as a caption text

## pylablib.core.gui.widgets.combo\_box module

**class** `pylablib.core.gui.widgets.combo_box.ComboBox` (*parent*)

Bases: `object`

Expanded combo box.

Maintains internally stored consistent value (which can be, e.g., accessed from different threads). Allows setting values which are reported via *value\_changed* signal instead of simple indices.

**wheelEvent** (*event*)

**set\_out\_of\_range** (*action='error'*)

Set behavior when out-of-range value is applied.

Can be "error" (raise error), "reset" (reset to no-value position), "reset\_start" (reset to the first position) or "ignore" (keep current value).

**set\_direct\_index\_action** (*action='error'*)

Set behavior when index values are specified, but direct indexing is used.

Can be "ignore" (do not allow direct indexing and treat any value as index value), "value\_default" (allow direct indexing, but prioritize index values with the same value), or "index\_default" (allow direct indexing and prioritize it if index value with the same value exists).

**index\_to\_value** (*idx*)

Turn numerical index into value

**value\_to\_index** (*value*)

Turn value into a numerical index

**set\_index\_values** (*index\_values*, *value=None*, *index=None*)

Set a list of values corresponding to combo box indices.

Can be either a list of values, whose length must be equal to the number of options, or `None` (simply use indices). Note: if the number of combo box options changed (e.g., using `addItem` or `insertItem` methods), the index values need to be manually updated; otherwise, the errors might arise if the index is larger than the number of values. If *value* is specified, set as the new values. If *index* is specified, use it as the index of a new value; if both *value* and *index* are specified, the *value* takes priority.

**get\_index\_values()**  
Return the list of values corresponding to combo box indices

**get\_options()**  
Return the list of labels corresponding to combo box indices

**get\_options\_dict()**  
Return the dictionary {value: label} of the option labels

**set\_options(options, index\_values=None, value=None, index=None)**  
Set new set of options.  
  
If *index\_values* is not None, set these as the new index values; otherwise, index values are reset. If *options* is a dictionary, interpret it as a mapping {option: index\_value}. If *value* is specified, set as the new values. If *index* is specified, use it as the index of a new value; if both *value* and *index* are specified, the *value* takes priority.

**insert\_option(option, index\_value=None, index=None)**  
Insert or append a new option to the list  
  
Insertion (i.e., index is not None) only works for index-valued combo boxes.

**value\_changed = <Mock name='mock.QtCore.pyqtSignal()' id='140096333774352'>**  
Signal emitted when value is changed

**get\_value()**  
Get current numerical value

**set\_value(value, notify\_value\_change=True)**  
Set current value.  
  
If notify\_value\_change==True, emit the *value\_changed* signal; otherwise, change value silently.

**repr\_value(value)**  
Return representation of *value* as a combo box text

## pylablib.core.gui.widgets.container module

**class** pylablib.core.gui.widgets.container.**TTimer**(name, period, timer)  
Bases: tuple

**count()**  
Return number of occurrences of value.

**index()**  
Return first index of value.  
  
Raises ValueError if the value is not present.

**name**

**period**

**timer**

**class** pylablib.core.gui.widgets.container.**TTimerEvent**(start, loop, stop, timer)  
Bases: tuple

**count()**  
Return number of occurrences of value.

**index()**  
Return first index of value.  
Raises ValueError if the value is not present.

**loop**

**start**

**stop**

**timer**

**class** pylablib.core.gui.widgets.container.**TChild**(*name, widget, gui\_values\_path*)  
Bases: `tuple`

**count()**  
Return number of occurrences of value.

**gui\_values\_path**

**index()**  
Return first index of value.  
Raises ValueError if the value is not present.

**name**

**widget**

**class** pylablib.core.gui.widgets.container.**IQContainer**(*\*args, name=None, \*\*kwargs*)

Bases: `object`

Basic controller object which combines and controls several other widget.

Can either corresponds to a widget (e.g., a frame or a group box), or simply be an organizing entity.

**Parameters** **name** – entity name (used by default when adding this object to a values table)

Abstract mix-in class, which needs to be added to a class inheriting from `QObject`. Alternatively, one can directly use `QContainer`, which already inherits from `QObject`.

**TimerUIDGenerator** = `<pylablib.core.utils.general.NamedUIDGenerator object>`

**contained\_value\_changed** = `<Mock name='mock.QtCore.pyqtSignal()' id='140096333774352'>`

**setup\_name**(*name*)  
Set the object's name

**setup**(*name=None*)  
Setup the container by initializing its GUI values and setting the `ctl` attribute

**add\_timer**(*name, period, autostart=True*)  
Add a periodic timer with the given *name* and *period*.

Rarely needs to be called explicitly (one is created automatically if timer event is created). If `autostart==True` and the container has been started (by calling `start()` method), start the timer as well.

**start\_timer**(*name*)  
Start the timer with the given name (also called automatically on `start()` method)

**stop\_timer**(*name*)  
Stop the timer with the given name (also called automatically on `stop()` method)

**is\_timer\_running** (*name*)

Check if the timer with the given name is running

**add\_timer\_event** (*name*, *loop*=None, *start*=None, *stop*=None, *period*=None, *timer*=None, *autostart*=True)

Add timer event with the given *name*.

Add an event which should be called periodically (e.g., a GUI update). Internally implemented through Qt timers. *loop*, *start* and *stop* are the functions called, correspondingly, on timer (periodically), when timer is start, and when it's finished. One can either specify the timer by name (*timer* parameter), or create a new one with the given *period*. If *autostart*==True and the container has been started (by calling *start*() method), start the timer as well.

**add\_child\_values** (*name*, *widget*, *path*, *add\_change\_event*=True)

Add child's values to the container's table.

If *widget* is a container and *path*==" " or ends in "/\*" (e.g., "subpath/\*"), use its *setup\_gui\_values* to make it share the same GUI values; otherwise, simply add it to the GUI values under the given path. if *add\_change\_event*==True, changing of the widget's value emits the container's *contained\_value\_changed* event

**add\_child** (*name*, *widget*, *gui\_values\_path*=True, *add\_change\_event*=True)

Add a contained child widget.

If *gui\_values\_path* is False or None, do not add it to the GUI values table; if it is True, add it under the same root (*path*==" ") if it's a container, and under *name* if it's not; otherwise, *gui\_values\_path* specifies the path under which the widget values are stored. if *add\_change\_event*==True, changing of the widget's value emits the container's *contained\_value\_changed* event

**get\_child** (*name*)

Get the child widget with the given name

**remove\_child** (*name*)

Remove child from the container and clear it

**add\_virtual\_element** (*name*, *value*=None, *multivalued*=False, *add\_indicator*=True)

Add a virtual value element.

Doesn't correspond to any actual widget, but behaves very similarly from the application point of view (its value can be set or read, it has on-change events, it can have indicator). The element value is simply stored on set and retrieved on get. If *multivalued*==True, the internal value is assumed to be complex, so it is forced to be a *Dictionary* every time it is set. If *add\_indicator*==True, add default indicator handler as well.

**add\_property\_element** (*name*, *getter*=None, *setter*=None, *add\_indicator*=True)

Add a property value element.

Doesn't correspond to any actual widget, but behaves very similarly from the application point of view; each time the value is set or get, the corresponding setter and getter methods are called. If *add\_indicator*==True, add default (stored value) indicator handler as well.

**start** ()

Start the container.

Starts all the internal timers, and calls *start* method for all the contained widgets.

**stop** ()

Stop the container.

Stops all the internal timers, and calls *stop* method for all the contained widgets.

**is\_running** ()

Check if the container is running (started and not yet stopped)

**is\_stopping** ()

Check if the container is stopping (stopping initialized and not yet done)

**clear** ()

Clear the container.

Stop all timers and widgets, and call `clear` methods of all contained widgets, remove all widgets from the values table, remove all widgets from the table.

**get\_handler** (*name*)

Get value handler of a widget with the given name

**get\_widget** (*name*)

Get a widget corresponding to a value with the given name

**get\_value** (*name=None*)

Get value of a widget with the given name (None means all values)

**get\_all\_values** ()

Get values of all widget in the container

**set\_value** (*name, value*)

Set value of a widget with the given name (None means all values)

**set\_all\_values** (*value*)

Set values of all widgets in the container

**get\_value\_changed\_signal** (*name*)

Get a value-changed signal for a widget with the given name

**update\_value** (*name=None*)

Send update signal for a handler with a given name or list of names.

Emit a value changed signal with the current value to notify the subscribed slots. If *name* is None, emit for all values in the table.

**get\_indicator** (*name=None*)

Get indicator value for a widget with the given name (None means all indicators)

**get\_all\_indicators** ()

Get indicator values of all widget in the container

**set\_indicator** (*name, value, ignore\_missing=False*)

Set indicator value for a widget or a branch with the given name

**set\_all\_indicators** (*value, ignore\_missing=True*)

**update\_indicators** ()

Update all indicators to represent current values

**class** `pylablib.core.gui.widgets.container.QContainer` (*\*args, name=None, \*\*kwargs*)

Bases: `pylablib.core.gui.widgets.container.IQContainer`, `object`

Basic controller object which combines and controls several other widget.

Can either corresponds to a widget (e.g., a frame or a group box), or simply be an organizing entity.

**Parameters** *name* – entity name (used by default when adding this object to a values table)

Simply a combination of `IQContainer` and `QObject`.

**TimerUIDGenerator** = `<pylablib.core.utils.general.NamedUIDGenerator object>`

**add\_child** (*name, widget, gui\_values\_path=True, add\_change\_event=True*)

Add a contained child widget.

If `gui_values_path` is `False` or `None`, do not add it to the GUI values table; if it is `True`, add it under the same root (`path==""`) if it's a container, and under `name` if it's not; otherwise, `gui_values_path` specifies the path under which the widget values are stored. if `add_change_event==True`, changing of the widget's value emits the container's `contained_value_changed` event

**add\_child\_values** (*name, widget, path, add\_change\_event=True*)

Add child's values to the container's table.

If `widget` is a container and `path==""` or ends in `"/"` (e.g., `"subpath/"`), use its `setup_gui_values` to make it share the same GUI values; otherwise, simply add it to the GUI values under the given path. if `add_change_event==True`, changing of the widget's value emits the container's `contained_value_changed` event

**add\_property\_element** (*name, getter=None, setter=None, add\_indicator=True*)

Add a property value element.

Doesn't correspond to any actual widget, but behaves very similarly from the application point of view; each time the value is set or get, the corresponding setter and getter methods are called. If `add_indicator==True`, add default (stored value) indicator handler as well.

**add\_timer** (*name, period, autostart=True*)

Add a periodic timer with the given *name* and *period*.

Rarely needs to be called explicitly (one is created automatically if timer event is created). If `autostart==True` and the container has been started (by calling `start()` method), start the timer as well.

**add\_timer\_event** (*name, loop=None, start=None, stop=None, period=None, timer=None, autostart=True*)

Add timer event with the given *name*.

Add an event which should be called periodically (e.g., a GUI update). Internally implemented through Qt timers. *loop*, *start* and *stop* are the functions called, correspondingly, on timer (periodically), when timer is start, and when it's finished. One can either specify the timer by name (*timer* parameter), or create a new one with the given *period*. If `autostart==True` and the container has been started (by calling `start()` method), start the timer as well.

**add\_virtual\_element** (*name, value=None, multivalued=False, add\_indicator=True*)

Add a virtual value element.

Doesn't correspond to any actual widget, but behaves very similarly from the application point of view (its value can be set or read, it has on-change events, it can have indicator). The element value is simply stored on set and retrieved on get. If `multivalued==True`, the internal value is assumed to be complex, so it is forced to be a *Dictionary* every time it is set. If `add_indicator==True`, add default indicator handler as well.

**clear** ()

Clear the container.

Stop all timers and widgets, and call `clear` methods of all contained widgets, remove all widgets from the values table, remove all widgets from the table.

**contained\_value\_changed** = `<Mock name='mock.QtCore.pyqtSignal()' id='140096333774352'>`

**get\_all\_indicators** ()

Get indicator values of all widget in the container

**get\_all\_values** ()

Get values of all widget in the container

**get\_child** (*name*)

Get the child widget with the given name

**get\_handler** (*name*)  
Get value handler of a widget with the given name

**get\_indicator** (*name=None*)  
Get indicator value for a widget with the given name (*None* means all indicators)

**get\_value** (*name=None*)  
Get value of a widget with the given name (*None* means all values)

**get\_value\_changed\_signal** (*name*)  
Get a value-changed signal for a widget with the given name

**get\_widget** (*name*)  
Get a widget corresponding to a value with the given name

**is\_running** ()  
Check if the container is running (started and not yet stopped)

**is\_stopping** ()  
Check if the container is stopping (stopping initialized and not yet done)

**is\_timer\_running** (*name*)  
Check if the timer with the given name is running

**remove\_child** (*name*)  
Remove child from the container and clear it

**set\_all\_indicators** (*value, ignore\_missing=True*)

**set\_all\_values** (*value*)  
Set values of all widgets in the container

**set\_indicator** (*name, value, ignore\_missing=False*)  
Set indicator value for a widget or a branch with the given name

**set\_value** (*name, value*)  
Set value of a widget with the given name (*None* means all values)

**setup** (*name=None*)  
Setup the container by initializing its GUI values and setting the `ctl` attribute

**setup\_name** (*name*)  
Set the object's name

**start** ()  
Start the container.  
  
Starts all the internal timers, and calls `start` method for all the contained widgets.

**start\_timer** (*name*)  
Start the timer with the given name (also called automatically on `start()` method)

**stop** ()  
Stop the container.  
  
Stops all the internal timers, and calls `stop` method for all the contained widgets.

**stop\_timer** (*name*)  
Stop the timer with the given name (also called automatically on `stop()` method)

**update\_indicators** ()  
Update all indicators to represent current values



**update\_value** (*name=None*)

Send update signal for a handler with a given name or list of names.

Emit a value changed signal with the current value to notify the subscribed slots. If *name* is *None*, emit for all values in the table.

**class** `pylablib.core.gui.widgets.container.IQWidgetContainer(*args, **kwargs)`

Bases: `pylablib.core.gui.widgets.layout_manager.IQLayoutManagedWidget`,  
`pylablib.core.gui.widgets.container.IQContainer`

Generic widget container.

Combines `IQContainer` management of GUI values and timers with `IQLayoutManagedWidget` management of the contained widget's layout.

Typically, adding widget adds them both to the container values and to the layout; however, this can be skipped by either using `IQLayoutManagedWidget.add_to_layout()` (only add to the layout), or specifying `location="skip"` in `add_child()` (only add to the container).

Abstract mix-in class, which needs to be added to a class inheriting from `QWidget`. Alternatively, one can directly use `QWidgetContainer`, which already inherits from `QWidget`.

**setup** (*layout='vbox', no\_margins=False, name=None*)

Setup the layout.

#### Parameters

- **layout** – layout kind; can be "grid", "vbox" (vertical single-column box), or "hbox" (horizontal single-row box).
- **no\_margins** – if *True*, set all layout margins to zero (useful when the widget is in the middle of layout hierarchy)

**add\_child** (*name, widget, location=None, gui\_values\_path=True*)

Add a contained child widget.

*name* specifies the child storage name; if *name==False*, only add the widget to they layout, but not to the container. *location* specifies the layout location to which the widget is added; if *location=="skip"*, skip adding it to the layout (can be manually added later). Note that if the widget is added to the layout, it will be completely deleted when `clear` or `remove_child` methods are called; otherwise, simply its `clear` method will be called, and its GUI values will be deleted.

If *gui\_values\_path* is *False* or *None*, do not add it to the GUI values table; if it is *True*, add it under the same root (*path=="*) if it's a container, and under *name* if it's not; otherwise, *gui\_values\_path* specifies the path under which the widget values are stored.

**remove\_child** (*name*)

Remove widget from the container and the layout, clear it, and remove it

**add\_frame** (*name, layout='vbox', location=None, gui\_values\_path=True, no\_margins=True*)

Add a new frame container to the layout.

*layout* specifies the layout ("vbox", "hbox", or "grid") of the new frame, and *location* specifies its location within the container layout. If *no\_margins==True*, the frame will have no inner layout margins. The other parameters are the same as in `add_child()` method.

**add\_group\_box** (*name, caption, layout='vbox', location=None, gui\_values\_path=True, no\_margins=True*)

Add a new group box container with the given *caption* to the layout.

*layout* specifies the layout ("vbox", "hbox", or "grid") of the new frame, and *location* specifies its location within the container layout. If *no\_margins==True*, the frame will have no inner layout margins. The other parameters are the same as in `add_child()` method.

**clear()**

Clear the container.

All the timers are stopped, all the contained widgets are cleared and removed.

**TimerUIDGenerator** = <pylablib.core.utils.general.NamedUIDGenerator object>

**add\_child\_values** (*name, widget, path, add\_change\_event=True*)

Add child's values to the container's table.

If *widget* is a container and *path*="" or ends in "/"\* (e.g., "subpath/\*"), use its *setup\_gui\_values* to make it share the same GUI values; otherwise, simply add it to the GUI values under the given path. if *add\_change\_event*==True, changing of the widget's value emits the container's *contained\_value\_changed* event

**add\_decoration\_label** (*text, location='next'*)

Add a decoration text label with the given text

**add\_padding** (*kind='auto', location='next', stretch=0*)

Add a padding (expandable spacer) of the given kind to the given location.

*kind* can be "vertical", "horizontal", "auto" (vertical for *grid* and *vbox* layouts, horizontal for *hbox*), or "both" (stretches in both directions). If *stretch* is not None, it specifies stretch of the spacer the corresponding direction (applied to the upper row and leftmost column for multi-cell spacer); can also be a tuple with two stretches along vertical and horizontal directions.

**add\_property\_element** (*name, getter=None, setter=None, add\_indicator=True*)

Add a property value element.

Doesn't correspond to any actual widget, but behaves very similarly from the application point of view; each time the value is set or get, the corresponding setter and getter methods are called. If *add\_indicator*==True, add default (stored value) indicator handler as well.

**add\_spacer** (*height=0, width=0, stretch\_height=False, stretch\_width=False, stretch=0, location='next'*)

Add a spacer with the given width and height to the given location.

If *stretch\_height*==True or *stretch\_width*==True, the widget will stretch in these directions; otherwise, the widget size is fixed. If *stretch* is not None, it specifies stretch of the spacer the corresponding direction (applied to the upper row and leftmost column for multi-cell spacer); if *kind*=="both", it can also be a tuple with two stretches along vertical and horizontal directions.

**add\_sublayout** (*name, kind='grid', location=None*)

Add a sublayout to the given location.

*name* specifies the sublayout name, which can be used to refer to it in specifying locations later. *kind* can be "grid", "vbox" (vertical single-column box), or "hbox" (horizontal single-row box).

**add\_timer** (*name, period, autostart=True*)

Add a periodic timer with the given *name* and *period*.

Rarely needs to be called explicitly (one is created automatically if timer event is created). If *autostart*==True and the container has been started (by calling *start()* method), start the timer as well.

**add\_timer\_event** (*name, loop=None, start=None, stop=None, period=None, timer=None, autostart=True*)

Add timer event with the given *name*.

Add an event which should be called periodically (e.g., a GUI update). Internally implemented through Qt timers. *loop*, *start* and *stop* are the functions called, correspondingly, on timer (periodically), when timer is start, and when it's finished. One can either specify the timer by name (*timer* parameter), or create

a new one with the given *period*. If `autostart==True` and the container has been started (by calling `start()` method), start the timer as well.

**add\_to\_layout** (*element*, *location=None*, *kind='widget'*)

Add an existing *element* to the layout at the given *location*.

*kind* can be "widget" for widgets, "layout" for other layouts, or "item" for layout items (spacers).

**add\_virtual\_element** (*name*, *value=None*, *multivalued=False*, *add\_indicator=True*)

Add a virtual value element.

Doesn't correspond to any actual widget, but behaves very similarly from the application point of view (its value can be set or read, it has on-change events, it can have indicator). The element value is simply stored on set and retrieved on get. If `multivalued==True`, the internal value is assumed to be complex, so it is forced to be a *Dictionary* every time it is set. If `add_indicator==True`, add default indicator handler as well.

**contained\_value\_changed** = `<Mock name='mock.QtCore.pyqtSignal()' id='140096333774352'>`

**get\_all\_indicators** ()

Get indicator values of all widget in the container

**get\_all\_values** ()

Get values of all widget in the container

**get\_child** (*name*)

Get the child widget with the given name

**get\_element\_position** (*element*)

Get the sublayout and the position of the given widget.

Return tuple (*sublayout*, *location*), where *sublayout* is the sublayout name ("name" for the main layout), and *location* is a tuple (*row*, *column*, *rowspan*, *colspan*). If the given widget is not in this layout, return *None*.

**get\_handler** (*name*)

Get value handler of a widget with the given name

**get\_indicator** (*name=None*)

Get indicator value for a widget with the given name (*None* means all indicators)

**get\_layout\_shape** (*name=None*)

Get shape (*rows*, *cols*) of the current layout

**get\_sublayout** (*name=None*)

Get the previously added sublayout

**get\_sublayout\_kind** (*name=None*)

Get the kind of the previously added sublayout

**get\_value** (*name=None*)

Get value of a widget with the given name (*None* means all values)

**get\_value\_changed\_signal** (*name*)

Get a value-changed signal for a widget with the given name

**get\_widget** (*name*)

Get a widget corresponding to a value with the given name

**insert\_column** (*col*, *sublayout=None*, *stretch=0*)

Insert a new column at the given location in the grid layout

**insert\_row** (*row*, *sublayout=None*, *stretch=0*)

Insert a new row at the given location in the grid layout

**is\_running()**  
Check if the container is running (started and not yet stopped)

**is\_stopping()**  
Check if the container is stopping (stopping initialized and not yet done)

**is\_timer\_running(name)**  
Check if the timer with the given name is running

**iter\_sublayout\_items(name=None, include=('widget', ), nested=False)**  
Iterate over items contained in a given sublayout.  
  
*include* is a tuple which contains items to iterate over; can include "widget" or "layout". If *nested==True*, iterate over items in contained layouts as well.

**remove\_layout\_element(element)**  
Remove a previously added layout element

**set\_all\_indicators(value, ignore\_missing=True)**

**set\_all\_values(value)**  
Set values of all widgets in the container

**set\_column\_stretch(\*args, layout=None)**  
Set column stretch for a given layout.  
  
Takes either two arguments *index* and *stretch*, or a single list of stretches for all columns.

**set\_indicator(name, value, ignore\_missing=False)**  
Set indicator value for a widget or a branch with the given name

**set\_row\_stretch(\*args, layout=None)**  
Set row stretch for a given layout.  
  
Takes either two arguments *index* and *stretch*, or a single list of stretches for all rows.

**set\_value(name, value)**  
Set value of a widget with the given name (None means all values)

**setup\_name(name)**  
Set the object's name

**start()**  
Start the container.  
  
Starts all the internal timers, and calls `start` method for all the contained widgets.

**start\_timer(name)**  
Start the timer with the given name (also called automatically on `start()` method)

**stop()**  
Stop the container.  
  
Stops all the internal timers, and calls `stop` method for all the contained widgets.

**stop\_timer(name)**  
Stop the timer with the given name (also called automatically on `stop()` method)

**update\_indicators()**  
Update all indicators to represent current values

**update\_value(name=None)**  
Send update signal for a handler with a given name or list of names.

Emit a value changed signal with the current value to notify the subscribed slots. If *name* is *None*, emit for all values in the table.

**using\_layout** (*name*)

Use a different sublayout as default inside the *with* block

**using\_new\_sublayout** (*name*, *kind*='grid', *location*=*None*)

Create a different sublayout and use it as default inside the *with* block.

*kind* can be "grid", "vbox" (vertical single-column box), or "hbox" (horizontal single-row box).

**class** pylablib.core.gui.widgets.container.QWidgetContainer (\*args, \*\*kwargs)

Bases: *pylablib.core.gui.widgets.container.IQWidgetContainer*, *object*

Generic widget container.

Combines *IQContainer* management of GUI values and timers with *IQLayoutManagedWidget* management of the contained widget's layout.

Typically, adding widget adds them both to the container values and to the layout; however, this can be skipped by either using *QLayoutManagedWidget.add\_to\_layout()* (only add to the layout), or specifying *location*="skip" in *add\_child()* (only add to the container).

Simply a combination of *IQWidgetContainer* and *QWidget*.

**TimerUIDGenerator** = <pylablib.core.utils.general.NamedUIDGenerator object>

**add\_child** (*name*, *widget*, *location*=*None*, *gui\_values\_path*=*True*)

Add a contained child widget.

*name* specifies the child storage name; if *name*==*False*, only add the widget to the layout, but not to the container. *location* specifies the layout location to which the widget is added; if *location*=="skip", skip adding it to the layout (can be manually added later). Note that if the widget is added to the layout, it will be completely deleted when *clear* or *remove\_child* methods are called; otherwise, simply its *clear* method will be called, and its GUI values will be deleted.

If *gui\_values\_path* is *False* or *None*, do not add it to the GUI values table; if it is *True*, add it under the same root (*path*=="") if it's a container, and under *name* if it's not; otherwise, *gui\_values\_path* specifies the path under which the widget values are stored.

**add\_child\_values** (*name*, *widget*, *path*, *add\_change\_event*=*True*)

Add child's values to the container's table.

If *widget* is a container and *path*==" or ends in "/"\* (e.g., "subpath/\*"), use its *setup\_gui\_values* to make it share the same GUI values; otherwise, simply add it to the GUI values under the given path. if *add\_change\_event*==*True*, changing of the widget's value emits the container's *contained\_value\_changed* event

**add\_decoration\_label** (*text*, *location*='next')

Add a decoration text label with the given text

**add\_frame** (*name*, *layout*='vbox', *location*=*None*, *gui\_values\_path*=*True*, *no\_margins*=*True*)

Add a new frame container to the layout.

*layout* specifies the layout ("vbox", "hbox", or "grid") of the new frame, and *location* specifies its location within the container layout. If *no\_margins*==*True*, the frame will have no inner layout margins. The other parameters are the same as in *add\_child()* method.

**add\_group\_box** (*name*, *caption*, *layout*='vbox', *location*=*None*, *gui\_values\_path*=*True*, *no\_margins*=*True*)

Add a new group box container with the given *caption* to the layout.

*layout* specifies the layout ("vbox", "hbox", or "grid") of the new frame, and *location* specifies its location within the container layout. If *no\_margins==True*, the frame will have no inner layout margins. The other parameters are the same as in *add\_child()* method.

**add\_padding** (*kind='auto', location='next', stretch=0*)

Add a padding (expandable spacer) of the given kind to the given location.

*kind* can be "vertical", "horizontal", "auto" (vertical for grid and vbox layouts, horizontal for hbox), or "both" (stretches in both directions). If *stretch* is not None, it specifies stretch of the spacer the corresponding direction (applied to the upper row and leftmost column for multi-cell spacer); can also be a tuple with two stretches along vertical and horizontal directions.

**add\_property\_element** (*name, getter=None, setter=None, add\_indicator=True*)

Add a property value element.

Doesn't correspond to any actual widget, but behaves very similarly from the application point of view; each time the value is set or get, the corresponding setter and getter methods are called. If *add\_indicator==True*, add default (stored value) indicator handler as well.

**add\_spacer** (*height=0, width=0, stretch\_height=False, stretch\_width=False, stretch=0, location='next'*)

Add a spacer with the given width and height to the given location.

If *stretch\_height==True* or *stretch\_width==True*, the widget will stretch in these directions; otherwise, the widget size is fixed. If *stretch* is not None, it specifies stretch of the spacer the corresponding direction (applied to the upper row and leftmost column for multi-cell spacer); if *kind=="both"*, it can also be a tuple with two stretches along vertical and horizontal directions.

**add\_sublayout** (*name, kind='grid', location=None*)

Add a sublayout to the given location.

*name* specifies the sublayout name, which can be used to refer to it in specifying locations later. *kind* can be "grid", "vbox" (vertical single-column box), or "hbox" (horizontal single-row box).

**add\_timer** (*name, period, autostart=True*)

Add a periodic timer with the given *name* and *period*.

Rarely needs to be called explicitly (one is created automatically if timer event is created). If *autostart==True* and the container has been started (by calling *start()* method), start the timer as well.

**add\_timer\_event** (*name, loop=None, start=None, stop=None, period=None, timer=None, autostart=True*)

Add timer event with the given *name*.

Add an event which should be called periodically (e.g., a GUI update). Internally implemented through Qt timers. *loop*, *start* and *stop* are the functions called, correspondingly, on timer (periodically), when timer is start, and when it's finished. One can either specify the timer by name (*timer* parameter), or create a new one with the given *period*. If *autostart==True* and the container has been started (by calling *start()* method), start the timer as well.

**add\_to\_layout** (*element, location=None, kind='widget'*)

Add an existing *element* to the layout at the given *location*.

*kind* can be "widget" for widgets, "layout" for other layouts, or "item" for layout items (spacers).

**add\_virtual\_element** (*name, value=None, multivalued=False, add\_indicator=True*)

Add a virtual value element.

Doesn't correspond to any actual widget, but behaves very similarly from the application point of view (its value can be set or read, it has on-change events, it can have indicator). The element value is simply stored on set and retrieved on get. If *multivalued==True*, the internal value is assumed to be complex, so it

is forced to be a *Dictionary* every time it is set. If `add_indicator==True`, add default indicator handler as well.

**clear()**

Clear the container.

All the timers are stopped, all the contained widgets are cleared and removed.

**contained\_value\_changed** = `<Mock name='mock.QtCore.pyqtSignal()' id='140096333774352'>`

**get\_all\_indicators()**

Get indicator values of all widget in the container

**get\_all\_values()**

Get values of all widget in the container

**get\_child(name)**

Get the child widget with the given name

**get\_element\_position(element)**

Get the sublayout and the position of the given widget.

Return tuple (sublayout, location), where sublayout is the sublayout name ("name" for the main layout), and location is a tuple (row, column, rowspan, colspan). If the given widget is not in this layout, return None.

**get\_handler(name)**

Get value handler of a widget with the given name

**get\_indicator(name=None)**

Get indicator value for a widget with the given name (None means all indicators)

**get\_layout\_shape(name=None)**

Get shape (rows, cols) of the current layout

**get\_sublayout(name=None)**

Get the previously added sublayout

**get\_sublayout\_kind(name=None)**

Get the kind of the previously added sublayout

**get\_value(name=None)**

Get value of a widget with the given name (None means all values)

**get\_value\_changed\_signal(name)**

Get a value-changed signal for a widget with the given name

**get\_widget(name)**

Get a widget corresponding to a value with the given name

**insert\_column(col, sublayout=None, stretch=0)**

Insert a new column at the given location in the grid layout

**insert\_row(row, sublayout=None, stretch=0)**

Insert a new row at the given location in the grid layout

**is\_running()**

Check if the container is running (started and not yet stopped)

**is\_stopping()**

Check if the container is stopping (stopping initialized and not yet done)

**is\_timer\_running(name)**

Check if the timer with the given name is running



**iter\_sublayout\_items** (*name=None, include=('widget', ), nested=False*)

Iterate over items contained in a given sublayout.

*include* is a tuple which contains items to iterate over; can include "widget" or "layout". If *nested==True*, iterate over items in contained layouts as well.

**remove\_child** (*name*)

Remove widget from the container and the layout, clear it, and remove it

**remove\_layout\_element** (*element*)

Remove a previously added layout element

**set\_all\_indicators** (*value, ignore\_missing=True*)

**set\_all\_values** (*value*)

Set values of all widgets in the container

**set\_column\_stretch** (*\*args, layout=None*)

Set column stretch for a given layout.

Takes either two arguments *index* and *stretch*, or a single list of stretches for all columns.

**set\_indicator** (*name, value, ignore\_missing=False*)

Set indicator value for a widget or a branch with the given name

**set\_row\_stretch** (*\*args, layout=None*)

Set row stretch for a given layout.

Takes either two arguments *index* and *stretch*, or a single list of stretches for all rows.

**set\_value** (*name, value*)

Set value of a widget with the given name (None means all values)

**setup** (*layout='vbox', no\_margins=False, name=None*)

Setup the layout.

#### Parameters

- **layout** – layout kind; can be "grid", "vbox" (vertical single-column box), or "hbox" (horizontal single-row box).
- **no\_margins** – if *True*, set all layout margins to zero (useful when the widget is in the middle of layout hierarchy)

**setup\_name** (*name*)

Set the object's name

**start** ()

Start the container.

Starts all the internal timers, and calls *start* method for all the contained widgets.

**start\_timer** (*name*)

Start the timer with the given name (also called automatically on *start()* method)

**stop** ()

Stop the container.

Stops all the internal timers, and calls *stop* method for all the contained widgets.

**stop\_timer** (*name*)

Stop the timer with the given name (also called automatically on *stop()* method)

**update\_indicators** ()

Update all indicators to represent current values



**update\_value** (*name=None*)

Send update signal for a handler with a given name or list of names.

Emit a value changed signal with the current value to notify the subscribed slots. If *name* is *None*, emit for all values in the table.

**using\_layout** (*name*)

Use a different sublayout as default inside the *with* block

**using\_new\_sublayout** (*name*, *kind='grid'*, *location=None*)

Create a different sublayout and use it as default inside the *with* block.

*kind* can be "grid", "vbox" (vertical single-column box), or "hbox" (horizontal single-row box).

**class** pylablib.core.gui.widgets.container.QFrameContainer (\*args, \*\*kwargs)

Bases: [pylablib.core.gui.widgets.container.IQWidgetContainer](#), [object](#)

An extension of [IQWidgetContainer](#) for a QFrame Qt base class

**TimerUIDGenerator** = <pylablib.core.utils.general.NamedUIDGenerator object>

**add\_child** (*name*, *widget*, *location=None*, *gui\_values\_path=True*)

Add a contained child widget.

*name* specifies the child storage name; if *name==False*, only add the widget to the layout, but not to the container. *location* specifies the layout location to which the widget is added; if *location=="skip"*, skip adding it to the layout (can be manually added later). Note that if the widget is added to the layout, it will be completely deleted when *clear* or *remove\_child* methods are called; otherwise, simply its *clear* method will be called, and its GUI values will be deleted.

If *gui\_values\_path* is *False* or *None*, do not add it to the GUI values table; if it is *True*, add it under the same root (*path==""*) if it's a container, and under *name* if it's not; otherwise, *gui\_values\_path* specifies the path under which the widget values are stored.

**add\_child\_values** (*name*, *widget*, *path*, *add\_change\_event=True*)

Add child's values to the container's table.

If *widget* is a container and *path==""* or ends in *"/"* (e.g., *"subpath/"*), use its *setup\_gui\_values* to make it share the same GUI values; otherwise, simply add it to the GUI values under the given path. if *add\_change\_event==True*, changing of the widget's value emits the container's *contained\_value\_changed* event

**add\_decoration\_label** (*text*, *location='next'*)

Add a decoration text label with the given text

**add\_frame** (*name*, *layout='vbox'*, *location=None*, *gui\_values\_path=True*, *no\_margins=True*)

Add a new frame container to the layout.

*layout* specifies the layout ("vbox", "hbox", or "grid") of the new frame, and *location* specifies its location within the container layout. If *no\_margins==True*, the frame will have no inner layout margins. The other parameters are the same as in [add\\_child\(\)](#) method.

**add\_group\_box** (*name*, *caption*, *layout='vbox'*, *location=None*, *gui\_values\_path=True*, *no\_margins=True*)

Add a new group box container with the given *caption* to the layout.

*layout* specifies the layout ("vbox", "hbox", or "grid") of the new frame, and *location* specifies its location within the container layout. If *no\_margins==True*, the frame will have no inner layout margins. The other parameters are the same as in [add\\_child\(\)](#) method.

**add\_padding** (*kind='auto'*, *location='next'*, *stretch=0*)

Add a padding (expandable spacer) of the given kind to the given location.

*kind* can be "vertical", "horizontal", "auto" (vertical for grid and vbox layouts, horizontal for hbox), or "both" (stretches in both directions). If *stretch* is not None, it specifies stretch of the spacer the corresponding direction (applied to the upper row and leftmost column for multi-cell spacer); can also be a tuple with two stretches along vertical and horizontal directions.

**add\_property\_element** (*name*, *getter=None*, *setter=None*, *add\_indicator=True*)

Add a property value element.

Doesn't correspond to any actual widget, but behaves very similarly from the application point of view; each time the value is set or get, the corresponding setter and getter methods are called. If *add\_indicator==True*, add default (stored value) indicator handler as well.

**add\_spacer** (*height=0*, *width=0*, *stretch\_height=False*, *stretch\_width=False*, *stretch=0*, *location='next'*)

Add a spacer with the given width and height to the given location.

If *stretch\_height==True* or *stretch\_width==True*, the widget will stretch in these directions; otherwise, the widget size is fixed. If *stretch* is not None, it specifies stretch of the spacer the corresponding direction (applied to the upper row and leftmost column for multi-cell spacer); if *kind=="both"*, it can also be a tuple with two stretches along vertical and horizontal directions.

**add\_sublayout** (*name*, *kind='grid'*, *location=None*)

Add a sublayout to the given location.

*name* specifies the sublayout name, which can be used to refer to it in specifying locations later. *kind* can be "grid", "vbox" (vertical single-column box), or "hbox" (horizontal single-row box).

**add\_timer** (*name*, *period*, *autostart=True*)

Add a periodic timer with the given *name* and *period*.

Rarely needs to be called explicitly (one is created automatically if timer event is created). If *autostart==True* and the container has been started (by calling *start()* method), start the timer as well.

**add\_timer\_event** (*name*, *loop=None*, *start=None*, *stop=None*, *period=None*, *timer=None*, *autostart=True*)

Add timer event with the given *name*.

Add an event which should be called periodically (e.g., a GUI update). Internally implemented through Qt timers. *loop*, *start* and *stop* are the functions called, correspondingly, on timer (periodically), when timer is start, and when it's finished. One can either specify the timer by name (*timer* parameter), or create a new one with the given *period*. If *autostart==True* and the container has been started (by calling *start()* method), start the timer as well.

**add\_to\_layout** (*element*, *location=None*, *kind='widget'*)

Add an existing *element* to the layout at the given *location*.

*kind* can be "widget" for widgets, "layout" for other layouts, or "item" for layout items (spacers).

**add\_virtual\_element** (*name*, *value=None*, *multivalued=False*, *add\_indicator=True*)

Add a virtual value element.

Doesn't correspond to any actual widget, but behaves very similarly from the application point of view (its value can be set or read, it has on-change events, it can have indicator). The element value is simply stored on set and retrieved on get. If *multivalued==True*, the internal value is assumed to be complex, so it is forced to be a *Dictionary* every time it is set. If *add\_indicator==True*, add default indicator handler as well.

**clear()**

Clear the container.

All the timers are stopped, all the contained widgets are cleared and removed.

**contained\_value\_changed** = <Mock name='mock.QtCore.pyqtSignal()' id='140096333774352'>

**get\_all\_indicators**()  
Get indicator values of all widget in the container

**get\_all\_values**()  
Get values of all widget in the container

**get\_child**(*name*)  
Get the child widget with the given name

**get\_element\_position**(*element*)  
Get the sublayout and the position of the given widget.  
  
Return tuple (*sublayout*, *location*), where *sublayout* is the sublayout name ("name" for the main layout), and *location* is a tuple (*row*, *column*, *rowspan*, *colspan*). If the given widget is not in this layout, return None.

**get\_handler**(*name*)  
Get value handler of a widget with the given name

**get\_indicator**(*name=None*)  
Get indicator value for a widget with the given name (None means all indicators)

**get\_layout\_shape**(*name=None*)  
Get shape (*rows*, *cols*) of the current layout

**get\_sublayout**(*name=None*)  
Get the previously added sublayout

**get\_sublayout\_kind**(*name=None*)  
Get the kind of the previously added sublayout

**get\_value**(*name=None*)  
Get value of a widget with the given name (None means all values)

**get\_value\_changed\_signal**(*name*)  
Get a value-changed signal for a widget with the given name

**get\_widget**(*name*)  
Get a widget corresponding to a value with the given name

**insert\_column**(*col*, *sublayout=None*, *stretch=0*)  
Insert a new column at the given location in the grid layout

**insert\_row**(*row*, *sublayout=None*, *stretch=0*)  
Insert a new row at the given location in the grid layout

**is\_running**()  
Check if the container is running (started and not yet stopped)

**is\_stopping**()  
Check if the container is stopping (stopping initialized and not yet done)

**is\_timer\_running**(*name*)  
Check if the timer with the given name is running

**iter\_sublayout\_items**(*name=None*, *include=('widget',)*, *nested=False*)  
Iterate over items contained in a given sublayout.  
  
*include* is a tuple which contains items to iterate over; can include "widget" or "layout". If *nested==True*, iterate over items in contained layouts as well.

**remove\_child** (*name*)

Remove widget from the container and the layout, clear it, and remove it

**remove\_layout\_element** (*element*)

Remove a previously added layout element

**set\_all\_indicators** (*value, ignore\_missing=True*)

**set\_all\_values** (*value*)

Set values of all widgets in the container

**set\_column\_stretch** (*\*args, layout=None*)

Set column stretch for a given layout.

Takes either two arguments *index* and *stretch*, or a single list of stretches for all columns.

**set\_indicator** (*name, value, ignore\_missing=False*)

Set indicator value for a widget or a branch with the given name

**set\_row\_stretch** (*\*args, layout=None*)

Set row stretch for a given layout.

Takes either two arguments *index* and *stretch*, or a single list of stretches for all rows.

**set\_value** (*name, value*)

Set value of a widget with the given name (None means all values)

**setup** (*layout='vbox', no\_margins=False, name=None*)

Setup the layout.

#### Parameters

- **layout** – layout kind; can be "grid", "vbox" (vertical single-column box), or "hbox" (horizontal single-row box).
- **no\_margins** – if True, set all layout margins to zero (useful when the widget is in the middle of layout hierarchy)

**setup\_name** (*name*)

Set the object's name

**start** ()

Start the container.

Starts all the internal timers, and calls *start* method for all the contained widgets.

**start\_timer** (*name*)

Start the timer with the given name (also called automatically on *start()* method)

**stop** ()

Stop the container.

Stops all the internal timers, and calls *stop* method for all the contained widgets.

**stop\_timer** (*name*)

Stop the timer with the given name (also called automatically on *stop()* method)

**update\_indicators** ()

Update all indicators to represent current values

**update\_value** (*name=None*)

Send update signal for a handler with a given name or list of names.

Emit a value changed signal with the current value to notify the subscribed slots. If *name* is None, emit for all values in the table.

**using\_layout** (*name*)

Use a different sublayout as default inside the with block

**using\_new\_sublayout** (*name*, *kind*='grid', *location*=None)

Create a different sublayout and use it as default inside the with block.

*kind* can be "grid", "vbox" (vertical single-column box), or "hbox" (horizontal single-row box).

**class** pylablib.core.gui.widgets.container.QDialogContainer (\*args, \*\*kwargs)

Bases: [pylablib.core.gui.widgets.container.IQWidgetContainer](#), object

An extension of [IQWidgetContainer](#) for a QDialog Qt base class

**TimerUIDGenerator** = <pylablib.core.utils.general.NamedUIDGenerator object>

**add\_child** (*name*, *widget*, *location*=None, *gui\_values\_path*=True)

Add a contained child widget.

*name* specifies the child storage name; if *name*==False, only add the widget to they layout, but not to the container. *location* specifies the layout location to which the widget is added; if *location*=="skip", skip adding it to the layout (can be manually added later). Note that if the widget is added to the layout, it will be completely deleted when `clear` or `remove_child` methods are called; otherwise, simply its `clear` method will be called, and its GUI values will be deleted.

If *gui\_values\_path* is False or None, do not add it to the GUI values table; if it is True, add it under the same root (*path*=="") if it's a container, and under *name* if it's not; otherwise, *gui\_values\_path* specifies the path under which the widget values are stored.

**add\_child\_values** (*name*, *widget*, *path*, *add\_change\_event*=True)

Add child's values to the container's table.

If *widget* is a container and *path*==" " or ends in "/"\* (e.g., "subpath/\*"), use its `setup_gui_values` to make it share the same GUI values; otherwise, simply add it to the GUI values under the given path. if *add\_change\_event*==True, changing of the widget's value emits the container's `contained_value_changed` event

**add\_decoration\_label** (*text*, *location*='next')

Add a decoration text label with the given text

**add\_frame** (*name*, *layout*='vbox', *location*=None, *gui\_values\_path*=True, *no\_margins*=True)

Add a new frame container to the layout.

*layout* specifies the layout ("vbox", "hbox", or "grid") of the new frame, and *location* specifies its location within the container layout. If *no\_margins*==True, the frame will have no inner layout margins. The other parameters are the same as in [add\\_child\(\)](#) method.

**add\_group\_box** (*name*, *caption*, *layout*='vbox', *location*=None, *gui\_values\_path*=True, *no\_margins*=True)

Add a new group box container with the given *caption* to the layout.

*layout* specifies the layout ("vbox", "hbox", or "grid") of the new frame, and *location* specifies its location within the container layout. If *no\_margins*==True, the frame will have no inner layout margins. The other parameters are the same as in [add\\_child\(\)](#) method.

**add\_padding** (*kind*='auto', *location*='next', *stretch*=0)

Add a padding (expandable spacer) of the given kind to the given location.

*kind* can be "vertical", "horizontal", "auto" (vertical for grid and vbox layouts, horizontal for hbox), or "both" (stretches in both directions). If *stretch* is not None, it specifies stretch of the spacer the corresponding direction (applied to the upper row and leftmost column for multi-cell spacer); can also be a tuple with two stretches along vertical and horizontal directions.

**add\_property\_element** (*name*, *getter=None*, *setter=None*, *add\_indicator=True*)

Add a property value element.

Doesn't correspond to any actual widget, but behaves very similarly from the application point of view; each time the value is set or get, the corresponding setter and getter methods are called. If *add\_indicator==True*, add default (stored value) indicator handler as well.

**add\_spacer** (*height=0*, *width=0*, *stretch\_height=False*, *stretch\_width=False*, *stretch=0*, *location='next'*)

Add a spacer with the given width and height to the given location.

If *stretch\_height==True* or *stretch\_width==True*, the widget will stretch in these directions; otherwise, the widget size is fixed. If *stretch* is not *None*, it specifies stretch of the spacer the corresponding direction (applied to the upper row and leftmost column for multi-cell spacer); if *kind=="both"*, it can also be a tuple with two stretches along vertical and horizontal directions.

**add\_sublayout** (*name*, *kind='grid'*, *location=None*)

Add a sublayout to the given location.

*name* specifies the sublayout name, which can be used to refer to it in specifying locations later. *kind* can be "grid", "vbox" (vertical single-column box), or "hbox" (horizontal single-row box).

**add\_timer** (*name*, *period*, *autostart=True*)

Add a periodic timer with the given *name* and *period*.

Rarely needs to be called explicitly (one is created automatically if timer event is created). If *autostart==True* and the container has been started (by calling *start()* method), start the timer as well.

**add\_timer\_event** (*name*, *loop=None*, *start=None*, *stop=None*, *period=None*, *timer=None*, *autostart=True*)

Add timer event with the given *name*.

Add an event which should be called periodically (e.g., a GUI update). Internally implemented through Qt timers. *loop*, *start* and *stop* are the functions called, correspondingly, on timer (periodically), when timer is start, and when it's finished. One can either specify the timer by name (*timer* parameter), or create a new one with the given *period*. If *autostart==True* and the container has been started (by calling *start()* method), start the timer as well.

**add\_to\_layout** (*element*, *location=None*, *kind='widget'*)

Add an existing *element* to the layout at the given *location*.

*kind* can be "widget" for widgets, "layout" for other layouts, or "item" for layout items (spacers).

**add\_virtual\_element** (*name*, *value=None*, *multivalued=False*, *add\_indicator=True*)

Add a virtual value element.

Doesn't correspond to any actual widget, but behaves very similarly from the application point of view (its value can be set or read, it has on-change events, it can have indicator). The element value is simply stored on set and retrieved on get. If *multivalued==True*, the internal value is assumed to be complex, so it is forced to be a *Dictionary* every time it is set. If *add\_indicator==True*, add default indicator handler as well.

**clear** ()

Clear the container.

All the timers are stopped, all the contained widgets are cleared and removed.

**contained\_value\_changed** = <Mock name='mock.QtCore.pyqtSignal()' id='140096333774352'>

**get\_all\_indicators** ()

Get indicator values of all widget in the container

**get\_all\_values** ()  
Get values of all widget in the container

**get\_child** (*name*)  
Get the child widget with the given name

**get\_element\_position** (*element*)  
Get the sublayout and the position of the given widget.  
  
Return tuple (*sublayout*, *location*), where *sublayout* is the sublayout name ("name" for the main layout), and *location* is a tuple (*row*, *column*, *rowspan*, *colspan*). If the given widget is not in this layout, return *None*.

**get\_handler** (*name*)  
Get value handler of a widget with the given name

**get\_indicator** (*name=None*)  
Get indicator value for a widget with the given name (*None* means all indicators)

**get\_layout\_shape** (*name=None*)  
Get shape (*rows*, *cols*) of the current layout

**get\_sublayout** (*name=None*)  
Get the previously added sublayout

**get\_sublayout\_kind** (*name=None*)  
Get the kind of the previously added sublayout

**get\_value** (*name=None*)  
Get value of a widget with the given name (*None* means all values)

**get\_value\_changed\_signal** (*name*)  
Get a value-changed signal for a widget with the given name

**get\_widget** (*name*)  
Get a widget corresponding to a value with the given name

**insert\_column** (*col*, *sublayout=None*, *stretch=0*)  
Insert a new column at the given location in the grid layout

**insert\_row** (*row*, *sublayout=None*, *stretch=0*)  
Insert a new row at the given location in the grid layout

**is\_running** ()  
Check if the container is running (started and not yet stopped)

**is\_stopping** ()  
Check if the container is stopping (stopping initialized and not yet done)

**is\_timer\_running** (*name*)  
Check if the timer with the given name is running

**iter\_sublayout\_items** (*name=None*, *include=('widget',)*, *nested=False*)  
Iterate over items contained in a given sublayout.  
  
*include* is a tuple which contains items to iterate over; can include "widget" or "layout". If *nested==True*, iterate over items in contained layouts as well.

**remove\_child** (*name*)  
Remove widget from the container and the layout, clear it, and remove it

**remove\_layout\_element** (*element*)  
Remove a previously added layout element

**set\_all\_indicators** (*value, ignore\_missing=True*)

**set\_all\_values** (*value*)

Set values of all widgets in the container

**set\_column\_stretch** (*\*args, layout=None*)

Set column stretch for a given layout.

Takes either two arguments `index` and `stretch`, or a single list of stretches for all columns.

**set\_indicator** (*name, value, ignore\_missing=False*)

Set indicator value for a widget or a branch with the given name

**set\_row\_stretch** (*\*args, layout=None*)

Set row stretch for a given layout.

Takes either two arguments `index` and `stretch`, or a single list of stretches for all rows.

**set\_value** (*name, value*)

Set value of a widget with the given name (None means all values)

**setup** (*layout='vbox', no\_margins=False, name=None*)

Setup the layout.

#### Parameters

- **layout** – layout kind; can be "grid", "vbox" (vertical single-column box), or "hbox" (horizontal single-row box).
- **no\_margins** – if True, set all layout margins to zero (useful when the widget is in the middle of layout hierarchy)

**setup\_name** (*name*)

Set the object's name

**start** ()

Start the container.

Starts all the internal timers, and calls `start` method for all the contained widgets.

**start\_timer** (*name*)

Start the timer with the given name (also called automatically on `start()` method)

**stop** ()

Stop the container.

Stops all the internal timers, and calls `stop` method for all the contained widgets.

**stop\_timer** (*name*)

Stop the timer with the given name (also called automatically on `stop()` method)

**update\_indicators** ()

Update all indicators to represent current values

**update\_value** (*name=None*)

Send update signal for a handler with a given name or list of names.

Emit a value changed signal with the current value to notify the subscribed slots. If *name* is None, emit for all values in the table.

**using\_layout** (*name*)

Use a different sublayout as default inside the `with` block

**using\_new\_sublayout** (*name, kind='grid', location=None*)

Create a different sublayout and use it as default inside the `with` block.



*kind* can be "grid", "vbox" (vertical single-column box), or "hbox" (horizontal single-row box).

**class** pylablib.core.gui.widgets.container.QGroupBoxContainer (\*args, \*\*kwargs)  
 Bases: [pylablib.core.gui.widgets.container.IQWidgetContainer](#), object

An extension of [IQWidgetContainer](#) for a QGroupBox Qt base class

**setup** (caption=None, layout='vbox', no\_margins=False, name=None)  
 Setup the layout.

#### Parameters

- **layout** – layout kind; can be "grid", "vbox" (vertical single-column box), or "hbox" (horizontal single-row box).
- **no\_margins** – if True, set all layout margins to zero (useful when the widget is in the middle of layout hierarchy)

**TimerUIDGenerator** = <pylablib.core.utils.general.NamedUIDGenerator object>

**add\_child** (name, widget, location=None, gui\_values\_path=True)  
 Add a contained child widget.

*name* specifies the child storage name; if *name*==False, only add the widget to the layout, but not to the container. *location* specifies the layout location to which the widget is added; if *location*=="skip", skip adding it to the layout (can be manually added later). Note that if the widget is added to the layout, it will be completely deleted when `clear` or `remove_child` methods are called; otherwise, simply its `clear` method will be called, and its GUI values will be deleted.

If *gui\_values\_path* is False or None, do not add it to the GUI values table; if it is True, add it under the same root (*path*=="") if it's a container, and under *name* if it's not; otherwise, *gui\_values\_path* specifies the path under which the widget values are stored.

**add\_child\_values** (name, widget, path, add\_change\_event=True)  
 Add child's values to the container's table.

If *widget* is a container and *path*==" " or ends in "/"\* (e.g., "subpath/\*"), use its `setup_gui_values` to make it share the same GUI values; otherwise, simply add it to the GUI values under the given path. If *add\_change\_event*==True, changing of the widget's value emits the container's `contained_value_changed` event

**add\_decoration\_label** (text, location='next')  
 Add a decoration text label with the given text

**add\_frame** (name, layout='vbox', location=None, gui\_values\_path=True, no\_margins=True)  
 Add a new frame container to the layout.

*layout* specifies the layout ("vbox", "hbox", or "grid") of the new frame, and *location* specifies its location within the container layout. If *no\_margins*==True, the frame will have no inner layout margins. The other parameters are the same as in [add\\_child\(\)](#) method.

**add\_group\_box** (name, caption, layout='vbox', location=None, gui\_values\_path=True, no\_margins=True)  
 Add a new group box container with the given *caption* to the layout.

*layout* specifies the layout ("vbox", "hbox", or "grid") of the new frame, and *location* specifies its location within the container layout. If *no\_margins*==True, the frame will have no inner layout margins. The other parameters are the same as in [add\\_child\(\)](#) method.

**add\_padding** (kind='auto', location='next', stretch=0)  
 Add a padding (expandable spacer) of the given kind to the given location.

*kind* can be "vertical", "horizontal", "auto" (vertical for grid and vbox layouts, horizontal for hbox), or "both" (stretches in both directions). If *stretch* is not None, it specifies stretch of the

spacer the corresponding direction (applied to the upper row and leftmost column for multi-cell spacer); can also be a tuple with two stretches along vertical and horizontal directions.

**add\_property\_element** (*name*, *getter=None*, *setter=None*, *add\_indicator=True*)

Add a property value element.

Doesn't correspond to any actual widget, but behaves very similarly from the application point of view; each time the value is set or get, the corresponding setter and getter methods are called. If *add\_indicator==True*, add default (stored value) indicator handler as well.

**add\_spacer** (*height=0*, *width=0*, *stretch\_height=False*, *stretch\_width=False*, *stretch=0*, *location='next'*)

Add a spacer with the given width and height to the given location.

If *stretch\_height==True* or *stretch\_width==True*, the widget will stretch in these directions; otherwise, the widget size is fixed. If *stretch* is not *None*, it specifies stretch of the spacer the corresponding direction (applied to the upper row and leftmost column for multi-cell spacer); if *kind=="both"*, it can also be a tuple with two stretches along vertical and horizontal directions.

**add\_sublayout** (*name*, *kind='grid'*, *location=None*)

Add a sublayout to the given location.

*name* specifies the sublayout name, which can be used to refer to it in specifying locations later. *kind* can be "grid", "vbox" (vertical single-column box), or "hbox" (horizontal single-row box).

**add\_timer** (*name*, *period*, *autostart=True*)

Add a periodic timer with the given *name* and *period*.

Rarely needs to be called explicitly (one is created automatically if timer event is created). If *autostart==True* and the container has been started (by calling *start()* method), start the timer as well.

**add\_timer\_event** (*name*, *loop=None*, *start=None*, *stop=None*, *period=None*, *timer=None*, *autostart=True*)

Add timer event with the given *name*.

Add an event which should be called periodically (e.g., a GUI update). Internally implemented through Qt timers. *loop*, *start* and *stop* are the functions called, correspondingly, on timer (periodically), when timer is start, and when it's finished. One can either specify the timer by name (*timer* parameter), or create a new one with the given *period*. If *autostart==True* and the container has been started (by calling *start()* method), start the timer as well.

**add\_to\_layout** (*element*, *location=None*, *kind='widget'*)

Add an existing *element* to the layout at the given *location*.

*kind* can be "widget" for widgets, "layout" for other layouts, or "item" for layout items (spacers).

**add\_virtual\_element** (*name*, *value=None*, *multivalued=False*, *add\_indicator=True*)

Add a virtual value element.

Doesn't correspond to any actual widget, but behaves very similarly from the application point of view (its value can be set or read, it has on-change events, it can have indicator). The element value is simply stored on set and retrieved on get. If *multivalued==True*, the internal value is assumed to be complex, so it is forced to be a *Dictionary* every time it is set. If *add\_indicator==True*, add default indicator handler as well.

**clear()**

Clear the container.

All the timers are stopped, all the contained widgets are cleared and removed.

**contained\_value\_changed** = <Mock name='mock.QtCore.pyqtSignal()' id='140096333774352'>

**get\_all\_indicators()**  
Get indicator values of all widget in the container

**get\_all\_values()**  
Get values of all widget in the container

**get\_child(*name*)**  
Get the child widget with the given name

**get\_element\_position(*element*)**  
Get the sublayout and the position of the given widget.  
  
Return tuple (*sublayout*, *location*), where *sublayout* is the sublayout name ("name" for the main layout), and *location* is a tuple (*row*, *column*, *rowspan*, *colspan*). If the given widget is not in this layout, return None.

**get\_handler(*name*)**  
Get value handler of a widget with the given name

**get\_indicator(*name=None*)**  
Get indicator value for a widget with the given name (None means all indicators)

**get\_layout\_shape(*name=None*)**  
Get shape (*rows*, *cols*) of the current layout

**get\_sublayout(*name=None*)**  
Get the previously added sublayout

**get\_sublayout\_kind(*name=None*)**  
Get the kind of the previously added sublayout

**get\_value(*name=None*)**  
Get value of a widget with the given name (None means all values)

**get\_value\_changed\_signal(*name*)**  
Get a value-changed signal for a widget with the given name

**get\_widget(*name*)**  
Get a widget corresponding to a value with the given name

**insert\_column(*col*, *sublayout=None*, *stretch=0*)**  
Insert a new column at the given location in the grid layout

**insert\_row(*row*, *sublayout=None*, *stretch=0*)**  
Insert a new row at the given location in the grid layout

**is\_running()**  
Check if the container is running (started and not yet stopped)

**is\_stopping()**  
Check if the container is stopping (stopping initialized and not yet done)

**is\_timer\_running(*name*)**  
Check if the timer with the given name is running

**iter\_sublayout\_items(*name=None*, *include=('widget',)*, *nested=False*)**  
Iterate over items contained in a given sublayout.  
  
*include* is a tuple which contains items to iterate over; can include "widget" or "layout". If *nested==True*, iterate over items in contained layouts as well.

**remove\_child(*name*)**  
Remove widget from the container and the layout, clear it, and remove it

**remove\_layout\_element** (*element*)

Remove a previously added layout element

**set\_all\_indicators** (*value, ignore\_missing=True*)

**set\_all\_values** (*value*)

Set values of all widgets in the container

**set\_column\_stretch** (*\*args, layout=None*)

Set column stretch for a given layout.

Takes either two arguments `index` and `stretch`, or a single list of stretches for all columns.

**set\_indicator** (*name, value, ignore\_missing=False*)

Set indicator value for a widget or a branch with the given name

**set\_row\_stretch** (*\*args, layout=None*)

Set row stretch for a given layout.

Takes either two arguments `index` and `stretch`, or a single list of stretches for all rows.

**set\_value** (*name, value*)

Set value of a widget with the given name (None means all values)

**setup\_name** (*name*)

Set the object's name

**start** ()

Start the container.

Starts all the internal timers, and calls `start` method for all the contained widgets.

**start\_timer** (*name*)

Start the timer with the given name (also called automatically on `start()` method)

**stop** ()

Stop the container.

Stops all the internal timers, and calls `stop` method for all the contained widgets.

**stop\_timer** (*name*)

Stop the timer with the given name (also called automatically on `stop()` method)

**update\_indicators** ()

Update all indicators to represent current values

**update\_value** (*name=None*)

Send update signal for a handler with a given name or list of names.

Emit a value changed signal with the current value to notify the subscribed slots. If *name* is `None`, emit for all values in the table.

**using\_layout** (*name*)

Use a different sublayout as default inside the `with` block

**using\_new\_sublayout** (*name, kind='grid', location=None*)

Create a different sublayout and use it as default inside the `with` block.

*kind* can be "grid", "vbox" (vertical single-column box), or "hbox" (horizontal single-row box).

```
class pylablib.core.gui.widgets.container.QScrollAreaContainer (*args,  
                                                                name=None,  
                                                                **kwargs)
```

Bases: `pylablib.core.gui.widgets.container.IQContainer, object`

An extension of `IQWidgetContainer` for a `QScrollArea` Qt base class.

Due to Qt organization, this container is “intermediate”: it contains only a single *QWidgetContainer* widget (named “widget”), which in turn has all of the standard container traits: layout, multiple widgets, etc. Hence, when dealing with any container methods (adding children, changing layout, etc.), this widget (accessible with `.widget()` method) should be used.

```
class QContainedWidget (*args, **kwargs)
    Bases: pylablib.core.gui.widgets.container.QWidgetContainer

    resizeEvent (event)

    TimerUIDGenerator = <pylablib.core.utils.general.NamedUIDGenerator object>

    add_child (name, widget, location=None, gui_values_path=True)
        Add a contained child widget.

        name specifies the child storage name; if name==False, only add the widget to they layout,
        but not to the container. location specifies the layout location to which the widget is added; if
        location=="skip", skip adding it to the layout (can be manually added later). Note that if
        the widget is added to the layout, it will be completely deleted when clear or remove_child
        methods are called; otherwise, simply its clear method will be called, and its GUI values will be
        deleted.

        If gui_values_path is False or None, do not add it to the GUI values table; if it is True, add
        it under the same root (path=="") if it's a container, and under name if it's not; otherwise,
        gui_values_path specifies the path under which the widget values are stored.

    add_child_values (name, widget, path, add_change_event=True)
        Add child's values to the container's table.

        If widget is a container and path==" or ends in /* (e.g., "subpath/*"), use its
        setup_gui_values to make it share the same GUI values; otherwise, simply add it to the GUI
        values under the given path. if add_change_event==True, changing of the widget's value
        emits the container's contained_value_changed event

    add_decoration_label (text, location='next')
        Add a decoration text label with the given text

    add_frame (name, layout='vbox', location=None, gui_values_path=True, no_margins=True)
        Add a new frame container to the layout.

        layout specifies the layout ("vbox", "hbox", or "grid") of the new frame, and location specifies
        its location within the container layout. If no_margins==True, the frame will have no inner
        layout margins. The other parameters are the same as in add_child() method.

    add_group_box (name, caption, layout='vbox', location=None, gui_values_path=True,
                    no_margins=True)
        Add a new group box container with the given caption to the layout.

        layout specifies the layout ("vbox", "hbox", or "grid") of the new frame, and location specifies
        its location within the container layout. If no_margins==True, the frame will have no inner
        layout margins. The other parameters are the same as in add_child() method.

    add_padding (kind='auto', location='next', stretch=0)
        Add a padding (expandable spacer) of the given kind to the given location.

        kind can be "vertical", "horizontal", "auto" (vertical for grid and vbox layouts, hor-
        izontal for hbox), or "both" (stretches in both directions). If stretch is not None, it specifies
        stretch of the spacer the corresponding direction (applied to the upper row and leftmost column for
        multi-cell spacer); can also be a tuple with two stretches along vertical and horizontal directions.

    add_property_element (name, getter=None, setter=None, add_indicator=True)
        Add a property value element.
```

Doesn't correspond to any actual widget, but behaves very similarly from the application point of view; each time the value is set or get, the corresponding setter and getter methods are called. If `add_indicator==True`, add default (stored value) indicator handler as well.

**add\_spacer** (*height=0, width=0, stretch\_height=False, stretch\_width=False, stretch=0, location='next'*)

Add a spacer with the given width and height to the given location.

If `stretch_height==True` or `stretch_width==True`, the widget will stretch in these directions; otherwise, the widget size is fixed. If *stretch* is not `None`, it specifies stretch of the spacer the corresponding direction (applied to the upper row and leftmost column for multi-cell spacer); if *kind*==`"both"`, it can also be a tuple with two stretches along vertical and horizontal directions.

**add\_sublayout** (*name, kind='grid', location=None*)

Add a sublayout to the given location.

*name* specifies the sublayout name, which can be used to refer to it in specifying locations later. *kind* can be `"grid"`, `"vbox"` (vertical single-column box), or `"hbox"` (horizontal single-row box).

**add\_timer** (*name, period, autostart=True*)

Add a periodic timer with the given *name* and *period*.

Rarely needs to be called explicitly (one is created automatically if timer event is created). If `autostart==True` and the container has been started (by calling `start()` method), start the timer as well.

**add\_timer\_event** (*name, loop=None, start=None, stop=None, period=None, timer=None, autostart=True*)

Add timer event with the given *name*.

Add an event which should be called periodically (e.g., a GUI update). Internally implemented through Qt timers. *loop*, *start* and *stop* are the functions called, correspondingly, on timer (periodically), when timer is start, and when it's finished. One can either specify the timer by name (*timer* parameter), or create a new one with the given *period*. If `autostart==True` and the container has been started (by calling `start()` method), start the timer as well.

**add\_to\_layout** (*element, location=None, kind='widget'*)

Add an existing *element* to the layout at the given *location*.

*kind* can be `"widget"` for widgets, `"layout"` for other layouts, or `"item"` for layout items (spacers).

**add\_virtual\_element** (*name, value=None, multivalued=False, add\_indicator=True*)

Add a virtual value element.

Doesn't correspond to any actual widget, but behaves very similarly from the application point of view (its value can be set or read, it has on-change events, it can have indicator). The element value is simply stored on set and retrieved on get. If `multivalued==True`, the internal value is assumed to be complex, so it is forced to be a *Dictionary* every time it is set. If `add_indicator==True`, add default indicator handler as well.

**clear** ()

Clear the container.

All the timers are stopped, all the contained widgets are cleared and removed.

**contained\_value\_changed** = `<Mock name='mock.QtCore.pyqtSignal()' id='14009633377435'`

**get\_all\_indicators** ()

Get indicator values of all widget in the container

**get\_all\_values()**  
Get values of all widget in the container

**get\_child(name)**  
Get the child widget with the given name

**get\_element\_position(element)**  
Get the sublayout and the position of the given widget.  
  
Return tuple (sublayout, location), where sublayout is the sublayout name ("name" for the main layout), and location is a tuple (row, column, rowspan, colspan). If the given widget is not in this layout, return None.

**get\_handler(name)**  
Get value handler of a widget with the given name

**get\_indicator(name=None)**  
Get indicator value for a widget with the given name (None means all indicators)

**get\_layout\_shape(name=None)**  
Get shape (rows, cols) of the current layout

**get\_sublayout(name=None)**  
Get the previously added sublayout

**get\_sublayout\_kind(name=None)**  
Get the kind of the previously added sublayout

**get\_value(name=None)**  
Get value of a widget with the given name (None means all values)

**get\_value\_changed\_signal(name)**  
Get a value-changed signal for a widget with the given name

**get\_widget(name)**  
Get a widget corresponding to a value with the given name

**insert\_column(col, sublayout=None, stretch=0)**  
Insert a new column at the given location in the grid layout

**insert\_row(row, sublayout=None, stretch=0)**  
Insert a new row at the given location in the grid layout

**is\_running()**  
Check if the container is running (started and not yet stopped)

**is\_stopping()**  
Check if the container is stopping (stopping initialized and not yet done)

**is\_timer\_running(name)**  
Check if the timer with the given name is running

**iter\_sublayout\_items(name=None, include=('widget', ), nested=False)**  
Iterate over items contained in a given sublayout.  
  
*include* is a tuple which contains items to iterate over; can include "widget" or "layout". If *nested==True*, iterate over items in contained layouts as well.

**remove\_child(name)**  
Remove widget from the container and the layout, clear it, and remove it

**remove\_layout\_element(element)**  
Remove a previously added layout element

**set\_all\_indicators** (*value, ignore\_missing=True*)

**set\_all\_values** (*value*)

Set values of all widgets in the container

**set\_column\_stretch** (*\*args, layout=None*)

Set column stretch for a given layout.

Takes either two arguments *index* and *stretch*, or a single list of stretches for all columns.

**set\_indicator** (*name, value, ignore\_missing=False*)

Set indicator value for a widget or a branch with the given name

**set\_row\_stretch** (*\*args, layout=None*)

Set row stretch for a given layout.

Takes either two arguments *index* and *stretch*, or a single list of stretches for all rows.

**set\_value** (*name, value*)

Set value of a widget with the given name (None means all values)

**setup** (*layout='vbox', no\_margins=False, name=None*)

Setup the layout.

#### Parameters

- **layout** – layout kind; can be "grid", "vbox" (vertical single-column box), or "hbox" (horizontal single-row box).
- **no\_margins** – if True, set all layout margins to zero (useful when the widget is in the middle of layout hierarchy)

**setup\_name** (*name*)

Set the object's name

**start** ()

Start the container.

Starts all the internal timers, and calls `start` method for all the contained widgets.

**start\_timer** (*name*)

Start the timer with the given name (also called automatically on `start()` method)

**stop** ()

Stop the container.

Stops all the internal timers, and calls `stop` method for all the contained widgets.

**stop\_timer** (*name*)

Stop the timer with the given name (also called automatically on `stop()` method)

**update\_indicators** ()

Update all indicators to represent current values

**update\_value** (*name=None*)

Send update signal for a handler with a given name or list of names.

Emit a value changed signal with the current value to notify the subscribed slots. If *name* is None, emit for all values in the table.

**using\_layout** (*name*)

Use a different sublayout as default inside the `with` block

**using\_new\_sublayout** (*name, kind='grid', location=None*)

Create a different sublayout and use it as default inside the `with` block.



*kind* can be "grid", "vbox" (vertical single-column box), or "hbox" (horizontal single-row box).

**setup** (*layout='vbox', no\_margins=False, name=None, fix\_width=True, fix\_height=False*)

Setup the container.

*layout* specifies the container layout, *no\_margins* determines whether margins within the container are removed, *name* specifies the widget name (if not specified yet). *fix\_width* and *fix\_height* determine whether the corresponding direction behaves as a scroll window (i.e., the size is fixed when the content changes), or as a standard widget container (the size is determined by the content).

**clear** ()

Clear the container.

Stop all timers and widgets, and call `clear` methods of all contained widgets, remove all widgets from the values table, remove all widgets from the table.

**TimerUIDGenerator** = <pylablib.core.utils.general.NamedUIDGenerator object>

**add\_child** (*name, widget, gui\_values\_path=True, add\_change\_event=True*)

Add a contained child widget.

If *gui\_values\_path* is `False` or `None`, do not add it to the GUI values table; if it is `True`, add it under the same root (*path*=="") if it's a container, and under *name* if it's not; otherwise, *gui\_values\_path* specifies the path under which the widget values are stored. if *add\_change\_event*==`True`, changing of the widget's value emits the container's `contained_value_changed` event

**add\_child\_values** (*name, widget, path, add\_change\_event=True*)

Add child's values to the container's table.

If *widget* is a container and *path*==" " or ends in `"/*` (e.g., `"subpath/*"`), use its `setup_gui_values` to make it share the same GUI values; otherwise, simply add it to the GUI values under the given path. if *add\_change\_event*==`True`, changing of the widget's value emits the container's `contained_value_changed` event

**add\_property\_element** (*name, getter=None, setter=None, add\_indicator=True*)

Add a property value element.

Doesn't correspond to any actual widget, but behaves very similarly from the application point of view; each time the value is set or get, the corresponding setter and getter methods are called. If *add\_indicator*==`True`, add default (stored value) indicator handler as well.

**add\_timer** (*name, period, autostart=True*)

Add a periodic timer with the given *name* and *period*.

Rarely needs to be called explicitly (one is created automatically if timer event is created). If *autostart*==`True` and the container has been started (by calling `start()` method), start the timer as well.

**add\_timer\_event** (*name, loop=None, start=None, stop=None, period=None, timer=None, autostart=True*)

Add timer event with the given *name*.

Add an event which should be called periodically (e.g., a GUI update). Internally implemented through Qt timers. *loop*, *start* and *stop* are the functions called, correspondingly, on timer (periodically), when timer is start, and when it's finished. One can either specify the timer by name (*timer* parameter), or create a new one with the given *period*. If *autostart*==`True` and the container has been started (by calling `start()` method), start the timer as well.

**add\_virtual\_element** (*name, value=None, multivalued=False, add\_indicator=True*)

Add a virtual value element.

Doesn't correspond to any actual widget, but behaves very similarly from the application point of view (its value can be set or read, it has on-change events, it can have indicator). The element value is simply stored on set and retrieved on get. If `multivalued==True`, the internal value is assumed to be complex, so it is forced to be a *Dictionary* every time it is set. If `add_indicator==True`, add default indicator handler as well.

```
contained_value_changed = <Mock name='mock.QtCore.pyqtSignal()' id='140096333774352'>
get_all_indicators()
    Get indicator values of all widget in the container
get_all_values()
    Get values of all widget in the container
get_child(name)
    Get the child widget with the given name
get_handler(name)
    Get value handler of a widget with the given name
get_indicator(name=None)
    Get indicator value for a widget with the given name (None means all indicators)
get_value(name=None)
    Get value of a widget with the given name (None means all values)
get_value_changed_signal(name)
    Get a value-changed signal for a widget with the given name
get_widget(name)
    Get a widget corresponding to a value with the given name
is_running()
    Check if the container is running (started and not yet stopped)
is_stopping()
    Check if the container is stopping (stopping initialized and not yet done)
is_timer_running(name)
    Check if the timer with the given name is running
remove_child(name)
    Remove child from the container and clear it
set_all_indicators(value, ignore_missing=True)
set_all_values(value)
    Set values of all widgets in the container
set_indicator(name, value, ignore_missing=False)
    Set indicator value for a widget or a branch with the given name
set_value(name, value)
    Set value of a widget with the given name (None means all values)
setup_name(name)
    Set the object's name
start()
    Start the container.

    Starts all the internal timers, and calls start method for all the contained widgets.
```

**start\_timer** (*name*)

Start the timer with the given name (also called automatically on *start()* method)

**stop** ()

Stop the container.

Stops all the internal timers, and calls *stop* method for all the contained widgets.

**stop\_timer** (*name*)

Stop the timer with the given name (also called automatically on *stop()* method)

**update\_indicators** ()

Update all indicators to represent current values

**update\_value** (*name=None*)

Send update signal for a handler with a given name or list of names.

Emit a value changed signal with the current value to notify the subscribed slots. If *name* is *None*, emit for all values in the table.

**class** `pylablib.core.gui.widgets.container.QTabContainer` (*\*args, \*\*kwargs*)

Bases: `pylablib.core.gui.widgets.container.IQContainer`, `object`

Container which manages tab widget.

Does not have its own layout, but can add or remove tabs, which are represented as `QFrameContainer` widgets.

**add\_tab** (*name*, *caption*, *index=None*, *widget=None*, *layout='vbox'*, *gui\_values\_path=True*, *no\_margins=True*)

Add a new tab container with the given *caption* to the widget.

*index* specifies the new tab's index (*None* means adding to the end, negative values count from the end). If *widget* is *None*, create a new frame widget using the given *layout* ("vbox", "hbox", or "grid") and *no\_margins* (specifies whether the frame has inner margins) arguments; otherwise, use the supplied widget. The other parameters are the same as in *add\_child()* method.

**remove\_tab** (*name*)

Remove a tab with the given name.

Clear it, remove its GUI values, and delete it and all contained widgets.

**clear** ()

Clear the container.

Stop all timers and widgets, and call *clear* methods of all contained widgets, remove all widgets from the values table, remove all widgets from the table.

**get\_current\_name** ()

Get current tab name

**set\_by\_name** (*name*)

Set tab by name

**TimerUIDGenerator** = `<pylablib.core.utils.general.NamedUIDGenerator object>`

**add\_child** (*name*, *widget*, *gui\_values\_path=True*, *add\_change\_event=True*)

Add a contained child widget.

If *gui\_values\_path* is *False* or *None*, do not add it to the GUI values table; if it is *True*, add it under the same root (*path=""*) if it's a container, and under *name* if it's not; otherwise, *gui\_values\_path* specifies the path under which the widget values are stored. if *add\_change\_event==True*, changing of the widget's value emits the container's *contained\_value\_changed* event

**add\_child\_values** (*name, widget, path, add\_change\_event=True*)

Add child's values to the container's table.

If *widget* is a container and *path*=="" or ends in "/"\* (e.g., "subpath/\*"), use its `setup_gui_values` to make it share the same GUI values; otherwise, simply add it to the GUI values under the given path. If *add\_change\_event*==True, changing of the widget's value emits the container's `contained_value_changed` event

**add\_property\_element** (*name, getter=None, setter=None, add\_indicator=True*)

Add a property value element.

Doesn't correspond to any actual widget, but behaves very similarly from the application point of view; each time the value is set or get, the corresponding setter and getter methods are called. If *add\_indicator*==True, add default (stored value) indicator handler as well.

**add\_timer** (*name, period, autostart=True*)

Add a periodic timer with the given *name* and *period*.

Rarely needs to be called explicitly (one is created automatically if timer event is created). If *autostart*==True and the container has been started (by calling `start()` method), start the timer as well.

**add\_timer\_event** (*name, loop=None, start=None, stop=None, period=None, timer=None, autostart=True*)

Add timer event with the given *name*.

Add an event which should be called periodically (e.g., a GUI update). Internally implemented through Qt timers. *loop*, *start* and *stop* are the functions called, correspondingly, on timer (periodically), when timer is start, and when it's finished. One can either specify the timer by name (*timer* parameter), or create a new one with the given *period*. If *autostart*==True and the container has been started (by calling `start()` method), start the timer as well.

**add\_virtual\_element** (*name, value=None, multivalued=False, add\_indicator=True*)

Add a virtual value element.

Doesn't correspond to any actual widget, but behaves very similarly from the application point of view (its value can be set or read, it has on-change events, it can have indicator). The element value is simply stored on set and retrieved on get. If *multivalued*==True, the internal value is assumed to be complex, so it is forced to be a *Dictionary* every time it is set. If *add\_indicator*==True, add default indicator handler as well.

**contained\_value\_changed** = <Mock name='mock.QtCore.pyqtSignal()' id='140096333774352'>

**get\_all\_indicators** ()

Get indicator values of all widget in the container

**get\_all\_values** ()

Get values of all widget in the container

**get\_child** (*name*)

Get the child widget with the given name

**get\_handler** (*name*)

Get value handler of a widget with the given name

**get\_indicator** (*name=None*)

Get indicator value for a widget with the given name (None means all indicators)

**get\_value** (*name=None*)

Get value of a widget with the given name (None means all values)

**get\_value\_changed\_signal** (*name*)

Get a value-changed signal for a widget with the given name

**get\_widget** (*name*)  
Get a widget corresponding to a value with the given name

**is\_running** ()  
Check if the container is running (started and not yet stopped)

**is\_stopping** ()  
Check if the container is stopping (stopping initialized and not yet done)

**is\_timer\_running** (*name*)  
Check if the timer with the given name is running

**remove\_child** (*name*)  
Remove child from the container and clear it

**set\_all\_indicators** (*value*, *ignore\_missing=True*)

**set\_all\_values** (*value*)  
Set values of all widgets in the container

**set\_indicator** (*name*, *value*, *ignore\_missing=False*)  
Set indicator value for a widget or a branch with the given name

**set\_value** (*name*, *value*)  
Set value of a widget with the given name (None means all values)

**setup** (*name=None*)  
Setup the container by initializing its GUI values and setting the `ctl` attribute

**setup\_name** (*name*)  
Set the object's name

**start** ()  
Start the container.  
  
Starts all the internal timers, and calls `start` method for all the contained widgets.

**start\_timer** (*name*)  
Start the timer with the given name (also called automatically on `start()` method)

**stop** ()  
Stop the container.  
  
Stops all the internal timers, and calls `stop` method for all the contained widgets.

**stop\_timer** (*name*)  
Stop the timer with the given name (also called automatically on `stop()` method)

**update\_indicators** ()  
Update all indicators to represent current values

**update\_value** (*name=None*)  
Send update signal for a handler with a given name or list of names.  
  
Emit a value changed signal with the current value to notify the subscribed slots. If *name* is None, emit for all values in the table.

## pylablib.core.gui.widgets.edit module

**class** pylablib.core.gui.widgets.edit.**TextEdit** (*parent*, *value=None*)  
Bases: `object`  
Expanded text edit.

Maintains internally stored consistent value (which can be, e.g., accessed from different threads).

**keyPressEvent** (*event*)

**set\_expandable** (*left=0, right=0, top=0, bottom=0*)

Make text edit expandable.

If it is expandable, the edit size is expanded by the given size into the corresponding directions. If all are zero, the widget behaves as normal.

**focusInEvent** (*evt*)

**focusOutEvent** (*evt*)

**value\_entered** = <Mock name='mock.QtCore.pyqtSignal()' id='140096333774352'>

Signal emitted when value is entered (regardless of whether it stayed the same)

**value\_changed** = <Mock name='mock.QtCore.pyqtSignal()' id='140096333774352'>

Signal emitted when value is changed

**get\_value** ()

Get current text value

**show\_value** (*interrupt\_edit=False*)

Display currently stored text value

If *interrupt\_edit==True* and the edit is currently being modified by the user, don't update the display.

**set\_value** (*value, notify\_value\_change=True, interrupt\_edit=False*)

Set current text value.

If *notify\_value\_change==True*, emit the *value\_changed* signal; otherwise, change value silently. If *interrupt\_edit==True* and the edit is currently being modified by the user, don't update the display (but still update the internally stored value).

**class** `pylablib.core.gui.widgets.edit.NumEdit` (*parent, value=None, limiter=None, formatter=None, custom\_steps=None*)

Bases: `object`

Labview-style numerical edit.

Maintains internally stored consistent value (which can be, e.g., accessed from different threads). Supports different number representations, metric prefixes (in input or output), keyboard shortcuts (up/down for changing number, escape for cancelling).

#### Parameters

- **parent** – parent widget
- **value** – initial value (*None* means no value is set)
- **limiter** – number limiter (for details, see `set_limiter()`)
- **formatter** – number formatter (for details, see `set_formatter()`)
- **custom\_steps** – if not *None*, can specify custom fixed value steps when up/down keys are pressed with a modifier key (Control, Alt, or Shift) specifies a dictionary `{'ctrl':ctrl_step, 'alt':alt_step, 'shift':shift_step}` with the corresponding steps (missing elements mean that the modifier key is ignored)

**keyPressEvent** (*event*)

**set\_limiter** (*limiter, new\_value=None*)

Change current numerical limiter.

Limiter can be a callable object which takes a single value and either returns a limited value, or raises `limiter.LimitError` if it should be ignored; or it can be a tuple (`lower`, `upper`, `action`, `value_type`), where `lower` and `upper` are the limits (`None` means no limits), `action` defines out-of-limit action (either "ignore" to ignore entered value, or "coerce" to truncate to the nearest limit), and `value_type` can be `None` (keep value as is), "float" (cast value to float), "int" (cast value to int). If the tuple is shorter, the missing parts are filled by default values (`None`, `None`, "ignore", `None`).

**set\_formatter** (*formatter*)

Change current numerical formatter.

Formatter can be a callable object turning value into a string, a string ("float", "int", or a format string, e.g., ".5f"), or a tuple starting with "float" which contains arguments to the `formatter.FloatFormatter`.

**set\_float\_formatter** (*output\_format='auto', digits=9, add\_trailing\_zeros=True, leading\_zeros=0, explicit\_sign=False*)

Set up float formatter.

Has the same functionality as `set_formatter()` (i.e., `set_float_formatter(*args)` is equivalent to `set_formatter(("float",)+args)`), but explicitly lists the arguments.

#### Parameters

- **output\_format** (*str*) – can be "auto" (use standard Python conversion), "SI" (use SI prefixes if possible), or "sci" (scientific "E" notation).
- **digits** (*int*) – if `add_trailing_zeros==False`, determines the number of significant digits; otherwise, determines precision (number of digits after decimal point).
- **add\_trailing\_zeros** (*bool*) – if `True`, always show fixed number of digits after the decimal point, with zero padding if necessary.
- **leading\_zeros** (*bool*) – determines the minimal size of the integer part (before the decimal point) of the number; pads with zeros if necessary.
- **explicit\_sign** (*bool*) – if `True`, always add explicit plus sign.

**set\_custom\_steps** (*custom\_steps=None*)

Specify custom fixed value steps when up/down keys are pressed with a modifier key (Control, Alt, or Shift).

`custom_steps` is a dictionary {'ctrl':ctrl\_step, 'alt':alt\_step, 'shift':shift\_step} with the corresponding steps (missing elements mean that the modifier key is ignored).

**get\_cursor\_order** ()

Get a decimal order of the text cursor

**set\_cursor\_order** (*order*)

Move text cursor to a given decimal order

**repr\_value** (*value*)

Return representation of *value* according to the current numerical format

**value\_entered** = <Mock name='mock.QtCore.pyqtSignal()' id='140096333774352'>  
Signal emitted when value is entered (regardless of whether it stayed the same)

**value\_changed** = <Mock name='mock.QtCore.pyqtSignal()' id='140096333774352'>  
Signal emitted when value is changed

**get\_value()**

Get current numerical value

**show\_value** (*interrupt\_edit=False, preserve\_cursor\_order=True*)

Display currently stored numerical value

If `interrupt_edit==False` and the edit is currently being modified by the user, don't update the display. If `preserve_cursor_order==True` and the display value is being edited, keep the decimal order of the cursor position after change.

**set\_value** (*value, notify\_value\_change=True, interrupt\_edit=False, preserve\_cursor\_order=True*)

Set and display current numerical value.

If `notify_value_change==True`, emit the `value_changed` signal; otherwise, change value silently. If `interrupt_edit==False` and the edit is currently being modified by the user, don't update the display (but still update the internally stored value). If `preserve_cursor_order==True` and the display value is being edited, keep the decimal order of the cursor position after change.

## pylablib.core.gui.widgets.label module

**class** `pylablib.core.gui.widgets.label.TextLabel` (*parent, value=None*)

Bases: `object`

Labview-style text label.

The main difference from the standard `QLabel` is the changed event.

**clicked** = `<Mock name='mock.QtCore.pyqtSignal()' id='140096333774352'>`

**mousePressEvent** (*ev*)

**value\_changed** = `<Mock name='mock.QtCore.pyqtSignal()' id='140096333774352'>`

Signal emitted when value is changed

**get\_value()**

Get current numerical value

**set\_value** (*value*)

Set and display current text value

**class** `pylablib.core.gui.widgets.label.EnumLabel` (*parent, options, value=None, prep=None*)

Bases: `object`

Labview-style label for enumerated values.

Can automatically convert input enum values into corresponding text labels based on the *options* dictionary. Can also specify a function which takes a single value argument and converts into a enum value before checking *options*; useful for “fuzzy” options (e.g., when 0 and `False` mean the same thing)

**clicked** = `<Mock name='mock.QtCore.pyqtSignal()' id='140096333774352'>`

**mousePressEvent** (*ev*)

**set\_out\_of\_range** (*action='error'*)

Set behavior when out-of-range value is applied.

Can be `"error"` (raise error), `"text"` (turn value into text and display it), or `"ignore"` (keep current value).

**set\_options** (*options, value=None, index=None*)

Set new set of options.



If *index\_values* is not `None`, set these as the new index values; otherwise, index values are reset. If *options* is a dictionary, interpret it as a mapping {*option*: *index\_value*}. If *value* is specified, set as the new values. If *index* is specified, use it as the index of a new value; if both *value* and *index* are specified, the *value* takes priority.

**value\_changed** = `<Mock name='mock.QtCore.pyqtSignal()' id='140096333774352'>`  
Signal emitted when value is changed

**get\_value()**  
Get current numerical value

**set\_value(value)**  
Set and display current text value

**repr\_value(value)**  
Return representation of *value* as a combo box text

**class** `pylablib.core.gui.widgets.label.NumLabel` (*parent*, *value=None*, *limiter=None*, *formatter=None*, *allow\_text=True*)

Bases: `object`

Labview-style numerical label.

Supports different number representations and metric prefixes.

#### Parameters

- **parent** – parent widget
- **value** – initial value (`None` means no value is set)
- **limiter** – number limiter (for details, see `set_limiter()`)
- **formatter** – number formatter (for details, see `set_formatter()`)
- **allow\_text** – if `True`, can also take text values (which are displayed as is); otherwise, raise an error.

**clicked** = `<Mock name='mock.QtCore.pyqtSignal()' id='140096333774352'>`

**mousePressEvent** (*ev*)

**set\_limiter** (*limiter*, *new\_value=None*)  
Change current numerical limiter.

Limiter can be a callable object which takes a single value and either returns a limited value, or raises `limiter.LimitError` if it should be ignored; or it can be a tuple (*lower*, *upper*, *action*, *value\_type*), where *lower* and *upper* are the limits (`None` means no limits), *action* defines out-of-limit action (either "ignore" to ignore entered value, or "coerce" to truncate to the nearest limit), and *value\_type* can be `None` (keep value as is), "float" (cast value to float), "int" (cast value to int). If the tuple is shorter, the missing parts are filled by default values (`None`, `None`, "ignore", `None`).

**set\_formatter** (*formatter*)  
Change current numerical formatter.

Formatter can be a callable object turning value into a string, a string ("float", "int", or a format string, e.g., "%.5f"), or a tuple starting with "float" which contains arguments to the `formatter.FloatFormatter`.

**set\_float\_formatter** (*output\_format='auto'*, *digits=9*, *add\_trailing\_zeros=True*, *leading\_zeros=0*, *explicit\_sign=False*)  
Set up float formatter.

Has the same functionality as `set_formatter()` (i.e., `set_float_formatter(*args)` is equivalent to `set_formatter("float",)+args)`), but explicitly lists the arguments.

#### Parameters

- **output\_format** (*str*) – can be "auto" (use standard Python conversion), "SI" (use SI prefixes if possible), or "sci" (scientific “E” notation).
- **digits** (*int*) – if `add_trailing_zeros==False`, determines the number of significant digits; otherwise, determines precision (number of digits after decimal point).
- **add\_trailing\_zeros** (*bool*) – if `True`, always show fixed number of digits after the decimal point, with zero padding if necessary.
- **leading\_zeros** (*bool*) – determines the minimal size of the integer part (before the decimal point) of the number; pads with zeros if necessary.
- **explicit\_sign** (*bool*) – if `True`, always add explicit plus sign.

**repr\_value** (*value*)

Return representation of *value* according to the current numerical format

**value\_changed** = <Mock name='mock.QtCore.pyqtSignal()' id='140096333774352'>

Signal emitted when value is changed

**get\_value** ()

Get current numerical value

**set\_value** (*value*)

Set and display current numerical value

### pylablib.core.gui.widgets.layout\_manager module

**class** pylablib.core.gui.widgets.layout\_manager.**QLayoutManagedWidget** (\*args,  
\*\*kwargs)

Bases: `object`

GUI widget which can manage layouts.

Typically, first it is set up using `setup()` method to specify the master layout kind; afterwards, widgets and sublayout can be added using `add_to_layout()`. In addition, it can directly add named sublayouts using `add_sublayout()` method.

Abstract mix-in class, which needs to be added to a class inheriting from `QWidget`. Alternatively, one can directly use `QLayoutManagedWidget`, which already inherits from `QWidget`.

**setup** (*layout='grid', no\_margins=False*)

Setup the layout.

#### Parameters

- **layout** – layout kind; can be "grid", "vbox" (vertical single-column box), or "hbox" (horizontal single-row box).
- **no\_margins** – if `True`, set all layout margins to zero (useful when the widget is in the middle of layout hierarchy)

**using\_layout** (*name*)

Use a different sublayout as default inside the `with` block

**add\_to\_layout** (*element*, *location=None*, *kind='widget'*)

Add an existing *element* to the layout at the given *location*.

*kind* can be "widget" for widgets, "layout" for other layouts, or "item" for layout items (spacers).

**remove\_layout\_element** (*element*)

Remove a previously added layout element

**get\_element\_position** (*element*)

Get the sublayout and the position of the given widget.

Return tuple (*sublayout*, *location*), where *sublayout* is the sublayout name ("name" for the main layout), and *location* is a tuple (*row*, *column*, *rowspan*, *colspan*). If the given widget is not in this layout, return None.

**add\_sublayout** (*name*, *kind='grid'*, *location=None*)

Add a sublayout to the given location.

*name* specifies the sublayout name, which can be used to refer to it in specifying locations later. *kind* can be "grid", "vbox" (vertical single-column box), or "hbox" (horizontal single-row box).

**using\_new\_sublayout** (*name*, *kind='grid'*, *location=None*)

Create a different sublayout and use it as default inside the *with* block.

*kind* can be "grid", "vbox" (vertical single-column box), or "hbox" (horizontal single-row box).

**get\_sublayout** (*name=None*)

Get the previously added sublayout

**iter\_sublayout\_items** (*name=None*, *include=('widget',)*, *nested=False*)

Iterate over items contained in a given sublayout.

*include* is a tuple which contains items to iterate over; can include "widget" or "layout". If *nested==True*, iterate over items in contained layouts as well.

**get\_sublayout\_kind** (*name=None*)

Get the kind of the previously added sublayout

**get\_layout\_shape** (*name=None*)

Get shape (*rows*, *cols*) of the current layout

**add\_spacer** (*height=0*, *width=0*, *stretch\_height=False*, *stretch\_width=False*, *stretch=0*, *location='next'*)

Add a spacer with the given width and height to the given location.

If *stretch\_height==True* or *stretch\_width==True*, the widget will stretch in these directions; otherwise, the widget size is fixed. If *stretch* is not None, it specifies stretch of the spacer the corresponding direction (applied to the upper row and leftmost column for multi-cell spacer); if *kind=="both"*, it can also be a tuple with two stretches along vertical and horizontal directions.

**add\_padding** (*kind='auto'*, *location='next'*, *stretch=0*)

Add a padding (expandable spacer) of the given kind to the given location.

*kind* can be "vertical", "horizontal", "auto" (vertical for *grid* and *vbox* layouts, horizontal for *hbox*), or "both" (stretches in both directions). If *stretch* is not None, it specifies stretch of the spacer the corresponding direction (applied to the upper row and leftmost column for multi-cell spacer); can also be a tuple with two stretches along vertical and horizontal directions.

**set\_row\_stretch** (*\*args*, *layout=None*)

Set row stretch for a given layout.

Takes either two arguments *index* and *stretch*, or a single list of stretches for all rows.

**set\_column\_stretch** (\*args, layout=None)

Set column stretch for a given layout.

Takes either two arguments `index` and `stretch`, or a single list of stretches for all columns.

**add\_decoration\_label** (text, location='next')

Add a decoration text label with the given text

**insert\_row** (row, sublayout=None, stretch=0)

Insert a new row at the given location in the grid layout

**insert\_column** (col, sublayout=None, stretch=0)

Insert a new column at the given location in the grid layout

**clear** ()

Clear the layout and remove all the added elements

**class** pylablib.core.gui.widgets.layout\_manager.QLayoutManagedWidget (\*args,  
\*\*kwargs)

Bases: `pylablib.core.gui.widgets.layout_manager.IQLayoutManagedWidget`, `object`

GUI widget which can manage layouts.

Typically, first it is set up using `setup()` method to specify the master layout kind; afterwards, widgets and sublayout can be added using `add_to_layout()`. In addition, it can directly add named sublayouts using `add_sublayout()` method.

Simply a combination of `IQLayoutManagedWidget` and `QWidget`.

**add\_decoration\_label** (text, location='next')

Add a decoration text label with the given text

**add\_padding** (kind='auto', location='next', stretch=0)

Add a padding (expandable spacer) of the given kind to the given location.

`kind` can be "vertical", "horizontal", "auto" (vertical for grid and vbox layouts, horizontal for hbox), or "both" (stretches in both directions). If `stretch` is not None, it specifies stretch of the spacer the corresponding direction (applied to the upper row and leftmost column for multi-cell spacer); can also be a tuple with two stretches along vertical and horizontal directions.

**add\_spacer** (height=0, width=0, stretch\_height=False, stretch\_width=False, stretch=0, location='next')

Add a spacer with the given width and height to the given location.

If `stretch_height==True` or `stretch_width==True`, the widget will stretch in these directions; otherwise, the widget size is fixed. If `stretch` is not None, it specifies stretch of the spacer the corresponding direction (applied to the upper row and leftmost column for multi-cell spacer); if `kind=="both"`, it can also be a tuple with two stretches along vertical and horizontal directions.

**add\_sublayout** (name, kind='grid', location=None)

Add a sublayout to the given location.

`name` specifies the sublayout name, which can be used to refer to it in specifying locations later. `kind` can be "grid", "vbox" (vertical single-column box), or "hbox" (horizontal single-row box).

**add\_to\_layout** (element, location=None, kind='widget')

Add an existing `element` to the layout at the given `location`.

`kind` can be "widget" for widgets, "layout" for other layouts, or "item" for layout items (spacers).

**clear** ()

Clear the layout and remove all the added elements

**get\_element\_position** (*element*)

Get the sublayout and the position of the given widget.

Return tuple (*sublayout*, *location*), where *sublayout* is the sublayout name ("name" for the main layout), and *location* is a tuple (*row*, *column*, *rowspan*, *colspan*). If the given widget is not in this layout, return None.

**get\_layout\_shape** (*name=None*)

Get shape (*rows*, *cols*) of the current layout

**get\_sublayout** (*name=None*)

Get the previously added sublayout

**get\_sublayout\_kind** (*name=None*)

Get the kind of the previously added sublayout

**insert\_column** (*col*, *sublayout=None*, *stretch=0*)

Insert a new column at the given location in the grid layout

**insert\_row** (*row*, *sublayout=None*, *stretch=0*)

Insert a new row at the given location in the grid layout

**iter\_sublayout\_items** (*name=None*, *include=('widget',)*, *nested=False*)

Iterate over items contained in a given sublayout.

*include* is a tuple which contains items to iterate over; can include "widget" or "layout". If *nested==True*, iterate over items in contained layouts as well.

**remove\_layout\_element** (*element*)

Remove a previously added layout element

**set\_column\_stretch** (*\*args*, *layout=None*)

Set column stretch for a given layout.

Takes either two arguments *index* and *stretch*, or a single list of stretches for all columns.

**set\_row\_stretch** (*\*args*, *layout=None*)

Set row stretch for a given layout.

Takes either two arguments *index* and *stretch*, or a single list of stretches for all rows.

**setup** (*layout='grid'*, *no\_margins=False*)

Setup the layout.

#### Parameters

- **layout** – layout kind; can be "grid", "vbox" (vertical single-column box), or "hbox" (horizontal single-row box).
- **no\_margins** – if True, set all layout margins to zero (useful when the widget is in the middle of layout hierarchy)

**using\_layout** (*name*)

Use a different sublayout as default inside the `with` block

**using\_new\_sublayout** (*name*, *kind='grid'*, *location=None*)

Create a different sublayout and use it as default inside the `with` block.

*kind* can be "grid", "vbox" (vertical single-column box), or "hbox" (horizontal single-row box).

**pylablib.core.gui.widgets.param\_table module**

**class** pylablib.core.gui.widgets.param\_table.**ParamTable** (parent=None, name=None)

Bases: *pylablib.core.gui.widgets.container.QWidgetContainer*

GUI parameter table.

Simplifies creating code-generated controls and displays table layouts.

Has methods for adding various kinds of controls (labels, edit boxes, combo boxes, check boxes), automatically creates values table for easy settings/getting. By default supports 2-column (label-control) and 3-column (label-control-indicator) layout, depending on the parameters given to *setup()*.

Similar to *GUIValues*, has three container-like accessor: *.h* for getting the value handler (i.e., *self.get\_handler(name)* is equivalent to *self.h[name]*), *.w* for getting the underlying widget (i.e., *self.get\_widget(name)* is equivalent to *self.w[name]*), *.v* for settings/getting values using the default getting method (equivalent to *.wv* if *cache\_values=False* in *setup()*, and to *.cv* otherwise), *.wv* for settings/getting current current widget values without caching (i.e., *self.get\_value(name)* is equivalent to *self.v[name]*, and *self.set\_value(name, value)* is equivalent to *self.v[name]=value*), *.cv* for settings/getting values using cached value's table for getting (i.e., *self.current\_values[name]* is equivalent to *self.cv[name]*, and *self.set\_value(name, value)* is equivalent to *self.cv[name]=value*), (i.e., *self.get\_value(name)* is equivalent to *self.v[name]*, and *self.set\_value(name, value)* is equivalent to *self.v[name]=value*), *.i* for settings/getting indicator values (i.e., *self.get\_indicator(name)* is equivalent to *self.i[name]*, and *self.set\_indicator(name, value)* is equivalent to *self.i[name]=value*) *.vs* for getting the value changed Qt signal (i.e., *self.get\_value\_changed\_signal(name)* is equivalent to *self.s[name]*),

Like most widgets, requires calling *setup()* to set up before usage.

**Parameters** *parent* – parent widget

**setup** (name=None, add\_indicator=True, gui\_thread\_safe=False, cache\_values=False, change\_focused\_control=False)  
Setup the table.

**Parameters**

- **name** (*str*) – table widget name
- **add\_indicator** (*bool*) – if True, add indicators for all added widgets by default.
- **gui\_thread\_safe** (*bool*) – if True, all value-access and indicator-access calls (*get/set\_value*, *get/set\_all\_values*, *get/set\_indicator*, and *update\_indicators*) are automatically called in the GUI thread.
- **cache\_values** (*bool*) – if True or "update\_one", store a dictionary with all the current values and update it every time a GUI value is changed; provides a thread-safe way to check current parameters without lag (unlike *get\_value()* or *get\_all\_values()* with *gui\_thread\_safe==True*, which re-route calls to a GUI thread and may cause up to 100ms delay) can also be set to "update\_all", in which case change of any value will cause value update of all variables; otherwise, change of a value will only cause update of that same value (might potentially miss some value updates for custom controls).
- **change\_focused\_control** (*bool*) – if False and *set\_value()* method is called while the widget has user focus, ignore the value; note that *set\_all\_values()* will still set the widget value.

**add\_sublayout** (name, kind='grid', location=('next', 0, 1, 'end'))

Add a sublayout to the given location.

*name* specifies the sublayout name, which can be used to refer to it in specifying locations later. *kind* can be "grid", "vbox" (vertical single-column box), or "hbox" (horizontal single-row box).

**using\_new\_sublayout** (*name*, *kind*='grid', *location*=('next', 0, 1, 'end'))

Create a different sublayout and use it as default inside the `with` block.

*kind* can be "grid", "vbox" (vertical single-column box), or "hbox" (horizontal single-row box).

**pad\_borders** (*kind*='both', *stretch*=0)

Add expandable paddings on the bottom and/or right border.

*kind* can be "bottom", "right", "both", or "none" (do nothing). Note that if more elements are added, they will be placed after the padding, so the table will be padded in the middle.

**add\_frame** (*name*, *layout*='vbox', *location*=('next', 0, 1, 'end'), *gui\_values\_path*=True, *no\_margins*=True)

Add a new frame container to the layout.

*layout* specifies the layout ("vbox", "hbox", or "grid") of the new frame, and *location* specifies its location within the container layout. If *no\_margins*==True, the frame will have no inner layout margins. The other parameters are the same as in `add_child()` method.

**add\_group\_box** (*name*, *caption*, *layout*='vbox', *location*=('next', 0, 1, 'end'), *gui\_values\_path*=True, *no\_margins*=True)

Add a new group box container with the given *caption* to the layout.

*layout* specifies the layout ("vbox", "hbox", or "grid") of the new frame, and *location* specifies its location within the container layout. If *no\_margins*==True, the frame will have no inner layout margins. The other parameters are the same as in `add_child()` method.

**class ParamRow** (*widget*, *label*, *indicator*, *value\_handler*, *indicator\_handler*)

Bases: `tuple`

**count** ()

Return number of occurrences of value.

**index** ()

Return first index of value.

Raises `ValueError` if the value is not present.

**indicator**

**indicator\_handler**

**label**

**value\_handler**

**widget**

**add\_simple\_widget** (*name*, *widget*, *label*=None, *value\_handler*=None, *add\_indicator*=None, *location*=None, *tooltip*=None, *add\_change\_event*=True)

Add a 'simple' (single-spaced, single-valued) widget to the table.

#### Parameters

- **name** (*str*) – widget name (used to reference its value in the values table)
- **widget** – widget to add
- **label** (*str*) – if not None, specifies label to put in front of the widget in the layout
- **value\_handler** – value handler of the widget; by default, use auto-detected value handler (works for many simple built-in or custom widgets)

- **add\_indicator** – if `True`, add an indicator label in the third column and a corresponding indicator handler in the built-in values table; by default, use the default value supplied to `setup()`
- **location** (*tuple*) – tuple (`row`, `column`) specifying location of the widget (or widget label, if it is specified); by default, add to a new row in the end and into the first column can also be a string `"skip"`, which means that the widget is added to some other location manually later (this option only works if `label=None`, and doesn't add any indicator)
- **tooltip** – widget tooltip (mouseover text)
- **add\_change\_event** (*bool*) – if `True`, changing of the widget's value emits the table's `contained_value_changed` event

Return the widget's value handler

**add\_custom\_widget** (*name*, *widget*, *value\_handler=None*, *indicator\_handler=None*, *location=None*, *tooltip=None*, *add\_change\_event=True*)

Add a 'custom' (multi-spaced, possibly complex-valued) widget to the table.

#### Parameters

- **name** (*str*) – widget name (used to reference its value in the values table)
- **widget** – widget to add
- **value\_handler** – value handler of the widget; by default, use auto-detected value handler (works for many simple built-in or custom widgets)
- **indicator\_handler** – indicator handler of the widget; by default, use auto-detected indicator handler (use `set/get_indicator` methods if present, or no indicator otherwise)
- **location** (*tuple*) – tuple (`row`, `column`, `rowspan`, `colspan`) specifying location of the widget; by default, add to a new row in the end and into the first column, span one row and all table columns can also be a string `"skip"`, which means that the widget is added to some other location manually later
- **add\_change\_event** (*bool*) – if `True`, changing of the widget's value emits the table's `contained_value_changed` event

Return the widget's value handler

**remove\_widget** (*name*)

Remove the widget and, if applicable, its indicator and label

**add\_virtual\_element** (*name*, *value=None*, *multivalued=False*, *add\_indicator=None*)

Add a virtual table element.

Doesn't correspond to any actual widget, but behaves very similarly from the application point of view (its value can be set or read, it has on-change events, it can have indicator). The element value is simply stored on set and retrieved on get. If `multivalued==True`, the internal value is assumed to be complex, so it is forced to be a *Dictionary* every time it is set. If `add_indicator==True`, add default indicator handler as well.

**add\_property\_element** (*name*, *getter=None*, *setter=None*, *add\_indicator=True*)

Add a property value element.

Doesn't correspond to any actual widget, but behaves very similarly from the application point of view; each time the value is set or get, the corresponding setter and getter methods are called. If `add_indicator==True`, add default (stored value) indicator handler as well.



**add\_button** (*name*, *caption*, *label=None*, *add\_indicator=None*, *location=None*, *tooltip=None*, *add\_change\_event=True*, *virtual=False*)

Add a button to the table.

#### Parameters

- **name** (*str*) – widget name (used to reference its value in the values table)
- **caption** (*str*) – text on the button
- **virtual** (*bool*) – if `True`, the widget is not added, and a virtual handler is added instead

Rest of the arguments and the return value are the same as `add_simple_widget()`.

**add\_toggle\_button** (*name*, *caption*, *value=False*, *label=None*, *add\_indicator=None*, *location=None*, *tooltip=None*, *add\_change\_event=True*, *virtual=False*)

Add a toggle button to the table.

#### Parameters

- **name** (*str*) – widget name (used to reference its value in the values table)
- **caption** (*str* or *list*) – text on the button; can be a single string, or a list of two strings which specifies the caption for off and on states
- **value** (*bool*) – specifies initial value
- **virtual** (*bool*) – if `True`, the widget is not added, and a virtual handler is added instead

Rest of the arguments and the return value are the same as `add_simple_widget()`.

**add\_dropdown\_button** (*name*, *caption*, *options=None*, *index\_values=None*, *label=None*, *add\_indicator=None*, *location=None*, *tooltip=None*, *add\_change\_event=True*, *virtual=False*)

Add a button which shows a dropdown menu when clicked.

Similar in behavior to a regular button, but its changed event provides a single argument which is the name of the selected item.

#### Parameters

- **name** (*str*) – widget name (used to reference its value in the values table)
- **caption** (*str* or *list*) – text on the button
- **options** (*list*) – list of strings specifying menu options
- **index\_values** (*list*) – list of values corresponding to menu options; if supplied, these values are used when setting/getting values or sending signals.
- **virtual** (*bool*) – if `True`, the widget is not added, and a virtual handler is added instead

Rest of the arguments and the return value are the same as `add_simple_widget()`.

**add\_check\_box** (*name*, *caption*, *value=False*, *label=None*, *add\_indicator=None*, *location=None*, *tooltip=None*, *add\_change\_event=True*, *virtual=False*)

Add a checkbox to the table.

#### Parameters

- **name** (*str*) – widget name (used to reference its value in the values table)
- **caption** (*str*) – text on the checkbox
- **value** (*bool*) – specifies initial value

- **virtual** (*bool*) – if `True`, the widget is not added, and a virtual handler is added instead

Rest of the arguments and the return value are the same as `add_simple_widget()`.

**add\_text\_label** (*name*, *value*=", *label*=None, *location*=None, *tooltip*=None, *add\_change\_event*=False, *virtual*=False)

Add a text label to the table.

#### Parameters

- **name** (*str*) – widget name (used to reference its value in the values table)
- **value** (*bool*) – specifies initial value
- **virtual** (*bool*) – if `True`, the widget is not added, and a virtual handler is added instead

Rest of the arguments and the return value are the same as `add_simple_widget()`.

**add\_enum\_label** (*name*, *options*, *value*=None, *out\_of\_range*='error', *prep*=None, *label*=None, *location*=None, *tooltip*=None, *add\_change\_event*=False, *virtual*=False)

Add a text label to the table.

#### Parameters

- **name** (*str*) – widget name (used to reference its value in the values table)
- **options** (*list*) – dictionary {*option*: *index\_value*} which converts values into text
- **out\_of\_range** (*str*) – behavior when out-of-range value is applied; can be "error" (raise error), "text" (convert value into text), or "ignore" (keep current value).
- **prep** – a function which takes a single value argument and converts into an option; useful for “fuzzy” options (e.g., when 0 and `False` mean the same thing)
- **virtual** (*bool*) – if `True`, the widget is not added, and a virtual handler is added instead

Rest of the arguments and the return value are the same as `add_simple_widget()`.

**add\_num\_label** (*name*, *value*=0, *limiter*=None, *formatter*=None, *label*=None, *tooltip*=None, *location*=None, *add\_change\_event*=False, *virtual*=False)

Add a numerical label to the table.

#### Parameters

- **name** (*str*) – widget name (used to reference its value in the values table)
- **value** (*float*) – specifies initial value
- **limiter** (*tuple*) – tuple (*upper\_limit*, *lower\_limit*, *action*, *value\_type*) specifying value limits; see `limiter.as_limiter()` for details
- **formatter** (*tuple*) – either "int" (for integer values), or tuple specifying floating value format; see `formatter.as_formatter()` for details
- **virtual** (*bool*) – if `True`, the widget is not added, and a virtual handler is added instead

Rest of the arguments and the return value are the same as `add_simple_widget()`.

**add\_text\_edit** (*name*, *value*=", *label*=None, *multiline*=False, *add\_indicator*=None, *location*=None, *tooltip*=None, *add\_change\_event*=True, *virtual*=False)

Add a text edit to the table.

**Parameters**

- **name** (*str*) – widget name (used to reference its value in the values table)
- **value** (*bool*) – specifies initial value
- **multiline** (*bool*) – if `True`, use multi-line text edit widget; otherwise, use a standard single-line edit
- **virtual** (*bool*) – if `True`, the widget is not added, and a virtual handler is added instead

Rest of the arguments and the return value are the same as `add_simple_widget()`.

**add\_num\_edit** (*name*, *value=None*, *limiter=None*, *formatter=None*, *custom\_steps=None*, *label=None*, *add\_indicator=None*, *location=None*, *tooltip=None*, *add\_change\_event=True*, *virtual=False*)

Add a numerical edit to the table.

**Parameters**

- **name** (*str*) – widget name (used to reference its value in the values table)
- **value** (*bool*) – specifies initial value
- **limiter** (*tuple*) – tuple (*upper\_limit*, *lower\_limit*, *action*, *value\_type*) specifying value limits; see `NumEdit.set_limiter()` for details
- **formatter** (*tuple*) – either `"int"` (for integer values), or tuple specifying floating value format; see `NumEdit.set_formatter()` for details
- **custom\_steps** – if not `None`, can specify custom fixed value steps when up/down keys are pressed with a modifier key (`Control`, `Alt`, or `Shift`) specifies a dictionary `{'ctrl':ctrl_step, 'alt':alt_step, 'shift':shift_step}` with the corresponding steps (missing elements mean that the modifier key is ignored)
- **virtual** (*bool*) – if `True`, the widget is not added, and a virtual handler is added instead

Rest of the arguments and the return value are the same as `add_simple_widget()`.

**add\_progress\_bar** (*name*, *value=None*, *label=None*, *location=None*, *tooltip=None*, *add\_change\_event=True*, *virtual=False*)

Add a progress bar to the table.

**Parameters**

- **name** (*str*) – widget name (used to reference its value in the values table)
- **value** (*bool*) – specifies initial value
- **virtual** (*bool*) – if `True`, the widget is not added, and a virtual handler is added instead

Rest of the arguments and the return value are the same as `add_simple_widget()`.

**add\_combo\_box** (*name*, *value=None*, *options=None*, *index\_values=None*, *out\_of\_range='reset'*, *label=None*, *add\_indicator=None*, *location=None*, *tooltip=None*, *add\_change\_event=True*, *virtual=False*)

Add a combo box to the table.

**Parameters**

- **name** (*str*) – widget name (used to reference its value in the values table)

- **value** – specifies initial value
- **options** (*list*) – list of strings specifying box options or a dictionary {option: index\_value}
- **index\_values** (*list*) – list of values corresponding to box options; if supplied, these values are used when setting/getting values or sending signals; if *options* is a dictionary, this parameter is ignored
- **out\_of\_range** (*str*) – behavior when out-of-range value is applied; can be "error" (raise error), "reset" (reset to no-value position), or "ignore" (keep current value).
- **virtual** (*bool*) – if True, the widget is not added, and a virtual handler is added instead

Rest of the arguments and the return value are the same as `add_simple_widget()`.

**set\_enabled** (*names=None, enabled=True, include\_indicator=True, include\_label=True*)

Enable or disable widgets with the given names (by default, all widgets)

**set\_visible** (*names=None, visible=True, include\_indicator=True, include\_label=True*)

Show or hide widgets with the given names (by default, all widgets)

**get\_value** (*name=None*)

**get\_all\_values** ()

**set\_value** (*name, value, force=True*)

Set value of a widget with the given name.

If *force==True*, force widget value (e.g., ignoring restriction on not changing values of focused widgets)

**set\_all\_values** (*value, force=True*)

Set values of all widgets in the table.

If *force==True*, force widget value (e.g., ignoring restriction on not changing values of focused widgets)

**update\_value** (*name=None*)

**get\_widget** (*name*)

Get a widget corresponding to a value with the given name

**get\_indicator\_widget** (*name*)

Get indicator widget for a parameter with the given name, or `None` if this parameter has no indicator label

**get\_label\_widget** (*name*)

Get label widget for a parameter with the given name, or `None` if this parameter has no label

**get\_child** (*name*)

Get the child widget with the given name

**remove\_child** (*name*)

Remove widget from the container and the layout, clear it, and remove it

**get\_indicator** (*name=None*)

**get\_all\_indicators** ()

**set\_indicator** (*name, value, ignore\_missing=False*)

**set\_all\_indicators** (*value, ignore\_missing=True*)

**update\_indicators** ()

**clear** (*disconnect=False*)

Clear the table (remove all widgets)

If *disconnect==True*, also disconnect all slots connected to the `contained_value_changed` signal.

**TimerUIDGenerator** = `<pylablib.core.utils.general.NamedUIDGenerator object>`

**add\_child** (*name, widget, location=None, gui\_values\_path=True*)

Add a contained child widget.

*name* specifies the child storage name; if *name==False*, only add the widget to the layout, but not to the container. *location* specifies the layout location to which the widget is added; if *location=="skip"*, skip adding it to the layout (can be manually added later). Note that if the widget is added to the layout, it will be completely deleted when `clear` or `remove_child` methods are called; otherwise, simply its `clear` method will be called, and its GUI values will be deleted.

If *gui\_values\_path* is `False` or `None`, do not add it to the GUI values table; if it is `True`, add it under the same root (*path==""*) if it's a container, and under *name* if it's not; otherwise, *gui\_values\_path* specifies the path under which the widget values are stored.

**add\_child\_values** (*name, widget, path, add\_change\_event=True*)

Add child's values to the container's table.

If *widget* is a container and *path==""* or ends in `"/*"` (e.g., `"subpath/*"`), use its `setup_gui_values` to make it share the same GUI values; otherwise, simply add it to the GUI values under the given path. If *add\_change\_event==True*, changing of the widget's value emits the container's `contained_value_changed` event

**add\_decoration\_label** (*text, location='next'*)

Add a decoration text label with the given text

**add\_padding** (*kind='auto', location='next', stretch=0*)

Add a padding (expandable spacer) of the given kind to the given location.

*kind* can be `"vertical"`, `"horizontal"`, `"auto"` (vertical for grid and vbox layouts, horizontal for hbox), or `"both"` (stretches in both directions). If *stretch* is not `None`, it specifies stretch of the spacer the corresponding direction (applied to the upper row and leftmost column for multi-cell spacer); can also be a tuple with two stretches along vertical and horizontal directions.

**add\_spacer** (*height=0, width=0, stretch\_height=False, stretch\_width=False, stretch=0, location='next'*)

Add a spacer with the given width and height to the given location.

If *stretch\_height==True* or *stretch\_width==True*, the widget will stretch in these directions; otherwise, the widget size is fixed. If *stretch* is not `None`, it specifies stretch of the spacer the corresponding direction (applied to the upper row and leftmost column for multi-cell spacer); if *kind=="both"*, it can also be a tuple with two stretches along vertical and horizontal directions.

**add\_timer** (*name, period, autostart=True*)

Add a periodic timer with the given *name* and *period*.

Rarely needs to be called explicitly (one is created automatically if timer event is created). If *autostart==True* and the container has been started (by calling `start()` method), start the timer as well.

**add\_timer\_event** (*name, loop=None, start=None, stop=None, period=None, timer=None, autostart=True*)

Add timer event with the given *name*.

Add an event which should be called periodically (e.g., a GUI update). Internally implemented through Qt timers. *loop*, *start* and *stop* are the functions called, correspondingly, on timer (periodically), when

timer is start, and when it's finished. One can either specify the timer by name (*timer* parameter), or create a new one with the given *period*. If `autostart==True` and the container has been started (by calling `start()` method), start the timer as well.

**add\_to\_layout** (*element*, *location=None*, *kind='widget'*)

Add an existing *element* to the layout at the given *location*.

*kind* can be "widget" for widgets, "layout" for other layouts, or "item" for layout items (spacers).

**contained\_value\_changed** = `<Mock name='mock.QtCore.pyqtSignal()' id='140096333774352'>`

**get\_element\_position** (*element*)

Get the sublayout and the position of the given widget.

Return tuple (*sublayout*, *location*), where *sublayout* is the sublayout name ("name" for the main layout), and *location* is a tuple (row, column, rowspan, colspan). If the given widget is not in this layout, return None.

**get\_handler** (*name*)

Get value handler of a widget with the given name

**get\_layout\_shape** (*name=None*)

Get shape (rows, cols) of the current layout

**get\_sublayout** (*name=None*)

Get the previously added sublayout

**get\_sublayout\_kind** (*name=None*)

Get the kind of the previously added sublayout

**get\_value\_changed\_signal** (*name*)

Get a value-changed signal for a widget with the given name

**insert\_column** (*col*, *sublayout=None*, *stretch=0*)

Insert a new column at the given location in the grid layout

**insert\_row** (*row*, *sublayout=None*, *stretch=0*)

Insert a new row at the given location in the grid layout

**is\_running** ()

Check if the container is running (started and not yet stopped)

**is\_stopping** ()

Check if the container is stopping (stopping initialized and not yet done)

**is\_timer\_running** (*name*)

Check if the timer with the given name is running

**iter\_sublayout\_items** (*name=None*, *include=('widget',)*, *nested=False*)

Iterate over items contained in a given sublayout.

*include* is a tuple which contains items to iterate over; can include "widget" or "layout". If *nested==True*, iterate over items in contained layouts as well.

**remove\_layout\_element** (*element*)

Remove a previously added layout element

**set\_column\_stretch** (*\*args*, *layout=None*)

Set column stretch for a given layout.

Takes either two arguments *index* and *stretch*, or a single list of stretches for all columns.

**set\_row\_stretch** (*\*args*, *layout=None*)

Set row stretch for a given layout.

Takes either two arguments `index` and `stretch`, or a single list of stretches for all rows.

**setup\_name** (*name*)  
Set the object's name

**start** ()  
Start the container.

Starts all the internal timers, and calls `start` method for all the contained widgets.

**start\_timer** (*name*)  
Start the timer with the given name (also called automatically on `start()` method)

**stop** ()  
Stop the container.

Stops all the internal timers, and calls `stop` method for all the contained widgets.

**stop\_timer** (*name*)  
Stop the timer with the given name (also called automatically on `stop()` method)

**using\_layout** (*name*)  
Use a different sublayout as default inside the `with` block

**class** `pylablib.core.gui.widgets.param_table.StatusTable` (*parent=None*,  
*name=None*)  
Bases: `pylablib.core.gui.widgets.param_table.ParamTable`

Expansion of `ParamTable` which adds status lines, which automatically subscribe to signals and update values.

**setup** (*name=None*, *add\_indicator=True*, *gui\_thread\_safe=False*, *cache\_values=False*,  
*change\_focused\_control=False*)  
Setup the table.

#### Parameters

- **name** (*str*) – table widget name
- **add\_indicator** (*bool*) – if `True`, add indicators for all added widgets by default.
- **gui\_thread\_safe** (*bool*) – if `True`, all value-access and indicator-access calls (`get/set_value`, `get/set_all_values`, `get/set_indicator`, and `update_indicators`) are automatically called in the GUI thread.
- **cache\_values** (*bool*) – if `True` or `"update_one"`, store a dictionary with all the current values and update it every time a GUI value is changed; provides a thread-safe way to check current parameters without lag (unlike `get_value()` or `get_all_values()` with `gui_thread_safe==True`, which re-route calls to a GUI thread and may cause up to 100ms delay) can also be set to `"update_all"`, in which case change of any value will cause value update of all variables; otherwise, change of a value will only cause update of that same value (might potentially miss some value updates for custom controls).
- **change\_focused\_control** (*bool*) – if `False` and `set_value()` method is called while the widget has user focus, ignore the value; note that `set_all_values()` will still set the widget value.

**add\_status\_line** (*name*, *label=None*, *srcs=None*, *tags=None*, *filt=None*, *fnt=None*)  
Add a status line to the table:

#### Parameters

- **name** (*str*) – widget name (used to reference its value in the values table)

- **label** (*str*) – if not `None`, specifies label to put in front of the status line
- **srcs** (*list*) – status signal sources
- **tags** (*list*) – status signal tags
- **filt** (*list*) – filter function for the signals
- **fnt** – if not `None`, specifies a function which takes 3 arguments (signal source, tag, and value) and generates a status line text.

**update\_status\_line** (*name*, *thread=None*, *path=None*)

Update status line to the variable with the given *path* from the thread with the given *thread* name.

If *thread* is `None`, use `srcs` name provided upon creation. If *path* is `None`, use `tags` name provided upon creation.

**class ParamRow** (*widget*, *label*, *indicator*, *value\_handler*, *indicator\_handler*)

Bases: `tuple`

**count** ()

Return number of occurrences of value.

**index** ()

Return first index of value.

Raises `ValueError` if the value is not present.

**indicator**

**indicator\_handler**

**label**

**value\_handler**

**widget**

**TimerUIDGenerator** = `<pylablib.core.utils.general.NamedUIDGenerator object>`

**add\_button** (*name*, *caption*, *label=None*, *add\_indicator=None*, *location=None*, *tooltip=None*,  
*add\_change\_event=True*, *virtual=False*)

Add a button to the table.

#### Parameters

- **name** (*str*) – widget name (used to reference its value in the values table)
- **caption** (*str*) – text on the button
- **virtual** (*bool*) – if `True`, the widget is not added, and a virtual handler is added instead

Rest of the arguments and the return value are the same as `add_simple_widget()`.

**add\_check\_box** (*name*, *caption*, *value=False*, *label=None*, *add\_indicator=None*, *location=None*,  
*tooltip=None*, *add\_change\_event=True*, *virtual=False*)

Add a checkbox to the table.

#### Parameters

- **name** (*str*) – widget name (used to reference its value in the values table)
- **caption** (*str*) – text on the checkbox
- **value** (*bool*) – specifies initial value



- **virtual** (*bool*) – if `True`, the widget is not added, and a virtual handler is added instead

Rest of the arguments and the return value are the same as `add_simple_widget()`.

**add\_child** (*name, widget, location=None, gui\_values\_path=True*)

Add a contained child widget.

*name* specifies the child storage name; if *name*==`False`, only add the widget to they layout, but not to the container. *location* specifies the layout location to which the widget is added; if *location*==`"skip"`, skip adding it to the layout (can be manually added later). Note that if the widget is added to the layout, it will be completely deleted when `clear` or `remove_child` methods are called; otherwise, simply its `clear` method will be called, and its GUI values will be deleted.

If *gui\_values\_path* is `False` or `None`, do not add it to the GUI values table; if it is `True`, add it under the same root (*path*==`" "`) if it's a container, and under *name* if it's not; otherwise, *gui\_values\_path* specifies the path under which the widget values are stored.

**add\_child\_values** (*name, widget, path, add\_change\_event=True*)

Add child's values to the container's table.

If *widget* is a container and *path*==`" "` or ends in `"/"` (e.g., `"subpath/"`), use its `setup_gui_values` to make it share the same GUI values; otherwise, simply add it to the GUI values under the given path. if *add\_change\_event*==`True`, changing of the widget's value emits the container's `contained_value_changed` event

**add\_combo\_box** (*name, value=None, options=None, index\_values=None, out\_of\_range='reset', label=None, add\_indicator=None, location=None, tooltip=None, add\_change\_event=True, virtual=False*)

Add a combo box to the table.

#### Parameters

- **name** (*str*) – widget name (used to reference its value in the values table)
- **value** – specifies initial value
- **options** (*list*) – list of strings specifying box options or a dictionary {*option*: *index\_value*}
- **index\_values** (*list*) – list of values corresponding to box options; if supplied, these values are used when setting/getting values or sending signals; if *options* is a dictionary, this parameter is ignored
- **out\_of\_range** (*str*) – behavior when out-of-range value is applied; can be `"error"` (raise error), `"reset"` (reset to no-value position), or `"ignore"` (keep current value).
- **virtual** (*bool*) – if `True`, the widget is not added, and a virtual handler is added instead

Rest of the arguments and the return value are the same as `add_simple_widget()`.

**add\_custom\_widget** (*name, widget, value\_handler=None, indicator\_handler=None, location=None, tooltip=None, add\_change\_event=True*)

Add a 'custom' (multi-spaced, possibly complex-valued) widget to the table.

#### Parameters

- **name** (*str*) – widget name (used to reference its value in the values table)
- **widget** – widget to add

- **value\_handler** – value handler of the widget; by default, use auto-detected value handler (works for many simple built-in or custom widgets)
- **indicator\_handler** – indicator handler of the widget; by default, use auto-detected indicator handler (use `set/get_indicator` methods if present, or no indicator otherwise)
- **location** (*tuple*) – tuple (row, column, rowspan, colspan) specifying location of the widget; by default, add to a new row in the end and into the first column, span one row and all table columns can also be a string "skip", which means that the widget is added to some other location manually later
- **add\_change\_event** (*bool*) – if `True`, changing of the widget's value emits the table's `contained_value_changed` event

Return the widget's value handler

**add\_decoration\_label** (*text*, *location*='next')

Add a decoration text label with the given text

**add\_dropdown\_button** (*name*, *caption*, *options*=None, *index\_values*=None, *label*=None, *add\_indicator*=None, *location*=None, *tooltip*=None, *add\_change\_event*=True, *virtual*=False)

Add a button which shows a dropdown menu when clicked.

Similar in behavior to a regular button, but its changed event provides a single argument which is the name of the selected item.

#### Parameters

- **name** (*str*) – widget name (used to reference its value in the values table)
- **caption** (*str* or *list*) – text on the button
- **options** (*list*) – list of strings specifying menu options
- **index\_values** (*list*) – list of values corresponding to menu options; if supplied, these values are used when setting/getting values or sending signals.
- **virtual** (*bool*) – if `True`, the widget is not added, and a virtual handler is added instead

Rest of the arguments and the return value are the same as `add_simple_widget()`.

**add\_enum\_label** (*name*, *options*, *value*=None, *out\_of\_range*='error', *prep*=None, *label*=None, *location*=None, *tooltip*=None, *add\_change\_event*=False, *virtual*=False)

Add a text label to the table.

#### Parameters

- **name** (*str*) – widget name (used to reference its value in the values table)
- **options** (*list*) – dictionary {option: index\_value} which converts values into text
- **out\_of\_range** (*str*) – behavior when out-of-range value is applied; can be "error" (raise error), "text" (convert value into text), or "ignore" (keep current value).
- **prep** – a function which takes a single value argument and converts into an option; useful for "fuzzy" options (e.g., when 0 and False mean the same thing)
- **virtual** (*bool*) – if `True`, the widget is not added, and a virtual handler is added instead

Rest of the arguments and the return value are the same as `add_simple_widget()`.

**add\_frame** (*name*, *layout*='vbox', *location*=('next', 0, 1, 'end'), *gui\_values\_path*=True, *no\_margins*=True)

Add a new frame container to the layout.

*layout* specifies the layout ("vbox", "hbox", or "grid") of the new frame, and *location* specifies its location within the container layout. If *no\_margins*==True, the frame will have no inner layout margins. The other parameters are the same as in `add_child()` method.

**add\_group\_box** (*name*, *caption*, *layout*='vbox', *location*=('next', 0, 1, 'end'), *gui\_values\_path*=True, *no\_margins*=True)

Add a new group box container with the given *caption* to the layout.

*layout* specifies the layout ("vbox", "hbox", or "grid") of the new frame, and *location* specifies its location within the container layout. If *no\_margins*==True, the frame will have no inner layout margins. The other parameters are the same as in `add_child()` method.

**add\_num\_edit** (*name*, *value*=None, *limiter*=None, *formatter*=None, *custom\_steps*=None, *label*=None, *add\_indicator*=None, *location*=None, *tooltip*=None, *add\_change\_event*=True, *virtual*=False)

Add a numerical edit to the table.

#### Parameters

- **name** (*str*) – widget name (used to reference its value in the values table)
- **value** (*bool*) – specifies initial value
- **limiter** (*tuple*) – tuple (*upper\_limit*, *lower\_limit*, *action*, *value\_type*) specifying value limits; see `NumEdit.set_limiter()` for details
- **formatter** (*tuple*) – either "int" (for integer values), or tuple specifying floating value format; see `NumEdit.set_formatter()` for details
- **custom\_steps** – if not None, can specify custom fixed value steps when up/down keys are pressed with a modifier key (Control, Alt, or Shift) specifies a dictionary {'ctrl':ctrl\_step, 'alt':alt\_step, 'shift':shift\_step} with the corresponding steps (missing elements mean that the modifier key is ignored)
- **virtual** (*bool*) – if True, the widget is not added, and a virtual handler is added instead

Rest of the arguments and the return value are the same as `add_simple_widget()`.

**add\_num\_label** (*name*, *value*=0, *limiter*=None, *formatter*=None, *label*=None, *tooltip*=None, *location*=None, *add\_change\_event*=False, *virtual*=False)

Add a numerical label to the table.

#### Parameters

- **name** (*str*) – widget name (used to reference its value in the values table)
- **value** (*float*) – specifies initial value
- **limiter** (*tuple*) – tuple (*upper\_limit*, *lower\_limit*, *action*, *value\_type*) specifying value limits; see `limiter.as_limiter()` for details
- **formatter** (*tuple*) – either "int" (for integer values), or tuple specifying floating value format; see `formatter.as_formatter()` for details
- **virtual** (*bool*) – if True, the widget is not added, and a virtual handler is added instead

Rest of the arguments and the return value are the same as `add_simple_widget()`.

**add\_padding** (*kind*='auto', *location*='next', *stretch*=0)

Add a padding (expandable spacer) of the given kind to the given location.

*kind* can be "vertical", "horizontal", "auto" (vertical for grid and vbox layouts, horizontal for hbox), or "both" (stretches in both directions). If *stretch* is not None, it specifies stretch of the spacer the corresponding direction (applied to the upper row and leftmost column for multi-cell spacer); can also be a tuple with two stretches along vertical and horizontal directions.

**add\_progress\_bar** (*name*, *value*=None, *label*=None, *location*=None, *tooltip*=None, *add\_change\_event*=True, *virtual*=False)

Add a progress bar to the table.

#### Parameters

- **name** (*str*) – widget name (used to reference its value in the values table)
- **value** (*bool*) – specifies initial value
- **virtual** (*bool*) – if True, the widget is not added, and a virtual handler is added instead

Rest of the arguments and the return value are the same as `add_simple_widget()`.

**add\_property\_element** (*name*, *getter*=None, *setter*=None, *add\_indicator*=True)

Add a property value element.

Doesn't correspond to any actual widget, but behaves very similarly from the application point of view; each time the value is set or get, the corresponding setter and getter methods are called. If `add_indicator==True`, add default (stored value) indicator handler as well.

**add\_simple\_widget** (*name*, *widget*, *label*=None, *value\_handler*=None, *add\_indicator*=None, *location*=None, *tooltip*=None, *add\_change\_event*=True)

Add a 'simple' (single-spaced, single-valued) widget to the table.

#### Parameters

- **name** (*str*) – widget name (used to reference its value in the values table)
- **widget** – widget to add
- **label** (*str*) – if not None, specifies label to put in front of the widget in the layout
- **value\_handler** – value handler of the widget; by default, use auto-detected value handler (works for many simple built-in or custom widgets)
- **add\_indicator** – if True, add an indicator label in the third column and a corresponding indicator handler in the built-in values table; by default, use the default value supplied to `setup()`
- **location** (*tuple*) – tuple (*row*, *column*) specifying location of the widget (or widget label, if it is specified); by default, add to a new row in the end and into the first column can also be a string "skip", which means that the widget is added to some other location manually later (this option only works if `label=None`, and doesn't add any indicator)
- **tooltip** – widget tooltip (mouseover text)
- **add\_change\_event** (*bool*) – if True, changing of the widget's value emits the table's `contained_value_changed` event

Return the widget's value handler

**add\_spacer** (*height=0, width=0, stretch\_height=False, stretch\_width=False, stretch=0, location='next'*)

Add a spacer with the given width and height to the given location.

If *stretch\_height*==True or *stretch\_width*==True, the widget will stretch in these directions; otherwise, the widget size is fixed. If *stretch* is not None, it specifies stretch of the spacer the corresponding direction (applied to the upper row and leftmost column for multi-cell spacer); if *kind*=="both", it can also be a tuple with two stretches along vertical and horizontal directions.

**add\_sublayout** (*name, kind='grid', location=('next', 0, 1, 'end')*)

Add a sublayout to the given location.

*name* specifies the sublayout name, which can be used to refer to it in specifying locations later. *kind* can be "grid", "vbox" (vertical single-column box), or "hbox" (horizontal single-row box).

**add\_text\_edit** (*name, value="", label=None, multiline=False, add\_indicator=None, location=None, tooltip=None, add\_change\_event=True, virtual=False*)

Add a text edit to the table.

#### Parameters

- **name** (*str*) – widget name (used to reference its value in the values table)
- **value** (*bool*) – specifies initial value
- **multiline** (*bool*) – if True, use multi-line text edit widget; otherwise, use a standard single-line edit
- **virtual** (*bool*) – if True, the widget is not added, and a virtual handler is added instead

Rest of the arguments and the return value are the same as `add_simple_widget()`.

**add\_text\_label** (*name, value="", label=None, location=None, tooltip=None, add\_change\_event=False, virtual=False*)

Add a text label to the table.

#### Parameters

- **name** (*str*) – widget name (used to reference its value in the values table)
- **value** (*bool*) – specifies initial value
- **virtual** (*bool*) – if True, the widget is not added, and a virtual handler is added instead

Rest of the arguments and the return value are the same as `add_simple_widget()`.

**add\_timer** (*name, period, autostart=True*)

Add a periodic timer with the given *name* and *period*.

Rarely needs to be called explicitly (one is created automatically if timer event is created). If *autostart*==True and the container has been started (by calling `start()` method), start the timer as well.

**add\_timer\_event** (*name, loop=None, start=None, stop=None, period=None, timer=None, autostart=True*)

Add timer event with the given *name*.

Add an event which should be called periodically (e.g., a GUI update). Internally implemented through Qt timers. *loop*, *start* and *stop* are the functions called, correspondingly, on timer (periodically), when timer is start, and when it's finished. One can either specify the timer by name (*timer* parameter), or create a new one with the given *period*. If *autostart*==True and the container has been started (by calling `start()` method), start the timer as well.

**add\_to\_layout** (*element*, *location*=None, *kind*='widget')

Add an existing *element* to the layout at the given *location*.

*kind* can be "widget" for widgets, "layout" for other layouts, or "item" for layout items (spacers).

**add\_toggle\_button** (*name*, *caption*, *value*=False, *label*=None, *add\_indicator*=None, *location*=None, *tooltip*=None, *add\_change\_event*=True, *virtual*=False)

Add a toggle button to the table.

#### Parameters

- **name** (*str*) – widget name (used to reference its value in the values table)
- **caption** (*str* or *list*) – text on the button; can be a single string, or a list of two strings which specifies the caption for off and on states
- **value** (*bool*) – specifies initial value
- **virtual** (*bool*) – if True, the widget is not added, and a virtual handler is added instead

Rest of the arguments and the return value are the same as `add_simple_widget()`.

**add\_virtual\_element** (*name*, *value*=None, *multivalued*=False, *add\_indicator*=None)

Add a virtual table element.

Doesn't correspond to any actual widget, but behaves very similarly from the application point of view (its value can be set or read, it has on-change events, it can have indicator). The element value is simply stored on set and retrieved on get. If `multivalued==True`, the internal value is assumed to be complex, so it is forced to be a *Dictionary* every time it is set. If `add_indicator==True`, add default indicator handler as well.

**clear** (*disconnect*=False)

Clear the table (remove all widgets)

If `disconnect==True`, also disconnect all slots connected to the `contained_value_changed` signal.

**contained\_value\_changed** = <Mock name='mock.QtCore.pyqtSignal()' id='140096333774352'>

**get\_all\_indicators** ()

**get\_all\_values** ()

**get\_child** (*name*)

Get the child widget with the given name

**get\_element\_position** (*element*)

Get the sublayout and the position of the given widget.

Return tuple (*sublayout*, *location*), where *sublayout* is the sublayout name ("name" for the main layout), and *location* is a tuple (*row*, *column*, *rowspan*, *colspan*). If the given widget is not in this layout, return None.

**get\_handler** (*name*)

Get value handler of a widget with the given name

**get\_indicator** (*name*=None)

**get\_indicator\_widget** (*name*)

Get indicator widget for a parameter with the given name, or None if this parameter has no indicator label

**get\_label\_widget** (*name*)

Get label widget for a parameter with the given name, or None if this parameter has no label

**get\_layout\_shape** (*name=None*)  
Get shape (*rows*, *cols*) of the current layout

**get\_sublayout** (*name=None*)  
Get the previously added sublayout

**get\_sublayout\_kind** (*name=None*)  
Get the kind of the previously added sublayout

**get\_value** (*name=None*)

**get\_value\_changed\_signal** (*name*)  
Get a value-changed signal for a widget with the given name

**get\_widget** (*name*)  
Get a widget corresponding to a value with the given name

**insert\_column** (*col*, *sublayout=None*, *stretch=0*)  
Insert a new column at the given location in the grid layout

**insert\_row** (*row*, *sublayout=None*, *stretch=0*)  
Insert a new row at the given location in the grid layout

**is\_running** ()  
Check if the container is running (started and not yet stopped)

**is\_stopping** ()  
Check if the container is stopping (stopping initialized and not yet done)

**is\_timer\_running** (*name*)  
Check if the timer with the given name is running

**iter\_sublayout\_items** (*name=None*, *include=('widget',)*, *nested=False*)  
Iterate over items contained in a given sublayout.  
  
*include* is a tuple which contains items to iterate over; can include "widget" or "layout". If *nested==True*, iterate over items in contained layouts as well.

**pad\_borders** (*kind='both'*, *stretch=0*)  
Add expandable paddings on the bottom and/or right border.  
  
*kind* can be "bottom", "right", "both", or "none" (do nothing). Note that if more elements are added, they will be placed after the padding, so the table will be padded in the middle.

**remove\_child** (*name*)  
Remove widget from the container and the layout, clear it, and remove it

**remove\_layout\_element** (*element*)  
Remove a previously added layout element

**remove\_widget** (*name*)  
Remove the widget and, if applicable, its indicator and label

**set\_all\_indicators** (*value*, *ignore\_missing=True*)

**set\_all\_values** (*value*, *force=True*)  
Set values of all widgets in the table.  
  
If *force==True*, force widget value (e.g., ignoring restriction on not changing values of focused widgets)

**set\_column\_stretch** (*\*args*, *layout=None*)  
Set column stretch for a given layout.  
  
Takes either two arguments *index* and *stretch*, or a single list of stretches for all columns.

**set\_enabled** (*names=None, enabled=True, include\_indicator=True, include\_label=True*)

Enable or disable widgets with the given names (by default, all widgets)

**set\_indicator** (*name, value, ignore\_missing=False*)

**set\_row\_stretch** (*\*args, layout=None*)

Set row stretch for a given layout.

Takes either two arguments `index` and `stretch`, or a single list of stretches for all rows.

**set\_value** (*name, value, force=True*)

Set value of a widget with the given name.

If `force==True`, force widget value (e.g., ignoring restriction on not changing values of focused widgets)

**set\_visible** (*names=None, visible=True, include\_indicator=True, include\_label=True*)

Show or hide widgets with the given names (by default, all widgets)

**setup\_name** (*name*)

Set the object's name

**start** ()

Start the container.

Starts all the internal timers, and calls `start` method for all the contained widgets.

**start\_timer** (*name*)

Start the timer with the given name (also called automatically on `start()` method)

**stop** ()

Stop the container.

Stops all the internal timers, and calls `stop` method for all the contained widgets.

**stop\_timer** (*name*)

Stop the timer with the given name (also called automatically on `stop()` method)

**update\_indicators** ()

**update\_value** (*name=None*)

**using\_layout** (*name*)

Use a different sublayout as default inside the `with` block

**using\_new\_sublayout** (*name, kind='grid', location=('next', 0, 1, 'end')*)

Create a different sublayout and use it as default inside the `with` block.

*kind* can be "grid", "vbox" (vertical single-column box), or "hbox" (horizontal single-row box).

## Module contents

### Submodules

#### pylablib.core.gui.formatter module

`pylablib.core.gui.formatter.parse_float` (*s*)

Parse string as a float, with metric prefixes recognition.

Return tuple (*sign, integer, dot, fractional, exponent, prefix*), where each entry has structure (*begin, end, text*). Return `None` if string is unrecognizable.



`pylablib.core.gui.formatter.pos_to_order(s, pos)`

For a given string representation of a float and position in the string, get the decimal order for this position.

Return None if string is un-parsable or position is out of range (not in mantissa section of the number).

`pylablib.core.gui.formatter.order_to_pos(s, order)`

For a given string representation of float and decimal order, get the position in the string corresponding to this order.

If order is out of range for a given representation, truncates to most/least significant digit position. Return None if string is un-parsable.

`pylablib.core.gui.formatter.str_to_float(s)`

Return float value of a string, with metric prefixes recognition.

Raise `ValueError` if string is unrecognizable.

`pylablib.core.gui.formatter.is_integer(n, tolerance=0.0)`

Check if *n* is less than *tolerance* away from the nearest integer.

`pylablib.core.gui.formatter.float_to_str_SI(n, digits=3, trailing_zeros=False)`

Represent float using SI metric prefixes.

For orders  $\geq 27$  and  $\leq -24$  use usual scientific notation with order being multiple of 3. If `trailing_zeros==True`, then digits define precision, rather than number significant digits

```
class pylablib.core.gui.formatter.FloatFormatter(output_format='auto',      digits=3,
                                                add_trailing_zeros=True,    leading_zeros=0, explicit_sign=False)
```

Bases: `object`

Floating point number formatter.

Callable object with takes a number as an argument and returns its string representation.

#### Parameters

- **output\_format** (*str*) – can be "auto" (use standard Python conversion), "SI" (use SI prefixes if possible), or "sci" (scientific “E” notation).
- **digits** (*int*) – if `add_trailing_zeros==False`, determines the number of significant digits; otherwise, determines precision (number of digits after decimal point).
- **add\_trailing\_zeros** (*bool*) – if `True`, always show fixed number of digits after the decimal point, with zero padding if necessary.
- **leading\_zeros** (*bool*) – determines the minimal size of the integer part (before the decimal point) of the number; pads with zeros if necessary.
- **explicit\_sign** (*bool*) – if `True`, always add explicit plus sign.

```
class pylablib.core.gui.formatter.IntegerFormatter
```

Bases: `object`

Simple integer number formatter.

Callable object with takes a number as an argument and returns its string representation.

For more flexibility (e.g., adding leading zeros) it is possible to use `FloatFormatter` with `digits=0` and `add_trailing_zeros=True`.

```
class pylablib.core.gui.formatter.FmtStringFormatter(fmt)
```

Bases: `object`

Formatter based on format string.

Callable object with takes a number as an argument and returns is string representation.

`pylablib.core.gui.formatter.as_formatter(formatter)`

Turn an object into a formatter.

Can be a callable object turning value into a string, a string ("float", "int", or a format string, e.g., ".5f"), or a tuple starting with "float" which contains arguments to the `FloatFormatter`.

## pylablib.core.gui.limiter module

**exception** `pylablib.core.gui.limiter.LimitError`(*value*, *lower\_limit=None*, *upper\_limit=None*)

Bases: `ArithmeticError`

Error raised when the value is out of limits and can't be coerced

**args**

**with\_traceback()**

Exception.with\_traceback(tb) – set self.\_\_traceback\_\_ to tb and return self.

**class** `pylablib.core.gui.limiter.NumberLimit`(*lower\_limit=None*, *upper\_limit=None*, *action='coerce'*, *value\_type=None*)

Bases: `object`

Number limiter, which checks validity of user inputs.

Callable object with takes a number as an argument and either returns its coerced version (or the number itself, if it is within limits), or raises `LimitError` if it should be ignored.

### Parameters

- **lower\_limit** – lower limit (inclusive), or `None` if there is no limit.
- **upper\_limit** – upper limit (inclusive), or `None` if there is no limit.
- **action** (*str*) – action taken if the number is out of limits; either "coerce" (return the closest valid value), or "ignore" (raise `LimitError`).
- **value\_type** (*str*) – determines value type coercion; can be `None` (do nothing, only check limits), "float" (cast to float), or "int" (cast to integer).

**cast** (*value*)

`pylablib.core.gui.limiter.filter_limiter(pred)`

Turn a predicate into a limiter.

Returns a function that raises `LimitError` if the predicate is false.

`pylablib.core.gui.limiter.as_limiter(limiter)`

Turn an object into a limiter.

Limiter can be a callable object which takes a single value and either returns a limited value, or raises `LimitError` if it should be ignored; or it can be a tuple (*lower*, *upper*, *action*, *value\_type*), where *lower* and *upper* are the limits (`None` means no limits), *action* defines out-of-limit action (either "ignore" to ignore entered value, or "coerce" to truncate to the nearest limit), and *value\_type* can be `None` (keep value as is), "float" (cast value to float), "int" (cast value to int). If the tuple is shorter, the missing parts are filled by default values (`None`, `None`, "ignore", `None`).

**pylablib.core.gui.utils module**

`pylablib.core.gui.utils.get_top_parent (widget)`

Find the top-level parent (parent which does not have further parents)

`pylablib.core.gui.utils.find_layout_element (layout, element)`

Find a layout element.

Can be a widget, a sublayout, or a layout element Return item index within the layout. If layout is empty or item is not present, return None

`pylablib.core.gui.utils.delete_layout_item (layout, idx)`

Remove and item with the given index (completely delete it)

`pylablib.core.gui.utils.clean_layout (layout, delete_layout=False)`

Delete all items from the layout.

If `delete_layout==True`, delete the layout as well.

`pylablib.core.gui.utils.get_layout_container (widget, top=None, kind='widget')`

Find a container widget or layout which contains the given widget.

Note that the container widget does not necessarily correspond to the element parent. If no container could be found, return None. If *kind* can be either "widget" (return the containing widget), or "layout" (return the containing layout, which is a layout or sublayout of the containing widget).

This method works by traversing the whole layout tree, so it can be relatively slow. *top* can specify the top container (widget or layout) which definitely contains the given widget; if not specified, use the top-level parent found by `get_top_parent()`.

`pylablib.core.gui.utils.get_all_layout_containers (widget, top=None, kind='widget')`

Get a list of all widgets or layouts containing the current widget.

The list is arranged from the bottom of the hierarchy (starting from *widget*) to the *top*. Note that the container widget does not necessarily correspond to the element parent. If no containers could be found, return None. If *kind* can be either "widget" (return the containing widgets), or "layout" (return the containing layouts, which are layouts or sublayouts of the containing widgets).

This method works by traversing the whole layout tree, so it can be relatively slow. *top* can specify the top container (widget or layout) which definitely contains the given widget; if not specified, use the top-level parent found by `get_top_parent()`.

`pylablib.core.gui.utils.delete_widget (widget)`

Remove widget from its layout container and delete it

**class** `pylablib.core.gui.utils.TWidgetLocation (layout, position)`

Bases: `tuple`

**count** ()

Return number of occurrences of value.

**index** ()

Return first index of value.

Raises ValueError if the value is not present.

**layout**

**position**

`pylablib.core.gui.utils.get_widget_location (widget, layout=None)`

Get location of a widget within the given layout.

Return tuple (layout, position), where layout is the layout object, and position is either a single position number (for box layouts), or a tuple (row, col, rowspan, colspan) for a grid layout. If layout is not specified, autodetect it.

`pylablib.core.gui.utils.place_widget_at_location(widget, location)`

Insert a widget within the given layout location.

*location* is a tuple tuple (layout, position), where layout is the layout object, and position is either a single position number (for box layouts), or a tuple (row, col, rowspan, colspan) for a grid layout. The tuple has the same format as returned by `get_widget_location()`.

`pylablib.core.gui.utils.is_layout_row_empty(layout, row)`

Check if the given row in a grid layout is empty

`pylablib.core.gui.utils.get_last_filled_row(layout, start_row=0)`

Find the last non-empty row in a grid layout after *start\_row* (inclusive).

If all rows after (and including) *start\_row* are empty, return None .

`pylablib.core.gui.utils.get_first_empty_row(layout, start_row=0)`

Find the first completely empty row in a grid layout after *start\_row* (inclusive)

`pylablib.core.gui.utils.insert_layout_row(layout, row, stretch=0, compress=False)`

Insert row in a grid layout at a given index.

Any multi-column item spanning over the row (i.e., starting at least one row before *row* and ending at least one row after *row*) gets stretched. Anything else either stays in place (if it's above *row*), or gets moved one row down. *stretch* determines the stretch factor of the new row. If `compress==True`, try to find an empty row below the inserted position and shift it to the new row's place; otherwise, add a completely new row.

`pylablib.core.gui.utils.is_layout_column_empty(layout, col)`

Check if the given column in a grid layout is empty

`pylablib.core.gui.utils.get_last_filled_column(layout, start_col=0)`

Find the last non-empty column in a grid layout after *start\_col* (inclusive).

If all rows after (and including) *start\_col* are empty, return None .

`pylablib.core.gui.utils.get_first_empty_column(layout, start_col=0)`

Find the first completely empty column in a grid layout after *start\_col* (inclusive)

`pylablib.core.gui.utils.insert_layout_column(layout, col, stretch=0, compress=False)`

Insert column in a grid layout at a given index.

Any multi-row item spanning over the column (i.e., starting at least one column before *col* and ending at least one column after *col*) gets stretched. Anything else either stays in place (if it's above *col*), or gets moved one column to the right. *stretch* determines the stretch factor of the new column. If `compress==True`, try to find an empty column below the inserted position and shift it to the new column's place; otherwise, add a completely new column.

`pylablib.core.gui.utils.compress_grid_layout(layout)`

Find all empty rows in a grid layout and shift them to the bottom

`pylablib.core.gui.utils.get_relative_position(widget, origin=None)`

Get widget's position relative to the origin (top-level parent if None)

`pylablib.core.gui.utils.get_relative_rectangle(widget, origin=None, border=0, trim_border=True)`

Get widget rectangle area relative to the origin (top-level parent if None).

If *border* is non-zero, it specifies a border (integer or 2-tuple) around the widget to add to the rectangle. If `trim_border==True`, the resulting rectangle is trimmed to lie within the origin area. Return `QRect` object.

`pylablib.core.gui.utils.get_screenshot` (*window=None, rect=None, widget=None, border=0, include\_titlebar=True*)

Take a screenshot of a given window or a given widget.

Either *window* or *widget* must be defined. If *rect* (type `QRect`) or *widget* are defined, they specify the area to include into screenshot; in this case, *border* can define an additional border to add to the rectangle. If rectangle is not defined, then *include\_titlebar* specifies whether the window titlebar is included.

## pylablib.core.gui.value\_handling module

Uniform representation of values from different widgets: numerical and text edits and labels, combo and check boxes, buttons.

`pylablib.core.gui.value_handling.build_children_tree` (*root, types\_include, is\_atomic=None, is\_excluded=None, self\_node='#'*)

`pylablib.core.gui.value_handling.has_methods` (*widget, methods\_sets*)

Check if the widget has methods from given set.

*methods\_sets* is a list of method sets. The function returns `True` if the widget has at least one method from each of the sets.

`pylablib.core.gui.value_handling.get_method_kind` (*method, add\_args=0*)

Determine whether the method takes name as its argument

*add\_args* specifies number of additional required arguments. Return "named" if the method has at least *add\_args*+1 arguments, and the first one is called "name". Otherwise, return "simple".

**exception** `pylablib.core.gui.value_handling.NoParameterError`

Bases: `KeyError`

Error raised by some handlers to indicate that the parameter is missing

**args**

**with\_traceback** ()

Exception.with\_traceback(tb) – set self.\_\_traceback\_\_ to tb and return self.

**class** `pylablib.core.gui.value_handling.IValueHandler` (*widget*)

Bases: `object`

Generic handler of a widget value.

Has method to get and set the value (or all values, if the widget has internal value structure), representing values as strings, and value changed signal.

**Parameters** *widget* – handled widget.

**get\_value** (*name=None*)

Get widget value.

If *name* is not `None`, it specifies the name of the value parameter inside the widget (for complex widgets).

**set\_value** (*value, name=None*)

Set widget value.

If *name* is not `None`, it specifies the name of the value parameter inside the widget (for complex widgets).

**repr\_value** (*value, name=None*)

Return textual representation of the value.

If `name` is not `None`, it specifies the name of the value parameter inside the widget (for complex widgets).

**get\_handler** (*name=None*)

Get handler of a contained widget (or same widget, if `name==None`)

**get\_value\_changed\_signal** ()

Get the Qt signal emitted when the value is changed

**connect\_value\_changed\_handler** (*handler, only\_signal=True*)

Connect value changed signal.

If `only_signal==True`, equivalent to connecting a handler function to `get_value_changed_signal()` signal; however, if `only_signal==False`, it also works for some objects (e.g., `QLabel`) don't have built-in on-changed signals by calling the handler explicitly every time the value is changed.

Note that the connection is always direct (i.e., it doesn't deal with message queues and different threads, but rather directly calls the handler function). If you need to connect a handler to a signal using some other connection method, you can use `get_value_changed_signal()` directly.

**can\_set\_value** (*allow\_focus=True*)

Check if setting value from the code is allowed.

**Parameters** **focus** – if `False`, indicates that settings of focused widgets isn't allowed, with some exceptions (buttons, check boxes, combo boxes)

**class** `pylablib.core.gui.value_handling.VirtualValueHandler` (*value=None, multivalued=False*)

Bases: `pylablib.core.gui.value_handling.IValueHandler`

Virtual value handler (to simulate controls which are not present in the GUI).

#### Parameters

- **value** – initial value
- **multivalued** (*bool*) – if `True`, the internal value is assumed to be complex, so it is forced to be a *Dictionary* every time it is set.

**get\_value** (*name=None*)

Get widget value.

If `name` is not `None`, it specifies the name of the value parameter inside the widget (for complex widgets).

**set\_value** (*value, name=None*)

Set widget value.

If `name` is not `None`, it specifies the name of the value parameter inside the widget (for complex widgets).

**can\_set\_value** (*allow\_focus=True*)

Check if setting value from the code is allowed.

**Parameters** **focus** – if `False`, indicates that settings of focused widgets isn't allowed, with some exceptions (buttons, check boxes, combo boxes)

**connect\_value\_changed\_handler** (*handler, only\_signal=True*)

Connect value changed signal.

If `only_signal==True`, equivalent to connecting a handler function to `get_value_changed_signal()` signal; however, if `only_signal==False`, it also works for some objects (e.g., `QLabel`) don't have built-in on-changed signals by calling the handler explicitly every time the value is changed.

Note that the connection is always direct (i.e., it doesn't deal with message queues and different threads, but rather directly calls the handler function). If you need to connect a handler to a signal using some other connection method, you can use `get_value_changed_signal()` directly.

**get\_handler** (*name=None*)

Get handler of a contained widget (or same widget, if *name==None*)

**get\_value\_changed\_signal** ()

Get the Qt signal emitted when the value is changed

**repr\_value** (*value, name=None*)

Return textual representation of the value.

If *name* is not *None*, it specifies the name of the value parameter inside the widget (for complex widgets).

```
class pylablib.core.gui.value_handling.PropertyValueHandler (getter=None, set-
ter=None, de-
fault_name=None)
```

Bases: `pylablib.core.gui.value_handling.IValueHandler`

Virtual value handler which uses custom getter/setter methods to simulate a value.

If getter or setter are not supplied but are called, they raise `NoParameterError`; this means that they are ignored in `GUIValues.get_all_values()` and `GUIValues.set_all_values()` methods, but raise an error when access directly (e.g., using `GUIValues.get_value()`).

#### Parameters

- **getter** – value getter method; takes 0 or 1 (*name*) arguments and returns the value
- **setter** – value setter method; takes 1 (*value*) or 2 (*name* and *value*) arguments and sets the value
- **default\_name** (*str*) – default name to be supplied to getter and setter methods if they require a name argument

**get\_value** (*name=None*)

Get widget value.

If *name* is not *None*, it specifies the name of the value parameter inside the widget (for complex widgets).

**set\_value** (*value, name=None*)

Set widget value.

If *name* is not *None*, it specifies the name of the value parameter inside the widget (for complex widgets).

**can\_set\_value** (*allow\_focus=True*)

Check if setting value from the code is allowed.

**Parameters** **focus** – if *False*, indicates that settings of focused widgets isn't allowed, with some exceptions (buttons, check boxes, combo boxes)

**connect\_value\_changed\_handler** (*handler, only\_signal=True*)

Connect value changed signal.

If *only\_signal==True*, equivalent to connecting a handler function to `get_value_changed_signal()` signal; however, if *only\_signal==False*, it also works for some objects (e.g., `QLabel`) don't have built-in on-changed signals by calling the handler explicitly every time the value is changed.

Note that the connection is always direct (i.e., it doesn't deal with message queues and different threads, but rather directly calls the handler function). If you need to connect a handler to a signal using some other connection method, you can use `get_value_changed_signal()` directly.



**get\_handler** (*name=None*)

Get handler of a contained widget (or same widget, if *name==None*)

**get\_value\_changed\_signal** ()

Get the Qt signal emitted when the value is changed

**repr\_value** (*value, name=None*)

Return textual representation of the value.

If *name* is not *None*, it specifies the name of the value parameter inside the widget (for complex widgets).

**class** `pylablib.core.gui.value_handling.StandardValueHandler` (*widget, de-  
fault\_name=None*)

Bases: `pylablib.core.gui.value_handling.IValueHandler`

Standard value handler, typically used for custom widgets.

To implement getting and setting values, looks for `get/set_value` and `get/set_all_values` methods for the widget and uses them accordingly. To implement value representing, looks for `repr_value` method (if not defined, use simple string conversion). To implement value change signal, looks for `value_changed` widget signal.

#### Parameters

- **widget** – handled widget
- **default\_name** (*str*) – default name to be supplied to `get/set_value` and `get/set_all_values` methods if they require a name argument.

**get\_value** (*name=None*)

Get widget value.

If *name* is not *None*, it specifies the name of the value parameter inside the widget (for complex widgets).

**set\_value** (*value, name=None*)

Set widget value.

If *name* is not *None*, it specifies the name of the value parameter inside the widget (for complex widgets).

**repr\_value** (*value, name=None*)

Return textual representation of the value.

If *name* is not *None*, it specifies the name of the value parameter inside the widget (for complex widgets).

**get\_handler** (*name=None*)

Get handler of a contained widget (or same widget, if *name==None*)

**can\_set\_value** (*allow\_focus=True*)

Check if setting value from the code is allowed.

**Parameters** **focus** – if *False*, indicates that settings of focused widgets isn't allowed, with some exceptions (buttons, check boxes, combo boxes)

**connect\_value\_changed\_handler** (*handler, only\_signal=True*)

Connect value changed signal.

If *only\_signal==True*, equivalent to connecting a handler function to `get_value_changed_signal()` signal; however, if *only\_signal==False*, it also works for some objects (e.g., `QLabel`) don't have built-in on-changed signals by calling the handler explicitly every time the value is changed.

Note that the connection is always direct (i.e., it doesn't deal with message queues and different threads, but rather directly calls the handler function). If you need to connect a handler to a signal using some other connection method, you can use `get_value_changed_signal()` directly.



**get\_value\_changed\_signal()**

Get the Qt signal emitted when the value is changed

**class** pylablib.core.gui.value\_handling.**ISingleValueHandler** (*widget*)

Bases: *pylablib.core.gui.value\_handling.IValueHandler*

Base class for single-value widget handler, typically used for built-in Qt widgets.

Defines new functions `get/set_single_value` which don't take a name argument; raises an error if the name is supplied to any of the standard functions.

**Parameters** **widget** – handled widget

**get\_single\_value()**

Get the widget value

**get\_value** (*name=None*)

Get widget value.

If name is not None raise an error.

**set\_single\_value** (*value*)

Set the widget value

**set\_value** (*value, name=None*)

Set widget value.

If name is not None raise an error.

**repr\_single\_value** (*value*)

Represent the widget value as a string

**repr\_value** (*value, name=None*)

Return textual representation of the value.

If name is not None raise an error.

**can\_set\_value** (*allow\_focus=True*)

Check if setting value from the code is allowed.

**Parameters** **focus** – if False, indicates that settings of focused widgets isn't allowed, with some exceptions (buttons, check boxes, combo boxes)

**connect\_value\_changed\_handler** (*handler, only\_signal=True*)

Connect value changed signal.

If `only_signal==True`, equivalent to connecting a handler function to `get_value_changed_signal()` signal; however, if `only_signal==False`, it also works for some objects (e.g., `QLabel`) don't have built-in on-changed signals by calling the handler explicitly every time the value is changed.

Note that the connection is always direct (i.e., it doesn't deal with message queues and different threads, but rather directly calls the handler function). If you need to connect a handler to a signal using some other connection method, you can use `get_value_changed_signal()` directly.

**get\_handler** (*name=None*)

Get handler of a contained widget (or same widget, if `name==None`)

**get\_value\_changed\_signal()**

Get the Qt signal emitted when the value is changed

**class** pylablib.core.gui.value\_handling.**LineEditValueHandler** (*widget*)

Bases: *pylablib.core.gui.value\_handling.ISingleValueHandler*

Value handler for `QLineEdit` widget

**get\_single\_value()**  
Get the widget value

**set\_single\_value(value)**  
Set the widget value

**get\_value\_changed\_signal()**  
Get the Qt signal emitted when the value is changed

**can\_set\_value(allow\_focus=True)**  
Check if setting value from the code is allowed.

**Parameters focus** – if `False`, indicates that settings of focused widgets isn't allowed, with some exceptions (buttons, check boxes, combo boxes)

**connect\_value\_changed\_handler(handler, only\_signal=True)**  
Connect value changed signal.

If `only_signal==True`, equivalent to connecting a handler function to `get_value_changed_signal()` signal; however, if `only_signal==False`, it also works for some objects (e.g., `QLabel`) don't have built-in on-changed signals by calling the handler explicitly every time the value is changed.

Note that the connection is always direct (i.e., it doesn't deal with message queues and different threads, but rather directly calls the handler function). If you need to connect a handler to a signal using some other connection method, you can use `get_value_changed_signal()` directly.

**get\_handler(name=None)**  
Get handler of a contained widget (or same widget, if `name==None`)

**get\_value(name=None)**  
Get widget value.

If `name` is not `None` raise an error.

**repr\_single\_value(value)**  
Represent the widget value as a string

**repr\_value(value, name=None)**  
Return textual representation of the value.  
If `name` is not `None` raise an error.

**set\_value(value, name=None)**  
Set widget value.

If `name` is not `None` raise an error.

**class** `pylablib.core.gui.value_handling.LabelValueHandler(widget)`  
Bases: `pylablib.core.gui.value_handling.ISingleValueHandler`

Value handler for `QLabel` widget

**get\_single\_value()**  
Get the widget value

**set\_single\_value(value)**  
Set the widget value

**can\_set\_value(allow\_focus=True)**  
Check if setting value from the code is allowed.

**Parameters focus** – if `False`, indicates that settings of focused widgets isn't allowed, with some exceptions (buttons, check boxes, combo boxes)

**connect\_value\_changed\_handler** (*handler, only\_signal=True*)

Connect value changed signal.

If `only_signal==True`, equivalent to connecting a handler function to `get_value_changed_signal()` signal; however, if `only_signal==False`, it also works for some objects (e.g., `QLabel`) don't have built-in on-changed signals by calling the handler explicitly every time the value is changed.

Note that the connection is always direct (i.e., it doesn't deal with message queues and different threads, but rather directly calls the handler function). If you need to connect a handler to a signal using some other connection method, you can use `get_value_changed_signal()` directly.

**get\_handler** (*name=None*)

Get handler of a contained widget (or same widget, if `name==None`)

**get\_value** (*name=None*)

Get widget value.

If `name` is not `None` raise an error.

**get\_value\_changed\_signal** ()

Get the Qt signal emitted when the value is changed

**repr\_single\_value** (*value*)

Represent the widget value as a string

**repr\_value** (*value, name=None*)

Return textual representation of the value.

If `name` is not `None` raise an error.

**set\_value** (*value, name=None*)

Set widget value.

If `name` is not `None` raise an error.

**class** `pylablib.core.gui.value_handling.IBoolValueHandler` (*widget, labels=('Off', 'On')*)

Bases: `pylablib.core.gui.value_handling.ISingleValueHandler`

Generic value handler for widgets with boolean values

**repr\_single\_value** (*value*)

Represent the widget value as a string

**can\_set\_value** (*allow\_focus=True*)

Check if setting value from the code is allowed.

**Parameters** `focus` – if `False`, indicates that settings of focused widgets isn't allowed, with some exceptions (buttons, check boxes, combo boxes)

**connect\_value\_changed\_handler** (*handler, only\_signal=True*)

Connect value changed signal.

If `only_signal==True`, equivalent to connecting a handler function to `get_value_changed_signal()` signal; however, if `only_signal==False`, it also works for some objects (e.g., `QLabel`) don't have built-in on-changed signals by calling the handler explicitly every time the value is changed.

Note that the connection is always direct (i.e., it doesn't deal with message queues and different threads, but rather directly calls the handler function). If you need to connect a handler to a signal using some other connection method, you can use `get_value_changed_signal()` directly.

**get\_handler** (*name=None*)  
Get handler of a contained widget (or same widget, if *name==None*)

**get\_single\_value** ()  
Get the widget value

**get\_value** (*name=None*)  
Get widget value.  
  
If *name* is not *None* raise an error.

**get\_value\_changed\_signal** ()  
Get the Qt signal emitted when the value is changed

**repr\_value** (*value, name=None*)  
Return textual representation of the value.  
  
If *name* is not *None* raise an error.

**set\_single\_value** (*value*)  
Set the widget value

**set\_value** (*value, name=None*)  
Set widget value.  
  
If *name* is not *None* raise an error.

**class** `pylablib.core.gui.value_handling.CheckboxValueHandler` (*widget, labels=('Off', 'On')*)

Bases: `pylablib.core.gui.value_handling.IBoolValueHandler`

Value handler for `QCheckBox` widget

**get\_single\_value** ()  
Get the widget value

**set\_single\_value** (*value*)  
Set the widget value

**repr\_single\_value** (*value*)  
Represent the widget value as a string

**get\_value\_changed\_signal** ()  
Get the Qt signal emitted when the value is changed

**can\_set\_value** (*allow\_focus=True*)  
Check if setting value from the code is allowed.

**Parameters** *focus* – if *False*, indicates that settings of focused widgets isn't allowed, with some exceptions (buttons, check boxes, combo boxes)

**connect\_value\_changed\_handler** (*handler, only\_signal=True*)  
Connect value changed signal.

If *only\_signal==True*, equivalent to connecting a handler function to `get_value_changed_signal()` signal; however, if *only\_signal==False*, it also works for some objects (e.g., `QLabel`) don't have built-in on-changed signals by calling the handler explicitly every time the value is changed.

Note that the connection is always direct (i.e., it doesn't deal with message queues and different threads, but rather directly calls the handler function). If you need to connect a handler to a signal using some other connection method, you can use `get_value_changed_signal()` directly.

**get\_handler** (*name=None*)  
Get handler of a contained widget (or same widget, if *name==None*)

**get\_value** (*name=None*)

Get widget value.

If name is not None raise an error.

**repr\_value** (*value, name=None*)

Return textual representation of the value.

If name is not None raise an error.

**set\_value** (*value, name=None*)

Set widget value.

If name is not None raise an error.

**class** `pylablib.core.gui.value_handling.PushButtonValueHandler` (*widget, labels=('Off', 'On')*)

Bases: `pylablib.core.gui.value_handling.IBoolValueHandler`

Value handler for QPushButton widget

**get\_single\_value** ()

Get the widget value

**set\_single\_value** (*value*)

Set the widget value

**get\_value\_changed\_signal** ()

Get the Qt signal emitted when the value is changed

**repr\_single\_value** (*value*)

Represent the widget value as a string

**can\_set\_value** (*allow\_focus=True*)

Check if setting value from the code is allowed.

**Parameters** **focus** – if False, indicates that settings of focused widgets isn't allowed, with some exceptions (buttons, check boxes, combo boxes)

**connect\_value\_changed\_handler** (*handler, only\_signal=True*)

Connect value changed signal.

If `only_signal==True`, equivalent to connecting a handler function to `get_value_changed_signal()` signal; however, if `only_signal==False`, it also works for some objects (e.g., `QLabel`) don't have built-in on-changed signals by calling the handler explicitly every time the value is changed.

Note that the connection is always direct (i.e., it doesn't deal with message queues and different threads, but rather directly calls the handler function). If you need to connect a handler to a signal using some other connection method, you can use `get_value_changed_signal()` directly.

**get\_handler** (*name=None*)

Get handler of a contained widget (or same widget, if `name==None`)

**get\_value** (*name=None*)

Get widget value.

If name is not None raise an error.

**repr\_value** (*value, name=None*)

Return textual representation of the value.

If name is not None raise an error.

**set\_value** (*value*, *name=None*)

Set widget value.

If name is not None raise an error.

**class** `pylablib.core.gui.value_handling.ToolButtonValueHandler` (*widget*, *labels=('Off', 'On')*, *la-*

Bases: `pylablib.core.gui.value_handling.IBoolValueHandler`

Value handler for `QToolButton` widget

**get\_single\_value** ()

Get the widget value

**set\_single\_value** (*value*)

Set the widget value

**get\_value\_changed\_signal** ()

Get the Qt signal emitted when the value is changed

**repr\_single\_value** (*value*)

Represent the widget value as a string

**can\_set\_value** (*allow\_focus=True*)

Check if setting value from the code is allowed.

**Parameters** **focus** – if False, indicates that settings of focused widgets isn't allowed, with some exceptions (buttons, check boxes, combo boxes)

**connect\_value\_changed\_handler** (*handler*, *only\_signal=True*)

Connect value changed signal.

If `only_signal==True`, equivalent to connecting a handler function to `get_value_changed_signal()` signal; however, if `only_signal==False`, it also works for some objects (e.g., `QLabel`) don't have built-in on-changed signals by calling the handler explicitly every time the value is changed.

Note that the connection is always direct (i.e., it doesn't deal with message queues and different threads, but rather directly calls the handler function). If you need to connect a handler to a signal using some other connection method, you can use `get_value_changed_signal()` directly.

**get\_handler** (*name=None*)

Get handler of a contained widget (or same widget, if `name==None`)

**get\_value** (*name=None*)

Get widget value.

If name is not None raise an error.

**repr\_value** (*value*, *name=None*)

Return textual representation of the value.

If name is not None raise an error.

**set\_value** (*value*, *name=None*)

Set widget value.

If name is not None raise an error.

**class** `pylablib.core.gui.value_handling.ComboBoxValueHandler` (*widget*)

Bases: `pylablib.core.gui.value_handling.ISingleValueHandler`

Value handler for `QComboBox` widget

**get\_single\_value()**  
Get the widget value

**set\_single\_value(value)**  
Set the widget value

**get\_value\_changed\_signal()**  
Get the Qt signal emitted when the value is changed

**repr\_single\_value(value)**  
Represent the widget value as a string

**can\_set\_value(allow\_focus=True)**  
Check if setting value from the code is allowed.

**Parameters focus** – if `False`, indicates that settings of focused widgets isn't allowed, with some exceptions (buttons, check boxes, combo boxes)

**connect\_value\_changed\_handler(handler, only\_signal=True)**  
Connect value changed signal.

If `only_signal==True`, equivalent to connecting a handler function to `get_value_changed_signal()` signal; however, if `only_signal==False`, it also works for some objects (e.g., `QLabel`) don't have built-in on-changed signals by calling the handler explicitly every time the value is changed.

Note that the connection is always direct (i.e., it doesn't deal with message queues and different threads, but rather directly calls the handler function). If you need to connect a handler to a signal using some other connection method, you can use `get_value_changed_signal()` directly.

**get\_handler(name=None)**  
Get handler of a contained widget (or same widget, if `name==None`)

**get\_value(name=None)**  
Get widget value.

If `name` is not `None` raise an error.

**repr\_value(value, name=None)**  
Return textual representation of the value.

If `name` is not `None` raise an error.

**set\_value(value, name=None)**  
Set widget value.

If `name` is not `None` raise an error.

**class** `pylablib.core.gui.value_handling.ProgressBarValueHandler(widget)`  
Bases: `pylablib.core.gui.value_handling.ISingleValueHandler`

Value handler for `QProgressBar` widget

**get\_single\_value()**  
Get the widget value

**set\_single\_value(value)**  
Set the widget value

**can\_set\_value(allow\_focus=True)**  
Check if setting value from the code is allowed.

**Parameters focus** – if `False`, indicates that settings of focused widgets isn't allowed, with some exceptions (buttons, check boxes, combo boxes)

**connect\_value\_changed\_handler** (*handler*, *only\_signal=True*)

Connect value changed signal.

If `only_signal==True`, equivalent to connecting a handler function to `get_value_changed_signal()` signal; however, if `only_signal==False`, it also works for some objects (e.g., `QLabel`) don't have built-in on-changed signals by calling the handler explicitly every time the value is changed.

Note that the connection is always direct (i.e., it doesn't deal with message queues and different threads, but rather directly calls the handler function). If you need to connect a handler to a signal using some other connection method, you can use `get_value_changed_signal()` directly.

**get\_handler** (*name=None*)

Get handler of a contained widget (or same widget, if `name==None`)

**get\_value** (*name=None*)

Get widget value.

If `name` is not `None` raise an error.

**get\_value\_changed\_signal** ()

Get the Qt signal emitted when the value is changed

**repr\_single\_value** (*value*)

Represent the widget value as a string

**repr\_value** (*value*, *name=None*)

Return textual representation of the value.

If `name` is not `None` raise an error.

**set\_value** (*value*, *name=None*)

Set widget value.

If `name` is not `None` raise an error.

`pylablib.core.gui.value_handling.is_handled_widget` (*widget*)

Check if the widget can be handles by `StandardValueHandler`

`pylablib.core.gui.value_handling.create_value_handler` (*widget*)

Autodetect value handler for the given widget

**class** `pylablib.core.gui.value_handling.IIndicatorHandler`

Bases: `object`

Generic handler of an indicator.

Has methods to get and set the indicator value.

**get\_value** (*name=None*)

Get indicator value.

If `name` is not `None`, it specifies the name of the indicator parameter inside the widget (for complex widgets).

**set\_value** (*value*, *name=None*)

Set indicator value.

If `name` is not `None`, it specifies the name of the indicator parameter inside the widget (for complex widgets).

`pylablib.core.gui.value_handling.VirtualIndicatorHandler`

alias of `pylablib.core.gui.value_handling.VirtualValueHandler`



```
class pylablib.core.gui.value_handling.StandardIndicatorHandler (widget, de-  
fault_name=None)
```

Bases: `pylablib.core.gui.value_handling.IIndicatorHandler`

Default indicator handler, typically used for custom widgets.

To implement getting and setting values, looks for `get/set_indicator` and `get/set_all_indicators` methods for the widget and uses them accordingly.

#### Parameters

- **widget** – handled widget
- **default\_name** (*str*) – default name to be supplied to `get/set_indicator` methods if they require a name argument.

```
get_value (name=None)  
Get indicator value.
```

If name is not None, it specifies the name of the indicator parameter inside the widget (for complex widgets).

```
set_value (value, name=None)  
Set indicator value.
```

If name is not None, it specifies the name of the indicator parameter inside the widget (for complex widgets).

```
class pylablib.core.gui.value_handling.LabelIndicatorHandler (label, for-  
matter=None, repr_value_name=None)
```

Bases: `pylablib.core.gui.value_handling.IIndicatorHandler`

Indicator handler which uses a label to show the value.

Can takes optional widget or widget handler which converts values into strings using its `repr_value` method (by default, use the standard string conversion).

#### Parameters

- **label** – widget or value handler used to represent the value (takes string values)
- **formatter** – specifies a way to turn values into string representation; can be a widget handler or a widget (its `repr_func` method is used to represent its value), a function (it takes either a single value argument or two arguments `name` and `value` and returns string value), or None (use simple string conversion)
- **repr\_value\_name** (*str*) – default name to be supplied to `repr_value` if it requires a name argument and name is not supplied

```
get_value (name=None)  
Get indicator value.
```

If name is not None, it specifies the name of the indicator parameter inside the widget (for complex widgets).

```
repr_value (value, name=None)  
Represent a value with a given name
```

```
set_value (value, name=None)  
Set indicator value.
```

If name is not None, it specifies the name of the indicator parameter inside the widget (for complex widgets).

`pylablib.core.gui.value_handling.create_indicator_handler` (*widget*, *label=None*, *require\_setter=False*)

Autodetect indicator handler for the given widget and optional indicator label

**exception** `pylablib.core.gui.value_handling.MissingGUIHandlerError`

Bases: `KeyError`

Missing GUI handler

**args**

**with\_traceback()**

Exception.with\_traceback(tb) – set self.\_\_traceback\_\_ to tb and return self.

**class** `pylablib.core.gui.value_handling.GUIValues` (*gui\_thread\_safe=True*)

Bases: `object`

A collection of values which can be used to manipulate many value handlers at once and represent them as a hierarchical structure.

Has four container-like accessor: `.h` for getting/adding/removing the value handler (i.e., `self.get_handler(name)` is equivalent to `self.h[name]`, and `self.add_handler(name, handler)` is equivalent to `self.h[name]=handler`, and `self.remove_handler(name)` is equivalent to `del self.h[name]`), `.w` for getting the underlying widget (i.e., `self.get_widget(name)` is equivalent to `self.w[name]`), `.v` for settings/getting values (i.e., `self.get_value(name)` is equivalent to `self.v[name]`, and `self.set_value(name, value)` is equivalent to `self.v[name]=value`), `.i` for settings/getting indicator values (i.e., `self.get_indicator(name)` is equivalent to `self.i[name]`, and `self.set_indicator(name, value)` is equivalent to `self.i[name]=value`) `.vs` for getting the value changed Qt signal (i.e., `self.get_value_changed_signal(name)` is equivalent to `self.s[name]`),

**Parameters** `gui_thread_safe` (*bool*) – if True, all value-access and indicator-access calls (`get/set_value`, `get/set_all_values`, `get/set_indicator`, `get/set_all_indicators`, and `update_indicators`) are automatically called in the GUI thread.

**add\_handler** (*name*, *handler*)

Add a value handler under a given name

**remove\_handler** (*name*, *remove\_indicator=True*, *disconnect=False*)

Remove the value handler with a given name.

If `remove_indicator==True`, also try to remove the indicator widget. If `disconnect==True`, also disconnect all slots connected to the `value_changed` signal. Unlike most methods (e.g., `get_value()` or `get_handler()`), does not recursively query the children, so it only works if the handler is contained in this table.

**get\_handler** (*name*)

Get the value handler with the given name

**add\_widget** (*name*, *widget*, *add\_indicator=True*)

Add a widget under a given name (value handler type is auto-detected)

**get\_widget** (*name*)

Get the widget corresponding to the handler under the given name

**add\_nested** (*name*, *gui\_values*, *add\_indicator=True*)

Add a nested `GUIValues` under a given name

**add\_virtual\_element** (*name*, *value=None*, *multivalued=False*, *add\_indicator=True*)

Add a virtual value element.

Doesn't correspond to any actual widget, but behaves very similarly from the application point of view (its value can be set or read, it has on-change events, it can have indicator). The element value is simply stored on set and retrieved on get. If `multivalued==True`, the internal value is assumed to be complex, so it is forced to be a *Dictionary* every time it is set. If `add_indicator==True`, add default indicator handler as well.

**add\_property\_element** (*name*, *getter=None*, *setter=None*, *add\_indicator=True*)

Add a property value element.

Doesn't correspond to any actual widget, but behaves very similarly from the application point of view; each time the value is set or get, the corresponding setter and getter methods are called. If `add_indicator==True`, add default (stored value) indicator handler as well.

**add\_all\_children** (*root*, *root\_name=None*, *types\_include=None*, *types\_exclude=()*,  
*names\_exclude=None*)

Add a widget and all its children to the values set.

The result is organized as a tree using parent-child relations (note that it implies that only children widgets correspond to tree nodes, i.e., only their values can be get/set).

#### Parameters

- **root** – root widget
- **root\_name** – path to the sub-branch where the values will be placed
- **types\_include** – if not `None`, specifies list of widget classes (e.g., `QCheckBox`) to include
- **types\_exclude** – specifies list of widget classes to exclude
- **names\_exclude** – if not `None`, specifies list of widget names to exclude

**class IndicatorsSet** (*ind*)

Bases: `tuple`

**count** ()

Return number of occurrences of value.

**ind**

**index** ()

Return first index of value.

Raises `ValueError` if the value is not present.

**add\_indicator\_handler** (*name*, *handler*, *ind\_name='\_\_default\_\_'*)

Add indicator handler with a given name.

*ind\_name* can distinguish different sub-indicators with the same name, if the same value has multiple indicators.

**remove\_indicator\_handler** (*name*, *ind\_name=None*)

Remove indicator handler with a given name.

*ind\_name* can distinguish different sub-indicators with the same name, if the same value has multiple indicators. By default, remove all indicators with this name

**add\_widget\_indicator** (*name*, *widget*, *label=None*, *ind\_name='\_\_default\_\_'*)

Add widget-based indicator with a given name.

If *label* is `None`, use *widget*'s `get/set_indicator` or `get/set_all_indicators` functions to indicate the value. Otherwise, use the given *label* to indicate the value (*label* is used to show the value, *widget* is used to represent it). *ind\_name* can distinguish different sub-indicators with the same name, if the same value has multiple indicators.

**add\_label\_indicator** (*name, label, formatter=None, ind\_name='\_\_default\_\_'*)

Add label-based indicator with a given name.

*formatter* specifies a way to turn values into string representation; can be a widget handler or a widget (its `repr_func` method is used to represent its value), a function (it takes either a single value argument or two arguments *name* and *value* and returns string value), or `None` (use simple string conversion) *ind\_name* can distinguish different sub-indicators with the same name, if the same value has multiple indicators.

**get\_value** (*name=None*)

Get a value or a set of values in a subtree under a given name (all values by default).

Automatically handles complex widgets and sub-names. If *name* refers to a branch, return a *Dictionary* object containing tree structure of the names. If supplied, *include* and *exclude* are containers specifying included and excluded names (relative to the root); by default, include everything and exclude nothing.

**get\_all\_values** (*root=None*)

Get all values in the given sub-branch.

Same as `get_value()`, but returns an empty dictionary if the *name* is missing.

**set\_value** (*name, value*)

Set value under a given name.

Automatically handles complex widgets and sub-names

**set\_all\_values** (*value, root=None*)

**get\_indicator** (*name=None, ind\_name='\_\_default\_\_'*)

Get indicator value with a given name.

*ind\_name* can distinguish different sub-indicators with the same name, if the same value has multiple indicators. If supplied, *include* and *exclude* are containers specifying included and excluded names (relative to the root); by default, include everything and exclude nothing.

**get\_all\_indicators** (*root=None, ind\_name='\_\_default\_\_'*)

Get all indicator values in the given sub-branch.

Same as `get_indicator()`, but returns an empty dictionary if the *root* is missing.

**set\_indicator** (*name, value, ind\_name=None, ignore\_missing=False*)

Set indicator value with a given name.

*ind\_name* can distinguish different sub-indicators with the same name, if the same value has multiple indicators. By default, set all sub-indicators to the given value. If supplied, *include* and *exclude* are containers specifying included and excluded names (relative to the root); by default, include everything and exclude nothing. If *ignore\_missing==True* and the given indicator and sub-indicator names are missing, raise an error; otherwise, do nothing.

**set\_all\_indicators** (*value, root="", ind\_name=None, ignore\_missing=True*)

**update\_indicators** (*root=""*)

Update all indicators in a subtree with the given root (all values by default) to represent current values.

If supplied, *include* and *exclude* are containers specifying included and excluded names (relative to the root); by default, include everything and exclude nothing.

**repr\_value** (*name, value*)

Get a textual representation of a value under a given name.

Automatically handles complex widgets and sub-names.

**get\_value\_changed\_signal** (*name*)

Get changed events for a value under a given name

**update\_value** (*name=None*)

Send update signal for a handler with a given name or list of names.

Emit a value changed signal with the current value to notify the subscribed slots. If *name* is `None`, emit for all values in the table.

`pylablib.core.gui.value_handling.get_gui_values` (*gui\_values=None*,  
*gui\_values\_path=""*)

Get new or existing `GUIValues` object and the sub-branch path inside it based on the supplied arguments.

If *gui\_values* is `None` or `"new"`, create a new object and set empty root path. If *gui\_values* itself has `gui_values` attribute, get this attribute, and prepend object's `gui_values_path` attribute to the given path. Otherwise, assume that *gui\_values* is `GUIValues` object, and use the supplied root.

`pylablib.core.gui.value_handling.virtual_gui_values` (\*\**kwargs*)

Create a gui values set with all virtual values.

*kwargs* define element names and default values.

## Module contents

### pylablib.core.thread package

#### Submodules

#### pylablib.core.thread.callsync module

**class** `pylablib.core.thread.callsync.QCallResultSynchronizer` (*skippable=True*)

Bases: `pylablib.core.thread.synchronizing.QThreadNotifier`

**get\_progress** ()

Get the progress of the call execution.

Can be `"waiting"` (call is not done executing), `"done"` (call done successfully), `"fail"` (call failed, probably due to thread being stopped), `"skip"` (call was skipped), or `"exception"` (call raised an exception).

**is\_call\_done** ()

Check if the call is done

**skipped** ()

Check if the call was skipped

**failed** ()

Check if the call failed

**get\_value\_sync** (*timeout=None*, *default=None*, *error\_on\_fail=True*, *error\_on\_skip=True*,  
*pass\_exception=True*)

Wait (with the given *timeout*) for the value passed by the notifier

If `error_on_fail==True` and the controlled thread notifies of a fail (usually, if it's stopped before it executed the call), raise `threadprop.NoControllerThreadError`; otherwise, return *default*. If `error_on_skip==True` and the call was skipped (e.g., due to full call queue), raise `threadprop.SkippedCallError`; otherwise, return *default*. If `pass_exception==True` and the returned value represents exception, re-raise it in the caller thread; otherwise, return *default*.

**done\_notify()**

Check if notifying is done

**done\_wait()**

Check if waiting is done

**get\_value()**

Get the value passed by the notifier (doesn't check if it has been passed already)

**notify(\*args, \*\*kwargs)**

Notify the waiting process.

Can only be called once per notifier lifetime. If the notifier allows skipping, and this method is called before `wait()`, return immediately.

**notifying\_state()**

**success\_wait()**

Check if waiting is done successfully

**wait(\*args, \*\*kwargs)**

Wait for the notification.

Can only be called once per notifier lifetime. If the notifier allows skipping, and this method is called after `notify()`, return immediately.

**waiting()**

Check if waiting is in progress

**waiting\_state()**

**class** `pylablib.core.thread.callsync.QDummyResultSynchronizer`

Bases: `object`

Dummy result synchronizer for call which don't require result synchronization (e.g., multicasts)

**notify(value)**

**class** `pylablib.core.thread.callsync.QDirectResultSynchronizer(value)`

Bases: `object`

Result "synchronizer" for direct calls.

Behaves as a regular result synchronizer with an already executed call.

**get\_progress()**

Get the progress of the call execution (always return "done")

**is\_call\_done()**

Check if the call is done (always return True)

**skipped()**

Check if the call was skipped (always return False)

**failed()**

Check if the call failed (always return False)

**get\_value()**

Return stored value

**get\_value\_sync**(*timeout=None, default=None, error\_on\_fail=True, error\_on\_skip=True, pass\_exception=True*)

Return stored value.

Parameters are only for compatibility with `QCallResultSynchronizer`.

```

wait (*args, **kwargs)
    Do nothing (present only for compatibility with QCallResultSynchronizer)

notify (*args, **kwargs)
    Do nothing (present only for compatibility with QCallResultSynchronizer)

waiting()
    Check if waiting is in progress (always return False)

done_wait()
    Check if waiting is done (always return True)

success_wait()
    Check if waiting is done successfully (always return True)

done_notify()
    Check if notifying is done (always return True)

waiting_state()

notifying_state()

class pylablib.core.thread.callsync.QScheduledCall (func,                args=None,
                                                    kwargs=None,        silent=False,
                                                    result_synchronizer=None)

Bases: object

Object representing a scheduled remote call.

Can be executed, skipped, or failed in the target thread, in which case it notifies the result synchronizer (if
supplied).

Parameters

    • func – callable to be invoked in the destination thread

    • args – arguments to be passed to func

    • kwargs – keyword arguments to be passed to func

    • silent – if True, silence the exception in the execution thread and simply pass it to
      the caller thread; otherwise, the exception is raised in both threads

    • result_synchronizer – result synchronizer object; can be None (create
      new QCallResultSynchronizer), "async" (no result synchronization), or a
      QCallResultSynchronizer object.

class Callback (func, pass_result, call_on_exception, call_on_unschedule)
    Bases: tuple

    call_on_exception

    call_on_unschedule

    count()
        Return number of occurrences of value.

    func

    index()
        Return first index of value.

        Raises ValueError if the value is not present.

    pass_result

```

**execute** (*silent=None*)

Execute the call and notify the result synchronizer (invoked by the destination thread)

**add\_callback** (*callback*, *pass\_result=True*, *call\_on\_exception=False*, *call\_on\_unschedule=False*,  
*front=False*)

Set the callback to be executed after the main call is done.

If *pass\_result==True*, pass function result to the callback (or *None* if call failed); otherwise, pass no arguments. If *call\_on\_exception==True*, call it even if the original call raised an exception. If *call\_on\_unschedule==True*, call it for any call unscheduling event, including using *skip()* or *fail()* methods (this effectively ignores *call\_on\_exception*, since the callback is called regardless of the exception). If *front==True*, add the callback in the front of the line (executes first).

**fail** ()

Notify that the call is failed (invoked by the destination thread)

**skip** ()

Notify that the call is skipped (invoked by the destination thread)

**class** `pylablib.core.thread.callsync.TDefaultCallInfo` (*call\_time*)

Bases: `tuple`

**call\_time**

**count** ()

Return number of occurrences of value.

**index** ()

Return first index of value.

Raises `ValueError` if the value is not present.

**class** `pylablib.core.thread.callsync.QScheduler` (*call\_info\_argname=None*)

Bases: `object`

Generic call scheduler.

Two methods are used by the external scheduling routines: *build\_call()* to create a *QScheduledCall* with appropriate parameters, and *schedule()*, which takes a call and schedules it. The *schedule()* method should return *True* if the scheduling was successful (at least, for now), and *False* otherwise.

**Parameters** *call\_info\_argname* – if not *None*, supplies a name of a keyword argument via which call info (generated by *build\_call\_info()*) is passed on function call

**build\_call\_info** ()

Build call info tuple which can be passed to scheduled calls

**build\_call** (*func*, *args=None*, *kwargs=None*, *callback=None*, *pass\_result=True*, *call-  
back\_on\_exception=True*, *sync\_result=True*)

Build *QScheduledCall* for subsequent scheduling.

#### Parameters

- **func** – function to be called
- **args** – arguments to be passed to *func*
- **kwargs** – keyword arguments to be passed to *func*
- **callback** – optional callback to be called when *func* is done
- **pass\_result** (*bool*) – if *True*, pass *func* result as a single argument to the call-back; otherwise, give no arguments



- **callback\_on\_exception** (*bool*) – if True, execute the callback on call fail or skip (if it requires an argument, None is supplied); otherwise, only execute it if the call was successful
- **sync\_result** – if True, the call has a default result synchronizer; otherwise, no synchronization is made.

**schedule** (*call*)

Schedule the call

**clear** ()

Clear the scheduler

**class** `pylablib.core.thread.callsync.QDirectCallScheduler` (*call\_info\_argname=None*)

Bases: `pylablib.core.thread.callsync.QScheduler`

Simplest call scheduler: directly executes the calls on scheduling in the scheduling thread.

**Parameters** **call\_info\_argname** – if not None, supplies a name of a keyword argument via which call info (generated by `QScheduler.build_call_info()`) is passed on function call

**build\_call** (*func*, *args=None*, *kwargs=None*, *callback=None*, *pass\_result=True*, *callback\_on\_exception=True*, *sync\_result=False*)

Build `QScheduledCall` for subsequent scheduling.

**Parameters**

- **func** – function to be called
- **args** – arguments to be passed to *func*
- **kwargs** – keyword arguments to be passed to *func*
- **callback** – optional callback to be called when *func* is done
- **pass\_result** (*bool*) – if True, pass *func* result as a single argument to the callback; otherwise, give no arguments
- **callback\_on\_exception** (*bool*) – if True, execute the callback on call fail or skip (if it requires an argument, None is supplied); otherwise, only execute it if the call was successful
- **sync\_result** – if True, the call has a default result synchronizer; otherwise, no synchronization is made.

**schedule** (*call*)

Schedule the call

**build\_call\_info** ()

Build call info tuple which can be passed to scheduled calls

**clear** ()

Clear the scheduler

**class** `pylablib.core.thread.callsync.QQueueScheduler` (*on\_full\_queue='skip\_current'*, *call\_info\_argname=None*)

Bases: `pylablib.core.thread.callsync.QScheduler`

Call scheduler with a builtin call queue.

Supports placing the calls and retrieving them (from the destination thread). Has ability to skip some calls if, e.g., the queue is too full. Whether the call should be skipped is determined by `can_schedule()` (should be overloaded in subclasses). Used as a default command scheduler.

### Parameters

- **on\_full\_queue** – action to be taken if the call can't be scheduled (i.e., `can_schedule()` returns `False`); can be "skip\_current" (skip the call which is being scheduled), "skip\_newest" (skip the most recent call; place the current), "skip\_oldest" (skip the oldest call in the queue; place the current), "call\_current" (execute the call which is being scheduled immediately in the caller thread), "call\_newest" (execute the most recent call immediately in the caller thread), "call\_oldest" (execute the oldest call in the queue immediately in the caller thread), or "wait" (wait until the call can be scheduled, which is checked after every call removal from the queue; place the call)
- **call\_info\_argname** – if not `None`, supplies a name of a keyword argument via which call info (generated by `QScheduler.build_call_info()`) is passed on function call

### Methods to overload:

- `can_schedule()`: check if the call can be scheduled
- `call_added()`: called when a new call has been added to the queue
- `call_popped()`: called when a call has been removed from the queue (either for execution, or for skipping)

**can\_schedule** (*call*)

Check if the call can be scheduled

**call\_added** (*call*)

Called whenever *call* has been added to the queue

**call\_popped** (*call*, *idx*)

Called whenever *call* has been removed from the queue

*idx* determines the call position within the queue.

**schedule** (*call*)

Schedule a call

**pop\_call** ()

Pop the call from the queue head.

If the queue is empty, return `None`

**unschedule** (*call*)

Unschedule a given call.

Designed for joint queue operation, so the call is not notified (assume that it has been already notified elsewhere).

**has\_calls** ()

Check if there are queued calls

**clear** (*close=True*)

Clear the call queue.

If *close*==`True`, mark the queue as closed (any attempt to schedule more calls fails automatically) and fail all calls in the queue; otherwise, skip all calls currently in the queue.

**build\_call** (*func*, *args=None*, *kwargs=None*, *callback=None*, *pass\_result=True*, *call-back\_on\_exception=True*, *sync\_result=True*)

Build `QScheduledCall` for subsequent scheduling.

### Parameters

- **func** – function to be called
- **args** – arguments to be passed to *func*
- **kwargs** – keyword arguments to be passed to *func*
- **callback** – optional callback to be called when *func* is done
- **pass\_result** (*bool*) – if *True*, pass *func* result as a single argument to the callback; otherwise, give no arguments
- **callback\_on\_exception** (*bool*) – if *True*, execute the callback on call fail or skip (if it requires an argument, *None* is supplied); otherwise, only execute it if the call was successful
- **sync\_result** – if *True*, the call has a default result synchronizer; otherwise, no synchronization is made.

**build\_call\_info()**

Build call info tuple which can be passed to scheduled calls

```
class pylablib.core.thread.callsync.QQueueLengthLimitScheduler (max_len=1,
                                                                on_full_queue='skip_current',
                                                                call_info_argname=None)
```

Bases: `pylablib.core.thread.callsync.QQueueScheduler`

Queued call scheduler with a length limit.

### Parameters

- **max\_len** – maximal queue length; non-positive values are interpreted as no limit can also be a tuple (*arg\_name*, *max\_len*), in which case the length is calculated separately for every value of the parameter *arg\_name* supplied to the method
- **on\_full\_queue** – action to be taken if the call can't be scheduled (the queue is full); can be "skip\_current" (skip the call which is being scheduled), "skip\_newest" (skip the most recent call; place the current), "skip\_oldest" (skip the oldest call in the queue; place the current), "call\_current" (execute the call which is being scheduled immediately in the caller thread), "call\_newest" (execute the most recent call immediately in the caller thread), "call\_oldest" (execute the oldest call in the queue immediately in the caller thread), or "wait" (wait until the call can be scheduled, which is checked after every call removal from the queue; place the call)
- **call\_info\_argname** – if not *None*, supplies a name of a keyword argument via which call info (generated by `QScheduler.build_call_info()`) is passed on function call

**change\_max\_len(max\_len)**

Change maximal length of the call queue (doesn't affect already scheduled calls)

**get\_current\_len()**

Get current number of calls in the queue

**call\_added(call)**

Called whenever *call* has been added to the queue

**call\_popped(call, idx)**

Called whenever *call* has been removed from the queue

*idx* determines the call position within the queue.

**can\_schedule** (*call*)

Check if the call can be scheduled

**build\_call** (*func*, *args=None*, *kwargs=None*, *callback=None*, *pass\_result=True*, *call-back\_on\_exception=True*, *sync\_result=True*)

Build *QScheduledCall* for subsequent scheduling.

#### Parameters

- **func** – function to be called
- **args** – arguments to be passed to *func*
- **kwargs** – keyword arguments to be passed to *func*
- **callback** – optional callback to be called when *func* is done
- **pass\_result** (*bool*) – if *True*, pass *func* result as a single argument to the call-back; otherwise, give no arguments
- **callback\_on\_exception** (*bool*) – if *True*, execute the callback on call fail or skip (if it requires an argument, *None* is supplied); otherwise, only execute it if the call was successful
- **sync\_result** – if *True*, the call has a default result synchronizer; otherwise, no synchronization is made.

**build\_call\_info** ()

Build call info tuple which can be passed to scheduled calls

**clear** (*close=True*)

Clear the call queue.

If *close==True*, mark the queue as closed (any attempt to schedule more calls fails automatically) and fail all calls in the queue; otherwise, skip all calls currently in the queue.

**has\_calls** ()

Check if there are queued calls

**pop\_call** ()

Pop the call from the queue head.

If the queue is empty, return *None*

**schedule** (*call*)

Schedule a call

**unschedule** (*call*)

Unschedule a given call.

Designed for joint queue operation, so the call is not notified (assume that it has been already notified elsewhere).

```
class pylablib.core.thread.callsync.QQueueSizeLimitScheduler (max_size=1,  
                                                             size_calc=None,  
                                                             on_full_queue='skip_current',  
                                                             call_info_argname=None)
```

Bases: *pylablib.core.thread.callsync.QQueueScheduler*

Queued call scheduler with a generic size limit; similar to *QQueueLengthLimitScheduler*, but more flexible and can implement more restrictions (e.g., queue length and arguments RAM size).

#### Parameters

- **max\_size** – maximal total size of the arguments; can be either a single number, or a tuple (if several different size metrics are involved); non-positive values are interpreted as no limit
- **size\_calc** – function that takes a single argument (call to be placed) and returns its size; can be either a single number, or a tuple (if several different size metrics are involved); by default, simply returns 1, which makes the scheduler behavior identical to [QQueueLengthLimitScheduler](#)
- **on\_full\_queue** – action to be taken if the call can't be scheduled (the queue is full); can be "skip\_current" (skip the call which is being scheduled), "skip\_newest" (skip the most recent call; place the current), "skip\_oldest" (skip the oldest call in the queue; place the current), "call\_current" (execute the call which is being scheduled immediately in the caller thread), "call\_newest" (execute the most recent call immediately in the caller thread), "call\_oldest" (execute the oldest call in the queue immediately in the caller thread), or "wait" (wait until the call can be scheduled, which is checked after every call removal from the queue; place the call)
- **call\_info\_argname** – if not None, supplies a name of a keyword argument via which call info (generated by [QScheduler.build\\_call\\_info\(\)](#)) is passed on function call

**change\_max\_size** (*max\_size*)

Change size restrictions

**get\_current\_size** ()

Get current size metrics

**call\_added** (*call*)

Called whenever *call* has been added to the queue

**call\_popped** (*call*, *idx*)

Called whenever *call* has been removed from the queue

*idx* determines the call position within the queue.

**can\_schedule** (*call*)

Check if the call can be scheduled

**build\_call** (*func*, *args=None*, *kwargs=None*, *callback=None*, *pass\_result=True*, *call-back\_on\_exception=True*, *sync\_result=True*)

Build [QScheduledCall](#) for subsequent scheduling.

#### Parameters

- **func** – function to be called
- **args** – arguments to be passed to *func*
- **kwargs** – keyword arguments to be passed to *func*
- **callback** – optional callback to be called when *func* is done
- **pass\_result** (*bool*) – if True, pass *func* result as a single argument to the callback; otherwise, give no arguments
- **callback\_on\_exception** (*bool*) – if True, execute the callback on call fail or skip (if it requires an argument, None is supplied); otherwise, only execute it if the call was successful
- **sync\_result** – if True, the call has a default result synchronizer; otherwise, no synchronization is made.

**build\_call\_info()**

Build call info tuple which can be passed to scheduled calls

**clear** (*close=True*)

Clear the call queue.

If *close==True*, mark the queue as closed (any attempt to schedule more calls fails automatically) and fail all calls in the queue; otherwise, skip all calls currently in the queue.

**has\_calls()**

Check if there are queued calls

**pop\_call()**

Pop the call from the queue head.

If the queue is empty, return *None*

**schedule** (*call*)

Schedule a call

**unschedule** (*call*)

Unschedule a given call.

Designed for joint queue operation, so the call is not notified (assume that it has been already notified elsewhere).

`pylablib.core.thread.callsync.schedule_multiple_queues(call, queues)`

Schedule the call simultaneously in several queues.

Go through queues in the given order and schedule call in every one of them. If one of the schedules failed or the call has been executed there, unschedule it from all the previous queues and return *False*; otherwise, return *True*.

**class** `pylablib.core.thread.callsync.QMultiQueueScheduler` (*schedulers*, *notifiers*)

Bases: `object`

Wrapper around `schedule_multiple_queues()` which acts as a single scheduler.

Support additional notifiers, which are called if the scheduling is successful (e.g., to notify and wake up the destination thread).

**build\_call** (*\*args*, *\*\*kwargs*)

**schedule** (*call*)

**class** `pylablib.core.thread.callsync.QThreadCallScheduler` (*thread=None*,  
*tag=None*, *priority=0*, *interrupt=True*,  
*call\_info\_argname=None*)

Bases: `pylablib.core.thread.callsync.QScheduler`

Call scheduler via thread calls (`QThreadController.call_in_thread_callback()`)

#### Parameters

- **thread** – destination thread (by default, thread which creates the scheduler)
- **tag** – if supplied, send the call in a message with the given tag; otherwise, use the interrupt call (generally, higher priority method).
- **priority** – message priority (only when *tag* is not *None*)
- **interrupt** – whether the call is an interrupt (call inside any loop, e.g., during waiting or sleeping), or it should be called in the main event loop

- **call\_info\_argname** – if not `None`, supplies a name of a keyword argument via which call info (generated by `QScheduler.build_call_info()`) is passed on function call

**schedule** (*call*)

Schedule the call

**build\_call** (*func*, *args=None*, *kwargs=None*, *callback=None*, *pass\_result=True*, *callback\_on\_exception=True*, *sync\_result=True*)

Build `QScheduledCall` for subsequent scheduling.

#### Parameters

- **func** – function to be called
- **args** – arguments to be passed to *func*
- **kwargs** – keyword arguments to be passed to *func*
- **callback** – optional callback to be called when *func* is done
- **pass\_result** (*bool*) – if `True`, pass *func* result as a single argument to the callback; otherwise, give no arguments
- **callback\_on\_exception** (*bool*) – if `True`, execute the callback on call fail or skip (if it requires an argument, `None` is supplied); otherwise, only execute it if the call was successful
- **sync\_result** – if `True`, the call has a default result synchronizer; otherwise, no synchronization is made.

**build\_call\_info** ()

Build call info tuple which can be passed to scheduled calls

**clear** ()

Clear the scheduler

```
class pylablib.core.thread.callsync.QMulticastThreadCallScheduler (thread=None,
                                                                    limit_queue=1,
                                                                    tag=None,
                                                                    priority=0,
                                                                    interrupt=True,
                                                                    call_info_argname=None)
```

Bases: `pylablib.core.thread.callsync.QThreadCallScheduler`

Extended call scheduler via thread calls, which can limit number of queued calls.

#### Parameters

- **thread** – destination thread (by default, thread which creates the scheduler)
- **limit\_queue** – call queue limit (non-positive numbers are interpreted as no limit)
- **tag** – if supplied, send the call in a message with the given tag; otherwise, use the interrupt call (generally, higher priority method).
- **priority** – message priority (only when *tag* is not `None`)
- **interrupt** – whether the call is an interrupt (call inside any loop, e.g., during waiting or sleeping), or it should be called in the main event loop
- **call\_info\_argname** – if not `None`, supplies a name of a keyword argument via which call info (generated by `QScheduler.build_call_info()`) is passed on function call

**schedule** (*call*)

Schedule the call

**build\_call** (*func*, *args=None*, *kwargs=None*, *callback=None*, *pass\_result=True*, *callback\_on\_exception=True*, *sync\_result=True*)

Build *QScheduledCall* for subsequent scheduling.

**Parameters**

- **func** – function to be called
- **args** – arguments to be passed to *func*
- **kwargs** – keyword arguments to be passed to *func*
- **callback** – optional callback to be called when *func* is done
- **pass\_result** (*bool*) – if *True*, pass *func* result as a single argument to the callback; otherwise, give no arguments
- **callback\_on\_exception** (*bool*) – if *True*, execute the callback on call fail or skip (if it requires an argument, *None* is supplied); otherwise, only execute it if the call was successful
- **sync\_result** – if *True*, the call has a default result synchronizer; otherwise, no synchronization is made.

**build\_call\_info** ()

Build call info tuple which can be passed to scheduled calls

**clear** ()

Clear the scheduler

## pylablib.core.thread.controller module

`pylablib.core.thread.controller.exint (error_msg_template='{ }:',  
pass_stop_exception=False)`

Context that intercepts exceptions and stops the execution in a controlled manner (quitting the main thread)

`pylablib.core.thread.controller.add_exception_hook (name, func, single_call=False)`

Add an exception hook, which is called whenever exception is caught via *exint* () wrapper.

If *single\_call==True*, the hook is removed from the set when it is called.

`pylablib.core.thread.controller.remove_exception_hook (name)`

Remove the exception hook with the given name

`pylablib.core.thread.controller.exsafe (func)`

Decorator that intercepts exceptions raised by *func* and stops the execution in a controlled manner (quitting the main thread)

`pylablib.core.thread.controller.exsafeSlot (*slargs, **slkwargs)`

Wrapper around Qt slot which intercepts exceptions and stops the execution in a controlled manner

`pylablib.core.thread.controller.toploopSlot (*slargs, **slkwargs)`

Wrapper around Qt slot which intercepts exceptions and stops the execution in a controlled manner

**class** `pylablib.core.thread.controller.QThreadControllerThread (controller)`

Bases: `object`

`finalized = <Mock name='mock.QtCore.pyqtSignal()' id='140096333774352'>`

`run ()`



**quit\_sync()**

`pylablib.core.thread.controller.remote_call(func)`

Decorator that turns a controller method into a remote call (call from a different thread is passed synchronously)

`pylablib.core.thread.controller.call_in_thread(thread_name, interrupt=True,  
pass_exception=True, silent=False,  
sync=True)`

Decorator that turns any function into a remote call in a thread with a given name (call from a different thread is passed synchronously)

`pylablib.core.thread.controller.call_in_gui_thread(func=None, interrupt=True,  
pass_exception=True,  
silent=False, sync=True)`

Decorator that turns any function into a remote call in a GUI thread (call from a different thread is passed synchronously)

`pylablib.core.thread.controller.gui_thread_method(func)`

Decorator for an object's method that checks if the object's `gui_thread_safe` attribute is true, in which case the call is routed to the GUI thread

**class** `pylablib.core.thread.controller.QThreadController(name=None, kind='loop',  
multicast_pool=None)`

Bases: `object`

Generic Qt thread controller.

Responsible for all inter-thread synchronization. There is one controller per thread, and

#### Parameters

- **name** (*str*) – thread name (by default, generate a new unique name); this name can be used to obtain thread controller via `get_controller()`
- **kind** (*str*) – thread kind; can be "loop" (thread is running in the Qt message loop; behavior is implemented in `process_message()` and remote calls), "run" (thread executes `run()` method and quits after it is complete), or "main" (can only be created in the main GUI thread)
- **multicast\_pool** – `MulticastPool` for this thread (by default, use the default common pool)

#### Methods to overload:

- `on_start()`: executed on the thread startup (between synchronization points "start" and "run")
- `on_finish()`: executed on thread cleanup (attempts to execute in any case, including exceptions)
- `run()`: executed once per thread; thread is stopped afterwards (only if `kind=="run"`)
- `process_message()`: function that takes 2 arguments (tag and value) of the message and processes it; return (in which case it is stored and can be recovered via `wait_for_message()/pop_message()`); by default, always return False
- `process_interrupt()`: function that takes 2 arguments (tag and value) of the interrupt message (message by default, assumes that any value with tag "execute" is a function and executes it)

#### Signals:

- started: emitted on thread start (after `on_start()` is executed)
- finished: emitted on thread finish (before `on_finish()` is executed)

**started** = `<Mock name='mock.QtCore.pyqtSignal()' id='140096333774352'>`

This signal is emitted after the thread has started (after the setup code has been executed, before its lifetime state is changed)

**finished** = `<Mock name='mock.QtCore.pyqtSignal()' id='140096333774352'>`

This signal is emitted before the thread has finished (before the cleanup code has been executed, after its lifetime state is changed)

**allowing\_toploop** (*depth=1*)

Context manager which temporarily treats the current loop level and several deeper levels as a top loop.

All event loops which lie up to *depth* below this one are treated as top loops.

**blocking\_control\_signals** (*kinds='all', ignore=None*)

Context manager which temporarily blocks external control signals.

After leaving the wrapped code segment, all of the blocked but not ignored calls are executed. *kind* determines the kind of calls to block; it is a collection of elements among "message", "stop", and "call" and blocks, correspondingly, messages, stop signals, and any `call_in_thread`-related requests; can be also be "all", which includes all of these categories. *ignore* specifies kinds which are completely ignored if sent during the blocking interval; can also be "all", which includes all of the *kinds* categories. Useful to temporarily “suspend” the thread communication with other threads, especially for the main GUI thread (e.g., to show a blocking message box). Local call method.

**wait\_for\_message** (*tag, timeout=None, top\_loop=False*)

Wait for a single message with a given tag.

Return value of a received message with this tag. If timeout is passed, raise `threadprop.TimeoutThreadError`. If `top_loop==True`, treat the waiting as the top message loop (i.e., any top loop message or signal can be executed here). Local call method.

**new\_messages\_number** (*tag*)

Get the number of queued messages with a given tag.

Local call method.

**pop\_message** (*tag*)

Pop the latest message with the given tag.

Select the message with the highest priority, and among those the oldest one. If no messages are available, raise `threadprop.NoMessageThreadError`. Local call method.

**wait\_for\_sync** (*tag, uid, timeout=None*)

Wait for synchronization signal with the given tag and UID.

This method is rarely invoked directly, and is usually used by synchronizers code. If timeout is passed, raise `threadprop.TimeoutThreadError`. Local call method.

**wait\_for\_any\_message** (*timeout=None, top\_loop=False*)

Wait for any message (including synchronization messages or pokes).

If timeout is passed, raise `threadprop.TimeoutThreadError`. If `top_loop==True`, treat the waiting as the top message loop (i.e., any top loop message or signal can be executed here). Local call method.

**wait\_until** (*check, timeout=None, top\_loop=False*)

Wait until a given condition is true.

Condition is given by the *check* function, which is called after every new received message and should return True if the condition is met. If `top_loop==True`, treat the waiting as the top message loop (i.e., any top loop message or signal can be executed here). If timeout is passed, raise `threadprop.TimeoutThreadError`. Local call method.

**check\_messages** (*top\_loop=False*)

Receive new messages.

Runs the underlying message loop to process newly received message and signals (and place them in corresponding queues if necessary). This method is rarely invoked, and only should be used periodically during long computations to not ‘freeze’ the thread. If `top_loop==True`, treat the waiting as the top message loop (i.e., any top loop message or signal can be executed here). Local call method.

**sleep** (*timeout, wake\_on\_message=False, top\_loop=False*)

Sleep for a given time (in seconds).

Unlike `time.sleep()`, constantly checks the event loop for new messages (e.g., if stop or interrupt commands are issued). In addition, if `wake_on_message==True`, wake up if any message has been received; in this case, return `True` if the wait has been completed, and `False` if it has been interrupted by a message. If `top_loop==True`, treat the waiting as the top message loop (i.e., any top loop message or signal can be executed here). If `timeout` is `None`, wait forever (usually, until the application is closed, or some interrupt message raises an error). Local call method.

**no\_stopping** ()

Context manager, which temporarily suspends stop requests (`InterruptExceptionStop` exceptions).

If the stop request has been made within this block, raise the exception on exit. Note that `stop()` method and, correspondingly, `stop_controller()` still work, when called from the controlled thread.

**process\_interrupt** (*tag, value*)

Process a new interrupt.

If the function returns `False`, the interrupt is put in the corresponding queue. Otherwise, the message is interrupt to be already, and it gets ‘absorbed’. Local call method, called automatically.

**process\_message** (*tag, value*)

Process a new message.

If the function returns `False`, the message is put in the corresponding queue. Otherwise, the message is considered to be already, and it gets ‘absorbed’. Local call method, called automatically.

**on\_start** ()

Method invoked on the start of the thread.

Local call method, called automatically.

**on\_finish** ()

Method invoked in the end of the thread.

Called regardless of the stopping reason (normal finishing, exception, application finishing). Local call method, called automatically.

**run** ()

Method called to run the main thread code (only for "run" thread kind).

Local call method, called automatically.

**subscribe\_sync** (*callback, srcs='any', tags=None, dsts='any', filt=None, subscription\_priority=0, limit\_queue=None, call\_interrupt=True, add\_call\_info=False, return\_result=False, sid=None*)

Subscribe a synchronous callback to a multicast.

If a multicast is sent, `callback` is called from the `dest_controller` thread (by default, thread which is calling this function) via the thread call mechanism (`QThreadController.call_in_thread_callback()`). In Qt, analogous to making a signal connection with a queued call. By default, the subscribed destination is the thread’s name. Local call method.

### Parameters

- **callback** – callback function, which takes 3 arguments: source, tag, and value.
- **srcs** (*str* or [*str*]) – multicast source name or list of source names to filter the subscription; can be "any" (any source) or "all" (only multicasts specifically having "all" as a source).
- **tags** – multicast tag or list of tags to filter the subscription (any tag by default); can also contain Unix shell style pattern ("\*" matches everything, "?" matches one symbol, etc.)
- **dsts** (*str* or [*str*]) – multicast destination name or list of destination names to filter the subscription; can be "any" (any destination) or "all" (only source specifically having "all" as a destination).
- **filt** (*callable*) – additional filter function which takes 4 arguments: source, destination, tag, and value, and checks whether multicast passes the requirements.
- **subscription\_priority** (*int*) – subscription priority (higher priority subscribers are called first).
- **limit\_queue** (*int*) – limits the maximal number of scheduled calls (if the multicast is sent while at least *limit\_queue* callbacks are already in queue to be executed, ignore it) 0 or negative value means no limit (not recommended, as it can increase the queue indefinitely if the multicast rate is high enough)
- **call\_interrupt** – whether the call is an interrupt (call inside any loop, e.g., during waiting or sleeping), or it should be called in the main event loop
- **add\_call\_info** (*bool*) – if True, add a fourth argument containing a call information (tuple with a single element, a timestamps of the call).
- **return\_result** – if True, use a result synchronizer to return the result of the subscribed call; otherwise, ignore the result
- **sid** (*int*) – subscription ID (by default, generate a new unique name).

**subscribe\_direct** (*callback*, *srcs*='any', *tags*=None, *dsts*='any', *filt*=None, *subscription\_priority*=0, *scheduler*=None, *return\_result*=False, *sid*=None)

Subscribe asynchronous callback to a multicast.

If a multicast is sent, *callback* is called from the sending thread (not subscribed thread). Therefore, should be used with care. In Qt, analogous to making a signal connection with a direct call. By default, the subscribed destination is the thread's name. Local call method.

### Parameters

- **callback** – callback function, which takes 3 arguments: source, tag, and value.
- **srcs** (*str* or [*str*]) – multicast source name or list of source names to filter the subscription; can be "any" (any source) or "all" (only multicasts specifically having "all" as a source).
- **tags** – multicast tag or list of tags to filter the subscription (any tag by default); can also contain Unix shell style pattern ("\*" matches everything, "?" matches one symbol, etc.)
- **dsts** (*str* or [*str*]) – multicast destination name or list of destination names to filter the subscription; can be "any" (any destination) or "all" (only source specifically having "all" as a destination).
- **filt** (*callable*) – additional filter function which takes 4 arguments: source, destination, tag, and value, and checks whether multicast passes the requirements.

- **subscription\_priority** (*int*) – subscription priority (higher priority subscribers are called first).
- **scheduler** – if defined, multicast call gets scheduled using this scheduler instead of being called directly (which is the default behavior)
- **return\_result** – if `True`, use a result synchronizer to return the result of the subscribed call; otherwise, ignore the result
- **sid** (*int*) – subscription ID (by default, generate a new unique id and return it).

#### **unsubscribe** (*sid*)

Unsubscribe from a subscription with a given ID.

Note that multicasts which are already emitted but not processed will remain in the queue; if they need to be ignored, it should be handled explicitly. Local call method.

#### **send\_multicast** (*dst='any', tag=None, value=None, src=None, filter\_results=True*)

Send a multicast to the multicast pool.

By default, the multicast source is the thread's name. Return result synchronizers for all executed subscribed methods. Local call method.

##### **Parameters**

- **dst** (*str*) – multicast destination; can be a name, "all" (will pass all subscribers' destination filters), or "any" (will only be passed to subscribers specifically subscribed to multicast with "any" destination).
- **tag** (*str*) – multicast tag.
- **value** – multicast value.
- **src** (*str*) – multicast source; can be `None` (current thread name), a specific name, "all" (will pass all subscribers' source filters), or "any" (will only be passed to subscribers specifically subscribed to multicast with "any" source).
- **filter\_results** – if `True`, filter the results to exclude dummy synchronizers, which correspond to calls which do not return anything

#### **send\_multicast\_sync** (*dst='any', tag=None, value=None, src=None, timeout=None, default\_result=None, pass\_exception=True*)

Send a multicast to the multicast pool and synchronize the results, if available.

By default, the multicast source is the thread's name. Results are collected and synchronized only from the subscriptions which return them (i.e., set `return_result=True`). Local call method.

##### **Parameters**

- **dst** (*str*) – multicast destination; can be a name, "all" (will pass all subscribers' destination filters), or "any" (will only be passed to subscribers specifically subscribed to multicast with "any" destination).
- **tag** (*str*) – multicast tag.
- **value** – multicast value.
- **src** (*str*) – multicast source; can be `None` (current thread name), a specific name, "all" (will pass all subscribers' source filters), or "any" (will only be passed to subscribers specifically subscribed to multicast with "any" source).
- **timeout** – synchronization timeout (`None` means waiting forever)
- **default\_result** – default result value if synchronization failed (timed out, thread stopped, etc.)

- **pass\_exception** – if `True` and the signal processor raised an exception, raise it in this thread as well. If `pass_exception==True` and the returned value represents exception, re-raise it in the caller thread; otherwise, return *default*.

**set\_variable** (*name, value, update=False, notify=False, notify\_tag='changed/\*', simple=False*)  
Set thread variable.

Can be called in any thread (controlled or external). If `notify==True`, send an multicast with the given *notify\_tag* (where "\*" symbol is replaced by the variable name). If `update==True` and the value is a dictionary, update the branch rather than overwrite it. If `simple==True`, assume that the result is a single atomic variable, in which case the lock is not used; note that in this case the threads waiting on this variable (or branches containing it) will not be notified. Local call method.

**delete\_variable** (*name, missing\_error=False*)  
Delete thread variable.

If `missing_error==False` and no variable exists, do nothing; otherwise, raise an error. Local call method.

**set\_func\_variable** (*name, func, use\_lock=True*)  
Set a 'function' variable.

Acts as a thread variable to the external user, but instead of reading a stored value, it executed a function instead. Note, that the function is executed in the caller thread (i.e., the thread which tries to access the variable), so use of synchronization methods (commands, signals, locks) is highly advised.

If `use_lock==True`, then the function call will be wrapped into the usual variable lock, i.e., it won't run concurrently with other variable access. Local call method.

**add\_thread\_method** (*name, method, interrupt=True*)  
Add a thread method.

Adds a named method to the thread, which can be called later using `call_thread_method()`. This method will be called in this thread.

Useful for GUI thread to set up some global access methods, which other threads can safely use. For `QTaskThread` threads it's a better idea to set up a command instead. Local call method.

**delete\_thread\_method** (*name*)  
Delete a thread method.

Local call method.

**call\_thread\_method** (*name, \*args, \*\*kwargs*)  
Call a thread method.

Method needs to be set up beforehand using `add_thread_method()`. It is always executed in the current thread. Local call method.

**send\_message** (*tag, value, priority=0*)  
Send a message to the thread with a given tag, value and priority.  
External call method.

**send\_interrupt** (*tag, value, priority=0*)  
Send an interrupt message to the thread with a given tag, value and priority.  
External call method.

**send\_sync** (*tag, uid*)  
Send a synchronization signal with the given tag and UID.

This method is rarely invoked directly, and is usually used by synchronizers code (e.g., `QThreadNotifier`). External call method.

**get\_variable** (*name*, *default=None*, *copy\_branch=True*, *missing\_error=False*, *simple=False*)

Get thread variable.

If *missing\_error==False* and no variable exists, return *default*; otherwise, raise an error. If *copy\_branch==True* and the variable is a *Dictionary* branch, return its copy to ensure that it stays unaffected on possible further variable assignments. If *simple==True*, assume that the result is a single atomic variable, in which case the lock is not used; this only works with actual variables and not function variables. Universal call method.

**sync\_variable** (*name*, *pred*, *timeout=None*)

Wait until thread variable with the given *name* satisfies the condition given by *pred*.

*pred* can be a variable value, a container (list, set, tuple) of possible values, or a function which takes one argument (variable value) and returns whether the condition is satisfied. It is executed in the caller thread. External call method.

**start** ()

Start the thread.

External call method.

**request\_stop** ()

Request thread stop (send a stop command).

External call method.

**stop** (*code=0*, *sync=False*)

Stop the thread.

If called from the thread, stop immediately by raising a *threadprop.InterruptExceptionStop* exception. Otherwise, schedule thread stop. If the thread kind is "main", stop the whole application with the given exit code. Otherwise, stop the thread. If *sync==True* and the thread is not main or current, wait until it is completely stopped. Universal call method.

**sync\_stop** ()

Wait until the controller and the thread are stopped.

External call method.

**poke** ()

Send a dummy message to the thread.

A cheap way to notify the thread that something happened (useful for, e.g., making thread leave *wait\_for\_any\_message()* method). External call method.

**running** ()

Check if the thread is running

**finishing** ()

Check if the thread is finishing

**notify\_exec\_point** (*point*)

Mark the given execution point as passed.

Automatically invoked points include "start" (thread starting), "run" (thread setup and ready to run), "cleanup" (thread stopping is invoked, starting to clean up) and "stop" (thread finished). Can be extended for arbitrary points. Local call method.

**fail\_exec\_point** (*point*)

Mark the given execution point as failed.



Automatically invoked for "run" (thread setup and ready to run) if the startup raised an error before the thread properly started ("start", "cleanup", and "stop" are notified in any case) Can be extended for arbitrary points. Local call method.

**get\_exec\_counter** (*point*)

Get the counter (number of notifications) for the given point.

See `sync_exec_point()` for details. External call.

**sync\_exec\_point** (*point*, *timeout=None*, *counter=1*)

Wait for the given execution point.

Automatically invoked points include "start" (thread starting), "run" (thread setup and ready to run), "cleanup" (thread stopping is invoked, starting to clean up) and "stop" (thread finished). If *timeout* is passed, raise `threadprop.TimeoutThreadError`. *counter* specifies the minimal number of pre-requisite `notify_exec_point()` calls to finish the waiting (by default, a single call is enough). Return actual number of notifier calls up to date. External call method.

**add\_stop\_notifier** (*func*, *call\_if\_stopped=True*)

Add stop notifier: a function which is called when the thread is about to be stopped (left the main message loop).

The supplied function is called in the controlled thread close to its shutdown, so it should be short, non-blocking, and thread-safe. If the thread is already stopped and *call\_if\_stopped==True*, call *func* immediately (from the caller's thread). Return `True` if the thread is still running and the notifier is added, and `False` otherwise. Local call method.

**remove\_stop\_notifier** (*func*)

Remove the stop notifier from this controller.

Return `True` if the notifier was in this thread and is now removed, and `False` otherwise. Local call method.

**is\_in\_controlled** ()

Check if the thread executing this code is controlled by this controller

**call\_in\_thread\_callback** (*func*, *args=None*, *kwargs=None*, *callback=None*, *tag=None*, *priority=0*, *interrupt=True*)

Call a function in this thread with the given arguments.

If *callback* is supplied, call it with the result as a single argument (call happens in the controller thread). If *tag* is supplied, send the call in a message with the given tag; otherwise, use the interrupt call (generally, higher priority method). If *interrupt==True*, method can be called inside any control loop (either main loop, or during waiting); otherwise, only call it in the top loop. Universal call method.

**call\_in\_thread\_sync** (*func*, *args=None*, *kwargs=None*, *sync=True*, *callback=None*, *timeout=None*, *default\_result=None*, *pass\_exception=True*, *silent=False*, *tag=None*, *priority=0*, *interrupt=True*, *error\_on\_stopped=True*, *same\_thread\_shortcut=True*)

Call a function in this thread with the given arguments.

If *sync==True*, calling thread is blocked until the controlled thread executes the function, and the function result is returned (in essence, the fact that the function executes in a different thread is transparent). Otherwise, exit call immediately, and return a synchronizer object (`QCallResultSynchronizer`), which can be used to check if the call is done (method `is_done`) and obtain the result (method `QCallResultSynchronizer.get_value_sync()`). If *callback* is not `None`, call it after the function is successfully executed (from the target thread), with a single parameter being function result. If *pass\_exception==True* and *func* raises an exception, re-raise it in the caller thread (applies only if *sync==True*). If *silent==True* and *func* raises an exception, silence it in the execution thread and only re-raise it in the caller thread; note that if *pass\_exception==False* and *silent==True*, the exception is ignored in both threads. If *tag* is supplied, send the call in a message with the given tag and



priority; otherwise, use the interrupt call (generally, higher priority method). If `interrupt==True`, method can be called inside any control loop (either main loop, or during waiting); otherwise, only call it in the top loop. If `error_on_stopped==True` and the controlled thread is stopped before it executed the call, raise `threadprop.NoControllerThreadError`; otherwise, return `default_result`. If `same_thread_shortcut==True` (default), the call is synchronous, and the caller thread is the same as the controlled thread, call the function directly. Universal call method.

```
class pylablib.core.thread.controller.QTaskThread(name=None, args=None,  
                                                kwargs=None, multi-  
                                                cast_pool=None)
```

Bases: `pylablib.core.thread.controller.QThreadController`

Thread which allows to set up and run jobs and batch jobs with a certain time period, and execute commands in the meantime.

#### Parameters

- **name** (*str*) – thread name (by default, generate a new unique name)
- **args** – args supplied to `setup_task()` method
- **kwargs** – keyword args supplied to `setup_task()` method
- **multicast\_pool** – `MulticastPool` for this thread (by default, use the default common pool)

#### ca

asynchronous command accessor, which makes calls more function-like; `ctl.ca.comm(*args, **kwarg)` is equivalent to `ctl.call_command("comm", args, kwargs, sync=False)`

#### cai

asynchronous command accessor which ignores and silences any exceptions (including missing /stopped controller) useful for sending queries during thread finalizing / application shutdown, when it's not guaranteed that the command recipient is running `ctl.cai.comm(*args, **kwarg)` is equivalent to `ctl.call_command("comm", args, kwargs, sync=False, ignore_errors=True)`

#### cad

asynchronous command accessor returning a result synchronizer, which makes calls more function-like; `ctl.cad.comm(*args, **kwarg)` is equivalent to `ctl.call_command("comm", args, kwargs, sync="delayed")`

#### cs

synchronous command accessor, which makes calls more function-like; `ctl.cs.comm(*args, **kwarg)` is equivalent to `ctl.call_command("comm", args, kwargs, sync=True)`

#### css

synchronous command accessor which is made 'exception-safe' via `exsafe()` wrapper (i.e., safe to directly connect to slots) `ctl.css.comm(*args, **kwarg)` is equivalent to `with exint():`  
`ctl.call_command("comm", args, kwargs, sync=True)`

#### csi

synchronous command accessor which ignores and silences any exceptions (including missing /stopped controller) useful for sending queries during thread finalizing / application shutdown, when it's not guaranteed that the command recipient is running

#### m

method accessor; directly calls the method corresponding to the command; `ctl.m.comm(*args, **kwarg)` is equivalent to `ctl.call_command("comm", *args, **kwargs)`, which is often also equivalent to `ctl.comm(*args, **kwargs)`; for most practical purposes it's the same as directly invoking the class method, but it makes intent more explicit (as command methods are usually

not called directly from other threads), and it doesn't invoke warning about calling method instead of command from another thread.

**Methods to overload:**

- `setup_task()`: executed on the thread startup (between synchronization points "start" and "run")
- `finalize_task()`: executed on thread cleanup (attempts to execute in any case, including exceptions)

**class** `TBatchJob` (*job, cleanup, min\_run\_time, priority*)

Bases: `tuple`

**cleanup**

**count** ()

Return number of occurrences of value.

**index** ()

Return first index of value.

Raises `ValueError` if the value is not present.

**job**

**min\_run\_time**

**priority**

**class** `TCommand` (*command, scheduler, priority*)

Bases: `tuple`

**command**

**count** ()

Return number of occurrences of value.

**index** ()

Return first index of value.

Raises `ValueError` if the value is not present.

**priority**

**scheduler**

**class** `Job` (*job, period, queue, jobs\_order*)

Bases: `object`

A single job loop.

Deals with scheduling, time counting, pausing, and cleanup.

**Parameters**

- **job** – job function
- **period** – job period
- **queue** – thread controller's scheduling queue, to which the job must be added
- **jobs\_order** – thread controller's job queue which determines the jobs scheduling order

**schedule()**  
Schedule the job

**mark\_unscheduled()**  
Mark the job as unscheduled.  
Called automatically on job completion.

**unschedule()**  
Manually unschedule the job (e.g., when paused or removed)

**clear()**  
Clear the job and remove it from the jobs list

**change\_period(*period*)**  
Change the job period

**pause(*paused=True, unschedule=True*)**  
Pause or resume the job.  
If pausing and `unschedule==True`, remove already scheduled job from the queue.

**time\_left(*t=None*)**  
Get the amount of time left till the next call, or `None` if the job is paused

**add\_job(*name, job, period, initial\_call=True, priority=-10*)**  
Add a recurrent *job* which is called every *period* seconds.  
  
The job starts running automatically when the main thread loop start executing. If `initial_call==True`, call *job* once immediately after adding. *priority* specifies the call priority in the scheduling queue; by default, it is lower than the command and multicasts (0). Local call method.

**change\_job\_period(*name, period*)**  
Change the period of the job *name*.  
  
Local call method.

**remove\_job(*name*)**  
Remove the job *name* from the job list.  
  
Local call method.

**add\_batch\_job(*name, job, cleanup=None, min\_runtime=0, priority=-10*)**  
Add a batch *job* which is executed once, but with continuations.  
  
After this call the job is just created, but is not running. To start it, call `start_batch_job()`. If specified, *cleanup* is a finalizing function which is called both when the job terminates normally, and when it is forcibly stopped (including thread termination). *min\_runtime* specifies minimal expected runtime of a job; if a job executes faster than this time, it is repeated again unless at least *min\_runtime* seconds passed; useful for high-throughput jobs, as it reduces overhead from the job scheduling mechanism (repeating within *min\_runtime* time window is fast)  
  
Unlike the usual recurrent jobs, here *job* is a generator (usually defined by a function with `yield` statement). When the job is running, the generator is periodically called until it raises `StopIteration` exception, which signifies that the job is finished. From generator function point of view, after the job is started, the function is executed normally, but every time `yield` statement is encountered, the execution is suspended for *period* seconds (specified in `start_batch_job()`). *priority* specifies the call priority in the scheduling queue; by default, it is lower than the command and multicasts (0). Local call method.

**change\_batch\_job\_parameters(*name, job='keep', cleanup='keep', min\_runtime='keep', priority='keep', stop=False, restart=False*)**  
Change parameters (main body, cleanup function, and minimal runtime) of the batch job.

The parameters are the same as for `add_batch_job()`. If any of them are "keep", don't change them. If `stop==True`, stop the job before changing the parameters; otherwise the job is continued with the previous parameters (including cleanup) until it is stopped and restarted. If `restart==True`, restart the job after changing the parameters. Local call method.

**remove\_batch\_job** (*name*)

Remove the batch job *name*, stopping it if necessary.

Local call method.

**start\_batch\_job** (*name*, *period*, *\*args*, *start\_immediate=True*, *\*\*kwargs*)

Start the batch job with the given name.

*period* specifies suspension period. Optional arguments are passed to the job and the cleanup functions. If `start_immediate==True`, start the job (i.e., run the first iteration) immediately during the call; otherwise, start it only when it is scheduled, after the currently running call is complete. Local call method.

**is\_batch\_job\_running** (*name*)

Check if a given batch job running.

Local call method.

**stop\_batch\_job** (*name*, *stop\_immediate=True*, *error\_on\_stopped=False*)

Stop a given batch job.

If `error_on_stopped==True` and the job is not currently running, raise an error. Otherwise, do nothing. If `stop_immediate==True`, stop the job (i.e., unschedule it and run the cleanup code) immediately during the call; otherwise, stop it when its next iteration is called. Local call method.

**restart\_batch\_job** (*name*, *start\_immediate=True*, *error\_on\_stopped=False*)

Restart the running batch job with its current arguments.

If `error_on_stopped==True` and the job is not currently running, raise an error. Otherwise, do nothing. Local call method.

**run\_as\_batch\_job** (*job*, *period*, *cleanup=None*, *name=None*, *priority=-10*, *start\_immediate=True*, *args=None*, *kwargs=None*)

Create a temporarily batch job and immediately run it.

If *name* is `None`, generate a new unique name. The job is removed after it is complete (i.e., after cleanup). Note that this implies, that it can not be restarted using `restart_batch_job()`, as it will be removed after the stopping before the restart. All the parameters are the same as for `add_batch_job()` and `start_batch_job()`. Return the batch job name (either supplied or newly generated).

**run** ()

Method called to run the main thread code (only for "run" thread kind).

Local call method, called automatically.

**on\_start** ()

Method invoked on the start of the thread.

Local call method, called automatically.

**on\_finish** ()

Method invoked in the end of the thread.

Called regardless of the stopping reason (normal finishing, exception, application finishing). Local call method, called automatically.

**setup\_task** (*\*args*, *\*\*kwargs*)

Setup the thread (called before the main task loop).

Local call method, called automatically.

**finalize\_task()**

Finalize the thread (always called on thread termination, regardless of the reason).

Local call method, called automatically.

**update\_status**(*kind, status, text=None, notify=True*)

Update status represented in thread variables.

*kind* is the status kind and *status* is its value. Status variable name is "status/"+*kind*. If *text* is not *None*, it specifies new status text stored in "status/"+*kind*+"\_text". If *notify*==*True*, send an multicast about the status change. Local call method.

**add\_command**(*name, command=None, scheduler=None, limit\_queue=None, on\_full\_queue='skip\_current', priority=0*)

Add a new command to the command set.

Return scheduler, which can be used for adding another command (if the same queue should be used for several commands). Local call method.

#### Parameters

- **name** – command name
- **command** – command function; if *None*, look for the method with the given *name*.
- **scheduler** – a command scheduler; by default, it is a *QQueueLengthLimitScheduler*, which maintains a call queue with the given length limit and full queue behavior; can also be a name of a different command, with which it will share a single queue with the same limitations; if supplied, *limit\_queue* and *on\_full\_queue* parameters are ignored
- **limit\_queue** – command call queue limit; *None* means no limit
- **on\_full\_queue** – action to be taken if the call can't be scheduled (the queue is full); can be "skip\_current" (skip the call which is being scheduled), "skip\_newest" (skip the most recent call; place the current), "skip\_oldest" (skip the oldest call in the queue; place the current), "call\_current" (execute the call which is being scheduled immediately in the caller thread), "call\_newest" (execute the most recent call immediately in the caller thread), "call\_oldest" (execute the oldest call in the queue immediately in the caller thread), or "wait" (wait until the call can be scheduled, which is checked after every call removal from the queue; place the call)
- **priority** – command priority; higher-priority multicasts and commands are always executed before the lower-priority ones.

**add\_direct\_call\_command**(*name, command=None, error\_on\_async=True*)

Add a direct method call which appears as a command.

Unlike regular commands, the call is executed directly in the caller thread (i.e., it is identical to the direct method call). Useful for lightweight and/or lock-wrapped methods, which can be called in a thread-safe way, but which still use command interface for consistency. Note that this kind of commands doesn't have the same level of synchronization as regular commands (e.g., it can be executed during execution of another command, or commsync multicast method). Local call method.

#### Parameters

- **name** – command name
- **command** – command function; if *None*, look for the method with the given *name*.

- **error\_on\_async** – if `True` and the command is called asynchronously, raise an error; otherwise, substitute for a synchronous call

**subscribe\_commsync** (*callback*, *srcs*='any', *tags*=None, *dsts*='any', *filt*=None, *subscription\_priority*=0, *scheduler*=None, *limit\_queue*=None, *on\_full\_queue*='skip\_current', *priority*=0, *add\_call\_info*=False, *return\_result*=False, *sid*=None)

Subscribe a callback to a multicast which is synchronized with commands and jobs execution.

Unlike the standard `QThreadController.subscribe_sync()` method, the subscribed callback will only be executed between jobs or commands, not during one of these. Local call method.

#### Parameters

- **callback** – callback function, which takes 3 arguments: source, tag, and value.
- **srcs** (*str* or [*str*]) – multicast source name or list of source names to filter the subscription; can be "any" (any source) or "all" (only multicasts specifically having "all" as a source).
- **tags** – multicast tag or list of tags to filter the subscription (any tag by default); can also contain Unix shell style pattern ("\*" matches everything, "?" matches one symbol, etc.)
- **dsts** (*str* or [*str*]) – multicast destination name or list of destination names to filter the subscription; can be "any" (any destination) or "all" (only source specifically having "all" as a destination).
- **filt** (*callable*) – additional filter function which takes 4 arguments: source, destination, tag, and value, and checks whether multicast passes the requirements.
- **subscription\_priority** (*int*) – subscription priority (higher priority subscribers are called first).
- **scheduler** – if defined, multicast call gets scheduled using this scheduler; by default, create a new call queue scheduler with the given *limit\_queue*, *on\_full\_queue* and *add\_call\_info* arguments.
- **limit\_queue** (*int*) – limits the maximal number of scheduled calls (if the multicast is sent while at least *limit\_queue* callbacks are already in queue to be executed, ignore it) 0 or negative value means no limit (not recommended, as it can increase the queue indefinitely if the multicast rate is high enough)
- **on\_full\_queue** – action to be taken if the call can't be scheduled (the queue is full); can be "skip\_current" (skip the call which is being scheduled), "skip\_newest" (skip the most recent call; place the current), "skip\_oldest" (skip the oldest call in the queue; place the current), "call\_current" (execute the call which is being scheduled immediately in the caller thread), "call\_newest" (execute the most recent call immediately in the caller thread), "call\_oldest" (execute the oldest call in the queue immediately in the caller thread), or "wait" (wait until the call can be scheduled, which is checked after every call removal from the queue; place the call)
- **add\_call\_info** (*bool*) – if `True`, add a fourth argument containing a call information (tuple with a single element, a timestamps of the call).
- **return\_result** – if `True`, use a result synchronizer to return the result of the subscribed call; otherwise, ignore the result
- **sid** (*int*) – subscription ID (by default, generate a new unique id and return it).

**call\_command\_direct** (*name*, *args=None*, *kwargs=None*)

Invoke a command directly and immediately in the current thread.

Universal call method.

**call\_command** (*name*, *args=None*, *kwargs=None*, *sync=False*, *callback=None*, *timeout=None*, *ignore\_errors=False*)

Invoke command call with the given name and arguments

If *callback* is not *None*, call it after the command is successfully executed (from the target thread), with a single parameter being the command result. If *sync==True*, pause caller thread execution (for at most *timeout* seconds) until the command has been executed by the target thread, and then return the command result. If *sync=="delayed"*, return *QCallResultSynchronizer* object which can be used to wait for and read the command result; otherwise, return *None*. In the *sync==True* case, if *ignore\_errors==True*, ignore all possible problems with the call (controller stopped, call raised an exception, call was skipped) and return *None* instead; otherwise, these problems raise exceptions in the caller thread. Universal call method.

**call\_in\_thread\_commsync** (*func*, *args=None*, *kwargs=None*, *sync=True*, *timeout=None*, *priority=0*, *ignore\_errors=False*, *same\_thread\_shortcut=True*)

Call a function in this thread such that it is synchronous with other commands, and jobs.

Mostly equivalent to calling a command, only the command function is supplied instead of its name, and the advanced scheduling (maximal schedule size, sharing with different commands, etc.) is not used. *args* and *kwargs* specify the function arguments. If *sync==True*, pause caller thread execution (for at most *timeout* seconds) until the command has been executed by the target thread, and then return the command result. If *sync=="delayed"*, return *QCallResultSynchronizer* object which can be used to wait for and read the command result; otherwise, return *None*. *priority* sets the call priority (by default, the same as the standard commands). In the *sync==True* case, if *ignore\_errors==True*, ignore all possible problems with the call (controller stopped, call raised an exception, call was skipped) and return *None* instead; otherwise, these problems raise exceptions in the caller thread. If *same\_thread\_shortcut==True* (default) and the caller thread is the same as the controlled thread, call the function directly. Universal call method.

**comm\_paused** ()

Context manager, which allows to temporarily pause all calls (commands, jobs, etc.)

**class CommandAccess** (*parent*, *sync*, *direct=False*, *timeout=None*, *safe=False*, *ignore\_errors=False*)

Bases: *object*

Accessor object designed to simplify command syntax.

Automatically created by the thread, so doesn't need to be invoked externally.

**add\_stop\_notifier** (*func*, *call\_if\_stopped=True*)

Add stop notifier: a function which is called when the thread is about to be stopped (left the main message loop).

The supplied function is called in the controlled thread close to its shutdown, so it should be short, non-blocking, and thread-safe. If the thread is already stopped and *call\_if\_stopped==True*, call *func* immediately (from the caller's thread). Return *True* if the thread is still running and the notifier is added, and *False* otherwise. Local call method.

**add\_thread\_method** (*name*, *method*, *interrupt=True*)

Add a thread method.

Adds a named method to the thread, which can be called later using *call\_thread\_method()*. This method will be called in this thread.

Useful for GUI thread to set up some global access methods, which other threads can safely use. For *QTaskThread* threads it's a better idea to set up a command instead. Local call method.



**allowing\_toploop** (*depth=1*)

Context manager which temporarily treats the current loop level and several deeper levels as a top loop.

All event loops which lie up to *depth* below this one are treated as top loops.

**blocking\_control\_signals** (*kinds='all', ignore=None*)

Context manager which temporarily blocks external control signals.

After leaving the wrapped code segment, all of the blocked but not ignored calls are executed. *kind* determines the kind of calls to block; it is a collection of elements among "message", "stop", and "call" and blocks, correspondingly, messages, stop signals, and any `call_in_thread`-related requests; can be also be "all", which includes all of these categories. *ignore* specifies kinds which are completely ignored if sent during the blocking interval; can also be "all", which includes all of the *kinds* categories. Useful to temporarily “suspend” the thread communication with other threads, especially for the main GUI thread (e.g., to show a blocking message box). Local call method.

**call\_in\_thread\_callback** (*func, args=None, kwargs=None, callback=None, tag=None, priority=0, interrupt=True*)

Call a function in this thread with the given arguments.

If *callback* is supplied, call it with the result as a single argument (call happens in the controller thread). If *tag* is supplied, send the call in a message with the given tag; otherwise, use the interrupt call (generally, higher priority method). If `interrupt==True`, method can be called inside any control loop (either main loop, or during waiting); otherwise, only call it in the top loop. Universal call method.

**call\_in\_thread\_sync** (*func, args=None, kwargs=None, sync=True, callback=None, timeout=None, default\_result=None, pass\_exception=True, silent=False, tag=None, priority=0, interrupt=True, error\_on\_stopped=True, same\_thread\_shortcut=True*)

Call a function in this thread with the given arguments.

If `sync==True`, calling thread is blocked until the controlled thread executes the function, and the function result is returned (in essence, the fact that the function executes in a different thread is transparent). Otherwise, exit call immediately, and return a synchronizer object (`QCallResultSynchronizer`), which can be used to check if the call is done (method `is_done`) and obtain the result (method `QCallResultSynchronizer.get_value_sync()`). If *callback* is not `None`, call it after the function is successfully executed (from the target thread), with a single parameter being function result. If `pass_exception==True` and *func* raises an exception, re-raise it in the caller thread (applies only if `sync==True`). If `silent==True` and *func* raises an exception, silence it in the execution thread and only re-raise it in the caller thread; note that if `pass_exception==False` and `silent==True`, the exception is ignored in both threads. If *tag* is supplied, send the call in a message with the given tag and priority; otherwise, use the interrupt call (generally, higher priority method). If `interrupt==True`, method can be called inside any control loop (either main loop, or during waiting); otherwise, only call it in the top loop. If `error_on_stopped==True` and the controlled thread is stopped before it executed the call, raise `threadprop.NoControllerThreadError`; otherwise, return *default\_result*. If `same_thread_shortcut==True` (default), the call is synchronous, and the caller thread is the same as the controlled thread, call the function directly. Universal call method.

**call\_thread\_method** (*name, \*args, \*\*kwargs*)

Call a thread method.

Method needs to be set up beforehand using `add_thread_method()`. It is always executed in the current thread. Local call method.

**check\_messages** (*top\_loop=False*)

Receive new messages.

Runs the underlying message loop to process newly received message and signals (and place them in corresponding queues if necessary). This method is rarely invoked, and only should be used periodically



during long computations to not ‘freeze’ the thread. If `top_loop==True`, treat the waiting as the top message loop (i.e., any top loop message or signal can be executed here). Local call method.

**delete\_thread\_method** (*name*)

Delete a thread method.

Local call method.

**delete\_variable** (*name*, *missing\_error=False*)

Delete thread variable.

If `missing_error==False` and no variable exists, do nothing; otherwise, raise and error. Local call method.

**fail\_exec\_point** (*point*)

Mark the given execution point as failed.

Automatically invoked for "run" (thread setup and ready to run) if the startup raised an error before the thread properly started ("start", "cleanup", and "stop" are notified in any case) Can be extended for arbitrary points. Local call method.

**finished** = `<Mock name='mock.QtCore.pyqtSignal()' id='140096333774352'>`

**finishing** ()

Check if the thread is finishing

**get\_exec\_counter** (*point*)

Get the counter (number of notifications) for the given point.

See `sync_exec_point()` for details. External call.

**get\_variable** (*name*, *default=None*, *copy\_branch=True*, *missing\_error=False*, *simple=False*)

Get thread variable.

If `missing_error==False` and no variable exists, return *default*; otherwise, raise and error. If `copy_branch==True` and the variable is a *Dictionary* branch, return its copy to ensure that it stays unaffected on possible further variable assignments. If `simple==True`, assume that the result is a single atomic variable, in which case the lock is not used; this only works with actual variables and not function variables. Universal call method.

**is\_in\_controlled** ()

Check if the thread executing this code is controlled by this controller

**new\_messages\_number** (*tag*)

Get the number of queued messages with a given tag.

Local call method.

**no\_stopping** ()

Context manager, which temporarily suspends stop requests (*InterruptExceptionStop* exceptions).

If the stop request has been made within this block, raise the exception on exit. Note that `stop()` method and, correspondingly, `stop_controller()` still work, when called from the controlled thread.

**notify\_exec\_point** (*point*)

Mark the given execution point as passed.

Automatically invoked points include "start" (thread starting), "run" (thread setup and ready to run), "cleanup" (thread stopping is invoked, starting to clean up) and "stop" (thread finished). Can be extended for arbitrary points. Local call method.

**poke** ()

Send a dummy message to the thread.

A cheap way to notify the thread that something happened (useful for, e.g., making thread leave `wait_for_any_message()` method). External call method.

**pop\_message** (*tag*)

Pop the latest message with the given tag.

Select the message with the highest priority, and among those the oldest one. If no messages are available, raise `threadprop.NoMessageThreadError`. Local call method.

**process\_interrupt** (*tag, value*)

Process a new interrupt.

If the function returns `False`, the interrupt is put in the corresponding queue. Otherwise, the the message is interrupt to be already, and it gets ‘absorbed’. Local call method, called automatically.

**process\_message** (*tag, value*)

Process a new message.

If the function returns `False`, the message is put in the corresponding queue. Otherwise, the the message is considered to be already, and it gets ‘absorbed’. Local call method, called automatically.

**remove\_stop\_notifier** (*func*)

Remove the stop notifier from this controller.

Return `True` if the notifier was in this thread and is now removed, and `False` otherwise. Local call method.

**request\_stop** ()

Request thread stop (send a stop command).

External call method.

**running** ()

Check if the thread is running

**send\_interrupt** (*tag, value, priority=0*)

Send an interrupt message to the thread with a given tag, value and priority.

External call method.

**send\_message** (*tag, value, priority=0*)

Send a message to the thread with a given tag, value and priority.

External call method.

**send\_multicast** (*dst='any', tag=None, value=None, src=None, filter\_results=True*)

Send a multicast to the multicast pool.

By default, the multicast source is the thread’s name. Return result synchronizers for all executed subscribed methods. Local call method.

#### Parameters

- **dst** (*str*) – multicast destination; can be a name, "all" (will pass all subscribers’ destination filters), or "any" (will only be passed to subscribers specifically subscribed to multicast with "any" destination).
- **tag** (*str*) – multicast tag.
- **value** – multicast value.
- **src** (*str*) – multicast source; can be `None` (current thread name), a specific name, "all" (will pass all subscribers’ source filters), or "any" (will only be passed to subscribers specifically subscribed to multicast with "any" source).

- **filter\_results** – if `True`, filter the results to exclude dummy synchronizers, which correspond to calls which do not return anything

**send\_multicast\_sync** (*dst='any', tag=None, value=None, src=None, timeout=None, default\_result=None, pass\_exception=True*)

Send a multicast to the multicast pool and synchronize the results, if available.

By default, the multicast source is the thread's name. Results are collected and synchronized only from the subscriptions which return them (i.e., set `return_result=True`). Local call method.

#### Parameters

- **dst** (*str*) – multicast destination; can be a name, "all" (will pass all subscribers' destination filters), or "any" (will only be passed to subscribers specifically subscribed to multicast with "any" destination).
- **tag** (*str*) – multicast tag.
- **value** – multicast value.
- **src** (*str*) – multicast source; can be `None` (current thread name), a specific name, "all" (will pass all subscribers' source filters), or "any" (will only be passed to subscribers specifically subscribed to multicast with "any" source).
- **timeout** – synchronization timeout (`None` means waiting forever)
- **default\_result** – default result value if synchronization failed (timed out, thread stopped, etc.)
- **pass\_exception** – if `True` and the signal processor raised an exception, raise it in this thread as well. If `pass_exception==True` and the returned value represents exception, re-raise it in the caller thread; otherwise, return *default*.

**send\_sync** (*tag, uid*)

Send a synchronization signal with the given tag and UID.

This method is rarely invoked directly, and is usually used by synchronizers code (e.g., `QThreadNotifier`). External call method.

**set\_func\_variable** (*name, func, use\_lock=True*)

Set a 'function' variable.

Acts as a thread variable to the external user, but instead of reading a stored value, it executed a function instead. Note, that the function is executed in the caller thread (i.e., the thread which tries to access the variable), so use of synchronization methods (commands, signals, locks) is highly advised.

If `use_lock==True`, then the function call will be wrapped into the usual variable lock, i.e., it won't run concurrently with other variable access. Local call method.

**set\_variable** (*name, value, update=False, notify=False, notify\_tag='changed/\*', simple=False*)

Set thread variable.

Can be called in any thread (controlled or external). If `notify==True`, send an multicast with the given *notify\_tag* (where "\*" symbol is replaced by the variable name). If `update==True` and the value is a dictionary, update the branch rather than overwrite it. If `simple==True`, assume that the result is a single atomic variable, in which case the lock is not used; note that in this case the threads waiting on this variable (or branches containing it) will not be notified. Local call method.

**sleep** (*timeout, wake\_on\_message=False, top\_loop=False*)

Sleep for a given time (in seconds).

Unlike `time.sleep()`, constantly checks the event loop for new messages (e.g., if stop or interrupt commands are issued). In addition, if `wake_on_message==True`, wake up if any message has been received; in this case, return `True` if the wait has been completed, and `False` if it has been interrupted by

a message. If `top_loop==True`, treat the waiting as the top message loop (i.e., any top loop message or signal can be executed here). If `timeout` is `None`, wait forever (usually, until the application is closed, or some interrupt message raises an error). Local call method.

**start** ()

Start the thread.

External call method.

**started** = `<Mock name='mock.QtCore.pyqtSignal()' id='140096333774352'>`

**stop** (*code=0, sync=False*)

Stop the thread.

If called from the thread, stop immediately by raising a `threadprop.InterruptExceptionStop` exception. Otherwise, schedule thread stop. If the thread kind is "main", stop the whole application with the given exit code. Otherwise, stop the thread. If `sync==True` and the thread is not main or current, wait until it is completely stopped. Universal call method.

**subscribe\_direct** (*callback, srcs='any', tags=None, dsts='any', filt=None, subscription\_priority=0, scheduler=None, return\_result=False, sid=None*)

Subscribe asynchronous callback to a multicast.

If a multicast is sent, *callback* is called from the sending thread (not subscribed thread). Therefore, should be used with care. In Qt, analogous to making a signal connection with a direct call. By default, the subscribed destination is the thread's name. Local call method.

#### Parameters

- **callback** – callback function, which takes 3 arguments: source, tag, and value.
- **srcs** (*str or [str]*) – multicast source name or list of source names to filter the subscription; can be "any" (any source) or "all" (only multicasts specifically having "all" as a source).
- **tags** – multicast tag or list of tags to filter the subscription (any tag by default); can also contain Unix shell style pattern ("\*" matches everything, "?" matches one symbol, etc.)
- **dsts** (*str or [str]*) – multicast destination name or list of destination names to filter the subscription; can be "any" (any destination) or "all" (only source specifically having "all" as a destination).
- **filt** (*callable*) – additional filter function which takes 4 arguments: source, destination, tag, and value, and checks whether multicast passes the requirements.
- **subscription\_priority** (*int*) – subscription priority (higher priority subscribers are called first).
- **scheduler** – if defined, multicast call gets scheduled using this scheduler instead of being called directly (which is the default behavior)
- **return\_result** – if `True`, use a result synchronizer to return the result of the subscribed call; otherwise, ignore the result
- **sid** (*int*) – subscription ID (by default, generate a new unique id and return it).

**subscribe\_sync** (*callback, srcs='any', tags=None, dsts='any', filt=None, subscription\_priority=0, limit\_queue=None, call\_interrupt=True, add\_call\_info=False, return\_result=False, sid=None*)

Subscribe a synchronous callback to a multicast.

If a multicast is sent, *callback* is called from the *dest\_controller* thread (by default, thread which is calling this function) via the thread call mechanism (*QThreadController.call\_in\_thread\_callback()*). In Qt, analogous to making a signal connection with a queued call. By default, the subscribed destination is the thread's name. Local call method.

#### Parameters

- **callback** – callback function, which takes 3 arguments: source, tag, and value.
- **srcs** (*str* or [*str*]) – multicast source name or list of source names to filter the subscription; can be "any" (any source) or "all" (only multicasts specifically having "all" as a source).
- **tags** – multicast tag or list of tags to filter the subscription (any tag by default); can also contain Unix shell style pattern ("\*" matches everything, "?" matches one symbol, etc.)
- **dsts** (*str* or [*str*]) – multicast destination name or list of destination names to filter the subscription; can be "any" (any destination) or "all" (only source specifically having "all" as a destination).
- **filt** (*callable*) – additional filter function which takes 4 arguments: source, destination, tag, and value, and checks whether multicast passes the requirements.
- **subscription\_priority** (*int*) – subscription priority (higher priority subscribers are called first).
- **limit\_queue** (*int*) – limits the maximal number of scheduled calls (if the multicast is sent while at least *limit\_queue* callbacks are already in queue to be executed, ignore it) 0 or negative value means no limit (not recommended, as it can increase the queue indefinitely if the multicast rate is high enough)
- **call\_interrupt** – whether the call is an interrupt (call inside any loop, e.g., during waiting or sleeping), or it should be called in the main event loop
- **add\_call\_info** (*bool*) – if True, add a fourth argument containing a call information (tuple with a single element, a timestamps of the call).
- **return\_result** – if True, use a result synchronizer to return the result of the subscribed call; otherwise, ignore the result
- **sid** (*int*) – subscription ID (by default, generate a new unique name).

**sync\_exec\_point** (*point*, *timeout=None*, *counter=1*)

Wait for the given execution point.

Automatically invoked points include "start" (thread starting), "run" (thread setup and ready to run), "cleanup" (thread stopping is invoked, starting to clean up) and "stop" (thread finished). If *timeout* is passed, raise *threadprop.TimeoutThreadError*. *counter* specifies the minimal number of pre-requisite *notify\_exec\_point()* calls to finish the waiting (by default, a single call is enough). Return actual number of notifier calls up to date. External call method.

**sync\_stop** ()

Wait until the controller and the thread are stopped.

External call method.

**sync\_variable** (*name*, *pred*, *timeout=None*)

Wait until thread variable with the given *name* satisfies the condition given by *pred*.

*pred* can be a variable values, a container (list, set, tuple) of possible values, or a function which takes one argument (variable value) and returns whether the condition is satisfied. It is executed in the caller thread. External call method.

**unsubscribe** (*sid*)

Unsubscribe from a subscription with a given ID.

Note that multicasts which are already emitted but not processed will remain in the queue; if they need to be ignored, it should be handled explicitly. Local call method.

**wait\_for\_any\_message** (*timeout=None, top\_loop=False*)

Wait for any message (including synchronization messages or pokes).

If timeout is passed, raise `threadprop.TimeoutThreadError`. If `top_loop==True`, treat the waiting as the top message loop (i.e., any top loop message or signal can be executed here). Local call method.

**wait\_for\_message** (*tag, timeout=None, top\_loop=False*)

Wait for a single message with a given tag.

Return value of a received message with this tag. If timeout is passed, raise `threadprop.TimeoutThreadError`. If `top_loop==True`, treat the waiting as the top message loop (i.e., any top loop message or signal can be executed here). Local call method.

**wait\_for\_sync** (*tag, uid, timeout=None*)

Wait for synchronization signal with the given tag and UID.

This method is rarely invoked directly, and is usually used by synchronizers code. If timeout is passed, raise `threadprop.TimeoutThreadError`. Local call method.

**wait\_until** (*check, timeout=None, top\_loop=False*)

Wait until a given condition is true.

Condition is given by the *check* function, which is called after every new received message and should return `True` if the condition is met. If `top_loop==True`, treat the waiting as the top message loop (i.e., any top loop message or signal can be executed here). If timeout is passed, raise `threadprop.TimeoutThreadError`. Local call method.

`pylablib.core.thread.controller.get_controller` (*name=None, sync=True, timeout=None, sync\_point=None*)

Find a controller with a given name.

If *name* is not supplied, yield current controller instead. If *name* is of `int` type, interpret it as a thread id. If the controller is not present and `sync==True`, wait (with the given timeout) until the controller is running; otherwise, raise error if the controller is not running. If *sync\_point* is not `None`, synchronize to the thread *sync\_point* point (by default, "run", i.e., after the setup is done) before returning.

`pylablib.core.thread.controller.sync_controller` (*name, sync\_point='run', timeout=None*)

Find a controller with a given name and synchronize to the given point.

If the controller is not present and `sync==True`, wait (with the given timeout) until the controller is running; otherwise, raise error if the controller is not running. Analogous to `get_controller(name, sync=True, timeout=timeout, sync_point=sync_point)`.

`pylablib.core.thread.controller.get_gui_controller` (*sync=False, timeout=None, create\_if\_missing=True*)

Get GUI thread controller.

If the controller is not present and `sync==True`, wait (with the given timeout) until the controller is running. If the controller is still not present and `create_if_missing==True`, initialize the standard GUI controller.

`pylablib.core.thread.controller.stop_controller` (*name=None, code=0, sync=True, require\_controller=False*)

Stop a controller with a given name (current controller by default).

*code* specifies controller exit code (only applies to the main thread controller). If `require_controller==True` and the controller is not present, raise an error; otherwise, do nothing. If `sync==True`, wait until the controller is stopped.

```
pylablib.core.thread.controller.stop_all_controllers(sync=True, concurrent=True,
                                                    stop_self=True)
```

Stop all running threads.

If `sync==True`, wait until all of the controllers are stopped. If `sync==True` and `concurrent==True` stop threads in concurrent manner (first issue stop messages to all of them, then wait until all are stopped). If `sync==True` and `concurrent==False` stop threads in consecutive manner (wait for each thread to stop before stopping the next one). If `stop_self==True` stop current thread after stopping all other threads.

```
pylablib.core.thread.controller.stop_app(code=0, sync=False)
Initialize stopping the application.
```

Do this either by stopping the GUI controller (if it exists), or by stopping all controllers. If `sync` is `True` and the thread is not the main one, wait at this point until the process is stopped during the app shutdown; otherwise, the execution will continue as normal, and the thread will be stopped at a later time during the app shutdown.

```
pylablib.core.thread.controller.restart_app(code=0, sync=False)
Restart the application.
```

Equivalent to `stop_app()` followed by the script restart. If `sync` is `True` and the thread is not the main one, wait at this point until the process is stopped during the app shutdown; otherwise, the execution will continue as normal, and the thread will be stopped at a later time during the app shutdown.

## pylablib.core.thread.multicast\_pool module

```
class pylablib.core.thread.multicast_pool.TMulticast(src, tag, value)
Bases: tuple
```

```
count()
Return number of occurrences of value.
```

```
index()
Return first index of value.

Raises ValueError if the value is not present.
```

```
src
```

```
tag
```

```
value
```

```
class pylablib.core.thread.multicast_pool.MulticastPool
Bases: object
```

Multicast dispatcher (somewhat similar in functionality to Qt signals).

Manages dispatching multicasts between sources and destinations (callback functions). Each multicast has defined source, destination (both can also be "all" or "any", see methods descriptions for details), tag and value. Any thread can send a multicast or subscribe for a multicast with given filters (source, destination, tag, additional filters). If a multicast is emitted, it is checked against filters for all subscribers, and the passing ones are then called.

```
subscribe_direct(callback, srcs='any', dsts='any', tags=None, filt=None, priority=0, scheduler=None, return_result=False, sid=None)
Subscribe an asynchronous callback to a multicast.
```



If a multicast is sent, *callback* is called from the sending thread (not subscribed thread). Therefore, should be used with care. In Qt, analogous to making a signal connection with a direct call.

#### Parameters

- **callback** – callback function, which takes 3 arguments: source, tag, and value.
- **srcs** (*str* or [*str*]) – multicast source name or list of source names to filter the subscription; can be "any" (any source) or "all" (only multicasts specifically having "all" as a source).
- **dsts** (*str* or [*str*]) – multicast destination name or list of destination names to filter the subscription; can be "any" (any destination) or "all" (only source specifically having "all" as a destination).
- **tags** – multicast tag or list of tags to filter the subscription (any tag by default); can also contain Unix shell style pattern ("\*" matches everything, "?" matches one symbol, etc.)
- **filt** (*callable*) – additional filter function which takes 4 arguments: source, destination, tag, and value, and checks whether multicast passes the requirements.
- **priority** (*int*) – subscription priority (higher priority subscribers are called first).
- **scheduler** – if defined, multicast call gets scheduled using this scheduler instead of being called directly (which is the default behavior)
- **return\_result** – if `True`, use a result synchronizer to return the result of the subscribed call; otherwise, ignore the result
- **sid** (*int*) – subscription ID (by default, generate a new unique name).

**Returns** subscription ID, which can be used to unsubscribe later.

**unsubscribe** (*sid*)

Unsubscribe from a subscription with a given ID

**send** (*src*, *dst*=*'any'*, *tag*=*None*, *value*=*None*)

Send a multicast.

#### Parameters

- **src** (*str*) – multicast source; can be a name, "all" (will pass all subscribers' source filters), or "any" (will only be passed to subscribers specifically subscribed to multicasts with "any" source).
- **dst** (*str*) – multicast destination; can be a name, "all" (will pass all subscribers' destination filters), or "any" (will only be passed to subscribers specifically subscribed to multicasts with "any" destination).
- **tag** (*str*) – multicast tag.
- **value** – multicast value.

### pylablib.core.thread.notifier module

**class** pylablib.core.thread.notifier.**ISkipableNotifier** (*skippable*=*False*)

Bases: `object`

Generic skippable notifier.

The main methods are `wait()` (wait until the event happened) and `notify()` (notify that the event happened). Only calls underlying waiting and notifying methods once, duplicate calls are ignored.



Parameters **skippable** (*bool*) – if `True`, allows for skippable wait events (if `notify()` is called before `wait()`, neither methods are actually called).

**wait** (\*args, \*\*kwargs)

Wait for the notification.

Can only be called once per notifier lifetime. If the notifier allows skipping, and this method is called after `notify()`, return immediately.

**notify** (\*args, \*\*kwargs)

Notify the waiting process.

Can only be called once per notifier lifetime. If the notifier allows skipping, and this method is called before `wait()`, return immediately.

**waiting** ()

Check if waiting is in progress

**done\_wait** ()

Check if waiting is done

**success\_wait** ()

Check if waiting is done successfully

**done\_notify** ()

Check if notifying is done

**waiting\_state** ()

**notifying\_state** ()

## pylablib.core.thread.profile module

`pylablib.core.thread.profile.start(reset=True)`

Start yappi profile logging.

If `reset==True`, reset the stats.

`pylablib.core.thread.profile.reset()`

Reset yappi profiling stats

`pylablib.core.thread.profile.stop()`

Stop yappi profiling

`pylablib.core.thread.profile.get_stats()`

Get yappi profiling stats.

Return tuple ((`ttime`, `wtime`), (`threads`, `ctls`)). Here `ttime` and `wtime` are total execution time (sum of all thread times) and the wall time (since the last reset) respectively. `threads` are yappi-generated stats, and `ctls` is the list [(`name`, `ctl`)] with the controller names and thread controllers, which are ordered in the same way as `threads` (for any non-controlled or stopped thread these are set to `None`).

`pylablib.core.thread.profile.print_stats(nfunc=None, ntotfunc=None, min_func_frac=0.001)`

Print yappi profiling stats.

`nfunc` is the number of top (most expensive) functions to print per each thread, `ntotfunc` is the number of global top function to print; `None` for either means that they are not printed. `min_func_frac` specifies the minimal fraction of the total time for which the function stats are still printed (to prevent lost of printouts for “cheap” threads).

**pylablib.core.thread.synchronizing module**

**class** pylablib.core.thread.synchronizing.QThreadNotifier(*skippable=True*)

Bases: *pylablib.core.thread.notifier.ISkippableNotifier*

Wait-notify thread synchronizer for controlled Qt threads based on *notifier.ISkippableNotifier*.

Like *notifier.ISkippableNotifier*, the main functions are *ISkippableNotifier.wait()* (wait in a message loop until notified or until timeout expires) and *ISkippableNotifier.notify()* (notify the waiting thread). Both of these can only be called once and will raise an error on repeating calls. Along with notifying a variable can be passed, which can be accessed using *get\_value()* and *get\_value\_sync()*.

**Parameters** *skippable* (*bool*) – if True, allows for skippable wait events (if *ISkippableNotifier.notify()* is called before *ISkippableNotifier.wait()*, neither methods are actually called).

**get\_value()**

Get the value passed by the notifier (doesn't check if it has been passed already)

**get\_value\_sync** (*timeout=None*)

Wait (with the given *timeout*) for the value passed by the notifier

**done\_notify()**

Check if notifying is done

**done\_wait()**

Check if waiting is done

**notify** (*\*args, \*\*kwargs*)

Notify the waiting process.

Can only be called once per notifier lifetime. If the notifier allows skipping, and this method is called before *wait()*, return immediately.

**notifying\_state()**

**success\_wait()**

Check if waiting is done successfully

**wait** (*\*args, \*\*kwargs*)

Wait for the notification.

Can only be called once per notifier lifetime. If the notifier allows skipping, and this method is called after *notify()*, return immediately.

**waiting()**

Check if waiting is in progress

**waiting\_state()**

**class** pylablib.core.thread.synchronizing.QMultiThreadNotifier

Bases: *object*

Wait-notify thread synchronizer that can be used for multiple threads and called multiple times.

Performs similar function to conditional variables. The synchronizer has an internal counter which is increased by 1 every time it is notified. The wait functions have an option to wait until the counter reaches the specific counter value (usually, 1 above the last wait call).

**wait** (*state=1, timeout=None*)

Wait until notifier counter is equal to at least *state*

Return current counter state plus 1, which is the next smallest value resulting in waiting.

**wait\_until** (*condition*, *timeout=None*)

Wait until *condition* is met.

*condition* is a function which is called (in the waiting thread) every time the synchronizer is notified. If it return non-False, the waiting is complete and its result is returned.

**notify** ()

Notify all waiting threads

**fail** ()

Mark notifier as fails

Fails all waiting notifiers. All subsequent wait calls raise an error

**class** `pylablib.core.thread.synchronizing.QLockNotifier`

Bases: `object`

Resource lock.

Behaves similarly to the regular lock, but waiting is done in the message loop, which still allows interrupts.

**acquire** (*timeout=None*)

**release** ()

## pylablib.core.thread.threadprop module

**exception** `pylablib.core.thread.threadprop.ThreadError` (*msg=None*)

Bases: `RuntimeError`

Generic thread error

**args**

**with\_traceback** ()

Exception.with\_traceback(tb) – set self.\_\_traceback\_\_ to tb and return self.

**exception** `pylablib.core.thread.threadprop.NoControllerThreadError` (*msg=None*)

Bases: `pylablib.core.thread.threadprop.ThreadError`

Thread error for a case of thread having no controllers

**args**

**with\_traceback** ()

Exception.with\_traceback(tb) – set self.\_\_traceback\_\_ to tb and return self.

**exception** `pylablib.core.thread.threadprop.DuplicateControllerThreadError` (*msg=None*)

Bases: `pylablib.core.thread.threadprop.ThreadError`

Thread error for a case of a duplicate thread controller

**args**

**with\_traceback** ()

Exception.with\_traceback(tb) – set self.\_\_traceback\_\_ to tb and return self.

**exception** `pylablib.core.thread.threadprop.TimeoutThreadError` (*msg=None*)

Bases: `pylablib.core.thread.threadprop.ThreadError`, `TimeoutError`

Thread error for a case of a wait timeout

**args**

**characters\_written**

**errno**  
POSIX exception code

**filename**  
exception filename

**filename2**  
second exception filename

**strerror**  
exception strerror

**with\_traceback()**  
Exception.with\_traceback(tb) – set self.\_\_traceback\_\_ to tb and return self.

**exception** `pylablib.core.thread.threadprop.NoMessageThreadError(msg=None)`  
Bases: `pylablib.core.thread.threadprop.ThreadError`  
Thread error for a case of trying to get a non-existing message

**args**  
**with\_traceback()**  
Exception.with\_traceback(tb) – set self.\_\_traceback\_\_ to tb and return self.

**exception** `pylablib.core.thread.threadprop.SkippedCallError(msg=None)`  
Bases: `pylablib.core.thread.threadprop.ThreadError`  
Thread error for a case of external call getting skipped (unscheduled)

**args**  
**with\_traceback()**  
Exception.with\_traceback(tb) – set self.\_\_traceback\_\_ to tb and return self.

**exception** `pylablib.core.thread.threadprop.InterruptException(msg=None)`  
Bases: `Exception`  
Generic interrupt exception (raised by some function to signal interrupts from other threads)

**args**  
**with\_traceback()**  
Exception.with\_traceback(tb) – set self.\_\_traceback\_\_ to tb and return self.

**exception** `pylablib.core.thread.threadprop.InterruptExceptionStop(msg=None)`  
Bases: `pylablib.core.thread.threadprop.InterruptException`  
Interrupt exception denoting thread stop request

**args**  
**with\_traceback()**  
Exception.with\_traceback(tb) – set self.\_\_traceback\_\_ to tb and return self.

`pylablib.core.thread.threadprop.get_app()`  
Get current application instance

`pylablib.core.thread.threadprop.get_gui_thread()`  
Get main (GUI) thread, or None if application is not running

`pylablib.core.thread.threadprop.is_gui_running()`  
Check if GUI is running

`pylablib.core.thread.threadprop.is_gui_thread()`  
Check if the current thread is the one running the GUI loop

`pylablib.core.thread.threadprop.current_controller` (*require\_controller=True*)

Get controller of the current thread.

If the current thread has not controller *and* '*require\_controller==True*', raise an error; otherwise, return `None`.

## pylablib.core.thread.utils module

**class** `pylablib.core.thread.utils.ReadChangeLock`

Bases: `object`

Lock based on condition variables which handles a state which can be read or changed.

Any number of threads can read simultaneously, but changing is incompatible with other reading or changing.

**can\_read** ()

Check if the state can be read

**can\_change** ()

Check if the state can be changed

**reading** ()

Context manager denoting reading event

**changing** ()

Context manager denoting changing event

## Module contents

### pylablib.core.utils package

#### Submodules

### pylablib.core.utils.array\_utils module

`pylablib.core.utils.array_utils.as_array` (*data*, *force\_copy=False*, *try\_object=True*)

Turn *data* into a numpy array.

If *force\_copy==True*, copy the data if it's already a numpy array. If *try\_object==False*, only try to convert to numerical numpy arrays; otherwise, generic numpy arrays (with *dtype=="object"*) are acceptable.

`pylablib.core.utils.array_utils.get_shape` (*data*, *strict=False*)

Get the data shape.

If the data is a nested list and *strict==True*, raise an error unless all sublists have the same length (i.e., the data is rectangular).

### pylablib.core.utils.ctypes\_wrap module

`pylablib.core.utils.ctypes_wrap.get_value` (*rval*)

Get value of a ctypes variable

`pylablib.core.utils.ctypes_wrap.setup_func` (*func*, *argtypes*, *restype=None*, *er-  
rcheck=None*)

Setup a ctypes function.

Assign argtypes (list of argument types), restype (return value type) and errcheck (error checking function called for the return value).

```
class pylablib.core.utils.ctypes_wrap.CFunctionWrapper (restype=None,          er-
                                                         rcheck=None,          tu-
                                                         ple_single_retval=False,
                                                         return_res='auto',      de-
                                                         fault_rvals='rest',
                                                         pointer_byref=False)
```

Bases: `object`

Wrapper object for ctypes function.

The main methods are `wrap_annotated()` and `wrap_bare()`, which wrap a ctypes function and returns a Python function with a proper signature. These methods can also handle some standard use cases such as passing parameters by reference, or setting up the function arguments, or parsing the results. These methods can also be invoked when the wrapper is used as a callable; in this case, the exact method is determined by the presence of `.argtypes` attribute in the supplied function.

#### Parameters

- **restype** – default type of the function return value when calling `wrap_bare()` and `restype` is not supplied there explicitly (defaults to `ctypes.int`)
- **errcheck** – default error-checking function which is automatically called for the return value; can also be overridden explicitly when calling wrapping methods if `None`, no error checking method
- **tuple\_single\_retval** (*bool*) – determines if a single return values gets turned into a single-element tuple
- **return\_res** (*bool*) – determined if the function result gets returned; only used when list of return arguments (`rvals`) to wrapping functions is not explicitly supplied; can also be set to "auto" (default), which means that function returns its return value when no other `rvals` are found, and omits it otherwise.
- **default\_rvals** – default value for `rvals` in `wrap_annotated()` and `wrap_bare()`, if it is specified as `None` (default for those methods).
- **pointer\_byref** (*bool*) – if `True`, use explicit pointer creation instead of `byref` (in rare cases use of `byref` crashes the call).

**byref** (*value*)

**wrap\_bare** (*func*, *argtypes*, *argnames=None*, *restype=None*, *args='nonrval'*, *rvals='default'*, *argprep=None*, *rconv=None*, *byref='all'*, *errcheck=None*)

Annotate and wrap bare C function in a Python call.

Same as `wrap_annotated()`, but annotates the function first.

#### Parameters

- **func** – C function
- **argtypes** – list of ctypes types corresponding to function arguments; gets assigned as `func.argtypes`
- **argnames** – list of argument names; if not supplied, generated automatically as "arg1", "arg2", etc. Same for names which are defined as `None`.
- **restype** – type of the function return value; if `None`, use the value supplied to the wrapper constructor (defaults to `ctypes.int`)

- **args** – names of Python function arguments; can also be "all" (all C function arguments in that order), or "nonrval" (same, but with return value arguments excluded) by default, use "nonrval"
- **rvals** – names of return value arguments; can include either a C function argument name, or None (which means the function return value); can also be "rest" (lists all the arguments not included into args; if args=="nonrval", assume that there are no rvals), "pointer" (assume that all pointer arguments are rvals; this does not include `c_void_p`, `c_char_p`, or `c_wchar_p`); by default, use the value supplied on the wrapper creation ("rest" by default)
- **argprep** – dictionary {name: prep} of ways to prepare of C function arguments; each prep can be a value (which is assumed to be default argument value), or a callable, which is given values of Python function arguments
- **rconv** – dictionary {name: conv} of converters of the return values; each conv is a function which takes 3 arguments: unconverted ctypes value, dictionary of all C function arguments, and dictionary of all Python function arguments if conv takes less than 3 argument, then the arguments list is trimmed (e.g., if it takes only one argument, it will be an unconverted value) conv can also be "ctypes" (return raw ctypes value), or "raw" (return raw value for buffers).
- **byref** – list of all argument names which should be passed by reference; by default, it includes all arguments listed in rvals
- **errcheck** – error-checking function which is automatically called for the return value; if None, use the value supplied to the wrapper constructor (none by default)

**wrap\_annotated** (*func*, *args*='nonrval', *rvals*='default', *alias*=None, *argprep*=None, *rconv*=None, *byref*='all', *errcheck*=None)  
 Wrap annotated C function in a Python call.

Assumes that the functions has defined `.argtypes` (list of argument types) and `.argnames` (list of argument names) attributes.

#### Parameters

- **func** – C function
- **args** – names of Python function arguments; can also be "all" (all C function arguments in that order), or "nonrval" (same, but with return value arguments excluded); by default, use "nonrval"
- **rvals** – names of return value arguments; can include either a C function argument name, or None (which means the function return value); can also be "rest" (lists all the arguments not included into args; if args=="nonrval", assume that there are no rvals), "pointer" (assume that all pointer arguments are rvals; this does not include `c_void_p`, `c_char_p`, or `c_wchar_p`); by default, use the value supplied on the wrapper creation ("rest" by default)
- **alias** – either a list of argument names which replace `.argnames`, or a dictionary {argname: alias} which transforms names; all names in all other parameters (rvals, argprep, rconv, and byref) take aliased names
- **argprep** – dictionary {name: prep} of ways to prepare of C function arguments; each prep can be a value (which is assumed to be default argument value), or a callable, which is given values of Python function arguments
- **rconv** – dictionary {name: conv} of converters of the return values; each conv is a function which takes 3 arguments: unconverted ctypes value, dictionary of all C function arguments, and dictionary of all Python function arguments if conv takes

less than 3 argument, then the arguments list is trimmed (e.g., if it takes only one argument, it will be an unconverted value)

- **byref** – list of all argument names which should be passed by reference; by default, it includes all arguments listed in `rvals`
- **errcheck** – error-checking function which is automatically called for the return value; if `None`, use the value supplied to the wrapper constructor (none by default)

`pylablib.core.utils.ctypes_wrap.strprep(l, ctype=None, unicode=False)`

Make a string preparation function.

Return a function which creates a string with a fixed length of `l` bytes and returns a pointer to it. `ctype` can specify the type of the result (by default, `ctypes.c_char_p`).

`pylablib.core.utils.ctypes_wrap.buffprep(size_arg_pos, dtype)`

Make a buffer preparation function.

Return a function which creates a string with a variable size (specified by an argument at a position `size_arg_pos`). The buffer size is given in elements. `dtype` specifies the datatype of the buffer, whose size is used to determine buffer size in bytes.

`pylablib.core.utils.ctypes_wrap.buffconv(size_arg_pos, dtype)`

Make a buffer conversion function.

Return a function which converts a pointer of a variable size (specified by an argument at a position `size_arg_pos`) into a numpy array. The buffer size is given in elements. `dtype` specifies the datatype of the resulting array.

**class** `pylablib.core.utils.ctypes_wrap.CStructWrapper(struct=None)`

Bases: `object`

Wrapper around a ctypes structure, which allows for easier creation of parsing of these structures.

When created, all structure fields can be accessed/modified as attributes of the wrapper object. It can also be converted into tuple using `tup()` method, or back into C structure using `to_struct()` method.

Class variable `_struct` should be set to the ctypes structure which is being wrapped. Several other class variables determine the behavior when generating and parsing:

- `_prep`: dictionary {name: prep} of methods to prepare individual structure parameters; can be either a value or a function (which takes as ordered arguments all structure fields as ctypes values)
- `_conv`: dictionary {name: conv} of methods to convert individual structure parameters when parsing a C structure; can be either a function (which takes ctypes value of the field as a single argument) or a value; also can be used as a source of default values on wrapper creation
- `_tup`: dictionary {name: conv} of functions to convert structure values when generating a tuple
- `_tup_exc`: list of values to exclude from the resulting tuple
- `_tup_inc`: list of values to include in the resulting tuple (if `None`, include all)
- `_tup_add`: list of values to add to the resulting tuple (these values must then exist either as attributes, or as entries in `_tup` dictionary)
- `_tup_order`: order of fields in the returned tuple (by default, same as structure order)

Also specifies two overloaded methods for a more flexible preparation/conversion of structures. `conv()` takes no arguments and is called in the end of wrapper creation to finish setting up attributes. `prep()` takes a single argument (C structure) and is called when converting into a C structure to finish setting up the fields (e.g., size field).

**Parameters** `struct` – C structure to wrap (if `None`, create a new ‘blank’ structure).



**to\_struct()**  
Convert wrapper into a C structure

**prep(struct)**  
Prepare C structure after creation (by default, do nothing)

**conv()**  
Prepare wrapper after setting up the fields from the wrapped structure

**tup()**  
Convert wrapper into a named tuple

**classmethod prep\_struct(\*args, \*\*kwargs)**  
Prepare a blank C structure

**classmethod prep\_struct\_args(\*args, \*\*kwargs)**  
Prepare a C structure with the given supplied fields

**classmethod tup\_struct(struct, \*args, \*\*kwargs)**  
Convert C structure into a named tuple

`pylablib.core.utils.ctypes_wrap.class_tuple_to_dict(val, norm_strings=True, expand_lists=False)`

Convert a named tuple (usually, a tuple returned by `CStructWrapper.tup()`) into a dictionary.

Iterate recursively over all named tuple elements as well. If `norm_strings==True`, automatically translate byte strings into regular ones. If `expand_lists==True`, iterate recursively over lists members.

## pylablib.core.utils.dictionary module

Tree-like multi-level dictionary with advanced indexing options.

`pylablib.core.utils.dictionary.split_path(path, omit_empty=True, sep=None)`

Split generic path into individual path entries.

### Parameters

- **path** – Generic path. Lists and tuples (possible nested) are flattened; strings are split according to separators; non-strings are converted into strings first.
- **omit\_empty** (*bool*) – Determines if empty entries are skipped.
- **sep** (*str*) – If not `None`, defines regex for path separators; default separator is `'/'`.

**Returns** A list of individual entries.

**Return type** `list`

`pylablib.core.utils.dictionary.normalize_path_entry(entry,`

`case_normalization=None)`

Normalize the case of the entry if it's not case-sensitive. Normalization is either `None` (no normalization, names are case-sensitive), `'lower'` or `'upper'`

`pylablib.core.utils.dictionary.normalize_path(path, omit_empty=True, case_normalization=None, sep=None, force=False)`

Split and normalize generic path into individual path entries.

### Parameters

- **path** – Generic path. Lists and tuples (possible nested) are flattened; strings are split according to separators; non-strings are converted into strings first.
- **omit\_empty** (*bool*) – Determines if empty entries are skipped.

- **case\_normalization** (*str*) – Case normalization rules; can be `None` (no normalization, names are case-sensitive), `'lower'` or `'upper'`.
- **sep** (*str*) – If not `None`, defines regex for path separators; default separator is `'/'`.
- **force** (*bool*) – If `False`, treat lists as if they're already normalized.

**Returns** A list of individual normalized entries.

**Return type** `list`

`pylablib.core.utils.dictionary.is_dictionary(obj, generic=False)`

Determine if the object is a dictionary.

**Parameters**

- **obj** – object
- **generic** (*bool*) – if `False`, passes only `Dictionary` (or subclasses) objects; otherwise, passes any dictionary-like object.

**Returns** `bool`

`pylablib.core.utils.dictionary.as_dictionary(obj, case_normalization=None)`

Convert object into `Dictionary` with the given parameters.

If object is already a `Dictionary` (or its subclass), return unchanged, even if its parameters are different.

`pylablib.core.utils.dictionary.as_dict(obj, style='nested', copy=True)`

Convert object into standard `dict` with the given parameters.

If object is already a `dict`, return unchanged, even if the parameters are different.

```
class pylablib.core.utils.dictionary.Dictionary(root=None,  
                                              case_normalization=None,  
                                              copy=True)
```

Bases: `object`

Multi-level dictionary.

Access is done by path (all path elements are converted into strings and concatenated to form a single string path). If dictionary is not case-sensitive, all inserted and accessed paths are normalized to lower or upper case.

**Parameters**

- **root** (*dict* or `Dictionary`) – Initial value.
- **case\_normalization** (*str*) – Case normalization rules; can be `None` (no normalization, names are case-sensitive), `'lower'` or `'upper'`.
- **copy** (*bool*) – If `True`, make copy of the supplied data; otherwise, just make it the root.

**Warning:** If `copy==False`, the root data is already assumed to be normalized. If it isn't, the behavior might be incorrect.

**static** `is_dictionary(obj, generic=True)`

Determine if the object is a dictionary.

**Parameters**

- **obj** –

- **generic** (*bool*) – if *False*, passes only *Dictionary* (or subclasses) objects; otherwise, passes any dictionary-like object.

**Returns** *bool*

**static as\_dictionary** (*obj*, *case\_normalization=None*)

Convert object into *Dictionary* with the given parameters.

If object is already a *Dictionary* (or its subclass), return unchanged, even if its parameters are different.

**add\_entry** (*path*, *value*, *force=False*, *branch\_option='normalize'*)

Add value to a given path (overwrite leaf value if necessary).

Doesn't replace leaves with branches and vice-versa if *force==False*.

#### Parameters

- **path** –
- **value** –
- **force** (*bool*) – If *True*, change leaf into a branch and vice-versa; otherwise, raises *ValueError* if the conversion is necessary.
- **branch\_option** (*str*) –

**Decides what to do if the value is dictionary-like:**

- 'attach' – just attach the root,
- 'copy' – copy and attach,
- 'normalize' – copy while normalizing all the keys according to the current rules.

**get\_entry** (*path*, *as\_pointer=False*)

Get entry at a given path

#### Parameters

- **path** –
- **as\_pointer** (*bool*) – If *True* and entry is not a leaf, return *DictionaryPointer*; otherwise, return *Dictionary*

**has\_entry** (*path*, *kind='all'*)

Determine if the path is in the dictionary.

*kind* determines which kind of path to consider and can be 'leaf', 'branch' or 'all'.

**is\_leaf\_path** (*path*)

Determine if the path is in the dictionary and points to a leaf

**is\_branch\_path** (*path*)

Determine if the path is in the dictionary and points to a branch

**get\_max\_prefix** (*path*, *kind='all'*)

Find the longest prefix of *path* contained in the dictionary.

Return tuple (*prefix*, *rest*), where both path entries are normalized according to the dictionary rules (i.e., these are lists representing normalized paths). *kind* determines which kind of path to consider and can be 'leaf', 'branch' or 'all'. If the longest prefix is of a different kind, return (*None*, *None*).

**del\_entry** (*path*)

Delete entry from the dictionary.

Return `True` if the path was present. Note that it never raises *KeyError*.

**size** ()

Return the total size of the dictionary (number of nodes)

**get** (*path*, *default=None*)

Analog of `dict.get(): D.get(k, d) -> D[k] if k in D else d`

**pop** (*path*, *default=None*)

Analog of `dict.pop():` remove value at *path* and return it if *path* in *D*, otherwise return *default*

Note that it never raises *KeyError*.

**setdefault** (*path*, *default=None*)

Analog of `dict.setdefault(): D.setdefault(k, d) -> D.get(k, d)`, also sets `D[k]=d` if *k* not in *D*.

**items** (*ordered=False*, *leafs=False*, *path\_kind='split'*, *wrap\_branches=True*)

Analog of `dict.items()`, by default iterating only over the immediate children of the root.

#### Parameters

- **ordered** (*bool*) – If `True`, loop over keys in alphabetic order.
- **leafs** (*bool*) – If `True`, loop over leaf nodes (i.e., behave as ‘flat’ dictionary); otherwise, loop over immediate children (i.e., behave as ‘nested’ dictionary)
- **path\_kind** (*str*) – either `"split"` (each path is a tuple of individual keys), or `"joined"` (each path is a single string)
- **wrap\_branches** (*bool*) – if `True`, wrap sub-branches into *DictionaryPointer* objects; otherwise, return them as nested built-in dictionaries

**iteritems** (*ordered=False*, *leafs=False*, *path\_kind='split'*, *wrap\_branches=True*)

Analog of `dict.items()`, by default iterating only over the immediate children of the root.

#### Parameters

- **ordered** (*bool*) – If `True`, loop over keys in alphabetic order.
- **leafs** (*bool*) – If `True`, loop over leaf nodes (i.e., behave as ‘flat’ dictionary); otherwise, loop over immediate children (i.e., behave as ‘nested’ dictionary)
- **path\_kind** (*str*) – either `"split"` (each path is a tuple of individual keys), or `"joined"` (each path is a single string)
- **wrap\_branches** (*bool*) – if `True`, wrap sub-branches into *DictionaryPointer* objects; otherwise, return them as nested built-in dictionaries

**viewitems** (*ordered=False*, *leafs=False*, *path\_kind='split'*, *wrap\_branches=True*)

Analog of `dict.items()`, by default iterating only over the immediate children of the root.

#### Parameters

- **ordered** (*bool*) – If `True`, loop over keys in alphabetic order.
- **leafs** (*bool*) – If `True`, loop over leaf nodes (i.e., behave as ‘flat’ dictionary); otherwise, loop over immediate children (i.e., behave as ‘nested’ dictionary)

- **path\_kind** (*str*) – either "split" (each path is a tuple of individual keys), or "joined" (each path is a single string)
- **wrap\_branches** (*bool*) – if True, wrap sub-branches into *DictionaryPointer* objects; otherwise, return them as nested built-in dictionaries

**values** (*ordered=False, leafs=False, wrap\_branches=True*)

Analog of `dict.values()`, iterating only over the immediate children of the root.

#### Parameters

- **ordered** (*bool*) – If True, loop over keys in alphabetic order.
- **leafs** (*bool*) – If True, loop over leaf nodes (i.e., behave as 'flat' dictionary); otherwise, loop over immediate children (i.e., behave as 'nested' dictionary)
- **wrap\_branches** (*bool*) – if True, wrap sub-branches into *DictionaryPointer* objects; otherwise, return them as nested built-in dictionaries

**viewvalues** (*ordered=False, leafs=False, wrap\_branches=True*)

Analog of `dict.values()`, iterating only over the immediate children of the root.

#### Parameters

- **ordered** (*bool*) – If True, loop over keys in alphabetic order.
- **leafs** (*bool*) – If True, loop over leaf nodes (i.e., behave as 'flat' dictionary); otherwise, loop over immediate children (i.e., behave as 'nested' dictionary)
- **wrap\_branches** (*bool*) – if True, wrap sub-branches into *DictionaryPointer* objects; otherwise, return them as nested built-in dictionaries

**itervalues** (*ordered=False, leafs=False, wrap\_branches=True*)

Analog of `dict.values()`, iterating only over the immediate children of the root.

#### Parameters

- **ordered** (*bool*) – If True, loop over keys in alphabetic order.
- **leafs** (*bool*) – If True, loop over leaf nodes (i.e., behave as 'flat' dictionary); otherwise, loop over immediate children (i.e., behave as 'nested' dictionary)
- **wrap\_branches** (*bool*) – if True, wrap sub-branches into *DictionaryPointer* objects; otherwise, return them as nested built-in dictionaries

**keys** (*ordered=False, leafs=False, path\_kind='split'*)

Analog of `dict.keys()`, iterating only over the immediate children of the root.

#### Parameters

- **ordered** (*bool*) – If True, loop over keys in alphabetic order.
- **leafs** (*bool*) – If True, loop over leaf nodes (i.e., behave as 'flat' dictionary); otherwise, loop over immediate children (i.e., behave as 'nested' dictionary)
- **path\_kind** (*str*) – either "split" (each path is a tuple of individual keys), or "joined" (each path is a single string)

**viewkeys** (*ordered=False, leafs=False, path\_kind='split'*)

Analog of `dict.keys()`, iterating only over the immediate children of the root.

**Parameters**

- **ordered** (*bool*) – If `True`, loop over keys in alphabetic order.
- **leafs** (*bool*) – If `True`, loop over leaf nodes (i.e., behave as ‘flat’ dictionary); otherwise, loop over immediate children (i.e., behave as ‘nested’ dictionary)
- **path\_kind** (*str*) – either `"split"` (each path is a tuple of individual keys), or `"joined"` (each path is a single string)

**iterkeys** (*ordered=False, leafs=False, path\_kind='split'*)

Analog of `dict.keys()`, iterating only over the immediate children of the root.

**Parameters**

- **ordered** (*bool*) – If `True`, loop over keys in alphabetic order.
- **leafs** (*bool*) – If `True`, loop over leaf nodes (i.e., behave as ‘flat’ dictionary); otherwise, loop over immediate children (i.e., behave as ‘nested’ dictionary)
- **path\_kind** (*str*) – either `"split"` (each path is a tuple of individual keys), or `"joined"` (each path is a single string)

**paths** (*ordered=False, topdown=False, path\_kind='split'*)

Return list of all paths (leafs and nodes).

**Parameters**

- **ordered** (*bool*) – If `True`, loop over paths in alphabetic order.
- **topdown** (*bool*) – If `True`, return node’s leafs before its subtrees leafs.
- **path\_kind** (*str*) – either `"split"` (each path is a tuple of individual keys), or `"joined"` (each path is a single string)

**internodes** (*to\_visit='leafs', ordered=False, include\_path=False, topdown=False*)

Iterate over nodes.

**Parameters**

- **to\_visit** (*str*) – Can be `'leafs'`, `'branches'` or `'all'` and determines which parts of the dictionary are visited.
- **ordered** (*bool*) – If `True`, loop over paths in alphabetic order.
- **include\_path** (*bool*) – Include in the return value.
- **topdown** (*bool*) – If `True`, visit node and its leafs before its subtrees leafs.

**Yields** Values for leafs and *DictionaryPointer* for branches. If `include_path==True`, yields tuple (*path*, *value*), where *path* is in the form of a normalized list.

**nodes** (*to\_visit='leafs', ordered=False, include\_path=False, topdown=False*)

Iterate over nodes.

**Parameters**

- **to\_visit** (*str*) – Can be `'leafs'`, `'branches'` or `'all'` and determines which parts of the dictionary are visited.
- **ordered** (*bool*) – If `True`, loop over paths in alphabetic order.
- **include\_path** (*bool*) – Include in the return value.
- **topdown** (*bool*) – If `True`, visit node and its leafs before its subtrees leafs.

**Yields** Values for leafs and *DictionaryPointer* for branches. If `include_path==True`, yields tuple (*path*, *value*), where *path* is in the form of a normalized list.

**merge** (*source*, *path*=", *overwrite*=True, *normalize\_paths*=True)

Attach source (*dict* or other *Dictionary*) to a given branch; source is automatically deep-copied.

If *source* is not a dictionary, simply assign it (i.e., `D.merge(v, p)` is equivalent to `D.add_entry(p, v, force=True)` in this case). Compared to *add\_entry()*, merges two branches instead of removing the old branch completely.

#### Parameters

- **source** (*dict* or *Dictionary*) –
- **branch** (*tuple* or *str*) – Destination path.
- **overwrite** (*bool*) – If True, replaces the old entries with the new ones (it only matters for leaf assignments).
- **normalize\_paths** (*bool*) – If True and the dictionary isn't case sensitive, perform normalization if the *source*.

**update** (*source*, *path*=", *overwrite*=True, *normalize\_paths*=True)

Attach source (*dict* or other *Dictionary*) to a given branch; source is automatically deep-copied.

If *source* is not a dictionary, simply assign it (i.e., `D.merge(v, p)` is equivalent to `D.add_entry(p, v, force=True)` in this case). Compared to *add\_entry()*, merges two branches instead of removing the old branch completely.

#### Parameters

- **source** (*dict* or *Dictionary*) –
- **branch** (*tuple* or *str*) – Destination path.
- **overwrite** (*bool*) – If True, replaces the old entries with the new ones (it only matters for leaf assignments).
- **normalize\_paths** (*bool*) – If True and the dictionary isn't case sensitive, perform normalization if the *source*.

**detach** (*path*)

Remove a branch or a leaf from the current dictionary.

Branch is returned as a separate *Dictionary*. If *path* is missing, raise a *KeyError*.

**collect** (*paths*, *detach*=False, *ignore\_missing*=True)

Collect a set of subpaths into a separate dictionary.

#### Parameters

- **paths** – list or set of paths
- **detach** – if True, added branches are removed from this dictionary
- **ignore\_missing** – if True, ignore paths from the list which are not present in this dictionary; otherwise, raise a *KeyError*.

**branch\_copy** (*branch*=")

Get a copy of the branch as a *Dictionary*

**copy** ()

Get a full copy the dictionary

**updated** (*source*, *path=""*, *overwrite=True*, *normalize\_paths=True*)

Get a copy of the dictionary and attach a new branch to it.

Parameters are the same as in the `Dictionary.merge()`.

**as\_dict** (*style='nested'*, *copy=True*)

Convert into a `dict` object.

#### Parameters

- **style** (*str*) –

Determines style of the result:

- 'nested' – subtrees are turned into nested dictionaries,
- 'flat' – single dictionary is formed with full paths as keys.

- **copy** (*bool*) – If False and `style=='nested'`, return the root dictionary.

**asdict** (*style='nested'*, *copy=True*)

Convert into a `dict` object.

#### Parameters

- **style** (*str*) –

Determines style of the result:

- 'nested' – subtrees are turned into nested dictionaries,
- 'flat' – single dictionary is formed with full paths as keys.

- **copy** (*bool*) – If False and `style=='nested'`, return the root dictionary.

**as\_json** (*style='nested'*)

Convert into a JSON string.

**Parameters** **style** (*str*) – Determines style of the result: - 'nested' – subtrees are turned into nested dictionaries, - 'flat' – single dictionary is formed with full paths as keys.

**classmethod from\_json** (*data*, *case\_normalization=None*)

Convert JSON representations of a dictionary into a `Dictionary` object

**as\_pandas** (*index\_key=True*, *as\_series=True*)

Convert into a pandas DataFrame or Series object.

#### Parameters

- **index\_key** (*bool*) – If False, create a 2-column table with the first column ("key") containing string path and the second column ("value") containing value; otherwise, move key to the table index.
- **as\_series** (*bool*) – If `index_key==True` and `as_series==True`, convert the resulting DataFrame into 1D Series (the key is the index); otherwise, keep it as a single-column table

**get\_path** ()

**branch\_pointer** (*branch=""*)

Get a `DictionaryPointer` of a given branch

**map\_self** (*func*, *to\_visit='leaves'*, *pass\_path=False*, *topdown=False*, *branch\_option='normalize'*)

Apply *func* to the nodes in the dictionary.

Note that any pointers to the replaced branches or their sub-branches will become invalid.



**Parameters**

- **func** (*callable*) – Mapping function. Leafs are passed by value, branches (if visited) are passed as *DictionaryPointer*.
- **to\_visit** (*str*) – Can be 'leaves', 'branches' or 'all' and determines which parts of the dictionary passed to the map function.
- **pass\_path** (*bool*) – If True, pass the node path (in the form of a normalized list) as a first argument to *func*.
- **topdown** (*bool*) – If True, visit node and its leaves before its subtrees leaves.
- **branch\_option** (*str*) – If the function returns a dict-like object, determines how to incorporate into the dictionary; can be "normalize" (make a copy with normalized paths and insert that), "copy" (make a copy without normalization), or "attach" (simply replace the value without copying and normalization)

**filter\_self** (*pred*, *to\_visit*='leaves', *pass\_path*=False, *topdown*=False)

Remove all the nodes from the dictionary for which *pred* returns False.

**Parameters**

- **pred** (*callable*) – Filter function. Leafs are passed to *pred* by value, branches (if visited) are passed as *DictionaryPointer*.
- **to\_visit** (*str*) – Can be 'leaves', 'branches' or 'all' and determines which parts of the dictionary passed to the predicate.
- **pass\_path** (*bool*) – If True, pass the node path (in the form of a normalized list) as a first argument to *pred*.
- **topdown** (*bool*) – If True, visit node and its leaves before its subtrees leaves.

**diff** (*other*)

Perform an element-wise comparison to another Dictionary.

If the other Dictionary has a different case sensitivity, raise *ValueError*.

**Returns** *DictionaryDiff*

**static diff\_flatdict** (*first*, *second*)

Find the difference between flat *dict* objects.

**Returns** *DictionaryDiff*

**static find\_intersection** (*dicts*, *use\_flatten*=False)

Find intersection of multiple dictionaries.

**Parameters**

- **dicts** (*[Dictionary]*) –
- **use\_flatten** (*bool*) – If True flatten all dictionaries before comparison (works faster for a large number of dictionaries).

**Returns** *DictionaryIntersection*

**get\_matching\_paths** (*pattern*, *wildkey*='\*', *wildpath*='\*\*', *only\_leaves*=True)

Get all paths in the tree that match the provided pattern.

**Parameters**

- **wildkey** (*str*) – Pattern symbol that matches any key.
- **wildpath** (*str*) – Pattern symbol that matches any subpath (possibly empty).

- **only\_leaves** (*bool*) – If `True`, only check leaf paths; otherwise, check subtree paths (i.e., incomplete leaf paths) as well. Basically, `only_leaves=False` is analogous to adding `wildpath` at the end of the pattern.

**get\_matching\_subtree** (*pattern*, *wildkey*='\*', *wildpath*='\*\*', *only\_leaves*=*True*)

Get a subtree containing nodes with paths matching the provided pattern.

#### Parameters

- **wildkey** (*str*) – Pattern symbol that matches any key.
- **wildpath** (*str*) – Pattern symbol that matches any subpath (possibly empty).
- **only\_leaves** (*bool*) – If `True`, only check leaf paths; otherwise, check subtree paths (i.e., incomplete leaf paths) as well. Basically, `only_leaves=False` is analogous to adding `wildpath` at the end of the pattern.

**class** `pylablib.core.utils.dictionary.DictionaryDiff`

Bases: `pylablib.core.utils.dictionary.DictionaryDiff`

Describes a difference between the two dictionaries.

#### **same**

Contains the leafs which is the same.

Type `Dictionary`

#### **changed\_from**

Contains the leafs from the first dictionary which have different values in the second dictionary.

Type `Dictionary`

#### **changed\_to**

Contains the leafs from the second dictionary which have different values in the first dictionary.

Type `Dictionary`

#### **removed**

Contains the leafs from the first dictionary which are absent in the second dictionary.

Type `Dictionary`

#### **added**

Contains the leafs from the second dictionary which are absent in the first dictionary.

Type `Dictionary`

#### **added**

#### **changed\_from**

#### **changed\_to**

#### **count** ()

Return number of occurrences of value.

#### **index** ()

Return first index of value.

Raises `ValueError` if the value is not present.

#### **removed**

#### **same**

**class** `pylablib.core.utils.dictionary.DictionaryIntersection`  
 Bases: `pylablib.core.utils.dictionary.DictionaryIntersection`

Describes the result of finding intersection of multiple dictionaries.

**common**

Contains the intersection of all dictionaries.

Type `Dictionary`

**individual**

Contains list of difference from intersection for all dictionaries.

Type `[Dictionary]`

**common**

**count()**

Return number of occurrences of value.

**index()**

Return first index of value.

Raises `ValueError` if the value is not present.

**individual**

**class** `pylablib.core.utils.dictionary.DictionaryPointer` (*root=None, pointer=None, case\_normalization=None, copy=True*)

Bases: `pylablib.core.utils.dictionary.Dictionary`

Similar to `Dictionary`, but can point at one of the branches instead of the full dictionary.

Effect is mostly equivalent to prepending some path to all queries.

**Parameters**

- **root** (*dict or Dictionary*) – Complete tree.
- **pointer** – Path to the pointer location.
- **case\_normalization** (*str*) – Case normalization rules; can be `None` (no normalization, names are case-sensitive), `'lower'` or `'upper'`.
- **copy** (*bool*) – If `True`, make copy of the supplied data; otherwise, just make it the root.

**Warning:** If `copy==False`, the root data is already assumed to be normalized. If it isn't, the behavior might be incorrect.

**get\_path()**

Return pointer path in the whole dictionary.

**move\_to** (*path="", absolute=True*)

Move the pointer to a new path.

**Parameters**

- **path** –
- **absolute** (*bool*) – If `True`, path is specified with respect to the root; otherwise, it's specified with respect to the current position (and can only go deeper).

**move\_up** (*levels*, *strict=True*)

Move the pointer by the given number of levels up.

If *strict==True* and there are not enough levels above, raise an error. Otherwise, stop at the top dictionary level.

**branch\_pointer** (*branch=""*)

Get a *DictionaryPointer* of a given branch.

**add\_entry** (*path*, *value*, *force=False*, *branch\_option='normalize'*)

Add value to a given path (overwrite leaf value if necessary).

Doesn't replace leaves with branches and vice-versa if *force==False*.

#### Parameters

- **path** –
- **value** –
- **force** (*bool*) – If *True*, change leaf into a branch and vice-versa; otherwise, raises *ValueError* if the conversion is necessary.
- **branch\_option** (*str*) –

**Decides what to do if the value is dictionary-like:**

- 'attach' – just attach the root,
- 'copy' – copy and attach,
- 'normalize' – copy while normalizing all the keys according to the current rules.

**as\_dict** (*style='nested'*, *copy=True*)

Convert into a *dict* object.

#### Parameters

- **style** (*str*) –

**Determines style of the result:**

- 'nested' – subtrees are turned into nested dictionaries,
- 'flat' – single dictionary is formed with full paths as keys.

- **copy** (*bool*) – If *False* and *style=='nested'*, return the root dictionary.

**static as\_dictionary** (*obj*, *case\_normalization=None*)

Convert object into *Dictionary* with the given parameters.

If object is already a *Dictionary* (or its subclass), return unchanged, even if its parameters are different.

**as\_json** (*style='nested'*)

Convert into a JSON string.

**Parameters style** (*str*) – Determines style of the result: - 'nested' – subtrees are turned into nested dictionaries, - 'flat' – single dictionary is formed with full paths as keys.

**as\_pandas** (*index\_key=True*, *as\_series=True*)

Convert into a pandas DataFrame or Series object.

#### Parameters

- **index\_key** (*bool*) – If `False`, create a 2-column table with the first column ("key") containing string path and the second column ("value") containing value; otherwise, move key to the table index.
- **as\_series** (*bool*) – If `index_key==True` and `as_series==True`, convert the resulting `DataFrame` into 1D Series (the key is the index); otherwise, keep it as a single-column table

**asdict** (*style='nested', copy=True*)  
Convert into a `dict` object.

#### Parameters

- **style** (*str*) –

**Determines style of the result:**

- 'nested' – subtrees are turned into nested dictionaries,
- 'flat' – single dictionary is formed with full paths as keys.

- **copy** (*bool*) – If `False` and `style=='nested'`, return the root dictionary.

**branch\_copy** (*branch=""*)  
Get a copy of the branch as a `Dictionary`

**collect** (*paths, detach=False, ignore\_missing=True*)  
Collect a set of subpaths into a separate dictionary.

#### Parameters

- **paths** – list or set of paths
- **detach** – if `True`, added branches are removed from this dictionary
- **ignore\_missing** – if `True`, ignore paths from the list which are not present in this dictionary; otherwise, raise a `KeyError`.

**copy** ()  
Get a full copy the dictionary

**del\_entry** (*path*)  
Delete entry from the dictionary.

Return `True` if the path was present. Note that it never raises `KeyError`.

**detach** (*path*)  
Remove a branch or a leaf from the current dictionary.  
Branch is returned as a separate `Dictionary`. If *path* is missing, raise a `KeyError`.

**diff** (*other*)  
Perform an element-wise comparison to another `Dictionary`.  
If the other `Dictionary` has a different case sensitivity, raise `ValueError`.

**Returns** `DictionaryDiff`

**static diff\_flatdict** (*first, second*)  
Find the difference between flat `dict` objects.

**Returns** `DictionaryDiff`

**filter\_self** (*pred, to\_visit='leaves', pass\_path=False, topdown=False*)  
Remove all the nodes from the dictionary for which *pred* returns `False`.

#### Parameters

- **pred** (*callable*) – Filter function. Leafs are passed to *pred* by value, branches (if visited) are passed as *DictionaryPointer*.
- **to\_visit** (*str*) – Can be 'leafs', 'branches' or 'all' and determines which parts of the dictionary passed to the predicate.
- **pass\_path** (*bool*) – If True, pass the node path (in the form of a normalized list) as a first argument to *pred*.
- **topdown** (*bool*) – If True, visit node and its leafs before its subtrees leafs.

**static find\_intersection** (*dicts*, *use\_flatten=False*)

Find intersection of multiple dictionaries.

#### Parameters

- **dicts** (*[Dictionary]*) –
- **use\_flatten** (*bool*) – If True flatten all dictionaries before comparison (works faster for a large number of dictionaries).

Returns *DictionaryIntersection*

**classmethod from\_json** (*data*, *case\_normalization=None*)

Convert JSON representations of a dictionary into a *Dictionary* object

**get** (*path*, *default=None*)

Analog of `dict.get(): D.get(k, d) -> D[k] if k in D else d`

**get\_entry** (*path*, *as\_pointer=False*)

Get entry at a given path

#### Parameters

- **path** –
- **as\_pointer** (*bool*) – If True and entry is not a leaf, return *DictionaryPointer*; otherwise, return *Dictionary*

**get\_matching\_paths** (*pattern*, *wildkey='\**', *wildpath='\*\*'*, *only\_leaves=True*)

Get all paths in the tree that match the provided pattern.

#### Parameters

- **wildkey** (*str*) – Pattern symbol that matches any key.
- **wildpath** (*str*) – Pattern symbol that matches any subpath (possibly empty).
- **only\_leaves** (*bool*) – If True, only check leaf paths; otherwise, check subtree paths (i.e., incomplete leaf paths) as well. Basically, *only\_leaves=False* is analogous to adding *wildpath* at the end of the pattern.

**get\_matching\_subtree** (*pattern*, *wildkey='\**', *wildpath='\*\*'*, *only\_leaves=True*)

Get a subtree containing nodes with paths matching the provided pattern.

#### Parameters

- **wildkey** (*str*) – Pattern symbol that matches any key.
- **wildpath** (*str*) – Pattern symbol that matches any subpath (possibly empty).
- **only\_leaves** (*bool*) – If True, only check leaf paths; otherwise, check subtree paths (i.e., incomplete leaf paths) as well. Basically, *only\_leaves=False* is analogous to adding *wildpath* at the end of the pattern.

**get\_max\_prefix** (*path*, *kind*='all')

Find the longest prefix of *path* contained in the dictionary.

Return tuple (*prefix*, *rest*), where both path entries are normalized according to the dictionary rules (i.e., these are lists representing normalized paths). *kind* determines which kind of path to consider and can be 'leaf', 'branch' or 'all'. If the longest prefix is of a different kind, return (None, None).

**has\_entry** (*path*, *kind*='all')

Determine if the path is in the dictionary.

*kind* determines which kind of path to consider and can be 'leaf', 'branch' or 'all'.

**is\_branch\_path** (*path*)

Determine if the path is in the dictionary and points to a branch

**static is\_dictionary** (*obj*, *generic*=True)

Determine if the object is a dictionary.

#### Parameters

- **obj** –
- **generic** (*bool*) – if False, passes only *Dictionary* (or subclasses) objects; otherwise, passes any dictionary-like object.

**Returns** bool

**is\_leaf\_path** (*path*)

Determine if the path is in the dictionary and points to a leaf

**items** (*ordered*=False, *leafs*=False, *path\_kind*='split', *wrap\_branches*=True)

Analog of `dict.items()`, by default iterating only over the immediate children of the root.

#### Parameters

- **ordered** (*bool*) – If True, loop over keys in alphabetic order.
- **leafs** (*bool*) – If True, loop over leaf nodes (i.e., behave as ‘flat’ dictionary); otherwise, loop over immediate children (i.e., behave as ‘nested’ dictionary)
- **path\_kind** (*str*) – either "split" (each path is a tuple of individual keys), or "joined" (each path is a single string)
- **wrap\_branches** (*bool*) – if True, wrap sub-branches into *DictionaryPointer* objects; otherwise, return them as nested built-in dictionaries

**iteritems** (*ordered*=False, *leafs*=False, *path\_kind*='split', *wrap\_branches*=True)

Analog of `dict.items()`, by default iterating only over the immediate children of the root.

#### Parameters

- **ordered** (*bool*) – If True, loop over keys in alphabetic order.
- **leafs** (*bool*) – If True, loop over leaf nodes (i.e., behave as ‘flat’ dictionary); otherwise, loop over immediate children (i.e., behave as ‘nested’ dictionary)
- **path\_kind** (*str*) – either "split" (each path is a tuple of individual keys), or "joined" (each path is a single string)
- **wrap\_branches** (*bool*) – if True, wrap sub-branches into *DictionaryPointer* objects; otherwise, return them as nested built-in dictionaries

**iterkeys** (*ordered=False, leafs=False, path\_kind='split'*)

Analog of `dict.keys()`, iterating only over the immediate children of the root.

#### Parameters

- **ordered** (*bool*) – If `True`, loop over keys in alphabetic order.
- **leafs** (*bool*) – If `True`, loop over leaf nodes (i.e., behave as ‘flat’ dictionary); otherwise, loop over immediate children (i.e., behave as ‘nested’ dictionary)
- **path\_kind** (*str*) – either `"split"` (each path is a tuple of individual keys), or `"joined"` (each path is a single string)

**iternodes** (*to\_visit='leafs', ordered=False, include\_path=False, topdown=False*)

Iterate over nodes.

#### Parameters

- **to\_visit** (*str*) – Can be `'leafs'`, `'branches'` or `'all'` and determines which parts of the dictionary are visited.
- **ordered** (*bool*) – If `True`, loop over paths in alphabetic order.
- **include\_path** (*bool*) – Include in the return value.
- **topdown** (*bool*) – If `True`, visit node and its leafs before its subtrees leafs.

**Yields** Values for leafs and `DictionaryPointer` for branches. If `include_path==True`, yields tuple (`path`, `value`), where `path` is in the form of a normalized list.

**intervalues** (*ordered=False, leafs=False, wrap\_branches=True*)

Analog of `dict.values()`, iterating only over the immediate children of the root.

#### Parameters

- **ordered** (*bool*) – If `True`, loop over keys in alphabetic order.
- **leafs** (*bool*) – If `True`, loop over leaf nodes (i.e., behave as ‘flat’ dictionary); otherwise, loop over immediate children (i.e., behave as ‘nested’ dictionary)
- **wrap\_branches** (*bool*) – if `True`, wrap sub-branches into `DictionaryPointer` objects; otherwise, return them as nested built-in dictionaries

**keys** (*ordered=False, leafs=False, path\_kind='split'*)

Analog of `dict.keys()`, iterating only over the immediate children of the root.

#### Parameters

- **ordered** (*bool*) – If `True`, loop over keys in alphabetic order.
- **leafs** (*bool*) – If `True`, loop over leaf nodes (i.e., behave as ‘flat’ dictionary); otherwise, loop over immediate children (i.e., behave as ‘nested’ dictionary)
- **path\_kind** (*str*) – either `"split"` (each path is a tuple of individual keys), or `"joined"` (each path is a single string)

**map\_self** (*func, to\_visit='leafs', pass\_path=False, topdown=False, branch\_option='normalize'*)

Apply `func` to the nodes in the dictionary.

Note that any pointers to the replaced branches or their sub-branches will become invalid.

#### Parameters



- **func** (*callable*) – Mapping function. Leafs are passed by value, branches (if visited) are passed as *DictionaryPointer*.
- **to\_visit** (*str*) – Can be 'leafs', 'branches' or 'all' and determines which parts of the dictionary passed to the map function.
- **pass\_path** (*bool*) – If True, pass the node path (in the form of a normalized list) as a first argument to *func*.
- **topdown** (*bool*) – If True, visit node and its leafs before its subtrees leafs.
- **branch\_option** (*str*) – If the function returns a dict-like object, determines how to incorporate into the dictionary; can be "normalize" (make a copy with normalized paths and insert that), "copy" (make a copy without normalization), or "attach" (simply replace the value without copying and normalization)

**merge** (*source*, *path*='', *overwrite*=True, *normalize\_paths*=True)

Attach source (*dict* or other *Dictionary*) to a given branch; source is automatically deep-copied.

If *source* is not a dictionary, simply assign it (i.e., `D.merge(v, p)` is equivalent to `D.add_entry(p, v, force=True)` in this case). Compared to *add\_entry()*, merges two branches instead of removing the old branch completely.

#### Parameters

- **source** (*dict* or *Dictionary*) –
- **branch** (*tuple* or *str*) – Destination path.
- **overwrite** (*bool*) – If True, replaces the old entries with the new ones (it only matters for leaf assignments).
- **normalize\_paths** (*bool*) – If True and the dictionary isn't case sensitive, perform normalization if the *source*.

**nodes** (*to\_visit*='leafs', *ordered*=False, *include\_path*=False, *topdown*=False)

Iterate over nodes.

#### Parameters

- **to\_visit** (*str*) – Can be 'leafs', 'branches' or 'all' and determines which parts of the dictionary are visited.
- **ordered** (*bool*) – If True, loop over paths in alphabetic order.
- **include\_path** (*bool*) – Include in the return value.
- **topdown** (*bool*) – If True, visit node and its leafs before its subtrees leafs.

**Yields** Values for leafs and *DictionaryPointer* for branches. If *include\_path*==True, yields tuple (*path*, *value*), where *path* is in the form of a normalized list.

**paths** (*ordered*=False, *topdown*=False, *path\_kind*='split')

Return list of all paths (leafs and nodes).

#### Parameters

- **ordered** (*bool*) – If True, loop over paths in alphabetic order.
- **topdown** (*bool*) – If True, return node's leafs before its subtrees leafs.
- **path\_kind** (*str*) – either "split" (each path is a tuple of individual keys), or "joined" (each path is a single string)

**pop** (*path*, *default=None*)

Analog of `dict.pop()`: remove value at *path* and return it if *path* in *D*, otherwise return *default*

Note that it never raises `KeyError`.

**setdefault** (*path*, *default=None*)

Analog of `dict.setdefault()`: `D.setdefault(k,d) -> D.get(k,d)`, also sets `D[k]=d` if *k* not in *D*.

**size** ()

Return the total size of the dictionary (number of nodes)

**update** (*source*, *path="*, *overwrite=True*, *normalize\_paths=True*)

Attach *source* (`dict` or other *Dictionary*) to a given branch; *source* is automatically deep-copied.

If *source* is not a dictionary, simply assign it (i.e., `D.merge(v,p)` is equivalent to `D.add_entry(p,v,force=True)` in this case). Compared to `add_entry()`, merges two branches instead of removing the old branch completely.

#### Parameters

- **source** (*dict* or *Dictionary*) –
- **branch** (*tuple* or *str*) – Destination path.
- **overwrite** (*bool*) – If `True`, replaces the old entries with the new ones (it only matters for leaf assignments).
- **normalize\_paths** (*bool*) – If `True` and the dictionary isn't case sensitive, perform normalization if the *source*.

**updated** (*source*, *path="*, *overwrite=True*, *normalize\_paths=True*)

Get a copy of the dictionary and attach a new branch to it.

Parameters are the same as in the `Dictionary.merge()`.

**values** (*ordered=False*, *leafs=False*, *wrap\_branches=True*)

Analog of `dict.values()`, iterating only over the immediate children of the root.

#### Parameters

- **ordered** (*bool*) – If `True`, loop over keys in alphabetic order.
- **leafs** (*bool*) – If `True`, loop over leaf nodes (i.e., behave as 'flat' dictionary); otherwise, loop over immediate children (i.e., behave as 'nested' dictionary)
- **wrap\_branches** (*bool*) – if `True`, wrap sub-branches into *DictionaryPointer* objects; otherwise, return them as nested built-in dictionaries

**viewitems** (*ordered=False*, *leafs=False*, *path\_kind='split'*, *wrap\_branches=True*)

Analog of `dict.items()`, by default iterating only over the immediate children of the root.

#### Parameters

- **ordered** (*bool*) – If `True`, loop over keys in alphabetic order.
- **leafs** (*bool*) – If `True`, loop over leaf nodes (i.e., behave as 'flat' dictionary); otherwise, loop over immediate children (i.e., behave as 'nested' dictionary)
- **path\_kind** (*str*) – either `"split"` (each path is a tuple of individual keys), or `"joined"` (each path is a single string)

- **wrap\_branches** (*bool*) – if `True`, wrap sub-branches into `DictionaryPointer` objects; otherwise, return them as nested built-in dictionaries

**viewkeys** (*ordered=False, leafs=False, path\_kind='split'*)

Analog of `dict.keys()`, iterating only over the immediate children of the root.

#### Parameters

- **ordered** (*bool*) – If `True`, loop over keys in alphabetic order.
- **leafs** (*bool*) – If `True`, loop over leaf nodes (i.e., behave as ‘flat’ dictionary); otherwise, loop over immediate children (i.e., behave as ‘nested’ dictionary)
- **path\_kind** (*str*) – either `"split"` (each path is a tuple of individual keys), or `"joined"` (each path is a single string)

**viewvalues** (*ordered=False, leafs=False, wrap\_branches=True*)

Analog of `dict.values()`, iterating only over the immediate children of the root.

#### Parameters

- **ordered** (*bool*) – If `True`, loop over keys in alphabetic order.
- **leafs** (*bool*) – If `True`, loop over leaf nodes (i.e., behave as ‘flat’ dictionary); otherwise, loop over immediate children (i.e., behave as ‘nested’ dictionary)
- **wrap\_branches** (*bool*) – if `True`, wrap sub-branches into `DictionaryPointer` objects; otherwise, return them as nested built-in dictionaries

`pylablib.core.utils.dictionary.combine_dictionaries` (*dicts, func, select='all', pass\_missing=False*)

Combine several dictionaries element-wise (only for leafs) using a given function.

#### Parameters

- **dicts** (*list or tuple*) – list of dictionaries (`Dictionary` or `dict`) to be combined
- **func** (*callable*) – combination function. Takes a single argument, which is a list of elements to be combined.
- **select** (*str*) – determines which keys are selected for the resulting dictionary. Can be either `"all"` (only keep keys which are present in all the dictionaries), or `"any"` (keep keys which are present in at least one dictionary). Only keys that point to leafs count; if a key points to a non-leaf branch in some dictionary, it is considered absent from this dictionary.
- **pass\_missing** (*bool*) – if `select=="any"`, this parameter determines whether missing elements will be passed to `func` as `None`, or omitted entirely.

**class** `pylablib.core.utils.dictionary.PrefixTree` (*root=None, case\_normalization=None, wildcard='\*', matchcard='.', copy=True*)

Bases: `pylablib.core.utils.dictionary.Dictionary`

Expansion of a `Dictionary` designed to store data related to prefixes.

Each branch node can have a leaf with a name given by wildcard (`'*'` by default) or matchcard (`'.'` by default). Wildcard assumes that the branch node path is a prefix; matchcard assumes exact match. These leafs are inspected when specific prefix tree functions (`find_largest_prefix()` and `find_all_prefixes()`) are used.

**Parameters**

- **root** (*dict* or *Dictionary*) – Complete tree.
- **case\_normalization** (*str*) – Case normalization rules; can be `None` (no normalization, names are case-sensitive), `'lower'` or `'upper'`.
- **wildcard** (*str*) – Symbol for a wildcard entry.
- **matchcard** (*str*) – Symbol for a matchcard entry.
- **copy** (*bool*) – If `True`, make copy of the supplied data; otherwise, just make it the root.

**Warning:** If `copy==False`, the root data is already assumed to be normalized. If it isn't, the behavior might be incorrect.

**copy()**

Get a full copy the prefix tree

**find\_largest\_prefix** (*path*, *default=None*, *allow\_nomatch\_exact=True*, *return\_path=False*, *return\_subpath=False*)

Find the entry which is the largest prefix of a given path.

**Parameters**

- **path** –
- **default** – Default value if the path isn't found.
- **allow\_nomatch\_exact** (*bool*) – If `True`, just element with the given path can be returned; otherwise, only elements stored under wildcards and matchcards are considered.
- **return\_path** (*bool*) – If `True`, return path to the element (i.e., the largest prefix) instead of the element itself.
- **return\_subpath** (*bool*) – If `True`, return tuple with a second element being part of the *path* left after subtraction of the prefix.

**find\_all\_prefixes** (*path*, *allow\_nomatch\_exact=True*, *return\_path=True*, *return\_subpath=False*)

Find list of all the entries which are prefixes of a given path.

**Parameters**

- **path** –
- **default** – Default value if the path isn't found.
- **allow\_nomatch\_exact** (*bool*) – If `True`, just element with the given path can be returned; otherwise, only elements stored under wildcards and matchcards are considered.
- **return\_path** (*bool*) – If `True`, return path to the element (i.e., the largest prefix) instead of the element itself.
- **return\_subpath** (*bool*) – If `True`, return tuple with a second element being part of the *path* left after subtraction of the prefix.

**add\_entry** (*path*, *value*, *force=False*, *branch\_option='normalize'*)

Add value to a given path (overwrite leaf value if necessary).

Doesn't replace leaves with branches and vice-versa if `force==False`.

#### Parameters

- **path** –
- **value** –
- **force** (*bool*) – If `True`, change leaf into a branch and vice-versa; otherwise, raises `ValueError` if the conversion is necessary.
- **branch\_option** (*str*) –

Decides what to do if the value is dictionary-like:

- `'attach'` – just attach the root,
- `'copy'` – copy and attach,
- `'normalize'` – copy while normalizing all the keys according to the current rules.

**as\_dict** (*style='nested', copy=True*)  
Convert into a `dict` object.

#### Parameters

- **style** (*str*) –

Determines style of the result:

- `'nested'` – subtrees are turned into nested dictionaries,
- `'flat'` – single dictionary is formed with full paths as keys.

- **copy** (*bool*) – If `False` and `style=='nested'`, return the root dictionary.

**static as\_dictionary** (*obj, case\_normalization=None*)  
Convert object into `Dictionary` with the given parameters.

If object is already a `Dictionary` (or its subclass), return unchanged, even if its parameters are different.

**as\_json** (*style='nested'*)  
Convert into a JSON string.

**Parameters style** (*str*) – Determines style of the result: - `'nested'` – subtrees are turned into nested dictionaries, - `'flat'` – single dictionary is formed with full paths as keys.

**as\_pandas** (*index\_key=True, as\_series=True*)  
Convert into a pandas `DataFrame` or `Series` object.

#### Parameters

- **index\_key** (*bool*) – If `False`, create a 2-column table with the first column (`"key"`) containing string path and the second column (`"value"`) containing value; otherwise, move key to the table index.
- **as\_series** (*bool*) – If `index_key==True` and `as_series==True`, convert the resulting `DataFrame` into 1D `Series` (the key is the index); otherwise, keep it as a single-column table

**asdict** (*style='nested', copy=True*)  
Convert into a `dict` object.

#### Parameters

- **style** (*str*) –  
**Determines style of the result:**
  - 'nested' – subtrees are turned into nested dictionaries,
  - 'flat' – single dictionary is formed with full paths as keys.
- **copy** (*bool*) – If `False` and `style=='nested'`, return the root dictionary.

**branch\_copy** (*branch=""*)

Get a copy of the branch as a *Dictionary*

**branch\_pointer** (*branch=""*)

Get a *DictionaryPointer* of a given branch

**collect** (*paths, detach=False, ignore\_missing=True*)

Collect a set of subpaths into a separate dictionary.

#### Parameters

- **paths** – list or set of paths
- **detach** – if `True`, added branches are removed from this dictionary
- **ignore\_missing** – if `True`, ignore paths from the list which are not present in this dictionary; otherwise, raise a *KeyError*.

**del\_entry** (*path*)

Delete entry from the dictionary.

Return `True` if the path was present. Note that it never raises *KeyError*.

**detach** (*path*)

Remove a branch or a leaf from the current dictionary.

Branch is returned as a separate *Dictionary*. If *path* is missing, raise a *KeyError*.

**diff** (*other*)

Perform an element-wise comparison to another Dictionary.

If the other Dictionary has a different case sensitivity, raise *ValueError*.

**Returns** *DictionaryDiff*

**static diff\_flatdict** (*first, second*)

Find the difference between flat *dict* objects.

**Returns** *DictionaryDiff*

**filter\_self** (*pred, to\_visit='leafs', pass\_path=False, topdown=False*)

Remove all the nodes from the dictionary for which *pred* returns `False`.

#### Parameters

- **pred** (*callable*) – Filter function. Leafs are passed to *pred* by value, branches (if visited) are passed as *DictionaryPointer*.
- **to\_visit** (*str*) – Can be 'leafs', 'branches' or 'all' and determines which parts of the dictionary passed to the predicate.
- **pass\_path** (*bool*) – If `True`, pass the node path (in the form of a normalized list) as a first argument to *pred*.
- **topdown** (*bool*) – If `True`, visit node and its leafs before its subtrees leafs.

**static find\_intersection** (*dicts*, *use\_flatten=False*)

Find intersection of multiple dictionaries.

#### Parameters

- **dicts** (*Dictionary*) –
- **use\_flatten** (*bool*) – If `True` flatten all dictionaries before comparison (works faster for a large number of dictionaries).

**Returns** *DictionaryIntersection*

**classmethod from\_json** (*data*, *case\_normalization=None*)

Convert JSON representations of a dictionary into a *Dictionary* object

**get** (*path*, *default=None*)

Analog of `dict.get(): D.get(k, d) -> D[k] if k in D else d`

**get\_entry** (*path*, *as\_pointer=False*)

Get entry at a given path

#### Parameters

- **path** –
- **as\_pointer** (*bool*) – If `True` and entry is not a leaf, return *DictionaryPointer*; otherwise, return *Dictionary*

**get\_matching\_paths** (*pattern*, *wildkey='\**', *wildpath='\*\*'*, *only\_leaves=True*)

Get all paths in the tree that match the provided pattern.

#### Parameters

- **wildkey** (*str*) – Pattern symbol that matches any key.
- **wildpath** (*str*) – Pattern symbol that matches any subpath (possibly empty).
- **only\_leaves** (*bool*) – If `True`, only check leaf paths; otherwise, check subtree paths (i.e., incomplete leaf paths) as well. Basically, `only_leaves=False` is analogous to adding `wildpath` at the end of the pattern.

**get\_matching\_subtree** (*pattern*, *wildkey='\**', *wildpath='\*\*'*, *only\_leaves=True*)

Get a subtree containing nodes with paths matching the provided pattern.

#### Parameters

- **wildkey** (*str*) – Pattern symbol that matches any key.
- **wildpath** (*str*) – Pattern symbol that matches any subpath (possibly empty).
- **only\_leaves** (*bool*) – If `True`, only check leaf paths; otherwise, check subtree paths (i.e., incomplete leaf paths) as well. Basically, `only_leaves=False` is analogous to adding `wildpath` at the end of the pattern.

**get\_max\_prefix** (*path*, *kind='all'*)

Find the longest prefix of *path* contained in the dictionary.

Return tuple (*prefix*, *rest*), where both path entries are normalized according to the dictionary rules (i.e., these are lists representing normalized paths). *kind* determines which kind of path to consider and can be 'leaf', 'branch' or 'all'. If the longest prefix is of a different kind, return (*None*, *None*).

**get\_path** ()

**has\_entry** (*path*, *kind*='all')

Determine if the path is in the dictionary.

*kind* determines which kind of path to consider and can be 'leaf', 'branch' or 'all'.

**is\_branch\_path** (*path*)

Determine if the path is in the dictionary and points to a branch

**static is\_dictionary** (*obj*, *generic*=True)

Determine if the object is a dictionary.

#### Parameters

- **obj** –
- **generic** (*bool*) – if False, passes only *Dictionary* (or subclasses) objects; otherwise, passes any dictionary-like object.

**Returns** bool

**is\_leaf\_path** (*path*)

Determine if the path is in the dictionary and points to a leaf

**items** (*ordered*=False, *leafs*=False, *path\_kind*='split', *wrap\_branches*=True)

Analog of `dict.items()`, by default iterating only over the immediate children of the root.

#### Parameters

- **ordered** (*bool*) – If True, loop over keys in alphabetic order.
- **leafs** (*bool*) – If True, loop over leaf nodes (i.e., behave as 'flat' dictionary); otherwise, loop over immediate children (i.e., behave as 'nested' dictionary)
- **path\_kind** (*str*) – either "split" (each path is a tuple of individual keys), or "joined" (each path is a single string)
- **wrap\_branches** (*bool*) – if True, wrap sub-branches into *DictionaryPointer* objects; otherwise, return them as nested built-in dictionaries

**iteritems** (*ordered*=False, *leafs*=False, *path\_kind*='split', *wrap\_branches*=True)

Analog of `dict.items()`, by default iterating only over the immediate children of the root.

#### Parameters

- **ordered** (*bool*) – If True, loop over keys in alphabetic order.
- **leafs** (*bool*) – If True, loop over leaf nodes (i.e., behave as 'flat' dictionary); otherwise, loop over immediate children (i.e., behave as 'nested' dictionary)
- **path\_kind** (*str*) – either "split" (each path is a tuple of individual keys), or "joined" (each path is a single string)
- **wrap\_branches** (*bool*) – if True, wrap sub-branches into *DictionaryPointer* objects; otherwise, return them as nested built-in dictionaries

**iterkeys** (*ordered*=False, *leafs*=False, *path\_kind*='split')

Analog of `dict.keys()`, iterating only over the immediate children of the root.

#### Parameters

- **ordered** (*bool*) – If True, loop over keys in alphabetic order.
- **leafs** (*bool*) – If True, loop over leaf nodes (i.e., behave as 'flat' dictionary); otherwise, loop over immediate children (i.e., behave as 'nested' dictionary)



- **path\_kind** (*str*) – either "split" (each path is a tuple of individual keys), or "joined" (each path is a single string)

**internodes** (*to\_visit='leaves', ordered=False, include\_path=False, topdown=False*)

Iterate over nodes.

#### Parameters

- **to\_visit** (*str*) – Can be 'leaves', 'branches' or 'all' and determines which parts of the dictionary are visited.
- **ordered** (*bool*) – If True, loop over paths in alphabetic order.
- **include\_path** (*bool*) – Include in the return value.
- **topdown** (*bool*) – If True, visit node and its leaves before its subtrees leaves.

**Yields** Values for leaves and *DictionaryPointer* for branches. If `include_path==True`, yields tuple (*path*, *value*), where *path* is in the form of a normalized list.

**intervalues** (*ordered=False, leaves=False, wrap\_branches=True*)

Analog of `dict.values()`, iterating only over the immediate children of the root.

#### Parameters

- **ordered** (*bool*) – If True, loop over keys in alphabetic order.
- **leaves** (*bool*) – If True, loop over leaf nodes (i.e., behave as 'flat' dictionary); otherwise, loop over immediate children (i.e., behave as 'nested' dictionary)
- **wrap\_branches** (*bool*) – if True, wrap sub-branches into *DictionaryPointer* objects; otherwise, return them as nested built-in dictionaries

**keys** (*ordered=False, leaves=False, path\_kind='split'*)

Analog of `dict.keys()`, iterating only over the immediate children of the root.

#### Parameters

- **ordered** (*bool*) – If True, loop over keys in alphabetic order.
- **leaves** (*bool*) – If True, loop over leaf nodes (i.e., behave as 'flat' dictionary); otherwise, loop over immediate children (i.e., behave as 'nested' dictionary)
- **path\_kind** (*str*) – either "split" (each path is a tuple of individual keys), or "joined" (each path is a single string)

**map\_self** (*func, to\_visit='leaves', pass\_path=False, topdown=False, branch\_option='normalize'*)

Apply *func* to the nodes in the dictionary.

Note that any pointers to the replaced branches or their sub-branches will become invalid.

#### Parameters

- **func** (*callable*) – Mapping function. Leaves are passed by value, branches (if visited) are passed as *DictionaryPointer*.
- **to\_visit** (*str*) – Can be 'leaves', 'branches' or 'all' and determines which parts of the dictionary passed to the map function.
- **pass\_path** (*bool*) – If True, pass the node path (in the form of a normalized list) as a first argument to *func*.
- **topdown** (*bool*) – If True, visit node and its leaves before its subtrees leaves.

- **branch\_option** (*str*) – If the function returns a dict-like object, determines how to incorporate into the dictionary; can be "normalize" (make a copy with normalized paths and insert that), "copy" (make a copy without normalization), or "attach" (simply replace the value without copying and normalization)

**merge** (*source*, *path*=", *overwrite*=True, *normalize\_paths*=True)

Attach *source* (*dict* or other *Dictionary*) to a given branch; *source* is automatically deep-copied.

If *source* is not a dictionary, simply assign it (i.e., `D.merge(v, p)` is equivalent to `D.add_entry(p, v, force=True)` in this case). Compared to `add_entry()`, merges two branches instead of removing the old branch completely.

#### Parameters

- **source** (*dict* or *Dictionary*) –
- **branch** (*tuple* or *str*) – Destination path.
- **overwrite** (*bool*) – If True, replaces the old entries with the new ones (it only matters for leaf assignments).
- **normalize\_paths** (*bool*) – If True and the dictionary isn't case sensitive, perform normalization if the *source*.

**nodes** (*to\_visit*='leafs', *ordered*=False, *include\_path*=False, *topdown*=False)

Iterate over nodes.

#### Parameters

- **to\_visit** (*str*) – Can be 'leafs', 'branches' or 'all' and determines which parts of the dictionary are visited.
- **ordered** (*bool*) – If True, loop over paths in alphabetic order.
- **include\_path** (*bool*) – Include in the return value.
- **topdown** (*bool*) – If True, visit node and its leafs before its subtrees leafs.

**Yields** Values for leafs and *DictionaryPointer* for branches. If *include\_path*==True, yields tuple (*path*, *value*), where *path* is in the form of a normalized list.

**paths** (*ordered*=False, *topdown*=False, *path\_kind*='split')

Return list of all paths (leafs and nodes).

#### Parameters

- **ordered** (*bool*) – If True, loop over paths in alphabetic order.
- **topdown** (*bool*) – If True, return node's leafs before its subtrees leafs.
- **path\_kind** (*str*) – either "split" (each path is a tuple of individual keys), or "joined" (each path is a single string)

**pop** (*path*, *default*=None)

Analog of `dict.pop()`: remove value at *path* and return it if *path* in *D*, otherwise return *default*

Note that it never raises *KeyError*.

**setdefault** (*path*, *default*=None)

Analog of `dict.setdefault()`: `D.setdefault(k, d) -> D.get(k, d)`, also sets `D[k]=d` if *k* not in *D*.

**size** ()

Return the total size of the dictionary (number of nodes)

**update** (*source*, *path*="", *overwrite*=True, *normalize\_paths*=True)

Attach *source* (*dict* or other *Dictionary*) to a given branch; *source* is automatically deep-copied.

If *source* is not a dictionary, simply assign it (i.e., `D.merge(v, p)` is equivalent to `D.add_entry(p, v, force=True)` in this case). Compared to `add_entry()`, merges two branches instead of removing the old branch completely.

#### Parameters

- **source** (*dict* or *Dictionary*) –
- **branch** (*tuple* or *str*) – Destination path.
- **overwrite** (*bool*) – If True, replaces the old entries with the new ones (it only matters for leaf assignments).
- **normalize\_paths** (*bool*) – If True and the dictionary isn't case sensitive, perform normalization if the *source*.

**updated** (*source*, *path*="", *overwrite*=True, *normalize\_paths*=True)

Get a copy of the dictionary and attach a new branch to it.

Parameters are the same as in the `Dictionary.merge()`.

**values** (*ordered*=False, *leafs*=False, *wrap\_branches*=True)

Analog of `dict.values()`, iterating only over the immediate children of the root.

#### Parameters

- **ordered** (*bool*) – If True, loop over keys in alphabetic order.
- **leafs** (*bool*) – If True, loop over leaf nodes (i.e., behave as 'flat' dictionary); otherwise, loop over immediate children (i.e., behave as 'nested' dictionary)
- **wrap\_branches** (*bool*) – if True, wrap sub-branches into *DictionaryPointer* objects; otherwise, return them as nested built-in dictionaries

**viewitems** (*ordered*=False, *leafs*=False, *path\_kind*='split', *wrap\_branches*=True)

Analog of `dict.items()`, by default iterating only over the immediate children of the root.

#### Parameters

- **ordered** (*bool*) – If True, loop over keys in alphabetic order.
- **leafs** (*bool*) – If True, loop over leaf nodes (i.e., behave as 'flat' dictionary); otherwise, loop over immediate children (i.e., behave as 'nested' dictionary)
- **path\_kind** (*str*) – either "split" (each path is a tuple of individual keys), or "joined" (each path is a single string)
- **wrap\_branches** (*bool*) – if True, wrap sub-branches into *DictionaryPointer* objects; otherwise, return them as nested built-in dictionaries

**viewkeys** (*ordered*=False, *leafs*=False, *path\_kind*='split')

Analog of `dict.keys()`, iterating only over the immediate children of the root.

#### Parameters

- **ordered** (*bool*) – If True, loop over keys in alphabetic order.
- **leafs** (*bool*) – If True, loop over leaf nodes (i.e., behave as 'flat' dictionary); otherwise, loop over immediate children (i.e., behave as 'nested' dictionary)

- **path\_kind** (*str*) – either "split" (each path is a tuple of individual keys), or "joined" (each path is a single string)

**viewvalues** (*ordered=False, leafs=False, wrap\_branches=True*)

Analog of `dict.values()`, iterating only over the immediate children of the root.

#### Parameters

- **ordered** (*bool*) – If `True`, loop over keys in alphabetic order.
- **leafs** (*bool*) – If `True`, loop over leaf nodes (i.e., behave as 'flat' dictionary); otherwise, loop over immediate children (i.e., behave as 'nested' dictionary)
- **wrap\_branches** (*bool*) – if `True`, wrap sub-branches into `DictionaryPointer` objects; otherwise, return them as nested built-in dictionaries

```
class pylablib.core.utils.dictionary.FilterTree (root=None,  
                                              case_normalization=None, de-  
                                              fault=False, match_prefix=False,  
                                              copy=True)
```

Bases: `pylablib.core.utils.dictionary.Dictionary`

Expansion of a `Dictionary` designed to store hierarchical path filtering rules.

Store path templates and the corresponding values (usually `True` or `False` for a filter tree, but other values are possible). The `match()` method is then tested against this templates, and the value of the closest matching template (or default value, if none match) is returned. The templates can contain direct matches (e.g., "a/b/c/", which matches only "a/b/c/" path), "\*" path entries for a single level wildcard (e.g., "a/\*/c" matches "a/b/c" or "a/d/c", but not "a/c" or "a/b/d/c"), or "\*\*" path entries for a multi-level wildcard (e.g., "a/\*\*/c" matches "a/b/c", "a/c", or "a/b/d/c"). The paths are always tested first for direct match, then for "\*" match, then for "\*\*" match starting from the smallest subpath matching "\*\*".

#### Parameters

- **root** (*dict or Dictionary*) – A filter tree or a list of filter tree paths (which are all assumed to be have the `True` value)s
- **case\_normalization** (*str*) – Case normalization rules; can be `None` (no normalization, names are case-sensitive), 'lower' or 'upper'.
- **default** – Default value to return if no match is found.
- **match\_prefix** – if `True`, match the result even if only its prefix matches the tree content (same effect as adding "/\*" to every tree path)
- **copy** (*bool*) – If `True`, make copy of the supplied data; otherwise, just make it the root.

**Warning:** If `copy==False`, the root data is already assumed to be normalized. If it isn't, the behavior might be incorrect.

**copy** ()

Get a full copy the prefix tree

**match** (*path*)

Return the match result for the path

**add\_entry** (*path, value, force=False, branch\_option='normalize'*)

Add value to a given path (overwrite leaf value if necessary).

Doesn't replace leaves with branches and vice-versa if `force==False`.

#### Parameters

- **path** –
- **value** –
- **force** (*bool*) – If `True`, change leaf into a branch and vice-versa; otherwise, raises `ValueError` if the conversion is necessary.
- **branch\_option** (*str*) –

Decides what to do if the value is dictionary-like:

- 'attach' – just attach the root,
- 'copy' – copy and attach,
- 'normalize' – copy while normalizing all the keys according to the current rules.

**as\_dict** (*style='nested', copy=True*)  
Convert into a `dict` object.

#### Parameters

- **style** (*str*) –

Determines style of the result:

- 'nested' – subtrees are turned into nested dictionaries,
- 'flat' – single dictionary is formed with full paths as keys.

- **copy** (*bool*) – If `False` and `style=='nested'`, return the root dictionary.

**static as\_dictionary** (*obj, case\_normalization=None*)  
Convert object into `Dictionary` with the given parameters.

If object is already a `Dictionary` (or its subclass), return unchanged, even if its parameters are different.

**as\_json** (*style='nested'*)  
Convert into a JSON string.

**Parameters style** (*str*) – Determines style of the result: - 'nested' – subtrees are turned into nested dictionaries, - 'flat' – single dictionary is formed with full paths as keys.

**as\_pandas** (*index\_key=True, as\_series=True*)  
Convert into a pandas DataFrame or Series object.

#### Parameters

- **index\_key** (*bool*) – If `False`, create a 2-column table with the first column ("key") containing string path and the second column ("value") containing value; otherwise, move key to the table index.
- **as\_series** (*bool*) – If `index_key==True` and `as_series==True`, convert the resulting DataFrame into 1D Series (the key is the index); otherwise, keep it as a single-column table

**asdict** (*style='nested', copy=True*)  
Convert into a `dict` object.

#### Parameters

- **style** (*str*) –  
**Determines style of the result:**
  - 'nested' – subtrees are turned into nested dictionaries,
  - 'flat' – single dictionary is formed with full paths as keys.
- **copy** (*bool*) – If `False` and `style=='nested'`, return the root dictionary.

**branch\_copy** (*branch=""*)

Get a copy of the branch as a *Dictionary*

**branch\_pointer** (*branch=""*)

Get a *DictionaryPointer* of a given branch

**collect** (*paths, detach=False, ignore\_missing=True*)

Collect a set of subpaths into a separate dictionary.

#### Parameters

- **paths** – list or set of paths
- **detach** – if `True`, added branches are removed from this dictionary
- **ignore\_missing** – if `True`, ignore paths from the list which are not present in this dictionary; otherwise, raise a *KeyError*.

**del\_entry** (*path*)

Delete entry from the dictionary.

Return `True` if the path was present. Note that it never raises *KeyError*.

**detach** (*path*)

Remove a branch or a leaf from the current dictionary.

Branch is returned as a separate *Dictionary*. If *path* is missing, raise a *KeyError*.

**diff** (*other*)

Perform an element-wise comparison to another Dictionary.

If the other Dictionary has a different case sensitivity, raise *ValueError*.

**Returns** *DictionaryDiff*

**static diff\_flatdict** (*first, second*)

Find the difference between flat *dict* objects.

**Returns** *DictionaryDiff*

**filter\_self** (*pred, to\_visit='leafs', pass\_path=False, topdown=False*)

Remove all the nodes from the dictionary for which *pred* returns `False`.

#### Parameters

- **pred** (*callable*) – Filter function. Leafs are passed to *pred* by value, branches (if visited) are passed as *DictionaryPointer*.
- **to\_visit** (*str*) – Can be 'leafs', 'branches' or 'all' and determines which parts of the dictionary passed to the predicate.
- **pass\_path** (*bool*) – If `True`, pass the node path (in the form of a normalized list) as a first argument to *pred*.
- **topdown** (*bool*) – If `True`, visit node and its leafs before its subtrees leafs.

**static find\_intersection** (*dicts*, *use\_flatten=False*)

Find intersection of multiple dictionaries.

#### Parameters

- **dicts** (*Dictionary*) –
- **use\_flatten** (*bool*) – If `True` flatten all dictionaries before comparison (works faster for a large number of dictionaries).

**Returns** *DictionaryIntersection*

**classmethod from\_json** (*data*, *case\_normalization=None*)

Convert JSON representations of a dictionary into a *Dictionary* object

**get** (*path*, *default=None*)

Analog of `dict.get(): D.get(k, d) -> D[k] if k in D else d`

**get\_entry** (*path*, *as\_pointer=False*)

Get entry at a given path

#### Parameters

- **path** –
- **as\_pointer** (*bool*) – If `True` and entry is not a leaf, return *DictionaryPointer*; otherwise, return *Dictionary*

**get\_matching\_paths** (*pattern*, *wildkey='\**', *wildpath='\*\*'*, *only\_leaves=True*)

Get all paths in the tree that match the provided pattern.

#### Parameters

- **wildkey** (*str*) – Pattern symbol that matches any key.
- **wildpath** (*str*) – Pattern symbol that matches any subpath (possibly empty).
- **only\_leaves** (*bool*) – If `True`, only check leaf paths; otherwise, check subtree paths (i.e., incomplete leaf paths) as well. Basically, `only_leaves=False` is analogous to adding `wildpath` at the end of the pattern.

**get\_matching\_subtree** (*pattern*, *wildkey='\**', *wildpath='\*\*'*, *only\_leaves=True*)

Get a subtree containing nodes with paths matching the provided pattern.

#### Parameters

- **wildkey** (*str*) – Pattern symbol that matches any key.
- **wildpath** (*str*) – Pattern symbol that matches any subpath (possibly empty).
- **only\_leaves** (*bool*) – If `True`, only check leaf paths; otherwise, check subtree paths (i.e., incomplete leaf paths) as well. Basically, `only_leaves=False` is analogous to adding `wildpath` at the end of the pattern.

**get\_max\_prefix** (*path*, *kind='all'*)

Find the longest prefix of *path* contained in the dictionary.

Return tuple (*prefix*, *rest*), where both path entries are normalized according to the dictionary rules (i.e., these are lists representing normalized paths). *kind* determines which kind of path to consider and can be 'leaf', 'branch' or 'all'. If the longest prefix is of a different kind, return (*None*, *None*).

**get\_path** ()

**has\_entry** (*path*, *kind*='all')

Determine if the path is in the dictionary.

*kind* determines which kind of path to consider and can be 'leaf', 'branch' or 'all'.

**is\_branch\_path** (*path*)

Determine if the path is in the dictionary and points to a branch

**static is\_dictionary** (*obj*, *generic*=True)

Determine if the object is a dictionary.

**Parameters**

- **obj** –
- **generic** (*bool*) – if False, passes only *Dictionary* (or subclasses) objects; otherwise, passes any dictionary-like object.

**Returns** bool

**is\_leaf\_path** (*path*)

Determine if the path is in the dictionary and points to a leaf

**items** (*ordered*=False, *leafs*=False, *path\_kind*='split', *wrap\_branches*=True)

Analog of `dict.items()`, by default iterating only over the immediate children of the root.

**Parameters**

- **ordered** (*bool*) – If True, loop over keys in alphabetic order.
- **leafs** (*bool*) – If True, loop over leaf nodes (i.e., behave as 'flat' dictionary); otherwise, loop over immediate children (i.e., behave as 'nested' dictionary)
- **path\_kind** (*str*) – either "split" (each path is a tuple of individual keys), or "joined" (each path is a single string)
- **wrap\_branches** (*bool*) – if True, wrap sub-branches into *DictionaryPointer* objects; otherwise, return them as nested built-in dictionaries

**iteritems** (*ordered*=False, *leafs*=False, *path\_kind*='split', *wrap\_branches*=True)

Analog of `dict.items()`, by default iterating only over the immediate children of the root.

**Parameters**

- **ordered** (*bool*) – If True, loop over keys in alphabetic order.
- **leafs** (*bool*) – If True, loop over leaf nodes (i.e., behave as 'flat' dictionary); otherwise, loop over immediate children (i.e., behave as 'nested' dictionary)
- **path\_kind** (*str*) – either "split" (each path is a tuple of individual keys), or "joined" (each path is a single string)
- **wrap\_branches** (*bool*) – if True, wrap sub-branches into *DictionaryPointer* objects; otherwise, return them as nested built-in dictionaries

**iterkeys** (*ordered*=False, *leafs*=False, *path\_kind*='split')

Analog of `dict.keys()`, iterating only over the immediate children of the root.

**Parameters**

- **ordered** (*bool*) – If True, loop over keys in alphabetic order.
- **leafs** (*bool*) – If True, loop over leaf nodes (i.e., behave as 'flat' dictionary); otherwise, loop over immediate children (i.e., behave as 'nested' dictionary)



- **path\_kind** (*str*) – either "split" (each path is a tuple of individual keys), or "joined" (each path is a single string)

**internodes** (*to\_visit='leaves', ordered=False, include\_path=False, topdown=False*)

Iterate over nodes.

#### Parameters

- **to\_visit** (*str*) – Can be 'leaves', 'branches' or 'all' and determines which parts of the dictionary are visited.
- **ordered** (*bool*) – If True, loop over paths in alphabetic order.
- **include\_path** (*bool*) – Include in the return value.
- **topdown** (*bool*) – If True, visit node and its leaves before its subtrees leaves.

**Yields** Values for leaves and *DictionaryPointer* for branches. If `include_path==True`, yields tuple (*path*, *value*), where *path* is in the form of a normalized list.

**intervalues** (*ordered=False, leaves=False, wrap\_branches=True*)

Analog of `dict.values()`, iterating only over the immediate children of the root.

#### Parameters

- **ordered** (*bool*) – If True, loop over keys in alphabetic order.
- **leaves** (*bool*) – If True, loop over leaf nodes (i.e., behave as 'flat' dictionary); otherwise, loop over immediate children (i.e., behave as 'nested' dictionary)
- **wrap\_branches** (*bool*) – if True, wrap sub-branches into *DictionaryPointer* objects; otherwise, return them as nested built-in dictionaries

**keys** (*ordered=False, leaves=False, path\_kind='split'*)

Analog of `dict.keys()`, iterating only over the immediate children of the root.

#### Parameters

- **ordered** (*bool*) – If True, loop over keys in alphabetic order.
- **leaves** (*bool*) – If True, loop over leaf nodes (i.e., behave as 'flat' dictionary); otherwise, loop over immediate children (i.e., behave as 'nested' dictionary)
- **path\_kind** (*str*) – either "split" (each path is a tuple of individual keys), or "joined" (each path is a single string)

**map\_self** (*func, to\_visit='leaves', pass\_path=False, topdown=False, branch\_option='normalize'*)

Apply *func* to the nodes in the dictionary.

Note that any pointers to the replaced branches or their sub-branches will become invalid.

#### Parameters

- **func** (*callable*) – Mapping function. Leaves are passed by value, branches (if visited) are passed as *DictionaryPointer*.
- **to\_visit** (*str*) – Can be 'leaves', 'branches' or 'all' and determines which parts of the dictionary passed to the map function.
- **pass\_path** (*bool*) – If True, pass the node path (in the form of a normalized list) as a first argument to *func*.
- **topdown** (*bool*) – If True, visit node and its leaves before its subtrees leaves.

- **branch\_option** (*str*) – If the function returns a dict-like object, determines how to incorporate into the dictionary; can be "normalize" (make a copy with normalized paths and insert that), "copy" (make a copy without normalization), or "attach" (simply replace the value without copying and normalization)

**merge** (*source*, *path*=", *overwrite*=True, *normalize\_paths*=True)

Attach *source* (*dict* or other *Dictionary*) to a given branch; *source* is automatically deep-copied.

If *source* is not a dictionary, simply assign it (i.e., `D.merge(v, p)` is equivalent to `D.add_entry(p, v, force=True)` in this case). Compared to `add_entry()`, merges two branches instead of removing the old branch completely.

#### Parameters

- **source** (*dict* or *Dictionary*) –
- **branch** (*tuple* or *str*) – Destination path.
- **overwrite** (*bool*) – If True, replaces the old entries with the new ones (it only matters for leaf assignments).
- **normalize\_paths** (*bool*) – If True and the dictionary isn't case sensitive, perform normalization if the *source*.

**nodes** (*to\_visit*='leafs', *ordered*=False, *include\_path*=False, *topdown*=False)

Iterate over nodes.

#### Parameters

- **to\_visit** (*str*) – Can be 'leafs', 'branches' or 'all' and determines which parts of the dictionary are visited.
- **ordered** (*bool*) – If True, loop over paths in alphabetic order.
- **include\_path** (*bool*) – Include in the return value.
- **topdown** (*bool*) – If True, visit node and its leafs before its subtrees leafs.

**Yields** Values for leafs and *DictionaryPointer* for branches. If *include\_path*==True, yields tuple (*path*, *value*), where *path* is in the form of a normalized list.

**paths** (*ordered*=False, *topdown*=False, *path\_kind*='split')

Return list of all paths (leafs and nodes).

#### Parameters

- **ordered** (*bool*) – If True, loop over paths in alphabetic order.
- **topdown** (*bool*) – If True, return node's leafs before its subtrees leafs.
- **path\_kind** (*str*) – either "split" (each path is a tuple of individual keys), or "joined" (each path is a single string)

**pop** (*path*, *default*=None)

Analog of `dict.pop()`: remove value at *path* and return it if *path* in *D*, otherwise return *default*

Note that it never raises *KeyError*.

**setdefault** (*path*, *default*=None)

Analog of `dict.setdefault()`: `D.setdefault(k, d) -> D.get(k, d)`, also sets `D[k]=d` if *k* not in *D*.

**size** ()

Return the total size of the dictionary (number of nodes)

**update** (*source*, *path*=", *overwrite*=True, *normalize\_paths*=True)

Attach *source* (*dict* or other *Dictionary*) to a given branch; *source* is automatically deep-copied.

If *source* is not a dictionary, simply assign it (i.e., `D.merge(v, p)` is equivalent to `D.add_entry(p, v, force=True)` in this case). Compared to `add_entry()`, merges two branches instead of removing the old branch completely.

#### Parameters

- **source** (*dict* or *Dictionary*) –
- **branch** (*tuple* or *str*) – Destination path.
- **overwrite** (*bool*) – If True, replaces the old entries with the new ones (it only matters for leaf assignments).
- **normalize\_paths** (*bool*) – If True and the dictionary isn't case sensitive, perform normalization if the *source*.

**updated** (*source*, *path*=", *overwrite*=True, *normalize\_paths*=True)

Get a copy of the dictionary and attach a new branch to it.

Parameters are the same as in the *Dictionary.merge()*.

**values** (*ordered*=False, *leafs*=False, *wrap\_branches*=True)

Analog of `dict.values()`, iterating only over the immediate children of the root.

#### Parameters

- **ordered** (*bool*) – If True, loop over keys in alphabetic order.
- **leafs** (*bool*) – If True, loop over leaf nodes (i.e., behave as 'flat' dictionary); otherwise, loop over immediate children (i.e., behave as 'nested' dictionary)
- **wrap\_branches** (*bool*) – if True, wrap sub-branches into *DictionaryPointer* objects; otherwise, return them as nested built-in dictionaries

**viewitems** (*ordered*=False, *leafs*=False, *path\_kind*='split', *wrap\_branches*=True)

Analog of `dict.items()`, by default iterating only over the immediate children of the root.

#### Parameters

- **ordered** (*bool*) – If True, loop over keys in alphabetic order.
- **leafs** (*bool*) – If True, loop over leaf nodes (i.e., behave as 'flat' dictionary); otherwise, loop over immediate children (i.e., behave as 'nested' dictionary)
- **path\_kind** (*str*) – either "split" (each path is a tuple of individual keys), or "joined" (each path is a single string)
- **wrap\_branches** (*bool*) – if True, wrap sub-branches into *DictionaryPointer* objects; otherwise, return them as nested built-in dictionaries

**viewkeys** (*ordered*=False, *leafs*=False, *path\_kind*='split')

Analog of `dict.keys()`, iterating only over the immediate children of the root.

#### Parameters

- **ordered** (*bool*) – If True, loop over keys in alphabetic order.
- **leafs** (*bool*) – If True, loop over leaf nodes (i.e., behave as 'flat' dictionary); otherwise, loop over immediate children (i.e., behave as 'nested' dictionary)

- **path\_kind** (*str*) – either "split" (each path is a tuple of individual keys), or "joined" (each path is a single string)

**viewvalues** (*ordered=False, leafs=False, wrap\_branches=True*)

Analog of `dict.values()`, iterating only over the immediate children of the root.

#### Parameters

- **ordered** (*bool*) – If `True`, loop over keys in alphabetic order.
- **leafs** (*bool*) – If `True`, loop over leaf nodes (i.e., behave as 'flat' dictionary); otherwise, loop over immediate children (i.e., behave as 'nested' dictionary)
- **wrap\_branches** (*bool*) – if `True`, wrap sub-branches into `DictionaryPointer` objects; otherwise, return them as nested built-in dictionaries

**class** `pylablib.core.utils.dictionary.PrefixShortcutTree` (*shortcuts=None*)

Bases: `object`

Convenient storage for dictionary path shortcuts.

**Parameters** **shortcuts** (*dict*) – Dictionary of shortcuts {`shortcut`: `full_path`}.

**copy** ()

Return full copy

**add\_shortcut** (*source, dest, exact=False*)

Add a single shortcut.

#### Parameters

- **source** – Shortcut path.
- **dest** – expanded path corresponding to the shortcut.
- **exact** (*bool*) – If `True`, the shortcut works only for the exact path; otherwise, it works for any path with 'source' as a prefix.

**add\_shortcuts** (*shortcuts, exact=False*)

Add a dictionary of shortcuts {`shortcut`: `full_path`}.

Arguments are the same as in `PrefixShortcutTree.add_shortcut()`.

**remove\_shortcut** (*source*)

Remove a shortcut from the tree

**updated** (*shortcuts, exact=False*)

Make a copy and add additional shortcuts.

Arguments are the same as in `PrefixShortcutTree.add_shortcuts()`.

**class** `pylablib.core.utils.dictionary.DictionaryNode` (*\*\*vargs*)

Bases: `object`

`pylablib.core.utils.dictionary.dict_to_object_local` (*data, name=None, object\_generator=<function \_default\_object\_generator>*)

**class** `pylablib.core.utils.dictionary.ItemAccessor` (*getter=None, setter=None, deleter=None, contains\_checker='auto', normalize\_names=True, path\_separator=None, missing\_error=None*)

Bases: `object`

Simple wrapper which implements array interface using supplied methods.

Also has an option to normalize requested paths (enabled by default)

#### Parameters

- **getter** – method for getting values (`None` means none is supplied, so getting raises an error)
- **setter** – method for setting values (`None` means none is supplied, so setting raises an error)
- **deleter** – method for deleting values (`None` means none is supplied, so deleting raises an error)
- **contains\_checker** – method for checking if variable is present (`None` means none is supplied, so checking containment raises an error; "auto" means that getter raising `KeyError` is used for checking)
- **normalize\_names** – if `True`, normalize a supplied path using the standard `Dictionary` rules and join it into a single string using the supplied separator
- **path\_separator** – path separator regex used for splitting and joining the supplied paths (by default, the standard `" / "` separator)
- **missing\_error** – if not `None`, specifies the error raised on the missing value; used in `__contains__`, `get()` and `setdefault()` to determine if the value is missing

**get** (*name*, *default=None*)

**setdefault** (*name*, *default=None*)

### pylablib.core.utils.files module

Utilities for working with the file system: creating/removing/listing folders, comparing folders and files, working with zip archives.

`pylablib.core.utils.files.eof` (*f*, *strict=False*)  
Standard EOF function.

Return `True` if the the marker is at the end of the file. If `strict==True`, only return `True` if the marker is exactly at the end of file; otherwise, return `True` if it's at the end of further.

`pylablib.core.utils.files.get_file_creation_time` (*path*, *timestamp=True*)  
Try to find a file creation time. Return current time if an error occurs.

If `timestamp==True`, return UNIX timestamp; otherwise, return `datetime.datetime`.

`pylablib.core.utils.files.get_file_modification_time` (*path*, *timestamp=True*)  
Try to find a file modification time. Return current time if an error occurs.

If `timestamp==True`, return UNIX timestamp; otherwise, return `datetime.datetime`

`pylablib.core.utils.files.touch` (*fname*, *times=None*)  
Update file access and modification times.

**Parameters** *times* (*tuple*) – Access and modification times; if *times* is `None`, use current time.

`pylablib.core.utils.files.generate_indexed_filename` (*name\_format*, *idx\_start=0*, *folder=""*)

Generate an unused indexed filename in *folder*.

The name has *name\_format* (using standard Python `format()` rules, e.g., `"data_{:03d}.dat"`), and the index starts with *idx\_start*.

```
pylablib.core.utils.files.generate_prefixed_filename (prefix="",          suffix="",
                                                    idx_start=None,  idx_fmt='d',
                                                    folder=None)
```

Generate an unused filename with the given *prefix* and *suffix* in the given *folder*.

By default, the format is `prefix_{:d}_suffix`, where the parameter is the index starting with *idx\_start*. If *idx\_start* is `None`, first check simply `prefix+suffix` name before using numbered indices.

```
pylablib.core.utils.files.generate_temp_filename (prefix='__tmp__',      idx_start=0,
                                                  idx_template='d', folder="")
```

Generate a temporary filename with a given prefix.

*idx\_template* is the number index format (only the parameter itself, not the whole string).

```
pylablib.core.utils.files.fullsplit (path, ignore_empty=True)
```

Split path into a list.

If `ignore_empty==True`, exclude empty folder names.

```
pylablib.core.utils.files.normalize_path (p)
```

Normalize filesystem path (case and origin). If two paths are identical, they should be equal when normalized

```
pylablib.core.utils.files.case_sensitive_path ()
```

Check if OS path names are case-sensitive (e.g., Linux)

```
pylablib.core.utils.files.paths_equal (a, b)
```

Determine if the two paths are equal (can be local or have different case)

```
pylablib.core.utils.files.relative_path (a, b, check_paths=True)
```

Determine return path *a* as seen from *b*.

If `check_paths==True`, check if *a* is contained in *b* and raise the `OSError` if it isn't.

```
pylablib.core.utils.files.is_path_valid (p)
```

Check if the string is a valid path.

Not guaranteed to have complete success rate, but catches most likely errors (invalid characters, reserved file names, too long, etc.) Does not check if the path actually exists or if it can be written into.

```
class pylablib.core.utils.files.TempFile (folder="",          name=None,      mode='w',
                                          wait_time=None, rep_time=None)
```

Bases: `object`

Temporary file context manager.

Upon creation, generate an unused temporary filename. Upon entry, create the file using supplied mode and return self. Upon exit, close and remove the file.

Can be mostly substituted by `tempfile.TemporaryFile()`, but generates file locally, and with specified/determined name. Preserved largely for legacy reasons.

#### Parameters

- **folder** (*str*) – Containing folder.
- **name** (*str*) – File name. If `None`, generate new temporary name.
- **mode** (*str*) – File opening mode.
- **wait\_time** (*float*) – Waiting time between attempts to create the file if the first try fails.
- **rep\_time** (*int*) – Number of attempts to create the file if the first try fails.

**f**  
File object.

**name**  
File name.

**Type** `str`

**full\_name**  
File name including containing folder.

**Type** `str`

`pylablib.core.utils.files.copy_file(source, dest, overwrite=True, cmp_on_overwrite=True, preserve_metadata=True)`

Copy file, creating a containing folder if necessary. Return `True` if the operation was performed.

#### Parameters

- **overwrite** (`bool`) – If `True`, overwrite existing file.
- **cmp\_on\_overwrite** (`bool`) – If `True` and the two files are compared to be the same, don't perform overwrite.
- **preserve\_metadata** (`bool`) – If `True`, preserve file metadata (such as creation time) by using `shutil.copy2()`; otherwise, use `shutil.copy()`

`pylablib.core.utils.files.move_file(source, dest, overwrite=True, cmp_on_overwrite=True, preserve_if_not_move=False)`

Move file, creating a containing folder if necessary. Returns `True` if the operation was performed.

#### Parameters

- **overwrite** (`bool`) – If `True`, overwrite existing file (if the existing file isn't overwritten, preserve the original).
- **cmp\_on\_overwrite** (`bool`) – If `True` and the two files are compared to be the same, don't perform overwrite.
- **preserve\_if\_not\_move** (`bool`) – If `True` and the files are identical, preserve the original.

`pylablib.core.utils.files.ensure_dir_singlelevel(path, error_on_file=True)`

`pylablib.core.utils.files.ensure_dir(path, error_on_file=True)`

Ensure that the folder exists (create a new one if necessary).

If `error_on_file==True`, raise `OSError` if there's a file with the same name.

`pylablib.core.utils.files.remove_dir(path, error_on_file=True)`

Remove the folder recursively if it exists.

If `error_on_file==True`, raise `OSError` if there's a file with the same name.

`pylablib.core.utils.files.remove_dir_if_empty(path, error_on_file=True)`

Remove the folder only if it's empty.

If `error_on_file==True`, raise `OSError` if there's a file with the same name.

`pylablib.core.utils.files.clean_dir(path, error_on_file=True)`

Remove the folder and then recreate it.

If `error_on_file==True`, raise `OSError` if there's a file with the same name.

**class** `pylablib.core.utils.files.FolderList`

Bases: `pylablib.core.utils.files.FolderList`

Describes folder content

**count()**

Return number of occurrences of value.

**files**

**folders**

**index()**

Return first index of value.

Raises ValueError if the value is not present.

`pylablib.core.utils.files.list_dir(folder="", folder_filter=None, file_filter=None, separate_kinds=True, error_on_file=True)`

Return folder content filtered by *folder\_filter* and *file\_filter*.

#### Parameters

- **folder** (*str*) – Path to the folder.
- **folder\_filter** – Folder filter function (more description at `string.get_string_filter()`).
- **file\_filter** – File filter function (more description at `string.get_string_filter()`).
- **separate\_kinds** (*bool*) – if True, return `FolderList` with files and folder separate; otherwise, return a single list (works much faster).
- **error\_on\_file** (*bool*) – if True, raise `OSError` if there's a file with the same name as the target folder.

`pylablib.core.utils.files.dir_empty(folder, folder_filter=None, file_filter=None, level='single', error_on_file=True)`

Check if the folder is empty (only checks content filtered by *folder\_filter* and *file\_filter*).

#### Parameters

- **folder** (*str*) – Path to the folder.
- **folder\_filter** – Folder filter function (more description at `string.get_string_filter()`).
- **file\_filter** – File filter function (more description at `string.get_string_filter()`).
- **level** (*str*) – if 'single', check only immediate folder content; if 'recursive', follow recursively in all folders passing *folder\_filter*.
- **error\_on\_file** (*bool*) – if True, raise `OSError` if there's a file with the same name as the target folder.

`pylablib.core.utils.files.walk_dir(folder, folder_filter=None, file_filter=None, rel_path=True, topdown=True, visit_folder_filter=None, max_depth=None)`

Modification of `os.walk()` function.

Acts in a similar way, but *followlinks* is always False and errors of `os.listdir()` are always passed.

#### Parameters

- **folder** (*str*) – Path to the folder.



- **folder\_filter** – Folder filter function (more description at `string.get_string_filter()`).
- **file\_filter** – File filter function (more description at `string.get_string_filter()`).
- **rel\_path** (*bool*) – If True, the returned folder path is specified relative to the initial path.
- **topdown** (*bool*) – If True, return folder before its subfolders.
- **visit\_folder\_filter** – Filter for visiting folders (more description at `string.get_string_filter()`). If not None, specifies filter for visiting folders which is different from `folder_filter` (filter for returned folders).
- **max\_depth** (*int*) – If not None, limits the recursion depth.

#### Yields

For each folder (including the original) yields a tuple (**folder\_path**, **folders**, **files**), where *folder\_path* is the containing folder name and *folders* and *files* are its content (similar to `list_dir()`).

```
pylablib.core.utils.files.list_dir_recursive(folder, folder_filter=None,
                                             file_filter=None, topdown=True,
                                             visit_folder_filter=None,
                                             max_depth=None)
```

Recursive walk analog of `list_dir()`.

Parameters are the same as `walk_dir()`.

#### Returns *FolderList*

```
pylablib.core.utils.files.copy_dir(source, dest, folder_filter=None, file_filter=None,
                                   overwrite=True, cmp_on_overwrite=True, pre-
                                   serve_metadata=True)
```

Copy files satisfying the filtering conditions.

#### Parameters

- **source** (*str*) – Source path.
- **dest** (*str*) – Destination path.
- **folder\_filter** – Folder filter function (more description at `string.get_string_filter()`).
- **file\_filter** – File filter function (more description at `string.get_string_filter()`).
- **overwrite** (*bool*) – If True, overwrite existing files.
- **cmp\_on\_overwrite** (*bool*) – If True and the two files are compared to be the same, don't perform overwrite.
- **preserve\_metadata** (*bool*) – If True, preserve file metadata (such as creation time) by using `shutil.copy2()`; otherwise, use `shutil.copy()`

```
pylablib.core.utils.files.move_dir(source, dest, folder_filter=None, file_filter=None,
                                   overwrite=True, cmp_on_overwrite=True, pre-
                                   serve_if_not_move=False)
```

Move files satisfying the filtering conditions.

#### Parameters

- **source** (*str*) – Source path.
- **dest** (*str*) – Destination path.
- **folder\_filter** – Folder filter function (more description at [string.get\\_string\\_filter\(\)](#)).
- **file\_filter** – File filter function (more description at [string.get\\_string\\_filter\(\)](#)).
- **overwrite** (*bool*) – If True, overwrite existing files (if the existing file isn't overwritten, preserve the original).
- **cmp\_on\_overwrite** (*bool*) – If True and the two files are compared to be the same, don't perform overwrite.
- **preserve\_if\_not\_move** (*bool*) – If True and the files are identical, preserve the original.

`pylablib.core.utils.files.combine_diff(d1, d2)`

`pylablib.core.utils.files.cmp_dirs(a, b, folder_filter=None, file_filter=None, shallow=True, return_difference=False)`

Compare the folders based on the content filtered by *folder\_filter* and *file\_filter*.

#### Parameters

- **a** (*str*) – First folder path
- **b** (*str*) – Second folder path
- **folder\_filter** – Folder filter function (more description at [string.get\\_string\\_filter\(\)](#)).
- **file\_filter** – File filter function (more description at [string.get\\_string\\_filter\(\)](#)).
- **shallow** – If True, do shallow comparison of the files (see [filecmp.cmp\(\)](#)).
- **return\_difference** – If False, simply return *bool*; otherwise, return difference type ('=', '+', '-' or '\*').

`pylablib.core.utils.files.retry_copy(source, dest, overwrite=True, cmp_on_overwrite=True, preserve_metadata=True, try_times=5, delay=0.3)`

Retrying version of [copy\\_file\(\)](#).

If the operation raises error, wait for *delay* (in seconds) and call it again. Try total of *try\_times* times.

`pylablib.core.utils.files.retry_move(source, dest, overwrite=True, cmp_on_overwrite=True, preserve_if_not_move=False, try_times=5, delay=0.3)`

Retrying version of [move\\_file\(\)](#) (see [retry\\_copy\(\)](#) for details on retrying).

`pylablib.core.utils.files.retry_remove(path, try_times=5, delay=0.3)`

Retrying version of `os.remove()` (see [retry\\_copy\(\)](#) for details on retrying).

`pylablib.core.utils.files.retry_ensure_dir(path, error_on_file=True, try_times=5, delay=0.3)`

Retrying version of [ensure\\_dir\(\)](#) (see [retry\\_copy\(\)](#) for details on retrying).

`pylablib.core.utils.files.retry_copy_dir(source, dest, folder_filter=None, file_filter=None, overwrite=True, cmp_on_overwrite=True, preserve_metadata=True, try_times=5, delay=0.3)`

Retrying version of [copy\\_dir\(\)](#) (see [retry\\_copy\(\)](#) for details on retrying).

```
pylablib.core.utils.files.retry_move_dir(source, dest, folder_filter=None, file_filter=None,
                                         overwrite=True, cmp_on_overwrite=True,
                                         preserve_if_not_move=False, try_times=5,
                                         delay=0.3)
```

Retrying version of `move_dir()` (see `retry_copy()` for details on retrying).

```
pylablib.core.utils.files.retry_remove_dir(path, error_on_file=True, try_times=5, delay=0.3)
```

Retrying version of `remove_dir()` (see `retry_copy()` for details on retrying).

```
pylablib.core.utils.files.retry_remove_dir_if_empty(path, error_on_file=True,
                                                    try_times=5, delay=0.3)
```

Retrying version of `remove_dir_if_empty()` (see `retry_copy()` for details on retrying).

```
pylablib.core.utils.files.retry_clean_dir(path, error_on_file=True, try_times=5, delay=0.3)
```

Retrying version of `clean_dir()` (see `retry_copy()` for details on retrying).

```
pylablib.core.utils.files.zip_folder(zip_path, source_path, inside_path="",
                                     folder_filter=None, file_filter=None, mode='a',
                                     compression=8, compresslevel=None)
```

Add a folder into a zip archive.

#### Parameters

- **zip\_path** (*str*) – Path to the .zip file.
- **source\_path** (*str*) – Path to the source folder.
- **inside\_path** (*str*) – Destination path inside the zip archive.
- **folder\_filter** – Folder filter function (more description at `string.get_string_filter()`).
- **file\_filter** – File filter function (more description at `string.get_string_filter()`).
- **mode** (*str*) – Zip archive adding mode (see `zipfile.ZipFile`).
- **compression** – Zip archive compression (see `zipfile.ZipFile`).
- **compresslevel** – Zip archive compression level (see `zipfile.ZipFile`); ignored for Python version below 3.7.

```
pylablib.core.utils.files.zip_file(zip_path, source_path, inside_name=None, mode='a',
                                   compression=8, compresslevel=None)
```

Add a file into a zip archive.

#### Parameters

- **zip\_path** (*str*) – Path to the .zip file.
- **source\_path** (*str*) – Path to the source file.
- **inside\_name** (*str*) – Destination file name inside the zip archive (source name on the top level by default).
- **mode** (*str*) – Zip archive adding mode (see `zipfile.ZipFile`).
- **compression** – Zip archive compression (see `zipfile.ZipFile`).
- **compresslevel** – Zip archive compression level (see `zipfile.ZipFile`); ignored for Python version below 3.7.

```
pylablib.core.utils.files.zip_multiple_files(zip_path, source_paths, in-  
side_names=None, mode='a', com-  
pression=8, compresslevel=None)
```

Add a multiple files into a zip archive.

#### Parameters

- **zip\_path** (*str*) – Path to the .zip file.
- **source\_paths** (*[str]*) – List of path to the source files.
- **inside\_names** (*[str]* or *None*) – List of destination file names inside the zip archive (source name on the top level by default).
- **mode** (*str*) – Zip archive adding mode (see `zipfile.ZipFile`).
- **compression** – Zip archive compression (see `zipfile.ZipFile`).
- **compresslevel** – Zip archive compression level (see `zipfile.ZipFile`); ignored for Python version below 3.7.

```
pylablib.core.utils.files.unzip_folder(zip_path, dest_path, inside_path="",  
folder_filter=None, file_filter=None)
```

Extract a folder from a zip archive (create containing folder if necessary).

#### Parameters

- **zip\_path** (*str*) – Path to the .zip file.
- **dest\_path** (*str*) – Path to the destination folder.
- **inside\_path** (*str*) – Source path inside the zip archive; extracted data paths are relative (i.e., they don't include *inside\_path*).
- **folder\_filter** – Folder filter function (more description at `string.get_string_filter()`).
- **file\_filter** – File filter function (more description at `string.get_string_filter()`).

```
pylablib.core.utils.files.unzip_file(zip_path, dest_path, inside_path)
```

Extract a file from a zip archive (create containing folder if necessary).

#### Parameters

- **zip\_path** (*str*) – Path to the .zip file.
- **dest\_path** (*str*) – Destination file path.
- **inside\_path** (*str*) – Source path inside the zip archive.

## pylablib.core.utils.funcargparse module

Contains routines for checking arguments passed into a function for better flexibility.

```
pylablib.core.utils.funcargparse.parameter_value_error(par_val, par_name,  
message=None, er-  
ror_type=None)
```

Raise parameter value error (`ValueError` by default).

```
pylablib.core.utils.funcargparse.parameter_range_error(par_val,      par_name,
                                                         par_set=None,    mes-
                                                         sage=None,      er-
                                                         ror_type=None)
```

Raise parameter range error (`ValueError` by default).

```
pylablib.core.utils.funcargparse.check_parameter_range(par_val,      par_name,
                                                         par_set,    message=None,
                                                         error_type=None)
```

Raise error if *par\_val* is not in in the *par\_set* (*par\_name* is used in the error message).

```
pylablib.core.utils.funcargparse.getdefault(value,      default_value,    unas-
                                              signed_value=None,    con-
                                              flict_action='ignore',    message=None,
                                              error_type=None)
```

Analog of dict's `getdefault`.

If *value* is *unassigned\_value*, return *default\_value* instead. If *conflict\_action*=='error' and *value* != *default\_value*, raise value error using *message* and *error\_type*.

```
pylablib.core.utils.funcargparse.is_sequence(value, sequence_type='builtin;nostring')
```

Check if *value* is a sequence.

***sequence\_type* is semicolon separated list of possible sequence types:**

- 'builtin' - list, tuple or str
- 'nostring' - str is not allows
- 'array' - list, tuple or `numpy.ndarray`
- 'indexable' - anything which can be indexed
- 'haslength' - anything with length property

```
pylablib.core.utils.funcargparse.make_sequence(element,      length=1,      se-
                                              quence_type='list')
```

Turn element into a sequence of *sequence\_type* ('list' or 'tuple') repeated *length* times.

```
pylablib.core.utils.funcargparse.as_sequence(value,      multiply_length=1,    al-
                                              lowed_type='builtin;nostring',
                                              wrapping_type='list',
                                              length_conflict_action='ignore',    mes-
                                              sage=None, error_type=None)
```

Ensure that *value* is a sequence.

If *value* is not a sequence of *allowed\_type* (as checked by `is_sequence()`), turn it into a sequence specified by *wrapping\_type* and *multiply\_length*.

If *value* is a sequence and *length\_conflict\_action*=='error', raise error with *error\_type* and *error\_message* if the length doesn't match *multiply\_length*. Otherwise, return value unchanged.

## pylablib.core.utils.functions module

Utilities for dealing with function, methods and function signatures.

```
class pylablib.core.utils.functions.FunctionSignature (arg_names=None,  
                                                    defaults=None,  
                                                    varg_name=None,  
                                                    kwarg_name=None,  
                                                    kwonly_arg_names=None,  
                                                    cls=None,      obj=None,  
                                                    name=None, doc=None)
```

Bases: `object`

Description of a function signature, including name, argument names, default values, names of varg and kwarg arguments, class and object (for methods) and docstring.

#### Parameters

- **arg\_names** (*list*) – Names of the arguments.
- **default** (*dict*) – Dictionary {name: value} of default values.
- **varg\_name** (*str*) – Name of \*varg parameter (None means no such parameter).
- **kwarg\_name** (*str*) – Name of \*\*kwarg parameter (None means no such parameter).
- **cls** – Caller class, for methods.
- **obj** – Caller object, for methods.
- **name** (*str*) – Function name.
- **doc** (*str*) – Function docstring.

**get\_defaults\_list** ()

Get list of default values for arguments in the order specified in the signature.

**signature** (*pass\_order=None*)

Get string containing a signature (arguments list) of the function (call or definition), including \*vargs and \*\*kwargs.

If *pass\_order* is not None, it specifies the order in which the arguments are passed.

**wrap\_function** (*func, pass\_order=None*)

Wrap a function *func* into a containing function with this signature.

Sets function name, argument names, default values, object and class (for methods) and docstring. If *pass\_order* is not None, it determines the order in which the positional arguments are passed to the wrapped function.

**as\_kwargs** (*args, kwargs, add\_defaults=False, exclude=None*)

Turn *args* and *kwargs* into a single *kwargs* dictionary using the names of positional arguments.

If *add\_defaults*==True, add all the non-specified default arguments as well. If the function takes \**args* argument and some of the supplied arguments go there, place them into a list under "\*" key in the result. If *exclude* is not None is specifies arguments which should be excluded.

**arg\_value** (*argname, args=None, kwargs=None*)

Get the value of the argument with the given name for given *args* and *kwargs*

**mandatory\_args\_num** ()

Get minimal number of arguments which have to be passed to the function.

The mandatory arguments are the ones which are not bound to caller object (i.e., not `self`) and don't have default values.

**max\_args\_num** (*include\_positional=True, include\_keywords=True*)

Get maximal number of arguments which can be passed to the function.

### Parameters

- **include\_positional** (*bool*) – If True and function accepts \*vargs, return None (unlimited number of arguments).
- **include\_keywords** (*bool*) – If True and function accepts \*\*kwargs, return None (unlimited number of arguments).

**static from\_function** (*func*, *follow\_wrapped=True*)

Get signature of the given function or method.

If *follow\_wrapped==True*, follow `__wrapped__` attributes until the innermost function (useful for getting signatures of functions wrapped using `functools` methods).

**copy** ()

Return a copy

**as\_simple\_func** ()

Turn the signature into a simple function (as opposed to a bound method).

If the signature corresponds to a bound method, get rid of the first argument in the signature (`self`) and the bound object. Otherwise, return unchanged.

**static merge** (*inner*, *outer*, *add\_place='front'*, *merge\_duplicates=True*, *overwrite=None*, *hide\_outer\_obj=False*)

Merge two signatures (used for wrapping functions).

The signature describes the function would take arguments according to the *outer* signature and pass them according to the *inner* signature.

#### The arguments are combined:

- if *add\_place=='front'*, the outer arguments are placed in the beginning, followed by inner arguments not already listed;
- if *add\_place=='back'*, the inner arguments are placed in the beginning, followed by outer arguments not already listed.

The default values are joined, with the outer values superseding the inner values.

*overwrite* is a set or a list specifying which inner parameters are overwritten by the outer. It includes 'name', 'doc', 'cls', 'obj', 'varg\_name' and 'kwarg\_name'; the default value is all parameters.

If the inner signature is a bound method and *hide\_inner\_obj==True*, treat it as a function (with `self` argument missing). In this case, the wrapped signature `.obj` field will be None.

### Returns

(signature, pass\_order)

*pass\_order* is the order in which the arguments of the combined signature may be passed to the inner signature; it may be different from the signature order if *add\_place=='front'*. If *merge\_duplicates==True*, duplicate entries in *pass\_order* are omitted; otherwise, they're repeated.

### Return type tuple

`pylablib.core.utils.functions.funcsig` (*func*, *follow\_wrapped=True*)

Return a function signature object

`pylablib.core.utils.functions.getargsfrom` (*source*, *\*\*merge\_params*)

Decorator factory.

Returns decorator that conforms function signature to the source function. `**merge_params` are passed to the `FunctionSignature.merge()` method merging wrapped and source signature.

The default behavior (conforming parameter names, default values args and kwargs names) is useful for wrapping universal functions like `g(*args, **kwargs)`.

Example:

```
def f(x, y=2):
    return x+y

@getargsfrom(f)
def g(*args): # Now g has the same signature as f, including parameter names and
    ↪ default values.
    return prod(args)
```

`pylablib.core.utils.functions.call_cut_args(func, *args, **kwargs)`

Call `func` with the given arguments, omitting the ones that don't fit its signature.

`pylablib.core.utils.functions.getattr_call(obj, attr_name, *args, **kwargs)`

Call the getter for the attribute `attr_name` of `obj`.

If the attribute is a property, pass `*args` and `**kwargs` to the getter (`fget`); otherwise, ignore them.

`pylablib.core.utils.functions.setattr_call(obj, attr_name, *args, **kwargs)`

Call the setter for the attribute `attr_name` of `obj`.

If the attribute is a property, pass `*args` and `**kwargs` to the setter (`fset`); otherwise, the set value is assumed to be either the first argument, or the keyword argument with the name `'value'`.

`pylablib.core.utils.functions.delattr_call(obj, attr_name, *args, **kwargs)`

Call the deleter for the attribute `attr_name` of `obj`.

If the attribute is a property, pass `*args` and `**kwargs` to the deleter (`fdel`); otherwise, ignore them.

**class** `pylablib.core.utils.functions.IObjectCall`

Bases: `object`

Universal interface for object method call (makes methods, attributes and properties look like methods).

Should be called with an object as a first argument.

**class** `pylablib.core.utils.functions.MethodObjectCall(method)`

Bases: `pylablib.core.utils.functions.IObjectCall`

Object call created from an object method.

**Parameters** `method` – Either a method object or a method name which is used for the call.

**class** `pylablib.core.utils.functions.AttrObjectCall(name, as_getter)`

Bases: `pylablib.core.utils.functions.IObjectCall`

Object call created from an object attribute (makes attributes and properties look like methods).

**Parameters**

- **name** (`str`) – Attribute name.
- **as\_getter** (`bool`) – If `True`, call the getter when invoked; otherwise, call the setter.

If an attribute is a simple attribute, than getter gets no arguments and setter gets one argument (either the first argument, or the keyword argument named `'value'`). If it's a property, pass all the parameters to the property call.



```
class pylablib.core.utils.functions.IObjectProperty
```

Bases: `object`

Universal interface for an object property (makes methods, attributes and properties look like properties).

Can be used to get, set or remove a property.

```
get (obj, params=None)
```

```
set (obj, value)
```

```
rem (obj, params=None)
```

```
class pylablib.core.utils.functions.MethodObjectProperty (getter=None,           set-  
                                                         ter=None,           re-  
                                                         mover=None,        ex-  
                                                         pand_tuple=True)
```

Bases: `pylablib.core.utils.functions.IObjectProperty`

Object property created from object methods (makes methods look like properties).

#### Parameters

- **getter** (*callable*) – Method invoked on `get()`. If `None`, raise `RuntimeError` when called.
- **setter** (*callable*) – Method invoked on `set()`. If `None`, raise `RuntimeError` when called.
- **remover** (*callable*) – Method invoked on `rem()`. If `None`, raise `RuntimeError` when called.
- **expand\_tuple** (*bool*) – If `True` and if the first argument in the method call is a tuple, expand it as an argument list for the underlying function call.

```
get (obj, params=None)
```

```
set (obj, value)
```

```
rem (obj, params=None)
```

```
class pylablib.core.utils.functions.AttrObjectProperty (name,           use_getter=True,  
                                                         use_setter=True,  
                                                         use_remover=True,        ex-  
                                                         pand_tuple=True)
```

Bases: `pylablib.core.utils.functions.IObjectProperty`

Object property created from object attribute. Works with attributes or properties.

#### Parameters

- **name** (*str*) – Attribute name.
- **use\_getter** (*bool*) – If `False`, raise `RuntimeError` when calling `get` method.
- **use\_setter** (*bool*) – If `False`, raise `RuntimeError` when calling `set` method.
- **use\_remover** (*bool*) – If `False`, raise `RuntimeError` when calling `rem` method.
- **expand\_tuple** (*bool*) – If `True` and if the first argument in the method call is a tuple, expand it as an argument list for the underlying function call.

```
get (obj, params=None)
```

```
set (obj, value)
```

```
rem (obj, params=None)
```

`pylablib.core.utils.functions.empty_object_property(value=None)`

Dummy property which does nothing and returns *value* on *get* (None by default).

`pylablib.core.utils.functions.obj_prop(*args, **kwargs)`

Build an object property wrapper.

If no arguments (or a single None argument) are supplied, return a dummy property. If one argument is supplied, return `AttrObjectProperty` for a property with a given name. Otherwise, return `MethodObjectProperty` property.

`pylablib.core.utils.functions.as_obj_prop(value)`

Turn value into an object property using `obj_prop()` function.

If it's already `IOObjectProperty`, return unchanged. If *value* is a tuple, expand as an argument list.

`pylablib.core.utils.functions.delaydef(gen)`

Wrapper for a delayed definition of a function inside of a module.

Useful if defining a function is computationally costly. The wrapped function should be a generator of the target function rather than the function itself.

On the first call the generator is executed to define the target function, which is then substituted for all subsequent calls.

## pylablib.core.utils.general module

Collection of small utilities.

`pylablib.core.utils.general.set_props(obj, prop_names, props)`

Set multiple attributes of *obj*.

Names are given by *prop\_names* list and values are given by *props* list.

`pylablib.core.utils.general.get_props(obj, prop_names)`

Get multiple attributes of *obj*.

Names are given by *prop\_names* list.

`pylablib.core.utils.general.getattr_multivar(obj, attrs, **kwargs)`

Try to get an attribute of *obj* given a list *attrs* of its potential names.

If no attributes are found and default keyword argument is supplied, return this default value; otherwise, raise `AttributeError`.

`pylablib.core.utils.general.using_method(func, method_name=None, inherit_signature=True)`

Decorator that makes the function attempt to call the first argument's method instead of *func*.

Before calling the function, try and call a method of the first argument named *method\_name* (*func* name by default). If the method exists, call it instead of the wrapped function. If *inherit\_signature*==True, completely copy the signature of the wrapped method (name, args list, docstring, etc.).

`pylablib.core.utils.general.to_predicate(x)`

Turn *x* into a predicate.

If *x* is callable, it will be called with a single argument and returned value determines if the argument passes. If *x* is a container, an argument passes if it's contained in *x*.

`pylablib.core.utils.general.map_container(value, func)`

Map values in the container.

*value* can be a tuple, a list or a dict (mapping is applied to the values) raises `ValueError` if it's something else.

`pylablib.core.utils.general.as_container(val, t)`

Turn iterable into a container of type *t*.

Can handle named tuples, which have different constructor signature.

`pylablib.core.utils.general.recursive_map(value, func)`

Map container recursively.

*value* can be a tuple, a list or a dict (mapping is applied to the values).

`pylablib.core.utils.general.make_flat_namedtuple(nt, fields=None, name=None, subfield_fmt='{field:}_{subfield:}')`

Turn a nested structure of named tuples into a single flat namedtuple.

#### Parameters

- **nt** – toplevel namedtuple class to be flattened
- **fields** – a dictionary {name: desc} of the fields, where name is the named tuple name, and desc is either a nested namedtuple class, or a list of arguments which are passed to the recursive call to this function (e.g., [TTuple, {"field": TNestedTuple}]). Any tuple field which is present in this dictionary gets recursively flattened, and the field names of the corresponding returned tuple are added to the full list of fields
- **name** – name of the resulting tuple
- **subfield\_fmt** – format string, which describes how the combined field name is built out of the original field name and the subtuple field name; by default, connect with "\_", i.e., `t.field.subfield` turns into `t.field_subfield`.

**Returns** a new namedtuple class, which describes the flattened structure

`pylablib.core.utils.general.any_item(d)`

Return arbitrary tuple (key, value) contained in the dictionary (works both in Python 2 and 3)

`pylablib.core.utils.general.merge_dicts(*dicts)`

Combine multiple dict objects together.

If multiple dictionaries have the same keys, later arguments have higher priority.

`pylablib.core.utils.general.filter_dict(pred, d, exclude=False)`

Filter dictionary based on a predicate.

*pred* can be a callable or a container (in which case the predicate is true if a value is in the container). If *exclude==True*, the predicate is inverted.

`pylablib.core.utils.general.map_dict_keys(func, d)`

Map dictionary keys with *func*

`pylablib.core.utils.general.map_dict_values(func, d)`

Map dictionary values with *func*

`pylablib.core.utils.general.to_dict(d, default=None)`

Convert a dict or a list of pairs or single keys (or mixed) into a dict.

If a list element is single, *default* value is used.

`pylablib.core.utils.general.to_pairs_list(d, default=None)`

Convert a dict or a list of pairs or single keys (or mixed) into a list of pairs.

If a list element is single, *default* value is used. When converting list into list, the order is preserved.

`pylablib.core.utils.general.invert_dict(d, kmap=None)`

Invert dictionary (switch keys and values).

If *kmap* is supplied, it's a function mapping dictionary values into inverted dictionary keys (identity by default).

`pylablib.core.utils.general.flatten_list(l)`  
Flatten nested list/tuple structure into a single list.

`pylablib.core.utils.general.partition_list(pred, l)`  
Split the list *l* into two parts based on the predicate.

`pylablib.core.utils.general.split_in_groups(key_func, l, continuous=True, max_group_size=None)`  
Split the list *l* into groups according to the *key\_func*.

Go over the list and group the elements with the same key value together. If `continuous==False`, groups all elements with the same key together regardless of where they are in the list. otherwise, group only continuous sequences of the elements with the same key together (element with different key in the middle will result in two groups). If `continuous==True` and *max\_group\_size* is not None, it determines the maximal size of a group; larger groups are split into separate groups.

`pylablib.core.utils.general.sort_set_by_list(s, l, keep_duplicates=True)`  
Convert the set *s* into a list ordered by a list *l*.

Elements in *s* which are not in *l* are omitted. If `keep_duplicates==True`, keep duplicate occurrences in *l* in the result; otherwise, only keep the first occurrence.

`pylablib.core.utils.general.compare_lists(l1, l2, sort_lists=False, keep_duplicates=True)`  
Return three lists (*l1* and *l2*, *l1-l2*, *l2-l1*).

If `sort_lists==True`, sort the first two lists by *l1*, and the last one by *l2*; otherwise, the order is undefined. If `sort_lists==True`, *keep\_duplicated* determines if duplicate elements show up in the result.

`pylablib.core.utils.general.topological_order(graph, visit_order=None)`  
Get a topological order of a graph.

Return a list of nodes where each node is listed after its children. If *visit\_order* is not None, it is a list specifying nodes visiting order (nodes earlier in the list are visited first). Otherwise, the visit order is undefined. *graph* is a dictionary {node: [children]}. If graph contains loops, raise `ValueError`.

**class** `pylablib.core.utils.general.DummyResource`  
Bases: `object`

Object that acts as a resource (has `__enter__` and `__exit__` methods), but doesn't do anything.

Analog of:

```
@contextlib.contextmanager
def dummy_resource():
    yield
```

**class** `pylablib.core.utils.general.RetryOnException` (*tries=None, exceptions=None*)  
Bases: `object`

Wrapper for repeating the same block of code several time if an exception occurs

Useful for filesystem or communication operations, where retrying a failed operation is a valid option.

#### Parameters

- **tries** (*int*) – Determines how many time will the chunk of code execute before re-raising the exception; None (default) means no limit
- **exceptions** (*Exception or list*) – A single exception class or a list of exception classes which are going to be silenced.

Example:

```
for t in RetryOnException(tries, exceptions):
    with t:
        ... do stuff ...
```

is analogue of:

```
for i in range(tries):
    try:
        ... do stuff ...
    except exceptions:
        if i==tries-1:
            raise
```

```
class ExceptionCatcher(retrier, try_number)
    Bases: object

    reraise()
```

`pylablib.core.utils.general.retry_wait` (*func*, *try\_times=1*, *delay=0.0*, *exceptions=None*)

Try calling function (with no arguments) at most *try\_times* as long as it keeps raising exception.

If *exceptions* is not *None*, it specifies which exception types should be silenced. If an exception has been raised, wait *delay* seconds before retrying.

```
class pylablib.core.utils.general.SilenceException(exceptions=None,
                                                    on_exception=None,
                                                    reraise=False)

    Bases: object
```

Context which silences exceptions raised in a block of code.

#### Parameters

- **exceptions** (*Exception* or *list*) – A single exception class or a list of exception classes which are going to be silenced.
- **on\_exception** (*callable*) – A callback to be invoked if an exception occurs.
- **reraise** (*bool*) – Defines if the exception is re-raised after the callback has been invoked.

A simple bit of syntax sugar. The code:

```
with SilenceException(exceptions, on_exception, reraise):
    ... do stuff ...
```

is exactly analogous to:

```
try:
    ... do stuff ...
except exceptions:
    on_exception()
    if reraise:
        raise
```

`pylablib.core.utils.general.full_exit` (*code*=<*Signals.SIGTERM: 15*>)

Terminate the current process and all of its threads.

Doesn't perform any cleanup or resource release; should only be used if the process is irrevocably damaged.

```
class pylablib.core.utils.general.UIDGenerator(thread_safe=False)
    Bases: object
```

Generator of unique numeric IDs.

**Parameters** **thread\_safe** (*bool*) – If True, using lock to ensure that simultaneous calls from different threads are handled properly.

**reset** (*value=0*)

Reset the generator to the given value

```
class pylablib.core.utils.general.NamedUIDGenerator (name_template='{0}{1:03d}',
                                                    thread_safe=False)
```

Bases: `object`

Generator of unique string IDs based on a name.

#### Parameters

- **name\_template** (*str*) – Format string with two parameters (name and numeric ID) used to generate string IDs.
- **thread\_safe** (*bool*) – If True, using lock to ensure that simultaneous calls from different threads are handled properly.

```
pylablib.core.utils.general.call_limit (func, period=1, cooldown=0.0, limit=None, default=None)
```

Wrap *func* such that calls to it are forwarded only under certain conditions.

If *period*>1, then *func* is called after at least *period* calls to the wrapped function. If *cooldown*>0, then *func* is called after at least *cooldown* seconds passed since the last call. if *limit* is not None, then *func* is called only first *limit* times. If several conditions are specified, they should be satisfied simultaneously. *default* specifies return value if *func* wasn't called. Returned function also has an added method *reset*, which resets the internal call and time counters.

```
pylablib.core.utils.general.doc_inherit (parent)
```

Wrapper for inheriting docstrings from parent classes.

Takes parent class as an argument and replaces the docstring of the wrapped function by the docstring of the same-named function from the parent class (if available).

```
class pylablib.core.utils.general.Countdown (timeout, start=True)
```

Bases: `object`

Object for convenient handling of timeouts and countdowns with interrupts.

#### Parameters

- **timeout** (*float*) – Countdown timeout; if None, assumed to be infinite.
- **start** (*bool*) – if True, automatically start the countdown; otherwise, wait until *trigger()* is called explicitly

**reset** (*start=True*)

Restart the countdown from the current moment

**trigger** (*restart=True*)

Trigger the countdown.

If *restart*==True, restart the countdown if it's running; otherwise, do nothing in that situation.

**running** ()

Check if the countdown is running

**stop** ()

Stop the timer if currently running

**time\_left** (*t=None, bound\_below=True*)

Return the amount of time left. For infinite timeout, return `None`.

If `bound_below==True`, instead of negative time return zero. If *t* is supplied, it indicates the current time; otherwise, use `time.time()`.

**add\_time** (*dt, t=None, bound\_below=True*)

Add a given amount of time (positive or negative) to the start time (timeout stays the same).

If `bound_below==True`, do not let the end time (start time plus timeout) to get below the current time. If *t* is supplied, it indicates the current time; otherwise, use `time.time()`.

**set\_timeout** (*timeout*)

Change the timer timeout

**time\_passed** ()

Return the amount of time passed since the countdown start/reset, or `None` if it is not started

**passed** ()

Check if the timeout has passed

**class** `pylablib.core.utils.general.Timer` (*period, skip\_first=False*)

Bases: `object`

Object for keeping time of repeating tasks.

**Parameters** *period* (*float*) – Timer period.

**change\_period** (*period, method='current'*)

Change the timer period.

*method* specifies the changing method. Could be "current" (change the period of the ongoing tick), "next" (change the period starting from the next tick), "reset\_skip" (reset the timer and skip the first tick) or "reset\_noskip" (reset the timer and don't skip the first tick).

**reset** (*skip\_first=False*)

Reset the timer.

If `skip_first==False`, timer ticks immediately; otherwise, it starts ticking only after one period.

**time\_left** (*t=None, bound\_below=True*)

Return the amount of time left before the next tick.

If `bound_below==True`, instead of negative time return zero.

**passed** (*t=None*)

Return the number of ticks passed.

If timer period is zero, always return 1.

**acknowledge** (*n=None, nmin=0*)

Acknowledge the timer tick.

*n* specifies the number of tick to acknowledge (by default, all passed). Return number of actually acknowledged ticks (0 if the timer hasn't ticked since the last acknowledgement).

**class** `pylablib.core.utils.general.StreamFileLogger` (*path, stream=None, lock=None, autoflush=False*)

Bases: `object`

Stream logger that replaces standard output stream (usually stdout or stderr) and logs them into a file.

**Parameters**

- **path** – path to the destination logfile. The file is always appended.

- **stream** – an optional output stream into which the output will be duplicated; usually, the original stream which is being replaced
- **lock** – a thread lock object, which is used for any file writing operation; necessary if replacing standard streams (such as `sys.stdout` or `sys.stderr`) in a multithreading environment.
- **autoflush** – if `True`, flush after any write operation into *stream*

It is also possible to subclass the file and overload `write_header()` method to write a header before the first file write operation during the execution.

The intended use is to log stdout or stderr streams:

```
import sys, threading
sys.stderr = StreamFileLogger("error_log.txt", stream=sys.stderr, lock=threading.
    ↳ Lock())
```

**write\_header**(*f*)  
Write header to file stream *f*

**add\_path**(*path*)  
Add another logging path to the list

**remove\_path**(*path*)  
Remove logging path to the list

**write**(*s*)

**flush**()

`pylablib.core.utils.general.setbp()`

`pylablib.core.utils.general.timing`(*n=1, name=None, profile=False*)  
Context manager for timing a piece of code.

Measures the time it takes to execute the wrapped code and prints the result.

#### Parameters

- **n** – can specify the number of repetitions, which is used to show time per single repetition.
- **name** – name which is printed alongside the time
- **profile** – if `True`, use `cProfile` and print its output instead of a simple timing

**class** `pylablib.core.utils.general.AccessIterator`(*obj, access\_function=None*)  
Bases: `object`

Simple sequential access iterator with customizable access function (by default it's 1D indexing).

Determines end of iterations by `IndexError`.

#### Parameters

- **obj** – Container to be iterated over.
- **access\_function** (*callable*) – A function which takes two parameters *obj* and *idx* and either returns the element or raises `IndexError`. By default, a simple `__getitem__` operation.

**next**()

`pylablib.core.utils.general.muxcall`(*argname, special\_args=None, mux\_argnames=None, return\_kind='list', allow\_partial=False*)

Wrap a function such that it can become multiplexable over a given argument.



## Parameters

- **argname** – name of the argument to loop over
- **special\_args** – if not `None`, defines a dictionary `{arg: func}` for special values of the argument (e.g., "all", `None`, etc.), where `arg` is its value, and `func` is the method taking the same arguments as the called function and returning the substitute argument (e.g., a list of all arguments)
- **mux\_argnames** – names of additional arguments which, when supplied list or dict values, and when the `argname` value is a list, specify different values for different calls
- **return\_kind** – method to combined multiple returned values; can be "list", "dict" (return dict `{arg: result}`), or "none" (simply return `None`)
- **allow\_partial** – if `True` and some of `mux_argnames` argument do not specify value for the full range of `argname` value, do not call the function for those unspecified values; otherwise (`allow_partial` is `True`), the error will be raised

`pylablib.core.utils.general.wait_for_keypress` (*message='Waiting...'*)

`pylablib.core.utils.general.restart` ()

Restart the script.

Execution will not resume after this call. Note: due to Windows limitations, this function does not replace the current process with a new one, but rather calls a new process and makes the current one wait for its execution. Hence, each nested call adds an additional loaded application into the memory. Therefore, nesting restart calls (i.e., calling several restarts in a row) should be avoided.

## pylablib.core.utils.indexing module

Processing and normalization of different indexing styles.

`pylablib.core.utils.indexing.string_list_idx` (*names\_to\_find,* *names\_list,*  
*only\_exact=False*)

Index through a list of strings in *names\_list*.

Return corresponding numerical indices. Case sensitive; first look for exact matching, then for prefix matching (unless `only_exact=True`).

`pylablib.core.utils.indexing.is_slice` (*idx*)

Check if *idx* is slice.

`pylablib.core.utils.indexing.is_range` (*idx*)

Check if *idx* is iterable (list, numpy array, or *builtins.range*).

`pylablib.core.utils.indexing.is_bool_array` (*idx*)

Check if *idx* is a boolean array.

`pylablib.core.utils.indexing.to_range` (*idx, length*)

Turn list, array, *builtins.range*, slice into an iterable.

`pylablib.core.utils.indexing.covers_all` (*idx, length, strict=False, ordered=True*)

Check if *idx* covers all of the elements (indices from 0 to *length*).

If `strict==True`, strictly checks the condition; otherwise may return `False` even if *idx* actually covers everything, but takes less time (i.e., can be used for optimization). If `ordered==True`, only returns `True` when indices follow in order.

**class** `pylablib.core.utils.indexing.IIndex`

Bases: `object`

A generic index object.

Used to transform a variety of indexes into a subset applicable for specific objects (numpy arrays or lists).

**Allowed input index types:**

- scalar: integer, string
- vector: integer lists or numpy arrays, bool lists or numpy arrays, string lists or numpy arrays, builtin.ranges, slices and string slices

**tup()**

Represent index as a tuple for easy unpacking.

**class** pylablib.core.utils.indexing.**NumpyIndex**(*idx, ndim=None*)

Bases: [pylablib.core.utils.indexing.IIndex](#)

NumPy compatible index: allows for integers, slices, numpy integer or boolean arrays, integer lists or builtin.ranges.

**Parameters**

- **idx** – raw index
- **ndim** – index dimensionality (either 0 or 1); if supplied, assume that *idx* is already normalized

**tup()**

Represent index as a tuple for easy unpacking.

**class** pylablib.core.utils.indexing.**ListIndex**(*idx, names=None, ndim=None*)

Bases: [pylablib.core.utils.indexing.IIndex](#)

List compatible index: allows for integers, slices, numpy integer arrays, integer lists or builtin.ranges.

**Parameters**

- **idx** – raw index
- **names** – list of allowed index string values, which is used to convert them into integers
- **ndim** – index dimensionality (either 0 or 1); if supplied, assume that *idx* is already normalized

**tup()**

Represent index as a tuple for easy unpacking.

**class** pylablib.core.utils.indexing.**ListIndexNoSlice**(*idx, names=None, length=None, ndim=None*)

Bases: [pylablib.core.utils.indexing.ListIndex](#)

List compatible index with slice unwrapped into builtin.range: allows for integers, numpy integer arrays, integer lists or builtin.ranges.

**Parameters**

- **idx** – raw index
- **names** – list of allowed index string values, which is used to convert them into integers
- **length** – length of the list (used to expand slice indices)
- **ndim** – index dimensionality (either 0 or 1); if supplied, assume that *idx* is already normalized

**tup()**

Represent index as a tuple for easy unpacking.

`pylablib.core.utils.indexing.to_double_index (idx, names)`

Convert double index into a pair of indexes.

Assume that one index is purely numerical, while the other can take names (out of the supplied list).

#### Parameters

- **idx** – raw double index
- **names** – list of allowed index string values, which is used to convert them into integers

## pylablib.core.utils.ipc module

Universal interface for inter-process communication.

Focus on higher throughput for large numpy arrays via shared memory.

**class** `pylablib.core.utils.ipc.IIPCChannel`

Bases: `object`

Generic IPC channel interface

**send** (*data*)

Send data

**recv** (*timeout=None*)

Receive data

**send\_numpy** (*data*)

Send numpy array

**recv\_numpy** (*timeout=None*)

Receive numpy array

**get\_peer\_args** ()

Get arguments required to create a peer connection

**classmethod from\_args** (\**args*)

Create a peer connection from the supplied arguments

**class** `pylablib.core.utils.ipc.TPipeMsg (id, data)`

Bases: `tuple`

**count** ()

Return number of occurrences of value.

**data**

**id**

**index** ()

Return first index of value.

Raises ValueError if the value is not present.

**class** `pylablib.core.utils.ipc.PipeIPCChannel (pipe_conn=None)`

Bases: `pylablib.core.utils.ipc.IIPCChannel`

Generic IPC channel interface using pipe.

**get\_peer\_args** ()

Get arguments required to create a peer connection

**send** (*data*)  
Send data

**recv** (*timeout=None*)  
Receive data

**classmethod from\_args** (*\*args*)  
Create a peer connection from the supplied arguments

**recv\_numpy** (*timeout=None*)  
Receive numpy array

**send\_numpy** (*data*)  
Send numpy array

**class** `pylablib.core.utils.ipc.SharedMemIPCChannel` (*pipe\_conn=None, arr=None, arr\_size=None*)

Bases: `pylablib.core.utils.ipc.PipeIPCChannel`

Generic IPC channel interface using pipe and shared memory for large arrays.

**get\_peer\_args** ()  
Get arguments required to create a peer connection

**send\_numpy** (*data, method='auto', timeout=None*)  
Send numpy array

**recv\_numpy** (*timeout=None*)  
Receive numpy array

**classmethod from\_args** (*\*args*)  
Create a peer connection from the supplied arguments

**recv** (*timeout=None*)  
Receive data

**send** (*data*)  
Send data

**class** `pylablib.core.utils.ipc.TShmemVarDesc` (*offset, size, kind, fixed\_size*)

Bases: `tuple`

**count** ()  
Return number of occurrences of value.

**fixed\_size**

**index** ()  
Return first index of value.  
  
Raises ValueError if the value is not present.

**kind**

**offset**

**size**

**class** `pylablib.core.utils.ipc.SharedMemIPCTable` (*pipe\_conn=None, arr=None, arr\_size=None, lock=True*)

Bases: `object`

Shared memory table for exchanging shared variables between processes.

Can be used instead of channels for variables which are rarely changed but frequently checked (e.g., status), or when synchronization of sending and receiving might be difficult

**add\_variable** (*name*, *size*, *kind*=*'pickle'*)

Add a variable with a given name.

The variable info is also communicated to the other endpoint. *size* determines maximal variable size in bytes. If the actual size ever exceeds it, an exception will be raised. *kind* determines the way to convert variable into bytes; can be *"pickle"* (universal, but large size overhead), *"nps\_###"* (where *###* can be any numpy scalar dtype description, e.g., *"float"* or *"<u2"*) for numpy scalars, or *"npa\_###"* (where *###* means the same as for *nps*) for numpy arrays (in this case the array size and shape need to be communicated separately).

**set\_variable** (*name*, *value*)

Set a variable with a given name.

If the variable is missing, raise an exception.

**get\_variable** (*name*, *default*=*None*)

Get a variable with a given name.

If the variable is missing, return *default*.

**is\_peer\_connected** ()

Check if the peer is connected (i.e., the other side of the pipe is initialized)

**close\_connection** ()

Mark the connection as closed

**is\_peer\_closed** ()

Check if the peer is closed

**get\_peer\_args** ()

Get arguments required to create a peer connection

**classmethod from\_args** (\**args*)

Create a peer connection from the supplied arguments

## pylablib.core.utils.library\_parameters module

Storage for global library parameters

`pylablib.core.utils.library_parameters.temp_library_parameters` (*restore*=*None*)

Context manager, which restores library parameters upon exit.

If *rester* is not *None*, it can specify a list of parameters to be restored (by default, all parameters).

## pylablib.core.utils.module module

Library for dealing with python module properties.

`pylablib.core.utils.module.get_package_version` (*pkg*)

Get the version of the package.

If the package version is unavailable, return *None*.

`pylablib.core.utils.module.cmp_versions` (*ver1*, *ver2*)

Compare two package versions.

Return *'<'* if the first version is older (smaller), *'>'* if it's younger (larger) or *'='* if it's the same.

`pylablib.core.utils.module.cmp_package_version(pkg, ver)`

Compare current package version to *ver*.

*ver* should be a name of the package (rather than the module). Return '<' if current version is older (smaller), '>' if it's younger (larger) or '=' if it's the same. If the package version is unavailable, return None.

`pylablib.core.utils.module.expand_relative_path(module_name, rel_path)`

Turn a relative module path into an absolute one.

*module\_name* is the absolute name of the reference module, *rel\_path* is the path relative to this module.

`pylablib.core.utils.module.get_loaded_package_modules(pkg_name)`

Get all modules in the package *pkg\_name*.

Returns a dict {name: module}.

`pylablib.core.utils.module.get_imported_modules(module, explicit=False)`

Get modules imported within a given module.

If *explicit*==True, take into account only toplevel objects which are modules (corresponds to `import module` or `from package import module` statements) If *explicit*==False, also include all modules containing toplevel objects (corresponds to `from module import Class` or `from package import function` statements). Return a dictionary {name: module} (modules with the same name are considered to be the same).

`pylablib.core.utils.module.get_reload_order(modules)`

Find reload order for modules which respects dependencies (a module is loaded before its dependents).

*modules* is a dict {name: module}.

The module dependencies (i.e., the modules which the current module depends on) are determined based on imported modules and modules containing toplevel module objects.

`pylablib.core.utils.module.reload_package_modules(pkg_name, ignore_errors=False)`

Reload package *pkg\_name*, while respecting dependencies of its submodules.

If *ignore\_errors*=True, ignore `ImportError` exceptions during the reloading process.

`pylablib.core.utils.module.unload_package_modules(pkg_name, ignore_errors=False)`

Reload package *pkg\_name*, while respecting dependencies of its submodules.

If *ignore\_errors*=True, ignore `ImportError` exceptions during the reloading process.

`pylablib.core.utils.module.get_library_path()`

Get a filesystem path for the pyLabLib library (the one containing current the module).

`pylablib.core.utils.module.get_library_name()`

Get the name for the pyLabLib library (the one containing current the module).

`pylablib.core.utils.module.get_executable(console=False)`

Get Python executable.

If *console*==True and the current executable is windowed (i.e., "pythonw.exe"), return the corresponding "python.exe" instead.

`pylablib.core.utils.module.get_python_folder()`

Return Python interpreter folder (the folder containing the python executable)

`pylablib.core.utils.module.pip_install(pkg, upgrade=False)`

Call `pip install` for a given package.

If *upgrade*==True, call with `--upgrade` key (upgrade current version if it is already installed).

`pylablib.core.utils.module.install_if_older(pkg, min_ver="")`

Install *pkg* from the default PyPI repository if its version is lower than *min\_ver*

If *min\_ver* is None, upgrade to the newest version regardless; if *min\_ver*="", install only if no version is installed. Return True if the package was installed.

### pylablib.core.utils.nbtools module

`pylablib.core.utils.nbtools.c_array`

Generate a numba C-ordered array type with the given element type, number of dimensions, and read-only flag

`pylablib.core.utils.nbtools.au1(x, off)`

Extract a little-endian 1-byte unsigned integer from a numpy byte array at the given offset

`pylablib.core.utils.nbtools.au2(x, off)`

Extract a little-endian 2-byte unsigned integer from a numpy byte array at the given offset

`pylablib.core.utils.nbtools.au4(x, off)`

Extract a little-endian 4-byte unsigned integer from a numpy byte array at the given offset

`pylablib.core.utils.nbtools.au8(x, off)`

Extract a little-endian 8-byte unsigned integer from a numpy byte array at the given offset

`pylablib.core.utils.nbtools.ai1(x, off)`

Extract a little-endian 1-byte unsigned integer from a numpy byte array at the given offset

`pylablib.core.utils.nbtools.ai2(x, off)`

Extract a little-endian 2-byte unsigned integer from a numpy byte array at the given offset

`pylablib.core.utils.nbtools.ai4(x, off)`

Extract a little-endian 4-byte unsigned integer from a numpy byte array at the given offset

`pylablib.core.utils.nbtools.ai8(x, off)`

Extract a little-endian 8-byte unsigned integer from a numpy byte array at the given offset

`pylablib.core.utils.nbtools.copy_array_chunks`

Generate and compile a numba function for copying an array in chunks. *base* specifies the base array type (by default, unsigned byte); if *par*==True, generate a parallelized implementation. if *nogil*==True, use the *nogil* numba option to release GIL during the execution.

The returned function takes 4 arguments: source array, destination array, number of chunks, and size (in elements) of each chunk.

`pylablib.core.utils.nbtools.copy_array_strided`

Generate and compile a numba function for copying an array in chunks with an arbitrary stride. *base* specifies the base array type (by default, unsigned byte); if *par*==True, generate a parallelized implementation. if *nogil*==True, use the *nogil* numba option to release GIL during the execution.

The returned function takes 6 arguments: source array, destination array, number of chunks, size (in elements) of each chunk, chunks stride (in elements) in the source array, and offset (in elements) from the beginning of the first array. If size is the same as stride and the offset is zero, this function would mimic the one generated by `copy_array_chunks()`.

### pylablib.core.utils.net module

A wrapper for built-in TCP/IP routines.

**exception** `pylablib.core.utils.net.SocketError`

Bases: `OSError`

Base socket error class.

**args**

**characters\_written**

**errno**

POSIX exception code

**filename**

exception filename

**filename2**

second exception filename

**strerror**

exception strerror

**with\_traceback()**

Exception.with\_traceback(tb) – set self.\_\_traceback\_\_ to tb and return self.

**exception** `pylablib.core.utils.net.SocketTimeout`

Bases: `pylablib.core.utils.net.SocketError`

Socket timeout error.

**args**

**characters\_written**

**errno**

POSIX exception code

**filename**

exception filename

**filename2**

second exception filename

**strerror**

exception strerror

**with\_traceback()**

Exception.with\_traceback(tb) – set self.\_\_traceback\_\_ to tb and return self.

`pylablib.core.utils.net.get_local_addr()`

Get local IP address

`pylablib.core.utils.net.get_all_local_addr()`

Get a list of all local IP addresses

`pylablib.core.utils.net.get_local_hostname(full=True)`

Get a local host name

`pylablib.core.utils.net.get_all_remote_addr(hostname)`

Get a list of all remote addresses of a remote host by name

`pylablib.core.utils.net.get_remote_hostname(addr, error_on_missing=False)`

Get a remote host name by its address

`pylablib.core.utils.net.as_addr_port(addr, port)`

Parse the given address and port combination.

*addr* can be a host address, a tuple (*addr*, *port*), or a string "*addr:port*"; in the first case the given *port* is used, while in the other two it is ignore. Return tuple (*addr*, *port*).



```
class pylablib.core.utils.net.ClientSocket (sock=None, timeout=None,  

                                             wait_callback=None, send_method='declen',  

                                             recv_method='declen', datatype='auto',  

                                             nodelay=False)
```

Bases: `object`

A client socket (used to connect to a server socket).

#### Parameters

- **sock** (*socket.socket*) – Socket to wrap; if `None` create a new one.
- **timeout** (*float*) – The timeout used for connecting and sending/receiving (`None` means no timeout).
- **wait\_callback** (*callable*) – Called periodically (every 100ms by default) while waiting for connecting or sending/receiving.
- **send\_method** (*str*) – Default sending method.
- **recv\_method** (*str*) – Default receiving method.
- **datatype** (*str*) – Type of the returned data; can be `"bytes"` (return *bytes* object), `"str"` (return *str* object), or `"auto"` (default Python result: *str* in Python 2 and *bytes* in Python 3)
- **nodelay** (*bool*) – Whether to enable `TCP_NODELAY`.

#### Possible sending/receiving methods are:

- `'fixedlen'`: data is sent as is, and receiving requires to know the length of the message;
- `'declen'`: data is prepended by a length, and receiving reads this length and doesn't need pre-determined length info.

#### **sock**

Corresponding Python socket.

**Type** `socket.socket`

#### **declen\_bo**

Byteorder of the prepended length for `'declen'` sending method. Can be either `'>'` (big-endian, default) or `'<'`.

**Type** `str`

#### **declen\_ll**

Length of the prepended length for `'declen'` sending method; default is 4 bytes (corresponding to maximum of 4Gb per single length-prepended message)

**Type** `int`

#### **set\_wait\_callback** (*wait\_callback=None*)

Set callback function for waiting during connecting or sending/receiving

#### **set\_timeout** (*timeout=None*)

Set timeout for connecting or sending/receiving

#### **get\_timeout** ()

Get timeout for connecting or sending/receiving

#### **using\_timeout** (*timeout=None*)

Context manager for usage of a different timeout inside a block

**connect** (*host, port*)

Connect to a remote host

**close** ()

Close the connection

**is\_connected** ()

Check if the connection is opened

**get\_local\_name** ()

Return IP address and port of this socket

**get\_peer\_name** ()

Return IP address and port of the peer socket

**recv\_fixedlen** (*l*)

Receive fixed-length message of length *l*

**recv\_delimiter** (*delim, lmax=None, chunk\_l=1024, strict=False*)

Receive a single message ending with a delimiter *delim* (can be several characters, or list several possible delimiter strings).

*lmax* specifies the maximal received length (*None* means no limit). *chunk\_l* specifies the size of data chunk to be read in one try. If *strict==False*, keep receiving as much data as possible until a delimiter is found in the end (only works properly if a single line is expected); otherwise, receive the data byte-by-byte and stop as soon as a delimiter is found (equivalent to setting *chunk\_l=1*).

**recv\_decllen** ()

Receive variable-length message (prepended by its length).

Length format is described by *declen\_bo* and *declen\_ll* attributes.

**recv** (*l=None*)

Receive a message using the default method.

**recv\_all** (*chunk\_l=1024*)

Receive all of the data currently in the socket.

*chunk\_l* specifies the size of data chunk to be read in one try. For technical reasons, use 1ms timeout (i.e., this operation takes 1ms).

**recv\_ack** (*l=None*)

Receive a message using the default method and send an acknowledgement (message length)

**send\_fixedlen** (*msg*)

Send a message as is

**send\_decllen** (*msg*)

Send a message as a variable-length (prepending its length in the sent message).

Length format is described by *declen\_bo* and *declen\_ll* attributes.

**send\_delimiter** (*msg, delimiter*)

Send a message with a delimiter *delim* (can be several characters)

**send** (*msg*)

Send a message using the default method.

**send\_ack** (*msg*)

Send a message using default method and wait for acknowledgement (message length).

If the acknowledgement message length doesn't agree, raise *SocketError*.

```
pylablib.core.utils.net.recv_JSON(sock, chunk_l=1024, strict=True)
```

Receive a complete JSON token from the socket.

*chunk\_l* specifies the size of data chunk to be read in one try. If *strict*==False, keep receiving as much data as possible until the received data forms a complete JSON token. otherwise, receive the data byte-by-byte and stop as soon as a token is formed (equivalent to setting *chunk\_l*=1).

```
pylablib.core.utils.net.listen(host, port, conn_func, port_func=None, wait_callback=None,
                               timeout=None, backlog=10, wrap_socket=True, connections_number=None, socket_kwargs=None)
```

Run a server socket at the given host and port.

#### Parameters

- **host** (*str*) – Server host address. If None, use the local host defined by `socket.gethostname()`.
- **port** (*int*) – Server port. If 0, generate an arbitrary free port.
- **conn\_func** (*callable*) – Called with the client socket as a single argument every time a connection is established.
- **port\_func** (*callable*) – Called with the port as a single argument when the listening starts (useful with *port*=0).
- **wait\_callback** (*callable*) – A callback function which is called periodically (every 100ms by default) while awaiting for connections.
- **timeout** (*float*) – Timeout for waiting for the connections (None is no timeout).
- **backlog** (*int*) – Backlog length for the socket (see `socket.socket.listen()`).
- **wrap\_socket** (*bool*) – If True, wrap the client socket of the connection into `ClientSocket` class; otherwise, return `socket.socket` object.
- **connections\_number** (*int*) – Specifies maximal number of connections before the listening function returns (by default, the number is unlimited).
- **socket\_kwargs** (*dict*) – additional keyword arguments passed to `ClientSocket` constructor.

Checking for connections is paused until *conn\_func* returns. If multiple simultaneous connections are expected, *conn\_func* should spawn a separate processing thread and return. If *connections\_number* is None (i.e., there's no limit on the number of connections before closing), this function never returns.

### pylablib.core.utils.numerical module

Numerical functions that don't deal with sequences.

```
pylablib.core.utils.numerical.gcd(*numbers)
```

Euclid's algorithm for GCD. Arguments are cast to integer

```
pylablib.core.utils.numerical.integer_distance(x)
```

Get distance to the closes integer

```
pylablib.core.utils.numerical.gcd_approx(a, b, min_fraction=1e-08, tolerance=1e-05)
```

Approximate Euclid's algorithm for possible non-integer values.

Try to find a number *d* such that *a*/*d* and *b*/*d* are less than *tolerance* away from a closest integer. If GCD becomes less than *min\_fraction* \* min(*a*, *b*), raise `ArithmeticError`.

```
pylablib.core.utils.numerical.round_significant(x, n)
```

Rounds *x* to *n* significant digits (not the same as *n* decimal places!).

`pylablib.core.utils.numerical.limit_to_range(x, min_val=None, max_val=None, default=0)`

Confine  $x$  to the given limit.

Default limit values are `None`, which means no limit. *default* specifies returned value if both  $x$ , *min\_val* and *max\_val* are `None`.

**class** `pylablib.core.utils.numerical.infinite_list(start=0, step=1)`

Bases: `object`

Mimics the behavior of the usual list, but is infinite and immutable.

Supports accessing elements, slicing (including slices giving infinite lists) and iterating. Iterating over it naturally leads to an infinite loop, so it should only be used either for finite slices or for loops with break condition.

#### Parameters

- **start** – The first element of the list.
- **step** – List step.

**class** `counter(lst)`

Bases: `object`

**next()**

`pylablib.core.utils.numerical.unity()`

Return a unity function

`pylablib.core.utils.numerical.constant(c)`

Return a function which returns a constant  $c$ .

$c$  can only be either a scalar, or an array-like object with the shape matching the expected argument.

`pylablib.core.utils.numerical.polynomial(coeffs)`

Return a polynomial function which with coefficients *coeffs*.

Coefficients are list lowest-order first, so that `coeffs[i]` is the coefficient in front of  $x^{**i}$ .

## pylablib.core.utils.observer\_pool module

A simple observer pool (notification pool) implementation.

**class** `pylablib.core.utils.observer_pool.ObserverPool(expand_tuple=True)`

Bases: `object`

An observer pool.

Stores notification functions (callbacks), and calls them whenever `notify()` is called. The callbacks can have priority (higher priority ones are called first) and filter (observer is only called if the filter function passes the notification tag).

**Parameters** **expand\_tuple** (*bool*) – if `True` and the notification value is a tuple, treat it as an argument list for the callback functions.

**class** `Observer(filt, callback, priority, attr, cacheable)`

Bases: `tuple`

**attr**

**cacheable**

**callback**

**count** ()  
Return number of occurrences of value.

**filt**

**index** ()  
Return first index of value.  
  
Raises ValueError if the value is not present.

**priority**

**add\_observer** (*callback*, *name=None*, *filt=None*, *priority=0*, *attr=None*, *cacheable=False*)  
Add the observer callback.

#### Parameters

- **callback** (*callable*) – callback function; takes at least one argument (notification tag), and possible more depending on the notification value.
- **name** (*str*) – stored callback name; by default, a unique name is auto-generated
- **filt** (*callable or None*) – a filter function for this observer (the observer is called only if the *notify()* function tag and value pass the filter); by default, all tags are accepted
- **priority** (*int*) – callback priority; higher priority callback are invoked first.
- **attr** – additional observer attributes (can be used by *ObserverPool* subclasses to change their behavior).
- **cacheable** (*bool*) – if True, assumes that the filter function only depends on the tag, so its calls can be cached.

**Returns** callback name (equal to *name* if supplied; an automatically generated name otherwise).

**remove\_observer** (*name*)  
Remove the observer callback with the given name

**find\_observers** (*tag*, *value*)

**notify** (*tag*, *value=()*)  
Notify the observers by calling their callbacks.

Return a dictionary of the callback results. By default the value is an empty tuple: for *expand\_tuple==True* this means that only one argument (*tag*) is passed to the callbacks.

## pylablib.core.utils.py3 module

Dealing with Python2 / Python3 compatibility.

**pylablib.core.utils.py3.as\_str** (*data*)  
Convert a string into a text string

**pylablib.core.utils.py3.as\_bytes** (*data*)  
Convert a string into bytes

**pylablib.core.utils.py3.as\_builtin\_bytes** (*data*)  
Convert a string into bytes

`pylablib.core.utils.py3.as_datatype (data, datatype)`

Convert a string into a given datatypes.

*datatype* can be "str" (text string), "bytes" (byte string), or "auto" (no conversion).

## pylablib.core.utils.rpyc\_utils module

Routines and classes related to RPyC package

`pylablib.core.utils.rpyc_utils.obtain (proxy, serv=None, deep=False, direct=False)`

Obtain a remote netref object by value (i.e., copy it to the local Python instance).

Wrapper around `rpyc.utils.classic.obtain()` with some special cases handling. *serv* specifies the current remote service. If it is of type `SocketTunnelService`, use its socket tunnel for faster transfer. If *deep*==True and *proxy* is a container (tuple, list, or dict), run the function recursively for all its sub-elements. If *direct*==True, directly use RPyC obtain method; otherwise use the custom method, which works better with large numpy arrays, but worse with composite types (e.g., lists).

`pylablib.core.utils.rpyc_utils.transfer (obj, serv)`

Send a local object to the remote PC by value (i.e., copy it to the remote Python instance).

A 'reversed' version of `obtain()`.

**class** `pylablib.core.utils.rpyc_utils.SocketTunnelService (server=False)`

Bases: `sphinx.ext.autodoc.importer._MockObject`

Extension of the standard `rpyc.core.service.SlaveService` with built-in network socket tunnel for faster data transfer.

In order for the tunnel to work, services on both ends need to be subclasses of `SocketTunnelService`. Because of the initial setup protocol, the two services are asymmetric: one should be 'server' (corresponding to the listening server), and one should be 'client' (external connection). The roles are decided by the *server* constructor parameter.

**tunnel\_send** (*obj*, *packer*=None)

Send data through the socket tunnel.

If *packer* is not None, it defines a function to convert *obj* to a bytes string.

**tunnel\_recv** (*unpacker*=None)

Receive data sent through the socket tunnel.

If *unpacker* is not None, it defines a function to convert the received bytes string into an object.

**obtain** (*proxy*)

Execute `obtain()` on the local instance

**transfer** (*obj*)

Execute `transfer()` on the local instance

**on\_connect** (*conn*)

**on\_disconnect** (*conn*)

**class** `pylablib.core.utils.rpyc_utils.DeviceService (verbose=False)`

Bases: `pylablib.core.utils.rpyc_utils.SocketTunnelService`

Device RPyC service.

Expands on `SocketTunnelService` by adding `get_device()` method, which opens local devices, tracks them, and closes them automatically on disconnect.

**on\_connect** (*conn*)

**on\_disconnect** (*conn*)

**get\_device\_class** (*cls*)  
Get remote device class.

*cls* is the full class name, including the module within `pylablib.devices` (e.g., `Attocube.ANC300`).

**get\_device** (*cls*, *\*args*, *\*\*kwargs*)  
Connect to a device.

*cls* is the full class name, including the module within `pylablib.devices` (e.g., `Attocube.ANC300`). Stores reference to the connected device and closes it automatically on disconnect.

**obtain** (*proxy*)  
Execute `obtain()` on the local instance

**transfer** (*obj*)  
Execute `transfer()` on the local instance

**tunnel\_recv** (*unpacker=None*)  
Receive data sent through the socket tunnel.

If *unpacker* is not `None`, it defines a function to convert the received bytes string into an object.

**tunnel\_send** (*obj*, *packer=None*)  
Send data through the socket tunnel.

If *packer* is not `None`, it defines a function to convert *obj* to a bytes string.

`pylablib.core.utils.rpyc_utils.run_device_service` (*port=18812*, *verbose=False*)  
Start `DeviceService` at the given port

`pylablib.core.utils.rpyc_utils.connect_device_service` (*addr*, *port=18812*, *timeout=3*, *attempts=2*, *error\_on\_fail=True*, *config=None*)

Connect to the `DeviceService` running at the given address and port

*timeout* and *attempts* define respectively timeout of a single connection attempt, and the number of attempts (RPyC default is 3 seconds timeout and 6 attempts). If *error\_on\_fail==True*, raise error if the connection failed; otherwise, return `None`

## pylablib.core.utils.strdump module

Utils for converting variables into standard python objects (lists, dictionaries, strings, etc.) and back (e.g., for a more predictable LAN transfer). Provides an extension for pickle for more customized classes (numpy arrays, Dictionary).

**class** `pylablib.core.utils.strdump.StrDumper`  
Bases: `object`

Class for dumping and loading an object.

Stores procedures for dumping and loading, i.e., conversion from complex classes (such as `Dictionary`) to simple built-in classes (such as `dict` or `str`).

**add\_class** (*cls*, *dumpf=None*, *loadf=None*, *name=None*, *allow\_subclass=True*, *recursive=False*)  
Add a rule for dumping/loading an object of class *cls*.

### Parameters

- *cls* –

- **dumpf** (*callable*) – Function for dumping an object of the class; None means identity function.
- **loadf** (*callable*) – Function for loading an object of the class; None means identity function.
- **name** (*str*) – Name of class, which is stored in the packed data (`cls.__name__` by default).
- **allow\_subclass** (*bool*) – If True, this rule is also used for subclasses of this class.
- **recursive** (*bool*) – If True, the functions are given a second argument, which is a dumping/loading function for their sub-elements.

**dump** (*obj*)

Convert an object into a dumped value

**load** (*obj*)

Convert a dumped value into an object

**loads** (*s*)

Convert a pickled string of a damped object into an object

**dumps** (*obj*)

Dump an object into a pickled string

`pylablib.core.utils.strdump.dumper = <pylablib.core.utils.strdump.StrDumper object>`

Default dumper for converting into standard Python classes and pickling.

Converts `numpy.ndarray` and `Dictionary` objects (these conversion routines are defined when corresponding modules are imported). The converted values include non-printable characters (conversion uses `numpy.load()` and `numpy.ndarray.dump()`), so they can't be saved into text files. However, they're suited for pickling.

`pylablib.core.utils.strdump.dump(obj)`

Convert obj into standard Python classes using the default dumper

`pylablib.core.utils.strdump.load(s)`

Convert standard Python class representation *s* into an object using the default dumper

`pylablib.core.utils.strdump.dumps(obj)`

Convert obj into a pickled string using the default dumper

`pylablib.core.utils.strdump.loads(s)`

Convert a pickled string into an object using the default dumper

## pylablib.core.utils.string module

String search, manipulation and conversion routines.

`pylablib.core.utils.string.string_equal(name1, name2, case_sensitive=True, as_prefix=False)`

Determine if *name1* and *name2* are equal with taking special rules (*case\_sensitive* and *as\_prefix*) into account.

If *as\_prefix*==True, strings match even if *name1* is just a prefix of *name2* (not the other way around).

`pylablib.core.utils.string.find_list_string(name, str_list, case_sensitive=True, as_prefix=False, first_matched=False)`

Find *name* in the string list.



Comparison parameters are defined in `string_equal()`. If `first_matched==True`, stop at the first match; otherwise if multiple occurrences happen, raise `ValueError`.

**Returns** tuple (index, value).

```
pylablib.core.utils.string.find_dict_string(name, str_dict, case_sensitive=True,
                                             as_prefix=False)
```

Find name in the string dictionary.

Comparison parameters are defined in `string_equal()`. If multiple occurrences happen, raise `ValueError`.

**Returns** tuple (key, value).

```
pylablib.core.utils.string.find_first_entry(line, elements, start=0, not_found_value=-1)
```

Find the index of the earliest position inside the *line* of any of the strings in *elements*, starting from *start*.

If none are found, return *not\_found\_value*.

```
pylablib.core.utils.string.find_all_first_locations(line, elements, start=0,
                                                    not_found_value=-1,
                                                    known_locations=None)
```

Find the indices of the earliest position inside the *line* of all of the strings in *elements*, starting from *start*.

Return dict {element: pos}, where pos is either position in the string, or *not\_found\_value* if no entries are present. *known\_locations* can specify a dictionary of already known locations of some of the elements. In this case, only missing elements or elements located before *start* will be re-evaluated.

```
pylablib.core.utils.string.translate_string_filter(filt, syntax, match_case=True,
                                                    default=False)
```

Turns *filt* into a matching function.

The matching function takes single *str* argument, returns *bool* value.

*filt* can be

- *None*: function always returns default,
- *bool*: function always returns this value,
- *str*: pattern, determined by *syntax*,
- anything else: returned as is (assumed to already be a callable).

*syntax* can be 're' (re), 'glob' (glob) or 'pred' (simply matching predicate). *match\_case* determines whether the filter cares about the string case when matching.

```
class pylablib.core.utils.string.StringFilter(include=None, exclude=None, syntax='re', match_case=False)
```

Bases: `object`

String filter function.

Matches string if it matches include (matches all strings by default) and doesn't match exclude (matches nothing by default).

**Parameters**

- **include** – Inclusion filter (translated by `translate_string_filter()` with syntax specified by *syntax*); include all by default.
- **exclude** – Exclusion filter (translated by `translate_string_filter()` with syntax specified by *syntax*); exclude none by default.

- **syntax** – Default syntax for pattern filters. Can be 're' (re), 'glob' (glob) or 'pred' (simply matching predicate).
- **match\_case** (*bool*) – Determines whether filter ignores case when matching.

`pylablib.core.utils.string.get_string_filter` (*include=None, exclude=None, syntax='re', match\_case=False*)

Generate *StringFilter* with the given parameters.

If the first argument is already *StringFilter*, return as is. If it's a tuple, expand as argument list.

`pylablib.core.utils.string.sfglob` (*include=None, exclude=None*)

Return string filter based on *glob* syntax

`pylablib.core.utils.string.sfregex` (*include=None, exclude=None, match\_case=False*)

Return string filter based on *re* syntax

`pylablib.core.utils.string.filter_string_list` (*l, filt*)

Filter string list based on the filter

`pylablib.core.utils.string.escape_string` (*value, location='element', escape\_convertible=True, quote\_type=''*)

Escape string.

Escaping can be partially skipped depending on *location*:

- **"parameter"**: escape only if it contains hard delimiters ("*\n\t\v\r*") anywhere or *\_border\_escaped* ("*'*" or space) on the sides (suited for parameters taking the full string);
- **"entry"**: same as above, plus containing soft delimiters (*,* or space) anywhere (suited for entries of a table);
- **"element"**: always escaped

If *escape\_convertible==True*, escape strings which can be misinterpreted as other values, such as "*1*" or "*[]*"; otherwise, escape only strings which contain special characters.

If *quote\_type* is not *None*, automatically put the string into the specified quotation marks; if *quote\_type* is *None*, all quotation marks are escaped; if it's not *None*, only *quote\_type* marks are escaped.

**class** `pylablib.core.utils.string.TConversionClass` (*label, cls, rep, conv*)

Bases: *tuple*

**cls**

**conv**

**count** ()

Return number of occurrences of value.

**index** ()

Return first index of value.

Raises *ValueError* if the value is not present.

**label**

**rep**

`pylablib.core.utils.string.add_conversion_class` (*label, cls, rep, conv*)

Add a string conversion class.

Some values (e.g., numpy arrays or named tuples) lose some of their associated information when converted into strings. With this function is possible to define custom conversion rules for such classes.

**Parameters**

- **label** (*str*) – class label (e.g., "array")
- **cls** – class which is used to determine if the value should use this conversion functions (e.g., `np.ndarray`)
- **rep** – function which takes a single argument (object of class *cls*) and returns its representations; can return a string or an object which is easier to convert to a string (e.g., a list or a tuple)
- **conv** – function which takes one or several arguments (converted values of the class representation) and returns the corresponding object; if *rep* returns a tuple, treat it as a list of several arguments, which are passed to *conv* separately; otherwise, *conv* gets a single argument which is the result of *rep*

When converting to string, if an object of class *cls* is encountered, it is converted in a string `label(str_rep)` (e.g., `"array([0, 1, 2])"`), where *str\_rep* is the result of calling *rep* (if this result is a tuple, avoid double parentheses, e.g., if the result is a tuple `(1, 2)`, the string becomes `"label(1, 2)"` instead of `"label((1, 2))"`). When converting from string, the values inside the parentheses are passed as arguments to *conv* function to get the resulting value.

`pylablib.core.utils.string.add_namedtuple_class(cls)`

Add conversion class for a given named tuple class.

For details, see `add_conversion_class()`.

`pylablib.core.utils.string.to_string(value, location='element', value_formats=None, parenthesis_rules='text', use_classes=False)`

Convert value to string with an option of modifying format string.

#### Parameters

- **value** –
- **location** (*str*) – Used for converting strings (see `escape_string()`).
- **value\_formats** (*dict*) – dictionary {*value\_type*: *fmt*}, where *value\_type* can be `int`, `float` or `complex` and *fmt* is a format string used to represent value of this type (e.g., `"5.3f"`); default formats are {`float`: `".12E"`, `complex`: `".12E"`, `int`: `"d"`}.
- **parenthesis\_rules** (*str*) – determine how to deal with single-element tuples and complex numbers can be `"text"` (single-element tuples are represented with simple parentheses, e.g., `"(1)"`; complex number are represented without parentheses, e.g., `"1+2j"`) or `"python"` (single-element tuples are represented with a comma in the end, e.g., `"(1,)"`; complex number are represented with parentheses, e.g., `"(1+2j)"`)
- **use\_classes** (*bool*) – if `True`, use additional representation classes for special objects (e.g., numpy arrays will be represented as `"array([1, 2, 3])"` instead of just `"[1, 2, 3]"`). This improves conversion fidelity, but makes result harder to parse (e.g., by external string parsers). See `add_conversion_class()` for more explanation.

`pylablib.core.utils.string.is_convertible(value)`

Check if the value can be converted to a string using standard `to_string()` function.

`pylablib.core.utils.string.extract_escaped_string(line, start=0)`

Extract escaped string in quotation marks from the *line*, starting from *start*.

*line[start]* should be a quotation mark (`'` or `"`) or `r` or `b` followed by a quotation mark (for raw or binary strings).

**Returns** tuple (end position, un-escaped string).

```
pylablib.core.utils.string.unescape_string(value)
```

Un-escape string.

Only attempt if the string starts a quotation mark " or '. Otherwise (including strings like 'r"' or 'b"'), return the string as is. Raise an error if the string starts with a quotation mark, but does not correspond to a proper escaped string (e.g., '"abc or '"abc"def).

```
pylablib.core.utils.string.to_range(range_tuple)
```

```
pylablib.core.utils.string.from_string(value, case_sensitive=True, parenthesis_rules='text', use_classes=True)
```

Parse a string.

Recognizes integers, floats, complex numbers (with i or j for complex part), strings (in quotation marks), dicts, sets, list and tuples, booleans and None. If item is unrecognizable, assumed to be a string.

#### Parameters

- **case\_sensitive** (*bool*) – applied when compared to None, True or False.
- **parenthesis\_rules** (*str*) – determines how to deal with empty entries (e.g., [1, 3]) and complex number representation ("1+2j" vs. "(1+2j)"):
  - 'text': any empty entries are translated into empty\_string (i.e., [, ] -> [empty\_string, empty\_string]), except for completely empty structures ([ ] or ( )); complex numbers are represented without parentheses, so that "(1+2j)" will be interpreted as a single-element tuple (1+2j, ).
  - 'python': empty entries in the middle are not allowed; empty entries at the end are ignored (i.e., [2, ] -> [2]) (single-element tuple can still be expressed in two ways: (e, ) or (e)); complex numbers are by default represented with parentheses, so that "(1+2j)" will be interpreted as a complex number, and only (1+2j, ), ((1+2j)) or ((1+2j), ) as a single-element tuple.
- **use\_classes** (*bool*) – if True, use additional representation classes for special objects (e.g., "array([1, 2, 3])" will be converted into a numpy array instead of raising an error). See [add\\_conversion\\_class\(\)](#) for more explanation.

```
pylablib.core.utils.string.from_string_partial(value, delimiters=re.compile('\s*,\s*\|s+'), case_sensitive=True, parenthesis_rules='text', use_classes=True, return_string=False)
```

Convert the first part of the supplied string (bounded by *delimiters*) into a value.

*delimiters* is a string or a regexp (default is "\s\*,\s\*\|s+", i.e., comma or spaces). If *return\_string*==False, convert the value string and return tuple (end\_position, converted\_value); otherwise, return tuple (end\_position, value\_string).

The rest of the parameters is the same as in [from\\_string\(\)](#).

```
pylablib.core.utils.string.from_row_string(value, delimiters=re.compile('\s*,\s*\|s+'), case_sensitive=True, parenthesis_rules='text', use_classes=True, return_string=False)
```

Convert the row string into a list of values, separated by delimiters.

If *return\_string*==False, return list of converted objects; otherwise, return list of unconverted strings.

The rest of the parameters is the same as in [from\\_string\\_partial\(\)](#).

## pylablib.core.utils.strpack module

Utilities for packing values into bitstrings. Small extension of the struct module.

`pylablib.core.utils.strpack.int2bytes (val, l, bo='>')`

Convert integer into a list of bytes of length *l*.

*bo* determines byte order: '>' is big-endian (MSB first), '<' is little-endian (LSB first).

`pylablib.core.utils.strpack.bytes2int (val, bo='>')`

Convert a list of bytes into an integer.

*bo* determines byte order: '>' is big-endian (MSB first), '<' is little-endian (LSB first).

`pylablib.core.utils.strpack.int2bits (val, l, bo='>')`

Convert integer into a list of bits of length *l*.

*bo* determines byte (and bit) order: '>' is big-endian (MSB first), '<' is little-endian (LSB first).

`pylablib.core.utils.strpack.bits2int (val, bo='>')`

Convert a list of bits into an integer.

*bo* determines byte (and bit) order: '>' is big-endian (MSB first), '<' is little-endian (LSB first).

`pylablib.core.utils.strpack.pack_uint (val, l, bo='>')`

Convert an unsigned integer into a bytestring of length *l*.

Return bytes object. *bo* determines byte order: '>' is big-endian (MSB first), '<' is little-endian (LSB first).

`pylablib.core.utils.strpack.pack_int (val, l, bo='>')`

Convert a signed integer into a bytestring of length *l*.

Return bytes object. *bo* determines byte order: '>' is big-endian (MSB first), '<' is little-endian (LSB first).

`pylablib.core.utils.strpack.unpack_uint (msg, bo='>')`

Convert a bytestring into an unsigned integer.

*bo* determines byte order: '>' is big-endian (MSB first), '<' is little-endian (LSB first).

`pylablib.core.utils.strpack.unpack_int (msg, bo='>')`

Convert a bytestring into a signed integer.

*bo* determines byte order: '>' is big-endian (MSB first), '<' is little-endian (LSB first).

`pylablib.core.utils.strpack.unpack_numpy_ul2bit (buffer, byteorder='<', count=-1)`

## pylablib.core.utils.units module

Routines for conversion of physical units.

`pylablib.core.utils.units.split_units (value)`

Split string value with a dimension.

Return tuple (*val*, *unit*), where *val* is the float part of the value, and *unit* is the string representing units.

`pylablib.core.utils.units.convert_length_units (value, value_unit='m', result_unit='m', case_sensitive=True)`

Convert *value* from *value\_unit* to *result\_unit*.

The possible length units are 'm', 'mm', 'um', 'nm', 'pm', 'fm'. If *case\_sensitive*==True, matching units is case sensitive.

```
pylablib.core.utils.units.convert_time_units(value, value_unit='s', result_unit='s',
                                              case_sensitive=True)
```

Convert *value* from *value\_unit* to *result\_unit*.

The possible time units are 's', 'ms', 'us', 'ns', 'ps', 'fs', 'as'. If *case\_sensitive*==True, matching units is case sensitive.

```
pylablib.core.utils.units.convert_frequency_units(value, value_unit='Hz',
                                                  result_unit='Hz',
                                                  case_sensitive=True)
```

Convert *value* from *value\_unit* to *result\_unit*.

The possible frequency units are 'Hz', 'kHz', 'MHz', 'GHz'. If *case\_sensitive*==True, matching units is case sensitive.

```
pylablib.core.utils.units.convert_power_units(value, value_unit='dBm', re-
                                              sult_unit='dBm', case_sensitive=True,
                                              impedance=50.0)
```

Convert *value* from *value\_unit* to *result\_unit*.

For conversion between voltage and power, assume RMS voltage and the given *impedance*. The possible power units are 'dBm', 'dBmV', 'dBuV', 'W', 'mW', 'uW', 'nW', 'mV', 'nV'. If *case\_sensitive*==True, matching units is case sensitive.

## Module contents

### Module contents

#### pylablib.devices package

#### Subpackages

#### pylablib.devices.AWG package

#### Submodules

#### pylablib.devices.AWG.generic module

**exception** pylablib.devices.AWG.generic.GenericAWGError

Bases: `pylablib.core.devio.base.DeviceError`

Generic AWG error

**args**

**with\_traceback()**

Exception.with\_traceback(tb) – set self.\_\_traceback\_\_ to tb and return self.

**exception** pylablib.devices.AWG.generic.GenericAWGBackendError(*exc*)

Bases: `pylablib.devices.AWG.generic.GenericAWGError`, `pylablib.core.devio.comm_backend.DeviceBackendError`

AWG backend communication error

**args**

**with\_traceback()**

Exception.with\_traceback(tb) – set self.\_\_traceback\_\_ to tb and return self.

```

class pylablib.devices.AWG.generic.GenericAWG(addr)
    Bases: pylablib.core.devio.SCPISCPIDevice

    Generic arbitrary wave generator, based on Agilent 33500.

    With slight modifications works for many other AWGs using largely the same syntax.

    Error
        alias of GenericAWGError

    ReraiseError
        alias of GenericAWGBackendError

    get_channels_number()
        Get the number of channels

    get_current_channel()
        Get current channel

    select_current_channel(channel)
        Select current default channel

    is_output_enabled(channel=None)
        Check if the output is enabled

    enable_output(enabled=True, channel=None)
        Turn the output on or off

    get_output_polarity(channel=None)
        Get output polarity.

        Can be either "norm" or "inv".

    set_output_polarity(polarity='norm', channel=None)
        Set output polarity.

        Can be either "norm" or "inv".

    is_sync_output_enabled(channel=None)
        Check if SYNC output is enabled

    enable_sync_output(enabled=True, channel=None)
        Enable or disable SYNC output

    get_load(channel=None)
        Get the output load

    set_load(load=None, channel=None)
        Set the output load (None means High-Z)

    get_function(channel=None)
        Get output function.

        Can be one of the following: "sine", "square", "ramp", "pulse", "noise", "prbs", "DC",
        "user", "arb". Not all functions can be available, depending on the particular model of the generator.

    set_function(func, channel=None)
        Set output function.

        Can be one of the following: "sine", "square", "ramp", "pulse", "noise", "prbs", "DC",
        "user", "arb". Not all functions can be available, depending on the particular model of the generator.

    get_amplitude(channel=None)
        Get output amplitude (i.e., half of the span)

```

**set\_amplitude** (*amplitude*, *channel=None*)  
Set output amplitude (i.e., half of the span)

**get\_offset** (*channel=None*)  
Get output offset

**set\_offset** (*offset*, *channel=None*)  
Set output offset

**get\_output\_range** (*channel=None*)  
Get output voltage range.  
  
Return tuple (*vmin*, *vmax*) with the low and high voltage values (i.e., *offset-amplitude* and *offset+amplitude*).

**set\_output\_range** (*rng*, *channel=None*)  
Set output voltage range.  
  
If span is less than  $1\text{E}-4$ , automatically switch to DC mode.

**get\_frequency** (*channel=None*)  
Get output frequency

**set\_frequency** (*frequency*, *channel=None*)  
Set output frequency

**get\_phase** (*channel=None*)  
Get output phase (in degrees)

**set\_phase** (*phase*, *channel=None*)  
Set output phase (in degrees)

**sync\_phase** ()  
Synchronize phase between two channels

**get\_duty\_cycle** (*channel=None*)  
Get output duty cycle (in percent).  
  
Only applies to "square" output function.

**set\_duty\_cycle** (*dcycle*, *channel=None*)  
Set output duty cycle (in percent).  
  
Only applies to "square" output function.

**get\_ramp\_symmetry** (*channel=None*)  
Get output ramp symmetry (in percent).  
  
Only applies to "ramp" output function.

**set\_ramp\_symmetry** (*rsymm*, *channel=None*)  
Set output ramp symmetry (in percent).  
  
Only applies to "ramp" output function.

**get\_pulse\_width** (*channel=None*)  
Get output pulse width (in seconds).  
  
Only applies to "pulse" output function.

**set\_pulse\_width** (*width*, *channel=None*)  
Set output pulse width (in seconds).  
  
Only applies to "pulse" output function.



**is\_burst\_enabled** (*channel=None*)  
Check if the burst mode is enabled

**enable\_burst** (*enabled=True, channel=None*)  
Enable burst mode

**get\_burst\_mode** (*channel=None*)  
Get burst mode.  
  
Can be either "trig" or "gate".

**set\_burst\_mode** (*mode, channel=None*)  
Set burst mode.  
  
Can be either "trig" or "gate".

**get\_burst\_ncycles** (*channel=None*)  
Get burst mode ncycles.  
  
Infinite corresponds to a large value (>1E37).

**set\_burst\_ncycles** (*ncycles=1, channel=None*)  
Set burst mode ncycles.  
  
Infinite corresponds to None

**get\_gate\_polarity** (*channel=None*)  
Get burst gate polarity.  
  
Can be either "norm" or "inv".

**set\_gate\_polarity** (*polarity='norm', channel=None*)  
Set burst gate polarity.  
  
Can be either "norm" or "inv".

**get\_trigger\_source** (*channel=None*)  
Get trigger source.  
  
Can be either "imm", "ext", or "bus".

**set\_trigger\_source** (*src, channel=None*)  
Set trigger source.  
  
Can be either "imm", "ext", or "bus".

**get\_trigger\_slope** (*channel=None*)  
Get trigger slope.  
  
Can be either "pos", or "neg".

**set\_trigger\_slope** (*slope, channel=None*)  
Set trigger slope.  
  
Can be either "pos", or "neg".

**is\_trigger\_output\_enabled** (*channel=None*)  
Check if the trigger output is enabled

**enable\_trigger\_output** (*enabled=True, channel=None*)  
Enable trigger output

**get\_output\_trigger\_slope** (*channel=None*)  
Get output trigger slope.  
  
Can be either "pos", or "neg".

**BackendError**

alias of `pylablib.core.devio.comm_backend.DeviceBackendError`

**class NoParameterCaller** (*device, kind*)

Bases: `object`

Class to simplify calling functions without a parameter

**apply\_settings** (*settings*)

Apply the settings.

*settings* is a dict {name: value} of the available device settings. Non-applicable settings are ignored.

**ask** (*msg, data\_type='string', delay=0.0, timeout=None, read\_echo=False*)

Write a message and read a reply.

*msg* is the query message, *delay* is the delay between write and read. Other parameters are the same as in `read()`. If `read_echo==True`, assume that the device first echoes the input and skip it.

**close** ()

Close the backend

**flush** (*one\_line=False*)

Flush the read buffer (read all the available data and return the number of bytes read).

If `one_line==True`, read only a single line.

**static get\_arg\_type** (*arg*)

Autodetect argument type

**get\_device\_variable** (*key*)

Get the value of a settings, status, or full info parameter

**get\_esr** (*timeout=None*)

Get the device status register (by default, "\*ESR?" command)

**get\_full\_info** (*include=0*)

Get dict {name: value} containing full device information (including status and settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting `include=-10` queries all available variables, which is equivalent to `include="all"`.

**get\_full\_status** (*include=0*)

Get dict {name: value} containing the device status (including settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting `include=-10` queries all available variables, which is equivalent to `include="all"`.

**get\_id** (*timeout=None*)

Get the device IDN. (query SCPI "\*IDN?" command)

**get\_settings** (*include=0*)

Get dict {name: value} containing all the device settings.

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting `include=-10` queries all available variables, which is equivalent to `include="all"`.

**is\_opened()**  
Check if the device is connected

**lock** (*timeout=None*)  
Lock the access to the device from other threads/processes (isn't necessarily implemented)

**locking** (*timeout=None*)  
Context manager for lock & unlock

**open()**  
Open the backend

**static parse\_array\_data** (*data, fmt, include\_header=False*)  
Parse the data returned by the device. *fmt* is *DataFormat* description in numpy format (e.g., "<u2").  
  
If *include\_header==True*, the data is assumed to be in a (somewhat) standard SCPI format: b'#', then a single digit *s* denoting length of the size block, then *s* digits denoting length of the data (in bytes) followed by the actual data. Otherwise (*include\_header==False*), assume that the header is already removed.

**read** (*data\_type='string', timeout=None*)  
Read data from the device.  
  
*data\_type* determines the type of the data. Can be 'raw' (just raw data), 'string' (with trailing and leading spaces stripped), 'int', 'float', 'bool' (interprets 0 or 'off' as False, anything else as True), 'value' (returns tuple (value, unit), where *value* is float), a callable (return the result of this callable applied to the string value), a dictionary (return the stored value corresponding to the string value, or to the value converted into integer if the string value is not present), or a list of data types (the result is treated as a list of values with the given types separated by commas). *timeout* overrides the default value.

**read\_binary\_array\_data** (*include\_header=False, timeout=None, flush\_term=True*)  
Read a binary data in the from the device.  
  
The data assumes the standard binary transfer header consisting of "#" symbol, then a single digit with the size of the length string, then the length string containing the length of the binary data (in bytes). If *include\_header==True*, return the data with the header; otherwise, return only the content. If *flush\_term==True*, flush the following line to skip terminator characters after the binary data, which are added by some devices. *timeout* overrides the default value.

**reconnect** (*new\_instrument=True, ignore\_error=True*)  
Remake the connection.  
  
If *new\_instrument==True*, create a new backend instance. If *ignore\_error==True*, ignore errors on closing.

**reset()**  
Reset the device (by default, "\*RST" command)

**set\_device\_variable** (*key, value*)  
Set the value of a settings parameter

**set\_output\_trigger\_slope** (*slope, channel=None*)  
Set output trigger slope.  
  
Can be either "pos", or "neg".

**sleep** (*delay*)  
Wait for *delay* seconds

**unlock()**  
Unlock the access to the device from other threads/processes (isn't necessarily implemented)

**using\_write\_buffer()**

Context manager for using a write buffer.

While it's active, all the consecutive `write()` operations are bundled together with `;` delimiter. The actual write is performed at the `read()/ask()` operation or at the end of the block.

**wait** (*wait\_type='sync', timeout=None, wait\_callback=None*)

Pause execution until device overlapped commands are complete.

*wait\_type* is either `'sync'` (perform `wait_sync()`), `'dev'` (perform `wait_dev()`) or `'none'` (do nothing).

**wait\_dev()**

Pause execution of the device commands until device overlapped commands (e.g., taking sweeps) are complete.

Note that the code execution is not paused.

**wait\_sync** (*timeout=None, wait\_callback=None*)

Pause execution of the script until device overlapped commands (e.g., taking sweeps) are complete.

*timeout* and *wait\_callback* override default constructor parameters.

**write** (*msg, arg=None, arg\_type=None, unit=None, bool\_selector=None, wait\_sync=None, read\_echo=False, read\_echo\_delay=0.0*)

Send a command.

**Parameters**

- **msg** (*str*) – Text message.
- **arg** – Optional argument to append in the end. If a list of arguments is supplied, the result is joined with `,`.
- **arg\_type** (*str*) – Argument type. Can be `'raw'` (in which case data is sent raw), `'string'`, `'int'`, `'float'`, `'bool'`, a format string (such as `'{: .3f}'`) or a list of argument types (for an iterable argument); if format string is used and the argument is a list or a tuple, then it is expanded as a list of arguments (e.g., `arg_type='{0};{1}'` with `arg=[1, 2]` will produce a string `'1; 2'`); if a list of types is used, each element of *arg* is converted using the corresponding type, and the result is joined with `,`.
- **unit** (*str*) – If not `None`, use it as a unit to append after the value.
- **bool\_selector** (*tuple*) – A tuple (*false\_value*, *true\_value*) of two strings to represent bool argument; by default, use `._bool_selector` attribute.
- **wait\_sync** – if `True`, append the sync command (specified as `._wait_sync_comm` attribute, `"*OPC?"` by default) after the message and pause the execution command is complete; useful in long set operations, where the device might ignore later inputs until the current command is complete; if `None`, use the class default `._default_write_sync` attribute (`False` by default).
- **read\_echo** (*bool*) – If `True`, read a single line after write.
- **read\_echo\_delay** (*float*) – The delay between write and read if `read_echo==True`.

**pylablib.devices.AWG.specific module**

```
class pylablib.devices.AWG.specific.Agilent33500(addr, channels_number='auto')
```

Bases: `pylablib.devices.AWG.generic.GenericAWG`

Agilent 33500 AWG.

**Parameters** `channels_number` – number of channels; if "auto", try to determine automatically (by certain commands causing errors)

**BackendError**

alias of `pylablib.core.devio.comm_backend.DeviceBackendError`

**Error**

alias of `pylablib.devices.AWG.generic.GenericAWGError`

**class NoParameterCaller** (*device, kind*)

Bases: `object`

Class to simplify calling functions without a parameter

**ReraiseError**

alias of `pylablib.devices.AWG.generic.GenericAWGBackendError`

**apply\_settings** (*settings*)

Apply the settings.

*settings* is a dict {name: value} of the available device settings. Non-applicable settings are ignored.

**ask** (*msg, data\_type='string', delay=0.0, timeout=None, read\_echo=False*)

Write a message and read a reply.

*msg* is the query message, *delay* is the delay between write and read. Other parameters are the same as in `read()`. If `read_echo==True`, assume that the device first echoes the input and skip it.

**close** ()

Close the backend

**enable\_burst** (*enabled=True, channel=None*)

Enable burst mode

**enable\_output** (*enabled=True, channel=None*)

Turn the output on or off

**enable\_sync\_output** (*enabled=True, channel=None*)

Enable or disable SYNC output

**enable\_trigger\_output** (*enabled=True, channel=None*)

Enable trigger output

**flush** (*one\_line=False*)

Flush the read buffer (read all the available data and return the number of bytes read).

If `one_line==True`, read only a single line.

**get\_amplitude** (*channel=None*)

Get output amplitude (i.e., half of the span)

**static get\_arg\_type** (*arg*)

Autodetect argument type

**get\_burst\_mode** (*channel=None*)

Get burst mode.

Can be either "trig" or "gate".

**get\_burst\_ncycles** (*channel=None*)

Get burst mode ncycles.

Infinite corresponds to a large value (>1E37).

**get\_channels\_number** ()

Get the number of channels

**get\_current\_channel** ()

Get current channel

**get\_device\_variable** (*key*)

Get the value of a settings, status, or full info parameter

**get\_duty\_cycle** (*channel=None*)

Get output duty cycle (in percent).

Only applies to "square" output function.

**get\_esr** (*timeout=None*)

Get the device status register (by default, "\*ESR?" command)

**get\_frequency** (*channel=None*)

Get output frequency

**get\_full\_info** (*include=0*)

Get dict {name: value} containing full device information (including status and settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_full\_status** (*include=0*)

Get dict {name: value} containing the device status (including settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_function** (*channel=None*)

Get output function.

Can be one of the following: "sine", "square", "ramp", "pulse", "noise", "prbs", "DC", "user", "arb". Not all functions can be available, depending on the particular model of the generator.

**get\_gate\_polarity** (*channel=None*)

Get burst gate polarity.

Can be either "norm" or "inv".

**get\_id** (*timeout=None*)

Get the device IDN. (query SCPI '\*IDN?' command)

**get\_load** (*channel=None*)

Get the output load

**get\_offset** (*channel=None*)

Get output offset

**get\_output\_polarity** (*channel=None*)

Get output polarity.

Can be either "norm" or "inv".

**get\_output\_range** (*channel=None*)

Get output voltage range.

Return tuple (*vmin*, *vmax*) with the low and high voltage values (i.e., *offset-amplitude* and *offset+amplitude*).

**get\_output\_trigger\_slope** (*channel=None*)

Get output trigger slope.

Can be either "pos", or "neg".

**get\_phase** (*channel=None*)

Get output phase (in degrees)

**get\_pulse\_width** (*channel=None*)

Get output pulse width (in seconds).

Only applies to "pulse" output function.

**get\_ramp\_symmetry** (*channel=None*)

Get output ramp symmetry (in percent).

Only applies to "ramp" output function.

**get\_settings** (*include=0*)

Get dict {*name*: *value*} containing all the device settings.

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_trigger\_slope** (*channel=None*)

Get trigger slope.

Can be either "pos", or "neg".

**get\_trigger\_source** (*channel=None*)

Get trigger source.

Can be either "imm", "ext", or "bus".

**is\_burst\_enabled** (*channel=None*)

Check if the burst mode is enabled

**is\_opened** ()

Check if the device is connected

**is\_output\_enabled** (*channel=None*)

Check if the output is enabled

**is\_sync\_output\_enabled** (*channel=None*)

Check if SYNC output is enabled

**is\_trigger\_output\_enabled** (*channel=None*)

Check if the trigger output is enabled

**lock** (*timeout=None*)

Lock the access to the device from other threads/processes (isn't necessarily implemented)

**locking** (*timeout=None*)

Context manager for lock & unlock

**open** ()

Open the backend

**static parse\_array\_data** (*data, fmt, include\_header=False*)

Parse the data returned by the device. *fmt* is *DataFormat* description in numpy format (e.g., "<u2").

If *include\_header==True*, the data is assumed to be in a (somewhat) standard SCPI format: *b' # '*, then a single digit *s* denoting length of the size block, then *s* digits denoting length of the data (in bytes) followed by the actual data. Otherwise (*include\_header==False*), assume that the header is already removed.

**read** (*data\_type='string', timeout=None*)

Read data from the device.

*data\_type* determines the type of the data. Can be 'raw' (just raw data), 'string' (with trailing and leading spaces stripped), 'int', 'float', 'bool' (interprets 0 or 'off' as False, anything else as True), 'value' (returns tuple (value, unit), where *value* is float), a callable (return the result of this callable applied to the string value), a dictionary (return the stored value corresponding to the string value, or to the value converted into integer if the string value is not present), or a list of data types (the result is treated as a list of values with the given types separated by commas). *timeout* overrides the default value.

**read\_binary\_array\_data** (*include\_header=False, timeout=None, flush\_term=True*)

Read a binary data in the from the device.

The data assumes the standard binary transfer header consisting of "# " symbol, then a single digit with the size of the length string, then the length string containing the length of the binary data (in bytes). If *include\_header==True*, return the data with the header; otherwise, return only the content. If *flush\_term==True*, flush the following line to skip terminator characters after the binary data, which are added by some devices. *timeout* overrides the default value.

**reconnect** (*new\_instrument=True, ignore\_error=True*)

Remake the connection.

If *new\_instrument==True*, create a new backend instance. If *ignore\_error==True*, ignore errors on closing.

**reset** ()

Reset the device (by default, "\*RST" command)

**select\_current\_channel** (*channel*)

Select current default channel

**set\_amplitude** (*amplitude, channel=None*)

Set output amplitude (i.e., half of the span)

**set\_burst\_mode** (*mode, channel=None*)

Set burst mode.

Can be either "trig" or "gate".

**set\_burst\_ncycles** (*ncycles=1, channel=None*)

Set burst mode ncycles.

Infinite corresponds to None

**set\_device\_variable** (*key, value*)

Set the value of a settings parameter

**set\_duty\_cycle** (*dcycle, channel=None*)

Set output duty cycle (in percent).

Only applies to "square" output function.

**set\_frequency** (*frequency, channel=None*)

Set output frequency



**set\_function** (*func, channel=None*)

Set output function.

Can be one of the following: "sine", "square", "ramp", "pulse", "noise", "prbs", "DC", "user", "arb". Not all functions can be available, depending on the particular model of the generator.

**set\_gate\_polarity** (*polarity='norm', channel=None*)

Set burst gate polarity.

Can be either "norm" or "inv".

**set\_load** (*load=None, channel=None*)

Set the output load (None means High-Z)

**set\_offset** (*offset, channel=None*)

Set output offset

**set\_output\_polarity** (*polarity='norm', channel=None*)

Set output polarity.

Can be either "norm" or "inv".

**set\_output\_range** (*rng, channel=None*)

Set output voltage range.

If span is less than  $1\text{E}-4$ , automatically switch to DC mode.

**set\_output\_trigger\_slope** (*slope, channel=None*)

Set output trigger slope.

Can be either "pos", or "neg".

**set\_phase** (*phase, channel=None*)

Set output phase (in degrees)

**set\_pulse\_width** (*width, channel=None*)

Set output pulse width (in seconds).

Only applies to "pulse" output function.

**set\_ramp\_symmetry** (*rsymm, channel=None*)

Set output ramp symmetry (in percent).

Only applies to "ramp" output function.

**set\_trigger\_slope** (*slope, channel=None*)

Set trigger slope.

Can be either "pos", or "neg".

**set\_trigger\_source** (*src, channel=None*)

Set trigger source.

Can be either "imm", "ext", or "bus".

**sleep** (*delay*)

Wait for *delay* seconds

**sync\_phase** ()

Synchronize phase between two channels

**unlock** ()

Unlock the access to the device from other threads/processes (isn't necessarily implemented)

**using\_write\_buffer()**

Context manager for using a write buffer.

While it's active, all the consecutive `write()` operations are bundled together with `;` delimiter. The actual write is performed at the `read()/ask()` operation or at the end of the block.

**wait** (*wait\_type='sync', timeout=None, wait\_callback=None*)

Pause execution until device overlapped commands are complete.

*wait\_type* is either `'sync'` (perform `wait_sync()`), `'dev'` (perform `wait_dev()`) or `'none'` (do nothing).

**wait\_dev()**

Pause execution of the device commands until device overlapped commands (e.g., taking sweeps) are complete.

Note that the code execution is not paused.

**wait\_sync** (*timeout=None, wait\_callback=None*)

Pause execution of the script until device overlapped commands (e.g., taking sweeps) are complete.

*timeout* and *wait\_callback* override default constructor parameters.

**write** (*msg, arg=None, arg\_type=None, unit=None, bool\_selector=None, wait\_sync=None, read\_echo=False, read\_echo\_delay=0.0*)

Send a command.

#### Parameters

- **msg** (*str*) – Text message.
- **arg** – Optional argument to append in the end. If a list of arguments is supplied, the result is joined with `" , "`.
- **arg\_type** (*str*) – Argument type. Can be `'raw'` (in which case data is sent raw), `'string'`, `'int'`, `'float'`, `'bool'`, a format string (such as `'{: .3f}'`) or a list of argument types (for an iterable argument); if format string is used and the argument is a list or a tuple, then it is expanded as a list of arguments (e.g., `arg_type='{0};{1}'` with `arg=[1, 2]` will produce a string `'1; 2'`); if a list of types is used, each element of *arg* is converted using the corresponding type, and the result is joined with `" , "`.
- **unit** (*str*) – If not `None`, use it as a unit to append after the value.
- **bool\_selector** (*tuple*) – A tuple (*false\_value*, *true\_value*) of two strings to represent bool argument; by default, use `._bool_selector` attribute.
- **wait\_sync** – if `True`, append the sync command (specified as `._wait_sync_comm` attribute, `"*OPC?"` by default) after the message and pause the execution command is complete; useful in long set operations, where the device might ignore later inputs until the current command is complete; if `None`, use the class default `._default_write_sync` attribute (`False` by default).
- **read\_echo** (*bool*) – If `True`, read a single line after write.
- **read\_echo\_delay** (*float*) – The delay between write and read if `read_echo==True`.

**class** `pylablib.devices.AWG.specific.Agilent33220A(addr)`

Bases: `pylablib.devices.AWG.generic.GenericAWG`

Agilent 33220A AWG.

**BackendError**

alias of `pylablib.core.devio.comm_backend.DeviceBackendError`

**Error**

alias of `pylablib.devices.AWG.generic.GenericAWGError`

**class NoParameterCaller (device, kind)**

Bases: `object`

Class to simplify calling functions without a parameter

**ReraiseError**

alias of `pylablib.devices.AWG.generic.GenericAWGBackendError`

**apply\_settings (settings)**

Apply the settings.

*settings* is a dict {name: value} of the available device settings. Non-applicable settings are ignored.

**ask (msg, data\_type='string', delay=0.0, timeout=None, read\_echo=False)**

Write a message and read a reply.

*msg* is the query message, *delay* is the delay between write and read. Other parameters are the same as in `read()`. If `read_echo==True`, assume that the device first echoes the input and skip it.

**close ()**

Close the backend

**enable\_burst (enabled=True, channel=None)**

Enable burst mode

**enable\_output (enabled=True, channel=None)**

Turn the output on or off

**enable\_sync\_output (enabled=True, channel=None)**

Enable or disable SYNC output

**enable\_trigger\_output (enabled=True, channel=None)**

Enable trigger output

**flush (one\_line=False)**

Flush the read buffer (read all the available data and return the number of bytes read).

If `one_line==True`, read only a single line.

**get\_amplitude (channel=None)**

Get output amplitude (i.e., half of the span)

**static get\_arg\_type (arg)**

Autodetect argument type

**get\_burst\_mode (channel=None)**

Get burst mode.

Can be either "trig" or "gate".

**get\_burst\_ncycles (channel=None)**

Get burst mode ncycles.

Infinite corresponds to a large value (>1E37).

**get\_channels\_number ()**

Get the number of channels

**get\_current\_channel** ()

Get current channel

**get\_device\_variable** (*key*)

Get the value of a settings, status, or full info parameter

**get\_duty\_cycle** (*channel=None*)

Get output duty cycle (in percent).

Only applies to "square" output function.

**get\_esr** (*timeout=None*)

Get the device status register (by default, "\*ESR?" command)

**get\_frequency** (*channel=None*)

Get output frequency

**get\_full\_info** (*include=0*)

Get dict {name: value} containing full device information (including status and settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_full\_status** (*include=0*)

Get dict {name: value} containing the device status (including settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_function** (*channel=None*)

Get output function.

Can be one of the following: "sine", "square", "ramp", "pulse", "noise", "prbs", "DC", "user", "arb". Not all functions can be available, depending on the particular model of the generator.

**get\_gate\_polarity** (*channel=None*)

Get burst gate polarity.

Can be either "norm" or "inv".

**get\_id** (*timeout=None*)

Get the device IDN. (query SCPI '\*IDN?' command)

**get\_load** (*channel=None*)

Get the output load

**get\_offset** (*channel=None*)

Get output offset

**get\_output\_polarity** (*channel=None*)

Get output polarity.

Can be either "norm" or "inv".

**get\_output\_range** (*channel=None*)

Get output voltage range.

Return tuple (*vmin*, *vmax*) with the low and high voltage values (i.e., offset-amplitude and offset+amplitude).

**get\_output\_trigger\_slope** (*channel=None*)

Get output trigger slope.

Can be either "pos", or "neg".

**get\_phase** (*channel=None*)

Get output phase (in degrees)

**get\_pulse\_width** (*channel=None*)

Get output pulse width (in seconds).

Only applies to "pulse" output function.

**get\_ramp\_symmetry** (*channel=None*)

Get output ramp symmetry (in percent).

Only applies to "ramp" output function.

**get\_settings** (*include=0*)

Get dict {name: value} containing all the device settings.

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_trigger\_slope** (*channel=None*)

Get trigger slope.

Can be either "pos", or "neg".

**get\_trigger\_source** (*channel=None*)

Get trigger source.

Can be either "imm", "ext", or "bus".

**is\_burst\_enabled** (*channel=None*)

Check if the burst mode is enabled

**is\_opened** ()

Check if the device is connected

**is\_output\_enabled** (*channel=None*)

Check if the output is enabled

**is\_sync\_output\_enabled** (*channel=None*)

Check if SYNC output is enabled

**is\_trigger\_output\_enabled** (*channel=None*)

Check if the trigger output is enabled

**lock** (*timeout=None*)

Lock the access to the device from other threads/processes (isn't necessarily implemented)

**locking** (*timeout=None*)

Context manager for lock & unlock

**open** ()

Open the backend

**static parse\_array\_data** (*data, fmt, include\_header=False*)

Parse the data returned by the device. *fmt* is [DataFormat](#) description in numpy format (e.g., "<u2").

If *include\_header==True*, the data is assumed to be in a (somewhat) standard SCPI format: b'#', then a single digit *s* denoting length of the size block, then *s* digits denoting length of the data (in

bytes) followed by the actual data. Otherwise (`include_header==False`), assume that the header is already removed.

**read** (*data\_type='string', timeout=None*)

Read data from the device.

*data\_type* determines the type of the data. Can be 'raw' (just raw data), 'string' (with trailing and leading spaces stripped), 'int', 'float', 'bool' (interprets 0 or 'off' as False, anything else as True), 'value' (returns tuple (value, unit), where *value* is float), a callable (return the result of this callable applied to the string value), a dictionary (return the stored value corresponding to the string value, or to the value converted into integer if the string value is not present), or a list of data types (the result is treated as a list of values with the given types separated by commas). *timeout* overrides the default value.

**read\_binary\_array\_data** (*include\_header=False, timeout=None, flush\_term=True*)

Read a binary data in the from the device.

The data assumes the standard binary transfer header consisting of "#" symbol, then a single digit with the size of the length string, then the length string containing the length of the binary data (in bytes). If `include_header==True`, return the data with the header; otherwise, return only the content. If `flush_term==True`, flush the following line to skip terminator characters after the binary data, which are added by some devices. *timeout* overrides the default value.

**reconnect** (*new\_instrument=True, ignore\_error=True*)

Remake the connection.

If `new_instrument==True`, create a new backend instance. If `ignore_error==True`, ignore errors on closing.

**reset** ()

Reset the device (by default, "\*RST" command)

**select\_current\_channel** (*channel*)

Select current default channel

**set\_amplitude** (*amplitude, channel=None*)

Set output amplitude (i.e., half of the span)

**set\_burst\_mode** (*mode, channel=None*)

Set burst mode.

Can be either "trig" or "gate".

**set\_burst\_ncycles** (*ncycles=1, channel=None*)

Set burst mode ncycles.

Infinite corresponds to None

**set\_device\_variable** (*key, value*)

Set the value of a settings parameter

**set\_duty\_cycle** (*dcycle, channel=None*)

Set output duty cycle (in percent).

Only applies to "square" output function.

**set\_frequency** (*frequency, channel=None*)

Set output frequency

**set\_function** (*func, channel=None*)

Set output function.

Can be one of the following: "sine", "square", "ramp", "pulse", "noise", "prbs", "DC", "user", "arb". Not all functions can be available, depending on the particular model of the generator.

**set\_gate\_polarity** (*polarity='norm', channel=None*)

Set burst gate polarity.

Can be either "norm" or "inv".

**set\_load** (*load=None, channel=None*)

Set the output load (None means High-Z)

**set\_offset** (*offset, channel=None*)

Set output offset

**set\_output\_polarity** (*polarity='norm', channel=None*)

Set output polarity.

Can be either "norm" or "inv".

**set\_output\_range** (*rng, channel=None*)

Set output voltage range.

If span is less than 1E-4, automatically switch to DC mode.

**set\_output\_trigger\_slope** (*slope, channel=None*)

Set output trigger slope.

Can be either "pos", or "neg".

**set\_phase** (*phase, channel=None*)

Set output phase (in degrees)

**set\_pulse\_width** (*width, channel=None*)

Set output pulse width (in seconds).

Only applies to "pulse" output function.

**set\_ramp\_symmetry** (*rsymm, channel=None*)

Set output ramp symmetry (in percent).

Only applies to "ramp" output function.

**set\_trigger\_slope** (*slope, channel=None*)

Set trigger slope.

Can be either "pos", or "neg".

**set\_trigger\_source** (*src, channel=None*)

Set trigger source.

Can be either "imm", "ext", or "bus".

**sleep** (*delay*)

Wait for *delay* seconds

**sync\_phase** ()

Synchronize phase between two channels

**unlock** ()

Unlock the access to the device from other threads/processes (isn't necessarily implemented)

**using\_write\_buffer** ()

Context manager for using a write buffer.

While it's active, all the consecutive `write()` operations are bundled together with ; delimiter. The actual write is performed at the `read()/ask()` operation or at the end of the block.

**wait** (*wait\_type='sync', timeout=None, wait\_callback=None*)

Pause execution until device overlapped commands are complete.

*wait\_type* is either 'sync' (perform `wait_sync()`), 'dev' (perform `wait_dev()`) or 'none' (do nothing).

**wait\_dev** ()

Pause execution of the device commands until device overlapped commands (e.g., taking sweeps) are complete.

Note that the code execution is not paused.

**wait\_sync** (*timeout=None, wait\_callback=None*)

Pause execution of the script until device overlapped commands (e.g., taking sweeps) are complete.

*timeout* and *wait\_callback* override default constructor parameters.

**write** (*msg, arg=None, arg\_type=None, unit=None, bool\_selector=None, wait\_sync=None, read\_echo=False, read\_echo\_delay=0.0*)

Send a command.

#### Parameters

- **msg** (*str*) – Text message.
- **arg** – Optional argument to append in the end. If a list of arguments is supplied, the result is joined with ", ".
- **arg\_type** (*str*) – Argument type. Can be 'raw' (in which case data is sent raw), 'string', 'int', 'float', 'bool', a format string (such as '{:.3f}') or a list of argument types (for an iterable argument); if format string is used and the argument is a list or a tuple, then it is expanded as a list of arguments (e.g., `arg_type='{0};{1}'` with `arg=[1, 2]` will produce a string '1; 2'); if a list of types is used, each element of *arg* is converted using the corresponding type, and the result is joined with ", ".
- **unit** (*str*) – If not None, use it as a unit to append after the value.
- **bool\_selector** (*tuple*) – A tuple (*false\_value*, *true\_value*) of two strings to represent bool argument; by default, use `._bool_selector` attribute.
- **wait\_sync** – if True, append the sync command (specified as `._wait_sync_comm` attribute, `"*OPC?"` by default) after the message and pause the execution command is complete; useful in long set operations, where the device might ignore later inputs until the current command is complete; if None, use the class default `._default_write_sync` attribute (False by default).
- **read\_echo** (*bool*) – If True, read a single line after write.
- **read\_echo\_delay** (*float*) – The delay between write and read if `read_echo==True`.

**class** `pylablib.devices.AWG.specific.InstekAFG2225` (*addr*)

Bases: `pylablib.devices.AWG.generic.GenericAWG`

Instek AFG2225 AWG.

Compared to 2000/2100 series, has one extra channel and a bit more capabilities (burst trigger, pulse function)

**get\_offset** (*channel=None*)

Get output offset

**set\_offset** (*offset, channel=None*)

Set output offset



**get\_amplitude** (*channel=None*)  
Get output amplitude (i.e., half of the span)

**set\_amplitude** (*amplitude, channel=None*)  
Set output amplitude (i.e., half of the span)

**BackendError**  
alias of `pylablib.core.devio.comm_backend.DeviceBackendError`

**Error**  
alias of `pylablib.devices.AWG.generic.GenericAWGError`

**class NoParameterCaller** (*device, kind*)  
Bases: `object`  
Class to simplify calling functions without a parameter

**ReraiseError**  
alias of `pylablib.devices.AWG.generic.GenericAWGBackendError`

**apply\_settings** (*settings*)  
Apply the settings.  
*settings* is a dict {*name*: *value*} of the available device settings. Non-applicable settings are ignored.

**ask** (*msg, data\_type='string', delay=0.0, timeout=None, read\_echo=False*)  
Write a message and read a reply.  
*msg* is the query message, *delay* is the delay between write and read. Other parameters are the same as in `read()`. If `read_echo==True`, assume that the device first echoes the input and skip it.

**close** ()  
Close the backend

**enable\_burst** (*enabled=True, channel=None*)  
Enable burst mode

**enable\_output** (*enabled=True, channel=None*)  
Turn the output on or off

**enable\_sync\_output** (*enabled=True, channel=None*)  
Enable or disable SYNC output

**enable\_trigger\_output** (*enabled=True, channel=None*)  
Enable trigger output

**flush** (*one\_line=False*)  
Flush the read buffer (read all the available data and return the number of bytes read).  
If `one_line==True`, read only a single line.

**static get\_arg\_type** (*arg*)  
Autodetect argument type

**get\_burst\_mode** (*channel=None*)  
Get burst mode.  
Can be either "trig" or "gate".

**get\_burst\_ncycles** (*channel=None*)  
Get burst mode ncycles.  
Infinite corresponds to a large value (>1E37).

**get\_channels\_number** ()

Get the number of channels

**get\_current\_channel** ()

Get current channel

**get\_device\_variable** (*key*)

Get the value of a settings, status, or full info parameter

**get\_duty\_cycle** (*channel=None*)

Get output duty cycle (in percent).

Only applies to "square" output function.

**get\_esr** (*timeout=None*)

Get the device status register (by default, "\*ESR?" command)

**get\_frequency** (*channel=None*)

Get output frequency

**get\_full\_info** (*include=0*)

Get dict {name: value} containing full device information (including status and settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_full\_status** (*include=0*)

Get dict {name: value} containing the device status (including settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_function** (*channel=None*)

Get output function.

Can be one of the following: "sine", "square", "ramp", "pulse", "noise", "prbs", "DC", "user", "arb". Not all functions can be available, depending on the particular model of the generator.

**get\_gate\_polarity** (*channel=None*)

Get burst gate polarity.

Can be either "norm" or "inv".

**get\_id** (*timeout=None*)

Get the device IDN. (query SCPI '\*IDN?' command)

**get\_load** (*channel=None*)

Get the output load

**get\_output\_polarity** (*channel=None*)

Get output polarity.

Can be either "norm" or "inv".

**get\_output\_range** (*channel=None*)

Get output voltage range.

Return tuple (*vmin*, *vmax*) with the low and high voltage values (i.e., offset-amplitude and offset+amplitude).

**get\_output\_trigger\_slope** (*channel=None*)

Get output trigger slope.

Can be either "pos", or "neg".

**get\_phase** (*channel=None*)

Get output phase (in degrees)

**get\_pulse\_width** (*channel=None*)

Get output pulse width (in seconds).

Only applies to "pulse" output function.

**get\_ramp\_symmetry** (*channel=None*)

Get output ramp symmetry (in percent).

Only applies to "ramp" output function.

**get\_settings** (*include=0*)

Get dict {name: value} containing all the device settings.

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_trigger\_slope** (*channel=None*)

Get trigger slope.

Can be either "pos", or "neg".

**get\_trigger\_source** (*channel=None*)

Get trigger source.

Can be either "imm", "ext", or "bus".

**is\_burst\_enabled** (*channel=None*)

Check if the burst mode is enabled

**is\_opened** ()

Check if the device is connected

**is\_output\_enabled** (*channel=None*)

Check if the output is enabled

**is\_sync\_output\_enabled** (*channel=None*)

Check if SYNC output is enabled

**is\_trigger\_output\_enabled** (*channel=None*)

Check if the trigger output is enabled

**lock** (*timeout=None*)

Lock the access to the device from other threads/processes (isn't necessarily implemented)

**locking** (*timeout=None*)

Context manager for lock & unlock

**open** ()

Open the backend

**static parse\_array\_data** (*data, fmt, include\_header=False*)

Parse the data returned by the device. *fmt* is [DataFormat](#) description in numpy format (e.g., "<u2").

If *include\_header==True*, the data is assumed to be in a (somewhat) standard SCPI format: b'#', then a single digit *s* denoting length of the size block, then *s* digits denoting length of the data (in

bytes) followed by the actual data. Otherwise (`include_header==False`), assume that the header is already removed.

**read** (*data\_type='string', timeout=None*)

Read data from the device.

*data\_type* determines the type of the data. Can be 'raw' (just raw data), 'string' (with trailing and leading spaces stripped), 'int', 'float', 'bool' (interprets 0 or 'off' as False, anything else as True), 'value' (returns tuple (value, unit), where *value* is float), a callable (return the result of this callable applied to the string value), a dictionary (return the stored value corresponding to the string value, or to the value converted into integer if the string value is not present), or a list of data types (the result is treated as a list of values with the given types separated by commas). *timeout* overrides the default value.

**read\_binary\_array\_data** (*include\_header=False, timeout=None, flush\_term=True*)

Read a binary data in the from the device.

The data assumes the standard binary transfer header consisting of "#" symbol, then a single digit with the size of the length string, then the length string containing the length of the binary data (in bytes). If `include_header==True`, return the data with the header; otherwise, return only the content. If `flush_term==True`, flush the following line to skip terminator characters after the binary data, which are added by some devices. *timeout* overrides the default value.

**reconnect** (*new\_instrument=True, ignore\_error=True*)

Remake the connection.

If `new_instrument==True`, create a new backend instance. If `ignore_error==True`, ignore errors on closing.

**reset** ()

Reset the device (by default, "\*RST" command)

**select\_current\_channel** (*channel*)

Select current default channel

**set\_burst\_mode** (*mode, channel=None*)

Set burst mode.

Can be either "trig" or "gate".

**set\_burst\_ncycles** (*ncycles=1, channel=None*)

Set burst mode ncycles.

Infinite corresponds to None

**set\_device\_variable** (*key, value*)

Set the value of a settings parameter

**set\_duty\_cycle** (*dcycle, channel=None*)

Set output duty cycle (in percent).

Only applies to "square" output function.

**set\_frequency** (*frequency, channel=None*)

Set output frequency

**set\_function** (*func, channel=None*)

Set output function.

Can be one of the following: "sine", "square", "ramp", "pulse", "noise", "prbs", "DC", "user", "arb". Not all functions can be available, depending on the particular model of the generator.

**set\_gate\_polarity** (*polarity='norm', channel=None*)

Set burst gate polarity.

Can be either "norm" or "inv".

**set\_load** (*load=None, channel=None*)

Set the output load (None means High-Z)

**set\_output\_polarity** (*polarity='norm', channel=None*)

Set output polarity.

Can be either "norm" or "inv".

**set\_output\_range** (*rng, channel=None*)

Set output voltage range.

If span is less than 1E-4, automatically switch to DC mode.

**set\_output\_trigger\_slope** (*slope, channel=None*)

Set output trigger slope.

Can be either "pos", or "neg".

**set\_phase** (*phase, channel=None*)

Set output phase (in degrees)

**set\_pulse\_width** (*width, channel=None*)

Set output pulse width (in seconds).

Only applies to "pulse" output function.

**set\_ramp\_symmetry** (*rsymm, channel=None*)

Set output ramp symmetry (in percent).

Only applies to "ramp" output function.

**set\_trigger\_slope** (*slope, channel=None*)

Set trigger slope.

Can be either "pos", or "neg".

**set\_trigger\_source** (*src, channel=None*)

Set trigger source.

Can be either "imm", "ext", or "bus".

**sleep** (*delay*)

Wait for *delay* seconds

**sync\_phase** ()

Synchronize phase between two channels

**unlock** ()

Unlock the access to the device from other threads/processes (isn't necessarily implemented)

**using\_write\_buffer** ()

Context manager for using a write buffer.

While it's active, all the consecutive `write()` operations are bundled together with ; delimiter. The actual write is performed at the `read()/ask()` operation or at the end of the block.

**wait** (*wait\_type='sync', timeout=None, wait\_callback=None*)

Pause execution until device overlapped commands are complete.

*wait\_type* is either 'sync' (perform `wait_sync()`), 'dev' (perform `wait_dev()`) or 'none' (do nothing).

**wait\_dev()**

Pause execution of the device commands until device overlapped commands (e.g., taking sweeps) are complete.

Note that the code execution is not paused.

**wait\_sync** (*timeout=None, wait\_callback=None*)

Pause execution of the script until device overlapped commands (e.g., taking sweeps) are complete.

*timeout* and *wait\_callback* override default constructor parameters.

**write** (*msg, arg=None, arg\_type=None, unit=None, bool\_selector=None, wait\_sync=None, read\_echo=False, read\_echo\_delay=0.0*)

Send a command.

#### Parameters

- **msg** (*str*) – Text message.
- **arg** – Optional argument to append in the end. If a list of arguments is supplied, the result is joined with ", ".
- **arg\_type** (*str*) – Argument type. Can be 'raw' (in which case data is sent raw), 'string', 'int', 'float', 'bool', a format string (such as '{:.3f}') or a list of argument types (for an iterable argument); if format string is used and the argument is a list or a tuple, then it is expanded as a list of arguments (e.g., *arg\_type*='{0};{1}' with *arg*=[1, 2] will produce a string '1; 2'); if a list of types is used, each element of *arg* is converted using the corresponding type, and the result is joined with ", ".
- **unit** (*str*) – If not None, use it as a unit to append after the value.
- **bool\_selector** (*tuple*) – A tuple (*false\_value, true\_value*) of two strings to represent bool argument; by default, use *.\_bool\_selector* attribute.
- **wait\_sync** – if True, append the sync command (specified as *.\_wait\_sync\_comm* attribute, "\*OPC?" by default) after the message and pause the execution command is complete; useful in long set operations, where the device might ignore later inputs until the current command is complete; if None, use the class default *.\_default\_write\_sync* attribute (False by default).
- **read\_echo** (*bool*) – If True, read a single line after write.
- **read\_echo\_delay** (*float*) – The delay between write and read if *read\_echo*==True.

**class** `pylablib.devices.AWG.specific.InstekAFG2000(addr)`

Bases: `pylablib.devices.AWG.specific.InstekAFG2225`

Instek AFG2000/2100 series AWG.

Compared to AFG2225, has only one channel and fewer capabilities.

#### BackendError

alias of `pylablib.core.devio.comm_backend.DeviceBackendError`

#### Error

alias of `pylablib.devices.AWG.generic.GenericAWGError`

**class** `NoParameterCaller(device, kind)`

Bases: `object`

Class to simplify calling functions without a parameter

**ReraiseError**

alias of `pylablib.devices.AWG.generic.GenericAWGBackendError`

**apply\_settings** (*settings*)

Apply the settings.

*settings* is a dict {name: value} of the available device settings. Non-applicable settings are ignored.

**ask** (*msg*, *data\_type*='string', *delay*=0.0, *timeout*=None, *read\_echo*=False)

Write a message and read a reply.

*msg* is the query message, *delay* is the delay between write and read. Other parameters are the same as in `read()`. If *read\_echo*==True, assume that the device first echoes the input and skip it.

**close** ()

Close the backend

**enable\_burst** (*enabled*=True, *channel*=None)

Enable burst mode

**enable\_output** (*enabled*=True, *channel*=None)

Turn the output on or off

**enable\_sync\_output** (*enabled*=True, *channel*=None)

Enable or disable SYNC output

**enable\_trigger\_output** (*enabled*=True, *channel*=None)

Enable trigger output

**flush** (*one\_line*=False)

Flush the read buffer (read all the available data and return the number of bytes read).

If *one\_line*==True, read only a single line.

**get\_amplitude** (*channel*=None)

Get output amplitude (i.e., half of the span)

**static get\_arg\_type** (*arg*)

Autodetect argument type

**get\_burst\_mode** (*channel*=None)

Get burst mode.

Can be either "trig" or "gate".

**get\_burst\_ncycles** (*channel*=None)

Get burst mode ncycles.

Infinite corresponds to a large value (>1E37).

**get\_channels\_number** ()

Get the number of channels

**get\_current\_channel** ()

Get current channel

**get\_device\_variable** (*key*)

Get the value of a settings, status, or full info parameter

**get\_duty\_cycle** (*channel*=None)

Get output duty cycle (in percent).

Only applies to "square" output function.

**get\_esr** (*timeout=None*)

Get the device status register (by default, "\*ESR?" command)

**get\_frequency** (*channel=None*)

Get output frequency

**get\_full\_info** (*include=0*)

Get dict {name: value} containing full device information (including status and settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_full\_status** (*include=0*)

Get dict {name: value} containing the device status (including settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_function** (*channel=None*)

Get output function.

Can be one of the following: "sine", "square", "ramp", "pulse", "noise", "prbs", "DC", "user", "arb". Not all functions can be available, depending on the particular model of the generator.

**get\_gate\_polarity** (*channel=None*)

Get burst gate polarity.

Can be either "norm" or "inv".

**get\_id** (*timeout=None*)

Get the device IDN. (query SCPI '\*IDN?' command)

**get\_load** (*channel=None*)

Get the output load

**get\_offset** (*channel=None*)

Get output offset

**get\_output\_polarity** (*channel=None*)

Get output polarity.

Can be either "norm" or "inv".

**get\_output\_range** (*channel=None*)

Get output voltage range.

Return tuple (vmin, vmax) with the low and high voltage values (i.e., offset-amplitude and offset+amplitude).

**get\_output\_trigger\_slope** (*channel=None*)

Get output trigger slope.

Can be either "pos", or "neg".

**get\_phase** (*channel=None*)

Get output phase (in degrees)

**get\_pulse\_width** (*channel=None*)

Get output pulse width (in seconds).



Only applies to "pulse" output function.

**get\_ramp\_symmetry** (*channel=None*)

Get output ramp symmetry (in percent).

Only applies to "ramp" output function.

**get\_settings** (*include=0*)

Get dict {name: value} containing all the device settings.

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_trigger\_slope** (*channel=None*)

Get trigger slope.

Can be either "pos", or "neg".

**get\_trigger\_source** (*channel=None*)

Get trigger source.

Can be either "imm", "ext", or "bus".

**is\_burst\_enabled** (*channel=None*)

Check if the burst mode is enabled

**is\_opened** ()

Check if the device is connected

**is\_output\_enabled** (*channel=None*)

Check if the output is enabled

**is\_sync\_output\_enabled** (*channel=None*)

Check if SYNC output is enabled

**is\_trigger\_output\_enabled** (*channel=None*)

Check if the trigger output is enabled

**lock** (*timeout=None*)

Lock the access to the device from other threads/processes (isn't necessarily implemented)

**locking** (*timeout=None*)

Context manager for lock & unlock

**open** ()

Open the backend

**static parse\_array\_data** (*data, fmt, include\_header=False*)

Parse the data returned by the device. *fmt* is *DataFormat* description in numpy format (e.g., "<u2").

If *include\_header==True*, the data is assumed to be in a (somewhat) standard SCPI format: b'#', then a single digit *s* denoting length of the size block, then *s* digits denoting length of the data (in bytes) followed by the actual data. Otherwise (*include\_header==False*), assume that the header is already removed.

**read** (*data\_type='string', timeout=None*)

Read data from the device.

*data\_type* determines the type of the data. Can be 'raw' (just raw data), 'string' (with trailing and leading spaces stripped), 'int', 'float', 'bool' (interprets 0 or 'off' as False, anything else as True), 'value' (returns tuple (value, unit), where *value* is float), a callable (return the result of this callable applied to the string value), a dictionary (return the stored value corresponding to the

string value, or to the value converted into integer if the string value is not present), or a list of data types (the result is treated as a list of values with the given types separated by commas). *timeout* overrides the default value.

**read\_binary\_array\_data** (*include\_header=False, timeout=None, flush\_term=True*)

Read a binary data in the from the device.

The data assumes the standard binary transfer header consisting of "#" symbol, then a single digit with the size of the length string, then the length string containing the length of the binary data (in bytes). If *include\_header==True*, return the data with the header; otherwise, return only the content. If *flush\_term==True*, flush the following line to skip terminator characters after the binary data, which are added by some devices. *timeout* overrides the default value.

**reconnect** (*new\_instrument=True, ignore\_error=True*)

Remake the connection.

If *new\_instrument==True*, create a new backend instance. If *ignore\_error==True*, ignore errors on closing.

**reset** ()

Reset the device (by default, "\*RST" command)

**select\_current\_channel** (*channel*)

Select current default channel

**set\_amplitude** (*amplitude, channel=None*)

Set output amplitude (i.e., half of the span)

**set\_burst\_mode** (*mode, channel=None*)

Set burst mode.

Can be either "trig" or "gate".

**set\_burst\_ncycles** (*ncycles=1, channel=None*)

Set burst mode ncycles.

Infinite corresponds to None

**set\_device\_variable** (*key, value*)

Set the value of a settings parameter

**set\_duty\_cycle** (*dcycle, channel=None*)

Set output duty cycle (in percent).

Only applies to "square" output function.

**set\_frequency** (*frequency, channel=None*)

Set output frequency

**set\_function** (*func, channel=None*)

Set output function.

Can be one of the following: "sine", "square", "ramp", "pulse", "noise", "prbs", "DC", "user", "arb". Not all functions can be available, depending on the particular model of the generator.

**set\_gate\_polarity** (*polarity='norm', channel=None*)

Set burst gate polarity.

Can be either "norm" or "inv".

**set\_load** (*load=None, channel=None*)

Set the output load (None means High-Z)

**set\_offset** (*offset*, *channel=None*)  
Set output offset

**set\_output\_polarity** (*polarity='norm'*, *channel=None*)  
Set output polarity.  
Can be either "norm" or "inv".

**set\_output\_range** (*rng*, *channel=None*)  
Set output voltage range.  
If span is less than  $1\text{E}-4$ , automatically switch to DC mode.

**set\_output\_trigger\_slope** (*slope*, *channel=None*)  
Set output trigger slope.  
Can be either "pos", or "neg".

**set\_phase** (*phase*, *channel=None*)  
Set output phase (in degrees)

**set\_pulse\_width** (*width*, *channel=None*)  
Set output pulse width (in seconds).  
Only applies to "pulse" output function.

**set\_ramp\_symmetry** (*rsymm*, *channel=None*)  
Set output ramp symmetry (in percent).  
Only applies to "ramp" output function.

**set\_trigger\_slope** (*slope*, *channel=None*)  
Set trigger slope.  
Can be either "pos", or "neg".

**set\_trigger\_source** (*src*, *channel=None*)  
Set trigger source.  
Can be either "imm", "ext", or "bus".

**sleep** (*delay*)  
Wait for *delay* seconds

**sync\_phase** ()  
Synchronize phase between two channels

**unlock** ()  
Unlock the access to the device from other threads/processes (isn't necessarily implemented)

**using\_write\_buffer** ()  
Context manager for using a write buffer.

While it's active, all the consecutive `write()` operations are bundled together with ; delimiter. The actual write is performed at the `read()/ask()` operation or at the end of the block.

**wait** (*wait\_type='sync'*, *timeout=None*, *wait\_callback=None*)  
Pause execution until device overlapped commands are complete.

*wait\_type* is either 'sync' (perform `wait_sync()`), 'dev' (perform `wait_dev()`) or 'none' (do nothing).

**wait\_dev** ()  
Pause execution of the device commands until device overlapped commands (e.g., taking sweeps) are complete.

Note that the code execution is not paused.

**wait\_sync** (*timeout=None, wait\_callback=None*)

Pause execution of the script until device overlapped commands (e.g., taking sweeps) are complete.

*timeout* and *wait\_callback* override default constructor parameters.

**write** (*msg, arg=None, arg\_type=None, unit=None, bool\_selector=None, wait\_sync=None, read\_echo=False, read\_echo\_delay=0.0*)

Send a command.

#### Parameters

- **msg** (*str*) – Text message.
- **arg** – Optional argument to append in the end. If a list of arguments is supplied, the result is joined with ", ".
- **arg\_type** (*str*) – Argument type. Can be 'raw' (in which case data is sent raw), 'string', 'int', 'float', 'bool', a format string (such as '{:.3f}') or a list of argument types (for an iterable argument); if format string is used and the argument is a list or a tuple, then it is expanded as a list of arguments (e.g., *arg\_type*='{0};{1}' with *arg*=[1, 2] will produce a string '1; 2'); if a list of types is used, each element of *arg* is converted using the corresponding type, and the result is joined with ", ".
- **unit** (*str*) – If not None, use it as a unit to append after the value.
- **bool\_selector** (*tuple*) – A tuple (*false\_value, true\_value*) of two strings to represent bool argument; by default, use *.\_bool\_selector* attribute.
- **wait\_sync** – if True, append the sync command (specified as *.\_wait\_sync\_comm* attribute, "\*OPC?" by default) after the message and pause the execution command is complete; useful in long set operations, where the device might ignore later inputs until the current command is complete; if None, use the class default *.\_default\_write\_sync* attribute (False by default).
- **read\_echo** (*bool*) – If True, read a single line after write.
- **read\_echo\_delay** (*float*) – The delay between write and read if *read\_echo*==True.

**class** `pylablib.devices.AWG.specific.RSInstekAFG21000` (*addr*)

Bases: `pylablib.devices.AWG.specific.InstekAFG2000`

RS Instek AFG21000 series AWG.

Compared to Instek AFG2000, it takes care of the amplitude output bug.

**get\_offset** (*channel=None*)

Get output offset

**get\_amplitude** (*channel=None*)

Get output amplitude (i.e., half of the span)

#### BackendError

alias of `pylablib.core.devio.comm_backend.DeviceBackendError`

#### Error

alias of `pylablib.devices.AWG.generic.GenericAWGError`

**class** `NoParameterCaller` (*device, kind*)

Bases: `object`

Class to simplify calling functions without a parameter

**ReraiseError**

alias of `pylablib.devices.AWG.generic.GenericAWGBackendError`

**apply\_settings** (*settings*)

Apply the settings.

*settings* is a dict {name: value} of the available device settings. Non-applicable settings are ignored.

**ask** (*msg*, *data\_type*='string', *delay*=0.0, *timeout*=None, *read\_echo*=False)

Write a message and read a reply.

*msg* is the query message, *delay* is the delay between write and read. Other parameters are the same as in `read()`. If *read\_echo*==True, assume that the device first echoes the input and skip it.

**close** ()

Close the backend

**enable\_burst** (*enabled*=True, *channel*=None)

Enable burst mode

**enable\_output** (*enabled*=True, *channel*=None)

Turn the output on or off

**enable\_sync\_output** (*enabled*=True, *channel*=None)

Enable or disable SYNC output

**enable\_trigger\_output** (*enabled*=True, *channel*=None)

Enable trigger output

**flush** (*one\_line*=False)

Flush the read buffer (read all the available data and return the number of bytes read).

If *one\_line*==True, read only a single line.

**static get\_arg\_type** (*arg*)

Autodetect argument type

**get\_burst\_mode** (*channel*=None)

Get burst mode.

Can be either "trig" or "gate".

**get\_burst\_ncycles** (*channel*=None)

Get burst mode ncycles.

Infinite corresponds to a large value (>1E37).

**get\_channels\_number** ()

Get the number of channels

**get\_current\_channel** ()

Get current channel

**get\_device\_variable** (*key*)

Get the value of a settings, status, or full info parameter

**get\_duty\_cycle** (*channel*=None)

Get output duty cycle (in percent).

Only applies to "square" output function.

**get\_esr** (*timeout*=None)

Get the device status register (by default, "\*ESR?" command)

**get\_frequency** (*channel=None*)

Get output frequency

**get\_full\_info** (*include=0*)

Get dict {name: value} containing full device information (including status and settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_full\_status** (*include=0*)

Get dict {name: value} containing the device status (including settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_function** (*channel=None*)

Get output function.

Can be one of the following: "sine", "square", "ramp", "pulse", "noise", "prbs", "DC", "user", "arb". Not all functions can be available, depending on the particular model of the generator.

**get\_gate\_polarity** (*channel=None*)

Get burst gate polarity.

Can be either "norm" or "inv".

**get\_id** (*timeout=None*)

Get the device IDN. (query SCPI '\*IDN?' command)

**get\_load** (*channel=None*)

Get the output load

**get\_output\_polarity** (*channel=None*)

Get output polarity.

Can be either "norm" or "inv".

**get\_output\_range** (*channel=None*)

Get output voltage range.

Return tuple (*vmin*, *vmax*) with the low and high voltage values (i.e., offset-amplitude and offset+amplitude).

**get\_output\_trigger\_slope** (*channel=None*)

Get output trigger slope.

Can be either "pos", or "neg".

**get\_phase** (*channel=None*)

Get output phase (in degrees)

**get\_pulse\_width** (*channel=None*)

Get output pulse width (in seconds).

Only applies to "pulse" output function.

**get\_ramp\_symmetry** (*channel=None*)

Get output ramp symmetry (in percent).

Only applies to "ramp" output function.

**get\_settings** (*include=0*)

Get dict {name: value} containing all the device settings.

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_trigger\_slope** (*channel=None*)

Get trigger slope.

Can be either "pos", or "neg".

**get\_trigger\_source** (*channel=None*)

Get trigger source.

Can be either "imm", "ext", or "bus".

**is\_burst\_enabled** (*channel=None*)

Check if the burst mode is enabled

**is\_opened** ()

Check if the device is connected

**is\_output\_enabled** (*channel=None*)

Check if the output is enabled

**is\_sync\_output\_enabled** (*channel=None*)

Check if SYNC output is enabled

**is\_trigger\_output\_enabled** (*channel=None*)

Check if the trigger output is enabled

**lock** (*timeout=None*)

Lock the access to the device from other threads/processes (isn't necessarily implemented)

**locking** (*timeout=None*)

Context manager for lock & unlock

**open** ()

Open the backend

**static parse\_array\_data** (*data, fmt, include\_header=False*)

Parse the data returned by the device. *fmt* is *DataFormat* description in numpy format (e.g., "<u2").

If *include\_header==True*, the data is assumed to be in a (somewhat) standard SCPI format: b'#', then a single digit *s* denoting length of the size block, then *s* digits denoting length of the data (in bytes) followed by the actual data. Otherwise (*include\_header==False*), assume that the header is already removed.

**read** (*data\_type='string', timeout=None*)

Read data from the device.

*data\_type* determines the type of the data. Can be 'raw' (just raw data), 'string' (with trailing and leading spaces stripped), 'int', 'float', 'bool' (interprets 0 or 'off' as False, anything else as True), 'value' (returns tuple (value, unit), where *value* is float), a callable (return the result of this callable applied to the string value), a dictionary (return the stored value corresponding to the string value, or to the value converted into integer if the string value is not present), or a list of data types (the result is treated as a list of values with the given types separated by commas). *timeout* overrides the default value.

**read\_binary\_array\_data** (*include\_header=False, timeout=None, flush\_term=True*)

Read a binary data in the from the device.

The data assumes the standard binary transfer header consisting of "#" symbol, then a single digit with the size of the length string, then the length string containing the length of the binary data (in bytes). If `include_header==True`, return the data with the header; otherwise, return only the content. If `flush_term==True`, flush the following line to skip terminator characters after the binary data, which are added by some devices. `timeout` overrides the default value.

**reconnect** (*new\_instrument=True, ignore\_error=True*)

Remake the connection.

If `new_instrument==True`, create a new backend instance. If `ignore_error==True`, ignore errors on closing.

**reset** ()

Reset the device (by default, "\*RST" command)

**select\_current\_channel** (*channel*)

Select current default channel

**set\_amplitude** (*amplitude, channel=None*)

Set output amplitude (i.e., half of the span)

**set\_burst\_mode** (*mode, channel=None*)

Set burst mode.

Can be either "trig" or "gate".

**set\_burst\_ncycles** (*ncycles=1, channel=None*)

Set burst mode ncycles.

Infinite corresponds to None

**set\_device\_variable** (*key, value*)

Set the value of a settings parameter

**set\_duty\_cycle** (*dcycle, channel=None*)

Set output duty cycle (in percent).

Only applies to "square" output function.

**set\_frequency** (*frequency, channel=None*)

Set output frequency

**set\_function** (*func, channel=None*)

Set output function.

Can be one of the following: "sine", "square", "ramp", "pulse", "noise", "prbs", "DC", "user", "arb". Not all functions can be available, depending on the particular model of the generator.

**set\_gate\_polarity** (*polarity='norm', channel=None*)

Set burst gate polarity.

Can be either "norm" or "inv".

**set\_load** (*load=None, channel=None*)

Set the output load (None means High-Z)

**set\_offset** (*offset, channel=None*)

Set output offset

**set\_output\_polarity** (*polarity='norm', channel=None*)

Set output polarity.

Can be either "norm" or "inv".



**set\_output\_range** (*rng, channel=None*)

Set output voltage range.

If span is less than  $1\text{E}-4$ , automatically switch to DC mode.

**set\_output\_trigger\_slope** (*slope, channel=None*)

Set output trigger slope.

Can be either "pos", or "neg".

**set\_phase** (*phase, channel=None*)

Set output phase (in degrees)

**set\_pulse\_width** (*width, channel=None*)

Set output pulse width (in seconds).

Only applies to "pulse" output function.

**set\_ramp\_symmetry** (*rsymm, channel=None*)

Set output ramp symmetry (in percent).

Only applies to "ramp" output function.

**set\_trigger\_slope** (*slope, channel=None*)

Set trigger slope.

Can be either "pos", or "neg".

**set\_trigger\_source** (*src, channel=None*)

Set trigger source.

Can be either "imm", "ext", or "bus".

**sleep** (*delay*)

Wait for *delay* seconds

**sync\_phase** ()

Synchronize phase between two channels

**unlock** ()

Unlock the access to the device from other threads/processes (isn't necessarily implemented)

**using\_write\_buffer** ()

Context manager for using a write buffer.

While it's active, all the consecutive `write()` operations are bundled together with ; delimiter. The actual write is performed at the `read()/ask()` operation or at the end of the block.

**wait** (*wait\_type='sync', timeout=None, wait\_callback=None*)

Pause execution until device overlapped commands are complete.

*wait\_type* is either 'sync' (perform `wait_sync()`), 'dev' (perform `wait_dev()`) or 'none' (do nothing).

**wait\_dev** ()

Pause execution of the device commands until device overlapped commands (e.g., taking sweeps) are complete.

Note that the code execution is not paused.

**wait\_sync** (*timeout=None, wait\_callback=None*)

Pause execution of the script until device overlapped commands (e.g., taking sweeps) are complete.

*timeout* and *wait\_callback* override default constructor parameters.

```
write(msg, arg=None, arg_type=None, unit=None, bool_selector=None, wait_sync=None,
      read_echo=False, read_echo_delay=0.0)
    Send a command.
```

#### Parameters

- **msg** (*str*) – Text message.
- **arg** – Optional argument to append in the end. If a list of arguments is supplied, the result is joined with ", ".
- **arg\_type** (*str*) – Argument type. Can be 'raw' (in which case data is sent raw), 'string', 'int', 'float', 'bool', a format string (such as '{:.3f}') or a list of argument types (for an iterable argument); if format string is used and the argument is a list or a tuple, then it is expanded as a list of arguments (e.g., `arg_type='{0};{1}'` with `arg=[1, 2]` will produce a string `'1; 2'`); if a list of types is used, each element of `arg` is converted using the corresponding type, and the result is joined with ", ".
- **unit** (*str*) – If not None, use it as a unit to append after the value.
- **bool\_selector** (*tuple*) – A tuple (`false_value`, `true_value`) of two strings to represent bool argument; by default, use `._bool_selector` attribute.
- **wait\_sync** – if True, append the sync command (specified as `._wait_sync_comm` attribute, `"*OPC?"` by default) after the message and pause the execution command is complete; useful in long set operations, where the device might ignore later inputs until the current command is complete; if None, use the class default `._default_write_sync` attribute (False by default).
- **read\_echo** (*bool*) – If True, read a single line after write.
- **read\_echo\_delay** (*float*) – The delay between write and read if `read_echo==True`.

```
class pylablib.devices.AWG.specific.TektronixAFG1000(addr, channels_number='auto')
```

Bases: `pylablib.devices.AWG.generic.GenericAWG`

```
get_pulse_width(channel=None)
    Get output pulse width (in seconds).
```

Only applies to "pulse" output function.

```
set_pulse_width(width, channel=None)
    Set output pulse width (in seconds).
```

Only applies to "pulse" output function.

#### BackendError

alias of `pylablib.core.devio.comm_backend.DeviceBackendError`

#### Error

alias of `pylablib.devices.AWG.generic.GenericAWGError`

```
class NoParameterCaller(device, kind)
```

Bases: `object`

Class to simplify calling functions without a parameter

#### ReraiseError

alias of `pylablib.devices.AWG.generic.GenericAWGBackendError`

**apply\_settings** (*settings*)

Apply the settings.

*settings* is a dict {name: value} of the available device settings. Non-applicable settings are ignored.

**ask** (*msg*, *data\_type='string'*, *delay=0.0*, *timeout=None*, *read\_echo=False*)

Write a message and read a reply.

*msg* is the query message, *delay* is the delay between write and read. Other parameters are the same as in [read\(\)](#). If *read\_echo==True*, assume that the device first echoes the input and skip it.

**close** ()

Close the backend

**enable\_burst** (*enabled=True*, *channel=None*)

Enable burst mode

**enable\_output** (*enabled=True*, *channel=None*)

Turn the output on or off

**enable\_sync\_output** (*enabled=True*, *channel=None*)

Enable or disable SYNC output

**enable\_trigger\_output** (*enabled=True*, *channel=None*)

Enable trigger output

**flush** (*one\_line=False*)

Flush the read buffer (read all the available data and return the number of bytes read).

If *one\_line==True*, read only a single line.

**get\_amplitude** (*channel=None*)

Get output amplitude (i.e., half of the span)

**static get\_arg\_type** (*arg*)

Autodetect argument type

**get\_burst\_mode** (*channel=None*)

Get burst mode.

Can be either "trig" or "gate".

**get\_burst\_ncycles** (*channel=None*)

Get burst mode ncycles.

Infinite corresponds to a large value (>1E37).

**get\_channels\_number** ()

Get the number of channels

**get\_current\_channel** ()

Get current channel

**get\_device\_variable** (*key*)

Get the value of a settings, status, or full info parameter

**get\_duty\_cycle** (*channel=None*)

Get output duty cycle (in percent).

Only applies to "square" output function.

**get\_esr** (*timeout=None*)

Get the device status register (by default, "\*ESR?" command)

**get\_frequency** (*channel=None*)

Get output frequency

**get\_full\_info** (*include=0*)

Get dict {name: value} containing full device information (including status and settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_full\_status** (*include=0*)

Get dict {name: value} containing the device status (including settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_function** (*channel=None*)

Get output function.

Can be one of the following: "sine", "square", "ramp", "pulse", "noise", "prbs", "DC", "user", "arb". Not all functions can be available, depending on the particular model of the generator.

**get\_gate\_polarity** (*channel=None*)

Get burst gate polarity.

Can be either "norm" or "inv".

**get\_id** (*timeout=None*)

Get the device IDN. (query SCPI '\*IDN?' command)

**get\_load** (*channel=None*)

Get the output load

**get\_offset** (*channel=None*)

Get output offset

**get\_output\_polarity** (*channel=None*)

Get output polarity.

Can be either "norm" or "inv".

**get\_output\_range** (*channel=None*)

Get output voltage range.

Return tuple (vmin, vmax) with the low and high voltage values (i.e., offset+amplitude and offset+amplitude).

**get\_output\_trigger\_slope** (*channel=None*)

Get output trigger slope.

Can be either "pos", or "neg".

**get\_phase** (*channel=None*)

Get output phase (in degrees)

**get\_ramp\_symmetry** (*channel=None*)

Get output ramp symmetry (in percent).

Only applies to "ramp" output function.

**get\_settings** (*include=0*)

Get dict {name: value} containing all the device settings.

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_trigger\_slope** (*channel=None*)

Get trigger slope.

Can be either "pos", or "neg".

**get\_trigger\_source** (*channel=None*)

Get trigger source.

Can be either "imm", "ext", or "bus".

**is\_burst\_enabled** (*channel=None*)

Check if the burst mode is enabled

**is\_opened** ()

Check if the device is connected

**is\_output\_enabled** (*channel=None*)

Check if the output is enabled

**is\_sync\_output\_enabled** (*channel=None*)

Check if SYNC output is enabled

**is\_trigger\_output\_enabled** (*channel=None*)

Check if the trigger output is enabled

**lock** (*timeout=None*)

Lock the access to the device from other threads/processes (isn't necessarily implemented)

**locking** (*timeout=None*)

Context manager for lock & unlock

**open** ()

Open the backend

**static parse\_array\_data** (*data, fmt, include\_header=False*)

Parse the data returned by the device. *fmt* is *DataFormat* description in numpy format (e.g., "<u2").

If *include\_header==True*, the data is assumed to be in a (somewhat) standard SCPI format: b'#', then a single digit *s* denoting length of the size block, then *s* digits denoting length of the data (in bytes) followed by the actual data. Otherwise (*include\_header==False*), assume that the header is already removed.

**read** (*data\_type='string', timeout=None*)

Read data from the device.

*data\_type* determines the type of the data. Can be 'raw' (just raw data), 'string' (with trailing and leading spaces stripped), 'int', 'float', 'bool' (interprets 0 or 'off' as False, anything else as True), 'value' (returns tuple (value, unit), where *value* is float), a callable (return the result of this callable applied to the string value), a dictionary (return the stored value corresponding to the string value, or to the value converted into integer if the string value is not present), or a list of data types (the result is treated as a list of values with the given types separated by commas). *timeout* overrides the default value.

**read\_binary\_array\_data** (*include\_header=False, timeout=None, flush\_term=True*)

Read a binary data in the from the device.

The data assumes the standard binary transfer header consisting of "#" symbol, then a single digit with the size of the length string, then the length string containing the length of the binary data (in bytes). If `include_header==True`, return the data with the header; otherwise, return only the content. If `flush_term==True`, flush the following line to skip terminator characters after the binary data, which are added by some devices. `timeout` overrides the default value.

**reconnect** (*new\_instrument=True, ignore\_error=True*)

Remake the connection.

If `new_instrument==True`, create a new backend instance. If `ignore_error==True`, ignore errors on closing.

**reset** ()

Reset the device (by default, "\*RST" command)

**select\_current\_channel** (*channel*)

Select current default channel

**set\_amplitude** (*amplitude, channel=None*)

Set output amplitude (i.e., half of the span)

**set\_burst\_mode** (*mode, channel=None*)

Set burst mode.

Can be either "trig" or "gate".

**set\_burst\_ncycles** (*ncycles=1, channel=None*)

Set burst mode ncycles.

Infinite corresponds to None

**set\_device\_variable** (*key, value*)

Set the value of a settings parameter

**set\_duty\_cycle** (*dcycle, channel=None*)

Set output duty cycle (in percent).

Only applies to "square" output function.

**set\_frequency** (*frequency, channel=None*)

Set output frequency

**set\_function** (*func, channel=None*)

Set output function.

Can be one of the following: "sine", "square", "ramp", "pulse", "noise", "prbs", "DC", "user", "arb". Not all functions can be available, depending on the particular model of the generator.

**set\_gate\_polarity** (*polarity='norm', channel=None*)

Set burst gate polarity.

Can be either "norm" or "inv".

**set\_load** (*load=None, channel=None*)

Set the output load (None means High-Z)

**set\_offset** (*offset, channel=None*)

Set output offset

**set\_output\_polarity** (*polarity='norm', channel=None*)

Set output polarity.

Can be either "norm" or "inv".

**set\_output\_range** (*rng, channel=None*)

Set output voltage range.

If span is less than  $1\text{E}-4$ , automatically switch to DC mode.

**set\_output\_trigger\_slope** (*slope, channel=None*)

Set output trigger slope.

Can be either "pos", or "neg".

**set\_phase** (*phase, channel=None*)

Set output phase (in degrees)

**set\_ramp\_symmetry** (*rsymm, channel=None*)

Set output ramp symmetry (in percent).

Only applies to "ramp" output function.

**set\_trigger\_slope** (*slope, channel=None*)

Set trigger slope.

Can be either "pos", or "neg".

**set\_trigger\_source** (*src, channel=None*)

Set trigger source.

Can be either "imm", "ext", or "bus".

**sleep** (*delay*)

Wait for *delay* seconds

**sync\_phase** ()

Synchronize phase between two channels

**unlock** ()

Unlock the access to the device from other threads/processes (isn't necessarily implemented)

**using\_write\_buffer** ()

Context manager for using a write buffer.

While it's active, all the consecutive `write()` operations are bundled together with ; delimiter. The actual write is performed at the `read()/ask()` operation or at the end of the block.

**wait** (*wait\_type='sync', timeout=None, wait\_callback=None*)

Pause execution until device overlapped commands are complete.

*wait\_type* is either 'sync' (perform `wait_sync()`), 'dev' (perform `wait_dev()`) or 'none' (do nothing).

**wait\_dev** ()

Pause execution of the device commands until device overlapped commands (e.g., taking sweeps) are complete.

Note that the code execution is not paused.

**wait\_sync** (*timeout=None, wait\_callback=None*)

Pause execution of the script until device overlapped commands (e.g., taking sweeps) are complete.

*timeout* and *wait\_callback* override default constructor parameters.

**write** (*msg, arg=None, arg\_type=None, unit=None, bool\_selector=None, wait\_sync=None, read\_echo=False, read\_echo\_delay=0.0*)

Send a command.

#### Parameters

- **msg** (*str*) – Text message.
- **arg** – Optional argument to append in the end. If a list of arguments is supplied, the result is joined with ", ".
- **arg\_type** (*str*) – Argument type. Can be 'raw' (in which case data is sent raw), 'string', 'int', 'float', 'bool', a format string (such as '{:.3f}') or a list of argument types (for an iterable argument); if format string is used and the argument is a list or a tuple, then it is expanded as a list of arguments (e.g., `arg_type='{0};{1}'` with `arg=[1,2]` will produce a string '1;2'); if a list of types is used, each element of *arg* is converted using the corresponding type, and the result is joined with ", ".
- **unit** (*str*) – If not None, use it as a unit to append after the value.
- **bool\_selector** (*tuple*) – A tuple (*false\_value*, *true\_value*) of two strings to represent bool argument; by default, use `._bool_selector` attribute.
- **wait\_sync** – if True, append the sync command (specified as `._wait_sync_comm` attribute, `"*OPC?"` by default) after the message and pause the execution command is complete; useful in long set operations, where the device might ignore later inputs until the current command is complete; if None, use the class default `._default_write_sync` attribute (False by default).
- **read\_echo** (*bool*) – If True, read a single line after write.
- **read\_echo\_delay** (*float*) – The delay between write and read if `read_echo==True`.

**class** `pylablib.devices.AWG.specific.RigolDG1000(addr)`

Bases: `pylablib.devices.AWG.generic.GenericAWG`

Rigol DG1000 AWG.

**sync\_phase** ()

Synchronize phase between two channels

**BackendError**

alias of `pylablib.core.devio.comm_backend.DeviceBackendError`

**Error**

alias of `pylablib.devices.AWG.generic.GenericAWGError`

**class** `NoParameterCaller(device, kind)`

Bases: `object`

Class to simplify calling functions without a parameter

**ReraiseError**

alias of `pylablib.devices.AWG.generic.GenericAWGBackendError`

**apply\_settings** (*settings*)

Apply the settings.

*settings* is a dict {name: value} of the available device settings. Non-applicable settings are ignored.

**ask** (*msg*, *data\_type*='string', *delay*=0.0, *timeout*=None, *read\_echo*=False)

Write a message and read a reply.

*msg* is the query message, *delay* is the delay between write and read. Other parameters are the same as in `read()`. If `read_echo==True`, assume that the device first echoes the input and skip it.

**close** ()

Close the backend



**enable\_burst** (*enabled=True, channel=None*)

Enable burst mode

**enable\_output** (*enabled=True, channel=None*)

Turn the output on or off

**enable\_sync\_output** (*enabled=True, channel=None*)

Enable or disable SYNC output

**enable\_trigger\_output** (*enabled=True, channel=None*)

Enable trigger output

**flush** (*one\_line=False*)

Flush the read buffer (read all the available data and return the number of bytes read).

If *one\_line==True*, read only a single line.

**get\_amplitude** (*channel=None*)

Get output amplitude (i.e., half of the span)

**static get\_arg\_type** (*arg*)

Autodetect argument type

**get\_burst\_mode** (*channel=None*)

Get burst mode.

Can be either "trig" or "gate".

**get\_burst\_ncycles** (*channel=None*)

Get burst mode ncycles.

Infinite corresponds to a large value (>1E37).

**get\_channels\_number** ()

Get the number of channels

**get\_current\_channel** ()

Get current channel

**get\_device\_variable** (*key*)

Get the value of a settings, status, or full info parameter

**get\_duty\_cycle** (*channel=None*)

Get output duty cycle (in percent).

Only applies to "square" output function.

**get\_esr** (*timeout=None*)

Get the device status register (by default, "\*ESR?" command)

**get\_frequency** (*channel=None*)

Get output frequency

**get\_full\_info** (*include=0*)

Get dict {name: value} containing full device information (including status and settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_full\_status** (*include=0*)

Get dict {name: value} containing the device status (including settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting `include=-10` queries all available variables, which is equivalent to `include="all"`.

**get\_function** (*channel=None*)

Get output function.

Can be one of the following: "sine", "square", "ramp", "pulse", "noise", "prbs", "DC", "user", "arb". Not all functions can be available, depending on the particular model of the generator.

**get\_gate\_polarity** (*channel=None*)

Get burst gate polarity.

Can be either "norm" or "inv".

**get\_id** (*timeout=None*)

Get the device IDN. (query SCPI '\*IDN?' command)

**get\_load** (*channel=None*)

Get the output load

**get\_offset** (*channel=None*)

Get output offset

**get\_output\_polarity** (*channel=None*)

Get output polarity.

Can be either "norm" or "inv".

**get\_output\_range** (*channel=None*)

Get output voltage range.

Return tuple (*vmin*, *vmax*) with the low and high voltage values (i.e., offset-amplitude and offset+amplitude).

**get\_output\_trigger\_slope** (*channel=None*)

Get output trigger slope.

Can be either "pos", or "neg".

**get\_phase** (*channel=None*)

Get output phase (in degrees)

**get\_pulse\_width** (*channel=None*)

Get output pulse width (in seconds).

Only applies to "pulse" output function.

**get\_ramp\_symmetry** (*channel=None*)

Get output ramp symmetry (in percent).

Only applies to "ramp" output function.

**get\_settings** (*include=0*)

Get dict {*name*: *value*} containing all the device settings.

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting `include=-10` queries all available variables, which is equivalent to `include="all"`.

**get\_trigger\_slope** (*channel=None*)

Get trigger slope.

Can be either "pos", or "neg".

**get\_trigger\_source** (*channel=None*)

Get trigger source.

Can be either "imm", "ext", or "bus".

**is\_burst\_enabled** (*channel=None*)

Check if the burst mode is enabled

**is\_opened** ()

Check if the device is connected

**is\_output\_enabled** (*channel=None*)

Check if the output is enabled

**is\_sync\_output\_enabled** (*channel=None*)

Check if SYNC output is enabled

**is\_trigger\_output\_enabled** (*channel=None*)

Check if the trigger output is enabled

**lock** (*timeout=None*)

Lock the access to the device from other threads/processes (isn't necessarily implemented)

**locking** (*timeout=None*)

Context manager for lock & unlock

**open** ()

Open the backend

**static parse\_array\_data** (*data, fmt, include\_header=False*)

Parse the data returned by the device. *fmt* is *DataFormat* description in numpy format (e.g., "<u2").

If *include\_header==True*, the data is assumed to be in a (somewhat) standard SCPI format: b'#', then a single digit *s* denoting length of the size block, then *s* digits denoting length of the data (in bytes) followed by the actual data. Otherwise (*include\_header==False*), assume that the header is already removed.

**read** (*data\_type='string', timeout=None*)

Read data from the device.

*data\_type* determines the type of the data. Can be 'raw' (just raw data), 'string' (with trailing and leading spaces stripped), 'int', 'float', 'bool' (interprets 0 or 'off' as False, anything else as True), 'value' (returns tuple (value, unit), where *value* is float), a callable (return the result of this callable applied to the string value), a dictionary (return the stored value corresponding to the string value, or to the value converted into integer if the string value is not present), or a list of data types (the result is treated as a list of values with the given types separated by commas). *timeout* overrides the default value.

**read\_binary\_array\_data** (*include\_header=False, timeout=None, flush\_term=True*)

Read a binary data in the from the device.

The data assumes the standard binary transfer header consisting of "#" symbol, then a single digit with the size of the length string, then the length string containing the length of the binary data (in bytes). If *include\_header==True*, return the data with the header; otherwise, return only the content. If *flush\_term==True*, flush the following line to skip terminator characters after the binary data, which are added by some devices. *timeout* overrides the default value.

**reconnect** (*new\_instrument=True, ignore\_error=True*)

Remake the connection.

If `new_instrument==True`, create a new backend instance. If `ignore_error==True`, ignore errors on closing.

**reset** ()

Reset the device (by default, "\*RST" command)

**select\_current\_channel** (*channel*)

Select current default channel

**set\_amplitude** (*amplitude*, *channel=None*)

Set output amplitude (i.e., half of the span)

**set\_burst\_mode** (*mode*, *channel=None*)

Set burst mode.

Can be either "trig" or "gate".

**set\_burst\_ncycles** (*ncycles=1*, *channel=None*)

Set burst mode ncycles.

Infinite corresponds to None

**set\_device\_variable** (*key*, *value*)

Set the value of a settings parameter

**set\_duty\_cycle** (*dcycle*, *channel=None*)

Set output duty cycle (in percent).

Only applies to "square" output function.

**set\_frequency** (*frequency*, *channel=None*)

Set output frequency

**set\_function** (*func*, *channel=None*)

Set output function.

Can be one of the following: "sine", "square", "ramp", "pulse", "noise", "prbs", "DC", "user", "arb". Not all functions can be available, depending on the particular model of the generator.

**set\_gate\_polarity** (*polarity='norm'*, *channel=None*)

Set burst gate polarity.

Can be either "norm" or "inv".

**set\_load** (*load=None*, *channel=None*)

Set the output load (None means High-Z)

**set\_offset** (*offset*, *channel=None*)

Set output offset

**set\_output\_polarity** (*polarity='norm'*, *channel=None*)

Set output polarity.

Can be either "norm" or "inv".

**set\_output\_range** (*rng*, *channel=None*)

Set output voltage range.

If span is less than  $1\text{E}-4$ , automatically switch to DC mode.

**set\_output\_trigger\_slope** (*slope*, *channel=None*)

Set output trigger slope.

Can be either "pos", or "neg".

**set\_phase** (*phase*, *channel=None*)

Set output phase (in degrees)

**set\_pulse\_width** (*width*, *channel=None*)

Set output pulse width (in seconds).

Only applies to "pulse" output function.

**set\_ramp\_symmetry** (*rsymm*, *channel=None*)

Set output ramp symmetry (in percent).

Only applies to "ramp" output function.

**set\_trigger\_slope** (*slope*, *channel=None*)

Set trigger slope.

Can be either "pos", or "neg".

**set\_trigger\_source** (*src*, *channel=None*)

Set trigger source.

Can be either "imm", "ext", or "bus".

**sleep** (*delay*)

Wait for *delay* seconds

**unlock** ()

Unlock the access to the device from other threads/processes (isn't necessarily implemented)

**using\_write\_buffer** ()

Context manager for using a write buffer.

While it's active, all the consecutive `write()` operations are bundled together with ; delimiter. The actual write is performed at the `read()/ask()` operation or at the end of the block.

**wait** (*wait\_type='sync'*, *timeout=None*, *wait\_callback=None*)

Pause execution until device overlapped commands are complete.

*wait\_type* is either 'sync' (perform `wait_sync()`), 'dev' (perform `wait_dev()`) or 'none' (do nothing).

**wait\_dev** ()

Pause execution of the device commands until device overlapped commands (e.g., taking sweeps) are complete.

Note that the code execution is not paused.

**wait\_sync** (*timeout=None*, *wait\_callback=None*)

Pause execution of the script until device overlapped commands (e.g., taking sweeps) are complete.

*timeout* and *wait\_callback* override default constructor parameters.

**write** (*msg*, *arg=None*, *arg\_type=None*, *unit=None*, *bool\_selector=None*, *wait\_sync=None*, *read\_echo=False*, *read\_echo\_delay=0.0*)

Send a command.

#### Parameters

- **msg** (*str*) – Text message.
- **arg** – Optional argument to append in the end. If a list of arguments is supplied, the result is joined with ", ".
- **arg\_type** (*str*) – Argument type. Can be 'raw' (in which case data is sent raw), 'string', 'int', 'float', 'bool', a format string (such as '{:.3f}'), or a list of argument types (for an iterable argument); if format string is used and

the argument is a list or a tuple, then it is expanded as a list of arguments (e.g., `arg_type='{0};{1}'` with `arg=[1, 2]` will produce a string `'1; 2'`); if a list of types is used, each element of `arg` is converted using the corresponding type, and the result is joined with `", "`.

- **unit** (*str*) – If not `None`, use it as a unit to append after the value.
- **bool\_selector** (*tuple*) – A tuple (`false_value`, `true_value`) of two strings to represent bool argument; by default, use `._bool_selector` attribute.
- **wait\_sync** – if `True`, append the sync command (specified as `._wait_sync_comm` attribute, `"*OPC?"` by default) after the message and pause the execution command is complete; useful in long set operations, where the device might ignore later inputs until the current command is complete; if `None`, use the class default `._default_write_sync` attribute (`False` by default).
- **read\_echo** (*bool*) – If `True`, read a single line after write.
- **read\_echo\_delay** (*float*) – The delay between write and read if `read_echo==True`.

## Module contents

### pylablib.devices.AlliedVision package

#### Submodules

#### pylablib.devices.AlliedVision.Bonito module

**exception** `pylablib.devices.AlliedVision.Bonito.BonitoError`

Bases: `pylablib.core.device.base.DeviceError`

Generic AlliedVision Bonito error

**args**

**with\_traceback()**

Exception.with\_traceback(tb) – set self.\_\_traceback\_\_ to tb and return self.

**class** `pylablib.devices.AlliedVision.Bonito.TDeviceInfo` (*version*, *serial\_number*,  
*grabber\_info*)

Bases: `tuple`

**count()**

Return number of occurrences of value.

**grabber\_info**

**index()**

Return first index of value.

Raises `ValueError` if the value is not present.

**serial\_number**

**version**

**class** `pylablib.devices.AlliedVision.Bonito.IBonitoCamera` (*\*\*kwargs*)

Bases: `pylablib.devices.interface.camera.ICamera`

**Error**

alias of `pylablib.core.devio.base.DeviceError`

**GrabberClass = None**

**open()**

Open the connection

**serial\_query** (*query*, *timeout=3.0*)

**get\_serial\_parameter** (*comm*, *kind='int'*, *timeout=3.0*)

**set\_serial\_parameter** (*comm*, *value*)

**get\_device\_info()**

Get camera model data.

Return tuple (*model*, *serial\_number*, *grabber\_info*).

**get\_detector\_size()**

Get camera detector size (in pixels) as a tuple (*width*, *height*); as the camera does not provide this information, use the frame grabber parameters

**get\_roi()**

Get current ROI.

Return tuple (*hstart*, *hend*, *vstart*, *vend*).

**set\_roi** (*hstart=0*, *hend=None*, *vstart=0*, *vend=None*)

Setup camera ROI.

By default, all non-supplied parameters take extreme values.

**get\_roi\_limits** (*hbin=1*, *vbin=1*)

**setup\_acquisition** (*mode='sequence'*, *nframes=100*)

Setup acquisition.

Any non-specified acquisition parameters are assumed to be the same as previously set (or default, if not explicitly set before). Return the new acquisition parameters.

**get\_exposure\_control\_mode()**

Get the exposure control mode.

Return tuple (*timing\_mode*, *feature\_mode*), where *timing\_mode* determines how the exposure and frame period are timed (continuous, external trigger control, internal control, etc.), and *feature\_mode* controls additional features (permanent exposure, enhanced full well mode). See documentation for details.

**set\_exposure\_control\_mode** (*timing\_mode=None*, *feature\_mode=None*)

Set the exposure control mode.

*timing\_mode* determines how the exposure and frame period are timed (continuous, external trigger control, internal control, etc.), and *feature\_mode* controls additional features (permanent exposure, enhanced full well mode). See documentation for details.

**get\_exposure()**

Get current exposure.

Note that the actual exposure might be different, depending on the exposure control mode.

**set\_exposure** (*exposure*, *setup\_mode=True*)

Set current exposure.

Note that the actual exposure might be different, depending on the exposure control mode. If `setup_mode==True`, automatically set the exposure mode to take the given exposure value into account.

**get\_frame\_period()**

Get frame period (time between two consecutive frames in the internal trigger mode).

Note that the actual frame period might be different, depending on the exposure control mode.

**set\_frame\_period**(*frame\_period*, *setup\_mode=True*)

Set frame period (time between two consecutive frames in the internal trigger mode).

Note that the actual frame period might be different, depending on the exposure control mode. If `setup_mode==True`, automatically set the exposure mode to take the given exposure value into account.

**get\_frame\_timings()**

Get acquisition timing.

Return tuple (*exposure*, *frame\_period*).

**is\_status\_line\_enabled()**

Check if the status line is on

**enable\_status\_line**(*enabled=True*)

Enable or disable status line

**get\_black\_level\_offset()**

Get the black level offset

**set\_black\_level\_offset**(*offset*)

Set the black level offset

**get\_digital\_gain()**

Get the digital gain (0 for 1x, 1 for 2x, 2 for 4x)

**set\_digital\_gain**(*gain*)

Get the digital gain (0 for 1x, 1 for 2x, 2 for 4x)

**FrameTransferError**

alias of `pylablib.devices.interface.camera.DefaultFrameTransferError`

**class NoParameterCaller**(*device*, *kind*)

Bases: `object`

Class to simplify calling functions without a parameter

**TimeoutError**

alias of `pylablib.core.devio.base.DeviceError`

**acquisition\_in\_progress()**

Check if acquisition is in progress

**apply\_settings**(*settings*)

Apply the settings.

*settings* is the dict {*name*: *value*} of the device available settings. Non-applicable settings are ignored.

**clear\_acquisition()**

Clear acquisition settings

**close()**

Close the connection



**get\_acquisition\_parameters()**

Get acquisition parameters.

Return dictionary {name: value}

**get\_data\_dimensions()**

Get readout data dimensions (in pixels) as a tuple (width, height); take indexing mode into account

**get\_device\_variable(key)**

Get the value of a settings, status, or full info parameter

**get\_frame\_format()**

Get format for the returned images.

Can be "list" (list of 2D arrays), "array" (a single 3D array), or "chunks" (list of 3D "chunk" arrays; supported for some cameras and provides the best performance).

**get\_frame\_info\_fields()**

Get the names of frame info fields.

Applicable when frame info format (set by `set_frame_info_format()`) is "list" or "array".

**get\_frame\_info\_format()**

Get format of the frame info.

Can be "namedtuple" (potentially nested named tuples; convenient to get particular values), "list" (flat list of values, with field names are given by `get_frame_info_fields()`; convenient for building a table), "array" (same as "list", but with a numpy array, which is easier to use for "chunks" frame format), or "dict" (flat dictionary with the same fields as the "list" format; more resilient to future format changes)

**get\_frame\_info\_period()**

Get period of frame info acquisition.

Frame info might be skipped (set to None) except for frames which indices are divisible by *period*. Useful for certain cameras where acquiring frame info takes a lot of time and can reduce performance at higher frame rates. Note that this parameter can still be ignored (i.e., always set to 1) if the performance is not an issue for a given camera class.

**get\_frames\_status()**

Get acquisition and buffer status.

Return tuple (acquired, unread, skipped, size), where acquired is the total number of acquired frames, unread is the number of acquired but not read frames, skipped is the number of skipped (not read and then written over) frames, and `buffer_size` is the total buffer size (in frames).

**get\_full\_info(include=0)**

Get dict {name: value} containing full device information (including status and settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting `include=-10` queries all available variables, which is equivalent to `include="all"`.

**get\_full\_status(include=0)**

Get dict {name: value} containing the device status (including settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting `include=-10` queries all available variables, which is equivalent to `include="all"`.

**get\_image\_indexing()**

Get indexing for the returned images.

Can be "rct" (first index row, second index column, rows counted from the top), "rcb" (same as "rc", rows counted from the bottom), "xyt" (first index column, second index row, rows counted from the top), or "xyb" (same as "xyt", rows counted from the bottom)

**get\_new\_images\_range()**

Get the range of the new images.

Return tuple (*first*, *last*) with images range (first inclusive). If no images are available, return None. If some images were in the buffer were overwritten, exclude them from the range.

**get\_settings(include=0)**

Get dict {*name*: *value*} containing all the device settings.

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include*=-10 queries all available variables, which is equivalent to *include*="all".

**grab(nframes=1, frame\_timeout=5.0, missing\_frame='skip', return\_info=False, buff\_size=None)**

Snap *nframes* images (with preset image read mode parameters)

*buff\_size* determines buffer size (if None, use the default size). Timeout is specified for a single-frame acquisition, not for the whole acquisition time. *missing\_frame* determines what to do with frames which have been lost: can be "none" (replacing them with None), "zero" (replacing them with zero-filled frame), or "skip" (skipping them, while still keeping total returned frames number to *n*). If *return\_info*==True, return tuple (*frames*, *infos*), where *infos* is a list of frame info tuples (camera-dependent); if some frames are missing and *missing\_frame*!="skip", the corresponding frame info is None.

**is\_acquisition\_setup()**

Check if acquisition is set up.

If the camera does not support separate acquisition setup, always return True.

**is\_opened()**

Check if the device is connected

**pausing\_acquisition(clear=None, stop=True, setup\_after=None, start\_after=True, combine\_nested=True)**

Context manager which temporarily pauses acquisition during execution of with block.

Useful for applying certain settings which can't be changed during the acquisition. If *clear*==True, clear acquisition in addition to stopping (by default, use the class default specified as *\_clear\_pausing\_acquisition* attribute). If *stop*==True, stop the acquisition (if *clear*==True, stop regardless). If *setup\_after*==True, setup the acquisition after pause if necessary (None means setup only if clearing was required). If *start\_after*==True, start the acquisition after pause if necessary (None means start only if stopping was required). If *combine\_nested*==True, then any nested *pausing\_acquisition* calls will stop/clear acquisition as necessary, but won't setup/start it again until this *pausing\_acquisition* call is complete.

Yields tuple (*acq\_in\_progress*, *acq\_params*), which indicates whether acquisition is currently in progress, and what are the current acquisition parameters.

**read\_multiple\_images(rng=None, peek=False, missing\_frame='skip', return\_info=False, return\_rng=False)**

Read multiple images specified by *rng* (by default, all un-read images).

If *rng* is specified, it is a tuple (*first*, *last*) with images range (first inclusive). If no new frames are available, return an empty list; if no acquisition is running, return None. If *peek*==True, return images

but not mark them as read. *missing\_frame* determines what to do with frames which are out of range (missing or lost): can be "none" (replacing them with None), "zero" (replacing them with zero-filled frame), or "skip" (skipping them). If *return\_info*==True, return tuple (frames, infos), where infos is a list of frame info tuples (camera-dependent, by default, only the frame index); if some frames are missing and *missing\_frame*!="skip", the corresponding frame info is None. if *return\_rng*==True, return the range covered resulting frames; if *missing\_frame*=="skip", the range can be smaller than the supplied *rng* if some frames are skipped.

**read\_newest\_image** (*peek*=False, *return\_info*=False)

Read the newest un-read image.

If no un-read frames are available, return None. If *peek*==True, return the image but not mark it as read. If *return\_info*==True, return tuple (frame, info), where info is an info tuples (camera-dependent, see [read\\_multiple\\_images\(\)](#)).

**read\_oldest\_image** (*peek*=False, *return\_info*=False)

Read the oldest un-read image.

If no un-read frames are available, return None. If *peek*==True, return the image but not mark it as read. If *return\_info*==True, return tuple (frame, info), where info is an info tuples (camera-dependent, see [read\\_multiple\\_images\(\)](#)).

**set\_device\_variable** (*key*, *value*)

Set the value of a settings parameter

**set\_frame\_format** (*fmt*)

Set format for the returned images.

Can be "list" (list of 2D arrays), "array" (a single 3D array), "chunks" (list of 3D "chunk" arrays; supported for some cameras and provides the best performance), or "try\_chunks" (same as "chunks", but if chunks are not supported, set to "list" instead). If format is "chunks" and chunks are not supported by the camera, it results in one frame per chunk. Note that if the format is set to "array" or "chunks", the frame info format is also automatically set to "array". If the format is set to "chunks", then the image info is also returned in chunks form (list of 2D info arrays with the same length as the corresponding frame chunks).

**set\_frame\_info\_format** (*fmt*, *include\_fields*=None)

Set format of the frame info.

Can be "namedtuple" (potentially nested named tuples; convenient to get particular values), "list" (flat list of values, with field names are given by [get\\_frame\\_info\\_fields\(\)](#); convenient for building a table), "array" (same as "list", but with a numpy array, which is easier to use for "chunks" frame format), or "dict" (flat dictionary with the same fields as the "list" format; more resilient to future format changes) If *include\_fields* is not None, it specifies the fields included for non-"tuple" formats; note that order or *include\_fields* is ignored, and the resulting fields are always ordered same as in the original.

**set\_frame\_info\_period** (*period*=1)

Set period of frame info acquisition.

Frame info might be skipped (set to None) except for frames which indices are divisible by *period*. Useful for certain cameras where acquiring frame info takes a lot of time and can reduce performance at higher frame rates. Note that this parameter can still be ignored (i.e., always set to 1) if the performance is not an issue for a given camera class.

**set\_image\_indexing** (*indexing*)

Set up indexing for the returned images.

Can be "rct" (first index row, second index column, rows counted from the top), "rcb" (same as "rc", rows counted from the bottom), "xyt" (first index column, second index row, rows counted from

the top), or "xyb" (same as "xyt", rows counted from the bottom)

**snap** (*timeout=5.0, return\_info=False*)  
Snap a single frame

**start\_acquisition** (*\*args, \*\*kwargs*)  
Start acquisition.

Can take the same keyword parameters as :meth:“*setup\_acquisition*. If the acquisition is not set up yet, set it up using the supplied parameters (use default of *setup\_acquisition()*, if the parameter is None). Otherwise, if any supplied parameters are different from the current ones, change them and reset the acquisition.

**stop\_acquisition** ()  
Stop acquisition

**wait\_for\_frame** (*since='lastread', nframes=1, timeout=20.0, error\_on\_stopped=False*)  
Wait for one or several new camera frames.

*since* specifies the reference point for waiting to acquire *nframes* frames; can be “lastread” (from the last read frame), “lastwait” (wait for the last successful *wait\_for\_frame()* call), “now” (from the start of the current call), or “start” (from the acquisition start, i.e., wait until *nframes* frames have been acquired). *timeout* can be either a number, None (infinite timeout), or a tuple (*timeout, frame\_timeout*), in which case the call times out if the total time exceeds *timeout*, or a single frame wait exceeds *frame\_timeout*. If the call times out, raise *TimeoutError*. If *error\_on\_stopped==True* and the acquisition is not running, raise *Error*; otherwise, simply return False without waiting.

**class** `pylablib.devices.AlliedVision.Bonito.BonitoIMAQCamera` (*imaq\_name='img0'*)  
Bases: `pylablib.devices.AlliedVision.Bonito.IBonitoCamera`, `pylablib.devices.IMAQ.IMAQFrameGrabber`

IMAQ+PFCam interface to a AlliedVision Bonito camera.

Parameters **imaq\_name** – IMAQ interface name (can be learned by `IMAQ.list_cameras()`; usually, but not always, starts with “img”)

**Error**  
alias of `pylablib.core.devio.base.DeviceError`

**GrabberClass**  
alias of `pylablib.devices.IMAQ.IMAQFrameGrabber`

**FrameTransferError**  
alias of `pylablib.devices.interface.camera.DefaultFrameTransferError`

**class** `NoParameterCaller` (*device, kind*)  
Bases: `object`

Class to simplify calling functions without a parameter

**TimeoutError** = `<Mock spec='str' id='140096274107088'>`

**acquisition\_in\_progress** ()  
Check if acquisition is in progress

**apply\_settings** (*settings*)  
Apply the settings.

*settings* is the dict {name: value} of the device available settings. Non-applicable settings are ignored.

**clear\_acquisition()**

Clear all acquisition details and free all buffers

**clear\_all\_triggers** (*reset\_acquisition=True*)

Disable all triggers of the session

If the input triggers configuration has been changed, acquisition needs to be restart; if `reset_acquisition==True`, perform it automatically.

**close()**

Close connection to the camera

**configure\_trigger\_in** (*trig\_type*, *trig\_line=0*, *trig\_pol='high'*, *trig\_action='none'*, *timeout=None*, *reset\_acquisition=True*)

Configure input trigger.

#### Parameters

- **trig\_type** (*str*) – trigger source type; can be "ext", "rtsi", "iso\_in", or "software"
- **trig\_line** (*int*) – trigger line number
- **trig\_pol** (*str*) – trigger polarity; can be "high" or "low"
- **trig\_action** (*str*) – trigger action; can be "none" (disable trigger), "capture" (start capturing), "stop" (stop capturing), "buffer" (capture a single frame), or "bufflist" (capture the whole buffer list once)
- **timeout** (*float*) – timeout in seconds; None means not timeout.
- **reset\_acquisition** (*bool*) – if the input triggers configuration has been changed, acquisition needs to be restart; if True, perform it automatically

**configure\_trigger\_out** (*trig\_type*, *trig\_line=0*, *trig\_pol='high'*, *trig\_drive='disable'*)

Configure trigger output.

#### Parameters

- **trig\_type** (*str*) – trigger drive destination type; can be "ext", "rtsi", or "iso\_out"
- **trig\_line** (*int*) – trigger line number
- **trig\_pol** (*str*) – trigger polarity; can be "high" or "low"
- **trig\_drive** (*str*) – trigger output signal; can be "disable" (disable drive), "acq\_in\_progress" (asserted when acquisition is started), "acq\_done" (asserted when acquisition is done), "unasserted" (force unasserted level), "asserted" (force asserted level), "hsync" (asserted on start of a single line start), "vsync" (asserted on start of a frame scan), "frame\_start" (asserted when a single frame is captured), or "frame\_done" (asserted when a single frame is done)

**enable\_status\_line** (*enabled=True*)

Enable or disable status line

**get\_acquisition\_parameters()**

Get acquisition parameters.

Return dictionary {name: value}

**get\_all\_grabber\_attribute\_values()**

Get a dictionary of all readable attributes.

The attributes types are autodetected, and some of the types of uncommon attributes may be misrepresented.

**get\_black\_level\_offset()**

Get the black level offset

**get\_data\_dimensions()**

Get readout data dimensions (in pixels) as a tuple (width, height); take indexing mode into account

**get\_detector\_size()**

Get camera detector size (in pixels) as a tuple (width, height); as the camera does not provide this information, use the frame grabber parameters

**get\_device\_info()**

Get camera model data.

Return tuple (model, serial\_number, grabber\_info).

**get\_device\_variable(key)**

Get the value of a settings, status, or full info parameter

**get\_digital\_gain()**

Get the digital gain (0 for 1x, 1 for 2x, 2 for 4x)

**get\_exposure()**

Get current exposure.

Note that the actual exposure might be different, depending on the exposure control mode.

**get\_exposure\_control\_mode()**

Get the exposure control mode.

Return tuple (timing\_mode, feature\_mode), where timing\_mode determines how the exposure and frame period are timed (continuous, external trigger control, internal control, etc.), and feature\_mode controls additional features (permanent exposure, enhanced full well mode). See documentation for details.

**get\_frame\_format()**

Get format for the returned images.

Can be "list" (list of 2D arrays), "array" (a single 3D array), or "chunks" (list of 3D "chunk" arrays; supported for some cameras and provides the best performance).

**get\_frame\_info\_fields()**

Get the names of frame info fields.

Applicable when frame info format (set by `set_frame_info_format()`) is "list" or "array".

**get\_frame\_info\_format()**

Get format of the frame info.

Can be "namedtuple" (potentially nested named tuples; convenient to get particular values), "list" (flat list of values, with field names are given by `get_frame_info_fields()`; convenient for building a table), "array" (same as "list", but with a numpy array, which is easier to use for "chunks" frame format), or "dict" (flat dictionary with the same fields as the "list" format; more resilient to future format changes)

**get\_frame\_info\_period()**

Get period of frame info acquisition.

Frame info might be skipped (set to None) except for frames which indices are divisible by *period*. Useful for certain cameras where acquiring frame info takes a lot of time and can reduce performance at higher

frame rates. Note that this parameter can still be ignored (i.e., always set to 1) if the performance is not an issue for a given camera class.

#### **get\_frame\_period()**

Get frame period (time between two consecutive frames in the internal trigger mode).

Note that the actual frame period might be different, depending on the exposure control mode.

#### **get\_frame\_timings()**

Get acquisition timing.

Return tuple (exposure, frame\_period).

#### **get\_frames\_status()**

Get acquisition and buffer status.

Return tuple (acquired, unread, skipped, size), where *acquired* is the total number of acquired frames, *unread* is the number of acquired but not read frames, *skipped* is the number of skipped (not read and then written over) frames, and *buffer\_size* is the total buffer size (in frames).

#### **get\_full\_info(include=0)**

Get dict {name: value} containing full device information (including status and settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

#### **get\_full\_status(include=0)**

Get dict {name: value} containing the device status (including settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

#### **get\_grabber\_attribute\_value(attr, default=None, kind='auto')**

Get value of an attribute with a given name or index.

If *default* is not None, return *default* if the attribute is not supported; otherwise, raise an error. *kind* is the attribute kind, and it can be "uint32", "uint64", "double", or "auto" (autodetect based on the stored list of attribute kinds).

#### **get\_grabber\_detector\_size()**

Get camera detector size (in pixels) as a tuple (width, height)

#### **get\_grabber\_roi()**

Get current ROI.

Return tuple (hstart, hend, vstart, vend). *hstart* and *hend* specify horizontal image extent, *vstart* and *vend* specify vertical image extent (start is inclusive, stop is exclusive, starting from 0).

#### **get\_grabber\_roi\_limits(hbin=1, vbin=1)**

Get the minimal and maximal ROI parameters.

Return tuple (hlim, vlim), where each element is in turn a limit 5-tuple (min, max, pstep, sstep, maxbin) with, correspondingly, minimal and maximal size, position and size step, and the maximal binning (fixed to 1 if not binning is allowed). In some cameras, the step and the minimal size depend on the binning, which can be supplied.

#### **get\_image\_indexing()**

Get indexing for the returned images.

Can be "rct" (first index row, second index column, rows counted from the top), "rcb" (same as "rc", rows counted from the bottom), "xyt" (first index column, second index row, rows counted from the top), or "xyb" (same as "xyt", rows counted from the bottom)

**get\_new\_images\_range()**

Get the range of the new images.

Return tuple (first, last) with images range (first inclusive). If no images are available, return None. If some images were in the buffer were overwritten, exclude them from the range.

**get\_roi()**

Get current ROI.

Return tuple (hstart, hend, vstart, vend).

**get\_roi\_limits(hbin=1, vbin=1)**

**get\_serial\_parameter(comm, kind='int', timeout=3.0)**

**get\_serial\_params()**

Return serial parameters as a tuple (write\_term, datatype)

**get\_settings(include=0)**

Get dict {name: value} containing all the device settings.

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include*=-10 queries all available variables, which is equivalent to *include*="all".

**grab(nframes=1, frame\_timeout=5.0, missing\_frame='skip', return\_info=False, buff\_size=None)**

Snap *nframes* images (with preset image read mode parameters)

*buff\_size* determines buffer size (if None, use the default size). Timeout is specified for a single-frame acquisition, not for the whole acquisition time. *missing\_frame* determines what to do with frames which have been lost: can be "none" (replacing them with None), "zero" (replacing them with zero-filled frame), or "skip" (skipping them, while still keeping total returned frames number to *n*). If *return\_info*==True, return tuple (frames, infos), where infos is a list of frame info tuples (camera-dependent); if some frames are missing and *missing\_frame*!="skip", the corresponding frame info is None.

**is\_acquisition\_setup()**

Check if acquisition is set up.

If the camera does not support separate acquisition setup, always return True.

**is\_opened()**

Check if the device is connected

**is\_status\_line\_enabled()**

Check if the status line is on

**open()**

Open the connection

**pausing\_acquisition(clear=None, stop=True, setup\_after=None, start\_after=True, combine\_nested=True)**

Context manager which temporarily pauses acquisition during execution of with block.

Useful for applying certain settings which can't be changed during the acquisition. If *clear*==True, clear acquisition in addition to stopping (by default, use the class default specified as *\_clear\_pausing\_acquisition* attribute). If *stop*==True, stop the acquisition (if *clear*==True, stop regardless). If *setup\_after*==True, setup the acquisition after pause if



necessary (None means setup only if clearing was required). If `start_after==True`, start the acquisition after pause if necessary (None means start only if stopping was required). If `combine_nested==True`, then any nested `pausing_acquisition` calls will stop/clear acquisition as necessary, but won't setup/start it again until this `pausing_acquisition` call is complete.

Yields tuple (`acq_in_progress`, `acq_params`), which indicates whether acquisition is currently in progress, and what are the current acquisition parameters.

**`read_multiple_images`** (*rng=None, peek=False, missing\_frame='skip', return\_info=False, return\_rng=False*)

Read multiple images specified by *rng* (by default, all un-read images).

If *rng* is specified, it is a tuple (*first*, *last*) with images range (first inclusive). If no new frames are available, return an empty list; if no acquisition is running, return None. If `peek==True`, return images but not mark them as read. *missing\_frame* determines what to do with frames which are out of range (missing or lost): can be "none" (replacing them with None), "zero" (replacing them with zero-filled frame), or "skip" (skipping them). If `return_info==True`, return tuple (*frames*, *infos*), where *infos* is a list of frame info tuples (camera-dependent, by default, only the frame index); if some frames are missing and `missing_frame!="skip"`, the corresponding frame info is None. If `return_rng==True`, return the range covered resulting frames; if `missing_frame=="skip"`, the range can be smaller than the supplied *rng* if some frames are skipped.

**`read_newest_image`** (*peek=False, return\_info=False*)

Read the newest un-read image.

If no un-read frames are available, return None. If `peek==True`, return the image but not mark it as read. If `return_info==True`, return tuple (*frame*, *info*), where *info* is an info tuples (camera-dependent, see `read_multiple_images()`).

**`read_oldest_image`** (*peek=False, return\_info=False*)

Read the oldest un-read image.

If no un-read frames are available, return None. If `peek==True`, return the image but not mark it as read. If `return_info==True`, return tuple (*frame*, *info*), where *info* is an info tuples (camera-dependent, see `read_multiple_images()`).

**`read_trigger`** (*trig\_type, trig\_line=0, trig\_pol='high'*)

Read current value of a trigger (input or output).

#### Parameters

- **`trig_type`** (*str*) – trigger drive destination type; can be "ext", "rtsi", "iso\_in", or "iso\_out"
- **`trig_line`** (*int*) – trigger line number
- **`trig_pol`** (*str*) – trigger polarity; can be "high" or "low"

**`reset`** ()

Reset connection to the camera

**`send_software_trigger`** ()

Send software trigger signal

**`serial_flush`** ()

Flush CameraLink serial port

**`serial_query`** (*query, timeout=3.0*)

**`serial_read`** (*n, timeout=3.0, datatype=None*)

Read specified number of bytes from CameraLink serial port.

#### Parameters

- **n** – number of bytes to read
- **timeout** – operation timeout (in seconds)
- **datatype** – return datatype; can be "bytes" (return raw bytes), or "str" (convert into UTF-8 string) if None, use the value set up using `setup_serial_params()` (by default, "bytes")

**serial\_readline** (*timeout=3.0, datatype=None, maxn=1024*)

Read bytes from CameraLink serial port until the termination character (defined in camera file) is encountered.

#### Parameters

- **timeout** – operation timeout (in seconds)
- **datatype** – return datatype; can be "bytes" (return raw bytes), or "str" (convert into UTF-8 string) if None, use the value set up using `setup_serial_params()` (by default, "bytes")
- **maxn** – maximal number of bytes to read

**serial\_write** (*msg, timeout=3.0, term=None*)

Write message into CameraLink serial port.

#### Parameters

- **timeout** – operation timeout (in seconds)
- **term** – additional write terminator character to add to the message; if None, use the value set up using `setup_serial_params()` (by default, no additional terminator)

**set\_black\_level\_offset** (*offset*)

Set the black level offset

**set\_device\_variable** (*key, value*)

Set the value of a settings parameter

**set\_digital\_gain** (*gain*)

Get the digital gain (0 for 1x, 1 for 2x, 2 for 4x)

**set\_exposure** (*exposure, setup\_mode=True*)

Set current exposure.

Note that the actual exposure might be different, depending on the exposure control mode. If `setup_mode==True`, automatically set the exposure mode to take the given exposure value into account.

**set\_exposure\_control\_mode** (*timing\_mode=None, feature\_mode=None*)

Set the exposure control mode.

*timing\_mode* determines how the exposure and frame period are timed (continuous, external trigger control, internal control, etc.), and *feature\_mode* controls additional features (permanent exposure, enhanced full well mode). See documentation for details.

**set\_frame\_format** (*fmt*)

Set format for the returned images.

Can be "list" (list of 2D arrays), "array" (a single 3D array), "chunks" (list of 3D "chunk" arrays; supported for some cameras and provides the best performance), or "try\_chunks" (same as "chunks", but if chunks are not supported, set to "list" instead). If format is "chunks" and chunks are not supported by the camera, it results in one frame per chunk. Note that if the format is set to "array" or "chunks", the frame info format is also automatically set to "array". If the format

is set to "chunks", then the image info is also returned in chunks form (list of 2D info arrays with the same length as the corresponding frame chunks).

**set\_frame\_info\_format** (*fmt*, *include\_fields=None*)

Set format of the frame info.

Can be "namedtuple" (potentially nested named tuples; convenient to get particular values), "list" (flat list of values, with field names are given by `get_frame_info_fields()`; convenient for building a table), "array" (same as "list", but with a numpy array, which is easier to use for "chunks" frame format), or "dict" (flat dictionary with the same fields as the "list" format; more resilient to future format changes) If *include\_fields* is not None, it specifies the fields included for non-"tuple" formats; note that order of *include\_fields* is ignored, and the resulting fields are always ordered same as in the original.

**set\_frame\_info\_period** (*period=1*)

Set period of frame info acquisition.

Frame info might be skipped (set to None) except for frames which indices are divisible by *period*. Useful for certain cameras where acquiring frame info takes a lot of time and can reduce performance at higher frame rates. Note that this parameter can still be ignored (i.e., always set to 1) if the performance is not an issue for a given camera class.

**set\_frame\_period** (*frame\_period*, *setup\_mode=True*)

Set frame period (time between two consecutive frames in the internal trigger mode).

Note that the actual frame period might be different, depending on the exposure control mode. If *setup\_mode==True*, automatically set the exposure mode to take the given exposure value into account.

**set\_grabber\_attribute\_value** (*attr*, *value*, *kind='int32'*)

Set value of an attribute with a given name or index.

*kind* is the attribute kind, and it can be "uint32", "uint64", "double", or "auto" (autodetect based on the stored list of attribute kinds).

**set\_grabber\_roi** (*hstart=0*, *hend=None*, *vstart=0*, *vend=None*)

Setup camera ROI.

*hstart* and *hend* specify horizontal image extent, *vstart* and *vend* specify vertical image extent (start is inclusive, stop is exclusive, starting from 0). By default, all non-supplied parameters take extreme values (0 for start, maximal for end).

**set\_image\_indexing** (*indexing*)

Set up indexing for the returned images.

Can be "rct" (first index row, second index column, rows counted from the top), "rcb" (same as "rc", rows counted from the bottom), "xyt" (first index column, second index row, rows counted from the top), or "xyb" (same as "xyt", rows counted from the bottom)

**set\_roi** (*hstart=0*, *hend=None*, *vstart=0*, *vend=None*)

Setup camera ROI.

By default, all non-supplied parameters take extreme values.

**set\_serial\_parameter** (*comm*, *value*)

**setup\_acquisition** (*mode='sequence'*, *nframes=100*)

Setup acquisition.

Any non-specified acquisition parameters are assumed to be the same as previously set (or default, if not explicitly set before). Return the new acquisition parameters.

**setup\_serial\_params** (*write\_term*=", *datatype*='bytes')

Setup default serial communication parameters.

#### Parameters

- **write\_term** – default terminator character to be added to the sent messages
- **datatype** – type of the result of read commands; can be "bytes" (return raw bytes), or "str" (convert into UTF-8 string)

**snap** (*timeout*=5.0, *return\_info*=False)

Snap a single frame

**start\_acquisition** (*\*args*, *\*\*kwargs*)

Start acquisition.

Can take the same keyword parameters as `:meth:“setup_acquisition`. If the acquisition is not set up yet, set it up using the supplied parameters (use default of `setup_acquisition()`, if the parameter is None). Otherwise, if any supplied parameters are different from the current ones, change them and reset the acquisition.

**stop\_acquisition** ()

Stop acquisition

**wait\_for\_frame** (*since*='lastread', *nframes*=1, *timeout*=20.0, *error\_on\_stopped*=False)

Wait for one or several new camera frames.

*since* specifies the reference point for waiting to acquire *nframes* frames; can be "lastread" (from the last read frame), "lastwait" (wait for the last successful `wait_for_frame()` call), "now" (from the start of the current call), or "start" (from the acquisition start, i.e., wait until *nframes* frames have been acquired). *timeout* can be either a number, None (infinite timeout), or a tuple (*timeout*, *frame\_timeout*), in which case the call times out if the total time exceeds *timeout*, or a single frame wait exceeds *frame\_timeout*. If the call times out, raise `TimeoutError`. If *error\_on\_stopped*==True and the acquisition is not running, raise `Error`; otherwise, simply return False without waiting.

`pylablib.devices.AlliedVision.Bonito.check_grabber_association` (*cam*)

Check if the given IMAQ frame grabber corresponds to Bonito camera.

*cam* should be an opened instance of `IMAQCamera`.

**class** `pylablib.devices.AlliedVision.Bonito.TStatusLine` (*framestamp*)

Bases: `tuple`

**count** ()

Return number of occurrences of value.

**framestamp**

**index** ()

Return first index of value.

Raises `ValueError` if the value is not present.

`pylablib.devices.AlliedVision.Bonito.get_status_lines` (*frames*)

Get frame info from the binary status line.

*frames* can be 2D array (one frame), 3D array (stack of frames, first index is frame number), or list of 1D or 2D arrays. Assume that the status line is present; if it isn't, the returned frame info will be a random noise. Return a 1D or 2D numpy array, where the first axis (if present) is the frame number, and the last is the status line entry.

**class** `pylablib.devices.AlliedVision.Bonito.BonitoStatusLineChecker`

Bases: `pylablib.devices.interface.camera.StatusLineChecker`

**get\_framestamp** (*frames*)  
Get framestamps from status lines in the given frames

**check\_indices** (*indices, step=1*)  
Check if indices are consistent with the given step

## Module contents

### pylablib.devices.Andor package

#### Submodules

#### pylablib.devices.Andor.AndorSDK2 module

```
class pylablib.devices.Andor.AndorSDK2.LibraryController (lib)
    Bases: pylablib.devices.utils.load_lib.LibraryController

    close (opid)
        Mark device closing.

        Return tuple (close_result, uninit_result) with the results of the closing and the shutdown.
        If library does not need to be shut down yet, set uninit_result=None

    open ()
        Mark device opening.

        Return tuple (init_result, open_result, opid) with the results of the initialization and the
        opening, and the opening ID which should afterwards be used for closing. If library is already initialized,
        set init_result=None

    preinit ()
        Pre-initialize the library, if it hasn't been done already

    shutdown ()
        Close all opened connections and shutdown the library

    temp_open ()
        Context for temporarily opening a new device connection

pylablib.devices.Andor.AndorSDK2.restart_lib ()
pylablib.devices.Andor.AndorSDK2.get_cameras_number ()
    Get number of connected Andor cameras

class pylablib.devices.Andor.AndorSDK2.TDeviceInfo (controller_model, head_model, se-
                                                rial_number)
    Bases: tuple

    controller_model

    count ()
        Return number of occurrences of value.

    head_model

    index ()
        Return first index of value.

        Raises ValueError if the value is not present.

    serial_number
```

```
class pylablib.devices.Andor.AndorSDK2.TCycleTimings (exposure, accum_cycle_time,  
                                                    kinetic_cycle_time)
```

Bases: `tuple`

**accum\_cycle\_time**

**count** ()

Return number of occurrences of value.

**exposure**

**index** ()

Return first index of value.

Raises `ValueError` if the value is not present.

**kinetic\_cycle\_time**

```
class pylablib.devices.Andor.AndorSDK2.TAcqProgress (frames_done, cycles_done)
```

Bases: `tuple`

**count** ()

Return number of occurrences of value.

**cycles\_done**

**frames\_done**

**index** ()

Return first index of value.

Raises `ValueError` if the value is not present.

```
class pylablib.devices.Andor.AndorSDK2.AndorSDK2Camera (idx=0, ini_path="",  
                                                    temperature=None,  
                                                    fan_mode='off',  
                                                    amp_mode=None)  
  
Bases: pylablib.devices.interface.camera.IBinROICamera, pylablib.devices.  
interface.camera.IExposureCamera
```

Andor SDK2 camera.

Due to the library features, the camera needs to set up all of the parameters to some default values upon connection. Most of these parameters are chosen as reasonable defaults: full ROI, minimal exposure time, closed shutter, internal trigger, fastest recommended verticals shift speed, no EMCCD gain. However, some should be supplied during the connection: temperature setpoint (where appropriate), fan mode, and amplifier mode; while there is still a possibility to have default values of these parameters, they might not be appropriate in some settings, and frequently need to be changed.

Caution: the manufacturer DLL is designed such that if the camera is not closed on the program termination, the allocated resources are never released. If this happens, these resources are blocked until the complete OS restart.

#### Parameters

- **idx** (*int*) – camera index (use `get_cameras_number()` to get the total number of connected cameras)
- **ini\_path** (*str*) – path to .ini file, if required by the camera
- **temperature** – initial temperature setpoint (in C); can also be `None` (select the bottom 20% of the whole range), or "off" (turn the cooler off and set the maximal of the whole range)
- **fan\_mode** – initial fan mode

- **amp\_mode** – initial amplifier mode (a tuple like the one returned by `get_amp_mode()`); can also be `None`, which selects the slowest, smallest gain mode

**Error**

alias of `pylablib.devices.Andor.base.AndorError`

**TimeoutError**

alias of `pylablib.devices.Andor.base.AndorTimeoutError`

**open()**

Open connection to the camera

**close()**

Close connection to the camera

**is\_opened()**

Check if the device is connected

**get\_device\_info()**

Get camera device info.

Return tuple (controller\_mode, head\_model, serial\_number).

**get\_status()**

Get camera status.

Return either "idle" (no acquisition), "acquiring" (acquisition in progress) or "temp\_cycle" (temperature cycle in progress).

**acquisition\_in\_progress()**

Check if acquisition is in progress

**get\_capabilities()**

Get camera capabilities.

For description of the structure, see Andor SDK manual.

**get\_pixel\_size()**

Get camera pixel size (in m)

**is\_cooler\_on()**

Check if the cooler is on

**set\_cooler (on=True)**

Set the cooler on or off

**get\_temperature\_status()**

Get temperature status.

Can return "off" (cooler off), "not\_reached" (cooling in progress), "not\_stabilized" (reached but not stabilized yet), "stabilized" (completely stabilized) or "drifted".

**get\_temperature()**

Get the current camera temperature (in C)

**set\_temperature (temperature, enable\_cooler=True)**

Change the temperature setpoint (in C).

If enable\_cooler==True, turn the cooler on automatically.

**get\_temperature\_setpoint()**

Get the temperature setpoint (in C)

**get\_temperature\_range()**

Return the available range of temperatures (in C)

**get\_all\_amp\_modes()**

Get all available preamp modes.

Each preamp mode is characterized by an AD channel index, amplifier index, channel speed (horizontal scan speed) index and preamp gain index. Return list of tuples (*channel*, *channel\_bitdepth*, *oamp*, *oamp\_kind*, *hsspeed*, *hsspeed\_MHz*, *preamp*, *preamp\_gain*), where *channel*, *oamp*, *hsspeed* and *preamp* are indices, while *channel\_bitdepth*, *oamp\_kind*, *hsspeed\_MHz* and *preamp\_gain* are descriptions.

**get\_max\_vsspeed()**

Get maximal recommended vertical scan speed

**set\_amp\_mode** (*channel=None*, *oamp=None*, *hsspeed=None*, *preamp=None*)

Setup preamp mode.

Can specify AD channel index, amplifier index, channel speed (horizontal scan speed) index and preamp gain index. *None* (default) means leaving the current value.

**get\_amp\_mode** (*full=True*)

Return the current amplifier mode.

If *full==True*, return a full description (e.g., actual preamp gain or channel name); otherwise, return just the essential indices information (enough to set the mode for this camera, but no explanations).

**set\_vsspeed** (*vsspeed*)

Set vertical scan speed index

**get\_channel()**

Get current channel index

**get\_channel\_bitdepth** (*channel=None*)

Get channel bit depth corresponding to the given channel index (current by default)

**get\_oamp()**

Get current output amplifier index

**get\_oamp\_desc** (*oamp=None*)

Get output amplifier kind corresponding to the given oamp index (current by default)

**get\_hsspeed()**

Get current horizontal speed index

**get\_hsspeed\_frequency** (*hsspeed=None*)

Get horizontal scan frequency (in Hz) corresponding to the given hsspeed index (current by default)

**get\_preamp()**

Get current preamp index

**get\_preamp\_gain** (*preamp=None*)

Get preamp gain corresponding to the given preamp index (current by default)

**get\_vsspeed()**

Get current vertical speed index

**get\_vsspeed\_period** (*vsspeed=None*)

Get vertical scan period corresponding to the given vsspeed index (current by default)

**get\_EMCCD\_gain()**

Get current EMCCD gain.

Return tuple (*gain*, *advanced*).



**set\_EMCCD\_gain** (*gain, advanced=None*)

Set EMCCD gain.

Gain goes up to 300 if `advanced==False` or higher if `advanced==True` (in this mode the sensor can be permanently damaged by strong light).

**init\_amp\_mode** (*mode=None*)

Initialize the camera channel, frequencies and amp settings to some default mode.

If *mode* is supplied, use this mode; otherwise, use the slowest, lowest gain mode (the first one returned by `get_all_amp_modes()`). Also set the maximal recommended vertical shift speed and no EMCCD gain.

**get\_min\_shutter\_times** ()

Get minimal shutter opening and closing times

**setup\_shutter** (*mode, ttl\_mode=0, open\_time=None, close\_time=None*)

Setup shutter.

*mode* can be "auto", "open" or "closed", *ttl\_mode* can be 0 (low is open) or 1 (high is open), *open\_time* and *close\_time* specify opening and closing times (required to calculate the minimal exposure times). By default, these time are minimal allowed times.

**get\_shutter\_parameters** ()

Return shutter parameters as a tuple (*mode, ttl\_mode, open\_time, close\_time*)

**get\_shutter** ()

Get shutter state ("auto", "open", or "closed")

**set\_fan\_mode** (*mode*)

Set fan mode.

Can be "full", "low" or "off".

**get\_fan\_mode** ()

Return fan mode ("full", "low", or "off")

**read\_in\_aux\_port** (*port*)

Get state at a given auxiliary port

**set\_out\_aux\_port** (*port, state*)

Set state at a given auxiliary port

**set\_trigger\_mode** (*mode*)

Set trigger mode.

Can be "int" (internal), "ext" (external), "ext\_start" (external start), "ext\_exp" (external exposure), "ext\_fvb\_em" (external FVB EM), "software" (software trigger) or "ext\_charge\_shift" (external charge shifting).

For description, see Andor SDK manual.

**get\_trigger\_mode** ()

Return trigger mode

**get\_trigger\_level\_limits** ()

Get limits on the trigger level

**setup\_ext\_trigger** (*level=None, invert=None, term\_highZ=None*)

Setup external trigger (level, inversion, and high-Z termination).

Any None values are not changed. If any returned values are None, it means that this option is not supported.

**get\_ext\_trigger\_parameters()**

Return external trigger parameters (level, inversion, high-Z termination).

If any returned values are None, it means that this option is not supported.

**send\_software\_trigger()**

Send software trigger signal

**set\_acquisition\_mode(mode, setup\_params=True)**

Set the acquisition mode.

Can be "single", "accum", "kinetic", "fast\_kinetic" or "cont" (continuous). If setup\_params==True, make sure that the last specified parameters for this mode are set up. For description of each mode, see Andor SDK manual and corresponding setup\_\*\_mode functions.

**get\_acquisition\_mode()**

Get the current acquisition mode

**setup\_accum\_mode(num\_acc, cycle\_time\_acc=0)**

Switch into the accum acquisition mode and set up its parameters.

*num\_acc* is the number of accumulated frames, *cycle\_time\_acc* is the acquisition period (by default the minimal possible based on exposure and transfer time).

**get\_accum\_mode\_parameters()**

Return accum acquisition mode parameters (num\_acc, cycle\_time\_acc)

**setup\_kinetic\_mode(num\_cycle, cycle\_time=0.0, num\_acc=1, cycle\_time\_acc=0, num\_prescan=0)**

Switch into the kinetic acquisition mode and set up its parameters.

*num\_cycle* is the number of kinetic cycles frames, *cycle\_time* is the acquisition period between accum frames, *num\_accum* is the number of accumulated frames, *cycle\_time\_acc* is the accum acquisition period, *num\_prescan* is the number of prescans.

**get\_kinetic\_mode\_parameters()**

Return kinetic acquisition mode parameters (num\_cycle, cycle\_time, num\_acc, cycle\_time\_acc, num\_prescan)

**setup\_fast\_kinetic\_mode(num\_acc, cycle\_time\_acc=0.0)**

Switch into the fast kinetic acquisition mode and set up its parameters.

*num\_acc* is the number of accumulated frames, *cycle\_time\_acc* is the acquisition period (by default the minimal possible based on exposure and transfer time).

**get\_fast\_kinetic\_mode\_parameters()**

Return fast kinetic acquisition mode parameters (num\_acc, cycle\_time\_acc)

**setup\_cont\_mode(cycle\_time=0)**

Switch into the continuous acquisition mode and set up its parameters.

*cycle\_time* is the acquisition period (by default the minimal possible based on exposure and transfer time).

**get\_cont\_mode\_parameters()**

Return continuous acquisition mode parameters *cycle\_time*

**set\_exposure(exposure)**

Set camera exposure

**get\_exposure()**

Get current exposure

**set\_frame\_period(frame\_period)**

Set frame acquisition period for the continuous mode

**enable\_frame\_transfer\_mode** (*enable=True*)

Enable frame transfer mode.

For description, see Andor SDK manual.

**is\_frame\_transfer\_enabled** ()

Return whether the frame transfer mode is enabled

**get\_cycle\_timings** ()

Get acquisition timing.

Return tuple (*exposure*, *accum\_cycle\_time*, *kinetic\_cycle\_time*). In continuous mode, the relevant cycle time is *kinetic\_cycle\_time*.

**get\_frame\_timings** ()

Get acquisition timing.

Return tuple (*exposure*, *frame\_period*). Frame period is the rate of frame generation, not of internal frame acquisition (e.g., in accumulator or kinetic mode this is the rate of generating a single accumulated frame, which is *num\_acc* times larger than the internal frame period).

**get\_readout\_time** ()

Get frame readout time

**get\_keepclean\_time** ()

Get sensor keep-clean time

**set\_read\_mode** (*mode*)

Set camera read mode.

Can be "fvb" (average all image vertically and return it as one row), "single\_track" (read a single row or several rows averaged together), "multi\_track" (read multiple rows or averaged sets of rows), "random\_track" (read several arbitrary lines), or "image" (read a whole image or its rectangular part).

**get\_read\_mode** ()

Get the current read mode

**setup\_single\_track\_mode** (*center=0*, *width=1*)

Switch into the singe-track read mode and set up its parameters.

*center* and *width* specify selection of the rows to be averaged together.

**get\_single\_track\_mode\_parameters** ()

Return singe-track read mode parameters (*center*, *width*)

**setup\_multi\_track\_mode** (*number=1*, *height=1*, *offset=1*)

Switch into the multi-track read mode and set up its parameters.

*number* is the number of rows (or row sets) to read, *height* is number of one row set (1 for a single row), *offset* is the distance between the row sets. Return a tuple (*number*, *height*, *offset*, *top*, *gap*), where *top* is the offset of the first row from the top, and *gap* is the gap between the tracks.

**get\_multi\_track\_mode\_parameters** ()

Return multi-track read mode parameters (*number*, *height*, *offset*)

**setup\_random\_track\_mode** (*tracks=None*)

Switch into the random-track read mode and set up its parameters.

*tracks* is a list of tuples (*start*, *stop*) specifying track span (start are inclusive, stop are exclusive, starting from 0). Note that it does not affect the current read mode, which should be set using [`set\_read\_mode\(\)`](#).

**get\_random\_track\_mode\_parameters()**

Return random-track read mode parameters, i.e., the list of track positions

**setup\_image\_mode** (*hstart=0, hend=None, vstart=0, vend=None, hbin=1, vbin=1*)

Switch into the image read mode and set up its parameters.

*hstart* and *hend* specify horizontal image extent, *vstart* and *vend* specify vertical image extent (start are inclusive, stop are exclusive, starting from 0), *hbin* and *vbin* specify binning. By default, all non-supplied parameters take extreme values.

**get\_image\_mode\_parameters()**

Return image read mode parameters, (*hstart, hend, vstart, vend, hbin, vbin*)

**get\_detector\_size()**

Get camera detector size (in pixels) as a tuple (*width, height*)

**get\_roi()**

Get current ROI.

Return tuple (*hstart, hend, vstart, vend, hbin, vbin*). *hstart* and *hend* specify horizontal image extent, *vstart* and *vend* specify vertical image extent (start is inclusive, stop is exclusive, starting from 0), *hbin* and *vbin* specify binning.

**set\_roi** (*hstart=0, hend=None, vstart=0, vend=None, hbin=1, vbin=1*)

Setup camera ROI.

*hstart* and *hend* specify horizontal image extent, *vstart* and *vend* specify vertical image extent (start is inclusive, stop is exclusive, starting from 0), *hbin* and *vbin* specify binning. By default, all non-supplied parameters take extreme values (0 for start, maximal for end, 1 for binning).

**get\_roi\_limits** (*hbin=1, vbin=1*)

Get the minimal and maximal ROI parameters.

Return tuple (*hlim, vlim*), where each element is in turn a limit 5-tuple (*min, max, pstep, sstep, maxbin*) with, correspondingly, minimal and maximal size, position and size step, and the maximal binning. In some cameras, the step and the minimal size depend on the binning, which can be supplied.

**setup\_acquisition** (*mode=None, nframes=None*)

Setup acquisition.

Any non-specified acquisition parameters are assumed to be the same as previously set (or default, if not explicitly set before). Return the new acquisition parameters.

**clear\_acquisition()**

Clear acquisition settings

**start\_acquisition** (*\*args, \*\*kwargs*)

Start acquisition.

Can take the same keyword parameters as :meth:“*setup\_acquisition*”. If the acquisition is not set up yet, set it up using the supplied parameters (use default of *setup\_acquisition()*, if the parameter is None). Otherwise, if any supplied parameters are different from the current ones, change them and reset the acquisition.

**stop\_acquisition()**

Stop acquisition

**get\_acquisition\_progress()**

Get acquisition progress.

Return tuple (*frames\_done, acc\_done*) with the number of full transferred frames and the number of acquired sub-frames in the current accumulation cycle.

**get\_buffer\_size()**

Get the size of the image ring buffer

**FrameTransferError**

alias of `pylablib.devices.interface.camera.DefaultFrameTransferError`

**class NoParameterCaller** (*device, kind*)

Bases: `object`

Class to simplify calling functions without a parameter

**apply\_settings** (*settings*)

Apply the settings.

*settings* is the dict {name: value} of the device available settings. Non-applicable settings are ignored.

**get\_acquisition\_parameters()**

Get acquisition parameters.

Return dictionary {name: value}

**get\_data\_dimensions()**

Get readout data dimensions (in pixels) as a tuple (width, height); take indexing mode into account

**get\_device\_variable** (*key*)

Get the value of a settings, status, or full info parameter

**get\_frame\_format()**

Get format for the returned images.

Can be "list" (list of 2D arrays), "array" (a single 3D array), or "chunks" (list of 3D "chunk" arrays; supported for some cameras and provides the best performance).

**get\_frame\_info\_fields()**

Get the names of frame info fields.

Applicable when frame info format (set by `set_frame_info_format()`) is "list" or "array".

**get\_frame\_info\_format()**

Get format of the frame info.

Can be "namedtuple" (potentially nested named tuples; convenient to get particular values), "list" (flat list of values, with field names are given by `get_frame_info_fields()`; convenient for building a table), "array" (same as "list", but with a numpy array, which is easier to use for "chunks" frame format), or "dict" (flat dictionary with the same fields as the "list" format; more resilient to future format changes)

**get\_frame\_info\_period()**

Get period of frame info acquisition.

Frame info might be skipped (set to None) except for frames which indices are divisible by *period*. Useful for certain cameras where acquiring frame info takes a lot of time and can reduce performance at higher frame rates. Note that this parameter can still be ignored (i.e., always set to 1) if the performance is not an issue for a given camera class.

**get\_frame\_period()**

Get frame period (time between two consecutive frames in the internal trigger mode)

**get\_frames\_status()**

Get acquisition and buffer status.

Return tuple (acquired, unread, skipped, size), where acquired is the total number of acquired frames, unread is the number of acquired but not read frames, skipped is the number of skipped (not read and then written over) frames, and buffer\_size is the total buffer size (in frames).

**get\_full\_info** (include=0)

Get dict {name: value} containing full device information (including status and settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting include=-10 queries all available variables, which is equivalent to include="all".

**get\_full\_status** (include=0)

Get dict {name: value} containing the device status (including settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting include=-10 queries all available variables, which is equivalent to include="all".

**get\_image\_indexing** ()

Get indexing for the returned images.

Can be "rct" (first index row, second index column, rows counted from the top), "rcb" (same as "rc", rows counted from the bottom), "xyt" (first index column, second index row, rows counted from the top), or "xyb" (same as "xyt", rows counted from the bottom)

**get\_new\_images\_range** ()

Get the range of the new images.

Return tuple (first, last) with images range (first inclusive). If no images are available, return None. If some images were in the buffer were overwritten, exclude them from the range.

**get\_settings** (include=0)

Get dict {name: value} containing all the device settings.

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting include=-10 queries all available variables, which is equivalent to include="all".

**grab** (nframes=1, frame\_timeout=5.0, missing\_frame='skip', return\_info=False, buff\_size=None)

Snap *nframes* images (with preset image read mode parameters)

*buff\_size* determines buffer size (if None, use the default size). Timeout is specified for a single-frame acquisition, not for the whole acquisition time. *missing\_frame* determines what to do with frames which have been lost: can be "none" (replacing them with None), "zero" (replacing them with zero-filled frame), or "skip" (skipping them, while still keeping total returned frames number to *n*). If return\_info==True, return tuple (frames, infos), where infos is a list of frame info tuples (camera-dependent); if some frames are missing and missing\_frame!="skip", the corresponding frame info is None.

**is\_acquisition\_setup** ()

Check if acquisition is set up.

If the camera does not support separate acquisition setup, always return True.

**pausing\_acquisition** (clear=None, stop=True, setup\_after=None, start\_after=True, combine\_nested=True)

Context manager which temporarily pauses acquisition during execution of with block.

Useful for applying certain settings which can't be changed during the acquisition. If `clear==True`, clear acquisition in addition to stopping (by default, use the class default specified as `_clear_pausing_acquisition` attribute). If `stop==True`, stop the acquisition (if `clear==True`, stop regardless). If `setup_after==True`, setup the acquisition after pause if necessary (None means setup only if clearing was required). If `start_after==True`, start the acquisition after pause if necessary (None means start only if stopping was required). If `combine_nested==True`, then any nested `pausing_acquisition` calls will stop/clear acquisition as necessary, but won't setup/start it again until this `pausing_acquisition` call is complete.

Yields tuple (`acq_in_progress`, `acq_params`), which indicates whether acquisition is currently in progress, and what are the current acquisition parameters.

**read\_multiple\_images** (*rng=None*, *peek=False*, *missing\_frame='skip'*, *return\_info=False*, *return\_rng=False*)

Read multiple images specified by *rng* (by default, all un-read images).

If *rng* is specified, it is a tuple (*first*, *last*) with images range (*first* inclusive). If no new frames are available, return an empty list; if no acquisition is running, return None. If *peek==True*, return images but not mark them as read. *missing\_frame* determines what to do with frames which are out of range (missing or lost): can be "none" (replacing them with None), "zero" (replacing them with zero-filled frame), or "skip" (skipping them). If *return\_info==True*, return tuple (*frames*, *infos*), where *infos* is a list of frame info tuples (camera-dependent, by default, only the frame index); if some frames are missing and *missing\_frame!="skip"*, the corresponding frame info is None. If *return\_rng==True*, return the range covered resulting frames; if *missing\_frame=="skip"*, the range can be smaller than the supplied *rng* if some frames are skipped.

**read\_newest\_image** (*peek=False*, *return\_info=False*)

Read the newest un-read image.

If no un-read frames are available, return None. If *peek==True*, return the image but not mark it as read. If *return\_info==True*, return tuple (*frame*, *info*), where *info* is an info tuples (camera-dependent, see `read_multiple_images()`).

**read\_oldest\_image** (*peek=False*, *return\_info=False*)

Read the oldest un-read image.

If no un-read frames are available, return None. If *peek==True*, return the image but not mark it as read. If *return\_info==True*, return tuple (*frame*, *info*), where *info* is an info tuples (camera-dependent, see `read_multiple_images()`).

**set\_device\_variable** (*key*, *value*)

Set the value of a settings parameter

**set\_frame\_format** (*fnt*)

Set format for the returned images.

Can be "list" (list of 2D arrays), "array" (a single 3D array), "chunks" (list of 3D "chunk" arrays; supported for some cameras and provides the best performance), or "try\_chunks" (same as "chunks", but if chunks are not supported, set to "list" instead). If format is "chunks" and chunks are not supported by the camera, it results in one frame per chunk. Note that if the format is set to "array" or "chunks", the frame info format is also automatically set to "array". If the format is set to "chunks", then the image info is also returned in chunks form (list of 2D info arrays with the same length as the corresponding frame chunks).

**set\_frame\_info\_format** (*fnt*, *include\_fields=None*)

Set format of the frame info.

Can be "namedtuple" (potentially nested named tuples; convenient to get particular values), "list" (flat list of values, with field names are given by `get_frame_info_fields()`; convenient for building a table), "array" (same as "list", but with a numpy array, which is easier to use for "chunks"

frame format), or "dict" (flat dictionary with the same fields as the "list" format; more resilient to future format changes) If *include\_fields* is not None, it specifies the fields included for non-"tuple" formats; note that order of *include\_fields* is ignored, and the resulting fields are always ordered same as in the original.

**set\_frame\_info\_period** (*period=1*)

Set period of frame info acquisition.

Frame info might be skipped (set to None) except for frames which indices are divisible by *period*. Useful for certain cameras where acquiring frame info takes a lot of time and can reduce performance at higher frame rates. Note that this parameter can still be ignored (i.e., always set to 1) if the performance is not an issue for a given camera class.

**set\_image\_indexing** (*indexing*)

Set up indexing for the returned images.

Can be "rct" (first index row, second index column, rows counted from the top), "rcb" (same as "rc", rows counted from the bottom), "xyt" (first index column, second index row, rows counted from the top), or "xyb" (same as "xyt", rows counted from the bottom)

**snap** (*timeout=5.0, return\_info=False*)

Snap a single frame

**wait\_for\_frame** (*since='lastread', nframes=1, timeout=20.0, error\_on\_stopped=False*)

Wait for one or several new camera frames.

*since* specifies the reference point for waiting to acquire *nframes* frames; can be "lastread" (from the last read frame), "lastwait" (wait for the last successful *wait\_for\_frame()* call), "now" (from the start of the current call), or "start" (from the acquisition start, i.e., wait until *nframes* frames have been acquired). *timeout* can be either a number, None (infinite timeout), or a tuple (timeout, frame\_timeout), in which case the call times out if the total time exceeds timeout, or a single frame wait exceeds frame\_timeout. If the call times out, raise `TimeoutError`. If *error\_on\_stopped==True* and the acquisition is not running, raise `Error`; otherwise, simply return False without waiting.

## pylablib.devices.Andor.AndorSDK3 module

**class** `pylablib.devices.Andor.AndorSDK3.LibraryController` (*lib*)

Bases: `pylablib.devices.utils.load_lib.LibraryController`

**close** (*opid*)

Mark device closing.

Return tuple (close\_result, uninit\_result) with the results of the closing and the shutdown. If library does not need to be shut down yet, set uninit\_result=None

**open** ()

Mark device opening.

Return tuple (init\_result, open\_result, opid) with the results of the initialization and the opening, and the opening ID which should afterwards be used for closing. If library is already initialized, set init\_result=None

**preinit** ()

Pre-initialize the library, if it hasn't been done already

**shutdown** ()

Close all opened connections and shutdown the library



**temp\_open()**

Context for temporarily opening a new device connection

`pylablib.devices.Andor.AndorSDK3.restart_lib()`

`pylablib.devices.Andor.AndorSDK3.get_cameras_number()`

Get number of connected Andor cameras

**class** `pylablib.devices.Andor.AndorSDK3.AndorSDK3Attribute` (*handle*, *name*,  
*kind='auto'*)

Bases: `object`

Andor SDK3 camera attribute.

Allows to query and set values and get additional information. Usually created automatically by a Andor SDK3 camera instance, but could also be created manually.

#### Parameters

- **handle** – Andor SDK3 camera handle
- **pid** – attribute id
- **kind** – attribute kind; can be "float", "int", "str", "bool", "enum", or "comm" (command); can also be "auto" (default), in which case it is obtained from the stored feature table; newer features might be missing, in which case kind needs to be supplied explicitly, or it raises an error

**name**

attribute name

**kind**

attribute kind; can be "float", "int", "str", "bool", "enum", or "comm" (command)

**implemented**

whether attribute is implemented

**Type** `bool`

**readable**

whether attribute is readable

**Type** `bool`

**writable**

whether attribute is writable

**Type** `bool`

**min**

minimal attribute value (if applicable)

**Type** `float` or `int`

**max**

maximal attribute value (if applicable)

**Type** `float` or `int`

**ivalues**

list of possible integer values for enum attributes

**values**

list of possible text values for enum attributes

**labels**

dict {label: index} which shows all possible values of an enumerated attribute and their corresponding numerical values

**ilabels**

dict {index: label} which shows labels corresponding to numerical values of an enumerated attribute

**is\_command**

whether attribute is a command (same as kind=="comm")

Type `bool`

**update\_properties()**

Update all attribute properties: implemented, readable, writable, limits

**get\_value** (*enum\_as\_str=True, not\_implemented\_error=True, default=None*)

Get current value.

If `enum_as_str==True`, return enum values as strings; otherwise, return as indices. If `not_implemented_error==True` and the feature is not implemented, raise `AndorError`; otherwise, return *default* if it is not implemented.

**set\_value** (*value, not\_implemented\_error=True*)

Set current value.

If `not_implemented_error==True` and the feature is not implemented, raise `AndorError`; otherwise, do nothing.

**call\_command()**

Execute the given command

**get\_range** (*enum\_as\_str=True*)

Get allowed range of the given value.

For "int" or "float" values return tuple (min, max) (inclusive); for "enum" return list of possible values (if `enum_as_str==True`, return list of string values, otherwise return list of indices). For all other value kinds return None.

**update\_limits()**

Update minimal and maximal attribute limits and return tuple (min, max)

**truncate\_value** (*value*)

Limit value to lie within the allowed range

```
class pylablib.devices.Andor.AndorSDK3.TDeviceInfo(camera_name, camera_model,  
                                                serial_number, firmware_version,  
                                                software_version)
```

Bases: `tuple`

**camera\_model****camera\_name****count()**

Return number of occurrences of value.

**firmware\_version****index()**

Return first index of value.

Raises `ValueError` if the value is not present.

**serial\_number**

**software\_version**

**class** pylablib.devices.Andor.AndorSDK3.**TMissedFramesStatus** (*skipped, overflows*)  
 Bases: `tuple`

**count**()  
 Return number of occurrences of value.

**index**()  
 Return first index of value.  
 Raises ValueError if the value is not present.

**overflows**

**skipped**

**class** pylablib.devices.Andor.AndorSDK3.**TFrameInfo** (*frame\_index, timestamp\_dev, size, pixeltype, stride*)

Bases: `tuple`

**count**()  
 Return number of occurrences of value.

**frame\_index**

**index**()  
 Return first index of value.  
 Raises ValueError if the value is not present.

**pixeltype**

**size**

**stride**

**timestamp\_dev**

**class** pylablib.devices.Andor.AndorSDK3.**AndorSDK3Camera** (*idx=0*)  
 Bases: `pylablib.devices.interface.camera.IBinROICamera`, `pylablib.devices.interface.camera.IExposureCamera`, `pylablib.devices.interface.camera.IAttributeCamera`

Andor SDK3 camera.

**Parameters** *idx* (*int*) – camera index (use `get_cameras_number()` to get the total number of connected cameras)

**Error**  
 alias of `pylablib.devices.Andor.base.AndorError`

**TimeoutError**  
 alias of `pylablib.devices.Andor.base.AndorTimeoutError`

**FrameTransferError**  
 alias of `pylablib.devices.Andor.base.AndorFrameTransferError`

**open**()  
 Open connection to the camera

**close**()  
 Close connection to the camera

**is\_opened**()  
 Check if the device is connected

**add\_attribute** (*name*, *kind*)

Add a new attribute which is not currently present in the dictionary.

*kind* can be "float", "int", "str", "bool", "enum", or "comm" (command).

**get\_attribute** (*name*, *update\_properties=False*, *error\_on\_missing=True*)

Get the camera attribute with the given name.

If *update\_properties==True*, automatically update all attribute properties.

**get\_attribute\_value** (*name*, *enum\_as\_str=True*, *update\_properties=False*, *error\_on\_missing=True*, *default=None*)

Get value of an attribute with the given name.

If *update\_properties==True*, automatically update all attribute properties before settings. If the value doesn't exist or can not be read and *error\_on\_missing==True*, raise error; otherwise, return *default*. If *default* is not None, assume that *error\_on\_missing==False*.

**set\_attribute\_value** (*name*, *value*, *update\_properties=True*, *error\_on\_missing=True*)

Set value of an attribute with the given name.

If the value doesn't exist or can not be written and *error\_on\_missing==True*, raise error; otherwise, do nothing. If *update\_properties==True*, automatically update all attribute properties before settings.

**get\_all\_attribute\_values** (*root="*, *enum\_as\_str=True*, *update\_properties=False*)

Get values of all attributes.

If *update\_properties==True*, automatically update all attribute properties before settings.

**set\_all\_attribute\_values** (*settings*, *update\_properties=True*)

Set values of all attribute in the given dictionary.

If *update\_properties==True*, automatically update all attribute properties before settings.

**call\_command** (*name*)

Execute the given command

**get\_device\_info** ()

Get camera info.

Return tuple (camera\_name, camera\_model, serial\_number, firmware\_version, software\_version).

**get\_trigger\_mode** ()

Get trigger mode.

Can be "int" (internal), "ext" (external), "software" (software trigger), "ext\_start" (external start), or "ext\_exp" (external exposure).

**set\_trigger\_mode** (*mode*)

Set trigger mode.

Can be "int" (internal), "ext" (external), or "software" (software trigger).

**get\_shutter** ()

Get current shutter mode

**set\_shutter** (*mode*)

Set trigger mode.

Can be "open", "closed", or "auto".

**is\_cooler\_on** ()

Check if the cooler is on

**set\_cooler** (*on=True*)  
Set the cooler on or off

**get\_temperature** ()  
Get the current camera temperature

**get\_temperature\_setpoint** ()  
Get current temperature setpoint

**set\_temperature** (*temperature, enable\_cooler=True*)  
Change the temperature setpoint.  
  
If *enable\_cooler==True*, turn the cooler on automatically.

**get\_exposure** ()  
Get current exposure

**set\_exposure** (*exposure*)  
Set camera exposure

**get\_frame\_period** ()  
Get frame period (time between two consecutive frames in the internal trigger mode)

**set\_frame\_period** (*frame\_period*)  
Set frame period (time between two consecutive frames in the internal trigger mode)

**get\_frame\_timings** ()  
Get acquisition timing.  
  
Return tuple (*exposure, frame\_period*).

**is\_metadata\_enabled** ()  
Check if the metadata enabled

**enable\_metadata** (*enable=True*)  
Enable or disable metadata streaming

**class BufferManager** (*cam*)  
Bases: `object`  
  
Buffer manager: stores, constantly reads and re-schedules buffers, keeps track of acquired frames and buffer overflow events

**allocate\_buffers** (*nbuff, size, queued\_buffers=None*)  
Allocate and queue buffers.  
  
*queued\_buffers* specifies number of allocated buffers to keep queued at a given time (by default, all of them)

**queue\_buffer** (*queued\_buffers*)  
Queue extra buffers (only before any buffers are read)

**deallocate\_buffers** ()  
Deallocated buffers (flushing should be done manually)

**reset** ()  
Reset counter (on frame acquisition)

**read** (*idx*)  
Return the oldest available acquired but not read buffer, and mark it as read

**start\_loop** ()  
Start buffer scheduling loop

**stop\_loop()**  
Stop buffer scheduling loop

**wait\_for\_frame** (*idx=None, timeout=None*)  
Wait for a new frame acquisition

**on\_overflow()**  
Process buffer overflow event

**new\_overflow()**

**get\_status()**  
Get counter status: tuple (acquired, total\_length)

**setup\_acquisition** (*mode='sequence', nframes=100*)  
Setup acquisition.

*mode* can be either "snap" (single frame or a fixed number of frames) or "sequence" (continuous acquisition). *nframes* determines number of frames to acquire in the single mode, or size of the ring buffer in the "sequence" mode (by default, 100).

**clear\_acquisition()**  
Clear acquisition settings

**start\_acquisition** (*\*args, \*\*kwargs*)  
Start acquisition.

Can take the same keyword parameters as :meth:“*setup\_acquisition*. If the acquisition is not set up yet, set it up using the supplied parameters (use default of *setup\_acquisition()*, if the parameter is None). Otherwise, if any supplied parameters are different from the current ones, change them and reset the acquisition.

**stop\_acquisition()**  
Stop acquisition

**acquisition\_in\_progress()**  
Check if acquisition is in progress

**get\_missed\_frames\_status()**  
Get missed frames status.

Return tuple (*skipped, overflows*) with the number skipped frames (sent from camera to the PC, but not read and overwritten) and number of buffer overflows (events when the frame rate is too for the data transfer, so some unknown number of frames is skipped).

**reset\_overflows\_counter()**  
Reset buffer overflows counter

**set\_overflow\_behavior** (*behavior*)  
Choose the camera behavior if buffer overflow is encountered when waiting for a new frame.

Can be "error" (raise *AndorFrameTransferError*), "restart" (restart the acquisition), or "ignore" (ignore the overflow, which will cause the wait to time out).

**get\_detector\_size()**  
Get camera detector size (in pixels) as a tuple (*width, height*)

**get\_roi()**  
Get current ROI.

Return tuple (*hstart, hend, vstart, vend, hbin, vbin*).

**set\_roi** (*hstart=0, hend=None, vstart=0, vend=None, hbin=1, vbin=1*)  
Set current ROI.

By default, all non-supplied parameters take extreme values. Binning is the same for both axes.

**get\_roi\_limits** (*hbin=1, vbin=1*)

Get the minimal and maximal ROI parameters.

Return tuple (*hlim, vlim*), where each element is in turn a limit 5-tuple (*min, max, pstep, sstep, maxbin*) with, correspondingly, minimal and maximal size, position and size step, and the maximal binning.

Note that the minimal ROI size depends on the current (not just supplied) binning settings. For more accurate results, is it only after setting up the binning.

**class NoParameterCaller** (*device, kind*)

Bases: `object`

Class to simplify calling functions without a parameter

**apply\_settings** (*settings*)

Apply the settings.

*settings* is the dict {*name*: *value*} of the device available settings. Non-applicable settings are ignored.

**get\_acquisition\_parameters** ()

Get acquisition parameters.

Return dictionary {*name*: *value*}

**get\_all\_attributes** (*copy=False*)

Return a dictionary of all available attributes.

If *copy==True*, copy the dictionary; otherwise, return the internal dictionary structure (should not be modified).

**get\_data\_dimensions** ()

Get readout data dimensions (in pixels) as a tuple (*width, height*); take indexing mode into account

**get\_device\_variable** (*key*)

Get the value of a settings, status, or full info parameter

**get\_frame\_format** ()

Get format for the returned images.

Can be "list" (list of 2D arrays), "array" (a single 3D array), or "chunks" (list of 3D "chunk" arrays; supported for some cameras and provides the best performance).

**get\_frame\_info\_fields** ()

Get the names of frame info fields.

Applicable when frame info format (set by `set_frame_info_format()`) is "list" or "array".

**get\_frame\_info\_format** ()

Get format of the frame info.

Can be "namedtuple" (potentially nested named tuples; convenient to get particular values), "list" (flat list of values, with field names are given by `get_frame_info_fields()`; convenient for building a table), "array" (same as "list", but with a numpy array, which is easier to use for "chunks" frame format), or "dict" (flat dictionary with the same fields as the "list" format; more resilient to future format changes)

**get\_frame\_info\_period** ()

Get period of frame info acquisition.

Frame info might be skipped (set to `None`) except for frames which indices are divisible by *period*. Useful for certain cameras where acquiring frame info takes a lot of time and can reduce performance at higher frame rates. Note that this parameter can still be ignored (i.e., always set to 1) if the performance is not an issue for a given camera class.

**get\_frames\_status()**

Get acquisition and buffer status.

Return tuple (*acquired*, *unread*, *skipped*, *size*), where *acquired* is the total number of acquired frames, *unread* is the number of acquired but not read frames, *skipped* is the number of skipped (not read and then written over) frames, and *buffer\_size* is the total buffer size (in frames).

**get\_full\_info(include=0)**

Get dict {*name*: *value*} containing full device information (including status and settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include*=-10 queries all available variables, which is equivalent to *include*="all".

**get\_full\_status(include=0)**

Get dict {*name*: *value*} containing the device status (including settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include*=-10 queries all available variables, which is equivalent to *include*="all".

**get\_image\_indexing()**

Get indexing for the returned images.

Can be "rct" (first index row, second index column, rows counted from the top), "rcb" (same as "rc", rows counted from the bottom), "xyt" (first index column, second index row, rows counted from the top), or "xyb" (same as "xyt", rows counted from the bottom)

**get\_new\_images\_range()**

Get the range of the new images.

Return tuple (*first*, *last*) with images range (first inclusive). If no images are available, return `None`. If some images were in the buffer were overwritten, exclude them from the range.

**get\_settings(include=0)**

Get dict {*name*: *value*} containing all the device settings.

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include*=-10 queries all available variables, which is equivalent to *include*="all".

**grab(nframes=1, frame\_timeout=5.0, missing\_frame='skip', return\_info=False, buff\_size=None)**

Snap *nframes* images (with preset image read mode parameters)

*buff\_size* determines buffer size (if `None`, use the default size). Timeout is specified for a single-frame acquisition, not for the whole acquisition time. *missing\_frame* determines what to do with frames which have been lost: can be "none" (replacing them with `None`), "zero" (replacing them with zero-filled frame), or "skip" (skipping them, while still keeping total returned frames number to *n*). If *return\_info*==`True`, return tuple (*frames*, *infos*), where *infos* is a list of frame info tuples (camera-dependent); if some frames are missing and *missing\_frame*!="skip", the corresponding frame info is `None`.



**is\_acquisition\_setup()**

Check if acquisition is set up.

If the camera does not support separate acquisition setup, always return `True`.

**pausing\_acquisition**(*clear=None, stop=True, setup\_after=None, start\_after=True, combine\_nested=True*)

Context manager which temporarily pauses acquisition during execution of `with` block.

Useful for applying certain settings which can't be changed during the acquisition. If `clear==True`, clear acquisition in addition to stopping (by default, use the class default specified as `_clear_pausing_acquisition` attribute). If `stop==True`, stop the acquisition (if `clear==True`, stop regardless). If `setup_after==True`, setup the acquisition after pause if necessary (None means setup only if clearing was required). If `start_after==True`, start the acquisition after pause if necessary (None means start only if stopping was required). If `combine_nested==True`, then any nested `pausing_acquisition` calls will stop/clear acquisition as necessary, but won't setup/start it again until this `pausing_acquisition` call is complete.

Yields tuple (`acq_in_progress, acq_params`), which indicates whether acquisition is currently in progress, and what are the current acquisition parameters.

**read\_newest\_image**(*peek=False, return\_info=False*)

Read the newest un-read image.

If no un-read frames are available, return `None`. If `peek==True`, return the image but not mark it as read. If `return_info==True`, return tuple (`frame, info`), where `info` is an info tuples (camera-dependent, see `read_multiple_images()`).

**read\_oldest\_image**(*peek=False, return\_info=False*)

Read the oldest un-read image.

If no un-read frames are available, return `None`. If `peek==True`, return the image but not mark it as read. If `return_info==True`, return tuple (`frame, info`), where `info` is an info tuples (camera-dependent, see `read_multiple_images()`).

**set\_device\_variable**(*key, value*)

Set the value of a settings parameter

**set\_frame\_format**(*fnt*)

Set format for the returned images.

Can be "list" (list of 2D arrays), "array" (a single 3D array), "chunks" (list of 3D "chunk" arrays; supported for some cameras and provides the best performance), or "try\_chunks" (same as "chunks", but if chunks are not supported, set to "list" instead). If format is "chunks" and chunks are not supported by the camera, it results in one frame per chunk. Note that if the format is set to "array" or "chunks", the frame info format is also automatically set to "array". If the format is set to "chunks", then the image info is also returned in chunks form (list of 2D info arrays with the same length as the corresponding frame chunks).

**set\_frame\_info\_format**(*fnt, include\_fields=None*)

Set format of the frame info.

Can be "namedtuple" (potentially nested named tuples; convenient to get particular values), "list" (flat list of values, with field names are given by `get_frame_info_fields()`; convenient for building a table), "array" (same as "list", but with a numpy array, which is easier to use for "chunks" frame format), or "dict" (flat dictionary with the same fields as the "list" format; more resilient to future format changes) If `include_fields` is not `None`, it specifies the fields included for non-"tuple" formats; note that order or `include_fields` is ignored, and the resulting fields are always ordered same as in the original.

**set\_frame\_info\_period** (*period=1*)

Set period of frame info acquisition.

Frame info might be skipped (set to `None`) except for frames which indices are divisible by *period*. Useful for certain cameras where acquiring frame info takes a lot of time and can reduce performance at higher frame rates. Note that this parameter can still be ignored (i.e., always set to 1) if the performance is not an issue for a given camera class.

**set\_image\_indexing** (*indexing*)

Set up indexing for the returned images.

Can be "rct" (first index row, second index column, rows counted from the top), "rcb" (same as "rc", rows counted from the bottom), "xyt" (first index column, second index row, rows counted from the top), or "xyb" (same as "xyt", rows counted from the bottom)

**snap** (*timeout=5.0, return\_info=False*)

Snap a single frame

**wait\_for\_frame** (*since='lastread', nframes=1, timeout=20.0, error\_on\_stopped=False*)

Wait for one or several new camera frames.

*since* specifies the reference point for waiting to acquire *nframes* frames; can be "lastread" (from the last read frame), "lastwait" (wait for the last successful `wait_for_frame()` call), "now" (from the start of the current call), or "start" (from the acquisition start, i.e., wait until *nframes* frames have been acquired). *timeout* can be either a number, `None` (infinite timeout), or a tuple (*timeout*, *frame\_timeout*), in which case the call times out if the total time exceeds *timeout*, or a single frame wait exceeds *frame\_timeout*. If the call times out, raise `TimeoutError`. If *error\_on\_stopped*==`True` and the acquisition is not running, raise `Error`; otherwise, simply return `False` without waiting.

**read\_multiple\_images** (*rng=None, peek=False, missing\_frame='skip', return\_info=False, return\_rng=False*)

Read multiple images specified by *rng* (by default, all un-read images).

If *rng* is specified, it is a tuple (*first*, *last*) with images range (first inclusive). If no new frames are available, return an empty list; if no acquisition is running, return `None`. If *peek*==`True`, return images but not mark them as read. *missing\_frame* determines what to do with frames which are out of range (missing or lost): can be "none" (replacing them with `None`), "zero" (replacing them with zero-filled frame), or "skip" (skipping them). If *return\_info*==`True`, return tuple (*frames*, *infos*), where *infos* is a list of `TFrameInfo` instances describing frame index and frame metadata, which contains timestamp, image size, pixel format, and row stride; if some frames are missing and *missing\_frame*!="skip", the corresponding frame info is `None`. if *return\_rng*==`True`, return the range covered resulting frames; if *missing\_frame*=="skip", the range can be smaller than the supplied *rng* if some frames are skipped.

## pylablib.devices.Andor.atcore\_features module

## pylablib.devices.Andor.base module

**exception** `pylablib.devices.Andor.base.AndorError`

Bases: `pylablib.core.devio.base.DeviceError`

Generic Andor error

**args**

**with\_traceback** ()

Exception.with\_traceback(tb) – set self.\_\_traceback\_\_ to tb and return self.

**exception** `pylablib.devices.Andor.base.AndorTimeoutError`

Bases: `pylablib.devices.Andor.base.AndorError`

Andor timeout error

**args**

**with\_traceback()**

Exception.with\_traceback(tb) – set self.\_\_traceback\_\_ to tb and return self.

**exception** `pylablib.devices.Andor.base.AndorFrameTransferError`

Bases: `pylablib.devices.Andor.base.AndorError`

Andor frame transfer error

**args**

**with\_traceback()**

Exception.with\_traceback(tb) – set self.\_\_traceback\_\_ to tb and return self.

**exception** `pylablib.devices.Andor.base.AndorNotSupportedError`

Bases: `pylablib.devices.Andor.base.AndorError`

Option not supported error

**args**

**with\_traceback()**

Exception.with\_traceback(tb) – set self.\_\_traceback\_\_ to tb and return self.

## Module contents

### `pylablib.devices.Arcus` package

#### Submodules

#### `pylablib.devices.Arcus.base` module

**exception** `pylablib.devices.Arcus.base.ArcusError`

Bases: `pylablib.core.devio.base.DeviceError`

Generic Arcus error

**args**

**with\_traceback()**

Exception.with\_traceback(tb) – set self.\_\_traceback\_\_ to tb and return self.

**exception** `pylablib.devices.Arcus.base.ArcusBackendError` (*exc*)

Bases: `pylablib.devices.Arcus.base.ArcusError`, `pylablib.core.devio.comm_backend.DeviceBackendError`

Generic Arcus backend communication error

**args**

**with\_traceback()**

Exception.with\_traceback(tb) – set self.\_\_traceback\_\_ to tb and return self.

**pylablib.devices.Arcus.performax module**

`pylablib.devices.Arcus.performax.get_usb_device_info (devid)`

Get info for the given device index (starting from 0).

Return tuple (index, serial, model, desc, vid, pid).

`pylablib.devices.Arcus.performax.list_usb_performax_devices ()`

List all performax devices.

Return list of tuples (index, serial, model, desc, vid, pid), one per device.

**class** `pylablib.devices.Arcus.performax.GenericPerformaxStage (idx=0, conn=None)`

Bases: `pylablib.devices.interface.stage.IMultiaxisStage`

Generic Arcus Performax translation stage.

**Parameters**

- **idx** (*int*) – stage index; if using a USB connection, specifies a USB device index; if using RS485 connection, specifies device index on the bus
- **conn** – if not `None`, defines a connection to RS485 connection. Usually (e.g., for USB-to-RS485 adapters) this is a serial connection, which either a name (e.g., "COM1"), or a tuple (name, baudrate) (e.g., ("COM1", 9600)); if *conn* is `None`, assume direct USB connection and use the manufacturer-provided DLL

**Error**

alias of `pylablib.devices.Arcus.base.ArcusError`

**open ()**

Open the connection to the stage

**close ()**

Close the connection to the stage

**is\_opened ()**

Check if the device is connected

**get\_device\_info ()**

Get the device info

**query (comm)**

Send a query to the stage and return the reply

**get\_device\_number ()**

Get the device number used in RS-485 communications.

Usually it is a string with the format similar to "4EX00".

**set\_device\_number (number, store=True)**

Get the device number used in RS-485 communications.

*number* can be either a full device id (e.g., "4EX00"), or a single number between 0 and 99. In order for the change to take effect, the device needs to be power-cycled. If *store*==`True`, automatically store settings to the memory; otherwise, the settings will be lost unless `store_defaults ()` is called at some point before the power-cycle.

**store\_defaults ()**

Store some of the settings to the memory as defaults.

Applies to device number, baudrate, limit error behavior, polarity, and some other settings.

**class NoParameterCaller** (*device, kind*)

Bases: `object`

Class to simplify calling functions without a parameter

**apply\_settings** (*settings*)

Apply the settings.

*settings* is the dict {*name*: *value*} of the device available settings. Non-applicable settings are ignored.

**get\_all\_axes** ()

Get the list of all available axes (taking mapping into account)

**get\_device\_variable** (*key*)

Get the value of a settings, status, or full info parameter

**get\_full\_info** (*include=0*)

Get dict {*name*: *value*} containing full device information (including status and settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_full\_status** (*include=0*)

Get dict {*name*: *value*} containing the device status (including settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_settings** (*include=0*)

Get dict {*name*: *value*} containing all the device settings.

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**remap\_axes** (*mapping, accept\_original=True*)

Rename axes to the new labels.

*mapping* is the new axes mapping, which can be a list of new axes name (corresponding to the old axes in order returned by `get_all_axes()`), or a dictionary {*alias*: *original*} of the new axes aliases.

**set\_device\_variable** (*key, value*)

Set the value of a settings parameter

**class** `pylablib.devices.Arcus.performax.Performax4EXStage` (*idx=0, conn=None, enable=True*)

Bases: `pylablib.devices.Arcus.performax.GenericPerformaxStage`

Arcus Performax 4EX/4ET translation stage.

#### Parameters

- **idx** (*int*) – stage index; if using a USB connection, specifies a USB device index; if using RS485 connection, specifies device index on the bus
- **conn** – if not `None`, defines a connection to RS485 connection. Usually (e.g., for USB-to-RS485 adapters) this is a serial connection, which either a name (e.g., "COM1"), or

a tuple (name, baudrate) (e.g., ("COM1", 9600)); if *conn* is None, assume direct USB connection and use the manufacturer-provided DLL

- **enable** – if True, enable all axes on startup

**get\_baudrate()**

Get current baud rate

**set\_baudrate**(*baudrate*, *store=True*)

Set current baud rate.

Acceptable values are 9600 (default), 19200, 38400, 57600, and 115200. In order for the change to take effect, the device needs to be power-cycled. If *store==True*, automatically store settings to the memory; otherwise, the settings will be lost unless *store\_defaults()* is called at some point before the power-cycle.

**enable\_absolute\_mode**(*enable=True*)

Set absolute motion mode

**enable\_limit\_errors**(*enable=True*, *autoclear=True*)

Enable limit errors.

If on, reaching limit switch on an axis puts it into an error state, which immediately stops this and all other axes; any further motion command on this axis will raise an error (it is still possible to restart motion on other axes); the axis motion can only be resumed by calling *clear\_limit\_error()*. If off, the limited axis still stops, but the other axes are unaffected. If *autoclear==True* and *enable==False*, also clear the current limit errors on all axes.

**limit\_errors\_enabled()**

Check if global limit errors are enabled.

If on, reaching limit switch on an axis puts it into an error state, which immediately stops this and all other axes; any further motion command on this axis will raise an error (it is still possible to restart motion on other axes); the axis motion can only be resumed by calling *clear\_limit\_error()*. If off, the limited axis still stops, but the other axes are unaffected.

**is\_enabled**(*axis='all'*)

Check if the axis output is enabled

**enable\_axis**(*axis='all'*, *enable=True*)

Enable axis output.

If the output is disabled, the steps are generated by the controller, but not sent to the motors.

**get\_position**(*axis='all'*)

Get the current axis pulse position

**set\_position\_reference**(*axis*, *position=0*)

Set the current axis pulse position as a reference.

Re-calibrate the pulse position counter so that the current position is set as *position* (0 by default).

**get\_encoder**(*axis='all'*)

Get the current axis encoder value

**set\_encoder\_reference**(*axis*, *position=0*)

Set the current axis encoder value as a reference.

Re-calibrate the encoder counter so that the current position is set as *position* (0 by default).

**move\_to**(*axis*, *position*)

Move a given axis to a given position

**move\_by** (*axis*, *steps=1*)  
Move a given axis for a given number of steps

**jog** (*axis*, *direction*)  
Jog a given axis in a given direction.  
*direction* can be either "-" (negative) or "+" (positive). The motion continues until it is explicitly stopped, or until a limit is hit.

**stop** (*axis='all'*, *immediate=False*)  
Stop motion of a given axis.  
If *immediate==True* make an abrupt stop; otherwise, slow down gradually.

**home** (*axis*, *direction*, *home\_mode*)  
Home the given axis using a given home mode.  
*direction* can be "+" or "-". The mode can be "only\_home\_input", "only\_home\_input\_lowspeed", "only\_limit\_input", "only\_zidx\_input", or "home\_and\_zidx\_input". For meaning, see Arcus PMX manual.

**get\_global\_speed** ()  
Get the global speed setting (in Hz); overridden by a non-zero axis speed

**get\_axis\_speed** (*axis='all'*)  
Get the individual axis speed setting (in Hz); 0 means that the global speed is used

**set\_global\_speed** (*speed*)  
Set the global speed setting (in Hz); overridden by a non-zero axis speed

**set\_axis\_speed** (*axis*, *speed*)  
Set the individual axis speed setting (in Hz); 0 means that the global speed is used

**get\_current\_axis\_speed** (*axis='all'*)  
Get the instantaneous speed (in Hz)

**get\_status\_n** (*axis='all'*)  
Get the axis status as an integer

**get\_status** (*axis='all'*)  
Get the axis status as a set of string descriptors

**is\_moving** (*axis='all'*)  
Check if a given axis is moving

**wait\_move** (*axis*, *timeout=None*, *period=0.05*)  
Wait until motion is done

**check\_limit\_error** (*axis='all'*)  
Check if the axis hit limit errors.  
Return "" (not errors), "+" (positive limit error) or "-" (negative limit error).

**clear\_limit\_error** (*axis='all'*)  
Clear axis limit errors

**get\_analog\_input** (*channel*)  
Get voltage (in V) at a given input (starting with 1)

**get\_digital\_input** (*channel*)  
Get value (0 or 1) at a given digital input (1 through 8)

**get\_digital\_input\_register** ()  
Get all 8 digital inputs as a single 8-bit integer

**get\_digital\_output** (*channel*)

Get value (0 or 1) at a given digital output (1 through 8)

**get\_digital\_output\_register** ()

Get all 8 digital inputs as a single 8-bit integer

**set\_digital\_output** (*channel*, *value*)

Set value (0 or 1) at a given digital output (1 through 8)

**set\_digital\_output\_register** (*value*)

Set all 8 digital inputs as a single 8-bit integer

**Error**

alias of `pylablib.devices.Arcus.base.ArcusError`

**class NoParameterCaller** (*device*, *kind*)

Bases: `object`

Class to simplify calling functions without a parameter

**apply\_settings** (*settings*)

Apply the settings.

*settings* is the dict {*name*: *value*} of the device available settings. Non-applicable settings are ignored.

**close** ()

Close the connection to the stage

**get\_all\_axes** ()

Get the list of all available axes (taking mapping into account)

**get\_device\_info** ()

Get the device info

**get\_device\_number** ()

Get the device number used in RS-485 communications.

Usually it is a string with the format similar to "4EX00".

**get\_device\_variable** (*key*)

Get the value of a settings, status, or full info parameter

**get\_full\_info** (*include=0*)

Get dict {*name*: *value*} containing full device information (including status and settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include*=-10 queries all available variables, which is equivalent to *include*="all".

**get\_full\_status** (*include=0*)

Get dict {*name*: *value*} containing the device status (including settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include*=-10 queries all available variables, which is equivalent to *include*="all".

**get\_settings** (*include=0*)

Get dict {*name*: *value*} containing all the device settings.

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the



lowest priority is -10, setting `include=-10` queries all available variables, which is equivalent to `include="all"`.

**is\_opened()**

Check if the device is connected

**open()**

Open the connection to the stage

**query(comm)**

Send a query to the stage and return the reply

**remap\_axes(mapping, accept\_original=True)**

Rename axes to the new labels.

*mapping* is the new axes mapping, which can be a list of new axes name (corresponding to the old axes in order returned by `get_all_axes()`), or a dictionary {alias: original} of the new axes aliases.

**set\_device\_number(number, store=True)**

Get the device number used in RS-485 communications.

*number* can be either a full device id (e.g., "4EX00"), or a single number between 0 and 99. In order for the change to take effect, the device needs to be power-cycled. If *store*==True, automatically store settings to the memory; otherwise, the settings will be lost unless `store_defaults()` is called at some point before the power-cycle.

**set\_device\_variable(key, value)**

Set the value of a settings parameter

**store\_defaults()**

Store some of the settings to the memory as defaults.

Applies to device number, baudrate, limit error behavior, polarity, and some other settings.

**class** `pylablib.devices.Arcus.performax.Performax2EXStage` (*idx=0, conn=None, enable=True*)

Bases: `pylablib.devices.Arcus.performax.Performax4EXStage`

Arcus Performax 2EX/2ED translation stage.

#### Parameters

- **idx** (*int*) – stage index; if using a USB connection, specifies a USB device index; if using RS485 connection, specifies device index on the bus
- **conn** – if not `None`, defines a connection to RS485 connection. Usually (e.g., for USB-to-RS485 adapters) this is a serial connection, which either a name (e.g., "COM1"), or a tuple (name, baudrate) (e.g., ("COM1", 9600)); if *conn* is `None`, assume direct USB connection and use the manufacturer-provided DLL
- **enable** – if `True`, enable all axes on startup

#### Error

alias of `pylablib.devices.Arcus.base.ArcusError`

**class** `NoParameterCaller` (*device, kind*)

Bases: `object`

Class to simplify calling functions without a parameter

**apply\_settings(settings)**

Apply the settings.

*settings* is the dict {name: value} of the device available settings. Non-applicable settings are ignored.

**check\_limit\_error** (*axis='all'*)

Check if the axis hit limit errors.

Return "" (not errors), "+" (positive limit error) or "-" (negative limit error).

**clear\_limit\_error** (*axis='all'*)

Clear axis limit errors

**close** ()

Close the connection to the stage

**enable\_absolute\_mode** (*enable=True*)

Set absolute motion mode

**enable\_axis** (*axis='all', enable=True*)

Enable axis output.

If the output is disabled, the steps are generated by the controller, but not sent to the motors.

**enable\_limit\_errors** (*enable=True, autoclear=True*)

Enable limit errors.

If on, reaching limit switch on an axis puts it into an error state, which immediately stops this and all other axes; any further motion command on this axis will raise an error (it is still possible to restart motion on other axes); the axis motion can only be resumed by calling *clear\_limit\_error()*. If off, the limited axis still stops, but the other axes are unaffected. If *autoclear==True* and *enable==False*, also clear the current limit errors on all axes.

**get\_all\_axes** ()

Get the list of all available axes (taking mapping into account)

**get\_analog\_input** (*channel*)

Get voltage (in V) at a given input (starting with 1)

**get\_axis\_speed** (*axis='all'*)

Get the individual axis speed setting (in Hz); 0 means that the global speed is used

**get\_baudrate** ()

Get current baud rate

**get\_current\_axis\_speed** (*axis='all'*)

Get the instantaneous speed (in Hz)

**get\_device\_info** ()

Get the device info

**get\_device\_number** ()

Get the device number used in RS-485 communications.

Usually it is a string with the format similar to "4EX00".

**get\_device\_variable** (*key*)

Get the value of a settings, status, or full info parameter

**get\_digital\_input** (*channel*)

Get value (0 or 1) at a given digital input (1 through 8)

**get\_digital\_input\_register** ()

Get all 8 digital inputs as a single 8-bit integer

**get\_digital\_output** (*channel*)  
Get value (0 or 1) at a given digital output (1 through 8)

**get\_digital\_output\_register** ()  
Get all 8 digital inputs as a single 8-bit integer

**get\_encoder** (*axis='all'*)  
Get the current axis encoder value

**get\_full\_info** (*include=0*)  
Get dict {name: value} containing full device information (including status and settings).  
*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_full\_status** (*include=0*)  
Get dict {name: value} containing the device status (including settings).  
*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_global\_speed** ()  
Get the global speed setting (in Hz); overridden by a non-zero axis speed

**get\_position** (*axis='all'*)  
Get the current axis pulse position

**get\_settings** (*include=0*)  
Get dict {name: value} containing all the device settings.  
*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_status** (*axis='all'*)  
Get the axis status as a set of string descriptors

**get\_status\_n** (*axis='all'*)  
Get the axis status as an integer

**home** (*axis, direction, home\_mode*)  
Home the given axis using a given home mode.  
*direction* can be "+" or "-". The mode can be "only\_home\_input", "only\_home\_input\_lowspeed", "only\_limit\_input", "only\_zidx\_input", or "home\_and\_zidx\_input". For meaning, see Arcus PMX manual.

**is\_enabled** (*axis='all'*)  
Check if the axis output is enabled

**is\_moving** (*axis='all'*)  
Check if a given axis is moving

**is\_opened** ()  
Check if the device is connected

**jog** (*axis, direction*)  
Jog a given axis in a given direction.

*direction* can be either "-" (negative) or "+" (positive). The motion continues until it is explicitly stopped, or until a limit is hit.

**limit\_errors\_enabled()**

Check if global limit errors are enabled.

If on, reaching limit switch on an axis puts it into an error state, which immediately stops this and all other axes; any further motion command on this axis will raise an error (it is still possible to restart motion on other axes); the axis motion can only be resumed by calling `clear_limit_error()`. If off, the limited axis still stops, but the other axes are unaffected.

**move\_by**(*axis, steps=1*)

Move a given axis for a given number of steps

**move\_to**(*axis, position*)

Move a given axis to a given position

**open**()

Open the connection to the stage

**query**(*comm*)

Send a query to the stage and return the reply

**remap\_axes**(*mapping, accept\_original=True*)

Rename axes to the new labels.

*mapping* is the new axes mapping, which can be a list of new axes name (corresponding to the old axes in order returned by `get_all_axes()`), or a dictionary {*alias*: *original*} of the new axes aliases.

**set\_axis\_speed**(*axis, speed*)

Set the individual axis speed setting (in Hz); 0 means that the global speed is used

**set\_baudrate**(*baudrate, store=True*)

Set current baud rate.

Acceptable values are 9600 (default), 19200, 38400, 57600, and 115200. In order for the change to take effect, the device needs to be power-cycled. If *store==True*, automatically store settings to the memory; otherwise, the settings will be lost unless `store_defaults()` is called at some point before the power-cycle.

**set\_device\_number**(*number, store=True*)

Get the device number used in RS-485 communications.

*number* can be either a full device id (e.g., "4EX00"), or a single number between 0 and 99. In order for the change to take effect, the device needs to be power-cycled. If *store==True*, automatically store settings to the memory; otherwise, the settings will be lost unless `store_defaults()` is called at some point before the power-cycle.

**set\_device\_variable**(*key, value*)

Set the value of a settings parameter

**set\_digital\_output**(*channel, value*)

Set value (0 or 1) at a given digital output (1 through 8)

**set\_digital\_output\_register**(*value*)

Set all 8 digital inputs as a single 8-bit integer

**set\_encoder\_reference**(*axis, position=0*)

Set the current axis encoder value as a reference.

Re-calibrate the encoder counter so that the current position is set as *position* (0 by default).

**set\_global\_speed** (*speed*)

Set the global speed setting (in Hz); overridden by a non-zero axis speed

**set\_position\_reference** (*axis*, *position=0*)

Set the current axis pulse position as a reference.

Re-calibrate the pulse position counter so that the current position is set as *position* (0 by default).

**stop** (*axis='all'*, *immediate=False*)

Stop motion of a given axis.

If *immediate==True* make an abrupt stop; otherwise, slow down gradually.

**store\_defaults** ()

Store some of the settings to the memory as defaults.

Applies to device number, baudrate, limit error behavior, polarity, and some other settings.

**wait\_move** (*axis*, *timeout=None*, *period=0.05*)

Wait until motion is done

```
class pylablib.devices.Arcus.performax.PerformaxDMXJSASTage (idx=0, conn=None,  
                                                         enable=True, auto-  
                                                         clear=True)
```

Bases: `pylablib.devices.Arcus.performax.GenericPerformaxStage`

Arcus Performax DMX-J-SA translation stage.

#### Parameters

- **idx** (*int*) – stage index; if using a USB connection, specifies a USB device index; if using RS485 connection, specifies device index on the bus
- **conn** – if not `None`, defines a connection to RS485 connection. Usually (e.g., for USB-to-RS485 adapters) this is a serial connection, which either a name (e.g., "COM1"), or a tuple (*name*, *baudrate*) (e.g., ("COM1", 9600)); if *conn* is `None`, assume direct USB connection and use the manufacturer-provided DLL
- **enable** – if `True`, enable all axes on startup
- **autoclear** – if `True`, automatically clear limit error before the motion start

**enable\_absolute\_mode** (*enable=True*)

Set absolute motion mode

**is\_enabled** ()

Check if the output is enabled

**enable\_axis** (*enable=True*)

Enable output.

If the output is disabled, the steps are generated by the controller, but not sent to the motors.

**get\_position** ()

Get the current pulse position

**set\_position\_reference** (*position=0*)

Set the current pulse position as a reference.

Re-calibrate the pulse position counter so that the current position is set as *position* (0 by default).

**move\_to** (*position*)

Move to a given position

**move\_by** (*steps=1*)

Move for a given number of steps

**jog** (*direction*)

Jog in a given direction.

*direction* can be either "-" (negative) or "+" (positive). The motion continues until it is explicitly stopped, or until a limit is hit.

**stop** (*immediate=False*)

Stop motion.

If *immediate==True* make an abrupt stop; otherwise, slow down gradually.

**home** (*direction, home\_mode*)

Home using a given home mode.

*direction* can be "+" or "-". The mode can be "only\_home\_input", "only\_home\_input\_lowspeed", or "only\_limit\_input". For meaning, see Arcus PMX manual.

**get\_axis\_speed** ()

Get the speed setting (in Hz)

**set\_axis\_speed** (*speed*)

Set the speed setting (in Hz)

**get\_status\_n** ()

Get the status as an integer

**get\_status** ()

Get the status as a set of string descriptors

**is\_moving** ()

Check if motor is moving

**wait\_move** (*timeout=None, period=0.05*)

Wait until motion is done

**check\_limit\_error** ()

Check if the motor hit limit errors.

Return "" (not errors), "+" (positive limit error) or "-" (negative limit error).

**clear\_limit\_error** ()

Clear limit error

**get\_digital\_input** (*channel*)

Get value (0 or 1) at a given digital input (1 through 5)

**get\_digital\_input\_register** ()

Get all 5 digital inputs as a single 5-bit integer

**get\_digital\_output** (*channel*)

Get value (0 or 1) at a given digital output (1 through 2)

**get\_digital\_output\_register** ()

Get all 2 digital outputs as a single 2-bit integer

**Error**

alias of `pylablib.devices.Arcus.base.ArcusError`

**class NoParameterCaller** (*device, kind*)

Bases: `object`

Class to simplify calling functions without a parameter

**apply\_settings** (*settings*)

Apply the settings.

*settings* is the dict {name: value} of the device available settings. Non-applicable settings are ignored.

**close** ()

Close the connection to the stage

**get\_all\_axes** ()

Get the list of all available axes (taking mapping into account)

**get\_device\_info** ()

Get the device info

**get\_device\_number** ()

Get the device number used in RS-485 communications.

Usually it is a string with the format similar to "4EX00".

**get\_device\_variable** (*key*)

Get the value of a settings, status, or full info parameter

**get\_full\_info** (*include=0*)

Get dict {name: value} containing full device information (including status and settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_full\_status** (*include=0*)

Get dict {name: value} containing the device status (including settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_settings** (*include=0*)

Get dict {name: value} containing all the device settings.

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**is\_opened** ()

Check if the device is connected

**open** ()

Open the connection to the stage

**query** (*comm*)

Send a query to the stage and return the reply

**remap\_axes** (*mapping, accept\_original=True*)

Rename axes to the new labels.

*mapping* is the new axes mapping, which can be a list of new axes name (corresponding to the old axes in order returned by `get_all_axes()`), or a dictionary {alias: original} of the new axes aliases.

**set\_device\_number** (*number*, *store=True*)

Get the device number used in RS-485 communications.

*number* can be either a full device id (e.g., "4EX00"), or a single number between 0 and 99. In order for the change to take effect, the device needs to be power-cycled. If *store==True*, automatically store settings to the memory; otherwise, the settings will be lost unless *store\_defaults()* is called at some point before the power-cycle.

**set\_device\_variable** (*key*, *value*)

Set the value of a settings parameter

**set\_digital\_output** (*channel*, *value*)

Set value (0 or 1) at a given digital output (1 through 2)

**store\_defaults** ()

Store some of the settings to the memory as defaults.

Applies to device number, baudrate, limit error behavior, polarity, and some other settings.

**set\_digital\_output\_register** (*value*)

Set all 2 digital inputs as a single 2-bit integer

## Module contents

### pylablib.devices.Arduino package

#### Submodules

#### pylablib.devices.Arduino.base module

**exception** pylablib.devices.Arduino.base.**ArduinoError**

Bases: *pylablib.core.devio.base.DeviceError*

Generic Arduino devices error

**args**

**with\_traceback** ()

Exception.with\_traceback(tb) – set self.\_\_traceback\_\_ to tb and return self.

**exception** pylablib.devices.Arduino.base.**ArduinoBackendError** (*exc*)

Bases: *pylablib.devices.Arduino.base.ArduinoError*, *pylablib.core.devio.comm\_backend.DeviceBackendError*

Generic Arduino backend communication error

**args**

**with\_traceback** ()

Exception.with\_traceback(tb) – set self.\_\_traceback\_\_ to tb and return self.

**class** pylablib.devices.Arduino.base.**IArduinoDevice** (*port*, *rate=9600*, *timeout=10.0*,  
*term\_write='n'*, *term\_read='n'*,  
*flush\_before\_op=True*)

Bases: *pylablib.core.devio.comm\_backend.ICommBackendWrapper*

Generic Arduino device.

#### Parameters

- **port** – serial port name



- **rate** – baud rate
- **timeout** – default communication timeout
- **term\_write** – default write terminating character (automatically appended on every sent message)
- **term\_read** – default read terminating character (used to determine when the incoming message is received completely)
- **flush\_before\_op** – if `True` (default), automatically flush input buffer on comm/query

**Error**

alias of `ArduinoError`

**reopen()**

Close and reopen the device connection

**reset\_board()**

Reset the board by pulsing the DTR and RTS lines

**comm(comm, timeout=None, flush=False, flush\_delay=0.02)**

Send a device command.

If *timeout* is not `None`, it specifies a custom timeout for the operation. If *flush*==`True`, then wait for *flush\_delay* seconds after the write and read everything returned by the device.

**query(query, timeout=None, query\_delay=0, flush=False, flush\_delay=0.02)**

Send a device query and return the reply.

If *timeout* is not `None`, it specifies a custom timeout for the reply read operation. If *query\_delay*>0, it specifies the delay between write and subsequent read attempt. If *flush*==`True`, then wait for *flush\_delay* seconds after the write and read everything returned by the device.

**class NoParameterCaller(device, kind)**

Bases: `object`

Class to simplify calling functions without a parameter

**apply\_settings(settings)**

Apply the settings.

*settings* is the dict {name: value} of the device available settings. Non-applicable settings are ignored.

**close()**

Close the backend

**get\_device\_variable(key)**

Get the value of a settings, status, or full info parameter

**get\_full\_info(include=0)**

Get dict {name: value} containing full device information (including status and settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include*=-10 queries all available variables, which is equivalent to *include*="all".

**get\_full\_status(include=0)**

Get dict {name: value} containing the device status (including settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting `include=-10` queries all available variables, which is equivalent to `include="all"`.

**get\_settings** (*include=0*)

Get dict {name: value} containing all the device settings.

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting `include=-10` queries all available variables, which is equivalent to `include="all"`.

**is\_opened** ()

Check if the device is connected

**lock** (*timeout=None*)

Lock the access to the device from other threads/processes (isn't necessarily implemented)

**locking** (*timeout=None*)

Context manager for lock & unlock

**open** ()

Open the backend

**set\_device\_variable** (*key, value*)

Set the value of a settings parameter

**unlock** ()

Unlock the access to the device from other threads/processes (isn't necessarily implemented)

## Module contents

### pylablib.devices.Attocube package

#### Submodules

#### pylablib.devices.Attocube.anc300 module

`pylablib.devices.Attocube.anc300.muxaxis` (\*args, \*\*kwargs)

**class** `pylablib.devices.Attocube.anc300.TDeviceInfo` (*serial, version*)

Bases: `tuple`

**count** ()

Return number of occurrences of value.

**index** ()

Return first index of value.

Raises `ValueError` if the value is not present.

**serial**

**version**

**class** `pylablib.devices.Attocube.anc300.ANC300` (*conn, backend='auto', pwd='123456'*)

Bases: `pylablib.core.devio.comm_backend.ICommBackendWrapper`, `pylablib.devices.interface.stage.IMultiaxisStage`

Attocube ANC300 controller.

### Parameters

- **conn** – connection parameters; for Ethernet connection is a tuple (*addr*, *port*), a string "*addr:port*", or a string "*addr*" (default port 7240 us assumed)
- **backend** (*str*) – communication backend; by default, try to determine from the communication parameters
- **pwd** (*str*) – connection password for Ethernet connection (default is "123456")

### Error

alias of `pylablib.devices.Attocube.base.AttocubeError`

### **open** ()

Open the connection to the stage

### **query** (*msg*)

Send a query to the stage and return the reply

### **update\_available\_axes** ()

Update the list of available axes.

Need to call only if the hardware configuration of the ANC module has changed.

### **get\_device\_info** ()

Get the device info of the controller board: (*serial*, *version*)

### **get\_axis\_serial** (*axis*='all')

Get serial number of the controller board

### **set\_mode** (*axis*='all', *mode*='stp')

Set axis mode.

*axis* is either an axis index (starting from 1), or "all" (all axes). *mode* can be "gnd" (ground), "stp" (step), "cap" (measure capacitance, then ground), "offs" (offset only, no stepping), "stp+" (offset with added stepping waveform), "stp-" (offset with subtracted stepping). Note that not all modes are supported by all modules: ANM150 doesn't support offset voltage ("offs", "stp+", "stp-" modes), ANM200 doesn't support stepping ("stp", "stp+", "stp-" modes).

### **get\_mode** (*axis*='all')

Get axis mode.

*axis* is either an axis index (starting from 1), or "all" (all axes). See `set_mode()` for the description of the modes.

### **is\_enabled** (*axis*='all')

Check if the axis is enabled

### **enable\_axis** (*axis*='all', *mode*='stp')

Enable specific axis (set to step mode)

### **disable\_axis** (*axis*='all')

Disable specific axis (set to ground mode)

### **measure\_capacitance** (*axis*='all', *wait*=True)

Measure axis capacitance; finish in the GND mode.

If *wait*==True, wait until the capacitance measurement is finished (takes about a second per axis).

### **get\_voltage** (*axis*='all')

Get axis step amplitude in Volts

**set\_voltage** (*axis, voltage*)

Set axis step amplitude in Volts

**get\_offset** (*axis='all'*)

Get axis offset voltage in Volts

**set\_offset** (*axis, voltage*)

Set axis offset voltage in Volts

**get\_output** (*axis='all'*)

Get axis current output voltage in Volts

**get\_frequency** (*axis='all'*)

Get axis step frequency in Hz

**set\_frequency** (*axis, freq*)

Set axis step frequency in Hz

**get\_capacitance** (*axis='all', measure=False*)

Get capacitance measurement on the axis.

If `measure==True`, re-measure axis capacitance (takes about a second); otherwise, get the last measurement value.

**get\_voltage\_pattern** (*axis, kind*)

Get axis voltage pattern.

*kind* be either "up" for up pattern or "down" for down pattern. The pattern is a numpy array of 256 numbers from 0 to 255 corresponding to the output voltage from 0 to the axis voltage. This pattern is output (repeatedly) for each step. The default is a simple linear ramp.

**set\_voltage\_pattern** (*axis, kind, pattern=None*)

Set axis voltage pattern.

*kind* be either "up" for up pattern or "down" for down pattern. The pattern is an array of 256 numbers from 0 to 255 corresponding to the output voltage from 0 to the axis voltage. This pattern is output (repeatedly) for each step. The default is a simple linear ramp, which is set if `pattern` is `None`.

**get\_trigger\_input** (*axis='all'*)

Get trigger input lines for the given axis.

Return tuple (`up`, `down`) with values for up and down step triggers, which can be either integer with the trigger line number, or "off" if the trigger is off.

**set\_trigger\_input** (*axis, up=None, down=None*)

Set trigger input lines for the given axis.

*up* and *down* are can be integer with the trigger line number, "off" if the trigger is off, or `None` (keep the value unchanged).

**get\_external\_input\_modes** (*axis='all'*)

Get external BNC input modes.

Return tuple (`acin`, `dcin`) indicating whether AC-IN and DC-IN channels are enabled.

**set\_external\_input\_modes** (*axis, acin=None, dcin=None*)

Enable or disable external BNC inputs.

*acin* and *dcin* are can be boolean indicating if the corresponding input is enabled, or `None` (keep the value unchanged).

**get\_axis\_correction** (*axis*)

Get axis correction factor.

The factor is automatically applied when the motion is in the negative direction.

**set\_axis\_correction** (*axis*, *factor=1.0*)

Set axis correction factor.

The factor is automatically applied when the motion is in the negative direction.

**jog** (*axis*, *direction*)

Jog continuously in the given direction ("+" or "-").

The motion will continue until another move or stop command is called.

**move\_by** (*axis*, *steps=1*)

Move a given axis for a given number of steps

**wait\_move** (*axis*, *timeout=30.0*)

Wait for a given axis to stop moving.

If the motion is not finished after *timeout* seconds, raise a backend error.

**is\_moving** (*axis*)

Check if a given axis is moving

**stop** (*axis='all'*)

Stop motion of a given axis

**class NoParameterCaller** (*device*, *kind*)

Bases: `object`

Class to simplify calling functions without a parameter

**apply\_settings** (*settings*)

Apply the settings.

*settings* is the dict {name: value} of the device available settings. Non-applicable settings are ignored.

**close** ()

Close the backend

**get\_all\_axes** ()

Get the list of all available axes (taking mapping into account)

**get\_device\_variable** (*key*)

Get the value of a settings, status, or full info parameter

**get\_full\_info** (*include=0*)

Get dict {name: value} containing full device information (including status and settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_full\_status** (*include=0*)

Get dict {name: value} containing the device status (including settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_settings** (*include=0*)

Get dict {name: value} containing all the device settings.

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting `include=-10` queries all available variables, which is equivalent to `include="all"`.

**is\_opened()**

Check if the device is connected

**lock** (*timeout=None*)

Lock the access to the device from other threads/processes (isn't necessarily implemented)

**locking** (*timeout=None*)

Context manager for lock & unlock

**remap\_axes** (*mapping, accept\_original=True*)

Rename axes to the new labels.

*mapping* is the new axes mapping, which can be a list of new axes name (corresponding to the old axes in order returned by `get_all_axes()`), or a dictionary {*alias*: *original*} of the new axes aliases.

**set\_device\_variable** (*key, value*)

Set the value of a settings parameter

**unlock()**

Unlock the access to the device from other threads/processes (isn't necessarily implemented)

## pylablib.devices.Attocube.anc350 module

`pylablib.devices.Attocube.anc350.get_usb_devices_number()`

Get the number of controllers connected via USB

**class** `pylablib.devices.Attocube.anc350.ANC350` (*conn=0, timeout=5.0*)

Bases: `pylablib.core.devio.comm_backend.ICommBackendWrapper`, `pylablib.devices.interface.stage.IMultiaxisStage`

Attocube ANC350 controller.

### Parameters

- **conn** – connection parameters - index of the Attocube ANC350 in the system (for a single controller leave 0)
- **timeout** (*float*) – default operation timeout

### Error

alias of `pylablib.devices.Attocube.base.AttocubeError`

**class** `Telegram` (*opcode, address, index, data, corr\_number*)

Bases: `tuple`

**address**

**corr\_number**

**count()**

Return number of occurrences of value.

**data**

**index**

**opcode**

```

class Reply (address, index, reason, data)
    Bases: tuple

    address

    count ()
        Return number of occurrences of value.

    data

    index

    reason

check_tell (timeout=0.01)
    Check for queued TELL (periodic value update) commands

set_value (address, index, value, ack=False)
    Set device value at the given address and index.

    If ack==True, request ACK responds and return its value; otherwise, return immediately after set.

get_value (address, index, as_int=True)
    Get device value at the given address and index.

    If as_int==True, convert the result into a signed integer; otherwise return raw byte string.

enable_updates (enabled=True)
    Enable or disable periodic TELL updates

get_hardware_id ()
    Return device HWID (by default -1)

set_hardware_id (hwid, persist=False)
    Set device HWID (can be used to identify different devices).

    If persist==True, the value persists after power cycling.

is_connected (axis='all')
    Check if axis is connected

is_enabled (axis='all')
    Check if axis is enabled

enable_axis (axis='all', enabled=True)
    Enable a specific axis or all axes

disable_axis (axis='all')
    Disable a specific axis or all axes

is_moving (axis='all')
    Move a given axis for a given number of steps

check_limit (axis='all')
    Check if the end of travel has been reached.

    Return None if no limits are reached, "fwd" if forward limit is reached, "bwd" if backward limit is
    reached, or "both" if both are reached together (normally shouldn't happen).

get_status_n (axis='all')
    Get numerical status of the axis.

    For details, see ANC350 protocol.

status_bits = [(1, 'running'), (2, 'limit'), (256, 'sens_err'), (1024, 'sens_disconn')]

```

**get\_status** (*axis='all'*)

Get device status.

Return list of status strings, which can include "running" (axis is moving), "limit" (one of the limits is reached), "sens\_err" (sensor error), "sens\_disconn" (sensor disconnected), or "ref\_valid" (reference is valid).

**get\_target\_position** (*axis='all'*)

Get the target position for the given axis (the position towards which it is moving)

**get\_precision** (*axis='all'*)

Get the axis precision in m (used for checking if the target is reached)

**set\_precision** (*axis='all', precision=1e-06*)

Set the axis precision in m (used for checking if the target is reached)

**is\_target\_reached** (*axis='all', precision=None*)

Check if the target position is reached.

If *precision* is not None, it sets final position tolerance (in m).

**get\_sensor\_voltage** ()

Get position sensor voltage in Volts

**set\_sensor\_voltage** (*voltage*)

Set position sensor voltage in Volts

**get\_voltage** (*axis='all'*)

Get axis step voltage in Volts

**set\_voltage** (*axis, voltage*)

Set axis step voltage in Volts

**get\_offset** (*axis='all'*)

Get axis offset voltage in Volts

**set\_offset** (*axis, voltage*)

Set axis offset voltage in Volts

**get\_frequency** (*axis='all'*)

Get axis step frequency in Hz

**set\_frequency** (*axis, freq*)

Set axis step frequency in Hz

**get\_capacitance** (*axis='all', measure=False, delay=0.5*)

Get axis capacitance in F.

If *measure==True*, initialize the measurement and get the result after the measurement *delay*. Otherwise, return the last measured value.

**get\_position** (*axis='all'*)

Get axis position (in m)

**move\_to** (*axis, position, precision=None*)

Move to target position (in m).

If *precision* is not None, it sets final position tolerance.

**move\_by** (*axis, dist*)

Move along a given axis by a given distance (in m)

**move\_by\_steps** (*axis, steps=1*)

Move along a given axis by a given number of steps



**wait\_move** (*axis*, *precision=1e-06*, *timeout=10.0*, *period=0.01*)

Wait for a given axis to stop moving or to reach target position.

If the motion is not finished after *timeout* seconds, raise a backend error. Precision sets the final positioning precision (in m).

**stop** (*axis='all'*)

Stop motion of a given axis

**jog** (*axis*, *direction*)

Jog a given axis in a given direction.

*direction* can be either "-" (negative) or "+" (positive). The motion continues until it is explicitly stopped, or until a limit is hit.

**class NoParameterCaller** (*device*, *kind*)

Bases: `object`

Class to simplify calling functions without a parameter

**apply\_settings** (*settings*)

Apply the settings.

*settings* is the dict {*name*: *value*} of the device available settings. Non-applicable settings are ignored.

**close** ()

Close the backend

**get\_all\_axes** ()

Get the list of all available axes (taking mapping into account)

**get\_device\_variable** (*key*)

Get the value of a settings, status, or full info parameter

**get\_full\_info** (*include=0*)

Get dict {*name*: *value*} containing full device information (including status and settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_full\_status** (*include=0*)

Get dict {*name*: *value*} containing the device status (including settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_settings** (*include=0*)

Get dict {*name*: *value*} containing all the device settings.

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**is\_opened** ()

Check if the device is connected

**lock** (*timeout=None*)

Lock the access to the device from other threads/processes (isn't necessarily implemented)

**locking** (*timeout=None*)

Context manager for lock & unlock

**open** ()

Open the backend

**remap\_axes** (*mapping, accept\_original=True*)

Rename axes to the new labels.

*mapping* is the new axes mapping, which can be a list of new axes name (corresponding to the old axes in order returned by `get_all_axes()`), or a dictionary {*alias*: *original*} of the new axes aliases.

**set\_device\_variable** (*key, value*)

Set the value of a settings parameter

**unlock** ()

Unlock the access to the device from other threads/processes (isn't necessarily implemented)

## pylablib.devices.Attocube.base module

**exception** `pylablib.devices.Attocube.base.AttocubeError`

Bases: `pylablib.core.devio.base.DeviceError`

Generic Attocube error

**args**

**with\_traceback** ()

Exception.with\_traceback(tb) – set self.\_\_traceback\_\_ to tb and return self.

**exception** `pylablib.devices.Attocube.base.AttocubeBackendError` (*exc*)

Bases: `pylablib.devices.Attocube.base.AttocubeError`, `pylablib.core.devio.comm_backend.DeviceBackendError`

Attocube backend communication error

**args**

**with\_traceback** ()

Exception.with\_traceback(tb) – set self.\_\_traceback\_\_ to tb and return self.

## Module contents

### pylablib.devices.Basler package

#### Submodules

### pylablib.devices.Basler.pylon module

**class** `pylablib.devices.Basler.pylon.LibraryController` (*lib*)

Bases: `pylablib.devices.utils.load_lib.LibraryController`

**close** (*opid*)

Mark device closing.

Return tuple (*close\_result*, *uninit\_result*) with the results of the closing and the shutdown. If library does not need to be shut down yet, set *uninit\_result*=None

**open** ()

Mark device opening.

Return tuple (*init\_result*, *open\_result*, *opid*) with the results of the initialization and the opening, and the opening ID which should afterwards be used for closing. If library is already initialized, set *init\_result*=None

**preinit** ()

Pre-initialize the library, if it hasn't been done already

**shutdown** ()

Close all opened connections and shutdown the library

**temp\_open** ()

Context for temporarily opening a new device connection

`pylablib.devices.Basler.pylon.restart_lib()`

**class** `pylablib.devices.Basler.pylon.TCameraInfo` (*name*, *model*, *serial*, *devclass*,  
*devversion*, *vendor*, *friendly\_name*,  
*user\_name*, *props*)

Bases: `tuple`

**count** ()

Return number of occurrences of value.

**devclass**

**devversion**

**friendly\_name**

**index** ()

Return first index of value.

Raises ValueError if the value is not present.

**model**

**name**

**props**

**serial**

**user\_name**

**vendor**

`pylablib.devices.Basler.pylon.get_device_info` (*index*)

Get Pylon camera info for a camera with the given index

`pylablib.devices.Basler.pylon.list_cameras` (*desc*=True)

List all cameras available through Basler Pylon interface

If *desc*==True, return complete camera descriptions; otherwise, simply return the names.

`pylablib.devices.Basler.pylon.get_cameras_number` ()

Get number of connected dx cameras

```
class pylablib.devices.Basler.pylon.BaslerPylonAttribute (node, full_name=None)
```

Bases: `object`

Object representing an Pylon GenAPI attribute.

Allows to query and set values and get additional information. Usually created automatically by an *BaslerPylonCamera* instance.

#### Parameters

- **node** – pylon GenApi node handler
- **full\_name** – if supplied, attribute's full name, including the tree structure

**name**  
attribute name

**kind**  
attribute kind; can be "int", "float", "bool", "enum", "str", "command", "category", or "unknown"

**display\_name**  
attribute display name (short description name)

**tooltip**  
longer attribute description

**description**  
full attribute description (usually, same as *tooltip*)

**visibility**  
attribute visibility; can be "simple", "intermediate", "advanced", "invisible", or "unknown"

**access**  
attribute access level; can be "read\_only", "write\_only", "read\_write", "na" (not applicable, e.g., command), or "not\_implemented"

**readable**  
whether attribute is readable

Type `bool`

**writable**  
whether attribute is writable

Type `bool`

**implemented**  
whether the attribute is implemented in the given camera (normally always `True`)

Type `bool`

**available**  
whether the attribute can be changed or called

Type `bool`

**min**  
minimal attribute value (if applicable)

Type `float` or `int`

**max**  
maximal attribute value (if applicable)

Type float or int

**inc**

minimal attribute increment value (if applicable)

Type float or int

**units**

attribute units (if applicable)

**repr**

shows what a numerical unit represents; can be "lin", "log", "bool", "pure", "hex", or "unknown"

**ivalues**

list of possible integer values for enum attributes

**values**

list of possible text values for enum attributes

**labels**

dict {label: index} which shows all possible values of an enumerated attribute and their corresponding numerical values

**ilabels**

dict {index: label} which shows labels corresponding to numerical values of an enumerated attribute

**update\_limits()**

Update minimal and maximal attribute limits and return tuple (min, max, inc)

**truncate\_value(value)**

Truncate value to lie within attribute limits

**get\_value(enum\_as\_str=True)**

Get attribute value.

If enum\_as\_str==True, return enum-style values as strings; otherwise, return corresponding integer values.

**set\_value(value, truncate=True)**

Set attribute value.

If truncate==True, automatically truncate value to lie within allowed range.

**call\_command()**

Execute the given command

```
class pylablib.devices.Basler.pylon.TDeviceInfo(name, model, serial, devclass,
                                              devversion, vendor, friendly_name,
                                              user_name, props)
```

Bases: tuple

**count()**

Return number of occurrences of value.

**devclass**

**devversion**

**friendly\_name**

**index()**

Return first index of value.

Raises ValueError if the value is not present.

**model**

**name**

**props**

**serial**

**user\_name**

**vendor**

```
class pylablib.devices.Basler.pylon.BaslerPylonCamera (idx=0, name=None)
    Bases:      pylablib.devices.interface.camera.IROICamera,      pylablib.devices.
                interface.camera.IAttributeCamera,                  pylablib.devices.interface.camera.
                IExposureCamera
```

Generic Basler pylon camera interface.

#### Parameters

- **idx** – camera index among the cameras listed using `list_cameras()`
- **name** – camera name; if specified, then *idx* is ignored and the camera is determined based on the provided name

```
Error = <Mock name='mock.BaslerError' id='140096294766992'>
```

```
TimeoutError = <Mock spec='str' id='140096274462032'>
```

```
open ()
```

Open connection to the camera

```
close ()
```

Close connection to the camera

```
is_opened ()
```

Check if the device is connected

```
post_open ()
```

Additional setup after camera opening

```
get_attribute_value (name, error_on_missing=True, default=None, enum_as_str=True)
```

Get value of an attribute with the given name.

If the value doesn't exist or can not be read and `error_on_missing==True`, raise error; otherwise, return *default*. If *default* is not None, assume that `error_on_missing==False`. If *name* points at a dictionary branch, return a dictionary with all values in this branch. If `enum_as_str==True`, return enum-style values as strings; otherwise, return corresponding integer values.

```
set_attribute_value (name, value, truncate=True, error_on_missing=True)
```

Set value of an attribute with the given name.

If the value doesn't exist or can not be written and `error_on_missing==True`, raise error; otherwise, do nothing. If *name* points at a dictionary branch, set all values in this branch (in this case *value* must be a dictionary). If `truncate==True`, truncate value to lie within attribute range.

```
call_command (name)
```

Execute the given command

```
get_all_attribute_values (root="", enum_as_str=True, ignore_errors=True)
```

Get values of all attributes with the given *root*

**set\_all\_attribute\_values** (*settings*, *root=""*, *truncate=True*)

Set values of all attributes with the given *root*.

If *truncate==True*, truncate value to lie within attribute range.

**get\_device\_info** ()

Get camera information.

Return tuple (name, model, serial, devclass, devversion, vendor, friendly\_name, user\_name, props).

**get\_detector\_size** ()

Get camera detector size (in pixels) as a tuple (width, height)

**get\_roi** ()

Get current ROI.

Return tuple (*hstart*, *hend*, *vstart*, *vend*). *hstart* and *hend* specify horizontal image extent, *vstart* and *vend* specify vertical image extent (start is inclusive, stop is exclusive, starting from 0).

**set\_roi** (*hstart=0*, *hend=None*, *vstart=0*, *vend=None*)

Setup camera ROI.

*hstart* and *hend* specify horizontal image extent, *vstart* and *vend* specify vertical image extent (start is inclusive, stop is exclusive, starting from 0). By default, all non-supplied parameters take extreme values (0 for start, maximal for end).

**get\_roi\_limits** (*hbin=1*, *vbin=1*)

Get the minimal and maximal ROI parameters.

Return tuple (*hlim*, *vlim*), where each element is in turn a limit 5-tuple (*min*, *max*, *pstep*, *sstep*, *maxbin*) with, correspondingly, minimal and maximal size, position and size step, and the maximal binning (fixed to 1 if not binning is allowed). In some cameras, the step and the minimal size depend on the binning, which can be supplied.

**get\_exposure** ()

Get current exposure

**set\_exposure** (*exposure*)

Set camera exposure

**get\_frame\_period** ()

Get frame period (time between two consecutive frames in the internal trigger mode)

**set\_frame\_period** (*frame\_period*)

Set frame period (time between two consecutive frames in the internal trigger mode)

**get\_frame\_timings** ()

Get acquisition timing.

Return tuple (*exposure*, *frame\_period*).

**class BufferManager** (*strm*, *size*, *nbuff*)

Bases: `object`

Buffer manager, which deals with buffer memory allocation, registering and deregistering, and retrieving the result and the leftovers

**register** ()

Register buffers

**deregister** ()

Deregister buffers

**get\_buffer** (*fidx*)  
Get buffer corresponding to the given frame index

**get\_handle** (*fidx*)  
Get buffer handle corresponding to the given frame index

**get\_all\_handles** ()  
Get all buffer handles as a numpy array

**queue** (*fidx=None*)  
Queue a buffer with the given index or all buffers

**retrieve** ()  
Retrieve the next buffer and return its info and whether it is ready

**flush** ()  
Retrieve all leftover buffers

**class ScheduleLooper**  
Bases: `object`  
Schedule loop handler, which starts and stop the separate schedule loop thread

**start\_loop** (*buff\_mgr*)  
Start loop serving the given buffer manager

**stop\_loop** ()  
Stop the loop thread

**is\_looping** ()  
Check if the loop is running

**get\_status** ()  
Get the current loop status, which is the tuple (*acquired*,)

**setup\_acquisition** (*mode='sequence', nframes=100*)  
Setup acquisition mode.  
  
*mode* can be either "snap" (single frame or a fixed number of frames) or "sequence" (continuous acquisition). *nframes* sets up number of frame buffers.

**clear\_acquisition** ()  
Clear acquisition settings

**start\_acquisition** (\**args*, \*\**kwargs*)  
Start acquisition.  
  
Can take the same keyword parameters as :meth:“*setup\_acquisition*. If the acquisition is not set up yet, set it up using the supplied parameters (use default of `setup_acquisition()`, if the parameter is `None`). Otherwise, if any supplied parameters are different from the current ones, change them and reset the acquisition.

**stop\_acquisition** ()  
Stop acquisition

**acquisition\_in\_progress** ()  
Check if acquisition is in progress

**enable\_raw\_readout** (*enable='rows'*)  
Enable raw frame transfer.  
  
Should be used if the camera uses unsupported pixel format. Can be "frame" (return the whole frame as a 1D "u1" numpy array), "rows" (return a 2D array, where each row corresponds to a single image row), or `False` (convert to image data, or raise an error if the format is not supported; default)



**FrameTransferError**

alias of `pylablib.devices.interface.camera.DefaultFrameTransferError`

**class NoParameterCaller** (*device, kind*)

Bases: `object`

Class to simplify calling functions without a parameter

**apply\_settings** (*settings*)

Apply the settings.

*settings* is the dict {name: value} of the device available settings. Non-applicable settings are ignored.

**get\_acquisition\_parameters** ()

Get acquisition parameters.

Return dictionary {name: value}

**get\_all\_attributes** (*copy=False*)

Return a dictionary of all available attributes.

If *copy*==True, copy the dictionary; otherwise, return the internal dictionary structure (should not be modified).

**get\_attribute** (*name, error\_on\_missing=True*)

Get the camera attribute with the given name

**get\_data\_dimensions** ()

Get readout data dimensions (in pixels) as a tuple (width, height); take indexing mode into account

**get\_device\_variable** (*key*)

Get the value of a settings, status, or full info parameter

**get\_frame\_format** ()

Get format for the returned images.

Can be "list" (list of 2D arrays), "array" (a single 3D array), or "chunks" (list of 3D "chunk" arrays; supported for some cameras and provides the best performance).

**get\_frame\_info\_fields** ()

Get the names of frame info fields.

Applicable when frame info format (set by `set_frame_info_format()`) is "list" or "array".

**get\_frame\_info\_format** ()

Get format of the frame info.

Can be "namedtuple" (potentially nested named tuples; convenient to get particular values), "list" (flat list of values, with field names are given by `get_frame_info_fields()`; convenient for building a table), "array" (same as "list", but with a numpy array, which is easier to use for "chunks" frame format), or "dict" (flat dictionary with the same fields as the "list" format; more resilient to future format changes)

**get\_frame\_info\_period** ()

Get period of frame info acquisition.

Frame info might be skipped (set to None) except for frames which indices are divisible by *period*. Useful for certain cameras where acquiring frame info takes a lot of time and can reduce performance at higher frame rates. Note that this parameter can still be ignored (i.e., always set to 1) if the performance is not an issue for a given camera class.

**get\_frames\_status** ()

Get acquisition and buffer status.

Return tuple (acquired, unread, skipped, size), where acquired is the total number of acquired frames, unread is the number of acquired but not read frames, skipped is the number of skipped (not read and then written over) frames, and buffer\_size is the total buffer size (in frames).

**get\_full\_info** (include=0)

Get dict {name: value} containing full device information (including status and settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include*=-10 queries all available variables, which is equivalent to *include*="all".

**get\_full\_status** (include=0)

Get dict {name: value} containing the device status (including settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include*=-10 queries all available variables, which is equivalent to *include*="all".

**get\_image\_indexing** ()

Get indexing for the returned images.

Can be "rct" (first index row, second index column, rows counted from the top), "rcb" (same as "rc", rows counted from the bottom), "xyt" (first index column, second index row, rows counted from the top), or "xyb" (same as "xyt", rows counted from the bottom)

**get\_new\_images\_range** ()

Get the range of the new images.

Return tuple (first, last) with images range (first inclusive). If no images are available, return None. If some images were in the buffer were overwritten, exclude them from the range.

**get\_settings** (include=0)

Get dict {name: value} containing all the device settings.

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include*=-10 queries all available variables, which is equivalent to *include*="all".

**grab** (nframes=1, frame\_timeout=5.0, missing\_frame='skip', return\_info=False, buff\_size=None)

Snap *nframes* images (with preset image read mode parameters)

*buff\_size* determines buffer size (if None, use the default size). Timeout is specified for a single-frame acquisition, not for the whole acquisition time. *missing\_frame* determines what to do with frames which have been lost: can be "none" (replacing them with None), "zero" (replacing them with zero-filled frame), or "skip" (skipping them, while still keeping total returned frames number to *n*). If *return\_info*==True, return tuple (frames, infos), where *infos* is a list of frame info tuples (camera-dependent); if some frames are missing and *missing\_frame*!="skip", the corresponding frame info is None.

**is\_acquisition\_setup** ()

Check if acquisition is set up.

If the camera does not support separate acquisition setup, always return True.

**pausing\_acquisition** (clear=None, stop=True, setup\_after=None, start\_after=True, combine\_nested=True)

Context manager which temporarily pauses acquisition during execution of with block.

Useful for applying certain settings which can't be changed during the acquisition. If `clear==True`, clear acquisition in addition to stopping (by default, use the class default specified as `_clear_pausing_acquisition` attribute). If `stop==True`, stop the acquisition (if `clear==True`, stop regardless). If `setup_after==True`, setup the acquisition after pause if necessary (None means setup only if clearing was required). If `start_after==True`, start the acquisition after pause if necessary (None means start only if stopping was required). If `combine_nested==True`, then any nested `pausing_acquisition` calls will stop/clear acquisition as necessary, but won't setup/start it again until this `pausing_acquisition` call is complete.

Yields tuple (`acq_in_progress`, `acq_params`), which indicates whether acquisition is currently in progress, and what are the current acquisition parameters.

**read\_multiple\_images** (*rng=None, peek=False, missing\_frame='skip', return\_info=False, return\_rng=False*)

Read multiple images specified by *rng* (by default, all un-read images).

If *rng* is specified, it is a tuple (*first*, *last*) with images range (*first* inclusive). If no new frames are available, return an empty list; if no acquisition is running, return None. If `peek==True`, return images but not mark them as read. *missing\_frame* determines what to do with frames which are out of range (missing or lost): can be "none" (replacing them with None), "zero" (replacing them with zero-filled frame), or "skip" (skipping them). If `return_info==True`, return tuple (*frames*, *infos*), where *infos* is a list of frame info tuples (camera-dependent, by default, only the frame index); if some frames are missing and `missing_frame!="skip"`, the corresponding frame info is None. If `return_rng==True`, return the range covered resulting frames; if `missing_frame=="skip"`, the range can be smaller than the supplied *rng* if some frames are skipped.

**read\_newest\_image** (*peek=False, return\_info=False*)

Read the newest un-read image.

If no un-read frames are available, return None. If `peek==True`, return the image but not mark it as read. If `return_info==True`, return tuple (*frame*, *info*), where *info* is an info tuples (camera-dependent, see `read_multiple_images()`).

**read\_oldest\_image** (*peek=False, return\_info=False*)

Read the oldest un-read image.

If no un-read frames are available, return None. If `peek==True`, return the image but not mark it as read. If `return_info==True`, return tuple (*frame*, *info*), where *info* is an info tuples (camera-dependent, see `read_multiple_images()`).

**set\_device\_variable** (*key, value*)

Set the value of a settings parameter

**set\_frame\_format** (*fnt*)

Set format for the returned images.

Can be "list" (list of 2D arrays), "array" (a single 3D array), "chunks" (list of 3D "chunk" arrays; supported for some cameras and provides the best performance), or "try\_chunks" (same as "chunks", but if chunks are not supported, set to "list" instead). If format is "chunks" and chunks are not supported by the camera, it results in one frame per chunk. Note that if the format is set to "array" or "chunks", the frame info format is also automatically set to "array". If the format is set to "chunks", then the image info is also returned in chunks form (list of 2D info arrays with the same length as the corresponding frame chunks).

**set\_frame\_info\_format** (*fnt, include\_fields=None*)

Set format of the frame info.

Can be "namedtuple" (potentially nested named tuples; convenient to get particular values), "list" (flat list of values, with field names are given by `get_frame_info_fields()`; convenient for building a table), "array" (same as "list", but with a numpy array, which is easier to use for "chunks"

frame format), or "dict" (flat dictionary with the same fields as the "list" format; more resilient to future format changes) If *include\_fields* is not None, it specifies the fields included for non-"tuple" formats; note that order of *include\_fields* is ignored, and the resulting fields are always ordered same as in the original.

**set\_frame\_info\_period** (*period=1*)

Set period of frame info acquisition.

Frame info might be skipped (set to None) except for frames which indices are divisible by *period*. Useful for certain cameras where acquiring frame info takes a lot of time and can reduce performance at higher frame rates. Note that this parameter can still be ignored (i.e., always set to 1) if the performance is not an issue for a given camera class.

**set\_image\_indexing** (*indexing*)

Set up indexing for the returned images.

Can be "rct" (first index row, second index column, rows counted from the top), "rcb" (same as "rc", rows counted from the bottom), "xyt" (first index column, second index row, rows counted from the top), or "xyb" (same as "xyt", rows counted from the bottom)

**snap** (*timeout=5.0, return\_info=False*)

Snap a single frame

**wait\_for\_frame** (*since='lastread', nframes=1, timeout=20.0, error\_on\_stopped=False*)

Wait for one or several new camera frames.

*since* specifies the reference point for waiting to acquire *nframes* frames; can be "lastread" (from the last read frame), "lastwait" (wait for the last successful `wait_for_frame()` call), "now" (from the start of the current call), or "start" (from the acquisition start, i.e., wait until *nframes* frames have been acquired). *timeout* can be either a number, None (infinite timeout), or a tuple (timeout, frame\_timeout), in which case the call times out if the total time exceeds timeout, or a single frame wait exceeds frame\_timeout. If the call times out, raise `TimeoutError`. If *error\_on\_stopped==True* and the acquisition is not running, raise `Error`; otherwise, simply return False without waiting.

## Module contents

### pylablib.devices.BitFlow package

#### Submodules

#### pylablib.devices.BitFlow.BitFlow module

**exception** pylablib.devices.BitFlow.BitFlow.**BitFlowError**

Bases: `pylablib.core.devio.base.DeviceError`

Generic BitFlow devices error

**args**

**with\_traceback** ()

Exception.with\_traceback(tb) – set self.\_\_traceback\_\_ to tb and return self.

**exception** pylablib.devices.BitFlow.BitFlow.**BitFlowTimeoutError**

Bases: `pylablib.devices.BitFlow.BitFlow.BitFlowError`

BitFlow frame timeout error

**args**

**with\_traceback()**

Exception.with\_traceback(tb) – set self.\_\_traceback\_\_ to tb and return self.

**class** pylablib.devices.BitFlow.BitFlow.**TDeviceInfo**(*idx, model, idreg*)

Bases: [tuple](#)

**count()**

Return number of occurrences of value.

**idreg**

**idx**

**index()**

Return first index of value.

Raises ValueError if the value is not present.

**model**

pylablib.devices.BitFlow.BitFlow.**list\_cameras()**

List all cameras available through BitFlow interface

pylablib.devices.BitFlow.BitFlow.**get\_cameras\_number()**

Get number of connected BitFlow cameras

**class** pylablib.devices.BitFlow.BitFlow.**BitFlowFrameGrabber**(*bitflow\_idx=0, bitflow\_camfile=None, do\_open=True, \*\*kwargs*)

Bases: [pylablib.devices.interface.camera.IROICamera](#)

Generic BitFlow frame grabber interface.

Compared to [BitFlowCamera](#), has more permissive initialization arguments, which simplifies its use as a base class for expanded cameras.

#### Parameters

- **bitflow\_idx** – board index, starting from 0
- **bitflow\_camfile** – if not None, a path to a valid camera file used for this frame grabber and camera combination; in this case, a temporary camera file is generated based on the provided one and used to change some otherwise unavailable camera parameters such as ROI and pixel bit depth (they are otherwise fixed to whatever is specified in the default camera file)
- **do\_open** – if False, skip the last step of opening the device (should be opened in a subclass)

#### Error

alias of [BitFlowError](#)

#### TimeoutError

alias of [BitFlowTimeoutError](#)

**open()**

Open connection to the camera

**close()**

Close connection to the camera

**is\_opened()**

Check if the device is connected

**get\_device\_info()**

Get camera model data.

Return tuple (*idx*, *model*, *idreg*) with the board index, model number and the setting of the ID switch on the board

**get\_detector\_size()**

Get camera detector size (in pixels) as a tuple (*width*, *height*)

**get\_grabber\_detector\_size()**

Get camera detector size (in pixels) as a tuple (*width*, *height*)

**get\_roi()**

Get current ROI.

Return tuple (*hstart*, *hend*, *vstart*, *vend*). *hstart* and *hend* specify horizontal image extent, *vstart* and *vend* specify vertical image extent (start is inclusive, stop is exclusive, starting from 0).

**get\_grabber\_roi()**

Get current ROI.

Return tuple (*hstart*, *hend*, *vstart*, *vend*). *hstart* and *hend* specify horizontal image extent, *vstart* and *vend* specify vertical image extent (start is inclusive, stop is exclusive, starting from 0).

**set\_roi(*hstart=0, hend=None, vstart=0, vend=None*)**

Setup camera ROI.

*hstart* and *hend* specify horizontal image extent, *vstart* and *vend* specify vertical image extent (start is inclusive, stop is exclusive, starting from 0). By default, all non-supplied parameters take extreme values (0 for start, maximal for end).

**set\_grabber\_roi(*hstart=0, hend=None, vstart=0, vend=None*)**

Setup camera ROI.

*hstart* and *hend* specify horizontal image extent, *vstart* and *vend* specify vertical image extent (start is inclusive, stop is exclusive, starting from 0). By default, all non-supplied parameters take extreme values (0 for start, maximal for end).

**get\_roi\_limits(*hbin=1, vbin=1*)**

Get the minimal and maximal ROI parameters.

Return tuple (*hlim*, *vlim*), where each element is in turn a limit 5-tuple (*min*, *max*, *pstep*, *sstep*, *maxbin*) with, correspondingly, minimal and maximal size, position and size step, and the maximal binning (fixed to 1 if not binning is allowed). In some cameras, the step and the minimal size depend on the binning, which can be supplied.

**get\_grabber\_roi\_limits(*hbin=1, vbin=1*)**

Get the minimal and maximal ROI parameters.

Return tuple (*hlim*, *vlim*), where each element is in turn a limit 5-tuple (*min*, *max*, *pstep*, *sstep*, *maxbin*) with, correspondingly, minimal and maximal size, position and size step, and the maximal binning (fixed to 1 if not binning is allowed). In some cameras, the step and the minimal size depend on the binning, which can be supplied.

**class BufferManager(*cam*)**

Bases: `object`

Buffer manager: stores, constantly reads and re-schedules buffers, keeps track of acquired frames and buffer overflow events

```

reset ()
    Reset counter (on frame acquisition)

start_loop ()
    Start buffer scheduling loop

stop_loop ()
    Stop buffer scheduling loop

is_running ()
    Check if the buffer loop is running

get_status ()
    Get counter status: tuple (acquired,)

setup_acquisition (mode='sequence', nframes=100, frame_merge=1)
    Setup acquisition mode.

    mode can be either "snap" (single frame or a fixed number of frames) or "sequence" (continuous
    acquisition). nframes sets up number of frames in the acquisition buffer. frame_merge specifies the
    number of frames to merge together to form one buffer; if it is larger than 1, several camera frames will
    be merged into a single frame grabber "super-frame" for acquisition, to lower the effective frame rate
    (which is capped at 2-4kFPS due to the necessity of Python loops). This is done transparently for the user,
    so the only visible change is the fact that the number of acquired frames is always updated in quanta of
    frame_merge.

clear_acquisition ()
    Clear all acquisition details and free all buffers

start_acquisition (*args, **kwargs)
    Start acquisition.

    Can take the same keyword parameters as :meth:~setup_acquisition. If the acquisition is not set up yet,
    set it up using the supplied parameters (use default of setup_acquisition(), if the parameter is
    None). Otherwise, if any supplied parameters are different from the current ones, change them and reset
    the acquisition.

stop_acquisition ()
    Stop acquisition

acquisition_in_progress ()
    Check if acquisition is in progress

FrameTransferError
    alias of pylablib.devices.interface.camera.DefaultFrameTransferError

class NoParameterCaller (device, kind)
    Bases: object

    Class to simplify calling functions without a parameter

apply_settings (settings)
    Apply the settings.

    settings is the dict {name: value} of the device available settings. Non-applicable settings are ig-
    nored.

get_acquisition_parameters ()
    Get acquisition parameters.

    Return dictionary {name: value}

```

**get\_data\_dimensions()**

Get readout data dimensions (in pixels) as a tuple (width, height); take indexing mode into account

**get\_device\_variable(key)**

Get the value of a settings, status, or full info parameter

**get\_frame\_format()**

Get format for the returned images.

Can be "list" (list of 2D arrays), "array" (a single 3D array), or "chunks" (list of 3D "chunk" arrays; supported for some cameras and provides the best performance).

**get\_frame\_info\_fields()**

Get the names of frame info fields.

Applicable when frame info format (set by `set_frame_info_format()`) is "list" or "array".

**get\_frame\_info\_format()**

Get format of the frame info.

Can be "namedtuple" (potentially nested named tuples; convenient to get particular values), "list" (flat list of values, with field names are given by `get_frame_info_fields()`; convenient for building a table), "array" (same as "list", but with a numpy array, which is easier to use for "chunks" frame format), or "dict" (flat dictionary with the same fields as the "list" format; more resilient to future format changes)

**get\_frame\_info\_period()**

Get period of frame info acquisition.

Frame info might be skipped (set to None) except for frames which indices are divisible by *period*. Useful for certain cameras where acquiring frame info takes a lot of time and can reduce performance at higher frame rates. Note that this parameter can still be ignored (i.e., always set to 1) if the performance is not an issue for a given camera class.

**get\_frames\_status()**

Get acquisition and buffer status.

Return tuple (acquired, unread, skipped, size), where acquired is the total number of acquired frames, unread is the number of acquired but not read frames, skipped is the number of skipped (not read and then written over) frames, and buffer\_size is the total buffer size (in frames).

**get\_full\_info(include=0)**

Get dict {name: value} containing full device information (including status and settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include*=-10 queries all available variables, which is equivalent to *include*="all".

**get\_full\_status(include=0)**

Get dict {name: value} containing the device status (including settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include*=-10 queries all available variables, which is equivalent to *include*="all".

**get\_image\_indexing()**

Get indexing for the returned images.

Can be "rct" (first index row, second index column, rows counted from the top), "rcb" (same as "rc", rows counted from the bottom), "xyt" (first index column, second index row, rows counted from



the top), or "xyb" (same as "xyt", rows counted from the bottom)

**get\_new\_images\_range()**

Get the range of the new images.

Return tuple (*first*, *last*) with images range (first inclusive). If no images are available, return None. If some images were in the buffer were overwritten, exclude them from the range.

**get\_settings(include=0)**

Get dict {*name*: *value*} containing all the device settings.

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include*=-10 queries all available variables, which is equivalent to *include*="all".

**grab(nframes=1, frame\_timeout=5.0, missing\_frame='skip', return\_info=False, buff\_size=None)**

Snap *nframes* images (with preset image read mode parameters)

*buff\_size* determines buffer size (if None, use the default size). Timeout is specified for a single-frame acquisition, not for the whole acquisition time. *missing\_frame* determines what to do with frames which have been lost: can be "none" (replacing them with None), "zero" (replacing them with zero-filled frame), or "skip" (skipping them, while still keeping total returned frames number to *n*). If *return\_info*==True, return tuple (*frames*, *infos*), where *infos* is a list of frame info tuples (camera-dependent); if some frames are missing and *missing\_frame*!="skip", the corresponding frame info is None.

**is\_acquisition\_setup()**

Check if acquisition is set up.

If the camera does not support separate acquisition setup, always return True.

**pausing\_acquisition(clear=None, stop=True, setup\_after=None, start\_after=True, combine\_nested=True)**

Context manager which temporarily pauses acquisition during execution of *with* block.

Useful for applying certain settings which can't be changed during the acquisition. If *clear*==True, clear acquisition in addition to stopping (by default, use the class default specified as *\_clear\_pausing\_acquisition* attribute). If *stop*==True, stop the acquisition (if *clear*==True, stop regardless). If *setup\_after*==True, setup the acquisition after pause if necessary (None means setup only if clearing was required). If *start\_after*==True, start the acquisition after pause if necessary (None means start only if stopping was required). If *combine\_nested*==True, then any nested *pausing\_acquisition* calls will stop/clear acquisition as necessary, but won't setup/start it again until this *pausing\_acquisition* call is complete.

Yields tuple (*acq\_in\_progress*, *acq\_params*), which indicates whether acquisition is currently in progress, and what are the current acquisition parameters.

**read\_multiple\_images(rng=None, peek=False, missing\_frame='skip', return\_info=False, return\_rng=False)**

Read multiple images specified by *rng* (by default, all un-read images).

If *rng* is specified, it is a tuple (*first*, *last*) with images range (first inclusive). If no new frames are available, return an empty list; if no acquisition is running, return None. If *peek*==True, return images but not mark them as read. *missing\_frame* determines what to do with frames which are out of range (missing or lost): can be "none" (replacing them with None), "zero" (replacing them with zero-filled frame), or "skip" (skipping them). If *return\_info*==True, return tuple (*frames*, *infos*), where *infos* is a list of frame info tuples (camera-dependent, by default, only the frame index); if some frames are missing and *missing\_frame*!="skip", the corresponding frame info is None. if *return\_rng*==True, return the range covered resulting frames; if *missing\_frame*=="skip", the range can be smaller than the supplied *rng* if some frames are skipped.

**read\_newest\_image** (*peek=False, return\_info=False*)

Read the newest un-read image.

If no un-read frames are available, return None. If `peek==True`, return the image but not mark it as read. If `return_info==True`, return tuple (frame, info), where info is an info tuples (camera-dependent, see [read\\_multiple\\_images\(\)](#)).

**read\_oldest\_image** (*peek=False, return\_info=False*)

Read the oldest un-read image.

If no un-read frames are available, return None. If `peek==True`, return the image but not mark it as read. If `return_info==True`, return tuple (frame, info), where info is an info tuples (camera-dependent, see [read\\_multiple\\_images\(\)](#)).

**set\_device\_variable** (*key, value*)

Set the value of a settings parameter

**set\_frame\_format** (*fmt*)

Set format for the returned images.

Can be "list" (list of 2D arrays), "array" (a single 3D array), "chunks" (list of 3D "chunk" arrays; supported for some cameras and provides the best performance), or "try\_chunks" (same as "chunks", but if chunks are not supported, set to "list" instead). If format is "chunks" and chunks are not supported by the camera, it results in one frame per chunk. Note that if the format is set to "array" or "chunks", the frame info format is also automatically set to "array". If the format is set to "chunks", then the image info is also returned in chunks form (list of 2D info arrays with the same length as the corresponding frame chunks).

**set\_frame\_info\_format** (*fmt, include\_fields=None*)

Set format of the frame info.

Can be "namedtuple" (potentially nested named tuples; convenient to get particular values), "list" (flat list of values, with field names are given by [get\\_frame\\_info\\_fields\(\)](#); convenient for building a table), "array" (same as "list", but with a numpy array, which is easier to use for "chunks" frame format), or "dict" (flat dictionary with the same fields as the "list" format; more resilient to future format changes) If `include_fields` is not None, it specifies the fields included for non-"tuple" formats; note that order of `include_fields` is ignored, and the resulting fields are always ordered same as in the original.

**set\_frame\_info\_period** (*period=1*)

Set period of frame info acquisition.

Frame info might be skipped (set to None) except for frames which indices are divisible by *period*. Useful for certain cameras where acquiring frame info takes a lot of time and can reduce performance at higher frame rates. Note that this parameter can still be ignored (i.e., always set to 1) if the performance is not an issue for a given camera class.

**set\_image\_indexing** (*indexing*)

Set up indexing for the returned images.

Can be "rct" (first index row, second index column, rows counted from the top), "rcb" (same as "rc", rows counted from the bottom), "xyt" (first index column, second index row, rows counted from the top), or "xyb" (same as "xyt", rows counted from the bottom)

**snap** (*timeout=5.0, return\_info=False*)

Snap a single frame

**wait\_for\_frame** (*since='lastread', nframes=1, timeout=20.0, error\_on\_stopped=False*)

Wait for one or several new camera frames.

*since* specifies the reference point for waiting to acquire *nframes* frames; can be "lastread" (from the last read frame), "lastwait" (wait for the last successful `wait_for_frame()` call), "now" (from the start of the current call), or "start" (from the acquisition start, i.e., wait until *nframes* frames have been acquired). *timeout* can be either a number, None (infinite timeout), or a tuple (timeout, frame\_timeout), in which case the call times out if the total time exceeds timeout, or a single frame wait exceeds frame\_timeout. If the call times out, raise `TimeoutError`. If `error_on_stopped==True` and the acquisition is not running, raise `Error`; otherwise, simply return False without waiting.

**class** `pylablib.devices.BitFlow.BitFlow.BitFlowCamera (idx=0, camfile=None)`

Bases: `pylablib.devices.BitFlow.BitFlow.BitFlowFrameGrabber`

Generic BitFlow camera interface.

**Parameters** *idx* – board index, starting from 0

**class** `BufferManager (cam)`

Bases: `object`

Buffer manager: stores, constantly reads and re-schedules buffers, keeps track of acquired frames and buffer overflow events

**get\_status ()**

Get counter status: tuple (acquired,)

**is\_running ()**

Check if the buffer loop is running

**reset ()**

Reset counter (on frame acquisition)

**start\_loop ()**

Start buffer scheduling loop

**stop\_loop ()**

Stop buffer scheduling loop

**Error**

alias of `BitFlowError`

**FrameTransferError**

alias of `pylablib.devices.interface.camera.DefaultFrameTransferError`

**class** `NoParameterCaller (device, kind)`

Bases: `object`

Class to simplify calling functions without a parameter

**TimeoutError**

alias of `BitFlowTimeoutError`

**acquisition\_in\_progress ()**

Check if acquisition is in progress

**apply\_settings (settings)**

Apply the settings.

*settings* is the dict {name: value} of the device available settings. Non-applicable settings are ignored.

**clear\_acquisition ()**

Clear all acquisition details and free all buffers

**close()**  
Close connection to the camera

**get\_acquisition\_parameters()**  
Get acquisition parameters.  
  
Return dictionary {name: value}

**get\_data\_dimensions()**  
Get readout data dimensions (in pixels) as a tuple (width, height); take indexing mode into account

**get\_detector\_size()**  
Get camera detector size (in pixels) as a tuple (width, height)

**get\_device\_info()**  
Get camera model data.  
  
Return tuple (idx, model, idreg) with the board index, model number and the setting of the ID switch on the board

**get\_device\_variable(key)**  
Get the value of a settings, status, or full info parameter

**get\_frame\_format()**  
Get format for the returned images.  
  
Can be "list" (list of 2D arrays), "array" (a single 3D array), or "chunks" (list of 3D "chunk" arrays; supported for some cameras and provides the best performance).

**get\_frame\_info\_fields()**  
Get the names of frame info fields.  
  
Applicable when frame info format (set by `set_frame_info_format()`) is "list" or "array".

**get\_frame\_info\_format()**  
Get format of the frame info.  
  
Can be "namedtuple" (potentially nested named tuples; convenient to get particular values), "list" (flat list of values, with field names are given by `get_frame_info_fields()`; convenient for building a table), "array" (same as "list", but with a numpy array, which is easier to use for "chunks" frame format), or "dict" (flat dictionary with the same fields as the "list" format; more resilient to future format changes)

**get\_frame\_info\_period()**  
Get period of frame info acquisition.  
  
Frame info might be skipped (set to None) except for frames which indices are divisible by *period*. Useful for certain cameras where acquiring frame info takes a lot of time and can reduce performance at higher frame rates. Note that this parameter can still be ignored (i.e., always set to 1) if the performance is not an issue for a given camera class.

**get\_frames\_status()**  
Get acquisition and buffer status.  
  
Return tuple (acquired, unread, skipped, size), where acquired is the total number of acquired frames, unread is the number of acquired but not read frames, skipped is the number of skipped (not read and then written over) frames, and buffer\_size is the total buffer size (in frames).

**get\_full\_info(include=0)**  
Get dict {name: value} containing full device information (including status and settings).  
  
*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the

lowest priority is -10, setting `include=-10` queries all available variables, which is equivalent to `include="all"`.

**get\_full\_status** (*include=0*)

Get dict {name: value} containing the device status (including settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting `include=-10` queries all available variables, which is equivalent to `include="all"`.

**get\_grabber\_detector\_size** ()

Get camera detector size (in pixels) as a tuple (width, height)

**get\_grabber\_roi** ()

Get current ROI.

Return tuple (*hstart*, *hend*, *vstart*, *vend*). *hstart* and *hend* specify horizontal image extent, *vstart* and *vend* specify vertical image extent (start is inclusive, stop is exclusive, starting from 0).

**get\_grabber\_roi\_limits** (*hbin=1*, *vbin=1*)

Get the minimal and maximal ROI parameters.

Return tuple (*hlim*, *vlim*), where each element is in turn a limit 5-tuple (*min*, *max*, *pstep*, *sstep*, *maxbin*) with, correspondingly, minimal and maximal size, position and size step, and the maximal binning (fixed to 1 if not binning is allowed). In some cameras, the step and the minimal size depend on the binning, which can be supplied.

**get\_image\_indexing** ()

Get indexing for the returned images.

Can be "rct" (first index row, second index column, rows counted from the top), "rcb" (same as "rc", rows counted from the bottom), "xyt" (first index column, second index row, rows counted from the top), or "xyb" (same as "xyt", rows counted from the bottom)

**get\_new\_images\_range** ()

Get the range of the new images.

Return tuple (*first*, *last*) with images range (first inclusive). If no images are available, return None. If some images were in the buffer were overwritten, exclude them from the range.

**get\_roi** ()

Get current ROI.

Return tuple (*hstart*, *hend*, *vstart*, *vend*). *hstart* and *hend* specify horizontal image extent, *vstart* and *vend* specify vertical image extent (start is inclusive, stop is exclusive, starting from 0).

**get\_roi\_limits** (*hbin=1*, *vbin=1*)

Get the minimal and maximal ROI parameters.

Return tuple (*hlim*, *vlim*), where each element is in turn a limit 5-tuple (*min*, *max*, *pstep*, *sstep*, *maxbin*) with, correspondingly, minimal and maximal size, position and size step, and the maximal binning (fixed to 1 if not binning is allowed). In some cameras, the step and the minimal size depend on the binning, which can be supplied.

**get\_settings** (*include=0*)

Get dict {name: value} containing all the device settings.

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting `include=-10` queries all available variables, which is equivalent to `include="all"`.

**grab** (*nframes=1, frame\_timeout=5.0, missing\_frame='skip', return\_info=False, buff\_size=None*)  
Snap *nframes* images (with preset image read mode parameters)

*buff\_size* determines buffer size (if *None*, use the default size). Timeout is specified for a single-frame acquisition, not for the whole acquisition time. *missing\_frame* determines what to do with frames which have been lost: can be "none" (replacing them with *None*), "zero" (replacing them with zero-filled frame), or "skip" (skipping them, while still keeping total returned frames number to *n*). If *return\_info==True*, return tuple (*frames, infos*), where *infos* is a list of frame info tuples (camera-dependent); if some frames are missing and *missing\_frame!="skip"*, the corresponding frame info is *None*.

**is\_acquisition\_setup** ()  
Check if acquisition is set up.

If the camera does not support separate acquisition setup, always return *True*.

**is\_opened** ()  
Check if the device is connected

**open** ()  
Open connection to the camera

**pausing\_acquisition** (*clear=None, stop=True, setup\_after=None, start\_after=True, combine\_nested=True*)  
Context manager which temporarily pauses acquisition during execution of *with* block.

Useful for applying certain settings which can't be changed during the acquisition. If *clear==True*, clear acquisition in addition to stopping (by default, use the class default specified as *\_clear\_pausing\_acquisition* attribute). If *stop==True*, stop the acquisition (if *clear==True*, stop regardless). If *setup\_after==True*, setup the acquisition after pause if necessary (*None* means setup only if clearing was required). If *start\_after==True*, start the acquisition after pause if necessary (*None* means start only if stopping was required). If *combine\_nested==True*, then any nested *pausing\_acquisition* calls will stop/clear acquisition as necessary, but won't setup/start it again until this *pausing\_acquisition* call is complete.

Yields tuple (*acq\_in\_progress, acq\_params*), which indicates whether acquisition is currently in progress, and what are the current acquisition parameters.

**read\_multiple\_images** (*rng=None, peek=False, missing\_frame='skip', return\_info=False, return\_rng=False*)  
Read multiple images specified by *rng* (by default, all un-read images).

If *rng* is specified, it is a tuple (*first, last*) with images range (first inclusive). If no new frames are available, return an empty list; if no acquisition is running, return *None*. If *peek==True*, return images but not mark them as read. *missing\_frame* determines what to do with frames which are out of range (missing or lost): can be "none" (replacing them with *None*), "zero" (replacing them with zero-filled frame), or "skip" (skipping them). If *return\_info==True*, return tuple (*frames, infos*), where *infos* is a list of frame info tuples (camera-dependent, by default, only the frame index); if some frames are missing and *missing\_frame!="skip"*, the corresponding frame info is *None*. if *return\_rng==True*, return the range covered resulting frames; if *missing\_frame=="skip"*, the range can be smaller than the supplied *rng* if some frames are skipped.

**read\_newest\_image** (*peek=False, return\_info=False*)  
Read the newest un-read image.

If no un-read frames are available, return *None*. If *peek==True*, return the image but not mark it as read. If *return\_info==True*, return tuple (*frame, info*), where *info* is an info tuples (camera-dependent, see [read\\_multiple\\_images\(\)](#)).

**read\_oldest\_image** (*peek=False, return\_info=False*)  
Read the oldest un-read image.

If no un-read frames are available, return None. If `peek==True`, return the image but not mark it as read. If `return_info==True`, return tuple (frame, info), where info is an info tuples (camera-dependent, see `read_multiple_images()`).

**set\_device\_variable** (*key, value*)  
Set the value of a settings parameter

**set\_frame\_format** (*fmt*)  
Set format for the returned images.

Can be "list" (list of 2D arrays), "array" (a single 3D array), "chunks" (list of 3D "chunk" arrays; supported for some cameras and provides the best performance), or "try\_chunks" (same as "chunks", but if chunks are not supported, set to "list" instead). If format is "chunks" and chunks are not supported by the camera, it results in one frame per chunk. Note that if the format is set to "array" or "chunks", the frame info format is also automatically set to "array". If the format is set to "chunks", then the image info is also returned in chunks form (list of 2D info arrays with the same length as the corresponding frame chunks).

**set\_frame\_info\_format** (*fmt, include\_fields=None*)  
Set format of the frame info.

Can be "namedtuple" (potentially nested named tuples; convenient to get particular values), "list" (flat list of values, with field names are given by `get_frame_info_fields()`; convenient for building a table), "array" (same as "list", but with a numpy array, which is easier to use for "chunks" frame format), or "dict" (flat dictionary with the same fields as the "list" format; more resilient to future format changes) If *include\_fields* is not None, it specifies the fields included for non-"tuple" formats; note that order or *include\_fields* is ignored, and the resulting fields are always ordered same as in the original.

**set\_frame\_info\_period** (*period=1*)  
Set period of frame info acquisition.

Frame info might be skipped (set to None) except for frames which indices are divisible by *period*. Useful for certain cameras where acquiring frame info takes a lot of time and can reduce performance at higher frame rates. Note that this parameter can still be ignored (i.e., always set to 1) if the performance is not an issue for a given camera class.

**set\_grabber\_roi** (*hstart=0, hend=None, vstart=0, vend=None*)  
Setup camera ROI.

*hstart* and *hend* specify horizontal image extent, *vstart* and *vend* specify vertical image extent (start is inclusive, stop is exclusive, starting from 0). By default, all non-supplied parameters take extreme values (0 for start, maximal for end).

**set\_image\_indexing** (*indexing*)  
Set up indexing for the returned images.

Can be "rct" (first index row, second index column, rows counted from the top), "rcb" (same as "rc", rows counted from the bottom), "xyt" (first index column, second index row, rows counted from the top), or "xyb" (same as "xyt", rows counted from the bottom)

**set\_roi** (*hstart=0, hend=None, vstart=0, vend=None*)  
Setup camera ROI.

*hstart* and *hend* specify horizontal image extent, *vstart* and *vend* specify vertical image extent (start is inclusive, stop is exclusive, starting from 0). By default, all non-supplied parameters take extreme values (0 for start, maximal for end).

**setup\_acquisition** (*mode='sequence', nframes=100, frame\_merge=1*)  
Setup acquisition mode.

*mode* can be either "snap" (single frame or a fixed number of frames) or "sequence" (continuous acquisition). *nframes* sets up number of frames in the acquisition buffer. *frame\_merge* specifies the number of frames to merge together to form one buffer; if it is larger than 1, several camera frames will be merged into a single frame grabber "super-frame" for acquisition, to lower the effective frame rate (which is capped at 2-4kFPS due to the necessity of Python loops). This is done transparently for the user, so the only visible change is the fact that the number of acquired frames is always updated in quanta of *frame\_merge*.

**snap** (*timeout=5.0, return\_info=False*)  
Snap a single frame

**start\_acquisition** (*\*args, \*\*kwargs*)  
Start acquisition.

Can take the same keyword parameters as `:meth:~setup_acquisition`. If the acquisition is not set up yet, set it up using the supplied parameters (use default of `setup_acquisition()`, if the parameter is None). Otherwise, if any supplied parameters are different from the current ones, change them and reset the acquisition.

**stop\_acquisition** ()  
Stop acquisition

**wait\_for\_frame** (*since='lastread', nframes=1, timeout=20.0, error\_on\_stopped=False*)  
Wait for one or several new camera frames.

*since* specifies the reference point for waiting to acquire *nframes* frames; can be "lastread" (from the last read frame), "lastwait" (wait for the last successful `wait_for_frame()` call), "now" (from the start of the current call), or "start" (from the acquisition start, i.e., wait until *nframes* frames have been acquired). *timeout* can be either a number, None (infinite timeout), or a tuple (*timeout, frame\_timeout*), in which case the call times out if the total time exceeds *timeout*, or a single frame wait exceeds *frame\_timeout*. If the call times out, raise `TimeoutError`. If `error_on_stopped==True` and the acquisition is not running, raise `Error`; otherwise, simply return False without waiting.

**class** `pylablib.devices.BitFlow.BitFlow.CameraFileEditor`

Bases: `object`

Camera file editor based on XML ElementTree parser.

Provides methods for loading and saving the tree, and to change basic parameters in the default operational mode.

**load** (*path, clean=True*)  
Load file from the given path and optionally check the structure remove the non-default modes

**save** (*path*)  
Save file to the given path

**clean\_modes** ()  
Check the loaded tree structure and remove non-default operational modes

**get\_mode\_parameters** ()  
Get default operational mode parameters.

Return tuple (*size, fmt, bpp*) with the acquisition size (*xsize, ysize*), format (e.g., "1X2-1Y") and the number of bits per pixel. If the tree is not loaded or mode is not present, return None.

**set\_mode\_parameters** (*size=None, fmt=None, bpp=None*)  
Get default operational mode parameters.



*size* is the acquisition size (*xsize*, *ysize*), *fmt* is the tap format (e.g., "1X2-1Y"), and *bpp* is the number of bits per pixel. Parameters set to `None` stay unchanged. Return `True` if any parameters have changed their values and `False` otherwise.

## Module contents

### pylablib.devices.Conrad package

#### Submodules

#### pylablib.devices.Conrad.base module

**exception** `pylablib.devices.Conrad.base.ConradError`

Bases: `pylablib.core.devio.base.DeviceError`

Generic Conrad devices error

**args**

**with\_traceback()**

Exception.with\_traceback(tb) – set self.\_\_traceback\_\_ to tb and return self.

**exception** `pylablib.devices.Conrad.base.ConradBackendError` (*exc*)

Bases: `pylablib.devices.Conrad.base.ConradError`, `pylablib.core.devio.comm_backend.DeviceBackendError`

Generic Conrad backend communication error

**args**

**with\_traceback()**

Exception.with\_traceback(tb) – set self.\_\_traceback\_\_ to tb and return self.

**class** `pylablib.devices.Conrad.base.RelayBoard` (*conn*, *start\_addr=1*)

Bases: `pylablib.core.devio.comm_backend.ICommBackendWrapper`

Conrad relay board controller

#### Parameters

- **conn** – serial connection parameters (usually port or a tuple containing port and baudrate)
- **start\_addr** – address which is assigned to the first board in the chain upon initialization; all following boards increase the address by 1

#### Error

alias of `ConradError`

**open()**

Open the connection to the board

**class** `TMessage` (*comm*, *addr*, *data*)

Bases: `tuple`

**addr**

**comm**

**count()**

Return number of occurrences of value.

**data**

**index()**  
Return first index of value.  
Raises ValueError if the value is not present.

**query** (*comm, addr=1, data=0, multi\_result=False*)  
Send a query with the given command, address and data.  
If *multi\_result==False*, read a single reply frame; otherwise, keep reading until reply with the same command as sent is received (used in initialization and broadcast queries).

**get\_all\_relays** (*addr=1*)  
Get all relay states.  
If *addr* is not 0, return dictionary {*relay*:*value*}, where *relay* is the relay index on the board (between 1 and 8 inclusive). If *addr==0* (broadcast), return dictionary {*addr*:*board\_state*}, where *board\_state* is in turn a state dictionary is described above.

**set\_all\_relays** (*values, addr=1*)  
Set all relay states.  
*values* can be a list (listing relay states from lowest to highest), or a dictionary {*relay*:*value*}, where relays are numbered from 1 to 8. Relays without values are kept unchanged. If *addr==0*, broadcast to all boards

**get\_relay** (*relay, addr=1*)  
Get the state at a given relay (indexed from 1 to 8 inclusive)

**set\_relay** (*relay, enable=True, addr=1*)  
Get the state at a given relay (indexed from 1 to 8 inclusive)

**class NoParameterCaller** (*device, kind*)  
Bases: `object`  
Class to simplify calling functions without a parameter

**apply\_settings** (*settings*)  
Apply the settings.  
*settings* is the dict {*name*: *value*} of the device available settings. Non-applicable settings are ignored.

**close()**  
Close the backend

**get\_device\_variable** (*key*)  
Get the value of a settings, status, or full info parameter

**get\_full\_info** (*include=0*)  
Get dict {*name*: *value*} containing full device information (including status and settings).  
*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_full\_status** (*include=0*)  
Get dict {*name*: *value*} containing the device status (including settings).  
*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the

lowest priority is -10, setting `include=-10` queries all available variables, which is equivalent to `include="all"`.

**get\_settings** (*include=0*)

Get dict {name: value} containing all the device settings.

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting `include=-10` queries all available variables, which is equivalent to `include="all"`.

**is\_opened** ()

Check if the device is connected

**lock** (*timeout=None*)

Lock the access to the device from other threads/processes (isn't necessarily implemented)

**locking** (*timeout=None*)

Context manager for lock & unlock

**set\_device\_variable** (*key, value*)

Set the value of a settings parameter

**unlock** ()

Unlock the access to the device from other threads/processes (isn't necessarily implemented)

## Module contents

### pylablib.devices.Cryocon package

#### Submodules

#### pylablib.devices.Cryocon.base module

**exception** `pylablib.devices.Cryocon.base.CryoconError`

Bases: `pylablib.core.devio.base.DeviceError`

Generic Cryocon devices error

**args**

**with\_traceback** ()

Exception.with\_traceback(tb) – set self.\_\_traceback\_\_ to tb and return self.

**exception** `pylablib.devices.Cryocon.base.CryoconBackendError` (*exc*)

Bases: `pylablib.devices.Cryocon.base.CryoconError`, `pylablib.core.devio.comm_backend.DeviceBackendError`

Generic Lakeshore backend communication error

**args**

**with\_traceback** ()

Exception.with\_traceback(tb) – set self.\_\_traceback\_\_ to tb and return self.

**class** `pylablib.devices.Cryocon.base.Cryocon1x` (*conn, nchannels='auto'*)

Bases: `pylablib.core.devio.SCPISCPIDevice`

Cryocon 1x series (12C, 14C, 18C) temperature controller.

**Parameters** `conn` – serial connection parameters (usually port or a tuple containing port and baudrate)

**Error**

alias of *CryoconError*

**ReraiseError**

alias of *CryoconBackendError*

**get\_number\_of\_channels** ()

Return total number of channels in the device (2, 4, or 8)

**get\_display\_units** (*channel*)

**set\_display\_units** (*channel*, *units*)

**get\_temperature** (*channel*, *display\_units=False*)

Get a reading on a given channel.

If `display_units==True`, return reading in the display units; otherwise, return reading in Kelvin. If in this case the display units are "S" (sensor), set them to Kelvin to get the reading. If sensor is disconnected, return `None`.

**get\_all\_temperatures** (*display\_units=False*)

Get readings on all channels.

If `display_units==True`, return reading in the display units; otherwise, return reading in Kelvin. If in this case the display units are "S" (sensor), set them to Kelvin to get the reading. If sensor is disconnected, return `None`.

**get\_sensor\_reading** (*channel*)

Get readings (in sensor units) on a given channel (1 to 8)

**get\_all\_sensor\_readings** ()

Get readings (in sensor units) on all channels

**get\_sensor\_kind** (*channel*)

Get sensor kind of a given channel (1 to 8)

**get\_all\_sensor\_kinds** ()

Get readings (in sensor units) on all channels

**set\_sensor\_kind** (*channel*, *kind*)

Set sensor kind of a given channel (1 to 8).

Can be an integer using internal classification (see manual), or one of "none", "S900", "DT670", "DT470", "S950", "SI410", "Pt100", "Pt1k", "Pt10k", "ThFe", "RO105", "RO600". Setting kind to "none" disables the sensor.

**BackendError**

alias of *pylablib.core.devio.comm\_backend.DeviceBackendError*

**class NoParameterCaller** (*device*, *kind*)

Bases: *object*

Class to simplify calling functions without a parameter

**apply\_settings** (*settings*)

Apply the settings.

*settings* is a dict {`name`: `value`} of the available device settings. Non-applicable settings are ignored.

**ask** (*msg*, *data\_type='string'*, *delay=0.0*, *timeout=None*, *read\_echo=False*)

Write a message and read a reply.

*msg* is the query message, *delay* is the delay between write and read. Other parameters are the same as in `read()`. If `read_echo==True`, assume that the device first echoes the input and skip it.

**close()**

Close the backend

**flush** (*one\_line=False*)

Flush the read buffer (read all the available data and return the number of bytes read).

If `one_line==True`, read only a single line.

**static get\_arg\_type** (*arg*)

Autodetect argument type

**get\_device\_variable** (*key*)

Get the value of a settings, status, or full info parameter

**get\_esr** (*timeout=None*)

Get the device status register (by default, "\*ESR?" command)

**get\_full\_info** (*include=0*)

Get dict {name: value} containing full device information (including status and settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting `include=-10` queries all available variables, which is equivalent to `include="all"`.

**get\_full\_status** (*include=0*)

Get dict {name: value} containing the device status (including settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting `include=-10` queries all available variables, which is equivalent to `include="all"`.

**get\_id** (*timeout=None*)

Get the device IDN. (query SCPI '\*IDN?' command)

**get\_settings** (*include=0*)

Get dict {name: value} containing all the device settings.

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting `include=-10` queries all available variables, which is equivalent to `include="all"`.

**is\_opened()**

Check if the device is connected

**lock** (*timeout=None*)

Lock the access to the device from other threads/processes (isn't necessarily implemented)

**locking** (*timeout=None*)

Context manager for lock & unlock

**open()**

Open the backend

**static parse\_array\_data** (*data, fmt, include\_header=False*)

Parse the data returned by the device. *fmt* is `DataFormat` description in numpy format (e.g., "<u2").

If `include_header==True`, the data is assumed to be in a (somewhat) standard SCPI format: `b'#'`, then a single digit `s` denoting length of the size block, then `s` digits denoting length of the data (in bytes) followed by the actual data. Otherwise (`include_header==False`), assume that the header is already removed.

**read** (*data\_type='string', timeout=None*)

Read data from the device.

*data\_type* determines the type of the data. Can be `'raw'` (just raw data), `'string'` (with trailing and leading spaces stripped), `'int'`, `'float'`, `'bool'` (interprets 0 or `'off'` as `False`, anything else as `True`), `'value'` (returns tuple (*value*, *unit*), where *value* is float), a callable (return the result of this callable applied to the string value), a dictionary (return the stored value corresponding to the string value, or to the value converted into integer if the string value is not present), or a list of data types (the result is treated as a list of values with the given types separated by commas). *timeout* overrides the default value.

**read\_binary\_array\_data** (*include\_header=False, timeout=None, flush\_term=True*)

Read a binary data in the from the device.

The data assumes the standard binary transfer header consisting of `"#"` symbol, then a single digit with the size of the length string, then the length string containing the length of the binary data (in bytes). If `include_header==True`, return the data with the header; otherwise, return only the content. If `flush_term==True`, flush the following line to skip terminator characters after the binary data, which are added by some devices. *timeout* overrides the default value.

**reconnect** (*new\_instrument=True, ignore\_error=True*)

Remake the connection.

If `new_instrument==True`, create a new backend instance. If `ignore_error==True`, ignore errors on closing.

**reset** ()

Reset the device (by default, `"*RST"` command)

**set\_device\_variable** (*key, value*)

Set the value of a settings parameter

**sleep** (*delay*)

Wait for *delay* seconds

**unlock** ()

Unlock the access to the device from other threads/processes (isn't necessarily implemented)

**using\_write\_buffer** ()

Context manager for using a write buffer.

While it's active, all the consecutive `write()` operations are bundled together with `;` delimiter. The actual write is performed at the `read()/ask()` operation or at the end of the block.

**wait** (*wait\_type='sync', timeout=None, wait\_callback=None*)

Pause execution until device overlapped commands are complete.

*wait\_type* is either `'sync'` (perform `wait_sync()`), `'dev'` (perform `wait_dev()`) or `'none'` (do nothing).

**wait\_dev** ()

Pause execution of the device commands until device overlapped commands (e.g., taking sweeps) are complete.

Note that the code execution is not paused.

**wait\_sync** (*timeout=None, wait\_callback=None*)

Pause execution of the script until device overlapped commands (e.g., taking sweeps) are complete.

*timeout* and *wait\_callback* override default constructor parameters.

**write** (*msg, arg=None, arg\_type=None, unit=None, bool\_selector=None, wait\_sync=None, read\_echo=False, read\_echo\_delay=0.0*)

Send a command.

#### Parameters

- **msg** (*str*) – Text message.
- **arg** – Optional argument to append in the end. If a list of arguments is supplied, the result is joined with ", ".
- **arg\_type** (*str*) – Argument type. Can be 'raw' (in which case data is sent raw), 'string', 'int', 'float', 'bool', a format string (such as '{:.3f}') or a list of argument types (for an iterable argument); if format string is used and the argument is a list or a tuple, then it is expanded as a list of arguments (e.g., `arg_type='{0};{1}'` with `arg=[1, 2]` will produce a string '1; 2'); if a list of types is used, each element of *arg* is converted using the corresponding type, and the result is joined with ", ".
- **unit** (*str*) – If not None, use it as a unit to append after the value.
- **bool\_selector** (*tuple*) – A tuple (*false\_value*, *true\_value*) of two strings to represent bool argument; by default, use `._bool_selector` attribute.
- **wait\_sync** – if True, append the sync command (specified as `._wait_sync_comm` attribute, `"*OPC?"` by default) after the message and pause the execution command is complete; useful in long set operations, where the device might ignore later inputs until the current command is complete; if None, use the class default `._default_write_sync` attribute (False by default).
- **read\_echo** (*bool*) – If True, read a single line after write.
- **read\_echo\_delay** (*float*) – The delay between write and read if `read_echo==True`.

## Module contents

### pylablib.devices.Cryomagnetics package

#### Submodules

### pylablib.devices.Cryomagnetics.base module

**exception** `pylablib.devices.Cryomagnetics.base.CryomagneticsError`

Bases: `pylablib.core.devio.base.DeviceError`

Generic Cryomagnetics devices error

**args**

**with\_traceback** ()

Exception.with\_traceback(tb) – set self.\_\_traceback\_\_ to tb and return self.

**exception** `pylablib.devices.Cryomagnetics.base.CryomagneticsBackendError` (*exc*)

Bases: `pylablib.devices.Cryomagnetics.base.CryomagneticsError`, `pylablib.core.devio.comm_backend.DeviceBackendError`

Generic Cryomagnetics backend communication error

**args**

**with\_traceback** ()

Exception.with\_traceback(tb) – set self.\_\_traceback\_\_ to tb and return self.

**class** `pylablib.devices.Cryomagnetics.base.LM500` (*conn*)

Bases: `pylablib.core.devio.SCPIDevice`

Cryomagnetics LM500/510 level monitor.

Channels are enumerated from 1. To abort filling or reset a timeout, call `SCPIDevice.reset()` method.

**Parameters** **conn** – serial connection parameters (usually port or a tuple containing port and baudrate)

**Error**

alias of `CryomagneticsError`

**ReraiseError**

alias of `CryomagneticsBackendError`

**close** ()

Close connection to the device

**get\_channel** ()

Get current measurement channel

**select\_channel** (*channel=1*)

Select the current measurement channel

**get\_type** (*channel=None*)

Get type of a given channel ("lhc" or "ln")

**get\_mode** (*channel=None*)

Get measurement mode at the given channel (None for the currently selected channel).

Can be either 'sample\_hold', or 'continuous'.

**set\_mode** (*mode, channel=None*)

Set measurement mode at the given channel (None for the current channel).

Can be either 'sample\_hold', or 'continuous'.

**get\_interval** (*channel=None*)

Get measurement interval (in seconds) in sample/hold mode at the given channel (None for the current channel)

**set\_interval** (*intvl, channel=None*)

Set measurement interval (in seconds) in sample/hold mode at the given channel (None for the current channel)

**start\_measurement** (*channel=None*)

Initialize measurement on a given channel

**wait\_for\_measurement** (*channel=None, timeout=None*)

Wait for the measurement on a given channel to finish

**get\_level** (*channel=None*)

Get level reading on a given channel



**measure\_level** (*channel=None*)  
 Measure the level (perform the measurement and return the result) on a given channel

**start\_fill** (*channel=None*)  
 Initialize filling at a given channel (*None* for the current channel)

**get\_fill\_status** (*channel=None*)  
 Get filling status at a given channels (*None* for the current channel).  
 Return either "off" (filling is off), "timeout" (filling timed out) or a float (time since filling started, in seconds).

**get\_low\_level** (*channel=None*)  
 Get low level (automated refill start) setting on a given channel (*None* for the current channel)

**set\_low\_level** (*level, channel=None*)  
 Set low level (automated refill start) setting on a given channel (*None* for the current channel)

**get\_high\_level** (*channel=None*)  
 Get high level (automated refill stop) setting on a given channel (*None* for the current channel)

**set\_high\_level** (*level, channel=None*)  
 Set high level (automated refill stop) setting on a given channel (*None* for the current channel)

**BackendError**  
 alias of `pylablib.core.devio.comm_backend.DeviceBackendError`

**class NoParameterCaller** (*device, kind*)  
 Bases: `object`  
 Class to simplify calling functions without a parameter

**apply\_settings** (*settings*)  
 Apply the settings.  
*settings* is a dict {name: value} of the available device settings. Non-applicable settings are ignored.

**ask** (*msg, data\_type='string', delay=0.0, timeout=None, read\_echo=False*)  
 Write a message and read a reply.  
*msg* is the query message, *delay* is the delay between write and read. Other parameters are the same as in `read()`. If `read_echo==True`, assume that the device first echoes the input and skip it.

**flush** (*one\_line=False*)  
 Flush the read buffer (read all the available data and return the number of bytes read).  
 If `one_line==True`, read only a single line.

**static get\_arg\_type** (*arg*)  
 Autodetect argument type

**get\_device\_variable** (*key*)  
 Get the value of a settings, status, or full info parameter

**get\_esr** (*timeout=None*)  
 Get the device status register (by default, "\*ESR?" command)

**get\_full\_info** (*include=0*)  
 Get dict {name: value} containing full device information (including status and settings).  
*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting `include=-10` queries all available variables, which is equivalent to `include="all"`.

**get\_full\_status** (*include=0*)

Get dict {name: value} containing the device status (including settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_id** (*timeout=None*)

Get the device IDN. (query SCPI '\*IDN?' command)

**get\_settings** (*include=0*)

Get dict {name: value} containing all the device settings.

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**is\_opened** ()

Check if the device is connected

**lock** (*timeout=None*)

Lock the access to the device from other threads/processes (isn't necessarily implemented)

**locking** (*timeout=None*)

Context manager for lock & unlock

**open** ()

Open the backend

**static parse\_array\_data** (*data, fmt, include\_header=False*)

Parse the data returned by the device. *fmt* is *DataFormat* description in numpy format (e.g., "<u2").

If *include\_header==True*, the data is assumed to be in a (somewhat) standard SCPI format: b'#', then a single digit *s* denoting length of the size block, then *s* digits denoting length of the data (in bytes) followed by the actual data. Otherwise (*include\_header==False*), assume that the header is already removed.

**read** (*data\_type='string', timeout=None*)

Read data from the device.

*data\_type* determines the type of the data. Can be 'raw' (just raw data), 'string' (with trailing and leading spaces stripped), 'int', 'float', 'bool' (interprets 0 or 'off' as False, anything else as True), 'value' (returns tuple (value, unit), where *value* is float), a callable (return the result of this callable applied to the string value), a dictionary (return the stored value corresponding to the string value, or to the value converted into integer if the string value is not present), or a list of data types (the result is treated as a list of values with the given types separated by commas). *timeout* overrides the default value.

**read\_binary\_array\_data** (*include\_header=False, timeout=None, flush\_term=True*)

Read a binary data in the from the device.

The data assumes the standard binary transfer header consisting of "#" symbol, then a single digit with the size of the length string, then the length string containing the length of the binary data (in bytes). If *include\_header==True*, return the data with the header; otherwise, return only the content. If *flush\_term==True*, flush the following line to skip terminator characters after the binary data, which are added by some devices. *timeout* overrides the default value.

**reconnect** (*new\_instrument=True, ignore\_error=True*)

Remake the connection.

If `new_instrument==True`, create a new backend instance. If `ignore_error==True`, ignore errors on closing.

**reset()**

Reset the device (by default, "`*RST`" command)

**set\_device\_variable**(*key, value*)

Set the value of a settings parameter

**sleep**(*delay*)

Wait for *delay* seconds

**unlock()**

Unlock the access to the device from other threads/processes (isn't necessarily implemented)

**using\_write\_buffer()**

Context manager for using a write buffer.

While it's active, all the consecutive `write()` operations are bundled together with `;` delimiter. The actual write is performed at the `read()/ask()` operation or at the end of the block.

**wait**(*wait\_type='sync', timeout=None, wait\_callback=None*)

Pause execution until device overlapped commands are complete.

*wait\_type* is either `'sync'` (perform `wait_sync()`), `'dev'` (perform `wait_dev()`) or `'none'` (do nothing).

**wait\_dev()**

Pause execution of the device commands until device overlapped commands (e.g., taking sweeps) are complete.

Note that the code execution is not paused.

**wait\_sync**(*timeout=None, wait\_callback=None*)

Pause execution of the script until device overlapped commands (e.g., taking sweeps) are complete.

*timeout* and *wait\_callback* override default constructor parameters.

**write**(*msg, arg=None, arg\_type=None, unit=None, bool\_selector=None, wait\_sync=None, read\_echo=False, read\_echo\_delay=0.0*)

Send a command.

#### Parameters

- **msg**(*str*) – Text message.
- **arg** – Optional argument to append in the end. If a list of arguments is supplied, the result is joined with `" "`.
- **arg\_type**(*str*) – Argument type. Can be `'raw'` (in which case data is sent raw), `'string'`, `'int'`, `'float'`, `'bool'`, a format string (such as `'{: .3f}'`) or a list of argument types (for an iterable argument); if format string is used and the argument is a list or a tuple, then it is expanded as a list of arguments (e.g., `arg_type='{0};{1}'` with `arg=[1, 2]` will produce a string `'1; 2'`); if a list of types is used, each element of *arg* is converted using the corresponding type, and the result is joined with `" "`.
- **unit**(*str*) – If not `None`, use it as a unit to append after the value.
- **bool\_selector**(*tuple*) – A tuple (*false\_value*, *true\_value*) of two strings to represent bool argument; by default, use `._bool_selector` attribute.
- **wait\_sync** – if `True`, append the sync command (specified as `._wait_sync_comm` attribute, `"*OPC?"` by default) after the message and

pause the execution command is complete; useful in long set operations, where the device might ignore later inputs until the current command is complete; if `None`, use the class default `._default_write_sync` attribute (`False` by default).

- **read\_echo** (*bool*) – If `True`, read a single line after write.
- **read\_echo\_delay** (*float*) – The delay between write and read if `read_echo==True`.

**class** `pylablib.devices.Cryomagnetics.base.LM510` (*conn*)

Bases: `pylablib.devices.Cryomagnetics.base.LM500`

Cryomagnetics 510 level monitor.

Compared to `LM500`, adds additional specific methods to enable/disable automatic refill.

Channels are enumerated from 1. To abort filling or reset a timeout, call `SCPIDevice.reset()` method.

**Parameters** `conn` – serial connection parameters (usually port or a tuple containing port and baudrate)

**set\_control\_mode** (*mode*, *channel=None*)

Set automated refill mode on a given channel (`None` for the current channel); can be "off" or "auto"

**BackendError**

alias of `pylablib.core.devio.comm_backend.DeviceBackendError`

**Error**

alias of `CryomagneticsError`

**class** `NoParameterCaller` (*device*, *kind*)

Bases: `object`

Class to simplify calling functions without a parameter

**ReraiseError**

alias of `CryomagneticsBackendError`

**apply\_settings** (*settings*)

Apply the settings.

*settings* is a dict {`name`: `value`} of the available device settings. Non-applicable settings are ignored.

**ask** (*msg*, *data\_type='string'*, *delay=0.0*, *timeout=None*, *read\_echo=False*)

Write a message and read a reply.

*msg* is the query message, *delay* is the delay between write and read. Other parameters are the same as in `read()`. If `read_echo==True`, assume that the device first echoes the input and skip it.

**close** ()

Close connection to the device

**flush** (*one\_line=False*)

Flush the read buffer (read all the available data and return the number of bytes read).

If `one_line==True`, read only a single line.

**static get\_arg\_type** (*arg*)

Autodetect argument type

**get\_channel** ()

Get current measurement channel

**get\_device\_variable** (*key*)

Get the value of a settings, status, or full info parameter

**get\_esr** (*timeout=None*)

Get the device status register (by default, "\*ESR?" command)

**get\_fill\_status** (*channel=None*)

Get filling status at a given channels (None for the current channel).

Return either "off" (filling is off), "timeout" (filling timed out) or a float (time since filling started, in seconds).

**get\_full\_info** (*include=0*)

Get dict {name: value} containing full device information (including status and settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_full\_status** (*include=0*)

Get dict {name: value} containing the device status (including settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_high\_level** (*channel=None*)

Get high level (automated refill stop) setting on a given channel (None for the current channel)

**get\_id** (*timeout=None*)

Get the device IDN. (query SCPI '\*IDN?' command)

**get\_interval** (*channel=None*)

Get measurement interval (in seconds) in sample/hold mode at the given channel (None for the current channel)

**get\_level** (*channel=None*)

Get level reading on a given channel

**get\_low\_level** (*channel=None*)

Get low level (automated refill start) setting on a given channel (None for the current channel)

**get\_mode** (*channel=None*)

Get measurement mode at the given channel (None for the currently selected channel).

Can be either 'sample\_hold', or 'continuous'.

**get\_settings** (*include=0*)

Get dict {name: value} containing all the device settings.

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_type** (*channel=None*)

Get type of a given channel ("lhc" or "ln")

**is\_opened** ()

Check if the device is connected

**lock** (*timeout=None*)

Lock the access to the device from other threads/processes (isn't necessarily implemented)

**locking** (*timeout=None*)

Context manager for lock & unlock

**measure\_level** (*channel=None*)

Measure the level (perform the measurement and return the result) on a given channel

**open** ()

Open the backend

**static parse\_array\_data** (*data, fmt, include\_header=False*)

Parse the data returned by the device. *fmt* is *DataFormat* description in numpy format (e.g., "<u2").

If *include\_header==True*, the data is assumed to be in a (somewhat) standard SCPI format: b'#', then a single digit *s* denoting length of the size block, then *s* digits denoting length of the data (in bytes) followed by the actual data. Otherwise (*include\_header==False*), assume that the header is already removed.

**read** (*data\_type='string', timeout=None*)

Read data from the device.

*data\_type* determines the type of the data. Can be 'raw' (just raw data), 'string' (with trailing and leading spaces stripped), 'int', 'float', 'bool' (interprets 0 or 'off' as False, anything else as True), 'value' (returns tuple (value, unit), where *value* is float), a callable (return the result of this callable applied to the string value), a dictionary (return the stored value corresponding to the string value, or to the value converted into integer if the string value is not present), or a list of data types (the result is treated as a list of values with the given types separated by commas). *timeout* overrides the default value.

**read\_binary\_array\_data** (*include\_header=False, timeout=None, flush\_term=True*)

Read a binary data in the from the device.

The data assumes the standard binary transfer header consisting of "#" symbol, then a single digit with the size of the length string, then the length string containing the length of the binary data (in bytes). If *include\_header==True*, return the data with the header; otherwise, return only the content. If *flush\_term==True*, flush the following line to skip terminator characters after the binary data, which are added by some devices. *timeout* overrides the default value.

**reconnect** (*new\_instrument=True, ignore\_error=True*)

Remake the connection.

If *new\_instrument==True*, create a new backend instance. If *ignore\_error==True*, ignore errors on closing.

**reset** ()

Reset the device (by default, "\*RST" command)

**select\_channel** (*channel=1*)

Select the current measurement channel

**set\_device\_variable** (*key, value*)

Set the value of a settings parameter

**set\_high\_level** (*level, channel=None*)

Set high level (automated refill stop) setting on a given channel (None for the current channel)

**set\_interval** (*intvl, channel=None*)

Set measurement interval (in seconds) in sample/hold mode at the given channel (None for the current channel)

**set\_low\_level** (*level, channel=None*)

Set low level (automated refill start) setting on a given channel (None for the current channel)

**set\_mode** (*mode*, *channel=None*)

Set measurement mode at the given channel (None for the current channel).

Can be either 'sample\_hold', or 'continuous'.

**sleep** (*delay*)

Wait for *delay* seconds

**start\_fill** (*channel=None*)

Initialize filling at a given channel (None for the current channel)

**start\_measurement** (*channel=None*)

Initialize measurement on a given channel

**unlock** ()

Unlock the access to the device from other threads/processes (isn't necessarily implemented)

**using\_write\_buffer** ()

Context manager for using a write buffer.

While it's active, all the consecutive `write()` operations are bundled together with ; delimiter. The actual write is performed at the `read()/ask()` operation or at the end of the block.

**wait** (*wait\_type='sync'*, *timeout=None*, *wait\_callback=None*)

Pause execution until device overlapped commands are complete.

*wait\_type* is either 'sync' (perform `wait_sync()`), 'dev' (perform `wait_dev()`) or 'none' (do nothing).

**wait\_dev** ()

Pause execution of the device commands until device overlapped commands (e.g., taking sweeps) are complete.

Note that the code execution is not paused.

**wait\_for\_measurement** (*channel=None*, *timeout=None*)

Wait for the measurement on a given channel to finish

**wait\_sync** (*timeout=None*, *wait\_callback=None*)

Pause execution of the script until device overlapped commands (e.g., taking sweeps) are complete.

*timeout* and *wait\_callback* override default constructor parameters.

**write** (*msg*, *arg=None*, *arg\_type=None*, *unit=None*, *bool\_selector=None*, *wait\_sync=None*, *read\_echo=False*, *read\_echo\_delay=0.0*)

Send a command.

#### Parameters

- **msg** (*str*) – Text message.
- **arg** – Optional argument to append in the end. If a list of arguments is supplied, the result is joined with ", ".
- **arg\_type** (*str*) – Argument type. Can be 'raw' (in which case data is sent raw), 'string', 'int', 'float', 'bool', a format string (such as '{:.3f}') or a list of argument types (for an iterable argument); if format string is used and the argument is a list or a tuple, then it is expanded as a list of arguments (e.g., `arg_type='{0};{1}'` with `arg=[1, 2]` will produce a string '1; 2'); if a list of types is used, each element of *arg* is converted using the corresponding type, and the result is joined with ", ".
- **unit** (*str*) – If not None, use it as a unit to append after the value.

- **bool\_selector** (*tuple*) – A tuple (*false\_value*, *true\_value*) of two strings to represent bool argument; by default, use `._bool_selector` attribute.
- **wait\_sync** – if `True`, append the sync command (specified as `._wait_sync_comm` attribute, `"*OPC?"` by default) after the message and pause the execution command is complete; useful in long set operations, where the device might ignore later inputs until the current command is complete; if `None`, use the class default `._default_write_sync` attribute (`False` by default).
- **read\_echo** (*bool*) – If `True`, read a single line after write.
- **read\_echo\_delay** (*float*) – The delay between write and read if `read_echo==True`.

## Module contents

### pylablib.devices.DCAM package

#### Submodules

#### pylablib.devices.DCAM.DCAM module

**class** `pylablib.devices.DCAM.DCAM.LibraryController` (*lib*)

Bases: `pylablib.devices.utils.load_lib.LibraryController`

**close** (*opid*)

Mark device closing.

Return tuple (*close\_result*, *uninit\_result*) with the results of the closing and the shutdown. If library does not need to be shut down yet, set *uninit\_result*=None

**open** ()

Mark device opening.

Return tuple (*init\_result*, *open\_result*, *opid*) with the results of the initialization and the opening, and the opening ID which should afterwards be used for closing. If library is already initialized, set *init\_result*=None

**preinit** ()

Pre-initialize the library, if it hasn't been done already

**shutdown** ()

Close all opened connections and shutdown the library

**temp\_open** ()

Context for temporarily opening a new device connection

`pylablib.devices.DCAM.DCAM.restart_lib` ()

`pylablib.devices.DCAM.DCAM.get_cameras_number` ()

Get number of connected DCAM cameras

**class** `pylablib.devices.DCAM.DCAM.DCAMAttribute` (*handle*, *pid*)

Bases: `object`

DCAM camera attribute.

Allows to query and set values and get additional information. Usually created automatically by a DCAM camera instance, but could also be created manually.



**Parameters**

- **handle** – DCAM camera handle
- **pid** – attribute id

**name**

attribute name

**kind**

attribute kind; can be "int", "float", "enum", or "none" (can't determine)

**Type** `str`

**readable**

whether attribute is readable

**Type** `bool`

**writable**

whether attribute is writable

**Type** `bool`

**min**

minimal attribute value (if applicable)

**Type** `float`

**max**

maximal attribute value (if applicable)

**Type** `float`

**step**

attribute value step (if applicable)

**Type** `float`

**unit**

attribute units (index value)

**Type** `int`

**ivalues**

list of possible integer values for enum attributes

**values**

list of possible text values for enum attributes

**labels**

dict {label: index} which shows all possible values of an enumerated attribute and their corresponding numerical values

**ilabels**

dict {index: label} which shows labels corresponding to numerical values of an enumerated attribute

**as\_text** (*value=None*)

Get the given attribute value as text (by default, current value)

**update\_limits** ()

Update minimal and maximal attribute limits and return tuple (min, max)

**get\_value** (*enum\_as\_str=False*)

Get current attribute value.

If `enum_as_str==True`, try to represent enums as their string values; otherwise, return their integer values (only integers can be used for setting).

**set\_value** (*value*)

Set attribute value

**class** `pylablib.devices.DCAM.DCAM.TDeviceInfo` (*vendor, model, serial\_number, camera\_version*)

Bases: `tuple`

**camera\_version**

**count** ()

Return number of occurrences of value.

**index** ()

Return first index of value.

Raises `ValueError` if the value is not present.

**model**

**serial\_number**

**vendor**

**class** `pylablib.devices.DCAM.DCAM.TFrameInfo` (*frame\_index, framestamp, timestamp\_us, camerastamp, position, pixeltype*)

Bases: `tuple`

**camerastamp**

**count** ()

Return number of occurrences of value.

**frame\_index**

**framestamp**

**index** ()

Return first index of value.

Raises `ValueError` if the value is not present.

**pixeltype**

**position**

**timestamp\_us**

**class** `pylablib.devices.DCAM.DCAM.DCAMCamera` (*idx=0*)

Bases: `pylablib.devices.interface.camera.IBinROICamera`, `pylablib.devices.interface.camera.IExposureCamera`, `pylablib.devices.interface.camera.IAttributeCamera`

**Error** = `<Mock name='mock.DCAMError' id='140096283900304'>`

**TimeoutError** = `<Mock spec='str' id='140096270275280'>`

**open** ()

Open connection to the camera

**close** ()

Close connection to the camera

**is\_opened()**  
Check if the device is connected

**get\_device\_info()**  
Get camera model data.  
  
Return tuple (vendor, model, serial\_number, camera\_version).

**get\_attribute\_value** (name, enum\_as\_str=False, error\_on\_missing=True, default=None)  
Get value of an attribute with the given name.  
  
If the value doesn't exist or can not be read and error\_on\_missing==True, raise error; otherwise, return *default*. If *default* is not None, assume that error\_on\_missing==False. If enum\_as\_str==True, try to represent enums as their string values; otherwise, return their integer values (only integers can be used for setting).

**set\_attribute\_value** (name, value, error\_on\_missing=True)  
Set value of an attribute with the given name.  
  
If the value doesn't exist or can not be written and error\_on\_missing==True, raise error; otherwise, do nothing.

**get\_all\_attribute\_values** (root="", enum\_as\_str=False)  
Get values of all attributes.  
  
If enum\_as\_str==True, try to represent enums as their string values; otherwise, return their integer values (only integers can be used for setting).

**set\_all\_attribute\_values** (settings)  
Set values of all attribute in the given dictionary

**set\_trigger\_mode** (mode)  
Set trigger mode.  
  
Can be "int" (internal), "ext" (external), or "software" (software trigger).

**get\_trigger\_mode** ()  
Get trigger mode.  
  
Can be "int" (internal), "ext" (external), or "software" (software trigger).

**get\_all\_trigger\_modes** ()  
Return the list of all available trigger modes

**setup\_ext\_trigger** (invert=False, delay=0.0)  
Setup external trigger (inversion and delay)

**get\_ext\_trigger\_parameters** ()  
Return external trigger parameters (inversion and delay)

**send\_software\_trigger** ()  
Send software trigger signal

**set\_exposure** (exposure)  
Set camera exposure

**get\_exposure** ()  
Set current exposure

**set\_readout\_speed** (speed='fast')  
Set readout speed (can be "fast" or "slow")

**get\_readout\_speed** ()  
Set current readout speed

**get\_all\_readout\_speeds()**

Return the list of all available readout speeds

**get\_frame\_readout\_time()**

Set current frame readout time

**get\_frame\_timings()**

Get acquisition timing.

Return tuple (exposure, frame\_period).

**get\_defect\_correct\_mode()**

Check if the defect pixel correction mode is on

**set\_defect\_correct\_mode(enabled=True)**

Enable or disable the defect pixel correction mode

**get\_detector\_size()**

Get camera detector size (in pixels) as a tuple (width, height)

**get\_roi()**

Get current ROI.

Return tuple (hstart, hend, vstart, vend, hbin, vbin).

**set\_roi(hstart=0, hend=None, vstart=0, vend=None, hbin=1, vbin=1)**

Set current ROI.

By default, all non-supplied parameters take extreme values. Binning is the same for both axes, so value of *vbin* is ignored (it is left for compatibility).

**get\_roi\_limits(hbin=1, vbin=1)**

Get the minimal and maximal ROI parameters.

Return tuple (hlim, vlim), where each element is in turn a limit 5-tuple (min, max, pstep, sstep, maxbin) with, correspondingly, minimal and maximal size, position and size step, and the maximal binning. In some cameras, the step and the minimal size depend on the binning, which can be supplied.

**setup\_acquisition(mode='sequence', nframes=100)**

Setup acquisition.

*mode* can be either "snap" (single frame or a fixed number of frames) or "sequence" (continuous acquisition). *nframes* determines number of frames to acquire in the single mode, or size of the ring buffer in the "sequence" mode (by default, 100).

**clear\_acquisition()**

Clear acquisition settings

**start\_acquisition(\*args, \*\*kwargs)**

Start acquisition.

Can take the same keyword parameters as :meth:“*setup\_acquisition*”. If the acquisition is not set up yet, set it up using the supplied parameters (use default of *setup\_acquisition()*, if the parameter is None). Otherwise, if any supplied parameters are different from the current ones, change them and reset the acquisition.

**stop\_acquisition()**

Stop acquisition

**get\_status()**

Get acquisition status.

Can be "busy" (capturing in progress), "ready" (ready for capturing), "stable" (not prepared for capturing), "unstable" (can't be prepared for capturing), or "error" (some other error).

**acquisition\_in\_progress()**

Check if acquisition is in progress

**get\_transfer\_info()**

Get frame transfer info.

Return tuple (last\_buff, frame\_count), where last\_buff is the index of the last filled buffer, and frame\_count is the total number of acquired frames.

**FrameTransferError**

alias of `pylablib.devices.interface.camera.DefaultFrameTransferError`

**class NoParameterCaller(device, kind)**

Bases: `object`

Class to simplify calling functions without a parameter

**apply\_settings(settings)**

Apply the settings.

settings is the dict {name: value} of the device available settings. Non-applicable settings are ignored.

**get\_acquisition\_parameters()**

Get acquisition parameters.

Return dictionary {name: value}

**get\_all\_attributes(copy=False)**

Return a dictionary of all available attributes.

If copy==True, copy the dictionary; otherwise, return the internal dictionary structure (should not be modified).

**get\_attribute(name, error\_on\_missing=True)**

Get the camera attribute with the given name

**get\_data\_dimensions()**

Get readout data dimensions (in pixels) as a tuple (width, height); take indexing mode into account

**get\_device\_variable(key)**

Get the value of a settings, status, or full info parameter

**get\_frame\_format()**

Get format for the returned images.

Can be "list" (list of 2D arrays), "array" (a single 3D array), or "chunks" (list of 3D "chunk" arrays; supported for some cameras and provides the best performance).

**get\_frame\_info\_fields()**

Get the names of frame info fields.

Applicable when frame info format (set by `set_frame_info_format()`) is "list" or "array".

**get\_frame\_info\_format()**

Get format of the frame info.

Can be "namedtuple" (potentially nested named tuples; convenient to get particular values), "list" (flat list of values, with field names are given by `get_frame_info_fields()`; convenient for building a table), "array" (same as "list", but with a numpy array, which is easier to use for "chunks"

frame format), or "dict" (flat dictionary with the same fields as the "list" format; more resilient to future format changes)

**get\_frame\_info\_period()**

Get period of frame info acquisition.

Frame info might be skipped (set to `None`) except for frames which indices are divisible by *period*. Useful for certain cameras where acquiring frame info takes a lot of time and can reduce performance at higher frame rates. Note that this parameter can still be ignored (i.e., always set to 1) if the performance is not an issue for a given camera class.

**get\_frame\_period()**

Get frame period (time between two consecutive frames in the internal trigger mode)

**get\_frames\_status()**

Get acquisition and buffer status.

Return tuple (*acquired*, *unread*, *skipped*, *size*), where *acquired* is the total number of acquired frames, *unread* is the number of acquired but not read frames, *skipped* is the number of skipped (not read and then written over) frames, and *buffer\_size* is the total buffer size (in frames).

**get\_full\_info(include=0)**

Get dict {name: value} containing full device information (including status and settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include*=-10 queries all available variables, which is equivalent to *include*="all".

**get\_full\_status(include=0)**

Get dict {name: value} containing the device status (including settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include*=-10 queries all available variables, which is equivalent to *include*="all".

**get\_image\_indexing()**

Get indexing for the returned images.

Can be "rct" (first index row, second index column, rows counted from the top), "rcb" (same as "rc", rows counted from the bottom), "xyt" (first index column, second index row, rows counted from the top), or "xyb" (same as "xyt", rows counted from the bottom)

**get\_new\_images\_range()**

Get the range of the new images.

Return tuple (*first*, *last*) with images range (first inclusive). If no images are available, return `None`. If some images were in the buffer were overwritten, exclude them from the range.

**get\_settings(include=0)**

Get dict {name: value} containing all the device settings.

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include*=-10 queries all available variables, which is equivalent to *include*="all".

**grab(nframes=1, frame\_timeout=5.0, missing\_frame='skip', return\_info=False, buff\_size=None)**

Snap *nframes* images (with preset image read mode parameters)

*buff\_size* determines buffer size (if `None`, use the default size). Timeout is specified for a single-frame acquisition, not for the whole acquisition time. *missing\_frame* determines what to do with frames which have been lost: can be "none" (replacing them with `None`), "zero" (replacing them with zero-filled frame), or "skip" (skipping them, while still keeping total returned frames number to *n*). If `return_info==True`, return tuple (*frames*, *infos*), where *infos* is a list of frame info tuples (camera-dependent); if some frames are missing and `missing_frame!="skip"`, the corresponding frame info is `None`.

**is\_acquisition\_setup()**

Check if acquisition is set up.

If the camera does not support separate acquisition setup, always return `True`.

**pausing\_acquisition**(*clear=None*, *stop=True*, *setup\_after=None*, *start\_after=True*, *combine\_nested=True*)

Context manager which temporarily pauses acquisition during execution of `with` block.

Useful for applying certain settings which can't be changed during the acquisition. If `clear==True`, clear acquisition in addition to stopping (by default, use the class default specified as `_clear_pausing_acquisition` attribute). If `stop==True`, stop the acquisition (if `clear==True`, stop regardless). If `setup_after==True`, setup the acquisition after pause if necessary (`None` means setup only if clearing was required). If `start_after==True`, start the acquisition after pause if necessary (`None` means start only if stopping was required). If `combine_nested==True`, then any nested `pausing_acquisition` calls will stop/clear acquisition as necessary, but won't setup/start it again until this `pausing_acquisition` call is complete.

Yields tuple (*acq\_in\_progress*, *acq\_params*), which indicates whether acquisition is currently in progress, and what are the current acquisition parameters.

**read\_newest\_image**(*peek=False*, *return\_info=False*)

Read the newest un-read image.

If no un-read frames are available, return `None`. If `peek==True`, return the image but not mark it as read. If `return_info==True`, return tuple (*frame*, *info*), where *info* is an info tuples (camera-dependent, see `read_multiple_images()`).

**read\_oldest\_image**(*peek=False*, *return\_info=False*)

Read the oldest un-read image.

If no un-read frames are available, return `None`. If `peek==True`, return the image but not mark it as read. If `return_info==True`, return tuple (*frame*, *info*), where *info* is an info tuples (camera-dependent, see `read_multiple_images()`).

**set\_device\_variable**(*key*, *value*)

Set the value of a settings parameter

**set\_frame\_format**(*fnt*)

Set format for the returned images.

Can be "list" (list of 2D arrays), "array" (a single 3D array), "chunks" (list of 3D "chunk" arrays; supported for some cameras and provides the best performance), or "try\_chunks" (same as "chunks", but if chunks are not supported, set to "list" instead). If format is "chunks" and chunks are not supported by the camera, it results in one frame per chunk. Note that if the format is set to "array" or "chunks", the frame info format is also automatically set to "array". If the format is set to "chunks", then the image info is also returned in chunks form (list of 2D info arrays with the same length as the corresponding frame chunks).

**set\_frame\_info\_format**(*fnt*, *include\_fields=None*)

Set format of the frame info.

Can be "namedtuple" (potentially nested named tuples; convenient to get particular values), "list" (flat list of values, with field names are given by `get_frame_info_fields()`; convenient for building a table), "array" (same as "list", but with a numpy array, which is easier to use for "chunks" frame format), or "dict" (flat dictionary with the same fields as the "list" format; more resilient to future format changes) If `include_fields` is not None, it specifies the fields included for non-"tuple" formats; note that order or `include_fields` is ignored, and the resulting fields are always ordered same as in the original.

**set\_frame\_info\_period** (*period=1*)

Set period of frame info acquisition.

Frame info might be skipped (set to None) except for frames which indices are divisible by *period*. Useful for certain cameras where acquiring frame info takes a lot of time and can reduce performance at higher frame rates. Note that this parameter can still be ignored (i.e., always set to 1) if the performance is not an issue for a given camera class.

**set\_image\_indexing** (*indexing*)

Set up indexing for the returned images.

Can be "rct" (first index row, second index column, rows counted from the top), "rcb" (same as "rc", rows counted from the bottom), "xyt" (first index column, second index row, rows counted from the top), or "xyb" (same as "xyt", rows counted from the bottom)

**snap** (*timeout=5.0, return\_info=False*)

Snap a single frame

**wait\_for\_frame** (*since='lastread', nframes=1, timeout=20.0, error\_on\_stopped=False*)

Wait for one or several new camera frames.

*since* specifies the reference point for waiting to acquire *nframes* frames; can be "lastread" (from the last read frame), "lastwait" (wait for the last successful `wait_for_frame()` call), "now" (from the start of the current call), or "start" (from the acquisition start, i.e., wait until *nframes* frames have been acquired). *timeout* can be either a number, None (infinite timeout), or a tuple (*timeout*, *frame\_timeout*), in which case the call times out if the total time exceeds *timeout*, or a single frame wait exceeds *frame\_timeout*. If the call times out, raise `TimeoutError`. If `error_on_stopped==True` and the acquisition is not running, raise `Error`; otherwise, simply return False without waiting.

**read\_multiple\_images** (*rng=None, peek=False, missing\_frame='skip', return\_info=False, return\_rng=False*)

Read multiple images specified by *rng* (by default, all un-read images).

If *rng* is specified, it is a tuple (*first*, *last*) with images range (first inclusive). If no new frames are available, return an empty list; if no acquisition is running, return None. If `peek==True`, return images but not mark them as read. *missing\_frame* determines what to do with frames which are out of range (missing or lost): can be "none" (replacing them with None), "zero" (replacing them with zero-filled frame), or "skip" (skipping them). If `return_info==True`, return tuple (*frames*, *infos*), where *infos* is a list of `TFrameInfo` instances describing frame index, framestamp and timestamp, camera stamp, frame location on the sensor, and pixel type; if some frames are missing and `missing_frame!="skip"`, the corresponding frame info is None. if `return_rng==True`, return the range covered resulting frames; if `missing_frame=="skip"`, the range can be smaller than the supplied *rng* if some frames are skipped.

## Module contents

### pylablib.devices.HighFinesse package



## Submodules

### pylablib.devices.HighFinesse.wlm module

`pylablib.devices.HighFinesse.wlm.muxchannel(*args, **kwargs)`

Multiplex the function over its channel argument

**class** `pylablib.devices.HighFinesse.wlm.TDeviceInfo(model, serial_number, revision_number, compilation_number)`

Bases: `tuple`

**compilation\_number**

**count()**

Return number of occurrences of value.

**index()**

Return first index of value.

Raises `ValueError` if the value is not present.

**model**

**revision\_number**

**serial\_number**

**class** `pylablib.devices.HighFinesse.wlm.WLM(version=None, dll_path=None, app_path=None, autostart=True)`

Bases: `pylablib.core.devio.interface.IDevice`

Generic HighFinesse wavemeter.

#### Parameters

- **version** (*int*) – wavemeter version; if `None`, use any available version
- **dll\_path** – path to `wlmData.dll`; if `None`, use standard locations or search based on the version
- **app\_path** – path to the wavemeter server application (looks like `wlm_ws.exe` or `wlm_ws7.exe`); if `None`, try to autodetect, or rely on the server already running
- **autostart** – if `True`, start measurements automatically (if the wavemeter server app is not running, it will launch with the measurements stopped).

**Error** = `<Mock name='mock.HighFinesseError' id='140096293235792'>`

**open()**

Open the connection to the wavemeter

**close()**

Close the connection to the wavemeter

**is\_opened()**

Check if the device is connected

**get\_device\_info()**

Get the wavemeter info.

Return tuple (model, serial\_number, revision\_number, compilation\_number).

**start\_measurement()**

Start wavemeter measurement

**stop\_measurement()**

Stop wavemeter measurement

**is\_measurement\_running()**

Check if the measurement is running

**set\_read\_mode(mode)**

Set value read mode, which applies to *get\_frequency()* and *get\_wavelength()*.

Can be "latest" (always return the latest measurement result; default), or "single" (if there's no new measurement since the last call, the result is "noval" which, depending on the arguments, causes wait, is returned as is, or raises an error).

**get\_read\_mode()**

Get value read mode, which applies to *get\_frequency()* and *get\_wavelength()*.

Can be "latest" (always return the latest measurement result; default), or "single" (if there's no new measurement since the last call, the result is "noval" which, depending on the arguments, causes wait, is returned as is, or raises an error).

**get\_channels\_number(refresh=True)**

Get number of channels in the wavemeter

**get\_default\_channel()**

Get the default channel (starting from 1) which is used for querying

**set\_default\_channel(channel)**

Set the default channel (starting from 1) which is used for querying

**get\_frequency(channel=None, error\_on\_invalid=True, wait=True, timeout=5.0)**

Get the wavemeter readings (in Hz) on a given channel.

*channel* is the measurement channel (starting from 1); if *None*, use the default channel. If *error\_on\_invalid==True*, raise an error if the measurement is invalid (e.g., over- or underexposure); otherwise, the method can return "under" if the meter is underexposed or "over" if it is overexposed, "badsig" if there is no calculable signal, "noval" if there are no values acquired yet, or "nosig" if there is no signal. If *wait==True* and the result is "noval" (e.g., if the read mode is "single" and no new value has been acquired since the last call), wait for at most *timeout* until a new value appears; if the *timeout* has passed, use the default behavior (error or "noval" result).

**get\_wavelength(channel=None, error\_on\_invalid=True, wait=True, timeout=5.0)**

Get the wavemeter readings (in m, and in vacuum).

*channel* is the measurement channel (starting from 1); if *None*, use the default channel. If *error\_on\_invalid==True*, raise an error if the measurement is invalid (e.g., over- or underexposure); otherwise, the method can return "under" if the meter is underexposed or "over" if it is overexposed, "badsig" if there is no calculable signal, "noval" if there are no values acquired yet, or "nosig" if there is no signal. If *wait==True* and the result is "noval" (e.g., if the read mode is "single" and no new value has been acquired since the last call), wait for at most *timeout* until a new value appears; if the *timeout* has passed, use the default behavior (error or "noval" result).

**get\_exposure\_mode(channel=None)**

Get the exposure mode ("manual" or "auto") at the given channel

**set\_exposure\_mode(mode='auto', channel=None)**

Set the exposure mode ("manual" or "auto") at the given channel

**get\_exposure(sensor=1, channel=None)**

Get the exposure for a given channel and sensor (starting from 1)

**set\_exposure(exposure, sensor=1, channel=None)**

Manually set the exposure for a given channel and sensor (starting from 1)

**get\_switcher\_mode()**  
Get the switcher mode ("off" for manual switching or "on" for cycling mode)

**set\_switcher\_mode(mode='on')**  
Set the switcher mode ("off" for manual switching or "on" for cycling mode)

**get\_active\_channel()**  
Get the current active channel

**set\_active\_channel(channel, automode=True)**  
Set the current switcher channel.  
  
Only makes sense in the manual ("off") switcher mode. If automode==True, switch to this mode automatically.

**is\_switcher\_channel\_enabled(channel, automode=True)**  
Check whether the switcher channel enabled.  
  
Only works in the cycling ("on") switcher mode. If automode==True, switch to this mode automatically.

**is\_switcher\_channel\_shown(channel, automode=True)**  
Check whether the switcher channel is shown in the wavemeter control application.  
  
Only works in the cycling ("on") switcher mode. If automode==True, switch to this mode automatically.

**enable\_switcher\_channel(channel, enable=True, show=None, automode=True)**  
Enable or disable the current switcher channel in the switch mode.  
  
Only works in the cycling ("on") switcher mode. If automode==True, switch to this mode automatically.

**get\_pulse\_mode()**  
Get the current pulse mode.  
  
Can be "cw" (CW laser mode), "int" (standard single-laser internally triggered mode), "ext" (single- or double-laser mode with external TTL trigger), or "opt" (double-laser mode with optical triggering).

**set\_pulse\_mode(mode)**  
Set the current pulse mode.  
  
Can be "cw" (CW laser mode), "int" (standard single-laser internally triggered mode), "ext" (single- or double-laser mode with external TTL trigger), or "opt" (double-laser mode with optical triggering).

**get\_precision\_mode()**  
Set the current precision mode ("fine", "wide", or "grating")

**set\_precision\_mode(mode)**  
Set the current precision mode ("fine", "wide", or "grating")

**get\_measurement\_interval()**  
Set measurement interval (per channel), or None if the interval mode is off

**set\_measurement\_interval(interval=None)**  
Set measurement interval (per channel).  
  
None means that the interval mode is off.

**calibrate(source\_type, source\_frequency, channel=None)**  
Initialize the calibration.  
  
*source\_type* is the calibration source type, which can be "hene\_633" (HeNe 633nm laser), "hene\_1152" (HeNe 1152nm laser), "hene\_free" (free-running HeNe laser), "nel" (Ne lamp),

or "other" (other source). *source\_frequency* is the exact source frequency (in Hz) sent through the given *channel*.

**get\_autocalibration\_parameters()**

Get up the automatic calibration parameters.

Return tuple (*enable*, *unit*, *period*), where *enable* determines if it is enabled, and *unit* and *period* together specify the calibration period. *unit* can be "start" (once on the measurement start; *period* is irrelevant here), "meas" (once every *period* frequency measurements), "min" (once every *period* minutes), "hours", or "days".

**setup\_autocalibration** (*enable=True*, *unit=None*, *period=None*)

Set up the automatic calibration parameters.

*enable* determines if it is enabled. *unit* and *period* together specify the calibration period. *unit* can be "start" (once on the measurement start; *period* is irrelevant here), "meas" (once every *period* frequency measurements), "min" (once every *period* minutes), "hours", or "days". Any None parameters are kept at the present value.

**class NoParameterCaller** (*device*, *kind*)

Bases: `object`

Class to simplify calling functions without a parameter

**apply\_settings** (*settings*)

Apply the settings.

*settings* is the dict {*name*: *value*} of the device available settings. Non-applicable settings are ignored.

**get\_device\_variable** (*key*)

Get the value of a settings, status, or full info parameter

**get\_full\_info** (*include=0*)

Get dict {*name*: *value*} containing full device information (including status and settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_full\_status** (*include=0*)

Get dict {*name*: *value*} containing the device status (including settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_settings** (*include=0*)

Get dict {*name*: *value*} containing all the device settings.

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**set\_device\_variable** (*key*, *value*)

Set the value of a settings parameter

## Module contents

### pylablib.devices.IMAQ package

#### Submodules

#### pylablib.devices.IMAQ.IMAQ module

`pylablib.devices.IMAQ.IMAQ.list_cameras()`

List all cameras available through IMAQ interface

`pylablib.devices.IMAQ.IMAQ.get_cameras_number()`

Get number of connected IMAQ cameras

**class** `pylablib.devices.IMAQ.IMAQ.TDeviceInfo(serial_number, interface)`

Bases: `tuple`

**count()**

Return number of occurrences of value.

**index()**

Return first index of value.

Raises `ValueError` if the value is not present.

**interface**

**serial\_number**

**class** `pylablib.devices.IMAQ.IMAQ.IMAQFrameGrabber(imaq_name='img0', do_open=True, **kwargs)`

Bases: `pylablib.devices.interface.camera.IROICamera`

Generic IMAQ frame grabber interface.

Compared to `IMAQCamera`, has more permissive initialization arguments, which simplifies its use as a base class for expanded cameras.

#### Parameters

- **imaq\_name** – interface name (can be learned by `list_cameras()`; usually, but not always, starts with "cam" or "img")
- **do\_open** – if `False`, skip the last step of opening the device (should be opened in a subclass)

**Error** = `<Mock name='mock.IMAQError' id='140096294147280'>`

**TimeoutError** = `<Mock spec='str' id='140096274107088'>`

**open()**

Open connection to the camera

**close()**

Close connection to the camera

**reset()**

Reset connection to the camera

**is\_opened()**

Check if the device is connected

**get\_grabber\_attribute\_value** (*attr*, *default=None*, *kind='auto'*)

Get value of an attribute with a given name or index.

If *default* is not *None*, return *default* if the attribute is not supported; otherwise, raise an error. *kind* is the attribute kind, and it can be "uint32", "uint64", "double", or "auto" (autodetect based on the stored list of attribute kinds).

**set\_grabber\_attribute\_value** (*attr*, *value*, *kind='int32'*)

Set value of an attribute with a given name or index.

*kind* is the attribute kind, and it can be "uint32", "uint64", "double", or "auto" (autodetect based on the stored list of attribute kinds).

**get\_all\_grabber\_attribute\_values** ()

Get a dictionary of all readable attributes.

The attributes types are autodetected, and some of the types of uncommon attributes may be misrepresented.

**get\_device\_info** ()

Get camera model data.

Return tuple (*serial*, *interface*) with the board serial number and an the interface type (e.g., "1430" for NI PCIe-1430)

**get\_detector\_size** ()

Get camera detector size (in pixels) as a tuple (*width*, *height*)

**get\_grabber\_detector\_size** ()

Get camera detector size (in pixels) as a tuple (*width*, *height*)

**get\_roi** ()

Get current ROI.

Return tuple (*hstart*, *hend*, *vstart*, *vend*). *hstart* and *hend* specify horizontal image extent, *vstart* and *vend* specify vertical image extent (start is inclusive, stop is exclusive, starting from 0).

**get\_grabber\_roi** ()

Get current ROI.

Return tuple (*hstart*, *hend*, *vstart*, *vend*). *hstart* and *hend* specify horizontal image extent, *vstart* and *vend* specify vertical image extent (start is inclusive, stop is exclusive, starting from 0).

**set\_roi** (*hstart=0*, *hend=None*, *vstart=0*, *vend=None*)

Setup camera ROI.

*hstart* and *hend* specify horizontal image extent, *vstart* and *vend* specify vertical image extent (start is inclusive, stop is exclusive, starting from 0). By default, all non-supplied parameters take extreme values (0 for start, maximal for end).

**set\_grabber\_roi** (*hstart=0*, *hend=None*, *vstart=0*, *vend=None*)

Setup camera ROI.

*hstart* and *hend* specify horizontal image extent, *vstart* and *vend* specify vertical image extent (start is inclusive, stop is exclusive, starting from 0). By default, all non-supplied parameters take extreme values (0 for start, maximal for end).

**get\_roi\_limits** (*hbin=1*, *vbin=1*)

Get the minimal and maximal ROI parameters.

Return tuple (*hlim*, *vlim*), where each element is in turn a limit 5-tuple (*min*, *max*, *pstep*, *sstep*, *maxbin*) with, correspondingly, minimal and maximal size, position and size step, and the

maximal binning (fixed to 1 if not binning is allowed). In some cameras, the step and the minimal size depend on the binning, which can be supplied.

**get\_grabber\_roi\_limits** (*hbin=1, vbin=1*)

Get the minimal and maximal ROI parameters.

Return tuple (*hlim, vlim*), where each element is in turn a limit 5-tuple (*min, max, pstep, sstep, maxbin*) with, correspondingly, minimal and maximal size, position and size step, and the maximal binning (fixed to 1 if not binning is allowed). In some cameras, the step and the minimal size depend on the binning, which can be supplied.

**configure\_trigger\_in** (*trig\_type, trig\_line=0, trig\_pol='high', trig\_action='none', time-out=None, reset\_acquisition=True*)

Configure input trigger.

#### Parameters

- **trig\_type** (*str*) – trigger source type; can be "ext", "rtsi", "iso\_in", or "software"
- **trig\_line** (*int*) – trigger line number
- **trig\_pol** (*str*) – trigger polarity; can be "high" or "low"
- **trig\_action** (*str*) – trigger action; can be "none" (disable trigger), "capture" (start capturing), "stop" (stop capturing), "buffer" (capture a single frame), or "bufflist" (capture the whole buffer list once)
- **timeout** (*float*) – timeout in seconds; None means not timeout.
- **reset\_acquisition** (*bool*) – if the input triggers configuration has been changed, acquisition needs to be restart; if True, perform it automatically

**send\_software\_trigger** ()

Send software trigger signal

**configure\_trigger\_out** (*trig\_type, trig\_line=0, trig\_pol='high', trig\_drive='disable'*)

Configure trigger output.

#### Parameters

- **trig\_type** (*str*) – trigger drive destination type; can be "ext", "rtsi", or "iso\_out"
- **trig\_line** (*int*) – trigger line number
- **trig\_pol** (*str*) – trigger polarity; can be "high" or "low"
- **trig\_drive** (*str*) – trigger output signal; can be "disable" (disable drive), "acq\_in\_progress" (asserted when acquisition is started), "acq\_done" (asserted when acquisition is done), "unasserted" (force unasserted level), "asserted" (force asserted level), "hsync" (asserted on start of a single line start), "vsync" (asserted on start of a frame scan), "frame\_start" (asserted when a single frame is captured), or "frame\_done" (asserted when a single frame is done)

**read\_trigger** (*trig\_type, trig\_line=0, trig\_pol='high'*)

Read current value of a trigger (input or output).

#### Parameters

- **trig\_type** (*str*) – trigger drive destination type; can be "ext", "rtsi", "iso\_in", or "iso\_out"
- **trig\_line** (*int*) – trigger line number

- **trig\_pol** (*str*) – trigger polarity; can be "high" or "low"

**clear\_all\_triggers** (*reset\_acquisition=True*)

Disable all triggers of the session

If the input triggers configuration has been changed, acquisition needs to be restart; if `reset_acquisition==True`, perform it automatically.

**setup\_serial\_params** (*write\_term="", datatype='bytes'*)

Setup default serial communication parameters.

#### Parameters

- **write\_term** – default terminator character to be added to the sent messages
- **datatype** – type of the result of read commands; can be "bytes" (return raw bytes), or "str" (convert into UTF-8 string)

**get\_serial\_params** ()

Return serial parameters as a tuple (*write\_term*, *datatype*)

**serial\_write** (*msg, timeout=3.0, term=None*)

Write message into CameraLink serial port.

#### Parameters

- **timeout** – operation timeout (in seconds)
- **term** – additional write terminator character to add to the message; if `None`, use the value set up using `setup_serial_params()` (by default, no additional terminator)

**serial\_read** (*n, timeout=3.0, datatype=None*)

Read specified number of bytes from CameraLink serial port.

#### Parameters

- **n** – number of bytes to read
- **timeout** – operation timeout (in seconds)
- **datatype** – return datatype; can be "bytes" (return raw bytes), or "str" (convert into UTF-8 string) if `None`, use the value set up using `setup_serial_params()` (by default, "bytes")

**serial\_readline** (*timeout=3.0, datatype=None, maxn=1024*)

Read bytes from CameraLink serial port until the termination character (defined in camera file) is encountered.

#### Parameters

- **timeout** – operation timeout (in seconds)
- **datatype** – return datatype; can be "bytes" (return raw bytes), or "str" (convert into UTF-8 string) if `None`, use the value set up using `setup_serial_params()` (by default, "bytes")
- **maxn** – maximal number of bytes to read

**serial\_flush** ()

Flush CameraLink serial port

**setup\_acquisition** (*mode='sequence', nframes=100*)

Setup acquisition mode.



*mode* can be either "snap" (single frame or a fixed number of frames) or "sequence" (continuous acquisition). (note that `IMAQCamera.acquisition_in_progress()` would still return True in this case, even though new frames are no longer acquired). *nframes* sets up number of frame buffers.

**clear\_acquisition()**

Clear all acquisition details and free all buffers

**start\_acquisition(\*args, \*\*kwargs)**

Start acquisition.

Can take the same keyword parameters as `:meth:‘‘setup_acquisition`. If the acquisition is not set up yet, set it up using the supplied parameters (use default of `setup_acquisition()`, if the parameter is None). Otherwise, if any supplied parameters are different from the current ones, change them and reset the acquisition.

**stop\_acquisition()**

Stop acquisition

**acquisition\_in\_progress()**

Check if acquisition is in progress

**FrameTransferError**

alias of `pylablib.devices.interface.camera.DefaultFrameTransferError`

**class NoParameterCaller(device, kind)**

Bases: `object`

Class to simplify calling functions without a parameter

**apply\_settings(settings)**

Apply the settings.

*settings* is the dict {name: value} of the device available settings. Non-applicable settings are ignored.

**get\_acquisition\_parameters()**

Get acquisition parameters.

Return dictionary {name: value}

**get\_data\_dimensions()**

Get readout data dimensions (in pixels) as a tuple (width, height); take indexing mode into account

**get\_device\_variable(key)**

Get the value of a settings, status, or full info parameter

**get\_frame\_format()**

Get format for the returned images.

Can be "list" (list of 2D arrays), "array" (a single 3D array), or "chunks" (list of 3D "chunk" arrays; supported for some cameras and provides the best performance).

**get\_frame\_info\_fields()**

Get the names of frame info fields.

Applicable when frame info format (set by `set_frame_info_format()`) is "list" or "array".

**get\_frame\_info\_format()**

Get format of the frame info.

Can be "namedtuple" (potentially nested named tuples; convenient to get particular values), "list" (flat list of values, with field names are given by `get_frame_info_fields()`; convenient for building a table), "array" (same as "list", but with a numpy array, which is easier to use for "chunks"

frame format), or "dict" (flat dictionary with the same fields as the "list" format; more resilient to future format changes)

**get\_frame\_info\_period()**

Get period of frame info acquisition.

Frame info might be skipped (set to `None`) except for frames which indices are divisible by *period*. Useful for certain cameras where acquiring frame info takes a lot of time and can reduce performance at higher frame rates. Note that this parameter can still be ignored (i.e., always set to 1) if the performance is not an issue for a given camera class.

**get\_frames\_status()**

Get acquisition and buffer status.

Return tuple (*acquired*, *unread*, *skipped*, *size*), where *acquired* is the total number of acquired frames, *unread* is the number of acquired but not read frames, *skipped* is the number of skipped (not read and then written over) frames, and *buffer\_size* is the total buffer size (in frames).

**get\_full\_info(include=0)**

Get dict {*name*: *value*} containing full device information (including status and settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include*=-10 queries all available variables, which is equivalent to *include*="all".

**get\_full\_status(include=0)**

Get dict {*name*: *value*} containing the device status (including settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include*=-10 queries all available variables, which is equivalent to *include*="all".

**get\_image\_indexing()**

Get indexing for the returned images.

Can be "rct" (first index row, second index column, rows counted from the top), "rcb" (same as "rc", rows counted from the bottom), "xyt" (first index column, second index row, rows counted from the top), or "xyb" (same as "xyt", rows counted from the bottom)

**get\_new\_images\_range()**

Get the range of the new images.

Return tuple (*first*, *last*) with images range (first inclusive). If no images are available, return `None`. If some images were in the buffer were overwritten, exclude them from the range.

**get\_settings(include=0)**

Get dict {*name*: *value*} containing all the device settings.

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include*=-10 queries all available variables, which is equivalent to *include*="all".

**grab(nframes=1, frame\_timeout=5.0, missing\_frame='skip', return\_info=False, buff\_size=None)**

Snap *nframes* images (with preset image read mode parameters)

*buff\_size* determines buffer size (if `None`, use the default size). Timeout is specified for a single-frame acquisition, not for the whole acquisition time. *missing\_frame* determines what to do with frames which have been lost: can be "none" (replacing them with `None`), "zero" (replacing them with zero-filled frame), or "skip" (skipping them, while still keeping total returned frames number to *n*). If

`return_info==True`, return tuple (frames, infos), where infos is a list of frame info tuples (camera-dependent); if some frames are missing and `missing_frame!="skip"`, the corresponding frame info is None.

**is\_acquisition\_setup()**

Check if acquisition is set up.

If the camera does not support separate acquisition setup, always return True.

**pausing\_acquisition**(clear=None, stop=True, setup\_after=None, start\_after=True, combine\_nested=True)

Context manager which temporarily pauses acquisition during execution of with block.

Useful for applying certain settings which can't be changed during the acquisition. If `clear==True`, clear acquisition in addition to stopping (by default, use the class default specified as `_clear_pausing_acquisition` attribute). If `stop==True`, stop the acquisition (if `clear==True`, stop regardless). If `setup_after==True`, setup the acquisition after pause if necessary (None means setup only if clearing was required). If `start_after==True`, start the acquisition after pause if necessary (None means start only if stopping was required). If `combine_nested==True`, then any nested `pausing_acquisition` calls will stop/clear acquisition as necessary, but won't setup/start it again until this `pausing_acquisition` call is complete.

Yields tuple (acq\_in\_progress, acq\_params), which indicates whether acquisition is currently in progress, and what are the current acquisition parameters.

**read\_multiple\_images**(rng=None, peek=False, missing\_frame='skip', return\_info=False, return\_rng=False)

Read multiple images specified by rng (by default, all un-read images).

If rng is specified, it is a tuple (first, last) with images range (first inclusive). If no new frames are available, return an empty list; if no acquisition is running, return None. If `peek==True`, return images but not mark them as read. `missing_frame` determines what to do with frames which are out of range (missing or lost): can be "none" (replacing them with None), "zero" (replacing them with zero-filled frame), or "skip" (skipping them). If `return_info==True`, return tuple (frames, infos), where infos is a list of frame info tuples (camera-dependent, by default, only the frame index); if some frames are missing and `missing_frame!="skip"`, the corresponding frame info is None. if `return_rng==True`, return the range covered resulting frames; if `missing_frame=="skip"`, the range can be smaller than the supplied rng if some frames are skipped.

**read\_newest\_image**(peek=False, return\_info=False)

Read the newest un-read image.

If no un-read frames are available, return None. If `peek==True`, return the image but not mark it as read. If `return_info==True`, return tuple (frame, info), where info is an info tuples (camera-dependent, see `read_multiple_images()`).

**read\_oldest\_image**(peek=False, return\_info=False)

Read the oldest un-read image.

If no un-read frames are available, return None. If `peek==True`, return the image but not mark it as read. If `return_info==True`, return tuple (frame, info), where info is an info tuples (camera-dependent, see `read_multiple_images()`).

**set\_device\_variable**(key, value)

Set the value of a settings parameter

**set\_frame\_format**(fmt)

Set format for the returned images.

Can be "list" (list of 2D arrays), "array" (a single 3D array), "chunks" (list of 3D "chunk" arrays; supported for some cameras and provides the best performance), or "try\_chunks" (same as

"chunks", but if chunks are not supported, set to "list" instead). If format is "chunks" and chunks are not supported by the camera, it results in one frame per chunk. Note that if the format is set to "array" or "chunks", the frame info format is also automatically set to "array". If the format is set to "chunks", then the image info is also returned in chunks form (list of 2D info arrays with the same length as the corresponding frame chunks).

**set\_frame\_info\_format** (*fmt*, *include\_fields=None*)

Set format of the frame info.

Can be "namedtuple" (potentially nested named tuples; convenient to get particular values), "list" (flat list of values, with field names are given by `get_frame_info_fields()`; convenient for building a table), "array" (same as "list", but with a numpy array, which is easier to use for "chunks" frame format), or "dict" (flat dictionary with the same fields as the "list" format; more resilient to future format changes) If *include\_fields* is not None, it specifies the fields included for non-"tuple" formats; note that order of *include\_fields* is ignored, and the resulting fields are always ordered same as in the original.

**set\_frame\_info\_period** (*period=1*)

Set period of frame info acquisition.

Frame info might be skipped (set to None) except for frames which indices are divisible by *period*. Useful for certain cameras where acquiring frame info takes a lot of time and can reduce performance at higher frame rates. Note that this parameter can still be ignored (i.e., always set to 1) if the performance is not an issue for a given camera class.

**set\_image\_indexing** (*indexing*)

Set up indexing for the returned images.

Can be "rct" (first index row, second index column, rows counted from the top), "rcb" (same as "rc", rows counted from the bottom), "xyt" (first index column, second index row, rows counted from the top), or "xyb" (same as "xyt", rows counted from the bottom)

**snap** (*timeout=5.0*, *return\_info=False*)

Snap a single frame

**wait\_for\_frame** (*since='lastread'*, *nframes=1*, *timeout=20.0*, *error\_on\_stopped=False*)

Wait for one or several new camera frames.

*since* specifies the reference point for waiting to acquire *nframes* frames; can be "lastread" (from the last read frame), "lastwait" (wait for the last successful `wait_for_frame()` call), "now" (from the start of the current call), or "start" (from the acquisition start, i.e., wait until *nframes* frames have been acquired). *timeout* can be either a number, None (infinite timeout), or a tuple (timeout, frame\_timeout), in which case the call times out if the total time exceeds timeout, or a single frame wait exceeds frame\_timeout. If the call times out, raise `TimeoutError`. If *error\_on\_stopped==True* and the acquisition is not running, raise `Error`; otherwise, simply return False without waiting.

**class** `pylablib.devices.IMAQ.IMAQ.IMAQCamera` (*name='img0'*)

Bases: `pylablib.devices.IMAQ.IMAQ.IMAQFrameGrabber`

Generic IMAQ camera interface.

**Parameters** *name* – interface name (can be learned by `list_cameras()`; usually, but not always, starts with "cam" or "img")

**Error** = `<Mock name='mock.IMAQError' id='140096294147280'>`

**FrameTransferError**

alias of `pylablib.devices.interface.camera.DefaultFrameTransferError`

```

class NoParameterCaller (device, kind)
    Bases: object

    Class to simplify calling functions without a parameter

TimeoutError = <Mock spec='str' id='140096274107088'>

acquisition_in_progress ()
    Check if acquisition is in progress

apply_settings (settings)
    Apply the settings.

    settings is the dict {name: value} of the device available settings. Non-applicable settings are ignored.

clear_acquisition ()
    Clear all acquisition details and free all buffers

clear_all_triggers (reset_acquisition=True)
    Disable all triggers of the session

    If the input triggers configuration has been changed, acquisition needs to be restart; if reset_acquisition==True, perform it automatically.

close ()
    Close connection to the camera

configure_trigger_in (trig_type, trig_line=0, trig_pol='high', trig_action='none', timeout=None, reset_acquisition=True)
    Configure input trigger.

```

#### Parameters

- **trig\_type** (*str*) – trigger source type; can be "ext", "rtsi", "iso\_in", or "software"
- **trig\_line** (*int*) – trigger line number
- **trig\_pol** (*str*) – trigger polarity; can be "high" or "low"
- **trig\_action** (*str*) – trigger action; can be "none" (disable trigger), "capture" (start capturing), "stop" (stop capturing), "buffer" (capture a single frame), or "bufflist" (capture the whole buffer list once)
- **timeout** (*float*) – timeout in seconds; None means not timeout.
- **reset\_acquisition** (*bool*) – if the input triggers configuration has been changed, acquisition needs to be restart; if True, perform it automatically

```

configure_trigger_out (trig_type, trig_line=0, trig_pol='high', trig_drive='disable')
    Configure trigger output.

```

#### Parameters

- **trig\_type** (*str*) – trigger drive destination type; can be "ext", "rtsi", or "iso\_out"
- **trig\_line** (*int*) – trigger line number
- **trig\_pol** (*str*) – trigger polarity; can be "high" or "low"
- **trig\_drive** (*str*) – trigger output signal; can be "disable" (disable drive), "acq\_in\_progress" (asserted when acquisition is started), "acq\_done" (asserted when acquisition is done), "unasserted" (force unasserted level), "asserted" (force asserted level), "hsync" (asserted on start of a single line)

start), "vsync" (asserted on start of a frame scan), "frame\_start" (asserted when a single frame is captured), or "frame\_done" (asserted when a single frame is done)

**get\_acquisition\_parameters()**

Get acquisition parameters.

Return dictionary {name: value}

**get\_all\_grabber\_attribute\_values()**

Get a dictionary of all readable attributes.

The attributes types are autodetected, and some of the types of uncommon attributes may be misrepresented.

**get\_data\_dimensions()**

Get readout data dimensions (in pixels) as a tuple (width, height); take indexing mode into account

**get\_detector\_size()**

Get camera detector size (in pixels) as a tuple (width, height)

**get\_device\_info()**

Get camera model data.

Return tuple (serial, interface) with the board serial number and an the interface type (e.g., "1430" for NI PCIe-1430)

**get\_device\_variable(key)**

Get the value of a settings, status, or full info parameter

**get\_frame\_format()**

Get format for the returned images.

Can be "list" (list of 2D arrays), "array" (a single 3D array), or "chunks" (list of 3D "chunk" arrays; supported for some cameras and provides the best performance).

**get\_frame\_info\_fields()**

Get the names of frame info fields.

Applicable when frame info format (set by `set_frame_info_format()`) is "list" or "array".

**get\_frame\_info\_format()**

Get format of the frame info.

Can be "namedtuple" (potentially nested named tuples; convenient to get particular values), "list" (flat list of values, with field names are given by `get_frame_info_fields()`; convenient for building a table), "array" (same as "list", but with a numpy array, which is easier to use for "chunks" frame format), or "dict" (flat dictionary with the same fields as the "list" format; more resilient to future format changes)

**get\_frame\_info\_period()**

Get period of frame info acquisition.

Frame info might be skipped (set to None) except for frames which indices are divisible by *period*. Useful for certain cameras where acquiring frame info takes a lot of time and can reduce performance at higher frame rates. Note that this parameter can still be ignored (i.e., always set to 1) if the performance is not an issue for a given camera class.

**get\_frames\_status()**

Get acquisition and buffer status.

Return tuple (*acquired*, *unread*, *skipped*, *size*), where *acquired* is the total number of acquired frames, *unread* is the number of acquired but not read frames, *skipped* is the number of skipped (not read and then written over) frames, and *buffer\_size* is the total buffer size (in frames).

**get\_full\_info** (*include=0*)

Get dict {*name*: *value*} containing full device information (including status and settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_full\_status** (*include=0*)

Get dict {*name*: *value*} containing the device status (including settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_grabber\_attribute\_value** (*attr*, *default=None*, *kind='auto'*)

Get value of an attribute with a given name or index.

If *default* is not *None*, return *default* if the attribute is not supported; otherwise, raise an error. *kind* is the attribute kind, and it can be "uint32", "uint64", "double", or "auto" (autodetect based on the stored list of attribute kinds).

**get\_grabber\_detector\_size** ()

Get camera detector size (in pixels) as a tuple (*width*, *height*)

**get\_grabber\_roi** ()

Get current ROI.

Return tuple (*hstart*, *hend*, *vstart*, *vend*). *hstart* and *hend* specify horizontal image extent, *vstart* and *vend* specify vertical image extent (start is inclusive, stop is exclusive, starting from 0).

**get\_grabber\_roi\_limits** (*hbin=1*, *vbin=1*)

Get the minimal and maximal ROI parameters.

Return tuple (*hlim*, *vlim*), where each element is in turn a limit 5-tuple (*min*, *max*, *pstep*, *sstep*, *maxbin*) with, correspondingly, minimal and maximal size, position and size step, and the maximal binning (fixed to 1 if not binning is allowed). In some cameras, the step and the minimal size depend on the binning, which can be supplied.

**get\_image\_indexing** ()

Get indexing for the returned images.

Can be "rct" (first index row, second index column, rows counted from the top), "rcb" (same as "rc", rows counted from the bottom), "xyt" (first index column, second index row, rows counted from the top), or "xyb" (same as "xyt", rows counted from the bottom)

**get\_new\_images\_range** ()

Get the range of the new images.

Return tuple (*first*, *last*) with images range (first inclusive). If no images are available, return *None*. If some images were in the buffer were overwritten, exclude them from the range.

**get\_roi** ()

Get current ROI.

Return tuple (*hstart*, *hend*, *vstart*, *vend*). *hstart* and *hend* specify horizontal image extent, *vstart* and *vend* specify vertical image extent (start is inclusive, stop is exclusive, starting from 0).

**get\_roi\_limits** (*hbin=1, vbin=1*)

Get the minimal and maximal ROI parameters.

Return tuple (*hlim, vlim*), where each element is in turn a limit 5-tuple (*min, max, pstep, sstep, maxbin*) with, correspondingly, minimal and maximal size, position and size step, and the maximal binning (fixed to 1 if not binning is allowed). In some cameras, the step and the minimal size depend on the binning, which can be supplied.

**get\_serial\_params** ()

Return serial parameters as a tuple (*write\_term, datatype*)

**get\_settings** (*include=0*)

Get dict {*name: value*} containing all the device settings.

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**grab** (*nframes=1, frame\_timeout=5.0, missing\_frame='skip', return\_info=False, buff\_size=None*)

Snap *nframes* images (with preset image read mode parameters)

*buff\_size* determines buffer size (if *None*, use the default size). Timeout is specified for a single-frame acquisition, not for the whole acquisition time. *missing\_frame* determines what to do with frames which have been lost: can be "none" (replacing them with *None*), "zero" (replacing them with zero-filled frame), or "skip" (skipping them, while still keeping total returned frames number to *n*). If *return\_info==True*, return tuple (*frames, infos*), where *infos* is a list of frame info tuples (camera-dependent); if some frames are missing and *missing\_frame!="skip"*, the corresponding frame info is *None*.

**is\_acquisition\_setup** ()

Check if acquisition is set up.

If the camera does not support separate acquisition setup, always return *True*.

**is\_opened** ()

Check if the device is connected

**open** ()

Open connection to the camera

**pausing\_acquisition** (*clear=None, stop=True, setup\_after=None, start\_after=True, combine\_nested=True*)

Context manager which temporarily pauses acquisition during execution of *with* block.

Useful for applying certain settings which can't be changed during the acquisition. If *clear==True*, clear acquisition in addition to stopping (by default, use the class default specified as *\_clear\_pausing\_acquisition* attribute). If *stop==True*, stop the acquisition (if *clear==True*, stop regardless). If *setup\_after==True*, setup the acquisition after pause if necessary (*None* means setup only if clearing was required). If *start\_after==True*, start the acquisition after pause if necessary (*None* means start only if stopping was required). If *combine\_nested==True*, then any nested *pausing\_acquisition* calls will stop/clear acquisition as necessary, but won't setup/start it again until this *pausing\_acquisition* call is complete.

Yields tuple (*acq\_in\_progress, acq\_params*), which indicates whether acquisition is currently in progress, and what are the current acquisition parameters.

**read\_multiple\_images** (*rng=None, peek=False, missing\_frame='skip', return\_info=False, return\_rng=False*)

Read multiple images specified by *rng* (by default, all un-read images).



If *rng* is specified, it is a tuple (*first*, *last*) with images range (first inclusive). If no new frames are available, return an empty list; if no acquisition is running, return *None*. If *peek*==True, return images but not mark them as read. *missing\_frame* determines what to do with frames which are out of range (missing or lost): can be "none" (replacing them with *None*), "zero" (replacing them with zero-filled frame), or "skip" (skipping them). If *return\_info*==True, return tuple (*frames*, *infos*), where *infos* is a list of frame info tuples (camera-dependent, by default, only the frame index); if some frames are missing and *missing\_frame*!="skip", the corresponding frame info is *None*. if *return\_rng*==True, return the range covered resulting frames; if *missing\_frame*=="skip", the range can be smaller than the supplied *rng* if some frames are skipped.

**read\_newest\_image** (*peek=False*, *return\_info=False*)

Read the newest un-read image.

If no un-read frames are available, return *None*. If *peek*==True, return the image but not mark it as read. If *return\_info*==True, return tuple (*frame*, *info*), where *info* is an info tuples (camera-dependent, see [read\\_multiple\\_images\(\)](#)).

**read\_oldest\_image** (*peek=False*, *return\_info=False*)

Read the oldest un-read image.

If no un-read frames are available, return *None*. If *peek*==True, return the image but not mark it as read. If *return\_info*==True, return tuple (*frame*, *info*), where *info* is an info tuples (camera-dependent, see [read\\_multiple\\_images\(\)](#)).

**read\_trigger** (*trig\_type*, *trig\_line=0*, *trig\_pol='high'*)

Read current value of a trigger (input or output).

#### Parameters

- **trig\_type** (*str*) – trigger drive destination type; can be "ext", "rtsi", "iso\_in", or "iso\_out"
- **trig\_line** (*int*) – trigger line number
- **trig\_pol** (*str*) – trigger polarity; can be "high" or "low"

**reset** ()

Reset connection to the camera

**send\_software\_trigger** ()

Send software trigger signal

**serial\_flush** ()

Flush CameraLink serial port

**serial\_read** (*n*, *timeout=3.0*, *datatype=None*)

Read specified number of bytes from CameraLink serial port.

#### Parameters

- **n** – number of bytes to read
- **timeout** – operation timeout (in seconds)
- **datatype** – return datatype; can be "bytes" (return raw bytes), or "str" (convert into UTF-8 string) if *None*, use the value set up using [setup\\_serial\\_params\(\)](#) (by default, "bytes")

**serial\_readline** (*timeout=3.0*, *datatype=None*, *maxn=1024*)

Read bytes from CameraLink serial port until the termination character (defined in camera file) is encountered.

#### Parameters

- **timeout** – operation timeout (in seconds)
- **datatype** – return datatype; can be "bytes" (return raw bytes), or "str" (convert into UTF-8 string) if None, use the value set up using `setup_serial_params()` (by default, "bytes")
- **maxn** – maximal number of bytes to read

**serial\_write** (*msg, timeout=3.0, term=None*)

Write message into CameraLink serial port.

#### Parameters

- **timeout** – operation timeout (in seconds)
- **term** – additional write terminator character to add to the message; if None, use the value set up using `setup_serial_params()` (by default, no additional terminator)

**set\_device\_variable** (*key, value*)

Set the value of a settings parameter

**set\_frame\_format** (*fmt*)

Set format for the returned images.

Can be "list" (list of 2D arrays), "array" (a single 3D array), "chunks" (list of 3D "chunk" arrays; supported for some cameras and provides the best performance), or "try\_chunks" (same as "chunks", but if chunks are not supported, set to "list" instead). If format is "chunks" and chunks are not supported by the camera, it results in one frame per chunk. Note that if the format is set to "array" or "chunks", the frame info format is also automatically set to "array". If the format is set to "chunks", then the image info is also returned in chunks form (list of 2D info arrays with the same length as the corresponding frame chunks).

**set\_frame\_info\_format** (*fmt, include\_fields=None*)

Set format of the frame info.

Can be "namedtuple" (potentially nested named tuples; convenient to get particular values), "list" (flat list of values, with field names are given by `get_frame_info_fields()`; convenient for building a table), "array" (same as "list", but with a numpy array, which is easier to use for "chunks" frame format), or "dict" (flat dictionary with the same fields as the "list" format; more resilient to future format changes) If `include_fields` is not None, it specifies the fields included for non-"tuple" formats; note that order or `include_fields` is ignored, and the resulting fields are always ordered same as in the original.

**set\_frame\_info\_period** (*period=1*)

Set period of frame info acquisition.

Frame info might be skipped (set to None) except for frames which indices are divisible by *period*. Useful for certain cameras where acquiring frame info takes a lot of time and can reduce performance at higher frame rates. Note that this parameter can still be ignored (i.e., always set to 1) if the performance is not an issue for a given camera class.

**set\_grabber\_attribute\_value** (*attr, value, kind='int32'*)

Set value of an attribute with a given name or index.

*kind* is the attribute kind, and it can be "uint32", "uint64", "double", or "auto" (autodetect based on the stored list of attribute kinds).

**set\_grabber\_roi** (*hstart=0, hend=None, vstart=0, vend=None*)

Setup camera ROI.

*hstart* and *hend* specify horizontal image extent, *vstart* and *vend* specify vertical image extent (start is inclusive, stop is exclusive, starting from 0). By default, all non-supplied parameters take extreme values (0 for start, maximal for end).

**set\_image\_indexing** (*indexing*)

Set up indexing for the returned images.

Can be "rct" (first index row, second index column, rows counted from the top), "rcb" (same as "rc", rows counted from the bottom), "xyt" (first index column, second index row, rows counted from the top), or "xyb" (same as "xyt", rows counted from the bottom)

**set\_roi** (*hstart=0, hend=None, vstart=0, vend=None*)

Setup camera ROI.

*hstart* and *hend* specify horizontal image extent, *vstart* and *vend* specify vertical image extent (start is inclusive, stop is exclusive, starting from 0). By default, all non-supplied parameters take extreme values (0 for start, maximal for end).

**setup\_acquisition** (*mode='sequence', nframes=100*)

Setup acquisition mode.

*mode* can be either "snap" (single frame or a fixed number of frames) or "sequence" (continuous acquisition). (note that `IMAQCamera.acquisition_in_progress()` would still return True in this case, even though new frames are no longer acquired). *nframes* sets up number of frame buffers.

**setup\_serial\_params** (*write\_term=", datatype='bytes'*)

Setup default serial communication parameters.

#### Parameters

- **write\_term** – default terminator character to be added to the sent messages
- **datatype** – type of the result of read commands; can be "bytes" (return raw bytes), or "str" (convert into UTF-8 string)

**snap** (*timeout=5.0, return\_info=False*)

Snap a single frame

**start\_acquisition** (*\*args, \*\*kwargs*)

Start acquisition.

Can take the same keyword parameters as `:meth:'setup_acquisition`. If the acquisition is not set up yet, set it up using the supplied parameters (use default of `setup_acquisition()`, if the parameter is None). Otherwise, if any supplied parameters are different from the current ones, change them and reset the acquisition.

**stop\_acquisition** ()

Stop acquisition

**wait\_for\_frame** (*since='lastread', nframes=1, timeout=20.0, error\_on\_stopped=False*)

Wait for one or several new camera frames.

*since* specifies the reference point for waiting to acquire *nframes* frames; can be "lastread" (from the last read frame), "lastwait" (wait for the last successful `wait_for_frame()` call), "now" (from the start of the current call), or "start" (from the acquisition start, i.e., wait until *nframes* frames have been acquired). *timeout* can be either a number, None (infinite timeout), or a tuple (*timeout*, *frame\_timeout*), in which case the call times out if the total time exceeds *timeout*, or a single frame wait exceeds *frame\_timeout*. If the call times out, raise `TimeoutError`. If `error_on_stopped==True` and the acquisition is not running, raise `Error`; otherwise, simply return False without waiting.

## pylablib.devices.IMAQ.niimaq\_attrtypes module

### Module contents

## pylablib.devices.IMAQdx package

### Submodules

## pylablib.devices.IMAQdx.IMAQdx module

**class** pylablib.devices.IMAQdx.IMAQdx.**TCameraInfo** (*name, type, version, flags, serial\_number, bus, vendor, model, camera\_file, attr\_url*)

Bases: `tuple`

**attr\_url**

**bus**

**camera\_file**

**count** ()

Return number of occurrences of value.

**flags**

**index** ()

Return first index of value.

Raises ValueError if the value is not present.

**model**

**name**

**serial\_number**

**type**

**vendor**

**version**

pylablib.devices.IMAQdx.IMAQdx.**list\_cameras** (*connected=True, desc=True*)

List all cameras available through IMAQdx interface

If desc==True, return complete camera descriptions; otherwise, simply return the names.

pylablib.devices.IMAQdx.IMAQdx.**get\_cameras\_number** ()

Get number of connected dx cameras

**class** pylablib.devices.IMAQdx.IMAQdx.**IMAQdxAttribute** (*sid, name*)

Bases: `object`

Object representing an IMAQdx camera parameter.

Allows to query and set values and get additional information. Usually created automatically by an *IMAQdxCamera* instance, but could be created manually.

#### Parameters

- **sid** – camera session ID

- **name** – attribute text name

**name**  
attribute name

**kind**  
attribute kind; can be "u32", "i64", "f64", "str", "enum", "bool", "command", or "blob"

**display\_name**  
attribute display name (short description name)

**tooltip**  
longer attribute description

**description**  
full attribute description (usually, same as *tooltip*)

**units**  
attribute units (if applicable)

**visibility**  
attribute visibility ("simple", "intermediate", or "advanced")

**readable**  
whether attribute is readable  
**Type** `bool`

**writable**  
whether attribute is writable  
**Type** `bool`

**min**  
minimal attribute value (if applicable)  
**Type** `float` or `int`

**max**  
maximal attribute value (if applicable)  
**Type** `float` or `int`

**inc**  
minimal attribute increment value (if applicable)  
**Type** `float` or `int`

**ivalues**  
list of possible integer values for enum attributes

**values**  
list of possible text values for enum attributes

**labels**  
dict {label: index} which shows all possible values of an enumerated attribute and their corresponding numerical values

**ilabels**  
dict {index: label} which shows labels corresponding to numerical values of an enumerated attribute

**update\_limits()**  
Update minimal and maximal attribute limits and return tuple (min, max, inc)

**truncate\_value** (*value*)

Truncate value to lie within attribute limits

**get\_value** (*enum\_as\_str=True*)

Get attribute value.

If *enum\_as\_str*==True, return enum-style values as strings; otherwise, return corresponding integer values.

**set\_value** (*value, truncate=True*)

Get attribute value.

If *truncate*==True, automatically truncate value to lie within allowed range.

```
class pylablib.devices.IMAQdx.IMAQdx.TDeviceInfo (vendor, model, serial_number,  
                                                bus_type)
```

Bases: `tuple`

**bus\_type**

**count** ()

Return number of occurrences of value.

**index** ()

Return first index of value.

Raises ValueError if the value is not present.

**model**

**serial\_number**

**vendor**

```
class pylablib.devices.IMAQdx.IMAQdx.IMAQdxCamera (name='cam0', mode='controller',  
                                                  visibility='advanced')  
Bases: pylablib.devices.interface.camera.IROICamera, pylablib.devices.  
interface.camera.IAttributeCamera
```

Generic IMAQdx camera interface.

#### Parameters

- **name** – interface name (can be learned by `list_cameras()`; usually, but not always, starts with "cam")
- **mode** – connection mode; can be "controller" (full control) or "listener" (only reading)
- **visibility** – attribute visibility when listing attributes; can be "simple", "intermediate" or "advanced" (higher mode exposes more attributes).

```
Error = <Mock name='mock.IMAQdxError' id='140096275815696'>
```

```
TimeoutError = <Mock spec='str' id='140096275815376'>
```

**open** ()

Open connection to the camera

**close** ()

Close connection to the camera

**reset** ()

Reset connection to the camera

**is\_opened()**

Check if the device is connected

**post\_open()**

Additional setup after camera opening

**get\_attribute\_value** (*name*, *error\_on\_missing=True*, *default=None*, *enum\_as\_str=True*)

Get value of an attribute with the given name.

If the value doesn't exist or can not be read and *error\_on\_missing==True*, raise error; otherwise, return *default*. If *default* is not *None*, assume that *error\_on\_missing==False*. If *name* points at a dictionary branch, return a dictionary with all values in this branch. If *enum\_as\_str==True*, return enum-style values as strings; otherwise, return corresponding integer values.

**set\_attribute\_value** (*name*, *value*, *truncate=True*, *error\_on\_missing=True*)

Set value of an attribute with the given name.

If the value doesn't exist or can not be written and *error\_on\_missing==True*, raise error; otherwise, do nothing. If *name* points at a dictionary branch, set all values in this branch (in this case *value* must be a dictionary). If *truncate==True*, truncate value to lie within attribute range.

**get\_all\_attribute\_values** (*root="*, *enum\_as\_str=True*, *ignore\_errors=True*)

Get values of all attributes with the given *root*

**set\_all\_attribute\_values** (*settings*, *root="*, *truncate=True*)

Set values of all attributes with the given *root*.

If *truncate==True*, truncate value to lie within attribute range.

**get\_device\_info()**

Get camera information.

Return tuple (*vendor*, *model*, *serial\_number*, *bus\_type*).

**get\_detector\_size()**

Get camera detector size (in pixels) as a tuple (*width*, *height*)

**get\_roi()**

Get current ROI.

Return tuple (*hstart*, *hend*, *vstart*, *vend*). *hstart* and *hend* specify horizontal image extent, *vstart* and *vend* specify vertical image extent (start is inclusive, stop is exclusive, starting from 0).

**set\_roi** (*hstart=0*, *hend=None*, *vstart=0*, *vend=None*)

Setup camera ROI.

*hstart* and *hend* specify horizontal image extent, *vstart* and *vend* specify vertical image extent (start is inclusive, stop is exclusive, starting from 0). By default, all non-supplied parameters take extreme values (0 for start, maximal for end).

**get\_roi\_limits** (*hbin=1*, *vbin=1*)

Get the minimal and maximal ROI parameters.

Return tuple (*hlim*, *vlim*), where each element is in turn a limit 5-tuple (*min*, *max*, *pstep*, *sstep*, *maxbin*) with, correspondingly, minimal and maximal size, position and size step, and the maximal binning (fixed to 1 if not binning is allowed). In some cameras, the step and the minimal size depend on the binning, which can be supplied.

**setup\_acquisition** (*mode='sequence'*, *nframes=100*)

Setup acquisition mode.

*mode* can be either "snap" (single frame or a fixed number of frames) or "sequence" (continuous acquisition). (note that `IMAQdxCamera.acquisition_in_progress()` would still return True in this case, even though new frames are no longer acquired). *nframes* sets up number of frame buffers.

**clear\_acquisition()**

Clear acquisition settings

**start\_acquisition(\*args, \*\*kwargs)**

Start acquisition.

Can take the same keyword parameters as `:meth:~setup_acquisition`. If the acquisition is not set up yet, set it up using the supplied parameters (use default of `setup_acquisition()`, if the parameter is None). Otherwise, if any supplied parameters are different from the current ones, change them and reset the acquisition.

**stop\_acquisition()**

Stop acquisition

**acquisition\_in\_progress()**

Check if acquisition is in progress

**refresh\_acquisition(delay=0.005)**

Stop and restart the acquisition, waiting *delay* seconds in between

**enable\_raw\_readout(enable='rows')**

Enable raw frame transfer.

Should be used if the camera uses unsupported pixel format. Can be "frame" (return the whole frame as a 1D "u1" numpy array), "rows" (return a 2D array, where each row corresponds to a single image row), or False (convert to image data, or raise an error if the format is not supported; default)

**FrameTransferError**

alias of `pylablib.devices.interface.camera.DefaultFrameTransferError`

**class NoParameterCaller(device, kind)**

Bases: `object`

Class to simplify calling functions without a parameter

**apply\_settings(settings)**

Apply the settings.

*settings* is the dict {name: value} of the device available settings. Non-applicable settings are ignored.

**get\_acquisition\_parameters()**

Get acquisition parameters.

Return dictionary {name: value}

**get\_all\_attributes(copy=False)**

Return a dictionary of all available attributes.

If *copy*==True, copy the dictionary; otherwise, return the internal dictionary structure (should not be modified).

**get\_attribute(name, error\_on\_missing=True)**

Get the camera attribute with the given name

**get\_data\_dimensions()**

Get readout data dimensions (in pixels) as a tuple (width, height); take indexing mode into account

**get\_device\_variable(key)**

Get the value of a settings, status, or full info parameter



**get\_frame\_format()**

Get format for the returned images.

Can be "list" (list of 2D arrays), "array" (a single 3D array), or "chunks" (list of 3D "chunk" arrays; supported for some cameras and provides the best performance).

**get\_frame\_info\_fields()**

Get the names of frame info fields.

Applicable when frame info format (set by `set_frame_info_format()`) is "list" or "array".

**get\_frame\_info\_format()**

Get format of the frame info.

Can be "namedtuple" (potentially nested named tuples; convenient to get particular values), "list" (flat list of values, with field names are given by `get_frame_info_fields()`; convenient for building a table), "array" (same as "list", but with a numpy array, which is easier to use for "chunks" frame format), or "dict" (flat dictionary with the same fields as the "list" format; more resilient to future format changes)

**get\_frame\_info\_period()**

Get period of frame info acquisition.

Frame info might be skipped (set to `None`) except for frames which indices are divisible by *period*. Useful for certain cameras where acquiring frame info takes a lot of time and can reduce performance at higher frame rates. Note that this parameter can still be ignored (i.e., always set to 1) if the performance is not an issue for a given camera class.

**get\_frames\_status()**

Get acquisition and buffer status.

Return tuple (acquired, unread, skipped, size), where acquired is the total number of acquired frames, unread is the number of acquired but not read frames, skipped is the number of skipped (not read and then written over) frames, and `buffer_size` is the total buffer size (in frames).

**get\_full\_info(include=0)**

Get dict {name: value} containing full device information (including status and settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting `include=-10` queries all available variables, which is equivalent to `include="all"`.

**get\_full\_status(include=0)**

Get dict {name: value} containing the device status (including settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting `include=-10` queries all available variables, which is equivalent to `include="all"`.

**get\_image\_indexing()**

Get indexing for the returned images.

Can be "rct" (first index row, second index column, rows counted from the top), "rcb" (same as "rc", rows counted from the bottom), "xyt" (first index column, second index row, rows counted from the top), or "xyb" (same as "xyt", rows counted from the bottom)

**get\_new\_images\_range()**

Get the range of the new images.

Return tuple (*first*, *last*) with images range (first inclusive). If no images are available, return None. If some images were in the buffer were overwritten, exclude them from the range.

**get\_settings** (*include=0*)

Get dict {*name*: *value*} containing all the device settings.

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**grab** (*nframes=1*, *frame\_timeout=5.0*, *missing\_frame='skip'*, *return\_info=False*, *buff\_size=None*)

Snap *nframes* images (with preset image read mode parameters)

*buff\_size* determines buffer size (if None, use the default size). Timeout is specified for a single-frame acquisition, not for the whole acquisition time. *missing\_frame* determines what to do with frames which have been lost: can be "none" (replacing them with None), "zero" (replacing them with zero-filled frame), or "skip" (skipping them, while still keeping total returned frames number to *n*). If *return\_info==True*, return tuple (*frames*, *infos*), where *infos* is a list of frame info tuples (camera-dependent); if some frames are missing and *missing\_frame!="skip"*, the corresponding frame info is None.

**is\_acquisition\_setup** ()

Check if acquisition is set up.

If the camera does not support separate acquisition setup, always return True.

**pausing\_acquisition** (*clear=None*, *stop=True*, *setup\_after=None*, *start\_after=True*, *combine\_nested=True*)

Context manager which temporarily pauses acquisition during execution of *with* block.

Useful for applying certain settings which can't be changed during the acquisition. If *clear==True*, clear acquisition in addition to stopping (by default, use the class default specified as *\_clear\_pausing\_acquisition* attribute). If *stop==True*, stop the acquisition (if *clear==True*, stop regardless). If *setup\_after==True*, setup the acquisition after pause if necessary (None means setup only if clearing was required). If *start\_after==True*, start the acquisition after pause if necessary (None means start only if stopping was required). If *combine\_nested==True*, then any nested *pausing\_acquisition* calls will stop/clear acquisition as necessary, but won't setup/start it again until this *pausing\_acquisition* call is complete.

Yields tuple (*acq\_in\_progress*, *acq\_params*), which indicates whether acquisition is currently in progress, and what are the current acquisition parameters.

**read\_multiple\_images** (*rng=None*, *peek=False*, *missing\_frame='skip'*, *return\_info=False*, *return\_rng=False*)

Read multiple images specified by *rng* (by default, all un-read images).

If *rng* is specified, it is a tuple (*first*, *last*) with images range (first inclusive). If no new frames are available, return an empty list; if no acquisition is running, return None. If *peek==True*, return images but not mark them as read. *missing\_frame* determines what to do with frames which are out of range (missing or lost): can be "none" (replacing them with None), "zero" (replacing them with zero-filled frame), or "skip" (skipping them). If *return\_info==True*, return tuple (*frames*, *infos*), where *infos* is a list of frame info tuples (camera-dependent, by default, only the frame index); if some frames are missing and *missing\_frame!="skip"*, the corresponding frame info is None. if *return\_rng==True*, return the range covered resulting frames; if *missing\_frame=="skip"*, the range can be smaller than the supplied *rng* if some frames are skipped.

**read\_newest\_image** (*peek=False*, *return\_info=False*)

Read the newest un-read image.

If no un-read frames are available, return None. If `peek==True`, return the image but not mark it as read. If `return_info==True`, return tuple (frame, info), where info is an info tuples (camera-dependent, see [read\\_multiple\\_images\(\)](#)).

**read\_oldest\_image** (*peek=False, return\_info=False*)

Read the oldest un-read image.

If no un-read frames are available, return None. If `peek==True`, return the image but not mark it as read. If `return_info==True`, return tuple (frame, info), where info is an info tuples (camera-dependent, see [read\\_multiple\\_images\(\)](#)).

**set\_device\_variable** (*key, value*)

Set the value of a settings parameter

**set\_frame\_format** (*fmt*)

Set format for the returned images.

Can be "list" (list of 2D arrays), "array" (a single 3D array), "chunks" (list of 3D "chunk" arrays; supported for some cameras and provides the best performance), or "try\_chunks" (same as "chunks", but if chunks are not supported, set to "list" instead). If format is "chunks" and chunks are not supported by the camera, it results in one frame per chunk. Note that if the format is set to "array" or "chunks", the frame info format is also automatically set to "array". If the format is set to "chunks", then the image info is also returned in chunks form (list of 2D info arrays with the same length as the corresponding frame chunks).

**set\_frame\_info\_format** (*fmt, include\_fields=None*)

Set format of the frame info.

Can be "namedtuple" (potentially nested named tuples; convenient to get particular values), "list" (flat list of values, with field names are given by [get\\_frame\\_info\\_fields\(\)](#); convenient for building a table), "array" (same as "list", but with a numpy array, which is easier to use for "chunks" frame format), or "dict" (flat dictionary with the same fields as the "list" format; more resilient to future format changes) If *include\_fields* is not None, it specifies the fields included for non-"tuple" formats; note that order or *include\_fields* is ignored, and the resulting fields are always ordered same as in the original.

**set\_frame\_info\_period** (*period=1*)

Set period of frame info acquisition.

Frame info might be skipped (set to None) except for frames which indices are divisible by *period*. Useful for certain cameras where acquiring frame info takes a lot of time and can reduce performance at higher frame rates. Note that this parameter can still be ignored (i.e., always set to 1) if the performance is not an issue for a given camera class.

**set\_image\_indexing** (*indexing*)

Set up indexing for the returned images.

Can be "rct" (first index row, second index column, rows counted from the top), "rcb" (same as "rc", rows counted from the bottom), "xyt" (first index column, second index row, rows counted from the top), or "xyb" (same as "xyt", rows counted from the bottom)

**snap** (*timeout=5.0, return\_info=False*)

Snap a single frame

**wait\_for\_frame** (*since='lastread', nframes=1, timeout=20.0, error\_on\_stopped=False*)

Wait for one or several new camera frames.

*since* specifies the reference point for waiting to acquire *nframes* frames; can be "lastread" (from the last read frame), "lastwait" (wait for the last successful [wait\\_for\\_frame\(\)](#) call), "now" (from the start of the current call), or "start" (from the acquisition start, i.e., wait until *nframes* frames have been acquired). *timeout* can be either a number, None (infinite timeout), or a tuple

(timeout, frame\_timeout), in which case the call times out if the total time exceeds timeout, or a single frame wait exceeds frame\_timeout. If the call times out, raise `TimeoutError`. If `error_on_stopped==True` and the acquisition is not running, raise `Error`; otherwise, simply return `False` without waiting.

```
class pylablib.devices.IMAQdx.IMAQdx.EthernetIMAQdxCamera (name='cam0',
                                                         mode='controller',
                                                         visibility='advanced',
                                                         small_packet=False)
```

Bases: `pylablib.devices.IMAQdx.IMAQdx.IMAQdxCamera`

LAN-controlled IMAQdx camera.

Compared to the standard camera, has an option of automatically switching to a smaller TCP/IP packet size (can be useful if the PC network adapter can't handle jumbo packets).

#### Parameters

- **name** – interface name (can be learned by `list_cameras()`; usually, but not always, starts with "cam")
- **mode** – connection mode; can be "controller" (full control) or "listener" (only reading)
- **visibility** – default attribute visibility when listing attributes; can be "simple", "intermediate" or "advanced" (higher mode exposes more attributes).
- **small\_packet** – if `True`, automatically set small packet size (1500 bytes).

**post\_open()**

Additional setup after camera opening

**Error** = <Mock name='mock.IMAQdxError' id='140096275815696'>

**FrameTransferError**

alias of `pylablib.devices.interface.camera.DefaultFrameTransferError`

**class NoParameterCaller** (device, kind)

Bases: `object`

Class to simplify calling functions without a parameter

**TimeoutError** = <Mock spec='str' id='140096275815376'>

**acquisition\_in\_progress()**

Check if acquisition is in progress

**apply\_settings** (settings)

Apply the settings.

*settings* is the dict {name: value} of the device available settings. Non-applicable settings are ignored.

**clear\_acquisition()**

Clear acquisition settings

**close()**

Close connection to the camera

**enable\_raw\_readout** (enable='rows')

Enable raw frame transfer.

Should be used if the camera uses unsupported pixel format. Can be "frame" (return the whole frame as a 1D "u1" numpy array), "rows" (return a 2D array, where each row corresponds to a single image row), or `False` (convert to image data, or raise an error if the format is not supported; default)

**get\_acquisition\_parameters()**

Get acquisition parameters.

Return dictionary {name: value}

**get\_all\_attribute\_values** (*root=""*, *enum\_as\_str=True*, *ignore\_errors=True*)

Get values of all attributes with the given *root*

**get\_all\_attributes** (*copy=False*)

Return a dictionary of all available attributes.

If *copy==True*, copy the dictionary; otherwise, return the internal dictionary structure (should not be modified).

**get\_attribute** (*name*, *error\_on\_missing=True*)

Get the camera attribute with the given name

**get\_attribute\_value** (*name*, *error\_on\_missing=True*, *default=None*, *enum\_as\_str=True*)

Get value of an attribute with the given name.

If the value doesn't exist or can not be read and *error\_on\_missing==True*, raise error; otherwise, return *default*. If *default* is not *None*, assume that *error\_on\_missing==False*. If *name* points at a dictionary branch, return a dictionary with all values in this branch. If *enum\_as\_str==True*, return enum-style values as strings; otherwise, return corresponding integer values.

**get\_data\_dimensions()**

Get readout data dimensions (in pixels) as a tuple (*width*, *height*); take indexing mode into account

**get\_detector\_size()**

Get camera detector size (in pixels) as a tuple (*width*, *height*)

**get\_device\_info()**

Get camera information.

Return tuple (*vendor*, *model*, *serial\_number*, *bus\_type*).

**get\_device\_variable** (*key*)

Get the value of a settings, status, or full info parameter

**get\_frame\_format()**

Get format for the returned images.

Can be "list" (list of 2D arrays), "array" (a single 3D array), or "chunks" (list of 3D "chunk" arrays; supported for some cameras and provides the best performance).

**get\_frame\_info\_fields()**

Get the names of frame info fields.

Applicable when frame info format (set by `set_frame_info_format()`) is "list" or "array".

**get\_frame\_info\_format()**

Get format of the frame info.

Can be "namedtuple" (potentially nested named tuples; convenient to get particular values), "list" (flat list of values, with field names are given by `get_frame_info_fields()`; convenient for building a table), "array" (same as "list", but with a numpy array, which is easier to use for "chunks" frame format), or "dict" (flat dictionary with the same fields as the "list" format; more resilient to future format changes)

**get\_frame\_info\_period()**

Get period of frame info acquisition.

Frame info might be skipped (set to *None*) except for frames which indices are divisible by *period*. Useful for certain cameras where acquiring frame info takes a lot of time and can reduce performance at higher

frame rates. Note that this parameter can still be ignored (i.e., always set to 1) if the performance is not an issue for a given camera class.

**get\_frames\_status()**

Get acquisition and buffer status.

Return tuple (*acquired*, *unread*, *skipped*, *size*), where *acquired* is the total number of acquired frames, *unread* is the number of acquired but not read frames, *skipped* is the number of skipped (not read and then written over) frames, and *buffer\_size* is the total buffer size (in frames).

**get\_full\_info(include=0)**

Get dict {*name*: *value*} containing full device information (including status and settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include*=-10 queries all available variables, which is equivalent to *include*="all".

**get\_full\_status(include=0)**

Get dict {*name*: *value*} containing the device status (including settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include*=-10 queries all available variables, which is equivalent to *include*="all".

**get\_image\_indexing()**

Get indexing for the returned images.

Can be "rct" (first index row, second index column, rows counted from the top), "rcb" (same as "rc", rows counted from the bottom), "xyt" (first index column, second index row, rows counted from the top), or "xyb" (same as "xyt", rows counted from the bottom)

**get\_new\_images\_range()**

Get the range of the new images.

Return tuple (*first*, *last*) with images range (first inclusive). If no images are available, return None. If some images were in the buffer were overwritten, exclude them from the range.

**get\_roi()**

Get current ROI.

Return tuple (*hstart*, *hend*, *vstart*, *vend*). *hstart* and *hend* specify horizontal image extent, *vstart* and *vend* specify vertical image extent (start is inclusive, stop is exclusive, starting from 0).

**get\_roi\_limits(hbin=1, vbin=1)**

Get the minimal and maximal ROI parameters.

Return tuple (*hlim*, *vlim*), where each element is in turn a limit 5-tuple (*min*, *max*, *pstep*, *sstep*, *maxbin*) with, correspondingly, minimal and maximal size, position and size step, and the maximal binning (fixed to 1 if not binning is allowed). In some cameras, the step and the minimal size depend on the binning, which can be supplied.

**get\_settings(include=0)**

Get dict {*name*: *value*} containing all the device settings.

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include*=-10 queries all available variables, which is equivalent to *include*="all".

**grab** (*nframes=1, frame\_timeout=5.0, missing\_frame='skip', return\_info=False, buff\_size=None*)  
 Snap *nframes* images (with preset image read mode parameters)

*buff\_size* determines buffer size (if *None*, use the default size). Timeout is specified for a single-frame acquisition, not for the whole acquisition time. *missing\_frame* determines what to do with frames which have been lost: can be "none" (replacing them with *None*), "zero" (replacing them with zero-filled frame), or "skip" (skipping them, while still keeping total returned frames number to *n*). If *return\_info==True*, return tuple (*frames*, *infos*), where *infos* is a list of frame info tuples (camera-dependent); if some frames are missing and *missing\_frame!="skip"*, the corresponding frame info is *None*.

**is\_acquisition\_setup** ()  
 Check if acquisition is set up.

If the camera does not support separate acquisition setup, always return *True*.

**is\_opened** ()  
 Check if the device is connected

**open** ()  
 Open connection to the camera

**pausing\_acquisition** (*clear=None, stop=True, setup\_after=None, start\_after=True, combine\_nested=True*)  
 Context manager which temporarily pauses acquisition during execution of *with* block.

Useful for applying certain settings which can't be changed during the acquisition. If *clear==True*, clear acquisition in addition to stopping (by default, use the class default specified as *\_clear\_pausing\_acquisition* attribute). If *stop==True*, stop the acquisition (if *clear==True*, stop regardless). If *setup\_after==True*, setup the acquisition after pause if necessary (*None* means setup only if clearing was required). If *start\_after==True*, start the acquisition after pause if necessary (*None* means start only if stopping was required). If *combine\_nested==True*, then any nested *pausing\_acquisition* calls will stop/clear acquisition as necessary, but won't setup/start it again until this *pausing\_acquisition* call is complete.

Yields tuple (*acq\_in\_progress*, *acq\_params*), which indicates whether acquisition is currently in progress, and what are the current acquisition parameters.

**read\_multiple\_images** (*rng=None, peek=False, missing\_frame='skip', return\_info=False, return\_rng=False*)  
 Read multiple images specified by *rng* (by default, all un-read images).

If *rng* is specified, it is a tuple (*first*, *last*) with images range (first inclusive). If no new frames are available, return an empty list; if no acquisition is running, return *None*. If *peek==True*, return images but not mark them as read. *missing\_frame* determines what to do with frames which are out of range (missing or lost): can be "none" (replacing them with *None*), "zero" (replacing them with zero-filled frame), or "skip" (skipping them). If *return\_info==True*, return tuple (*frames*, *infos*), where *infos* is a list of frame info tuples (camera-dependent, by default, only the frame index); if some frames are missing and *missing\_frame!="skip"*, the corresponding frame info is *None*. If *return\_rng==True*, return the range covered resulting frames; if *missing\_frame=="skip"*, the range can be smaller than the supplied *rng* if some frames are skipped.

**read\_newest\_image** (*peek=False, return\_info=False*)  
 Read the newest un-read image.

If no un-read frames are available, return *None*. If *peek==True*, return the image but not mark it as read. If *return\_info==True*, return tuple (*frame*, *info*), where *info* is an info tuples (camera-dependent, see [read\\_multiple\\_images\(\)](#)).

**read\_oldest\_image** (*peek=False, return\_info=False*)  
 Read the oldest un-read image.

If no un-read frames are available, return None. If `peek==True`, return the image but not mark it as read. If `return_info==True`, return tuple (frame, info), where info is an info tuples (camera-dependent, see [read\\_multiple\\_images\(\)](#)).

**refresh\_acquisition** (*delay=0.005*)

Stop and restart the acquisition, waiting *delay* seconds in between

**reset** ()

Reset connection to the camera

**set\_all\_attribute\_values** (*settings, root="", truncate=True*)

Set values of all attributes with the given *root*.

If `truncate==True`, truncate value to lie within attribute range.

**set\_attribute\_value** (*name, value, truncate=True, error\_on\_missing=True*)

Set value of an attribute with the given name.

If the value doesn't exist or can not be written and `error_on_missing==True`, raise error; otherwise, do nothing. If *name* points at a dictionary branch, set all values in this branch (in this case *value* must be a dictionary). If `truncate==True`, truncate value to lie within attribute range.

**set\_device\_variable** (*key, value*)

Set the value of a settings parameter

**set\_frame\_format** (*fmt*)

Set format for the returned images.

Can be "list" (list of 2D arrays), "array" (a single 3D array), "chunks" (list of 3D "chunk" arrays; supported for some cameras and provides the best performance), or "try\_chunks" (same as "chunks", but if chunks are not supported, set to "list" instead). If format is "chunks" and chunks are not supported by the camera, it results in one frame per chunk. Note that if the format is set to "array" or "chunks", the frame info format is also automatically set to "array". If the format is set to "chunks", then the image info is also returned in chunks form (list of 2D info arrays with the same length as the corresponding frame chunks).

**set\_frame\_info\_format** (*fmt, include\_fields=None*)

Set format of the frame info.

Can be "namedtuple" (potentially nested named tuples; convenient to get particular values), "list" (flat list of values, with field names are given by [get\\_frame\\_info\\_fields\(\)](#); convenient for building a table), "array" (same as "list", but with a numpy array, which is easier to use for "chunks" frame format), or "dict" (flat dictionary with the same fields as the "list" format; more resilient to future format changes) If *include\_fields* is not None, it specifies the fields included for non-"tuple" formats; note that order or *include\_fields* is ignored, and the resulting fields are always ordered same as in the original.

**set\_frame\_info\_period** (*period=1*)

Set period of frame info acquisition.

Frame info might be skipped (set to None) except for frames which indices are divisible by *period*. Useful for certain cameras where acquiring frame info takes a lot of time and can reduce performance at higher frame rates. Note that this parameter can still be ignored (i.e., always set to 1) if the performance is not an issue for a given camera class.

**set\_image\_indexing** (*indexing*)

Set up indexing for the returned images.

Can be "rct" (first index row, second index column, rows counted from the top), "rcb" (same as "rc", rows counted from the bottom), "xyt" (first index column, second index row, rows counted from the top), or "xyb" (same as "xyt", rows counted from the bottom)



**set\_roi** (*hstart=0, hend=None, vstart=0, vend=None*)

Setup camera ROI.

*hstart* and *hend* specify horizontal image extent, *vstart* and *vend* specify vertical image extent (start is inclusive, stop is exclusive, starting from 0). By default, all non-supplied parameters take extreme values (0 for start, maximal for end).

**setup\_acquisition** (*mode='sequence', nframes=100*)

Setup acquisition mode.

*mode* can be either "snap" (single frame or a fixed number of frames) or "sequence" (continuous acquisition). (note that `IMAQdxCamera.acquisition_in_progress()` would still return True in this case, even though new frames are no longer acquired). *nframes* sets up number of frame buffers.

**snap** (*timeout=5.0, return\_info=False*)

Snap a single frame

**start\_acquisition** (*\*args, \*\*kwargs*)

Start acquisition.

Can take the same keyword parameters as `:meth:'setup_acquisition`. If the acquisition is not set up yet, set it up using the supplied parameters (use default of `setup_acquisition()`, if the parameter is None). Otherwise, if any supplied parameters are different from the current ones, change them and reset the acquisition.

**stop\_acquisition** ()

Stop acquisition

**wait\_for\_frame** (*since='lastread', nframes=1, timeout=20.0, error\_on\_stopped=False*)

Wait for one or several new camera frames.

*since* specifies the reference point for waiting to acquire *nframes* frames; can be "lastread" (from the last read frame), "lastwait" (wait for the last successful `wait_for_frame()` call), "now" (from the start of the current call), or "start" (from the acquisition start, i.e., wait until *nframes* frames have been acquired). *timeout* can be either a number, None (infinite timeout), or a tuple (*timeout*, *frame\_timeout*), in which case the call times out if the total time exceeds *timeout*, or a single frame wait exceeds *frame\_timeout*. If the call times out, raise `TimeoutError`. If `error_on_stopped==True` and the acquisition is not running, raise `Error`; otherwise, simply return False without waiting.

## Module contents

### pylablib.devices.KJL package

#### Submodules

#### pylablib.devices.KJL.base module

**exception** `pylablib.devices.KJL.base.KJLError`

Bases: `pylablib.core.devio.base.DeviceError`

Generic KJL device error

**args**

**with\_traceback** ()

Exception.with\_traceback(tb) – set self.\_\_traceback\_\_ to tb and return self.

**exception** `pylablib.devices.KJL.base.KJLBackendError` (*exc*)

Bases: `pylablib.devices.KJL.base.KJLError`, `pylablib.core.devio.comm_backend.DeviceBackendError`

Generic KJL backend communication error

**args**

**with\_traceback** ()

Exception.with\_traceback(tb) – set self.\_\_traceback\_\_ to tb and return self.

**class** `pylablib.devices.KJL.base.TKJL300DeviceInfo` (*swver*)

Bases: `tuple`

**count** ()

Return number of occurrences of value.

**index** ()

Return first index of value.

Raises ValueError if the value is not present.

**swver**

**class** `pylablib.devices.KJL.base.KJL300` (*conn*, *addr=1*)

Bases: `pylablib.core.devio.comm_backend.ICommBackendWrapper`

KJL300 series pressure gauge.

**Parameters**

- **conn** – serial connection parameters (usually port or a tuple containing port and baudrate)
- **addr** – RS485 address (required both for RS-485 and for RS-232 communication; factory default is 1)

**Error**

alias of `KJLError`

**comm** (*msg*)

Send a command to the device

**query** (*msg*)

**get\_device\_info** ()

Get device info (a tuple (*swver*))

**reset** (*confirm\_addr=False*)

Reset the controller.

If *confirm\_addr*==True, set current RS485 address again (required for resetting after some commands).

**get\_pressure** ()

Get current pressure in Pa

**get\_relay\_setpoints** (*relay=1*)

Get relay setpoints (in Pa).

*relay* is the relay index (either 1 or 2). Return tuple (*on*, *off*) for on-below and off-above pressures (*on* is always smaller than *off*)

**set\_relay\_setpoints** (*relay=1*, *on=None*, *off=None*, *reset=True*)

Set relay setpoints (in Pa).

*relay* is the relay index (either 1 or 2). *on* and *off* are on-below and off-above pressures (*on* is always smaller than *off*). If *reset==True*, reset the device after changing the setpoints (required to take effect). None values are left unchanged.

**set\_zero** (*pressure=0*)

Set vacuum calibration point (in Pa)

**set\_span** (*pressure=100000.0*)

Set atmosphere calibration point (in Pa)

**class NoParameterCaller** (*device, kind*)

Bases: `object`

Class to simplify calling functions without a parameter

**apply\_settings** (*settings*)

Apply the settings.

*settings* is the dict {*name*: *value*} of the device available settings. Non-applicable settings are ignored.

**close** ()

Close the backend

**get\_device\_variable** (*key*)

Get the value of a settings, status, or full info parameter

**get\_full\_info** (*include=0*)

Get dict {*name*: *value*} containing full device information (including status and settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_full\_status** (*include=0*)

Get dict {*name*: *value*} containing the device status (including settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_settings** (*include=0*)

Get dict {*name*: *value*} containing all the device settings.

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**is\_opened** ()

Check if the device is connected

**lock** (*timeout=None*)

Lock the access to the device from other threads/processes (isn't necessarily implemented)

**locking** (*timeout=None*)

Context manager for lock & unlock

**open** ()

Open the backend

**set\_device\_variable** (*key, value*)

Set the value of a settings parameter

**unlock** ()

Unlock the access to the device from other threads/processes (isn't necessarily implemented)

## Module contents

### pylablib.devices.Lakeshore package

#### Submodules

#### pylablib.devices.Lakeshore.base module

**exception** pylablib.devices.Lakeshore.base.**LakeshoreError**

Bases: *pylablib.core.devio.base.DeviceError*

Generic Lakeshore devices error

**args**

**with\_traceback** ()

Exception.with\_traceback(tb) – set self.\_\_traceback\_\_ to tb and return self.

**exception** pylablib.devices.Lakeshore.base.**LakeshoreBackendError** (*exc*)

Bases: *pylablib.devices.Lakeshore.base.LakeshoreError, pylablib.core.devio.comm\_backend.DeviceBackendError*

Generic Lakeshore backend communication error

**args**

**with\_traceback** ()

Exception.with\_traceback(tb) – set self.\_\_traceback\_\_ to tb and return self.

**class** pylablib.devices.Lakeshore.base.**TLakeshore218AnalogSettings** (*bipolar,*  
*mode,*  
*channel,*  
*source,*  
*high\_value,*  
*low\_value,*  
*man\_value*)

Bases: *tuple*

**bipolar**

**channel**

**count** ()

Return number of occurrences of value.

**high\_value**

**index** ()

Return first index of value.

Raises ValueError if the value is not present.

**low\_value**

**man\_value**

**mode****source**

```
class pylablib.devices.Lakeshore.base.TLakeshore218FilterSettings (enabled,  

                                                                    points,  

                                                                    window)
```

Bases: `tuple`**count** ()

Return number of occurrences of value.

**enabled****index** ()

Return first index of value.

Raises ValueError if the value is not present.

**points****window**

```
class pylablib.devices.Lakeshore.base.TLakeshore218CurveHeader (name,      serial,  

                                                                    fmt,      limit,  

                                                                    coeff)
```

Bases: `tuple`**coeff****count** ()

Return number of occurrences of value.

**fmt****index** ()

Return first index of value.

Raises ValueError if the value is not present.

**limit****name****serial**

```
class pylablib.devices.Lakeshore.base.Lakeshore218 (conn)
```

Bases: `pylablib.core.devio.SCPIDevice`

Lakeshore 218 temperature controller.

The channels are enumerated from 1 to 8 and are split into 2 groups: "A" for 1-4 and "B" for 5-8.

**Parameters** **conn** – serial connection parameters (usually port or a tuple containing port and baudrate)**Error**alias of `LakeshoreError`**ReraiseError**alias of `LakeshoreBackendError`**is\_enabled** (*channel*)

Check if a given channel is enabled

**set\_enabled** (*channel, enabled=True*)

Enable or disable a given channel

**get\_sensor\_type** (*group*)

Get sensor type for a given group ("A" for sensors 1-4 or "B" for sensors 5-8).

For types, see INTYPE command description in the Lakeshore 218 programming manual.

**set\_sensor\_type** (*group*, *sensor\_type*)

Set sensor type for a given group ("A" for sensors 1-4 or "B" for sensors 5-8).

For types, see INTYPE command description in the Lakeshore 218 programming manual.

**get\_sensor\_curve\_index** (*channel*)

Get sensor curve index for a given channel (1 to 8).

For curve descriptions, see INCRV command description in the Lakeshore 218 programming manual.

**set\_sensor\_curve\_index** (*channel*, *index*)

Get sensor curve index for a given channel (1 to 8).

For curve descriptions, see INCRV command description in the Lakeshore 218 programming manual.

**get\_curve\_header** (*index*)

Get header of a given curve (1-9 or 21-28).

Return tuple (*name*, *serial*, *fmt*, *limit*, *coeff*). For values descriptions, see CRVHDR command description in the Lakeshore 218 programming manual.

**set\_curve\_header** (*index*, *name=None*, *serial=None*, *fmt=None*, *limit=None*, *coeff=None*)

Set header of a given user curve (21-28).

For values descriptions, see CRVHDR command description in the Lakeshore 218 programming manual.

**get\_curve** (*index*, *trim\_zeros=True*)

Get values of a given curve (1-9 or 21-28).

Return 2-column numpy array with up to 200 points, where the first column is sensor reading, and the second is temperature; for associated sensor units, see [get\\_curve\\_header\(\)](#). If *trim\_zeros==True*, trim the trailing zero-valued points. Note, that it takes about 10 seconds to complete.

**set\_curve** (*index*, *curve*)

Set values of a given user curve (21-28).

*curve* is a 2-column numpy array with up to 200 points, where the first column is sensor reading, and the second is temperature; for associated sensor units, see [get\\_curve\\_header\(\)](#). Note, that it takes about 20 seconds to complete.

**get\_temperature** (*channel*)

Get readings (in Kelvin) on a given channel (1 to 8)

**get\_all\_temperatures** ()

Get readings (in Kelvin) on all channels

**get\_sensor\_reading** (*channel*)

Get readings (in sensor units) on a given channel (1 to 8)

**get\_all\_sensor\_readings** ()

Get readings (in sensor units) on all channels

**get\_analog\_output\_settings** (*output*)

Get analog output settings for a given output (1 or 2).

For parameters, see [setup\\_analog\\_output\(\)](#) and ANALOG command description in the Lakeshore 218 programming manual.

**setup\_analog\_output** (*output*, *bipolar=None*, *mode=None*, *channel=None*, *source=None*, *high\_value=None*, *low\_value=None*, *man\_value=None*)

Setup analog output settings for a given output (1 or 2).

For parameters, see ANALOG command description in the Lakeshore 218 programming manual. Value of None means keeping the current parameter value.

**set\_analog\_output\_value** (*output*, *value*, *bipolar=False*, *enabled=True*)

Set manual analog output value.

A simplified version of `setup_analog_output()`.

**get\_analog\_output** (*output*)

Get value (in percents of the total range) at a given output (1 or 2)

**get\_filter\_settings** (*channel*)

Get input filter settings for a given channel (1 to 8).

For parameters, see `setup_filter()` and FILTER command description in the Lakeshore 218 programming manual.

**setup\_filter** (*channel*, *enabled=None*, *points=None*, *window=None*)

Setup input filter settings for a given channel (1 to 8).

For parameters, see FILTER command description in the Lakeshore 218 programming manual. Value of None means keeping the current parameter value.

**BackendError**

alias of `pylablib.core.devio.comm_backend.DeviceBackendError`

**class NoParameterCaller** (*device*, *kind*)

Bases: `object`

Class to simplify calling functions without a parameter

**apply\_settings** (*settings*)

Apply the settings.

*settings* is a dict {name: value} of the available device settings. Non-applicable settings are ignored.

**ask** (*msg*, *data\_type='string'*, *delay=0.0*, *timeout=None*, *read\_echo=False*)

Write a message and read a reply.

*msg* is the query message, *delay* is the delay between write and read. Other parameters are the same as in `read()`. If *read\_echo==True*, assume that the device first echoes the input and skip it.

**close** ()

Close the backend

**flush** (*one\_line=False*)

Flush the read buffer (read all the available data and return the number of bytes read).

If *one\_line==True*, read only a single line.

**static get\_arg\_type** (*arg*)

Autodetect argument type

**get\_device\_variable** (*key*)

Get the value of a settings, status, or full info parameter

**get\_esr** (*timeout=None*)

Get the device status register (by default, "\*ESR?" command)

**get\_full\_info** (*include=0*)

Get dict {name: value} containing full device information (including status and settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include*=-10 queries all available variables, which is equivalent to *include*="all".

**get\_full\_status** (*include*=0)

Get dict {name: value} containing the device status (including settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include*=-10 queries all available variables, which is equivalent to *include*="all".

**get\_id** (*timeout*=None)

Get the device IDN. (query SCPI '\*IDN?' command)

**get\_settings** (*include*=0)

Get dict {name: value} containing all the device settings.

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include*=-10 queries all available variables, which is equivalent to *include*="all".

**is\_opened** ()

Check if the device is connected

**lock** (*timeout*=None)

Lock the access to the device from other threads/processes (isn't necessarily implemented)

**locking** (*timeout*=None)

Context manager for lock & unlock

**open** ()

Open the backend

**static parse\_array\_data** (*data*, *fmt*, *include\_header*=False)

Parse the data returned by the device. *fmt* is *DataFormat* description in numpy format (e.g., "<u2").

If *include\_header*==True, the data is assumed to be in a (somewhat) standard SCPI format: b'# ', then a single digit *s* denoting length of the size block, then *s* digits denoting length of the data (in bytes) followed by the actual data. Otherwise (*include\_header*==False), assume that the header is already removed.

**read** (*data\_type*='string', *timeout*=None)

Read data from the device.

*data\_type* determines the type of the data. Can be 'raw' (just raw data), 'string' (with trailing and leading spaces stripped), 'int', 'float', 'bool' (interprets 0 or 'off' as False, anything else as True), 'value' (returns tuple (value, unit), where *value* is float), a callable (return the result of this callable applied to the string value), a dictionary (return the stored value corresponding to the string value, or to the value converted into integer if the string value is not present), or a list of data types (the result is treated as a list of values with the given types separated by commas). *timeout* overrides the default value.

**read\_binary\_array\_data** (*include\_header*=False, *timeout*=None, *flush\_term*=True)

Read a binary data in the from the device.

The data assumes the standard binary transfer header consisting of "#" symbol, then a single digit with the size of the length string, then the length string containing the length of the binary data (in bytes). If *include\_header*==True, return the data with the header; otherwise, return only the content. If



`flush_term==True`, flush the following line to skip terminator characters after the binary data, which are added by some devices. `timeout` overrides the default value.

**reconnect** (*new\_instrument=True, ignore\_error=True*)

Remake the connection.

If `new_instrument==True`, create a new backend instance. If `ignore_error==True`, ignore errors on closing.

**reset** ()

Reset the device (by default, "`*RST`" command)

**set\_device\_variable** (*key, value*)

Set the value of a settings parameter

**sleep** (*delay*)

Wait for *delay* seconds

**unlock** ()

Unlock the access to the device from other threads/processes (isn't necessarily implemented)

**using\_write\_buffer** ()

Context manager for using a write buffer.

While it's active, all the consecutive `write()` operations are bundled together with `;` delimiter. The actual write is performed at the `read()/ask()` operation or at the end of the block.

**wait** (*wait\_type='sync', timeout=None, wait\_callback=None*)

Pause execution until device overlapped commands are complete.

*wait\_type* is either `'sync'` (perform `wait_sync()`), `'dev'` (perform `wait_dev()`) or `'none'` (do nothing).

**wait\_dev** ()

Pause execution of the device commands until device overlapped commands (e.g., taking sweeps) are complete.

Note that the code execution is not paused.

**wait\_sync** (*timeout=None, wait\_callback=None*)

Pause execution of the script until device overlapped commands (e.g., taking sweeps) are complete.

*timeout* and *wait\_callback* override default constructor parameters.

**write** (*msg, arg=None, arg\_type=None, unit=None, bool\_selector=None, wait\_sync=None, read\_echo=False, read\_echo\_delay=0.0*)

Send a command.

#### Parameters

- **msg** (*str*) – Text message.
- **arg** – Optional argument to append in the end. If a list of arguments is supplied, the result is joined with `" , "`.
- **arg\_type** (*str*) – Argument type. Can be `'raw'` (in which case data is sent raw), `'string'`, `'int'`, `'float'`, `'bool'`, a format string (such as `'{: .3f}'`) or a list of argument types (for an iterable argument); if format string is used and the argument is a list or a tuple, then it is expanded as a list of arguments (e.g., `arg_type='{0};{1}'` with `arg=[1, 2]` will produce a string `'1; 2'`); if a list of types is used, each element of *arg* is converted using the corresponding type, and the result is joined with `" , "`.
- **unit** (*str*) – If not `None`, use it as a unit to append after the value.

- **bool\_selector** (*tuple*) – A tuple (*false\_value*, *true\_value*) of two strings to represent bool argument; by default, use `._bool_selector` attribute.
- **wait\_sync** – if `True`, append the sync command (specified as `._wait_sync_comm` attribute, `"*OPC?"` by default) after the message and pause the execution command is complete; useful in long set operations, where the device might ignore later inputs until the current command is complete; if `None`, use the class default `._default_write_sync` attribute (`False` by default).
- **read\_echo** (*bool*) – If `True`, read a single line after write.
- **read\_echo\_delay** (*float*) – The delay between write and read if `read_echo==True`.

```
class pylablib.devices.Lakeshore.base.TLakeshore370RangeSettings(exc_mode,  
                                                                exc_range,  
                                                                res_range,  
                                                                autorange,  
                                                                enable)
```

Bases: `tuple`

**autorange**

**count** ()

Return number of occurrences of value.

**enable**

**exc\_mode**

**exc\_range**

**index** ()

Return first index of value.

Raises `ValueError` if the value is not present.

**res\_range**

```
class pylablib.devices.Lakeshore.base.TLakeshore370AnalogSettings(bipolar,  
                                                                    mode,  
                                                                    channel,  
                                                                    source,  
                                                                    high_value,  
                                                                    low_value,  
                                                                    man_value)
```

Bases: `tuple`

**bipolar**

**channel**

**count** ()

Return number of occurrences of value.

**high\_value**

**index** ()

Return first index of value.

Raises `ValueError` if the value is not present.

**low\_value**

**man\_value**

**mode****source**

**class** `pylablib.devices.Lakeshore.base.TLakeshore370FilterSettings` (*enabled*,  
*settle\_time*,  
*window*)

Bases: `tuple`**count** ()

Return number of occurrences of value.

**enabled****index** ()

Return first index of value.

Raises `ValueError` if the value is not present.**settle\_time****window**

**class** `pylablib.devices.Lakeshore.base.Lakeshore370` (*conn*)

Bases: `pylablib.core.devio.SCPI.SCPIDevice`

Lakeshore 370 resistance bridge / temperature controller.

All channels are enumerated from 0.

**Parameters** *conn* – serial connection parameters (usually port or a tuple containing port and baudrate)

**Error**alias of `LakeshoreError`**ReraiseError**alias of `LakeshoreBackendError`**get\_temperature** (*channel*)

Get temperature readings (in K) on a given channel

**get\_resistance** (*channel*)

Get resistance readings (in Ohm) on a given channel

**get\_sensor\_power** (*channel*)

Get dissipated power (in W) on a given channel

**select\_channel** (*channel*)

Select measurement channel

**get\_channel** ()

Get current measurement channel

**get\_channel\_range\_settings** (*channel*)

Setup the current measurement channel range parameters.

For parameters, see `setup_channel_range()` and RDGRNG command description in the Lakeshore 370 programming manual.

**setup\_channel\_range** (*channel=None*, *exc\_mode='v'*, *exc\_range=1*, *res\_range=22*, *autorange=True*, *enable=True*)

Setup the measurement channel range (all channels by default).

*exc\_mode* is the excitation mode ("i" or "v"), *exc\_range* is the excitation range (1 is smallest), *res\_range* is the resistance range (1 is smallest). For range descriptions, see Lakeshore 370 programming manual.

**get\_analog\_output\_settings** (*output*)

Get analog output settings for a given output (1 or 2).

For parameters, see [`setup\_analog\_output\(\)`](#) and ANALOG command description in the Lakeshore 370 programming manual.

**setup\_analog\_output** (*output*, *bipolar=None*, *mode=None*, *channel=None*, *source=None*,  
*high\_value=None*, *low\_value=None*, *man\_value=None*)

Setup analog output settings for a given output (1 or 2).

For parameters, see ANALOG command description in the Lakeshore 370 programming manual. Value of None means keeping the current parameter value.

**set\_analog\_output\_value** (*output*, *value*, *bipolar=False*, *enabled=True*)

Set manual analog output value.

A simplified version of [`setup\_analog\_output\(\)`](#).

**get\_analog\_output** (*output*)

Get value (in percents of the total range) at a given output (1 or 2)

**get\_filter\_settings** (*channel*)

Get input filter settings for a given channel (1 to 16).

For parameters, see [`setup\_filter\(\)`](#) and FILTER command description in the Lakeshore 370 programming manual.

**setup\_filter** (*channel*, *enabled=None*, *settle\_time=None*, *window=None*)

Setup input filter settings for a given channel (1 to 16).

For parameters, see FILTER command description in the Lakeshore 370 programming manual. Value of None means keeping the current parameter value.

**BackendError**

alias of [`pylablib.core.devio.comm\_backend.DeviceBackendError`](#)

**class NoParameterCaller** (*device*, *kind*)

Bases: [`object`](#)

Class to simplify calling functions without a parameter

**apply\_settings** (*settings*)

Apply the settings.

*settings* is a dict {name: value} of the available device settings. Non-applicable settings are ignored.

**ask** (*msg*, *data\_type='string'*, *delay=0.0*, *timeout=None*, *read\_echo=False*)

Write a message and read a reply.

*msg* is the query message, *delay* is the delay between write and read. Other parameters are the same as in [`read\(\)`](#). If *read\_echo==True*, assume that the device first echoes the input and skip it.

**close** ()

Close the backend

**flush** (*one\_line=False*)

Flush the read buffer (read all the available data and return the number of bytes read).

If *one\_line==True*, read only a single line.

**static get\_arg\_type** (*arg*)

Autodetect argument type

**get\_device\_variable** (*key*)

Get the value of a settings, status, or full info parameter

**get\_esr** (*timeout=None*)

Get the device status register (by default, "\*ESR?" command)

**get\_full\_info** (*include=0*)

Get dict {name: value} containing full device information (including status and settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_full\_status** (*include=0*)

Get dict {name: value} containing the device status (including settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_id** (*timeout=None*)

Get the device IDN. (query SCPI '\*IDN?' command)

**get\_settings** (*include=0*)

Get dict {name: value} containing all the device settings.

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**is\_opened** ()

Check if the device is connected

**lock** (*timeout=None*)

Lock the access to the device from other threads/processes (isn't necessarily implemented)

**locking** (*timeout=None*)

Context manager for lock & unlock

**open** ()

Open the backend

**static parse\_array\_data** (*data, fmt, include\_header=False*)

Parse the data returned by the device. *fmt* is *DataFormat* description in numpy format (e.g., "<u2").

If *include\_header==True*, the data is assumed to be in a (somewhat) standard SCPI format: b'#', then a single digit *s* denoting length of the size block, then *s* digits denoting length of the data (in bytes) followed by the actual data. Otherwise (*include\_header==False*), assume that the header is already removed.

**read** (*data\_type='string', timeout=None*)

Read data from the device.

*data\_type* determines the type of the data. Can be 'raw' (just raw data), 'string' (with trailing and leading spaces stripped), 'int', 'float', 'bool' (interprets 0 or 'off' as False, anything else as True), 'value' (returns tuple (value, unit), where *value* is float), a callable (return the result of this callable applied to the string value), a dictionary (return the stored value corresponding to the string value, or to the value converted into integer if the string value is not present), or a list of data types (the result is treated as a list of values with the given types separated by commas). *timeout* overrides the default value.

**read\_binary\_array\_data** (*include\_header=False, timeout=None, flush\_term=True*)

Read a binary data in the from the device.

The data assumes the standard binary transfer header consisting of "#" symbol, then a single digit with the size of the length string, then the length string containing the length of the binary data (in bytes). If *include\_header==True*, return the data with the header; otherwise, return only the content. If *flush\_term==True*, flush the following line to skip terminator characters after the binary data, which are added by some devices. *timeout* overrides the default value.

**reconnect** (*new\_instrument=True, ignore\_error=True*)

Remake the connection.

If *new\_instrument==True*, create a new backend instance. If *ignore\_error==True*, ignore errors on closing.

**reset** ()

Reset the device (by default, "\*RST" command)

**set\_device\_variable** (*key, value*)

Set the value of a settings parameter

**sleep** (*delay*)

Wait for *delay* seconds

**unlock** ()

Unlock the access to the device from other threads/processes (isn't necessarily implemented)

**using\_write\_buffer** ()

Context manager for using a write buffer.

While it's active, all the consecutive *write()* operations are bundled together with ; delimiter. The actual write is performed at the *read()/ask()* operation or at the end of the block.

**wait** (*wait\_type='sync', timeout=None, wait\_callback=None*)

Pause execution until device overlapped commands are complete.

*wait\_type* is either 'sync' (perform *wait\_sync()*), 'dev' (perform *wait\_dev()*) or 'none' (do nothing).

**wait\_dev** ()

Pause execution of the device commands until device overlapped commands (e.g., taking sweeps) are complete.

Note that the code execution is not paused.

**wait\_sync** (*timeout=None, wait\_callback=None*)

Pause execution of the script until device overlapped commands (e.g., taking sweeps) are complete.

*timeout* and *wait\_callback* override default constructor parameters.

**write** (*msg, arg=None, arg\_type=None, unit=None, bool\_selector=None, wait\_sync=None, read\_echo=False, read\_echo\_delay=0.0*)

Send a command.

#### Parameters

- **msg** (*str*) – Text message.
- **arg** – Optional argument to append in the end. If a list of arguments is supplied, the result is joined with ", ".
- **arg\_type** (*str*) – Argument type. Can be 'raw' (in which case data is sent raw), 'string', 'int', 'float', 'bool', a format string (such as '{:.3f}'), or a list of argument types (for an iterable argument); if format string is used and

the argument is a list or a tuple, then it is expanded as a list of arguments (e.g., `arg_type='{0};{1}'` with `arg=[1,2]` will produce a string `'1;2'`); if a list of types is used, each element of `arg` is converted using the corresponding type, and the result is joined with `", "`.

- **unit** (*str*) – If not `None`, use it as a unit to append after the value.
- **bool\_selector** (*tuple*) – A tuple (`false_value`, `true_value`) of two strings to represent bool argument; by default, use `._bool_selector` attribute.
- **wait\_sync** – if `True`, append the sync command (specified as `._wait_sync_comm` attribute, `"*OPC?"` by default) after the message and pause the execution command is complete; useful in long set operations, where the device might ignore later inputs until the current command is complete; if `None`, use the class default `._default_write_sync` attribute (`False` by default).
- **read\_echo** (*bool*) – If `True`, read a single line after write.
- **read\_echo\_delay** (*float*) – The delay between write and read if `read_echo==True`.

## Module contents

### pylablib.devices.LaserQuantum package

#### Submodules

#### pylablib.devices.LaserQuantum.base module

**exception** `pylablib.devices.LaserQuantum.base.LaserQuantumError`

Bases: `pylablib.core.devio.base.DeviceError`

Generic Laser Quantum devices error

**args**

**with\_traceback()**

Exception.with\_traceback(tb) – set self.\_\_traceback\_\_ to tb and return self.

**exception** `pylablib.devices.LaserQuantum.base.LaserQuantumBackendError` (*exc*)

Bases: `pylablib.devices.LaserQuantum.base.LaserQuantumError`, `pylablib.core.devio.comm_backend.DeviceBackendError`

Generic Laser Quantum backend communication error

**args**

**with\_traceback()**

Exception.with\_traceback(tb) – set self.\_\_traceback\_\_ to tb and return self.

**class** `pylablib.devices.LaserQuantum.base.TDeviceInfo` (*serial*, *software\_version*, *cal\_date*)

Bases: `tuple`

**cal\_date**

**count()**

Return number of occurrences of value.

```
index()
    Return first index of value.

    Raises ValueError if the value is not present.

serial

software_version

class pylablib.devices.LaserQuantum.base.TWorkHours (psu, laser_enabled,
                                                    laser_threshold)
    Bases: tuple

count()
    Return number of occurrences of value.

index()
    Return first index of value.

    Raises ValueError if the value is not present.

laser_enabled

laser_threshold

psu

class pylablib.devices.LaserQuantum.base.TTemperatures (head, psu)
    Bases: tuple

count()
    Return number of occurrences of value.

head

index()
    Return first index of value.

    Raises ValueError if the value is not present.

psu

class pylablib.devices.LaserQuantum.base.Finesse (conn)
    Bases: pylablib.core.devio.comm_backend.ICommBackendWrapper

    Laser Quantum Finesse pump laser.

    Parameters conn – serial connection parameters (usually port)

Error
    alias of LaserQuantumError

query (comm, reply_lines=1)
    Send a query to the device and read the reply.

    reply_lines specify the number of lines to read as a reply (almost all queries have only one line).

get_device_info()
    Get device information (serial, software_version, cal_date)

get_work_hours()
    Get the work hours (PSU run time, laser run time, laser above threshold time)

get_temperatures()
    Get device status, head temperature, and PSU temperature
```



**get\_output\_status()**  
Get output status.  
  
Can be "enabled" or "disabled".

**get\_interlock\_status()**  
Get manual interlock status

**get\_shutter\_status()**  
Get the shutter status

**is\_shutter\_opened()**  
Check if shutter is opened

**set\_shutter(*opened=True*)**  
Open or close the shutter

**is\_enabled()**  
Check if the output is enabled

**enable(*enabled=True*)**  
Turn the output on or off

**get\_output\_power()**  
Get the output power (in Watts)

**get\_output\_setpoint()**  
Get the output setpoint power (in Watts)

**set\_output\_power(*level*)**  
Set the output power setpoint (in Watts)

**get\_current()**  
Get the laser drive current (in %)

**class NoParameterCaller(*device, kind*)**  
Bases: `object`  
  
Class to simplify calling functions without a parameter

**apply\_settings(*settings*)**  
Apply the settings.  
  
*settings* is the dict {*name*: *value*} of the device available settings. Non-applicable settings are ignored.

**close()**  
Close the backend

**get\_device\_variable(*key*)**  
Get the value of a settings, status, or full info parameter

**get\_full\_info(*include=0*)**  
Get dict {*name*: *value*} containing full device information (including status and settings).  
  
*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_full\_status(*include=0*)**  
Get dict {*name*: *value*} containing the device status (including settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting `include=-10` queries all available variables, which is equivalent to `include="all"`.

**get\_settings** (*include=0*)

Get dict {name: value} containing all the device settings.

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting `include=-10` queries all available variables, which is equivalent to `include="all"`.

**is\_opened** ()

Check if the device is connected

**lock** (*timeout=None*)

Lock the access to the device from other threads/processes (isn't necessarily implemented)

**locking** (*timeout=None*)

Context manager for lock & unlock

**open** ()

Open the backend

**set\_device\_variable** (*key, value*)

Set the value of a settings parameter

**unlock** ()

Unlock the access to the device from other threads/processes (isn't necessarily implemented)

## Module contents

### pylablib.devices.Leybold package

#### Submodules

#### pylablib.devices.Leybold.base module

**exception** pylablib.devices.Leybold.base.LeyboldError

Bases: `pylablib.core.devio.base.DeviceError`

Generic Leybold device error

**args**

**with\_traceback** ()

Exception.with\_traceback(tb) – set self.\_\_traceback\_\_ to tb and return self.

**exception** pylablib.devices.Leybold.base.LeyboldBackendError (*exc*)

Bases: `pylablib.devices.Leybold.base.LeyboldError`, `pylablib.core.devio.comm_backend.DeviceBackendError`

Generic Leybold backend communication error

**args**

**with\_traceback** ()

Exception.with\_traceback(tb) – set self.\_\_traceback\_\_ to tb and return self.

```
class pylablib.devices.Leybold.base.TDeviceInfo (sensor, page, swver)
    Bases: tuple

    count ()
        Return number of occurrences of value.

    index ()
        Return first index of value.

        Raises ValueError if the value is not present.

    page

    sensor

    swver

class pylablib.devices.Leybold.base.TUpdateValue (value, display_units, status, error, de-
                                                    vice_info)
    Bases: tuple

    count ()
        Return number of occurrences of value.

    device_info

    display_units

    error

    index ()
        Return first index of value.

        Raises ValueError if the value is not present.

    status

    value

class pylablib.devices.Leybold.base.GenericITR (conn)
    Bases: pylablib.core.devio.comm_backend.ICommBackendWrapper

    Generic Leybold ITR pressure gauge.

        Parameters conn – serial connection parameters (usually port or a tuple containing port and bau-
            drate)

    Error
        alias of LeyboldError

    get_update (refresh=True)
        Get device state update.

        Return tuple (value, display_units, status, error, device_info), where value
        is the pressure in Pa, display_units are display units ("pa", "mbar", or "torr"), status is the
        devices status (e.g., emission status), error is the device error ("ok" if no errors), and device_info
        is a tuple (sensor, page, swver) with the sensor kind ID, data page, and software version.

        If refresh==True, get the latest update value; otherwise, get the latest read value.

    send_command (byte1, byte2, byte3)
        Send command to the device.

        Arguments represent the three command bytes. Values of these bytes for different commands are de-
        scribed in the manual.
```

**get\_device\_info()**

Get device info.

Return tuple (sensor, page, swver) with the sensor kind ID, data page, and software version.

**get\_units()**

Get device readout units ("mbar", "pa", or "torr")

**get\_pressure** (display\_units=False)

Get pressure.

If display\_units==False, return result in Pa; otherwise, use display units obtained using `get_units()`.

**class NoParameterCaller** (device, kind)

Bases: `object`

Class to simplify calling functions without a parameter

**apply\_settings** (settings)

Apply the settings.

*settings* is the dict {name: value} of the device available settings. Non-applicable settings are ignored.

**close()**

Close the backend

**get\_device\_variable** (key)

Get the value of a settings, status, or full info parameter

**get\_full\_info** (include=0)

Get dict {name: value} containing full device information (including status and settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include*=-10 queries all available variables, which is equivalent to *include*="all".

**get\_full\_status** (include=0)

Get dict {name: value} containing the device status (including settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include*=-10 queries all available variables, which is equivalent to *include*="all".

**get\_settings** (include=0)

Get dict {name: value} containing all the device settings.

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include*=-10 queries all available variables, which is equivalent to *include*="all".

**is\_opened()**

Check if the device is connected

**lock** (timeout=None)

Lock the access to the device from other threads/processes (isn't necessarily implemented)

**locking** (timeout=None)

Context manager for lock & unlock

```

open ()
    Open the backend

set_device_variable (key, value)
    Set the value of a settings parameter

unlock ()
    Unlock the access to the device from other threads/processes (isn't necessarily implemented)

class pylablib.devices.Leybold.base.TITR90Status (emission, atm_adj)
    Bases: tuple

    atm_adj

    count ()
        Return number of occurrences of value.

    emission

    index ()
        Return first index of value.

        Raises ValueError if the value is not present.

class pylablib.devices.Leybold.base.ITR90 (conn)
    Bases: pylablib.devices.Leybold.base.GenericITR

    Leybold ITR90 pressure gauge.

    Parameters conn – serial connection parameters (usually port or a tuple containing port and baudrate)

    set_units (units, store=True)
        Get device readout units ("mbar", "pa", or "torr").

        If store==True, store the value in the non-volatile power-independent memory.

    start_degas ()
        Start degas (turns off automatically after 3 minutes)

    stop_degas ()
        Stop degas

    Error
        alias of LeyboldError

class NoParameterCaller (device, kind)
    Bases: object

    Class to simplify calling functions without a parameter

    apply_settings (settings)
        Apply the settings.

        settings is the dict {name: value} of the device available settings. Non-applicable settings are ignored.

    close ()
        Close the backend

    get_device_info ()
        Get device info.

        Return tuple (sensor, page, swver) with the sensor kind ID, data page, and software version.

```

**get\_device\_variable** (*key*)

Get the value of a settings, status, or full info parameter

**get\_full\_info** (*include=0*)

Get dict {name: value} containing full device information (including status and settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_full\_status** (*include=0*)

Get dict {name: value} containing the device status (including settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_pressure** (*display\_units=False*)

Get pressure.

If *display\_units==False*, return result in Pa; otherwise, use display units obtained using *get\_units()*.

**get\_settings** (*include=0*)

Get dict {name: value} containing all the device settings.

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_units** ()

Get device readout units ("mbar", "pa", or "torr")

**get\_update** (*refresh=True*)

Get device state update.

Return tuple (value, display\_units, status, error, device\_info), where value is the pressure in Pa, display\_units are display units ("pa", "mbar", or "torr"), status is the devices status (e.g., emission status), error is the device error ("ok" if no errors), and device\_info is a tuple (sensor, page, swver) with the sensor kind ID, data page, and software version.

If *refresh==True*, get the latest update value; otherwise, get the latest read value.

**is\_opened** ()

Check if the device is connected

**lock** (*timeout=None*)

Lock the access to the device from other threads/processes (isn't necessarily implemented)

**locking** (*timeout=None*)

Context manager for lock & unlock

**open** ()

Open the backend

**send\_command** (*byte1, byte2, byte3*)

Send command to the device.

Arguments represent the three command bytes. Values of these bytes for different commands are described in the manual.

**set\_device\_variable** (*key, value*)

Set the value of a settings parameter

**unlock** ()

Unlock the access to the device from other threads/processes (isn't necessarily implemented)

## Module contents

### pylablib.devices.LighthousePhotonics package

#### Submodules

#### pylablib.devices.LighthousePhotonics.base module

**exception** pylablib.devices.LighthousePhotonics.base.**LighthousePhotonicsError**

Bases: [pylablib.core.devio.base.DeviceError](#)

Generic Lighthouse Photonics devices error

**args**

**with\_traceback** ()

Exception.with\_traceback(tb) – set self.\_\_traceback\_\_ to tb and return self.

**exception** pylablib.devices.LighthousePhotonics.base.**LighthousePhotonicsBackendError** (*exc*)

Bases: [pylablib.devices.LighthousePhotonics.base.LighthousePhotonicsError](#),  
[pylablib.core.devio.comm\\_backend.DeviceBackendError](#)

Generic Lighthouse Photonics backend communication error

**args**

**with\_traceback** ()

Exception.with\_traceback(tb) – set self.\_\_traceback\_\_ to tb and return self.

**class** pylablib.devices.LighthousePhotonics.base.**TDeviceInfo** (*product, version, serial, configuration*)

Bases: [tuple](#)

**configuration**

**count** ()

Return number of occurrences of value.

**index** ()

Return first index of value.

Raises ValueError if the value is not present.

**product**

**serial**

**version**

**class** pylablib.devices.LighthousePhotonics.base.**TWorkHours** (*controller, laser*)

Bases: [tuple](#)

**controller**

**count** ()

Return number of occurrences of value.

**index()**

Return first index of value.

Raises ValueError if the value is not present.

**laser**

**class** `pylablib.devices.LighthousePhotonics.base.SproutG(conn)`

Bases: `pylablib.core.devio.comm_backend.ICommBackendWrapper`

Lighthouse Photonics Sprout G laser.

**Parameters** `conn` – serial connection parameters (usually port)

**Error**

alias of `LighthousePhotonicsError`

**query** (`comm`, `allowed_replies=('0',)`)

Send a query to the device and parse the reply

**get\_device\_info()**

Get device information (product name, product version, serial number, configuration)

**get\_work\_hours()**

Return device operation hours (controller on) and run hours (laser on)

**get\_warning\_status()**

Get device warnings

**get\_interlock\_status()**

Get manual interlock status

**get\_shutter\_status()**

Get manual shutter status ("open" or "close")

**get\_output\_mode()**

Get output mode.

Can be "on", "off", "idle" (power standby mode), "calibrate", "interlock" (manual interlock is off), "warmup" (warmup mode), or "calibration" (calibration mode).

**set\_output\_mode** (`mode='on'`)

Set output mode.

`mode` can be "on", "off", "idle" (power standby mode), or "calibrate" (calibration mode).

**is\_enabled()**

Check if the output is on (idle or warmup don't count as on)

**enable** (`enabled=True`)

Turn the output on or off

**get\_output\_power()**

Set the actual output power (in Watts)

**get\_output\_setpoint()**

Get the output setpoint power (in Watts)

**set\_output\_power** (`level`)

Get the output power setpoint (in Watts)

**class** `NoParameterCaller(device, kind)`

Bases: `object`

Class to simplify calling functions without a parameter



**apply\_settings** (*settings*)

Apply the settings.

*settings* is the dict {name: value} of the device available settings. Non-applicable settings are ignored.

**close** ()

Close the backend

**get\_device\_variable** (*key*)

Get the value of a settings, status, or full info parameter

**get\_full\_info** (*include=0*)

Get dict {name: value} containing full device information (including status and settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_full\_status** (*include=0*)

Get dict {name: value} containing the device status (including settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_settings** (*include=0*)

Get dict {name: value} containing all the device settings.

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**is\_opened** ()

Check if the device is connected

**lock** (*timeout=None*)

Lock the access to the device from other threads/processes (isn't necessarily implemented)

**locking** (*timeout=None*)

Context manager for lock & unlock

**open** ()

Open the backend

**set\_device\_variable** (*key, value*)

Set the value of a settings parameter

**unlock** ()

Unlock the access to the device from other threads/processes (isn't necessarily implemented)

## Module contents

### pylablib.devices.M2 package

### Submodules

**pylablib.devices.M2.solstis module****exception** `pylablib.devices.M2.solstis.M2Error`Bases: `pylablib.core.devio.base.DeviceError`

Generic M2 error

**args****with\_traceback()**

Exception.with\_traceback(tb) – set self.\_\_traceback\_\_ to tb and return self.

**exception** `pylablib.devices.M2.solstis.M2CommunicationError` (*exc*)Bases: `pylablib.devices.M2.solstis.M2Error`, `pylablib.core.devio.comm_backend.DeviceBackendError`

M2 network communication error

**args****with\_traceback()**

Exception.with\_traceback(tb) – set self.\_\_traceback\_\_ to tb and return self.

**class** `pylablib.devices.M2.solstis.Solstis`(*addr*, *port*, *timeout=5.0*, *start\_link=True*,  
*use\_websocket='auto'*, *use\_cavity=True*)Bases: `pylablib.core.devio.interface.IDevice`

M2 Solstis Ice Bloc device.

**Parameters**

- **addr** (*str*) – IP address of the Ice Bloc device.
- **port** (*int*) – port of the Ice Bloc device.
- **timeout** (*float*) – default timeout of synchronous operations.
- **start\_link** (*bool*) – if True, initialize device link on creation.
- **use\_websocket** (*bool*) – if True, use websocket interface (same as used by the web interface) for additional functionality (wavemeter connection, etalon value, improved operation stopping); "auto" enables it if websocket package is installed, and disables otherwise
- **use\_cavity** – if False and any reference cavity methods are used, either ignore them, or use closest available methods instead

**Error**alias of `M2Error`**ReraiseError**alias of `M2CommunicationError`**BackendError**alias of `builtins.OSError`**open()**

Open the connection

**close()**

Close the connection

**is\_opened()**

Check if the device is connected

**set\_timeout** (*timeout*)  
Set timeout for connecting or sending/receiving

**flush** ()  
Flush read buffer

**query** (*op, params, reply\_op='auto', report=False*)  
Send a query using the standard device interface.  
  
*reply\_op* is the name of the reply operation (by default, its the operation name plus "\_reply"). If *report==True*, request completion report (does not apply to all operation).

**update\_reports** (*timeout=0.0*)  
Check for fresh operation reports

**get\_last\_report** (*op*)  
Get the latest report for the given operation

**check\_report** (*op*)  
Check and return the latest report for the given operation

**wait\_for\_report** (*op, error\_msg=None, timeout=None*)  
Wait for a report for the given operation  
  
*error\_msg* specifies the exception message if the report results in an error.

**start\_link** ()  
Initialize device link (called automatically on creation)

**connect\_wavemeter** (*sync=True*)  
Connect to the wavemeter (if *sync==True*, wait until the connection is established)

**disconnect\_wavemeter** (*sync=True*)  
Disconnect from the wavemeter (if *sync==True*, wait until the connection is broken)

**is\_wavemeter\_connected** ()  
Check if the wavemeter is connected

**get\_system\_status** ()  
Get the device system status

**get\_full\_web\_status** ()  
Get full websocket status.  
  
Return a large dictionary containing all the information available in the web interface.

**get\_full\_fine\_tuning\_status** ()  
Get full fine-tuning status (see M2 Solstis JSON protocol manual for "poll\_wave\_m" command)

**lock\_wavemeter** (*lock=True, sync=True, error\_on\_fail=True*)  
Lock or unlock the laser to the wavemeter (if *sync==True*, wait until the operation is complete)

**is\_wavemeter\_lock\_on** ()  
Check if the laser is locked to the wavemeter

**fine\_tune\_wavelength** (*wavelength, sync=True, timeout=None*)  
Fine-tune the wavelength.  
  
Only works if the wavemeter is connected. If *sync==True*, wait until the operation is complete (might take from several seconds up to several minutes).

**check\_fine\_tuning\_report** ()  
Check wavelength fine-tuning report  
  
Return "success" or "fail" if the operation is complete, or None if it is still in progress.

**wait\_for\_fine\_tuning** (*timeout=None*)

Wait until wavelength fine-tuning is complete

**get\_fine\_tuning\_status** ()

Get fine-tuning status.

Return either "idle" (no tuning or locking), "nolink" (no wavemeter link), "tuning" (tuning in progress), or "locked" (tuned and locked to the wavemeter).

**get\_fine\_wavelength** ()

Get fine-tuned wavelength.

Only works if the wavemeter is connected.

**stop\_fine\_tuning** ()

Stop fine wavelength tuning

**coarse\_tune\_wavelength** (*wavelength, sync=True*)

Coarse-tune the wavelength.

Only works if the wavemeter is disconnected. If *sync==True*, wait until the operation is complete.

**get\_full\_coarse\_tuning\_status** ()

Get full coarse-tuning status (see M2 M2 Solstis JSON protocol manual for "poll\_move\_wave\_t" command)

**get\_coarse\_tuning\_status** ()

Get coarse-tuning status.

Return either "done" (tuning is done), "tuning" (tuning in progress), or "fail" (tuning failed).

**get\_coarse\_wavelength** ()

Get course-tuned wavelength.

Only works if the wavemeter is disconnected.

**stop\_coarse\_tuning** ()

Stop coarse wavelength tuning

**tune\_etalon** (*value, sync=True*)

Tune the etalon to *value* percent.

Only works if the wavemeter is disconnected. If *sync==True*, wait until the operation is complete.

**lock\_etalon** (*sync=True*)

Lock the etalon.

If *sync==True*, wait until the operation is complete.

**unlock\_etalon** (*sync=True*)

Unlock the etalon .

If *sync==True*, wait until the operation is complete. Automatically unlock the reference cavity first (otherwise the operation fails).

**get\_etalon\_lock\_status** ()

Get etalon lock status.

Return either "off" (lock is off), "on" (lock is on), "debug" (lock in debug condition), "error" (lock had an error), "search" (lock is searching), or "low" (lock is off due to low output).

**tune\_laser\_resonator** (*value, fine=False, sync=True*)

Tune the laser cavity to *value* percent.

If `fine==True`, adjust fine tuning; otherwise, adjust coarse tuning. Only works if the wavemeter is disconnected. If `sync==True`, wait until the operation is complete.

**tune\_reference\_cavity** (*value*, *fine=False*, *sync=True*)

Tune the reference cavity to *value* percent.

If `fine==True`, adjust fine tuning; otherwise, adjust coarse tuning. Only works if the wavemeter is disconnected. If `sync==True`, wait until the operation is complete. If reference cavity is disabled by setting `use_cavity=False` on creation, do nothing.

**lock\_reference\_cavity** (*sync=True*)

Lock the laser to the reference cavity.

Automatically lock etalon first (otherwise the operation fails). If `sync==True`, wait until the operation is complete. If reference cavity is disabled by setting `use_cavity=False` on creation, do nothing.

**unlock\_reference\_cavity** (*sync=True*)

Unlock the laser from the reference cavity.

If `sync==True`, wait until the operation is complete. If reference cavity is disabled by setting `use_cavity=False` on creation, do nothing.

**get\_reference\_cavity\_lock\_status** ()

Get the reference cavity lock status.

Return either "off" (lock is off), "on" (lock is on), "debug" (lock in debug condition), "error" (lock had an error), "search" (lock is searching), "low" (lock is off due to low output), or "disabled" (reference cavity is disabled by setting `use_cavity=False` on creation).

**setup\_terascan** (*scan\_type*, *scan\_range*, *rate*, *trunc\_rate=True*)

Setup terascan.

#### Parameters

- **scan\_type** (*str*) – scan type. Can be "medium" (BRF+etalon, rate from 100 GHz/s to 1 GHz/s), "fine" (all elements, rate from 20 GHz/s to 1 MHz/s), or "line" (all elements, rate from 20 GHz/s to 50 kHz/s).
- **scan\_range** (*tuple*) – tuple (*start*, *stop*) with the scan range (in Hz).
- **rate** (*float*) – scan rate (in Hz/s).
- **trunc\_rate** (*bool*) – if `True`, truncate the scan rate to the nearest available rate (otherwise, incorrect rate would raise an error).

If reference cavity is disabled by setting `use_cavity=False` on creation and *scan\_type* is "line", use "fine" instead.

**start\_terascan** (*scan\_type*, *sync=False*, *sync\_done=False*)

Start terascan.

Scan parameters are set up separately using `setup_terascan()`. Scan type can be "medium" (BRF+etalon, rate from 100 GHz/s to 1 GHz/s), "fine" (all elements, rate from 20 GHz/s to 1 MHz/s), or "line" (all elements, rate from 20 GHz/s to 50 kHz/s). If reference cavity is disabled by setting `use_cavity=False` on creation and *scan\_type* is "line", use "fine" instead. If `sync==True`, wait until the scan is set up (not until the whole scan is complete). If `sync_done==True`, wait until the whole scan is complete (not recommended, as it can take hours).

**enable\_terascan\_updates** (*enable=True*, *update\_period=0*)

Enable sending periodic terascan updates.

If enabled, laser will send updates in the beginning and in the end of every terascan segment. If `update_period!=0`, it will also send updates every `update_period` percents of the segment (this option is not currently supported by M2 firmware).

**check\_terascan\_update()**

Check the latest terascan update.

Return `None` if none are available, or a dictionary `{"wavelength":current_wavelength, "operation":op}`, where `op` is "scanning" (scanning in progress), "stitching" (stitching in progress), "finished" (scan is finished), or "repeat" (segment is repeated).

**wait\_for\_terascan\_update()**

Wait until a new terascan update is available

**check\_terascan\_start\_report()**

Check report on terascan start.

Return "success" or "fail" if the operation is complete, or `None` if it is still in progress.

**stop\_terascan(scan\_type, sync=False)**

Stop terascan of the given type.

If reference cavity is disabled by setting `use_cavity=False` on creation and `scan_type` is "line", use "fine" instead. If `sync==True`, wait until the operation is complete.

**get\_terascan\_status(scan\_type, web\_status=True)**

Get status of a terascan of a given type.

**Return a dictionary with 4 items:** "current": current laser frequency (or `None` if no scan is in progress) "range": tuple with the fill scan range (or `None` if no frequency is available) "status": can be "stopped" (scan is not in progress), "scanning" (scan is in progress), or "stitching" (scan is in progress, but currently stitching) "web": whether scan is running in web interface (some failure modes still report "scanning" through the usual interface); only available if the laser web connection is on and if `web_status==True`.

If reference cavity is disabled by setting `use_cavity=False` on creation and `scan_type` is "line", use "fine" instead.

**start\_fast\_scan(scan\_type, width, period, sync=False, setup\_locks=True)**

Setup and start fast scan.

#### Parameters

- **scan\_type** (*str*) – scan type. Can be "cavity\_continuous", "cavity\_single", "cavity\_triangular", "etalon\_continuous", "etalon\_single", "resonator\_continuous", "resonator\_single", "resonator\_ramp", "resonator\_triangular", "ecd\_continuous", "ecd\_ramp", or "fringe\_test" (see M2 Solstis JSON protocol manual for details)
- **width** (*float*) – scan width (in Hz).
- **period** (*float*) – scan time/period (in s).
- **sync** (*bool*) – if `True`, wait until the scan is set up (not until the whole scan is complete).
- **setup\_locks** (*bool*) – if `True`, automatically setup etalon and reference cavity locks in the appropriate states for etalon, cavity, or resonator scans.

If reference cavity is disabled by setting `use_cavity=False` on creation, use resonator scans instead of cavity scans.

**check\_fast\_scan\_start\_report ()**

Check fast scan start report.

Return "success" or "fail" if the operation is complete, or None if it is still in progress.

**stop\_fast\_scan** (*scan\_type*, *return\_to\_start=True*, *sync=False*)

Stop fast scan of the given type.

If reference cavity is disabled by setting `use_cavity=False` on creation, use resonator scans instead of cavity scans. If `return_to_start==True`, return to the center frequency after stopping; otherwise, stay at the current instantaneous frequency. If `sync==True`, wait until the operation is complete.

**get\_fast\_scan\_status** (*scan\_type*)

Get status of a fast scan of a given type.

**Return dictionary with 2 items:** "status": can be "stopped" (scan is not in progress), "scanning" (scan is in progress). "value": current tuner value (in percent); does not necessary correspond to the scan progress.

If reference cavity is disabled by setting `use_cavity=False` on creation, use resonator scans instead of cavity scans.

**stop\_scan\_web** (*scan\_type*)

Stop scan of the current type (terascan or fine scan) using web interface.

More reliable than native programming interface, but requires activated web interface. If reference cavity is disabled by setting `use_cavity=False` on creation, use resonator scans instead of cavity scans.

**stop\_all\_operation** (*repeated=True*, *attempt=0*)

Stop all laser operations (tuning and scanning).

More reliable than native programming interface, but requires activated web interface. If `repeated==True`, repeat trying to stop the operations until succeeded (more reliable, but takes more time). If `attempt>0`, it can supply the number of already tried attempts to stop (with `repeated=False`); the more attempts failed, the more drastic measures will be taken to stop (e.g., initialize short terascan or a fast scan, cycle wavemeter connection, etc.) Return `True` if the operation is success and `False` otherwise.

**class NoParameterCaller** (*device*, *kind*)

Bases: `object`

Class to simplify calling functions without a parameter

**apply\_settings** (*settings*)

Apply the settings.

*settings* is the dict {name: value} of the device available settings. Non-applicable settings are ignored.

**get\_device\_variable** (*key*)

Get the value of a settings, status, or full info parameter

**get\_full\_info** (*include=0*)

Get dict {name: value} containing full device information (including status and settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting `include=-10` queries all available variables, which is equivalent to `include="all"`.

**get\_full\_status** (*include=0*)

Get dict {name: value} containing the device status (including settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting `include=-10` queries all available variables, which is equivalent to `include="all"`.

**get\_settings** (*include=0*)

Get dict {name: value} containing all the device settings.

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting `include=-10` queries all available variables, which is equivalent to `include="all"`.

**set\_device\_variable** (*key, value*)

Set the value of a settings parameter

## Module contents

### pylablib.devices.NI package

#### Submodules

#### pylablib.devices.NI.daq module

**exception** `pylablib.devices.NI.daq.NIError`

Bases: `pylablib.core.devio.base.DeviceError`

Generic NI error

**args**

**with\_traceback** ()

Exception.with\_traceback(tb) – set self.\_\_traceback\_\_ to tb and return self.

**exception** `pylablib.devices.NI.daq.NIDAQmxError` (*exc*)

Bases: `pylablib.devices.NI.daq.NIError`, `pylablib.core.devio.comm_backend.DeviceBackendError`

NI DAQmx backend operation error

**args**

**with\_traceback** ()

Exception.with\_traceback(tb) – set self.\_\_traceback\_\_ to tb and return self.

**class** `pylablib.devices.NI.daq.TDeviceInfo` (*name, model, serial\_number*)

Bases: `tuple`

**count** ()

Return number of occurrences of value.

**index** ()

Return first index of value.

Raises `ValueError` if the value is not present.

**model**

**name**



**serial\_number**

```
class pylablib.devices.NI.daq.TVoltageOutputClockParameters (rate, sync_with_ai,
                                                             continuous,
                                                             samps_per_chan,
                                                             autoloop)
```

Bases: `tuple`

**autoloop**

**continuous**

**count()**

Return number of occurrences of value.

**index()**

Return first index of value.

Raises ValueError if the value is not present.

**rate**

**samps\_per\_chan**

**sync\_with\_ai**

```
pylablib.devices.NI.daq.get_device_info (name)
```

Get device info.

Return tuple (name, model, serial).

```
pylablib.devices.NI.daq.list_devices ()
```

List all connected NI DAQ devices

```
class pylablib.devices.NI.daq.NIDAQ (dev_name='dev0', rate=100.0, buffer_size=100000.0,
                                     reset=False)
```

Bases: `pylablib.core.devio.interface.IDevice`

National Instruments DAQ device interface (wrapper around nidaqmx library).

Simplified interface to NI DAQ devices. Supports voltage, digital, and counter inputs (all synchronized to the same clock), and digital and voltage outputs (asynchronous).

#### Parameters

- **dev\_name** (*str*) – root device name.
- **rate** (*float*) – analog input sampling rate (can be adjusted later).
- **buffer\_size** (*int*) – size of the input buffer.
- **reset** (*int*) – if True, reset the device upon connection.

#### Error

alias of `NIError`

#### ReraiseError

alias of `NIDAQmxError`

#### BackendError

Used by autodoc\_mock\_imports.

**open()**

Open the connection

**close()**

Close the connection

**is\_opened()**

Check if the device is connected

**reset()**

Reset the device. All channels will be removed

**get\_device\_info()**

Get device info.

Return tuple (name, model, serial).

**setup\_clock(rate, src=None)**

Setup analog input clock (which is the main system clock).

If `src==None`, use internal clock with the given rate; otherwise use `src` terminal as a clock source (in this case, `rate` should be higher than the expected source rate).

**get\_clock\_parameters()**

Get analog input clock configuration.

Return tuple (rate, src).

**export\_clock(terminal)**

Export system clock to the given terminal (None to disconnect all terminals)

Only terminal one can be active at a time.

**get\_export\_clock\_terminal()**

Return terminal which outputs system clock (None if none is connected)

**add\_voltage\_input(name, channel, rng=(-10, 10), terminal\_cfg='default')**

Add analog voltage input.

Readout is synchronized to the system clock.

#### Parameters

- **name** (*str*) – channel name to refer to it later.
- **channel** (*str*) – terminal name (e.g., "ai0").
- **rng** – voltage range
- **terminal\_cfg** – terminal configuration; can be "default", "rse" (single-ended, referenced to AI SENSE input), "nrse" (single-ended, referenced to AI GND), "diff" (differential), or "pseudodiff" (see NI DAQ manual for details).

**add\_counter\_input(name, counter, terminal, clk\_src='ai/SampleClock', output\_format='rate')**

Add counter input (value is related to the number of counts).

Readout is synchronized to the system clock.

#### Parameters

- **name** (*str*) – channel name.
- **counter** (*str*) – on-board counter name (e.g., "ctr0").
- **terminal** (*str*) – terminal name (e.g., "pfi0").
- **clk\_src** (*str*) – source of the counter sampling clock. By default it is the analog input clock, which requires at least one voltage input channel (could be a dummy channel) to be set up first.
- **output\_format** (*str*) – output format. Can be "acc" (return accumulated number of counts since the sampling start), "diff" (return number of counts passed

between the two consecutive sampling points; essentially, a derivative of "acc"), or "rate" (return count rate based on the "diff" samples).

**add\_clock\_period\_input** (*counter*, *clk\_src*='ai/SampleClock')

Add clock period counter.

Useful when using external sample clock with unknown period. The clock input can be returned during *read()* operation, and it is used to calculate counter inputs in "rate" mode. Readout is synchronized to the system clock.

#### Parameters

- **counter** (*str*) – on-board counter name (e.g., "ctr0") to be used for clock measure.
- **clk\_src** (*str*) – source of the counter sampling clock. By default it is the analog input clock, which requires at least one voltage input channel (could be dummy channel) to operate.

**add\_digital\_input** (*name*, *channel*)

Add digital input.

Readout is synchronized to the system clock. :param name: channel name. :type name: str :param channel: terminal name (e.g., "port0/line12"). :type channel: str

**get\_input\_channels** (*include*=('ai', 'ci', 'di'))

Get names of all input channels (voltage input and counter input).

*include* specifies which channel types to include into the list ("ai" for voltage inputs, "ci" for counter inputs, "di" for digital inputs, "cpi" for clock period channel). The channels order is always fixed: first voltage inputs, then counter inputs, then digital inputs.

**get\_voltage\_input\_parameters** ()

Get parameters (names, channels, output ranges, and terminal configurations) of all analog voltage input channels

**get\_counter\_input\_parameters** ()

Get parameters (names, counters, terminals, clock sources, and output formats) of all counter input channels

**get\_digital\_input\_parameters** ()

Get parameters (names and channels) of all digital input channels

**get\_clock\_period\_input\_parameters** ()

Get parameters (counter input) of the clock period input channel

**start** (*flush\_read*=0, *finite*=None)

Start the sampling and output task.

*flush\_read* specifies number of samples to read and discard after start. If *finite* is not None, it specifies finite number of sample to acquire before stopping.

If counter channels are used, the first sample is usually unreliable, so *flush\_read*=1 is recommended; however, if exactly *finite* pulses are required at the clock export channel, *flush\_read*=0 is needed (the total number of pulses is *flush\_read*+*finite*).

**stop** ()

Stop the sampling task

**is\_running** ()

Check if the task is running

**available\_samples()**

Get number of available samples to read (return 0 if the task is not running)

**get\_buffer\_size()**

Get the sampling buffer size

**wait\_for\_sample**(*num=1, timeout=10.0, wait\_time=0.001*)

Wait until at least *num* samples are available.

If they are not available immediately, loop while checking every *wait\_time* interval until enough samples are accumulated. Return the number of available samples if successful, or 0 if the execution timed out.

**read**(*n=1, flush\_read=0, timeout=10.0, include=('ai', 'ci', 'di')*)

Read *n* samples. If the task is not running, automatically start before reading and stop after.

#### Parameters

- **n** (*int*) – number of samples to read. If  $n \leq 0$ , read all available samples.
- **flush\_read** (*int*) – number of initial samples to skip if the task is currently stopped and needs to be started. If counter channels are used, the first sample is usually unreliable, so *flush\_read*=1 is recommended; however, if exactly *n* pulses are required at the clock export channel, *flush\_read*=0 is needed.
- **include** (*tuple*) – specifies which channel types to include into the list ("ai" for voltage inputs, "ci" for counter inputs, "di" for digital inputs, "cpi" for clock period channel).

**Returns** 2D numpy array of values arranged according to *get\_input\_channels()* order with the given *include* parameter.

**add\_digital\_output**(*name, channel*)

Add digital output.

#### Parameters

- **name** (*str*) – channel name.
- **channel** (*str*) – terminal name (e.g., "do0").

**get\_digital\_output\_channels()**

Get names of all digital output channels

**get\_digital\_output\_parameters()**

Get parameters (names and channels) of all digital output channels

**set\_digital\_outputs**(*names, values*)

Set values of one or several digital outputs.

#### Parameters

- **names** (*str* or [*str*]) – name or list of names of outputs.
- **values** – output value or list of values.

**get\_digital\_outputs**(*names=None*)

Get values of one or several digital outputs.

**Parameters** **names** (*str* or [*str*] or *None*) – name or list of names of outputs (None means all outputs).

Return list of values ordered by *names* (or by *get\_digital\_output\_channels()* if *names*==None).

**add\_voltage\_output** (*name, channel, rng=(-10, 10), initial\_value=0.0*)

Add analog voltage output.

#### Parameters

- **name** (*str*) – channel name.
- **channel** (*str*) – terminal name (e.g., "a00").
- **rng** – voltage range.
- **initial\_value** (*float*) – initial output value (has to be initialized).

**get\_voltage\_output\_channels** ()

Get names of all analog voltage output channels

**get\_voltage\_output\_parameters** ()

Get parameters (names, channels and output ranges) of all analog voltage output channels

**set\_voltage\_outputs** (*names, values, minsamp=1, force\_restart=True, single\_shot=0*)

Set values of one or several analog voltage outputs.

#### Parameters

- **names** (*str or [str]*) – name or list of names of outputs.
- **values** – output value or list values. These can be single numbers, or arrays if the output clock is setup (see `setup_voltage_output_clock()`). In the latter case it sets up the output waveforms; note that waveforms for all channels must have the same length (a single number signifying a constant output is also allowed) If the analog output is set up to the finite mode (`continuous==False`), the finite waveform output happens right away, with the number of samples determined by `samps_per_channel` parameter of `setup_voltage_output_clock()`. In this case, if the supplied waveform is shorter than the number of samples, it gets repeated; if it's longer, it gets cut off.
- **minsamp** – in non-autoloop mode, specifies the minimal number of samples to write to the output buffer; if the length of *values* is less than this number, than the waveform is repeated by a required integer number of times to produce at least *minsamp* samples
- **force\_restart** – if `True`, restart the output after writing to immediately start outputting the new waveforms; otherwise, add it to the end of the buffer; only applies in non-autoloop mode (autoloop mode always restarts)
- **single\_shot** – specifies some number of samples from the start as “single-shot”, so whenever the waveform is repeated (either to reach *minsamp* samples, or when `fill_voltage_output_buffer()` is called), this part is ignored, and only the rest is repeated

**get\_voltage\_output\_buffer\_fill** ()

Get the number of samples still in the output buffer.

Only applies to non-autoloop mode, and return `None` otherwise.

**fill\_voltage\_output\_buffer** (*minsamp=1*)

Add samples to the output buffer until there are at least *minsamp* samples there.

Only applies to non-autoloop mode, and does nothing otherwise. The added samples are determined based on the last data written by `set_voltage_outputs()` and the `single_shot` argument specified there.

**get\_voltage\_outputs** (*names=None*)

Get values of one or several analog voltage outputs.

**Parameters** *names* (*str* or [*str*] or *None*) – name or list of names of outputs (None means all outputs).

Return list of values ordered by *names* (or by `get_voltage_output_channels()` if *names*==None). For continuous waveforms, return the array containing a single repetition of the waveform. For finite waveforms, repeat the array containing the last outputted waveform.

**class** `NoParameterCaller` (*device*, *kind*)

Bases: `object`

Class to simplify calling functions without a parameter

**apply\_settings** (*settings*)

Apply the settings.

*settings* is the dict {name: value} of the device available settings. Non-applicable settings are ignored.

**get\_device\_variable** (*key*)

Get the value of a settings, status, or full info parameter

**get\_full\_info** (*include*=0)

Get dict {name: value} containing full device information (including status and settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include*=-10 queries all available variables, which is equivalent to *include*="all".

**get\_full\_status** (*include*=0)

Get dict {name: value} containing the device status (including settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include*=-10 queries all available variables, which is equivalent to *include*="all".

**get\_settings** (*include*=0)

Get dict {name: value} containing all the device settings.

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include*=-10 queries all available variables, which is equivalent to *include*="all".

**set\_device\_variable** (*key*, *value*)

Set the value of a settings parameter

**setup\_voltage\_output\_clock** (*rate*=0, *sync\_with\_ai*=False, *continuous*=True, *samps\_per\_chan*=1000, *autoloop*=True, *minsamp*=1)

Setup analog output clock configuration.

#### Parameters

- **rate** – clock rate; if 0, assume constant voltage output (default)
- **sync\_with\_ai** – if True, the clock is synchronized to the analog input clock (the main clock); note that in this case output changes only when the analog read task is running
- **continuous** – if True, any written waveform gets repeated continuously; otherwise, it outputs written waveform only once, and then latches the output on the last value

- **samps\_per\_chan** – if `continuous==False`, it determines number of samples to output before stopping; otherwise, it determines the size of the output buffer
- **autoloop** – if it is `True`, then the specified output waveforms are automatically repeated to create a periodic output signal (referred to as “regeneration mode” in NI DAQ terminology); otherwise, written output data is “exhausted” once sent to the output, so the application needs to continuously write output waveforms to avoid output buffer from running empty (which causes an error). This mode gives better control over the output and allows to seamlessly adjust it in real time, but it is more demanding on the application.
- **minsamp** – if the waveform has been specified before, this argument sets the minimal number of samples to write to the output buffer after the clock is set up and the output is restarted

**get\_voltage\_output\_clock\_parameters()**

Get analog output clock configuration.

Return tuple (rate, sync\_with\_ai, continuous, samps\_per\_chan, autoloop).

## Module contents

### pylablib.devices.Newport package

#### Submodules

#### pylablib.devices.Newport.base module

**exception** `pylablib.devices.Newport.base.NewportError`

Bases: `pylablib.core.devio.base.DeviceError`

Generic Newport device error

**args**

**with\_traceback()**

Exception.with\_traceback(tb) – set self.\_\_traceback\_\_ to tb and return self.

**exception** `pylablib.devices.Newport.base.NewportBackendError` (*exc*)

Bases: `pylablib.devices.Newport.base.NewportError`, `pylablib.core.devio.comm_backend.DeviceBackendError`

Newport backend communication error

**args**

**with\_traceback()**

Exception.with\_traceback(tb) – set self.\_\_traceback\_\_ to tb and return self.

#### pylablib.devices.Newport.picomotor module

`pylablib.devices.Newport.picomotor.get_usb_devices_number()`

Get the number of controllers connected via USB

`pylablib.devices.Newport.picomotor.muxaddr(*args, **kwargs)`

Multiplex the function over its addr argument

```
class pylablib.devices.Newport.picomotor.TDeviceInfo(id)
```

Bases: `tuple`

```
count()
```

Return number of occurrences of value.

```
id
```

```
index()
```

Return first index of value.

Raises `ValueError` if the value is not present.

```
class pylablib.devices.Newport.picomotor.Picomotor8742(conn=0, backend='auto',  
                                                    timeout=5.0, multi-  
                                                    addr=False, scan=True)
```

Bases: `pylablib.core.devio.comm_backend.ICommBackendWrapper`, `pylablib.devices.interface.stage.IMultiaxisStage`

Picomotor 8742 4-axis controller.

#### Parameters

- **conn** – connection parameters; can be an index (starting from 0) for USB devices, or an IP address (e.g., "192.168.0.2") or host name (e.g., "8742-12345") for Ethernet devices
- **backend** – communication backend; by default, try to determine from the communication parameters
- **timeout** (*float*) – default operation timeout
- **multiaddr** – if `True`, assume that there are several daisy-chained devices connected to the current one; in this case, `get_device_info` and related methods return dictionaries `{addr: value}` for all connected controllers instead of simply values for the given controller
- **scan** – if `True` and `multiaddr==True`, scan for all connected devices (call `scan_devices()`) upon connection

#### Error

alias of `pylablib.devices.Newport.base.NewportError`

```
query(comm, axis=None, addr=None, read_reply=None)
```

```
get_id(addr=None)
```

Get the device identification string

```
get_device_info(addr=None)
```

Get the device info of the controller board: (`id_string`,)

```
reset(addr=None)
```

Restart the device.

Reboots the CPU and restores all saved settings from the parameter memory.

```
save_parameters(addr=None)
```

Store current parameters to the non-volatile memory.

Affects axes speed and acceleration, motor types, and Ethernet parameters.

```
restore_parameters(src='memory', addr=None)
```

Restore parameters from the non-volatile memory (if `src=="memory"`) for factory parameters (if `src=="factory"`).



Affects axes speed and acceleration, motor types, and Ethernet parameters.

**scan\_devices** (*reassign='conflict', sync=True*)

Scan for devices connected to the current host device via RS-485 daisy-chaining.

*reassign* controls how device addresses are assigned during the scan; can be "none" (keep current values; can lead to conflicts if several devices have the same address), "conflict" (change conflicting addresses), or "all" (assigned all new addresses in sequence starting from the host)

If *sync==True*, wait until the scan is done (might take several seconds).

**get\_addr\_map** ()

Get address map for devices connected to the current host device via RS-485 daisy-chaining.

Return tuple (*addresses, conflict*), where *addresses* is the list of all device addresses, and *conflict==True* if there address conflicts (several devices having the same address).

**wait\_for\_scan** (*timeout=10.0*)

Wait for the device connection scan to finish

**get\_addr** (*addr=None*)

Get RS-485 address of the given device (host if *addr* is None)

**set\_addr** (*new\_addr, addr=None*)

Set RS-485 address of the given device (host if *addr* is None)

**get\_ethernet\_parameters** (*addr=None*)

Get Ethernet connection parameters.

Return tuple (*hostname, ipaddr, ipmode, gateway, netmask*).

**setup\_ethernet** (*hostname=None, ipmode=None, ipaddr=None, gateway=None, netmask=None, addr=None*)

Setup Ethernet connection parameters.

Any None value remains unchanged. Note that these settings only take effect after saving parameters to the memory (*save\_parameters()*) and restarting the device (*reset()*). If the connection is made through Ethernet, then it will likely be invalidated, in which case a new device object with the updated parameters should be created after reset.

**autodetect\_motors** (*addr=None*)

Autodetect connected motors.

The command involves sending single-step commands to the motors, so it requires all axes to be stopped, and it might slightly affect the current position. After the detection the types can be stored in the memory via *save\_parameters()*.

**get\_motor\_type** (*axis='all', addr=None*)

Get type of the given axis motor

**set\_motor\_type** (*axis='all', motor\_type='standard', addr=None*)

Manually set type of the given axis motor

**move\_to** (*axis, position, addr=None*)

Move to a given position

**move\_by** (*axis, steps=1, addr=None*)

Move by a given number of steps

**get\_position** (*axis='all', addr=None*)

Get the current axis position

**set\_position\_reference** (*axis, position=0, addr=None*)

Set the current axis position as a reference (the actual motor position stays the same)

**jog** (*axis*, *direction*, *addr*=0)

Jog a given axis in a given direction.

*direction* can be either "-" (negative) or "+" (positive). The motion continues until it is explicitly stopped.

**is\_moving** (*axis*='all', *addr*=None)

Check if the axis is moving

**wait\_move** (*axis*='all', *addr*=None)

Wait until axis motion is done

**stop** (*axis*='all', *immediate*=False, *addr*=None)

Stop motion of a given axis.

If *immediate*==True make an abrupt stop; otherwise, slow down gradually. Note that immediate stop has to stop all axes simultaneously, so it only takes *axis*=="all".

**get\_velocity\_parameters** (*axis*='all', *addr*=None)

Return velocity parameters (*speed*, *accel*) for the given axis and controller.

*speed* and *accel* denote, correspondingly, maximal (i.e., steady regime) moving speed and acceleration in steps/s and steps/s<sup>2</sup>.

**setup\_velocity** (*axis*='all', *speed*=None, *accel*=None, *addr*=None)

Setup velocity parameters (*speed*, *accel*) for the given axis and controller.

*speed* and *accel* denote, correspondingly, maximal (i.e., steady regime) moving speed and acceleration in steps/s and steps/s<sup>2</sup>. None values are left unchanged.

**class NoParameterCaller** (*device*, *kind*)

Bases: `object`

Class to simplify calling functions without a parameter

**apply\_settings** (*settings*)

Apply the settings.

*settings* is the dict {*name*: *value*} of the device available settings. Non-applicable settings are ignored.

**close** ()

Close the backend

**get\_all\_axes** ()

Get the list of all available axes (taking mapping into account)

**get\_device\_variable** (*key*)

Get the value of a settings, status, or full info parameter

**get\_full\_info** (*include*=0)

Get dict {*name*: *value*} containing full device information (including status and settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include*=-10 queries all available variables, which is equivalent to *include*="all".

**get\_full\_status** (*include*=0)

Get dict {*name*: *value*} containing the device status (including settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the

lowest priority is -10, setting `include=-10` queries all available variables, which is equivalent to `include="all"`.

**get\_settings** (*include=0*)

Get dict {name: value} containing all the device settings.

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting `include=-10` queries all available variables, which is equivalent to `include="all"`.

**is\_opened** ()

Check if the device is connected

**lock** (*timeout=None*)

Lock the access to the device from other threads/processes (isn't necessarily implemented)

**locking** (*timeout=None*)

Context manager for lock & unlock

**open** ()

Open the backend

**remap\_axes** (*mapping, accept\_original=True*)

Rename axes to the new labels.

*mapping* is the new axes mapping, which can be a list of new axes name (corresponding to the old axes in order returned by `get_all_axes()`), or a dictionary {alias: original} of the new axes aliases.

**set\_device\_variable** (*key, value*)

Set the value of a settings parameter

**unlock** ()

Unlock the access to the device from other threads/processes (isn't necessarily implemented)

## Module contents

### pylablib.devices.OZOptics package

#### Submodules

#### pylablib.devices.OZOptics.base module

**exception** `pylablib.devices.OZOptics.base.OZOpticsError`

Bases: `pylablib.core.devio.base.DeviceError`

Generic OZOptics devices error

**args**

**with\_traceback** ()

Exception.with\_traceback(tb) – set self.\_\_traceback\_\_ to tb and return self.

**exception** `pylablib.devices.OZOptics.base.OZOpticsBackendError` (*exc*)

Bases: `pylablib.devices.OZOptics.base.OZOpticsError`, `pylablib.core.devio.comm_backend.DeviceBackendError`

Generic OZOptics backend communication error

**args**

**with\_traceback()**

Exception.with\_traceback(tb) – set self.\_\_traceback\_\_ to tb and return self.

**class** pylablib.devices.OZOptics.base.OZOpticsDevice(*conn, timeout=20.0*)

Bases: *pylablib.core.devio.comm\_backend.ICommBackendWrapper*

Generic OZOptics device.

**Parameters** *conn* – serial connection parameters (usually port or a tuple containing port and baudrate)

**Error**

alias of *OZOpticsError*

**query**(*comm, prefix=None, prefix\_line=None, timeout=None*)

Query the device.

If *prefix* is not *None*, it can specify a string which should be at the beginning of the *prefix\_line* line of the reply. If it is present, it is removed and the rest of that line is returned; otherwise, an error is raised. If *prefix\_line* is *None*, return the first reply line beginning with the given prefix value (or raise an error if not such line is present).

**restart()**

Restart the device

**get\_config()**

Get device configuration

**class** NoParameterCaller(*device, kind*)

Bases: *object*

Class to simplify calling functions without a parameter

**apply\_settings**(*settings*)

Apply the settings.

*settings* is the dict {name: value} of the device available settings. Non-applicable settings are ignored.

**close()**

Close the backend

**get\_device\_variable**(*key*)

Get the value of a settings, status, or full info parameter

**get\_full\_info**(*include=0*)

Get dict {name: value} containing full device information (including status and settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_full\_status**(*include=0*)

Get dict {name: value} containing the device status (including settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_settings** (*include=0*)

Get dict {name: value} containing all the device settings.

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**is\_opened** ()

Check if the device is connected

**lock** (*timeout=None*)

Lock the access to the device from other threads/processes (isn't necessarily implemented)

**locking** (*timeout=None*)

Context manager for lock & unlock

**open** ()

Open the backend

**set\_device\_variable** (*key, value*)

Set the value of a settings parameter

**unlock** ()

Unlock the access to the device from other threads/processes (isn't necessarily implemented)

**class** `pylablib.devices.OZOptics.base.TF100` (*conn, timeout=20.0*)

Bases: `pylablib.devices.OZOptics.base.OZOpticsDevice`

OZOptics TF100 tunable filter.

**Parameters** *conn* – serial connection parameters (usually port or a tuple containing port and baudrate)

**get\_wavelength\_correction** ()

Get the current wavelength correction parameters (*shift, scale*).

The relation between the set/get wavelength and the wavelength set to the device is calculated as  $\text{device\_wavelength} = \text{set\_wavelength} * \text{scale} + \text{shift}$

**set\_wavelength\_correction** (*shift=0.0, scale=1.0*)

Set the wavelength correction parameters.

The relation between the set/get wavelength and the wavelength set to the device is calculated as  $\text{device\_wavelength} = \text{set\_wavelength} * \text{scale} + \text{shift}$

**home** ()

Home the motor (needs to be called first after startup)

**get\_wavelength** ()

Get the currently set wavelength (or None if unknown / not homed)

**set\_wavelength** (*wavelength*)

Set the current wavelength

**Error**

alias of `OZOpticsError`

**class** `NoParameterCaller` (*device, kind*)

Bases: `object`

Class to simplify calling functions without a parameter

**apply\_settings** (*settings*)

Apply the settings.

*settings* is the dict {name: value} of the device available settings. Non-applicable settings are ignored.

**close** ()

Close the backend

**get\_config** ()

Get device configuration

**get\_device\_variable** (*key*)

Get the value of a settings, status, or full info parameter

**get\_full\_info** (*include=0*)

Get dict {name: value} containing full device information (including status and settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_full\_status** (*include=0*)

Get dict {name: value} containing the device status (including settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_settings** (*include=0*)

Get dict {name: value} containing all the device settings.

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**is\_opened** ()

Check if the device is connected

**lock** (*timeout=None*)

Lock the access to the device from other threads/processes (isn't necessarily implemented)

**locking** (*timeout=None*)

Context manager for lock & unlock

**open** ()

Open the backend

**query** (*comm, prefix=None, prefix\_line=None, timeout=None*)

Query the device.

If *prefix* is not None, it can specify a string which should be at the beginning of the *prefix\_line* line of the reply. If it is present, it is removed and the rest of that line is returned; otherwise, an error is raised. If *prefix\_line* is None, return the first reply line beginning with the given prefix value (or raise an error if not such line is present).

**restart** ()

Restart the device

```

set_device_variable (key, value)
    Set the value of a settings parameter

unlock ()
    Unlock the access to the device from other threads/processes (isn't necessarily implemented)

class pylablib.devices.OZOptics.base.DD100 (conn, timeout=20.0)
    Bases: pylablib.devices.OZOptics.base.OZOpticsDevice
    OZOptics DD100 variable attenuator.

        Parameters conn – serial connection parameters (usually port or a tuple containing port and baudrate)

home ()
    Home the motor (needs to be called first after startup)

get_min_attenuation ()
    Get the minimal possible attenuation (i.e., insertion loss)

get_max_attenuation ()
    Get the maximal possible possible attenuation in dB

get_attenuation ()
    Get the current attenuation in dB

set_attenuation (att)
    Set the current attenuation in dB

Error
    alias of OZOpticsError

class NoParameterCaller (device, kind)
    Bases: object
    Class to simplify calling functions without a parameter

apply_settings (settings)
    Apply the settings.

    settings is the dict {name: value} of the device available settings. Non-applicable settings are ignored.

close ()
    Close the backend

get_config ()
    Get device configuration

get_device_variable (key)
    Get the value of a settings, status, or full info parameter

get_full_info (include=0)
    Get dict {name: value} containing full device information (including status and settings).

    include specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting include=-10 queries all available variables, which is equivalent to include="all".

get_full_status (include=0)
    Get dict {name: value} containing the device status (including settings).

```

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting `include=-10` queries all available variables, which is equivalent to `include="all"`.

**get\_settings** (*include=0*)

Get dict {name: value} containing all the device settings.

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting `include=-10` queries all available variables, which is equivalent to `include="all"`.

**is\_opened** ()

Check if the device is connected

**lock** (*timeout=None*)

Lock the access to the device from other threads/processes (isn't necessarily implemented)

**locking** (*timeout=None*)

Context manager for lock & unlock

**open** ()

Open the backend

**query** (*comm, prefix=None, prefix\_line=None, timeout=None*)

Query the device.

If *prefix* is not `None`, it can specify a string which should be at the beginning of the *prefix\_line* line of the reply. If it is present, it is removed and the rest of that line is returned; otherwise, an error is raised. If *prefix\_line* is `None`, return the first reply line beginning with the given prefix value (or raise an error if not such line is present).

**restart** ()

Restart the device

**set\_device\_variable** (*key, value*)

Set the value of a settings parameter

**unlock** ()

Unlock the access to the device from other threads/processes (isn't necessarily implemented)

**class** `pylablib.devices.OZOptics.base.EPC04` (*conn, timeout=20.0*)

Bases: `pylablib.core.devio.comm_backend.ICommBackendWrapper`

OZOptics EPC04 polarization controller.

**Parameters** *conn* – serial connection parameters (usually port or a tuple containing port and baudrate)

**Error**

alias of `OZOpticsError`

**query** (*comm*)

**get\_voltages** ()

Get all voltages

**set\_voltage** (*channel, voltage*)

Set voltage at a given channel (0 through 3)

**set\_all\_voltages** (*voltages*)

Set all channel voltages.



*voltages* is a list of size 4 containing the voltage values.

**step\_voltage** (*channel*, *step*)

Step voltage at the given channel by the given step

**get\_mode** ()

Get current operating mode.

Can be "dc" (constant voltage) or "ac" (scrambling).

**set\_mode** (*mode*='dc')

Set current operating mode.

Can be "dc" (constant voltage) or "ac" (scrambling).

**get\_frequencies** ()

Get all scrambling frequencies

**set\_frequency** (*channel*, *frequency*)

Set scrambling frequency a given channel (0 through 3)

**set\_all\_frequencies** (*frequencies*)

Set all channel scrambling frequencies.

*frequencies* is a list of size 4 containing the frequency values.

**get\_waveform** ()

Get current scrambling waveform.

Can be "sin" (sine wave) or "tri" (triangle wave).

**set\_waveform** (*waveform*)

Set current scrambling waveform.

Can be "sin" (sine wave) or "tri" (triangle wave).

**save\_preset** ()

Save current state as a power-up preset

**class NoParameterCaller** (*device*, *kind*)

Bases: `object`

Class to simplify calling functions without a parameter

**apply\_settings** (*settings*)

Apply the settings.

*settings* is the dict {name: value} of the device available settings. Non-applicable settings are ignored.

**close** ()

Close the backend

**get\_device\_variable** (*key*)

Get the value of a settings, status, or full info parameter

**get\_full\_info** (*include*=0)

Get dict {name: value} containing full device information (including status and settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include*=-10 queries all available variables, which is equivalent to *include*="all".

**get\_full\_status** (*include=0*)

Get dict {name: value} containing the device status (including settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_settings** (*include=0*)

Get dict {name: value} containing all the device settings.

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**is\_opened** ()

Check if the device is connected

**lock** (*timeout=None*)

Lock the access to the device from other threads/processes (isn't necessarily implemented)

**locking** (*timeout=None*)

Context manager for lock & unlock

**open** ()

Open the backend

**set\_device\_variable** (*key, value*)

Set the value of a settings parameter

**unlock** ()

Unlock the access to the device from other threads/processes (isn't necessarily implemented)

## Module contents

### pylablib.devices.Ophir package

#### Submodules

#### pylablib.devices.Ophir.base module

**exception** pylablib.devices.Ophir.base.OphirError

Bases: `pylablib.core.devio.base.DeviceError`

Generic Ophir device error

**args****with\_traceback** ()

Exception.with\_traceback(tb) – set self.\_\_traceback\_\_ to tb and return self.

**exception** pylablib.devices.Ophir.base.OphirBackendError (*exc*)

Bases: `pylablib.devices.Ophir.base.OphirError`, `pylablib.core.devio.comm_backend.DeviceBackendError`

Generic Ophir backend communication error

**args**

**with\_traceback()**

Exception.with\_traceback(tb) – set self.\_\_traceback\_\_ to tb and return self.

**class** pylablib.devices.Ophir.base.OphirDevice(*conn*)

Bases: *pylablib.core.devio.comm\_backend.ICommBackendWrapper*

Generic Ophir device.

**Parameters** *conn* – serial connection parameters (usually port or a tuple containing port and baudrate)

**Error**

alias of *OphirError*

**query**(*comm*)

Send a query to the device and parse the reply

**class** NoParameterCaller(*device, kind*)

Bases: *object*

Class to simplify calling functions without a parameter

**apply\_settings**(*settings*)

Apply the settings.

*settings* is the dict {name: value} of the device available settings. Non-applicable settings are ignored.

**close**()

Close the backend

**get\_device\_variable**(*key*)

Get the value of a settings, status, or full info parameter

**get\_full\_info**(*include=0*)

Get dict {name: value} containing full device information (including status and settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_full\_status**(*include=0*)

Get dict {name: value} containing the device status (including settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_settings**(*include=0*)

Get dict {name: value} containing all the device settings.

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**is\_opened**()

Check if the device is connected

**lock**(*timeout=None*)

Lock the access to the device from other threads/processes (isn't necessarily implemented)

```
locking (timeout=None)
    Context manager for lock & unlock

open ()
    Open the backend

set_device_variable (key, value)
    Set the value of a settings parameter

unlock ()
    Unlock the access to the device from other threads/processes (isn't necessarily implemented)

class pylablib.devices.Ophir.base.THeadInfo (type, serial, name, capabilities)
    Bases: tuple

    capabilities

    count ()
        Return number of occurrences of value.

    index ()
        Return first index of value.

        Raises ValueError if the value is not present.

    name

    serial

    type

class pylablib.devices.Ophir.base.TDeviceInfo (id, serial, name, rom_version)
    Bases: tuple

    count ()
        Return number of occurrences of value.

    id

    index ()
        Return first index of value.

        Raises ValueError if the value is not present.

    name

    rom_version

    serial

class pylablib.devices.Ophir.base.TWavelengthInfo (mode, rng, curr_idx, presets,
                                                    curr_wavelength)
    Bases: tuple

    count ()
        Return number of occurrences of value.

    curr_idx

    curr_wavelength

    index ()
        Return first index of value.

        Raises ValueError if the value is not present.

    mode
```

**presets**

**rng**

**class** `pylablib.devices.Ophir.base.TRangeInfo` (*curr\_idx, ranges, curr\_range*)

Bases: `tuple`

**count** ()

Return number of occurrences of value.

**curr\_idx**

**curr\_range**

**index** ()

Return first index of value.

Raises `ValueError` if the value is not present.

**ranges**

**class** `pylablib.devices.Ophir.base.VegaPowerMeter` (*conn*)

Bases: `pylablib.devices.Ophir.base.OphirDevice`

Ophir Vega power meter.

**Parameters** **conn** – serial connection parameters (usually port or a tuple containing port and baudrate)

**get\_head\_info** ()

Get head information.

Return tuple (*type, serial, name, capabilities*).

**get\_device\_info** ()

Get device information.

Return tuple (*id, serial, name, rom\_version*).

**reset** ()

Reset the device

**get\_power** ()

Get the current power readings.

Return either measured power, or "over", if the power is overrange.

**get\_energy** ()

Get the current energy readings.

Return either measured energy, or "over", if the energy is overrange.

**get\_frequency** ()

Get the current frequency readings.

Return either measured frequency, or "over", if the power is overrange.

**get\_units** ()

Get device reading units

**get\_wavelength\_info** ()

Get wavelength setting info.

Return tuple (*mode, rng, curr\_idx, presets, curr\_wavelength*), where *mode* is the measurement mode ("continuous" or "discrete"), *rng* is a 2-tuple with the full wavelength range

(in m) for continuous mode or a set of all wavelengths for discrete mode, *curr\_idx* is the current wavelength preset index, *presets* is the list of all preset wavelengths (in m) for continuous mode or a set of all wavelengths for discrete mode, and *curr\_wavelength* is the current measurement wavelength (in m) for continuous mode or the current wavelength name for discrete mode.

**get\_wavelength()**

Get current wavelength

**set\_wavelength(wavelength)**

Set current wavelength.

*wavelength* is either a wavelength (in m) for the continuous mode, or a wavelength preset (as a string) for a discrete mode.

**get\_range\_info()**

Get power range info.

Return tuple (*curr\_idx*, *ranges*, *curr\_range*), where *curr\_idx* is the current power range index, *ranges* is the list of ranges (in W) for all indices and *curr\_range* is the current range (in W).

**get\_range()**

Get current power range (maximal power in W)

**get\_range\_idx()**

Get current power range index

Index goes from 0 (highest) to maximal (lowest); auto-ranging is -1.

**set\_range\_idx(rng\_idx)**

Set current range index.

*rng\_idx* is the range index from 0 (highest) to maximal (lowest); auto-ranging is -1. The corresponding ranges are given by [\*get\\_range\\_info\(\)\*](#).

**get\_battery\_condition()**

Check if the batter is OK

**get\_baudrate()**

Get current baud rate

**get\_supported\_baudrates()**

Get a list of all supported baud rates

**set\_baudrate(baudrate)**

Set current baud rate.

If the baudrate is different from the current one, close the device connection. The device object will need to be re-created with the newly specified baud rate.

**is\_filter\_in()**

Check if the filter is set to be on at the power meter

**set\_filter(filter\_in=True)**

Change the filter setting at the power meter (on or off)

**is\_diffuser\_in()**

Check if the diffuser is set to be on at the power meter

**set\_diffuser(diffuser\_in=True)**

Change the diffuser setting at the power meter (on or off)

**Error**

alias of [\*OphirError\*](#)

**class NoParameterCaller** (*device, kind*)

Bases: `object`

Class to simplify calling functions without a parameter

**apply\_settings** (*settings*)

Apply the settings.

*settings* is the dict {*name*: *value*} of the device available settings. Non-applicable settings are ignored.

**close** ()

Close the backend

**get\_device\_variable** (*key*)

Get the value of a settings, status, or full info parameter

**get\_full\_info** (*include=0*)

Get dict {*name*: *value*} containing full device information (including status and settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_full\_status** (*include=0*)

Get dict {*name*: *value*} containing the device status (including settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_settings** (*include=0*)

Get dict {*name*: *value*} containing all the device settings.

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**is\_opened** ()

Check if the device is connected

**lock** (*timeout=None*)

Lock the access to the device from other threads/processes (isn't necessarily implemented)

**locking** (*timeout=None*)

Context manager for lock & unlock

**open** ()

Open the backend

**query** (*comm*)

Send a query to the device and parse the reply

**set\_device\_variable** (*key, value*)

Set the value of a settings parameter

**unlock** ()

Unlock the access to the device from other threads/processes (isn't necessarily implemented)

## Module contents

### pylablib.devices.PCO package

#### Submodules

#### pylablib.devices.PCO.SC2 module

`pylablib.devices.PCO.SC2.list_cameras` (*cam\_interface=None*)

List camera connections (interface kind and camera index).

If *cam\_interface* is supplied, it defines one of camera interfaces to check (e.g., "usb3" or "clhs"). Otherwise, check all interfaces.

`pylablib.devices.PCO.SC2.get_cameras_number` (*cam\_interface=None*)

Get the total number of connected PCOSC2 cameras.

If *cam\_interface* is supplied, it defines one of camera interfaces to check (e.g., "usb3" or "clhs"). Otherwise, check all interfaces.

`pylablib.devices.PCO.SC2.reset_api` ()

Reset API.

All cameras must be closed; otherwise, the prompt to reboot will appear.

**class** `pylablib.devices.PCO.SC2.TDeviceInfo` (*model, interface, sensor, serial\_number*)

Bases: `tuple`

**count** ()

Return number of occurrences of value.

**index** ()

Return first index of value.

Raises `ValueError` if the value is not present.

**interface**

**model**

**sensor**

**serial\_number**

**class** `pylablib.devices.PCO.SC2.TCameraStatus` (*status, warnings, errors*)

Bases: `tuple`

**count** ()

Return number of occurrences of value.

**errors**

**index** ()

Return first index of value.

Raises `ValueError` if the value is not present.

**status**

**warnings**



```

class pylablib.devices.PCO.SC2.TInternalBufferStatus (scheduled, scheduled_max,
                                                    overruns)
    Bases: tuple
    count ()
        Return number of occurrences of value.
    index ()
        Return first index of value.
        Raises ValueError if the value is not present.
    overruns
    scheduled
    scheduled_max

class pylablib.devices.PCO.SC2.TFrameInfo (frame_index)
    Bases: tuple
    count ()
        Return number of occurrences of value.
    frame_index
    index ()
        Return first index of value.
        Raises ValueError if the value is not present.

class pylablib.devices.PCO.SC2.PCOSC2Camera (idx=0, cam_interface=None, re-
                                           boot_on_fail=True)
    Bases: pylablib.devices.interface.camera.IBinROICamera, pylablib.devices.
           interface.camera.IExposureCamera
    PCO SC2 camera.

    Parameters
    • idx (int) – camera index (use get_cameras_number() to get the total number of
      connected cameras)
    • cam_interface – camera interface; if it is None, get the first available connected
      camera (in this case idx is ignored); if not, then value of idx is used to connect to a
      particular camera (interfaces and indices can be obtain from list_cameras())
    • reboot_on_fail (bool) – if True and the camera raised an error during initial-
      ization (but after opening), reboot the camera and try to connect again useful when the
      camera is in a broken state (e.g., wrong ROI or pixel clock settings)

    Error = <Mock name='mock.PCOSC2Error' id='140096280915600'>
    TimeoutError = <Mock spec='str' id='140096270971024'>

    open ()
        Open connection to the camera
    close ()
        Close connection to the camera
    is_opened ()
        Check if the device is connected

```

**reboot** (*wait=True*)

Reboot the camera.

If *wait==True*, wait for the recommended time (10 seconds) after reboot for the camera to fully restart; attempt to open the camera before that can lead to an error.

**get\_full\_camera\_data** ()

Get a dictionary the all camera data available through the SDK

**update\_full\_data** ()

Update internal full camera data settings.

Takes some time (about 50ms), so more specific function are preferable for specific parameters.

**get\_device\_info** ()

Get camera model data.

Return tuple (*model*, *interface*, *sensor*, *serial\_number*).

**get\_capabilities** ()

Get camera capabilities.

For description of the capabilities, see PCO SC2 manual.

**get\_camera\_status** (*full=False*)

Get camera status.

If *full==True*, return current camera status as a set of enabled status states; otherwise, return tuple (*status*, *warnings*, *errors*) with additional information about warnings and error.

**get\_temperature** ()

Get the current camera temperature

Return tuple (*CCD*, *cam*, *power*) with temperatures of the sensor, camera, and power supply respectively.

**get\_conversion\_factor** ()

Get camera conversion factor (electrons per pixel value)

**get\_trigger\_mode** ()

Get current trigger mode (see [\*set\\_trigger\\_mode\(\)\*](#) for description)

**set\_trigger\_mode** (*mode*)

Set trigger mode.

Can be "int" (internal), "software" (software), "ext" (external+software), "ext\_exp" (external exposure), "ext\_sync" (external PLL sync), "ext\_exp\_fast" (fast external exposure), "ext\_cds" (external CDS control), "ext\_exp\_slow" (slow external exposure), or "ext\_sync\_hdsdi" (external synchronized SD/HDI).

For description, see PCO SDK manual.

**send\_software\_trigger** ()

Send software trigger signal

**class ScheduleLooper**

Bases: [\*object\*](#)

Numba schedule loop manager.

Compiles the loop function and provides callback storage.

**reset** ()

**notify** ()

```
class BufferManager (nbuff, size, metadata_size=0)
```

Bases: `object`

Frame buffer managers.

Stores and accesses frame buffer and status arrays and buffer info.

```
get_buffer_ptr (n)
```

Get address of n'th frame buffer

```
get_internal_buffer_status ()
```

Get the status of the internal smaller API buffer, showing the number of scheduled frames there, and the maximal number that can be scheduled

```
set_exposure (exposure)
```

Set camera exposure

```
get_exposure ()
```

Get current exposure

```
set_frame_delay (frame_delay)
```

Set camera frame delay

```
get_frame_delay ()
```

Get current frame delay

```
set_frame_period (frame_time=0, adjust_exposure=False)
```

Set frame time (frame acquisition period).

If the time can't be achieved even with zero frame delay and `adjust_exposure==True`, try to reduce the exposure to get the desired frame time; otherwise, keep the exposure the same.

```
get_frame_period ()
```

Get current frame time (frame acquisition period)

```
get_frame_timings ()
```

Get acquisition timing.

Return tuple (`exposure, frame_period`).

```
get_pixel_rate ()
```

Get camera pixel rate (in Hz)

```
get_available_pixel_rates ()
```

Get all available pixel rates

```
set_pixel_rate (rate=None)
```

Set camera pixel rate (in Hz)

The rate is always rounded to the closest available. If `rate` is `None`, set the maximal possible rate.

```
setup_acquisition (nframes=100)
```

Setup acquisition.

`nframes` determines number of size of the ring buffer (by default, 100).

```
start_acquisition (*args, **kwargs)
```

Start acquisition.

Can take the same keyword parameters as `:meth:~pylablib.BufferManager.setup_acquisition`. If the acquisition is not set up yet, set it up using the supplied parameters (use default of `setup_acquisition()`, if the parameter is `None`). Otherwise, if any supplied parameters are different from the current ones, change them and reset the acquisition.

**stop\_acquisition()**

Stop acquisition.

Clears buffers as well, so any readout afterwards is impossible.

**acquisition\_in\_progress()**

Check if the acquisition is in progress

**clear\_acquisition()**

Clear acquisition settings

**get\_detector\_size()**

Get camera detector size (in pixels) as a tuple (width, height)

**get\_roi()**

Get current ROI.

Return tuple (hstart, hend, vstart, vend, hbin, vbin). *hstart* and *hend* specify horizontal image extent, *vstart* and *vend* specify vertical image extent (start is inclusive, stop is exclusive, starting from 0), *hbin* and *vbin* specify binning.

**set\_roi(hstart=0, hend=None, vstart=0, vend=None, hbin=1, vbin=1, symmetric=False)**

Setup camera ROI.

*hstart* and *hend* specify horizontal image extent, *vstart* and *vend* specify vertical image extent (start is inclusive, stop is exclusive, starting from 0), *hbin* and *vbin* specify binning. By default, all non-supplied parameters take extreme values (0 for start, maximal for end, 1 for binning). If *symmetric*==True and camera requires symmetric ROI (see [requires\\_symmetric\\_roi\(\)](#)), respect this symmetry in the resulting ROI; otherwise, try to use software ROI feature to set up the required ranges (note: while software ROI does affect the size of the read out frame, it does not change the readout time, which would be the same as with *symmetric*==True).

**requires\_symmetric\_roi()**

Check if the camera requires horizontally or vertically symmetric ROI.

Return a tuple (horizontal, vertical). If True, one might still set up an asymmetric ROI for some cameras using the software ROI feature, but it does not affect camera readout rate

**get\_roi\_limits(hbin=1, vbin=1)**

Get the minimal and maximal ROI parameters.

Return tuple (hlim, vlim), where each element is in turn a limit 5-tuple (min, max, pstep, sstep, maxbin) with, correspondingly, minimal and maximal size, position and size step, and the maximal binning. In some cameras, the step and the minimal size depend on the binning, which can be supplied.

**enable\_pixel\_correction(enable=True)**

Enable or disable hotpixel correction

**is\_pixel\_correction\_enabled()**

Check if hotpixel correction is enabled

**get\_noise\_filter\_mode()**

Get the noise filter mode (for details, see [set\\_noise\\_filter\\_mode\(\)](#))

**set\_noise\_filter\_mode(mode='on')**

Set the noise filter mode.

Can be "off", "on", or "on\_hpc" (on + hot pixel correction).

**set\_status\_line\_mode(binary=True, text=False)**

Set status line mode.

*binary* determines if the binary line is present (it occupies first 14 pixels of the image). *text* determines if the text line is present (it is plane text timestamp, which takes first 8 rows and about 300 columns).

It is recommended to always have *binary* option on, since it is used to determine frame index for checking if there are any missing frames.

**get\_status\_line\_mode()**

Get status line mode.

Return tuple (binary, text) (see [set\\_status\\_line\\_mode\(\)](#) for description)

**get\_bit\_alignment()**

Get data bit alignment

Can be "LSB" (normal alignment) or "MSB" (if camera data is less than 16 bit, it is padded with zeros on the right to match 16 bit).

**set\_bit\_alignment(mode)**

Get data bit alignment

Can be "LSB" (normal alignment) or "MSB" (if camera data is less than 16 bit, it is padded with zeros on the right to match 16 bit).

**set\_metadata\_mode(mode=True)**

Set metadata mode

**get\_metadata\_mode()**

Get metadata mode.

Return tuple (enabled, size, version)

**FrameTransferError**

alias of [pylablib.devices.interface.camera.DefaultFrameTransferError](#)

**class NoParameterCaller(device, kind)**

Bases: [object](#)

Class to simplify calling functions without a parameter

**apply\_settings(settings)**

Apply the settings.

*settings* is the dict {name: value} of the device available settings. Non-applicable settings are ignored.

**get\_acquisition\_parameters()**

Get acquisition parameters.

Return dictionary {name: value}

**get\_data\_dimensions()**

Get readout data dimensions (in pixels) as a tuple (width, height); take indexing mode into account

**get\_device\_variable(key)**

Get the value of a settings, status, or full info parameter

**get\_frame\_format()**

Get format for the returned images.

Can be "list" (list of 2D arrays), "array" (a single 3D array), or "chunks" (list of 3D "chunk" arrays; supported for some cameras and provides the best performance).

**get\_frame\_info\_fields()**

Get the names of frame info fields.

Applicable when frame info format (set by `set_frame_info_format()`) is "list" or "array".

**get\_frame\_info\_format()**

Get format of the frame info.

Can be "namedtuple" (potentially nested named tuples; convenient to get particular values), "list" (flat list of values, with field names are given by `get_frame_info_fields()`; convenient for building a table), "array" (same as "list", but with a numpy array, which is easier to use for "chunks" frame format), or "dict" (flat dictionary with the same fields as the "list" format; more resilient to future format changes)

**get\_frame\_info\_period()**

Get period of frame info acquisition.

Frame info might be skipped (set to None) except for frames which indices are divisible by *period*. Useful for certain cameras where acquiring frame info takes a lot of time and can reduce performance at higher frame rates. Note that this parameter can still be ignored (i.e., always set to 1) if the performance is not an issue for a given camera class.

**get\_frames\_status()**

Get acquisition and buffer status.

Return tuple (*acquired*, *unread*, *skipped*, *size*), where *acquired* is the total number of acquired frames, *unread* is the number of acquired but not read frames, *skipped* is the number of skipped (not read and then written over) frames, and *buffer\_size* is the total buffer size (in frames).

**get\_full\_info(include=0)**

Get dict {name: value} containing full device information (including status and settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include*=-10 queries all available variables, which is equivalent to *include*="all".

**get\_full\_status(include=0)**

Get dict {name: value} containing the device status (including settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include*=-10 queries all available variables, which is equivalent to *include*="all".

**get\_image\_indexing()**

Get indexing for the returned images.

Can be "rct" (first index row, second index column, rows counted from the top), "rcb" (same as "rc", rows counted from the bottom), "xyt" (first index column, second index row, rows counted from the top), or "xyb" (same as "xyt", rows counted from the bottom)

**get\_new\_images\_range()**

Get the range of the new images.

Return tuple (*first*, *last*) with images range (first inclusive). If no images are available, return None. If some images were in the buffer were overwritten, exclude them from the range.

**get\_settings(include=0)**

Get dict {name: value} containing all the device settings.

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include*=-10 queries all available variables, which is equivalent to *include*="all".

**grab** (*nframes=1, frame\_timeout=5.0, missing\_frame='skip', return\_info=False, buff\_size=None*)  
 Snap *nframes* images (with preset image read mode parameters)

*buff\_size* determines buffer size (if *None*, use the default size). Timeout is specified for a single-frame acquisition, not for the whole acquisition time. *missing\_frame* determines what to do with frames which have been lost: can be "none" (replacing them with *None*), "zero" (replacing them with zero-filled frame), or "skip" (skipping them, while still keeping total returned frames number to *n*). If *return\_info==True*, return tuple (*frames*, *infos*), where *infos* is a list of frame info tuples (camera-dependent); if some frames are missing and *missing\_frame!="skip"*, the corresponding frame info is *None*.

**is\_acquisition\_setup** ()  
 Check if acquisition is set up.

If the camera does not support separate acquisition setup, always return *True*.

**pausing\_acquisition** (*clear=None, stop=True, setup\_after=None, start\_after=True, combine\_nested=True*)  
 Context manager which temporarily pauses acquisition during execution of *with* block.

Useful for applying certain settings which can't be changed during the acquisition. If *clear==True*, clear acquisition in addition to stopping (by default, use the class default specified as *\_clear\_pausing\_acquisition* attribute). If *stop==True*, stop the acquisition (if *clear==True*, stop regardless). If *setup\_after==True*, setup the acquisition after pause if necessary (*None* means setup only if clearing was required). If *start\_after==True*, start the acquisition after pause if necessary (*None* means start only if stopping was required). If *combine\_nested==True*, then any nested *pausing\_acquisition* calls will stop/clear acquisition as necessary, but won't setup/start it again until this *pausing\_acquisition* call is complete.

Yields tuple (*acq\_in\_progress*, *acq\_params*), which indicates whether acquisition is currently in progress, and what are the current acquisition parameters.

**read\_multiple\_images** (*rng=None, peek=False, missing\_frame='skip', return\_info=False, return\_rng=False*)  
 Read multiple images specified by *rng* (by default, all un-read images).

If *rng* is specified, it is a tuple (*first*, *last*) with images range (first inclusive). If no new frames are available, return an empty list; if no acquisition is running, return *None*. If *peek==True*, return images but not mark them as read. *missing\_frame* determines what to do with frames which are out of range (missing or lost): can be "none" (replacing them with *None*), "zero" (replacing them with zero-filled frame), or "skip" (skipping them). If *return\_info==True*, return tuple (*frames*, *infos*), where *infos* is a list of frame info tuples (camera-dependent, by default, only the frame index); if some frames are missing and *missing\_frame!="skip"*, the corresponding frame info is *None*. if *return\_rng==True*, return the range covered resulting frames; if *missing\_frame=="skip"*, the range can be smaller than the supplied *rng* if some frames are skipped.

**read\_newest\_image** (*peek=False, return\_info=False*)  
 Read the newest un-read image.

If no un-read frames are available, return *None*. If *peek==True*, return the image but not mark it as read. If *return\_info==True*, return tuple (*frame*, *info*), where *info* is an info tuples (camera-dependent, see [read\\_multiple\\_images\(\)](#)).

**read\_oldest\_image** (*peek=False, return\_info=False*)  
 Read the oldest un-read image.

If no un-read frames are available, return *None*. If *peek==True*, return the image but not mark it as read. If *return\_info==True*, return tuple (*frame*, *info*), where *info* is an info tuples (camera-dependent, see [read\\_multiple\\_images\(\)](#)).

**set\_device\_variable** (*key, value*)  
Set the value of a settings parameter

**set\_frame\_format** (*fmt*)  
Set format for the returned images.

Can be "list" (list of 2D arrays), "array" (a single 3D array), "chunks" (list of 3D "chunk" arrays; supported for some cameras and provides the best performance), or "try\_chunks" (same as "chunks", but if chunks are not supported, set to "list" instead). If format is "chunks" and chunks are not supported by the camera, it results in one frame per chunk. Note that if the format is set to "array" or "chunks", the frame info format is also automatically set to "array". If the format is set to "chunks", then the image info is also returned in chunks form (list of 2D info arrays with the same length as the corresponding frame chunks).

**set\_frame\_info\_format** (*fmt, include\_fields=None*)  
Set format of the frame info.

Can be "namedtuple" (potentially nested named tuples; convenient to get particular values), "list" (flat list of values, with field names are given by `get_frame_info_fields()`; convenient for building a table), "array" (same as "list", but with a numpy array, which is easier to use for "chunks" frame format), or "dict" (flat dictionary with the same fields as the "list" format; more resilient to future format changes) If *include\_fields* is not None, it specifies the fields included for non-"tuple" formats; note that order of *include\_fields* is ignored, and the resulting fields are always ordered same as in the original.

**set\_frame\_info\_period** (*period=1*)  
Set period of frame info acquisition.

Frame info might be skipped (set to None) except for frames which indices are divisible by *period*. Useful for certain cameras where acquiring frame info takes a lot of time and can reduce performance at higher frame rates. Note that this parameter can still be ignored (i.e., always set to 1) if the performance is not an issue for a given camera class.

**set\_image\_indexing** (*indexing*)  
Set up indexing for the returned images.

Can be "rct" (first index row, second index column, rows counted from the top), "rcb" (same as "rc", rows counted from the bottom), "xyt" (first index column, second index row, rows counted from the top), or "xyb" (same as "xyt", rows counted from the bottom)

**snap** (*timeout=5.0, return\_info=False*)  
Snap a single frame

**wait\_for\_frame** (*since='lastread', nframes=1, timeout=20.0, error\_on\_stopped=False*)  
Wait for one or several new camera frames.

*since* specifies the reference point for waiting to acquire *nframes* frames; can be "lastread" (from the last read frame), "lastwait" (wait for the last successful `wait_for_frame()` call), "now" (from the start of the current call), or "start" (from the acquisition start, i.e., wait until *nframes* frames have been acquired). *timeout* can be either a number, None (infinite timeout), or a tuple (*timeout, frame\_timeout*), in which case the call times out if the total time exceeds *timeout*, or a single frame wait exceeds *frame\_timeout*. If the call times out, raise `TimeoutError`. If *error\_on\_stopped==True* and the acquisition is not running, raise `Error`; otherwise, simply return False without waiting.

**class** `pylablib.devices.PCO.SC2.TStatusLine` (*framestamp*)  
Bases: `tuple`

**count** ()  
Return number of occurrences of value.



**framestamp**

**index()**

Return first index of value.

Raises ValueError if the value is not present.

`pylablib.devices.PCO.SC2.get_status_line(frame)`

Get frame info from the binary status line.

Assume that the status line is present; if it isn't, the returned frame info will be a random noise.

`pylablib.devices.PCO.SC2.get_status_lines(frames)`

Get frame info from the binary status line.

*frames* can be 2D array (one frame), 3D array (stack of frames, first index is frame number), or list of 1D or 2D arrays. Assume that the status line is present; if it isn't, the returned frame info will be a random noise. Return a 1D or 2D numpy array, where the first axis (if present) is the frame number, and the last is the status line entry.

**class** `pylablib.devices.PCO.SC2.StatusLineChecker`

Bases: `pylablib.devices.interface.camera.StatusLineChecker`

**get\_framestamp**(*frames*)

Get framestamps from status lines in the given frames

**check\_indices**(*indices, step=1*)

Check if indices are consistent with the given step

## Module contents

### pylablib.devices.Pfeiffer package

#### Submodules

#### pylablib.devices.Pfeiffer.base module

**exception** `pylablib.devices.Pfeiffer.base.PfeifferError`

Bases: `pylablib.core.devio.base.DeviceError`

Generic Pfeiffer device error

**args**

**with\_traceback()**

Exception.with\_traceback(tb) – set self.\_\_traceback\_\_ to tb and return self.

**exception** `pylablib.devices.Pfeiffer.base.PfeifferBackendError`(*exc*)

Bases: `pylablib.devices.Pfeiffer.base.PfeifferError`, `pylablib.core.devio.comm_backend.DeviceBackendError`

Generic Pfeiffer backend communication error

**args**

**with\_traceback()**

Exception.with\_traceback(tb) – set self.\_\_traceback\_\_ to tb and return self.

**class** `pylablib.devices.Pfeiffer.base.TTPG260SwitchSettings`(*channel, low\_thresh, high\_thresh*)

Bases: `tuple`

**channel**

**count** ()

Return number of occurrences of value.

**high\_thresh**

**index** ()

Return first index of value.

Raises ValueError if the value is not present.

**low\_thresh**

**class** pylablib.devices.Pfeiffer.base.**TPG260GaugeControlSettings** (*activation\_control*,  
*deactivation\_control*,  
*on\_thresh*,  
*off\_thresh*)

Bases: *tuple*

**activation\_control**

**count** ()

Return number of occurrences of value.

**deactivation\_control**

**index** ()

Return first index of value.

Raises ValueError if the value is not present.

**off\_thresh**

**on\_thresh**

**class** pylablib.devices.Pfeiffer.base.**TPG260** (*conn*)

Bases: *pylablib.core.devio.comm\_backend.ICommBackendWrapper*

TPG260 series (TPG261/262) pressure gauge.

**Parameters** *conn* – serial connection parameters (usually port or a tuple containing port and baudrate)

**Error**

alias of *PfeifferError*

**comm** (*msg*)

Send a command to the device

**query** (*msg*, *data\_type*='str')

Send a query to the device and return the reply

**get\_units** ()

Get device units for indication/reading ("mbar", "torr", or "pa")

**set\_units** (*units*)

Set device units for indication/reading ("mbar", "torr", or "pa")

**to\_Pa** (*value*, *units*=None)

Convert value in the given units to Pa.

If *units* is None, use the current display units.

**from\_Pa** (*value*, *units=None*)  
 Convert value in the given units from Pa.  
 If *units* is *None*, use the current display units.

**get\_display\_channel** ()  
 Get controller display channel

**set\_display\_channel** (*channel=1*)  
 Set controller display channel

**get\_display\_resolution** ()  
 Get controller display resolution (number of digits)

**set\_display\_resolution** (*resolution=2*)  
 Set controller display resolution (number of digits)

**is\_enabled** (*channel=1*)  
 Check if the gauge at the given channel is enabled.  
 If the gauge cannot be turned on/off (e.g., not connected), return *None*.

**enable** (*enable=True*, *channel=1*)  
 Enable or disable the gauge at the given channel

**get\_channel\_status** (*channel=1*)  
 Get channel status.  
 Can be "ok", "under" (underrange), "over" (overrange), "sensor\_error", "sensor\_off", "no\_sensor", or "id\_error".

**get\_pressure** (*channel=1*, *display\_units=False*, *status\_error=True*)  
 Get pressure at a given channel.  
 If *display\_units==False*, return result in Pa; otherwise, use display units obtained using *get\_units()*. If *status\_error==True* and the channel status is not "ok", raise an error; otherwise, return *None*.

**get\_gauge\_kind** (*channel=1*)

**get\_measurement\_filter** (*channel=1*)  
 Get gauge measurement filter ("fast", "medium", or "slow")

**set\_measurement\_filter** (*meas\_filter*, *channel=1*)  
 Set gauge measurement filter ("fast", "medium", or "slow")

**get\_calibration\_factor** (*channel=1*)  
 Get gauge calibration factor

**set\_calibration\_factor** (*coefficient*, *channel=1*)  
 Set gauge calibration factor

**get\_switch\_settings** (*switch\_function*)  
 Get settings for the given switch function (between 1 and 4).  
 Return tuple (*channel*, *low\_thresh*, *high\_thresh*). The thresholds are given in Pa.

**setup\_switch** (*switch\_function*, *channel*, *low\_thresh*, *high\_thresh*)  
 Get settings for the given switch function (between 1 and 4).  
 Return tuple (*channel*, *low\_thresh*, *high\_thresh*). The thresholds are given in Pa.

**get\_switch\_status** ()  
 Return status of the 4 switch functions

**get\_gauge\_control\_settings** (*channel*)

Get settings for the gauge control on the given channel.

Return tuple (activation\_control, deactivation\_control, on\_thresh, off\_thresh). The thresholds are given in Pa.

**setup\_gauge\_control** (*channel, activation\_control, deactivation\_control, on\_thresh, off\_thresh*)

Setup gauge control on the given channel.

Return tuple (activation\_control, deactivation\_control, on\_thresh, off\_thresh). The thresholds are given in Pa.

**get\_current\_errors** ()

Get a list of all present error messages.

If there are no errors, return a single-element list ["no\_error"].

**reset\_error** ()

Cancel currently active errors and return to measurement mode.

Return the list of currently present errors. If there are no errors, return a single-element list ["no\_error"].

**class NoParameterCaller** (*device, kind*)

Bases: `object`

Class to simplify calling functions without a parameter

**apply\_settings** (*settings*)

Apply the settings.

*settings* is the dict {name: value} of the device available settings. Non-applicable settings are ignored.

**close** ()

Close the backend

**get\_device\_variable** (*key*)

Get the value of a settings, status, or full info parameter

**get\_full\_info** (*include=0*)

Get dict {name: value} containing full device information (including status and settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include*=-10 queries all available variables, which is equivalent to *include*="all".

**get\_full\_status** (*include=0*)

Get dict {name: value} containing the device status (including settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include*=-10 queries all available variables, which is equivalent to *include*="all".

**get\_settings** (*include=0*)

Get dict {name: value} containing all the device settings.

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include*=-10 queries all available variables, which is equivalent to *include*="all".

**is\_opened()**

Check if the device is connected

**lock** (*timeout=None*)

Lock the access to the device from other threads/processes (isn't necessarily implemented)

**locking** (*timeout=None*)

Context manager for lock & unlock

**open()**

Open the backend

**set\_device\_variable** (*key, value*)

Set the value of a settings parameter

**unlock()**

Unlock the access to the device from other threads/processes (isn't necessarily implemented)

**class** `pylablib.devices.Pfeiffer.base.DPG202` (*conn*)

Bases: `pylablib.core.devio.comm_backend.ICommBackendWrapper`

DPG202/TPG202 control unit.

**Parameters** *conn* – serial connection parameters (usually port or a tuple containing port and baudrate)

**Error**

alias of `PfeifferError`

**query** (*parameter, value='?', action=0, address=1, send\_type=None, recv\_type=None*)

Send a query to the device and parse the reply.

**Parameters**

- **parameter** – parameter number
- **value** – value to send ("=?" for a value request)
- **action** – request action (0 for value request, 1 for a command)
- **address** – unit address
- **send\_type** – data type for the sent value (ignored for value requests)
- **recv\_type** – data type for the received value (None means returning a raw string value)

**get\_value** (*parameter, data\_type, address=1*)

Send a data request to the device.

**Parameters**

- **parameter** – parameter number
- **data\_type** – data type for the received value
- **address** – unit address

**comm** (*parameter, value, data\_type, address=1*)

Send a control command to the device.

**Parameters**

- **parameter** – parameter number
- **value** – associated command value

- **data\_type** – data type for the sent value
- **address** – unit address

**get\_pressure** (*address=1*)

Get pressure at a given unit address

**get\_error\_code** (*address=1*)

Get the current error code at a given unit address

**get\_software\_version** (*address=1*)

Get the software version at a given unit address

**get\_device\_name** (*address=1*)

Get the name of the gauge at a given unit address

**class NoParameterCaller** (*device, kind*)

Bases: `object`

Class to simplify calling functions without a parameter

**apply\_settings** (*settings*)

Apply the settings.

*settings* is the dict {name: value} of the device available settings. Non-applicable settings are ignored.

**close** ()

Close the backend

**get\_device\_variable** (*key*)

Get the value of a settings, status, or full info parameter

**get\_full\_info** (*include=0*)

Get dict {name: value} containing full device information (including status and settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_full\_status** (*include=0*)

Get dict {name: value} containing the device status (including settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_settings** (*include=0*)

Get dict {name: value} containing all the device settings.

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**is\_opened** ()

Check if the device is connected

**lock** (*timeout=None*)

Lock the access to the device from other threads/processes (isn't necessarily implemented)

**locking** (*timeout=None*)  
Context manager for lock & unlock

**open** ()  
Open the backend

**set\_device\_variable** (*key, value*)  
Set the value of a settings parameter

**unlock** ()  
Unlock the access to the device from other threads/processes (isn't necessarily implemented)

## Module contents

### pylablib.devices.Photometrics package

#### Submodules

#### pylablib.devices.Photometrics.pvcam module

**class** pylablib.devices.Photometrics.pvcam.**LibraryController** (*lib*)  
Bases: [pylablib.devices.utils.load\\_lib.LibraryController](#)

**close** (*opid*)  
Mark device closing.  
  
Return tuple (*close\_result, uninit\_result*) with the results of the closing and the shutdown.  
If library does not need to be shut down yet, set *uninit\_result=None*

**open** ()  
Mark device opening.  
  
Return tuple (*init\_result, open\_result, opid*) with the results of the initialization and the opening, and the opening ID which should afterwards be used for closing. If library is already initialized, set *init\_result=None*

**preinit** ()  
Pre-initialize the library, if it hasn't been done already

**shutdown** ()  
Close all opened connections and shutdown the library

**temp\_open** ()  
Context for temporarily opening a new device connection

[pylablib.devices.Photometrics.pvcam.list\\_cameras](#) ()  
List all cameras available through Pvcam interface

[pylablib.devices.Photometrics.pvcam.get\\_cameras\\_number](#) ()  
Get number of connected Pvcam cameras

**class** pylablib.devices.Photometrics.pvcam.**PvcamAttribute** (*handle, pid, cam=None*)  
Bases: [object](#)  
  
Object representing an Pvcam camera parameter.  
  
Allows to query and set values and get additional information. Usually created automatically by an [PvcamCamera](#) instance, but could be created manually.

#### Parameters

- **handle** – camera handle
- **pid** – parameter id of the attribute

**name**

attribute name

**kind**

attribute kind; can be "INT8", "INT16", "INT32", "INT64", "UNS8", "UNS16", "UNS32", "UNS64", "FLT32", "FLT64", "ENUM", "BOOLEAN", or "CHAR\_PTR"

**kind**

attribute kind (e.g., "INT8" or "FLT32")

**available**

whether attribute is available on the current hardware

**Type** `bool`

**readable**

whether attribute is readable

**Type** `bool`

**writable**

whether attribute is writable

**Type** `bool`

**min**

minimal attribute value (if applicable)

**Type** `float` or `int`

**max**

maximal attribute value (if applicable)

**Type** `float` or `int`

**inc**

minimal attribute increment value (if applicable)

**Type** `float` or `int`

**ivalues**

list of possible integer values for enum attributes

**values**

list of possible text values for enum attributes

**labels**

dict {label: index} which shows all possible values of an enumerated attribute and their corresponding numerical values

**ilabels**

dict {index: label} which shows labels corresponding to numerical values of an enumerated attribute

**default**

default values of the attribute

**update\_limits()**

Update attribute constraints



**truncate\_value** (*value*)

Truncate value to lie within attribute limits

**get\_value** (*enum\_as\_str=True, error\_on\_noacq=True*)

Get attribute value.

If *enum\_as\_str*==True, return enum-style values as strings; otherwise, return corresponding integer values.

**set\_value** (*value, truncate=True*)

Set attribute value.

If *truncate*==True, automatically truncate value to lie within allowed range.

**class** pylablib.devices.Photometrics.pvcam.**TDeviceInfo** (*vendor, product, chip, system, part, serial*)

Bases: `tuple`

**chip**

**count** ()

Return number of occurrences of value.

**index** ()

Return first index of value.

Raises ValueError if the value is not present.

**part**

**product**

**serial**

**system**

**vendor**

**class** pylablib.devices.Photometrics.pvcam.**TFrameInfo** (*frame\_index, times-tamp\_start\_ns, times-tamp\_end\_ns, framestamp, flags, exposure\_ns*)

Bases: `tuple`

**count** ()

Return number of occurrences of value.

**exposure\_ns**

**flags**

**frame\_index**

**framestamp**

**index** ()

Return first index of value.

Raises ValueError if the value is not present.

**timestamp\_end\_ns**

**timestamp\_start\_ns**

```
class pylablib.devices.Photometrics.pvcam.TReadoutInfo (port_idx, port_name,  
                                                    speed_idx, speed_freq,  
                                                    gain_idx, gain_name)  
  
    Bases: tuple  
  
    count ()  
        Return number of occurrences of value.  
  
    gain_idx  
  
    gain_name  
  
    index ()  
        Return first index of value.  
  
        Raises ValueError if the value is not present.  
  
    port_idx  
  
    port_name  
  
    speed_freq  
  
    speed_idx  
  
class pylablib.devices.Photometrics.pvcam.PvcamCamera (name=None)  
    Bases: pylablib.devices.interface.camera.IBinROICamera, pylablib.devices.  
           interface.camera.IExposureCamera, pylablib.devices.interface.camera.  
           IAttributeCamera  
  
    Generic Pvcam camera interface.  
  
        Parameters serial_number – camera serial number; if None, connect to the first non-used  
            camera  
  
    Error = <Mock name='mock.PvcamError' id='140096282471056'>  
    TimeoutError = <Mock spec='str' id='140096282695184'>  
  
    open ()  
        Open connection to the camera  
  
    close ()  
        Close connection to the camera  
  
    is_opened ()  
        Check if the device is connected  
  
    get_attribute_value (name, error_on_missing=True, error_on_noacq=False, default=None,  
                        enum_as_str=True)  
        Get value of an attribute with the given name.  
  
        If the value doesn't exist or can not be read and error_on_missing==True, raise error; otherwise,  
        return default. If default is not None, assume that error_on_missing==False. If name points at  
        a dictionary branch, return a dictionary with all values in this branch. If enum_as_str==True, return  
        enum-style values as strings; otherwise, return corresponding integer values.  
  
    set_attribute_value (name, value, truncate=True, error_on_missing=True)  
        Set value of an attribute with the given name.  
  
        If the value doesn't exist or can not be written and error_on_missing==True, raise error; otherwise,  
        do nothing. If name points at a dictionary branch, set all values in this branch (in this case value must be  
        a dictionary). If truncate==True, truncate value to lie within attribute range.
```

**get\_all\_attribute\_values** (*root=""*, *enum\_as\_str=True*, *error\_on\_noacq=False*)

Get values of all attributes with the given *root*

**set\_all\_attribute\_values** (*settings*, *root=""*, *truncate=True*)

Set values of all attributes with the given *root*.

If *truncate==True*, truncate value to lie within attribute range.

**get\_attribute\_range** (*name*, *error\_on\_missing=True*, *default=None*, *parameter=None*)

Return attribute range.

For numerical attributes it is a tuple (*min*, *max*), while for enum attributes it is a dictionary {*index*: *name*}. If *parameter* is specified, it is a parameter class used to convert the index for a enum attribute.

**get\_all\_readout\_modes** ()

Get a list of all possible readout modes.

Return a list of tuples (*port\_idx*, *port\_name*, *speed\_idx*, *speed\_freq*, *gain\_idx*, *gain\_name*). The indices (*port*, *speed*, and *gain*) can be used to set up a particular mode using `set_readout_mode()`.

**get\_readout\_mode** (*full=True*)

Get current readout mode.

If *full==True*, return a full tuple (*port\_idx*, *port\_name*, *speed\_idx*, *speed\_freq*, *gain\_idx*, *gain\_name*) containing the descriptions; otherwise, return only indices (*port\_idx*, *speed\_idx*, *gain\_idx*).

**set\_readout\_mode** (*port\_idx=None*, *speed\_idx=None*, *gain\_idx=None*)

Set the readout mode.

Any *None* value stays unchanged.

**get\_device\_info** ()

Get camera information.

Return tuple (*vendor*, *product*, *chip*, *system*, *part*, *serial*).

**get\_pixel\_size** ()

Get camera pixel size (in m)

**get\_pixel\_distance** ()

Get camera pixel distance (in m)

**get\_temperature\_setpoint** ()

Get the temperature setpoint (in C)

**get\_temperature** ()

Get the current camera temperature (in C)

**set\_temperature** (*temp*)

Change the temperature setpoint (in C)

**get\_fan\_mode** ()

Get current fan mode

**set\_fan\_mode** (*fan\_mode='high'*)

Set current fan mode

**is\_metadata\_enabled** ()

Check if metadata is enabled

**enable\_metadata** (*enable=True*)

Enable or disable metadata

**get\_exposure()**  
Get current exposure

**set\_exposure(exposure)**  
Set camera exposure

**get\_clear\_mode()**  
Get sensor clear mode

**set\_clear\_mode(mode)**  
Set sensor clear mode

**get\_clear\_cycles()**  
Get sensor clear cycles

**set\_clear\_cycles(ncycles)**  
Set sensor clear cycles

**get\_clearing\_time()**  
Get sensor clearing time (regardless of the mode)

**get\_readout\_time(include\_clear=True)**  
Get frame readout time.  
  
If `include_clear==True` and the clear mode is per-exposure ("Pre-Exposure" or "Pre-Exposure and Post-Sequence"), include it into this time.

**get\_frame\_timings()**  
Get acquisition timing.  
  
Return tuple (exposure, frame\_period).

**get\_trigger\_mode()**  
Get trigger mode

**set\_trigger\_mode(mode, out\_mode=None)**  
Set trigger mode

**send\_software\_trigger()**  
Send software trigger signal and return whether it has been accepted

**get\_detector\_size()**  
Get camera detector size (in pixels) as a tuple (width, height)

**get\_roi()**  
Get current ROI.  
  
Return tuple (hstart, hend, vstart, vend, hbin, vbin). *hstart* and *hend* specify horizontal image extent, *vstart* and *vend* specify vertical image extent (start is inclusive, stop is exclusive, starting from 0), *hbin* and *vbin* specify binning.

**set\_roi(hstart=0, hend=None, vstart=0, vend=None, hbin=1, vbin=1)**  
Setup camera ROI.  
  
*hstart* and *hend* specify horizontal image extent, *vstart* and *vend* specify vertical image extent (start is inclusive, stop is exclusive, starting from 0), *hbin* and *vbin* specify binning. By default, all non-supplied parameters take extreme values (0 for start, maximal for end, 1 for binning).

**get\_roi\_limits(hbin=1, vbin=1)**  
Get the minimal and maximal ROI parameters.  
  
Return tuple (hlim, vlim), where each element is in turn a limit 5-tuple (min, max, pstep, sstep, maxbin) with, correspondingly, minimal and maximal size, position and size step, and the

maximal binning. In some cameras, the step and the minimal size depend on the binning, which can be supplied.

**get\_supported\_binning\_modes()**

Get all possible binning combinations as a list [(hbin, vbin)]

**setup\_acquisition(mode='sequence', nframes=100)**

Setup acquisition mode.

*mode* can be either "snap" (single frame or a fixed number of frames) or "sequence" (continuous acquisition). *nframes* sets up number of frame buffers.

**clear\_acquisition()**

Clear acquisition settings

**start\_acquisition(\*args, \*\*kwargs)**

Start acquisition.

Can take the same keyword parameters as :meth:“*setup\_acquisition*”. If the acquisition is not set up yet, set it up using the supplied parameters (use default of *setup\_acquisition()*, if the parameter is None). Otherwise, if any supplied parameters are different from the current ones, change them and reset the acquisition.

**stop\_acquisition()**

Stop acquisition

**acquisition\_in\_progress()**

Check if acquisition is in progress

**read\_multiple\_images(rng=None, peek=False, missing\_frame='skip', return\_info=False, return\_rng=False)**

Read multiple images specified by *rng* (by default, all un-read images).

If *rng* is specified, it is a tuple (first, last) with images range (first inclusive). If no new frames are available, return an empty list; if no acquisition is running, return None. If *peek*==True, return images but not mark them as read. *missing\_frame* determines what to do with frames which are out of range (missing or lost): can be "none" (replacing them with None), "zero" (replacing them with zero-filled frame), or "skip" (skipping them). If *return\_info*==True, return tuple (frames, infos), where *infos* is a list of *TFrameInfo* instances describing frame index and frame metadata, which contains start and stop timestamps, framestamp, frame flags, and exposure; if some frames are missing and *missing\_frame*!="skip", the corresponding frame info is None. if *return\_rng*==True, return the range covered resulting frames; if *missing\_frame*=="skip", the range can be smaller than the supplied *rng* if some frames are skipped.

**FrameTransferError**

alias of *pylablib.devices.interface.camera.DefaultFrameTransferError*

**class NoParameterCaller(device, kind)**

Bases: *object*

Class to simplify calling functions without a parameter

**apply\_settings(settings)**

Apply the settings.

*settings* is the dict {name: value} of the device available settings. Non-applicable settings are ignored.

**get\_acquisition\_parameters()**

Get acquisition parameters.

Return dictionary {name: value}

**get\_all\_attributes** (*copy=False*)

Return a dictionary of all available attributes.

If `copy==True`, copy the dictionary; otherwise, return the internal dictionary structure (should not be modified).

**get\_attribute** (*name, error\_on\_missing=True*)

Get the camera attribute with the given name

**get\_data\_dimensions** ()

Get readout data dimensions (in pixels) as a tuple (`width`, `height`); take indexing mode into account

**get\_device\_variable** (*key*)

Get the value of a settings, status, or full info parameter

**get\_frame\_format** ()

Get format for the returned images.

Can be "list" (list of 2D arrays), "array" (a single 3D array), or "chunks" (list of 3D "chunk" arrays; supported for some cameras and provides the best performance).

**get\_frame\_info\_fields** ()

Get the names of frame info fields.

Applicable when frame info format (set by `set_frame_info_format()`) is "list" or "array".

**get\_frame\_info\_format** ()

Get format of the frame info.

Can be "namedtuple" (potentially nested named tuples; convenient to get particular values), "list" (flat list of values, with field names are given by `get_frame_info_fields()`; convenient for building a table), "array" (same as "list", but with a numpy array, which is easier to use for "chunks" frame format), or "dict" (flat dictionary with the same fields as the "list" format; more resilient to future format changes)

**get\_frame\_info\_period** ()

Get period of frame info acquisition.

Frame info might be skipped (set to `None`) except for frames which indices are divisible by *period*. Useful for certain cameras where acquiring frame info takes a lot of time and can reduce performance at higher frame rates. Note that this parameter can still be ignored (i.e., always set to 1) if the performance is not an issue for a given camera class.

**get\_frame\_period** ()

Get frame period (time between two consecutive frames in the internal trigger mode)

**get\_frames\_status** ()

Get acquisition and buffer status.

Return tuple (`acquired`, `unread`, `skipped`, `size`), where `acquired` is the total number of acquired frames, `unread` is the number of acquired but not read frames, `skipped` is the number of skipped (not read and then written over) frames, and `buffer_size` is the total buffer size (in frames).

**get\_full\_info** (*include=0*)

Get dict {`name`: `value`} containing full device information (including status and settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting `include=-10` queries all available variables, which is equivalent to `include="all"`.

**get\_full\_status** (*include=0*)

Get dict {`name`: `value`} containing the device status (including settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include*=-10 queries all available variables, which is equivalent to *include*="all".

#### **get\_image\_indexing()**

Get indexing for the returned images.

Can be "rct" (first index row, second index column, rows counted from the top), "rcb" (same as "rc", rows counted from the bottom), "xyt" (first index column, second index row, rows counted from the top), or "xyb" (same as "xyt", rows counted from the bottom)

#### **get\_new\_images\_range()**

Get the range of the new images.

Return tuple (*first*, *last*) with images range (first inclusive). If no images are available, return None. If some images were in the buffer were overwritten, exclude them from the range.

#### **get\_settings(include=0)**

Get dict {*name*: *value*} containing all the device settings.

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include*=-10 queries all available variables, which is equivalent to *include*="all".

#### **grab(nframes=1, frame\_timeout=5.0, missing\_frame='skip', return\_info=False, buff\_size=None)**

Snap *nframes* images (with preset image read mode parameters)

*buff\_size* determines buffer size (if None, use the default size). Timeout is specified for a single-frame acquisition, not for the whole acquisition time. *missing\_frame* determines what to do with frames which have been lost: can be "none" (replacing them with None), "zero" (replacing them with zero-filled frame), or "skip" (skipping them, while still keeping total returned frames number to *n*). If *return\_info*==True, return tuple (*frames*, *infos*), where *infos* is a list of frame info tuples (camera-dependent); if some frames are missing and *missing\_frame*!="skip", the corresponding frame info is None.

#### **is\_acquisition\_setup()**

Check if acquisition is set up.

If the camera does not support separate acquisition setup, always return True.

#### **pausing\_acquisition(clear=None, stop=True, setup\_after=None, start\_after=True, combine\_nested=True)**

Context manager which temporarily pauses acquisition during execution of with block.

Useful for applying certain settings which can't be changed during the acquisition. If *clear*==True, clear acquisition in addition to stopping (by default, use the class default specified as *\_clear\_pausing\_acquisition* attribute). If *stop*==True, stop the acquisition (if *clear*==True, stop regardless). If *setup\_after*==True, setup the acquisition after pause if necessary (None means setup only if clearing was required). If *start\_after*==True, start the acquisition after pause if necessary (None means start only if stopping was required). If *combine\_nested*==True, then any nested *pausing\_acquisition* calls will stop/clear acquisition as necessary, but won't setup/start it again until this *pausing\_acquisition* call is complete.

Yields tuple (*acq\_in\_progress*, *acq\_params*), which indicates whether acquisition is currently in progress, and what are the current acquisition parameters.

#### **read\_newest\_image(peek=False, return\_info=False)**

Read the newest un-read image.

If no un-read frames are available, return `None`. If `peek==True`, return the image but not mark it as read. If `return_info==True`, return tuple (frame, info), where info is an info tuples (camera-dependent, see [read\\_multiple\\_images\(\)](#)).

**read\_oldest\_image** (*peek=False, return\_info=False*)

Read the oldest un-read image.

If no un-read frames are available, return `None`. If `peek==True`, return the image but not mark it as read. If `return_info==True`, return tuple (frame, info), where info is an info tuples (camera-dependent, see [read\\_multiple\\_images\(\)](#)).

**set\_device\_variable** (*key, value*)

Set the value of a settings parameter

**set\_frame\_format** (*fmt*)

Set format for the returned images.

Can be "list" (list of 2D arrays), "array" (a single 3D array), "chunks" (list of 3D "chunk" arrays; supported for some cameras and provides the best performance), or "try\_chunks" (same as "chunks", but if chunks are not supported, set to "list" instead). If format is "chunks" and chunks are not supported by the camera, it results in one frame per chunk. Note that if the format is set to "array" or "chunks", the frame info format is also automatically set to "array". If the format is set to "chunks", then the image info is also returned in chunks form (list of 2D info arrays with the same length as the corresponding frame chunks).

**set\_frame\_info\_format** (*fmt, include\_fields=None*)

Set format of the frame info.

Can be "namedtuple" (potentially nested named tuples; convenient to get particular values), "list" (flat list of values, with field names are given by [get\\_frame\\_info\\_fields\(\)](#); convenient for building a table), "array" (same as "list", but with a numpy array, which is easier to use for "chunks" frame format), or "dict" (flat dictionary with the same fields as the "list" format; more resilient to future format changes) If *include\_fields* is not `None`, it specifies the fields included for non-"tuple" formats; note that order or *include\_fields* is ignored, and the resulting fields are always ordered same as in the original.

**set\_frame\_info\_period** (*period=1*)

Set period of frame info acquisition.

Frame info might be skipped (set to `None`) except for frames which indices are divisible by *period*. Useful for certain cameras where acquiring frame info takes a lot of time and can reduce performance at higher frame rates. Note that this parameter can still be ignored (i.e., always set to 1) if the performance is not an issue for a given camera class.

**set\_image\_indexing** (*indexing*)

Set up indexing for the returned images.

Can be "rct" (first index row, second index column, rows counted from the top), "rcb" (same as "rc", rows counted from the bottom), "xyt" (first index column, second index row, rows counted from the top), or "xyb" (same as "xyt", rows counted from the bottom)

**snap** (*timeout=5.0, return\_info=False*)

Snap a single frame

**wait\_for\_frame** (*since='lastread', nframes=1, timeout=20.0, error\_on\_stopped=False*)

Wait for one or several new camera frames.

*since* specifies the reference point for waiting to acquire *nframes* frames; can be "lastread" (from the last read frame), "lastwait" (wait for the last successful [wait\\_for\\_frame\(\)](#) call), "now" (from the start of the current call), or "start" (from the acquisition start, i.e., wait until *nframes* frames have been acquired). *timeout* can be either a number, `None` (infinite timeout), or a tuple



(*timeout*, *frame\_timeout*), in which case the call times out if the total time exceeds *timeout*, or a single frame wait exceeds *frame\_timeout*. If the call times out, raise `TimeoutError`. If *error\_on\_stopped*==True and the acquisition is not running, raise `Error`; otherwise, simply return False without waiting.

`pylablib.devices.Photometrics.pvcam.get_roi_parameters(buffer)`

Extract ROI parameters from the buffer.

*buffer* is the buffer represented as bytes numpy byte array. Return numpy array with one row per ROI and 4 columns: data offset from the frame start, data bytes per pixel, ROI height, and ROI width.

`pylablib.devices.Photometrics.pvcam.parse_metainfo_v1(buffer, nframes, stride)`

Extract frames metainfo for frames with v1 or v2 header.

*buffer* is the buffer represented as bytes numpy byte array, *nframes* is the number of frames in it, and *stride* is the frame stride (in bytes).

Return a 2D array with *nframes* rows and 7 columns: *framestamp*, *timestampBOF*, *timestampEOF*, *timestampRes*, *exposure*, *exposureRes*, *flags*.

`pylablib.devices.Photometrics.pvcam.parse_metainfo_v3(buffer, nframes, stride)`

Extract frames metainfo for frames with v3 header.

*buffer* is the buffer represented as bytes numpy byte array, *nframes* is the number of frames in it, and *stride* is the frame stride (in bytes).

Return a 2D array with *nframes* rows and 5 columns: *framestamp*, *timestampBOF*, *timestampEOF*, *exposure*, *flags*.

## Module contents

### pylablib.devices.PhotonFocus package

#### Submodules

#### pylablib.devices.PhotonFocus.PhotonFocus module

**class** `pylablib.devices.PhotonFocus.PhotonFocus.LibraryController(lib)`

Bases: `pylablib.devices.utils.load_lib.LibraryController`

**close**(*opid*)

Mark device closing.

Return tuple (*close\_result*, *uninit\_result*) with the results of the closing and the shutdown. If library does not need to be shut down yet, set *uninit\_result*=None

**open**()

Mark device opening.

Return tuple (*init\_result*, *open\_result*, *opid*) with the results of the initialization and the opening, and the opening ID which should afterwards be used for closing. If library is already initialized, set *init\_result*=None

**preinit**()

Pre-initialize the library, if it hasn't been done already

**shutdown**()

Close all opened connections and shutdown the library

**temp\_open()**

Context for temporarily opening a new device connection

`pylablib.devices.PhotonFocus.PhotonFocus.query_camera_name(port)`

Query cameras name at a given port in PFCam interface

**class** `pylablib.devices.PhotonFocus.PhotonFocus.TCameraInfo` (*manufacturer, port, version, type*)

Bases: `tuple`

**count()**

Return number of occurrences of value.

**index()**

Return first index of value.

Raises `ValueError` if the value is not present.

**manufacturer**

**port**

**type**

**version**

`pylablib.devices.PhotonFocus.PhotonFocus.list_cameras(only_supported=True)`

List all cameras available through PFCam interface.

If `only_supported==True`, only return cameras which support PFCam protocol (this check only works if the camera is not currently accessed by some other software). Return a list `[(port, info)]`, where `port` is the pfcam port given to `IPhotonFocusCamera` and its subclasses, and `info` is the information returned by `query_camera_name()`.

`pylablib.devices.PhotonFocus.PhotonFocus.get_cameras_number(only_supported=True)`

Get the total number of connected PFCam cameras

`pylablib.devices.PhotonFocus.PhotonFocus.get_port_index(manufacturer, port)`

Find PhotonFocus port index based on the manufacturer and port

**class** `pylablib.devices.PhotonFocus.PhotonFocus.PFCamAttribute` (*port, name*)

Bases: `object`

PFCam camera attribute.

Allows to query and set values and get additional information. Usually created automatically by a PhotonFocus camera instance, but could also be created manually.

#### Parameters

- **sid** – camera session ID
- **name** – attribute text name

**name**

attribute name

**kind**

attribute kind; can be "INT", "UINT", "FLOAT", "BOOL", "MODE", "STRING", or "COMMAND"

**readable**

whether attribute is readable

**Type** `bool`

**writable**  
whether attribute is writable  
Type `bool`

**is\_command**  
whether attribute is a command  
Type `bool`

**min**  
minimal attribute value (if applicable)  
Type `float` or `int`

**max**  
maximal attribute value (if applicable)  
Type `float` or `int`

**ivalues**  
list of possible integer values for enum attributes

**values**  
list of possible text values for enum attributes

**labels**  
dict {label: index} which shows all possible values of an enumerated attribute and their corresponding numerical values

**ilabels**  
dict {index: label} which shows labels corresponding to numerical values of an enumerated attribute

**update\_limits()**  
Update minimal and maximal attribute limits and return tuple (min, max)

**truncate\_value(value)**  
Truncate value to lie within attribute limits

**get\_value(enum\_as\_str=True)**  
Get attribute value.  
  
If `enum_as_str==True`, return enum-style values as strings; otherwise, return corresponding integer values.

**set\_value(value, truncate=True)**  
Get attribute value.  
  
If `truncate==True`, automatically truncate value to lie within allowed range.

**call\_command(arg=0)**  
If attribute is a command, call it with a given argument; otherwise, raise an error

**class** `pylablib.devices.PhotonFocus.PhotonFocus.TDeviceInfo` (*model, serial\_number, grabber\_info*)

Bases: `tuple`

**count()**  
Return number of occurrences of value.

**grabber\_info**

**index()**  
Return first index of value.

Raises ValueError if the value is not present.

**model**

**serial\_number**

```
class pylablib.devices.PhotonFocus.PhotonFocus.IPhotonFocusCamera (pfcam_port=0,  
                                                                    **kwargs)
```

Bases: `pylablib.devices.interface.camera.IAttributeCamera`

Generic PFCam interface to a PhotonFocus camera. Does not handle frames acquisition, so needs to be mixed with a frame grabber class to be fully operational. In this mixing, the class attribute `GrabberClass` should be set to this frame grabber class.

#### Parameters

- **pfcam\_port** – port number for pfcam interface (can be learned by `list_cameras()`; port number is the first element of the camera data tuple) can also be a tuple (manufacturer, port), e.g., ("National Instruments", "port0").
- **kwargs** – keyword arguments passed to the frame grabber initialization

#### Error

alias of `pylablib.core.devio.base.DeviceError`

**GrabberClass** = None

**setup\_max\_baudrate**()

Setup the maximal available baudrate

**get\_baudrate**()

Get the current baud rate

**open**()

Open connection to the camera

**close**()

Close connection to the camera

**get\_attribute\_value**(name, enum\_as\_str=True, error\_on\_missing=True, default=None)

Get value of an attribute with the given name.

If the value doesn't exist or can not be read and `error_on_missing==True`, raise error; otherwise, return *default*. If *default* is not None, assume that `error_on_missing==False`. If `enum_as_str==True`, try to represent enums as their string values; If *name* points at a dictionary branch, return a dictionary with all values in this branch.

**set\_attribute\_value**(name, value, truncate=True, error\_on\_missing=True)

Set value of an attribute with the given name.

If the value doesn't exist or can not be written and `error_on_missing==True`, raise error; otherwise, do nothing. If *name* points at a dictionary branch, set all values in this branch (in this case *value* must be a dictionary). If `truncate==True`, truncate value to lie within attribute range.

**get\_all\_attribute\_values**(root="", enum\_as\_str=True)

Get values of all attributes with the given *root*

**set\_all\_attribute\_values**(settings, root="", truncate=True)

Set values of all attributes with the given *root*.

If `truncate==True`, truncate value to lie within attribute range.

**update\_attribute\_value** (*name, value, error\_on\_missing=True, truncate=True*)

Set value of the attribute with a given name, but only if it's different from the current value.

Can take less time on some version of PFRremote (where single attribute setting is about 50ms). Arguments are the same as `set_attribute_value()`.

**call\_command** (*name, arg=0, error\_on\_missing=True*)

Execute the given command with the given argument.

If the command doesn't exist and `error_on_missing==True`, raise error; otherwise, do nothing.

**get\_device\_info** ()

Get camera model data.

Return tuple (model, serial\_number, grabber\_info).

**get\_detector\_size** ()

Get camera detector size (in pixels) as a tuple (width, height)

**get\_roi** ()

Get current ROI.

Return tuple (hstart, hend, vstart, vend).

**fast\_shift\_roi** (*hstart=0, vstart=0*)

Shift ROI by only changing its origin, but keeping the shape the same.

Note that if the ROI is invalid, it won't be truncated (as is the standard behavior of `set_roi()`), which might lead to errors later.

**set\_roi** (*hstart=0, hend=None, vstart=0, vend=None*)

Setup camera ROI.

By default, all non-supplied parameters take extreme values.

**get\_roi\_limits** (*hbin=1, vbin=1*)

**get\_exposure** ()

Get current exposure

**set\_exposure** (*exposure*)

Set current exposure

**get\_frame\_period** ()

Get frame period (time between two consecutive frames in the internal trigger mode)

**set\_frame\_period** (*frame\_period*)

Set frame period (time between two consecutive frames in the internal trigger mode)

**get\_frame\_timings** ()

Get acquisition timing.

Return tuple (exposure, frame\_period).

**is\_CFR\_enabled** ()

Check if the constant frame rate mode is enabled

**enable\_CFR** (*enabled=True*)

Enable constant frame rate mode

**get\_trigger\_interleave** ()

Check if the trigger interleave is on

**set\_trigger\_interleave** (*enabled*)

Set the trigger interleave option on or off

**is\_status\_line\_enabled()**

Check if the status line is on

**enable\_status\_line** (*enabled=True*)

Enable or disable status line

**get\_black\_level\_offset()**

Get the black level offset

**set\_black\_level\_offset** (*offset*)

Set the black level offset

**FrameTransferError**

alias of *pylablib.devices.interface.camera.DefaultFrameTransferError*

**class NoParameterCaller** (*device, kind*)

Bases: *object*

Class to simplify calling functions without a parameter

**TimeoutError**

alias of *pylablib.core.devio.base.DeviceError*

**acquisition\_in\_progress()**

Check if acquisition is in progress

**apply\_settings** (*settings*)

Apply the settings.

*settings* is the dict {name: value} of the device available settings. Non-applicable settings are ignored.

**clear\_acquisition()**

Clear acquisition settings

**get\_acquisition\_parameters()**

Get acquisition parameters.

Return dictionary {name: value}

**get\_all\_attributes** (*copy=False*)

Return a dictionary of all available attributes.

If *copy==True*, copy the dictionary; otherwise, return the internal dictionary structure (should not be modified).

**get\_attribute** (*name, error\_on\_missing=True*)

Get the camera attribute with the given name

**get\_data\_dimensions()**

Get readout data dimensions (in pixels) as a tuple (width, height); take indexing mode into account

**get\_device\_variable** (*key*)

Get the value of a settings, status, or full info parameter

**get\_frame\_format()**

Get format for the returned images.

Can be "list" (list of 2D arrays), "array" (a single 3D array), or "chunks" (list of 3D "chunk" arrays; supported for some cameras and provides the best performance).

**get\_frame\_info\_fields()**

Get the names of frame info fields.

Applicable when frame info format (set by *set\_frame\_info\_format()*) is "list" or "array".

**get\_frame\_info\_format()**

Get format of the frame info.

Can be "namedtuple" (potentially nested named tuples; convenient to get particular values), "list" (flat list of values, with field names are given by `get_frame_info_fields()`; convenient for building a table), "array" (same as "list", but with a numpy array, which is easier to use for "chunks" frame format), or "dict" (flat dictionary with the same fields as the "list" format; more resilient to future format changes)

**get\_frame\_info\_period()**

Get period of frame info acquisition.

Frame info might be skipped (set to `None`) except for frames which indices are divisible by *period*. Useful for certain cameras where acquiring frame info takes a lot of time and can reduce performance at higher frame rates. Note that this parameter can still be ignored (i.e., always set to 1) if the performance is not an issue for a given camera class.

**get\_frames\_status()**

Get acquisition and buffer status.

Return tuple (*acquired*, *unread*, *skipped*, *size*), where *acquired* is the total number of acquired frames, *unread* is the number of acquired but not read frames, *skipped* is the number of skipped (not read and then written over) frames, and *buffer\_size* is the total buffer size (in frames).

**get\_full\_info(include=0)**

Get dict {*name*: *value*} containing full device information (including status and settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include*=-10 queries all available variables, which is equivalent to *include*="all".

**get\_full\_status(include=0)**

Get dict {*name*: *value*} containing the device status (including settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include*=-10 queries all available variables, which is equivalent to *include*="all".

**get\_image\_indexing()**

Get indexing for the returned images.

Can be "rct" (first index row, second index column, rows counted from the top), "rcb" (same as "rc", rows counted from the bottom), "xyt" (first index column, second index row, rows counted from the top), or "xyb" (same as "xyt", rows counted from the bottom)

**get\_new\_images\_range()**

Get the range of the new images.

Return tuple (*first*, *last*) with images range (first inclusive). If no images are available, return `None`. If some images were in the buffer were overwritten, exclude them from the range.

**get\_settings(include=0)**

Get dict {*name*: *value*} containing all the device settings.

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include*=-10 queries all available variables, which is equivalent to *include*="all".

**grab** (*nframes=1, frame\_timeout=5.0, missing\_frame='skip', return\_info=False, buff\_size=None*)  
Snap *nframes* images (with preset image read mode parameters)

*buff\_size* determines buffer size (if *None*, use the default size). Timeout is specified for a single-frame acquisition, not for the whole acquisition time. *missing\_frame* determines what to do with frames which have been lost: can be "none" (replacing them with *None*), "zero" (replacing them with zero-filled frame), or "skip" (skipping them, while still keeping total returned frames number to *n*). If *return\_info==True*, return tuple (*frames, infos*), where *infos* is a list of frame info tuples (camera-dependent); if some frames are missing and *missing\_frame!="skip"*, the corresponding frame info is *None*.

**is\_acquisition\_setup** ()  
Check if acquisition is set up.

If the camera does not support separate acquisition setup, always return *True*.

**is\_opened** ()  
Check if the device is connected

**pausing\_acquisition** (*clear=None, stop=True, setup\_after=None, start\_after=True, combine\_nested=True*)  
Context manager which temporarily pauses acquisition during execution of *with* block.

Useful for applying certain settings which can't be changed during the acquisition. If *clear==True*, clear acquisition in addition to stopping (by default, use the class default specified as *\_clear\_pausing\_acquisition* attribute). If *stop==True*, stop the acquisition (if *clear==True*, stop regardless). If *setup\_after==True*, setup the acquisition after pause if necessary (*None* means setup only if clearing was required). If *start\_after==True*, start the acquisition after pause if necessary (*None* means start only if stopping was required). If *combine\_nested==True*, then any nested *pausing\_acquisition* calls will stop/clear acquisition as necessary, but won't setup/start it again until this *pausing\_acquisition* call is complete.

Yields tuple (*acq\_in\_progress, acq\_params*), which indicates whether acquisition is currently in progress, and what are the current acquisition parameters.

**read\_multiple\_images** (*rng=None, peek=False, missing\_frame='skip', return\_info=False, return\_rng=False*)  
Read multiple images specified by *rng* (by default, all un-read images).

If *rng* is specified, it is a tuple (*first, last*) with images range (first inclusive). If no new frames are available, return an empty list; if no acquisition is running, return *None*. If *peek==True*, return images but not mark them as read. *missing\_frame* determines what to do with frames which are out of range (missing or lost): can be "none" (replacing them with *None*), "zero" (replacing them with zero-filled frame), or "skip" (skipping them). If *return\_info==True*, return tuple (*frames, infos*), where *infos* is a list of frame info tuples (camera-dependent, by default, only the frame index); if some frames are missing and *missing\_frame!="skip"*, the corresponding frame info is *None*. If *return\_rng==True*, return the range covered resulting frames; if *missing\_frame=="skip"*, the range can be smaller than the supplied *rng* if some frames are skipped.

**read\_newest\_image** (*peek=False, return\_info=False*)  
Read the newest un-read image.

If no un-read frames are available, return *None*. If *peek==True*, return the image but not mark it as read. If *return\_info==True*, return tuple (*frame, info*), where *info* is an info tuples (camera-dependent, see [read\\_multiple\\_images\(\)](#)).

**read\_oldest\_image** (*peek=False, return\_info=False*)  
Read the oldest un-read image.

If no un-read frames are available, return *None*. If *peek==True*, return the image but not mark it as read. If *return\_info==True*, return tuple (*frame, info*), where *info* is an info tuples



(camera-dependent, see `read_multiple_images()`).

**set\_device\_variable** (*key, value*)  
Set the value of a settings parameter

**set\_frame\_format** (*fmt*)  
Set format for the returned images.

Can be "list" (list of 2D arrays), "array" (a single 3D array), "chunks" (list of 3D "chunk" arrays; supported for some cameras and provides the best performance), or "try\_chunks" (same as "chunks", but if chunks are not supported, set to "list" instead). If format is "chunks" and chunks are not supported by the camera, it results in one frame per chunk. Note that if the format is set to "array" or "chunks", the frame info format is also automatically set to "array". If the format is set to "chunks", then the image info is also returned in chunks form (list of 2D info arrays with the same length as the corresponding frame chunks).

**set\_frame\_info\_format** (*fmt, include\_fields=None*)  
Set format of the frame info.

Can be "namedtuple" (potentially nested named tuples; convenient to get particular values), "list" (flat list of values, with field names are given by `get_frame_info_fields()`; convenient for building a table), "array" (same as "list", but with a numpy array, which is easier to use for "chunks" frame format), or "dict" (flat dictionary with the same fields as the "list" format; more resilient to future format changes) If *include\_fields* is not None, it specifies the fields included for non-"tuple" formats; note that order or *include\_fields* is ignored, and the resulting fields are always ordered same as in the original.

**set\_frame\_info\_period** (*period=1*)  
Set period of frame info acquisition.

Frame info might be skipped (set to None) except for frames which indices are divisible by *period*. Useful for certain cameras where acquiring frame info takes a lot of time and can reduce performance at higher frame rates. Note that this parameter can still be ignored (i.e., always set to 1) if the performance is not an issue for a given camera class.

**set\_image\_indexing** (*indexing*)  
Set up indexing for the returned images.

Can be "rct" (first index row, second index column, rows counted from the top), "rcb" (same as "rc", rows counted from the bottom), "xyt" (first index column, second index row, rows counted from the top), or "xyb" (same as "xyt", rows counted from the bottom)

**setup\_acquisition** (*\*\*kwargs*)  
Setup acquisition.

Any non-specified acquisition parameters are assumed to be the same as previously set (or default, if not explicitly set before). Return the new acquisition parameters.

**snap** (*timeout=5.0, return\_info=False*)  
Snap a single frame

**start\_acquisition** (*\*args, \*\*kwargs*)  
Start acquisition.

Can take the same keyword parameters as *meth: "setup\_acquisition"*. If the acquisition is not set up yet, set it up using the supplied parameters (use default of `setup_acquisition()`, if the parameter is None). Otherwise, if any supplied parameters are different from the current ones, change them and reset the acquisition.

**stop\_acquisition** ()  
Stop acquisition

**wait\_for\_frame** (*since='lastread', nframes=1, timeout=20.0, error\_on\_stopped=False*)

Wait for one or several new camera frames.

*since* specifies the reference point for waiting to acquire *nframes* frames; can be "lastread" (from the last read frame), "lastwait" (wait for the last successful `wait_for_frame()` call), "now" (from the start of the current call), or "start" (from the acquisition start, i.e., wait until *nframes* frames have been acquired). *timeout* can be either a number, None (infinite timeout), or a tuple (timeout, frame\_timeout), in which case the call times out if the total time exceeds timeout, or a single frame wait exceeds frame\_timeout. If the call times out, raise `TimeoutError`. If `error_on_stopped==True` and the acquisition is not running, raise `Error`; otherwise, simply return False without waiting.

**class** `pylablib.devices.PhotonFocus.PhotonFocus.PhotonFocusIMAQCamera` (*imaq\_name='img0', pf-  
cam\_port=0*)

Bases: `pylablib.devices.PhotonFocus.PhotonFocus.IPhotonFocusCamera`, `pylablib.devices.IMAQ.IMAQ.IMAQFrameGrabber`

IMAQ+PFCam interface to a PhotonFocus camera.

#### Parameters

- **imaq\_name** – IMAQ interface name (can be learned by `IMAQ.list_cameras()`; usually, but not always, starts with "img")
- **pfcam\_port** – port number for pfcam interface (can be learned by `list_cameras()`; port number is the first element of the camera data tuple) can also be a tuple (manufacturer, port), e.g., ("National Instruments", "port0").

#### Error

alias of `pylablib.core.devio.base.DeviceError`

#### GrabberClass

alias of `pylablib.devices.IMAQ.IMAQ.IMAQFrameGrabber`

#### open()

Open connection to the camera

#### FrameTransferError

alias of `pylablib.devices.interface.camera.DefaultFrameTransferError`

**class** `NoParameterCaller` (*device, kind*)

Bases: `object`

Class to simplify calling functions without a parameter

**TimeoutError** = `<Mock spec='str' id='140096274107088'>`

**acquisition\_in\_progress()**

Check if acquisition is in progress

**apply\_settings** (*settings*)

Apply the settings.

*settings* is the dict {name: value} of the device available settings. Non-applicable settings are ignored.

**call\_command** (*name, arg=0, error\_on\_missing=True*)

Execute the given command with the given argument.

If the command doesn't exist and `error_on_missing==True`, raise error; otherwise, do nothing.

**clear\_acquisition()**

Clear all acquisition details and free all buffers

**clear\_all\_triggers** (*reset\_acquisition=True*)

Disable all triggers of the session

If the input triggers configuration has been changed, acquisition needs to be restart; if `reset_acquisition==True`, perform it automatically.

**close()**

Close connection to the camera

**configure\_trigger\_in** (*trig\_type*, *trig\_line=0*, *trig\_pol='high'*, *trig\_action='none'*, *timeout=None*, *reset\_acquisition=True*)

Configure input trigger.

#### Parameters

- **trig\_type** (*str*) – trigger source type; can be "ext", "rtsi", "iso\_in", or "software"
- **trig\_line** (*int*) – trigger line number
- **trig\_pol** (*str*) – trigger polarity; can be "high" or "low"
- **trig\_action** (*str*) – trigger action; can be "none" (disable trigger), "capture" (start capturing), "stop" (stop capturing), "buffer" (capture a single frame), or "bufflist" (capture the whole buffer list once)
- **timeout** (*float*) – timeout in seconds; None means not timeout.
- **reset\_acquisition** (*bool*) – if the input triggers configuration has been changed, acquisition needs to be restart; if True, perform it automatically

**configure\_trigger\_out** (*trig\_type*, *trig\_line=0*, *trig\_pol='high'*, *trig\_drive='disable'*)

Configure trigger output.

#### Parameters

- **trig\_type** (*str*) – trigger drive destination type; can be "ext", "rtsi", or "iso\_out"
- **trig\_line** (*int*) – trigger line number
- **trig\_pol** (*str*) – trigger polarity; can be "high" or "low"
- **trig\_drive** (*str*) – trigger output signal; can be "disable" (disable drive), "acq\_in\_progress" (asserted when acquisition is started), "acq\_done" (asserted when acquisition is done), "unasserted" (force unasserted level), "asserted" (force asserted level), "hsync" (asserted on start of a single line start), "vsync" (asserted on start of a frame scan), "frame\_start" (asserted when a single frame is captured), or "frame\_done" (asserted when a single frame is done)

**enable\_CFR** (*enabled=True*)

Enable constant frame rate mode

**enable\_status\_line** (*enabled=True*)

Enable or disable status line

**fast\_shift\_roi** (*hstart=0*, *vstart=0*)

Shift ROI by only changing its origin, but keeping the shape the same.

Note that if the ROI is invalid, it won't be truncated (as is the standard behavior of `set_roi()`), which might lead to errors later.

**get\_acquisition\_parameters()**

Get acquisition parameters.

Return dictionary {name: value}

**get\_all\_attribute\_values** (*root=""*, *enum\_as\_str=True*)

Get values of all attributes with the given *root*

**get\_all\_attributes** (*copy=False*)

Return a dictionary of all available attributes.

If *copy==True*, copy the dictionary; otherwise, return the internal dictionary structure (should not be modified).

**get\_all\_grabber\_attribute\_values()**

Get a dictionary of all readable attributes.

The attributes types are autodetected, and some of the types of uncommon attributes may be misrepresented.

**get\_attribute** (*name*, *error\_on\_missing=True*)

Get the camera attribute with the given name

**get\_attribute\_value** (*name*, *enum\_as\_str=True*, *error\_on\_missing=True*, *default=None*)

Get value of an attribute with the given name.

If the value doesn't exist or can not be read and *error\_on\_missing==True*, raise error; otherwise, return *default*. If *default* is not None, assume that *error\_on\_missing==False*. If *enum\_as\_str==True*, try to represent enums as their string values; If *name* points at a dictionary branch, return a dictionary with all values in this branch.

**get\_baudrate()**

Get the current baud rate

**get\_black\_level\_offset()**

Get the black level offset

**get\_data\_dimensions()**

Get readout data dimensions (in pixels) as a tuple (*width*, *height*); take indexing mode into account

**get\_detector\_size()**

Get camera detector size (in pixels) as a tuple (*width*, *height*)

**get\_device\_info()**

Get camera model data.

Return tuple (*model*, *serial\_number*, *grabber\_info*).

**get\_device\_variable** (*key*)

Get the value of a settings, status, or full info parameter

**get\_exposure()**

Get current exposure

**get\_frame\_format()**

Get format for the returned images.

Can be "list" (list of 2D arrays), "array" (a single 3D array), or "chunks" (list of 3D "chunk" arrays; supported for some cameras and provides the best performance).

**get\_frame\_info\_fields()**

Get the names of frame info fields.

Applicable when frame info format (set by `set_frame_info_format()`) is "list" or "array".

**get\_frame\_info\_format()**

Get format of the frame info.

Can be "namedtuple" (potentially nested named tuples; convenient to get particular values), "list" (flat list of values, with field names are given by `get_frame_info_fields()`; convenient for building a table), "array" (same as "list", but with a numpy array, which is easier to use for "chunks" frame format), or "dict" (flat dictionary with the same fields as the "list" format; more resilient to future format changes)

**get\_frame\_info\_period()**

Get period of frame info acquisition.

Frame info might be skipped (set to None) except for frames which indices are divisible by *period*. Useful for certain cameras where acquiring frame info takes a lot of time and can reduce performance at higher frame rates. Note that this parameter can still be ignored (i.e., always set to 1) if the performance is not an issue for a given camera class.

**get\_frame\_period()**

Get frame period (time between two consecutive frames in the internal trigger mode)

**get\_frame\_timings()**

Get acquisition timing.

Return tuple (exposure, frame\_period).

**get\_frames\_status()**

Get acquisition and buffer status.

Return tuple (acquired, unread, skipped, size), where *acquired* is the total number of acquired frames, *unread* is the number of acquired but not read frames, *skipped* is the number of skipped (not read and then written over) frames, and *buffer\_size* is the total buffer size (in frames).

**get\_full\_info(include=0)**

Get dict {name: value} containing full device information (including status and settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include*=-10 queries all available variables, which is equivalent to *include*="all".

**get\_full\_status(include=0)**

Get dict {name: value} containing the device status (including settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include*=-10 queries all available variables, which is equivalent to *include*="all".

**get\_grabber\_attribute\_value(attr, default=None, kind='auto')**

Get value of an attribute with a given name or index.

If *default* is not None, return *default* if the attribute is not supported; otherwise, raise an error. *kind* is the attribute kind, and it can be "uint32", "uint64", "double", or "auto" (autodetect based on the stored list of attribute kinds).

**get\_grabber\_detector\_size()**

Get camera detector size (in pixels) as a tuple (width, height)

**get\_grabber\_roi()**

Get current ROI.

Return tuple (*hstart*, *hend*, *vstart*, *vend*). *hstart* and *hend* specify horizontal image extent, *vstart* and *vend* specify vertical image extent (start is inclusive, stop is exclusive, starting from 0).

**get\_grabber\_roi\_limits** (*hbin=1*, *vbin=1*)

Get the minimal and maximal ROI parameters.

Return tuple (*hlim*, *vlim*), where each element is in turn a limit 5-tuple (*min*, *max*, *pstep*, *sstep*, *maxbin*) with, correspondingly, minimal and maximal size, position and size step, and the maximal binning (fixed to 1 if not binning is allowed). In some cameras, the step and the minimal size depend on the binning, which can be supplied.

**get\_image\_indexing** ()

Get indexing for the returned images.

Can be "rct" (first index row, second index column, rows counted from the top), "rcb" (same as "rc", rows counted from the bottom), "xyt" (first index column, second index row, rows counted from the top), or "xyb" (same as "xyt", rows counted from the bottom)

**get\_new\_images\_range** ()

Get the range of the new images.

Return tuple (*first*, *last*) with images range (first inclusive). If no images are available, return None. If some images were in the buffer were overwritten, exclude them from the range.

**get\_roi** ()

Get current ROI.

Return tuple (*hstart*, *hend*, *vstart*, *vend*).

**get\_roi\_limits** (*hbin=1*, *vbin=1*)

**get\_serial\_params** ()

Return serial parameters as a tuple (*write\_term*, *datatype*)

**get\_settings** (*include=0*)

Get dict {*name*: *value*} containing all the device settings.

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_trigger\_interleave** ()

Check if the trigger interleave is on

**grab** (*nframes=1*, *frame\_timeout=5.0*, *missing\_frame='skip'*, *return\_info=False*, *buff\_size=None*)

Snap *nframes* images (with preset image read mode parameters)

*buff\_size* determines buffer size (if None, use the default size). Timeout is specified for a single-frame acquisition, not for the whole acquisition time. *missing\_frame* determines what to do with frames which have been lost: can be "none" (replacing them with None), "zero" (replacing them with zero-filled frame), or "skip" (skipping them, while still keeping total returned frames number to *n*). If *return\_info==True*, return tuple (*frames*, *infos*), where *infos* is a list of frame info tuples (camera-dependent); if some frames are missing and *missing\_frame!="skip"*, the corresponding frame info is None.

**is\_CFR\_enabled** ()

Check if the constant frame rate mode is enabled

**is\_acquisition\_setup** ()

Check if acquisition is set up.

If the camera does not support separate acquisition setup, always return True.

**is\_opened()**

Check if the device is connected

**is\_status\_line\_enabled()**

Check if the status line is on

**pausing\_acquisition** (*clear=None, stop=True, setup\_after=None, start\_after=True, combine\_nested=True*)

Context manager which temporarily pauses acquisition during execution of with block.

Useful for applying certain settings which can't be changed during the acquisition. If `clear==True`, clear acquisition in addition to stopping (by default, use the class default specified as `_clear_pausing_acquisition` attribute). If `stop==True`, stop the acquisition (if `clear==True`, stop regardless). If `setup_after==True`, setup the acquisition after pause if necessary (None means setup only if clearing was required). If `start_after==True`, start the acquisition after pause if necessary (None means start only if stopping was required). If `combine_nested==True`, then any nested `pausing_acquisition` calls will stop/clear acquisition as necessary, but won't setup/start it again until this `pausing_acquisition` call is complete.

Yields tuple (`acq_in_progress, acq_params`), which indicates whether acquisition is currently in progress, and what are the current acquisition parameters.

**read\_multiple\_images** (*rng=None, peek=False, missing\_frame='skip', return\_info=False, return\_rng=False*)

Read multiple images specified by *rng* (by default, all un-read images).

If *rng* is specified, it is a tuple (*first, last*) with images range (first inclusive). If no new frames are available, return an empty list; if no acquisition is running, return None. If `peek==True`, return images but not mark them as read. *missing\_frame* determines what to do with frames which are out of range (missing or lost): can be "none" (replacing them with None), "zero" (replacing them with zero-filled frame), or "skip" (skipping them). If `return_info==True`, return tuple (*frames, infos*), where *infos* is a list of frame info tuples (camera-dependent, by default, only the frame index); if some frames are missing and `missing_frame!="skip"`, the corresponding frame info is None. if `return_rng==True`, return the range covered resulting frames; if `missing_frame=="skip"`, the range can be smaller than the supplied *rng* if some frames are skipped.

**read\_newest\_image** (*peek=False, return\_info=False*)

Read the newest un-read image.

If no un-read frames are available, return None. If `peek==True`, return the image but not mark it as read. If `return_info==True`, return tuple (*frame, info*), where *info* is an info tuples (camera-dependent, see `read_multiple_images()`).

**read\_oldest\_image** (*peek=False, return\_info=False*)

Read the oldest un-read image.

If no un-read frames are available, return None. If `peek==True`, return the image but not mark it as read. If `return_info==True`, return tuple (*frame, info*), where *info* is an info tuples (camera-dependent, see `read_multiple_images()`).

**read\_trigger** (*trig\_type, trig\_line=0, trig\_pol='high'*)

Read current value of a trigger (input or output).

#### Parameters

- **trig\_type** (*str*) – trigger drive destination type; can be "ext", "rtsi", "iso\_in", or "iso\_out"
- **trig\_line** (*int*) – trigger line number
- **trig\_pol** (*str*) – trigger polarity; can be "high" or "low"

**reset()**

Reset connection to the camera

**send\_software\_trigger()**

Send software trigger signal

**serial\_flush()**

Flush CameraLink serial port

**serial\_read(*n*, *timeout*=3.0, *datatype*=None)**

Read specified number of bytes from CameraLink serial port.

#### Parameters

- **n** – number of bytes to read
- **timeout** – operation timeout (in seconds)
- **datatype** – return datatype; can be "bytes" (return raw bytes), or "str" (convert into UTF-8 string) if None, use the value set up using `setup_serial_params()` (by default, "bytes")

**serial\_readline(*timeout*=3.0, *datatype*=None, *maxn*=1024)**

Read bytes from CameraLink serial port until the termination character (defined in camera file) is encountered.

#### Parameters

- **timeout** – operation timeout (in seconds)
- **datatype** – return datatype; can be "bytes" (return raw bytes), or "str" (convert into UTF-8 string) if None, use the value set up using `setup_serial_params()` (by default, "bytes")
- **maxn** – maximal number of bytes to read

**serial\_write(*msg*, *timeout*=3.0, *term*=None)**

Write message into CameraLink serial port.

#### Parameters

- **timeout** – operation timeout (in seconds)
- **term** – additional write terminator character to add to the message; if None, use the value set up using `setup_serial_params()` (by default, no additional terminator)

**set\_all\_attribute\_values(*settings*, *root*=", *truncate*=True)**

Set values of all attributes with the given *root*.

If `truncate==True`, truncate value to lie within attribute range.

**set\_attribute\_value(*name*, *value*, *truncate*=True, *error\_on\_missing*=True)**

Set value of an attribute with the given name.

If the value doesn't exist or can not be written and `error_on_missing==True`, raise error; otherwise, do nothing. If *name* points at a dictionary branch, set all values in this branch (in this case *value* must be a dictionary). If `truncate==True`, truncate value to lie within attribute range.

**set\_black\_level\_offset(*offset*)**

Set the black level offset

**set\_device\_variable(*key*, *value*)**

Set the value of a settings parameter



**set\_exposure** (*exposure*)

Set current exposure

**set\_frame\_format** (*fmt*)

Set format for the returned images.

Can be "list" (list of 2D arrays), "array" (a single 3D array), "chunks" (list of 3D "chunk" arrays; supported for some cameras and provides the best performance), or "try\_chunks" (same as "chunks", but if chunks are not supported, set to "list" instead). If format is "chunks" and chunks are not supported by the camera, it results in one frame per chunk. Note that if the format is set to "array" or "chunks", the frame info format is also automatically set to "array". If the format is set to "chunks", then the image info is also returned in chunks form (list of 2D info arrays with the same length as the corresponding frame chunks).

**set\_frame\_info\_format** (*fmt, include\_fields=None*)

Set format of the frame info.

Can be "namedtuple" (potentially nested named tuples; convenient to get particular values), "list" (flat list of values, with field names are given by `get_frame_info_fields()`; convenient for building a table), "array" (same as "list", but with a numpy array, which is easier to use for "chunks" frame format), or "dict" (flat dictionary with the same fields as the "list" format; more resilient to future format changes) If *include\_fields* is not None, it specifies the fields included for non-"tuple" formats; note that order of *include\_fields* is ignored, and the resulting fields are always ordered same as in the original.

**set\_frame\_info\_period** (*period=1*)

Set period of frame info acquisition.

Frame info might be skipped (set to None) except for frames which indices are divisible by *period*. Useful for certain cameras where acquiring frame info takes a lot of time and can reduce performance at higher frame rates. Note that this parameter can still be ignored (i.e., always set to 1) if the performance is not an issue for a given camera class.

**set\_frame\_period** (*frame\_period*)

Set frame period (time between two consecutive frames in the internal trigger mode)

**set\_grabber\_attribute\_value** (*attr, value, kind='int32'*)

Set value of an attribute with a given name or index.

*kind* is the attribute kind, and it can be "uint32", "uint64", "double", or "auto" (autodetect based on the stored list of attribute kinds).

**set\_grabber\_roi** (*hstart=0, hend=None, vstart=0, vend=None*)

Setup camera ROI.

*hstart* and *hend* specify horizontal image extent, *vstart* and *vend* specify vertical image extent (start is inclusive, stop is exclusive, starting from 0). By default, all non-supplied parameters take extreme values (0 for start, maximal for end).

**set\_image\_indexing** (*indexing*)

Set up indexing for the returned images.

Can be "rct" (first index row, second index column, rows counted from the top), "rcb" (same as "rc", rows counted from the bottom), "xyt" (first index column, second index row, rows counted from the top), or "xyb" (same as "xyt", rows counted from the bottom)

**set\_roi** (*hstart=0, hend=None, vstart=0, vend=None*)

Setup camera ROI.

By default, all non-supplied parameters take extreme values.

**set\_trigger\_interleave** (*enabled*)

Set the trigger interleave option on or off

**setup\_acquisition** (*mode='sequence', nframes=100*)

Setup acquisition mode.

*mode* can be either "snap" (single frame or a fixed number of frames) or "sequence" (continuous acquisition). (note that *IMAQCamera.acquisition\_in\_progress()* would still return True in this case, even though new frames are no longer acquired). *nframes* sets up number of frame buffers.

**setup\_max\_baudrate** ()

Setup the maximal available baudrate

**setup\_serial\_params** (*write\_term=", datatype='bytes'*)

Setup default serial communication parameters.

#### Parameters

- **write\_term** – default terminator character to be added to the sent messages
- **datatype** – type of the result of read commands; can be "bytes" (return raw bytes), or "str" (convert into UTF-8 string)

**snap** (*timeout=5.0, return\_info=False*)

Snap a single frame

**start\_acquisition** (*\*args, \*\*kwargs*)

Start acquisition.

Can take the same keyword parameters as *meth: "setup\_acquisition"*. If the acquisition is not set up yet, set it up using the supplied parameters (use default of *setup\_acquisition()*, if the parameter is None). Otherwise, if any supplied parameters are different from the current ones, change them and reset the acquisition.

**stop\_acquisition** ()

Stop acquisition

**update\_attribute\_value** (*name, value, error\_on\_missing=True, truncate=True*)

Set value of the attribute with a given name, but only if it's different from the current value.

Can take less time on some version of PFRremote (where single attribute setting is about 50ms). Arguments are the same as *set\_attribute\_value()*.

**wait\_for\_frame** (*since='lastread', nframes=1, timeout=20.0, error\_on\_stopped=False*)

Wait for one or several new camera frames.

*since* specifies the reference point for waiting to acquire *nframes* frames; can be "lastread" (from the last read frame), "lastwait" (wait for the last successful *wait\_for\_frame()* call), "now" (from the start of the current call), or "start" (from the acquisition start, i.e., wait until *nframes* frames have been acquired). *timeout* can be either a number, None (infinite timeout), or a tuple (*timeout, frame\_timeout*), in which case the call times out if the total time exceeds *timeout*, or a single frame wait exceeds *frame\_timeout*. If the call times out, raise *TimeoutError*. If *error\_on\_stopped==True* and the acquisition is not running, raise *Error*; otherwise, simply return False without waiting.

```
class pylablib.devices.PhotonFocus.PhotonFocus.PhotonFocusSiSoCamera (siso_board,
                                                                    siso_applet,
                                                                    siso_port=0,
                                                                    pf-
                                                                    cam_port=0)
```

Bases: *pylablib.devices.PhotonFocus.PhotonFocus.IPhotonFocusCamera, pylablib.devices.SiliconSoftware.fgrab.SiliconSoftwareFrameGrabber*

IMAQ+PFCam interface to a PhotonFocus camera.

### Parameters

- **siso\_board** – Silicon Software board index, starting from 0; available boards can be learned by `fgrab.list_boards()`
- **siso\_applet** – Silicon Software applet name, which can be learned by `fgrab.list_applets()`; usually, a simple applet like "DualLineGray16" or "MediumLineGray16" are most appropriate; can be either an applet name, or a direct path to the applet DLL
- **siso\_port** – Silicon Software port number, if several ports are supported by the camera and the applet
- **pfcam\_port** – port number for pfcam interface (can be learned by `list_cameras()`; port number is the first element of the camera data tuple) can also be a tuple (manufacturer, port), e.g., ("National Instruments", "port0").

### Error

alias of `pylablib.core.devio.base.DeviceError`

### GrabberClass

alias of `pylablib.devices.SiliconSoftware.fgrab.SiliconSoftwareFrameGrabber`

### open()

Open connection to the camera

### FrameTransferError

alias of `pylablib.devices.interface.camera.DefaultFrameTransferError`

### class NoParameterCaller(device, kind)

Bases: `object`

Class to simplify calling functions without a parameter

**TimeoutError** = <Mock spec='str' id='140096289529296'>

### acquisition\_in\_progress()

Check if acquisition is in progress

### apply\_settings(settings)

Apply the settings.

*settings* is the dict {name: value} of the device available settings. Non-applicable settings are ignored.

### call\_command(name, arg=0, error\_on\_missing=True)

Execute the given command with the given argument.

If the command doesn't exist and `error_on_missing==True`, raise error; otherwise, do nothing.

### clear\_acquisition()

Clear all acquisition details and free all buffers

### close()

Close connection to the camera

### enable\_CFR(enabled=True)

Enable constant frame rate mode

### enable\_status\_line(enabled=True)

Enable or disable status line

**fast\_shift\_roi** (*hstart=0, vstart=0*)

Shift ROI by only changing its origin, but keeping the shape the same.

Note that if the ROI is invalid, it won't be truncated (as is the standard behavior of `set_roi()`), which might lead to errors later.

**get\_acquisition\_parameters** ()

Get acquisition parameters.

Return dictionary {name: value}

**get\_all\_attribute\_values** (*root=", enum\_as\_str=True*)

Get values of all attributes with the given *root*

**get\_all\_attributes** (*copy=False*)

Return a dictionary of all available attributes.

If *copy==True*, copy the dictionary; otherwise, return the internal dictionary structure (should not be modified).

**get\_all\_grabber\_attribute\_values** (*root=", \*\*kwargs*)

Get values of all frame grabber attributes with the given *root*.

Additional arguments are passed to `get_value` methods of individual attributes.

**get\_all\_grabber\_attributes** (*copy=False*)

Return a dictionary of all available frame grabber `grabber_attributes`.

If *copy==True*, copy the dictionary; otherwise, return the internal dictionary structure (should not be modified).

**get\_attribute** (*name, error\_on\_missing=True*)

Get the camera attribute with the given name

**get\_attribute\_value** (*name, enum\_as\_str=True, error\_on\_missing=True, default=None*)

Get value of an attribute with the given name.

If the value doesn't exist or can not be read and *error\_on\_missing==True*, raise error; otherwise, return *default*. If *default* is not None, assume that *error\_on\_missing==False*. If *enum\_as\_str==True*, try to represent enums as their string values; If *name* points at a dictionary branch, return a dictionary with all values in this branch.

**get\_available\_camlink\_pixel\_formats** ()

Get all available CamLink pixel formats and the output pixel formats as a tuple of 2 lists

**get\_baudrate** ()

Get the current baud rate

**get\_black\_level\_offset** ()

Get the black level offset

**get\_camlink\_pixel\_format** ()

Get CamLink pixel format and the output pixel format as a tuple

**get\_data\_dimensions** ()

Get readout data dimensions (in pixels) as a tuple (*width, height*); take indexing mode into account

**get\_detector\_size** ()

Get camera detector size (in pixels) as a tuple (*width, height*)

**get\_device\_info** ()

Get camera model data.

Return tuple (*model, serial\_number, grabber\_info*).

**get\_device\_variable** (*key*)

Get the value of a settings, status, or full info parameter

**get\_exposure** ()

Get current exposure

**get\_frame\_format** ()

Get format for the returned images.

Can be "list" (list of 2D arrays), "array" (a single 3D array), or "chunks" (list of 3D "chunk" arrays; supported for some cameras and provides the best performance).

**get\_frame\_info\_fields** ()

Get the names of frame info fields.

Applicable when frame info format (set by `set_frame_info_format()`) is "list" or "array".

**get\_frame\_info\_format** ()

Get format of the frame info.

Can be "namedtuple" (potentially nested named tuples; convenient to get particular values), "list" (flat list of values, with field names are given by `get_frame_info_fields()`; convenient for building a table), "array" (same as "list", but with a numpy array, which is easier to use for "chunks" frame format), or "dict" (flat dictionary with the same fields as the "list" format; more resilient to future format changes)

**get\_frame\_info\_period** ()

Get period of frame info acquisition.

Frame info might be skipped (set to `None`) except for frames which indices are divisible by *period*. Useful for certain cameras where acquiring frame info takes a lot of time and can reduce performance at higher frame rates. Note that this parameter can still be ignored (i.e., always set to 1) if the performance is not an issue for a given camera class.

**get\_frame\_period** ()

Get frame period (time between two consecutive frames in the internal trigger mode)

**get\_frame\_timings** ()

Get acquisition timing.

Return tuple (exposure, frame\_period).

**get\_frames\_status** ()

Get acquisition and buffer status.

Return tuple (acquired, unread, skipped, size), where acquired is the total number of acquired frames, unread is the number of acquired but not read frames, skipped is the number of skipped (not read and then written over) frames, and `buffer_size` is the total buffer size (in frames).

**get\_full\_info** (*include=0*)

Get dict {name: value} containing full device information (including status and settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_full\_status** (*include=0*)

Get dict {name: value} containing the device status (including settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the

lowest priority is -10, setting `include=-10` queries all available variables, which is equivalent to `include="all"`.

**get\_genicam\_info\_xml()**

Get description in Genicam-compatible XML format

**get\_grabber\_attribute(name, error\_on\_missing=True)**

Get the camera attribute with the given name

**get\_grabber\_attribute\_value(name, error\_on\_missing=True, default=None, \*\*kwargs)**

Get value of a frame grabber attribute with the given name.

If the value doesn't exist and `error_on_missing==True`, raise error; otherwise, return *default*. If *default* is not None, automatically assume that `error_on_missing==False`. If *name* points at a dictionary branch, return a dictionary with all values in this branch. Additional arguments are passed to `get_value` methods of the individual attribute.

**get\_grabber\_detector\_size()**

Get camera detector size (in pixels) as a tuple (width, height)

**get\_grabber\_roi()**

Get current ROI.

Return tuple (hstart, hend, vstart, vend). *hstart* and *hend* specify horizontal image extent, *vstart* and *vend* specify vertical image extent (start is inclusive, stop is exclusive, starting from 0).

**get\_grabber\_roi\_limits(hbin=1, vbin=1)**

Get the minimal and maximal ROI parameters.

Return tuple (hlim, vlim), where each element is in turn a limit 5-tuple (min, max, pstep, sstep, maxbin) with, correspondingly, minimal and maximal size, position and size step, and the maximal binning (fixed to 1 if not binning is allowed). In some cameras, the step and the minimal size depend on the binning, which can be supplied.

**get\_image\_indexing()**

Get indexing for the returned images.

Can be "rct" (first index row, second index column, rows counted from the top), "rcb" (same as "rc", rows counted from the bottom), "xyt" (first index column, second index row, rows counted from the top), or "xyb" (same as "xyt", rows counted from the bottom)

**get\_new\_images\_range()**

Get the range of the new images.

Return tuple (first, last) with images range (first inclusive). If no images are available, return None. If some images were in the buffer were overwritten, exclude them from the range.

**get\_roi()**

Get current ROI.

Return tuple (hstart, hend, vstart, vend).

**get\_roi\_limits(hbin=1, vbin=1)**

**get\_settings(include=0)**

Get dict {name: value} containing all the device settings.

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting `include=-10` queries all available variables, which is equivalent to `include="all"`.

**get\_system\_info()**

Get the dictionary with all system information parameters

**get\_trigger\_interleave()**

Check if the trigger interleave is on

**grab** (*nframes=1, frame\_timeout=5.0, missing\_frame='skip', return\_info=False, buff\_size=None*)

Snap *nframes* images (with preset image read mode parameters)

*buff\_size* determines buffer size (if *None*, use the default size). Timeout is specified for a single-frame acquisition, not for the whole acquisition time. *missing\_frame* determines what to do with frames which have been lost: can be "none" (replacing them with *None*), "zero" (replacing them with zero-filled frame), or "skip" (skipping them, while still keeping total returned frames number to *n*). If *return\_info==True*, return tuple (*frames*, *infos*), where *infos* is a list of frame info tuples (camera-dependent); if some frames are missing and *missing\_frame!="skip"*, the corresponding frame info is *None*.

**is\_CFR\_enabled()**

Check if the constant frame rate mode is enabled

**is\_acquisition\_setup()**

Check if acquisition is set up.

If the camera does not support separate acquisition setup, always return *True*.

**is\_opened()**

Check if the device is connected

**is\_status\_line\_enabled()**

Check if the status line is on

**pausing\_acquisition** (*clear=None, stop=True, setup\_after=None, start\_after=True, combine\_nested=True*)

Context manager which temporarily pauses acquisition during execution of *with* block.

Useful for applying certain settings which can't be changed during the acquisition. If *clear==True*, clear acquisition in addition to stopping (by default, use the class default specified as *\_clear\_pausing\_acquisition* attribute). If *stop==True*, stop the acquisition (if *clear==True*, stop regardless). If *setup\_after==True*, setup the acquisition after pause if necessary (*None* means setup only if clearing was required). If *start\_after==True*, start the acquisition after pause if necessary (*None* means start only if stopping was required). If *combine\_nested==True*, then any nested *pausing\_acquisition* calls will stop/clear acquisition as necessary, but won't setup/start it again until this *pausing\_acquisition* call is complete.

Yields tuple (*acq\_in\_progress*, *acq\_params*), which indicates whether acquisition is currently in progress, and what are the current acquisition parameters.

**read\_multiple\_images** (*rng=None, peek=False, missing\_frame='skip', return\_info=False, return\_rng=False*)

Read multiple images specified by *rng* (by default, all un-read images).

If *rng* is specified, it is a tuple (*first*, *last*) with images range (first inclusive). If no new frames are available, return an empty list; if no acquisition is running, return *None*. If *peek==True*, return images but not mark them as read. *missing\_frame* determines what to do with frames which are out of range (missing or lost): can be "none" (replacing them with *None*), "zero" (replacing them with zero-filled frame), or "skip" (skipping them). If *return\_info==True*, return tuple (*frames*, *infos*), where *infos* is a list of frame info tuples (camera-dependent, by default, only the frame index); if some frames are missing and *missing\_frame!="skip"*, the corresponding frame info is *None*. if *return\_rng==True*, return the range covered resulting frames; if *missing\_frame=="skip"*, the range can be smaller than the supplied *rng* if some frames are skipped.

**read\_newest\_image** (*peek=False, return\_info=False*)

Read the newest un-read image.

If no un-read frames are available, return `None`. If `peek==True`, return the image but not mark it as read. If `return_info==True`, return tuple (`frame`, `info`), where `info` is an info tuples (camera-dependent, see [read\\_multiple\\_images\(\)](#)).

**read\_oldest\_image** (*peek=False, return\_info=False*)

Read the oldest un-read image.

If no un-read frames are available, return `None`. If `peek==True`, return the image but not mark it as read. If `return_info==True`, return tuple (`frame`, `info`), where `info` is an info tuples (camera-dependent, see [read\\_multiple\\_images\(\)](#)).

**set\_all\_attribute\_values** (*settings, root=", truncate=True*)

Set values of all attributes with the given `root`.

If `truncate==True`, truncate value to lie within attribute range.

**set\_all\_grabber\_attribute\_values** (*settings, root=", \*\*kwargs*)

Set values of all frame grabber attributes with the given `root`.

Additional arguments are passed to `set_value` methods of individual attributes.

**set\_attribute\_value** (*name, value, truncate=True, error\_on\_missing=True*)

Set value of an attribute with the given name.

If the value doesn't exist or can not be written and `error_on_missing==True`, raise error; otherwise, do nothing. If `name` points at a dictionary branch, set all values in this branch (in this case `value` must be a dictionary). If `truncate==True`, truncate value to lie within attribute range.

**set\_black\_level\_offset** (*offset*)

Set the black level offset

**set\_device\_variable** (*key, value*)

Set the value of a settings parameter

**set\_exposure** (*exposure*)

Set current exposure

**set\_frame\_format** (*fmt*)

Set format for the returned images.

Can be "list" (list of 2D arrays), "array" (a single 3D array), "chunks" (list of 3D "chunk" arrays; supported for some cameras and provides the best performance), or "try\_chunks" (same as "chunks", but if chunks are not supported, set to "list" instead). If format is "chunks" and chunks are not supported by the camera, it results in one frame per chunk. Note that if the format is set to "array" or "chunks", the frame info format is also automatically set to "array". If the format is set to "chunks", then the image info is also returned in chunks form (list of 2D info arrays with the same length as the corresponding frame chunks).

**set\_frame\_info\_format** (*fmt, include\_fields=None*)

Set format of the frame info.

Can be "namedtuple" (potentially nested named tuples; convenient to get particular values), "list" (flat list of values, with field names are given by [get\\_frame\\_info\\_fields\(\)](#); convenient for building a table), "array" (same as "list", but with a numpy array, which is easier to use for "chunks" frame format), or "dict" (flat dictionary with the same fields as the "list" format; more resilient to future format changes) If `include_fields` is not `None`, it specifies the fields included for non-"tuple" formats; note that order of `include_fields` is ignored, and the resulting fields are always ordered same as in the original.



**set\_frame\_info\_period** (*period=1*)

Set period of frame info acquisition.

Frame info might be skipped (set to `None`) except for frames which indices are divisible by *period*. Useful for certain cameras where acquiring frame info takes a lot of time and can reduce performance at higher frame rates. Note that this parameter can still be ignored (i.e., always set to 1) if the performance is not an issue for a given camera class.

**set\_frame\_merge** (*frame\_merge=1*)

**set\_frame\_period** (*frame\_period*)

Set frame period (time between two consecutive frames in the internal trigger mode)

**set\_grabber\_attribute\_value** (*name, value, error\_on\_missing=True, \*\*kwargs*)

Set value of a frame grabber attribute with the given name.

If the value doesn't exist and `error_on_missing==True`, raise error; otherwise, do nothing. If *name* points at a dictionary branch, set all values in this branch (in this case *value* must be a dictionary). Additional arguments are passed to `set_value` methods of the individual attribute.

**set\_grabber\_roi** (*hstart=0, hend=None, vstart=0, vend=None*)

Setup camera ROI.

*hstart* and *hend* specify horizontal image extent, *vstart* and *vend* specify vertical image extent (start is inclusive, stop is exclusive, starting from 0). By default, all non-supplied parameters take extreme values (0 for start, maximal for end).

**set\_image\_indexing** (*indexing*)

Set up indexing for the returned images.

Can be "rct" (first index row, second index column, rows counted from the top), "rcb" (same as "rc", rows counted from the bottom), "xyt" (first index column, second index row, rows counted from the top), or "xyb" (same as "xyt", rows counted from the bottom)

**set\_roi** (*hstart=0, hend=None, vstart=0, vend=None*)

Setup camera ROI.

By default, all non-supplied parameters take extreme values.

**set\_trigger\_interleave** (*enabled*)

Set the trigger interleave option on or off

**setup\_acquisition** (*mode='sequence', nframes=100*)

Setup acquisition mode.

*mode* can be either "snap" (single frame or a fixed number of frames) or "sequence" (continuous acquisition). (note that `IMAQCamera.acquisition_in_progress()` would still return `True` in this case, even though new frames are no longer acquired). *nframes* sets up number of frame buffers.

**setup\_camlink\_pixel\_format** (*bits\_per\_pixel=8, taps=1, output\_fmt=None, fmt=None*)

Set up CameraLink pixel format.

If *fmt* is `None`, use supplied *bits\_per\_pixel* (8, 10, 12, 14, or 16) and *taps* (1 or 2) to figure out the format; otherwise, *fmt* should be a numerical (e.g., 210) or string (e.g., "FG\_CL\_MEDIUM\_10\_BIT") format. *output\_fmt* specifies the result frame format; if `None`, use grayscale with the given *bits\_per\_pixel* if *fmt* is `None`, or 16 bit grayscale otherwise.

**setup\_max\_baudrate** ()

Setup the maximal available baudrate

**snap** (*timeout=5.0, return\_info=False*)

Snap a single frame

**start\_acquisition** (\*args, \*\*kwargs)

Start acquisition.

Can take the same keyword parameters as :meth:“*setup\_acquisition*”. If the acquisition is not set up yet, set it up using the supplied parameters (use default of *setup\_acquisition()*, if the parameter is None). Otherwise, if any supplied parameters are different from the current ones, change them and reset the acquisition.

**stop\_acquisition** ()

Stop acquisition

**update\_attribute\_value** (name, value, error\_on\_missing=True, truncate=True)

Set value of the attribute with a given name, but only if it's different from the current value.

Can take less time on some version of PFRremote (where single attribute setting is about 50ms). Arguments are the same as *set\_attribute\_value()*.

**wait\_for\_frame** (since='lastread', nframes=1, timeout=20.0, error\_on\_stopped=False)

Wait for one or several new camera frames.

*since* specifies the reference point for waiting to acquire *nframes* frames; can be “lastread” (from the last read frame), “lastwait” (wait for the last successful *wait\_for\_frame()* call), “now” (from the start of the current call), or “start” (from the acquisition start, i.e., wait until *nframes* frames have been acquired). *timeout* can be either a number, None (infinite timeout), or a tuple (timeout, frame\_timeout), in which case the call times out if the total time exceeds timeout, or a single frame wait exceeds frame\_timeout. If the call times out, raise *TimeoutError*. If *error\_on\_stopped==True* and the acquisition is not running, raise *Error*; otherwise, simply return False without waiting.

```
class pylablib.devices.PhotonFocus.PhotonFocus.PhotonFocusBitFlowCamera (bitflow_idx=0,
                                                                           bit-
                                                                           flow_camfile=None,
                                                                           pf-
                                                                           cam_port=0)
```

Bases: *pylablib.devices.PhotonFocus.PhotonFocus.IPhotonFocusCamera*, *pylablib.devices.BitFlow.BitFlow.BitFlowFrameGrabber*

BitFlow+PFCam interface to a PhotonFocus camera.

#### Parameters

- **bitflow\_idx** – board index, starting from 0
- **bitflow\_camfile** – if not None, a path to a valid camera file used for this frame grabber and camera combination; in this case, a temporary camera file is generated based on the provided one and used to change some otherwise unavailable camera parameters such as ROI and pixel bit depth (they are otherwise fixed to whatever is specified in the default camera file)
- **pfcam\_port** – port number for pfcam interface (can be learned by *list\_cameras()*; port number is the first element of the camera data tuple) can also be a tuple (manufacturer, port), e.g., (“National Instruments”, “port0”).

#### Error

alias of *pylablib.core.devio.base.DeviceError*

#### GrabberClass

alias of *pylablib.devices.BitFlow.BitFlow.BitFlowFrameGrabber*

**open** ()

Open connection to the camera

**setup\_acquisition** (*mode='sequence', nframes=100, frame\_merge=None*)

Setup acquisition.

Any non-specified acquisition parameters are assumed to be the same as previously set (or default, if not explicitly set before). Return the new acquisition parameters.

**class BufferManager** (*cam*)

Bases: `object`

Buffer manager: stores, constantly reads and re-schedules buffers, keeps track of acquired frames and buffer overflow events

**get\_status** ()

Get counter status: tuple (acquired,)

**is\_running** ()

Check if the buffer loop is running

**reset** ()

Reset counter (on frame acquisition)

**start\_loop** ()

Start buffer scheduling loop

**stop\_loop** ()

Stop buffer scheduling loop

**FrameTransferError**

alias of `pylablib.devices.interface.camera.DefaultFrameTransferError`

**class NoParameterCaller** (*device, kind*)

Bases: `object`

Class to simplify calling functions without a parameter

**TimeoutError**

alias of `pylablib.devices.BitFlow.BitFlow.BitFlowTimeoutError`

**acquisition\_in\_progress** ()

Check if acquisition is in progress

**apply\_settings** (*settings*)

Apply the settings.

*settings* is the dict {name: value} of the device available settings. Non-applicable settings are ignored.

**call\_command** (*name, arg=0, error\_on\_missing=True*)

Execute the given command with the given argument.

If the command doesn't exist and `error_on_missing==True`, raise error; otherwise, do nothing.

**clear\_acquisition** ()

Clear all acquisition details and free all buffers

**close** ()

Close connection to the camera

**enable\_CFR** (*enabled=True*)

Enable constant frame rate mode

**enable\_status\_line** (*enabled=True*)

Enable or disable status line

**fast\_shift\_roi** (*hstart=0, vstart=0*)

Shift ROI by only changing its origin, but keeping the shape the same.

Note that if the ROI is invalid, it won't be truncated (as is the standard behavior of `set_roi()`), which might lead to errors later.

**get\_acquisition\_parameters** ()

Get acquisition parameters.

Return dictionary {name: value}

**get\_all\_attribute\_values** (*root=", enum\_as\_str=True*)

Get values of all attributes with the given *root*

**get\_all\_attributes** (*copy=False*)

Return a dictionary of all available attributes.

If *copy==True*, copy the dictionary; otherwise, return the internal dictionary structure (should not be modified).

**get\_attribute** (*name, error\_on\_missing=True*)

Get the camera attribute with the given name

**get\_attribute\_value** (*name, enum\_as\_str=True, error\_on\_missing=True, default=None*)

Get value of an attribute with the given name.

If the value doesn't exist or can not be read and *error\_on\_missing==True*, raise error; otherwise, return *default*. If *default* is not None, assume that *error\_on\_missing==False*. If *enum\_as\_str==True*, try to represent enums as their string values; If *name* points at a dictionary branch, return a dictionary with all values in this branch.

**get\_baudrate** ()

Get the current baud rate

**get\_black\_level\_offset** ()

Get the black level offset

**get\_data\_dimensions** ()

Get readout data dimensions (in pixels) as a tuple (*width, height*); take indexing mode into account

**get\_detector\_size** ()

Get camera detector size (in pixels) as a tuple (*width, height*)

**get\_device\_info** ()

Get camera model data.

Return tuple (*model, serial\_number, grabber\_info*).

**get\_device\_variable** (*key*)

Get the value of a settings, status, or full info parameter

**get\_exposure** ()

Get current exposure

**get\_frame\_format** ()

Get format for the returned images.

Can be "list" (list of 2D arrays), "array" (a single 3D array), or "chunks" (list of 3D "chunk" arrays; supported for some cameras and provides the best performance).

**get\_frame\_info\_fields** ()

Get the names of frame info fields.

Applicable when frame info format (set by `set_frame_info_format()`) is "list" or "array".

**get\_frame\_info\_format()**

Get format of the frame info.

Can be "namedtuple" (potentially nested named tuples; convenient to get particular values), "list" (flat list of values, with field names are given by `get_frame_info_fields()`; convenient for building a table), "array" (same as "list", but with a numpy array, which is easier to use for "chunks" frame format), or "dict" (flat dictionary with the same fields as the "list" format; more resilient to future format changes)

**get\_frame\_info\_period()**

Get period of frame info acquisition.

Frame info might be skipped (set to `None`) except for frames which indices are divisible by *period*. Useful for certain cameras where acquiring frame info takes a lot of time and can reduce performance at higher frame rates. Note that this parameter can still be ignored (i.e., always set to 1) if the performance is not an issue for a given camera class.

**get\_frame\_period()**

Get frame period (time between two consecutive frames in the internal trigger mode)

**get\_frame\_timings()**

Get acquisition timing.

Return tuple (*exposure*, *frame\_period*).

**get\_frames\_status()**

Get acquisition and buffer status.

Return tuple (*acquired*, *unread*, *skipped*, *size*), where *acquired* is the total number of acquired frames, *unread* is the number of acquired but not read frames, *skipped* is the number of skipped (not read and then written over) frames, and *buffer\_size* is the total buffer size (in frames).

**get\_full\_info(include=0)**

Get dict {*name*: *value*} containing full device information (including status and settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include*=-10 queries all available variables, which is equivalent to *include*="all".

**get\_full\_status(include=0)**

Get dict {*name*: *value*} containing the device status (including settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include*=-10 queries all available variables, which is equivalent to *include*="all".

**get\_grabber\_detector\_size()**

Get camera detector size (in pixels) as a tuple (*width*, *height*)

**get\_grabber\_roi()**

Get current ROI.

Return tuple (*hstart*, *hend*, *vstart*, *vend*). *hstart* and *hend* specify horizontal image extent, *vstart* and *vend* specify vertical image extent (start is inclusive, stop is exclusive, starting from 0).

**get\_grabber\_roi\_limits(hbin=1, vbin=1)**

Get the minimal and maximal ROI parameters.

Return tuple (*hlim*, *vlim*), where each element is in turn a limit 5-tuple (*min*, *max*, *pstep*, *sstep*, *maxbin*) with, correspondingly, minimal and maximal size, position and size step, and the

maximal binning (fixed to 1 if not binning is allowed). In some cameras, the step and the minimal size depend on the binning, which can be supplied.

**get\_image\_indexing()**

Get indexing for the returned images.

Can be "rct" (first index row, second index column, rows counted from the top), "rcb" (same as "rc", rows counted from the bottom), "xyt" (first index column, second index row, rows counted from the top), or "xyb" (same as "xyt", rows counted from the bottom)

**get\_new\_images\_range()**

Get the range of the new images.

Return tuple (*first*, *last*) with images range (first inclusive). If no images are available, return None. If some images were in the buffer were overwritten, exclude them from the range.

**get\_roi()**

Get current ROI.

Return tuple (*hstart*, *hend*, *vstart*, *vend*).

**get\_roi\_limits(hbin=1, vbin=1)**

**get\_settings(include=0)**

Get dict {*name*: *value*} containing all the device settings.

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include*=-10 queries all available variables, which is equivalent to *include*="all".

**get\_trigger\_interleave()**

Check if the trigger interleave is on

**grab(nframes=1, frame\_timeout=5.0, missing\_frame='skip', return\_info=False, buff\_size=None)**

Snap *nframes* images (with preset image read mode parameters)

*buff\_size* determines buffer size (if None, use the default size). Timeout is specified for a single-frame acquisition, not for the whole acquisition time. *missing\_frame* determines what to do with frames which have been lost: can be "none" (replacing them with None), "zero" (replacing them with zero-filled frame), or "skip" (skipping them, while still keeping total returned frames number to *n*). If *return\_info*==True, return tuple (*frames*, *infos*), where *infos* is a list of frame info tuples (camera-dependent); if some frames are missing and *missing\_frame*!="skip", the corresponding frame info is None.

**is\_CFR\_enabled()**

Check if the constant frame rate mode is enabled

**is\_acquisition\_setup()**

Check if acquisition is set up.

If the camera does not support separate acquisition setup, always return True.

**is\_opened()**

Check if the device is connected

**is\_status\_line\_enabled()**

Check if the status line is on

**pausing\_acquisition(clear=None, stop=True, setup\_after=None, start\_after=True, combine\_nested=True)**

Context manager which temporarily pauses acquisition during execution of `with` block.

Useful for applying certain settings which can't be changed during the acquisition. If `clear==True`, clear acquisition in addition to stopping (by default, use the class default specified as `_clear_pausing_acquisition` attribute). If `stop==True`, stop the acquisition (if `clear==True`, stop regardless). If `setup_after==True`, setup the acquisition after pause if necessary (None means setup only if clearing was required). If `start_after==True`, start the acquisition after pause if necessary (None means start only if stopping was required). If `combine_nested==True`, then any nested `pausing_acquisition` calls will stop/clear acquisition as necessary, but won't setup/start it again until this `pausing_acquisition` call is complete.

Yields tuple (`acq_in_progress`, `acq_params`), which indicates whether acquisition is currently in progress, and what are the current acquisition parameters.

**read\_multiple\_images** (*rng=None*, *peek=False*, *missing\_frame='skip'*, *return\_info=False*, *return\_rng=False*)

Read multiple images specified by *rng* (by default, all un-read images).

If *rng* is specified, it is a tuple (*first*, *last*) with images range (*first* inclusive). If no new frames are available, return an empty list; if no acquisition is running, return None. If *peek==True*, return images but not mark them as read. *missing\_frame* determines what to do with frames which are out of range (missing or lost): can be "none" (replacing them with None), "zero" (replacing them with zero-filled frame), or "skip" (skipping them). If *return\_info==True*, return tuple (*frames*, *infos*), where *infos* is a list of frame info tuples (camera-dependent, by default, only the frame index); if some frames are missing and *missing\_frame!="skip"*, the corresponding frame info is None. If *return\_rng==True*, return the range covered resulting frames; if *missing\_frame=="skip"*, the range can be smaller than the supplied *rng* if some frames are skipped.

**read\_newest\_image** (*peek=False*, *return\_info=False*)

Read the newest un-read image.

If no un-read frames are available, return None. If *peek==True*, return the image but not mark it as read. If *return\_info==True*, return tuple (*frame*, *info*), where *info* is an info tuples (camera-dependent, see `read_multiple_images()`).

**read\_oldest\_image** (*peek=False*, *return\_info=False*)

Read the oldest un-read image.

If no un-read frames are available, return None. If *peek==True*, return the image but not mark it as read. If *return\_info==True*, return tuple (*frame*, *info*), where *info* is an info tuples (camera-dependent, see `read_multiple_images()`).

**set\_all\_attribute\_values** (*settings*, *root=""*, *truncate=True*)

Set values of all attributes with the given *root*.

If *truncate==True*, truncate value to lie within attribute range.

**set\_attribute\_value** (*name*, *value*, *truncate=True*, *error\_on\_missing=True*)

Set value of an attribute with the given name.

If the value doesn't exist or can not be written and *error\_on\_missing==True*, raise error; otherwise, do nothing. If *name* points at a dictionary branch, set all values in this branch (in this case *value* must be a dictionary). If *truncate==True*, truncate value to lie within attribute range.

**set\_black\_level\_offset** (*offset*)

Set the black level offset

**set\_device\_variable** (*key*, *value*)

Set the value of a settings parameter

**set\_exposure** (*exposure*)

Set current exposure

**set\_frame\_format** (*fmt*)

Set format for the returned images.

Can be "list" (list of 2D arrays), "array" (a single 3D array), "chunks" (list of 3D "chunk" arrays; supported for some cameras and provides the best performance), or "try\_chunks" (same as "chunks", but if chunks are not supported, set to "list" instead). If format is "chunks" and chunks are not supported by the camera, it results in one frame per chunk. Note that if the format is set to "array" or "chunks", the frame info format is also automatically set to "array". If the format is set to "chunks", then the image info is also returned in chunks form (list of 2D info arrays with the same length as the corresponding frame chunks).

**set\_frame\_info\_format** (*fmt, include\_fields=None*)

Set format of the frame info.

Can be "namedtuple" (potentially nested named tuples; convenient to get particular values), "list" (flat list of values, with field names are given by `get_frame_info_fields()`; convenient for building a table), "array" (same as "list", but with a numpy array, which is easier to use for "chunks" frame format), or "dict" (flat dictionary with the same fields as the "list" format; more resilient to future format changes) If *include\_fields* is not None, it specifies the fields included for non-"tuple" formats; note that order or *include\_fields* is ignored, and the resulting fields are always ordered same as in the original.

**set\_frame\_info\_period** (*period=1*)

Set period of frame info acquisition.

Frame info might be skipped (set to None) except for frames which indices are divisible by *period*. Useful for certain cameras where acquiring frame info takes a lot of time and can reduce performance at higher frame rates. Note that this parameter can still be ignored (i.e., always set to 1) if the performance is not an issue for a given camera class.

**set\_frame\_period** (*frame\_period*)

Set frame period (time between two consecutive frames in the internal trigger mode)

**set\_grabber\_roi** (*hstart=0, hend=None, vstart=0, vend=None*)

Setup camera ROI.

*hstart* and *hend* specify horizontal image extent, *vstart* and *vend* specify vertical image extent (start is inclusive, stop is exclusive, starting from 0). By default, all non-supplied parameters take extreme values (0 for start, maximal for end).

**set\_image\_indexing** (*indexing*)

Set up indexing for the returned images.

Can be "rct" (first index row, second index column, rows counted from the top), "rcb" (same as "rc", rows counted from the bottom), "xyt" (first index column, second index row, rows counted from the top), or "xyb" (same as "xyt", rows counted from the bottom)

**set\_roi** (*hstart=0, hend=None, vstart=0, vend=None*)

Setup camera ROI.

By default, all non-supplied parameters take extreme values.

**set\_trigger\_interleave** (*enabled*)

Set the trigger interleave option on or off

**setup\_max\_baudrate** ()

Setup the maximal available baudrate

**snap** (*timeout=5.0, return\_info=False*)

Snap a single frame



**start\_acquisition** (\*args, \*\*kwargs)

Start acquisition.

Can take the same keyword parameters as :meth:“*setup\_acquisition*”. If the acquisition is not set up yet, set it up using the supplied parameters (use default of *setup\_acquisition()*, if the parameter is None). Otherwise, if any supplied parameters are different from the current ones, change them and reset the acquisition.

**stop\_acquisition** ()

Stop acquisition

**update\_attribute\_value** (name, value, error\_on\_missing=True, truncate=True)

Set value of the attribute with a given name, but only if it's different from the current value.

Can take less time on some version of PFRemote (where single attribute setting is about 50ms). Arguments are the same as *set\_attribute\_value()*.

**wait\_for\_frame** (since='lastread', nframes=1, timeout=20.0, error\_on\_stopped=False)

Wait for one or several new camera frames.

*since* specifies the reference point for waiting to acquire *nframes* frames; can be “lastread” (from the last read frame), “lastwait” (wait for the last successful *wait\_for\_frame()* call), “now” (from the start of the current call), or “start” (from the acquisition start, i.e., wait until *nframes* frames have been acquired). *timeout* can be either a number, None (infinite timeout), or a tuple (timeout, frame\_timeout), in which case the call times out if the total time exceeds timeout, or a single frame wait exceeds frame\_timeout. If the call times out, raise *TimeoutError*. If *error\_on\_stopped*==True and the acquisition is not running, raise *Error*; otherwise, simply return False without waiting.

`pylablib.devices.PhotonFocus.PhotonFocus.check_grabber_association(cam)`

Check if PhotonFocus camera has correct association between the frame grabber and the PFRemote interface.

*cam* should be an opened instance of *PhotonFocusIMAQCamera* or *PhotonFocusSiSoCamera*. Note that this function changes camera parameters such as exposure, frame period, ROI, trigger source, and status line.

`pylablib.devices.PhotonFocus.PhotonFocus.get_status_lines(frames,  
check_transposed=True,  
drop_magic=True)`

Extract status lines (up to first 6 entries) from the given frames.

*frames* can be 2D array (one frame), 3D array (stack of frames, first index is frame number), or list of 1D or 2D arrays. Automatically check if the status line is present; return None if it's not. If *check\_transposed*==True, check for the case where the image is transposed (i.e., line becomes a column). If *drop\_magic*==True, remove the first status line entry, which is simply a special number marking the status line presence. Return a 1D or 2D numpy array, where the first axis (if present) is the frame number, and the last is the status line entry. The entries after the magic are the frame index, timestamp (in us), missed trigger counters (up to 255), average frame value, and the integration time (in pixel clock cycles, which depend on the camera).

`pylablib.devices.PhotonFocus.PhotonFocus.get_status_line_position(frame,  
check_transposed=True)`

Check whether status line is present in the frame, and return its location.

Return tuple (row, transposed), where *row* is the status line row (can be -1 or -2) and *transposed* is True if the line is present in the transposed image. If no status line is found, return None. If *check\_transposed*==True, check for the case where the image is transposed (i.e., line becomes a column).

```
pylablib.devices.PhotonFocus.PhotonFocus.remove_status_line(frame,
                                                             sl_pos='calculate',
                                                             policy='duplicate',
                                                             copy=True)
```

Remove status line from the frame.

#### Parameters

- **frame** – a frame to process (2D or 3D numpy array; if 3D, the first axis is the frame number)
- **sl\_pos** – status line position (returned by `get_status_line_position()`); if equal to "calculate", calculate here; for a 3D array, assumed to be the same for all frames
- **policy** – determines way to deal with the status line; can be "keep" (keep as is), "cut" (cut off the status line row), "zero" (set it to zero), "median" (set it to the image median), or "duplicate" (set it equal to the previous row; default)
- **copy** – if True, make copy of the original frames; otherwise, attempt to remove the line in-place

```
pylablib.devices.PhotonFocus.PhotonFocus.find_skipped_frames(lines, step=1)
```

Check if there are skipped frames based on status line reading.

*step* specifies expected index step between neighboring frames.

Return list [ (*idx*, *skipped*) ], where *idx* is the index after which *skipped* frames were skipped.

```
class pylablib.devices.PhotonFocus.PhotonFocus.StatusLineChecker
```

Bases: `pylablib.devices.interface.camera.StatusLineChecker`

```
check_indices(indices, step=1)
```

Check if indices are consistent with the given step

```
get_framestamp(frames)
```

Get framestamps from status lines in the given frames

## Module contents

### pylablib.devices.PhysikInstrumente package

#### Submodules

#### pylablib.devices.PhysikInstrumente.base module

```
exception pylablib.devices.PhysikInstrumente.base.PhysikInstrumenteError
```

Bases: `pylablib.core.devio.base.DeviceError`

Generic Physik Instrumente error

```
args
```

```
with_traceback()
```

Exception.with\_traceback(tb) – set self.\_\_traceback\_\_ to tb and return self.

```
exception pylablib.devices.PhysikInstrumente.base.PhysikInstrumenteBackendError(exc)
```

Bases: `pylablib.devices.PhysikInstrumente.base.PhysikInstrumenteError`,  
`pylablib.core.devio.comm_backend.DeviceBackendError`

Generic Physik Instrumente backend communication error

**args**

**with\_traceback()**

Exception.with\_traceback(tb) – set self.\_\_traceback\_\_ to tb and return self.

```
class pylablib.devices.PhysikInstrumente.base.GenericPIController(conn,
                                                                auto_online=True)
Bases: pylablib.core.devio.comm_backend.ICommBackendWrapper, pylablib.
        devices.interface.stage.IMultiaxisStage
```

Generic Physik Instrumente controller.

#### Parameters

- **conn** – connection parameters (usually port or a tuple containing port and baudrate)
- **auto\_online** – if True, switch to the online mode upon connection; in this online mode controller parameters are controlled remotely instead of the front panel (including external voltages), while in the offline mode most of the parameters are still controlled manually, and the remote connection is mostly used for readout

#### Error

alias of *PhysikInstrumenteError*

**open()**

Open the backend

**query(comm, multiline=False, reply=True)**

Query a single command to the controller.

If *multiline*==True, expect a multi-line reply and return a list with separate reply lines; otherwise, expect a single-line reply and raise an error if multi-line reply is received.

If *reply*==False, expect no reply at all (used for, e.g., set commands).

**query\_axis(comm, axis=None, subidx=None, kind='str')**

Query the given command for the given axis.

*axis* can be a single axis name (e.g., "A"), a list of axes, or None, which queries all axes. If *axis* is a single axis, simply return the corresponding value; otherwise, return a dict {*axis*: value}. *kind* can specify value kind: "str" (return as is), "float", "int", or "bool".

**set\_axis(comm, value, axis=None, subidx=None, reply=False)**

Query the given value for the given axis.

*value* can be a single value (set the same for all specified axes), a list of values (one per axis), or a dict {*axis*: value}. *axis* can be a single axis name (e.g., "A"), a list of axes, or None, which queries all axes. If *reply*==False, expect no reply.

**get\_id()**

Get the device ID string

**get\_help()**

Get the help for all commands; might take a long time on low-speed serial connections

**is\_online\_enabled()**

Check if online mode is enabled

**enable\_online(enable=True)**

Enable or disable online mode

**get\_axis\_parameter(pid, axis=None, kind='str')**

Get value of the given parameter id for the given axis (all axes by default)

**set\_axis\_parameter** (*pid, value, axis=None, kind='str'*)

Get value of the given parameter id for the given axis (all axes by default)

**class NoParameterCaller** (*device, kind*)

Bases: `object`

Class to simplify calling functions without a parameter

**apply\_settings** (*settings*)

Apply the settings.

*settings* is the dict {name: value} of the device available settings. Non-applicable settings are ignored.

**close** ()

Close the backend

**get\_all\_axes** ()

Get the list of all available axes (taking mapping into account)

**get\_device\_variable** (*key*)

Get the value of a settings, status, or full info parameter

**get\_full\_info** (*include=0*)

Get dict {name: value} containing full device information (including status and settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_full\_status** (*include=0*)

Get dict {name: value} containing the device status (including settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_settings** (*include=0*)

Get dict {name: value} containing all the device settings.

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**is\_opened** ()

Check if the device is connected

**lock** (*timeout=None*)

Lock the access to the device from other threads/processes (isn't necessarily implemented)

**locking** (*timeout=None*)

Context manager for lock & unlock

**remap\_axes** (*mapping, accept\_original=True*)

Rename axes to the new labels.

*mapping* is the new axes mapping, which can be a list of new axes name (corresponding to the old axes in order returned by `get_all_axes()`), or a dictionary {alias: original} of the new axes aliases.

**set\_device\_variable** (*key, value*)  
Set the value of a settings parameter

**unlock** ()  
Unlock the access to the device from other threads/processes (isn't necessarily implemented)

**class** pylablib.devices.PhysikInstrumente.base.**PIE516** (*conn, auto\_online=True*)

Bases: *pylablib.devices.PhysikInstrumente.base.GenericPIController*

Physik Instrumente E-516 controller.

#### Parameters

- **conn** – serial connection parameters (usually port or a tuple containing port and baudrate)
- **auto\_online** – if `True`, switch to the online mode upon connection; in this online mode controller parameters such as voltages or positions are controlled remotely instead of the front panel (including external voltages), while in the offline mode most of the parameters are still controlled manually, and the remote connection is mostly used for readout

**is\_servo\_enabled** (*axis=None*)  
Check if the servo is enabled on the given axis (all axes by default)

**enable\_servo** (*enable=True, axis=None*)  
Enable or disable servo on the given axis (all axes by default)

**is\_drift\_compensation\_enabled** (*axis=None*)  
Check if the drift compensation is enabled on the given axis (all axes by default)

**enable\_drift\_compensation** (*enable=True, axis=None*)  
Enable or disable drift compensation on the given axis (all axes by default)

**is\_velocity\_control\_enabled** (*axis=None*)  
Check if the velocity control is enabled on the given axis (all axes by default)

**enable\_velocity\_control** (*enable=True, axis=None*)  
Enable or disable velocity control on the given axis (all axes by default)

**get\_voltage\_setpoint** (*axis=None*)  
Get the current voltage setpoint on the given axis (all axes by default)

**get\_voltage** (*axis=None*)  
Get the actual voltage value on the given axis (all axes by default)

**set\_voltage** (*voltage, axis=None*)  
Get the target voltage on the given axis (all axes by default)

**get\_voltage\_lower\_limit** (*axis=None*)  
Get the lower output voltage limit on the given axis (all axes by default)

**set\_voltage\_lower\_limit** (*voltage, axis=None*)  
Get the lower output voltage limit on the given axis (all axes by default)

**get\_voltage\_upper\_limit** (*axis=None*)  
Get the upper output voltage limit on the given axis (all axes by default)

**set\_voltage\_upper\_limit** (*voltage, axis=None*)  
Get the upper output voltage limit on the given axis (all axes by default)

**get\_velocity** (*axis=None*)  
Get velocity on the given axis (all axes by default)

**set\_velocity** (*velocity*, *axis=None*)

Set velocity on the given axis (all axes by default)

**get\_position** (*axis=None*)

Get the current position on the given axis

**get\_target\_position** (*axis=None*)

Get the target motion position on the given axis

**move\_to** (*position*, *axis=None*)

Move the given axis to the given position

**move\_by** (*distance*, *axis=None*)

Move the given axis by the given distance

**stop** (*axis=None*)

Stop motion on the given axis (all axes by default)

**get\_position\_lower\_limit** (*axis=None*)

Get the lower position limit on the given axis (all axes by default)

**set\_position\_lower\_limit** (*position*, *axis=None*)

Get the lower position limit on the given axis (all axes by default)

**get\_position\_upper\_limit** (*axis=None*)

Get the upper position limit on the given axis (all axes by default)

**set\_position\_upper\_limit** (*position*, *axis=None*)

Get the upper position limit on the given axis (all axes by default)

**Error**

alias of *PhysikInstrumenteError*

**class NoParameterCaller** (*device*, *kind*)

Bases: *object*

Class to simplify calling functions without a parameter

**apply\_settings** (*settings*)

Apply the settings.

*settings* is the dict {*name*: *value*} of the device available settings. Non-applicable settings are ignored.

**close** ()

Close the backend

**enable\_online** (*enable=True*)

Enable or disable online mode

**get\_all\_axes** ()

Get the list of all available axes (taking mapping into account)

**get\_axis\_parameter** (*pid*, *axis=None*, *kind='str'*)

Get value of the given parameter id for the given axis (all axes by default)

**get\_device\_variable** (*key*)

Get the value of a settings, status, or full info parameter

**get\_full\_info** (*include=0*)

Get dict {*name*: *value*} containing full device information (including status and settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the

lowest priority is -10, setting `include=-10` queries all available variables, which is equivalent to `include="all"`.

**get\_full\_status** (*include=0*)

Get dict {name: value} containing the device status (including settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting `include=-10` queries all available variables, which is equivalent to `include="all"`.

**get\_help** ()

Get the help for all commands; might take a long time on low-speed serial connections

**get\_id** ()

Get the device ID string

**get\_settings** (*include=0*)

Get dict {name: value} containing all the device settings.

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting `include=-10` queries all available variables, which is equivalent to `include="all"`.

**is\_online\_enabled** ()

Check if online mode is enabled

**is\_opened** ()

Check if the device is connected

**lock** (*timeout=None*)

Lock the access to the device from other threads/processes (isn't necessarily implemented)

**locking** (*timeout=None*)

Context manager for lock & unlock

**open** ()

Open the backend

**query** (*comm, multiline=False, reply=True*)

Query a single command to the controller.

If `multiline==True`, expect a multi-line reply and return a list with separate reply lines; otherwise, expect a single-line reply and raise an error if multi-line reply is received.

If `reply==False`, expect no reply at all (used for, e.g., set commands).

**query\_axis** (*comm, axis=None, subidx=None, kind='str'*)

Query the given command for the given axis.

*axis* can be a single axis name (e.g., "A"), a list of axes, or `None`, which queries all axes. If *axis* is a single axis, simply return the corresponding value; otherwise, return a dict {axis: value}. *kind* can specify value kind: "str" (return as is), "float", "int", or "bool".

**remap\_axes** (*mapping, accept\_original=True*)

Rename axes to the new labels.

*mapping* is the new axes mapping, which can be a list of new axes name (corresponding to the old axes in order returned by `get_all_axes()`), or a dictionary {alias: original} of the new axes aliases.

**set\_axis** (*comm, value, axis=None, subidx=None, reply=False*)

Query the given value for the given axis.

*value* can be a single value (set the same for all specified axes), a list of values (one per axis), or a dict {axis: value}. *axis* can be a single axis name (e.g., "A"), a list of axes, or None, which queries all axes. If *reply==False*, expect no reply.

**set\_axis\_parameter** (*pid, value, axis=None, kind='str'*)

Get value of the given parameter id for the given axis (all axes by default)

**set\_device\_variable** (*key, value*)

Set the value of a settings parameter

**unlock** ()

Unlock the access to the device from other threads/processes (isn't necessarily implemented)

## Module contents

### pylablib.devices.PrincetonInstruments package

#### Submodules

#### pylablib.devices.PrincetonInstruments.picam module

**class** pylablib.devices.PrincetonInstruments.picam.**LibraryController** (*lib*)

Bases: *pylablib.devices.utils.load\_lib.LibraryController*

**close** (*opid*)

Mark device closing.

Return tuple (*close\_result, uninit\_result*) with the results of the closing and the shutdown. If library does not need to be shut down yet, set *uninit\_result=None*

**open** ()

Mark device opening.

Return tuple (*init\_result, open\_result, opid*) with the results of the initialization and the opening, and the opening ID which should afterwards be used for closing. If library is already initialized, set *init\_result=None*

**preinit** ()

Pre-initialize the library, if it hasn't been done already

**shutdown** ()

Close all opened connections and shutdown the library

**temp\_open** ()

Context for temporarily opening a new device connection

**class** pylablib.devices.PrincetonInstruments.picam.**TCameraInfo** (*name, serial\_number, model, interface*)

Bases: *tuple*

**count** ()

Return number of occurrences of value.

**index** ()

Return first index of value.



Raises ValueError if the value is not present.

**interface**

**model**

**name**

**serial\_number**

`pylablib.devices.PrincetonInstruments.picam.list_cameras()`

List all cameras available through Picam interface

`pylablib.devices.PrincetonInstruments.picam.get_cameras_number()`

Get number of connected Picam cameras

**class** `pylablib.devices.PrincetonInstruments.picam.TROIConstraints` (*flags, nrois, xrng, wrng, xbins, yrng, hrng, ybins*)

Bases: `tuple`

**count()**

Return number of occurrences of value.

**flags**

**hrng**

**index()**

Return first index of value.

Raises ValueError if the value is not present.

**nrois**

**wrng**

**xbins**

**xrng**

**ybins**

**yrng**

**class** `pylablib.devices.PrincetonInstruments.picam.PicamAttribute` (*handle, pid*)

Bases: `object`

Object representing an Picam camera parameter.

Allows to query and set values and get additional information. Usually created automatically by an `PicamCamera` instance, but could be created manually.

#### Parameters

- **handle** – camera handle
- **pid** – parameter id of the attribute

**name**

attribute name

**kind**

attribute kind; can be "Integer", "Large Integer", "Floating Point", "Enumeration", "Boolean", or "Rois"

**exists**

whether attribute is available on the current hardware

Type `bool`

**relevant**

whether attribute value is applicable to the hardware

Type `bool`

**read\_directly**

whether value can be read directly from the device; if `True`, then `get_value()` will automatically use the appropriate method

Type `bool`

**value\_access**

value access kind, which shows whether value can be written

Type `str`

**writable**

whether value is read-only

Type `bool`

**default**

default parameter value (only for writable parameters)

**can\_set\_online**

whether value can be changed during acquisition

Type `bool`

**cons\_type**

constraint type, e.g., "Collection", "Range", or "None"

Type `str`

**cons\_permanent**

whether the constraint is permanent, or dependent on other parameters; if `False`, then use `update_limits()` to update the constraints

Type `bool`

**cons\_error**

whether setting the out-of-range parameter causes error or just warning

Type `bool`

**cons\_novalid**

whether no parameter value is valid

Type `bool`

**min**

minimal attribute value (if applicable)

Type `float` or `int`

**max**

maximal attribute value (if applicable)

Type `float` or `int`

**inc**  
minimal attribute increment value (if applicable)  
Type `float` or `int`

**cons\_excluded**  
list of special parameters which are within the range but are excluded

**cons\_included**  
list of special parameters which are outside the range but are included

**ivalues**  
list of possible integer values for enum attributes

**values**  
list of possible text values for enum attributes

**labels**  
dict {label: index} which shows all possible values of an enumerated attribute and their corresponding numerical values

**ilabels**  
dict {index: label} which shows labels corresponding to numerical values of an enumerated attribute

**update\_limits** (*force=False*)  
Update attribute constraints.  
  
If *force==False* and the constraints are permanent, skip the update.

**truncate\_value** (*value*)  
Truncate value to lie within attribute limits

**get\_value** (*enum\_as\_str=True*)  
Get attribute value.  
  
If *enum\_as\_str==True*, return enum-style values as strings; otherwise, return corresponding integer values.

**set\_value** (*value, truncate=True*)  
Get attribute value.  
  
If *truncate==True*, automatically truncate value to lie within allowed range.

**class** `pylablib.devices.PrincetonInstruments.picam.TDeviceInfo` (*name, serial\_number, model, interface*)

Bases: `tuple`

**count** ()  
Return number of occurrences of value.

**index** ()  
Return first index of value.  
  
Raises `ValueError` if the value is not present.

**interface**

**model**

**name**

**serial\_number**

```
class pylablib.devices.PrincetonInstruments.picam.TFrameInfo (frame_index, times-  
tamp_start, times-  
tamp_end, frames-  
tamp)  
  
Bases: tuple  
  
count ()  
    Return number of occurrences of value.  
  
frame_index  
  
framestamp  
  
index ()  
    Return first index of value.  
  
    Raises ValueError if the value is not present.  
  
timestamp_end  
  
timestamp_start  
  
class pylablib.devices.PrincetonInstruments.picam.PicamCamera (serial_number=None)  
Bases: pylablib.devices.interface.camera.IBinROICamera, pylablib.devices.  
interface.camera.IExposureCamera, pylablib.devices.interface.camera.  
IAttributeCamera  
  
Generic Picam camera interface.  
  
    Parameters serial_number – camera serial number; if None, connect to the first non-used  
        camera  
  
Error = <Mock name='mock.PicamError' id='140096294851280'>  
TimeoutError = <Mock spec='str' id='140096267056336'>  
  
open ()  
    Open connection to the camera  
  
close ()  
    Close connection to the camera  
  
is_opened ()  
    Check if the device is connected  
  
get_attribute_value (name, error_on_missing=True, default=None, enum_as_str=True)  
    Get value of an attribute with the given name.  
  
    If the value doesn't exist or can not be read and error_on_missing==True, raise error; otherwise,  
    return default. If default is not None, assume that error_on_missing==False. If name points at  
    a dictionary branch, return a dictionary with all values in this branch. If enum_as_str==True, return  
    enum-style values as strings; otherwise, return corresponding integer values.  
  
set_attribute_value (name, value, truncate=True, error_on_missing=True)  
    Set value of an attribute with the given name.  
  
    If the value doesn't exist or can not be written and error_on_missing==True, raise error; otherwise,  
    do nothing. If name points at a dictionary branch, set all values in this branch (in this case value must be  
    a dictionary). If truncate==True, truncate value to lie within attribute range.  
  
get_all_attribute_values (root=", enum_as_str=True)  
    Get values of all attributes with the given root
```

**set\_all\_attribute\_values** (*settings*, *root=""*, *truncate=True*)

Set values of all attributes with the given *root*.

If *truncate==True*, truncate value to lie within attribute range.

**get\_device\_info** ()

Get camera information.

Return tuple (*vendor*, *model*, *serial\_number*, *bus\_type*).

**get\_pixel\_size** ()

Get camera pixel size (in m)

**enable\_metadata** (*enable=True*)

Enable or disable metadata

**is\_metadata\_enabled** (*individual=False*)

Check if metadata is enabled.

If *individual==True*, return individual metadata info (*time\_stamp\_start*, *time\_stamp\_end*, *frame\_stamp*, *gate\_delay*, *modulation\_phase*). Otherwise, return simply True or False depending on whether the basic group (time- and frame-stamps) is enabled. In this case, if the value is inconsistent with either for the groups, fix this to be consistent.

**get\_exposure** ()

Get current exposure

**set\_exposure** (*exposure*)

Set camera exposure

**get\_frame\_period** (*per\_readout=False*)

Get frame period (time between two consecutive frames in the internal trigger mode)

If *per\_readout==True*, return time difference between readouts, which can contain more than one frame; otherwise, return average time per frame (keep in mind that the frames still come in single unbroken readout).

**get\_frame\_timings** (*per\_readout=False*)

Get acquisition timing.

Return tuple (*exposure*, *frame\_period*). If *per\_readout==True*, frame period difference between readouts, which can contain more than one frame; otherwise, it is the time per frame (keep in mind that the frames still come in single unbroken readout).

**get\_detector\_size** ()

Get camera detector size (in pixels) as a tuple (*width*, *height*)

**get\_roi** ()

Get current ROI.

Return tuple (*hstart*, *hend*, *vstart*, *vend*, *hbin*, *vbin*). *hstart* and *hend* specify horizontal image extent, *vstart* and *vend* specify vertical image extent (start is inclusive, stop is exclusive, starting from 0), *hbin* and *vbin* specify binning.

**set\_roi** (*hstart=0*, *hend=None*, *vstart=0*, *vend=None*, *hbin=1*, *vbin=1*)

Setup camera ROI.

*hstart* and *hend* specify horizontal image extent, *vstart* and *vend* specify vertical image extent (start is inclusive, stop is exclusive, starting from 0), *hbin* and *vbin* specify binning. By default, all non-supplied parameters take extreme values (0 for start, maximal for end, 1 for binning).

**get\_roi\_limits** (*hbin=1*, *vbin=1*)

Get the minimal and maximal ROI parameters.

Return tuple (hlim, vlim), where each element is in turn a limit 5-tuple (min, max, pstep, sstep, maxbin) with, correspondingly, minimal and maximal size, position and size step, and the maximal binning. In some cameras, the step and the minimal size depend on the binning, which can be supplied.

**setup\_acquisition** (mode='sequence', nframes=100)

Setup acquisition mode.

*mode* can be either "snap" (single frame or a fixed number of frames) or "sequence" (continuous acquisition). *nframes* sets up number of frame buffers. If there are multiple frames per readout, it still means the number of frames, and the number of readouts is set up to contain all required frames (e.g., 10 frames per readout and 15 frames result in 2 readouts).

**clear\_acquisition** ()

Clear acquisition settings

**start\_acquisition** (\*args, \*\*kwargs)

Start acquisition.

Can take the same keyword parameters as :meth:“*setup\_acquisition*. If the acquisition is not set up yet, set it up using the supplied parameters (use default of *setup\_acquisition* (), if the parameter is None). Otherwise, if any supplied parameters are different from the current ones, change them and reset the acquisition.

**stop\_acquisition** ()

Stop acquisition

**acquisition\_in\_progress** ()

Check if acquisition is in progress

**FrameTransferError**

alias of *pylablib.devices.interface.camera.DefaultFrameTransferError*

**class NoParameterCaller** (device, kind)

Bases: *object*

Class to simplify calling functions without a parameter

**apply\_settings** (settings)

Apply the settings.

*settings* is the dict {name: value} of the device available settings. Non-applicable settings are ignored.

**get\_acquisition\_parameters** ()

Get acquisition parameters.

Return dictionary {name: value}

**get\_all\_attributes** (copy=False)

Return a dictionary of all available attributes.

If copy==True, copy the dictionary; otherwise, return the internal dictionary structure (should not be modified).

**get\_attribute** (name, error\_on\_missing=True)

Get the camera attribute with the given name

**get\_data\_dimensions** ()

Get readout data dimensions (in pixels) as a tuple (width, height); take indexing mode into account

**get\_device\_variable** (key)

Get the value of a settings, status, or full info parameter

**get\_frame\_format()**

Get format for the returned images.

Can be "list" (list of 2D arrays), "array" (a single 3D array), or "chunks" (list of 3D "chunk" arrays; supported for some cameras and provides the best performance).

**get\_frame\_info\_fields()**

Get the names of frame info fields.

Applicable when frame info format (set by `set_frame_info_format()`) is "list" or "array".

**get\_frame\_info\_format()**

Get format of the frame info.

Can be "namedtuple" (potentially nested named tuples; convenient to get particular values), "list" (flat list of values, with field names are given by `get_frame_info_fields()`; convenient for building a table), "array" (same as "list", but with a numpy array, which is easier to use for "chunks" frame format), or "dict" (flat dictionary with the same fields as the "list" format; more resilient to future format changes)

**get\_frame\_info\_period()**

Get period of frame info acquisition.

Frame info might be skipped (set to `None`) except for frames which indices are divisible by *period*. Useful for certain cameras where acquiring frame info takes a lot of time and can reduce performance at higher frame rates. Note that this parameter can still be ignored (i.e., always set to 1) if the performance is not an issue for a given camera class.

**get\_frames\_status()**

Get acquisition and buffer status.

Return tuple (acquired, unread, skipped, size), where acquired is the total number of acquired frames, unread is the number of acquired but not read frames, skipped is the number of skipped (not read and then written over) frames, and `buffer_size` is the total buffer size (in frames).

**get\_full\_info(include=0)**

Get dict {name: value} containing full device information (including status and settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting `include=-10` queries all available variables, which is equivalent to `include="all"`.

**get\_full\_status(include=0)**

Get dict {name: value} containing the device status (including settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting `include=-10` queries all available variables, which is equivalent to `include="all"`.

**get\_image\_indexing()**

Get indexing for the returned images.

Can be "rct" (first index row, second index column, rows counted from the top), "rcb" (same as "rc", rows counted from the bottom), "xyt" (first index column, second index row, rows counted from the top), or "xyb" (same as "xyt", rows counted from the bottom)

**get\_new\_images\_range()**

Get the range of the new images.

Return tuple (*first*, *last*) with images range (first inclusive). If no images are available, return *None*. If some images were in the buffer were overwritten, exclude them from the range.

**get\_settings** (*include=0*)

Get dict {*name*: *value*} containing all the device settings.

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**grab** (*nframes=1*, *frame\_timeout=5.0*, *missing\_frame='skip'*, *return\_info=False*, *buff\_size=None*)

Snap *nframes* images (with preset image read mode parameters)

*buff\_size* determines buffer size (if *None*, use the default size). Timeout is specified for a single-frame acquisition, not for the whole acquisition time. *missing\_frame* determines what to do with frames which have been lost: can be "none" (replacing them with *None*), "zero" (replacing them with zero-filled frame), or "skip" (skipping them, while still keeping total returned frames number to *n*). If *return\_info==True*, return tuple (*frames*, *infos*), where *infos* is a list of frame info tuples (camera-dependent); if some frames are missing and *missing\_frame!="skip"*, the corresponding frame info is *None*.

**is\_acquisition\_setup** ()

Check if acquisition is set up.

If the camera does not support separate acquisition setup, always return *True*.

**pausing\_acquisition** (*clear=None*, *stop=True*, *setup\_after=None*, *start\_after=True*, *combine\_nested=True*)

Context manager which temporarily pauses acquisition during execution of *with* block.

Useful for applying certain settings which can't be changed during the acquisition. If *clear==True*, clear acquisition in addition to stopping (by default, use the class default specified as *\_clear\_pausing\_acquisition* attribute). If *stop==True*, stop the acquisition (if *clear==True*, stop regardless). If *setup\_after==True*, setup the acquisition after pause if necessary (*None* means setup only if clearing was required). If *start\_after==True*, start the acquisition after pause if necessary (*None* means start only if stopping was required). If *combine\_nested==True*, then any nested *pausing\_acquisition* calls will stop/clear acquisition as necessary, but won't setup/start it again until this *pausing\_acquisition* call is complete.

Yields tuple (*acq\_in\_progress*, *acq\_params*), which indicates whether acquisition is currently in progress, and what are the current acquisition parameters.

**read\_newest\_image** (*peek=False*, *return\_info=False*)

Read the newest un-read image.

If no un-read frames are available, return *None*. If *peek==True*, return the image but not mark it as read. If *return\_info==True*, return tuple (*frame*, *info*), where *info* is an info tuples (camera-dependent, see [read\\_multiple\\_images\(\)](#)).

**read\_oldest\_image** (*peek=False*, *return\_info=False*)

Read the oldest un-read image.

If no un-read frames are available, return *None*. If *peek==True*, return the image but not mark it as read. If *return\_info==True*, return tuple (*frame*, *info*), where *info* is an info tuples (camera-dependent, see [read\\_multiple\\_images\(\)](#)).

**set\_device\_variable** (*key*, *value*)

Set the value of a settings parameter

**set\_frame\_format** (*fnt*)

Set format for the returned images.



Can be "list" (list of 2D arrays), "array" (a single 3D array), "chunks" (list of 3D "chunk" arrays; supported for some cameras and provides the best performance), or "try\_chunks" (same as "chunks", but if chunks are not supported, set to "list" instead). If format is "chunks" and chunks are not supported by the camera, it results in one frame per chunk. Note that if the format is set to "array" or "chunks", the frame info format is also automatically set to "array". If the format is set to "chunks", then the image info is also returned in chunks form (list of 2D info arrays with the same length as the corresponding frame chunks).

**set\_frame\_info\_format** (*fmt, include\_fields=None*)

Set format of the frame info.

Can be "namedtuple" (potentially nested named tuples; convenient to get particular values), "list" (flat list of values, with field names are given by `get_frame_info_fields()`; convenient for building a table), "array" (same as "list", but with a numpy array, which is easier to use for "chunks" frame format), or "dict" (flat dictionary with the same fields as the "list" format; more resilient to future format changes) If *include\_fields* is not None, it specifies the fields included for non-"tuple" formats; note that order of *include\_fields* is ignored, and the resulting fields are always ordered same as in the original.

**set\_frame\_info\_period** (*period=1*)

Set period of frame info acquisition.

Frame info might be skipped (set to None) except for frames which indices are divisible by *period*. Useful for certain cameras where acquiring frame info takes a lot of time and can reduce performance at higher frame rates. Note that this parameter can still be ignored (i.e., always set to 1) if the performance is not an issue for a given camera class.

**set\_image\_indexing** (*indexing*)

Set up indexing for the returned images.

Can be "rct" (first index row, second index column, rows counted from the top), "rcb" (same as "rc", rows counted from the bottom), "xyt" (first index column, second index row, rows counted from the top), or "xyb" (same as "xyt", rows counted from the bottom)

**snap** (*timeout=5.0, return\_info=False*)

Snap a single frame

**wait\_for\_frame** (*since='lastread', nframes=1, timeout=20.0, error\_on\_stopped=False*)

Wait for one or several new camera frames.

*since* specifies the reference point for waiting to acquire *nframes* frames; can be "lastread" (from the last read frame), "lastwait" (wait for the last successful `wait_for_frame()` call), "now" (from the start of the current call), or "start" (from the acquisition start, i.e., wait until *nframes* frames have been acquired). *timeout* can be either a number, None (infinite timeout), or a tuple (*timeout, frame\_timeout*), in which case the call times out if the total time exceeds *timeout*, or a single frame wait exceeds *frame\_timeout*. If the call times out, raise `TimeoutError`. If *error\_on\_stopped*==True and the acquisition is not running, raise `Error`; otherwise, simply return False without waiting.

**read\_multiple\_images** (*rng=None, peek=False, missing\_frame='skip', return\_info=False, return\_rng=False*)

Read multiple images specified by *rng* (by default, all un-read images).

If *rng* is specified, it is a tuple (*first, last*) with images range (first inclusive). If no new frames are available, return an empty list; if no acquisition is running, return None. If *peek*==True, return images but not mark them as read. *missing\_frame* determines what to do with frames which are out of range (missing or lost): can be "none" (replacing them with None), "zero" (replacing them with zero-filled frame), or "skip" (skipping them). If *return\_info*==True, return tuple (*frames, infos*), where *infos* is a list of `TFrameInfo` instances describing frame index and frame metadata, which

contains start and stop timestamps, and framestamp; if some frames are missing and `missing_frame!="skip"`, the corresponding frame info is `None`. if `return_rng==True`, return the range covered resulting frames; if `missing_frame=="skip"`, the range can be smaller than the supplied *rng* if some frames are skipped.

## Module contents

### pylablib.devices.SiliconSoftware package

#### Submodules

#### pylablib.devices.SiliconSoftware.fgrab module

```
class pylablib.devices.SiliconSoftware.fgrab.TBoardInfo(name, full_name)
```

Bases: `tuple`

**count**()

Return number of occurrences of value.

**full\_name**

**index**()

Return first index of value.

Raises `ValueError` if the value is not present.

**name**

```
pylablib.devices.SiliconSoftware.fgrab.get_board_info(board)
```

Get board info for a given index (starting from 0)

```
pylablib.devices.SiliconSoftware.fgrab.list_boards()
```

List all boards available through Silicon Software interface

```
pylablib.devices.SiliconSoftware.fgrab.get_boards_number()
```

List number of connected Silicon Software boards

```
class pylablib.devices.SiliconSoftware.fgrab.TAppletInfo(name, file)
```

Bases: `tuple`

**count**()

Return number of occurrences of value.

**file**

**index**()

Return first index of value.

Raises `ValueError` if the value is not present.

**name**

```
class pylablib.devices.SiliconSoftware.fgrab.TFullAppletInfo(name, uid, desc,  
                                                             category, platform,  
                                                             tags, version, path,  
                                                             file, flags, info)
```

Bases: `tuple`

**category**

**count()**  
Return number of occurrences of value.

**desc**

**file**

**flags**

**index()**  
Return first index of value.  
Raises ValueError if the value is not present.

**info**

**name**

**path**

**platform**

**tags**

**uid**

**version**

`pylablib.devices.SiliconSoftware.fgrab.list_applets(board, full_desc=False, valid=True, on_board=False)`

List all applets available for this board.

*board* is the board index (starting from 0) given by its position in the list returned by `list_boards()`. If `full_desc==True`, return full description for each applet; otherwise, return only name and file name. If `valid==True`, list only valid and compatible applets; otherwise, list all applets. If `on_board==True`, list applets running on board; otherwise, list all applets contained in the system.

`pylablib.devices.SiliconSoftware.fgrab.get_applet_info(board, **kwargs)`

Return full information for an applet with the given parameters (e.g., name, or full path)

**class** `pylablib.devices.SiliconSoftware.fgrab.FGrabAttribute(fg, aid, port=0, system=False)`

Bases: `object`

Object representing an Silicon Software frame grabber parameter.

Allows to query and set values and get additional information. Usually created automatically by an `:class:` instance, but could be created manually.

#### Parameters

- **fg** – opened frame grabber handle
- **aid** – attribute ID
- **port** – camera port within the frame grabber
- **system** – if `True`, this is a system attribute; otherwise, it is a camera attribute

**name**  
attribute name

**kind**  
attribute kind; can be "i32", "i64", "u32", "u64", "f64", or "str"

**min**  
minimal attribute value (if applicable)

**Type** float or int

**max**

maximal attribute value (if applicable)

**Type** float or int

**inc**

minimal attribute increment value (if applicable)

**Type** float or int

**ivalues**

list of possible integer values for enum attributes

**values**

list of possible text values for enum attributes

**labels**

dict {label: index} which shows all possible values of an enumerated attribute and their corresponding numerical values

**ilabels**

dict {index: label} which shows labels corresponding to numerical values of an enumerated attribute

**update\_limits()**

Update minimal and maximal attribute limits and return tuple (min, max, inc)

**truncate\_value** (value)

Truncate value to lie within attribute limits

**get\_value** (enum\_as\_str=True)

Get attribute value.

If enum\_as\_str==True, return enum-style values as strings; otherwise, return corresponding integer values.

**set\_value** (value, truncate=True)

Get attribute value.

If truncate==True, automatically truncate value to lie within allowed range.

**class** pylablib.devices.SiliconSoftware.fgrab.**TDeviceInfo** (applet\_info, system\_info, software\_version)

Bases: tuple

**applet\_info**

**count** ()

Return number of occurrences of value.

**index** ()

Return first index of value.

Raises ValueError if the value is not present.

**software\_version**

**system\_info**

**class** pylablib.devices.SiliconSoftware.fgrab.**TFrameInfo** (frame\_index, frames-tamp, timestamp, times-tamp\_long)

Bases: tuple

**count()**  
Return number of occurrences of value.

**frame\_index**

**framestamp**

**index()**  
Return first index of value.

Raises ValueError if the value is not present.

**timestamp**

**timestamp\_long**

```
class pylablib.devices.SiliconSoftware.fgrab.SiliconSoftwareFrameGrabber (siso_board=0,
                                                                    siso_applet='DualAreaGr
                                                                    siso_port=0,
                                                                    siso_detector_size=None,
                                                                    do_open=True,
                                                                    **kwargs)
```

Bases: `pylablib.devices.interface.camera.IGrabberAttributeCamera`, `pylablib.devices.interface.camera.IROICamera`

Generic Silicon Software frame grabber interface.

Compared to `SiliconSoftwareCamera`, has more permissive initialization arguments, which simplifies its use as a base class for expanded cameras.

#### Parameters

- **siso\_board** – board index, starting from 0; available boards can be learned by `list_boards()`
- **siso\_applet** – applet name, which can be learned by `list_applets()`; usually, a simple applet like "DualLineGray16" or "MediumLineGray16" are most appropriate; can be either an applet name, or a direct path to the applet DLL
- **siso\_port** – port number, if several ports are supported by the camera and the applet
- **siso\_detector\_size** – if not None, can specify the maximal detector size; by default, use the maximal available for the frame grabber (usually, 16384x16384)

```
Error = <Mock name='mock.SiliconSoftwareError' id='140096270265936'>
```

```
TimeoutError = <Mock spec='str' id='140096289529296'>
```

**open()**  
Open connection to the camera

**close()**  
Close connection to the camera

**is\_opened()**  
Check if the device is connected

**get\_all\_grabber\_attribute\_values** (*root*="", \*\**kwargs*)  
Get values of all frame grabber attributes with the given *root*.  
Additional arguments are passed to `get_value` methods of individual attributes.

**set\_all\_grabber\_attribute\_values** (*settings*, *root*="", \*\**kwargs*)  
Set values of all frame grabber attributes with the given *root*.  
Additional arguments are passed to `set_value` methods of individual attributes.

**get\_system\_info()**

Get the dictionary with all system information parameters

**get\_genicam\_info\_xml()**

Get description in Genicam-compatible XML format

**get\_device\_info()**

Get camera model data.

Return tuple (applet\_info, system\_info, software\_version) with the board serial number and an the interface type (e.g., "1430" for NI PCIe-1430)

**set\_frame\_merge** (*frame\_merge=1*)

**get\_detector\_size()**

Get camera detector size (in pixels) as a tuple (width, height)

**get\_grabber\_detector\_size()**

Get camera detector size (in pixels) as a tuple (width, height)

**get\_roi()**

Get current ROI.

Return tuple (hstart, hend, vstart, vend). *hstart* and *hend* specify horizontal image extent, *vstart* and *vend* specify vertical image extent (start is inclusive, stop is exclusive, starting from 0).

**get\_grabber\_roi()**

Get current ROI.

Return tuple (hstart, hend, vstart, vend). *hstart* and *hend* specify horizontal image extent, *vstart* and *vend* specify vertical image extent (start is inclusive, stop is exclusive, starting from 0).

**set\_roi** (*hstart=0, hend=None, vstart=0, vend=None*)

Setup camera ROI.

*hstart* and *hend* specify horizontal image extent, *vstart* and *vend* specify vertical image extent (start is inclusive, stop is exclusive, starting from 0). By default, all non-supplied parameters take extreme values (0 for start, maximal for end).

**set\_grabber\_roi** (*hstart=0, hend=None, vstart=0, vend=None*)

Setup camera ROI.

*hstart* and *hend* specify horizontal image extent, *vstart* and *vend* specify vertical image extent (start is inclusive, stop is exclusive, starting from 0). By default, all non-supplied parameters take extreme values (0 for start, maximal for end).

**get\_roi\_limits** (*hbin=1, vbin=1*)

Get the minimal and maximal ROI parameters.

Return tuple (hlim, vlim), where each element is in turn a limit 5-tuple (min, max, pstep, sstep, maxbin) with, correspondingly, minimal and maximal size, position and size step, and the maximal binning (fixed to 1 if not binning is allowed). In some cameras, the step and the minimal size depend on the binning, which can be supplied.

**get\_grabber\_roi\_limits** (*hbin=1, vbin=1*)

Get the minimal and maximal ROI parameters.

Return tuple (hlim, vlim), where each element is in turn a limit 5-tuple (min, max, pstep, sstep, maxbin) with, correspondingly, minimal and maximal size, position and size step, and the maximal binning (fixed to 1 if not binning is allowed). In some cameras, the step and the minimal size depend on the binning, which can be supplied.

**setup\_camlink\_pixel\_format** (*bits\_per\_pixel=8, taps=1, output\_fmt=None, fmt=None*)

Set up CameraLink pixel format.

If *fmt* is *None*, use supplied *bits\_per\_pixel* (8, 10, 12, 14, or 16) and *taps* (1 or 2) to figure out the format; otherwise, *fmt* should be a numerical (e.g., 210) or string (e.g., "FG\_CL\_MEDIUM\_10\_BIT") format. *output\_fmt* specifies the result frame format; if *None*, use grayscale with the given *bits\_per\_pixel* if *fmt* is *None*, or 16 bit grayscale otherwise.

**get\_camlink\_pixel\_format** ()

Get CamLink pixel format and the output pixel format as a tuple

**get\_available\_camlink\_pixel\_formats** ()

Get all available CamLink pixel formats and the output pixel formats as a tuple of 2 lists

**setup\_acquisition** (*mode='sequence', nframes=100*)

Setup acquisition mode.

*mode* can be either "snap" (single frame or a fixed number of frames) or "sequence" (continuous acquisition). (note that *IMAQCamera.acquisition\_in\_progress()* would still return *True* in this case, even though new frames are no longer acquired). *nframes* sets up number of frame buffers.

**clear\_acquisition** ()

Clear all acquisition details and free all buffers

**start\_acquisition** (\*args, \*\*kwargs)

Start acquisition.

Can take the same keyword parameters as :meth:“*setup\_acquisition*. If the acquisition is not set up yet, set it up using the supplied parameters (use default of *setup\_acquisition()*, if the parameter is *None*). Otherwise, if any supplied parameters are different from the current ones, change them and reset the acquisition.

**stop\_acquisition** ()

Stop acquisition

**acquisition\_in\_progress** ()

Check if acquisition is in progress

**FrameTransferError**

alias of *pylablib.devices.interface.camera.DefaultFrameTransferError*

**class NoParameterCaller** (*device, kind*)

Bases: *object*

Class to simplify calling functions without a parameter

**apply\_settings** (*settings*)

Apply the settings.

*settings* is the dict {name: value} of the device available settings. Non-applicable settings are ignored.

**get\_acquisition\_parameters** ()

Get acquisition parameters.

Return dictionary {name: value}

**get\_all\_grabber\_attributes** (*copy=False*)

Return a dictionary of all available frame grabber grabber\_attributes.

If *copy==True*, copy the dictionary; otherwise, return the internal dictionary structure (should not be modified).

**get\_data\_dimensions()**

Get readout data dimensions (in pixels) as a tuple (width, height); take indexing mode into account

**get\_device\_variable(key)**

Get the value of a settings, status, or full info parameter

**get\_frame\_format()**

Get format for the returned images.

Can be "list" (list of 2D arrays), "array" (a single 3D array), or "chunks" (list of 3D "chunk" arrays; supported for some cameras and provides the best performance).

**get\_frame\_info\_fields()**

Get the names of frame info fields.

Applicable when frame info format (set by `set_frame_info_format()`) is "list" or "array".

**get\_frame\_info\_format()**

Get format of the frame info.

Can be "namedtuple" (potentially nested named tuples; convenient to get particular values), "list" (flat list of values, with field names are given by `get_frame_info_fields()`; convenient for building a table), "array" (same as "list", but with a numpy array, which is easier to use for "chunks" frame format), or "dict" (flat dictionary with the same fields as the "list" format; more resilient to future format changes)

**get\_frame\_info\_period()**

Get period of frame info acquisition.

Frame info might be skipped (set to None) except for frames which indices are divisible by *period*. Useful for certain cameras where acquiring frame info takes a lot of time and can reduce performance at higher frame rates. Note that this parameter can still be ignored (i.e., always set to 1) if the performance is not an issue for a given camera class.

**get\_frames\_status()**

Get acquisition and buffer status.

Return tuple (acquired, unread, skipped, size), where acquired is the total number of acquired frames, unread is the number of acquired but not read frames, skipped is the number of skipped (not read and then written over) frames, and buffer\_size is the total buffer size (in frames).

**get\_full\_info(include=0)**

Get dict {name: value} containing full device information (including status and settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include*=-10 queries all available variables, which is equivalent to *include*="all".

**get\_full\_status(include=0)**

Get dict {name: value} containing the device status (including settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include*=-10 queries all available variables, which is equivalent to *include*="all".

**get\_grabber\_attribute(name, error\_on\_missing=True)**

Get the camera attribute with the given name

**get\_grabber\_attribute\_value(name, error\_on\_missing=True, default=None, \*\*kwargs)**

Get value of a frame grabber attribute with the given name.



If the value doesn't exist and `error_on_missing==True`, raise error; otherwise, return *default*. If *default* is not `None`, automatically assume that `error_on_missing==False`. If *name* points at a dictionary branch, return a dictionary with all values in this branch. Additional arguments are passed to `get_value` methods of the individual attribute.

#### **get\_image\_indexing()**

Get indexing for the returned images.

Can be "rct" (first index row, second index column, rows counted from the top), "rcb" (same as "rc", rows counted from the bottom), "xyt" (first index column, second index row, rows counted from the top), or "xyb" (same as "xyt", rows counted from the bottom)

#### **get\_new\_images\_range()**

Get the range of the new images.

Return tuple (*first*, *last*) with images range (first inclusive). If no images are available, return `None`. If some images were in the buffer were overwritten, exclude them from the range.

#### **get\_settings(include=0)**

Get dict {*name*: *value*} containing all the device settings.

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include*=-10 queries all available variables, which is equivalent to *include*="all".

#### **grab(nframes=1, frame\_timeout=5.0, missing\_frame='skip', return\_info=False, buff\_size=None)**

Snap *nframes* images (with preset image read mode parameters)

*buff\_size* determines buffer size (if `None`, use the default size). Timeout is specified for a single-frame acquisition, not for the whole acquisition time. *missing\_frame* determines what to do with frames which have been lost: can be "none" (replacing them with `None`), "zero" (replacing them with zero-filled frame), or "skip" (skipping them, while still keeping total returned frames number to *n*). If *return\_info*==`True`, return tuple (*frames*, *infos*), where *infos* is a list of frame info tuples (camera-dependent); if some frames are missing and *missing\_frame*!="skip", the corresponding frame info is `None`.

#### **is\_acquisition\_setup()**

Check if acquisition is set up.

If the camera does not support separate acquisition setup, always return `True`.

#### **pausing\_acquisition(clear=None, stop=True, setup\_after=None, start\_after=True, combine\_nested=True)**

Context manager which temporarily pauses acquisition during execution of `with` block.

Useful for applying certain settings which can't be changed during the acquisition. If *clear*==`True`, clear acquisition in addition to stopping (by default, use the class default specified as `_clear_pausing_acquisition` attribute). If *stop*==`True`, stop the acquisition (if *clear*==`True`, stop regardless). If *setup\_after*==`True`, setup the acquisition after pause if necessary (`None` means setup only if clearing was required). If *start\_after*==`True`, start the acquisition after pause if necessary (`None` means start only if stopping was required). If *combine\_nested*==`True`, then any nested `pausing_acquisition` calls will stop/clear acquisition as necessary, but won't setup/start it again until this `pausing_acquisition` call is complete.

Yields tuple (*acq\_in\_progress*, *acq\_params*), which indicates whether acquisition is currently in progress, and what are the current acquisition parameters.

#### **read\_multiple\_images(rng=None, peek=False, missing\_frame='skip', return\_info=False, return\_rng=False)**

Read multiple images specified by *rng* (by default, all un-read images).

If *rng* is specified, it is a tuple (*first*, *last*) with images range (first inclusive). If no new frames are available, return an empty list; if no acquisition is running, return *None*. If *peek*==True, return images but not mark them as read. *missing\_frame* determines what to do with frames which are out of range (missing or lost): can be "none" (replacing them with *None*), "zero" (replacing them with zero-filled frame), or "skip" (skipping them). If *return\_info*==True, return tuple (*frames*, *infos*), where *infos* is a list of frame info tuples (camera-dependent, by default, only the frame index); if some frames are missing and *missing\_frame*!="skip", the corresponding frame info is *None*. if *return\_rng*==True, return the range covered resulting frames; if *missing\_frame*=="skip", the range can be smaller than the supplied *rng* if some frames are skipped.

**read\_newest\_image** (*peek=False*, *return\_info=False*)

Read the newest un-read image.

If no un-read frames are available, return *None*. If *peek*==True, return the image but not mark it as read. If *return\_info*==True, return tuple (*frame*, *info*), where *info* is an info tuples (camera-dependent, see [read\\_multiple\\_images\(\)](#)).

**read\_oldest\_image** (*peek=False*, *return\_info=False*)

Read the oldest un-read image.

If no un-read frames are available, return *None*. If *peek*==True, return the image but not mark it as read. If *return\_info*==True, return tuple (*frame*, *info*), where *info* is an info tuples (camera-dependent, see [read\\_multiple\\_images\(\)](#)).

**set\_device\_variable** (*key*, *value*)

Set the value of a settings parameter

**set\_frame\_format** (*fmt*)

Set format for the returned images.

Can be "list" (list of 2D arrays), "array" (a single 3D array), "chunks" (list of 3D "chunk" arrays; supported for some cameras and provides the best performance), or "try\_chunks" (same as "chunks", but if chunks are not supported, set to "list" instead). If format is "chunks" and chunks are not supported by the camera, it results in one frame per chunk. Note that if the format is set to "array" or "chunks", the frame info format is also automatically set to "array". If the format is set to "chunks", then the image info is also returned in chunks form (list of 2D info arrays with the same length as the corresponding frame chunks).

**set\_frame\_info\_format** (*fmt*, *include\_fields=None*)

Set format of the frame info.

Can be "namedtuple" (potentially nested named tuples; convenient to get particular values), "list" (flat list of values, with field names are given by [get\\_frame\\_info\\_fields\(\)](#); convenient for building a table), "array" (same as "list", but with a numpy array, which is easier to use for "chunks" frame format), or "dict" (flat dictionary with the same fields as the "list" format; more resilient to future format changes) If *include\_fields* is not *None*, it specifies the fields included for non-"tuple" formats; note that order or *include\_fields* is ignored, and the resulting fields are always ordered same as in the original.

**set\_frame\_info\_period** (*period=1*)

Set period of frame info acquisition.

Frame info might be skipped (set to *None*) except for frames which indices are divisible by *period*. Useful for certain cameras where acquiring frame info takes a lot of time and can reduce performance at higher frame rates. Note that this parameter can still be ignored (i.e., always set to 1) if the performance is not an issue for a given camera class.

**set\_grabber\_attribute\_value** (*name*, *value*, *error\_on\_missing=True*, *\*\*kwargs*)

Set value of a frame grabber attribute with the given name.

If the value doesn't exist and `error_on_missing==True`, raise error; otherwise, do nothing. If *name* points at a dictionary branch, set all values in this branch (in this case *value* must be a dictionary). Additional arguments are passed to `set_value` methods of the individual attribute.

**set\_image\_indexing** (*indexing*)

Set up indexing for the returned images.

Can be "rct" (first index row, second index column, rows counted from the top), "rcb" (same as "rc", rows counted from the bottom), "xyt" (first index column, second index row, rows counted from the top), or "xyb" (same as "xyt", rows counted from the bottom)

**snap** (*timeout=5.0, return\_info=False*)

Snap a single frame

**wait\_for\_frame** (*since='lastread', nframes=1, timeout=20.0, error\_on\_stopped=False*)

Wait for one or several new camera frames.

*since* specifies the reference point for waiting to acquire *nframes* frames; can be "lastread" (from the last read frame), "lastwait" (wait for the last successful `wait_for_frame()` call), "now" (from the start of the current call), or "start" (from the acquisition start, i.e., wait until *nframes* frames have been acquired). *timeout* can be either a number, None (infinite timeout), or a tuple (*timeout*, *frame\_timeout*), in which case the call times out if the total time exceeds *timeout*, or a single frame wait exceeds *frame\_timeout*. If the call times out, raise `TimeoutError`. If `error_on_stopped==True` and the acquisition is not running, raise `Error`; otherwise, simply return False without waiting.

```
class pylablib.devices.SiliconSoftware.fgrab.SiliconSoftwareCamera (board,  
                                                                applet,  
                                                                port=0,  
                                                                detec-  
                                                                tor_size=None)
```

Bases: `pylablib.devices.SiliconSoftware.fgrab.SiliconSoftwareFrameGrabber`

Generic Silicon Software frame grabber interface.

#### Parameters

- **board** – board index, starting from 0; available boards can be learned by `list_boards()`
- **applet** – applet name, which can be learned by `list_applets()`; usually, a simple applet like "DualLineGray16" or "MediumLineGray16" are most appropriate; can be either an applet name, or a direct path to the applet DLL
- **port** – port number, if several ports are supported by the camera and the applet
- **detector\_size** – if not None, can specify the maximal detector size; by default, use the maximal available for the frame grabber (usually, 16384x16384)

```
Error = <Mock name='mock.SiliconSoftwareError' id='140096270265936'>
```

```
FrameTransferError
```

alias of `pylablib.devices.interface.camera.DefaultFrameTransferError`

```
class NoParameterCaller (device, kind)
```

Bases: `object`

Class to simplify calling functions without a parameter

```
TimeoutError = <Mock spec='str' id='140096289529296'>
```

```
acquisition_in_progress ()
```

Check if acquisition is in progress

**apply\_settings** (*settings*)

Apply the settings.

*settings* is the dict {name: value} of the device available settings. Non-applicable settings are ignored.

**clear\_acquisition** ()

Clear all acquisition details and free all buffers

**close** ()

Close connection to the camera

**get\_acquisition\_parameters** ()

Get acquisition parameters.

Return dictionary {name: value}

**get\_all\_grabber\_attribute\_values** (*root*=", \*\**kwargs*)

Get values of all frame grabber attributes with the given *root*.

Additional arguments are passed to *get\_value* methods of individual attributes.

**get\_all\_grabber\_attributes** (*copy*=*False*)

Return a dictionary of all available frame grabber *grabber\_attributes*.

If *copy*==True, copy the dictionary; otherwise, return the internal dictionary structure (should not be modified).

**get\_available\_camlink\_pixel\_formats** ()

Get all available CamLink pixel formats and the output pixel formats as a tuple of 2 lists

**get\_camlink\_pixel\_format** ()

Get CamLink pixel format and the output pixel format as a tuple

**get\_data\_dimensions** ()

Get readout data dimensions (in pixels) as a tuple (*width*, *height*); take indexing mode into account

**get\_detector\_size** ()

Get camera detector size (in pixels) as a tuple (*width*, *height*)

**get\_device\_info** ()

Get camera model data.

Return tuple (*applet\_info*, *system\_info*, *software\_version*) with the board serial number and an the interface type (e.g., "1430" for NI PCIe-1430)

**get\_device\_variable** (*key*)

Get the value of a settings, status, or full info parameter

**get\_frame\_format** ()

Get format for the returned images.

Can be "list" (list of 2D arrays), "array" (a single 3D array), or "chunks" (list of 3D "chunk" arrays; supported for some cameras and provides the best performance).

**get\_frame\_info\_fields** ()

Get the names of frame info fields.

Applicable when frame info format (set by *set\_frame\_info\_format* ()) is "list" or "array".

**get\_frame\_info\_format** ()

Get format of the frame info.

Can be "namedtuple" (potentially nested named tuples; convenient to get particular values), "list" (flat list of values, with field names are given by `get_frame_info_fields()`; convenient for building a table), "array" (same as "list", but with a numpy array, which is easier to use for "chunks" frame format), or "dict" (flat dictionary with the same fields as the "list" format; more resilient to future format changes)

#### **get\_frame\_info\_period()**

Get period of frame info acquisition.

Frame info might be skipped (set to `None`) except for frames which indices are divisible by *period*. Useful for certain cameras where acquiring frame info takes a lot of time and can reduce performance at higher frame rates. Note that this parameter can still be ignored (i.e., always set to 1) if the performance is not an issue for a given camera class.

#### **get\_frames\_status()**

Get acquisition and buffer status.

Return tuple (*acquired*, *unread*, *skipped*, *size*), where *acquired* is the total number of acquired frames, *unread* is the number of acquired but not read frames, *skipped* is the number of skipped (not read and then written over) frames, and *buffer\_size* is the total buffer size (in frames).

#### **get\_full\_info(include=0)**

Get dict {*name*: *value*} containing full device information (including status and settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include*=-10 queries all available variables, which is equivalent to *include*="all".

#### **get\_full\_status(include=0)**

Get dict {*name*: *value*} containing the device status (including settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include*=-10 queries all available variables, which is equivalent to *include*="all".

#### **get\_genicam\_info\_xml()**

Get description in Genicam-compatible XML format

#### **get\_grabber\_attribute(name, error\_on\_missing=True)**

Get the camera attribute with the given name

#### **get\_grabber\_attribute\_value(name, error\_on\_missing=True, default=None, \*\*kwargs)**

Get value of a frame grabber attribute with the given name.

If the value doesn't exist and *error\_on\_missing*==True, raise error; otherwise, return *default*. If *default* is not None, automatically assume that *error\_on\_missing*==False. If *name* points at a dictionary branch, return a dictionary with all values in this branch. Additional arguments are passed to *get\_value* methods of the individual attribute.

#### **get\_grabber\_detector\_size()**

Get camera detector size (in pixels) as a tuple (*width*, *height*)

#### **get\_grabber\_roi()**

Get current ROI.

Return tuple (*hstart*, *hend*, *vstart*, *vend*). *hstart* and *hend* specify horizontal image extent, *vstart* and *vend* specify vertical image extent (start is inclusive, stop is exclusive, starting from 0).

#### **get\_grabber\_roi\_limits(hbin=1, vbin=1)**

Get the minimal and maximal ROI parameters.

Return tuple (*hlim*, *vlim*), where each element is in turn a limit 5-tuple (*min*, *max*, *pstep*, *sstep*, *maxbin*) with, correspondingly, minimal and maximal size, position and size step, and the maximal binning (fixed to 1 if not binning is allowed). In some cameras, the step and the minimal size depend on the binning, which can be supplied.

**get\_image\_indexing()**

Get indexing for the returned images.

Can be "rct" (first index row, second index column, rows counted from the top), "rcb" (same as "rc", rows counted from the bottom), "xyt" (first index column, second index row, rows counted from the top), or "xyb" (same as "xyt", rows counted from the bottom)

**get\_new\_images\_range()**

Get the range of the new images.

Return tuple (*first*, *last*) with images range (first inclusive). If no images are available, return None. If some images were in the buffer were overwritten, exclude them from the range.

**get\_roi()**

Get current ROI.

Return tuple (*hstart*, *hend*, *vstart*, *vend*). *hstart* and *hend* specify horizontal image extent, *vstart* and *vend* specify vertical image extent (start is inclusive, stop is exclusive, starting from 0).

**get\_roi\_limits(hbin=1, vbin=1)**

Get the minimal and maximal ROI parameters.

Return tuple (*hlim*, *vlim*), where each element is in turn a limit 5-tuple (*min*, *max*, *pstep*, *sstep*, *maxbin*) with, correspondingly, minimal and maximal size, position and size step, and the maximal binning (fixed to 1 if not binning is allowed). In some cameras, the step and the minimal size depend on the binning, which can be supplied.

**get\_settings(include=0)**

Get dict {*name*: *value*} containing all the device settings.

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include*=-10 queries all available variables, which is equivalent to *include*="all".

**get\_system\_info()**

Get the dictionary with all system information parameters

**grab(nframes=1, frame\_timeout=5.0, missing\_frame='skip', return\_info=False, buff\_size=None)**

Snap *nframes* images (with preset image read mode parameters)

*buff\_size* determines buffer size (if None, use the default size). Timeout is specified for a single-frame acquisition, not for the whole acquisition time. *missing\_frame* determines what to do with frames which have been lost: can be "none" (replacing them with None), "zero" (replacing them with zero-filled frame), or "skip" (skipping them, while still keeping total returned frames number to *n*). If *return\_info*==True, return tuple (*frames*, *infos*), where *infos* is a list of frame info tuples (camera-dependent); if some frames are missing and *missing\_frame*!="skip", the corresponding frame info is None.

**is\_acquisition\_setup()**

Check if acquisition is set up.

If the camera does not support separate acquisition setup, always return True.

**is\_opened()**

Check if the device is connected

**open()**

Open connection to the camera

**pausing\_acquisition**(*clear=None, stop=True, setup\_after=None, start\_after=True, combine\_nested=True*)

Context manager which temporarily pauses acquisition during execution of `with` block.

Useful for applying certain settings which can't be changed during the acquisition. If `clear==True`, clear acquisition in addition to stopping (by default, use the class default specified as `_clear_pausing_acquisition` attribute). If `stop==True`, stop the acquisition (if `clear==True`, stop regardless). If `setup_after==True`, setup the acquisition after pause if necessary (None means setup only if clearing was required). If `start_after==True`, start the acquisition after pause if necessary (None means start only if stopping was required). If `combine_nested==True`, then any nested `pausing_acquisition` calls will stop/clear acquisition as necessary, but won't setup/start it again until this `pausing_acquisition` call is complete.

Yields tuple (`acq_in_progress, acq_params`), which indicates whether acquisition is currently in progress, and what are the current acquisition parameters.

**read\_multiple\_images**(*rng=None, peek=False, missing\_frame='skip', return\_info=False, return\_rng=False*)

Read multiple images specified by `rng` (by default, all un-read images).

If `rng` is specified, it is a tuple (`first, last`) with images range (first inclusive). If no new frames are available, return an empty list; if no acquisition is running, return None. If `peek==True`, return images but not mark them as read. `missing_frame` determines what to do with frames which are out of range (missing or lost): can be "none" (replacing them with None), "zero" (replacing them with zero-filled frame), or "skip" (skipping them). If `return_info==True`, return tuple (`frames, infos`), where `infos` is a list of frame info tuples (camera-dependent, by default, only the frame index); if some frames are missing and `missing_frame!="skip"`, the corresponding frame info is None. If `return_rng==True`, return the range covered resulting frames; if `missing_frame=="skip"`, the range can be smaller than the supplied `rng` if some frames are skipped.

**read\_newest\_image**(*peek=False, return\_info=False*)

Read the newest un-read image.

If no un-read frames are available, return None. If `peek==True`, return the image but not mark it as read. If `return_info==True`, return tuple (`frame, info`), where `info` is an info tuples (camera-dependent, see `read_multiple_images()`).

**read\_oldest\_image**(*peek=False, return\_info=False*)

Read the oldest un-read image.

If no un-read frames are available, return None. If `peek==True`, return the image but not mark it as read. If `return_info==True`, return tuple (`frame, info`), where `info` is an info tuples (camera-dependent, see `read_multiple_images()`).

**set\_all\_grabber\_attribute\_values**(*settings, root="", \*\*kwargs*)

Set values of all frame grabber attributes with the given `root`.

Additional arguments are passed to `set_value` methods of individual attributes.

**set\_device\_variable**(*key, value*)

Set the value of a settings parameter

**set\_frame\_format**(*fmt*)

Set format for the returned images.

Can be "list" (list of 2D arrays), "array" (a single 3D array), "chunks" (list of 3D "chunk" arrays; supported for some cameras and provides the best performance), or "try\_chunks" (same as "chunks", but if chunks are not supported, set to "list" instead). If format is "chunks" and

chunks are not supported by the camera, it results in one frame per chunk. Note that if the format is set to "array" or "chunks", the frame info format is also automatically set to "array". If the format is set to "chunks", then the image info is also returned in chunks form (list of 2D info arrays with the same length as the corresponding frame chunks).

**set\_frame\_info\_format** (*fmt, include\_fields=None*)

Set format of the frame info.

Can be "namedtuple" (potentially nested named tuples; convenient to get particular values), "list" (flat list of values, with field names are given by `get_frame_info_fields()`; convenient for building a table), "array" (same as "list", but with a numpy array, which is easier to use for "chunks" frame format), or "dict" (flat dictionary with the same fields as the "list" format; more resilient to future format changes) If *include\_fields* is not None, it specifies the fields included for non-"tuple" formats; note that order of *include\_fields* is ignored, and the resulting fields are always ordered same as in the original.

**set\_frame\_info\_period** (*period=1*)

Set period of frame info acquisition.

Frame info might be skipped (set to None) except for frames which indices are divisible by *period*. Useful for certain cameras where acquiring frame info takes a lot of time and can reduce performance at higher frame rates. Note that this parameter can still be ignored (i.e., always set to 1) if the performance is not an issue for a given camera class.

**set\_frame\_merge** (*frame\_merge=1*)

**set\_grabber\_attribute\_value** (*name, value, error\_on\_missing=True, \*\*kwargs*)

Set value of a frame grabber attribute with the given name.

If the value doesn't exist and *error\_on\_missing==True*, raise error; otherwise, do nothing. If *name* points at a dictionary branch, set all values in this branch (in this case *value* must be a dictionary). Additional arguments are passed to *set\_value* methods of the individual attribute.

**set\_grabber\_roi** (*hstart=0, hend=None, vstart=0, vend=None*)

Setup camera ROI.

*hstart* and *hend* specify horizontal image extent, *vstart* and *vend* specify vertical image extent (start is inclusive, stop is exclusive, starting from 0). By default, all non-supplied parameters take extreme values (0 for start, maximal for end).

**set\_image\_indexing** (*indexing*)

Set up indexing for the returned images.

Can be "rct" (first index row, second index column, rows counted from the top), "rcb" (same as "rc", rows counted from the bottom), "xyt" (first index column, second index row, rows counted from the top), or "xyb" (same as "xyt", rows counted from the bottom)

**set\_roi** (*hstart=0, hend=None, vstart=0, vend=None*)

Setup camera ROI.

*hstart* and *hend* specify horizontal image extent, *vstart* and *vend* specify vertical image extent (start is inclusive, stop is exclusive, starting from 0). By default, all non-supplied parameters take extreme values (0 for start, maximal for end).

**setup\_acquisition** (*mode='sequence', nframes=100*)

Setup acquisition mode.

*mode* can be either "snap" (single frame or a fixed number of frames) or "sequence" (continuous acquisition). (note that `IMAQCamera.acquisition_in_progress()` would still return True in this case, even though new frames are no longer acquired). *nframes* sets up number of frame buffers.



**setup\_camlink\_pixel\_format** (*bits\_per\_pixel=8, taps=1, output\_fmt=None, fmt=None*)

Set up CameraLink pixel format.

If *fmt* is *None*, use supplied *bits\_per\_pixel* (8, 10, 12, 14, or 16) and *taps* (1 or 2) to figure out the format; otherwise, *fmt* should be a numerical (e.g., 210) or string (e.g., "FG\_CL\_MEDIUM\_10\_BIT") format. *output\_fmt* specifies the result frame format; if *None*, use grayscale with the given *bits\_per\_pixel* if *fmt* is *None*, or 16 bit grayscale otherwise.

**snap** (*timeout=5.0, return\_info=False*)

Snap a single frame

**start\_acquisition** (*\*args, \*\*kwargs*)

Start acquisition.

Can take the same keyword parameters as *meth:“setup\_acquisition*. If the acquisition is not set up yet, set it up using the supplied parameters (use default of *setup\_acquisition()*, if the parameter is *None*). Otherwise, if any supplied parameters are different from the current ones, change them and reset the acquisition.

**stop\_acquisition** ()

Stop acquisition

**wait\_for\_frame** (*since='lastread', nframes=1, timeout=20.0, error\_on\_stopped=False*)

Wait for one or several new camera frames.

*since* specifies the reference point for waiting to acquire *nframes* frames; can be “lastread” (from the last read frame), “lastwait” (wait for the last successful *wait\_for\_frame()* call), “now” (from the start of the current call), or “start” (from the acquisition start, i.e., wait until *nframes* frames have been acquired). *timeout* can be either a number, *None* (infinite timeout), or a tuple (*timeout, frame\_timeout*), in which case the call times out if the total time exceeds *timeout*, or a single frame wait exceeds *frame\_timeout*. If the call times out, raise *TimeoutError*. If *error\_on\_stopped==True* and the acquisition is not running, raise *Error*; otherwise, simply return *False* without waiting.

## Module contents

### pylablib.devices.Sirah package

#### Submodules

### pylablib.devices.Sirah.Matisse module

**exception** `pylablib.devices.Sirah.Matisse.GenericSirahError`

Bases: `pylablib.core.devio.base.DeviceError`

Generic Sirah error

**args**

**with\_traceback** ()

Exception.with\_traceback(tb) – set self.\_\_traceback\_\_ to tb and return self.

**exception** `pylablib.devices.Sirah.Matisse.GenericSirahBackendError` (*exc*)

Bases: `pylablib.devices.Sirah.Matisse.GenericSirahError, pylablib.core.devio.comm_backend.DeviceBackendError`

Sirah backend communication error

**args**

**with\_traceback()**

Exception.with\_traceback(tb) – set self.\_\_traceback\_\_ to tb and return self.

**class** pylablib.devices.Sirah.Matisse.**TThinnetCtlParameters** (*setpoint, P, I, avg*)

Bases: `tuple`

**I**

**P**

**avg**

**count()**

Return number of occurrences of value.

**index()**

Return first index of value.

Raises ValueError if the value is not present.

**setpoint**

**class** pylablib.devices.Sirah.Matisse.**TPiezoetDriveParameters** (*amplitude, rate, oversamp*)

Bases: `tuple`

**amplitude**

**count()**

Return number of occurrences of value.

**index()**

Return first index of value.

Raises ValueError if the value is not present.

**oversamp**

**rate**

**class** pylablib.devices.Sirah.Matisse.**TPiezoetFeedbackParameters** (*P, avg, phase*)

Bases: `tuple`

**P**

**avg**

**count()**

Return number of occurrences of value.

**index()**

Return first index of value.

Raises ValueError if the value is not present.

**phase**

**class** pylablib.devices.Sirah.Matisse.**TPiezoetFeedforwardParameters** (*ampl, phase*)

Bases: `tuple`

**ampl**

**count()**

Return number of occurrences of value.

**index()**  
Return first index of value.  
Raises ValueError if the value is not present.

**phase**

**class** pylablib.devices.Sirah.Matisse.**TSlowpiezoCtlParameters** (*setpoint, P, I, freeP*)

Bases: `tuple`

**I**

**P**

**count()**  
Return number of occurrences of value.

**freeP**

**index()**  
Return first index of value.  
Raises ValueError if the value is not present.

**setpoint**

**class** pylablib.devices.Sirah.Matisse.**TFastpiezoCtlParameters** (*setpoint, I, lockpoint*)

Bases: `tuple`

**I**

**count()**  
Return number of occurrences of value.

**index()**  
Return first index of value.  
Raises ValueError if the value is not present.

**lockpoint**

**setpoint**

**class** pylablib.devices.Sirah.Matisse.**TRefcellWaveformParameters** (*lower\_limit, upper\_limit, oversamp, mode*)

Bases: `tuple`

**count()**  
Return number of occurrences of value.

**index()**  
Return first index of value.  
Raises ValueError if the value is not present.

**lower\_limit**

**mode**

**oversamp**

**upper\_limit**

```
class pylablib.devices.Sirah.Matisse.TScanMode (falling, stop_lower, stop_upper)
    Bases: tuple

    count ()
        Return number of occurrences of value.

    falling

    index ()
        Return first index of value.

        Raises ValueError if the value is not present.

    stop_lower

    stop_upper

class pylablib.devices.Sirah.Matisse.TScanParameters (device, mode, lower_limit,
                                                    upper_limit, rise_speed,
                                                    fall_speed)

    Bases: tuple

    count ()
        Return number of occurrences of value.

    device

    fall_speed

    index ()
        Return first index of value.

        Raises ValueError if the value is not present.

    lower_limit

    mode

    rise_speed

    upper_limit

class pylablib.devices.Sirah.Matisse.SirahMatisse (addr)
    Bases: pylablib.core.devio.SCPISCPIDevice

    Sirah Matisse laser control.

        Parameters addr – device address (usually a VISA name).

    Error
        alias of GenericSirahError

    ReraiseError
        alias of GenericSirahBackendError

    get_diode_power ()
        Get the current laser resonator power

    get_diode_power_waveform ()
        Get the current laser resonator power waveform

    get_diode_power_lowlevel ()
        Get the low-level cutoff current laser resonator power

    set_diode_power_lowlevel (cutoff)
        Set the low-level cutoff current laser resonator power
```

**get\_thinet\_power()**  
Get the current thin etalon reflex power

**get\_refcell\_waveform()**  
Get the reference cell signal waveform

**bifi\_get\_position()**  
Get the current position of the birefringent filter motor

**bifi\_get\_range()**  
Get the maximum position of the birefringent filter motor

**bifi\_get\_status\_n()**  
Get the numerical status of the birefringent filter motor

**bifi\_get\_status()**  
Get the parsed status of the birefringent filter motor.  
  
Return tuple (code, bits) with, correspondingly, the general status/error code (e.g., "idle", "moving\_abs", or "position\_out\_of\_range"), and a set of active status bits (e.g., "moving", "error", or "limit\_sw1").

**bifi\_clear\_errors()**  
Clear the indicated errors of the birefringent filter motor

**bifi\_is\_moving()**  
Check if the birefringent filter is moving

**bifi\_wait\_move** (timeout=30.0)  
Wait until birefringent filter is done moving

**bifi\_move\_to** (position, wait=True, wait\_timeout=30.0)  
Move the birefringent filter to the current position

**bifi\_stop()**  
Stop the birefringent filter motor

**bifi\_home** (wait=True, wait\_timeout=30.0)  
Home the birefringent filter motor

**thinet\_get\_position()**  
Get the current position of the thin etalon motor

**thinet\_get\_range()**  
Get the maximum position of the thin etalon motor

**thinet\_get\_status\_n()**  
Get the numerical status of the thin etalon motor

**thinet\_get\_status()**  
Get the parsed status of the thin etalon motor.  
  
Return tuple (code, bits) with, correspondingly, the general status/error code (e.g., "idle", "moving\_abs", or "position\_out\_of\_range"), and a set of active status bits (e.g., "moving", "error", or "limit\_sw1").

**thinet\_clear\_errors()**  
Clear the indicated errors of the thin etalon motor

**thinet\_is\_moving()**  
Check if the thin etalon is moving

**thinet\_wait\_move** (timeout=30.0)  
Wait until thin etalon is done moving

**thinet\_move\_to** (*position*, *wait=True*, *wait\_timeout=30.0*)  
Move the thin etalon to the current position

**thinet\_stop** ()  
Stop the thin etalon motor

**thinet\_home** (*wait=True*, *wait\_timeout=30.0*)  
Home the thin etalon motor

**get\_thinet\_ctl\_status** ()  
Get thin etalon lock status ("run" or "stop")

**set\_thinet\_ctl\_status** (*status='run'*)  
Set thin etalon lock status ("run" or "stop")

**get\_thinet\_error\_signal** ()  
Get error signal of the thin etalon lock (emulated when not available on older firmware)

**get\_thinet\_ctl\_params** ()  
Get thin etalon lock control parameters.  
Return tuple (setpoint, P, I, avg).

**set\_thinet\_ctl\_params** (*setpoint=None*, *P=None*, *I=None*, *avg=None*)  
Set thin etalon lock control parameters.  
Any parameters which are None remain unchanged.

**get\_piezoet\_ctl\_status** ()  
Get piezo etalon lock status ("run" or "stop")

**set\_piezoet\_ctl\_status** (*status='run'*)  
Set piezo etalon lock status ("run" or "stop")

**get\_piezoet\_position** ()  
Get piezo etalon DC position

**set\_piezoet\_position** (*value*)  
Set piezo etalon lock DC position

**get\_piezoet\_drive\_params** ()  
Get piezo etalon drive parameters.  
Return tuple (amplitude, rate, oversamp).

**set\_piezoet\_drive\_params** (*amplitude=None*, *rate=None*, *oversamp=None*)  
Set piezo etalon drive parameters.  
*oversamp* should be between 8 and 32. *rate* can take values "8k", "32k", "48k", or "96k". Any parameters which are None remain unchanged.

**get\_piezoet\_feedback\_params** ()  
Get piezo etalon feedback parameters.  
Return tuple (P, avg, phase) (phase is integer between 0 and oversampling).

**set\_piezoet\_feedback\_params** (*P=None*, *avg=None*, *phase=None*)  
Set piezo etalon feedback parameters.  
Phase is integer between 0 and oversampling. Any parameters which are None remain unchanged.

**get\_piezoet\_feedforward\_params** ()  
Get piezo etalon feedforward parameters.  
Return tuple (amp, phase) (phase is integer between 0 and oversampling).

**set\_piezoet\_feedforward\_params** (*amp=None, phase=None*)

Set piezo etalon feedforward parameters.

Phase is integer between 0 and oversampling. Any parameters which are *None* remain unchanged.

**get\_slowpiezo\_ctl\_status** ()

Get slow piezo lock status ("run" or "stop")

**set\_slowpiezo\_ctl\_status** (*status='run'*)

Set slow piezo lock status ("run" or "stop")

**get\_slowpiezo\_position** ()

Get slow piezo DC position

**set\_slowpiezo\_position** (*value*)

Set slow piezo DC position

**get\_slowpiezo\_ctl\_params** ()

Get slow piezo lock control parameters.

Return tuple (setpoint, P, I, freeP).

**set\_slowpiezo\_ctl\_params** (*setpoint=None, P=None, I=None, freeP=None*)

Set slow piezo lock control parameters.

Any parameters which are *None* remain unchanged.

**get\_fastpiezo\_ctl\_status** ()

Get fast piezo lock status ("run" or "stop")

**set\_fastpiezo\_ctl\_status** (*status='run'*)

Set fast piezo lock status ("run" or "stop")

**is\_fastpiezo\_locked** ()

Check if the fast piezo is locked (output is between 5% and 95%)

**get\_fastpiezo\_position** ()

Get fast piezo DC position between 0 and 1

**set\_fastpiezo\_position** (*value*)

Set fast piezo DC position between 0 and 1

**get\_fastpiezo\_ctl\_params** ()

Get fast piezo lock control parameters.

Return tuple (setpoint, I, lockpoint).

**set\_fastpiezo\_ctl\_params** (*setpoint=None, I=None, lockpoint=None*)

Set fast piezo lock control parameters.

Any parameters which are *None* remain unchanged.

**get\_refcell\_position** ()

Get reference cell DC position between 0 and 1

**set\_refcell\_position** (*value*)

Set reference cell DC position between 0 and 1

**get\_refcell\_waveform\_params** ()

Get reference cell waveform parameters.

Return tuple (lower\_limit, upper\_limit, oversamp, mode). mode can be "none", "avg", "min", or "max".

**set\_refcell\_waveform\_params** (*lower\_limit=None, upper\_limit=None, oversamp=None, mode=None*)

Set reference cell waveform parameters.

Any parameters which are None remain unchanged. *mode* can be "none", "avg", "min", or "max". *oversamp* should be between 4 and 512.

**get\_scan\_status** ()

Get scan status ("run" or "stop")

**set\_scan\_status** (*status='run'*)

Set scan status ("run" or "stop")

**wait\_scan** (*timeout=None*)

Wait until scan is stopped

**get\_scan\_position** ()

Get scan position

**set\_scan\_position** (*value*)

Set scan position

**get\_scan\_params** ()

Get scan parameters.

Return tuple (device, mode, lower\_limit, upper\_limit, rise\_speed, fall\_speed). device can be "none", "slow\_piezo", or "ref\_cell". mode is a tuple (falling, stop\_lower, stop\_upper).

**set\_scan\_params** (*device=None, mode=None, lower\_limit=None, upper\_limit=None, rise\_speed=None, fall\_speed=None*)

Set slow piezo lock control parameters.

*device* can be "none", "slow\_piezo", or "ref\_cell". *mode* is a tuple (falling, stop\_lower, stop\_upper). Any parameters which are None remain unchanged.

**BackendError**

alias of `pylablib.core.devio.comm_backend.DeviceBackendError`

**class NoParameterCaller** (*device, kind*)

Bases: `object`

Class to simplify calling functions without a parameter

**apply\_settings** (*settings*)

Apply the settings.

*settings* is a dict {name: value} of the available device settings. Non-applicable settings are ignored.

**ask** (*msg, data\_type='string', delay=0.0, timeout=None, read\_echo=False*)

Write a message and read a reply.

*msg* is the query message, *delay* is the delay between write and read. Other parameters are the same as in `read()`. If *read\_echo*==True, assume that the device first echoes the input and skip it.

**close** ()

Close the backend

**flush** (*one\_line=False*)

Flush the read buffer (read all the available data and return the number of bytes read).

If *one\_line*==True, read only a single line.

**static get\_arg\_type** (*arg*)

Autodetect argument type



**get\_device\_variable** (*key*)

Get the value of a settings, status, or full info parameter

**get\_esr** (*timeout=None*)

Get the device status register (by default, "\*ESR?" command)

**get\_full\_info** (*include=0*)

Get dict {name: value} containing full device information (including status and settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_full\_status** (*include=0*)

Get dict {name: value} containing the device status (including settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_id** (*timeout=None*)

Get the device IDN. (query SCPI '\*IDN?' command)

**get\_settings** (*include=0*)

Get dict {name: value} containing all the device settings.

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**is\_opened** ()

Check if the device is connected

**lock** (*timeout=None*)

Lock the access to the device from other threads/processes (isn't necessarily implemented)

**locking** (*timeout=None*)

Context manager for lock & unlock

**open** ()

Open the backend

**static parse\_array\_data** (*data, fmt, include\_header=False*)

Parse the data returned by the device. *fmt* is [DataFormat](#) description in numpy format (e.g., "<u2").

If *include\_header==True*, the data is assumed to be in a (somewhat) standard SCPI format: b'#', then a single digit *s* denoting length of the size block, then *s* digits denoting length of the data (in bytes) followed by the actual data. Otherwise (*include\_header==False*), assume that the header is already removed.

**read** (*data\_type='string', timeout=None*)

Read data from the device.

*data\_type* determines the type of the data. Can be 'raw' (just raw data), 'string' (with trailing and leading spaces stripped), 'int', 'float', 'bool' (interprets 0 or 'off' as False, anything else as True), 'value' (returns tuple (value, unit), where *value* is float), a callable (return the result of this callable applied to the string value), a dictionary (return the stored value corresponding to the string value, or to the value converted into integer if the string value is not present), or a list of data types

(the result is treated as a list of values with the given types separated by commas). *timeout* overrides the default value.

**read\_binary\_array\_data** (*include\_header=False, timeout=None, flush\_term=True*)

Read a binary data in the from the device.

The data assumes the standard binary transfer header consisting of "#" symbol, then a single digit with the size of the length string, then the length string containing the length of the binary data (in bytes). If *include\_header==True*, return the data with the header; otherwise, return only the content. If *flush\_term==True*, flush the following line to skip terminator characters after the binary data, which are added by some devices. *timeout* overrides the default value.

**reconnect** (*new\_instrument=True, ignore\_error=True*)

Remake the connection.

If *new\_instrument==True*, create a new backend instance. If *ignore\_error==True*, ignore errors on closing.

**reset** ()

Reset the device (by default, "\*RST" command)

**set\_device\_variable** (*key, value*)

Set the value of a settings parameter

**sleep** (*delay*)

Wait for *delay* seconds

**unlock** ()

Unlock the access to the device from other threads/processes (isn't necessarily implemented)

**using\_write\_buffer** ()

Context manager for using a write buffer.

While it's active, all the consecutive *write()* operations are bundled together with ; delimiter. The actual write is performed at the *read()/ask()* operation or at the end of the block.

**wait** (*wait\_type='sync', timeout=None, wait\_callback=None*)

Pause execution until device overlapped commands are complete.

*wait\_type* is either 'sync' (perform *wait\_sync()*), 'dev' (perform *wait\_dev()*) or 'none' (do nothing).

**wait\_dev** ()

Pause execution of the device commands until device overlapped commands (e.g., taking sweeps) are complete.

Note that the code execution is not paused.

**wait\_sync** (*timeout=None, wait\_callback=None*)

Pause execution of the script until device overlapped commands (e.g., taking sweeps) are complete.

*timeout* and *wait\_callback* override default constructor parameters.

**write** (*msg, arg=None, arg\_type=None, unit=None, bool\_selector=None, wait\_sync=None, read\_echo=False, read\_echo\_delay=0.0*)

Send a command.

#### Parameters

- **msg** (*str*) – Text message.
- **arg** – Optional argument to append in the end. If a list of arguments is supplied, the result is joined with ", ".

- **arg\_type** (*str*) – Argument type. Can be 'raw' (in which case data is sent raw), 'string', 'int', 'float', 'bool', a format string (such as '{:.3f}') or a list of argument types (for an iterable argument); if format string is used and the argument is a list or a tuple, then it is expanded as a list of arguments (e.g., `arg_type='{0};{1}'` with `arg=[1,2]` will produce a string '1;2'); if a list of types is used, each element of *arg* is converted using the corresponding type, and the result is joined with ", ".
- **unit** (*str*) – If not None, use it as a unit to append after the value.
- **bool\_selector** (*tuple*) – A tuple (*false\_value*, *true\_value*) of two strings to represent bool argument; by default, use `._bool_selector` attribute.
- **wait\_sync** – if True, append the sync command (specified as `._wait_sync_comm` attribute, `"*OPC?"` by default) after the message and pause the execution command is complete; useful in long set operations, where the device might ignore later inputs until the current command is complete; if None, use the class default `._default_write_sync` attribute (False by default).
- **read\_echo** (*bool*) – If True, read a single line after write.
- **read\_echo\_delay** (*float*) – The delay between write and read if `read_echo==True`.

## pylablib.devices.Sirah.tuner module

**class** pylablib.devices.Sirah.tuner.**MatisseTuner** (*laser*, *wavemeter*)

Bases: `object`

Matisse tuner.

Helps to coordinate with an external wavemeter to perform more complicated tasks: motors calibration, fine frequency tuning, and stitching scans.

### Parameters

- **laser** – opened Matisse laser object
- **wavemeter** – opened wavemeter object (currently only HighFinesse wavemeters are supported)

**get\_frequency** (*timeout=1.0*)

Get current frequency reading.

The only method relying on the wavemeter. Can be extended or overloaded to support different wavemeters.

**scan\_steps** (*motor*, *start*, *stop*, *step*)

Scan the given motor ("bifi" or "thinet") in discrete steps within the given range with a given step.

Return a 4-column numpy array containing motor position, internal diode power, thin etalon reflection power, and wavemeter frequency.

**scan\_centered** (*motor*, *span*, *step*)

Scan the given motor ("bifi" or "thinet") in discrete steps in a given span around the current position.

After the scan, return the motor to the original position.

Return a 4-column numpy array containing motor position, internal diode power, thin etalon reflection power, and wavemeter frequency.

**scan\_quick** (*motor, start, stop, autodir=True*)

Do a quick continuous scan of the given motor ("bifi" or "thinet") within the given range.

Compared to `scan_steps()`, which does a series of discrete defined moves, this method does a single continuous move and records values in its progress. This is quicker, but does not allow for the step size control, and results in non-uniform recorded positions. If `autodir==False`, first initialize the motor to *start* and then move to *stop*; otherwise, initialize to whichever border is closer.

Return a 4-column numpy array containing motor position, internal diode power, thin etalon reflection power, and wavemeter frequency.

**scan\_quick\_centered** (*motor, span*)

Do a quick continuous scan of the given motor ("bifi" or "thinet") in a given span around the current position.

After the scan, return the motor to the original position.

Return a 4-column numpy array containing motor position, internal diode power, thin etalon reflection power, and wavemeter frequency.

**scan\_both\_motors** (*bifi\_rng, te\_rng, verbose=False*)

Perform a 2D grid scan changing positions of both birefringent filter and thin etalon motors.

*bifi\_rng* and *te\_rng* are both 3-tuples (*start, stop, step*) specifying the scan ranges. If `verbose==True`, print a message per every birefringent filter position indicating the scan progress.

Return a 5-column numpy array containing birefringent filter motor position, thin etalon motor position, internal diode power, thin etalon reflection power, and wavemeter frequency.

**scan\_both\_motors\_quick** (*bifi\_rng, te\_rng, verbose=False*)

Perform a quick 2D grid scan changing positions of both birefringent filter and thin etalon motors.

For each discrete position of a birefringent filter motor perform a quick scan of the thin etalon motor. *bifi\_rng* is a 3-tuple (*start, stop, step*), while *te\_rng* is a 2-tuple (*start, stop*) specifying the scan ranges. If `verbose==True`, print a message per every birefringent filter position indicating the scan progress.

Return a 5-column numpy array containing birefringent filter motor position, thin etalon motor position, internal diode power, thin etalon reflection power, and wavemeter frequency.

**unlock\_all** ()

Unlock all relevant locks (slow piezo, fast piezo, piezo etalon, thin etalon)

**slow\_piezo\_tune\_to** (*target*)

Fine tune the laser to the given target frequency using only slow piezo tuning

**tune\_to** (*target*)

Tune the laser to the given frequency using all elements (bifi, thin etalon, piezo etalon, slow piezo)

**stitched\_scan** (*full\_rng, single\_span, speed, overlap=0.1, freq\_step=None*)

Perform a stitched laser scan.

#### Parameters

- **full\_rng** – 2-tuple (*start, stop*) with the full frequency scan range.
- **single\_span** – magnitude of a single continuous scan segment given in the slow piezo scan units (between 0 and 1)
- **speed** – single segment scan speed
- **overlap** – overlap of consecutive segments, as a fraction of *single\_span*

- **freq\_step** – if `None`, the start of the next segment is calculated based on the end of the previous segment and *overlap*; otherwise, it specifies a fixed frequency step between segments.

## Module contents

### pylablib.devices.SmarAct package

#### Submodules

#### pylablib.devices.SmarAct.scu3d module

```

class pylablib.devices.SmarAct.scu3d.LibraryController(lib)
    Bases: pylablib.devices.utils.load_lib.LibraryController

    close(opid)
        Mark device closing.

        Return tuple (close_result, uninit_result) with the results of the closing and the shutdown.
        If library does not need to be shut down yet, set uninit_result=None

    open()
        Mark device opening.

        Return tuple (init_result, open_result, opid) with the results of the initialization and the
        opening, and the opening ID which should afterwards be used for closing. If library is already initialized,
        set init_result=None

    preinit()
        Pre-initialize the library, if it hasn't been done already

    shutdown()
        Close all opened connections and shutdown the library

    temp_open()
        Context for temporarily opening a new device connection

class pylablib.devices.SmarAct.scu3d.TDeviceInfo(device_id,          firmware_version,
                                                  dll_version)
    Bases: tuple

    count()
        Return number of occurrences of value.

    device_id

    dll_version

    firmware_version

    index()
        Return first index of value.

        Raises ValueError if the value is not present.

pylablib.devices.SmarAct.scu3d.get_device_info(idx)
    Get info of the devices with the given index.

    Return tuple (device_id, firmware_version, dll_version).

```

```
pylablib.devices.SmarAct.scu3d.list_devices()
```

List all connected devices

```
pylablib.devices.SmarAct.scu3d.get_devices_number()
```

Get number of connected SCU3D controller

```
class pylablib.devices.SmarAct.scu3d.SCU3D (idx=0, axis_dir='+++')
```

Bases: `pylablib.devices.interface.stage.IMultiaxisStage`

SmarAct SCU3D translation stage controller.

#### Parameters

- **idx** (*int*) – stage index
- **axis\_dir** (*str*) – 3-symbol string specifying default directions of the axes (each symbol be "+" or "-")

```
Error = <Mock name='mock.SmarActError' id='140096265476688'>
```

```
open()
```

Open the connection to the stage

```
close()
```

Close the connection to the stage

```
is_opened()
```

Check if the device is connected

```
get_device_info()
```

Get info of the devices with the given index.

Return tuple (device\_id, firmware\_version, dll\_version).

```
get_axis_dir()
```

Get axis direction convention (a string of 3 symbols which are either "+" or "-" determining if the axis direction is flipped)

```
set_axis_dir (axis_dir)
```

Set axis direction convention (a string of 3 symbols which are either "+" or "-" determining if the axis direction is flipped)

```
move_macrostep (axis, steps, voltage, frequency)
```

Move along a given axis by a single “macrostep”, which consists of several regular steps.

*voltage* (in Volts) and *frequency* (in Hz) specify the motion parameters. This simulates the controller operation, where one “step” at large step sizes consists of several small steps.

```
move_by (axis, steps=1, stepsize=10)
```

Move along a given axis with a given number of macrosteps using one of the predefined step size.

*stepsize* can range from 1 (smallest) to 20 (largest), and roughly corresponds to the handheld controller parameters.

```
get_status (axis='all')
```

Get the axis status.

Can be "stopped" (default state), "setting\_amplitude" (setting open-loop step amplitude), "moving" (open-loop movement), "targeting" (closed-loop movement), "holding" (closed-loop position holding), "calibrating" (sensor calibration), or "moving\_to\_reference" (calibrating position sensor).

```
wait_for_status (axis, status='stopped', timeout=30.0)
```

Wait until the axis reaches a given status.

By default wait for "stopped" status (i.e., wait until the motion is finished).

**wait\_move** (*axis*, *timeout=30.0*)

Wait for a given axis to stop moving

**is\_moving** (*axis='all'*)

Check if a given axis is moving

**stop** (*axis='all'*)

Stop motion at a given axis

**class NoParameterCaller** (*device*, *kind*)

Bases: `object`

Class to simplify calling functions without a parameter

**apply\_settings** (*settings*)

Apply the settings.

*settings* is the dict {*name*: *value*} of the device available settings. Non-applicable settings are ignored.

**get\_all\_axes** ()

Get the list of all available axes (taking mapping into account)

**get\_device\_variable** (*key*)

Get the value of a settings, status, or full info parameter

**get\_full\_info** (*include=0*)

Get dict {*name*: *value*} containing full device information (including status and settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_full\_status** (*include=0*)

Get dict {*name*: *value*} containing the device status (including settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_settings** (*include=0*)

Get dict {*name*: *value*} containing all the device settings.

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**remap\_axes** (*mapping*, *accept\_original=True*)

Rename axes to the new labels.

*mapping* is the new axes mapping, which can be a list of new axes name (corresponding to the old axes in order returned by `get_all_axes()`), or a dictionary {*alias*: *original*} of the new axes aliases.

**set\_device\_variable** (*key*, *value*)

Set the value of a settings parameter

## Module contents

### pylablib.devices.Tektronix package

#### Submodules

#### pylablib.devices.Tektronix.base module

**exception** pylablib.devices.Tektronix.base.**TektronixError**

Bases: *pylablib.core.devio.base.DeviceError*

Generic Tektronix devices error

**args**

**with\_traceback()**

Exception.with\_traceback(tb) – set self.\_\_traceback\_\_ to tb and return self.

**exception** pylablib.devices.Tektronix.base.**TektronixBackendError** (*exc*)

Bases: *pylablib.devices.Tektronix.base.TektronixError, pylablib.core.devio.comm\_backend.DeviceBackendError*

Generic Tektronix backend communication error

**args**

**with\_traceback()**

Exception.with\_traceback(tb) – set self.\_\_traceback\_\_ to tb and return self.

*pylablib.devices.Tektronix.base.muxchannel* (\*args, \*\*kwargs)

Multiplex the function over its channel argument

**class** pylablib.devices.Tektronix.base.**TTriggerParameters** (*source, level, coupling, slope*)

Bases: *tuple*

**count()**

Return number of occurrences of value.

**coupling**

**index()**

Return first index of value.

Raises ValueError if the value is not present.

**level**

**slope**

**source**

**class** pylablib.devices.Tektronix.base.**ITektronixScope** (*addr, nchannels='auto'*)

Bases: *pylablib.core.devio.SCPISCPIDevice*

Generic Tektronix oscilloscope.

**Parameters**

- **addr** – device address; usually a VISA address string such as "USB0::0x0699::0x0364::C000000::INSTR"



- **nchannels** – can specify number of channels on the oscilloscope; by default, autode-  
tect number of channels (might take several seconds on connection)

**Error**

alias of *TektronixError*

**ReraiseError**

alias of *TektronixBackendError*

**get\_channels\_number()**

Get the number of channels

**get\_channels()** (*only\_main=False*)

Get the list of all input channels (if *only\_main==True*) or all available channels (if *only\_main==False*)

**normalize\_channel\_name()** (*channel*)

Normalize channel name as represented by the oscilloscope

**grab\_single()** (*wait=True, software\_trigger=False, wait\_timeout=None*)

Set single waveform grabbing and wait for acquisition.

If *wait==True*, wait until the acquisition is complete; otherwise, return immediately. if *software\_trigger==True*, send the software trigger after setup (i.e., the device triggers immediately regardless of the input).

**wait\_for\_grabbing()** (*timeout=None*)

Wait until the acquisition is complete

**grab\_continuous()** (*enable=True*)

Start or stop continuous grabbing

**stop\_grabbing()**

Stop grabbing or waiting (equivalent to *self.grab\_continuous(False)*)

**is\_continuous()**

Check if grabbing is continuous or single

**is\_grabbing()**

Check if acquisition is in progress.

Return *True* if the oscilloscope is recording data, or if the trigger is armed/ready and waiting; return *False* if the acquisition is stopped. To check if the trigger has been triggered, use *get\_trigger\_state()*.

**get\_edge\_trigger\_source()**

Get edge trigger source.

Can be an integer indicating channel number or a name of a special channel.

**set\_edge\_trigger\_source()** (*channel*)

Get edge trigger source.

Can be an integer indicating channel number or a name of a special channel.

**get\_edge\_trigger\_coupling()**

Get edge trigger coupling ("ac" or "dc")

**set\_edge\_trigger\_coupling()** (*coupling*)

Set edge trigger coupling ("ac" or "dc")

**get\_edge\_trigger\_slope()**

Get edge trigger slope ("fall" or "rise")

**set\_edge\_trigger\_slope** (*slope*)  
Set edge trigger slope ("fall" or "rise")

**get\_trigger\_level** ()  
Get edge trigger level (in Volts)

**set\_trigger\_level** (*level*)  
Set edge trigger level (in Volts)

**setup\_edge\_trigger** (*source, level, coupling='dc', slope='rise'*)  
Setup edge trigger.  
Set source, level, coupling and slope (see corresponding methods for details).

**get\_trigger\_mode** ()  
Get trigger mode.  
Can be either "auto" or "norm".

**set\_trigger\_mode** (*trigger\_mode='auto'*)  
Set trigger mode.  
Can be either "auto" or "norm".

**get\_trigger\_state** ()  
Get trigger state.  
Can be "armed" (acquiring pretrigger), "ready" (pretrigger acquired, wait for trigger event), "trigger" (triggered, acquiring the rest of the waveform), "auto" ("auto" mode trigger is acquiring data in the absence of trigger), "save" (acquisition is stopped), or "scan" (oscilloscope in the scan mode)

**force\_trigger** ()  
Force trigger event

**get\_horizontal\_span** ()  
Get horizontal span (in seconds)

**set\_horizontal\_span** (*span*)  
Set horizontal span (in seconds)

**get\_horizontal\_offset** ()  
Get horizontal offset (position of the center of the sweep; in seconds)

**set\_horizontal\_offset** (*offset=0.0*)  
Set horizontal offset (position of the center of the sweep; in seconds)

**get\_vertical\_span** (*channel*)  
Get channel vertical span (in V)

**set\_vertical\_span** (*channel, span*)  
Set channel vertical span (in V)

**get\_vertical\_position** (*channel*)  
Get channel vertical position (offset of the zero volt line; in V)

**set\_vertical\_position** (*channel, position*)  
Set channel vertical position (offset of the zero volt line; in V)

**is\_channel\_enabled** (*channel*)  
Check if channel is enabled

**enable\_channel** (*channel, enabled=True*)  
Enable or disable given channel

**get\_selected\_channel()**

Get selected source channel.

Return number if it is a real channel, or a string name otherwise.

**select\_channel(channel)**

Select a channel to read data.

Doesn't need to be called explicitly, if `read_multiple_sweeps()` or `read_sweep()` are used.

**get\_coupling(channel)**

Get channel coupling.

Can be "ac", "dc", or "gnd".

**set\_coupling(channel, coupling='dc')**

Set channel coupling.

Can be "ac", "dc", or "gnd".

**get\_probe\_attenuation(channel)**

Get channel probe attenuation

**set\_probe\_attenuation(channel, attenuation)**

Set channel probe attenuation

**get\_points\_number(kind='send')**

Get number of datapoints in various context.

`kind` defines the context. It can be "acq" (number of points acquired), "trace" (number of points in the source of the read-out trace; can be lower than "acq" if the data resolution is reduced, or if the source is not a channel data), or "send" (number of points in the sent waveform; can be lower than "trace" if `get_data_pts_range()` is used to specify and incomplete range). Not all kinds are defined for all scope model (e.g., "trace" is not defined for TDS2000 series oscilloscopes).

For length of read-out trace, see also `get_data_pts_range()`.

**set\_points\_number(pts\_num, reset\_limits=True)**

Set number of datapoints to record when acquiring a trace.

If `reset_limits==True`, reset the datapoints range (`set_data_pts_range()`) to the full range. The actual set value (returned by this method) can be different from the requested value.

**get\_data\_pts\_range()**

Get range of data points to read.

The range is defined from 1 to the points number (returned by `get_points_number()`).

**set\_data\_pts\_range(rng=None)**

Set range of data points to read.

The range is defined from 1 to the points number (returned by `get_points_number()` with `kind="acq"`). If `rng` is `None`, set the full range.

**set\_data\_format(fmt='default')**

Set data transfer format.

`fmt` is a string describing the format; can be either "ascii", or a numpy-style format string (e.g., "<u2"). If "default", use the oscilloscope default format (usually binary with smallest appropriate byte size).

**get\_data\_format()**

Get data transfer format.

Return a string describing the format; can be either "ascii", or a numpy-style format string (e.g., "<u2").

**get\_wfmpre** (*channel=None, enable=True*)

Get preamble dictionary describing all scaling and format data for the given channel or a list of channels.

Can be acquired once and used in subsequent multiple reads to save time on re-requesting. If *channel* is *None*, use the currently selected channel. If *enable==True*, make sure that the requested channel is enabled; getting preamble for disabled channels raises an error.

**read\_raw\_data** (*channel=None, fmt=None, timeout=None*)

Request, read and parse raw data at a given channel.

*fmt* is data format (e.g., "i1", "<i2", or "ascii") or "default", which uses the default oscilloscope format (usually binary with smallest appropriate byte size). If *fmt* is *None*, use the current format. If *channel* is *None*, use the currently selected channel.

Returned data is raw (i.e., not scaled and without x axis).

**read\_multiple\_sweeps** (*channels, wfmpres=None, ensure\_fmt=False, timeout=None, return\_wfmpres=None*)

Read data from a multiple channels channel.

#### Parameters

- **channels** – list of channel indices or names
- **wfmpres** – optional list or dictionary of preambles (obtained using `get_wfmpre()`); if it is *None*, obtain during reading, which slows down the data acquisition a bit
- **ensure\_fmt** – if *True*, make sure that oscilloscope data format agrees with the one in *wfmpre*
- **timeout** – read timeout
- **return\_wfmpres** – if *True*, return tuple (*sweeps, wfmpres*), where *wfmpres* can be used for further sweep readouts.

**read\_sweep** (*channel, wfmpre=None, ensure\_fmt=True, timeout=None*)

Read data from a single channel.

#### Parameters

- **channel** – channel index or name
- **wfmpre** – optional preamble dictionary (obtained using `get_wfmpre()`); if it is *None*, obtain during reading, which slows down the data acquisition a bit
- **ensure\_fmt** – if *True*, make sure that oscilloscope data format agrees with the one in *wfmpre*
- **timeout** – read timeout

#### BackendError

alias of `pylablib.core.devio.comm_backend.DeviceBackendError`

**class NoParameterCaller** (*device, kind*)

Bases: `object`

Class to simplify calling functions without a parameter

**apply\_settings** (*settings*)

Apply the settings.

*settings* is a dict {*name*: *value*} of the available device settings. Non-applicable settings are ignored.

**ask** (*msg*, *data\_type*='string', *delay*=0.0, *timeout*=None, *read\_echo*=False)

Write a message and read a reply.

*msg* is the query message, *delay* is the delay between write and read. Other parameters are the same as in `read()`. If `read_echo==True`, assume that the device first echoes the input and skip it.

**close** ()

Close the backend

**flush** (*one\_line*=False)

Flush the read buffer (read all the available data and return the number of bytes read).

If `one_line==True`, read only a single line.

**static get\_arg\_type** (*arg*)

Autodetect argument type

**get\_device\_variable** (*key*)

Get the value of a settings, status, or full info parameter

**get\_esr** (*timeout*=None)

Get the device status register (by default, "\*ESR?" command)

**get\_full\_info** (*include*=0)

Get dict {name: value} containing full device information (including status and settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting `include=-10` queries all available variables, which is equivalent to `include="all"`.

**get\_full\_status** (*include*=0)

Get dict {name: value} containing the device status (including settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting `include=-10` queries all available variables, which is equivalent to `include="all"`.

**get\_id** (*timeout*=None)

Get the device IDN. (query SCPI '\*IDN?' command)

**get\_settings** (*include*=0)

Get dict {name: value} containing all the device settings.

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting `include=-10` queries all available variables, which is equivalent to `include="all"`.

**is\_opened** ()

Check if the device is connected

**lock** (*timeout*=None)

Lock the access to the device from other threads/processes (isn't necessarily implemented)

**locking** (*timeout*=None)

Context manager for lock & unlock

**open** ()

Open the backend

**static parse\_array\_data** (*data*, *fmt*, *include\_header=False*)

Parse the data returned by the device. *fmt* is *DataFormat* description in numpy format (e.g., "<u2").

If *include\_header==True*, the data is assumed to be in a (somewhat) standard SCPI format: *b' # '*, then a single digit *s* denoting length of the size block, then *s* digits denoting length of the data (in bytes) followed by the actual data. Otherwise (*include\_header==False*), assume that the header is already removed.

**read** (*data\_type='string'*, *timeout=None*)

Read data from the device.

*data\_type* determines the type of the data. Can be 'raw' (just raw data), 'string' (with trailing and leading spaces stripped), 'int', 'float', 'bool' (interprets 0 or 'off' as False, anything else as True), 'value' (returns tuple (value, unit), where *value* is float), a callable (return the result of this callable applied to the string value), a dictionary (return the stored value corresponding to the string value, or to the value converted into integer if the string value is not present), or a list of data types (the result is treated as a list of values with the given types separated by commas). *timeout* overrides the default value.

**read\_binary\_array\_data** (*include\_header=False*, *timeout=None*, *flush\_term=True*)

Read a binary data in the from the device.

The data assumes the standard binary transfer header consisting of "# " symbol, then a single digit with the size of the length string, then the length string containing the length of the binary data (in bytes). If *include\_header==True*, return the data with the header; otherwise, return only the content. If *flush\_term==True*, flush the following line to skip terminator characters after the binary data, which are added by some devices. *timeout* overrides the default value.

**reconnect** (*new\_instrument=True*, *ignore\_error=True*)

Remake the connection.

If *new\_instrument==True*, create a new backend instance. If *ignore\_error==True*, ignore errors on closing.

**reset** ()

Reset the device (by default, "\*RST" command)

**set\_device\_variable** (*key*, *value*)

Set the value of a settings parameter

**sleep** (*delay*)

Wait for *delay* seconds

**unlock** ()

Unlock the access to the device from other threads/processes (isn't necessarily implemented)

**using\_write\_buffer** ()

Context manager for using a write buffer.

While it's active, all the consecutive *write()* operations are bundled together with ; delimiter. The actual write is performed at the *read()/ask()* operation or at the end of the block.

**wait** (*wait\_type='sync'*, *timeout=None*, *wait\_callback=None*)

Pause execution until device overlapped commands are complete.

*wait\_type* is either 'sync' (perform *wait\_sync()*), 'dev' (perform *wait\_dev()*) or 'none' (do nothing).

**wait\_dev** ()

Pause execution of the device commands until device overlapped commands (e.g., taking sweeps) are complete.

Note that the code execution is not paused.

**wait\_sync** (*timeout=None, wait\_callback=None*)

Pause execution of the script until device overlapped commands (e.g., taking sweeps) are complete.

*timeout* and *wait\_callback* override default constructor parameters.

**write** (*msg, arg=None, arg\_type=None, unit=None, bool\_selector=None, wait\_sync=None, read\_echo=False, read\_echo\_delay=0.0*)

Send a command.

#### Parameters

- **msg** (*str*) – Text message.
- **arg** – Optional argument to append in the end. If a list of arguments is supplied, the result is joined with ", ".
- **arg\_type** (*str*) – Argument type. Can be 'raw' (in which case data is sent raw), 'string', 'int', 'float', 'bool', a format string (such as '{:.3f}') or a list of argument types (for an iterable argument); if format string is used and the argument is a list or a tuple, then it is expanded as a list of arguments (e.g., `arg_type='{0};{1}'` with `arg=[1, 2]` will produce a string '1; 2'); if a list of types is used, each element of *arg* is converted using the corresponding type, and the result is joined with ", ".
- **unit** (*str*) – If not None, use it as a unit to append after the value.
- **bool\_selector** (*tuple*) – A tuple (*false\_value*, *true\_value*) of two strings to represent bool argument; by default, use `._bool_selector` attribute.
- **wait\_sync** – if True, append the sync command (specified as `._wait_sync_comm` attribute, `"*OPC?"` by default) after the message and pause the execution command is complete; useful in long set operations, where the device might ignore later inputs until the current command is complete; if None, use the class default `._default_write_sync` attribute (False by default).
- **read\_echo** (*bool*) – If True, read a single line after write.
- **read\_echo\_delay** (*float*) – The delay between write and read if `read_echo==True`.

**class** `pylablib.devices.Tektronix.base.TDS2000` (*addr, nchannels='auto'*)

Bases: `pylablib.devices.Tektronix.base.ITektronixScope`

Tektronix TDS2000 series oscilloscope.

#### Parameters

- **addr** – device address; usually a VISA address string such as `"USB0::0x0699::0x0364::C000000::INSTR"`
- **nchannels** – can specify number of channels on the oscilloscope; by default, autodetect number of channels (might take several seconds on connection)

#### BackendError

alias of `pylablib.core.devio.comm_backend.DeviceBackendError`

#### Error

alias of `TektronixError`

**class** `NoParameterCaller` (*device, kind*)

Bases: `object`

Class to simplify calling functions without a parameter

**ReraiseError**

alias of *TektronixBackendError*

**apply\_settings** (*settings*)

Apply the settings.

*settings* is a dict {name: value} of the available device settings. Non-applicable settings are ignored.

**ask** (*msg*, *data\_type*='string', *delay*=0.0, *timeout*=None, *read\_echo*=False)

Write a message and read a reply.

*msg* is the query message, *delay* is the delay between write and read. Other parameters are the same as in *read()*. If *read\_echo*==True, assume that the device first echoes the input and skip it.

**close** ()

Close the backend

**enable\_channel** (*channel*, *enabled*=True)

Enable or disable given channel

**flush** (*one\_line*=False)

Flush the read buffer (read all the available data and return the number of bytes read).

If *one\_line*==True, read only a single line.

**force\_trigger** ()

Force trigger event

**static get\_arg\_type** (*arg*)

Autodetect argument type

**get\_channels** (*only\_main*=False)

Get the list of all input channels (if *only\_main*==True) or all available channels (if *only\_main*==False)

**get\_channels\_number** ()

Get the number of channels

**get\_coupling** (*channel*)

Get channel coupling.

Can be "ac", "dc", or "gnd".

**get\_data\_format** ()

Get data transfer format.

Return a string describing the format; can be either "ascii", or a numpy-style format string (e.g., "<u2").

**get\_data\_pts\_range** ()

Get range of data points to read.

The range is defined from 1 to the points number (returned by *get\_points\_number()*).

**get\_device\_variable** (*key*)

Get the value of a settings, status, or full info parameter

**get\_edge\_trigger\_coupling** ()

Get edge trigger coupling ("ac" or "dc")

**get\_edge\_trigger\_slope** ()

Get edge trigger slope ("fall" or "rise")



**get\_edge\_trigger\_source()**

Get edge trigger source.

Can be an integer indicating channel number or a name of a special channel.

**get\_esr** (*timeout=None*)

Get the device status register (by default, "\*ESR?" command)

**get\_full\_info** (*include=0*)

Get dict {name: value} containing full device information (including status and settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_full\_status** (*include=0*)

Get dict {name: value} containing the device status (including settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_horizontal\_offset** ()

Get horizontal offset (position of the center of the sweep; in seconds)

**get\_horizontal\_span** ()

Get horizontal span (in seconds)

**get\_id** (*timeout=None*)

Get the device IDN. (query SCPI '\*IDN?' command)

**get\_points\_number** (*kind='send'*)

Get number of datapoints in various context.

*kind* defines the context. It can be "acq" (number of points acquired), "trace" (number of points in the source of the read-out trace; can be lower than "acq" if the data resolution is reduced, or if the source is not a channel data), or "send" (number of points in the sent waveform; can be lower than "trace" if [get\\_data\\_pts\\_range\(\)](#) is used to specify and incomplete range). Not all kinds are defined for all scope model (e.g., "trace" is not defined for TDS2000 series oscilloscopes).

For length of read-out trace, see also [get\\_data\\_pts\\_range\(\)](#).

**get\_probe\_attenuation** (*channel*)

Get channel probe attenuation

**get\_selected\_channel** ()

Get selected source channel.

Return number if it is a real channel, or a string name otherwise.

**get\_settings** (*include=0*)

Get dict {name: value} containing all the device settings.

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_trigger\_level** ()

Get edge trigger level (in Volts)

**get\_trigger\_mode()**

Get trigger mode.

Can be either "auto" or "norm".

**get\_trigger\_state()**

Get trigger state.

Can be "armed" (acquiring pretrigger), "ready" (pretrigger acquired, wait for trigger event), "trigger" (triggered, acquiring the rest of the waveform), "auto" ("auto" mode trigger is acquiring data in the absence of trigger), "save" (acquisition is stopped), or "scan" (oscilloscope in the scan mode)

**get\_vertical\_position(channel)**

Get channel vertical position (offset of the zero volt line; in V)

**get\_vertical\_span(channel)**

Get channel vertical span (in V)

**get\_wfmpre(channel=None, enable=True)**

Get preamble dictionary describing all scaling and format data for the given channel or a list of channels.

Can be acquired once and used in subsequent multiple reads to save time on re-requesting. If *channel* is *None*, use the currently selected channel. If *enable==True*, make sure that the requested channel is enabled; getting preamble for disabled channels raises an error.

**grab\_continuous(enable=True)**

Start or stop continuous grabbing

**grab\_single(wait=True, software\_trigger=False, wait\_timeout=None)**

Set single waveform grabbing and wait for acquisition.

If *wait==True*, wait until the acquisition is complete; otherwise, return immediately. If *software\_trigger==True*, send the software trigger after setup (i.e., the device triggers immediately regardless of the input).

**is\_channel\_enabled(channel)**

Check if channel is enabled

**is\_continuous()**

Check if grabbing is continuous or single

**is\_grabbing()**

Check if acquisition is in progress.

Return *True* if the oscilloscope is recording data, or if the trigger is armed/ready and waiting; return *False* if the acquisition is stopped. To check if the trigger has been triggered, use `get_trigger_state()`.

**is\_opened()**

Check if the device is connected

**lock(timeout=None)**

Lock the access to the device from other threads/processes (isn't necessarily implemented)

**locking(timeout=None)**

Context manager for lock & unlock

**normalize\_channel\_name(channel)**

Normalize channel name as represented by the oscilloscope

**open()**

Open the backend

**static parse\_array\_data** (*data, fmt, include\_header=False*)

Parse the data returned by the device. *fmt* is *DataFormat* description in numpy format (e.g., "<u2").

If *include\_header==True*, the data is assumed to be in a (somewhat) standard SCPI format: *b' # '*, then a single digit *s* denoting length of the size block, then *s* digits denoting length of the data (in bytes) followed by the actual data. Otherwise (*include\_header==False*), assume that the header is already removed.

**read** (*data\_type='string', timeout=None*)

Read data from the device.

*data\_type* determines the type of the data. Can be 'raw' (just raw data), 'string' (with trailing and leading spaces stripped), 'int', 'float', 'bool' (interprets 0 or 'off' as False, anything else as True), 'value' (returns tuple (value, unit), where *value* is float), a callable (return the result of this callable applied to the string value), a dictionary (return the stored value corresponding to the string value, or to the value converted into integer if the string value is not present), or a list of data types (the result is treated as a list of values with the given types separated by commas). *timeout* overrides the default value.

**read\_binary\_array\_data** (*include\_header=False, timeout=None, flush\_term=True*)

Read a binary data in the from the device.

The data assumes the standard binary transfer header consisting of "# " symbol, then a single digit with the size of the length string, then the length string containing the length of the binary data (in bytes). If *include\_header==True*, return the data with the header; otherwise, return only the content. If *flush\_term==True*, flush the following line to skip terminator characters after the binary data, which are added by some devices. *timeout* overrides the default value.

**read\_multiple\_sweeps** (*channels, wfmpres=None, ensure\_fmt=False, timeout=None, return\_wfmpres=None*)

Read data from a multiple channels channel.

#### Parameters

- **channels** – list of channel indices or names
- **wfmpres** – optional list or dictionary of preambles (obtained using *get\_wfmpre()*); if it is None, obtain during reading, which slows down the data acquisition a bit
- **ensure\_fmt** – if True, make sure that oscilloscope data format agrees with the one in *wfmpre*
- **timeout** – read timeout
- **return\_wfmpres** – if True, return tuple (sweeps, wfmpres), where wfmpres can be used for further sweep readouts.

**read\_raw\_data** (*channel=None, fmt=None, timeout=None*)

Request, read and parse raw data at a given channel.

*fmt* is data format (e.g., "i1", "<i2", or "ascii") or "default", which uses the default oscilloscope format (usually binary with smallest appropriate byte size). If *fmt* is None, use the current format. If *channel* is None, use the currently selected channel.

Returned data is raw (i.e., not scaled and without x axis).

**read\_sweep** (*channel, wfmpre=None, ensure\_fmt=True, timeout=None*)

Read data from a single channel.

#### Parameters

- **channel** – channel index or name

- **wfmpre** – optional preamble dictionary (obtained using `get_wfmpre()`); if it is `None`, obtain during reading, which slows down the data acquisition a bit
- **ensure\_fmt** – if `True`, make sure that oscilloscope data format agrees with the one in `wfmpre`
- **timeout** – read timeout

**reconnect** (*new\_instrument=True, ignore\_error=True*)

Remake the connection.

If `new_instrument==True`, create a new backend instance. If `ignore_error==True`, ignore errors on closing.

**reset** ()

Reset the device (by default, "`*RST`" command)

**select\_channel** (*channel*)

Select a channel to read data.

Doesn't need to be called explicitly, if `read_multiple_sweeps()` or `read_sweep()` are used.

**set\_coupling** (*channel, coupling='dc'*)

Set channel coupling.

Can be "`ac`", "`dc`", or "`gnd`".

**set\_data\_format** (*fmt='default'*)

Set data transfer format.

*fmt* is a string describing the format; can be either "`ascii`", or a numpy-style format string (e.g., "`<u2`"). If "`default`", use the oscilloscope default format (usually binary with smallest appropriate byte size).

**set\_data\_pts\_range** (*rng=None*)

Set range of data points to read.

The range is defined from 1 to the points number (returned by `get_points_number()` with `kind="acq"`). If *rng* is `None`, set the full range.

**set\_device\_variable** (*key, value*)

Set the value of a settings parameter

**set\_edge\_trigger\_coupling** (*coupling*)

Set edge trigger coupling ("`ac`" or "`dc`")

**set\_edge\_trigger\_slope** (*slope*)

Set edge trigger slope ("`fall`" or "`rise`")

**set\_edge\_trigger\_source** (*channel*)

Get edge trigger source.

Can be an integer indicating channel number or a name of a special channel.

**set\_horizontal\_offset** (*offset=0.0*)

Set horizontal offset (position of the center of the sweep; in seconds)

**set\_horizontal\_span** (*span*)

Set horizontal span (in seconds)

**set\_points\_number** (*pts\_num, reset\_limits=True*)

Set number of datapoints to record when acquiring a trace.

If `reset_limits==True`, reset the datapoints range (`set_data_pts_range()`) to the full range. The actual set value (returned by this method) can be different from the requested value.

**set\_probe\_attenuation** (*channel*, *attenuation*)

Set channel probe attenuation

**set\_trigger\_level** (*level*)

Set edge trigger level (in Volts)

**set\_trigger\_mode** (*trigger\_mode*='auto')

Set trigger mode.

Can be either "auto" or "norm".

**set\_vertical\_position** (*channel*, *position*)

Set channel vertical position (offset of the zero volt line; in V)

**set\_vertical\_span** (*channel*, *span*)

Set channel vertical span (in V)

**setup\_edge\_trigger** (*source*, *level*, *coupling*='dc', *slope*='rise')

Setup edge trigger.

Set source, level, coupling and slope (see corresponding methods for details).

**sleep** (*delay*)

Wait for *delay* seconds

**stop\_grabbing** ()

Stop grabbing or waiting (equivalent to `self.grab_continuous(False)`)

**unlock** ()

Unlock the access to the device from other threads/processes (isn't necessarily implemented)

**using\_write\_buffer** ()

Context manager for using a write buffer.

While it's active, all the consecutive `write()` operations are bundled together with ; delimiter. The actual write is performed at the `read()/ask()` operation or at the end of the block.

**wait** (*wait\_type*='sync', *timeout*=None, *wait\_callback*=None)

Pause execution until device overlapped commands are complete.

*wait\_type* is either 'sync' (perform `wait_sync()`), 'dev' (perform `wait_dev()`) or 'none' (do nothing).

**wait\_dev** ()

Pause execution of the device commands until device overlapped commands (e.g., taking sweeps) are complete.

Note that the code execution is not paused.

**wait\_for\_grabbing** (*timeout*=None)

Wait until the acquisition is complete

**wait\_sync** (*timeout*=None, *wait\_callback*=None)

Pause execution of the script until device overlapped commands (e.g., taking sweeps) are complete.

*timeout* and *wait\_callback* override default constructor parameters.

**write** (*msg*, *arg*=None, *arg\_type*=None, *unit*=None, *bool\_selector*=None, *wait\_sync*=None,

*read\_echo*=False, *read\_echo\_delay*=0.0)

Send a command.

#### Parameters

- **msg** (*str*) – Text message.

- **arg** – Optional argument to append in the end. If a list of arguments is supplied, the result is joined with ", ".
- **arg\_type** (*str*) – Argument type. Can be 'raw' (in which case data is sent raw), 'string', 'int', 'float', 'bool', a format string (such as '{:.3f}') or a list of argument types (for an iterable argument); if format string is used and the argument is a list or a tuple, then it is expanded as a list of arguments (e.g., `arg_type='{0};{1}'` with `arg=[1,2]` will produce a string '1;2'); if a list of types is used, each element of *arg* is converted using the corresponding type, and the result is joined with ", ".
- **unit** (*str*) – If not None, use it as a unit to append after the value.
- **bool\_selector** (*tuple*) – A tuple (*false\_value*, *true\_value*) of two strings to represent bool argument; by default, use `._bool_selector` attribute.
- **wait\_sync** – if True, append the sync command (specified as `._wait_sync_comm` attribute, `"*OPC?"` by default) after the message and pause the execution command is complete; useful in long set operations, where the device might ignore later inputs until the current command is complete; if None, use the class default `._default_write_sync` attribute (False by default).
- **read\_echo** (*bool*) – If True, read a single line after write.
- **read\_echo\_delay** (*float*) – The delay between write and read if `read_echo==True`.

**class** `pylablib.devices.Tektronix.base.DPO2000` (*addr*, *nchannels*='auto')

Bases: `pylablib.devices.Tektronix.base.ITektronixScope`

Tektronix DPO2000 series oscilloscope.

#### Parameters

- **addr** – device address; usually a VISA address string such as `"USB0::0x0699::0x0364::C000000::INSTR"`
- **nchannels** – can specify number of channels on the oscilloscope; by default, autodetect number of channels (might take several seconds on connection)

#### BackendError

alias of `pylablib.core.devio.comm_backend.DeviceBackendError`

#### Error

alias of `TektronixError`

**class** `NoParameterCaller` (*device*, *kind*)

Bases: `object`

Class to simplify calling functions without a parameter

#### ReraiseError

alias of `TektronixBackendError`

**apply\_settings** (*settings*)

Apply the settings.

*settings* is a dict {*name*: *value*} of the available device settings. Non-applicable settings are ignored.

**ask** (*msg*, *data\_type*='string', *delay*=0.0, *timeout*=None, *read\_echo*=False)

Write a message and read a reply.

*msg* is the query message, *delay* is the delay between write and read. Other parameters are the same as in `read()`. If `read_echo==True`, assume that the device first echoes the input and skip it.

**close()**  
Close the backend

**enable\_channel** (*channel*, *enabled=True*)  
Enable or disable given channel

**flush** (*one\_line=False*)  
Flush the read buffer (read all the available data and return the number of bytes read).  
  
If *one\_line==True*, read only a single line.

**force\_trigger()**  
Force trigger event

**static get\_arg\_type** (*arg*)  
Autodetect argument type

**get\_channels** (*only\_main=False*)  
Get the list of all input channels (if *only\_main==True*) or all available channels (if *only\_main==False*)

**get\_channels\_number()**  
Get the number of channels

**get\_coupling** (*channel*)  
Get channel coupling.  
  
Can be "ac", "dc", or "gnd".

**get\_data\_format()**  
Get data transfer format.  
  
Return a string describing the format; can be either "ascii", or a numpy-style format string (e.g., "<u2").

**get\_data\_pts\_range()**  
Get range of data points to read.  
  
The range is defined from 1 to the points number (returned by *get\_points\_number()*).

**get\_device\_variable** (*key*)  
Get the value of a settings, status, or full info parameter

**get\_edge\_trigger\_coupling()**  
Get edge trigger coupling ("ac" or "dc")

**get\_edge\_trigger\_slope()**  
Get edge trigger slope ("fall" or "rise")

**get\_edge\_trigger\_source()**  
Get edge trigger source.  
  
Can be an integer indicating channel number or a name of a special channel.

**get\_esr** (*timeout=None*)  
Get the device status register (by default, "\*ESR?" command)

**get\_full\_info** (*include=0*)  
Get dict {name: value} containing full device information (including status and settings).  
  
*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_full\_status** (*include=0*)

Get dict {name: value} containing the device status (including settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_horizontal\_offset** ()

Get horizontal offset (position of the center of the sweep; in seconds)

**get\_horizontal\_span** ()

Get horizontal span (in seconds)

**get\_id** (*timeout=None*)

Get the device IDN. (query SCPI '\*IDN?' command)

**get\_points\_number** (*kind='send'*)

Get number of datapoints in various context.

*kind* defines the context. It can be "acq" (number of points acquired), "trace" (number of points in the source of the read-out trace; can be lower than "acq" if the data resolution is reduced, or if the source is not a channel data), or "send" (number of points in the sent waveform; can be lower than "trace" if [get\\_data\\_pts\\_range\(\)](#) is used to specify and incomplete range). Not all kinds are defined for all scope model (e.g., "trace" is not defined for TDS2000 series oscilloscopes).

For length of read-out trace, see also [get\\_data\\_pts\\_range\(\)](#).

**get\_probe\_attenuation** (*channel*)

Get channel probe attenuation

**get\_selected\_channel** ()

Get selected source channel.

Return number if it is a real channel, or a string name otherwise.

**get\_settings** (*include=0*)

Get dict {name: value} containing all the device settings.

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_trigger\_level** ()

Get edge trigger level (in Volts)

**get\_trigger\_mode** ()

Get trigger mode.

Can be either "auto" or "norm".

**get\_trigger\_state** ()

Get trigger state.

Can be "armed" (acquiring pretrigger), "ready" (pretrigger acquired, wait for trigger event), "trigger" (triggered, acquiring the rest of the waveform), "auto" ("auto" mode trigger is acquiring data in the absence of trigger), "save" (acquisition is stopped), or "scan" (oscilloscope in the scan mode)

**get\_vertical\_position** (*channel*)

Get channel vertical position (offset of the zero volt line; in V)



**get\_vertical\_span** (*channel*)  
Get channel vertical span (in V)

**get\_wfmpre** (*channel=None, enable=True*)  
Get preamble dictionary describing all scaling and format data for the given channel or a list of channels.  
  
Can be acquired once and used in subsequent multiple reads to save time on re-requesting. If *channel* is *None*, use the currently selected channel. If *enable==True*, make sure that the requested channel is enabled; getting preamble for disabled channels raises an error.

**grab\_continuous** (*enable=True*)  
Start or stop continuous grabbing

**grab\_single** (*wait=True, software\_trigger=False, wait\_timeout=None*)  
Set single waveform grabbing and wait for acquisition.  
  
If *wait==True*, wait until the acquisition is complete; otherwise, return immediately. if *software\_trigger==True*, send the software trigger after setup (i.e., the device triggers immediately regardless of the input).

**is\_channel\_enabled** (*channel*)  
Check if channel is enabled

**is\_continuous** ()  
Check if grabbing is continuous or single

**is\_grabbing** ()  
Check if acquisition is in progress.  
  
Return *True* if the oscilloscope is recording data, or if the trigger is armed/ready and waiting; return *False* if the acquisition is stopped. To check if the trigger has been triggered, use `get_trigger_state()`.

**is\_opened** ()  
Check if the device is connected

**lock** (*timeout=None*)  
Lock the access to the device from other threads/processes (isn't necessarily implemented)

**locking** (*timeout=None*)  
Context manager for lock & unlock

**normalize\_channel\_name** (*channel*)  
Normalize channel name as represented by the oscilloscope

**open** ()  
Open the backend

**static parse\_array\_data** (*data, fmt, include\_header=False*)  
Parse the data returned by the device. *fmt* is *DataFormat* description in numpy format (e.g., "<u2").  
  
If *include\_header==True*, the data is assumed to be in a (somewhat) standard SCPI format: b '# ', then a single digit *s* denoting length of the size block, then *s* digits denoting length of the data (in bytes) followed by the actual data. Otherwise (*include\_header==False*), assume that the header is already removed.

**read** (*data\_type='string', timeout=None*)  
Read data from the device.  
  
*data\_type* determines the type of the data. Can be 'raw' (just raw data), 'string' (with trailing and leading spaces stripped), 'int', 'float', 'bool' (interprets 0 or 'off' as *False*, anything else as *True*), 'value' (returns tuple (*value*, *unit*), where *value* is float), a callable (return the result of this callable applied to the string value), a dictionary (return the stored value corresponding to the

string value, or to the value converted into integer if the string value is not present), or a list of data types (the result is treated as a list of values with the given types separated by commas). *timeout* overrides the default value.

**read\_binary\_array\_data** (*include\_header=False, timeout=None, flush\_term=True*)

Read a binary data in the from the device.

The data assumes the standard binary transfer header consisting of "#" symbol, then a single digit with the size of the length string, then the length string containing the length of the binary data (in bytes). If *include\_header==True*, return the data with the header; otherwise, return only the content. If *flush\_term==True*, flush the following line to skip terminator characters after the binary data, which are added by some devices. *timeout* overrides the default value.

**read\_multiple\_sweeps** (*channels, wfmpres=None, ensure\_fmt=False, timeout=None, return\_wfmpres=None*)

Read data from a multiple channels channel.

#### Parameters

- **channels** – list of channel indices or names
- **wfmpres** – optional list or dictionary of preambles (obtained using `get_wfmpre()`); if it is *None*, obtain during reading, which slows down the data acquisition a bit
- **ensure\_fmt** – if *True*, make sure that oscilloscope data format agrees with the one in *wfmpre*
- **timeout** – read timeout
- **return\_wfmpres** – if *True*, return tuple (*sweeps*, *wfmpres*), where *wfmpres* can be used for further sweep readouts.

**read\_raw\_data** (*channel=None, fmt=None, timeout=None*)

Request, read and parse raw data at a given channel.

*fmt* is data format (e.g., "i1", "<i2", or "ascii") or "default", which uses the default oscilloscope format (usually binary with smallest appropriate byte size). If *fmt* is *None*, use the current format. If *channel* is *None*, use the currently selected channel.

Returned data is raw (i.e., not scaled and without x axis).

**read\_sweep** (*channel, wfmpre=None, ensure\_fmt=True, timeout=None*)

Read data from a single channel.

#### Parameters

- **channel** – channel index or name
- **wfmpre** – optional preamble dictionary (obtained using `get_wfmpre()`); if it is *None*, obtain during reading, which slows down the data acquisition a bit
- **ensure\_fmt** – if *True*, make sure that oscilloscope data format agrees with the one in *wfmpre*
- **timeout** – read timeout

**reconnect** (*new\_instrument=True, ignore\_error=True*)

Remake the connection.

If *new\_instrument==True*, create a new backend instance. If *ignore\_error==True*, ignore errors on closing.

**reset** ()

Reset the device (by default, "\*RST" command)

**select\_channel** (*channel*)

Select a channel to read data.

Doesn't need to be called explicitly, if `read_multiple_sweeps()` or `read_sweep()` are used.

**set\_coupling** (*channel*, *coupling*='dc')

Set channel coupling.

Can be "ac", "dc", or "gnd".

**set\_data\_format** (*fmt*='default')

Set data transfer format.

*fmt* is a string describing the format; can be either "ascii", or a numpy-style format string (e.g., "<u2"). If "default", use the oscilloscope default format (usually binary with smallest appropriate byte size).

**set\_data\_pts\_range** (*rng*=None)

Set range of data points to read.

The range is defined from 1 to the points number (returned by `get_points_number()` with *kind*="acq"). If *rng* is None, set the full range.

**set\_device\_variable** (*key*, *value*)

Set the value of a settings parameter

**set\_edge\_trigger\_coupling** (*coupling*)

Set edge trigger coupling ("ac" or "dc")

**set\_edge\_trigger\_slope** (*slope*)

Set edge trigger slope ("fall" or "rise")

**set\_edge\_trigger\_source** (*channel*)

Get edge trigger source.

Can be an integer indicating channel number or a name of a special channel.

**set\_horizontal\_offset** (*offset*=0.0)

Set horizontal offset (position of the center of the sweep; in seconds)

**set\_horizontal\_span** (*span*)

Set horizontal span (in seconds)

**set\_points\_number** (*pts\_num*, *reset\_limits*=True)

Set number of datapoints to record when acquiring a trace.

If *reset\_limits*==True, reset the datapoints range (`set_data_pts_range()`) to the full range. The actual set value (returned by this method) can be different from the requested value.

**set\_probe\_attenuation** (*channel*, *attenuation*)

Set channel probe attenuation

**set\_trigger\_level** (*level*)

Set edge trigger level (in Volts)

**set\_trigger\_mode** (*trigger\_mode*='auto')

Set trigger mode.

Can be either "auto" or "norm".

**set\_vertical\_position** (*channel*, *position*)

Set channel vertical position (offset of the zero volt line; in V)

**set\_vertical\_span** (*channel*, *span*)

Set channel vertical span (in V)

**setup\_edge\_trigger** (*source, level, coupling='dc', slope='rise'*)

Setup edge trigger.

Set source, level, coupling and slope (see corresponding methods for details).

**sleep** (*delay*)

Wait for *delay* seconds

**stop\_grabbing** ()

Stop grabbing or waiting (equivalent to `self.grab_continuous(False)`)

**unlock** ()

Unlock the access to the device from other threads/processes (isn't necessarily implemented)

**using\_write\_buffer** ()

Context manager for using a write buffer.

While it's active, all the consecutive `write()` operations are bundled together with `;` delimiter. The actual write is performed at the `read()/ask()` operation or at the end of the block.

**wait** (*wait\_type='sync', timeout=None, wait\_callback=None*)

Pause execution until device overlapped commands are complete.

*wait\_type* is either 'sync' (perform `wait_sync()`), 'dev' (perform `wait_dev()`) or 'none' (do nothing).

**wait\_dev** ()

Pause execution of the device commands until device overlapped commands (e.g., taking sweeps) are complete.

Note that the code execution is not paused.

**wait\_for\_grabbing** (*timeout=None*)

Wait until the acquisition is complete

**wait\_sync** (*timeout=None, wait\_callback=None*)

Pause execution of the script until device overlapped commands (e.g., taking sweeps) are complete.

*timeout* and *wait\_callback* override default constructor parameters.

**write** (*msg, arg=None, arg\_type=None, unit=None, bool\_selector=None, wait\_sync=None, read\_echo=False, read\_echo\_delay=0.0*)

Send a command.

#### Parameters

- **msg** (*str*) – Text message.
- **arg** – Optional argument to append in the end. If a list of arguments is supplied, the result is joined with `" , "`.
- **arg\_type** (*str*) – Argument type. Can be 'raw' (in which case data is sent raw), 'string', 'int', 'float', 'bool', a format string (such as `'{: .3f}'`) or a list of argument types (for an iterable argument); if format string is used and the argument is a list or a tuple, then it is expanded as a list of arguments (e.g., `arg_type='{0};{1}'` with `arg=[1, 2]` will produce a string `'1; 2'`); if a list of types is used, each element of *arg* is converted using the corresponding type, and the result is joined with `" , "`.
- **unit** (*str*) – If not None, use it as a unit to append after the value.
- **bool\_selector** (*tuple*) – A tuple (*false\_value*, *true\_value*) of two strings to represent bool argument; by default, use `._bool_selector` attribute.

- **wait\_sync** – if `True`, append the sync command (specified as `._wait_sync_comm` attribute, `"*OPC?"` by default) after the message and pause the execution command is complete; useful in long set operations, where the device might ignore later inputs until the current command is complete; if `None`, use the class default `._default_write_sync` attribute (`False` by default).
- **read\_echo** (*bool*) – If `True`, read a single line after write.
- **read\_echo\_delay** (*float*) – The delay between write and read if `read_echo==True`.

## Module contents

### pylablib.devices.Thorlabs package

#### Submodules

#### pylablib.devices.Thorlabs.TLCamera module

```
class pylablib.devices.Thorlabs.TLCamera.LibraryController (lib)
    Bases: pylablib.devices.utils.load_lib.LibraryController

    close (opid)
        Mark device closing.

        Return tuple (close_result, uninit_result) with the results of the closing and the shutdown.
        If library does not need to be shut down yet, set uninit_result=None

    open ()
        Mark device opening.

        Return tuple (init_result, open_result, opid) with the results of the initialization and the
        opening, and the opening ID which should afterwards be used for closing. If library is already initialized,
        set init_result=None

    preinit ()
        Pre-initialize the library, if it hasn't been done already

    shutdown ()
        Close all opened connections and shutdown the library

    temp_open ()
        Context for temporarily opening a new device connection

pylablib.devices.Thorlabs.TLCamera.list_cameras ()
    List connected TLCamera cameras

pylablib.devices.Thorlabs.TLCamera.get_cameras_number ()
    Get number of connected TLCamera cameras

class pylablib.devices.Thorlabs.TLCamera.TDeviceInfo (model, name, serial_number,
                                                    firmware_version)
    Bases: tuple

    count ()
        Return number of occurrences of value.

    firmware_version
```

**index()**  
Return first index of value.  
Raises ValueError if the value is not present.

**model**

**name**

**serial\_number**

**class** pylablib.devices.Thorlabs.TLCamera.**TFrameInfo** (*frame\_index, framestamp, pixel-clock, pixeltype, offset*)

Bases: `tuple`

**count()**  
Return number of occurrences of value.

**frame\_index**

**framestamp**

**index()**  
Return first index of value.  
Raises ValueError if the value is not present.

**offset**

**pixelclock**

**pixeltype**

**class** pylablib.devices.Thorlabs.TLCamera.**ThorlabsTLCamera** (*serial=None*)  
Bases: `pylablib.devices.interface.camera.IBinROICamera`, `pylablib.devices.interface.camera.IExposureCamera`

Thorlabs TSI camera.

**Parameters** **serial** (*str*) – camera serial number; can be either a string obtained using `list_cameras()` function, or None, which means connecting to the first available camera (not recommended unless only one camera is connected)

**Error** = <Mock name='mock.ThorlabsTLCameraError' id='140096289851216'>

**TimeoutError** = <Mock spec='str' id='140096289849936'>

**open()**  
Open connection to the camera

**close()**  
Close connection to the camera

**is\_opened()**  
Check if the device is connected

**get\_device\_info()**  
Get camera model data.  
Return tuple (model, name, serial\_number, firmware\_version).

**class** **RingBuffer**

Bases: `object`

Frames ring buffer.

Reacts to each new frame and stores it in the internal buffer.

**reset()**  
Reset buffer and internal counters

**setup(buffsize, frame\_dim)**  
Setup a new buffer with the given maximal number of frames and frame dimensions

**cleanup()**  
Cleanup the buffer

**new\_frame(handle, buffer, idx, metadata, metadata\_size, context)**  
Callback for receiving a new frame

**wait\_for\_frame(idx=None, timeout=None)**  
Wait for a new frame acquisition

**get\_frame(idx)**  
Get the frame with the given index (or None if it is outside the buffer range)

**get\_status()**  
Get buffer status (acquired, missed, stored)

**get\_frame\_timings()**  
Get acquisition timing.  
  
Return tuple (exposure, frame\_period).

**set\_exposure(exposure)**  
Set camera exposure

**get\_trigger\_mode()**  
Get trigger mode.  
  
Can be "int" (internal/software), "ext" (external/hardware), or "bulb" (bulb trigger).

**set\_trigger\_mode(mode)**  
Set trigger mode.  
  
Can be "int" (internal/software), "ext" (external/hardware), or "bulb" (bulb trigger).

**get\_ext\_trigger\_parameters()**  
Return external trigger polarity

**setup\_ext\_trigger(polarity)**  
Setup external trigger polarity ("rise" or "fall")

**send\_software\_trigger()**  
Send software trigger signal

**get\_pixel\_correction\_parameters()**  
Return pixel correction parameters (enabled, threshold)

**setup\_pixel\_correction(enable=True, threshold=None)**  
Enable or disable hotpixel correction and set its threshold (None means keep unchanged)

**get\_timestamp\_clock\_frequency()**  
Return frequency of the frame timestamp clock (in Hz)

**setup\_acquisition(nframes=100)**  
Setup acquisition.  
  
*nframes* determines number of size of the ring buffer (by default, 100).

**clear\_acquisition()**  
Clear acquisition settings

**start\_acquisition** (*frames\_per\_trigger='default', auto\_start=True, nframes=None*)

Start camera acquisition.

#### Parameters

- **frames\_per\_trigger** – number of frames to acquire per trigger (software of hardware); *None* means unlimited number; by default, set to *None* for software trigger (i.e., run until stopped), and 1 for hardware trigger (i.e., one frame per trigger pulse)
- **auto\_start** – if *True* and the trigger is set into software mode, automatically start recording; otherwise, only start recording when *send\_software\_trigger()* is called explicitly; this value is meaningless in the hardware or bulb trigger mode
- **nframes** – number of frames in the ring buffer

**stop\_acquisition** ()

Stop acquisition

**acquisition\_in\_progress** ()

Check if acquisition is in progress

**get\_detector\_size** ()

Get camera detector size (in pixels) as a tuple (*width, height*)

**get\_roi** ()

Get current ROI.

Return tuple (*hstart, hend, vstart, vend, hbin, vbin*). *hstart* and *hend* specify horizontal image extent, *vstart* and *vend* specify vertical image extent (start is inclusive, stop is exclusive, starting from 0), *hbin* and *vbin* specify binning.

**set\_roi** (*hstart=0, hend=None, vstart=0, vend=None, hbin=1, vbin=1*)

Setup camera ROI.

*hstart* and *hend* specify horizontal image extent, *vstart* and *vend* specify vertical image extent (start is inclusive, stop is exclusive, starting from 0), *hbin* and *vbin* specify binning. By default, all non-supplied parameters take extreme values (0 for start, maximal for end, 1 for binning).

**get\_roi\_limits** (*hbin=1, vbin=1*)

Get the minimal and maximal ROI parameters.

Return tuple (*hlim, vlim*), where each element is in turn a limit 5-tuple (*min, max, pstep, sstep, maxbin*) with, correspondingly, minimal and maximal size, position and size step, and the maximal binning. In some cameras, the step and the minimal size depend on the binning, which can be supplied.

**read\_multiple\_images** (*rng=None, peek=False, missing\_frame='skip', return\_info=False, return\_rng=False*)

Read multiple images specified by *rng* (by default, all un-read images).

If *rng* is specified, it is a tuple (*first, last*) with images range (first inclusive). If no new frames are available, return an empty list; if no acquisition is running, return *None*. If *peek==True*, return images but not mark them as read. *missing\_frame* determines what to do with frames which are out of range (missing or lost): can be "none" (replacing them with *None*), "zero" (replacing them with zero-filled frame), or "skip" (skipping them). If *return\_info==True*, return tuple (*frames, infos*), where *infos* is a list of *TFrameInfo* instances describing frame index and frame metadata, which contains framestamp, pixel clock, pixel format, and pixel offset; if some frames are missing and *missing\_frame!="skip"*, the corresponding frame info is *None*. if *return\_rng==True*, return the range covered resulting frames; if *missing\_frame=="skip"*, the range can be smaller than the supplied *rng* if some frames are skipped.



**FrameTransferError**

alias of `pylablib.devices.interface.camera.DefaultFrameTransferError`

**class NoParameterCaller** (*device, kind*)

Bases: `object`

Class to simplify calling functions without a parameter

**apply\_settings** (*settings*)

Apply the settings.

*settings* is the dict {name: value} of the device available settings. Non-applicable settings are ignored.

**get\_acquisition\_parameters** ()

Get acquisition parameters.

Return dictionary {name: value}

**get\_data\_dimensions** ()

Get readout data dimensions (in pixels) as a tuple (width, height); take indexing mode into account

**get\_device\_variable** (*key*)

Get the value of a settings, status, or full info parameter

**get\_exposure** ()

Get current exposure

**get\_frame\_format** ()

Get format for the returned images.

Can be "list" (list of 2D arrays), "array" (a single 3D array), or "chunks" (list of 3D "chunk" arrays; supported for some cameras and provides the best performance).

**get\_frame\_info\_fields** ()

Get the names of frame info fields.

Applicable when frame info format (set by `set_frame_info_format()`) is "list" or "array".

**get\_frame\_info\_format** ()

Get format of the frame info.

Can be "namedtuple" (potentially nested named tuples; convenient to get particular values), "list" (flat list of values, with field names are given by `get_frame_info_fields()`; convenient for building a table), "array" (same as "list", but with a numpy array, which is easier to use for "chunks" frame format), or "dict" (flat dictionary with the same fields as the "list" format; more resilient to future format changes)

**get\_frame\_info\_period** ()

Get period of frame info acquisition.

Frame info might be skipped (set to None) except for frames which indices are divisible by *period*. Useful for certain cameras where acquiring frame info takes a lot of time and can reduce performance at higher frame rates. Note that this parameter can still be ignored (i.e., always set to 1) if the performance is not an issue for a given camera class.

**get\_frame\_period** ()

Get frame period (time between two consecutive frames in the internal trigger mode)

**get\_frames\_status** ()

Get acquisition and buffer status.

Return tuple (acquired, unread, skipped, size), where acquired is the total number of acquired frames, unread is the number of acquired but not read frames, skipped is the number of skipped (not read and then written over) frames, and buffer\_size is the total buffer size (in frames).

**get\_full\_info** (include=0)

Get dict {name: value} containing full device information (including status and settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include*=-10 queries all available variables, which is equivalent to *include*="all".

**get\_full\_status** (include=0)

Get dict {name: value} containing the device status (including settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include*=-10 queries all available variables, which is equivalent to *include*="all".

**get\_image\_indexing** ()

Get indexing for the returned images.

Can be "rct" (first index row, second index column, rows counted from the top), "rcb" (same as "rc", rows counted from the bottom), "xyt" (first index column, second index row, rows counted from the top), or "xyb" (same as "xyt", rows counted from the bottom)

**get\_new\_images\_range** ()

Get the range of the new images.

Return tuple (first, last) with images range (first inclusive). If no images are available, return None. If some images were in the buffer were overwritten, exclude them from the range.

**get\_settings** (include=0)

Get dict {name: value} containing all the device settings.

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include*=-10 queries all available variables, which is equivalent to *include*="all".

**grab** (nframes=1, frame\_timeout=5.0, missing\_frame='skip', return\_info=False, buff\_size=None)

Snap *nframes* images (with preset image read mode parameters)

*buff\_size* determines buffer size (if None, use the default size). Timeout is specified for a single-frame acquisition, not for the whole acquisition time. *missing\_frame* determines what to do with frames which have been lost: can be "none" (replacing them with None), "zero" (replacing them with zero-filled frame), or "skip" (skipping them, while still keeping total returned frames number to *n*). If *return\_info*==True, return tuple (frames, infos), where *infos* is a list of frame info tuples (camera-dependent); if some frames are missing and *missing\_frame*!="skip", the corresponding frame info is None.

**is\_acquisition\_setup** ()

Check if acquisition is set up.

If the camera does not support separate acquisition setup, always return True.

**pausing\_acquisition** (clear=None, stop=True, setup\_after=None, start\_after=True, combine\_nested=True)

Context manager which temporarily pauses acquisition during execution of with block.

Useful for applying certain settings which can't be changed during the acquisition. If `clear==True`, clear acquisition in addition to stopping (by default, use the class default specified as `_clear_pausing_acquisition` attribute). If `stop==True`, stop the acquisition (if `clear==True`, stop regardless). If `setup_after==True`, setup the acquisition after pause if necessary (None means setup only if clearing was required). If `start_after==True`, start the acquisition after pause if necessary (None means start only if stopping was required). If `combine_nested==True`, then any nested `pausing_acquisition` calls will stop/clear acquisition as necessary, but won't setup/start it again until this `pausing_acquisition` call is complete.

Yields tuple (`acq_in_progress`, `acq_params`), which indicates whether acquisition is currently in progress, and what are the current acquisition parameters.

**read\_newest\_image** (*peek=False, return\_info=False*)

Read the newest un-read image.

If no un-read frames are available, return None. If `peek==True`, return the image but not mark it as read. If `return_info==True`, return tuple (`frame`, `info`), where `info` is an info tuples (camera-dependent, see `read_multiple_images()`).

**read\_oldest\_image** (*peek=False, return\_info=False*)

Read the oldest un-read image.

If no un-read frames are available, return None. If `peek==True`, return the image but not mark it as read. If `return_info==True`, return tuple (`frame`, `info`), where `info` is an info tuples (camera-dependent, see `read_multiple_images()`).

**set\_device\_variable** (*key, value*)

Set the value of a settings parameter

**set\_frame\_format** (*fmt*)

Set format for the returned images.

Can be "list" (list of 2D arrays), "array" (a single 3D array), "chunks" (list of 3D "chunk" arrays; supported for some cameras and provides the best performance), or "try\_chunks" (same as "chunks", but if chunks are not supported, set to "list" instead). If format is "chunks" and chunks are not supported by the camera, it results in one frame per chunk. Note that if the format is set to "array" or "chunks", the frame info format is also automatically set to "array". If the format is set to "chunks", then the image info is also returned in chunks form (list of 2D info arrays with the same length as the corresponding frame chunks).

**set\_frame\_info\_format** (*fmt, include\_fields=None*)

Set format of the frame info.

Can be "namedtuple" (potentially nested named tuples; convenient to get particular values), "list" (flat list of values, with field names are given by `get_frame_info_fields()`; convenient for building a table), "array" (same as "list", but with a numpy array, which is easier to use for "chunks" frame format), or "dict" (flat dictionary with the same fields as the "list" format; more resilient to future format changes) If `include_fields` is not None, it specifies the fields included for non-"tuple" formats; note that order or `include_fields` is ignored, and the resulting fields are always ordered same as in the original.

**set\_frame\_info\_period** (*period=1*)

Set period of frame info acquisition.

Frame info might be skipped (set to None) except for frames which indices are divisible by *period*. Useful for certain cameras where acquiring frame info takes a lot of time and can reduce performance at higher frame rates. Note that this parameter can still be ignored (i.e., always set to 1) if the performance is not an issue for a given camera class.

**set\_image\_indexing** (*indexing*)

Set up indexing for the returned images.

Can be "rct" (first index row, second index column, rows counted from the top), "rcb" (same as "rc", rows counted from the bottom), "xyt" (first index column, second index row, rows counted from the top), or "xyb" (same as "xyt", rows counted from the bottom)

**snap** (*timeout=5.0, return\_info=False*)

Snap a single frame

**wait\_for\_frame** (*since='lastread', nframes=1, timeout=20.0, error\_on\_stopped=False*)

Wait for one or several new camera frames.

*since* specifies the reference point for waiting to acquire *nframes* frames; can be "lastread" (from the last read frame), "lastwait" (wait for the last successful `wait_for_frame()` call), "now" (from the start of the current call), or "start" (from the acquisition start, i.e., wait until *nframes* frames have been acquired). *timeout* can be either a number, None (infinite timeout), or a tuple (*timeout*, *frame\_timeout*), in which case the call times out if the total time exceeds *timeout*, or a single frame wait exceeds *frame\_timeout*. If the call times out, raise `TimeoutError`. If *error\_on\_stopped*==True and the acquisition is not running, raise `Error`; otherwise, simply return False without waiting.

## pylablib.devices.Thorlabs.base module

**exception** `pylablib.devices.Thorlabs.base.ThorlabsError`

Bases: `pylablib.core.devio.base.DeviceError`

Generic Thorlabs error

**args**

**with\_traceback** ()

Exception.with\_traceback(tb) – set self.\_\_traceback\_\_ to tb and return self.

**exception** `pylablib.devices.Thorlabs.base.ThorlabsBackendError` (*exc*)

Bases: `pylablib.devices.Thorlabs.base.ThorlabsError`, `pylablib.core.devio.comm_backend.DeviceBackendError`

Thorlabs backend communication error

**args**

**with\_traceback** ()

Exception.with\_traceback(tb) – set self.\_\_traceback\_\_ to tb and return self.

**exception** `pylablib.devices.Thorlabs.base.ThorlabsTimeoutError`

Bases: `pylablib.devices.Thorlabs.base.ThorlabsError`

Thorlabs timeout error

**args**

**with\_traceback** ()

Exception.with\_traceback(tb) – set self.\_\_traceback\_\_ to tb and return self.

## pylablib.devices.Thorlabs.elliptec module

`pylablib.devices.Thorlabs.elliptec.muxaddr` (*\*args, argname='addr', \*\*kwargs*)

Multiplex the function over its *addr* argument

```
class pylablib.devices.Thorlabs.elliptec.TDeviceInfo(serial_no, model_no, year,  
                                                    fw_ver, hw_ver, travel, pulse)
```

Bases: `tuple`

**count**()

Return number of occurrences of value.

**fw\_ver**

**hw\_ver**

**index**()

Return first index of value.

Raises `ValueError` if the value is not present.

**model\_no**

**pulse**

**serial\_no**

**travel**

**year**

```
class pylablib.devices.Thorlabs.elliptec.TMotorInfo(loop, motor, current, ramp_up,  
                                                    ramp_down, fw_freq, bk_freq)
```

Bases: `tuple`

**bk\_freq**

**count**()

Return number of occurrences of value.

**current**

**fw\_freq**

**index**()

Return first index of value.

Raises `ValueError` if the value is not present.

**loop**

**motor**

**ramp\_down**

**ramp\_up**

```
class pylablib.devices.Thorlabs.elliptec.ElliptecMotor(conn,                addr='all',  
                                                    default_addr=None,  
                                                    scale='stage', timeout=3.0)
```

Bases: `pylablib.core.devio.comm_backend.ICommBackendWrapper`

Basic Elliptec stage device.

#### Parameters

- **conn** – serial connection parameters (usually port or a tuple containing port and baudrate)
- **addr** – list of device addresses (between 0 and 15) connected to this serial port; if "all", automatically detect all connected devices

- **default\_addr** – address used by default when not supplied; by default, use the first address among the connected
- **scale** – scale of the position units to the internals units; can be "stage" (use stage units such as mm or deg based on its internal calibration), "step" (directly use step units), or a number which multiplies user-supplied units to produce steps
- **timeout** – default communication timeout

**Error**

alias of `pylablib.devices.Thorlabs.base.ThorlabsError`

**get\_connected\_addrs()**

Get a list of all connected device addresses

**get\_default\_addr()**

Get the current default address

**set\_default\_addr(addr)**

Set the current default address

**using\_default\_addr(addr)**

Context manager which temporarily changes the default address

**send\_comm(comm, data="", addr=None)**

Send a message with the given data to the devices at a given address.

For details, see ELLx communications protocol.

**class CommData(comm, data, addr)**

Bases: `tuple`

**addr**

**comm**

**count()**

Return number of occurrences of value.

**data**

**index()**

Return first index of value.

Raises `ValueError` if the value is not present.

**recv\_comm(comm=None, addr=None, datalen='auto', timeout=None)**

Receive a message.

*comm*, *addr*, and *datalen* can specify the expected return command, address, or the length of the data field (if "auto", determine based on the return command). *timeout* specifies the waiting timeout (by default, same as supplied upon the device connection).

For details, see ELLx communications protocol.

**query(comm, data="", addr=None, reply\_comm=None, reply\_addr='auto', reply\_datalen='auto', timeout=None)**

Send a query to the device and receive the reply.

A combination of `send_comm()` and `recv_comm()`.

**add\_background\_comm(comm)**

Mark given *comm* as a 'background' message, which can be sent by the device at any point without prompt (e.g., some operation confirmation).

If it is received instead during `recv_comm` or `query` operations, it is ignored, and the corresponding counter is increased.

**check\_background\_comm** (*comm*, *addr=None*)

Return message counter and the last message value (`None` if not message received yet) of a given 'background' message received from the given address

**change\_addr** (*newaddr*, *addr=None*)

Change the device address to a new value (between 0 and 15)

**store\_parameters** (*addr=None*)

Store current device parameters (e.g., frequencies) to the energy-independent memory

**get\_device\_info** (*addr=None*)

Get device info.

Return tuple (*serial\_no*, *model\_no*, *year*, *fw\_ver*, *hw\_ver*, *travel*, *pulse*).

**get\_status** (*addr=None*)

Get device status

**get\_motor\_info** (*motor=1*, *addr=None*)

Get info for a given motor (between 1 and 3).

Return tuple (*loop\_ena*, *motor\_ena*, *current*, *ramp\_up*, *ramp\_down*, *fw\_freq*, *bk\_freq*).

**get\_scale** (*addr=None*)

Get scale parameter for the specified address.

Can be "stage", "step", or a proportionality coefficient.

**set\_scale** (*scale*, *addr=None*)

Set scale parameter for the specified address.

Can be "stage", "step", or a proportionality coefficient.

**home** (*home\_dir='cw'*, *paddles='all'*, *addr=None*)

Home the device.

The operation is synchronous, i.e., it will not finish until the homing is done. If the device is a rotary stage, then *home\_dir* specifies homing direction ("cw" or "ccw"). If the device is a paddle polarization controller, then *paddles* is a list of all paddle indices (1 to 3) to home ("all" is the same as [1, 2, 3]).

**get\_home\_offset** (*addr=None*)

Get homing offset

**set\_home\_offset** (*offset*, *addr=None*)

Set homing offset (note: the manufacturer advises against it)

**get\_velocity** (*addr=None*)

Get velocity as a percentage from the maximal velocity (0 to 100)

**set\_velocity** (*velocity=100*, *addr=None*)

Set velocity as a percentage from the maximal velocity (0 to 100)

**get\_frequency** (*motor=1*, *addr=None*)

Get frequencies at a given motor as a tuple (*fw\_freq*, *bk\_freq*)

**set\_frequency** (*fw\_freq=None*, *bk\_freq=None*, *motor=1*, *addr=None*)

Set frequencies at a given motor.

Values set as `None` stay the same.

**search\_frequency** (*motor=1, addr=None*)

Run the automated frequency search on a given motor.

The position might change slightly throughout the process.

**get\_position** (*addr=None*)

Get the current position

**move\_to** (*position, addr=None, timeout=30.0*)

Move to the given position.

The operation is synchronous, i.e., it will not finish until the motion is stopped. Return `True` if the position was reached successfully or `False` otherwise.

**move\_by** (*distance, addr=None, timeout=30.0*)

Move by the given distance.

The operation is synchronous, i.e., it will not finish until the motion is stopped. Return `True` if the position was reached successfully or `False` otherwise.

**class NoParameterCaller** (*device, kind*)

Bases: `object`

Class to simplify calling functions without a parameter

**apply\_settings** (*settings*)

Apply the settings.

*settings* is the dict `{name: value}` of the device available settings. Non-applicable settings are ignored.

**close** ()

Close the backend

**get\_device\_variable** (*key*)

Get the value of a settings, status, or full info parameter

**get\_full\_info** (*include=0*)

Get dict `{name: value}` containing full device information (including status and settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or `"all"` (all available variables). Since the lowest priority is -10, setting `include=-10` queries all available variables, which is equivalent to `include="all"`.

**get\_full\_status** (*include=0*)

Get dict `{name: value}` containing the device status (including settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or `"all"` (all available variables). Since the lowest priority is -10, setting `include=-10` queries all available variables, which is equivalent to `include="all"`.

**get\_settings** (*include=0*)

Get dict `{name: value}` containing all the device settings.

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or `"all"` (all available variables). Since the lowest priority is -10, setting `include=-10` queries all available variables, which is equivalent to `include="all"`.

**is\_opened** ()

Check if the device is connected



**lock** (*timeout=None*)  
Lock the access to the device from other threads/processes (isn't necessarily implemented)

**locking** (*timeout=None*)  
Context manager for lock & unlock

**open** ()  
Open the backend

**set\_device\_variable** (*key, value*)  
Set the value of a settings parameter

**unlock** ()  
Unlock the access to the device from other threads/processes (isn't necessarily implemented)

## pylablib.devices.Thorlabs.kinesis module

`pylablib.devices.Thorlabs.kinesis.list_kinesis_devices` (*filter\_ids=True*)  
List all Thorlabs APT/Kinesis devices connected to this PC.

**Return list of tuples (conn, description).** If `filter_ids==True`, only leave devices with Thorlabs-like IDs (8-digit numbers). Otherwise, show all devices (some of them might not be Thorlabs-related).

**class** `pylablib.devices.Thorlabs.kinesis.TDeviceInfo` (*serial\_no, model\_no, fw\_ver, hw\_type, hw\_ver, mod\_state, nchannels, notes*)

Bases: `tuple`

**count** ()  
Return number of occurrences of value.

**fw\_ver**

**hw\_type**

**hw\_ver**

**index** ()  
Return first index of value.  
Raises `ValueError` if the value is not present.

**mod\_state**

**model\_no**

**nchannels**

**notes**

**serial\_no**

**class** `pylablib.devices.Thorlabs.kinesis.BasicKinesisDevice` (*conn, timeout=3.0, is\_rack\_system=False*)

Bases: `pylablib.core.devio.comm_backend.ICommBackendWrapper`

Generic Kinesis device.

Implements FTDI chip connectivity via pyft232 (virtual serial interface).

### Parameters

- **conn** – serial connection parameters (usually an 8-digit device serial number).

- **is\_rack\_system** – specify whether the device is a rack system or a standalone USB device (default mode).

**Error**

alias of `pylablib.devices.Thorlabs.base.ThorlabsError`

**static list\_devices** (*filter\_ids=True*)

List all connected devices.

Return list of tuples (*conn*, *description*). If *filter\_ids==True*, only leave devices with Thorlabs-like IDs (8-digit numbers). Otherwise, show all devices (some of them might not be Thorlabs-related).

**send\_comm** (*messageID*, *param1=0*, *param2=0*, *source=1*, *dest='host'*)

Send a message with no associated data.

For details, see APT communications protocol.

**send\_comm\_data** (*messageID*, *data*, *source=1*, *dest='host'*)

Send a message with associated data.

For details, see APT communications protocol.

**class CommShort** (*messageID*, *param1*, *param2*, *source*, *dest*)

Bases: `tuple`

**count** ()

Return number of occurrences of value.

**dest****index** ()

Return first index of value.

Raises `ValueError` if the value is not present.

**messageID****param1****param2****source****class CommData** (*messageID*, *data*, *source*, *dest*)

Bases: `tuple`

**count** ()

Return number of occurrences of value.

**data****dest****index** ()

Return first index of value.

Raises `ValueError` if the value is not present.

**messageID****source****recv\_comm** (*expected\_id=None*)

Receive a message.

Return either *CommShort* or *CommData* depending on the message type (fixed length with two parameters, or variable length with associated data). If *expected\_id* is not *None* and the received message ID is different from *expected\_id*, raise an error. For details, see APT communications protocol.

**query** (*messageID*, *param1*=0, *param2*=0, *source*=1, *dest*='host', *replyID*=-1)

Send a query to the device and receive the reply.

A combination of *send\_comm()* and *recv\_comm()*. If *replyID* is not *None*, specifies the expected reply message ID; if -1 (default), set to be *messageID*+1 (the standard convention).

**add\_background\_comm** (*messageID*)

Mark given *messageID* as a 'background' message, which can be sent at any point without prompt (e.g., some operation confirmation).

If it is received instead during *recv\_comm\_* operations, it is ignored, and the corresponding counter is increased.

**check\_background\_comm** (*messageID*)

Return message counter and the last message value (*None* if not message received yet) of a given 'background' message

**get\_device\_info** (*dest*='host')

Get device info.

Return tuple (serial\_no, model\_no, fw\_ver, hw\_type, hw\_ver, mod\_state, nchannels, notes).

**get\_number\_of\_channels** ()

Get number of channels on the device

**blink** (*channel*=1, *dest*='host')

Identify the physical device (by, e.g., blinking status LED or screen)

**class NoParameterCaller** (*device*, *kind*)

Bases: *object*

Class to simplify calling functions without a parameter

**apply\_settings** (*settings*)

Apply the settings.

*settings* is the dict {name: value} of the device available settings. Non-applicable settings are ignored.

**close** ()

Close the backend

**get\_device\_variable** (*key*)

Get the value of a settings, status, or full info parameter

**get\_full\_info** (*include*=0)

Get dict {name: value} containing full device information (including status and settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include*=-10 queries all available variables, which is equivalent to *include*="all".

**get\_full\_status** (*include*=0)

Get dict {name: value} containing the device status (including settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the

lowest priority is -10, setting `include=-10` queries all available variables, which is equivalent to `include="all"`.

**get\_settings** (*include=0*)

Get dict {name: value} containing all the device settings.

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting `include=-10` queries all available variables, which is equivalent to `include="all"`.

**is\_opened** ()

Check if the device is connected

**lock** (*timeout=None*)

Lock the access to the device from other threads/processes (isn't necessarily implemented)

**locking** (*timeout=None*)

Context manager for lock & unlock

**open** ()

Open the backend

**set\_device\_variable** (*key, value*)

Set the value of a settings parameter

**unlock** ()

Unlock the access to the device from other threads/processes (isn't necessarily implemented)

**class** `pylablib.devices.Thorlabs.kinesis.TVelocityParams` (*min\_velocity, acceleration, max\_velocity*)

Bases: `tuple`

**acceleration**

**count** ()

Return number of occurrences of value.

**index** ()

Return first index of value.

Raises `ValueError` if the value is not present.

**max\_velocity**

**min\_velocity**

**class** `pylablib.devices.Thorlabs.kinesis.TJogParams` (*mode, step\_size, min\_velocity, acceleration, max\_velocity, stop\_mode*)

Bases: `tuple`

**acceleration**

**count** ()

Return number of occurrences of value.

**index** ()

Return first index of value.

Raises `ValueError` if the value is not present.

**max\_velocity**

**min\_velocity**

```
mode
step_size
stop_mode

class pylablib.devices.Thorlabs.kinesis.TGenMoveParams(backlash_distance)
    Bases: tuple

    backlash_distance

    count()
        Return number of occurrences of value.

    index()
        Return first index of value.

        Raises ValueError if the value is not present.

class pylablib.devices.Thorlabs.kinesis.THomeParams(home_direction, limit_switch, ve-
                                                    locity, offset_distance)
    Bases: tuple

    count()
        Return number of occurrences of value.

    home_direction

    index()
        Return first index of value.

        Raises ValueError if the value is not present.

    limit_switch

    offset_distance

    velocity

class pylablib.devices.Thorlabs.kinesis.TLimitSwitchParams(hw_kind_cw,
                                                            hw_kind_ccw,
                                                            hw_swapped,
                                                            sw_position_cw,
                                                            sw_position_ccw,
                                                            sw_kind)
    Bases: tuple

    count()
        Return number of occurrences of value.

    hw_kind_ccw

    hw_kind_cw

    hw_swapped

    index()
        Return first index of value.

        Raises ValueError if the value is not present.

    sw_kind

    sw_position_ccw

    sw_position_cw
```

```
class pylablib.devices.Thorlabs.kinesis.TKCubeTrigIOParams(trig1_mode,  
                                                         trig1_pol, trig2_mode,  
                                                         trig2_pol)  
  
    Bases: tuple  
  
    count()  
        Return number of occurrences of value.  
  
    index()  
        Return first index of value.  
  
        Raises ValueError if the value is not present.  
  
    trig1_mode  
    trig1_pol  
    trig2_mode  
    trig2_pol  
  
class pylablib.devices.Thorlabs.kinesis.TKCubeTrigPosParams(start_fw,    step_fw,  
                                                             num_fw,    start_bk,  
                                                             step_bk,    num_bk,  
                                                             width, ncycles)  
  
    Bases: tuple  
  
    count()  
        Return number of occurrences of value.  
  
    index()  
        Return first index of value.  
  
        Raises ValueError if the value is not present.  
  
    ncycles  
    num_bk  
    num_fw  
    start_bk  
    start_fw  
    step_bk  
    step_fw  
    width  
  
class pylablib.devices.Thorlabs.kinesis.TPZMotorDriveParams(max_voltage, veloc-  
                                                             ity, acceleration)  
  
    Bases: tuple  
  
    acceleration  
  
    count()  
        Return number of occurrences of value.  
  
    index()  
        Return first index of value.  
  
        Raises ValueError if the value is not present.  
  
    max_voltage
```

**velocity**

```
class pylablib.devices.Thorlabs.kinesis.TPZMotorJogParams(mode, step_size_fw,
                                                         step_size_bk, velocity,
                                                         acceleration)
```

Bases: `tuple`

**acceleration**

**count()**

Return number of occurrences of value.

**index()**

Return first index of value.

Raises ValueError if the value is not present.

**mode**

**step\_size\_bk**

**step\_size\_fw**

**velocity**

```
pylablib.devices.Thorlabs.kinesis.muxchannel(*args, **kwargs)
```

```
class pylablib.devices.Thorlabs.kinesis.KinesisDevice(conn, timeout=3.0,
                                                         default_channel=1,
                                                         is_rack_system=False)
```

Bases: `pylablib.devices.interface.stage.IMultiaxisStage`, `pylablib.devices.Thorlabs.kinesis.BasicKinesisDevice`

Overarching Kinesis class containing all of the necessary private methods.

Subclasses are expected to make some of them visible by renaming, and to define device variables and opening behavior accordingly.

#### Parameters

- **conn** – serial connection parameters (usually an 8-digit device serial number).
- **timeout** – device communication timeout.
- **default\_channel** – starting default channel used when no channel is supplied to a channel-level command (such as `move_to` or `get_position`).
- **is\_rack\_system** – specify whether the device is a rack system or a standalone USB device (default mode).

**get\_all\_channels()**

Get the list of all available channels; alias of `get_all_axes` method

**set\_default\_channel(channel)**

Set the default channel for all channel-related methods

**using\_channel(channel)**

Context manager for temporarily using a different default channel

```
status_bits = [(1, 'sw_bk_lim'), (2, 'sw_fw_lim'), (16, 'moving_bk'), (32, 'moving_fw')]
```

```
class CommData(messageID, data, source, dest)
```

Bases: `tuple`

**count()**

Return number of occurrences of value.

**data**

**dest**

**index()**  
Return first index of value.  
Raises ValueError if the value is not present.

**messageID**

**source**

**class CommShort** (*messageID, param1, param2, source, dest*)  
Bases: `tuple`

**count()**  
Return number of occurrences of value.

**dest**

**index()**  
Return first index of value.  
Raises ValueError if the value is not present.

**messageID**

**param1**

**param2**

**source**

**Error**  
alias of `pylablib.devices.Thorlabs.base.ThorlabsError`

**class NoParameterCaller** (*device, kind*)  
Bases: `object`  
Class to simplify calling functions without a parameter

**add\_background\_comm** (*messageID*)  
Mark given messageID as a ‘background’ message, which can be sent at any point without prompt (e.g., some operation confirmation).  
  
If it is received instead during `recv_comm_` operations, it is ignored, and the corresponding counter is increased.

**apply\_settings** (*settings*)  
Apply the settings.  
  
*settings* is the dict {name: value} of the device available settings. Non-applicable settings are ignored.

**blink** (*channel=1, dest='host'*)  
Identify the physical device (by, e.g., blinking status LED or screen)

**check\_background\_comm** (*messageID*)  
Return message counter and the last message value (None if not message received yet) of a given ‘background’ message

**close()**  
Close the backend



**get\_all\_axes()**  
Get the list of all available axes (taking mapping into account)

**get\_device\_info** (*dest='host'*)  
Get device info.  
  
Return tuple (serial\_no, model\_no, fw\_ver, hw\_type, hw\_ver, mod\_state, nchannels, notes).

**get\_device\_variable** (*key*)  
Get the value of a settings, status, or full info parameter

**get\_full\_info** (*include=0*)  
Get dict {name: value} containing full device information (including status and settings).  
  
*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_full\_status** (*include=0*)  
Get dict {name: value} containing the device status (including settings).  
  
*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_number\_of\_channels()**  
Get number of channels on the device

**get\_settings** (*include=0*)  
Get dict {name: value} containing all the device settings.  
  
*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**is\_opened()**  
Check if the device is connected

**static list\_devices** (*filter\_ids=True*)  
List all connected devices.  
  
Return list of tuples (conn, description). If *filter\_ids==True*, only leave devices with Thorlabs-like IDs (8-digit numbers). Otherwise, show all devices (some of them might not be Thorlabs-related).

**lock** (*timeout=None*)  
Lock the access to the device from other threads/processes (isn't necessarily implemented)

**locking** (*timeout=None*)  
Context manager for lock & unlock

**open()**  
Open the backend

**query** (*messageID, param1=0, param2=0, source=1, dest='host', replyID=-1*)  
Send a query to the device and receive the reply.  
  
A combination of `send_comm()` and `recv_comm()`. If *replyID* is not None, specifies the expected reply message ID; if -1 (default), set to be *messageID+1* (the standard convention).

**recv\_comm** (*expected\_id=None*)

Receive a message.

Return either *CommShort* or *CommData* depending on the message type (fixed length with two parameters, or variable length with associated data). If *expected\_id* is not *None* and the received message ID is different from *expected\_id*, raise an error. For details, see APT communications protocol.

**remap\_axes** (*mapping, accept\_original=True*)

Rename axes to the new labels.

*mapping* is the new axes mapping, which can be a list of new axes name (corresponding to the old axes in order returned by *get\_all\_axes()*), or a dictionary {*alias*: *original*} of the new axes aliases.

**send\_comm** (*messageID, param1=0, param2=0, source=1, dest='host'*)

Send a message with no associated data.

For details, see APT communications protocol.

**send\_comm\_data** (*messageID, data, source=1, dest='host'*)

Send a message with associated data.

For details, see APT communications protocol.

**set\_device\_variable** (*key, value*)

Set the value of a settings parameter

**unlock** ()

Unlock the access to the device from other threads/processes (isn't necessarily implemented)

```
class pylablib.devices.Thorlabs.kinesis.TFlipperParameters (transit_time,  
                                                         io1_oper_mode,  
                                                         io1_sig_mode,  
                                                         io1_pulse_width,  
                                                         io2_oper_mode,  
                                                         io2_sig_mode,  
                                                         io2_pulse_width)
```

Bases: *tuple*

**count** ()

Return number of occurrences of value.

**index** ()

Return first index of value.

Raises *ValueError* if the value is not present.

**io1\_oper\_mode**

**io1\_pulse\_width**

**io1\_sig\_mode**

**io2\_oper\_mode**

**io2\_pulse\_width**

**io2\_sig\_mode**

**transit\_time**

```
class pylablib.devices.Thorlabs.kinesis.MFF (conn)
```

Bases: *pylablib.devices.Thorlabs.kinesis.KinesisDevice*

MFF (Motorized Filter Flip Mount) device.

Implements FTDI chip connectivity via pyft232 (virtual serial interface).

**Parameters** `conn` – serial connection parameters (usually 8-digit device serial number).

**get\_status\_n** (*channel=None*)

Get numerical status of the device.

For details, see APT communications protocol.

**get\_status** (*channel=None*)

Get device status.

Return list of status strings, which can include "sw\_fw\_lim" (forward limit switch reached), "sw\_bk\_lim" (backward limit switch reached), "moving\_fw" (moving forward), "moving\_bk" (moving backward), "jogging\_fw" (jogging forward), "jogging\_bk" (jogging backward), "homing" (homing), "homed" (homing done), "tracking", "settled", "motion\_error" (excessive position error), "current\_limit" (motor current limit exceeded), or "enabled" (motor is enabled).

**wait\_for\_status** (*status, enabled, channel=None, timeout=None, period=0.05*)

Wait until the given status (or list of status bits) is in the desired state.

*status* is a string or a list of strings describing the status bits to monitor; for possible values, see [get\\_status\(\)](#). If *enabled==True*, wait until one of the given statuses is enabled; otherwise, wait until all given statuses are disabled. *period* specifies status checking period (in s).

**move\_to\_state** (*state, channel=None*)

Move to the given flip mount state (either 0 or 1)

**get\_state** (*channel=None*)

Get the flip mount state (either 0 or 1).

Return None if the mount is current moving (i.e., the state os undefined)

**get\_flipper\_parameters** (*channel=None*)

Get current flipper parameters (*transit\_time, io1\_oper\_mode, io1\_sig\_mode, io1\_pulse\_width, io2\_oper\_mode, io2\_sig\_mode, io2\_pulse\_width*)

*transit\_time* specifies transit time (in seconds between 0.3 and 2.8); *io\*\_oper\_mode* specifies operation mode (in vs. out and position vs. motion input/indication), *io\*\_sig\_mode* specifies signal mode (button input, voltage edge input, edge output or pulse output). *io\*\_pulse\_width* specifies output pulse width if the corresponding output mode is selected. For detailed mode description, see the flip mirror or APT manual.

**setup\_flipper** (*transit\_time=None, io1\_oper\_mode=None, io1\_sig\_mode=None, io1\_pulse\_width=None, io2\_oper\_mode=None, io2\_sig\_mode=None, io2\_pulse\_width=None, channel=None*)

Set flipper parameters.

*transit\_time* specifies transit time (in seconds between 0.3 and 2.8); *io\*\_oper\_mode* specifies operation mode (in vs. out and position vs. motion input/indication), *io\*\_sig\_mode* specifies signal mode (button input, voltage edge input, edge output or pulse output). *io\*\_pulse\_width* specifies output pulse width if the corresponding output mode is selected. If any parameter is None, use the current value. For detailed mode description, see the flip mirror or APT manual.

**class CommData** (*messageID, data, source, dest*)

Bases: `tuple`

**count** ()

Return number of occurrences of value.

**data**

**dest**

**index()**  
Return first index of value.  
Raises ValueError if the value is not present.

**messageID**

**source**

**class CommShort** (*messageID, param1, param2, source, dest*)  
Bases: `tuple`

**count()**  
Return number of occurrences of value.

**dest**

**index()**  
Return first index of value.  
Raises ValueError if the value is not present.

**messageID**

**param1**

**param2**

**source**

**Error**  
alias of `pylablib.devices.Thorlabs.base.ThorlabsError`

**class NoParameterCaller** (*device, kind*)  
Bases: `object`  
Class to simplify calling functions without a parameter

**add\_background\_comm** (*messageID*)  
Mark given messageID as a ‘background’ message, which can be sent at any point without prompt (e.g., some operation confirmation).  
  
If it is received instead during `recv_comm_` operations, it is ignored, and the corresponding counter is increased.

**apply\_settings** (*settings*)  
Apply the settings.  
  
*settings* is the dict {`name`: `value`} of the device available settings. Non-applicable settings are ignored.

**blink** (*channel=1, dest='host'*)  
Identify the physical device (by, e.g., blinking status LED or screen)

**check\_background\_comm** (*messageID*)  
Return message counter and the last message value (`None` if not message received yet) of a given ‘background’ message

**close()**  
Close the backend

**get\_all\_axes()**  
Get the list of all available axes (taking mapping into account)

**get\_all\_channels()**

Get the list of all available channels; alias of `get_all_axes` method

**get\_device\_info** (*dest='host'*)

Get device info.

Return tuple (serial\_no, model\_no, fw\_ver, hw\_type, hw\_ver, mod\_state, nchannels, notes).

**get\_device\_variable** (*key*)

Get the value of a settings, status, or full info parameter

**get\_full\_info** (*include=0*)

Get dict {name: value} containing full device information (including status and settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_full\_status** (*include=0*)

Get dict {name: value} containing the device status (including settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_number\_of\_channels()**

Get number of channels on the device

**get\_settings** (*include=0*)

Get dict {name: value} containing all the device settings.

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**is\_opened()**

Check if the device is connected

**static list\_devices** (*filter\_ids=True*)

List all connected devices.

Return list of tuples (conn, description). If *filter\_ids==True*, only leave devices with Thorlabs-like IDs (8-digit numbers). Otherwise, show all devices (some of them might not be Thorlabs-related).

**lock** (*timeout=None*)

Lock the access to the device from other threads/processes (isn't necessarily implemented)

**locking** (*timeout=None*)

Context manager for lock & unlock

**open()**

Open the backend

**query** (*messageID, param1=0, param2=0, source=1, dest='host', replyID=-1*)

Send a query to the device and receive the reply.

A combination of `send_comm()` and `recv_comm()`. If *replyID* is not None, specifies the expected reply message ID; if -1 (default), set to be *messageID+1* (the standard convention).

**recv\_comm** (*expected\_id=None*)

Receive a message.

Return either `CommShort` or `CommData` depending on the message type (fixed length with two parameters, or variable length with associated data). If *expected\_id* is not `None` and the received message ID is different from *expected\_id*, raise an error. For details, see APT communications protocol.

**remap\_axes** (*mapping, accept\_original=True*)

Rename axes to the new labels.

*mapping* is the new axes mapping, which can be a list of new axes name (corresponding to the old axes in order returned by `get_all_axes()`), or a dictionary {*alias*: *original*} of the new axes aliases.

**send\_comm** (*messageID, param1=0, param2=0, source=1, dest='host'*)

Send a message with no associated data.

For details, see APT communications protocol.

**send\_comm\_data** (*messageID, data, source=1, dest='host'*)

Send a message with associated data.

For details, see APT communications protocol.

**set\_default\_channel** (*channel*)

Set the default channel for all channel-related methods

**set\_device\_variable** (*key, value*)

Set the value of a settings parameter

**status\_bits** = [(1, 'sw\_bk\_lim'), (2, 'sw\_fw\_lim'), (16, 'moving\_bk'), (32, 'moving\_fw']

**unlock** ()

Unlock the access to the device from other threads/processes (isn't necessarily implemented)

**using\_channel** (*channel*)

Context manager for temporarily using a different default channel

```
class pylablib.devices.Thorlabs.kinesis.KinesisMotor(conn, scale='step',  
                                                    default_channel=1,  
                                                    is_rack_system=False)
```

Bases: `pylablib.devices.Thorlabs.kinesis.KinesisDevice`

Thorlabs motor controller.

Implements FTDI chip connectivity via pyft232 (virtual serial interface).

The physical units are encoder steps for position (ratio to m or degrees depends on the connected stage), steps/sec for velocity, and steps/sec<sup>2</sup> for acceleration.

#### Parameters

- **conn** (*str*) – serial connection parameters (usually an 8-digit device serial number).
- **scale** – scale of the position, velocity, and acceleration units to the internals units; can be "stage" (attempt to autodetect motor and stage parameters), a string with the name of the stage, e.g., "MTS50-Z8" or "DDR100" (use the stage name to extract the scale; determine velocity and acceleration from this scale and the motor model), "step" (use encoder/motor steps as units; determine velocity and acceleration from this scale and the motor model), a single number (use this as the ratio of internal steps to physical units; determine velocity and acceleration from this scale and the motor model), or a 3-tuple of numbers (*position\_scale*, *velocity\_scale*, *acceleration\_scale*) which gives the ratio of internal units to physical units (useful for new or unrecognized controllers or stages, as no autodetection is required); in the

case of unrecognized devices, use internal units (same as setting `scale=(1, 1, 1)`); if the scale can't be autodetected, it can be obtained from the APT manual knowing the device and the stage model

- **default\_channel** – starting default channel used when no channel is supplied to a channel-level command (such as `move_to` or `get_position`).
- **is\_rack\_system** – specify whether the device is a rack system or a standalone USB device (default mode).

#### **get\_scale()**

Get the scaling coefficients.

Return a tuple (`position_scale`, `velocity_scale`, `acceleration_scale`) for scaling of the physical units in terms of internal units. To get the coefficients source and physical units, use `get_scale_units()`.

#### **get\_scale\_units()**

Get units used for calculating scaled position, velocity and acceleration values.

Can be "deg" (autodetected rotational stage: deg, deg/s and deg/s<sup>2</sup>), "m" (autodetected translational stage: m, m/sec and m/sec<sup>2</sup>), "step" (autodetected driver but not detected step scale: steps, steps/sec and steps/sec<sup>2</sup>) "user\_step" (autodetected driver and user supplied step scale: user-supplied step scale for position, same units per sec or sec<sup>2</sup> for velocity and acceleration), "user" (all three scales are supplied by user), or "internal" (no scales are supplied or detected, use device internal units)

#### **get\_stage()**

Return the name of the stage which was supplied by the user or autodetected.

If the stage is unknown, return None

#### **set\_supported\_channels(channels=1)**

Set the channels in the device.

Can be either a list of channels, a single number defining the number of channels numbered from 1 to `channels` (inclusive).

#### **get\_status\_n(channel=None)**

Get numerical status of the device.

For details, see APT communications protocol.

#### **get\_status(channel=None)**

Get device status.

Return list of status strings, which can include "sw\_fw\_lim" (forward limit switch reached), "sw\_bk\_lim" (backward limit switch reached), "moving\_fw" (moving forward), "moving\_bk" (moving backward), "jogging\_fw" (jogging forward), "jogging\_bk" (jogging backward), "homing" (homing), "homed" (homing done), "tracking", "settled", "motion\_error" (excessive position error), "current\_limit" (motor current limit exceeded), or "enabled" (motor is enabled).

#### **wait\_for\_status(status, enabled, channel=None, timeout=None, period=0.05)**

Wait until the given status (or list of status bits) is in the desired state.

`status` is a string or a list of strings describing the status bits to monitor; for possible values, see `get_status()`. If `enabled==True`, wait until one of the given statuses is enabled; otherwise, wait until all given statuses are disabled. `period` specifies status checking period (in s).

#### **home(sync=True, force=False, channel=None, timeout=None)**

Home the device.

If `sync==True`, wait until homing is done (with the given timeout). If `force==False`, only home if the device isn't homed already.

**is\_homing** (*channel=None*)

Check if homing is in progress

**is\_homed** (*channel=None*)

Check if the device is homed

**wait\_for\_home** (*channel=None, timeout=None*)

Wait until the device is homed

**get\_position** (*channel=None, scale=True*)

Get current position.

If `scale==True`, return value in the physical units (see class description); otherwise, return it in the device internal units (steps).

**set\_position\_reference** (*position=0, channel=None, scale=True*)

Set position reference (actual motor position stays the same).

If `scale==True`, assume that the position is in the physical units (see class description); otherwise, assume it is in the device internal units (steps).

**move\_by** (*distance=1, channel=None, scale=True*)

Move by a given amount (positive or negative) from the current position.

If `scale==True`, assume that the distance is in the physical units (see class description); otherwise, assume it is in the device internal units (steps).

**move\_to** (*position, channel=None, scale=True*)

Move to *position* (positive or negative).

If `scale==True`, assume that the position is in the physical units (see class description); otherwise, assume it is in the device internal units (steps).

**jog** (*direction, channel=None, kind='continuous'*)

Jog in the given direction ("**+**" or "**-**").

If `kind=="continuous"`, simply start motion in the given direction at the maximal speed until either the motor is stopped explicitly, or the limit is reached (this uses `MOT_MOVE_VELOCITY` command). If `kind=="builtin"`, use the built-in `MOT_MOVE_JOG` command, whose parameters are specified by `get_jog_parameters()`.

**is\_moving** (*channel=None*)

Check if motion is in progress

**wait\_move** (*channel=None, timeout=None*)

Wait until motion command is done

**stop** (*immediate=False, sync=True, channel=None, timeout=None*)

Stop the motion.

If `immediate==True` make an abrupt stop; otherwise, slow down gradually. If `sync==True`, wait until the motion is stopped.

**wait\_for\_stop** (*channel=None, timeout=None*)

Wait until motion or homing is done

**get\_velocity\_parameters** (*channel=None, scale=True*)

Get current velocity parameters (*min\_velocity*, *acceleration*, *max\_velocity*)

If `scale==True`, return values in the physical units (see class description); otherwise, return it in the device internal units.



**setup\_velocity** (*min\_velocity=None, acceleration=None, max\_velocity=None, channel=None, scale=True*)  
Set velocity parameters.

If any parameter is *None*, use the current value. If *scale==True*, assume that the specified values are in the physical units (see class description); otherwise, assume it is in the device internal units.

**get\_jog\_parameters** (*channel=None, scale=True*)  
Get current jog parameters (*mode, step\_size, min\_velocity, acceleration, max\_velocity, stop\_mode*)

If *scale==True*, return values in the physical units (see class description); otherwise, return it in the device internal units.

**setup\_jog** (*mode=None, step\_size=None, min\_velocity=None, acceleration=None, max\_velocity=None, stop\_mode=None, channel=None, scale=True*)  
Set jog parameters.

If any parameter is *None*, use the current value. If *scale==True*, assume that the specified values are in the physical units (see class description); otherwise, assume it is in the device internal units.

**get\_homing\_parameters** (*channel=None, scale=True*)  
Get current homing parameters (*home\_direction, limit\_switch, velocity, offset\_distance*)

If *scale==True*, return values are in the physical units (see class description); otherwise, return it in the device internal units.

**setup\_homing** (*home\_direction=None, limit\_switch=None, velocity=None, offset\_distance=None, channel=None, scale=True*)  
Set homing parameters.

If any parameter is *None*, use the current value. If *scale==True*, assume that the specified values are in the physical units (see class description); otherwise, assume it is in the device internal units.

**get\_gen\_move\_parameters** (*channel=None, scale=True*)  
Get general move parameters (*backlash\_distance*)

If *scale==True*, return values in the physical units (see class description); otherwise, return it in the device internal units.

**setup\_gen\_move** (*backlash\_distance=None, channel=None, scale=True*)  
Set jog parameters.

If any parameter is *None*, use the current value. If *scale==True*, assume that the specified value is in the physical units (see class description); otherwise, assume it is in the device internal units.

**get\_limit\_switch\_parameters** (*channel=None, scale=True*)  
Get current limit switch parameters (*hw\_kind\_cw, hw\_kind\_ccw, hw\_flipped, sw\_position\_cw, sw\_position\_ccw, sw\_kind*)

If *scale==True*, return values in the physical units (see class description); otherwise, return it in the device internal units (steps).

**setup\_limit\_switch** (*hw\_kind\_cw=None, hw\_kind\_ccw=None, hw\_swapped=None, sw\_position\_cw=None, sw\_position\_ccw=None, sw\_kind=None, channel=None, scale=True*)  
Set limit switch parameters.

If any parameter is *None*, use the current value. If *scale==True*, assume that the specified values are in the physical units (see class description); otherwise, assume it is in the device internal units (Steps).

**get\_kcube\_trigio\_parameters** ()  
Get KCube trigger IO parameters (*trig1\_pol, trig1\_pol, trig2\_mode, trig2\_pol*)

**setup\_kcube\_trigio** (*trig1\_mode=None, trig1\_pol=None, trig2\_mode=None, trig2\_pol=None*)

Set KCube trigger IO parameters.

If any parameter is `None`, use the current value.

**get\_kcube\_trigpos\_parameters** (*scale=True*)

Get KCube trigger position parameters (*start\_fw, step\_fw, num\_fw, start\_bk, step\_bk, num\_bk, width, ncycles*).

If *scale==True*, return positions and steps in the physical units (see class description); otherwise, return it in the device internal units (steps). Pulse width is always defined in seconds.

**setup\_kcube\_trigpos** (*start\_fw=None, step\_fw=None, num\_fw=None, start\_bk=None, step\_bk=None, num\_bk=None, width=None, ncycles=None, scale=True*)

Set KCube trigger IO parameters.

If any parameter is `None`, use the current value.

If *scale==True*, return positions and steps in the physical units (see class description); otherwise, return it in the device internal units (steps). Pulse width is always defined in seconds.

**class CommData** (*messageID, data, source, dest*)

Bases: `tuple`

**count** ()

Return number of occurrences of value.

**data**

**dest**

**index** ()

Return first index of value.

Raises `ValueError` if the value is not present.

**messageID**

**source**

**class CommShort** (*messageID, param1, param2, source, dest*)

Bases: `tuple`

**count** ()

Return number of occurrences of value.

**dest**

**index** ()

Return first index of value.

Raises `ValueError` if the value is not present.

**messageID**

**param1**

**param2**

**source**

**Error**

alias of `pylablib.devices.Thorlabs.base.ThorlabsError`

**class NoParameterCaller** (*device, kind*)

Bases: `object`

Class to simplify calling functions without a parameter

**add\_background\_comm** (*messageID*)

Mark given messageID as a 'background' message, which can be sent at any point without prompt (e.g., some operation confirmation).

If it is received instead during `recv_comm_` operations, it is ignored, and the corresponding counter is increased.

**apply\_settings** (*settings*)

Apply the settings.

*settings* is the dict {name: value} of the device available settings. Non-applicable settings are ignored.

**blink** (*channel=1, dest='host'*)

Identify the physical device (by, e.g., blinking status LED or screen)

**check\_background\_comm** (*messageID*)

Return message counter and the last message value (None if not message received yet) of a given 'background' message

**close** ()

Close the backend

**get\_all\_axes** ()

Get the list of all available axes (taking mapping into account)

**get\_all\_channels** ()

Get the list of all available channels; alias of `get_all_axes` method

**get\_device\_info** (*dest='host'*)

Get device info.

Return tuple (serial\_no, model\_no, fw\_ver, hw\_type, hw\_ver, mod\_state, nchannels, notes).

**get\_device\_variable** (*key*)

Get the value of a settings, status, or full info parameter

**get\_full\_info** (*include=0*)

Get dict {name: value} containing full device information (including status and settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_full\_status** (*include=0*)

Get dict {name: value} containing the device status (including settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_number\_of\_channels** ()

Get number of channels on the device

**get\_settings** (*include=0*)

Get dict {name: value} containing all the device settings.

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting `include=-10` queries all available variables, which is equivalent to `include="all"`.

**is\_opened()**

Check if the device is connected

**static list\_devices** (*filter\_ids=True*)

List all connected devices.

Return list of tuples (*conn*, *description*). If *filter\_ids==True*, only leave devices with Thorlabs-like IDs (8-digit numbers). Otherwise, show all devices (some of them might not be Thorlabs-related).

**lock** (*timeout=None*)

Lock the access to the device from other threads/processes (isn't necessarily implemented)

**locking** (*timeout=None*)

Context manager for lock & unlock

**open()**

Open the backend

**query** (*messageID*, *param1=0*, *param2=0*, *source=1*, *dest='host'*, *replyID=-1*)

Send a query to the device and receive the reply.

A combination of `send_comm()` and `recv_comm()`. If *replyID* is not `None`, specifies the expected reply message ID; if -1 (default), set to be *messageID*+1 (the standard convention).

**recv\_comm** (*expected\_id=None*)

Receive a message.

Return either `CommShort` or `CommData` depending on the message type (fixed length with two parameters, or variable length with associated data). If *expected\_id* is not `None` and the received message ID is different from *expected\_id*, raise an error. For details, see APT communications protocol.

**remap\_axes** (*mapping*, *accept\_original=True*)

Rename axes to the new labels.

*mapping* is the new axes mapping, which can be a list of new axes name (corresponding to the old axes in order returned by `get_all_axes()`), or a dictionary {*alias*: *original*} of the new axes aliases.

**send\_comm** (*messageID*, *param1=0*, *param2=0*, *source=1*, *dest='host'*)

Send a message with no associated data.

For details, see APT communications protocol.

**send\_comm\_data** (*messageID*, *data*, *source=1*, *dest='host'*)

Send a message with associated data.

For details, see APT communications protocol.

**set\_default\_channel** (*channel*)

Set the default channel for all channel-related methods

**set\_device\_variable** (*key*, *value*)

Set the value of a settings parameter

**status\_bits** = [(1, 'sw\_bk\_lim'), (2, 'sw\_fw\_lim'), (16, 'moving\_bk'), (32, 'moving\_fw')]

**unlock()**

Unlock the access to the device from other threads/processes (isn't necessarily implemented)

**using\_channel** (*channel*)

Context manager for temporarily using a different default channel

**class** `pylablib.devices.Thorlabs.kinesis.KinesisPiezoMotor` (*conn*, *default\_channel*, *fault\_channel=1*)

Bases: `pylablib.devices.Thorlabs.kinesis.KinesisDevice`

Thorlabs piezo motor (TIM/KIM series) controller.

Implements FTDI chip connectivity via pyft232 (virtual serial interface).

**Parameters** *conn* (*str*) – serial connection parameters (usually an 8-digit device serial number).

**get\_enabled\_channels** ()

Check which specific piezo motor channels are enabled.

Can be `None` (none enabled), or a tuple with either one or two channels: `(1,)` to `(4,)`, `(1,2)` or `(3,4)`.

**enable\_channels** (*channel*)

Enable specific piezo motor channel.

Can be `None` (none enabled), and integer 1 to 4, or a tuple `(1,2)` or `(3,4)` (enable 2 channel simultaneously).

**get\_status\_n** (*channel=None*)

Get numerical status of the piezo motor.

For details, see APT communications protocol.

**get\_status** (*channel=None*)

Get piezo motor status.

Return list of status strings, which can include "sw\_fw\_lim" (forward limit switch reached), "sw\_bk\_lim" (backward limit switch reached), "moving\_fw" (moving forward), "moving\_bk" (moving backward), "jogging\_fw" (jogging forward), "jogging\_bk" (jogging backward), "homing" (homing), "homed" (homing done), "tracking", "settled", "motion\_error" (excessive position error), "current\_limit" (motor current limit exceeded), or "enabled" (motor is enabled).

**wait\_for\_status** (*status*, *enabled*, *timeout=None*, *period=0.05*, *channel=None*)

**get\_position** (*channel=None*)

Get current piezo-motor position

**set\_position\_reference** (*position=0*, *channel=None*)

Set piezo-motor position reference (actual position stays the same)

**move\_by** (*distance=1*, *auto\_enable=True*, *channel=None*)

Move piezo-motor by a given *distance* (positive or negative)

**move\_to** (*position*, *auto\_enable=True*, *channel=None*)

Move piezo-motor to *position* (positive or negative)

**jog** (*direction*, *kind='continuous'*, *auto\_enable=True*, *channel=None*)

Jog piezo motor in the given direction ("`+`" or "`-`").

If *kind*=="continuous", simply start motion in the given direction at the standard jog speed until either the motor is stopped explicitly, or the limit is reached. If *kind*=="builtin", use the built-in jog command, whose parameters are specified by `get_jog_parameters()`. Note that *kind*=="continuous" is still implemented through the builtin jog, so it changes its parameters; hence, afterwards the jog parameters need to be manually restored.

**is\_moving** (*channel=None*)

Check if motion is in progress

**wait\_move** (*channel=None, timeout=None*)

Wait until motion command is done

**stop** (*channel=None, sync=True*)

Stop the piezo motor motion

**get\_drive\_parameters** (*channel=None*)

Get current piezo-motor drive parameters (*max\_voltage, velocity, acceleration*)

Voltage is defined in volts, velocity in steps/s, and acceleration in steps/s<sup>2</sup>.

**setup\_drive** (*max\_voltage=None, velocity=None, acceleration=None, channel=None*)

Set piezo-motor drive parameters.

If any parameter is *None*, use the current value. Voltage is defined in volts, velocity in steps/s, and acceleration in steps/s<sup>2</sup>.

**get\_jog\_parameters** (*channel=None*)

Get current piezo-motor jog parameters (*mode, step\_size\_fw, step\_size\_bk, velocity, acceleration*)

Step size is defined in piezo steps, velocity in steps/s, and acceleration in step/s<sup>2</sup>.

**setup\_jog** (*mode=None, step\_size\_fw=None, step\_size\_bk=None, velocity=None, acceleration=None, channel=None*)

Set piezo-motor jog parameters.

If any parameter is *None*, use the current value. Step size is defined in piezo steps, velocity in steps/s, and acceleration in step/s<sup>2</sup>. TIM-series controllers do not support separate step size, so *step\_size\_bk* is ignored for them.

**class CommData** (*messageID, data, source, dest*)

Bases: `tuple`

**count** ()

Return number of occurrences of value.

**data**

**dest**

**index** ()

Return first index of value.

Raises `ValueError` if the value is not present.

**messageID**

**source**

**class CommShort** (*messageID, param1, param2, source, dest*)

Bases: `tuple`

**count** ()

Return number of occurrences of value.

**dest**

**index** ()

Return first index of value.

Raises `ValueError` if the value is not present.

**messageID**

**param1**

**param2**

**source**

#### **Error**

alias of `pylablib.devices.Thorlabs.base.ThorlabsError`

**class NoParameterCaller** (*device, kind*)

Bases: `object`

Class to simplify calling functions without a parameter

**add\_background\_comm** (*messageID*)

Mark given messageID as a ‘background’ message, which can be sent at any point without prompt (e.g., some operation confirmation).

If it is received instead during `recv_comm` operations, it is ignored, and the corresponding counter is increased.

**apply\_settings** (*settings*)

Apply the settings.

*settings* is the dict {name: value} of the device available settings. Non-applicable settings are ignored.

**blink** (*channel=1, dest='host'*)

Identify the physical device (by, e.g., blinking status LED or screen)

**check\_background\_comm** (*messageID*)

Return message counter and the last message value (None if not message received yet) of a given ‘background’ message

**close** ()

Close the backend

**get\_all\_axes** ()

Get the list of all available axes (taking mapping into account)

**get\_all\_channels** ()

Get the list of all available channels; alias of `get_all_axes` method

**get\_device\_info** (*dest='host'*)

Get device info.

Return tuple (serial\_no, model\_no, fw\_ver, hw\_type, hw\_ver, mod\_state, nchannels, notes).

**get\_device\_variable** (*key*)

Get the value of a settings, status, or full info parameter

**get\_full\_info** (*include=0*)

Get dict {name: value} containing full device information (including status and settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting `include=-10` queries all available variables, which is equivalent to `include="all"`.

**get\_full\_status** (*include=0*)

Get dict {name: value} containing the device status (including settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting `include=-10` queries all available variables, which is equivalent to `include="all"`.

**get\_number\_of\_channels()**

Get number of channels on the device

**get\_settings(include=0)**

Get dict {name: value} containing all the device settings.

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting `include=-10` queries all available variables, which is equivalent to `include="all"`.

**is\_opened()**

Check if the device is connected

**static list\_devices(filter\_ids=True)**

List all connected devices.

Return list of tuples (conn, description). If `filter_ids==True`, only leave devices with Thorlabs-like IDs (8-digit numbers). Otherwise, show all devices (some of them might not be Thorlabs-related).

**lock(timeout=None)**

Lock the access to the device from other threads/processes (isn't necessarily implemented)

**locking(timeout=None)**

Context manager for lock & unlock

**open()**

Open the backend

**query(messageID, param1=0, param2=0, source=1, dest='host', replyID=-1)**

Send a query to the device and receive the reply.

A combination of `send_comm()` and `recv_comm()`. If `replyID` is not `None`, specifies the expected reply message ID; if -1 (default), set to be `messageID+1` (the standard convention).

**recv\_comm(expected\_id=None)**

Receive a message.

Return either `CommShort` or `CommData` depending on the message type (fixed length with two parameters, or variable length with associated data). If `expected_id` is not `None` and the received message ID is different from `expected_id`, raise an error. For details, see APT communications protocol.

**remap\_axes(mapping, accept\_original=True)**

Rename axes to the new labels.

*mapping* is the new axes mapping, which can be a list of new axes name (corresponding to the old axes in order returned by `get_all_axes()`), or a dictionary {alias: original} of the new axes aliases.

**send\_comm(messageID, param1=0, param2=0, source=1, dest='host')**

Send a message with no associated data.

For details, see APT communications protocol.

**send\_comm\_data(messageID, data, source=1, dest='host')**

Send a message with associated data.



For details, see APT communications protocol.

**set\_default\_channel** (*channel*)

Set the default channel for all channel-related methods

**set\_device\_variable** (*key, value*)

Set the value of a settings parameter

**status\_bits** = [(1, 'sw\_bk\_lim'), (2, 'sw\_fw\_lim'), (16, 'moving\_bk'), (32, 'moving\_fw')]

**unlock** ()

Unlock the access to the device from other threads/processes (isn't necessarily implemented)

**using\_channel** (*channel*)

Context manager for temporarily using a different default channel

**class** pylablib.devices.Thorlabs.kinesis.TQuadDetectorPIDParams (*p, i, d*)

Bases: `tuple`

**count** ()

Return number of occurrences of value.

**d**

**i**

**index** ()

Return first index of value.

Raises ValueError if the value is not present.

**p**

**class** pylablib.devices.Thorlabs.kinesis.TQuadDetectorSetpoint (*xpos, ypos*)

Bases: `tuple`

**count** ()

Return number of occurrences of value.

**index** ()

Return first index of value.

Raises ValueError if the value is not present.

**xpos**

**ypos**

**class** pylablib.devices.Thorlabs.kinesis.TQuadDetectorReadings (*xdiff, ydiff, sum,*

*xpos, ypos*)

Bases: `tuple`

**count** ()

Return number of occurrences of value.

**index** ()

Return first index of value.

Raises ValueError if the value is not present.

**sum**

**xdiff**

**xpos**

**ydiff**

**ypos**

```
class pylablib.devices.Thorlabs.kinesis.TQuadDetectorOutputParams (xmin, xmax,  
                                                                    ymin, ymax,  
                                                                    xgain,  
                                                                    ygain,  
                                                                    route,  
                                                                    open_loop_out)
```

Bases: `tuple`

**count** ()  
Return number of occurrences of value.

**index** ()  
Return first index of value.  
  
Raises ValueError if the value is not present.

**open\_loop\_out**

**route**

**xgain**

**xmax**

**xmin**

**ygain**

**ymax**

**ymin**

```
class pylablib.devices.Thorlabs.kinesis.KinesisQuadDetector (conn, timeout=3.0)  
Bases: pylablib.devices.Thorlabs.kinesis.BasicKinesisDevice
```

Kinesis quadrature detectors: KPA101, TPA101, TQD001.

Implements FTDI chip connectivity via pyft232 (virtual serial interface).

**Parameters** **conn** (*str*) – serial connection parameters (usually an 8-digit device serial number).

**get\_pid\_parameters** ()  
Get current PID gain parameters (*p, i, d*)

**set\_pid\_parameters** (*p=None, i=None, d=None*)  
Set current PID gain parameters (*p, i, d*).  
  
If any parameter is *None*, use the current value.

**get\_manual\_output** ()  
Get current manual output values (*xpos, ypos*) (used in open loop mode)

**set\_manual\_output** (*xpos=None, ypos=None*)  
Set current manual output values (used in open loop mode).  
  
If any parameter is *None*, use the current value.

**get\_readings** ()  
Get current readings (*xdiff, ydiff, sum, xpos, ypos*)

**get\_operation\_mode** ()  
Get current operation mode: "monitor", "open\_loop", "closed\_loop", or "auto\_loop"

**set\_operation\_mode** (*mode*)  
Set current operation mode: "monitor", "open\_loop", "closed\_loop", or "auto\_loop"

```

get_output_parameters()
    Get current output parameters (xmin, xmax, ymin, ymax, xgain, ygain, route,
    open_loop_out)

set_output_parameters(xmin=None, xmax=None, ymin=None, ymax=None, xgain=None,
    ygain=None, route=None, open_loop_out=None)
    Set current PID gain parameters (xmin, xmax, ymin, ymax, xgain, ygain, route,
    open_loop_out).

    xmin, xmax, ymin, and ymax specify output limits, xgain and ygain specify additional separate gain (be-
    tween -1 and 1), route sets where output is routed in the closed loop mode (either "sma_only" or
    "sma_hub"), open_loop_out specifies the output source in the open loop mode (either "zero" or
    "fixed"). If any parameter is None, use the current value.

class CommData(messageID, data, source, dest)
    Bases: tuple

    count()
        Return number of occurrences of value.

    data

    dest

    index()
        Return first index of value.

        Raises ValueError if the value is not present.

    messageID

    source

class CommShort(messageID, param1, param2, source, dest)
    Bases: tuple

    count()
        Return number of occurrences of value.

    dest

    index()
        Return first index of value.

        Raises ValueError if the value is not present.

    messageID

    param1

    param2

    source

Error
    alias of pylablib.devices.Thorlabs.base.ThorlabsError

class NoParameterCaller(device, kind)
    Bases: object

    Class to simplify calling functions without a parameter

add_background_comm(messageID)
    Mark given messageID as a ‘background’ message, which can be sent at any point without prompt (e.g.,
    some operation confirmation).

```

If it is received instead during `recv_comm_` operations, it is ignored, and the corresponding counter is increased.

**apply\_settings** (*settings*)

Apply the settings.

*settings* is the dict {name: value} of the device available settings. Non-applicable settings are ignored.

**blink** (*channel=1, dest='host'*)

Identify the physical device (by, e.g., blinking status LED or screen)

**check\_background\_comm** (*messageID*)

Return message counter and the last message value (None if not message received yet) of a given 'background' message

**close** ()

Close the backend

**get\_device\_info** (*dest='host'*)

Get device info.

Return tuple (serial\_no, model\_no, fw\_ver, hw\_type, hw\_ver, mod\_state, nchannels, notes).

**get\_device\_variable** (*key*)

Get the value of a settings, status, or full info parameter

**get\_full\_info** (*include=0*)

Get dict {name: value} containing full device information (including status and settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_full\_status** (*include=0*)

Get dict {name: value} containing the device status (including settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_number\_of\_channels** ()

Get number of channels on the device

**get\_settings** (*include=0*)

Get dict {name: value} containing all the device settings.

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**is\_opened** ()

Check if the device is connected

**static list\_devices** (*filter\_ids=True*)

List all connected devices.

Return list of tuples (*conn*, *description*). If *filter\_ids==True*, only leave devices with Thorlabs-like IDs (8-digit numbers). Otherwise, show all devices (some of them might not be Thorlabs-related).

**lock** (*timeout=None*)

Lock the access to the device from other threads/processes (isn't necessarily implemented)

**locking** (*timeout=None*)

Context manager for lock & unlock

**open** ()

Open the backend

**query** (*messageID*, *param1=0*, *param2=0*, *source=1*, *dest='host'*, *replyID=-1*)

Send a query to the device and receive the reply.

A combination of `send_comm()` and `recv_comm()`. If *replyID* is not `None`, specifies the expected reply message ID; if -1 (default), set to be *messageID*+1 (the standard convention).

**recv\_comm** (*expected\_id=None*)

Receive a message.

Return either `CommShort` or `CommData` depending on the message type (fixed length with two parameters, or variable length with associated data). If *expected\_id* is not `None` and the received message ID is different from *expected\_id*, raise an error. For details, see APT communications protocol.

**send\_comm** (*messageID*, *param1=0*, *param2=0*, *source=1*, *dest='host'*)

Send a message with no associated data.

For details, see APT communications protocol.

**send\_comm\_data** (*messageID*, *data*, *source=1*, *dest='host'*)

Send a message with associated data.

For details, see APT communications protocol.

**set\_device\_variable** (*key*, *value*)

Set the value of a settings parameter

**unlock** ()

Unlock the access to the device from other threads/processes (isn't necessarily implemented)

## pylablib.devices.Thorlabs.misc module

## pylablib.devices.Thorlabs.serial module

**class** `pylablib.devices.Thorlabs.serial.ThorlabsSerialInterface` (*conn*)

Bases: `pylablib.core.devio.SCPI.SCPIDevice`

Generic Thorlabs device interface using Serial communication.

**Parameters** *conn* – serial connection parameters (usually port or a tuple containing port and baudrate)

**Error**

alias of `pylablib.devices.Thorlabs.base.ThorlabsError`

**ReraiseError**

alias of `pylablib.devices.Thorlabs.base.ThorlabsBackendError`

**open** ()

Open the backend

**BackendError**

alias of `pylablib.core.devio.comm_backend.DeviceBackendError`

**class NoParameterCaller** (*device, kind*)

Bases: `object`

Class to simplify calling functions without a parameter

**apply\_settings** (*settings*)

Apply the settings.

*settings* is a dict {name: value} of the available device settings. Non-applicable settings are ignored.

**ask** (*msg, data\_type='string', delay=0.0, timeout=None, read\_echo=False*)

Write a message and read a reply.

*msg* is the query message, *delay* is the delay between write and read. Other parameters are the same as in `read()`. If `read_echo==True`, assume that the device first echoes the input and skip it.

**close** ()

Close the backend

**flush** (*one\_line=False*)

Flush the read buffer (read all the available data and return the number of bytes read).

If `one_line==True`, read only a single line.

**static get\_arg\_type** (*arg*)

Autodetect argument type

**get\_device\_variable** (*key*)

Get the value of a settings, status, or full info parameter

**get\_esr** (*timeout=None*)

Get the device status register (by default, "\*ESR?" command)

**get\_full\_info** (*include=0*)

Get dict {name: value} containing full device information (including status and settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting `include=-10` queries all available variables, which is equivalent to `include="all"`.

**get\_full\_status** (*include=0*)

Get dict {name: value} containing the device status (including settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting `include=-10` queries all available variables, which is equivalent to `include="all"`.

**get\_id** (*timeout=None*)

Get the device IDN. (query SCPI "\*IDN?" command)

**get\_settings** (*include=0*)

Get dict {name: value} containing all the device settings.

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting `include=-10` queries all available variables, which is equivalent to `include="all"`.

**is\_opened()**

Check if the device is connected

**lock** (*timeout=None*)

Lock the access to the device from other threads/processes (isn't necessarily implemented)

**locking** (*timeout=None*)

Context manager for lock & unlock

**static parse\_array\_data** (*data, fmt, include\_header=False*)

Parse the data returned by the device. *fmt* is *DataFormat* description in numpy format (e.g., "<u2").

If *include\_header==True*, the data is assumed to be in a (somewhat) standard SCPI format: b'#', then a single digit *s* denoting length of the size block, then *s* digits denoting length of the data (in bytes) followed by the actual data. Otherwise (*include\_header==False*), assume that the header is already removed.

**read** (*data\_type='string', timeout=None*)

Read data from the device.

*data\_type* determines the type of the data. Can be 'raw' (just raw data), 'string' (with trailing and leading spaces stripped), 'int', 'float', 'bool' (interprets 0 or 'off' as False, anything else as True), 'value' (returns tuple (value, unit), where *value* is float), a callable (return the result of this callable applied to the string value), a dictionary (return the stored value corresponding to the string value, or to the value converted into integer if the string value is not present), or a list of data types (the result is treated as a list of values with the given types separated by commas). *timeout* overrides the default value.

**read\_binary\_array\_data** (*include\_header=False, timeout=None, flush\_term=True*)

Read a binary data in the from the device.

The data assumes the standard binary transfer header consisting of "#" symbol, then a single digit with the size of the length string, then the length string containing the length of the binary data (in bytes). If *include\_header==True*, return the data with the header; otherwise, return only the content. If *flush\_term==True*, flush the following line to skip terminator characters after the binary data, which are added by some devices. *timeout* overrides the default value.

**reconnect** (*new\_instrument=True, ignore\_error=True*)

Remake the connection.

If *new\_instrument==True*, create a new backend instance. If *ignore\_error==True*, ignore errors on closing.

**reset** ()

Reset the device (by default, "\*RST" command)

**set\_device\_variable** (*key, value*)

Set the value of a settings parameter

**sleep** (*delay*)

Wait for *delay* seconds

**unlock** ()

Unlock the access to the device from other threads/processes (isn't necessarily implemented)

**using\_write\_buffer** ()

Context manager for using a write buffer.

While it's active, all the consecutive *write()* operations are bundled together with ; delimiter. The actual write is performed at the *read()/ask()* operation or at the end of the block.

**wait** (*wait\_type*='sync', *timeout*=None, *wait\_callback*=None)

Pause execution until device overlapped commands are complete.

*wait\_type* is either 'sync' (perform `wait_sync()`), 'dev' (perform `wait_dev()`) or 'none' (do nothing).

**wait\_dev** ()

Pause execution of the device commands until device overlapped commands (e.g., taking sweeps) are complete.

Note that the code execution is not paused.

**wait\_sync** (*timeout*=None, *wait\_callback*=None)

Pause execution of the script until device overlapped commands (e.g., taking sweeps) are complete.

*timeout* and *wait\_callback* override default constructor parameters.

**write** (*msg*, *arg*=None, *arg\_type*=None, *unit*=None, *bool\_selector*=None, *wait\_sync*=None, *read\_echo*=False, *read\_echo\_delay*=0.0)

Send a command.

#### Parameters

- **msg** (*str*) – Text message.
- **arg** – Optional argument to append in the end. If a list of arguments is supplied, the result is joined with ", ".
- **arg\_type** (*str*) – Argument type. Can be 'raw' (in which case data is sent raw), 'string', 'int', 'float', 'bool', a format string (such as '{:.3f}') or a list of argument types (for an iterable argument); if format string is used and the argument is a list or a tuple, then it is expanded as a list of arguments (e.g., `arg_type='{0};{1}'` with `arg=[1, 2]` will produce a string '1; 2'); if a list of types is used, each element of *arg* is converted using the corresponding type, and the result is joined with ", ".
- **unit** (*str*) – If not None, use it as a unit to append after the value.
- **bool\_selector** (*tuple*) – A tuple (*false\_value*, *true\_value*) of two strings to represent bool argument; by default, use `._bool_selector` attribute.
- **wait\_sync** – if True, append the sync command (specified as `._wait_sync_comm` attribute, `"*OPC?"` by default) after the message and pause the execution command is complete; useful in long set operations, where the device might ignore later inputs until the current command is complete; if None, use the class default `._default_write_sync` attribute (False by default).
- **read\_echo** (*bool*) – If True, read a single line after write.
- **read\_echo\_delay** (*float*) – The delay between write and read if `read_echo==True`.

**class** `pylablib.devices.Thorlabs.serial.FW` (*conn*, *respect\_bound*=True)

Bases: `pylablib.devices.Thorlabs.serial.ThorlabsSerialInterface`

Thorlabs FW102/212 motorized filter wheels.

#### Parameters

- **conn** – serial connection parameters (usually port or a tuple containing port and baudrate)
- **respect\_bound** (*bool*) – if True, avoid crossing the boundary between the first and the last position in the wheel



**ask** (*msg*, *data\_type*='string', *delay*=0.0, *timeout*=None, *read\_echo*=False)

Write a message and read a reply.

*msg* is the query message, *delay* is the delay between write and read. Other parameters are the same as in `read()`. If `read_echo==True`, assume that the device first echoes the input and skip it.

**get\_position** ()

Get the wheel position (starting from 1)

**set\_position** (*pos*)

Set the wheel position (starting from 1)

**get\_pcount** ()

Get the number of wheel positions (6 or 12)

**set\_pcount** (*pcount*)

Set the number of wheel positions (6 or 12)

**get\_speed\_mode** ()

Get the motion speed mode ("low" or "high")

**set\_speed\_mode** (*speed\_mode*)

Set the motion speed mode ("low" or "high")

**get\_trigger\_mode** ()

Get the trigger mode ("in" to input external trigger, "out" to output trigger)

**set\_trigger\_mode** (*trigger\_mode*)

Set the trigger mode ("in" to input external trigger, "out" to output trigger)

**get\_sensor\_mode** ()

Get the sensor mode ("off" to turn off when idle to eliminate stray light, "on" to remain on)

**set\_sensor\_mode** (*sensor\_mode*)

Set the sensor mode ("off" to turn off when idle to eliminate stray light, "on" to remain on)

**store\_settings** ()

Store current settings as default

**BackendError**

alias of `pylablib.core.devio.comm_backend.DeviceBackendError`

**Error**

alias of `pylablib.devices.Thorlabs.base.ThorlabsError`

**class NoParameterCaller** (*device*, *kind*)

Bases: `object`

Class to simplify calling functions without a parameter

**ReraiseError**

alias of `pylablib.devices.Thorlabs.base.ThorlabsBackendError`

**apply\_settings** (*settings*)

Apply the settings.

*settings* is a dict {name: value} of the available device settings. Non-applicable settings are ignored.

**close** ()

Close the backend

**flush** (*one\_line*=False)

Flush the read buffer (read all the available data and return the number of bytes read).

If `one_line==True`, read only a single line.

**static get\_arg\_type** (*arg*)

Autodetect argument type

**get\_device\_variable** (*key*)

Get the value of a settings, status, or full info parameter

**get\_esr** (*timeout=None*)

Get the device status register (by default, "\*ESR?" command)

**get\_full\_info** (*include=0*)

Get dict {name: value} containing full device information (including status and settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_full\_status** (*include=0*)

Get dict {name: value} containing the device status (including settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_id** (*timeout=None*)

Get the device IDN. (query SCPI '\*IDN?' command)

**get\_settings** (*include=0*)

Get dict {name: value} containing all the device settings.

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**is\_opened** ()

Check if the device is connected

**lock** (*timeout=None*)

Lock the access to the device from other threads/processes (isn't necessarily implemented)

**locking** (*timeout=None*)

Context manager for lock & unlock

**open** ()

Open the backend

**static parse\_array\_data** (*data, fmt, include\_header=False*)

Parse the data returned by the device. *fmt* is *DataFormat* description in numpy format (e.g., "<u2").

If *include\_header==True*, the data is assumed to be in a (somewhat) standard SCPI format: b'#', then a single digit *s* denoting length of the size block, then *s* digits denoting length of the data (in bytes) followed by the actual data. Otherwise (*include\_header==False*), assume that the header is already removed.

**read** (*data\_type='string', timeout=None*)

Read data from the device.

*data\_type* determines the type of the data. Can be 'raw' (just raw data), 'string' (with trailing and leading spaces stripped), 'int', 'float', 'bool' (interprets 0 or 'off' as False, anything else as True), 'value' (returns tuple (value, unit), where *value* is float), a callable (return the result

of this callable applied to the string value), a dictionary (return the stored value corresponding to the string value, or to the value converted into integer if the string value is not present), or a list of data types (the result is treated as a list of values with the given types separated by commas). *timeout* overrides the default value.

**read\_binary\_array\_data** (*include\_header=False, timeout=None, flush\_term=True*)

Read a binary data in the from the device.

The data assumes the standard binary transfer header consisting of "#" symbol, then a single digit with the size of the length string, then the length string containing the length of the binary data (in bytes). If *include\_header==True*, return the data with the header; otherwise, return only the content. If *flush\_term==True*, flush the following line to skip terminator characters after the binary data, which are added by some devices. *timeout* overrides the default value.

**reconnect** (*new\_instrument=True, ignore\_error=True*)

Remake the connection.

If *new\_instrument==True*, create a new backend instance. If *ignore\_error==True*, ignore errors on closing.

**reset** ()

Reset the device (by default, "\*RST" command)

**set\_device\_variable** (*key, value*)

Set the value of a settings parameter

**sleep** (*delay*)

Wait for *delay* seconds

**unlock** ()

Unlock the access to the device from other threads/processes (isn't necessarily implemented)

**using\_write\_buffer** ()

Context manager for using a write buffer.

While it's active, all the consecutive *write()* operations are bundled together with ; delimiter. The actual write is performed at the *read()/ask()* operation or at the end of the block.

**wait** (*wait\_type='sync', timeout=None, wait\_callback=None*)

Pause execution until device overlapped commands are complete.

*wait\_type* is either 'sync' (perform *wait\_sync()*), 'dev' (perform *wait\_dev()*) or 'none' (do nothing).

**wait\_dev** ()

Pause execution of the device commands until device overlapped commands (e.g., taking sweeps) are complete.

Note that the code execution is not paused.

**wait\_sync** (*timeout=None, wait\_callback=None*)

Pause execution of the script until device overlapped commands (e.g., taking sweeps) are complete.

*timeout* and *wait\_callback* override default constructor parameters.

**write** (*msg, arg=None, arg\_type=None, unit=None, bool\_selector=None, wait\_sync=None, read\_echo=False, read\_echo\_delay=0.0*)

Send a command.

#### Parameters

- **msg** (*str*) – Text message.

- **arg** – Optional argument to append in the end. If a list of arguments is supplied, the result is joined with ", ".
- **arg\_type** (*str*) – Argument type. Can be 'raw' (in which case data is sent raw), 'string', 'int', 'float', 'bool', a format string (such as '{:.3f}') or a list of argument types (for an iterable argument); if format string is used and the argument is a list or a tuple, then it is expanded as a list of arguments (e.g., `arg_type='{0};{1}'` with `arg=[1,2]` will produce a string '1;2'); if a list of types is used, each element of *arg* is converted using the corresponding type, and the result is joined with ", ".
- **unit** (*str*) – If not None, use it as a unit to append after the value.
- **bool\_selector** (*tuple*) – A tuple (*false\_value*, *true\_value*) of two strings to represent bool argument; by default, use `._bool_selector` attribute.
- **wait\_sync** – if True, append the sync command (specified as `._wait_sync_comm` attribute, `"*OPC?"` by default) after the message and pause the execution command is complete; useful in long set operations, where the device might ignore later inputs until the current command is complete; if None, use the class default `._default_write_sync` attribute (False by default).
- **read\_echo** (*bool*) – If True, read a single line after write.
- **read\_echo\_delay** (*float*) – The delay between write and read if `read_echo==True`.

**class** `pylablib.devices.Thorlabs.serial.FWv1` (*conn*, *pcount*=6, *respect\_bound*=True)

Bases: `pylablib.devices.Thorlabs.serial.ThorlabsSerialInterface`

Thorlabs FW102/212 v1.0 (older version) motorized filter wheels.

#### Parameters

- **conn** – serial connection parameters (usually port or a tuple containing port and baudrate)
- **pcount** – number of positions in the wheel
- **respect\_bound** (*bool*) – if True, avoid crossing the boundary between the first and the last position in the wheel

**ask** (*msg*, *data\_type*='string', *delay*=0.0, *timeout*=None, *read\_echo*=False)

Write a message and read a reply.

*msg* is the query message, *delay* is the delay between write and read. Other parameters are the same as in `read()`. If `read_echo==True`, assume that the device first echoes the input and skip it.

**get\_position** ()

Get the wheel position (starting from 1)

**set\_position** (*pos*)

Set the wheel position (starting from 1)

**get\_pcount** ()

Get the number of wheel positions (6 or 12)

**get\_trigger\_mode** ()

Get the trigger mode ("in" to input external trigger, "out" to output trigger)

**set\_trigger\_mode** (*trigger\_mode*)

Set the trigger mode ("in" to input external trigger, "out" to output trigger)

**BackendError**

alias of `pylablib.core.devio.comm_backend.DeviceBackendError`

**Error**

alias of `pylablib.devices.Thorlabs.base.ThorlabsError`

**class NoParameterCaller** (*device, kind*)

Bases: `object`

Class to simplify calling functions without a parameter

**ReraiseError**

alias of `pylablib.devices.Thorlabs.base.ThorlabsBackendError`

**apply\_settings** (*settings*)

Apply the settings.

*settings* is a dict {name: value} of the available device settings. Non-applicable settings are ignored.

**close** ()

Close the backend

**flush** (*one\_line=False*)

Flush the read buffer (read all the available data and return the number of bytes read).

If *one\_line*==True, read only a single line.

**static get\_arg\_type** (*arg*)

Autodetect argument type

**get\_device\_variable** (*key*)

Get the value of a settings, status, or full info parameter

**get\_esr** (*timeout=None*)

Get the device status register (by default, "\*ESR?" command)

**get\_full\_info** (*include=0*)

Get dict {name: value} containing full device information (including status and settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include*=-10 queries all available variables, which is equivalent to *include*="all".

**get\_full\_status** (*include=0*)

Get dict {name: value} containing the device status (including settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include*=-10 queries all available variables, which is equivalent to *include*="all".

**get\_id** (*timeout=None*)

Get the device IDN. (query SCPI '\*IDN?' command)

**get\_settings** (*include=0*)

Get dict {name: value} containing all the device settings.

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include*=-10 queries all available variables, which is equivalent to *include*="all".

**is\_opened()**

Check if the device is connected

**lock** (*timeout=None*)

Lock the access to the device from other threads/processes (isn't necessarily implemented)

**locking** (*timeout=None*)

Context manager for lock & unlock

**open()**

Open the backend

**static parse\_array\_data** (*data, fmt, include\_header=False*)

Parse the data returned by the device. *fmt* is *DataFormat* description in numpy format (e.g., "<u2").

If *include\_header==True*, the data is assumed to be in a (somewhat) standard SCPI format: *b'#'*, then a single digit *s* denoting length of the size block, then *s* digits denoting length of the data (in bytes) followed by the actual data. Otherwise (*include\_header==False*), assume that the header is already removed.

**read** (*data\_type='string', timeout=None*)

Read data from the device.

*data\_type* determines the type of the data. Can be 'raw' (just raw data), 'string' (with trailing and leading spaces stripped), 'int', 'float', 'bool' (interprets 0 or 'off' as False, anything else as True), 'value' (returns tuple (*value*, *unit*), where *value* is float), a callable (return the result of this callable applied to the string value), a dictionary (return the stored value corresponding to the string value, or to the value converted into integer if the string value is not present), or a list of data types (the result is treated as a list of values with the given types separated by commas). *timeout* overrides the default value.

**read\_binary\_array\_data** (*include\_header=False, timeout=None, flush\_term=True*)

Read a binary data in the from the device.

The data assumes the standard binary transfer header consisting of *"#"* symbol, then a single digit with the size of the length string, then the length string containing the length of the binary data (in bytes). If *include\_header==True*, return the data with the header; otherwise, return only the content. If *flush\_term==True*, flush the following line to skip terminator characters after the binary data, which are added by some devices. *timeout* overrides the default value.

**reconnect** (*new\_instrument=True, ignore\_error=True*)

Remake the connection.

If *new\_instrument==True*, create a new backend instance. If *ignore\_error==True*, ignore errors on closing.

**reset()**

Reset the device (by default, *"\*RST"* command)

**set\_device\_variable** (*key, value*)

Set the value of a settings parameter

**sleep** (*delay*)

Wait for *delay* seconds

**unlock()**

Unlock the access to the device from other threads/processes (isn't necessarily implemented)

**using\_write\_buffer()**

Context manager for using a write buffer.

While it's active, all the consecutive `write()` operations are bundled together with `;` delimiter. The actual write is performed at the `read()/ask()` operation or at the end of the block.

**wait** (*wait\_type='sync', timeout=None, wait\_callback=None*)

Pause execution until device overlapped commands are complete.

*wait\_type* is either 'sync' (perform `wait_sync()`), 'dev' (perform `wait_dev()`) or 'none' (do nothing).

**wait\_dev** ()

Pause execution of the device commands until device overlapped commands (e.g., taking sweeps) are complete.

Note that the code execution is not paused.

**wait\_sync** (*timeout=None, wait\_callback=None*)

Pause execution of the script until device overlapped commands (e.g., taking sweeps) are complete.

*timeout* and *wait\_callback* override default constructor parameters.

**write** (*msg, arg=None, arg\_type=None, unit=None, bool\_selector=None, wait\_sync=None, read\_echo=False, read\_echo\_delay=0.0*)

Send a command.

#### Parameters

- **msg** (*str*) – Text message.
- **arg** – Optional argument to append in the end. If a list of arguments is supplied, the result is joined with `", "`.
- **arg\_type** (*str*) – Argument type. Can be 'raw' (in which case data is sent raw), 'string', 'int', 'float', 'bool', a format string (such as `'{: .3f}'`) or a list of argument types (for an iterable argument); if format string is used and the argument is a list or a tuple, then it is expanded as a list of arguments (e.g., `arg_type='{0};{1}'` with `arg=[1, 2]` will produce a string `'1; 2'`); if a list of types is used, each element of *arg* is converted using the corresponding type, and the result is joined with `", "`.
- **unit** (*str*) – If not None, use it as a unit to append after the value.
- **bool\_selector** (*tuple*) – A tuple (*false\_value*, *true\_value*) of two strings to represent bool argument; by default, use `._bool_selector` attribute.
- **wait\_sync** – if True, append the sync command (specified as `._wait_sync_comm` attribute, `"*OPC?"` by default) after the message and pause the execution command is complete; useful in long set operations, where the device might ignore later inputs until the current command is complete; if None, use the class default `._default_write_sync` attribute (False by default).
- **read\_echo** (*bool*) – If True, read a single line after write.
- **read\_echo\_delay** (*float*) – The delay between write and read if `read_echo==True`.

**class** `pylablib.devices.Thorlabs.serial.MDT69xA` (*conn*)

Bases: `pylablib.devices.Thorlabs.serial.ThorlabsSerialInterface`

Thorlabs MDT693A/4A high-voltage source.

Uses MDT693A program interface, so should be compatible with both A and B versions (though it doesn't support all functions of MDT693B/4B)

**Parameters** `conn` – serial connection parameters (usually port or a tuple containing port and baudrate)

**get\_voltage** (*channel*='x')

Get the output voltage in Volts at a given channel

**set\_voltage** (*voltage*, *channel*='x')

Set the output voltage in Volts at a given channel

**get\_voltage\_range** ()

Get the selected voltage range in Volts (75, 100 or 150)

**BackendError**

alias of `pylablib.core.devio.comm_backend.DeviceBackendError`

**Error**

alias of `pylablib.devices.Thorlabs.base.ThorlabsError`

**class NoParameterCaller** (*device*, *kind*)

Bases: `object`

Class to simplify calling functions without a parameter

**ReraiseError**

alias of `pylablib.devices.Thorlabs.base.ThorlabsBackendError`

**apply\_settings** (*settings*)

Apply the settings.

*settings* is a dict {`name`: `value`} of the available device settings. Non-applicable settings are ignored.

**ask** (*msg*, *data\_type*='string', *delay*=0.0, *timeout*=None, *read\_echo*=False)

Write a message and read a reply.

*msg* is the query message, *delay* is the delay between write and read. Other parameters are the same as in `read()`. If `read_echo==True`, assume that the device first echoes the input and skip it.

**close** ()

Close the backend

**flush** (*one\_line*=False)

Flush the read buffer (read all the available data and return the number of bytes read).

If `one_line==True`, read only a single line.

**static get\_arg\_type** (*arg*)

Autodetect argument type

**get\_device\_variable** (*key*)

Get the value of a settings, status, or full info parameter

**get\_esr** (*timeout*=None)

Get the device status register (by default, "`*ESR?`" command)

**get\_full\_info** (*include*=0)

Get dict {`name`: `value`} containing full device information (including status and settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "`all`" (all available variables). Since the lowest priority is -10, setting `include=-10` queries all available variables, which is equivalent to `include="all"`.

**get\_full\_status** (*include*=0)

Get dict {`name`: `value`} containing the device status (including settings).



*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include*=-10 queries all available variables, which is equivalent to *include*="all".

**get\_id** (*timeout=None*)

Get the device IDN. (query SCPI '\*IDN?' command)

**get\_settings** (*include=0*)

Get dict {name: value} containing all the device settings.

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include*=-10 queries all available variables, which is equivalent to *include*="all".

**is\_opened** ()

Check if the device is connected

**lock** (*timeout=None*)

Lock the access to the device from other threads/processes (isn't necessarily implemented)

**locking** (*timeout=None*)

Context manager for lock & unlock

**open** ()

Open the backend

**static parse\_array\_data** (*data, fmt, include\_header=False*)

Parse the data returned by the device. *fmt* is *DataFormat* description in numpy format (e.g., "<u2").

If *include\_header*==True, the data is assumed to be in a (somewhat) standard SCPI format: b'#', then a single digit *s* denoting length of the size block, then *s* digits denoting length of the data (in bytes) followed by the actual data. Otherwise (*include\_header*==False), assume that the header is already removed.

**read** (*data\_type='string', timeout=None*)

Read data from the device.

*data\_type* determines the type of the data. Can be 'raw' (just raw data), 'string' (with trailing and leading spaces stripped), 'int', 'float', 'bool' (interprets 0 or 'off' as False, anything else as True), 'value' (returns tuple (value, unit), where *value* is float), a callable (return the result of this callable applied to the string value), a dictionary (return the stored value corresponding to the string value, or to the value converted into integer if the string value is not present), or a list of data types (the result is treated as a list of values with the given types separated by commas). *timeout* overrides the default value.

**read\_binary\_array\_data** (*include\_header=False, timeout=None, flush\_term=True*)

Read a binary data in the from the device.

The data assumes the standard binary transfer header consisting of "#" symbol, then a single digit with the size of the length string, then the length string containing the length of the binary data (in bytes). If *include\_header*==True, return the data with the header; otherwise, return only the content. If *flush\_term*==True, flush the following line to skip terminator characters after the binary data, which are added by some devices. *timeout* overrides the default value.

**reconnect** (*new\_instrument=True, ignore\_error=True*)

Remake the connection.

If *new\_instrument*==True, create a new backend instance. If *ignore\_error*==True, ignore errors on closing.

**reset** ()

Reset the device (by default, "`*RST`" command)

**set\_device\_variable** (*key*, *value*)

Set the value of a settings parameter

**sleep** (*delay*)

Wait for *delay* seconds

**unlock** ()

Unlock the access to the device from other threads/processes (isn't necessarily implemented)

**using\_write\_buffer** ()

Context manager for using a write buffer.

While it's active, all the consecutive `write()` operations are bundled together with `;` delimiter. The actual write is performed at the `read()/ask()` operation or at the end of the block.

**wait** (*wait\_type*=`'sync'`, *timeout*=`None`, *wait\_callback*=`None`)

Pause execution until device overlapped commands are complete.

*wait\_type* is either `'sync'` (perform `wait_sync()`), `'dev'` (perform `wait_dev()`) or `'none'` (do nothing).

**wait\_dev** ()

Pause execution of the device commands until device overlapped commands (e.g., taking sweeps) are complete.

Note that the code execution is not paused.

**wait\_sync** (*timeout*=`None`, *wait\_callback*=`None`)

Pause execution of the script until device overlapped commands (e.g., taking sweeps) are complete.

*timeout* and *wait\_callback* override default constructor parameters.

**write** (*msg*, *arg*=`None`, *arg\_type*=`None`, *unit*=`None`, *bool\_selector*=`None`, *wait\_sync*=`None`, *read\_echo*=`False`, *read\_echo\_delay*=`0.0`)

Send a command.

#### Parameters

- **msg** (*str*) – Text message.
- **arg** – Optional argument to append in the end. If a list of arguments is supplied, the result is joined with `" , "`.
- **arg\_type** (*str*) – Argument type. Can be `'raw'` (in which case data is sent raw), `'string'`, `'int'`, `'float'`, `'bool'`, a format string (such as `'{: .3f}'`) or a list of argument types (for an iterable argument); if format string is used and the argument is a list or a tuple, then it is expanded as a list of arguments (e.g., `arg_type='{0};{1}'` with `arg=[1, 2]` will produce a string `'1; 2'`); if a list of types is used, each element of *arg* is converted using the corresponding type, and the result is joined with `" , "`.
- **unit** (*str*) – If not `None`, use it as a unit to append after the value.
- **bool\_selector** (*tuple*) – A tuple (*false\_value*, *true\_value*) of two strings to represent bool argument; by default, use `._bool_selector` attribute.
- **wait\_sync** – if `True`, append the sync command (specified as `._wait_sync_comm` attribute, `"*OPC?"` by default) after the message and pause the execution command is complete; useful in long set operations, where the device might ignore later inputs until the current command is complete; if `None`, use the class default `._default_write_sync` attribute (`False` by default).

- **read\_echo** (*bool*) – If True, read a single line after write.
- **read\_echo\_delay** (*float*) – The delay between write and read if read\_echo==True.

## Module contents

### pylablib.devices.Toptica package

#### Submodules

#### pylablib.devices.Toptica.base module

**exception** pylablib.devices.Toptica.base.**TopticaError**

Bases: `pylablib.core.devio.base.DeviceError`

Generic Toptica device error

**args**

**with\_traceback()**

Exception.with\_traceback(tb) – set self.\_\_traceback\_\_ to tb and return self.

**exception** pylablib.devices.Toptica.base.**TopticaBackendError** (*exc*)

Bases: `pylablib.devices.Toptica.base.TopticaError`, `pylablib.core.devio.comm_backend.DeviceBackendError`

Toptica backend communication error

**args**

**with\_traceback()**

Exception.with\_traceback(tb) – set self.\_\_traceback\_\_ to tb and return self.

#### pylablib.devices.Toptica.ibeam module

`pylablib.devices.Toptica.ibeam.muxchan(*args, **kwargs)`

Multiplex the function over its addr argument

**class** pylablib.devices.Toptica.ibeam.**TDeviceInfo** (*serial, version*)

Bases: `tuple`

**count()**

Return number of occurrences of value.

**index()**

Return first index of value.

Raises ValueError if the value is not present.

**serial**

**version**

**class** pylablib.devices.Toptica.ibeam.**TWorkHours** (*power\_up, laser\_on*)

Bases: `tuple`

**count()**

Return number of occurrences of value.

**index()**  
Return first index of value.  
Raises ValueError if the value is not present.

**laser\_on**

**power\_up**

**class** `pylablib.devices.Toptica.ibeam.TTemperatures` (*diode, baseplate*)

Bases: `tuple`

**baseplate**

**count()**  
Return number of occurrences of value.

**diode**

**index()**  
Return first index of value.  
Raises ValueError if the value is not present.

**class** `pylablib.devices.Toptica.ibeam.TopticaIBeam` (*conn='COM1'*)

Bases: `pylablib.core.devio.comm_backend.ICommBackendWrapper`

Toptica iBeam smart laser controller.

#### Parameters

- **conn** – connection parameters - index of the Attocube ANC350 in the system (for a single controller leave 0)
- **timeout** (*float*) – default operation timeout

#### Error

alias of `pylablib.devices.Toptica.base.TopticaError`

**open()**  
Open the backend

**query** (*comm, multiline=False, keep\_whitespace=False, check\_error='FEW', reply=True*)

**reboot()**  
Reboot the laser system

**get\_device\_info()**  
Get the device info of the laser system: (serial, version)

**get\_full\_data** (*formatted=False*)  
Return the comprehensive device data

**get\_work\_hours()**  
Get the work hours (power on time and laser on time)

**get\_channels\_number()**  
Get number of supported laser channels

**is\_enabled()**  
Check if the output is enabled

**enable** (*enabled=True*)  
Turn the output on or off

**is\_channel\_enabled** (*channel='all'*)  
 Check if the specific channel is enabled

**enable\_channel** (*channel, enabled=True*)  
 Turn the specific channel on or off

**get\_channel\_power** (*channel='all'*)  
 Get specified channel power (in W)

**set\_channel\_power** (*channel, power*)  
 Set channel power (in W)

**get\_output\_power** ()  
 Get current output power (in W)

**get\_drive\_current** ()  
 Get current diode drive current (in A)

**get\_current\_limits** ()  
 Get settings of all current limits (in A) as a dictionary

**get\_temperatures** ()  
 Get settings of all current limits (in A) as a dictionary

**class NoParameterCaller** (*device, kind*)  
 Bases: `object`  
 Class to simplify calling functions without a parameter

**apply\_settings** (*settings*)  
 Apply the settings.  
*settings* is the dict {name: value} of the device available settings. Non-applicable settings are ignored.

**close** ()  
 Close the backend

**get\_device\_variable** (*key*)  
 Get the value of a settings, status, or full info parameter

**get\_full\_info** (*include=0*)  
 Get dict {name: value} containing full device information (including status and settings).  
*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_full\_status** (*include=0*)  
 Get dict {name: value} containing the device status (including settings).  
*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_settings** (*include=0*)  
 Get dict {name: value} containing all the device settings.  
*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the

lowest priority is -10, setting `include=-10` queries all available variables, which is equivalent to `include="all"`.

**is\_opened()**

Check if the device is connected

**lock** (*timeout=None*)

Lock the access to the device from other threads/processes (isn't necessarily implemented)

**locking** (*timeout=None*)

Context manager for lock & unlock

**set\_device\_variable** (*key, value*)

Set the value of a settings parameter

**unlock()**

Unlock the access to the device from other threads/processes (isn't necessarily implemented)

## Module contents

### pylablib.devices.Trinamic package

#### Submodules

#### pylablib.devices.Trinamic.base module

**exception** `pylablib.devices.Trinamic.base.TrinamicError`

Bases: `pylablib.core.devio.base.DeviceError`

Generic Trinamic error

**args**

**with\_traceback()**

Exception.with\_traceback(tb) – set self.\_\_traceback\_\_ to tb and return self.

**exception** `pylablib.devices.Trinamic.base.TrinamicBackendError` (*exc*)

Bases: `pylablib.devices.Trinamic.base.TrinamicError`, `pylablib.core.devio.comm_backend.DeviceBackendError`

Generic Trinamic backend communication error

**args**

**with\_traceback()**

Exception.with\_traceback(tb) – set self.\_\_traceback\_\_ to tb and return self.

**class** `pylablib.devices.Trinamic.base.TMCM1110` (*conn*)

Bases: `pylablib.core.devio.comm_backend.ICommBackendWrapper`, `pylablib.devices.interface.stage.IStage`

Trinamic stepper motor controller TMCM-1110 controlled using TMCL Firmware.

**Parameters** `conn` – serial connection parameters (usually port or a tuple containing port and baudrate)

**Error**

alias of `TrinamicError`

```

open ()
    Open the backend

class ReplyData (comm, status, value, addr, module)
    Bases: tuple

    addr

    comm

    count ()
        Return number of occurrences of value.

    index ()
        Return first index of value.

        Raises ValueError if the value is not present.

    module

    status

    value

query (comm, comm_type, value, result_format='i', bank=0, addr=0)
    Send a query to the stage and return the reply.

    For details, see TMCM-1110 firmware manual.

get_axis_parameter (parameter, result_format='i', addr=0)
    Get a given axis parameter

set_axis_parameter (parameter, value, addr=0)
    Set a given axis parameter (volatile; resets on power cycling)

store_axis_parameter (parameter, value=None, addr=0)
    Store a given axis parameter in EEPROM (by default, value is the current value)

get_global_parameter (parameter, result_format='i', bank=0, addr=0)
    Get a given global parameter

set_global_parameter (parameter, value, bank=0, addr=0)
    Set a given global parameter

get_general_input (port=0, bank=0, addr=0)
    Get value of an input at a given bank (0-2) and port.

    Bank 0 is digital input (7 ports), bank 1 is analog input (1 port, value from 0 to 2**16-1), bank 2 is digital
    output (8 ports). For port assignments, see TMCM-1110 firmware manual.

set_general_output (value, port=0, bank=2, addr=0)
    Set value of a digital input at a given bank (only bank 2 is available) and port.

    For port assignments, see TMCM-1110 firmware manual.

move_to (position, addr=0)
    Move to a given position

move_by (steps=1, addr=0)
    Move by a given number of steps

get_position (addr=0)
    Get the current axis position

set_position_reference (pos=0, addr=0)
    Set the current axis position as a reference (the actual motor position stays the same)

```

**jog** (*direction*, *speed=None*, *addr=0*)

Jog in a given direction with a given speed.

*direction* can be either "-" (negative, left) or "+" (positive, right). The motion continues until it is explicitly stopped, or until a limit is hit. If *speed* is *None*, use the standard speed value.

**stop** (*addr=0*)

Stop motion

**get\_microstep\_resolution** (*addr=0*)

Get the number of microsteps per full step (always a power of 2)

**set\_microstep\_resolution** (*resolution*, *addr=0*)

Set the number of microsteps per full step (rounded to a nearest power of 2)

**get\_current\_parameters** (*addr=0*)

Return diving current parameter (*drive\_current*, *standby\_current*).

*drive\_current* is the maximal drive current, which is given as a fraction of the maximal generated current (which is either 1A or 2.8A depending on the hardware jumper). *standby\_current* is given as a fraction of *drive\_current*.

**setup\_current** (*drive\_current=None*, *standby\_current=None*, *addr=0*)

Set drive and standby currents.

WARNING: too high of a setting might damage the motor. *drive\_current* is the maximal drive current, which is given as a fraction of the maximal generated current (which is either 1A or 2.8A depending on the hardware jumper). *standby\_current* is given as a fraction of *drive\_current*. Any *None* parameters are left unchanged.

**get\_limit\_switches\_parameters** (*addr=0*)

Return limit switch parameters (*left\_enable*, *right\_enable*)

**setup\_limit\_switches** (*left\_enable=None*, *right\_enable=None*, *addr=0*)

Setup limit switch parameters

**get\_velocity\_parameters** (*addr=0*)

Return velocity parameters (*speed*, *accel*, *pulse\_divisor*, *ramp\_divisor*).

*speed* and *accel* denote, correspondingly, maximal (i.e., steady regime) moving speed and acceleration in *internal* units. *pulse\_divisor* is the driver pulse divisor, which defines how internal velocity units translate into microsteps/s (see [get\\_velocity\\_factor\(\)](#)); can only be a power of 2, higher values mean slower motion. *ramp\_divisor* is the driver ramp divisor, which, together with the pulse divisor, defines how internal acceleration units translate into microsteps/s<sup>2</sup> (see [get\\_acceleration\\_factor\(\)](#)); rounded to the nearest power of 2, higher values mean slower acceleration.

**setup\_velocity** (*speed=None*, *accel=None*, *pulse\_divisor=None*, *ramp\_divisor=None*, *addr=0*)

Setup velocity parameters (*speed*, *accel*, *pulse\_divisor*, *ramp\_divisor*).

*speed* and *accel* denote, correspondingly, maximal (i.e., steady regime) moving speed and acceleration in *internal* units. *pulse\_divisor* is the driver pulse divisor, which defines how internal velocity units translate into microsteps/s (see [get\\_velocity\\_factor\(\)](#)); rounded to the nearest power of 2, higher values mean slower motion. *ramp\_divisor* is the driver ramp divisor, which, together with the pulse divisor, defines how internal acceleration units translate into microsteps/s<sup>2</sup> (see [get\\_acceleration\\_factor\(\)](#)); rounded to the nearest power of 2, higher values mean slower acceleration. *None* values are left unchanged.

**get\_velocity\_factor** (*addr=0*)

Get the ratio between the real speed (in microsteps/s) and the internal units



**get\_acceleration\_factor** (*addr=0*)

Get the ratio between the real acceleration (in microsteps/s<sup>2</sup>) and the internal units

**get\_current\_speed** (*addr=0*)

Get the instantaneous speed in internal units

**is\_moving** (*addr=0*)

Check if the motor is moving

**wait\_move** (*addr=0*)

Wait until motion is done

**class NoParameterCaller** (*device, kind*)

Bases: `object`

Class to simplify calling functions without a parameter

**apply\_settings** (*settings*)

Apply the settings.

*settings* is the dict {name: value} of the device available settings. Non-applicable settings are ignored.

**close** ()

Close the backend

**get\_device\_variable** (*key*)

Get the value of a settings, status, or full info parameter

**get\_full\_info** (*include=0*)

Get dict {name: value} containing full device information (including status and settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_full\_status** (*include=0*)

Get dict {name: value} containing the device status (including settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_settings** (*include=0*)

Get dict {name: value} containing all the device settings.

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**is\_opened** ()

Check if the device is connected

**lock** (*timeout=None*)

Lock the access to the device from other threads/processes (isn't necessarily implemented)

**locking** (*timeout=None*)

Context manager for lock & unlock

**set\_device\_variable** (*key, value*)

Set the value of a settings parameter

**unlock** ()

Unlock the access to the device from other threads/processes (isn't necessarily implemented)

## Module contents

### pylablib.devices.interface package

#### Submodules

#### pylablib.devices.interface.camera module

**exception** pylablib.devices.interface.camera.**DefaultFrameTransferError**

Bases: *pylablib.core.devio.base.DeviceError*

Generic frame transfer error

**args**

**with\_traceback** ()

Exception.with\_traceback(tb) – set self.\_\_traceback\_\_ to tb and return self.

**class** pylablib.devices.interface.camera.**TFramesStatus** (*acquired, unread, skipped, buffer\_size*)

Bases: *tuple*

**acquired**

**buffer\_size**

**count** ()

Return number of occurrences of value.

**index** ()

Return first index of value.

Raises ValueError if the value is not present.

**skipped**

**unread**

**class** pylablib.devices.interface.camera.**TFrameSize** (*width, height*)

Bases: *tuple*

**count** ()

Return number of occurrences of value.

**height**

**index** ()

Return first index of value.

Raises ValueError if the value is not present.

**width**

**class** pylablib.devices.interface.camera.**TFramePosition** (*left, top*)

Bases: *tuple*

**count()**  
Return number of occurrences of value.

**index()**  
Return first index of value.  
  
Raises ValueError if the value is not present.

**left**

**top**

**class** `pylablib.devices.interface.camera.TFrameInfo` (*frame\_index*)  
Bases: `tuple`

**count()**  
Return number of occurrences of value.

**frame\_index**

**index()**  
Return first index of value.  
  
Raises ValueError if the value is not present.

**class** `pylablib.devices.interface.camera.ICamera` (*\*args, \*\*kwargs*)  
Bases: `pylablib.core.devio.interface.IDevice`

Generic camera class.

Provides a consistent common interface for the most frequently encountered camera functions.

**Error**  
alias of `pylablib.core.devio.base.DeviceError`

**TimeoutError**  
alias of `pylablib.core.devio.base.DeviceError`

**FrameTransferError**  
alias of `DefaultFrameTransferError`

**is\_acquisition\_setup()**  
Check if acquisition is set up.  
  
If the camera does not support separate acquisition setup, always return True.

**get\_acquisition\_parameters()**  
Get acquisition parameters.  
  
Return dictionary {name: value}

**setup\_acquisition** (*\*\*kwargs*)  
Setup acquisition.

Any non-specified acquisition parameters are assumed to be the same as previously set (or default, if not explicitly set before). Return the new acquisition parameters.

**clear\_acquisition()**  
Clear acquisition settings

**start\_acquisition** (*\*args, \*\*kwargs*)  
Start acquisition.

Can take the same keyword parameters as `:meth:~setup_acquisition`. If the acquisition is not set up yet, set it up using the supplied parameters (use default of `setup_acquisition()`, if the parameter is

None). Otherwise, if any supplied parameters are different from the current ones, change them and reset the acquisition.

**stop\_acquisition()**

Stop acquisition

**acquisition\_in\_progress()**

Check if acquisition is in progress

**pausing\_acquisition**(*clear=None, stop=True, setup\_after=None, start\_after=True, combine\_nested=True*)

Context manager which temporarily pauses acquisition during execution of `with` block.

Useful for applying certain settings which can't be changed during the acquisition. If `clear==True`, clear acquisition in addition to stopping (by default, use the class default specified as `_clear_pausing_acquisition` attribute). If `stop==True`, stop the acquisition (if `clear==True`, stop regardless). If `setup_after==True`, setup the acquisition after pause if necessary (None means setup only if clearing was required). If `start_after==True`, start the acquisition after pause if necessary (None means start only if stopping was required). If `combine_nested==True`, then any nested `pausing_acquisition` calls will stop/clear acquisition as necessary, but won't setup/start it again until this `pausing_acquisition` call is complete.

Yields tuple (`acq_in_progress, acq_params`), which indicates whether acquisition is currently in progress, and what are the current acquisition parameters.

**get\_detector\_size()**

Get camera detector size (in pixels) as a tuple (`width, height`)

**get\_frames\_status()**

Get acquisition and buffer status.

Return tuple (`acquired, unread, skipped, size`), where `acquired` is the total number of acquired frames, `unread` is the number of acquired but not read frames, `skipped` is the number of skipped (not read and then written over) frames, and `buffer_size` is the total buffer size (in frames).

**wait\_for\_frame**(*since='lastread', nframes=1, timeout=20.0, error\_on\_stopped=False*)

Wait for one or several new camera frames.

*since* specifies the reference point for waiting to acquire *nframes* frames; can be "lastread" (from the last read frame), "lastwait" (wait for the last successful `wait_for_frame()` call), "now" (from the start of the current call), or "start" (from the acquisition start, i.e., wait until *nframes* frames have been acquired). *timeout* can be either a number, None (infinite timeout), or a tuple (`timeout, frame_timeout`), in which case the call times out if the total time exceeds `timeout`, or a single frame wait exceeds `frame_timeout`. If the call times out, raise `TimeoutError`. If `error_on_stopped==True` and the acquisition is not running, raise `Error`; otherwise, simply return `False` without waiting.

**get\_image\_indexing()**

Get indexing for the returned images.

Can be "rct" (first index row, second index column, rows counted from the top), "rcb" (same as "rc", rows counted from the bottom), "xyt" (first index column, second index row, rows counted from the top), or "xyb" (same as "xyt", rows counted from the bottom)

**set\_image\_indexing**(*indexing*)

Set up indexing for the returned images.

Can be "rct" (first index row, second index column, rows counted from the top), "rcb" (same as "rc", rows counted from the bottom), "xyt" (first index column, second index row, rows counted from the top), or "xyb" (same as "xyt", rows counted from the bottom)

**get\_data\_dimensions()**

Get readout data dimensions (in pixels) as a tuple (width, height); take indexing mode into account

**get\_frame\_format()**

Get format for the returned images.

Can be "list" (list of 2D arrays), "array" (a single 3D array), or "chunks" (list of 3D "chunk" arrays; supported for some cameras and provides the best performance).

**set\_frame\_format(fmt)**

Set format for the returned images.

Can be "list" (list of 2D arrays), "array" (a single 3D array), "chunks" (list of 3D "chunk" arrays; supported for some cameras and provides the best performance), or "try\_chunks" (same as "chunks", but if chunks are not supported, set to "list" instead). If format is "chunks" and chunks are not supported by the camera, it results in one frame per chunk. Note that if the format is set to "array" or "chunks", the frame info format is also automatically set to "array". If the format is set to "chunks", then the image info is also returned in chunks form (list of 2D info arrays with the same length as the corresponding frame chunks).

**get\_frame\_info\_format()**

Get format of the frame info.

Can be "namedtuple" (potentially nested named tuples; convenient to get particular values), "list" (flat list of values, with field names are given by `get_frame_info_fields()`; convenient for building a table), "array" (same as "list", but with a numpy array, which is easier to use for "chunks" frame format), or "dict" (flat dictionary with the same fields as the "list" format; more resilient to future format changes)

**set\_frame\_info\_format(fmt, include\_fields=None)**

Set format of the frame info.

Can be "namedtuple" (potentially nested named tuples; convenient to get particular values), "list" (flat list of values, with field names are given by `get_frame_info_fields()`; convenient for building a table), "array" (same as "list", but with a numpy array, which is easier to use for "chunks" frame format), or "dict" (flat dictionary with the same fields as the "list" format; more resilient to future format changes) If `include_fields` is not None, it specifies the fields included for non-"tuple" formats; note that order of `include_fields` is ignored, and the resulting fields are always ordered same as in the original.

**get\_frame\_info\_period()**

Get period of frame info acquisition.

Frame info might be skipped (set to None) except for frames which indices are divisible by *period*. Useful for certain cameras where acquiring frame info takes a lot of time and can reduce performance at higher frame rates. Note that this parameter can still be ignored (i.e., always set to 1) if the performance is not an issue for a given camera class.

**set\_frame\_info\_period(period=1)**

Set period of frame info acquisition.

Frame info might be skipped (set to None) except for frames which indices are divisible by *period*. Useful for certain cameras where acquiring frame info takes a lot of time and can reduce performance at higher frame rates. Note that this parameter can still be ignored (i.e., always set to 1) if the performance is not an issue for a given camera class.

**get\_frame\_info\_fields()**

Get the names of frame info fields.

Applicable when frame info format (set by `set_frame_info_format()`) is "list" or "array".

**get\_new\_images\_range()**

Get the range of the new images.

Return tuple (*first*, *last*) with images range (first inclusive). If no images are available, return None. If some images were in the buffer were overwritten, exclude them from the range.

**read\_multiple\_images** (*rng=None*, *peek=False*, *missing\_frame='skip'*, *return\_info=False*, *return\_rng=False*)

Read multiple images specified by *rng* (by default, all un-read images).

If *rng* is specified, it is a tuple (*first*, *last*) with images range (first inclusive). If no new frames are available, return an empty list; if no acquisition is running, return None. If *peek==True*, return images but not mark them as read. *missing\_frame* determines what to do with frames which are out of range (missing or lost): can be "none" (replacing them with None), "zero" (replacing them with zero-filled frame), or "skip" (skipping them). If *return\_info==True*, return tuple (*frames*, *infos*), where *infos* is a list of frame info tuples (camera-dependent, by default, only the frame index); if some frames are missing and *missing\_frame!="skip"*, the corresponding frame info is None. if *return\_rng==True*, return the range covered resulting frames; if *missing\_frame=="skip"*, the range can be smaller than the supplied *rng* if some frames are skipped.

**read\_oldest\_image** (*peek=False*, *return\_info=False*)

Read the oldest un-read image.

If no un-read frames are available, return None. If *peek==True*, return the image but not mark it as read. If *return\_info==True*, return tuple (*frame*, *info*), where *info* is an info tuples (camera-dependent, see [read\\_multiple\\_images\(\)](#)).

**read\_newest\_image** (*peek=False*, *return\_info=False*)

Read the newest un-read image.

If no un-read frames are available, return None. If *peek==True*, return the image but not mark it as read. If *return\_info==True*, return tuple (*frame*, *info*), where *info* is an info tuples (camera-dependent, see [read\\_multiple\\_images\(\)](#)).

**grab** (*nframes=1*, *frame\_timeout=5.0*, *missing\_frame='skip'*, *return\_info=False*, *buff\_size=None*)

Snap *nframes* images (with preset image read mode parameters)

*buff\_size* determines buffer size (if None, use the default size). Timeout is specified for a single-frame acquisition, not for the whole acquisition time. *missing\_frame* determines what to do with frames which have been lost: can be "none" (replacing them with None), "zero" (replacing them with zero-filled frame), or "skip" (skipping them, while still keeping total returned frames number to *n*). If *return\_info==True*, return tuple (*frames*, *infos*), where *infos* is a list of frame info tuples (camera-dependent); if some frames are missing and *missing\_frame!="skip"*, the corresponding frame info is None.

**snap** (*timeout=5.0*, *return\_info=False*)

Snap a single frame

**class NoParameterCaller** (*device*, *kind*)

Bases: `object`

Class to simplify calling functions without a parameter

**apply\_settings** (*settings*)

Apply the settings.

*settings* is the dict {*name*: *value*} of the device available settings. Non-applicable settings are ignored.

**close()**

Close the connection

**get\_device\_variable** (*key*)

Get the value of a settings, status, or full info parameter

**get\_full\_info** (*include=0*)

Get dict {name: value} containing full device information (including status and settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_full\_status** (*include=0*)

Get dict {name: value} containing the device status (including settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_settings** (*include=0*)

Get dict {name: value} containing all the device settings.

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**is\_opened** ()

Check if the device is connected

**open** ()

Open the connection

**set\_device\_variable** (*key, value*)

Set the value of a settings parameter

`pylablib.devices.interface.camera.acqstopped (*args, **kwargs)`

Decorator which temporarily stops acquisition for the function call

`pylablib.devices.interface.camera.acqcleared (*args, **kwargs)`

Decorator which temporarily clears acquisition for the function call

`pylablib.devices.interface.camera.trim_frames (frames, l, info=None, chunks='auto')`

Trim frames in different formats to the desired length

**class** `pylablib.devices.interface.camera.FrameCounter`

Bases: `object`

Frame counter.

Keeps track of the buffer occupation, acquired/missed frames, last read and wait buffers, etc.

**reset** (*buffer\_size=None*)

Reset the counters.

If *buffer\_size* is None, assume the the buffer is deallocated. Otherwise, it specifies the frame buffer size (in frames).

**update\_acquired\_frames** (*acquired\_frames*)

Update the counter of acquired frames (needs to be called by the camera whenever necessary)

**wait\_start** (*acquired\_frames*)

Set up waiting routine (called in the beginning of `ICamera.wait_for_frame()`)

**is\_wait\_done** (*acquired\_frames=None, since='lastread', nframes=1*)

Check if the waiting condition is satisfied based on the counter values:

If not `None`, *acquired\_frames* specifies the most recent number of acquired frames (the internal counters is automatically updated). *since* and *nframes* have the same meaning as in `ICamera.wait_for_frame()`.

**wait\_done** ()

Clean up waiting routine (called in the end of `ICamera.wait_for_frame()`)

**get\_frames\_status** (*acquired\_frames=None*)

Get status of the internal counters.

Return tuple (*acquired, unread, skipped, buffer\_size*). If the buffer is not allocated, all counters are 0.

**get\_new\_frames\_range** (*acquired\_frames=None*)

Get the range of the new frames (acquired but not read)

**trim\_frames\_range** (*rng*)

Trim the given frames range to only contains frames which are still in the buffer (i.e., remove the frames which are too old and have been overwritten)

**advance\_read\_frames** (*rng*)

Mark the specified frames range as read and advance the last read counter

**set\_first\_valid\_frame** (*first\_valid\_frame*)

Set the first valid frame; all frames older than it are considered invalid when calculating skipped frames and trimming ranges

**class** `pylablib.devices.interface.camera.FrameNotifier` (*strict=False*)

Bases: `object`

Notifier for a new available frame.

Used when the camera runs a separate polling thread or a callback, which needs to notify the main thread that a new frame has been acquired.

**Parameters** **strict** – determines whether `wait()` waits for a specified frame index, or just for any new frame (which is checked later)

**reset** ()

Reset the internal frame counter

**inc** ()

Increment the internal frame counter, notify the waiting threads, and return the counter value

**wait** (*idx=None, timeout=None*)

Wait for a new frame with a given index (if `None`, for the next acquired frame)

**class** `pylablib.devices.interface.camera.ChunkBufferManager` (*chunk\_size=1048576*)

Bases: `object`

Buffer manager, which takes care of creating and removing the buffer chunks, and reading out some parts of them.

**Parameters** **chunk\_size** – the minimal size of a single buffer chunk (continuous memory segment potentially containing several frames).

**get\_ctypes\_frames\_list** (*ctype=<class 'ctypes.c\_char\_p'>*)

Get stored buffers as a ctypes array with pointer of the given type

**get\_frames\_data** (*idx, nframes=1*)

Get frames data starting from *idx* and spanning *nframes* frames.



Return a list of tuples (`nread`, `chunk_data`), where `nread` is the number of frames in the chunk, and `chunk_data` is the raw buffer pointer as a `ctypes.c_char_p` object.

**allocate** (*nframes*, *frame\_size*)

Allocate buffers for the given number of frames and frame size (in bytes)

**deallocate** ()

Deallocate the buffers

**class** `pylablib.devices.interface.camera.IAttributeCamera` (*\*args*, *\*\*kwargs*)

Bases: `pylablib.devices.interface.camera.ICamera`

Camera class which supports camera attributes.

The method `_list_attributes` must be defined in a subclass; it should produce a list of camera attributes, which have name attribute for placing them into a dictionary. Attributes can also have readable and writable attributes, which are used in `get_all_attribute_values()` and `set_all_attribute_values()` to determine if the attribute values should be collected or set. Method `_update_attributes` should be called on opening to populate the dictionary of available attributes.

One can also define `_normalize_attribute_name`, which normalizes the attribute name into a dictionary name (e.g., replaces separators, removes spaces, or normalizes case).

**get\_attribute** (*name*, *error\_on\_missing=True*)

Get the camera attribute with the given name

**get\_all\_attributes** (*copy=False*)

Return a dictionary of all available attributes.

If `copy==True`, copy the dictionary; otherwise, return the internal dictionary structure (should not be modified).

**get\_attribute\_value** (*name*, *error\_on\_missing=True*, *default=None*, *\*\*kwargs*)

Get value of an attribute with the given name.

If the value doesn't exist and `error_on_missing==True`, raise error; otherwise, return *default*. If *default* is not `None`, automatically assume that `error_on_missing==False`. If *name* points at a dictionary branch, return a dictionary with all values in this branch. Additional arguments are passed to `get_value` methods of the individual attribute.

**set\_attribute\_value** (*name*, *value*, *error\_on\_missing=True*, *\*\*kwargs*)

Set value of an attribute with the given name.

If the value doesn't exist and `error_on_missing==True`, raise error; otherwise, do nothing. If *name* points at a dictionary branch, set all values in this branch (in this case *value* must be a dictionary). Additional arguments are passed to `set_value` methods of the individual attribute.

**get\_all\_attribute\_values** (*root="*, *\*\*kwargs*)

Get values of all attributes with the given *root*.

Additional arguments are passed to `get_value` methods of individual attributes.

**set\_all\_attribute\_values** (*settings*, *root="*, *\*\*kwargs*)

Set values of all attributes with the given *root*.

Additional arguments are passed to `set_value` methods of individual attributes.

**Error**

alias of `pylablib.core.devio.base.DeviceError`

**FrameTransferError**

alias of `DefaultFrameTransferError`

**class NoParameterCaller** (*device, kind*)

Bases: `object`

Class to simplify calling functions without a parameter

**TimeoutError**

alias of `pylablib.core.devio.base.DeviceError`

**acquisition\_in\_progress** ()

Check if acquisition is in progress

**apply\_settings** (*settings*)

Apply the settings.

*settings* is the dict {name: value} of the device available settings. Non-applicable settings are ignored.

**clear\_acquisition** ()

Clear acquisition settings

**close** ()

Close the connection

**get\_acquisition\_parameters** ()

Get acquisition parameters.

Return dictionary {name: value}

**get\_data\_dimensions** ()

Get readout data dimensions (in pixels) as a tuple (width, height); take indexing mode into account

**get\_detector\_size** ()

Get camera detector size (in pixels) as a tuple (width, height)

**get\_device\_variable** (*key*)

Get the value of a settings, status, or full info parameter

**get\_frame\_format** ()

Get format for the returned images.

Can be "list" (list of 2D arrays), "array" (a single 3D array), or "chunks" (list of 3D "chunk" arrays; supported for some cameras and provides the best performance).

**get\_frame\_info\_fields** ()

Get the names of frame info fields.

Applicable when frame info format (set by `set_frame_info_format()`) is "list" or "array".

**get\_frame\_info\_format** ()

Get format of the frame info.

Can be "namedtuple" (potentially nested named tuples; convenient to get particular values), "list" (flat list of values, with field names are given by `get_frame_info_fields()`; convenient for building a table), "array" (same as "list", but with a numpy array, which is easier to use for "chunks" frame format), or "dict" (flat dictionary with the same fields as the "list" format; more resilient to future format changes)

**get\_frame\_info\_period** ()

Get period of frame info acquisition.

Frame info might be skipped (set to None) except for frames which indices are divisible by *period*. Useful for certain cameras where acquiring frame info takes a lot of time and can reduce performance at higher frame rates. Note that this parameter can still be ignored (i.e., always set to 1) if the performance is not an issue for a given camera class.

**get\_frames\_status()**

Get acquisition and buffer status.

Return tuple (*acquired*, *unread*, *skipped*, *size*), where *acquired* is the total number of acquired frames, *unread* is the number of acquired but not read frames, *skipped* is the number of skipped (not read and then written over) frames, and *buffer\_size* is the total buffer size (in frames).

**get\_full\_info(include=0)**

Get dict {*name*: *value*} containing full device information (including status and settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include*=-10 queries all available variables, which is equivalent to *include*="all".

**get\_full\_status(include=0)**

Get dict {*name*: *value*} containing the device status (including settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include*=-10 queries all available variables, which is equivalent to *include*="all".

**get\_image\_indexing()**

Get indexing for the returned images.

Can be "rct" (first index row, second index column, rows counted from the top), "rcb" (same as "rc", rows counted from the bottom), "xyt" (first index column, second index row, rows counted from the top), or "xyb" (same as "xyt", rows counted from the bottom)

**get\_new\_images\_range()**

Get the range of the new images.

Return tuple (*first*, *last*) with images range (first inclusive). If no images are available, return None. If some images were in the buffer were overwritten, exclude them from the range.

**get\_settings(include=0)**

Get dict {*name*: *value*} containing all the device settings.

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include*=-10 queries all available variables, which is equivalent to *include*="all".

**grab(nframes=1, frame\_timeout=5.0, missing\_frame='skip', return\_info=False, buff\_size=None)**

Snap *nframes* images (with preset image read mode parameters)

*buff\_size* determines buffer size (if None, use the default size). Timeout is specified for a single-frame acquisition, not for the whole acquisition time. *missing\_frame* determines what to do with frames which have been lost: can be "none" (replacing them with None), "zero" (replacing them with zero-filled frame), or "skip" (skipping them, while still keeping total returned frames number to *n*). If *return\_info*==True, return tuple (*frames*, *infos*), where *infos* is a list of frame info tuples (camera-dependent); if some frames are missing and *missing\_frame*!="skip", the corresponding frame info is None.

**is\_acquisition\_setup()**

Check if acquisition is set up.

If the camera does not support separate acquisition setup, always return True.

**is\_opened()**

Check if the device is connected

**open()**

Open the connection

**pausing\_acquisition**(*clear=None, stop=True, setup\_after=None, start\_after=True, combine\_nested=True*)

Context manager which temporarily pauses acquisition during execution of `with` block.

Useful for applying certain settings which can't be changed during the acquisition. If `clear==True`, clear acquisition in addition to stopping (by default, use the class default specified as `_clear_pausing_acquisition` attribute). If `stop==True`, stop the acquisition (if `clear==True`, stop regardless). If `setup_after==True`, setup the acquisition after pause if necessary (None means setup only if clearing was required). If `start_after==True`, start the acquisition after pause if necessary (None means start only if stopping was required). If `combine_nested==True`, then any nested `pausing_acquisition` calls will stop/clear acquisition as necessary, but won't setup/start it again until this `pausing_acquisition` call is complete.

Yields tuple (`acq_in_progress, acq_params`), which indicates whether acquisition is currently in progress, and what are the current acquisition parameters.

**read\_multiple\_images**(*rng=None, peek=False, missing\_frame='skip', return\_info=False, return\_rng=False*)

Read multiple images specified by `rng` (by default, all un-read images).

If `rng` is specified, it is a tuple (`first, last`) with images range (first inclusive). If no new frames are available, return an empty list; if no acquisition is running, return `None`. If `peek==True`, return images but not mark them as read. `missing_frame` determines what to do with frames which are out of range (missing or lost): can be "none" (replacing them with `None`), "zero" (replacing them with zero-filled frame), or "skip" (skipping them). If `return_info==True`, return tuple (`frames, infos`), where `infos` is a list of frame info tuples (camera-dependent, by default, only the frame index); if some frames are missing and `missing_frame!="skip"`, the corresponding frame info is `None`. if `return_rng==True`, return the range covered resulting frames; if `missing_frame=="skip"`, the range can be smaller than the supplied `rng` if some frames are skipped.

**read\_newest\_image**(*peek=False, return\_info=False*)

Read the newest un-read image.

If no un-read frames are available, return `None`. If `peek==True`, return the image but not mark it as read. If `return_info==True`, return tuple (`frame, info`), where `info` is an info tuples (camera-dependent, see `read_multiple_images()`).

**read\_oldest\_image**(*peek=False, return\_info=False*)

Read the oldest un-read image.

If no un-read frames are available, return `None`. If `peek==True`, return the image but not mark it as read. If `return_info==True`, return tuple (`frame, info`), where `info` is an info tuples (camera-dependent, see `read_multiple_images()`).

**set\_device\_variable**(*key, value*)

Set the value of a settings parameter

**set\_frame\_format**(*fmt*)

Set format for the returned images.

Can be "list" (list of 2D arrays), "array" (a single 3D array), "chunks" (list of 3D "chunk" arrays; supported for some cameras and provides the best performance), or "try\_chunks" (same as "chunks", but if chunks are not supported, set to "list" instead). If format is "chunks" and chunks are not supported by the camera, it results in one frame per chunk. Note that if the format is set to "array" or "chunks", the frame info format is also automatically set to "array". If the format is set to "chunks", then the image info is also returned in chunks form (list of 2D info arrays with the same length as the corresponding frame chunks).

**set\_frame\_info\_format** (*fmt*, *include\_fields=None*)

Set format of the frame info.

Can be "namedtuple" (potentially nested named tuples; convenient to get particular values), "list" (flat list of values, with field names are given by `get_frame_info_fields()`; convenient for building a table), "array" (same as "list", but with a numpy array, which is easier to use for "chunks" frame format), or "dict" (flat dictionary with the same fields as the "list" format; more resilient to future format changes) If *include\_fields* is not None, it specifies the fields included for non-"tuple" formats; note that order or *include\_fields* is ignored, and the resulting fields are always ordered same as in the original.

**set\_frame\_info\_period** (*period=1*)

Set period of frame info acquisition.

Frame info might be skipped (set to None) except for frames which indices are divisible by *period*. Useful for certain cameras where acquiring frame info takes a lot of time and can reduce performance at higher frame rates. Note that this parameter can still be ignored (i.e., always set to 1) if the performance is not an issue for a given camera class.

**set\_image\_indexing** (*indexing*)

Set up indexing for the returned images.

Can be "rct" (first index row, second index column, rows counted from the top), "rcb" (same as "rc", rows counted from the bottom), "xyt" (first index column, second index row, rows counted from the top), or "xyb" (same as "xyt", rows counted from the bottom)

**setup\_acquisition** (*\*\*kwargs*)

Setup acquisition.

Any non-specified acquisition parameters are assumed to be the same as previously set (or default, if not explicitly set before). Return the new acquisition parameters.

**snap** (*timeout=5.0*, *return\_info=False*)

Snap a single frame

**start\_acquisition** (*\*args*, *\*\*kwargs*)

Start acquisition.

Can take the same keyword parameters as *meth:“setup\_acquisition*. If the acquisition is not set up yet, set it up using the supplied parameters (use default of `setup_acquisition()`, if the parameter is None). Otherwise, if any supplied parameters are different from the current ones, change them and reset the acquisition.

**stop\_acquisition** ()

Stop acquisition

**wait\_for\_frame** (*since='lastread'*, *nframes=1*, *timeout=20.0*, *error\_on\_stopped=False*)

Wait for one or several new camera frames.

*since* specifies the reference point for waiting to acquire *nframes* frames; can be "lastread" (from the last read frame), "lastwait" (wait for the last successful `wait_for_frame()` call), "now" (from the start of the current call), or "start" (from the acquisition start, i.e., wait until *nframes* frames have been acquired). *timeout* can be either a number, None (infinite timeout), or a tuple (*timeout*, *frame\_timeout*), in which case the call times out if the total time exceeds *timeout*, or a single frame wait exceeds *frame\_timeout*. If the call times out, raise `TimeoutError`. If *error\_on\_stopped==True* and the acquisition is not running, raise `Error`; otherwise, simply return False without waiting.

**class** `pylablib.devices.interface.camera.IGrabberAttributeCamera` (*\*args*,  
*\*\*kwargs*)

Bases: `pylablib.devices.interface.camera.ICamera`

Camera class which supports frame grabber attributes.

Essentially the same as *IAtributeCamera*, but with relevant methods and attributes renamed to support both frame grabber and camera attributes handling simultaneously.

The method `_list_grabber_attributes` must be defined in a subclass; it should produce a list of camera attributes, which have name attribute for placing them into a dictionary. Attributes can also have readable and writable attributes, which are used in `get_all_grabber_attribute_values()` and `set_all_grabber_attribute_values()` to determine if the attribute values should be collected or set. Method `_update_grabber_attributes` should be called on opening to populate the dictionary of available attributes.

One can also define `_normalize_grabber_attribute_name`, which normalizes the attribute name into a dictionary name (e.g., replaces separators, removes spaces, or normalizes case).

**get\_grabber\_attribute** (*name*, *error\_on\_missing=True*)

Get the camera attribute with the given name

**get\_all\_grabber\_attributes** (*copy=False*)

Return a dictionary of all available frame grabber grabber\_attributes.

If *copy==True*, copy the dictionary; otherwise, return the internal dictionary structure (should not be modified).

**get\_grabber\_attribute\_value** (*name*, *error\_on\_missing=True*, *default=None*, *\*\*kwargs*)

Get value of a frame grabber attribute with the given name.

If the value doesn't exist and *error\_on\_missing==True*, raise error; otherwise, return *default*. If *default* is not None, automatically assume that *error\_on\_missing==False*. If *name* points at a dictionary branch, return a dictionary with all values in this branch. Additional arguments are passed to `get_value` methods of the individual attribute.

**set\_grabber\_attribute\_value** (*name*, *value*, *error\_on\_missing=True*, *\*\*kwargs*)

Set value of a frame grabber attribute with the given name.

If the value doesn't exist and *error\_on\_missing==True*, raise error; otherwise, do nothing. If *name* points at a dictionary branch, set all values in this branch (in this case *value* must be a dictionary). Additional arguments are passed to `set_value` methods of the individual attribute.

**get\_all\_grabber\_attribute\_values** (*root="*, *\*\*kwargs*)

Get values of all frame grabber attributes with the given *root*.

Additional arguments are passed to `get_value` methods of individual attributes.

**set\_all\_grabber\_attribute\_values** (*settings*, *root="*, *\*\*kwargs*)

Set values of all frame grabber attributes with the given *root*.

Additional arguments are passed to `set_value` methods of individual attributes.

**Error**

alias of `pylablib.core.devio.base.DeviceError`

**FrameTransferError**

alias of `DefaultFrameTransferError`

**class NoParameterCaller** (*device*, *kind*)

Bases: `object`

Class to simplify calling functions without a parameter

**TimeoutError**

alias of `pylablib.core.devio.base.DeviceError`

**acquisition\_in\_progress()**  
Check if acquisition is in progress

**apply\_settings(settings)**  
Apply the settings.  
  
*settings* is the dict {name: value} of the device available settings. Non-applicable settings are ignored.

**clear\_acquisition()**  
Clear acquisition settings

**close()**  
Close the connection

**get\_acquisition\_parameters()**  
Get acquisition parameters.  
  
Return dictionary {name: value}

**get\_data\_dimensions()**  
Get readout data dimensions (in pixels) as a tuple (width, height); take indexing mode into account

**get\_detector\_size()**  
Get camera detector size (in pixels) as a tuple (width, height)

**get\_device\_variable(key)**  
Get the value of a settings, status, or full info parameter

**get\_frame\_format()**  
Get format for the returned images.  
  
Can be "list" (list of 2D arrays), "array" (a single 3D array), or "chunks" (list of 3D "chunk" arrays; supported for some cameras and provides the best performance).

**get\_frame\_info\_fields()**  
Get the names of frame info fields.  
  
Applicable when frame info format (set by `set_frame_info_format()`) is "list" or "array".

**get\_frame\_info\_format()**  
Get format of the frame info.  
  
Can be "namedtuple" (potentially nested named tuples; convenient to get particular values), "list" (flat list of values, with field names are given by `get_frame_info_fields()`; convenient for building a table), "array" (same as "list", but with a numpy array, which is easier to use for "chunks" frame format), or "dict" (flat dictionary with the same fields as the "list" format; more resilient to future format changes)

**get\_frame\_info\_period()**  
Get period of frame info acquisition.  
  
Frame info might be skipped (set to None) except for frames which indices are divisible by *period*. Useful for certain cameras where acquiring frame info takes a lot of time and can reduce performance at higher frame rates. Note that this parameter can still be ignored (i.e., always set to 1) if the performance is not an issue for a given camera class.

**get\_frames\_status()**  
Get acquisition and buffer status.  
  
Return tuple (acquired, unread, skipped, size), where *acquired* is the total number of acquired frames, *unread* is the number of acquired but not read frames, *skipped* is the number of skipped (not read and then written over) frames, and *buffer\_size* is the total buffer size (in frames).

**get\_full\_info** (*include=0*)

Get dict {name: value} containing full device information (including status and settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_full\_status** (*include=0*)

Get dict {name: value} containing the device status (including settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_image\_indexing** ()

Get indexing for the returned images.

Can be "rct" (first index row, second index column, rows counted from the top), "rcb" (same as "rc", rows counted from the bottom), "xyt" (first index column, second index row, rows counted from the top), or "xyb" (same as "xyt", rows counted from the bottom)

**get\_new\_images\_range** ()

Get the range of the new images.

Return tuple (*first*, *last*) with images range (first inclusive). If no images are available, return None. If some images were in the buffer were overwritten, exclude them from the range.

**get\_settings** (*include=0*)

Get dict {name: value} containing all the device settings.

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**grab** (*nframes=1*, *frame\_timeout=5.0*, *missing\_frame='skip'*, *return\_info=False*, *buff\_size=None*)

Snap *nframes* images (with preset image read mode parameters)

*buff\_size* determines buffer size (if None, use the default size). Timeout is specified for a single-frame acquisition, not for the whole acquisition time. *missing\_frame* determines what to do with frames which have been lost: can be "none" (replacing them with None), "zero" (replacing them with zero-filled frame), or "skip" (skipping them, while still keeping total returned frames number to *n*). If *return\_info==True*, return tuple (*frames*, *infos*), where *infos* is a list of frame info tuples (camera-dependent); if some frames are missing and *missing\_frame!="skip"*, the corresponding frame info is None.

**is\_acquisition\_setup** ()

Check if acquisition is set up.

If the camera does not support separate acquisition setup, always return True.

**is\_opened** ()

Check if the device is connected

**open** ()

Open the connection

**pausing\_acquisition** (*clear=None*, *stop=True*, *setup\_after=None*, *start\_after=True*, *combine\_nested=True*)

Context manager which temporarily pauses acquisition during execution of `with` block.



Useful for applying certain settings which can't be changed during the acquisition. If `clear==True`, clear acquisition in addition to stopping (by default, use the class default specified as `_clear_pausing_acquisition` attribute). If `stop==True`, stop the acquisition (if `clear==True`, stop regardless). If `setup_after==True`, setup the acquisition after pause if necessary (None means setup only if clearing was required). If `start_after==True`, start the acquisition after pause if necessary (None means start only if stopping was required). If `combine_nested==True`, then any nested `pausing_acquisition` calls will stop/clear acquisition as necessary, but won't setup/start it again until this `pausing_acquisition` call is complete.

Yields tuple (`acq_in_progress`, `acq_params`), which indicates whether acquisition is currently in progress, and what are the current acquisition parameters.

**read\_multiple\_images** (*rng=None, peek=False, missing\_frame='skip', return\_info=False, return\_rng=False*)

Read multiple images specified by *rng* (by default, all un-read images).

If *rng* is specified, it is a tuple (*first*, *last*) with images range (*first* inclusive). If no new frames are available, return an empty list; if no acquisition is running, return None. If `peek==True`, return images but not mark them as read. *missing\_frame* determines what to do with frames which are out of range (missing or lost): can be "none" (replacing them with None), "zero" (replacing them with zero-filled frame), or "skip" (skipping them). If `return_info==True`, return tuple (*frames*, *infos*), where *infos* is a list of frame info tuples (camera-dependent, by default, only the frame index); if some frames are missing and `missing_frame!="skip"`, the corresponding frame info is None. If `return_rng==True`, return the range covered resulting frames; if `missing_frame=="skip"`, the range can be smaller than the supplied *rng* if some frames are skipped.

**read\_newest\_image** (*peek=False, return\_info=False*)

Read the newest un-read image.

If no un-read frames are available, return None. If `peek==True`, return the image but not mark it as read. If `return_info==True`, return tuple (*frame*, *info*), where *info* is an info tuples (camera-dependent, see `read_multiple_images()`).

**read\_oldest\_image** (*peek=False, return\_info=False*)

Read the oldest un-read image.

If no un-read frames are available, return None. If `peek==True`, return the image but not mark it as read. If `return_info==True`, return tuple (*frame*, *info*), where *info* is an info tuples (camera-dependent, see `read_multiple_images()`).

**set\_device\_variable** (*key, value*)

Set the value of a settings parameter

**set\_frame\_format** (*fnt*)

Set format for the returned images.

Can be "list" (list of 2D arrays), "array" (a single 3D array), "chunks" (list of 3D "chunk" arrays; supported for some cameras and provides the best performance), or "try\_chunks" (same as "chunks", but if chunks are not supported, set to "list" instead). If format is "chunks" and chunks are not supported by the camera, it results in one frame per chunk. Note that if the format is set to "array" or "chunks", the frame info format is also automatically set to "array". If the format is set to "chunks", then the image info is also returned in chunks form (list of 2D info arrays with the same length as the corresponding frame chunks).

**set\_frame\_info\_format** (*fnt, include\_fields=None*)

Set format of the frame info.

Can be "namedtuple" (potentially nested named tuples; convenient to get particular values), "list" (flat list of values, with field names are given by `get_frame_info_fields()`; convenient for building a table), "array" (same as "list", but with a numpy array, which is easier to use for "chunks"

frame format), or "dict" (flat dictionary with the same fields as the "list" format; more resilient to future format changes) If *include\_fields* is not None, it specifies the fields included for non-"tuple" formats; note that order or *include\_fields* is ignored, and the resulting fields are always ordered same as in the original.

**set\_frame\_info\_period** (*period=1*)

Set period of frame info acquisition.

Frame info might be skipped (set to None) except for frames which indices are divisible by *period*. Useful for certain cameras where acquiring frame info takes a lot of time and can reduce performance at higher frame rates. Note that this parameter can still be ignored (i.e., always set to 1) if the performance is not an issue for a given camera class.

**set\_image\_indexing** (*indexing*)

Set up indexing for the returned images.

Can be "rct" (first index row, second index column, rows counted from the top), "rcb" (same as "rc", rows counted from the bottom), "xyt" (first index column, second index row, rows counted from the top), or "xyb" (same as "xyt", rows counted from the bottom)

**setup\_acquisition** (*\*\*kwargs*)

Setup acquisition.

Any non-specified acquisition parameters are assumed to be the same as previously set (or default, if not explicitly set before). Return the new acquisition parameters.

**snap** (*timeout=5.0, return\_info=False*)

Snap a single frame

**start\_acquisition** (*\*args, \*\*kwargs*)

Start acquisition.

Can take the same keyword parameters as :meth:“*setup\_acquisition*”. If the acquisition is not set up yet, set it up using the supplied parameters (use default of *setup\_acquisition()*, if the parameter is None). Otherwise, if any supplied parameters are different from the current ones, change them and reset the acquisition.

**stop\_acquisition** ()

Stop acquisition

**wait\_for\_frame** (*since='lastread', nframes=1, timeout=20.0, error\_on\_stopped=False*)

Wait for one or several new camera frames.

*since* specifies the reference point for waiting to acquire *nframes* frames; can be “lastread” (from the last read frame), “lastwait” (wait for the last successful *wait\_for\_frame()* call), “now” (from the start of the current call), or “start” (from the acquisition start, i.e., wait until *nframes* frames have been acquired). *timeout* can be either a number, None (infinite timeout), or a tuple (*timeout, frame\_timeout*), in which case the call times out if the total time exceeds *timeout*, or a single frame wait exceeds *frame\_timeout*. If the call times out, raise *TimeoutError*. If *error\_on\_stopped==True* and the acquisition is not running, raise *Error*; otherwise, simply return False without waiting.

**class** pylablib.devices.interface.camera.**TAcqTimings** (*exposure, frame\_period*)

Bases: *tuple*

**count** ()

Return number of occurrences of value.

**exposure**

**frame\_period**

```

index()
    Return first index of value.

    Raises ValueError if the value is not present.

class pylablib.devices.interface.camera.IExposureCamera (*args, **kwargs)
    Bases: pylablib.devices.interface.camera.ICamera

    get_exposure()
        Get current exposure

    set_exposure(exposure)
        Set camera exposure

    get_frame_period()
        Get frame period (time between two consecutive frames in the internal trigger mode)

    get_frame_timings()
        Get acquisition timing.

        Return tuple (exposure, frame_period).

Error
    alias of pylablib.core.devio.base.DeviceError

FrameTransferError
    alias of DefaultFrameTransferError

class NoParameterCaller (device, kind)
    Bases: object

    Class to simplify calling functions without a parameter

TimeoutError
    alias of pylablib.core.devio.base.DeviceError

acquisition_in_progress()
    Check if acquisition is in progress

apply_settings(settings)
    Apply the settings.

    settings is the dict {name: value} of the device available settings. Non-applicable settings are ignored.

clear_acquisition()
    Clear acquisition settings

close()
    Close the connection

get_acquisition_parameters()
    Get acquisition parameters.

    Return dictionary {name: value}

get_data_dimensions()
    Get readout data dimensions (in pixels) as a tuple (width, height); take indexing mode into account

get_detector_size()
    Get camera detector size (in pixels) as a tuple (width, height)

get_device_variable(key)
    Get the value of a settings, status, or full info parameter

```

**get\_frame\_format()**

Get format for the returned images.

Can be "list" (list of 2D arrays), "array" (a single 3D array), or "chunks" (list of 3D "chunk" arrays; supported for some cameras and provides the best performance).

**get\_frame\_info\_fields()**

Get the names of frame info fields.

Applicable when frame info format (set by `set_frame_info_format()`) is "list" or "array".

**get\_frame\_info\_format()**

Get format of the frame info.

Can be "namedtuple" (potentially nested named tuples; convenient to get particular values), "list" (flat list of values, with field names are given by `get_frame_info_fields()`; convenient for building a table), "array" (same as "list", but with a numpy array, which is easier to use for "chunks" frame format), or "dict" (flat dictionary with the same fields as the "list" format; more resilient to future format changes)

**get\_frame\_info\_period()**

Get period of frame info acquisition.

Frame info might be skipped (set to `None`) except for frames which indices are divisible by *period*. Useful for certain cameras where acquiring frame info takes a lot of time and can reduce performance at higher frame rates. Note that this parameter can still be ignored (i.e., always set to 1) if the performance is not an issue for a given camera class.

**get\_frames\_status()**

Get acquisition and buffer status.

Return tuple (acquired, unread, skipped, size), where acquired is the total number of acquired frames, unread is the number of acquired but not read frames, skipped is the number of skipped (not read and then written over) frames, and `buffer_size` is the total buffer size (in frames).

**get\_full\_info(include=0)**

Get dict {name: value} containing full device information (including status and settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting `include=-10` queries all available variables, which is equivalent to `include="all"`.

**get\_full\_status(include=0)**

Get dict {name: value} containing the device status (including settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting `include=-10` queries all available variables, which is equivalent to `include="all"`.

**get\_image\_indexing()**

Get indexing for the returned images.

Can be "rct" (first index row, second index column, rows counted from the top), "rcb" (same as "rc", rows counted from the bottom), "xyt" (first index column, second index row, rows counted from the top), or "xyb" (same as "xyt", rows counted from the bottom)

**get\_new\_images\_range()**

Get the range of the new images.

Return tuple (*first*, *last*) with images range (first inclusive). If no images are available, return None. If some images were in the buffer were overwritten, exclude them from the range.

**get\_settings** (*include=0*)

Get dict {*name*: *value*} containing all the device settings.

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**grab** (*nframes=1*, *frame\_timeout=5.0*, *missing\_frame='skip'*, *return\_info=False*, *buff\_size=None*)

Snap *nframes* images (with preset image read mode parameters)

*buff\_size* determines buffer size (if None, use the default size). Timeout is specified for a single-frame acquisition, not for the whole acquisition time. *missing\_frame* determines what to do with frames which have been lost: can be "none" (replacing them with None), "zero" (replacing them with zero-filled frame), or "skip" (skipping them, while still keeping total returned frames number to *n*). If *return\_info==True*, return tuple (*frames*, *infos*), where *infos* is a list of frame info tuples (camera-dependent); if some frames are missing and *missing\_frame!="skip"*, the corresponding frame info is None.

**is\_acquisition\_setup** ()

Check if acquisition is set up.

If the camera does not support separate acquisition setup, always return True.

**is\_opened** ()

Check if the device is connected

**open** ()

Open the connection

**pausing\_acquisition** (*clear=None*, *stop=True*, *setup\_after=None*, *start\_after=True*, *combine\_nested=True*)

Context manager which temporarily pauses acquisition during execution of *with* block.

Useful for applying certain settings which can't be changed during the acquisition. If *clear==True*, clear acquisition in addition to stopping (by default, use the class default specified as *\_clear\_pausing\_acquisition* attribute). If *stop==True*, stop the acquisition (if *clear==True*, stop regardless). If *setup\_after==True*, setup the acquisition after pause if necessary (None means setup only if clearing was required). If *start\_after==True*, start the acquisition after pause if necessary (None means start only if stopping was required). If *combine\_nested==True*, then any nested *pausing\_acquisition* calls will stop/clear acquisition as necessary, but won't setup/start it again until this *pausing\_acquisition* call is complete.

Yields tuple (*acq\_in\_progress*, *acq\_params*), which indicates whether acquisition is currently in progress, and what are the current acquisition parameters.

**read\_multiple\_images** (*rng=None*, *peek=False*, *missing\_frame='skip'*, *return\_info=False*, *return\_rng=False*)

Read multiple images specified by *rng* (by default, all un-read images).

If *rng* is specified, it is a tuple (*first*, *last*) with images range (first inclusive). If no new frames are available, return an empty list; if no acquisition is running, return None. If *peek==True*, return images but not mark them as read. *missing\_frame* determines what to do with frames which are out of range (missing or lost): can be "none" (replacing them with None), "zero" (replacing them with zero-filled frame), or "skip" (skipping them). If *return\_info==True*, return tuple (*frames*, *infos*), where *infos* is a list of frame info tuples (camera-dependent, by default, only the frame index); if some frames are missing and *missing\_frame!="skip"*, the corresponding frame info is None. if

`return_rng==True`, return the range covered resulting frames; if `missing_frame=="skip"`, the range can be smaller than the supplied `rng` if some frames are skipped.

**read\_newest\_image** (*peek=False, return\_info=False*)

Read the newest un-read image.

If no un-read frames are available, return `None`. If `peek==True`, return the image but not mark it as read. If `return_info==True`, return tuple (`frame`, `info`), where `info` is an info tuples (camera-dependent, see [read\\_multiple\\_images\(\)](#)).

**read\_oldest\_image** (*peek=False, return\_info=False*)

Read the oldest un-read image.

If no un-read frames are available, return `None`. If `peek==True`, return the image but not mark it as read. If `return_info==True`, return tuple (`frame`, `info`), where `info` is an info tuples (camera-dependent, see [read\\_multiple\\_images\(\)](#)).

**set\_device\_variable** (*key, value*)

Set the value of a settings parameter

**set\_frame\_format** (*fmt*)

Set format for the returned images.

Can be `"list"` (list of 2D arrays), `"array"` (a single 3D array), `"chunks"` (list of 3D “chunk” arrays; supported for some cameras and provides the best performance), or `"try_chunks"` (same as `"chunks"`, but if chunks are not supported, set to `"list"` instead). If format is `"chunks"` and chunks are not supported by the camera, it results in one frame per chunk. Note that if the format is set to `"array"` or `"chunks"`, the frame info format is also automatically set to `"array"`. If the format is set to `"chunks"`, then the image info is also returned in chunks form (list of 2D info arrays with the same length as the corresponding frame chunks).

**set\_frame\_info\_format** (*fmt, include\_fields=None*)

Set format of the frame info.

Can be `"namedtuple"` (potentially nested named tuples; convenient to get particular values), `"list"` (flat list of values, with field names are given by [get\\_frame\\_info\\_fields\(\)](#); convenient for building a table), `"array"` (same as `"list"`, but with a numpy array, which is easier to use for `"chunks"` frame format), or `"dict"` (flat dictionary with the same fields as the `"list"` format; more resilient to future format changes) If `include_fields` is not `None`, it specifies the fields included for non-`"tuple"` formats; note that order or `include_fields` is ignored, and the resulting fields are always ordered same as in the original.

**set\_frame\_info\_period** (*period=1*)

Set period of frame info acquisition.

Frame info might be skipped (set to `None`) except for frames which indices are divisible by *period*. Useful for certain cameras where acquiring frame info takes a lot of time and can reduce performance at higher frame rates. Note that this parameter can still be ignored (i.e., always set to 1) if the performance is not an issue for a given camera class.

**set\_image\_indexing** (*indexing*)

Set up indexing for the returned images.

Can be `"rct"` (first index row, second index column, rows counted from the top), `"rcb"` (same as `"rc"`, rows counted from the bottom), `"xyt"` (first index column, second index row, rows counted from the top), or `"xyb"` (same as `"xyt"`, rows counted from the bottom)

**setup\_acquisition** (*\*\*kwargs*)

Setup acquisition.

Any non-specified acquisition parameters are assumed to be the same as previously set (or default, if not explicitly set before). Return the new acquisition parameters.

**snap** (*timeout=5.0, return\_info=False*)

Snap a single frame

**start\_acquisition** (*\*args, \*\*kwargs*)

Start acquisition.

Can take the same keyword parameters as :meth:“*setup\_acquisition*”. If the acquisition is not set up yet, set it up using the supplied parameters (use default of *setup\_acquisition()*, if the parameter is None). Otherwise, if any supplied parameters are different from the current ones, change them and reset the acquisition.

**stop\_acquisition** ()

Stop acquisition

**wait\_for\_frame** (*since='lastread', nframes=1, timeout=20.0, error\_on\_stopped=False*)

Wait for one or several new camera frames.

*since* specifies the reference point for waiting to acquire *nframes* frames; can be “lastread” (from the last read frame), “lastwait” (wait for the last successful *wait\_for\_frame()* call), “now” (from the start of the current call), or “start” (from the acquisition start, i.e., wait until *nframes* frames have been acquired). *timeout* can be either a number, None (infinite timeout), or a tuple (timeout, frame\_timeout), in which case the call times out if the total time exceeds timeout, or a single frame wait exceeds frame\_timeout. If the call times out, raise *TimeoutError*. If *error\_on\_stopped==True* and the acquisition is not running, raise *Error*; otherwise, simply return False without waiting.

**class** `pylablib.devices.interface.camera.TAxisROILimit` (*min, max, pstep, sstep, maxbin*)

Bases: `tuple`

**count** ()

Return number of occurrences of value.

**index** ()

Return first index of value.

Raises *ValueError* if the value is not present.

**max**

**maxbin**

**min**

**pstep**

**sstep**

`pylablib.devices.interface.camera.truncate_roi_axis` (*roi, lim, symmetric=False*)

Truncate ROI to conform to the given ROI limits.

*roi* is a tuple (start, stop, bin), and *lim* is a tuple (min, max, pstep, sstep, maxbin). Assume that pstep and sstep divide min and max, and that either pstep divides sstep or the other way around. If *symmetric==True*, then max should be even.

**class** `pylablib.devices.interface.camera.IROICamera` (*\*args, \*\*kwargs*)

Bases: `pylablib.devices.interface.camera.ICamera`

**get\_roi** ()

Get current ROI.

Return tuple (*hstart*, *hend*, *vstart*, *vend*). *hstart* and *hend* specify horizontal image extent, *vstart* and *vend* specify vertical image extent (start is inclusive, stop is exclusive, starting from 0).

**set\_roi** (*hstart=0*, *hend=None*, *vstart=0*, *vend=None*)

Setup camera ROI.

*hstart* and *hend* specify horizontal image extent, *vstart* and *vend* specify vertical image extent (start is inclusive, stop is exclusive, starting from 0). By default, all non-supplied parameters take extreme values (0 for start, maximal for end).

**get\_roi\_limits** (*hbin=1*, *vbin=1*)

Get the minimal and maximal ROI parameters.

Return tuple (*hlim*, *vlim*), where each element is in turn a limit 5-tuple (*min*, *max*, *pstep*, *sstep*, *maxbin*) with, correspondingly, minimal and maximal size, position and size step, and the maximal binning (fixed to 1 if not binning is allowed). In some cameras, the step and the minimal size depend on the binning, which can be supplied.

#### Error

alias of `pylablib.core.devio.base.DeviceError`

#### FrameTransferError

alias of `DefaultFrameTransferError`

**class NoParameterCaller** (*device*, *kind*)

Bases: `object`

Class to simplify calling functions without a parameter

#### TimeoutError

alias of `pylablib.core.devio.base.DeviceError`

**acquisition\_in\_progress** ()

Check if acquisition is in progress

**apply\_settings** (*settings*)

Apply the settings.

*settings* is the dict {*name*: *value*} of the device available settings. Non-applicable settings are ignored.

**clear\_acquisition** ()

Clear acquisition settings

**close** ()

Close the connection

**get\_acquisition\_parameters** ()

Get acquisition parameters.

Return dictionary {*name*: *value*}

**get\_data\_dimensions** ()

Get readout data dimensions (in pixels) as a tuple (*width*, *height*); take indexing mode into account

**get\_detector\_size** ()

Get camera detector size (in pixels) as a tuple (*width*, *height*)

**get\_device\_variable** (*key*)

Get the value of a settings, status, or full info parameter

**get\_frame\_format** ()

Get format for the returned images.



Can be "list" (list of 2D arrays), "array" (a single 3D array), or "chunks" (list of 3D "chunk" arrays; supported for some cameras and provides the best performance).

**get\_frame\_info\_fields()**

Get the names of frame info fields.

Applicable when frame info format (set by `set_frame_info_format()`) is "list" or "array".

**get\_frame\_info\_format()**

Get format of the frame info.

Can be "namedtuple" (potentially nested named tuples; convenient to get particular values), "list" (flat list of values, with field names are given by `get_frame_info_fields()`; convenient for building a table), "array" (same as "list", but with a numpy array, which is easier to use for "chunks" frame format), or "dict" (flat dictionary with the same fields as the "list" format; more resilient to future format changes)

**get\_frame\_info\_period()**

Get period of frame info acquisition.

Frame info might be skipped (set to None) except for frames which indices are divisible by *period*. Useful for certain cameras where acquiring frame info takes a lot of time and can reduce performance at higher frame rates. Note that this parameter can still be ignored (i.e., always set to 1) if the performance is not an issue for a given camera class.

**get\_frames\_status()**

Get acquisition and buffer status.

Return tuple (*acquired*, *unread*, *skipped*, *size*), where *acquired* is the total number of acquired frames, *unread* is the number of acquired but not read frames, *skipped* is the number of skipped (not read and then written over) frames, and *buffer\_size* is the total buffer size (in frames).

**get\_full\_info(include=0)**

Get dict {name: value} containing full device information (including status and settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include*=-10 queries all available variables, which is equivalent to *include*="all".

**get\_full\_status(include=0)**

Get dict {name: value} containing the device status (including settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include*=-10 queries all available variables, which is equivalent to *include*="all".

**get\_image\_indexing()**

Get indexing for the returned images.

Can be "rct" (first index row, second index column, rows counted from the top), "rcb" (same as "rc", rows counted from the bottom), "xyt" (first index column, second index row, rows counted from the top), or "xyb" (same as "xyt", rows counted from the bottom)

**get\_new\_images\_range()**

Get the range of the new images.

Return tuple (*first*, *last*) with images range (first inclusive). If no images are available, return None. If some images were in the buffer were overwritten, exclude them from the range.

**get\_settings** (*include=0*)

Get dict {name: value} containing all the device settings.

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**grab** (*nframes=1, frame\_timeout=5.0, missing\_frame='skip', return\_info=False, buff\_size=None*)

Snap *nframes* images (with preset image read mode parameters)

*buff\_size* determines buffer size (if *None*, use the default size). Timeout is specified for a single-frame acquisition, not for the whole acquisition time. *missing\_frame* determines what to do with frames which have been lost: can be "none" (replacing them with *None*), "zero" (replacing them with zero-filled frame), or "skip" (skipping them, while still keeping total returned frames number to *n*). If *return\_info==True*, return tuple (*frames, infos*), where *infos* is a list of frame info tuples (camera-dependent); if some frames are missing and *missing\_frame!="skip"*, the corresponding frame info is *None*.

**is\_acquisition\_setup** ()

Check if acquisition is set up.

If the camera does not support separate acquisition setup, always return *True*.

**is\_opened** ()

Check if the device is connected

**open** ()

Open the connection

**pausing\_acquisition** (*clear=None, stop=True, setup\_after=None, start\_after=True, combine\_nested=True*)

Context manager which temporarily pauses acquisition during execution of *with* block.

Useful for applying certain settings which can't be changed during the acquisition. If *clear==True*, clear acquisition in addition to stopping (by default, use the class default specified as *\_clear\_pausing\_acquisition* attribute). If *stop==True*, stop the acquisition (if *clear==True*, stop regardless). If *setup\_after==True*, setup the acquisition after pause if necessary (*None* means setup only if clearing was required). If *start\_after==True*, start the acquisition after pause if necessary (*None* means start only if stopping was required). If *combine\_nested==True*, then any nested *pausing\_acquisition* calls will stop/clear acquisition as necessary, but won't setup/start it again until this *pausing\_acquisition* call is complete.

Yields tuple (*acq\_in\_progress, acq\_params*), which indicates whether acquisition is currently in progress, and what are the current acquisition parameters.

**read\_multiple\_images** (*rng=None, peek=False, missing\_frame='skip', return\_info=False, return\_rng=False*)

Read multiple images specified by *rng* (by default, all un-read images).

If *rng* is specified, it is a tuple (*first, last*) with images range (first inclusive). If no new frames are available, return an empty list; if no acquisition is running, return *None*. If *peek==True*, return images but not mark them as read. *missing\_frame* determines what to do with frames which are out of range (missing or lost): can be "none" (replacing them with *None*), "zero" (replacing them with zero-filled frame), or "skip" (skipping them). If *return\_info==True*, return tuple (*frames, infos*), where *infos* is a list of frame info tuples (camera-dependent, by default, only the frame index); if some frames are missing and *missing\_frame!="skip"*, the corresponding frame info is *None*. if *return\_rng==True*, return the range covered resulting frames; if *missing\_frame=="skip"*, the range can be smaller than the supplied *rng* if some frames are skipped.

**read\_newest\_image** (*peek=False, return\_info=False*)

Read the newest un-read image.

If no un-read frames are available, return `None`. If `peek==True`, return the image but not mark it as read. If `return_info==True`, return tuple (`frame`, `info`), where `info` is an info tuples (camera-dependent, see [read\\_multiple\\_images\(\)](#)).

**read\_oldest\_image** (*peek=False, return\_info=False*)

Read the oldest un-read image.

If no un-read frames are available, return `None`. If `peek==True`, return the image but not mark it as read. If `return_info==True`, return tuple (`frame`, `info`), where `info` is an info tuples (camera-dependent, see [read\\_multiple\\_images\(\)](#)).

**set\_device\_variable** (*key, value*)

Set the value of a settings parameter

**set\_frame\_format** (*fmt*)

Set format for the returned images.

Can be "list" (list of 2D arrays), "array" (a single 3D array), "chunks" (list of 3D "chunk" arrays; supported for some cameras and provides the best performance), or "try\_chunks" (same as "chunks", but if chunks are not supported, set to "list" instead). If format is "chunks" and chunks are not supported by the camera, it results in one frame per chunk. Note that if the format is set to "array" or "chunks", the frame info format is also automatically set to "array". If the format is set to "chunks", then the image info is also returned in chunks form (list of 2D info arrays with the same length as the corresponding frame chunks).

**set\_frame\_info\_format** (*fmt, include\_fields=None*)

Set format of the frame info.

Can be "namedtuple" (potentially nested named tuples; convenient to get particular values), "list" (flat list of values, with field names are given by [get\\_frame\\_info\\_fields\(\)](#); convenient for building a table), "array" (same as "list", but with a numpy array, which is easier to use for "chunks" frame format), or "dict" (flat dictionary with the same fields as the "list" format; more resilient to future format changes) If `include_fields` is not `None`, it specifies the fields included for non-"tuple" formats; note that order of `include_fields` is ignored, and the resulting fields are always ordered same as in the original.

**set\_frame\_info\_period** (*period=1*)

Set period of frame info acquisition.

Frame info might be skipped (set to `None`) except for frames which indices are divisible by `period`. Useful for certain cameras where acquiring frame info takes a lot of time and can reduce performance at higher frame rates. Note that this parameter can still be ignored (i.e., always set to 1) if the performance is not an issue for a given camera class.

**set\_image\_indexing** (*indexing*)

Set up indexing for the returned images.

Can be "rct" (first index row, second index column, rows counted from the top), "rcb" (same as "rc", rows counted from the bottom), "xyt" (first index column, second index row, rows counted from the top), or "xyb" (same as "xyt", rows counted from the bottom)

**setup\_acquisition** (*\*\*kwargs*)

Setup acquisition.

Any non-specified acquisition parameters are assumed to be the same as previously set (or default, if not explicitly set before). Return the new acquisition parameters.

**snap** (*timeout=5.0, return\_info=False*)  
Snap a single frame

**start\_acquisition** (*\*args, \*\*kwargs*)  
Start acquisition.

Can take the same keyword parameters as *meth: "setup\_acquisition"*. If the acquisition is not set up yet, set it up using the supplied parameters (use default of *setup\_acquisition()*, if the parameter is None). Otherwise, if any supplied parameters are different from the current ones, change them and reset the acquisition.

**stop\_acquisition** ()  
Stop acquisition

**wait\_for\_frame** (*since='lastread', nframes=1, timeout=20.0, error\_on\_stopped=False*)  
Wait for one or several new camera frames.

*since* specifies the reference point for waiting to acquire *nframes* frames; can be "lastread" (from the last read frame), "lastwait" (wait for the last successful *wait\_for\_frame()* call), "now" (from the start of the current call), or "start" (from the acquisition start, i.e., wait until *nframes* frames have been acquired). *timeout* can be either a number, None (infinite timeout), or a tuple (*timeout, frame\_timeout*), in which case the call times out if the total time exceeds *timeout*, or a single frame wait exceeds *frame\_timeout*. If the call times out, raise *TimeoutError*. If *error\_on\_stopped==True* and the acquisition is not running, raise *Error*; otherwise, simply return False without waiting.

**class** `pylablib.devices.interface.camera.IBinROICamera` (*\*args, \*\*kwargs*)

Bases: `pylablib.devices.interface.camera.ICamera`

**get\_roi** ()  
Get current ROI.

Return tuple (*hstart, hend, vstart, vend, hbin, vbin*). *hstart* and *hend* specify horizontal image extent, *vstart* and *vend* specify vertical image extent (start is inclusive, stop is exclusive, starting from 0), *hbin* and *vbin* specify binning.

**set\_roi** (*hstart=0, hend=None, vstart=0, vend=None, hbin=1, vbin=1*)  
Setup camera ROI.

*hstart* and *hend* specify horizontal image extent, *vstart* and *vend* specify vertical image extent (start is inclusive, stop is exclusive, starting from 0), *hbin* and *vbin* specify binning. By default, all non-supplied parameters take extreme values (0 for start, maximal for end, 1 for binning).

**get\_roi\_limits** (*hbin=1, vbin=1*)  
Get the minimal and maximal ROI parameters.

Return tuple (*hlim, vlim*), where each element is in turn a limit 5-tuple (*min, max, pstep, sstep, maxbin*) with, correspondingly, minimal and maximal size, position and size step, and the maximal binning. In some cameras, the step and the minimal size depend on the binning, which can be supplied.

**Error**  
alias of `pylablib.core.devio.base.DeviceError`

**FrameTransferError**  
alias of `DefaultFrameTransferError`

**class** `NoParameterCaller` (*device, kind*)  
Bases: `object`  
Class to simplify calling functions without a parameter

**TimeoutError**

alias of `pylablib.core.devio.base.DeviceError`

**acquisition\_in\_progress()**

Check if acquisition is in progress

**apply\_settings(settings)**

Apply the settings.

*settings* is the dict {name: value} of the device available settings. Non-applicable settings are ignored.

**clear\_acquisition()**

Clear acquisition settings

**close()**

Close the connection

**get\_acquisition\_parameters()**

Get acquisition parameters.

Return dictionary {name: value}

**get\_data\_dimensions()**

Get readout data dimensions (in pixels) as a tuple (width, height); take indexing mode into account

**get\_detector\_size()**

Get camera detector size (in pixels) as a tuple (width, height)

**get\_device\_variable(key)**

Get the value of a settings, status, or full info parameter

**get\_frame\_format()**

Get format for the returned images.

Can be "list" (list of 2D arrays), "array" (a single 3D array), or "chunks" (list of 3D "chunk" arrays; supported for some cameras and provides the best performance).

**get\_frame\_info\_fields()**

Get the names of frame info fields.

Applicable when frame info format (set by `set_frame_info_format()`) is "list" or "array".

**get\_frame\_info\_format()**

Get format of the frame info.

Can be "namedtuple" (potentially nested named tuples; convenient to get particular values), "list" (flat list of values, with field names are given by `get_frame_info_fields()`; convenient for building a table), "array" (same as "list", but with a numpy array, which is easier to use for "chunks" frame format), or "dict" (flat dictionary with the same fields as the "list" format; more resilient to future format changes)

**get\_frame\_info\_period()**

Get period of frame info acquisition.

Frame info might be skipped (set to None) except for frames which indices are divisible by *period*. Useful for certain cameras where acquiring frame info takes a lot of time and can reduce performance at higher frame rates. Note that this parameter can still be ignored (i.e., always set to 1) if the performance is not an issue for a given camera class.

**get\_frames\_status()**

Get acquisition and buffer status.

Return tuple (acquired, unread, skipped, size), where acquired is the total number of acquired frames, unread is the number of acquired but not read frames, skipped is the number of skipped (not read and then written over) frames, and buffer\_size is the total buffer size (in frames).

**get\_full\_info** (include=0)

Get dict {name: value} containing full device information (including status and settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting include=-10 queries all available variables, which is equivalent to include="all".

**get\_full\_status** (include=0)

Get dict {name: value} containing the device status (including settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting include=-10 queries all available variables, which is equivalent to include="all".

**get\_image\_indexing** ()

Get indexing for the returned images.

Can be "rct" (first index row, second index column, rows counted from the top), "rcb" (same as "rc", rows counted from the bottom), "xyt" (first index column, second index row, rows counted from the top), or "xyb" (same as "xyt", rows counted from the bottom)

**get\_new\_images\_range** ()

Get the range of the new images.

Return tuple (first, last) with images range (first inclusive). If no images are available, return None. If some images were in the buffer were overwritten, exclude them from the range.

**get\_settings** (include=0)

Get dict {name: value} containing all the device settings.

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting include=-10 queries all available variables, which is equivalent to include="all".

**grab** (nframes=1, frame\_timeout=5.0, missing\_frame='skip', return\_info=False, buff\_size=None)

Snap *nframes* images (with preset image read mode parameters)

*buff\_size* determines buffer size (if None, use the default size). Timeout is specified for a single-frame acquisition, not for the whole acquisition time. *missing\_frame* determines what to do with frames which have been lost: can be "none" (replacing them with None), "zero" (replacing them with zero-filled frame), or "skip" (skipping them, while still keeping total returned frames number to *n*). If return\_info==True, return tuple (frames, infos), where infos is a list of frame info tuples (camera-dependent); if some frames are missing and missing\_frame!="skip", the corresponding frame info is None.

**is\_acquisition\_setup** ()

Check if acquisition is set up.

If the camera does not support separate acquisition setup, always return True.

**is\_opened** ()

Check if the device is connected

**open** ()

Open the connection

**pausing\_acquisition** (*clear=None, stop=True, setup\_after=None, start\_after=True, combine\_nested=True*)

Context manager which temporarily pauses acquisition during execution of with block.

Useful for applying certain settings which can't be changed during the acquisition. If `clear==True`, clear acquisition in addition to stopping (by default, use the class default specified as `_clear_pausing_acquisition` attribute). If `stop==True`, stop the acquisition (if `clear==True`, stop regardless). If `setup_after==True`, setup the acquisition after pause if necessary (None means setup only if clearing was required). If `start_after==True`, start the acquisition after pause if necessary (None means start only if stopping was required). If `combine_nested==True`, then any nested `pausing_acquisition` calls will stop/clear acquisition as necessary, but won't setup/start it again until this `pausing_acquisition` call is complete.

Yields tuple (`acq_in_progress, acq_params`), which indicates whether acquisition is currently in progress, and what are the current acquisition parameters.

**read\_multiple\_images** (*rng=None, peek=False, missing\_frame='skip', return\_info=False, return\_rng=False*)

Read multiple images specified by *rng* (by default, all un-read images).

If *rng* is specified, it is a tuple (*first, last*) with images range (first inclusive). If no new frames are available, return an empty list; if no acquisition is running, return None. If `peek==True`, return images but not mark them as read. *missing\_frame* determines what to do with frames which are out of range (missing or lost): can be "none" (replacing them with None), "zero" (replacing them with zero-filled frame), or "skip" (skipping them). If `return_info==True`, return tuple (*frames, infos*), where *infos* is a list of frame info tuples (camera-dependent, by default, only the frame index); if some frames are missing and `missing_frame!="skip"`, the corresponding frame info is None. if `return_rng==True`, return the range covered resulting frames; if `missing_frame=="skip"`, the range can be smaller than the supplied *rng* if some frames are skipped.

**read\_newest\_image** (*peek=False, return\_info=False*)

Read the newest un-read image.

If no un-read frames are available, return None. If `peek==True`, return the image but not mark it as read. If `return_info==True`, return tuple (*frame, info*), where *info* is an info tuples (camera-dependent, see `read_multiple_images()`).

**read\_oldest\_image** (*peek=False, return\_info=False*)

Read the oldest un-read image.

If no un-read frames are available, return None. If `peek==True`, return the image but not mark it as read. If `return_info==True`, return tuple (*frame, info*), where *info* is an info tuples (camera-dependent, see `read_multiple_images()`).

**set\_device\_variable** (*key, value*)

Set the value of a settings parameter

**set\_frame\_format** (*fmt*)

Set format for the returned images.

Can be "list" (list of 2D arrays), "array" (a single 3D array), "chunks" (list of 3D "chunk" arrays; supported for some cameras and provides the best performance), or "try\_chunks" (same as "chunks", but if chunks are not supported, set to "list" instead). If format is "chunks" and chunks are not supported by the camera, it results in one frame per chunk. Note that if the format is set to "array" or "chunks", the frame info format is also automatically set to "array". If the format is set to "chunks", then the image info is also returned in chunks form (list of 2D info arrays with the same length as the corresponding frame chunks).

**set\_frame\_info\_format** (*fmt, include\_fields=None*)

Set format of the frame info.

Can be "namedtuple" (potentially nested named tuples; convenient to get particular values), "list" (flat list of values, with field names are given by `get_frame_info_fields()`; convenient for building a table), "array" (same as "list", but with a numpy array, which is easier to use for "chunks" frame format), or "dict" (flat dictionary with the same fields as the "list" format; more resilient to future format changes) If `include_fields` is not None, it specifies the fields included for non-"tuple" formats; note that order or `include_fields` is ignored, and the resulting fields are always ordered same as in the original.

**set\_frame\_info\_period** (*period=1*)

Set period of frame info acquisition.

Frame info might be skipped (set to None) except for frames which indices are divisible by *period*. Useful for certain cameras where acquiring frame info takes a lot of time and can reduce performance at higher frame rates. Note that this parameter can still be ignored (i.e., always set to 1) if the performance is not an issue for a given camera class.

**set\_image\_indexing** (*indexing*)

Set up indexing for the returned images.

Can be "rct" (first index row, second index column, rows counted from the top), "rcb" (same as "rc", rows counted from the bottom), "xyt" (first index column, second index row, rows counted from the top), or "xyb" (same as "xyt", rows counted from the bottom)

**setup\_acquisition** (*\*\*kwargs*)

Setup acquisition.

Any non-specified acquisition parameters are assumed to be the same as previously set (or default, if not explicitly set before). Return the new acquisition parameters.

**snap** (*timeout=5.0, return\_info=False*)

Snap a single frame

**start\_acquisition** (*\*args, \*\*kwargs*)

Start acquisition.

Can take the same keyword parameters as `:meth:~.setup_acquisition`. If the acquisition is not set up yet, set it up using the supplied parameters (use default of `setup_acquisition()`, if the parameter is None). Otherwise, if any supplied parameters are different from the current ones, change them and reset the acquisition.

**stop\_acquisition** ()

Stop acquisition

**wait\_for\_frame** (*since='lastread', nframes=1, timeout=20.0, error\_on\_stopped=False*)

Wait for one or several new camera frames.

*since* specifies the reference point for waiting to acquire *nframes* frames; can be "lastread" (from the last read frame), "lastwait" (wait for the last successful `wait_for_frame()` call), "now" (from the start of the current call), or "start" (from the acquisition start, i.e., wait until *nframes* frames have been acquired). *timeout* can be either a number, None (infinite timeout), or a tuple (*timeout, frame\_timeout*), in which case the call times out if the total time exceeds *timeout*, or a single frame wait exceeds *frame\_timeout*. If the call times out, raise `TimeoutError`. If `error_on_stopped==True` and the acquisition is not running, raise `Error`; otherwise, simply return False without waiting.

```
class pylablib.devices.interface.camera.TStatusLineDescription(kind, roi,  
                                                             frames-  
                                                             tamp_checker)
```

Bases: tuple



**count()**  
Return number of occurrences of value.

**framestamp\_checker**

**index()**  
Return first index of value.  
  
Raises ValueError if the value is not present.

**kind**

**roi**

**class** pylablib.devices.interface.camera.**StatusLineChecker**

Bases: `object`

Class responsible for checking status line consistency

**get\_framestamp**(*frames*)  
Get framestamps from status lines in the given frames

**check\_indices**(*indices*, *step=1*)  
Check if indices are consistent with the given step

`pylablib.devices.interface.camera.remove_status_line`(*frame*, *status\_line*, *policy='duplicate'*, *copy=True*, *value=0*)

Remove status line, if present.

#### Parameters

- **frame** – a frame to process (2D or 3D numpy array; if 3D, the first axis is the frame number)
- **status\_line** – status line descriptor (from the frames message)
- **policy** – determines way to deal with the status line; can be "keep" (keep as is), "cut" (cut off the status-line-containing row/column), "zero" (set it to zero), "value" (set it to a given value), "median" (set it to the image median), or "duplicate" (set it equal to the previous row; default) "cut" is only possible if the status line is on the edge of the image.
- **copy** – if True, make copy of the original frames; otherwise, attempt to remove the line in-place

`pylablib.devices.interface.camera.extract_status_line`(*frame*, *status\_line*, *copy=True*)

Extract status line, if present.

#### Parameters

- **frame** – a frame to process (2D or 3D numpy array; if 3D, the first axis is the frame number)
- **status\_line** – status line descriptor (from the frames message)
- **copy** – if True, make copy of the original status line data.

`pylablib.devices.interface.camera.insert_status_line`(*frame*, *status\_line*, *value*, *copy=True*)

Insert status line, if present.

#### Parameters

- **frame** – a frame to process (2D or 3D numpy array; if 3D, the first axis is the frame number)
- **status\_line** – status line descriptor (from the frames message)
- **value** – status line value
- **copy** – if True, make copy of the original status line data.

`pylablib.devices.interface.camera.get_status_line_roi(frame, status_line)`

Return ROI taken by the status line in the given frame

## pylablib.devices.interface.stage module

**class** `pylablib.devices.interface.stage.IStage`

Bases: `pylablib.core.devio.interface.IDevice`

Generic stage class

**class** `NoParameterCaller(device, kind)`

Bases: `object`

Class to simplify calling functions without a parameter

**apply\_settings**(*settings*)

Apply the settings.

*settings* is the dict {name: value} of the device available settings. Non-applicable settings are ignored.

**close**()

Close the connection

**get\_device\_variable**(*key*)

Get the value of a settings, status, or full info parameter

**get\_full\_info**(*include=0*)

Get dict {name: value} containing full device information (including status and settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_full\_status**(*include=0*)

Get dict {name: value} containing the device status (including settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**get\_settings**(*include=0*)

Get dict {name: value} containing all the device settings.

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include=-10* queries all available variables, which is equivalent to *include="all"*.

**is\_opened**()

Check if the device is connected

**open()**  
Open the connection

**set\_device\_variable**(*key*, *value*)  
Set the value of a settings parameter

`pylablib.devices.interface.stage.muxaxis(*args, argname='axis', **kwargs)`  
Multiplex the function over its axis argument

**class** `pylablib.devices.interface.stage.IMultiaxisStage(*args, default_axis='all', **kwargs)`

Bases: `pylablib.devices.interface.stage.IStage`

Generic multiaxis stage class.

Has methods to assign and map axes and the axis device parameter.

**Parameters** **default\_axis** – default axis parameter value used when `axis=None` is provided

**get\_all\_axes()**  
Get the list of all available axes (taking mapping into account)

**remap\_axes**(*mapping*, *accept\_original=True*)  
Rename axes to the new labels.

*mapping* is the new axes mapping, which can be a list of new axes name (corresponding to the old axes in order returned by `get_all_axes()`), or a dictionary {*alias*: *original*} of the new axes aliases.

**class** `NoParameterCaller(device, kind)`  
Bases: `object`

Class to simplify calling functions without a parameter

**apply\_settings**(*settings*)  
Apply the settings.

*settings* is the dict {*name*: *value*} of the device available settings. Non-applicable settings are ignored.

**close()**  
Close the connection

**get\_device\_variable**(*key*)  
Get the value of a settings, status, or full info parameter

**get\_full\_info**(*include=0*)  
Get dict {*name*: *value*} containing full device information (including status and settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting `include=-10` queries all available variables, which is equivalent to `include="all"`.

**get\_full\_status**(*include=0*)  
Get dict {*name*: *value*} containing the device status (including settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting `include=-10` queries all available variables, which is equivalent to `include="all"`.

**get\_settings**(*include=0*)  
Get dict {*name*: *value*} containing all the device settings.

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting `include=-10` queries all available variables, which is equivalent to `include="all"`.

**is\_opened()**  
Check if the device is connected

**open()**  
Open the connection

**set\_device\_variable**(*key, value*)  
Set the value of a settings parameter

## Module contents

### pylablib.devices.uc480 package

#### Submodules

#### pylablib.devices.uc480.uc480 module

**class** pylablib.devices.uc480.uc480.**TCameraInfo**(*cam\_id, dev\_id, sens\_id, model, serial\_number, in\_use, status*)

Bases: `tuple`

**cam\_id**

**count()**  
Return number of occurrences of value.

**dev\_id**

**in\_use**

**index()**  
Return first index of value.  
  
Raises `ValueError` if the value is not present.

**model**

**sens\_id**

**serial\_number**

**status**

**pylablib.devices.uc480.uc480.list\_cameras**(*backend='uc480'*)  
List all uc480/ueye camera connections (interface kind and camera index).

*backend* is the camera DLL backend; can be either "uc480" for Thorlabs-associated cameras, or "ueye" for IDS uEye-associated cameras

**pylablib.devices.uc480.uc480.get\_cameras\_number**(*backend='uc480'*)  
Get the total number of connected uc480/ueye cameras.

*backend* is the camera DLL backend; can be either "uc480" for Thorlabs-associated cameras, or "ueye" for IDS uEye-associated cameras

```
pylablib.devices.uc480.uc480.find_by_serial(serial_number, backend='uc480')
```

Find device ID using its serial number.

*backend* is the camera DLL backend; can be either "uc480" for Thorlabs-associated cameras, or "ueye" for IDS uEye-associated cameras

```
class pylablib.devices.uc480.uc480.TDeviceInfo(cam_id, model, manufacturer, se-  
rial_number, usb_version, date,  
dll_version, camera_type)
```

Bases: `tuple`

**cam\_id**

**camera\_type**

**count**()

Return number of occurrences of value.

**date**

**dll\_version**

**index**()

Return first index of value.

Raises ValueError if the value is not present.

**manufacturer**

**model**

**serial\_number**

**usb\_version**

```
class pylablib.devices.uc480.uc480.TAcquiredFramesStatus(acquired, transfer_missed,  
frameskip_events)
```

Bases: `tuple`

**acquired**

**count**()

Return number of occurrences of value.

**frameskip\_events**

**index**()

Return first index of value.

Raises ValueError if the value is not present.

**transfer\_missed**

```
class pylablib.devices.uc480.uc480.TTimestamp(year, month, day, hour, minute, second,  
millisecond)
```

Bases: `tuple`

**count**()

Return number of occurrences of value.

**day**

**hour**

**index()**  
Return first index of value.  
Raises ValueError if the value is not present.

**millisecond**

**minute**

**month**

**second**

**year**

```
class pylablib.devices.uc480.uc480.TFrameInfo(frame_index, framestamp, timestamp,  
                                              timestamp_dev, size, io_status, flags)
```

Bases: `tuple`

**count()**  
Return number of occurrences of value.

**flags**

**frame\_index**

**framestamp**

**index()**  
Return first index of value.  
Raises ValueError if the value is not present.

**io\_status**

**size**

**timestamp**

**timestamp\_dev**

```
class pylablib.devices.uc480.uc480.UC480Camera(cam_id=0, roi_binning_mode='auto',  
                                              dev_id=None, backend='uc480')  
Bases: pylablib.devices.interface.camera.IBinROICamera, pylablib.devices.  
interface.camera.IExposureCamera
```

Thorlabs uc480 / IDS uEye camera.

#### Parameters

- **cam\_id** (*int*) – camera ID; use 0 to get the first available camera
- **roi\_binning\_mode** – determines whether binning in ROI refers to binning or sub-sampling; can be "bin", "subsample", or "auto" (since most cameras only support one, it will pick the one which has non-trivial value, or "bin" if both are available).
- **dev\_id** (*int*) – if None use *cam\_id* as a camera id (*cam\_id* field of the camera info returned by `list_cameras()`); otherwise, ignore value of *cam\_id* and use *dev\_id* as device id (*dev\_id* field of the camera info). The first method requires assigning camera IDs beforehand (otherwise IDs might overlap, in which case only one camera can be accessed), but the assigned IDs are permanent; the second method always has unique IDs, but they might change if the cameras are disconnected and reconnected. For a more reliable assignment, one can use `find_by_serial()` function to find device ID based on the camera serial number.

- **backend** – camera DLL backend; can be either "uc480" for Thorlabs-associated cameras, or "ueye" for IDS uEye-associated cameras

```

Error = <Mock name='mock.uc480Error' id='140096270832528'>
TimeoutError = <Mock spec='str' id='140096294936144'>
FrameTransferError = <Mock spec='str' id='140096294936336'>
static find_by_serial (serial_number, backend='uc480')

open ()
    Open connection to the camera

close ()
    Close connection to the camera

is_opened ()
    Check if the device is connected

get_device_info ()
    Get camera model data.

    Return tuple (model, manufacturer, serial_number, usb_version, date,
dll_version, camera_type).

get_camera_id ()
    Get the current camera id

set_camera_id (cam_id)
    Set the new camera id (stored in non-volatile memory, i.e., survives power cycling)

get_frame_timings ()
    Get acquisition timing.

    Return tuple (exposure, frame_period).

set_exposure (exposure)
    Set camera exposure

set_frame_period (frame_time)
    Set frame period (time between two consecutive frames in the internal trigger mode)

get_pixel_rate ()
    Get camera pixel rate (in Hz)

get_available_pixel_rates ()
    Get all available pixel rates (in Hz)

get_pixel_rates_range ()
    Get range of allowed pixel rates (in Hz).

    Return tuple (min, max, step) if minimal and maximal value, and a step.

set_pixel_rate (rate=None)
    Set camera pixel rate (in Hz)

    The rate is always rounded to the closest available. If rate is None, set the maximal possible rate.

get_all_color_modes ()
    Get a list of all available color modes

get_color_mode ()
    Get current color mode.

    For possible modes, see get\_all\_color\_modes\(\).

```

**set\_color\_mode** (*mode*)

Set current color mode.

For possible modes, see [`get\_all\_color\_modes\(\)`](#).

**get\_gains** ()

Get current gains.

Return tuple (*master*, *red*, *green*, *blue*) of corresponding gain factors.

**get\_max\_gains** ()

Get maximal gains.

Return tuple (*master*, *red*, *green*, *blue*) of corresponding maximal gain factors.

**set\_gains** (*master=None*, *red=None*, *green=None*, *blue=None*)

Set current gains.

If supplied value is *None*, keep it unchanged.

**get\_gain\_boost** ()

Check if gain boost is enabled

**set\_gain\_boost** (*enabled*)

Enable or disable gain boost

**setup\_acquisition** (*nframes=100*)

Setup acquisition.

*nframes* determines number of size of the ring buffer (by default, 100).

**clear\_acquisition** ()

Clear acquisition settings

**start\_acquisition** (\**args*, \*\**kwargs*)

Start acquisition.

Can take the same keyword parameters as :meth:“[`setup\_acquisition`](#). If the acquisition is not set up yet, set it up using the supplied parameters (use default of [`setup\_acquisition\(\)`](#), if the parameter is *None*). Otherwise, if any supplied parameters are different from the current ones, change them and reset the acquisition.

**stop\_acquisition** ()

Stop acquisition

**acquisition\_in\_progress** ()

Check if acquisition is in progress

**get\_frames\_status** ()

Get acquisition and buffer status.

Return tuple (*acquired*, *unread*, *skipped*, *size*), where *acquired* is the total number of acquired frames, *unread* is the number of acquired but not read frames, *skipped* is the number of skipped (not read and then written over) frames, and *buffer\_size* is the total buffer size (in frames).

**get\_acquired\_frame\_status** ()

**set\_frameskip\_behavior** (*behavior*)

Choose the camera behavior if frame skip event is encountered when waiting for a new frame, reading frames, getting buffer status, etc.

Can be "error" (raise `uc480FrameTransferError`), "ignore" (continue acquisition, ignore the gap), or "skip" (mark some number of frames as skipped, but keep the frame counters consistent).



**get\_supported\_subsampling\_modes()**

Get all supported subsampling modes.

Return tuple (horizontal, vertical) of lists with all possible supported subsampling factors.

**get\_subsampling()**

Get current subsampling

**set\_subsampling(hsub=1, vsub=1)**

Set subsampling.

If values are not supported, get the closest value below the requested. Automatically turns off binning.

**get\_supported\_binning\_modes()**

Get all supported binning modes.

Return tuple (horizontal, vertical) of lists with all possible supported binning factors.

**get\_binning()**

Get current binning

**set\_binning(hbin=1, vbin=1)**

Set binning.

If values are not supported, get the closest value below the requested. Automatically turns off subsampling.

**get\_detector\_size()**

Get camera detector size (in pixels) as a tuple (width, height)

**get\_roi()**

Get current ROI.

Return tuple (hstart, hend, vstart, vend, hbin, vbin).

**set\_roi(hstart=0, hend=None, vstart=0, vend=None, hbin=1, vbin=1)**

Setup camera ROI.

*hstart* and *hend* specify horizontal image extent, *vstart* and *vend* specify vertical image extent (start are inclusive, stop are exclusive, starting from 0), *hbin* and *vbin* specify binning. By default, all non-supplied parameters take extreme values.

**get\_roi\_limits(hbin=1, vbin=1)**

Get the minimal and maximal ROI parameters.

Return tuple (hlim, vlim), where each element is in turn a limit 5-tuple (min, max, pstep, sstep, maxbin) with, correspondingly, minimal and maximal size, position and size step, and the maximal binning. In some cameras, the step and the minimal size depend on the binning, which can be supplied.

**class NoParameterCaller(device, kind)**

Bases: `object`

Class to simplify calling functions without a parameter

**apply\_settings(settings)**

Apply the settings.

*settings* is the dict {name: value} of the device available settings. Non-applicable settings are ignored.

**get\_acquisition\_parameters()**

Get acquisition parameters.

Return dictionary {name: value}

**get\_data\_dimensions()**

Get readout data dimensions (in pixels) as a tuple (width, height); take indexing mode into account

**get\_device\_variable(key)**

Get the value of a settings, status, or full info parameter

**get\_exposure()**

Get current exposure

**get\_frame\_format()**

Get format for the returned images.

Can be "list" (list of 2D arrays), "array" (a single 3D array), or "chunks" (list of 3D "chunk" arrays; supported for some cameras and provides the best performance).

**get\_frame\_info\_fields()**

Get the names of frame info fields.

Applicable when frame info format (set by `set_frame_info_format()`) is "list" or "array".

**get\_frame\_info\_format()**

Get format of the frame info.

Can be "namedtuple" (potentially nested named tuples; convenient to get particular values), "list" (flat list of values, with field names are given by `get_frame_info_fields()`; convenient for building a table), "array" (same as "list", but with a numpy array, which is easier to use for "chunks" frame format), or "dict" (flat dictionary with the same fields as the "list" format; more resilient to future format changes)

**get\_frame\_info\_period()**

Get period of frame info acquisition.

Frame info might be skipped (set to None) except for frames which indices are divisible by *period*. Useful for certain cameras where acquiring frame info takes a lot of time and can reduce performance at higher frame rates. Note that this parameter can still be ignored (i.e., always set to 1) if the performance is not an issue for a given camera class.

**get\_frame\_period()**

Get frame period (time between two consecutive frames in the internal trigger mode)

**get\_full\_info(include=0)**

Get dict {name: value} containing full device information (including status and settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include*=-10 queries all available variables, which is equivalent to *include*="all".

**get\_full\_status(include=0)**

Get dict {name: value} containing the device status (including settings).

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include*=-10 queries all available variables, which is equivalent to *include*="all".

**get\_image\_indexing()**

Get indexing for the returned images.

Can be "rct" (first index row, second index column, rows counted from the top), "rcb" (same as "rc", rows counted from the bottom), "xyt" (first index column, second index row, rows counted from the top), or "xyb" (same as "xyt", rows counted from the bottom)

**get\_new\_images\_range()**

Get the range of the new images.

Return tuple (*first*, *last*) with images range (first inclusive). If no images are available, return None. If some images were in the buffer were overwritten, exclude them from the range.

**get\_settings(include=0)**

Get dict {*name*: *value*} containing all the device settings.

*include* specifies either a list of variables (only these variables are returned), a priority threshold (only values with the priority equal or higher are returned), or "all" (all available variables). Since the lowest priority is -10, setting *include*=-10 queries all available variables, which is equivalent to *include*="all".

**grab(nframes=1, frame\_timeout=5.0, missing\_frame='skip', return\_info=False, buff\_size=None)**

Snap *nframes* images (with preset image read mode parameters)

*buff\_size* determines buffer size (if None, use the default size). Timeout is specified for a single-frame acquisition, not for the whole acquisition time. *missing\_frame* determines what to do with frames which have been lost: can be "none" (replacing them with None), "zero" (replacing them with zero-filled frame), or "skip" (skipping them, while still keeping total returned frames number to *n*). If *return\_info*==True, return tuple (*frames*, *infos*), where *infos* is a list of frame info tuples (camera-dependent); if some frames are missing and *missing\_frame*!="skip", the corresponding frame info is None.

**is\_acquisition\_setup()**

Check if acquisition is set up.

If the camera does not support separate acquisition setup, always return True.

**pausing\_acquisition(clear=None, stop=True, setup\_after=None, start\_after=True, combine\_nested=True)**

Context manager which temporarily pauses acquisition during execution of with block.

Useful for applying certain settings which can't be changed during the acquisition. If *clear*==True, clear acquisition in addition to stopping (by default, use the class default specified as *\_clear\_pausing\_acquisition* attribute). If *stop*==True, stop the acquisition (if *clear*==True, stop regardless). If *setup\_after*==True, setup the acquisition after pause if necessary (None means setup only if clearing was required). If *start\_after*==True, start the acquisition after pause if necessary (None means start only if stopping was required). If *combine\_nested*==True, then any nested *pausing\_acquisition* calls will stop/clear acquisition as necessary, but won't setup/start it again until this *pausing\_acquisition* call is complete.

Yields tuple (*acq\_in\_progress*, *acq\_params*), which indicates whether acquisition is currently in progress, and what are the current acquisition parameters.

**read\_newest\_image(peek=False, return\_info=False)**

Read the newest un-read image.

If no un-read frames are available, return None. If *peek*==True, return the image but not mark it as read. If *return\_info*==True, return tuple (*frame*, *info*), where *info* is an info tuples (camera-dependent, see [read\\_multiple\\_images\(\)](#)).

**read\_oldest\_image(peek=False, return\_info=False)**

Read the oldest un-read image.

If no un-read frames are available, return None. If *peek*==True, return the image but not mark it as read. If *return\_info*==True, return tuple (*frame*, *info*), where *info* is an info tuples (camera-dependent, see [read\\_multiple\\_images\(\)](#)).

**set\_device\_variable(key, value)**

Set the value of a settings parameter

**set\_frame\_format** (*fmt*)

Set format for the returned images.

Can be "list" (list of 2D arrays), "array" (a single 3D array), "chunks" (list of 3D "chunk" arrays; supported for some cameras and provides the best performance), or "try\_chunks" (same as "chunks", but if chunks are not supported, set to "list" instead). If format is "chunks" and chunks are not supported by the camera, it results in one frame per chunk. Note that if the format is set to "array" or "chunks", the frame info format is also automatically set to "array". If the format is set to "chunks", then the image info is also returned in chunks form (list of 2D info arrays with the same length as the corresponding frame chunks).

**set\_frame\_info\_format** (*fmt*, *include\_fields=None*)

Set format of the frame info.

Can be "namedtuple" (potentially nested named tuples; convenient to get particular values), "list" (flat list of values, with field names are given by `get_frame_info_fields()`; convenient for building a table), "array" (same as "list", but with a numpy array, which is easier to use for "chunks" frame format), or "dict" (flat dictionary with the same fields as the "list" format; more resilient to future format changes) If *include\_fields* is not None, it specifies the fields included for non-"tuple" formats; note that order or *include\_fields* is ignored, and the resulting fields are always ordered same as in the original.

**set\_frame\_info\_period** (*period=1*)

Set period of frame info acquisition.

Frame info might be skipped (set to None) except for frames which indices are divisible by *period*. Useful for certain cameras where acquiring frame info takes a lot of time and can reduce performance at higher frame rates. Note that this parameter can still be ignored (i.e., always set to 1) if the performance is not an issue for a given camera class.

**set\_image\_indexing** (*indexing*)

Set up indexing for the returned images.

Can be "rct" (first index row, second index column, rows counted from the top), "rcb" (same as "rc", rows counted from the bottom), "xyt" (first index column, second index row, rows counted from the top), or "xyb" (same as "xyt", rows counted from the bottom)

**snap** (*timeout=5.0*, *return\_info=False*)

Snap a single frame

**wait\_for\_frame** (*since='lastread'*, *nframes=1*, *timeout=20.0*, *error\_on\_stopped=False*)

Wait for one or several new camera frames.

*since* specifies the reference point for waiting to acquire *nframes* frames; can be "lastread" (from the last read frame), "lastwait" (wait for the last successful `wait_for_frame()` call), "now" (from the start of the current call), or "start" (from the acquisition start, i.e., wait until *nframes* frames have been acquired). *timeout* can be either a number, None (infinite timeout), or a tuple (*timeout*, *frame\_timeout*), in which case the call times out if the total time exceeds *timeout*, or a single frame wait exceeds *frame\_timeout*. If the call times out, raise `TimeoutError`. If *error\_on\_stopped==True* and the acquisition is not running, raise `Error`; otherwise, simply return False without waiting.

**read\_multiple\_images** (*rng=None*, *peek=False*, *missing\_frame='skip'*, *return\_info=False*, *return\_rng=False*)

Read multiple images specified by *rng* (by default, all un-read images).

If *rng* is specified, it is a tuple (*first*, *last*) with images range (first inclusive). If no new frames are available, return an empty list; if no acquisition is running, return None. If *peek==True*, return images but not mark them as read. *missing\_frame* determines what to do with frames which are out of range (missing or lost): can be "none" (replacing them with None), "zero" (replacing them with

zero-filled frame), or "skip" (skipping them). If `return_info==True`, return tuple (frames, infos), where `infos` is a list of `TFrameInfo` instances describing frame index, framestamp, global timestamp (real time), device timestamp (time from camera restart, in 0.1us steps), frame size, digital input state, and additional flags; if some frames are missing and `missing_frame!="skip"`, the corresponding frame info is `None`. if `return_rng==True`, return the range covered resulting frames; if `missing_frame=="skip"`, the range can be smaller than the supplied `rng` if some frames are skipped. Note that obtaining frame info might take about 2ms, so at high frame rates it will become a limiting factor.

## Module contents

### pylablib.devices.utils package

#### Submodules

#### pylablib.devices.utils.load\_lib module

`pylablib.devices.utils.load_lib.get_os_lib_folder()`

Get default Windows DLL folder (System32 or SysWOW64, depending on Python and Windows bitness)

`pylablib.devices.utils.load_lib.get_program_files_folder(subfolder="", arch=None)`

Get default Windows Program Files folder or a subfolder within it.

If `arch` is `None`, use the current Python architecture to determine the folder; otherwise, it specifies the architecture ("32bit" for Program Files (x86), "64bit" for Program Files)

`pylablib.devices.utils.load_lib.load_lib(name, locations=('global', ), call_conv='cdecl', locally=False, depends=None, depends_required=True, error_message=None, check_order='location', return_location=False)`

Load DLL.

#### Parameters

- **name** – name or path of the library (can also be a list or a tuple with several names, which are tried in that order).
- **locations** – list or tuple of locations to search for a library; the function tries locations in order and returns the first successfully loaded library a location is a string which can be a path to the containing folder, "parameter/\*" (the remaining part is a subpath inside "devices/dlls" library parameters; if this parameter is defined, it names folder or file for the dll), or "global" (load path as is; also searches in the standard OS specified locations determined by PATH variable, e.g., System32 folder).
- **depends** – if specified, it is a list of dependency libraries which need to be loaded first before the main DLL; they are assumed to be in the same location as the main file
- **depends\_required** – if `False`, ignore errors during dependency loads
- **locally** (*bool*) – if `True`, prepend path to the DLL containing folder to the environment PATH folders; this is usually required, if the loaded DLL imports other DLLs in the same folder
- **call\_conv** (*str*) – DLL call convention; can be either "cdecl" (corresponds to `ctypes.cdll`) or "stdcall" (corresponds to `ctypes.windll`)

- **error\_message** (*str*) – error message to add in addition to the default error message shown when the DLL is not found
- **check\_order** (*str*) – determines the order in which possible combinations of names and locations are looped over; can be "location" (loop over locations, and for each location loop over names), "name" (loop over names, and for each name loop over locations), or a list of tuples [ (*loc*, *name*) ] specifying order of checking (in the latter case, *name* and *location* arguments are ignored, except for generating error message).
- **return\_location** (*bool*) – if True, return a tuple (*dll*, *location*, *folder*) instead of a single *dll*.

```
class pylablib.devices.utils.load_lib.TLibraryOpenResult (init_result, open_result,  
                                                         opid)
```

Bases: `tuple`

**count** ()

Return number of occurrences of value.

**index** ()

Return first index of value.

Raises ValueError if the value is not present.

**init\_result**

**open\_result**

**opid**

```
class pylablib.devices.utils.load_lib.TLibraryCloseResult (close_result,  
                                                         uninit_result)
```

Bases: `tuple`

**close\_result**

**count** ()

Return number of occurrences of value.

**index** ()

Return first index of value.

Raises ValueError if the value is not present.

**uninit\_result**

```
class pylablib.devices.utils.load_lib.LibraryController (lib)
```

Bases: `object`

Simple wrapper to control libraries which require initialization when a new device is opened or shutdown when all devices are closed.

**Parameters** *lib* – controlled library

**preinit** ()

Pre-initialize the library, if it hasn't been done already

**open** ()

Mark device opening.

Return tuple (*init\_result*, *open\_result*, *opid*) with the results of the initialization and the opening, and the opening ID which should afterwards be used for closing. If library is already initialized, set *init\_result*=None

**close** (*opid*)

Mark device closing.

Return tuple (*close\_result*, *uninit\_result*) with the results of the closing and the shutdown.

If library does not need to be shut down yet, set *uninit\_result*=None

**temp\_open** ()

Context for temporarily opening a new device connection

**shutdown** ()

Close all opened connections and shutdown the library

## Module contents

## Module contents

## Submodules

### pylablib.widgets module

## Module contents

`pylablib.reload_all` (*from\_load\_path=True*, *keep\_parameters=True*)

Reload all loaded modules.

If *keep\_parameters*==True, keep the current library parameters (`pylablib.par`); otherwise, reset them to default.

`pylablib.unload_all` ()

Reload all loaded modules.

`pylablib.load_par` (*path*)

Load library parameters from a file

`pylablib.setbp` ()





## CHAPTER 2

---

### Indices and tables

---

- `genindex`
- `modindex`
- `search`



### p

`pylablib`, 875  
`pylablib.core`, 402  
`pylablib.core.dataproc`, 138  
`pylablib.core.dataproc.callable`, 104  
`pylablib.core.dataproc.feature`, 108  
`pylablib.core.dataproc.filters`, 110  
`pylablib.core.dataproc.fitting`, 114  
`pylablib.core.dataproc.fourier`, 117  
`pylablib.core.dataproc.iir_transform`, 120  
`pylablib.core.dataproc.image`, 121  
`pylablib.core.dataproc.interpolate`, 122  
`pylablib.core.dataproc.specfunc`, 124  
`pylablib.core.dataproc.table_wrap`, 125  
`pylablib.core.dataproc.transform`, 134  
`pylablib.core.dataproc.utils`, 135  
`pylablib.core.devio`, 168  
`pylablib.core.devio.backend_logger`, 141  
`pylablib.core.devio.base`, 142  
`pylablib.core.devio.comm_backend`, 142  
`pylablib.core.devio.data_format`, 161  
`pylablib.core.devio.interface`, 162  
`pylablib.core.devio.SCPI`, 138  
`pylablib.core.fileio`, 197  
`pylablib.core.fileio.datafile`, 168  
`pylablib.core.fileio.dict_entry`, 169  
`pylablib.core.fileio.loadfile`, 176  
`pylablib.core.fileio.loadfile_utils`, 182  
`pylablib.core.fileio.location`, 183  
`pylablib.core.fileio.parse_csv`, 188  
`pylablib.core.fileio.savefile`, 190  
`pylablib.core.fileio.table_stream`, 196  
`pylablib.core.gui`, 281  
`pylablib.core.gui.formatter`, 260  
`pylablib.core.gui.limiter`, 262  
`pylablib.core.gui.utils`, 263  
`pylablib.core.gui.value_handling`, 265  
`pylablib.core.gui.widgets`, 260  
`pylablib.core.gui.widgets.button`, 197  
`pylablib.core.gui.widgets.combo_box`, 198  
`pylablib.core.gui.widgets.container`, 199  
`pylablib.core.gui.widgets.edit`, 233  
`pylablib.core.gui.widgets.label`, 236  
`pylablib.core.gui.widgets.layout_manager`, 238  
`pylablib.core.gui.widgets.param_table`, 242  
`pylablib.core.thread`, 321  
`pylablib.core.thread.callsync`, 281  
`pylablib.core.thread.controller`, 292  
`pylablib.core.thread.multicast_pool`, 315  
`pylablib.core.thread.notifier`, 316  
`pylablib.core.thread.profile`, 317  
`pylablib.core.thread.synchronizing`, 318  
`pylablib.core.thread.threadprop`, 319  
`pylablib.core.thread.utils`, 321  
`pylablib.core.utils`, 402  
`pylablib.core.utils.array_utils`, 321  
`pylablib.core.utils.ctypes_wrap`, 321  
`pylablib.core.utils.dictionary`, 325  
`pylablib.core.utils.files`, 361  
`pylablib.core.utils.funcargparse`, 368  
`pylablib.core.utils.functions`, 369  
`pylablib.core.utils.general`, 374  
`pylablib.core.utils.indexing`, 381  
`pylablib.core.utils.ipc`, 383  
`pylablib.core.utils.library_parameters`, 385  
`pylablib.core.utils.module`, 385  
`pylablib.core.utils.nbtools`, 387  
`pylablib.core.utils.net`, 387  
`pylablib.core.utils.numerical`, 391  
`pylablib.core.utils.observer_pool`, 392  
`pylablib.core.utils.py3`, 393  
`pylablib.core.utils.rpyc_utils`, 394  
`pylablib.core.utils.strdump`, 395  
`pylablib.core.utils.string`, 396  
`pylablib.core.utils.strpack`, 401

pylablib.core.utils.units, 401  
pylablib.devices, 875  
pylablib.devices.AlliedVision, 465  
pylablib.devices.AlliedVision.Bonito, 450  
pylablib.devices.Andor, 487  
pylablib.devices.Andor.AndorSDK2, 465  
pylablib.devices.Andor.AndorSDK3, 476  
pylablib.devices.Andor.atcore\_features, 486  
pylablib.devices.Andor.base, 486  
pylablib.devices.Arcus, 500  
pylablib.devices.Arcus.base, 487  
pylablib.devices.Arcus.performax, 488  
pylablib.devices.Arduino, 502  
pylablib.devices.Arduino.base, 500  
pylablib.devices.Attocube, 510  
pylablib.devices.Attocube.anc300, 502  
pylablib.devices.Attocube.anc350, 506  
pylablib.devices.Attocube.base, 510  
pylablib.devices.AWG, 450  
pylablib.devices.AWG.generic, 402  
pylablib.devices.AWG.specific, 408  
pylablib.devices.Basler, 520  
pylablib.devices.Basler.pylon, 510  
pylablib.devices.BitFlow, 533  
pylablib.devices.BitFlow.BitFlow, 520  
pylablib.devices.Conrad, 535  
pylablib.devices.Conrad.base, 533  
pylablib.devices.Cryocon, 539  
pylablib.devices.Cryocon.base, 535  
pylablib.devices.Cryomagnetics, 548  
pylablib.devices.Cryomagnetics.base, 539  
pylablib.devices.DCAM, 556  
pylablib.devices.DCAM.DCAM, 548  
pylablib.devices.HighFinesse, 561  
pylablib.devices.HighFinesse.wlm, 557  
pylablib.devices.IMAQ, 576  
pylablib.devices.IMAQ.IMAQ, 561  
pylablib.devices.IMAQ.niimaq\_attrtypes, 576  
pylablib.devices.IMAQdx, 589  
pylablib.devices.IMAQdx.IMAQdx, 576  
pylablib.devices.interface, 864  
pylablib.devices.interface.camera, 830  
pylablib.devices.interface.stage, 862  
pylablib.devices.KJL, 592  
pylablib.devices.KJL.base, 589  
pylablib.devices.Lakeshore, 603  
pylablib.devices.Lakeshore.base, 592  
pylablib.devices.LaserQuantum, 606  
pylablib.devices.LaserQuantum.base, 603  
pylablib.devices.Leybold, 611  
pylablib.devices.Leybold.base, 606  
pylablib.devices.LighthousePhotonics, 613  
pylablib.devices.LighthousePhotonics.base, 611  
pylablib.devices.M2, 620  
pylablib.devices.M2.solstis, 614  
pylablib.devices.Newport, 631  
pylablib.devices.Newport.base, 627  
pylablib.devices.Newport.picomotor, 627  
pylablib.devices.NI, 627  
pylablib.devices.NI.daq, 620  
pylablib.devices.Ophir, 644  
pylablib.devices.Ophir.base, 638  
pylablib.devices.OZOptics, 638  
pylablib.devices.OZOptics.base, 631  
pylablib.devices.PCO, 653  
pylablib.devices.PCO.SC2, 644  
pylablib.devices.Pfeiffer, 659  
pylablib.devices.Pfeiffer.base, 653  
pylablib.devices.Photometrics, 669  
pylablib.devices.Photometrics.pvcam, 659  
pylablib.devices.PhotonFocus, 702  
pylablib.devices.PhotonFocus.PhotonFocus, 669  
pylablib.devices.PhysikInstrumente, 708  
pylablib.devices.PhysikInstrumente.base, 702  
pylablib.devices.PrincetonInstruments, 718  
pylablib.devices.PrincetonInstruments.picam, 708  
pylablib.devices.SiliconSoftware, 733  
pylablib.devices.SiliconSoftware.fgrab, 718  
pylablib.devices.Sirah, 745  
pylablib.devices.Sirah.Matisse, 733  
pylablib.devices.Sirah.tuner, 743  
pylablib.devices.SmarAct, 748  
pylablib.devices.SmarAct.scu3d, 745  
pylablib.devices.Tektronix, 769  
pylablib.devices.Tektronix.base, 748  
pylablib.devices.Thorlabs, 823  
pylablib.devices.Thorlabs.base, 776  
pylablib.devices.Thorlabs.elliptec, 776  
pylablib.devices.Thorlabs.kinesis, 781  
pylablib.devices.Thorlabs.misc, 809  
pylablib.devices.Thorlabs.serial, 809  
pylablib.devices.Thorlabs.TLCAmer, 769  
pylablib.devices.Toptica, 826  
pylablib.devices.Toptica.base, 823  
pylablib.devices.Toptica.ibeam, 823  
pylablib.devices.Trinamic, 830  
pylablib.devices.Trinamic.base, 826  
pylablib.devices.uc480, 873

`pylablib.devices.uc480.uc480`, [864](#)  
`pylablib.devices.utils`, [875](#)  
`pylablib.devices.utils.load_lib`, [873](#)  
`pylablib.widgets`, [875](#)



## A

- `acceleration` (`pylablib.devices.Thorlabs.kinesis.TJogParams` attribute), 784
- `acceleration` (`pylablib.devices.Thorlabs.kinesis.TPZMotorDriveParams` attribute), 786
- `acceleration` (`pylablib.devices.Thorlabs.kinesis.TPZMotorJogParams` attribute), 787
- `acceleration` (`pylablib.devices.Thorlabs.kinesis.TVelocityParams` attribute), 784
- `access` (`pylablib.devices.Basler.pylon.BaslerPylonAttribute` attribute), 512
- `AccessIterator` (class in `pylablib.core.utils.general`), 380
- `accum_cycle_time` (`pylablib.devices.Andor.AndorSDK2.TCycleTimings` attribute), 466
- `acknowledge()` (`pylablib.core.utils.general.Timer` method), 379
- `acqcleared()` (in module `pylablib.devices.interface.camera`), 835
- `acqstopped()` (in module `pylablib.devices.interface.camera`), 835
- `acquire()` (`pylablib.core.thread.synchronizing.QLockNotifier` method), 319
- `acquired` (`pylablib.devices.interface.camera.TFramesStatus` attribute), 830
- `acquired` (`pylablib.devices.uc480.uc480.TAcquiredFramesStatus` attribute), 865
- `acquisition_in_progress()` (`pylablib.devices.AlliedVision.Bonito.BonitoIMAQCamera` method), 456
- `acquisition_in_progress()` (`pylablib.devices.AlliedVision.Bonito.IBonitoCamera` method), 452
- `acquisition_in_progress()` (`pylablib.devices.Andor.AndorSDK2.AndorSDK2Camera` method), 467
- `acquisition_in_progress()` (`pylablib.devices.Andor.AndorSDK3.AndorSDK3Camera` method), 482
- `acquisition_in_progress()` (`pylablib.devices.Basler.pylon.BaslerPylonCamera` method), 516
- `acquisition_in_progress()` (`pylablib.devices.BitFlow.BitFlow.BitFlowCamera` method), 527
- `acquisition_in_progress()` (`pylablib.devices.BitFlow.BitFlow.BitFlowFrameGrabber` method), 523
- `acquisition_in_progress()` (`pylablib.devices.DCAM.DCAM.DCAMCamera` method), 553
- `acquisition_in_progress()` (`pylablib.devices.IMAQ.IMAQ.IMAQCamera` method), 569
- `acquisition_in_progress()` (`pylablib.devices.IMAQ.IMAQ.IMAQFrameGrabber` method), 565
- `acquisition_in_progress()` (`pylablib.devices.IMAQdx.IMAQdx.EthernetIMAQdxCamera` method), 584
- `acquisition_in_progress()` (`pylablib.devices.IMAQdx.IMAQdx.IMAQdxCamera` method), 580
- `acquisition_in_progress()` (`pylablib.devices.interface.camera.IAttributeCamera` method), 838
- `acquisition_in_progress()` (`pylablib.devices.interface.camera.IBinROICamera` method), 857
- `acquisition_in_progress()` (`pylablib.devices.interface.camera.ICamera` method), 832
- `acquisition_in_progress()` (`pylablib.devices.interface.camera.IExposureCamera` method), 847
- `acquisition_in_progress()` (`pylablib.devices.interface.camera.IGrabberAttributeCamera` method), 842

`acquisition_in_progress()` (py- `method`), 783  
`lablib.devices.interface.camera.IROICamera` `add_background_comm()` (py-  
`method`), 852 `lablib.devices.Thorlabs.kinesis.KinesisDevice`  
`acquisition_in_progress()` (py- `method`), 788  
`lablib.devices.PCO.SC2.PCOS2Camera` `add_background_comm()` (py-  
`method`), 648 `lablib.devices.Thorlabs.kinesis.KinesisMotor`  
`acquisition_in_progress()` (py- `method`), 799  
`lablib.devices.Photometrics.pvcam.PvcamCamera` `add_background_comm()` (py-  
`method`), 665 `lablib.devices.Thorlabs.kinesis.KinesisPiezoMotor`  
`acquisition_in_progress()` (py- `method`), 803  
`lablib.devices.PhotonFocus.PhotonFocus.IPhotonFocusCamera` `add_background_comm()` (py-  
`method`), 674 `lablib.devices.Thorlabs.kinesis.KinesisQuadDetector`  
`acquisition_in_progress()` (py- `method`), 807  
`lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusBitFlowCamera` `add_background_comm()` (py-  
`method`), 695 `lablib.devices.Thorlabs.kinesis.MFF` `method`),  
`acquisition_in_progress()` (py- 792  
`lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusIMAQCamera` `method`), 678  
`lablib.core.thread.controller.QTaskThread`  
`acquisition_in_progress()` (py- `method`), 303  
`lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusSISOCamera` (pylablib.core.gui.widgets.param\_table.ParamTable  
`method`), 687 `method`), 244  
`acquisition_in_progress()` (py- `add_button()` (pylablib.core.gui.widgets.param\_table.StatusTable  
`lablib.devices.PrincetonInstruments.picam.PicamCamera` `method`), 252  
`method`), 714 `add_callback()` (py-  
`acquisition_in_progress()` (py- `lablib.core.thread.callsync.QScheduledCall`  
`lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareCamera` `method`), 284  
`method`), 727 `add_check_box()` (py-  
`acquisition_in_progress()` (py- `lablib.core.gui.widgets.param_table.ParamTable`  
`lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareFrameGrabber` `method`), 245  
`method`), 723 `add_check_box()` (py-  
`acquisition_in_progress()` (py- `lablib.core.gui.widgets.param_table.StatusTable`  
`lablib.devices.Thorlabs.TLCamera.ThorlabsTLCamera` `method`), 252  
`method`), 772 `add_child()` (pylablib.core.gui.widgets.container.IQContainer  
`acquisition_in_progress()` (py- `method`), 201  
`lablib.devices.uc480.uc480.UC480Camera` `add_child()` (pylablib.core.gui.widgets.container.IQWidgetContainer  
`method`), 868 `method`), 205  
`activation_control` (py- `add_child()` (pylablib.core.gui.widgets.container.QContainer  
`lablib.devices.Pfeiffer.base.TTPG260GaugeControlSettings` `method`), 202  
`attribute`), 654 `add_child()` (pylablib.core.gui.widgets.container.QDialogContainer  
`add()` (pylablib.core.dataproc.filters.RunningDebounceFilter `method`), 217  
`method`), 114 `add_child()` (pylablib.core.gui.widgets.container.QFrameContainer  
`add()` (pylablib.core.dataproc.filters.RunningDecimationFilter `method`), 213  
`method`), 114 `add_child()` (pylablib.core.gui.widgets.container.QGroupBoxContainer  
`add_all_children()` (py- `method`), 221  
`lablib.core.gui.value_handling.GUIValues` `add_child()` (pylablib.core.gui.widgets.container.QScrollAreaContainer  
`method`), 279 `method`), 229  
`add_attribute()` (py- `add_child()` (pylablib.core.gui.widgets.container.QScrollAreaContainer  
`lablib.devices.Andor.AndorSDK3.AndorSDK3Camera` `method`), 225  
`method`), 479 `add_child()` (pylablib.core.gui.widgets.container.QTabContainer  
`add_background_comm()` (py- `method`), 231  
`lablib.devices.Thorlabs.elliptec.ElliptecMotor` `add_child()` (pylablib.core.gui.widgets.container.QWidgetContainer  
`method`), 778 `method`), 209  
`add_background_comm()` (py- `add_child()` (pylablib.core.gui.widgets.param\_table.ParamTable  
`lablib.devices.Thorlabs.kinesis.BasicKinesisDevice` `method`), 249



`add_child()` (`pylablib.core.gui.widgets.param_table.StatusTable` method), 253  
`add_child_values()` (`pylablib.core.gui.widgets.container.IQContainer` method), 201  
`add_child_values()` (`pylablib.core.gui.widgets.container.IQWidgetContainer` method), 206  
`add_child_values()` (`pylablib.core.gui.widgets.container.QContainer` method), 203  
`add_child_values()` (`pylablib.core.gui.widgets.container.QDialogContainer` method), 217  
`add_child_values()` (`pylablib.core.gui.widgets.container.QFrameContainer` method), 213  
`add_child_values()` (`pylablib.core.gui.widgets.container.QGroupBoxContainer` method), 221  
`add_child_values()` (`pylablib.core.gui.widgets.container.QScrollAreaContainer` method), 229  
`add_child_values()` (`pylablib.core.gui.widgets.container.QScrollAreaContainer.QContainedWidget` method), 225  
`add_child_values()` (`pylablib.core.gui.widgets.container.QTabContainer` method), 231  
`add_child_values()` (`pylablib.core.gui.widgets.container.QWidgetContainer` method), 209  
`add_child_values()` (`pylablib.core.gui.widgets.param_table.ParamTable` method), 249  
`add_child_values()` (`pylablib.core.gui.widgets.param_table.StatusTable` method), 253  
`add_chunk()` (`pylablib.core.fileio.parse_csv.ChunksAccumulator` method), 189  
`add_class()` (`pylablib.core.utils.strdump.StrDumper` method), 395  
`add_clock_period_input()` (`pylablib.devices.NI.daq.NIDAQ` method), 623  
`add_columns()` (`pylablib.core.fileio.parse_csv.ChunksAccumulator` method), 189  
`add_combo_box()` (`pylablib.core.gui.widgets.param_table.ParamTable` method), 247  
`add_combo_box()` (`pylablib.core.gui.widgets.param_table.StatusTable` method), 253  
`add_command()` (`pylablib.core.thread.controller.QTaskThread` method), 305  
`add_conversion_class()` (in module `pylablib.core.utils.string`), 398  
`add_counter_input()` (`pylablib.devices.NI.daq.NIDAQ` method), 622  
`add_custom_widget()` (`pylablib.core.gui.widgets.param_table.ParamTable` method), 244  
`add_custom_widget()` (`pylablib.core.gui.widgets.param_table.StatusTable` method), 253  
`add_decoration_label()` (`pylablib.core.gui.widgets.container.IQWidgetContainer` method), 206  
`add_decoration_label()` (`pylablib.core.gui.widgets.container.QDialogContainer` method), 217  
`add_decoration_label()` (`pylablib.core.gui.widgets.container.QFrameContainer` method), 213  
`add_decoration_label()` (`pylablib.core.gui.widgets.container.QGroupBoxContainer` method), 221  
`add_decoration_label()` (`pylablib.core.gui.widgets.container.QScrollAreaContainer.QContainedWidget` method), 225  
`add_decoration_label()` (`pylablib.core.gui.widgets.container.QTabContainer` method), 231  
`add_decoration_label()` (`pylablib.core.gui.widgets.container.QWidgetContainer` method), 209  
`add_decoration_label()` (`pylablib.core.gui.widgets.layout_manager.IQLayoutManagedWidget` method), 240  
`add_decoration_label()` (`pylablib.core.gui.widgets.layout_manager.QLayoutManagedWidget` method), 240  
`add_decoration_label()` (`pylablib.core.gui.widgets.param_table.ParamTable` method), 249  
`add_decoration_label()` (`pylablib.core.gui.widgets.param_table.StatusTable` method), 254  
`add_dict_entry_builder()` (in module `pylablib.core.fileio.dict_entry`), 169  
`add_dict_entry_class()` (in module `pylablib.core.fileio.dict_entry`), 169  
`add_dict_entry_parser()` (in module `pylablib.core.fileio.dict_entry`), 169  
`add_digital_input()` (`pylablib.devices.NI.daq.NIDAQ` method), 623  
`add_digital_output()` (`pylablib.devices.NI.daq.NIDAQ` method), 624  
`add_direct_call_command()` (`pylablib.core.thread.controller.QTaskThread` method), 305

<code>method)</code> , 305	<code>add_group_box()</code> (py-lablib.core.gui.widgets.container.QGroupBoxContainer method), 221
<code>add_dropdown_button()</code> (py-lablib.core.gui.widgets.param_table.ParamTable method), 245	<code>add_group_box()</code> (py-lablib.core.gui.widgets.container.QScrollAreaContainer.QContainer method), 225
<code>add_dropdown_button()</code> (py-lablib.core.gui.widgets.param_table.StatusTable method), 254	<code>add_group_box()</code> (py-lablib.core.gui.widgets.container.QWidgetContainer method), 209
<code>add_entry()</code> (pylablib.core.utils.dictionary.Dictionary method), 327	<code>add_group_box()</code> (py-lablib.core.gui.widgets.param_table.ParamTable method), 243
<code>add_entry()</code> (pylablib.core.utils.dictionary.DictionaryPainter method), 336	<code>add_group_box()</code> (py-lablib.core.gui.widgets.param_table.StatusTable method), 255
<code>add_entry()</code> (pylablib.core.utils.dictionary.FilterTree method), 352	<code>add_handler()</code> (py-lablib.core.gui.value_handling.GUIValues method), 278
<code>add_entry()</code> (pylablib.core.utils.dictionary.PrefixTree method), 344	<code>add_indicator_handler()</code> (py-lablib.core.gui.value_handling.GUIValues method), 279
<code>add_enum_label()</code> (py-lablib.core.gui.widgets.param_table.ParamTable method), 246	<code>add_job()</code> (pylablib.core.thread.controller.QTaskThread method), 303
<code>add_enum_label()</code> (py-lablib.core.gui.widgets.param_table.StatusTable method), 254	<code>add_label_indicator()</code> (py-lablib.core.gui.value_handling.GUIValues method), 280
<code>add_exception_hook()</code> (in module py-lablib.core.thread.controller), 292	<code>add_namedtuple_class()</code> (in module py-lablib.core.utils.string), 399
<code>add_file_format()</code> (py-lablib.core.fileio.dict_entry.ExternalNumpyDictionaryEntry static method), 175	<code>add_nested()</code> (pylablib.core.gui.value_handling.GUIValues method), 278
<code>add_file_format()</code> (py-lablib.core.fileio.dict_entry.IExternalFileDictionaryEntry static method), 174	<code>add_num_edit()</code> (py-lablib.core.gui.widgets.param_table.ParamTable method), 247
<code>add_frame()</code> (pylablib.core.gui.widgets.container.IQWidgetContainer method), 205	<code>add_num_edit()</code> (py-lablib.core.gui.widgets.param_table.StatusTable method), 255
<code>add_frame()</code> (pylablib.core.gui.widgets.container.QDialogContainer method), 217	<code>add_num_label()</code> (py-lablib.core.gui.widgets.param_table.ParamTable method), 246
<code>add_frame()</code> (pylablib.core.gui.widgets.container.QFrameContainer method), 213	<code>add_observer()</code> (py-lablib.core.utils.observer_pool.ObserverPool method), 393
<code>add_frame()</code> (pylablib.core.gui.widgets.container.QGroupBoxContainer method), 221	<code>add_padding()</code> (py-lablib.core.gui.widgets.container.IQWidgetContainer method), 206
<code>add_frame()</code> (pylablib.core.gui.widgets.container.QScrollAreaContainer.QContainer method), 225	<code>add_padding()</code> (py-lablib.core.gui.widgets.container.QDialogContainer method), 217
<code>add_frame()</code> (pylablib.core.gui.widgets.container.QWidgetContainer method), 209	<code>add_padding()</code> (py-lablib.core.gui.widgets.container.QFrameContainer method), 213
<code>add_frame()</code> (pylablib.core.gui.widgets.param_table.ParamTable method), 243	
<code>add_frame()</code> (pylablib.core.gui.widgets.param_table.StatusTable method), 255	
<code>add_group_box()</code> (py-lablib.core.gui.widgets.container.IQWidgetContainer method), 205	
<code>add_group_box()</code> (py-lablib.core.gui.widgets.container.QDialogContainer method), 217	
<code>add_group_box()</code> (py-lablib.core.gui.widgets.container.QFrameContainer method), 213	

`add_padding()` (py-lablib.core.gui.widgets.container.QGroupBoxContainer method), 221  
`add_padding()` (py-lablib.core.gui.widgets.container.QScrollAreaContainer.QContainerWidget method), 225  
`add_padding()` (py-lablib.core.gui.widgets.container.QTabContainer method), 225  
`add_padding()` (py-lablib.core.gui.widgets.container.QWidgetContainer method), 210  
`add_padding()` (py-lablib.core.gui.widgets.layout\_manager.IQLayoutManagedWidget method), 239  
`add_padding()` (py-lablib.core.gui.widgets.layout\_manager.QLayoutManagedWidget method), 240  
`add_padding()` (py-lablib.core.gui.widgets.param\_table.ParamTable method), 249  
`add_padding()` (py-lablib.core.gui.widgets.param\_table.StatusTable method), 256  
`add_path()` (pylablib.core.utils.general.StreamFileLogger method), 380  
`add_progress_bar()` (py-lablib.core.gui.widgets.param\_table.ParamTable method), 247  
`add_progress_bar()` (py-lablib.core.gui.widgets.param\_table.StatusTable method), 256  
`add_property_element()` (py-lablib.core.gui.value\_handling.GUIValues method), 279  
`add_property_element()` (py-lablib.core.gui.widgets.container.IQContainer method), 201  
`add_property_element()` (py-lablib.core.gui.widgets.container.IQWidgetContainer method), 206  
`add_property_element()` (py-lablib.core.gui.widgets.container.QContainer method), 203  
`add_property_element()` (py-lablib.core.gui.widgets.container.QDialogContainer method), 217  
`add_property_element()` (py-lablib.core.gui.widgets.container.QFrameContainer method), 214  
`add_property_element()` (py-lablib.core.gui.widgets.container.QGroupBoxContainer method), 222  
`add_property_element()` (py-lablib.core.gui.widgets.container.QScrollAreaContainer method), 226  
`add_property_element()` (py-lablib.core.gui.widgets.container.QWidgetContainer method), 210  
`add_property_element()` (py-lablib.core.gui.widgets.layout\_manager.IQLayoutManager method), 239  
`add_property_element()` (py-lablib.core.gui.widgets.layout\_manager.QLayoutManager method), 240  
`add_property_element()` (py-lablib.core.gui.widgets.param\_table.ParamTable method), 249  
`add_property_element()` (py-lablib.core.gui.widgets.param\_table.StatusTable method), 256  
`add_status_line()` (py-lablib.core.gui.widgets.param\_table.StatusTable method), 251  
`add_stop_notifier()` (py-lablib.core.thread.controller.QTaskThread method), 307  
`add_stop_notifier()` (py-lablib.core.thread.controller.QThreadController method), 307

method), 300

add\_sublayout() (pylablib.core.gui.widgets.container.IQWidgetContainer method), 206

add\_sublayout() (pylablib.core.gui.widgets.container.QDialogContainer method), 218

add\_sublayout() (pylablib.core.gui.widgets.container.QFrameContainer method), 214

add\_sublayout() (pylablib.core.gui.widgets.container.QGroupBoxContainer method), 222

add\_sublayout() (pylablib.core.gui.widgets.container.QScrollAreaContainer method), 226

add\_sublayout() (pylablib.core.gui.widgets.container.QWidgetContainer method), 210

add\_sublayout() (pylablib.core.gui.widgets.layout\_manager.IQLayoutManager method), 239

add\_sublayout() (pylablib.core.gui.widgets.layout\_manager.QLayoutManager method), 240

add\_sublayout() (pylablib.core.gui.widgets.param\_table.ParamTable method), 242

add\_sublayout() (pylablib.core.gui.widgets.param\_table.StatusTable method), 257

add\_tab() (pylablib.core.gui.widgets.container.QTabContainer method), 231

add\_text\_edit() (pylablib.core.gui.widgets.param\_table.ParamTable method), 246

add\_text\_edit() (pylablib.core.gui.widgets.param\_table.StatusTable method), 257

add\_text\_label() (pylablib.core.gui.widgets.param\_table.ParamTable method), 246

add\_text\_label() (pylablib.core.gui.widgets.param\_table.StatusTable method), 257

add\_thread\_method() (pylablib.core.thread.controller.QTaskThread method), 307

add\_thread\_method() (pylablib.core.thread.controller.QThreadController method), 298

add\_time() (pylablib.core.utils.general.Countdown method), 379

add\_timer() (pylablib.core.gui.widgets.container.IQContainer method), 200

add\_timer() (pylablib.core.gui.widgets.container.IQWidgetContainer method), 206

add\_timer() (pylablib.core.gui.widgets.container.QContainer method), 203

add\_timer() (pylablib.core.gui.widgets.container.QDialogContainer method), 218

add\_timer() (pylablib.core.gui.widgets.container.QFrameContainer method), 214

add\_timer() (pylablib.core.gui.widgets.container.QGroupBoxContainer method), 222

add\_timer() (pylablib.core.gui.widgets.container.QScrollAreaContainer method), 229

add\_timer() (pylablib.core.gui.widgets.container.QScrollAreaContainer method), 226

add\_timer() (pylablib.core.gui.widgets.container.QTabContainer method), 232

add\_timer() (pylablib.core.gui.widgets.container.QWidgetContainer method), 210

add\_timer() (pylablib.core.gui.widgets.param\_table.ParamTable method), 249

add\_timer() (pylablib.core.gui.widgets.param\_table.StatusTable method), 257

add\_timer\_event() (pylablib.core.gui.widgets.container.IQContainer method), 201

add\_timer\_event() (pylablib.core.gui.widgets.container.IQWidgetContainer method), 206

add\_timer\_event() (pylablib.core.gui.widgets.container.QContainer method), 203

add\_timer\_event() (pylablib.core.gui.widgets.container.QDialogContainer method), 218

add\_timer\_event() (pylablib.core.gui.widgets.container.QFrameContainer method), 214

add\_timer\_event() (pylablib.core.gui.widgets.container.QGroupBoxContainer method), 222

add\_timer\_event() (pylablib.core.gui.widgets.container.QScrollAreaContainer method), 229

add\_timer\_event() (pylablib.core.gui.widgets.container.QScrollAreaContainer.QContainer method), 226

add\_timer\_event() (pylablib.core.gui.widgets.container.QTabContainer method), 232

add\_timer\_event() (pylablib.core.gui.widgets.container.QWidgetContainer method), 210

add\_timer\_event() (py-

<code>lablib.core.gui.widgets.param_table.ParamTable</code> <code>method</code> ), 249	<code>lablib.core.gui.widgets.container.QContainer</code> <code>method</code> ), 203
<code>add_timer_event()</code> (py- <code>lablib.core.gui.widgets.param_table.StatusTable</code> <code>method</code> ), 257	<code>add_virtual_element()</code> (py- <code>lablib.core.gui.widgets.container.QDialogContainer</code> <code>method</code> ), 218
<code>add_to_layout()</code> (py- <code>lablib.core.gui.widgets.container.IQWidgetContainer</code> <code>method</code> ), 207	<code>add_virtual_element()</code> (py- <code>lablib.core.gui.widgets.container.QFrameContainer</code> <code>method</code> ), 214
<code>add_to_layout()</code> (py- <code>lablib.core.gui.widgets.container.QDialogContainer</code> <code>method</code> ), 218	<code>add_virtual_element()</code> (py- <code>lablib.core.gui.widgets.container.QGroupBoxContainer</code> <code>method</code> ), 222
<code>add_to_layout()</code> (py- <code>lablib.core.gui.widgets.container.QFrameContainer</code> <code>method</code> ), 214	<code>add_virtual_element()</code> (py- <code>lablib.core.gui.widgets.container.QScrollAreaContainer</code> <code>method</code> ), 229
<code>add_to_layout()</code> (py- <code>lablib.core.gui.widgets.container.QGroupBoxContainer</code> <code>method</code> ), 222	<code>add_virtual_element()</code> (py- <code>lablib.core.gui.widgets.container.QScrollAreaContainer.QContain</code> <code>method</code> ), 226
<code>add_to_layout()</code> (py- <code>lablib.core.gui.widgets.container.QScrollAreaContainer.QContainerWidget</code> <code>method</code> ), 226	<code>add_virtual_element()</code> (py- <code>lablib.core.gui.widgets.container.QTabContainer</code> <code>method</code> ), 232
<code>add_to_layout()</code> (py- <code>lablib.core.gui.widgets.container.QWidgetContainer</code> <code>method</code> ), 210	<code>add_virtual_element()</code> (py- <code>lablib.core.gui.widgets.container.QWidgetContainer</code> <code>method</code> ), 210
<code>add_to_layout()</code> (py- <code>lablib.core.gui.widgets.layout_manager.IQLayoutManagedWidget</code> <code>method</code> ), 238	<code>add_virtual_element()</code> (py- <code>lablib.core.gui.widgets.param_table.ParamTable</code> <code>method</code> ), 244
<code>add_to_layout()</code> (py- <code>lablib.core.gui.widgets.layout_manager.QLayoutManagedWidget</code> <code>method</code> ), 240	<code>add_virtual_element()</code> (py- <code>lablib.core.gui.widgets.param_table.StatusTable</code> <code>method</code> ), 258
<code>add_to_layout()</code> (py- <code>lablib.core.gui.widgets.param_table.ParamTable</code> <code>method</code> ), 250	<code>add_voltage_input()</code> (py- <code>lablib.devices.NI.daq.NIDAQ</code> <code>method</code> ), 622
<code>add_to_layout()</code> (py- <code>lablib.core.gui.widgets.param_table.StatusTable</code> <code>method</code> ), 257	<code>add_voltage_output()</code> (py- <code>lablib.devices.NI.daq.NIDAQ</code> <code>method</code> ), 624
<code>add_toggle_button()</code> (py- <code>lablib.core.gui.widgets.param_table.ParamTable</code> <code>method</code> ), 245	<code>add_widget()</code> (pylablib.core.gui.value_handling.GUIValues <code>method</code> ), 278
<code>add_toggle_button()</code> (py- <code>lablib.core.gui.widgets.param_table.StatusTable</code> <code>method</code> ), 258	<code>add_widget_indicator()</code> (py- <code>lablib.core.gui.value_handling.GUIValues</code> <code>method</code> ), 279
<code>add_variable()</code> (py- <code>lablib.core.utils.ipc.SharedMemIPCTable</code> <code>method</code> ), 384	<code>added</code> (pylablib.core.utils.dictionary.DictionaryDiff <code>at-</code> <code>tribute</code> ), 334
<code>add_virtual_element()</code> (py- <code>lablib.core.gui.value_handling.GUIValues</code> <code>method</code> ), 278	<code>addr</code> (pylablib.devices.Conrad.base.RelayBoard.TMessage <code>attribute</code> ), 533
<code>add_virtual_element()</code> (py- <code>lablib.core.gui.widgets.container.IQContainer</code> <code>method</code> ), 201	<code>addr</code> (pylablib.devices.Thorlabs.elliptec.ElliptecMotor.CommData <code>attribute</code> ), 778
<code>add_virtual_element()</code> (py- <code>lablib.core.gui.widgets.container.IQWidgetContainer</code> <code>method</code> ), 207	<code>addr</code> (pylablib.devices.Trinamic.base.TMCM1110.ReplyData <code>attribute</code> ), 827
<code>add_virtual_element()</code> (py- <code>lablib.core.gui.widgets.container.IQWidgetContainer</code> <code>method</code> ), 207	<code>address</code> (pylablib.devices.Attocube.anc350.ANC350.Reply <code>attribute</code> ), 507
<code>add_virtual_element()</code> (py- <code>lablib.core.gui.widgets.container.IQWidgetContainer</code> <code>method</code> ), 207	<code>address</code> (pylablib.devices.Attocube.anc350.ANC350.Telegram <code>attribute</code> ), 506
<code>add_virtual_element()</code> (py- <code>lablib.core.gui.widgets.container.IQWidgetContainer</code> <code>method</code> ), 207	<code>advance_read_frames()</code> (py- <code>lablib.devices.interface.camera.FrameCounter</code> <code>method</code> ), 836
<code>add_virtual_element()</code> (py-	<code>Agilent33220A</code> (class in py-



- lablib.devices.AWG.specific*), 414
- Agilent33220A.NoParameterCaller (class in *pylablib.devices.AWG.specific*), 415
- Agilent33500 (class in *pylablib.devices.AWG.specific*), 408
- Agilent33500.NoParameterCaller (class in *pylablib.devices.AWG.specific*), 409
- ai1() (in module *pylablib.core.utils.nbtools*), 387
- ai2() (in module *pylablib.core.utils.nbtools*), 387
- ai4() (in module *pylablib.core.utils.nbtools*), 387
- ai8() (in module *pylablib.core.utils.nbtools*), 387
- allocate() (*pylablib.devices.interface.camera.ChunkBufferManager* attribute), 720
- allocate\_buffers() (*pylablib.devices.Andor.AndorSDK3.AndorSDK3Camera.BufferManager* method), 481
- allowing\_toploop() (*pylablib.core.thread.controller.QTaskThread* method), 308
- allowing\_toploop() (*pylablib.core.thread.controller.QThreadController* method), 294
- ampl (*pylablib.devices.Sirah.Matisse.TPiezoetFeedforwardParameter* attribute), 734
- amplitude (*pylablib.devices.Sirah.Matisse.TPiezoetDriverParameter* attribute), 734
- ANC300 (class in *pylablib.devices.Attocube.anc300*), 502
- ANC300.NoParameterCaller (class in *pylablib.devices.Attocube.anc300*), 505
- ANC350 (class in *pylablib.devices.Attocube.anc350*), 506
- ANC350.NoParameterCaller (class in *pylablib.devices.Attocube.anc350*), 509
- ANC350.Reply (class in *pylablib.devices.Attocube.anc350*), 506
- ANC350.Telegram (class in *pylablib.devices.Attocube.anc350*), 506
- AndorError, 486
- AndorFrameTransferError, 487
- AndorNotSupportedError, 487
- AndorSDK2Camera (class in *pylablib.devices.Andor.AndorSDK2*), 466
- AndorSDK2Camera.NoParameterCaller (class in *pylablib.devices.Andor.AndorSDK2*), 473
- AndorSDK3Attribute (class in *pylablib.devices.Andor.AndorSDK3*), 477
- AndorSDK3Camera (class in *pylablib.devices.Andor.AndorSDK3*), 479
- AndorSDK3Camera.BufferManager (class in *pylablib.devices.Andor.AndorSDK3*), 481
- AndorSDK3Camera.NoParameterCaller (class in *pylablib.devices.Andor.AndorSDK3*), 483
- AndorTimeoutError, 486
- any\_item() (in module *pylablib.core.utils.general*), 375
- append() (*pylablib.core.dataproc.table\_wrap.Array1DWrapper* method), 127
- append() (*pylablib.core.dataproc.table\_wrap.Array2DWrapper.ColumnA* method), 130
- append() (*pylablib.core.dataproc.table\_wrap.Array2DWrapper.RowA* method), 130
- append() (*pylablib.core.dataproc.table\_wrap.DataFrame2DWrapper.ColumnA* method), 132
- append() (*pylablib.core.dataproc.table\_wrap.DataFrame2DWrapper.RowA* method), 132
- applet\_info (*pylablib.devices.SiliconSoftware.fgrab.TDeviceInfo* attribute), 720
- apply\_settings() (*pylablib.core.devio.comm\_backend.ICommBackendWrapper* method), 160
- apply\_settings() (*pylablib.core.devio.interface.IDevice* method), 163
- apply\_settings() (*pylablib.core.devio.SCPi.SCPIDevice* method), 141
- apply\_settings() (*pylablib.devices.AlliedVision.Bonito.BonitoIMAQCamera* method), 456
- apply\_settings() (*pylablib.devices.AlliedVision.Bonito.IBonitoCamera* method), 452
- apply\_settings() (*pylablib.devices.Andor.AndorSDK2.AndorSDK2Camera* method), 473
- apply\_settings() (*pylablib.devices.Andor.AndorSDK3.AndorSDK3Camera* method), 483
- apply\_settings() (*pylablib.devices.Arcus.performax.GenericPerformaxStage* method), 489
- apply\_settings() (*pylablib.devices.Arcus.performax.Performax2EXStage* method), 493
- apply\_settings() (*pylablib.devices.Arcus.performax.Performax4EXStage* method), 492
- apply\_settings() (*pylablib.devices.Arcus.performax.PerformaxDMXJSStage* method), 498
- apply\_settings() (*pylablib.devices.Arduino.base.IArduinoDevice* method), 501
- apply\_settings() (*pylablib.devices.Attocube.anc300.ANC300* method), 505
- apply\_settings() (*pylablib.devices.Attocube.anc350.ANC350* method), 509

<code>apply_settings()</code> <i>lablib.devices.AWG.generic.GenericAWG</i> <i>method</i> ), 406	(py-	<code>apply_settings()</code> <i>lablib.devices.IMAQ.IMAQ.IMAQFrameGrabber</i> <i>method</i> ), 565	(py-
<code>apply_settings()</code> <i>lablib.devices.AWG.specific.Agilent33220A</i> <i>method</i> ), 415	(py-	<code>apply_settings()</code> <i>lablib.devices.IMAQdx.IMAQdx.EthernetIMAQdxCamera</i> <i>method</i> ), 584	(py-
<code>apply_settings()</code> <i>lablib.devices.AWG.specific.Agilent33500</i> <i>method</i> ), 409	(py-	<code>apply_settings()</code> <i>lablib.devices.IMAQdx.IMAQdx.IMAQdxCamera</i> <i>method</i> ), 580	(py-
<code>apply_settings()</code> <i>lablib.devices.AWG.specific.InstekAFG2000</i> <i>method</i> ), 427	(py-	<code>apply_settings()</code> <i>lablib.devices.interface.camera.IAttributeCamera</i> <i>method</i> ), 838	(py-
<code>apply_settings()</code> <i>lablib.devices.AWG.specific.InstekAFG2225</i> <i>method</i> ), 421	(py-	<code>apply_settings()</code> <i>lablib.devices.interface.camera.IBinROICamera</i> <i>method</i> ), 857	(py-
<code>apply_settings()</code> <i>lablib.devices.AWG.specific.RigolDG1000</i> <i>method</i> ), 444	(py-	<code>apply_settings()</code> <i>lablib.devices.interface.camera.ICamera</i> <i>method</i> ), 834	(py-
<code>apply_settings()</code> <i>lablib.devices.AWG.specific.RSInstekAFG21000</i> <i>method</i> ), 433	(py-	<code>apply_settings()</code> <i>lablib.devices.interface.camera.IExposureCamera</i> <i>method</i> ), 847	(py-
<code>apply_settings()</code> <i>lablib.devices.AWG.specific.TektronixAFG1000</i> <i>method</i> ), 438	(py-	<code>apply_settings()</code> <i>lablib.devices.interface.camera.IGrabberAttributeCamera</i> <i>method</i> ), 843	(py-
<code>apply_settings()</code> <i>lablib.devices.Basler.pylon.BaslerPylonCamera</i> <i>method</i> ), 517	(py-	<code>apply_settings()</code> <i>lablib.devices.interface.camera.IROICamera</i> <i>method</i> ), 852	(py-
<code>apply_settings()</code> <i>lablib.devices.BitFlow.BitFlow.BitFlowCamera</i> <i>method</i> ), 527	(py-	<code>apply_settings()</code> <i>lablib.devices.interface.stage.IMultiaxisStage</i> <i>method</i> ), 863	(py-
<code>apply_settings()</code> <i>lablib.devices.BitFlow.BitFlow.BitFlowFrameGrabber</i> <i>method</i> ), 523	(py-	<code>apply_settings()</code> <i>lablib.devices.interface.stage.IStage</i> <i>method</i> ), 862	(py-
<code>apply_settings()</code> <i>lablib.devices.Conrad.base.RelayBoard</i> <i>method</i> ), 534	(py-	<code>apply_settings()</code> <i>lablib.devices.KJL.base.KJL300</i> <i>method</i> ), 591	(py-
<code>apply_settings()</code> <i>lablib.devices.Cryocon.base.Cryocon1x</i> <i>method</i> ), 536	(py-	<code>apply_settings()</code> <i>lablib.devices.Lakeshore.base.Lakeshore218</i> <i>method</i> ), 595	(py-
<code>apply_settings()</code> <i>lablib.devices.Cryomagnetics.base.LM500</i> <i>method</i> ), 541	(py-	<code>apply_settings()</code> <i>lablib.devices.Lakeshore.base.Lakeshore370</i> <i>method</i> ), 600	(py-
<code>apply_settings()</code> <i>lablib.devices.Cryomagnetics.base.LM510</i> <i>method</i> ), 544	(py-	<code>apply_settings()</code> <i>lablib.devices.LaserQuantum.base.Finesse</i> <i>method</i> ), 605	(py-
<code>apply_settings()</code> <i>lablib.devices.DCAM.DCAM.DCAMCamera</i> <i>method</i> ), 553	(py-	<code>apply_settings()</code> <i>lablib.devices.Leybold.base.GenericITR</i> <i>method</i> ), 608	(py-
<code>apply_settings()</code> <i>lablib.devices.HighFinesse.wlm.WLM</i> <i>method</i> ), 560	(py-	<code>apply_settings()</code> <i>lablib.devices.Leybold.base.ITR90</i> <i>method</i> ), 609	(py-
<code>apply_settings()</code> <i>lablib.devices.IMAQ.IMAQ.IMAQCamera</i> <i>method</i> ), 569	(py-	<code>apply_settings()</code> <i>lablib.devices.LighthousePhotonics.base.SproutG</i> <i>method</i> ), 612	(py-

<code>apply_settings()</code> <i>lablib.devices.M2.solstis.Solstis</i> 619	(py- method),	<i>lablib.devices.PhysikInstrumente.base.PIE516</i> method), 706
<code>apply_settings()</code> <i>lablib.devices.Newport.picomotor.Picomotor8742</i> method), 630	(py- method),	<i>lablib.devices.PrincetonInstruments.picam.PicamCamera</i> method), 714
<code>apply_settings()</code> ( <i>pylablib.devices.NI.daq.NIDAQ</i> method), 626	(py- method),	<i>lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareCamera</i> method), 727
<code>apply_settings()</code> <i>lablib.devices.Ophir.base.OphirDevice</i> method), 639	(py- method),	<i>lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareFrameGrabb</i> method), 723
<code>apply_settings()</code> <i>lablib.devices.Ophir.base.VegaPowerMeter</i> method), 643	(py- method),	<i>lablib.devices.Sirah.Matisse.SirahMatisse</i> method), 740
<code>apply_settings()</code> <i>lablib.devices.OZOptics.base.DD100</i> method), 635	(py- method),	<i>lablib.devices.SmarAct.scu3d.SCU3D</i> method), 747
<code>apply_settings()</code> <i>lablib.devices.OZOptics.base.EPC04</i> method), 637	(py- method),	<i>lablib.devices.Tektronix.base.DPO2000</i> method), 762
<code>apply_settings()</code> <i>lablib.devices.OZOptics.base.OZOpticsDevice</i> method), 632	(py- method),	<i>lablib.devices.Tektronix.base.ITektronixScope</i> method), 752
<code>apply_settings()</code> <i>lablib.devices.OZOptics.base.TF100</i> method), 633	(py- method),	<i>lablib.devices.Tektronix.base.TDS2000</i> method), 756
<code>apply_settings()</code> <i>lablib.devices.PCO.SC2.PCOSC2Camera</i> method), 649	(py- method),	<i>lablib.devices.Thorlabs.elliptec.ElliptecMotor</i> method), 780
<code>apply_settings()</code> <i>lablib.devices.Pfeiffer.base.DPG202</i> method), 658	(py- method),	<i>lablib.devices.Thorlabs.kinesis.BasicKinesisDevice</i> method), 783
<code>apply_settings()</code> <i>lablib.devices.Pfeiffer.base.TPG260</i> method), 656	(py- method),	<i>lablib.devices.Thorlabs.kinesis.KinesisDevice</i> method), 788
<code>apply_settings()</code> <i>lablib.devices.Photometrics.pvcam.PvcamCamera</i> method), 665	(py- method),	<i>lablib.devices.Thorlabs.kinesis.KinesisMotor</i> method), 799
<code>apply_settings()</code> <i>lablib.devices.PhotonFocus.PhotonFocus.IPhotonFocusCamera</i> method), 674	(py- method),	<i>lablib.devices.Thorlabs.kinesis.KinesisPiezoMotor</i> method), 803
<code>apply_settings()</code> <i>lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusBitFluxCamera</i> method), 695	(py- method),	<i>lablib.devices.Thorlabs.kinesis.KinesisQuadDetector</i> method), 808
<code>apply_settings()</code> <i>lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusIMACamera</i> method), 678	(py- method),	<i>lablib.devices.Thorlabs.kinesis.MFF</i> method), 792
<code>apply_settings()</code> <i>lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusSiSoCamera</i> method), 687	(py- method),	<i>lablib.devices.Thorlabs.serial.FW</i> method), 813
<code>apply_settings()</code> <i>lablib.devices.PhysikInstrumente.base.GenericPIController</i> method), 704	(py- method),	<i>lablib.devices.Thorlabs.serial.FWv1</i> method), 817
<code>apply_settings()</code>	(py- method),	<i>lablib.devices.Thorlabs.serial.FWv1</i> method), 817



[lablib.devices.Thorlabs.serial.MDT69xA method](#)), 820  
[apply\\_settings\(\)](#) ([pylablib.devices.Thorlabs.serial.ThorlabsSerialInterface method](#)), 810  
[apply\\_settings\(\)](#) ([pylablib.devices.Thorlabs.TLCamera.ThorlabsTLCamera method](#)), 773  
[apply\\_settings\(\)](#) ([pylablib.devices.Topica.ibeam.TopicalBeam method](#)), 825  
[apply\\_settings\(\)](#) ([pylablib.devices.Trinamic.base.TMCM1110 method](#)), 829  
[apply\\_settings\(\)](#) ([pylablib.devices.uc480.uc480.UC480Camera method](#)), 869  
[apply\\_window\(\)](#) (in module [pylablib.core.dataproc.fourier](#)), 118  
[ArcusBackendError](#), 487  
[ArcusError](#), 487  
[ArduinoBackendError](#), 500  
[ArduinoError](#), 500  
[area\(\)](#) ([pylablib.core.dataproc.image.ROI method](#)), 121  
[arg\\_value\(\)](#) ([pylablib.core.utils.functions.FunctionSignature method](#)), 370  
[args](#) ([pylablib.core.devio.base.DeviceError attribute](#)), 142  
[args](#) ([pylablib.core.devio.comm\\_backend.DeviceBackendError attribute](#)), 142  
[args](#) ([pylablib.core.devio.comm\\_backend.DeviceFT232Error attribute](#)), 150  
[args](#) ([pylablib.core.devio.comm\\_backend.DeviceNetworkError attribute](#)), 152  
[args](#) ([pylablib.core.devio.comm\\_backend.DeviceRecordedError attribute](#)), 157  
[args](#) ([pylablib.core.devio.comm\\_backend.DeviceSerialError attribute](#)), 147  
[args](#) ([pylablib.core.devio.comm\\_backend.DeviceUSBError attribute](#)), 155  
[args](#) ([pylablib.core.devio.comm\\_backend.DeviceVisaError attribute](#)), 144  
[args](#) ([pylablib.core.gui.limiter.LimitError attribute](#)), 262  
[args](#) ([pylablib.core.gui.value\\_handling.MissingGUIHandlerError attribute](#)), 278  
[args](#) ([pylablib.core.gui.value\\_handling.NoParameterError attribute](#)), 265  
[args](#) ([pylablib.core.thread.threadprop.DuplicateControllerThreadError attribute](#)), 319  
[args](#) ([pylablib.core.thread.threadprop.InterruptException attribute](#)), 320  
[args](#) ([pylablib.core.thread.threadprop.InterruptExceptionStop attribute](#)), 320  
[args](#) ([pylablib.core.thread.threadprop.NoControllerThreadError attribute](#)), 319  
[args](#) ([pylablib.core.thread.threadprop.NoMessageThreadError attribute](#)), 320  
[args](#) ([pylablib.core.thread.threadprop.SkippedCallError attribute](#)), 320  
[args](#) ([pylablib.core.thread.threadprop.ThreadError attribute](#)), 319  
[args](#) ([pylablib.core.thread.threadprop.TimeoutThreadError attribute](#)), 319  
[args](#) ([pylablib.core.utils.net.SocketError attribute](#)), 388  
[args](#) ([pylablib.core.utils.net.SocketTimeout attribute](#)), 388  
[args](#) ([pylablib.devices.AlliedVision.Bonito.BonitoError attribute](#)), 450  
[args](#) ([pylablib.devices.Andor.base.AndorError attribute](#)), 486  
[args](#) ([pylablib.devices.Andor.base.AndorFrameTransferError attribute](#)), 487  
[args](#) ([pylablib.devices.Andor.base.AndorNotSupportedError attribute](#)), 487  
[args](#) ([pylablib.devices.Andor.base.AndorTimeoutError attribute](#)), 487  
[args](#) ([pylablib.devices.Arcus.base.ArcusBackendError attribute](#)), 487  
[args](#) ([pylablib.devices.Arcus.base.ArcusError attribute](#)), 487  
[args](#) ([pylablib.devices.Arduino.base.ArduinoBackendError attribute](#)), 500  
[args](#) ([pylablib.devices.Arduino.base.ArduinoError attribute](#)), 500  
[args](#) ([pylablib.devices.Attocube.base.AttocubeBackendError attribute](#)), 510  
[args](#) ([pylablib.devices.Attocube.base.AttocubeError attribute](#)), 510  
[args](#) ([pylablib.devices.AWG.generic.GenericAWGBackendError attribute](#)), 402  
[args](#) ([pylablib.devices.AWG.generic.GenericAWGError attribute](#)), 402  
[args](#) ([pylablib.devices.BitFlow.BitFlow.BitFlowError attribute](#)), 520  
[args](#) ([pylablib.devices.BitFlow.BitFlow.BitFlowTimeoutError attribute](#)), 520  
[args](#) ([pylablib.devices.Conrad.base.ConradBackendError attribute](#)), 533  
[args](#) ([pylablib.devices.Conrad.base.ConradError attribute](#)), 533  
[args](#) ([pylablib.devices.Cryocon.base.CryoconBackendError attribute](#)), 535  
[args](#) ([pylablib.devices.Cryocon.base.CryoconError attribute](#)), 535  
[args](#) ([pylablib.devices.Cryomagnetics.base.CryomagneticsBackendError attribute](#)), 540  
[args](#) ([pylablib.devices.Cryomagnetics.base.CryomagneticsError attribute](#)), 540

- [attribute](#)), 539
- [args \(pylablib.devices.interface.camera.DefaultFrameTransferError attribute\)](#)), 748
- [attribute](#)), 830
- [args \(pylablib.devices.KJL.base.KJLBackendError attribute\)](#)), 590
- [args \(pylablib.devices.KJL.base.KJLError attribute\)](#)), 589
- [args \(pylablib.devices.Lakeshore.base.LakeshoreBackendError attribute\)](#)), 592
- [args \(pylablib.devices.Lakeshore.base.LakeshoreError attribute\)](#)), 592
- [args \(pylablib.devices.LaserQuantum.base.LaserQuantumBackendError attribute\)](#)), 823
- [attribute](#)), 603
- [args \(pylablib.devices.LaserQuantum.base.LaserQuantumError attribute\)](#)), 603
- [args \(pylablib.devices.Leybold.base.LeyboldBackendError attribute\)](#)), 606
- [args \(pylablib.devices.Leybold.base.LeyboldError attribute\)](#)), 606
- [args \(pylablib.devices.LighthousePhotonics.base.LighthousePhotonicsBackendError attribute\)](#)), 611
- [args \(pylablib.devices.LighthousePhotonics.base.LighthousePhotonicsError attribute\)](#)), 611
- [args \(pylablib.devices.M2.solstis.M2CommunicationError attribute\)](#)), 614
- [args \(pylablib.devices.M2.solstis.M2Error attribute\)](#)), 614
- [args \(pylablib.devices.Newport.base.NewportBackendError attribute\)](#)), 627
- [args \(pylablib.devices.Newport.base.NewportError attribute\)](#)), 627
- [args \(pylablib.devices.NI.daq.NIDAQmxError attribute\)](#)), 620
- [args \(pylablib.devices.NI.daq.NIError attribute\)](#)), 620
- [args \(pylablib.devices.Ophir.base.OphirBackendError attribute\)](#)), 638
- [args \(pylablib.devices.Ophir.base.OphirError attribute\)](#)), 638
- [args \(pylablib.devices.OZOptics.base.OZOpticsBackendError attribute\)](#)), 631
- [args \(pylablib.devices.OZOptics.base.OZOpticsError attribute\)](#)), 631
- [args \(pylablib.devices.Pfeiffer.base.PfeifferBackendError attribute\)](#)), 653
- [args \(pylablib.devices.Pfeiffer.base.PfeifferError attribute\)](#)), 653
- [args \(pylablib.devices.PhysikInstrumente.base.PhysikInstrumenteBackendError attribute\)](#)), 703
- [args \(pylablib.devices.PhysikInstrumente.base.PhysikInstrumenteError attribute\)](#)), 702
- [args \(pylablib.devices.Sirah.Matisse.GenericSirahBackendError attribute\)](#)), 733
- [args \(pylablib.devices.Sirah.Matisse.GenericSirahError attribute\)](#)), 733
- [args \(pylablib.devices.Tektronix.base.TektronixBackendError attribute\)](#)), 748
- [args \(pylablib.devices.Tektronix.base.TektronixError attribute\)](#)), 748
- [args \(pylablib.devices.Thorlabs.base.ThorlabsBackendError attribute\)](#)), 776
- [args \(pylablib.devices.Thorlabs.base.ThorlabsError attribute\)](#)), 776
- [args \(pylablib.devices.Thorlabs.base.ThorlabsTimeoutError attribute\)](#)), 776
- [args \(pylablib.devices.Toptica.base.TopticaBackendError attribute\)](#)), 823
- [args \(pylablib.devices.Toptica.base.TopticaError attribute\)](#)), 823
- [args \(pylablib.devices.Trinamic.base.TrinamicBackendError attribute\)](#)), 826
- [args \(pylablib.devices.Trinamic.base.TrinamicError attribute\)](#)), 826
- [Array1DWrapper \(class in pylablib.core.dataproc.table\\_wrap\)](#)), 126
- [Array1DWrapper.Accessor \(class in pylablib.core.dataproc.table\\_wrap\)](#)), 127
- [Array2DWrapper \(class in pylablib.core.dataproc.table\\_wrap\)](#)), 129
- [Array2DWrapper.ColumnAccessor \(class in pylablib.core.dataproc.table\\_wrap\)](#)), 130
- [Array2DWrapper.RowAccessor \(class in pylablib.core.dataproc.table\\_wrap\)](#)), 130
- [Array2DWrapper.TableAccessor \(class in pylablib.core.dataproc.table\\_wrap\)](#)), 131
- [array\\_replaced\(\) \(pylablib.core.dataproc.table\\_wrap.Array1DWrapper method\)](#)), 127
- [array\\_replaced\(\) \(pylablib.core.dataproc.table\\_wrap.Array2DWrapper method\)](#)), 131
- [array\\_replaced\(\) \(pylablib.core.dataproc.table\\_wrap.DataFrame2DWrapper method\)](#)), 133
- [array\\_replaced\(\) \(pylablib.core.dataproc.table\\_wrap.I1DWrapper method\)](#)), 126
- [array\\_replaced\(\) \(pylablib.core.dataproc.table\\_wrap.I2DWrapper method\)](#)), 129
- [array\\_replaced\(\) \(pylablib.core.dataproc.table\\_wrap.Series1DWrapper method\)](#)), 128
- [array\\_wrap\(\) \(in module pylablib.core.utils.net\)](#)), 388
- [as\\_array\(\) \(in module pylablib.core.utils.array\\_utils\)](#)), 321
- [as\\_builtin\\_bytes\(\) \(in module pylablib.core.utils.py3\)](#)), 393

`as_bytes()` (in module `pylablib.core.utils.py3`), 393  
`as_container()` (in module `pylablib.core.utils.general`), 375  
`as_datatype()` (in module `pylablib.core.utils.py3`), 393  
`as_dict()` (in module `pylablib.core.utils.dictionary`), 326  
`as_dict()` (`pylablib.core.utils.dictionary.Dictionary` method), 332  
`as_dict()` (`pylablib.core.utils.dictionary.DictionaryPointer` method), 336  
`as_dict()` (`pylablib.core.utils.dictionary.FilterTree` method), 353  
`as_dict()` (`pylablib.core.utils.dictionary.PrefixTree` method), 345  
`as_dictionary()` (in module `pylablib.core.utils.dictionary`), 326  
`as_dictionary()` (`pylablib.core.utils.dictionary.Dictionary` static method), 327  
`as_dictionary()` (`pylablib.core.utils.dictionary.DictionaryPointer` static method), 336  
`as_dictionary()` (`pylablib.core.utils.dictionary.FilterTree` static method), 353  
`as_dictionary()` (`pylablib.core.utils.dictionary.PrefixTree` static method), 345  
`as_formatter()` (in module `pylablib.core.gui.formatter`), 262  
`as_json()` (`pylablib.core.utils.dictionary.Dictionary` method), 332  
`as_json()` (`pylablib.core.utils.dictionary.DictionaryPointer` method), 336  
`as_json()` (`pylablib.core.utils.dictionary.FilterTree` method), 353  
`as_json()` (`pylablib.core.utils.dictionary.PrefixTree` method), 345  
`as_kwargs()` (`pylablib.core.utils.functions.FunctionSignature` method), 370  
`as_limiter()` (in module `pylablib.core.gui.limiter`), 262  
`as_obj_prop()` (in module `pylablib.core.utils.functions`), 374  
`as_pandas()` (`pylablib.core.utils.dictionary.Dictionary` method), 332  
`as_pandas()` (`pylablib.core.utils.dictionary.DictionaryPointer` method), 336  
`as_pandas()` (`pylablib.core.utils.dictionary.FilterTree` method), 353  
`as_pandas()` (`pylablib.core.utils.dictionary.PrefixTree` method), 345  
`as_sequence()` (in module `pylablib.core.utils.funcargparse`), 369  
`as_simple_func()` (`pylablib.core.utils.functions.FunctionSignature` method), 371  
`as_str()` (in module `pylablib.core.utils.py3`), 393  
`as_text()` (`pylablib.devices.DCAM.DCAM.DCAMAttribute` method), 549  
`asdict()` (`pylablib.core.utils.dictionary.Dictionary` method), 332  
`asdict()` (`pylablib.core.utils.dictionary.DictionaryPointer` method), 337  
`asdict()` (`pylablib.core.utils.dictionary.FilterTree` method), 353  
`asdict()` (`pylablib.core.utils.dictionary.PrefixTree` method), 345  
`ask()` (`pylablib.core.devio.comm_backend.FT232DeviceBackend` method), 151  
`ask()` (`pylablib.core.devio.comm_backend.IDeviceCommBackend` method), 144  
`ask()` (`pylablib.core.devio.comm_backend.NetworkDeviceBackend` method), 154  
`ask()` (`pylablib.core.devio.comm_backend.PyUSBDeviceBackend` method), 156  
`ask()` (`pylablib.core.devio.comm_backend.RecordedDeviceBackend` method), 158  
`ask()` (`pylablib.core.devio.comm_backend.SerialDeviceBackend` method), 149  
`ask()` (`pylablib.core.devio.comm_backend.VisaDeviceBackend` method), 146  
`ask()` (`pylablib.core.devio.SCPi.SCPiDevice` method), 140  
`ask()` (`pylablib.devices.AWG.generic.GenericAWG` method), 406  
`ask()` (`pylablib.devices.AWG.specific.Agilent33220A` method), 415  
`ask()` (`pylablib.devices.AWG.specific.Agilent33500` method), 409  
`ask()` (`pylablib.devices.AWG.specific.InstekAFG2000` method), 427  
`ask()` (`pylablib.devices.AWG.specific.InstekAFG2225` method), 421  
`ask()` (`pylablib.devices.AWG.specific.RigolDG1000` method), 444  
`ask()` (`pylablib.devices.AWG.specific.RSInstekAFG21000` method), 433  
`ask()` (`pylablib.devices.AWG.specific.TektronixAFG1000` method), 439  
`ask()` (`pylablib.devices.Cryocon.base.Cryocon1x` method), 536  
`ask()` (`pylablib.devices.Cryomagnetics.base.LM500` method), 541  
`ask()` (`pylablib.devices.Cryomagnetics.base.LM510` method), 544  
`ask()` (`pylablib.devices.Lakeshore.base.Lakeshore218`

- method), 595
- ask () (pylablib.devices.Lakeshore.base.Lakeshore370 method), 600
- ask () (pylablib.devices.Sirah.Matisse.SirahMatisse method), 740
- ask () (pylablib.devices.Tektronix.base.DPO2000 method), 762
- ask () (pylablib.devices.Tektronix.base.ITektronixScope method), 753
- ask () (pylablib.devices.Tektronix.base.TDS2000 method), 756
- ask () (pylablib.devices.Thorlabs.serial.FW method), 812
- ask () (pylablib.devices.Thorlabs.serial.FWv1 method), 816
- ask () (pylablib.devices.Thorlabs.serial.MDT69xA method), 820
- ask () (pylablib.devices.Thorlabs.serial.ThorlabsSerialInterface method), 810
- atm\_adj (pylablib.devices.Leybold.base.TITR90Status attribute), 609
- AttocubeBackendError, 510
- AttocubeError, 510
- attr (pylablib.core.utils.observer\_pool.ObserverPool.Observer attribute), 392
- attr\_url (pylablib.devices.IMAQdx.IMAQdx.TCameraInfo attribute), 576
- AttrObjectCall (class in py-lablib.core.utils.functions), 372
- AttrObjectProperty (class in py-lablib.core.utils.functions), 373
- au1 () (in module pylablib.core.utils.nbtools), 387
- au2 () (in module pylablib.core.utils.nbtools), 387
- au4 () (in module pylablib.core.utils.nbtools), 387
- au8 () (in module pylablib.core.utils.nbtools), 387
- autodetect\_backend () (in module py-lablib.core.devio.comm\_backend), 159
- autodetect\_motors () (py-lablib.devices.Newport.picomotor.Picomotor method), 629
- autoloop (pylablib.devices.NI.daq.TVoltageOutputClockParameters attribute), 621
- autorange (pylablib.devices.Lakeshore.base.TLakeshore370RangeSettings attribute), 598
- available (pylablib.devices.Basler.pylon.BaslerPylonAttribute attribute), 512
- available (pylablib.devices.Photometrics.pvcam.PvcamAttribute attribute), 660
- available\_samples () (py-lablib.devices.NI.daq.NIDAQ method), 623
- average\_interpolate\_1D () (in module py-lablib.core.dataproc.interpolate), 123
- avg (pylablib.devices.Sirah.Matisse.TPiezoetFeedbackParameters attribute), 734
- avg (pylablib.devices.Sirah.Matisse.TThinnetCtlParameters attribute), 734
- ## B
- backend\_error () (in module py-lablib.core.devio.comm\_backend), 160
- BackendError (pylablib.core.devio.comm\_backend.FT232DeviceBackend attribute), 150
- BackendError (pylablib.core.devio.comm\_backend.IDeviceCommBackend attribute), 143
- BackendError (pylablib.core.devio.comm\_backend.NetworkDeviceBackend attribute), 153
- BackendError (pylablib.core.devio.comm\_backend.PyUSBDeviceBackend attribute), 155
- BackendError (pylablib.core.devio.comm\_backend.RecordedDeviceBackend attribute), 158
- BackendError (pylablib.core.devio.comm\_backend.SerialDeviceBackend attribute), 147
- BackendError (pylablib.core.devio.comm\_backend.VisaDeviceBackend attribute), 145
- BackendError (pylablib.core.devio.SCPIDevice attribute), 138
- BackendError (pylablib.devices.AWG.generic.GenericAWG attribute), 405
- BackendError (pylablib.devices.AWG.specific.Agilent33220A attribute), 414
- BackendError (pylablib.devices.AWG.specific.Agilent33500 attribute), 409
- BackendError (pylablib.devices.AWG.specific.InstekAFG2000 attribute), 426
- BackendError (pylablib.devices.AWG.specific.InstekAFG2225 attribute), 421
- BackendError (pylablib.devices.AWG.specific.RigolDG1000 attribute), 444
- BackendError (pylablib.devices.AWG.specific.RSInstekAFG21000 attribute), 432
- BackendError (pylablib.devices.AWG.specific.TektronixAFG1000 attribute), 438
- BackendError (pylablib.devices.Cryocon.base.Cryocon1x attribute), 536
- BackendError (pylablib.devices.Cryomagnetics.base.LM500 attribute), 541
- BackendError (pylablib.devices.Cryomagnetics.base.LM510 attribute), 544
- BackendError (pylablib.devices.Lakeshore.base.Lakeshore218 attribute), 595
- BackendError (pylablib.devices.Lakeshore.base.Lakeshore370 attribute), 600
- BackendError (pylablib.devices.M2.solstis.Solstis attribute), 614
- BackendError (pylablib.devices.NI.daq.NIDAQ attribute), 621
- BackendError (pylablib.devices.Sirah.Matisse.SirahMatisse attribute), 740



- BackendError (pylablib.devices.Tektronix.base.DPO2000 attribute), 762
- BackendError (pylablib.devices.Tektronix.base.ITektronixScope attribute), 752
- BackendError (pylablib.devices.Tektronix.base.TDS2000 attribute), 755
- BackendError (pylablib.devices.Thorlabs.serial.FW attribute), 813
- BackendError (pylablib.devices.Thorlabs.serial.FWv1 attribute), 816
- BackendError (pylablib.devices.Thorlabs.serial.MDT690 attribute), 820
- BackendError (pylablib.devices.Thorlabs.serial.ThorlabsSerialInterface attribute), 810
- BackendLogger (class in py-lablib.core.devio.backend\_logger), 141
- backlash\_distance (py-lablib.devices.Thorlabs.kinesis.TGenMoveParams attribute), 785
- Baseline (class in pylablib.core.dataproc.feature), 108
- baseplate (pylablib.devices.Toptica.ibeam.TTemperatures attribute), 824
- BasicKinesisDevice (class in py-lablib.devices.Thorlabs.kinesis), 781
- BasicKinesisDevice.CommData (class in py-lablib.devices.Thorlabs.kinesis), 782
- BasicKinesisDevice.CommShort (class in py-lablib.devices.Thorlabs.kinesis), 782
- BasicKinesisDevice.NoParameterCaller (class in pylablib.devices.Thorlabs.kinesis), 783
- BaslerPylonAttribute (class in py-lablib.devices.Basler.pylon), 511
- BaslerPylonCamera (class in py-lablib.devices.Basler.pylon), 514
- BaslerPylonCamera.BufferManager (class in pylablib.devices.Basler.pylon), 515
- BaslerPylonCamera.NoParameterCaller (class in pylablib.devices.Basler.pylon), 517
- BaslerPylonCamera.ScheduleLooper (class in pylablib.devices.Basler.pylon), 516
- bifi\_clear\_errors() (py-lablib.devices.Sirah.Matisse.SirahMatisse method), 737
- bifi\_get\_position() (py-lablib.devices.Sirah.Matisse.SirahMatisse method), 737
- bifi\_get\_range() (py-lablib.devices.Sirah.Matisse.SirahMatisse method), 737
- bifi\_get\_status() (py-lablib.devices.Sirah.Matisse.SirahMatisse method), 737
- bifi\_get\_status\_n() (py-lablib.devices.Sirah.Matisse.SirahMatisse method), 737
- bifi\_is\_moving() (py-lablib.devices.Sirah.Matisse.SirahMatisse method), 737
- bifi\_move\_to() (py-lablib.devices.Sirah.Matisse.SirahMatisse method), 737
- bifi\_stop() (pylablib.devices.Sirah.Matisse.SirahMatisse method), 737
- BinaryTableInputFileFormatter (class in py-lablib.core.fileio.loadfile), 178
- bind() (pylablib.core.dataproc.callable.FunctionCallable method), 107
- bind() (pylablib.core.dataproc.callable.ICallable method), 104
- bind() (pylablib.core.dataproc.callable.JoinedCallable method), 106
- bind() (pylablib.core.dataproc.callable.MethodCallable method), 107
- bind() (pylablib.core.dataproc.callable.MultiplexedCallable method), 105
- bind\_namelist() (py-lablib.core.dataproc.callable.FunctionCallable method), 107
- bind\_namelist() (py-lablib.core.dataproc.callable.ICallable method), 104
- bind\_namelist() (py-lablib.core.dataproc.callable.JoinedCallable method), 106
- bind\_namelist() (py-lablib.core.dataproc.callable.MethodCallable method), 108
- bind\_namelist() (py-lablib.core.dataproc.callable.MultiplexedCallable method), 105
- binning\_average() (in module py-lablib.core.dataproc.filters), 112
- bipolar (pylablib.devices.Lakeshore.base.TLakeshore218AnalogSettings attribute), 592
- bipolar (pylablib.devices.Lakeshore.base.TLakeshore370AnalogSettings attribute), 598
- BitFlowCamera (class in py-lablib.devices.BitFlow.BitFlow), 527
- BitFlowCamera.BufferManager (class in py-lablib.devices.BitFlow.BitFlow), 527
- BitFlowCamera.NoParameterCaller (class in pylablib.devices.BitFlow.BitFlow), 527

BitFlowError, 520  
 BitFlowFrameGrabber (class in py-  
     lablib.devices.BitFlow.BitFlow), 521  
 BitFlowFrameGrabber.BufferManager (class  
     in pylablib.devices.BitFlow.BitFlow), 522  
 BitFlowFrameGrabber.NoParameterCaller  
     (class in pylablib.devices.BitFlow.BitFlow),  
     523  
 BitFlowTimeoutError, 520  
 bits2int() (in module pylablib.core.utils.strpack),  
     401  
 bk\_freq (pylablib.devices.Thorlabs.elliptec.TMotorInfo  
     attribute), 777  
 blink() (pylablib.devices.Thorlabs.kinesis.BasicKinesisDevice  
     method), 783  
 blink() (pylablib.devices.Thorlabs.kinesis.KinesisDevice  
     method), 788  
 blink() (pylablib.devices.Thorlabs.kinesis.KinesisMotor  
     method), 799  
 blink() (pylablib.devices.Thorlabs.kinesis.KinesisPiezoMotor  
     method), 803  
 blink() (pylablib.devices.Thorlabs.kinesis.KinesisQuadDetector  
     method), 808  
 blink() (pylablib.devices.Thorlabs.kinesis.MFF  
     method), 792  
 blocking\_control\_signals() (py-  
     lablib.core.thread.controller.QTaskThread  
     method), 308  
 blocking\_control\_signals() (py-  
     lablib.core.thread.controller.QThreadController  
     method), 294  
 BonitoError, 450  
 BonitoIMAQCamera (class in py-  
     lablib.devices.AlliedVision.Bonito), 456  
 BonitoIMAQCamera.NoParameterCaller (class  
     in pylablib.devices.AlliedVision.Bonito), 456  
 BonitoStatusLineChecker (class in py-  
     lablib.devices.AlliedVision.Bonito), 464  
 branch\_copy() (py-  
     lablib.core.utils.dictionary.Dictionary method),  
     331  
 branch\_copy() (py-  
     lablib.core.utils.dictionary.DictionaryPointer  
     method), 337  
 branch\_copy() (py-  
     lablib.core.utils.dictionary.FilterTree method),  
     354  
 branch\_copy() (py-  
     lablib.core.utils.dictionary.PrefixTree method),  
     346  
 branch\_pointer() (py-  
     lablib.core.utils.dictionary.Dictionary method),  
     332  
 branch\_pointer() (py-  
     lablib.core.utils.dictionary.DictionaryPointer  
     method), 336  
 branch\_pointer() (py-  
     lablib.core.utils.dictionary.FilterTree method),  
     354  
 branch\_pointer() (py-  
     lablib.core.utils.dictionary.PrefixTree method),  
     346  
 buffconv() (in module py-  
     lablib.core.utils.ctypes\_wrap), 324  
 buffer\_size (pylablib.devices.interface.camera.TFramesStatus  
     attribute), 830  
 buffprep() (in module py-  
     lablib.core.utils.ctypes\_wrap), 324  
 build\_call() (pylablib.core.thread.callsync.QDirectCallScheduler  
     method), 285  
 build\_call() (pylablib.core.thread.callsync.QMulticastThreadCallScheduler  
     method), 292  
 build\_call() (pylablib.core.thread.callsync.QMultiQueueScheduler  
     method), 290  
 build\_call() (pylablib.core.thread.callsync.QQueueLengthLimitScheduler  
     method), 288  
 build\_call() (pylablib.core.thread.callsync.QQueueScheduler  
     method), 286  
 build\_call() (pylablib.core.thread.callsync.QQueueSizeLimitScheduler  
     method), 289  
 build\_call() (pylablib.core.thread.callsync.QScheduler  
     method), 284  
 build\_call() (pylablib.core.thread.callsync.QThreadCallScheduler  
     method), 291  
 build\_call\_info() (py-  
     lablib.core.thread.callsync.QDirectCallScheduler  
     method), 285  
 build\_call\_info() (py-  
     lablib.core.thread.callsync.QMulticastThreadCallScheduler  
     method), 292  
 build\_call\_info() (py-  
     lablib.core.thread.callsync.QQueueLengthLimitScheduler  
     method), 288  
 build\_call\_info() (py-  
     lablib.core.thread.callsync.QQueueScheduler  
     method), 287  
 build\_call\_info() (py-  
     lablib.core.thread.callsync.QQueueSizeLimitScheduler  
     method), 289  
 build\_call\_info() (py-  
     lablib.core.thread.callsync.QScheduler  
     method), 284  
 build\_call\_info() (py-  
     lablib.core.thread.callsync.QThreadCallScheduler  
     method), 291  
 build\_children\_tree() (in module py-  
     lablib.core.gui.value\_handling), 265  
 build\_file\_format() (in module py-

[lablib.core.fileio.loadfile](#)), 179  
 bus ([pylablib.devices.IMAQdx.IMAQdx.TCameraInfo](#) attribute), 576  
 bus\_type ([pylablib.devices.IMAQdx.IMAQdx.TDeviceInfo](#) attribute), 578  
 byref() ([pylablib.core.utils.ctypes\\_wrap.CFunctionWrapper](#) method), 322  
 bytes2int() (in module [pylablib.core.utils.strpack](#)), 401  
**C**  
 c2xy() (in module [pylablib.core.dataproc.utils](#)), 137  
 c\_array (in module [pylablib.core.utils.nbtools](#)), 387  
 ca ([pylablib.core.thread.controller.QTaskThread](#) attribute), 301  
 cacheable ([pylablib.core.utils.observer\\_pool.ObserverPool.ObserverPool](#) attribute), 392  
 cad ([pylablib.core.thread.controller.QTaskThread](#) attribute), 301  
 cai ([pylablib.core.thread.controller.QTaskThread](#) attribute), 301  
 cal\_date ([pylablib.devices.LaserQuantum.base.TDeviceInfo](#) attribute), 603  
 calibrate() ([pylablib.devices.HighFinesse.wlm.WLM](#) method), 559  
 call\_added() ([pylablib.core.thread.callsync.QQueueLengthLimitScheduler](#) method), 287  
 call\_added() ([pylablib.core.thread.callsync.QQueueScheduler](#) method), 286  
 call\_added() ([pylablib.core.thread.callsync.QQueueSizeLimitScheduler](#) method), 289  
 call\_command() (py-[lablib.core.thread.controller.QTaskThread](#) method), 307  
 call\_command() (py-[lablib.devices.Andor.AndorSDK3.AndorSDK3Attribute](#) method), 478  
 call\_command() (py-[lablib.devices.Andor.AndorSDK3.AndorSDK3Camera](#) method), 480  
 call\_command() (py-[lablib.devices.Basler.pylon.BaslerPylonAttribute](#) method), 513  
 call\_command() (py-[lablib.devices.Basler.pylon.BaslerPylonCamera](#) method), 514  
 call\_command() (py-[lablib.devices.PhotonFocus.PhotonFocus.IPhotonFocusCamera](#) method), 673  
 call\_command() (py-[lablib.devices.PhotonFocus.PhotonFocus.PFCamAttribute](#) method), 671  
 call\_command() (py-[lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusBitFlowController](#) method), 695  
 call\_command() (py-[lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusIMQCamera](#) method), 678  
 call\_command() (py-[lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusSiSoCamera](#) method), 687  
 call\_command\_direct() (py-[lablib.core.thread.controller.QTaskThread](#) method), 306  
 call\_cut\_args() (in module [pylablib.core.utils.functions](#)), 372  
 call\_in\_gui\_thread() (in module [pylablib.core.thread.controller](#)), 293  
 call\_in\_thread() (in module [pylablib.core.thread.controller](#)), 293  
 call\_in\_thread\_callback() (py-[lablib.core.thread.controller.QTaskThread](#) method), 308  
 call\_in\_thread\_callback() (py-[lablib.core.thread.controller.QThreadController](#) method), 300  
 call\_in\_thread\_commsync() (py-[lablib.core.thread.controller.QTaskThread](#) method), 307  
 call\_in\_thread\_sync() (py-[lablib.core.thread.controller.QTaskThread](#) method), 308  
 call\_in\_thread\_sync() (py-[lablib.core.thread.controller.QThreadController](#) method), 300  
 call\_limit() (in module [pylablib.core.utils.general](#)), 378  
 call\_on\_exception (py-[lablib.core.thread.callsync.QScheduledCall.Callback](#) attribute), 283  
 call\_on\_unschedule (py-[lablib.core.thread.callsync.QScheduledCall.Callback](#) attribute), 283  
 call\_popped() (py-[lablib.core.thread.callsync.QQueueLengthLimitScheduler](#) method), 287  
 call\_popped() (py-[lablib.core.thread.callsync.QQueueScheduler](#) method), 286  
 call\_popped() (py-[lablib.core.thread.callsync.QQueueSizeLimitScheduler](#) method), 289  
 call\_thread\_method() (py-[lablib.core.thread.controller.QTaskThread](#) method), 308  
 call\_thread\_method() (py-[lablib.core.thread.controller.QThreadController](#) method), 298

`call_time` (`pylablib.core.thread.callsync.TDefaultCallInfo` attribute), 284  
`callback` (`pylablib.core.utils.observer_pool.ObserverPool.Observer` method), 270  
`can_id` (`pylablib.devices.uc480.uc480.TCameraInfo` attribute), 864  
`can_id` (`pylablib.devices.uc480.uc480.TDeviceInfo` attribute), 865  
`camera_file` (`pylablib.devices.IMAQdx.IMAQdx.TCameraInfo` attribute), 576  
`camera_model` (`pylablib.devices.Andor.AndorSDK3.TDeviceInfo` attribute), 478  
`camera_name` (`pylablib.devices.Andor.AndorSDK3.TDeviceInfo` attribute), 478  
`camera_type` (`pylablib.devices.uc480.uc480.TDeviceInfo` attribute), 865  
`camera_version` (`pylablib.devices.DCAM.DCAM.TDeviceInfo` attribute), 550  
`CameraFileEditor` (class in `pylablib.devices.BitFlow.BitFlow`), 532  
`camerastamp` (`pylablib.devices.DCAM.DCAM.TFrameInfo` attribute), 550  
`can_change` () (`pylablib.core.thread.utils.ReadChangeLock` method), 321  
`can_read` () (`pylablib.core.thread.utils.ReadChangeLock` method), 321  
`can_schedule` () (`pylablib.core.thread.callsync.QQueueLengthLimitScheduler` method), 287  
`can_schedule` () (`pylablib.core.thread.callsync.QQueueScheduler` method), 286  
`can_schedule` () (`pylablib.core.thread.callsync.QQueueSizeLimitScheduler` method), 289  
`can_set_online` (`pylablib.devices.PrincetonInstruments.picam.PicamAttribute` attribute), 710  
`can_set_value` () (`pylablib.core.gui.value_handling.CheckboxValueHandler` method), 272  
`can_set_value` () (`pylablib.core.gui.value_handling.ComboBoxValueHandler` method), 275  
`can_set_value` () (`pylablib.core.gui.value_handling.IBoolValueHandler` method), 271  
`can_set_value` () (`pylablib.core.gui.value_handling.ISingleValueHandler` method), 269  
`can_set_value` () (`pylablib.core.gui.value_handling.IValueHandler` method), 266  
`can_set_value` () (`pylablib.core.gui.value_handling.LabelValueHandler` method), 270  
`can_set_value` () (`pylablib.core.gui.value_handling.LineEditValueHandler` method), 270  
`can_set_value` () (`pylablib.core.gui.value_handling.ProgressBarValueHandler` method), 275  
`can_set_value` () (`pylablib.core.gui.value_handling.PropertyValueHandler` method), 267  
`can_set_value` () (`pylablib.core.gui.value_handling.PushButtonValueHandler` method), 273  
`can_set_value` () (`pylablib.core.gui.value_handling.StandardValueHandler` method), 268  
`can_set_value` () (`pylablib.core.gui.value_handling.ToolButtonValueHandler` method), 274  
`capabilities` (`pylablib.devices.Ophir.base.THeadInfo` attribute), 640  
`case_sensitive_path` () (in module `pylablib.core.utils.files`), 362  
`catcher` (`pylablib.core.gui.limiter.NumberLimit` method), 262  
`category` (`pylablib.devices.SiliconSoftware.fgrab.TFullAppletInfo` attribute), 718  
`center` () (`pylablib.core.dataproc.image.ROI` method), 121  
`CallableFunctionWrapper` (class in `pylablib.core.utils.ctypes_wrap`), 322  
`change_addr` () (`pylablib.devices.Thorlabs.elliptec.ElliptecMotor` method), 779  
`change_batch_job_parameters` () (`pylablib.core.thread.controller.QTaskThread` method), 303  
`change_job_period` () (`pylablib.core.thread.controller.QTaskThread` method), 303  
`change_max_len` () (`pylablib.core.thread.callsync.QQueueLengthLimitScheduler` method), 287  
`change_max_size` () (`pylablib.core.thread.callsync.QQueueSizeLimitScheduler` method), 289  
`change_period` () (`pylablib.core.thread.controller.QTaskThread.Job` method), 303



`change_period()` (pylablib.core.utils.general.Timer method), 379  
`changed_from` (pylablib.core.utils.dictionary.DictionaryDiff attribute), 334  
`changed_to` (pylablib.core.utils.dictionary.DictionaryDiff attribute), 334  
`changing()` (pylablib.core.thread.utils.ReadChangeLock method), 321  
`channel` (pylablib.devices.Lakeshore.base.TLakeshore218AnalogSettings attribute), 592  
`channel` (pylablib.devices.Lakeshore.base.TLakeshore370AnalogSettings attribute), 598  
`channel` (pylablib.devices.Pfeiffer.base.TTPG260SwitchSettings attribute), 653  
`characters_written` (py-lablib.core.thread.threadprop.TimeoutThreadError attribute), 319  
`characters_written` (py-lablib.core.utils.net.SocketError attribute), 388  
`characters_written` (py-lablib.core.utils.net.SocketTimeout attribute), 388  
`check_alias()` (py-lablib.core.devio.interface.EnumParameterClass method), 166  
`check_alias()` (py-lablib.core.devio.interface.FunctionParameterClass method), 167  
`check_alias()` (py-lablib.core.devio.interface.ICheckingParameterClass method), 163  
`check_alias()` (py-lablib.core.devio.interface.IEnumParameterClass method), 165  
`check_alias()` (py-lablib.core.devio.interface.RangeParameterClass method), 164  
`check_background_comm()` (py-lablib.devices.Thorlabs.elliptec.ElliptecMotor method), 779  
`check_background_comm()` (py-lablib.devices.Thorlabs.kinesis.BasicKinesisDevice method), 783  
`check_background_comm()` (py-lablib.devices.Thorlabs.kinesis.KinesisDevice method), 788  
`check_background_comm()` (py-lablib.devices.Thorlabs.kinesis.KinesisMotor method), 799  
`check_background_comm()` (py-lablib.devices.Thorlabs.kinesis.KinesisPiezoMotor method), 803  
`check_background_comm()` (py-lablib.devices.Thorlabs.kinesis.KinesisQuadDetector method), 808  
`check_background_comm()` (py-lablib.devices.Thorlabs.kinesis.MFF method), 792  
`check_fast_scan_start_report()` (py-lablib.devices.M2.solstis.Solstis method), 618  
`check_settings_tuning_report()` (py-lablib.devices.M2.solstis.Solstis method), 615  
`check_grabber_association()` (in module py-lablib.devices.AlliedVision.Bonito), 464  
`check_grabber_association()` (in module py-lablib.devices.PhotonFocus.PhotonFocus), 701  
`check_indices()` (py-lablib.devices.AlliedVision.Bonito.BonitoStatusLineChecker method), 465  
`check_indices()` (py-lablib.devices.interface.camera.StatusLineChecker method), 861  
`check_indices()` (py-lablib.devices.PCO.SC2.StatusLineChecker method), 653  
`check_indices()` (py-lablib.devices.PhotonFocus.PhotonFocus.StatusLineChecker method), 702  
`check_limit()` (py-lablib.devices.Attocube.anc350.ANC350 method), 507  
`check_limit_error()` (py-lablib.devices.Arcus.performax.Performax2EXStage method), 494  
`check_limit_error()` (py-lablib.devices.Arcus.performax.Performax4EXStage method), 491  
`check_limit_error()` (py-lablib.devices.Arcus.performax.PerformaxDMXJSASStage method), 498  
`check_messages()` (py-lablib.core.thread.controller.QTaskThread method), 308  
`check_messages()` (py-lablib.core.thread.controller.QThreadController method), 294  
`check_parameter_range()` (in module py-lablib.core.utils.funcargparse), 369  
`check_report()` (pylablib.devices.M2.solstis.Solstis method), 615  
`check_tell()` (pylablib.devices.Attocube.anc350.ANC350 method), 507  
`check_terascan_start_report()` (py-lablib.devices.M2.solstis.Solstis method), 618

`check_terascan_update()` (pylablib.devices.M2.solstis.Solstis method), 618  
`check_value()` (pylablib.core.devio.interface.EnumParameterClass method), 166  
`check_value()` (pylablib.core.devio.interface.FunctionParameterClass method), 167  
`check_value()` (pylablib.core.devio.interface.ICheckingParameterClass method), 164  
`check_value()` (pylablib.core.devio.interface.IEnumParameterClass method), 165  
`check_value()` (pylablib.core.devio.interface.RangeParameterClass method), 164  
`CheckboxValueHandler` (class in pylablib.core.gui.value\_handling), 272  
`chip` (pylablib.devices.Photometrics.pvcam.TDDeviceInfo attribute), 661  
`ChunkBufferManager` (class in pylablib.devices.interface.camera), 836  
`ChunksAccumulator` (class in pylablib.core.fileio.parse\_csv), 188  
`class_tuple_to_dict()` (in module pylablib.core.utils.ctypes\_wrap), 325  
`clean_dir()` (in module pylablib.core.utils.files), 363  
`clean_layout()` (in module pylablib.core.gui.utils), 263  
`clean_modes()` (pylablib.devices.BitFlow.BitFlow.CameraFileEditor method), 532  
`cleanup` (pylablib.core.thread.controller.QTaskThread.TBAttribute attribute), 302  
`cleanup()` (pylablib.devices.Thorlabs.TLCAmra.ThorlabsTLCAmra.TBAttribute method), 771  
`clear()` (pylablib.core.gui.widgets.container.IQContainer method), 202  
`clear()` (pylablib.core.gui.widgets.container.IQWidgetContainer method), 205  
`clear()` (pylablib.core.gui.widgets.container.QContainer method), 203  
`clear()` (pylablib.core.gui.widgets.container.QDialogContainer method), 218  
`clear()` (pylablib.core.gui.widgets.container.QFrameContainer method), 214  
`clear()` (pylablib.core.gui.widgets.container.QGroupBoxContainer method), 222  
`clear()` (pylablib.core.gui.widgets.container.QScrollAreaContainer method), 229  
`clear()` (pylablib.core.gui.widgets.container.QScrollAreaContainer.QContainerWidget method), 226  
`clear()` (pylablib.core.gui.widgets.container.QTabContainer method), 231  
`clear()` (pylablib.core.gui.widgets.container.QWidgetContainer method), 211  
`clear()` (pylablib.core.gui.widgets.layout\_manager.IQLayoutManagedWidget method), 240  
`clear()` (pylablib.core.gui.widgets.layout\_manager.QLayoutManagedWidget method), 240  
`clear()` (pylablib.core.gui.widgets.param\_table.ParamTable method), 248  
`clear()` (pylablib.core.gui.widgets.param\_table.StatusTable method), 258  
`clear()` (pylablib.core.thread.callsync.QDirectCallScheduler method), 285  
`clear()` (pylablib.core.thread.callsync.QMulticastThreadCallScheduler method), 292  
`clear()` (pylablib.core.thread.callsync.QQueueLengthLimitScheduler method), 288  
`clear()` (pylablib.core.thread.callsync.QQueueScheduler method), 286  
`clear()` (pylablib.core.thread.callsync.QQueueSizeLimitScheduler method), 290  
`clear()` (pylablib.core.thread.callsync.QScheduler method), 285  
`clear()` (pylablib.core.thread.callsync.QThreadCallScheduler method), 291  
`clear()` (pylablib.core.thread.controller.QTaskThread.Job method), 303  
`clear_acquisition()` (pylablib.devices.AlliedVision.Bonito.BonitoIMAQCamera method), 456  
`clear_acquisition()` (pylablib.devices.AlliedVision.Bonito.IBonitoCamera method), 452  
`clear_acquisition()` (pylablib.devices.Andor.AndorSDK2.AndorSDK2Camera method), 482  
`clear_acquisition()` (pylablib.devices.Andor.AndorSDK3.AndorSDK3Camera method), 482  
`clear_acquisition()` (pylablib.devices.Basler.pylon.BaslerPylonCamera method), 516  
`clear_acquisition()` (pylablib.devices.BitFlow.BitFlow.BitFlowCamera method), 527  
`clear_acquisition()` (pylablib.devices.BitFlow.BitFlow.BitFlowFrameGrabber method), 523  
`clear_acquisition()` (pylablib.devices.DCAM.DCAM.DCAMCamera method), 552  
`clear_acquisition()` (pylablib.devices.IMAQ.IMAQ.IMAQCamera method), 552

method), 569

clear\_acquisition() (py-lablib.devices.IMAQ.IMAQ.IMAQFrameGrabber method), 565

clear\_acquisition() (py-lablib.devices.IMAQdx.IMAQdx.EthernetIMAQdxCamera method), 584

clear\_acquisition() (py-lablib.devices.IMAQdx.IMAQdx.IMAQdxCamera method), 580

clear\_acquisition() (py-lablib.devices.interface.camera.IAttributeCamera method), 838

clear\_acquisition() (py-lablib.devices.interface.camera.IBinROICamera method), 857

clear\_acquisition() (py-lablib.devices.interface.camera.ICamera method), 831

clear\_acquisition() (py-lablib.devices.interface.camera.IExposureCamera method), 847

clear\_acquisition() (py-lablib.devices.interface.camera.IGrabberAttributeCamera method), 843

clear\_acquisition() (py-lablib.devices.interface.camera.IROICamera method), 852

clear\_acquisition() (py-lablib.devices.PCO.SC2.PCOSC2Camera method), 648

clear\_acquisition() (py-lablib.devices.Photometrics.pvcam.PvcamCamera method), 665

clear\_acquisition() (py-lablib.devices.PhotonFocus.PhotonFocus.IPhotonFocusCamera method), 674

clear\_acquisition() (py-lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusBitFlowCamera method), 695

clear\_acquisition() (py-lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusIMAQCamera method), 678

clear\_acquisition() (py-lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusSiSoCamera method), 687

clear\_acquisition() (py-lablib.devices.PrincetonInstruments.picam.PicamCamera method), 714

clear\_acquisition() (py-lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareCamera method), 728

clear\_acquisition() (py-lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareFrameGrabber method), 723

clear\_acquisition() (py-lablib.devices.Thorlabs.TLCamera.ThorlabsTLCamera method), 771

clear\_acquisition() (py-lablib.devices.uc480.uc480.UC480Camera method), 868

clear\_all\_triggers() (py-lablib.devices.AlliedVision.Bonito.BonitoIMAQCamera method), 457

clear\_all\_triggers() (py-lablib.devices.IMAQ.IMAQ.IMAQCamera method), 569

clear\_all\_triggers() (py-lablib.devices.IMAQ.IMAQ.IMAQFrameGrabber method), 564

clear\_all\_triggers() (py-lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusIMAQCamera method), 679

clear\_limit\_error() (py-lablib.devices.Arcus.performax.Performax2EXStage method), 494

clear\_limit\_error() (py-lablib.devices.Arcus.performax.Performax4EXStage method), 491

clear\_limit\_error() (py-lablib.devices.Arcus.performax.PerformaxDMXJSASStage method), 498

clicked (pylablib.core.gui.widgets.label.EnumLabel attribute), 236

clicked (pylablib.core.gui.widgets.label.NumLabel attribute), 237

clicked (pylablib.core.gui.widgets.label.TextLabel attribute), 236

ClientSocket (class in pylablib.core.utils.net), 388

close() (pylablib.core.devio.comm\_backend.FT232DeviceBackend method), 150

close() (pylablib.core.devio.comm\_backend.ICommBackendWrapper method), 161

close() (pylablib.core.devio.comm\_backend.IDeviceCommBackend method), 143

close() (pylablib.core.devio.comm\_backend.NetworkDeviceBackend method), 153

close() (pylablib.core.devio.comm\_backend.PyUSBDeviceBackend method), 156

close() (pylablib.core.devio.comm\_backend.RecordedDeviceBackend method), 158

close() (pylablib.core.devio.comm\_backend.SerialDeviceBackend method), 148

close() (pylablib.core.devio.comm\_backend.VisaDeviceBackend method), 145

close() (pylablib.core.devio.interface.IDevice method), 162

close() (pylablib.core.devio.SCPISCPIDevice method), 162

method), 140

close() (pylablib.core.fileio.location.FolderFileSystemDataLocation method), 188

close() (pylablib.core.fileio.location.IDataLocation method), 185

close() (pylablib.core.fileio.location.IFileSystemDataLocation method), 185

close() (pylablib.core.fileio.location.LocationFile method), 184

close() (pylablib.core.fileio.location.OpenedFileLocation method), 185

close() (pylablib.core.fileio.location.PrefixedFileSystemDataLocation method), 187

close() (pylablib.core.fileio.location.SingleFileSystemDataLocation method), 186

close() (pylablib.core.utils.net.ClientSocket method), 390

close() (pylablib.devices.AlliedVision.Bonito.BonitoIMACamera method), 457

close() (pylablib.devices.AlliedVision.Bonito.IBonitoCamera method), 452

close() (pylablib.devices.Andor.AndorSDK2.AndorSDK2Camera method), 467

close() (pylablib.devices.Andor.AndorSDK2.LibraryController method), 465

close() (pylablib.devices.Andor.AndorSDK3.AndorSDK3Camera method), 479

close() (pylablib.devices.Andor.AndorSDK3.LibraryController method), 476

close() (pylablib.devices.Arcus.performax.GenericPerformanceStage method), 488

close() (pylablib.devices.Arcus.performax.Performax2EXStage method), 494

close() (pylablib.devices.Arcus.performax.Performax4EXStage method), 492

close() (pylablib.devices.Arcus.performax.PerformaxDMXISStage method), 499

close() (pylablib.devices.Arduino.base.IArduinoDevice method), 501

close() (pylablib.devices.Attocube.anc300.ANC300 method), 505

close() (pylablib.devices.Attocube.anc350.ANC350 method), 509

close() (pylablib.devices.AWG.generic.GenericAWG method), 406

close() (pylablib.devices.AWG.specific.Agilent33220A method), 415

close() (pylablib.devices.AWG.specific.Agilent33500 method), 409

close() (pylablib.devices.AWG.specific.InstekAFG2000 method), 427

close() (pylablib.devices.AWG.specific.InstekAFG2225 method), 421

close() (pylablib.devices.AWG.specific.RigolDG1000 method), 444

close() (pylablib.devices.AWG.specific.RSInstekAFG21000 method), 433

close() (pylablib.devices.AWG.specific.TektronixAFG1000 method), 439

close() (pylablib.devices.Basler.pylon.BaslerPylonCamera method), 514

close() (pylablib.devices.Basler.pylon.LibraryController method), 510

close() (pylablib.devices.BitFlow.BitFlow.BitFlowCamera method), 527

close() (pylablib.devices.BitFlow.BitFlow.BitFlowFrameGrabber method), 521

close() (pylablib.devices.Conrad.base.RelayBoard method), 534

close() (pylablib.devices.Cryocon.base.Cryocon1x method), 537

close() (pylablib.devices.Cryomagnetics.base.LM500 method), 540

close() (pylablib.devices.Cryomagnetics.base.LM510 method), 544

close() (pylablib.devices.DCAM.DCAM.DCAMCamera method), 550

close() (pylablib.devices.DCAM.DCAM.LibraryController method), 548

close() (pylablib.devices.HighFinesse.wlm.WLM method), 557

close() (pylablib.devices.IMAQ.IMAQ.IMAQCamera method), 569

close() (pylablib.devices.IMAQ.IMAQ.IMAQFrameGrabber method), 561

close() (pylablib.devices.IMAQdx.IMAQdx.EthernetIMAQdxCamera method), 584

close() (pylablib.devices.IMAQdx.IMAQdx.IMAQdxCamera method), 578

close() (pylablib.devices.interface.camera.IAttributeCamera method), 838

close() (pylablib.devices.interface.camera.IBinROICamera method), 857

close() (pylablib.devices.interface.camera.ICamera method), 834

close() (pylablib.devices.interface.camera.IExposureCamera method), 847

close() (pylablib.devices.interface.camera.IGrabberAttributeCamera method), 843

close() (pylablib.devices.interface.camera.IROICamera method), 852

close() (pylablib.devices.interface.stage.IMultiaxisStage method), 863

close() (pylablib.devices.interface.stage.IStage method), 862

close() (pylablib.devices.KJL.base.KJL300 method), 591

close() (pylablib.devices.Lakeshore.base.Lakeshore218 method), 444



- method), 595
- close() (pylablib.devices.Lakeshore.base.Lakeshore370 method), 600
- close() (pylablib.devices.LaserQuantum.base.Finesse method), 605
- close() (pylablib.devices.Leybold.base.GenericITR method), 608
- close() (pylablib.devices.Leybold.base.ITR90 method), 609
- close() (pylablib.devices.LighthousePhotonics.base.SproutG method), 613
- close() (pylablib.devices.M2.solstis.Solstis method), 614
- close() (pylablib.devices.Newport.picomotor.Picomotor8742 method), 630
- close() (pylablib.devices.NI.daq.NIDAQ method), 621
- close() (pylablib.devices.Ophir.base.OphirDevice method), 639
- close() (pylablib.devices.Ophir.base.VegaPowerMeter method), 643
- close() (pylablib.devices.OZOptics.base.DD100 method), 635
- close() (pylablib.devices.OZOptics.base.EPC04 method), 637
- close() (pylablib.devices.OZOptics.base.OZOpticsDevice method), 632
- close() (pylablib.devices.OZOptics.base.TF100 method), 634
- close() (pylablib.devices.PCO.SC2.PCOSC2Camera method), 645
- close() (pylablib.devices.Pfeiffer.base.DPG202 method), 658
- close() (pylablib.devices.Pfeiffer.base.TPG260 method), 656
- close() (pylablib.devices.Photometrics.pvcam.LibraryController method), 659
- close() (pylablib.devices.Photometrics.pvcam.PvcamCamera method), 662
- close() (pylablib.devices.PhotonFocus.PhotonFocus.IPhotonFocus method), 672
- close() (pylablib.devices.PhotonFocus.PhotonFocus.LibraryController method), 669
- close() (pylablib.devices.PhotonFocus.PhotonFocus.PhotonFocusBitFlyerG6 device method), 695
- close() (pylablib.devices.PhotonFocus.PhotonFocus.PhotonFocusIMAQ-Compact method), 679
- close() (pylablib.devices.PhotonFocus.PhotonFocus.PhotonFocusSiSo method), 687
- close() (pylablib.devices.PhysikInstrumente.base.GenericPIC method), 704
- close() (pylablib.devices.PhysikInstrumente.base.PIE516 method), 706
- close() (pylablib.devices.PrincetonInstruments.picam.LibraryController method), 708
- close() (pylablib.devices.PrincetonInstruments.picam.PicamCamera method), 712
- close() (pylablib.devices.SiliconSoftware.fgrab.SiliconSoftwareCamera method), 728
- close() (pylablib.devices.SiliconSoftware.fgrab.SiliconSoftwareFrameGrabber method), 721
- close() (pylablib.devices.Sirah.Matisse.SirahMatisse method), 740
- close() (pylablib.devices.SmarAct.scu3d.LibraryController method), 745
- close() (pylablib.devices.SmarAct.scu3d.SCU3D method), 746
- close() (pylablib.devices.Tektronix.base.DPO2000 method), 762
- close() (pylablib.devices.Tektronix.base.ITektronixScope method), 753
- close() (pylablib.devices.Tektronix.base.TDS2000 method), 756
- close() (pylablib.devices.Thorlabs.elliptec.ElliptecMotor method), 780
- close() (pylablib.devices.Thorlabs.kinesis.BasicKinesisDevice method), 783
- close() (pylablib.devices.Thorlabs.kinesis.KinesisDevice method), 788
- close() (pylablib.devices.Thorlabs.kinesis.KinesisMotor method), 799
- close() (pylablib.devices.Thorlabs.kinesis.KinesisPiezoMotor method), 803
- close() (pylablib.devices.Thorlabs.kinesis.KinesisQuadDetector method), 808
- close() (pylablib.devices.Thorlabs.kinesis.MFF method), 792
- close() (pylablib.devices.Thorlabs.serial.FW method), 813
- close() (pylablib.devices.Thorlabs.serial.FWv1 method), 817
- close() (pylablib.devices.Thorlabs.serial.MDT69xA method), 820
- close() (pylablib.devices.Thorlabs.serial.ThorlabsSerialInterface method), 810
- close() (pylablib.devices.Thorlabs.TLCamera.LibraryController method), 769
- close() (pylablib.devices.Thorlabs.TLCamera.ThorlabsTLCamera method), 770
- close() (pylablib.devices.Toptica.ibeam.TopticaIBeam method), 825
- close() (pylablib.devices.Trinamic.base.TMCM1110 method), 829
- close() (pylablib.devices.uc480.uc480.UC480Camera method), 867
- close() (pylablib.devices.utils.load\_lib.LibraryController method), 874
- close() (pylablib.devices.utils.load\_lib.Connection method), 874
- close() (pylablib.core.utils.ipc.SharedMemIPCTable method), 874

method), 385

close\_result (pylablib.devices.utils.load\_lib.TLibraryCloseResult attribute), 874

cls (pylablib.core.utils.string.TConversionClass attribute), 398

cmp\_dirs() (in module pylablib.core.utils.files), 366

cmp\_package\_version() (in module pylablib.core.utils.module), 385

cmp\_versions() (in module pylablib.core.utils.module), 385

coarse\_tune\_wavelength() (pylablib.devices.M2.solstis.Solstis method), 616

coeff (pylablib.devices.Lakeshore.base.TLakeshore218CurrentHeader attribute), 593

collect() (pylablib.core.utils.dictionary.Dictionary method), 331

collect() (pylablib.core.utils.dictionary.DictionaryPointer method), 337

collect() (pylablib.core.utils.dictionary.FilterTree method), 354

collect() (pylablib.core.utils.dictionary.PrefixTree method), 346

collect\_into\_bins() (in module pylablib.core.dataproc.filters), 113

column() (pylablib.core.dataproc.table\_wrap.Array2DWrapper method), 131

column() (pylablib.core.dataproc.table\_wrap.DataFrame2DWrapper method), 133

column() (pylablib.core.dataproc.table\_wrap.I2DWrapper method), 129

columns\_replaced() (pylablib.core.dataproc.table\_wrap.Array2DWrapper method), 131

columns\_replaced() (pylablib.core.dataproc.table\_wrap.DataFrame2DWrapper method), 133

columns\_replaced() (pylablib.core.dataproc.table\_wrap.I2DWrapper method), 129

columns\_to\_table() (in module pylablib.core.fileio.parse\_csv), 189

combine\_conn() (pylablib.core.devio.comm\_backend.FT232DeviceBackend class method), 151

combine\_conn() (pylablib.core.devio.comm\_backend.IDeviceCommBackend class method), 143

combine\_conn() (pylablib.core.devio.comm\_backend.NetworkDeviceBackend class method), 154

combine\_conn() (pylablib.core.devio.comm\_backend.PyUSBDeviceBackend class method), 157

combine\_conn() (pylablib.core.devio.comm\_backend.RecordedDeviceBackend class method), 158

combine\_conn() (pylablib.core.devio.comm\_backend.SerialDeviceBackend class method), 149

combine\_conn() (pylablib.core.devio.comm\_backend.VisaDeviceBackend class method), 146

combine\_dictionaries() (in module pylablib.core.utils.dictionary), 343

combine\_diff() (in module pylablib.core.utils.files), 366

ComboBoxHeaderParameterClass (class in pylablib.core.devio.interface), 167

ComboBox (class in pylablib.core.gui.widgets.combo\_box), 198

ComboBoxValueHandler (class in pylablib.core.gui.value\_handling), 274

comm (pylablib.devices.Conrad.base.RelayBoard.TMessage attribute), 533

comm (pylablib.devices.Thorlabs.elliptec.ElliptecMotor.CommData attribute), 778

comm (pylablib.devices.Trinamic.base.TMCM1110.ReplyData attribute), 827

comm() (pylablib.devices.Arduino.base.IArduinoDevice method), 501

comm() (pylablib.devices.KJL.base.KJL300 method), 590

comm() (pylablib.devices.Pfeiffer.base.DPG202 method), 657

comm() (pylablib.devices.Pfeiffer.base.TPG260 method), 654

comm\_paused() (pylablib.core.thread.controller.QTaskThread method), 307

command (pylablib.core.thread.controller.QTaskThread.TCommand attribute), 302

common (pylablib.core.utils.dictionary.DictionaryIntersection attribute), 335

compare\_lists() (in module pylablib.core.utils.general), 376

compilation\_number (pylablib.devices.HighFinesse.wlm.TDeviceInfo attribute), 557

complex\_lorentzian\_k() (in module pylablib.core.dataproc.specfunc), 124

compress\_grid\_layout() (in module pylablib.core.gui.utils), 264

configuration (pylablib.devices.LighthousePhotonics.base.TDeviceInfo attribute), 611

figure\_trigger\_in() (pylablib.devices.AlliedVision.Bonito.BonitoIMAQCamera

method), 457

configure\_trigger\_in() (py-lablib.devices.IMAQ.IMAQ.IMAQCamera method), 569

configure\_trigger\_in() (py-lablib.devices.IMAQ.IMAQ.IMAQFrameGrabber method), 563

configure\_trigger\_in() (py-lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusIMAQCamera method), 679

configure\_trigger\_out() (py-lablib.devices.AlliedVision.Bonito.BonitoIMAQCamera method), 457

configure\_trigger\_out() (py-lablib.devices.IMAQ.IMAQ.IMAQCamera method), 569

configure\_trigger\_out() (py-lablib.devices.IMAQ.IMAQ.IMAQFrameGrabber method), 563

configure\_trigger\_out() (py-lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusIMAQCamera method), 679

connect() (pylablib.core.utils.net.ClientSocket method), 389

connect\_device\_service() (in module py-lablib.core.utils.rpyc\_utils), 395

connect\_value\_changed\_handler() (py-lablib.core.gui.value\_handling.CheckboxValueHandler method), 272

connect\_value\_changed\_handler() (py-lablib.core.gui.value\_handling.ComboBoxValueHandler method), 275

connect\_value\_changed\_handler() (py-lablib.core.gui.value\_handling.IBoolValueHandler method), 271

connect\_value\_changed\_handler() (py-lablib.core.gui.value\_handling.ISingleValueHandler method), 269

connect\_value\_changed\_handler() (py-lablib.core.gui.value\_handling.IValueHandler method), 266

connect\_value\_changed\_handler() (py-lablib.core.gui.value\_handling.LabelValueHandler method), 270

connect\_value\_changed\_handler() (py-lablib.core.gui.value\_handling.LineEditValueHandler method), 270

connect\_value\_changed\_handler() (py-lablib.core.gui.value\_handling.ProgressBarValueHandler method), 275

connect\_value\_changed\_handler() (py-lablib.core.gui.value\_handling.PropertyValueHandler method), 267

connect\_value\_changed\_handler() (py-lablib.core.gui.value\_handling.PushButtonValueHandler method), 273

connect\_value\_changed\_handler() (py-lablib.core.gui.value\_handling.StandardValueHandler method), 268

connect\_value\_changed\_handler() (py-lablib.core.gui.value\_handling.ToolButtonValueHandler method), 274

connect\_value\_changed\_handler() (py-lablib.core.gui.value\_handling.VirtualValueHandler method), 266

connect\_wavemeter() (py-lablib.devices.M2.solstis.Solstis method), 615

ConradBackendError, 533

ConradError, 533

cons\_error (pylablib.devices.PrincetonInstruments.picam.PicamAttribute attribute), 710

cons\_excluded (py-lablib.devices.PrincetonInstruments.picam.PicamAttribute attribute), 711

cons\_included (py-lablib.devices.PrincetonInstruments.picam.PicamAttribute attribute), 711

cons\_novalid (pylablib.devices.PrincetonInstruments.picam.PicamAttribute attribute), 710

cons\_permanent (py-lablib.devices.PrincetonInstruments.picam.PicamAttribute attribute), 710

cons\_type (pylablib.devices.PrincetonInstruments.picam.PicamAttribute attribute), 710

constant() (in module pylablib.core.utils.numerical), 392

contained\_value\_changed (py-lablib.core.gui.widgets.container.IQContainer attribute), 200

contained\_value\_changed (py-lablib.core.gui.widgets.container.IQWidgetContainer attribute), 207

contained\_value\_changed (py-lablib.core.gui.widgets.container.QContainer attribute), 203

contained\_value\_changed (py-lablib.core.gui.widgets.container.QDialogContainer attribute), 218

contained\_value\_changed (py-lablib.core.gui.widgets.container.QFrameContainer attribute), 214

contained\_value\_changed (py-lablib.core.gui.widgets.container.QGroupBoxContainer attribute), 222

contained\_value\_changed (py-lablib.core.gui.widgets.container.QScrollAreaContainer attribute), 230

contained\_value\_changed (py- method), 151  
 lablib.core.gui.widgets.container.QScrollAreaContainer.QWidgetContainer (py- method), 143  
 attribute), 226  
 contained\_value\_changed (py- cooldown () (pylablib.core.devio.comm\_backend.IDeviceCommBackend method), 143  
 lablib.core.gui.widgets.container.QTabContainer method), 154  
 attribute), 232  
 contained\_value\_changed (py- cooldown () (pylablib.core.devio.comm\_backend.PyUSBDeviceBackend method), 157  
 lablib.core.gui.widgets.container.QWidgetContainer cooldown () (pylablib.core.devio.comm\_backend.RecordedDeviceBackend method), 159  
 attribute), 211  
 contained\_value\_changed (py- cooldown () (pylablib.core.devio.comm\_backend.SerialDeviceBackend method), 149  
 lablib.core.gui.widgets.param\_table.ParamTable method), 149  
 attribute), 250  
 contained\_value\_changed (py- cooldown () (pylablib.core.devio.comm\_backend.VisaDeviceBackend method), 146  
 lablib.core.gui.widgets.param\_table.StatusTable copy () (pylablib.core.dataproc.image.ROI method),  
 attribute), 258 121  
 contains () (pylablib.core.dataproc.utils.Range copy () (pylablib.core.dataproc.table\_wrap.Array1DWrapper method), 136  
 method), 136  
 continuous (pylablib.devices.NI.daq.TVoltageOutputClockPayoff (pylablib.core.dataproc.table\_wrap.Array2DWrapper method), 131  
 attribute), 621  
 controller (pylablib.devices.LighthousePhotonics.base.WorkHorse (pylablib.core.dataproc.table\_wrap.DataFrame2DWrapper method), 133  
 attribute), 611  
 controller\_model (py- copy () (pylablib.core.dataproc.table\_wrap.I1DWrapper method), 126  
 lablib.devices.Andor.AndorSDK2.TDeviceInfo copy () (pylablib.core.dataproc.table\_wrap.I2DWrapper method), 129  
 attribute), 465  
 conv (pylablib.core.utils.string.TConversionClass at- copy () (pylablib.core.dataproc.table\_wrap.IGenWrapper method), 125  
 tribute), 398  
 conv () (pylablib.core.utils.ctypes\_wrap.CStructWrapper copy () (pylablib.core.dataproc.table\_wrap.Series1DWrapper method), 128  
 method), 325  
 convert\_columns () (py- copy () (pylablib.core.fileio.location.LocationName method), 184  
 lablib.core.fileio.parse\_csv.ChunksAccumulator method), 189  
 convert\_frequency\_units () (in module py- copy () (pylablib.core.utils.dictionary.Dictionary method), 331  
 lablib.core.utils.units), 402  
 convert\_from\_str () (py- copy () (pylablib.core.utils.dictionary.DictionaryPointer method), 337  
 lablib.core.devio.data\_format.DataFormat method), 162  
 convert\_image\_indexing () (in module py- copy () (pylablib.core.utils.dictionary.FilterTree method), 352  
 lablib.core.dataproc.image), 121  
 convert\_length\_units () (in module py- copy () (pylablib.core.utils.dictionary.PrefixShortcutTree method), 360  
 lablib.core.utils.units), 401  
 convert\_power\_units () (in module py- copy () (pylablib.core.utils.dictionary.PrefixTree method), 344  
 lablib.core.utils.units), 402  
 convert\_shape\_indexing () (in module py- copy () (pylablib.core.utils.functions.FunctionSignature method), 371  
 lablib.core.dataproc.image), 121  
 convert\_time\_units () (in module py- copy\_array\_chunks (in module py-  
 lablib.core.utils.units), 401  
 convert\_to\_str () (py- lablib.core.utils.nbtools), 387  
 lablib.core.devio.data\_format.DataFormat copy\_array\_strided (in module py-  
 method), 162  
 lablib.core.utils.nbtools), 387  
 convolution\_filter () (in module py- copy\_dir () (in module pylablib.core.utils.files), 365  
 lablib.core.dataproc.filters), 110  
 copy\_file () (in module pylablib.core.utils.files), 363  
 convolve1d () (in module py- corr\_number (pylablib.devices.Attocube.anc350.ANC350.Telegram  
 lablib.core.dataproc.filters), 110  
 attribute), 506  
 corrupted\_number () (py-  
 lablib.core.fileio.parse\_csv.ChunksAccumulator  
 method), 189



- count () (pylablib.core.dataproc.feature.Baseline method), 108
- count () (pylablib.core.dataproc.feature.Peak method), 108
- count () (pylablib.core.devio.interface.TRawParameterValues method), 168
- count () (pylablib.core.gui.utils.TWidgetLocation method), 263
- count () (pylablib.core.gui.value\_handling.GUIValues.IndicatorSet method), 279
- count () (pylablib.core.gui.widgets.container.TChild method), 200
- count () (pylablib.core.gui.widgets.container.TTimer method), 199
- count () (pylablib.core.gui.widgets.container.TTimerEvent method), 199
- count () (pylablib.core.gui.widgets.param\_table.ParamTable.ParamTableRow method), 243
- count () (pylablib.core.gui.widgets.param\_table.StatusTable.ParamTableRow method), 252
- count () (pylablib.core.thread.callsync.QScheduledCall.Callback method), 283
- count () (pylablib.core.thread.callsync.TDefaultCallInfo method), 284
- count () (pylablib.core.thread.controller.QTaskThread.TBatchJob method), 302
- count () (pylablib.core.thread.controller.QTaskThread.TCommand method), 302
- count () (pylablib.core.thread.multicast\_pool.TMulticast method), 315
- count () (pylablib.core.utils.dictionary.DictionaryDiff method), 334
- count () (pylablib.core.utils.dictionary.DictionaryIntersection method), 335
- count () (pylablib.core.utils.files.FolderList method), 364
- count () (pylablib.core.utils.ipc.TPipeMsg method), 383
- count () (pylablib.core.utils.ipc.TShmemVarDesc method), 384
- count () (pylablib.core.utils.observer\_pool.ObserverPool.Observer method), 392
- count () (pylablib.core.utils.string.TConversionClass method), 398
- count () (pylablib.devices.AlliedVision.Bonito.TDeviceInfo method), 450
- count () (pylablib.devices.AlliedVision.Bonito.TStatusLine method), 464
- count () (pylablib.devices.Andor.AndorSDK2.TAcqProgress method), 466
- count () (pylablib.devices.Andor.AndorSDK2.TCycleTimings method), 466
- count () (pylablib.devices.Andor.AndorSDK2.TDeviceInfo method), 465
- count () (pylablib.devices.Andor.AndorSDK3.TDeviceInfo method), 478
- count () (pylablib.devices.Andor.AndorSDK3.TFrameInfo method), 479
- count () (pylablib.devices.Andor.AndorSDK3.TMissedFramesStatus method), 479
- count () (pylablib.devices.Attocube.anc300.TDeviceInfo method), 502
- count () (pylablib.devices.Attocube.anc350.ANC350.Reply method), 507
- count () (pylablib.devices.Attocube.anc350.ANC350.Telegram method), 506
- count () (pylablib.devices.Basler.pylon.TCameraInfo method), 511
- count () (pylablib.devices.Basler.pylon.TDeviceInfo method), 513
- count () (pylablib.devices.BitFlow.BitFlow.TDeviceInfo method), 521
- count () (pylablib.devices.Conrad.base.RelayBoard.TMessage method), 533
- count () (pylablib.devices.DCAM.DCAM.TDeviceInfo method), 550
- count () (pylablib.devices.DCAM.DCAM.TFrameInfo method), 550
- count () (pylablib.devices.HighFinesse.wlm.TDeviceInfo method), 557
- count () (pylablib.devices.IMAQ.IMAQ.TDeviceInfo method), 561
- count () (pylablib.devices.IMAQdx.IMAQdx.TCameraInfo method), 576
- count () (pylablib.devices.IMAQdx.IMAQdx.TDeviceInfo method), 578
- count () (pylablib.devices.interface.camera.TAcqTimings method), 846
- count () (pylablib.devices.interface.camera.TAxisROILimit method), 851
- count () (pylablib.devices.interface.camera.TFrameInfo method), 831
- count () (pylablib.devices.interface.camera.TFramePosition method), 830
- count () (pylablib.devices.interface.camera.TFrameSize method), 830
- count () (pylablib.devices.interface.camera.TFramesStatus method), 830
- count () (pylablib.devices.interface.camera.TStatusLineDescription method), 860
- count () (pylablib.devices.KJL.base.TKJL300DeviceInfo method), 590
- count () (pylablib.devices.Lakeshore.base.TLakeshore218AnalogSettings method), 592
- count () (pylablib.devices.Lakeshore.base.TLakeshore218CurveHeader method), 593
- count () (pylablib.devices.Lakeshore.base.TLakeshore218FilterSettings method), 593

- count () (*pylablib.devices.Lakeshore.base.TLakeshore370AnalogSettings* method), 598
- count () (*pylablib.devices.Lakeshore.base.TLakeshore370FilterSettings* method), 599
- count () (*pylablib.devices.Lakeshore.base.TLakeshore370RangeSettings* method), 598
- count () (*pylablib.devices.LaserQuantum.base.TDeviceInfo* method), 603
- count () (*pylablib.devices.LaserQuantum.base.TTemperature* method), 604
- count () (*pylablib.devices.LaserQuantum.base.TWorkHours* method), 604
- count () (*pylablib.devices.Leybold.base.TDeviceInfo* method), 607
- count () (*pylablib.devices.Leybold.base.TITR90Status* method), 609
- count () (*pylablib.devices.Leybold.base.TUpdateValue* method), 607
- count () (*pylablib.devices.LighthousePhotonics.base.TDeviceInfo* method), 611
- count () (*pylablib.devices.LighthousePhotonics.base.TWorkHours* method), 611
- count () (*pylablib.devices.Newport.picomotor.TDeviceInfo* method), 628
- count () (*pylablib.devices.NI.daq.TDeviceInfo* method), 620
- count () (*pylablib.devices.NI.daq.TVoltageOutputClockParameters* method), 621
- count () (*pylablib.devices.Ophir.base.TDeviceInfo* method), 640
- count () (*pylablib.devices.Ophir.base.THeadInfo* method), 640
- count () (*pylablib.devices.Ophir.base.TRangeInfo* method), 641
- count () (*pylablib.devices.Ophir.base.TWavelengthInfo* method), 640
- count () (*pylablib.devices.PCO.SC2.TCameraStatus* method), 644
- count () (*pylablib.devices.PCO.SC2.TDeviceInfo* method), 644
- count () (*pylablib.devices.PCO.SC2.TFrameInfo* method), 645
- count () (*pylablib.devices.PCO.SC2.TInternalBufferStatus* method), 645
- count () (*pylablib.devices.PCO.SC2.TStatusLine* method), 652
- count () (*pylablib.devices.Pfeiffer.base.TTPG260GaugeControls* method), 654
- count () (*pylablib.devices.Pfeiffer.base.TTPG260SwitchSettings* method), 654
- count () (*pylablib.devices.Photometrics.pvcam.TDeviceInfo* method), 661
- count () (*pylablib.devices.Photometrics.pvcam.TFrameInfo* method), 661
- count () (*pylablib.devices.Photometrics.pvcam.TReadoutInfo* method), 662
- count () (*pylablib.devices.PhotonFocus.PhotonFocus.TCameraInfo* method), 670
- count () (*pylablib.devices.PhotonFocus.PhotonFocus.TDeviceInfo* method), 671
- count () (*pylablib.devices.PrincetonInstruments.picam.TCameraInfo* method), 708
- count () (*pylablib.devices.PrincetonInstruments.picam.TDeviceInfo* method), 711
- count () (*pylablib.devices.PrincetonInstruments.picam.TFrameInfo* method), 712
- count () (*pylablib.devices.PrincetonInstruments.picam.TROIConstraints* method), 709
- count () (*pylablib.devices.SiliconSoftware.fgrab.TAppletInfo* method), 718
- count () (*pylablib.devices.SiliconSoftware.fgrab.TBoardInfo* method), 718
- count () (*pylablib.devices.SiliconSoftware.fgrab.TDeviceInfo* method), 720
- count () (*pylablib.devices.SiliconSoftware.fgrab.TFrameInfo* method), 720
- count () (*pylablib.devices.SiliconSoftware.fgrab.TFullAppletInfo* method), 718
- count () (*pylablib.devices.Sirah.Matisse.TFastpiezoCtlParameters* method), 735
- count () (*pylablib.devices.Sirah.Matisse.TPiezoetDriveParameters* method), 734
- count () (*pylablib.devices.Sirah.Matisse.TPiezoetFeedbackParameters* method), 734
- count () (*pylablib.devices.Sirah.Matisse.TPiezoetFeedforwardParameters* method), 734
- count () (*pylablib.devices.Sirah.Matisse.TRefcellWaveformParameters* method), 735
- count () (*pylablib.devices.Sirah.Matisse.TScanMode* method), 736
- count () (*pylablib.devices.Sirah.Matisse.TScanParameters* method), 736
- count () (*pylablib.devices.Sirah.Matisse.TSlowpiezoCtlParameters* method), 735
- count () (*pylablib.devices.Sirah.Matisse.TThinnetCtlParameters* method), 734
- count () (*pylablib.devices.SmarAct.scu3d.TDeviceInfo* method), 745
- count () (*pylablib.devices.Tektronix.base.TTriggerParameters* method), 748
- count () (*pylablib.devices.Thorlabs.elliptec.ElliptecMotor.CommData* method), 778
- count () (*pylablib.devices.Thorlabs.elliptec.TDeviceInfo* method), 777
- count () (*pylablib.devices.Thorlabs.elliptec.TMotorInfo* method), 777
- count () (*pylablib.devices.Thorlabs.kinesis.BasicKinesisDevice.CommData* method), 782

count () (pylablib.devices.Thorlabs.kinesis.BasicKinesisDevice.CommShort method), 782

count () (pylablib.devices.Thorlabs.kinesis.KinesisDevice.CommData (pylablib.devices.Toptica.ibeam.TDeviceInfo method), 787

count () (pylablib.devices.Thorlabs.kinesis.KinesisDevice.CommShort (pylablib.devices.Toptica.ibeam.TTemperatures method), 788

count () (pylablib.devices.Thorlabs.kinesis.KinesisMotor.CommData (pylablib.devices.Toptica.ibeam.TWorkHours method), 798

count () (pylablib.devices.Thorlabs.kinesis.KinesisMotor.CommShort (pylablib.devices.Trinamic.base.TMCM1110.ReplyData method), 798

count () (pylablib.devices.Thorlabs.kinesis.KinesisPiezoMotor.CommData (pylablib.devices.uc480.uc480.TAcquiredFramesStatus method), 802

count () (pylablib.devices.Thorlabs.kinesis.KinesisPiezoMotor.CommShort (pylablib.devices.uc480.uc480.TCameraInfo method), 802

count () (pylablib.devices.Thorlabs.kinesis.KinesisQuadDetector.CommData (pylablib.devices.uc480.uc480.TDeviceInfo method), 807

count () (pylablib.devices.Thorlabs.kinesis.KinesisQuadDetector.CommShort (pylablib.devices.uc480.uc480.TFrameInfo method), 807

count () (pylablib.devices.Thorlabs.kinesis.MFF.CommData count () (pylablib.devices.uc480.uc480.TTimestamp method), 791

count () (pylablib.devices.Thorlabs.kinesis.MFF.CommShort count () (pylablib.devices.utils.load\_lib.TLibraryCloseResult method), 792

count () (pylablib.devices.Thorlabs.kinesis.TDeviceInfo count () (pylablib.devices.utils.load\_lib.TLibraryOpenResult method), 781

count () (pylablib.devices.Thorlabs.kinesis.TFlipperParameters method), 790

count () (pylablib.devices.Thorlabs.kinesis.TGenMoveParams method), 785

count () (pylablib.devices.Thorlabs.kinesis.THomeParams method), 785

count () (pylablib.devices.Thorlabs.kinesis.TJogParams method), 784

count () (pylablib.devices.Thorlabs.kinesis.TKCubeTrigIOParams method), 786

count () (pylablib.devices.Thorlabs.kinesis.TKCubeTrigPosParams method), 786

count () (pylablib.devices.Thorlabs.kinesis.TLimitSwitchParams method), 785

count () (pylablib.devices.Thorlabs.kinesis.TPZMotorDriveParams method), 786

count () (pylablib.devices.Thorlabs.kinesis.TPZMotorJogParams method), 787

count () (pylablib.devices.Thorlabs.kinesis.TQuadDetectorOutputParameter method), 806

count () (pylablib.devices.Thorlabs.kinesis.TQuadDetectorPIDParameters method), 805

count () (pylablib.devices.Thorlabs.kinesis.TQuadDetectorReading attribute), 805

count () (pylablib.devices.Thorlabs.kinesis.TQuadDetectorSetpoint attribute), 805

count () (pylablib.devices.Thorlabs.kinesis.TVelocityParams method), 784

count () (pylablib.devices.Thorlabs.TLCamera.TDeviceInfo method), 769

countdown (class in pylablib.core.utils.general), 378

coupling (pylablib.devices.Tektronix.base.TTriggerParameters attribute), 748

covers\_all () (in module pylablib.core.utils.indexing), 381

create\_indicator\_handler () (in module pylablib.core.gui.value\_handling), 277

create\_value\_handler () (in module pylablib.core.gui.value\_handling), 276

Cryocon1x (class in pylablib.devices.Cryocon.base), 535

Cryocon1x.NoParameterCaller (class in pylablib.devices.Cryocon.base), 536

CryoconBackendError, 535

CryoconError, 535

CryomagneticsBackendError, 539

CryomagneticsError, 539

cs (pylablib.core.thread.controller.QTaskThread attribute), 301

csi (pylablib.core.thread.controller.QTaskThread attribute), 301

css (pylablib.core.thread.controller.QTaskThread attribute), 301

CStructWrapper (class in pylablib.core.utils.ctypes\_wrap), 324

CSVTableInputFileFormat (class in pylablib.core.fileio.loadfile), 177

CSVTableOutputFileFormat (class in pylablib.core.fileio.savefile), 191

curr\_idx (pylablib.devices.Ophir.base.TRangeInfo attribute), 301

tribute), 641

curr\_idx (pylablib.devices.Ophir.base.TWavelengthInfo attribute), 640

curr\_range (pylablib.devices.Ophir.base.TRangeInfo attribute), 641

curr\_wavelength (pylablib.devices.Ophir.base.TWavelengthInfo attribute), 640

current (pylablib.devices.Thorlabs.elliptec.TMotorInfo attribute), 777

current\_controller() (in module pylablib.core.thread.threadprop), 320

cut\_out\_regions() (in module pylablib.core.dataproc.utils), 137

cut\_to\_range() (in module pylablib.core.dataproc.utils), 137

cycles\_done (pylablib.devices.Andor.AndorSDK2.AcqProgress attribute), 466

**D**

d (pylablib.devices.Thorlabs.kinesis.TQuadDetectorPIDParams attribute), 805

data (pylablib.core.utils.ipc.TPipeMsg attribute), 383

data (pylablib.devices.Attocube.anc350.ANC350.Reply attribute), 507

data (pylablib.devices.Attocube.anc350.ANC350.Telegram attribute), 506

data (pylablib.devices.Conrad.base.RelayBoard.TMessage attribute), 533

data (pylablib.devices.Thorlabs.elliptec.ElliptecMotor.CommandData attribute), 778

data (pylablib.devices.Thorlabs.kinesis.BasicKinesisDevice.CommandData attribute), 782

data (pylablib.devices.Thorlabs.kinesis.KinesisDevice.CommandData attribute), 787

data (pylablib.devices.Thorlabs.kinesis.KinesisMotor.CommandData attribute), 798

data (pylablib.devices.Thorlabs.kinesis.KinesisPiezoMotor.CommandData attribute), 802

data (pylablib.devices.Thorlabs.kinesis.KinesisQuadDetector.CommandData attribute), 807

data (pylablib.devices.Thorlabs.kinesis.MFF.CommData attribute), 791

DataFile (class in pylablib.core.fileio.datafile), 168

DataFormat (class in pylablib.core.devio.data\_format), 161

DataFrame2DWrapper (class in pylablib.core.dataproc.table\_wrap), 131

DataFrame2DWrapper.ColumnAccessor (class in pylablib.core.dataproc.table\_wrap), 132

DataFrame2DWrapper.RowAccessor (class in pylablib.core.dataproc.table\_wrap), 131

DataFrame2DWrapper.TableAccessor (class in pylablib.core.dataproc.table\_wrap), 133

date (pylablib.devices.uc480.uc480.TDeviceInfo attribute), 865

day (pylablib.devices.uc480.uc480.TTimestamp attribute), 865

DCAMAttribute (class in pylablib.devices.DCAM.DCAM), 548

DCAMCamera (class in pylablib.devices.DCAM.DCAM), 550

DCAMCamera.NoParameterCaller (class in pylablib.devices.DCAM.DCAM), 553

DD100 (class in pylablib.devices.OZOptics.base), 635

DD100.NoParameterCaller (class in pylablib.devices.OZOptics.base), 635

deactivation\_control (pylablib.devices.Pfeiffer.base.TTPG260GaugeControlSettings attribute), 654

deallocate\_buffers() (pylablib.devices.interface.camera.ChunkBufferManager method), 837

deallocate\_buffers() (pylablib.devices.Andor.AndorSDK3.AndorSDK3Camera.BufferManager method), 481

decimate() (in module pylablib.core.dataproc.filters), 112

decimate\_datasets() (in module pylablib.core.dataproc.filters), 112

decimate\_full() (in module pylablib.core.dataproc.filters), 112

decilen\_bo (pylablib.core.utils.net.ClientSocket attribute), 389

decilen\_ll (pylablib.core.utils.net.ClientSocket attribute), 389

default (pylablib.devices.Photometrics.pvcam.PvcamAttribute attribute), 660

default (pylablib.devices.PrincetonInstruments.picam.PicamAttribute attribute), 710

DefaultFrameTransferError, 830

del\_entry() (pylablib.core.utils.dictionary.Dictionary method), 327

del\_entry() (pylablib.core.utils.dictionary.DictionaryPointer method), 337

del\_entry() (pylablib.core.utils.dictionary.FilterTree method), 354

del\_entry() (pylablib.core.utils.dictionary.PrefixTree method), 346

del\_fit\_parameters() (pylablib.core.dataproc.fitting.Fitter method), 115

del\_fixed\_parameters() (pylablib.core.dataproc.fitting.Fitter method), 115

delattr\_call() (in module pylablib.core.utils.functions), 372

delaydef() (in module pylablib.core.utils.functions), 374



`delete_layout_item()` (in module `pylablib.core.gui.utils`), 263  
`delete_thread_method()` (`pylablib.core.thread.controller.QTaskThread` method), 309  
`delete_thread_method()` (`pylablib.core.thread.controller.QThreadController` method), 298  
`delete_variable()` (`pylablib.core.thread.controller.QTaskThread` method), 309  
`delete_variable()` (`pylablib.core.thread.controller.QThreadController` method), 298  
`delete_widget()` (in module `pylablib.core.gui.utils`), 263  
`deregister()` (`pylablib.devices.Basler.pylon.BaslerPylonCamera.BufferManager` method), 515  
`desc` (`pylablib.devices.SiliconSoftware.fgrab.TFullAppletInfo` attribute), 719  
`description` (`pylablib.devices.Basler.pylon.BaslerPylonAttribute` attribute), 512  
`description` (`pylablib.devices.IMAQdx.IMAQdx.IMAQdxAttribute` attribute), 577  
`dest` (`pylablib.devices.Thorlabs.kinesis.BasicKinesisDevice.CommData` attribute), 782  
`dest` (`pylablib.devices.Thorlabs.kinesis.BasicKinesisDevice.CommShort` attribute), 782  
`dest` (`pylablib.devices.Thorlabs.kinesis.KinesisDevice.CommData` attribute), 788  
`dest` (`pylablib.devices.Thorlabs.kinesis.KinesisDevice.CommShort` attribute), 788  
`dest` (`pylablib.devices.Thorlabs.kinesis.KinesisMotor.CommData` attribute), 798  
`dest` (`pylablib.devices.Thorlabs.kinesis.KinesisMotor.CommShort` attribute), 798  
`dest` (`pylablib.devices.Thorlabs.kinesis.KinesisPiezoMotor.CommData` attribute), 802  
`dest` (`pylablib.devices.Thorlabs.kinesis.KinesisPiezoMotor.CommShort` attribute), 802  
`dest` (`pylablib.devices.Thorlabs.kinesis.KinesisQuadDetector.CommData` attribute), 807  
`dest` (`pylablib.devices.Thorlabs.kinesis.KinesisQuadDetector.CommShort` attribute), 807  
`dest` (`pylablib.devices.Thorlabs.kinesis.MFF.CommData` attribute), 791  
`dest` (`pylablib.devices.Thorlabs.kinesis.MFF.CommShort` attribute), 792  
`detach()` (`pylablib.core.utils.dictionary.Dictionary` method), 331  
`detach()` (`pylablib.core.utils.dictionary.DictionaryPointer` method), 337  
`detach()` (`pylablib.core.utils.dictionary.FilterTree` method), 354  
`detach()` (`pylablib.core.utils.dictionary.PrefixTree` method), 346  
`detect_binary_file()` (in module `pylablib.core.fileio.loadfile_utils`), 182  
`detect_file_format()` (`pylablib.core.fileio.loadfile.BinaryTableInputFileFormatter` static method), 179  
`detect_file_format()` (`pylablib.core.fileio.loadfile.CSVTableInputFileFormat` static method), 177  
`detect_file_format()` (`pylablib.core.fileio.loadfile.DictionaryInputFileFormat` static method), 178  
`detect_file_format()` (`pylablib.core.fileio.loadfile.IInputFileFormat` static method), 176  
`detect_file_format()` (`pylablib.core.fileio.loadfile.ITextInputFileFormat` static method), 177  
`detect_textfile_type()` (in module `pylablib.core.fileio.loadfile_utils`), 182  
`dev_id` (`pylablib.devices.uc480.uc480.TCameraInfo` attribute), 864  
`devclass` (`pylablib.devices.Basler.pylon.TCameraInfo` attribute), 511  
`devclass` (`pylablib.devices.Basler.pylon.TDeviceInfo` attribute), 513  
`device` (`pylablib.devices.Sirah.Matisse.TScanParameters` attribute), 736  
`device_id` (`pylablib.devices.SmarAct.scu3d.TDeviceInfo` attribute), 745  
`device_info` (`pylablib.devices.Leybold.base.TUpdateValue` attribute), 607  
`DeviceBackendError`, 142  
`DeviceNotFoundError`, 142  
`DeviceFT232Error`, 149  
`DeviceNetworkError`, 152  
`DeviceRecordedError`, 157  
`DeviceSerialError`, 147  
`DeviceService` (class in `pylablib.core.utils.rpyc_utils`), 394  
`DeviceUSBError`, 155  
`DeviceWriteError`, 144  
`devversion` (`pylablib.devices.Basler.pylon.TCameraInfo` attribute), 511  
`devversion` (`pylablib.devices.Basler.pylon.TDeviceInfo` attribute), 513  
`dict_to_object_local()` (in module `pylablib.core.utils.dictionary`), 360  
`DictEntryBuilder` (class in `pylablib.core.fileio.dict_entry`), 169  
`DictEntryParser` (class in `pylablib.core.fileio.dict_entry`), 169  
`Dictionary` (class in `pylablib.core.utils.dictionary`),

- 326
- DictionaryDiff (class in py-lablib.core.utils.dictionary), 334
- DictionaryInputFileFormat (class in py-lablib.core.fileio.loadfile), 177
- DictionaryIntersection (class in py-lablib.core.utils.dictionary), 334
- DictionaryNode (class in py-lablib.core.utils.dictionary), 360
- DictionaryOutputFileFormat (class in py-lablib.core.fileio.savefile), 192
- DictionaryPointer (class in py-lablib.core.utils.dictionary), 335
- diff() (pylablib.core.utils.dictionary.Dictionary method), 333
- diff() (pylablib.core.utils.dictionary.DictionaryPointer method), 337
- diff() (pylablib.core.utils.dictionary.FilterTree method), 354
- diff() (pylablib.core.utils.dictionary.PrefixTree method), 346
- diff\_flatdict() (py-lablib.core.utils.dictionary.Dictionary static method), 333
- diff\_flatdict() (py-lablib.core.utils.dictionary.DictionaryPointer static method), 337
- diff\_flatdict() (py-lablib.core.utils.dictionary.FilterTree static method), 354
- diff\_flatdict() (py-lablib.core.utils.dictionary.PrefixTree static method), 346
- differentiate() (in module py-lablib.core.dataproc.filters), 111
- diode (pylablib.devices.Toptica.ibeam.TTemperatures attribute), 824
- dir\_empty() (in module pylablib.core.utils.files), 364
- disable\_axis() (py-lablib.devices.Attocube.anc300.ANC300 method), 503
- disable\_axis() (py-lablib.devices.Attocube.anc350.ANC350 method), 507
- disconnect\_wavemeter() (py-lablib.devices.M2.solstis.Solstis method), 615
- display\_name (pylablib.devices.Basler.pylon.BaslerPylonAttribute attribute), 512
- display\_name (pylablib.devices.IMAQdx.IMAQdx.IMAQdxAttribute attribute), 577
- display\_units (py-lablib.devices.Leybold.base.TUpdateValue attribute), 607
- dll\_version (pylablib.devices.SmarAct.scu3d.TDeviceInfo attribute), 745
- dll\_version (pylablib.devices.uc480.uc480.TDeviceInfo attribute), 865
- doc\_inherit() (in module py-lablib.core.utils.general), 378
- docstring() (pylablib.core.devio.interface.CombinedParameterClass method), 167
- docstring() (pylablib.core.devio.interface.EnumParameterClass method), 166
- docstring() (pylablib.core.devio.interface.FunctionParameterClass method), 167
- docstring() (pylablib.core.devio.interface.ICheckingParameterClass method), 164
- docstring() (pylablib.core.devio.interface.IEnumParameterClass method), 165
- docstring() (pylablib.core.devio.interface.IParameterClass method), 163
- docstring() (pylablib.core.devio.interface.RangeParameterClass method), 164
- done\_notify() (py-lablib.core.thread.callsync.QCallResultSynchronizer method), 281
- done\_notify() (py-lablib.core.thread.callsync.QDirectResultSynchronizer method), 283
- done\_notify() (py-lablib.core.thread.notifier.ISkippableNotifier method), 317
- done\_notify() (py-lablib.core.thread.synchronizing.QThreadNotifier method), 318
- done\_wait() (pylablib.core.thread.callsync.QCallResultSynchronizer method), 282
- done\_wait() (pylablib.core.thread.callsync.QDirectResultSynchronizer method), 283
- done\_wait() (pylablib.core.thread.notifier.ISkippableNotifier method), 317
- done\_wait() (pylablib.core.thread.synchronizing.QThreadNotifier method), 318
- DPG202 (class in pylablib.devices.Pfeiffer.base), 657
- DPG202.NoParameterCaller (class in py-lablib.devices.Pfeiffer.base), 658
- DPO2000 (class in pylablib.devices.Tektronix.base), 762
- DPO2000.NoParameterCaller (class in py-lablib.devices.Tektronix.base), 762
- DummyResource (class in pylablib.core.utils.general), 378
- dump() (in module pylablib.core.utils.strdump), 396
- dump() (pylablib.core.utils.strdump.StrDumper method), 396
- dumper (in module pylablib.core.utils.strdump), 396
- dumps() (in module pylablib.core.utils.strdump), 396
- dumps() (pylablib.core.utils.strdump.StrDumper method), 396

*method*), 396  
 DuplicateControllerThreadError, 319

## E

ElliptecMotor (class in py-  
 lablib.devices.Thorlabs.elliptec), 777  
 ElliptecMotor.CommData (class in py-  
 lablib.devices.Thorlabs.elliptec), 778  
 ElliptecMotor.NoParameterCaller (class in  
 pylablib.devices.Thorlabs.elliptec), 780  
 emission (pylablib.devices.Leybold.base.TITR90Status  
 attribute), 609  
 empty\_object\_property() (in module py-  
 lablib.core.utils.functions), 373  
 enable (pylablib.devices.Lakeshore.base.TLakeshore370RangeSetting  
 attribute), 598  
 enable() (pylablib.devices.LaserQuantum.base.Finesse  
 method), 605  
 enable() (pylablib.devices.LighthousePhotonics.base.Sprangle  
 method), 612  
 enable() (pylablib.devices.Pfeiffer.base.TPG260  
 method), 655  
 enable() (pylablib.devices.Toptica.ibeam.TopticaIBeam  
 method), 824  
 enable\_absolute\_mode() (py-  
 lablib.devices.Arcus.performax.Performax2EXStage  
 method), 494  
 enable\_absolute\_mode() (py-  
 lablib.devices.Arcus.performax.Performax4EXStage  
 method), 490  
 enable\_absolute\_mode() (py-  
 lablib.devices.Arcus.performax.PerformaxDMXJSAStage  
 method), 497  
 enable\_axis() (py-  
 lablib.devices.Arcus.performax.Performax2EXStage  
 method), 494  
 enable\_axis() (py-  
 lablib.devices.Arcus.performax.Performax4EXStage  
 method), 490  
 enable\_axis() (py-  
 lablib.devices.Arcus.performax.PerformaxDMXJSAStage  
 method), 497  
 enable\_axis() (py-  
 lablib.devices.Attocube.anc300.ANC300  
 method), 503  
 enable\_axis() (py-  
 lablib.devices.Attocube.anc350.ANC350  
 method), 507  
 enable\_burst() (py-  
 lablib.devices.AWG.generic.GenericAWG  
 method), 405  
 enable\_burst() (py-  
 lablib.devices.AWG.specific.Agilent33220A  
 method), 415  
 enable\_burst() (py-  
 lablib.devices.AWG.specific.Agilent33500  
 method), 409  
 enable\_burst() (py-  
 lablib.devices.AWG.specific.InstekAFG2000  
 method), 427  
 enable\_burst() (py-  
 lablib.devices.AWG.specific.InstekAFG2225  
 method), 421  
 enable\_burst() (py-  
 lablib.devices.AWG.specific.RigolDG1000  
 method), 444  
 enable\_burst() (py-  
 lablib.devices.AWG.specific.RSInstekAFG21000  
 method), 433  
 enable\_burst() (py-  
 lablib.devices.AWG.specific.TektronixAFG1000  
 method), 439  
 enable\_CFR() (pylablib.devices.PhotonFocus.PhotonFocus.IPhotonFocus  
 method), 673  
 enable\_CFR() (pylablib.devices.PhotonFocus.PhotonFocus.PhotonFocus  
 method), 695  
 enable\_CFR() (pylablib.devices.PhotonFocus.PhotonFocus.PhotonFocus  
 method), 679  
 enable\_CFR() (pylablib.devices.PhotonFocus.PhotonFocus.PhotonFocus  
 method), 687  
 enable\_channel() (py-  
 lablib.devices.Tektronix.base.DPO2000  
 method), 763  
 enable\_channel() (py-  
 lablib.devices.Tektronix.base.ITektronixScope  
 method), 750  
 enable\_channel() (py-  
 lablib.devices.Tektronix.base.TDS2000  
 method), 756  
 enable\_channel() (py-  
 lablib.devices.Toptica.ibeam.TopticaIBeam  
 method), 825  
 enable\_channels() (py-  
 lablib.devices.Thorlabs.kinesis.KinesisPiezoMotor  
 method), 801  
 enable\_drift\_compensation() (py-  
 lablib.devices.PhysikInstrumente.base.PIE516  
 method), 705  
 enable\_frame\_transfer\_mode() (py-  
 lablib.devices.Andor.AndorSDK2.AndorSDK2Camera  
 method), 470  
 enable\_limit\_errors() (py-  
 lablib.devices.Arcus.performax.Performax2EXStage  
 method), 494  
 enable\_limit\_errors() (py-  
 lablib.devices.Arcus.performax.Performax4EXStage  
 method), 490  
 enable\_metadata() (py-

<code>lablib.devices.Andor.AndorSDK3.AndorSDK3Camera</code> <code>method)</code> , 481	<code>lablib.devices.AlliedVision.Bonito.BonitoIMAQCamera</code> <code>method)</code> , 457
<code>enable_metadata()</code> <code>lablib.devices.Photometrics.pvcam.PvcamCamera</code> <code>method)</code> , 663	<code>enable_status_line()</code> <code>lablib.devices.AlliedVision.Bonito.IBonitoCamera</code> <code>method)</code> , 452
<code>enable_metadata()</code> <code>lablib.devices.PrincetonInstruments.picam.PicamCamera</code> <code>method)</code> , 713	<code>enable_status_line()</code> <code>lablib.devices.PhotonFocus.PhotonFocus.IPhotonFocusCamera</code> <code>method)</code> , 674
<code>enable_online()</code> <code>lablib.devices.PhysikInstrumente.base.GenericPIController</code> <code>method)</code> , 703	<code>enable_status_line()</code> <code>lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusBitFlowCa</code> <code>method)</code> , 695
<code>enable_online()</code> <code>lablib.devices.PhysikInstrumente.base.PIE516</code> <code>method)</code> , 706	<code>enable_status_line()</code> <code>lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusIMAQCam</code> <code>method)</code> , 679
<code>enable_output()</code> <code>lablib.devices.AWG.generic.GenericAWG</code> <code>method)</code> , 403	<code>enable_status_line()</code> <code>lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusSiSoCamer</code> <code>method)</code> , 687
<code>enable_output()</code> <code>lablib.devices.AWG.specific.Agilent33220A</code> <code>method)</code> , 415	<code>enable_switcher_channel()</code> <code>lablib.devices.HighFinesse.wlm.WLM</code> <code>method)</code> , 559
<code>enable_output()</code> <code>lablib.devices.AWG.specific.Agilent33500</code> <code>method)</code> , 409	<code>enable_sync_output()</code> <code>lablib.devices.AWG.generic.GenericAWG</code> <code>method)</code> , 403
<code>enable_output()</code> <code>lablib.devices.AWG.specific.InstekAFG2000</code> <code>method)</code> , 427	<code>enable_sync_output()</code> <code>lablib.devices.AWG.specific.Agilent33220A</code> <code>method)</code> , 415
<code>enable_output()</code> <code>lablib.devices.AWG.specific.InstekAFG2225</code> <code>method)</code> , 421	<code>enable_sync_output()</code> <code>lablib.devices.AWG.specific.Agilent33500</code> <code>method)</code> , 409
<code>enable_output()</code> <code>lablib.devices.AWG.specific.RigolDG1000</code> <code>method)</code> , 445	<code>enable_sync_output()</code> <code>lablib.devices.AWG.specific.InstekAFG2000</code> <code>method)</code> , 427
<code>enable_output()</code> <code>lablib.devices.AWG.specific.RSInstekAFG21000</code> <code>method)</code> , 433	<code>enable_sync_output()</code> <code>lablib.devices.AWG.specific.InstekAFG2225</code> <code>method)</code> , 421
<code>enable_output()</code> <code>lablib.devices.AWG.specific.TektronixAFG1000</code> <code>method)</code> , 439	<code>enable_sync_output()</code> <code>lablib.devices.AWG.specific.RigolDG1000</code> <code>method)</code> , 445
<code>enable_pixel_correction()</code> <code>lablib.devices.PCO.SC2.PCOSC2Camera</code> <code>method)</code> , 648	<code>enable_sync_output()</code> <code>lablib.devices.AWG.specific.RSInstekAFG21000</code> <code>method)</code> , 433
<code>enable_raw_readout()</code> <code>lablib.devices.Basler.pylon.BaslerPylonCamera</code> <code>method)</code> , 516	<code>enable_sync_output()</code> <code>lablib.devices.AWG.specific.TektronixAFG1000</code> <code>method)</code> , 439
<code>enable_raw_readout()</code> <code>lablib.devices.IMAQdx.IMAQdx.EthernetIMAQdxCamera</code> <code>method)</code> , 584	<code>enable_terascan_updates()</code> <code>lablib.devices.M2.solstis.Solstis</code> <code>method)</code> , 617
<code>enable_raw_readout()</code> <code>lablib.devices.IMAQdx.IMAQdx.IMAQdxCamera</code> <code>method)</code> , 580	<code>enable_trigger_output()</code> <code>lablib.devices.AWG.generic.GenericAWG</code> <code>method)</code> , 405
<code>enable_servo()</code> <code>lablib.devices.PhysikInstrumente.base.PIE516</code> <code>method)</code> , 705	<code>enable_trigger_output()</code> <code>lablib.devices.AWG.specific.Agilent33220A</code> <code>method)</code> , 415
<code>enable_status_line()</code>	<code>enable_trigger_output()</code>



- `pylablib.devices.AWG.specific.Agilent33500`  
`method`), 409
- `enable_trigger_output()` (`py-`  
`pylablib.devices.AWG.specific.InstekAFG2000`  
`method`), 427
- `enable_trigger_output()` (`py-`  
`pylablib.devices.AWG.specific.InstekAFG2225`  
`method`), 421
- `enable_trigger_output()` (`py-`  
`pylablib.devices.AWG.specific.RigolDG1000`  
`method`), 445
- `enable_trigger_output()` (`py-`  
`pylablib.devices.AWG.specific.RSInstekAFG21000`  
`method`), 433
- `enable_trigger_output()` (`py-`  
`pylablib.devices.AWG.specific.TektronixAFG1000`  
`method`), 439
- `enable_updates()` (`py-`  
`pylablib.devices.Attocube.anc350.ANC350`  
`method`), 507
- `enable_velocity_control()` (`py-`  
`pylablib.devices.PhysikInstrumente.base.PIE516`  
`method`), 705
- `enabled` (`pylablib.devices.Lakeshore.base.TLakeshore218`  
`FilterSettings` attribute), 593
- `enabled` (`pylablib.devices.Lakeshore.base.TLakeshore370`  
`FilterSettings` attribute), 599
- `ensure_dir()` (in module `pylablib.core.utils.files`), 363
- `ensure_dir_singlelevel()` (in module `py-`  
`pylablib.core.utils.files`), 363
- `EnumLabel` (class in `pylablib.core.gui.widgets.label`), 236
- `EnumParameterClass` (class in `py-`  
`pylablib.core.devio.interface`), 166
- `eof()` (in module `pylablib.core.utils.files`), 361
- `EPC04` (class in `pylablib.devices.OZOptics.base`), 636
- `EPC04.NoParameterCaller` (class in `py-`  
`pylablib.devices.OZOptics.base`), 637
- `errno` (`pylablib.core.thread.threadprop.TimeoutThreadError`  
attribute), 319
- `errno` (`pylablib.core.utils.net.SocketError` attribute), 388
- `errno` (`pylablib.core.utils.net.SocketTimeout` attribute), 388
- `Error` (`pylablib.core.devio.comm_backend.FT232DeviceBackend`  
attribute), 150
- `Error` (`pylablib.core.devio.comm_backend.IDeviceCommBackend`  
attribute), 143
- `Error` (`pylablib.core.devio.comm_backend.NetworkDeviceBackend`  
attribute), 153
- `Error` (`pylablib.core.devio.comm_backend.PyUSBDeviceBackend`  
attribute), 155
- `Error` (`pylablib.core.devio.comm_backend.RecordedDeviceBackend`  
attribute), 158
- `Error` (`pylablib.core.devio.comm_backend.SerialDeviceBackend`  
attribute), 148
- `Error` (`pylablib.core.devio.comm_backend.VisaDeviceBackend`  
attribute), 145
- `Error` (`pylablib.core.devio.SCPI.SCPIDevice` attribute), 138
- `Error` (`pylablib.devices.AlliedVision.Bonito.BonitoIMAQCamera`  
attribute), 456
- `Error` (`pylablib.devices.AlliedVision.Bonito.IBonitoCamera`  
attribute), 450
- `Error` (`pylablib.devices.Andor.AndorSDK2.AndorSDK2Camera`  
attribute), 467
- `Error` (`pylablib.devices.Andor.AndorSDK3.AndorSDK3Camera`  
attribute), 479
- `Error` (`pylablib.devices.Arcus.performax.GenericPerformaxStage`  
attribute), 488
- `Error` (`pylablib.devices.Arcus.performax.Performax2EXStage`  
attribute), 493
- `Error` (`pylablib.devices.Arcus.performax.Performax4EXStage`  
attribute), 492
- `Error` (`pylablib.devices.Arcus.performax.PerformaxDMX.JSASStage`  
attribute), 498
- `filter_settings` (`pylablib.devices.Arduino.base.IArduinoDevice`  
attribute), 501
- `filter_settings` (`pylablib.devices.Attocube.anc300.ANC300` at-  
tribute), 503
- `filter_settings` (`pylablib.devices.Attocube.anc350.ANC350` at-  
tribute), 506
- `filter_settings` (`pylablib.devices.AWG.generic.GenericAWG` at-  
tribute), 403
- `filter_settings` (`pylablib.devices.AWG.specific.Agilent33220A`  
attribute), 415
- `filter_settings` (`pylablib.devices.AWG.specific.Agilent33500` at-  
tribute), 409
- `filter_settings` (`pylablib.devices.AWG.specific.InstekAFG2000`  
attribute), 426
- `filter_settings` (`pylablib.devices.AWG.specific.InstekAFG2225`  
attribute), 421
- `filter_settings` (`pylablib.devices.AWG.specific.RigolDG1000` at-  
tribute), 444
- `filter_settings` (`pylablib.devices.AWG.specific.RSInstekAFG21000`  
attribute), 432
- `filter_settings` (`pylablib.devices.AWG.specific.TektronixAFG1000`  
attribute), 438
- `filter_settings` (`pylablib.devices.Basler.pylon.BaslerPylonCamera`  
attribute), 514
- `filter_settings` (`pylablib.devices.BitFlow.BitFlow.BitFlowCamera`  
attribute), 527
- `filter_settings` (`pylablib.devices.BitFlow.BitFlow.BitFlowFrameGrabber`  
attribute), 521
- `filter_settings` (`pylablib.devices.Conrad.base.RelayBoard` at-  
tribute), 533
- `filter_settings` (`pylablib.devices.Cryocon.base.Cryocon1x` at-

- tribute*), 536
- Error (*pylablib.devices.Cryomagnetics.base.LM500 attribute*), 540
- Error (*pylablib.devices.Cryomagnetics.base.LM510 attribute*), 544
- Error (*pylablib.devices.DCAM.DCAM.DCAMCamera attribute*), 550
- Error (*pylablib.devices.HighFinesse.wlm.WLM attribute*), 557
- Error (*pylablib.devices.IMAQ.IMAQ.IMAQCamera attribute*), 568
- Error (*pylablib.devices.IMAQ.IMAQ.IMAQFrameGrabber attribute*), 561
- Error (*pylablib.devices.IMAQdx.IMAQdx.EthernetIMAQdxCamera attribute*), 584
- Error (*pylablib.devices.IMAQdx.IMAQdx.IMAQdxCamera attribute*), 578
- Error (*pylablib.devices.interface.camera.IAttributeCamera attribute*), 837
- Error (*pylablib.devices.interface.camera.IBinROICamera attribute*), 856
- Error (*pylablib.devices.interface.camera.ICamera attribute*), 831
- Error (*pylablib.devices.interface.camera.IExposureCamera attribute*), 847
- Error (*pylablib.devices.interface.camera.IGrabberAttributeCamera attribute*), 842
- Error (*pylablib.devices.interface.camera.IROICamera attribute*), 852
- Error (*pylablib.devices.KJL.base.KJL300 attribute*), 590
- Error (*pylablib.devices.Lakeshore.base.Lakeshore218 attribute*), 593
- Error (*pylablib.devices.Lakeshore.base.Lakeshore370 attribute*), 599
- Error (*pylablib.devices.LaserQuantum.base.Finesse attribute*), 604
- Error (*pylablib.devices.Leybold.base.GenericITR attribute*), 607
- Error (*pylablib.devices.Leybold.base.ITR90 attribute*), 609
- error (*pylablib.devices.Leybold.base.TUpdateValue attribute*), 607
- Error (*pylablib.devices.LighthousePhotonics.base.SproutG attribute*), 612
- Error (*pylablib.devices.M2.solstis.Solstis attribute*), 614
- Error (*pylablib.devices.Newport.picomotor.Picomotor8742 attribute*), 628
- Error (*pylablib.devices.NI.daq.NIDAQ attribute*), 621
- Error (*pylablib.devices.Ophir.base.OphirDevice attribute*), 639
- Error (*pylablib.devices.Ophir.base.VegaPowerMeter attribute*), 642
- Error (*pylablib.devices.OZOptics.base.DD100 attribute*), 635
- Error (*pylablib.devices.OZOptics.base.EPC04 attribute*), 636
- Error (*pylablib.devices.OZOptics.base.OZOpticsDevice attribute*), 632
- Error (*pylablib.devices.OZOptics.base.TF100 attribute*), 633
- Error (*pylablib.devices.PCO.SC2.PCOSC2Camera attribute*), 645
- Error (*pylablib.devices.Pfeiffer.base.DPG202 attribute*), 657
- Error (*pylablib.devices.Pfeiffer.base.TPG260 attribute*), 654
- Error (*pylablib.devices.Photometrics.pvcam.PvcamCamera attribute*), 662
- Error (*pylablib.devices.PhotonFocus.PhotonFocus.IPhotonFocusCamera attribute*), 672
- Error (*pylablib.devices.PhotonFocus.PhotonFocus.PhotonFocusBitFlowC attribute*), 694
- Error (*pylablib.devices.PhotonFocus.PhotonFocus.PhotonFocusIMAQCar attribute*), 678
- Error (*pylablib.devices.PhotonFocus.PhotonFocus.PhotonFocusSiSoCam attribute*), 687
- Error (*pylablib.devices.PhysikInstrumente.base.GenericPIController attribute*), 703
- Error (*pylablib.devices.PhysikInstrumente.base.PIE516 attribute*), 706
- Error (*pylablib.devices.PrincetonInstruments.picam.PicamCamera attribute*), 712
- Error (*pylablib.devices.SiliconSoftware.fgrab.SiliconSoftwareCamera attribute*), 727
- Error (*pylablib.devices.SiliconSoftware.fgrab.SiliconSoftwareFrameGrab attribute*), 721
- Error (*pylablib.devices.Sirah.Matisse.SirahMatisse attribute*), 736
- Error (*pylablib.devices.SmarAct.scu3d.SCU3D attribute*), 746
- Error (*pylablib.devices.Tektronix.base.DPO2000 attribute*), 762
- Error (*pylablib.devices.Tektronix.base.ITektronixScope attribute*), 749
- Error (*pylablib.devices.Tektronix.base.TDS2000 attribute*), 755
- Error (*pylablib.devices.Thorlabs.elliptec.ElliptecMotor attribute*), 778
- Error (*pylablib.devices.Thorlabs.kinesis.BasicKinesisDevice attribute*), 782
- Error (*pylablib.devices.Thorlabs.kinesis.KinesisDevice attribute*), 788
- Error (*pylablib.devices.Thorlabs.kinesis.KinesisMotor attribute*), 798
- Error (*pylablib.devices.Thorlabs.kinesis.KinesisPiezoMotor attribute*), 803

- Error (*pylablib.devices.Thorlabs.kinesis.KinesisQuadDetector* attribute), 807
- Error (*pylablib.devices.Thorlabs.kinesis.MFF* attribute), 792
- Error (*pylablib.devices.Thorlabs.serial.FW* attribute), 813
- Error (*pylablib.devices.Thorlabs.serial.FWv1* attribute), 817
- Error (*pylablib.devices.Thorlabs.serial.MDT69xA* attribute), 820
- Error (*pylablib.devices.Thorlabs.serial.ThorlabsSerialInterface* attribute), 809
- Error (*pylablib.devices.Thorlabs.TLCamera.ThorlabsTLCamera* attribute), 770
- Error (*pylablib.devices.Toptica.ibeam.TopticaIBeam* attribute), 824
- Error (*pylablib.devices.Trinamic.base.TMCM1110* attribute), 826
- Error (*pylablib.devices.uc480.uc480.UC480Camera* attribute), 867
- errors (*pylablib.devices.PCO.SC2.TCCameraStatus* attribute), 644
- escape\_string() (in module *pylablib.core.utils.string*), 398
- EthernetIMAQdxCamera (class in *pylablib.devices.IMAQdx.IMAQdx*), 584
- EthernetIMAQdxCamera.NoParameterCaller (class in *pylablib.devices.IMAQdx.IMAQdx*), 584
- exc\_mode (*pylablib.devices.Lakeshore.base.TLakeshore370RangeSettings* attribute), 598
- exc\_range (*pylablib.devices.Lakeshore.base.TLakeshore370RangeSettings* attribute), 598
- execute() (*pylablib.core.thread.callsync.QScheduledCall* method), 283
- exint() (in module *pylablib.core.thread.controller*), 292
- exists (*pylablib.devices.PrincetonInstruments.picam.PicamAttribute* attribute), 709
- exp\_decay\_k() (in module *pylablib.core.dataproc.specfunc*), 124
- expand\_relative\_path() (in module *pylablib.core.utils.module*), 386
- ExpandedContainerDictionaryEntry (class in *pylablib.core.fileio.dict\_entry*), 176
- export\_clock() (*pylablib.devices.NI.daq.NIDAQ* method), 622
- exposure (*pylablib.devices.Andor.AndorSDK2.TCcycleTimings* attribute), 466
- exposure (*pylablib.devices.interface.camera.TAcqTimings* attribute), 846
- exposure\_ns (*pylablib.devices.Photometrics.pvcam.TFrameInfo* attribute), 661
- exsafe() (in module *pylablib.core.thread.controller*), 292
- exsafeSlot() (in module *pylablib.core.thread.controller*), 292
- ExternalBinTableDictionaryEntry (class in *pylablib.core.fileio.dict\_entry*), 173
- ExternalNumpyDictionaryEntry (class in *pylablib.core.fileio.dict\_entry*), 175
- ExternalTextTableDictionaryEntry (class in *pylablib.core.fileio.dict\_entry*), 172
- extract\_escaped\_string() (in module *pylablib.core.utils.string*), 399
- extract\_status\_line() (in module *pylablib.devices.interface.camera*), 861
- ## F
- f (*pylablib.core.utils.files.TempFile* attribute), 362
- fail() (*pylablib.core.thread.callsync.QScheduledCall* method), 284
- fail() (*pylablib.core.thread.synchronizing.QMultiThreadNotifier* method), 319
- fail\_exec\_point() (*pylablib.core.thread.controller.QTaskThread* method), 309
- fail\_exec\_point() (*pylablib.core.thread.controller.QThreadController* method), 299
- failed() (*pylablib.core.thread.callsync.QCallResultSynchronizer* method), 281
- failed() (*pylablib.core.thread.callsync.QDirectResultSynchronizer* method), 282
- fall\_speed (*pylablib.devices.Sirah.Matisse.TSscanParameters* attribute), 736
- falling (*pylablib.devices.Sirah.Matisse.TSscanMode* attribute), 736
- fast\_shift\_roi() (*pylablib.devices.PhotonFocus.PhotonFocus.IPhotonFocusCamera* method), 673
- fast\_shift\_roi() (*pylablib.devices.PhotonFocus.PhotonFocus.PhotonFocusBitFlowCamera* method), 695
- fast\_shift\_roi() (*pylablib.devices.PhotonFocus.PhotonFocus.PhotonFocusIMAQCamera* method), 679
- fast\_shift\_roi() (*pylablib.devices.PhotonFocus.PhotonFocus.PhotonFocusSiSoCamera* method), 687
- FGrabAttribute (class in *pylablib.devices.SiliconSoftware.fgrab*), 719
- file (*pylablib.devices.SiliconSoftware.fgrab.TAppletInfo* attribute), 718
- file (*pylablib.devices.SiliconSoftware.fgrab.TFullAppletInfo* attribute), 719
- file\_format (*pylablib.core.fileio.dict\_entry.ExternalNumpyDictionaryEntry* attribute), 175

[file\\_format \(pylablib.core.fileio.dict\\_entry.IExternalFileDictEntry attribute\), 174](#)  
[filename \(pylablib.core.thread.threadprop.TimeoutThreadError attribute\), 320](#)  
[filename \(pylablib.core.utils.net.SocketError attribute\), 388](#)  
[filename \(pylablib.core.utils.net.SocketTimeout attribute\), 388](#)  
[filename2 \(pylablib.core.thread.threadprop.TimeoutThreadError attribute\), 320](#)  
[filename2 \(pylablib.core.utils.net.SocketError attribute\), 388](#)  
[filename2 \(pylablib.core.utils.net.SocketTimeout attribute\), 388](#)  
[files \(pylablib.core.utils.files.FolderList attribute\), 364](#)  
[fill\\_voltage\\_output\\_buffer \(\) \(pylablib.devices.NI.daq.NIDAQ method\), 625](#)  
[filt \(pylablib.core.utils.observer\\_pool.ObserverPool.ObserverPool attribute\), 393](#)  
[filter\\_args\\_dict \(\) \(pylablib.core.dataproc.callable.FunctionCallable method\), 107](#)  
[filter\\_args\\_dict \(\) \(pylablib.core.dataproc.callable.ICallable method\), 104](#)  
[filter\\_args\\_dict \(\) \(pylablib.core.dataproc.callable.JoinedCallable method\), 106](#)  
[filter\\_args\\_dict \(\) \(pylablib.core.dataproc.callable.MethodCallable method\), 108](#)  
[filter\\_args\\_dict \(\) \(pylablib.core.dataproc.callable.MultiplexedCallable method\), 105](#)  
[filter\\_by \(\) \(in module pylablib.core.dataproc.utils\), 136](#)  
[filter\\_dict \(\) \(in module pylablib.core.utils.general\), 375](#)  
[filter\\_limiter \(\) \(in module pylablib.core.gui.limiter\), 262](#)  
[filter\\_self \(\) \(pylablib.core.utils.dictionary.Dictionary method\), 333](#)  
[filter\\_self \(\) \(pylablib.core.utils.dictionary.DictionaryPointer method\), 337](#)  
[filter\\_self \(\) \(pylablib.core.utils.dictionary.FilterTree method\), 354](#)  
[filter\\_self \(\) \(pylablib.core.utils.dictionary.PrefixTree method\), 346](#)  
[filter\\_string\\_list \(\) \(in module pylablib.core.utils.string\), 398](#)  
[FileDictEntry \(class in pylablib.core.utils.dictionary\), 352](#)  
[finalize\\_task \(\) \(pylablib.core.thread.controller.QTaskThread method\), 305](#)  
[finalized \(pylablib.core.thread.controller.QThreadControllerThread attribute\), 292](#)  
[find\\_all\\_first\\_locations \(\) \(in module pylablib.core.utils.string\), 397](#)  
[find\\_all\\_prefixes \(\) \(pylablib.core.utils.dictionary.PrefixTree method\), 344](#)  
[find\\_by\\_serial \(\) \(in module pylablib.devices.uc480.uc480\), 864](#)  
[find\\_by\\_serial \(\) \(pylablib.devices.uc480.uc480.UC480Camera static method\), 867](#)  
[find\\_closest\\_arg \(\) \(in module pylablib.core.dataproc.utils\), 136](#)  
[find\\_closest\\_value \(\) \(in module pylablib.core.dataproc.utils\), 136](#)  
[find\\_columns\\_lines \(\) \(in module pylablib.core.fileio.loadfile\\_utils\), 182](#)  
[find\\_dict\\_string \(\) \(in module pylablib.core.utils.string\), 397](#)  
[find\\_discrete\\_step \(\) \(in module pylablib.core.dataproc.utils\), 137](#)  
[find\\_first\\_entry \(\) \(in module pylablib.core.utils.string\), 397](#)  
[find\\_intersection \(\) \(pylablib.core.utils.dictionary.Dictionary static method\), 333](#)  
[find\\_intersection \(\) \(pylablib.core.utils.dictionary.DictionaryPointer static method\), 338](#)  
[find\\_intersection \(\) \(pylablib.core.utils.dictionary.FilterTree static method\), 354](#)  
[find\\_intersection \(\) \(pylablib.core.utils.dictionary.PrefixTree static method\), 346](#)  
[find\\_largest\\_prefix \(\) \(pylablib.core.utils.dictionary.PrefixTree method\), 344](#)  
[find\\_layout\\_element \(\) \(in module pylablib.core.gui.utils\), 263](#)  
[find\\_list\\_string \(\) \(in module pylablib.core.utils.string\), 396](#)  
[find\\_local\\_extrema \(\) \(in module pylablib.core.dataproc.feature\), 110](#)  
[find\\_observers \(\) \(pylablib.core.utils.observer\\_pool.ObserverPool method\), 393](#)  
[find\\_peaks\\_cutoff \(\) \(in module py-](#)



`lablib.core.dataproc.feature`), 109  
`find_savetime_comment()` (in module `pylablib.core.fileio.loadfile_utils`), 182  
`find_skipped_frames()` (in module `pylablib.devices.PhotonFocus.PhotonFocus`), 702  
`fine_tune_wavelength()` (`pylablib.devices.M2.solstis.Solstis` method), 615  
`Finesse` (class in `pylablib.devices.LaserQuantum.base`), 604  
`Finesse.NoParameterCaller` (class in `pylablib.devices.LaserQuantum.base`), 605  
`finished` (`pylablib.core.thread.controller.QTaskThread` attribute), 309  
`finished` (`pylablib.core.thread.controller.QThreadController` attribute), 294  
`finishing()` (`pylablib.core.thread.controller.QTaskThread` method), 309  
`finishing()` (`pylablib.core.thread.controller.QThreadController` method), 299  
`firmware_version` (`pylablib.devices.Andor.AndorSDK3.TDeviceInfo` attribute), 478  
`firmware_version` (`pylablib.devices.SmarAct.scu3d.TDeviceInfo` attribute), 745  
`firmware_version` (`pylablib.devices.Thorlabs.TLCamera.TDeviceInfo` attribute), 769  
`fit()` (`pylablib.core.dataproc.fitting.Fitter` method), 115  
`Fitter` (class in `pylablib.core.dataproc.fitting`), 114  
`fixed_size` (`pylablib.core.utils.ipc.TShmemVarDesc` attribute), 384  
`flags` (`pylablib.devices.IMAQdx.IMAQdx.TCameraInfo` attribute), 576  
`flags` (`pylablib.devices.Photometrics.pvcam.TFrameInfo` attribute), 661  
`flags` (`pylablib.devices.PrincetonInstruments.picam.TROIConstraints` attribute), 709  
`flags` (`pylablib.devices.SiliconSoftware.fgrab.TFullAppletInfo` attribute), 719  
`flags` (`pylablib.devices.uc480.uc480.TFrameInfo` attribute), 866  
`flatten_list()` (in module `pylablib.core.utils.general`), 376  
`flip_byteorder()` (`pylablib.core.devio.data_format.DataFormat` method), 161  
`flip_fourier_transform()` (in module `pylablib.core.dataproc.fourier`), 119  
`float_to_str_SI()` (in module `pylablib.core.gui.formatter`), 261  
`FloatFormatter` (class in `pylablib.core.gui.formatter`), 261  
`flush()` (`pylablib.core.devio.SCPISCPIDevice` method), 141  
`flush()` (`pylablib.core.utils.general.StreamFileLogger` method), 380  
`flush()` (`pylablib.devices.AWG.generic.GenericAWG` method), 406  
`flush()` (`pylablib.devices.AWG.specific.Agilent33220A` method), 415  
`flush()` (`pylablib.devices.AWG.specific.Agilent33500` method), 409  
`flush()` (`pylablib.devices.AWG.specific.InstekAFG2000` method), 427  
`flush()` (`pylablib.devices.AWG.specific.InstekAFG2225` method), 421  
`flush()` (`pylablib.devices.AWG.specific.RigolDG1000` method), 445  
`flush()` (`pylablib.devices.AWG.specific.RSInstekAFG21000` method), 433  
`flush()` (`pylablib.devices.AWG.specific.TektronixAFG1000` method), 439  
`flush()` (`pylablib.devices.Basler.pylon.BaslerPylonCamera.BufferManager` method), 516  
`flush()` (`pylablib.devices.Cryocon.base.Cryocon1x` method), 537  
`flush()` (`pylablib.devices.Cryomagnetics.base.LM500` method), 541  
`flush()` (`pylablib.devices.Cryomagnetics.base.LM510` method), 544  
`flush()` (`pylablib.devices.Lakeshore.base.Lakeshore218` method), 595  
`flush()` (`pylablib.devices.Lakeshore.base.Lakeshore370` method), 600  
`flush()` (`pylablib.devices.M2.solstis.Solstis` method), 615  
`flush()` (`pylablib.devices.Sirah.Matisse.SirahMatisse` method), 740  
`flush()` (`pylablib.devices.Tektronix.base.DPO2000` method), 763  
`flush()` (`pylablib.devices.Tektronix.base.ITektronixScope` method), 753  
`flush()` (`pylablib.devices.Tektronix.base.TDS2000` method), 756  
`flush()` (`pylablib.devices.Thorlabs.serial.FW` method), 813  
`flush()` (`pylablib.devices.Thorlabs.serial.FWv1` method), 817  
`flush()` (`pylablib.devices.Thorlabs.serial.MDT69xA` method), 820  
`flush()` (`pylablib.devices.Thorlabs.serial.ThorlabsSerialInterface` method), 810  
`flush_read()` (`pylablib.core.devio.comm_backend.FT232DeviceBackend` method), 152

flush\_read() (pylablib.core.devio.comm\_backend.IDeviceCommBackend (pylablib.devices.interface.camera.TFrameInfo attribute), 144

flush\_read() (pylablib.core.devio.comm\_backend.NetworkDeviceBackend (pylablib.devices.PCO.SC2.TFrameInfo attribute), 154

flush\_read() (pylablib.core.devio.comm\_backend.PyUSBDeviceBackend (pylablib.devices.Photometrics.pvcam.TFrameInfo attribute), 157

flush\_read() (pylablib.core.devio.comm\_backend.RecordableDeviceBackend (pylablib.devices.PrincetonInstruments.picam.TFrameInfo attribute), 159

flush\_read() (pylablib.core.devio.comm\_backend.SerialDeviceBackend (pylablib.devices.SiliconSoftware.fgrab.TFrameInfo attribute), 149

flush\_read() (pylablib.core.devio.comm\_backend.VisaDeviceBackend (pylablib.devices.Thorlabs.TLCAmera.TFrameInfo attribute), 146

fmt (pylablib.devices.Lakeshore.base.TLakeshore218CurveHeader\_index (pylablib.devices.uc480.uc480.TFrameInfo attribute), 593

FmtStringFormatter (class in py-lablib.core.gui.formatter), 261

focusInEvent() (py-lablib.core.gui.widgets.edit.TextEdit method), 234

focusOutEvent() (py-lablib.core.gui.widgets.edit.TextEdit method), 234

FolderFileSystemDataLocation (class in py-lablib.core.fileio.location), 187

FolderList (class in pylablib.core.utils.files), 363

folders (pylablib.core.utils.files.FolderList attribute), 364

followed() (pylablib.core.dataproc.transform.Indexed2DTransform method), 134

followed() (pylablib.core.dataproc.transform.LinearTransform method), 134

force\_trigger() (py-lablib.devices.Tektronix.base.DPO2000 method), 763

force\_trigger() (py-lablib.devices.Tektronix.base.ITektronixScope method), 750

force\_trigger() (py-lablib.devices.Tektronix.base.TDS2000 method), 756

fourier\_filter() (in module py-lablib.core.dataproc.filters), 113

fourier\_filter\_bandpass() (in module py-lablib.core.dataproc.filters), 113

fourier\_filter\_bandstop() (in module py-lablib.core.dataproc.filters), 113

fourier\_make\_response\_real() (in module py-lablib.core.dataproc.filters), 113

fourier\_transform() (in module py-lablib.core.dataproc.fourier), 118

frame\_index (pylablib.devices.Andor.AndorSDK3.TFrameInfo attribute), 479

frame\_index (pylablib.devices.DCAM.DCAM.TFrameInfo attribute), 550

frame\_period (pylablib.devices.interface.camera.TAcqTimings attribute), 846

FrameCounter (class in py-lablib.devices.interface.camera), 835

FrameNotifier (class in py-lablib.devices.interface.camera), 836

frames\_done (pylablib.devices.Andor.AndorSDK2.TAcqProgress attribute), 466

frameskip\_events (py-lablib.devices.uc480.uc480.TAcquiredFramesStatus attribute), 865

framestamp (pylablib.devices.AlliedVision.Bonito.TStatusLine attribute), 464

framestamp (pylablib.devices.DCAM.DCAM.TFrameInfo attribute), 550

framestamp (pylablib.devices.PCO.SC2.TStatusLine attribute), 652

framestamp (pylablib.devices.Photometrics.pvcam.TFrameInfo attribute), 661

framestamp (pylablib.devices.PrincetonInstruments.picam.TFrameInfo attribute), 712

framestamp (pylablib.devices.SiliconSoftware.fgrab.TFrameInfo attribute), 721

framestamp (pylablib.devices.Thorlabs.TLCAmera.TFrameInfo attribute), 770

framestamp (pylablib.devices.uc480.uc480.TFrameInfo attribute), 866

framestamp\_checker (py-lablib.devices.interface.camera.TStatusLineDescription attribute), 861

FrameTransferError (py-lablib.devices.AlliedVision.Bonito.BonitoIMAQCamera attribute), 456

FrameTransferError (py-lablib.devices.AlliedVision.Bonito.IBonitoCamera attribute), 452

FrameTransferError (py-lablib.devices.Andor.AndorSDK2.AndorSDK2Camera attribute), 473

FrameTransferError (py-

923

from_columns()	(py- lablib.core.dataproc.table_wrap.Array1DWrapper class method), 127	from_json() (pylablib.core.utils.dictionary.Dictionary class method), 332
from_columns()	(py- lablib.core.dataproc.table_wrap.Array2DWrapper class method), 131	from_json() (pylablib.core.utils.dictionary.DictionaryPointer class method), 338
from_columns()	(py- lablib.core.dataproc.table_wrap.DataFrame2DWrapper class method), 133	from_json() (pylablib.core.utils.dictionary.FilterTree class method), 355
from_columns()	(py- lablib.core.dataproc.table_wrap.I1DWrapper class method), 126	from_json() (pylablib.core.utils.dictionary.PrefixTree class method), 347
from_columns()	(py- lablib.core.dataproc.table_wrap.I2DWrapper class method), 129	from_object() (py- lablib.core.fileio.location.LocationName static method), 184
from_columns()	(py- lablib.core.dataproc.table_wrap.Series1DWrapper class method), 128	from_Pa() (pylablib.devices.Pfeiffer.base.TPG260 method), 654
from_data() (in module py- lablib.core.fileio.dict_entry), 170		from_row_string() (in module py- lablib.core.utils.string), 400
from_data() (pylablib.core.fileio.dict_entry.DictEntryBuilder method), 169		from_string() (in module pylablib.core.utils.string), 400
from_desc() (pylablib.core.devio.data_format.DataFormat static method), 162		from_string() (py- lablib.core.fileio.location.LocationName static method), 184
from_desc_SCPI() (py- lablib.core.devio.data_format.DataFormat static method), 162		full_exit() (in module pylablib.core.utils.general), 377
from_dict() (in module py- lablib.core.fileio.dict_entry), 170		full_name (pylablib.core.utils.files.TempFile at- tribute), 363
from_dict() (pylablib.core.fileio.dict_entry.DictEntryParser method), 169		full_name (pylablib.devices.SiliconSoftware.fgrab.TBoardInfo attribute), 718
from_dict() (pylablib.core.fileio.dict_entry.ExpandedCommandParser class method), 176		fullsplit() (in module pylablib.core.utils.files), 362
from_dict() (pylablib.core.fileio.dict_entry.ExternalBinaryTableDictionaryEntry class method), 174		fullsplit() (in module pylablib.core.utils.functions), 371
from_dict() (pylablib.core.fileio.dict_entry.ExternalNumPyTableDictionaryEntry class method), 175		FullTableDictionaryEntry (class in py- lablib.core.dataproc.callable), 106
from_dict() (pylablib.core.fileio.dict_entry.ExternalTextTableDictionaryEntry class method), 173		FullTableDictionaryEntry.NamesBoundCall (class in pylablib.core.dataproc.callable), 106
from_dict() (pylablib.core.fileio.dict_entry.IDictionaryEntry class method), 170		FunctionParameterClass (class in py- lablib.core.devio.interface), 167
from_dict() (pylablib.core.fileio.dict_entry.IDictionaryEntry method), 170		FunctionSignature (class in py- lablib.core.utils.functions), 369
from_dict() (pylablib.core.fileio.dict_entry.IExternalFileDictionaryEntry class method), 174		FullTableDictionaryEntry (class in py- lablib.devices.Thorlabs.serial), 812
from_dict() (pylablib.core.fileio.dict_entry.IExternalTableDictionaryEntry class method), 172		FW.NoParameterCaller (class in py- lablib.devices.Thorlabs.serial), 813
from_dict() (pylablib.core.fileio.dict_entry.InlineTableDictionaryEntry class method), 172		fw_freq (pylablib.devices.Thorlabs.elliptec.TMotorInfo attribute), 777
from_dict() (pylablib.core.fileio.dict_entry.ITableDictionaryEntry class method), 171		fw_ver (pylablib.devices.Thorlabs.elliptec.TDeviceInfo attribute), 777
from_function() (py- lablib.core.utils.functions.FunctionSignature static method), 371		fw_ver (pylablib.devices.Thorlabs.kinesis.TDeviceInfo attribute), 781
		FWv1 (class in pylablib.devices.Thorlabs.serial), 816
		FWv1.NoParameterCaller (class in py-



*lablib.devices.Thorlabs.serial*), 817

## G

- gain\_idx* (*pylablib.devices.Photometrics.pvcam.TReadoutInfo* attribute), 662
- gain\_name* (*pylablib.devices.Photometrics.pvcam.TReadoutInfo* attribute), 662
- gaussian\_filter*() (in module *pylablib.core.dataproc.filters*), 111
- gaussian\_filter\_nd*() (in module *pylablib.core.dataproc.filters*), 111
- gaussian\_k*() (in module *pylablib.core.dataproc.specfunc*), 124
- gcd*() (in module *pylablib.core.utils.numerical*), 391
- gcd\_approx*() (in module *pylablib.core.utils.numerical*), 391
- gen\_hamming\_w*() (in module *pylablib.core.dataproc.specfunc*), 124
- gen\_hamming\_w\_ft*() (in module *pylablib.core.dataproc.specfunc*), 125
- generate\_indexed\_filename*() (in module *pylablib.core.utils.files*), 361
- generate\_new\_name*() (*pylablib.core.fileio.location.FolderFileSystemDataLocation* method), 188
- generate\_new\_name*() (*pylablib.core.fileio.location.IDataLocation* method), 184
- generate\_new\_name*() (*pylablib.core.fileio.location.IFileSystemDataLocation* method), 186
- generate\_new\_name*() (*pylablib.core.fileio.location.OpenedFileLocation* method), 185
- generate\_new\_name*() (*pylablib.core.fileio.location.PrefixedFileSystemDataLocation* method), 187
- generate\_new\_name*() (*pylablib.core.fileio.location.SingleFileSystemDataLocation* method), 186
- generate\_prefixed\_filename*() (in module *pylablib.core.utils.files*), 362
- generate\_temp\_filename*() (in module *pylablib.core.utils.files*), 362
- GenericAWG* (class in *pylablib.devices.AWG.generic*), 403
- GenericAWG.NoParameterCaller* (class in *pylablib.devices.AWG.generic*), 406
- GenericAWGBackendError*, 402
- GenericAWGError*, 402
- GenericITR* (class in *pylablib.devices.Leybold.base*), 607
- GenericITR.NoParameterCaller* (class in *pylablib.devices.Leybold.base*), 608
- GenericPerformaxStage* (class in *pylablib.devices.Arcus.performax*), 488
- GenericPerformaxStage.NoParameterCaller* (class in *pylablib.devices.Arcus.performax*), 488
- GenericPIController* (class in *pylablib.devices.PhysikInstrumente.base*), 703
- GenericPIController.NoParameterCaller* (class in *pylablib.devices.PhysikInstrumente.base*), 704
- GenericSirahBackendError*, 733
- GenericSirahError*, 733
- get*() (*pylablib.core.dataproc.filters.RunningDebounceFilter* method), 114
- get*() (*pylablib.core.dataproc.filters.RunningDecimationFilter* method), 114
- get*() (*pylablib.core.fileio.datafile.DataFile* method), 169
- get*() (*pylablib.core.utils.dictionary.Dictionary* method), 328
- get*() (*pylablib.core.utils.dictionary.DictionaryPointer* method), 338
- get*() (*pylablib.core.utils.dictionary.FilterTree* method), 355
- get*() (*pylablib.core.utils.dictionary.ItemAccessor* method), 361
- get*() (*pylablib.core.utils.dictionary.PrefixTree* method), 347
- get*() (*pylablib.core.utils.functions.AttrObjectProperty* method), 373
- get*() (*pylablib.core.utils.functions.IObjectProperty* method), 373
- get*() (*pylablib.core.utils.functions.MethodObjectProperty* method), 373
- get\_acceleration\_factor*() (*pylablib.devices.Trinamic.base.TMCM1110* method), 828
- get\_accum\_mode\_parameters*() (*pylablib.devices.Andor.AndorSDK2.AndorSDK2Camera* method), 470
- get\_acquired\_frame\_status*() (*pylablib.devices.uc480.uc480.UC480Camera* method), 868
- get\_acquisition\_mode*() (*pylablib.devices.Andor.AndorSDK2.AndorSDK2Camera* method), 470
- get\_acquisition\_parameters*() (*pylablib.devices.AlliedVision.Bonito.BonitoIMAQCamera* method), 457
- get\_acquisition\_parameters*() (*pylablib.devices.AlliedVision.Bonito.IBonitoCamera* method), 452
- get\_acquisition\_parameters*() (*pylablib.devices.Andor.AndorSDK2.AndorSDK2Camera* method), 470

method), 473

get\_acquisition\_parameters() (py-lablib.devices.Andor.AndorSDK3.AndorSDK3Camera method), 483

get\_acquisition\_parameters() (py-lablib.devices.Basler.pylon.BaslerPylonCamera method), 517

get\_acquisition\_parameters() (py-lablib.devices.BitFlow.BitFlow.BitFlowCamera method), 528

get\_acquisition\_parameters() (py-lablib.devices.BitFlow.BitFlow.BitFlowFrameGrabber method), 523

get\_acquisition\_parameters() (py-lablib.devices.DCAM.DCAM.DCAMCamera method), 553

get\_acquisition\_parameters() (py-lablib.devices.IMAQ.IMAQ.IMAQCamera method), 570

get\_acquisition\_parameters() (py-lablib.devices.IMAQ.IMAQ.IMAQFrameGrabber method), 565

get\_acquisition\_parameters() (py-lablib.devices.IMAQdx.IMAQdx.EthernetIMAQdxCamera method), 585

get\_acquisition\_parameters() (py-lablib.devices.IMAQdx.IMAQdx.IMAQdxCamera method), 580

get\_acquisition\_parameters() (py-lablib.devices.interface.camera.IAttributeCamera method), 838

get\_acquisition\_parameters() (py-lablib.devices.interface.camera.IBinROICamera method), 857

get\_acquisition\_parameters() (py-lablib.devices.interface.camera.ICamera method), 831

get\_acquisition\_parameters() (py-lablib.devices.interface.camera.IExposureCamera method), 847

get\_acquisition\_parameters() (py-lablib.devices.interface.camera.IGrabberAttributeCamera method), 843

get\_acquisition\_parameters() (py-lablib.devices.interface.camera.IROICamera method), 852

get\_acquisition\_parameters() (py-lablib.devices.PCO.SC2.PCOSC2Camera method), 649

get\_acquisition\_parameters() (py-lablib.devices.Photometrics.pvcam.PvcamCamera method), 665

get\_acquisition\_parameters() (py-lablib.devices.PhotonFocus.PhotonFocus.IPhotonFocusCamera method), 674

get\_acquisition\_parameters() (py-lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusBitFlowCamera method), 696

get\_acquisition\_parameters() (py-lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusIMAQCamera method), 680

get\_acquisition\_parameters() (py-lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusSiSoCamera method), 688

get\_acquisition\_parameters() (py-lablib.devices.PrincetonInstruments.picam.PicamCamera method), 714

get\_acquisition\_parameters() (py-lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareCamera method), 728

get\_acquisition\_parameters() (py-lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareFrameGrabber method), 723

get\_acquisition\_parameters() (py-lablib.devices.Thorlabs.TLCamera.ThorlabsTLCamera method), 773

get\_acquisition\_parameters() (py-lablib.devices.uc480.uc480.UC480Camera method), 869

get\_acquisition\_progress() (py-lablib.devices.Andor.AndorSDK2.AndorSDK2Camera method), 472

get\_active\_channel() (py-lablib.devices.HighFinesse.wlm.WLM method), 559

get\_addr() (pylablib.devices.Newport.picomotor.Picomotor8742 method), 629

get\_addr\_map() (py-lablib.devices.Newport.picomotor.Picomotor8742 method), 629

get\_all\_amp\_modes() (py-lablib.devices.Andor.AndorSDK2.AndorSDK2Camera method), 468

get\_all\_attribute\_values() (py-lablib.devices.Andor.AndorSDK3.AndorSDK3Camera method), 480

get\_all\_attribute\_values() (py-lablib.devices.Basler.pylon.BaslerPylonCamera method), 514

get\_all\_attribute\_values() (py-lablib.devices.DCAM.DCAM.DCAMCamera method), 551

get\_all\_attribute\_values() (py-lablib.devices.IMAQdx.IMAQdx.EthernetIMAQdxCamera method), 585

get\_all\_attribute\_values() (py-lablib.devices.IMAQdx.IMAQdx.IMAQdxCamera method), 579

<code>get_all_attribute_values()</code> lablib.devices.interface.camera.IAttributeCamera method), 837	(py- lablib.devices.PrincetonInstruments.picam.PicamCamera method), 714	<code>get_all_attributes()</code> (py- lablib.devices.Arcus.performax.GenericPerformaxStage method), 489
<code>get_all_attribute_values()</code> lablib.devices.Photometrics.pvcam.PvcamCamera method), 662	<code>get_all_axes()</code> lablib.devices.Arcus.performax.Performax2EXStage method), 494	<code>get_all_axes()</code> (py- lablib.devices.Arcus.performax.Performax4EXStage method), 492
<code>get_all_attribute_values()</code> lablib.devices.PhotonFocus.PhotonFocus.IPhotonFocusCamera method), 672	<code>get_all_axes()</code> lablib.devices.Arcus.performax.PerformaxDMXJSASStage method), 499	<code>get_all_axes()</code> (py- lablib.devices.Arcus.performax.PerformaxDMXJSASStage method), 499
<code>get_all_attribute_values()</code> lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusBitFluxCamera method), 696	<code>get_all_axes()</code> lablib.devices.Arcus.performax.PerformaxDMXJSASStage method), 499	<code>get_all_axes()</code> (py- lablib.devices.Arcus.performax.PerformaxDMXJSASStage method), 499
<code>get_all_attribute_values()</code> lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusIMAQdxCamera method), 680	<code>get_all_axes()</code> lablib.devices.Arcus.performax.PerformaxDMXJSASStage method), 499	<code>get_all_axes()</code> (py- lablib.devices.Arcus.performax.PerformaxDMXJSASStage method), 499
<code>get_all_attribute_values()</code> lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusSiSoCamera method), 688	<code>get_all_axes()</code> lablib.devices.Arcus.performax.PerformaxDMXJSASStage method), 499	<code>get_all_axes()</code> (py- lablib.devices.Arcus.performax.PerformaxDMXJSASStage method), 499
<code>get_all_attribute_values()</code> lablib.devices.PrincetonInstruments.picam.PicamCamera method), 712	<code>get_all_axes()</code> lablib.devices.Arcus.performax.PerformaxDMXJSASStage method), 499	<code>get_all_axes()</code> (py- lablib.devices.Arcus.performax.PerformaxDMXJSASStage method), 499
<code>get_all_attributes()</code> lablib.devices.Andor.AndorSDK3.AndorSDK3Camera method), 483	<code>get_all_axes()</code> lablib.devices.Arcus.performax.PerformaxDMXJSASStage method), 499	<code>get_all_axes()</code> (py- lablib.devices.Arcus.performax.PerformaxDMXJSASStage method), 499
<code>get_all_attributes()</code> lablib.devices.Basler.pylon.BaslerPylonCamera method), 517	<code>get_all_axes()</code> lablib.devices.Arcus.performax.PerformaxDMXJSASStage method), 499	<code>get_all_axes()</code> (py- lablib.devices.Arcus.performax.PerformaxDMXJSASStage method), 499
<code>get_all_attributes()</code> lablib.devices.DCAM.DCAM.DCAMCamera method), 553	<code>get_all_axes()</code> lablib.devices.Arcus.performax.PerformaxDMXJSASStage method), 499	<code>get_all_axes()</code> (py- lablib.devices.Arcus.performax.PerformaxDMXJSASStage method), 499
<code>get_all_attributes()</code> lablib.devices.IMAQdx.IMAQdx.EthernetIMAQdxCamera method), 585	<code>get_all_axes()</code> lablib.devices.Arcus.performax.PerformaxDMXJSASStage method), 499	<code>get_all_axes()</code> (py- lablib.devices.Arcus.performax.PerformaxDMXJSASStage method), 499
<code>get_all_attributes()</code> lablib.devices.IMAQdx.IMAQdx.IMAQdxCamera method), 580	<code>get_all_axes()</code> lablib.devices.Arcus.performax.PerformaxDMXJSASStage method), 499	<code>get_all_axes()</code> (py- lablib.devices.Arcus.performax.PerformaxDMXJSASStage method), 499
<code>get_all_attributes()</code> lablib.devices.interface.camera.IAttributeCamera method), 837	<code>get_all_axes()</code> lablib.devices.Arcus.performax.PerformaxDMXJSASStage method), 499	<code>get_all_axes()</code> (py- lablib.devices.Arcus.performax.PerformaxDMXJSASStage method), 499
<code>get_all_attributes()</code> lablib.devices.Photometrics.pvcam.PvcamCamera method), 665	<code>get_all_axes()</code> lablib.devices.Arcus.performax.PerformaxDMXJSASStage method), 499	<code>get_all_axes()</code> (py- lablib.devices.Arcus.performax.PerformaxDMXJSASStage method), 499
<code>get_all_attributes()</code> lablib.devices.PhotonFocus.PhotonFocus.IPhotonFocusCamera method), 674	<code>get_all_axes()</code> lablib.devices.Arcus.performax.PerformaxDMXJSASStage method), 499	<code>get_all_axes()</code> (py- lablib.devices.Arcus.performax.PerformaxDMXJSASStage method), 499
<code>get_all_attributes()</code> lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusBitFluxCamera method), 696	<code>get_all_axes()</code> lablib.devices.Arcus.performax.PerformaxDMXJSASStage method), 499	<code>get_all_axes()</code> (py- lablib.devices.Arcus.performax.PerformaxDMXJSASStage method), 499
<code>get_all_attributes()</code> lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusIMAQdxCamera method), 680	<code>get_all_axes()</code> lablib.devices.Arcus.performax.PerformaxDMXJSASStage method), 499	<code>get_all_axes()</code> (py- lablib.devices.Arcus.performax.PerformaxDMXJSASStage method), 499
<code>get_all_attributes()</code> lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusSiSoCamera method), 688	<code>get_all_axes()</code> lablib.devices.Arcus.performax.PerformaxDMXJSASStage method), 499	<code>get_all_axes()</code> (py- lablib.devices.Arcus.performax.PerformaxDMXJSASStage method), 499

get_all_channels()	(py- lablib.devices.Thorlabs.kinesis.KinesisPiezoMotor method), 803	get_all_indicators()	(py- lablib.core.gui.widgets.container.IQWidgetContainer method), 207
get_all_channels()	(py- lablib.devices.Thorlabs.kinesis.MFF method), 792	get_all_indicators()	(py- lablib.core.gui.widgets.container.QContainer method), 203
get_all_color_modes()	(py- lablib.devices.uc480.uc480.UC480Camera method), 867	get_all_indicators()	(py- lablib.core.gui.widgets.container.QDialogContainer method), 218
get_all_grabber_attribute_values()	(py- lablib.devices.AlliedVision.Bonito.BonitoIMAQCamera method), 457	get_all_indicators()	(py- lablib.core.gui.widgets.container.QFrameContainer method), 215
get_all_grabber_attribute_values()	(py- lablib.devices.IMAQ.IMAQ.IMAQCamera method), 570	get_all_indicators()	(py- lablib.core.gui.widgets.container.QGroupBoxContainer method), 222
get_all_grabber_attribute_values()	(py- lablib.devices.IMAQ.IMAQ.IMAQFrameGrabber method), 562	get_all_indicators()	(py- lablib.core.gui.widgets.container.QScrollAreaContainer method), 230
get_all_grabber_attribute_values()	(py- lablib.devices.interface.camera.IGrabberAttributeCamera method), 842	get_all_indicators()	(py- lablib.core.gui.widgets.container.QScrollAreaContainer.QContainer method), 226
get_all_grabber_attribute_values()	(py- lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusIMAQCamera method), 680	get_all_indicators()	(py- lablib.core.gui.widgets.container.QTabContainer method), 232
get_all_grabber_attribute_values()	(py- lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusSiSoCamera method), 688	get_all_indicators()	(py- lablib.core.gui.widgets.container.QWidgetContainer method), 211
get_all_grabber_attribute_values()	(py- lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareCamera method), 728	get_all_indicators()	(py- lablib.core.gui.widgets.param_table.ParamTable method), 248
get_all_grabber_attribute_values()	(py- lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareFrameGrabber method), 721	get_all_indicators()	(py- lablib.core.gui.widgets.param_table.StatusTable method), 258
get_all_grabber_attributes()	(py- lablib.devices.interface.camera.IGrabberAttributeCamera method), 842	get_all_layout_containers()	(in module py- lablib.core.gui.utils), 263
get_all_grabber_attributes()	(py- lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusSiSoCamera method), 688	get_all_local_addr()	(in module py- lablib.core.utils.net), 388
get_all_grabber_attributes()	(py- lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareCamera method), 728	get_all_readout_modes()	(py- lablib.devices.Photometrics.pvcam.PvcamCamera method), 663
get_all_grabber_attributes()	(py- lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareFrameGrabber method), 723	get_all_readout_speeds()	(py- lablib.devices.DCAM.DCAM.DCAMCamera method), 551
get_all_handles()	(py- lablib.devices.Basler.pylon.BaslerPylonCamera.BufferManager method), 516	get_all_relayboards()	(py- lablib.devices.Conrad.base.RelayBoard method), 534
get_all_indicators()	(py- lablib.core.gui.value_handling.GUIValues method), 280	get_all_remote_addr()	(in module py- lablib.core.utils.net), 388
get_all_indicators()	(py- lablib.core.gui.widgets.container.IQContainer method), 202	get_all_sensor_kinds()	(py- lablib.devices.Cryocon.base.Cryocon1x method), 536
		get_all_sensor_readings()	(py- lablib.devices.Cryocon.base.Cryocon1x method), 536

<code>get_all_sensor_readings()</code> <i>lablib.devices.Lakeshore.base.Lakeshore218</i> method), 594	(py-	<code>get_amplitude()</code> <i>lablib.devices.AWG.generic.GenericAWG</i> method), 403	(py-
<code>get_all_temperatures()</code> <i>lablib.devices.Cryocon.base.Cryocon1x</i> method), 536	(py-	<code>get_amplitude()</code> <i>lablib.devices.AWG.specific.Agilent33220A</i> method), 415	(py-
<code>get_all_temperatures()</code> <i>lablib.devices.Lakeshore.base.Lakeshore218</i> method), 594	(py-	<code>get_amplitude()</code> <i>lablib.devices.AWG.specific.Agilent33500</i> method), 409	(py-
<code>get_all_trigger_modes()</code> <i>lablib.devices.DCAM.DCAM.DCAMCamera</i> method), 551	(py-	<code>get_amplitude()</code> <i>lablib.devices.AWG.specific.InstekAFG2000</i> method), 427	(py-
<code>get_all_values()</code> <i>lablib.core.gui.value_handling.GUIValues</i> method), 280	(py-	<code>get_amplitude()</code> <i>lablib.devices.AWG.specific.InstekAFG2225</i> method), 420	(py-
<code>get_all_values()</code> <i>lablib.core.gui.widgets.container.IQContainer</i> method), 202	(py-	<code>get_amplitude()</code> <i>lablib.devices.AWG.specific.RigolDG1000</i> method), 445	(py-
<code>get_all_values()</code> <i>lablib.core.gui.widgets.container.IQWidgetContainer</i> method), 207	(py-	<code>get_amplitude()</code> <i>lablib.devices.AWG.specific.RSInstekAFG21000</i> method), 432	(py-
<code>get_all_values()</code> <i>lablib.core.gui.widgets.container.QContainer</i> method), 203	(py-	<code>get_amplitude()</code> <i>lablib.devices.AWG.specific.TektronixAFG1000</i> method), 439	(py-
<code>get_all_values()</code> <i>lablib.core.gui.widgets.container.QDialogContainer</i> method), 218	(py-	<code>get_analog_input()</code> <i>lablib.devices.Arcus.performax.Performax2EXStage</i> method), 494	(py-
<code>get_all_values()</code> <i>lablib.core.gui.widgets.container.QFrameContainer</i> method), 215	(py-	<code>get_analog_input()</code> <i>lablib.devices.Arcus.performax.Performax4EXStage</i> method), 491	(py-
<code>get_all_values()</code> <i>lablib.core.gui.widgets.container.QGroupBoxContainer</i> method), 223	(py-	<code>get_analog_output()</code> <i>lablib.devices.Lakeshore.base.Lakeshore218</i> method), 595	(py-
<code>get_all_values()</code> <i>lablib.core.gui.widgets.container.QScrollAreaContainer</i> method), 230	(py-	<code>get_analog_output()</code> <i>lablib.devices.Lakeshore.base.Lakeshore370</i> method), 600	(py-
<code>get_all_values()</code> <i>lablib.core.gui.widgets.container.QScrollAreaContainer.QComboBoxWidget</i> method), 226	(py-	<code>get_analog_output_settings()</code> <i>lablib.devices.Lakeshore.base.Lakeshore218</i> method), 594	(py-
<code>get_all_values()</code> <i>lablib.core.gui.widgets.container.QTabContainer</i> method), 232	(py-	<code>get_analog_output_settings()</code> <i>lablib.devices.Lakeshore.base.Lakeshore370</i> method), 599	(py-
<code>get_all_values()</code> <i>lablib.core.gui.widgets.container.QWidgetContainer</i> method), 211	(py-	<code>get_app()</code> (in module <i>pylablib.core.thread.threadprop</i> ), 320	(py-
<code>get_all_values()</code> <i>lablib.core.gui.widgets.param_table.ParamTable</i> method), 248	(py-	<code>get_appended()</code> <i>lablib.core.dataproc.table_wrap.Array1DWrapper</i> method), 127	(py-
<code>get_all_values()</code> <i>lablib.core.gui.widgets.param_table.StatusTable</i> method), 258	(py-	<code>get_appended()</code> <i>lablib.core.dataproc.table_wrap.Array2DWrapper.ColumnAccess</i> method), 130	(py-
<code>get_amp_mode()</code> <i>lablib.devices.Andor.AndorSDK2.AndorSDK2Camera</i> method), 468	(py-	<code>get_appended()</code> <i>lablib.core.dataproc.table_wrap.Array2DWrapper.RowAccessor</i> method), 130	(py-
		<code>get_appended()</code>	(py-



[lablib.core.dataproc.table\\_wrap.DataFrame2DWrapper.ColumnAccess](#) (py-  
 method), 132 [get\\_arg\\_type\(\)](#) (py-  
[lablib.devices.Cryomagnetics.base.LM500](#)  
 get\_appended() (py-[lablib.core.dataproc.table\\_wrap.DataFrame2DWrapper.RowsAccess](#) method), 541  
 method), 132 [get\\_arg\\_type\(\)](#) (py-  
 get\_appended() (py-[lablib.devices.Cryomagnetics.base.LM510](#)  
[lablib.core.dataproc.table\\_wrap.Series1DWrapper](#) static method), 544  
 method), 128 [get\\_arg\\_type\(\)](#) (py-  
 get\_applet\_info() (in module py-[lablib.devices.Lakeshore.base.Lakeshore218](#)  
[lablib.devices.SiliconSoftware.fgrab](#)), 719 static method), 595  
 get\_arg\_default() (py-[lablib.core.dataproc.callable.FunctionCallable](#)  
 method), 106 [get\\_arg\\_type\(\)](#) (py-  
[lablib.devices.Lakeshore.base.Lakeshore370](#)  
 static method), 600  
 get\_arg\_default() (py-[lablib.core.dataproc.callable.ICallable](#)  
 method), 104 [get\\_arg\\_type\(\)](#) (py-  
[lablib.devices.Sirah.Matisse.SirahMatisse](#)  
 static method), 740  
 get\_arg\_default() (py-[lablib.core.dataproc.callable.JoinedCallable](#)  
 method), 105 [get\\_arg\\_type\(\)](#) (py-  
[lablib.devices.Tektronix.base.DPO2000](#) static  
 method), 763  
 get\_arg\_default() (py-[lablib.core.dataproc.callable.MethodCallable](#)  
 method), 107 [get\\_arg\\_type\(\)](#) (py-  
[lablib.devices.Tektronix.base.ITektronixScope](#)  
 static method), 753  
 get\_arg\_default() (py-[lablib.core.dataproc.callable.MultiplexedCallable](#)  
 method), 105 [get\\_arg\\_type\(\)](#) (py-  
[lablib.devices.Tektronix.base.TDS2000](#) static  
 method), 756  
 get\_arg\_type() (py-[lablib.core.devio.SCPISCPIDevice](#) static  
 method), 139 [get\\_arg\\_type\(\)](#) (py-  
[lablib.devices.Thorlabs.serial.FW](#) static  
 method), 813  
 get\_arg\_type() (py-[lablib.devices.AWG.generic.GenericAWG](#)  
 static method), 406 [get\\_arg\\_type\(\)](#) (py-  
[lablib.devices.Thorlabs.serial.FWv1](#) static  
 method), 817  
 get\_arg\_type() (py-[lablib.devices.AWG.specific.Agilent33220A](#)  
 static method), 415 [get\\_arg\\_type\(\)](#) (py-  
[lablib.devices.Thorlabs.serial.MDT69xA](#)  
 static method), 820  
 get\_arg\_type() (py-[lablib.devices.AWG.specific.Agilent33500](#)  
 static method), 409 [get\\_arg\\_type\(\)](#) (py-  
[lablib.devices.Thorlabs.serial.ThorlabsSerialInterface](#)  
 static method), 810  
 get\_arg\_type() (py-[lablib.devices.AWG.specific.InstekAFG2000](#)  
 static method), 427 [get\\_attenuation\(\)](#) (py-  
[lablib.devices.OZOptics.base.DD100](#) method),  
 635  
 get\_arg\_type() (py-[lablib.devices.AWG.specific.InstekAFG2225](#)  
 static method), 421 [get\\_attribute\(\)](#) (py-  
[lablib.devices.Andor.AndorSDK3.AndorSDK3Camera](#)  
 method), 480  
 get\_arg\_type() (py-[lablib.devices.AWG.specific.RigolDG1000](#)  
 static method), 445 [get\\_attribute\(\)](#) (py-  
[lablib.devices.Basler.pylon.BaslerPylonCamera](#)  
 method), 517  
 get\_arg\_type() (py-[lablib.devices.AWG.specific.RSInstekAFG21000](#)  
 static method), 433 [get\\_attribute\(\)](#) (py-  
[lablib.devices.DCAM.DCAM.DCAMCamera](#)  
 method), 553  
 get\_arg\_type() (py-[lablib.devices.AWG.specific.TektronixAFG1000](#)  
 static method), 439 [get\\_attribute\(\)](#) (py-  
[lablib.devices.IMAQdx.IMAQdx.EthernetIMAQdxCamera](#)  
 method), 585  
 get\_arg\_type() (py-[lablib.devices.Cryocon.base.Cryocon1x](#) static  
 method), 537 [get\\_attribute\(\)](#) (py-  
[lablib.devices.IMAQdx.IMAQdx.IMAQdxCamera](#)

931

class method), 152  
 get\_backend\_name() (py-lablib.core.devio.comm\_backend.IDeviceCommBackend class method), 143  
 get\_backend\_name() (py-lablib.core.devio.comm\_backend.NetworkDeviceBackend class method), 154  
 get\_backend\_name() (py-lablib.core.devio.comm\_backend.PyUSBDeviceBackend class method), 157  
 get\_backend\_name() (py-lablib.core.devio.comm\_backend.RecordedDeviceBackend class method), 159  
 get\_backend\_name() (py-lablib.core.devio.comm\_backend.SerialDeviceBackend class method), 149  
 get\_backend\_name() (py-lablib.core.devio.comm\_backend.VisaDeviceBackend class method), 146  
 get\_baseline\_simple() (in module py-lablib.core.dataproc.feature), 108  
 get\_battery\_condition() (py-lablib.devices.Ophir.base.VegaPowerMeter method), 642  
 get\_baudrate() (py-lablib.devices.Arcus.performax.Performax2EXStage method), 494  
 get\_baudrate() (py-lablib.devices.Arcus.performax.Performax4EXStage method), 490  
 get\_baudrate() (py-lablib.devices.Ophir.base.VegaPowerMeter method), 642  
 get\_baudrate() (py-lablib.devices.PhotonFocus.PhotonFocus.IPhotonFocusCamera method), 672  
 get\_baudrate() (py-lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusBitFlowCamera method), 696  
 get\_baudrate() (py-lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusIMAQCamera method), 680  
 get\_baudrate() (py-lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusSiSoCamera method), 688  
 get\_best\_fit() (in module py-lablib.core.dataproc.fitting), 117  
 get\_binning() (py-lablib.devices.uc480.uc480.UC480Camera method), 869  
 get\_bit\_alignment() (py-lablib.devices.PCO.SC2.PCOSC2Camera method), 649  
 get\_black\_level\_offset() (py-lablib.devices.AlliedVision.Bonito.BonitoIMAQCamera method), 458  
 get\_black\_level\_offset() (py-lablib.devices.AlliedVision.Bonito.IBonitoCamera method), 452  
 get\_black\_level\_offset() (py-lablib.devices.PhotonFocus.PhotonFocus.IPhotonFocusCamera method), 674  
 get\_black\_level\_offset() (py-lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusBitFlowCamera method), 696  
 get\_black\_level\_offset() (py-lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusIMAQCamera method), 680  
 get\_black\_level\_offset() (py-lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusSiSoCamera method), 688  
 get\_board\_info() (in module py-lablib.devices.SiliconSoftware.fgrab), 718  
 get\_boards\_number() (in module py-lablib.devices.SiliconSoftware.fgrab), 718  
 get\_buffer() (py-lablib.devices.Basler.pylon.BaslerPylonCamera.BufferManager method), 515  
 get\_buffer\_ptr() (py-lablib.devices.PCO.SC2.PCOSC2Camera.BufferManager method), 647  
 get\_buffer\_size() (py-lablib.devices.Andor.AndorSDK2.AndorSDK2Camera method), 472  
 get\_buffer\_size() (py-lablib.devices.NI.daq.NIDAQ method), 624  
 get\_burst\_mode() (py-lablib.devices.AWG.generic.GenericAWG method), 405  
 get\_burst\_mode() (py-lablib.devices.AWG.specific.Agilent33220A method), 415  
 get\_burst\_mode() (py-lablib.devices.AWG.specific.Agilent33500 method), 409  
 get\_burst\_mode() (py-lablib.devices.AWG.specific.Agilent33500 method), 409  
 get\_burst\_mode() (py-lablib.devices.AWG.specific.InstekAFG2000 method), 427  
 get\_burst\_mode() (py-lablib.devices.AWG.specific.InstekAFG2225 method), 421  
 get\_burst\_mode() (py-lablib.devices.AWG.specific.RigolDG1000 method), 445  
 get\_burst\_mode() (py-lablib.devices.AWG.specific.RSInstekAFG21000 method), 433  
 get\_burst\_mode() (py-lablib.devices.AWG.specific.TektronixAFG1000 method), 433



- method*), 439
- `get_burst_ncycles()` (*py-lablib.devices.AWG.generic.GenericAWG method*), 405
- `get_burst_ncycles()` (*py-lablib.devices.AWG.specific.Agilent33220A method*), 415
- `get_burst_ncycles()` (*py-lablib.devices.AWG.specific.Agilent33500 method*), 409
- `get_burst_ncycles()` (*py-lablib.devices.AWG.specific.InstekAFG2000 method*), 427
- `get_burst_ncycles()` (*py-lablib.devices.AWG.specific.InstekAFG2225 method*), 421
- `get_burst_ncycles()` (*py-lablib.devices.AWG.specific.RigolDG1000 method*), 445
- `get_burst_ncycles()` (*py-lablib.devices.AWG.specific.RSInstekAFG21000 method*), 433
- `get_burst_ncycles()` (*py-lablib.devices.AWG.specific.TektronixAFG1000 method*), 439
- `get_calibration_factor()` (*py-lablib.devices.Pfeiffer.base.TPG260 method*), 655
- `get_camera_id()` (*py-lablib.devices.uc480.uc480.UC480Camera method*), 867
- `get_camera_status()` (*py-lablib.devices.PCO.SC2.PCOSC2Camera method*), 646
- `get_cameras_number()` (*in module py-lablib.devices.Andor.AndorSDK2*), 465
- `get_cameras_number()` (*in module py-lablib.devices.Andor.AndorSDK3*), 477
- `get_cameras_number()` (*in module py-lablib.devices.Basler.pylon*), 511
- `get_cameras_number()` (*in module py-lablib.devices.BitFlow.BitFlow*), 521
- `get_cameras_number()` (*in module py-lablib.devices.DCAM.DCAM*), 548
- `get_cameras_number()` (*in module py-lablib.devices.IMAQ.IMAQ*), 561
- `get_cameras_number()` (*in module py-lablib.devices.IMAQdx.IMAQdx*), 576
- `get_cameras_number()` (*in module py-lablib.devices.PCO.SC2*), 644
- `get_cameras_number()` (*in module py-lablib.devices.Photometrics.pvcam*), 659
- `get_cameras_number()` (*in module py-lablib.devices.PhotonFocus.PhotonFocus*), 670
- `get_cameras_number()` (*in module py-lablib.devices.PrincetonInstruments.picam*), 709
- `get_cameras_number()` (*in module py-lablib.devices.Thorlabs.TLCamera*), 769
- `get_cameras_number()` (*in module py-lablib.devices.uc480.uc480*), 864
- `get_camlink_pixel_format()` (*py-lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusSiSoCamera method*), 688
- `get_camlink_pixel_format()` (*py-lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareCamera method*), 728
- `get_camlink_pixel_format()` (*py-lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareFrameGrabber method*), 723
- `get_capabilities()` (*py-lablib.devices.Andor.AndorSDK2.AndorSDK2Camera method*), 467
- `get_capabilities()` (*py-lablib.devices.PCO.SC2.PCOSC2Camera method*), 646
- `get_capacitance()` (*py-lablib.devices.Attocube.anc300.ANC300 method*), 504
- `get_capacitance()` (*py-lablib.devices.Attocube.anc350.ANC350 method*), 508
- `get_channel()` (*py-lablib.devices.Andor.AndorSDK2.AndorSDK2Camera method*), 468
- `get_channel()` (*py-lablib.devices.Cryomagnetics.base.LM500 method*), 540
- `get_channel()` (*py-lablib.devices.Cryomagnetics.base.LM510 method*), 544
- `get_channel()` (*py-lablib.devices.Lakeshore.base.Lakeshore370 method*), 599
- `get_channel_bitdepth()` (*py-lablib.devices.Andor.AndorSDK2.AndorSDK2Camera method*), 468
- `get_channel_power()` (*py-lablib.devices.Toptica.ibeam.TopticaIBeam method*), 825
- `get_channel_range_settings()` (*py-lablib.devices.Lakeshore.base.Lakeshore370 method*), 599
- `get_channel_status()` (*py-lablib.devices.Pfeiffer.base.TPG260 method*), 655
- `get_channels()`

`lablib.devices.Tektronix.base.DPO2000`  
`method)`, 763  
`get_channels()` (py-  
`lablib.devices.Tektronix.base.ITektronixScope`  
`method)`, 749  
`get_channels()` (py-  
`lablib.devices.Tektronix.base.TDS2000`  
`method)`, 756  
`get_channels_number()` (py-  
`lablib.devices.AWG.generic.GenericAWG`  
`method)`, 403  
`get_channels_number()` (py-  
`lablib.devices.AWG.specific.Agilent33220A`  
`method)`, 415  
`get_channels_number()` (py-  
`lablib.devices.AWG.specific.Agilent33500`  
`method)`, 409  
`get_channels_number()` (py-  
`lablib.devices.AWG.specific.InstekAFG2000`  
`method)`, 427  
`get_channels_number()` (py-  
`lablib.devices.AWG.specific.InstekAFG2225`  
`method)`, 421  
`get_channels_number()` (py-  
`lablib.devices.AWG.specific.RigolDG1000`  
`method)`, 445  
`get_channels_number()` (py-  
`lablib.devices.AWG.specific.RSInstekAFG21000`  
`method)`, 433  
`get_channels_number()` (py-  
`lablib.devices.AWG.specific.TektronixAFG1000`  
`method)`, 439  
`get_channels_number()` (py-  
`lablib.devices.HighFinesse.wlm.WLM` `method)`,  
558  
`get_channels_number()` (py-  
`lablib.devices.Tektronix.base.DPO2000`  
`method)`, 763  
`get_channels_number()` (py-  
`lablib.devices.Tektronix.base.ITektronixScope`  
`method)`, 749  
`get_channels_number()` (py-  
`lablib.devices.Tektronix.base.TDS2000`  
`method)`, 756  
`get_channels_number()` (py-  
`lablib.devices.Toptica.ibeam.TopicalBeam`  
`method)`, 824  
`get_child()` (pylablib.core.gui.widgets.container.IQContainer  
`method)`, 201  
`get_child()` (pylablib.core.gui.widgets.container.IQWidgetContainer  
`method)`, 207  
`get_child()` (pylablib.core.gui.widgets.container.QContainer  
`method)`, 203  
`get_child()` (pylablib.core.gui.widgets.container.QDialogContainer  
`method)`, 219  
`get_child()` (pylablib.core.gui.widgets.container.QFrameContainer  
`method)`, 215  
`get_child()` (pylablib.core.gui.widgets.container.QGroupBoxContainer  
`method)`, 223  
`get_child()` (pylablib.core.gui.widgets.container.QScrollAreaContainer  
`method)`, 230  
`get_child()` (pylablib.core.gui.widgets.container.QScrollAreaContainer  
`method)`, 227  
`get_child()` (pylablib.core.gui.widgets.container.QTabContainer  
`method)`, 232  
`get_child()` (pylablib.core.gui.widgets.container.QWidgetContainer  
`method)`, 211  
`get_child()` (pylablib.core.gui.widgets.param\_table.ParamTable  
`method)`, 248  
`get_child()` (pylablib.core.gui.widgets.param\_table.StatusTable  
`method)`, 258  
`get_clear_cycles()` (py-  
`lablib.devices.Photometrics.pvcam.PvcamCamera`  
`method)`, 664  
`get_clear_mode()` (py-  
`lablib.devices.Photometrics.pvcam.PvcamCamera`  
`method)`, 664  
`get_clearing_time()` (py-  
`lablib.devices.Photometrics.pvcam.PvcamCamera`  
`method)`, 664  
`get_clock_parameters()` (py-  
`lablib.devices.NI.daq.NIDAQ` `method)`, 622  
`get_clock_period_input_parameters()` (py-  
`lablib.devices.NI.daq.NIDAQ` `method)`, 623  
`get_coarse_tuning_status()` (py-  
`lablib.devices.M2.solstis.Solstis` `method)`,  
616  
`get_coarse_wavelength()` (py-  
`lablib.devices.M2.solstis.Solstis` `method)`,  
616  
`get_color_mode()` (py-  
`lablib.devices.uc480.uc480.UC480Camera`  
`method)`, 867  
`get_column_index()` (py-  
`lablib.core.dataproc.table_wrap.Array2DWrapper.ColumnAccess`  
`method)`, 130  
`get_column_index()` (py-  
`lablib.core.dataproc.table_wrap.DataFrame2DWrapper.ColumnA`  
`method)`, 132  
`get_columns_line()` (py-  
`lablib.core.fileio.savefile.CSVTableOutputFileFormat`  
`method)`, 191  
`get_config()` (pylablib.devices.OZOptics.base.DD100  
`method)`, 635  
`get_config()` (pylablib.devices.OZOptics.base.OZOpticsDevice  
`method)`, 632  
`get_config()` (pylablib.devices.OZOptics.base.TF100  
`method)`, 634

<code>get_connected_addrs()</code> ( <i>pylablib.devices.Thorlabs.elliptec.ElliptecMotor</i> method), 778	<code>get_current_channel()</code> ( <i>pylablib.devices.AWG.specific.RSInstekAFG21000</i> method), 433
<code>get_cont_mode_parameters()</code> ( <i>pylablib.devices.Andor.AndorSDK2.AndorSDK2Camera</i> method), 470	<code>get_current_channel()</code> ( <i>pylablib.devices.AWG.specific.TektronixAFG1000</i> method), 439
<code>get_controller()</code> (in module <i>pylablib.core.thread.controller</i> ), 314	<code>get_current_errors()</code> ( <i>pylablib.devices.Pfeiffer.base.TPG260</i> method), 656
<code>get_conversion_factor()</code> ( <i>pylablib.devices.PCO.SC2.PCOS2Camera</i> method), 646	<code>get_current_len()</code> ( <i>pylablib.core.thread.callsync.QQueueLengthLimitScheduler</i> method), 287
<code>get_correlations_ft()</code> (in module <i>pylablib.core.dataproc.fourier</i> ), 120	<code>get_current_limits()</code> ( <i>pylablib.devices.Toptica.ibeam.TopticaIBeam</i> method), 825
<code>get_counter_input_parameters()</code> ( <i>pylablib.devices.NI.daq.NIDAQ</i> method), 623	<code>get_current_name()</code> ( <i>pylablib.core.gui.widgets.container.QTabContainer</i> method), 231
<code>get_coupling()</code> ( <i>pylablib.devices.Tektronix.base.DPO2000</i> method), 763	<code>get_current_parameters()</code> ( <i>pylablib.devices.Trinamic.base.TMCM1110</i> method), 828
<code>get_coupling()</code> ( <i>pylablib.devices.Tektronix.base.ITektronixScope</i> method), 751	<code>get_current_size()</code> ( <i>pylablib.core.thread.callsync.QQueueSizeLimitScheduler</i> method), 289
<code>get_coupling()</code> ( <i>pylablib.devices.Tektronix.base.TDS2000</i> method), 756	<code>get_current_speed()</code> ( <i>pylablib.devices.Trinamic.base.TMCM1110</i> method), 829
<code>get_ctypes_frames_list()</code> ( <i>pylablib.devices.interface.camera.ChunkBufferManager</i> method), 836	<code>get_cursor_order()</code> ( <i>pylablib.core.gui.widgets.edit.NumEdit</i> method), 235
<code>get_current()</code> ( <i>pylablib.devices.LaserQuantum.base.Finesse</i> method), 605	<code>get_curve()</code> ( <i>pylablib.devices.Lakeshore.base.Lakeshore218</i> method), 594
<code>get_current_axis_speed()</code> ( <i>pylablib.devices.Arcus.performax.Performax2EXStage</i> method), 494	<code>get_curve_header()</code> ( <i>pylablib.devices.Lakeshore.base.Lakeshore218</i> method), 594
<code>get_current_axis_speed()</code> ( <i>pylablib.devices.Arcus.performax.Performax4EXStage</i> method), 491	<code>get_cycle_timings()</code> ( <i>pylablib.devices.Andor.AndorSDK2.AndorSDK2Camera</i> method), 471
<code>get_current_channel()</code> ( <i>pylablib.devices.AWG.generic.GenericAWG</i> method), 403	<code>get_data_dimensions()</code> ( <i>pylablib.devices.AlliedVision.Bonito.BonitoIMACamera</i> method), 458
<code>get_current_channel()</code> ( <i>pylablib.devices.AWG.specific.Agilent33220A</i> method), 415	<code>get_data_dimensions()</code> ( <i>pylablib.devices.AlliedVision.Bonito.IBonitoCamera</i> method), 453
<code>get_current_channel()</code> ( <i>pylablib.devices.AWG.specific.Agilent33500</i> method), 410	<code>get_data_dimensions()</code> ( <i>pylablib.devices.Andor.AndorSDK2.AndorSDK2Camera</i> method), 473
<code>get_current_channel()</code> ( <i>pylablib.devices.AWG.specific.InstekAFG2000</i> method), 427	<code>get_data_dimensions()</code> ( <i>pylablib.devices.Andor.AndorSDK3.AndorSDK3Camera</i> method), 483
<code>get_current_channel()</code> ( <i>pylablib.devices.AWG.specific.InstekAFG2225</i> method), 422	<code>get_data_dimensions()</code> ( <i>pylablib.devices.Basler.pylon.BaslerPylonCamera</i> method), 517
<code>get_current_channel()</code> ( <i>pylablib.devices.AWG.specific.RigolDG1000</i> method), 445	<code>get_data_dimensions()</code> ( <i>py-</i>

<i>lablib.devices.BitFlow.BitFlow.BitFlowCamera</i> <i>method</i> ), 528	<i>lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusSiSoCamera</i> <i>method</i> ), 688
<i>get_data_dimensions()</i> <i>lablib.devices.BitFlow.BitFlow.BitFlowFrameGrabber</i> <i>method</i> ), 523	<i>get_data_dimensions()</i> <i>lablib.devices.PrincetonInstruments.picam.PicamCamera</i> <i>method</i> ), 714
<i>get_data_dimensions()</i> <i>lablib.devices.DCAM.DCAM.DCAMCamera</i> <i>method</i> ), 553	<i>get_data_dimensions()</i> <i>lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareCamera</i> <i>method</i> ), 728
<i>get_data_dimensions()</i> <i>lablib.devices.IMAQ.IMAQ.IMAQCamera</i> <i>method</i> ), 570	<i>get_data_dimensions()</i> <i>lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareFrameGrabber</i> <i>method</i> ), 723
<i>get_data_dimensions()</i> <i>lablib.devices.IMAQ.IMAQ.IMAQFrameGrabber</i> <i>method</i> ), 565	<i>get_data_dimensions()</i> <i>lablib.devices.Thorlabs.TLCamera.ThorlabsTLCamera</i> <i>method</i> ), 773
<i>get_data_dimensions()</i> <i>lablib.devices.IMAQdx.IMAQdx.EthernetIMAQdxCamera</i> <i>method</i> ), 585	<i>get_data_dimensions()</i> <i>lablib.devices.uc480.uc480.UC480Camera</i> <i>method</i> ), 869
<i>get_data_dimensions()</i> <i>lablib.devices.IMAQdx.IMAQdx.IMAQdxCamera</i> <i>method</i> ), 580	<i>get_data_format()</i> <i>lablib.devices.Tektronix.base.DPO2000</i> <i>method</i> ), 763
<i>get_data_dimensions()</i> <i>lablib.devices.interface.camera.IAttributeCamera</i> <i>method</i> ), 838	<i>get_data_format()</i> <i>lablib.devices.Tektronix.base.ITektronixScope</i> <i>method</i> ), 751
<i>get_data_dimensions()</i> <i>lablib.devices.interface.camera.IBinROICamera</i> <i>method</i> ), 857	<i>get_data_format()</i> <i>lablib.devices.Tektronix.base.TDS2000</i> <i>method</i> ), 756
<i>get_data_dimensions()</i> <i>lablib.devices.interface.camera.ICamera</i> <i>method</i> ), 832	<i>get_data_pts_range()</i> <i>lablib.devices.Tektronix.base.DPO2000</i> <i>method</i> ), 763
<i>get_data_dimensions()</i> <i>lablib.devices.interface.camera.IExposureCamera</i> <i>method</i> ), 847	<i>get_data_pts_range()</i> <i>lablib.devices.Tektronix.base.ITektronixScope</i> <i>method</i> ), 751
<i>get_data_dimensions()</i> <i>lablib.devices.interface.camera.IGrabberAttributeCamera</i> <i>method</i> ), 843	<i>get_data_pts_range()</i> <i>lablib.devices.Tektronix.base.TDS2000</i> <i>method</i> ), 756
<i>get_data_dimensions()</i> <i>lablib.devices.interface.camera.IROICamera</i> <i>method</i> ), 852	<i>get_default_addr()</i> <i>lablib.devices.Thorlabs.elliptec.ElliptecMotor</i> <i>method</i> ), 778
<i>get_data_dimensions()</i> <i>lablib.devices.PCO.SC2.PCOSC2Camera</i> <i>method</i> ), 649	<i>get_default_channel()</i> <i>lablib.devices.HighFinesse.wlm.WLM</i> <i>method</i> ), 558
<i>get_data_dimensions()</i> <i>lablib.devices.Photometrics.pvcam.PvcamCamera</i> <i>method</i> ), 666	<i>get_defaults_list()</i> <i>lablib.core.utils.functions.FunctionSignature</i> <i>method</i> ), 370
<i>get_data_dimensions()</i> <i>lablib.devices.PhotonFocus.PhotonFocus.IPhotonFocusCamera</i> <i>method</i> ), 674	<i>get_defect_correct_mode()</i> <i>lablib.devices.DCAM.DCAM.DCAMCamera</i> <i>method</i> ), 552
<i>get_data_dimensions()</i> <i>lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusBitFlowCamera</i> <i>method</i> ), 696	<i>get_deleted()</i> <i>lablib.devices.dataproc.table_wrap.Array1DWrapper</i> <i>method</i> ), 126
<i>get_data_dimensions()</i> <i>lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusIMAQdxCamera</i> <i>method</i> ), 680	<i>get_deleted()</i> <i>lablib.devices.dataproc.table_wrap.Array2DWrapper.ColumnAccess</i> <i>method</i> ), 130
<i>get_data_dimensions()</i>	<i>get_deleted()</i>

<code>lablib.core.dataproc.table_wrap.Array2DWrapper.RowAccess</code> <code>method), 130</code>	<code>lablib.devices.interface.camera.ICamera</code> <code>method), 832</code>
<code>get_deleted()</code> <code>lablib.core.dataproc.table_wrap.DataFrame2DWrapper.Collablib.devices.interface.camera.IExposureCamera</code> <code>method), 132</code>	<code>(py- get_detector_size()</code> <code>method), 847</code>
<code>get_deleted()</code> <code>lablib.core.dataproc.table_wrap.DataFrame2DWrapper.RowAccess</code> <code>method), 132</code>	<code>(py- get_detector_size()</code> <code>method), 843</code>
<code>get_deleted()</code> <code>lablib.core.dataproc.table_wrap.Series1DWrapper</code> <code>method), 128</code>	<code>(py- get_detector_size()</code> <code>lablib.devices.interface.camera.IGrabberAttributeCamera</code> <code>method), 852</code>
<code>get_detector_size()</code> <code>lablib.devices.AlliedVision.Bonito.BonitoIMAQCamera</code> <code>method), 458</code>	<code>(py- get_detector_size()</code> <code>lablib.devices.interface.camera.IROICamera</code> <code>method), 852</code>
<code>get_detector_size()</code> <code>lablib.devices.AlliedVision.Bonito.IBonitoCamera</code> <code>method), 451</code>	<code>(py- get_detector_size()</code> <code>lablib.devices.PCO.SC2.PCOSC2Camera</code> <code>method), 648</code>
<code>get_detector_size()</code> <code>lablib.devices.Andor.AndorSDK2.AndorSDK2Camera</code> <code>method), 472</code>	<code>(py- get_detector_size()</code> <code>lablib.devices.Photometrics.pvcam.PvcamCamera</code> <code>method), 664</code>
<code>get_detector_size()</code> <code>lablib.devices.Andor.AndorSDK3.AndorSDK3Camera</code> <code>method), 482</code>	<code>(py- get_detector_size()</code> <code>lablib.devices.PhotonFocus.PhotonFocus.IPhotonFocusCamera</code> <code>method), 673</code>
<code>get_detector_size()</code> <code>lablib.devices.Basler.pylon.BaslerPylonCamera</code> <code>method), 515</code>	<code>(py- get_detector_size()</code> <code>lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusBitFlowCa</code> <code>method), 696</code>
<code>get_detector_size()</code> <code>lablib.devices.BitFlow.BitFlow.BitFlowCamera</code> <code>method), 528</code>	<code>(py- get_detector_size()</code> <code>lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusIMAQCam</code> <code>method), 680</code>
<code>get_detector_size()</code> <code>lablib.devices.BitFlow.BitFlow.BitFlowFrameGrabber</code> <code>method), 522</code>	<code>(py- get_detector_size()</code> <code>lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusSiSoCamer</code> <code>method), 688</code>
<code>get_detector_size()</code> <code>lablib.devices.DCAM.DCAM.DCAMCamera</code> <code>method), 552</code>	<code>(py- get_detector_size()</code> <code>lablib.devices.PrincetonInstruments.picam.PicamCamera</code> <code>method), 713</code>
<code>get_detector_size()</code> <code>lablib.devices.IMAQ.IMAQ.IMAQCamera</code> <code>method), 570</code>	<code>(py- get_detector_size()</code> <code>lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareCamera</code> <code>method), 728</code>
<code>get_detector_size()</code> <code>lablib.devices.IMAQ.IMAQ.IMAQFrameGrabber</code> <code>method), 562</code>	<code>(py- get_detector_size()</code> <code>lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareFrameGrabb</code> <code>method), 722</code>
<code>get_detector_size()</code> <code>lablib.devices.IMAQdx.IMAQdx.EthernetIMAQdxCamera</code> <code>method), 585</code>	<code>(py- get_detector_size()</code> <code>lablib.devices.Thorlabs.TLCamera.ThorlabsTLCamera</code> <code>method), 772</code>
<code>get_detector_size()</code> <code>lablib.devices.IMAQdx.IMAQdx.IMAQdxCamera</code> <code>method), 579</code>	<code>(py- get_detector_size()</code> <code>lablib.devices.uc480.uc480.UC480Camera</code> <code>method), 869</code>
<code>get_detector_size()</code> <code>lablib.devices.interface.camera.IAttributeCamera</code> <code>method), 838</code>	<code>(py- get_device()</code> ( <code>pylablib.core.utils.rpyc_utils.DeviceService</code> <code>method), 395</code>
<code>get_detector_size()</code> <code>lablib.devices.interface.camera.IBinROICamera</code> <code>method), 857</code>	<code>get_device_class()</code> ( <code>py-</code> <code>method), 395</code>
<code>get_detector_size()</code>	<code>get_device_info()</code> ( <code>in module py-</code> <code>method), 511</code>
<code>get_detector_size()</code>	<code>get_device_info()</code> ( <code>in module py-</code> <code>method), 621</code>
<code>get_detector_size()</code>	<code>get_device_info()</code> ( <code>in module py-</code>



[lablib.devices.SmarAct.scu3d](#)), 745  
 get\_device\_info() (py- [lablib.devices.AlliedVision.Bonito.BonitoIMAQCamera](#) method), 458  
 get\_device\_info() (py- [lablib.devices.AlliedVision.Bonito.IBonitoCamera](#) method), 451  
 get\_device\_info() (py- [lablib.devices.Andor.AndorSDK2.AndorSDK2Camera](#) method), 467  
 get\_device\_info() (py- [lablib.devices.Andor.AndorSDK3.AndorSDK3Camera](#) method), 480  
 get\_device\_info() (py- [lablib.devices.Arcus.performax.GenericPerformaxStage](#) method), 488  
 get\_device\_info() (py- [lablib.devices.Arcus.performax.Performax2EXStage](#) method), 494  
 get\_device\_info() (py- [lablib.devices.Arcus.performax.Performax4EXStage](#) method), 492  
 get\_device\_info() (py- [lablib.devices.Arcus.performax.PerformaxDMXJSAStage](#) method), 499  
 get\_device\_info() (py- [lablib.devices.Attocube.anc300.ANC300](#) method), 503  
 get\_device\_info() (py- [lablib.devices.Basler.pylon.BaslerPylonCamera](#) method), 515  
 get\_device\_info() (py- [lablib.devices.BitFlow.BitFlow.BitFlowCamera](#) method), 528  
 get\_device\_info() (py- [lablib.devices.BitFlow.BitFlow.BitFlowFrameGrabber](#) method), 522  
 get\_device\_info() (py- [lablib.devices.DCAM.DCAM.DCAMCamera](#) method), 551  
 get\_device\_info() (py- [lablib.devices.HighFinesse.wlm.WLM](#) method), 557  
 get\_device\_info() (py- [lablib.devices.IMAQ.IMAQ.IMAQCamera](#) method), 570  
 get\_device\_info() (py- [lablib.devices.IMAQ.IMAQ.IMAQFrameGrabber](#) method), 562  
 get\_device\_info() (py- [lablib.devices.IMAQdx.IMAQdx.EthernetIMAQdxCamera](#) method), 585  
 get\_device\_info() (py- [lablib.devices.IMAQdx.IMAQdx.IMAQdxCamera](#) method), 746  
 method), 579  
 get\_device\_info() (py- [lablib.devices.KJL.base.KJL300](#) method), 590  
 get\_device\_info() (py- [lablib.devices.LaserQuantum.base.Finesse](#) method), 604  
 get\_device\_info() (py- [lablib.devices.Leybold.base.GenericITR](#) method), 607  
 get\_device\_info() (py- [lablib.devices.Leybold.base.ITR90](#) method), 609  
 get\_device\_info() (py- [lablib.devices.LighthousePhotonics.base.SproutG](#) method), 612  
 get\_device\_info() (py- [lablib.devices.Newport.picomotor.Picomotor8742](#) method), 628  
 get\_device\_info() (py- [lablib.devices.NI.daq.NIDAQ](#) method), 622  
 get\_device\_info() (py- [lablib.devices.Ophir.base.VegaPowerMeter](#) method), 641  
 get\_device\_info() (py- [lablib.devices.PCO.SC2.PCOSC2Camera](#) method), 646  
 get\_device\_info() (py- [lablib.devices.Photometrics.pvcam.PvcamCamera](#) method), 663  
 get\_device\_info() (py- [lablib.devices.PhotonFocus.PhotonFocus.IPhotonFocusCamera](#) method), 673  
 get\_device\_info() (py- [lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusBitFlowCa](#) method), 696  
 get\_device\_info() (py- [lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusIMAQCam](#) method), 680  
 get\_device\_info() (py- [lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusSiSoCamer](#) method), 688  
 get\_device\_info() (py- [lablib.devices.PrincetonInstruments.picam.PicamCamera](#) method), 713  
 get\_device\_info() (py- [lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareCamera](#) method), 728  
 get\_device\_info() (py- [lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareFrameGrabb](#) method), 722  
 get\_device\_info() (py- [lablib.devices.SmarAct.scu3d.SCU3D](#) method), 746

<code>get_device_info()</code> <i>lablib.devices.Thorlabs.elliptec.ElliptecMotor</i> <i>method</i> ), 779	(py-	<code>get_device_variable()</code> <i>lablib.devices.AlliedVision.Bonito.BonitoIMAQCamera</i> <i>method</i> ), 458	(py-
<code>get_device_info()</code> <i>lablib.devices.Thorlabs.kinesis.BasicKinesisDevice</i> <i>method</i> ), 783	(py-	<code>get_device_variable()</code> <i>lablib.devices.AlliedVision.Bonito.IBonitoCamera</i> <i>method</i> ), 453	(py-
<code>get_device_info()</code> <i>lablib.devices.Thorlabs.kinesis.KinesisDevice</i> <i>method</i> ), 789	(py-	<code>get_device_variable()</code> <i>lablib.devices.Andor.AndorSDK2.AndorSDK2Camera</i> <i>method</i> ), 473	(py-
<code>get_device_info()</code> <i>lablib.devices.Thorlabs.kinesis.KinesisMotor</i> <i>method</i> ), 799	(py-	<code>get_device_variable()</code> <i>lablib.devices.Andor.AndorSDK3.AndorSDK3Camera</i> <i>method</i> ), 483	(py-
<code>get_device_info()</code> <i>lablib.devices.Thorlabs.kinesis.KinesisPiezoMotor</i> <i>method</i> ), 803	(py-	<code>get_device_variable()</code> <i>lablib.devices.Arcus.performax.GenericPerformaxStage</i> <i>method</i> ), 489	(py-
<code>get_device_info()</code> <i>lablib.devices.Thorlabs.kinesis.KinesisQuadDetector</i> <i>method</i> ), 808	(py-	<code>get_device_variable()</code> <i>lablib.devices.Arcus.performax.Performax2EXStage</i> <i>method</i> ), 494	(py-
<code>get_device_info()</code> <i>lablib.devices.Thorlabs.kinesis.MFF</i> <i>method</i> ), 793	(py-	<code>get_device_variable()</code> <i>lablib.devices.Arcus.performax.Performax4EXStage</i> <i>method</i> ), 492	(py-
<code>get_device_info()</code> <i>lablib.devices.Thorlabs.TLCamera.ThorlabsTLCamera</i> <i>method</i> ), 770	(py-	<code>get_device_variable()</code> <i>lablib.devices.Arcus.performax.PerformaxDMXJSAStage</i> <i>method</i> ), 499	(py-
<code>get_device_info()</code> <i>lablib.devices.Topica.ibeam.TopicalBeam</i> <i>method</i> ), 824	(py-	<code>get_device_variable()</code> <i>lablib.devices.Arduino.base.IArduinoDevice</i> <i>method</i> ), 501	(py-
<code>get_device_info()</code> <i>lablib.devices.uc480.uc480.UC480Camera</i> <i>method</i> ), 867	(py-	<code>get_device_variable()</code> <i>lablib.devices.Attocube.anc300.ANC300</i> <i>method</i> ), 505	(py-
<code>get_device_name()</code> <i>lablib.devices.Pfeiffer.base.DPG202</i> <i>method</i> ), 658	(py-	<code>get_device_variable()</code> <i>lablib.devices.Attocube.anc350.ANC350</i> <i>method</i> ), 509	(py-
<code>get_device_number()</code> <i>lablib.devices.Arcus.performax.GenericPerformaxStage</i> <i>method</i> ), 488	(py-	<code>get_device_variable()</code> <i>lablib.devices.AWG.generic.GenericAWG</i> <i>method</i> ), 406	(py-
<code>get_device_number()</code> <i>lablib.devices.Arcus.performax.Performax2EXStage</i> <i>method</i> ), 494	(py-	<code>get_device_variable()</code> <i>lablib.devices.AWG.specific.Agilent33220A</i> <i>method</i> ), 416	(py-
<code>get_device_number()</code> <i>lablib.devices.Arcus.performax.Performax4EXStage</i> <i>method</i> ), 492	(py-	<code>get_device_variable()</code> <i>lablib.devices.AWG.specific.Agilent33500</i> <i>method</i> ), 410	(py-
<code>get_device_number()</code> <i>lablib.devices.Arcus.performax.PerformaxDMXJSAStage</i> <i>method</i> ), 499	(py-	<code>get_device_variable()</code> <i>lablib.devices.AWG.specific.InstekAFG2000</i> <i>method</i> ), 427	(py-
<code>get_device_variable()</code> <i>lablib.core.devio.comm_backend.ICommBackendWrapper</i> <i>method</i> ), 160	(py-	<code>get_device_variable()</code> <i>lablib.devices.AWG.specific.InstekAFG2225</i> <i>method</i> ), 422	(py-
<code>get_device_variable()</code> <i>lablib.core.devio.interface.IDevice</i> <i>method</i> ), 163	(py-	<code>get_device_variable()</code> <i>lablib.devices.AWG.specific.RigolDG1000</i> <i>method</i> ), 445	(py-
<code>get_device_variable()</code> <i>lablib.core.devio.SCPi.SCPIDevice</i> <i>method</i> ), 140	(py-	<code>get_device_variable()</code> <i>lablib.devices.AWG.specific.RSInstekAFG21000</i> <i>method</i> ), 433	(py-

<code>get_device_variable()</code> (py-lablib.devices.AWG.specific.TektronixAFG1000 method), 439	<code>get_device_variable()</code> (py-lablib.devices.interface.camera.IGrabberAttributeCamera method), 843
<code>get_device_variable()</code> (py-lablib.devices.Basler.pylon.BaslerPylonCamera method), 517	<code>get_device_variable()</code> (py-lablib.devices.interface.camera.IROICamera method), 852
<code>get_device_variable()</code> (py-lablib.devices.BitFlow.BitFlow.BitFlowCamera method), 528	<code>get_device_variable()</code> (py-lablib.devices.interface.stage.IMultiaxisStage method), 863
<code>get_device_variable()</code> (py-lablib.devices.BitFlow.BitFlow.BitFlowFrameGrabber method), 524	<code>get_device_variable()</code> (py-lablib.devices.interface.stage.IStage method), 862
<code>get_device_variable()</code> (py-lablib.devices.Conrad.base.RelayBoard method), 534	<code>get_device_variable()</code> (py-lablib.devices.KJL.base.KJL300 method), 591
<code>get_device_variable()</code> (py-lablib.devices.Cryocon.base.Cryocon1x method), 537	<code>get_device_variable()</code> (py-lablib.devices.Lakeshore.base.Lakeshore218 method), 595
<code>get_device_variable()</code> (py-lablib.devices.Cryomagnetics.base.LM500 method), 541	<code>get_device_variable()</code> (py-lablib.devices.Lakeshore.base.Lakeshore370 method), 600
<code>get_device_variable()</code> (py-lablib.devices.Cryomagnetics.base.LM510 method), 544	<code>get_device_variable()</code> (py-lablib.devices.LaserQuantum.base.Finesse method), 605
<code>get_device_variable()</code> (py-lablib.devices.DCAM.DCAM.DCAMCamera method), 553	<code>get_device_variable()</code> (py-lablib.devices.Leybold.base.GenericITR method), 608
<code>get_device_variable()</code> (py-lablib.devices.HighFinesse.wlm.WLM method), 560	<code>get_device_variable()</code> (py-lablib.devices.Leybold.base.ITR90 method), 609
<code>get_device_variable()</code> (py-lablib.devices.IMAQ.IMAQ.IMAQCamera method), 570	<code>get_device_variable()</code> (py-lablib.devices.LighthousePhotonics.base.SproutG method), 613
<code>get_device_variable()</code> (py-lablib.devices.IMAQ.IMAQ.IMAQFrameGrabber method), 565	<code>get_device_variable()</code> (py-lablib.devices.M2.solstis.Solstis method), 619
<code>get_device_variable()</code> (py-lablib.devices.IMAQdx.IMAQdx.EthernetIMAQdxCamera method), 585	<code>get_device_variable()</code> (py-lablib.devices.Newport.picomotor.Picomotor8742 method), 630
<code>get_device_variable()</code> (py-lablib.devices.IMAQdx.IMAQdx.IMAQdxCamera method), 580	<code>get_device_variable()</code> (py-lablib.devices.NI.daq.NIDAQ method), 626
<code>get_device_variable()</code> (py-lablib.devices.interface.camera.IAttributeCamera method), 838	<code>get_device_variable()</code> (py-lablib.devices.Ophir.base.OphirDevice method), 639
<code>get_device_variable()</code> (py-lablib.devices.interface.camera.IBinROICamera method), 857	<code>get_device_variable()</code> (py-lablib.devices.Ophir.base.VegaPowerMeter method), 643
<code>get_device_variable()</code> (py-lablib.devices.interface.camera.ICamera method), 834	<code>get_device_variable()</code> (py-lablib.devices.OZOptics.base.DD100 method), 635
<code>get_device_variable()</code> (py-lablib.devices.interface.camera.IExposureCamera method), 847	<code>get_device_variable()</code> (py-lablib.devices.OZOptics.base.EPC04 method), 637
	<code>get_device_variable()</code> (py-



<i>lablib.devices.OZOptics.base.OZOpticsDevice</i> <i>method</i> ), 632	<i>lablib.devices.Tektronix.base.ITektronixScope</i> <i>method</i> ), 753
<i>get_device_variable()</i> (py- <i>lablib.devices.OZOptics.base.TF100</i> <i>method</i> ), 634	<i>get_device_variable()</i> (py- <i>lablib.devices.Tektronix.base.TDS2000</i> <i>method</i> ), 756
<i>get_device_variable()</i> (py- <i>lablib.devices.PCO.SC2.PCOSC2Camera</i> <i>method</i> ), 649	<i>get_device_variable()</i> (py- <i>lablib.devices.Thorlabs.elliptec.ElliptecMotor</i> <i>method</i> ), 780
<i>get_device_variable()</i> (py- <i>lablib.devices.Pfeiffer.base.DPG202</i> <i>method</i> ), 658	<i>get_device_variable()</i> (py- <i>lablib.devices.Thorlabs.kinesis.BasicKinesisDevice</i> <i>method</i> ), 783
<i>get_device_variable()</i> (py- <i>lablib.devices.Pfeiffer.base.TPG260</i> <i>method</i> ), 656	<i>get_device_variable()</i> (py- <i>lablib.devices.Thorlabs.kinesis.KinesisDevice</i> <i>method</i> ), 789
<i>get_device_variable()</i> (py- <i>lablib.devices.Photometrics.pvcam.PvcamCamera</i> <i>method</i> ), 666	<i>get_device_variable()</i> (py- <i>lablib.devices.Thorlabs.kinesis.KinesisMotor</i> <i>method</i> ), 799
<i>get_device_variable()</i> (py- <i>lablib.devices.PhotonFocus.PhotonFocus.IPhotonFocusCamera</i> <i>method</i> ), 674	<i>get_device_variable()</i> (py- <i>lablib.devices.Thorlabs.kinesis.KinesisPiezoMotor</i> <i>method</i> ), 803
<i>get_device_variable()</i> (py- <i>lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusBitFluxCamera</i> <i>method</i> ), 696	<i>get_device_variable()</i> (py- <i>lablib.devices.Thorlabs.kinesis.KinesisQuadDetector</i> <i>method</i> ), 808
<i>get_device_variable()</i> (py- <i>lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusIMACamera</i> <i>method</i> ), 680	<i>get_device_variable()</i> (py- <i>lablib.devices.Thorlabs.kinesis.MFF</i> <i>method</i> ), 793
<i>get_device_variable()</i> (py- <i>lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusSiSoCamera</i> <i>method</i> ), 688	<i>get_device_variable()</i> (py- <i>lablib.devices.Thorlabs.serial.FW</i> <i>method</i> ), 814
<i>get_device_variable()</i> (py- <i>lablib.devices.PhysikInstrumente.base.GenericPIController</i> <i>method</i> ), 704	<i>get_device_variable()</i> (py- <i>lablib.devices.Thorlabs.serial.FWv1</i> <i>method</i> ), 817
<i>get_device_variable()</i> (py- <i>lablib.devices.PhysikInstrumente.base.PIE516</i> <i>method</i> ), 706	<i>get_device_variable()</i> (py- <i>lablib.devices.Thorlabs.serial.MDT69xA</i> <i>method</i> ), 820
<i>get_device_variable()</i> (py- <i>lablib.devices.PrincetonInstruments.picam.PicamCamera</i> <i>method</i> ), 714	<i>get_device_variable()</i> (py- <i>lablib.devices.Thorlabs.serial.ThorlabsSerialInterface</i> <i>method</i> ), 810
<i>get_device_variable()</i> (py- <i>lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareCamera</i> <i>method</i> ), 728	<i>get_device_variable()</i> (py- <i>lablib.devices.Thorlabs.TLCamera.ThorlabsTLCamera</i> <i>method</i> ), 773
<i>get_device_variable()</i> (py- <i>lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareFrameGrabber</i> <i>method</i> ), 724	<i>get_device_variable()</i> (py- <i>lablib.devices.Toptica.ibeam.TopticaIBeam</i> <i>method</i> ), 825
<i>get_device_variable()</i> (py- <i>lablib.devices.Sirah.Matisse.SirahMatisse</i> <i>method</i> ), 740	<i>get_device_variable()</i> (py- <i>lablib.devices.Trinamic.base.TMCM1110</i> <i>method</i> ), 829
<i>get_device_variable()</i> (py- <i>lablib.devices.SmarAct.scu3d.SC3D</i> <i>method</i> ), 747	<i>get_device_variable()</i> (py- <i>lablib.devices.uc480.uc480.UC480Camera</i> <i>method</i> ), 870
<i>get_device_variable()</i> (py- <i>lablib.devices.Tektronix.base.DPO2000</i> <i>method</i> ), 763	<i>get_devices_number()</i> (in module py- <i>lablib.devices.SmarAct.scu3d</i> ), 746
<i>get_device_variable()</i> (py-	<i>get_dictionary_line()</i> (py- <i>lablib.core.fileio.savefile.DictionaryOutputFileFormat</i>

<i>method</i> ), 193	<i>get_diode_power_lowlevel()</i> (py-
<i>get_digital_gain()</i> (py-	<i>lablib.devices.Sirah.Matisse.SirahMatisse</i>
<i>lablib.devices.AlliedVision.Bonito.BonitoIMAQCamera</i>	<i>method</i> ), 736
<i>method</i> ), 458	<i>get_diode_power_waveform()</i> (py-
<i>get_digital_gain()</i> (py-	<i>lablib.devices.Sirah.Matisse.SirahMatisse</i>
<i>lablib.devices.AlliedVision.Bonito.IBonitoCamera</i>	<i>method</i> ), 736
<i>method</i> ), 452	<i>get_display_channel()</i> (py-
<i>get_digital_input()</i> (py-	<i>lablib.devices.Pfeiffer.base.TPG260</i>
<i>lablib.devices.Arcus.performax.Performax2EXStage</i>	<i>method</i> ), 655
<i>method</i> ), 494	<i>get_display_resolution()</i> (py-
<i>get_digital_input()</i> (py-	<i>lablib.devices.Pfeiffer.base.TPG260</i>
<i>lablib.devices.Arcus.performax.Performax4EXStage</i>	<i>method</i> ), 655
<i>method</i> ), 491	<i>get_display_units()</i> (py-
<i>get_digital_input()</i> (py-	<i>lablib.devices.Cryocon.base.Cryocon1x</i>
<i>lablib.devices.Arcus.performax.PerformaxDMXJSAStage</i>	<i>method</i> ), 536
<i>method</i> ), 498	<i>get_drive_current()</i> (py-
<i>get_digital_input_parameters()</i> (py-	<i>lablib.devices.Toptica.ibeam.TopticaIBeam</i>
<i>lablib.devices.NI.daq.NIDAQ</i>	<i>method</i> ), 825
<i>method</i> ), 623	<i>get_drive_parameters()</i> (py-
<i>get_digital_input_register()</i> (py-	<i>lablib.devices.Thorlabs.kinesis.KinesisPiezoMotor</i>
<i>lablib.devices.Arcus.performax.Performax2EXStage</i>	<i>method</i> ), 802
<i>method</i> ), 494	<i>get_dtype()</i> (pylablib.core.fileio.savefile.TableBinaryOutputFileFormat
<i>get_digital_input_register()</i> (py-	<i>method</i> ), 193
<i>lablib.devices.Arcus.performax.Performax4EXStage</i>	<i>get_duty_cycle()</i> (py-
<i>method</i> ), 491	<i>lablib.devices.AWG.generic.GenericAWG</i>
<i>get_digital_input_register()</i> (py-	<i>method</i> ), 404
<i>lablib.devices.Arcus.performax.PerformaxDMXJSAStage</i>	<i>get_duty_cycle()</i> (py-
<i>method</i> ), 498	<i>lablib.devices.AWG.specific.Agilent33220A</i>
<i>get_digital_output()</i> (py-	<i>method</i> ), 416
<i>lablib.devices.Arcus.performax.Performax2EXStage</i>	<i>get_duty_cycle()</i> (py-
<i>method</i> ), 494	<i>lablib.devices.AWG.specific.Agilent33500</i>
<i>get_digital_output()</i> (py-	<i>method</i> ), 410
<i>lablib.devices.Arcus.performax.Performax4EXStage</i>	<i>get_duty_cycle()</i> (py-
<i>method</i> ), 491	<i>lablib.devices.AWG.specific.InstekAFG2000</i>
<i>get_digital_output()</i> (py-	<i>method</i> ), 427
<i>lablib.devices.Arcus.performax.PerformaxDMXJSAStage</i>	<i>get_duty_cycle()</i> (py-
<i>method</i> ), 498	<i>lablib.devices.AWG.specific.InstekAFG2225</i>
<i>get_digital_output_channels()</i> (py-	<i>method</i> ), 422
<i>lablib.devices.NI.daq.NIDAQ</i>	<i>get_duty_cycle()</i> (py-
<i>method</i> ), 624	<i>lablib.devices.AWG.specific.RigolDG1000</i>
<i>get_digital_output_parameters()</i> (py-	<i>method</i> ), 445
<i>lablib.devices.NI.daq.NIDAQ</i>	<i>get_duty_cycle()</i> (py-
<i>method</i> ), 624	<i>lablib.devices.AWG.specific.RSInstekAFG21000</i>
<i>get_digital_output_register()</i> (py-	<i>method</i> ), 433
<i>lablib.devices.Arcus.performax.Performax2EXStage</i>	<i>get_duty_cycle()</i> (py-
<i>method</i> ), 495	<i>lablib.devices.AWG.specific.TektronixAFG1000</i>
<i>get_digital_output_register()</i> (py-	<i>method</i> ), 439
<i>lablib.devices.Arcus.performax.Performax4EXStage</i>	<i>edge_trigger_coupling()</i> (py-
<i>method</i> ), 492	<i>lablib.devices.Tektronix.base.DPO2000</i>
<i>get_digital_output_register()</i> (py-	<i>method</i> ), 763
<i>lablib.devices.Arcus.performax.PerformaxDMXJSAStage</i>	<i>get_edge_trigger_coupling()</i> (py-
<i>method</i> ), 498	<i>lablib.devices.Tektronix.base.ITektronixScope</i>
<i>get_digital_outputs()</i> (py-	<i>method</i> ), 749
<i>lablib.devices.NI.daq.NIDAQ</i>	<i>get_edge_trigger_coupling()</i> (py-
<i>method</i> ), 624	
<i>get_diode_power()</i> (py-	
<i>lablib.devices.Sirah.Matisse.SirahMatisse</i>	
<i>method</i> ), 736	

<i>lablib.devices.Tektronix.base.TDS2000</i> method), 756		<i>lablib.devices.Thorlabs.kinesis.KinesisPiezoMotor</i> method), 801	
<i>get_edge_trigger_slope()</i> <i>lablib.devices.Tektronix.base.DPO2000</i> method), 763	(py-	<i>get_encoder()</i> <i>lablib.devices.Arcus.performax.Performax2EXStage</i> method), 495	(py-
<i>get_edge_trigger_slope()</i> <i>lablib.devices.Tektronix.base.ITektronixScope</i> method), 749	(py-	<i>get_encoder()</i> <i>lablib.devices.Arcus.performax.Performax4EXStage</i> method), 490	(py-
<i>get_edge_trigger_slope()</i> <i>lablib.devices.Tektronix.base.TDS2000</i> method), 756	(py-	<i>get_energy()</i> ( <i>pylablib.devices.Ophir.base.VegaPowerMeter</i> method), 641	
<i>get_edge_trigger_source()</i> <i>lablib.devices.Tektronix.base.DPO2000</i> method), 763	(py-	<i>get_entry()</i> ( <i>pylablib.core.utils.dictionary.Dictionary</i> method), 327	
<i>get_edge_trigger_source()</i> <i>lablib.devices.Tektronix.base.ITektronixScope</i> method), 749	(py-	<i>get_entry()</i> ( <i>pylablib.core.utils.dictionary.DictionaryPointer</i> method), 338	
<i>get_edge_trigger_source()</i> <i>lablib.devices.Tektronix.base.TDS2000</i> method), 756	(py-	<i>get_entry()</i> ( <i>pylablib.core.utils.dictionary.FilterTree</i> method), 355	
<i>get_element_position()</i> <i>lablib.core.gui.widgets.container.IQWidgetContainer</i> method), 207	(py-	<i>get_entry()</i> ( <i>pylablib.core.utils.dictionary.PrefixTree</i> method), 347	
<i>get_element_position()</i> <i>lablib.core.gui.widgets.container.QDialogContainer</i> method), 219	(py-	<i>get_error_code()</i> <i>lablib.devices.Pfeiffer.base.DPG202</i> method), 658	(py-
<i>get_element_position()</i> <i>lablib.core.gui.widgets.container.QFrameContainer</i> method), 215	(py-	<i>get_esr()</i> ( <i>pylablib.core.devio.SCPi.SCPIDevice</i> method), 139	
<i>get_element_position()</i> <i>lablib.core.gui.widgets.container.QGroupBoxContainer</i> method), 223	(py-	<i>get_esr()</i> ( <i>pylablib.devices.AWG.generic.GenericAWG</i> method), 406	
<i>get_element_position()</i> <i>lablib.core.gui.widgets.container.QScrollAreaContainer</i> method), 227	(py-	<i>get_esr()</i> ( <i>pylablib.devices.AWG.specific.Agilent33220A</i> method), 416	
<i>get_element_position()</i> <i>lablib.core.gui.widgets.container.QWidgetContainer</i> method), 211	(py-	<i>get_esr()</i> ( <i>pylablib.devices.AWG.specific.Agilent33500</i> method), 410	
<i>get_element_position()</i> <i>lablib.core.gui.widgets.layout_manager.IQLayoutManager</i> method), 239	(py-	<i>get_esr()</i> ( <i>pylablib.devices.AWG.specific.InstekAFG2000</i> method), 427	
<i>get_element_position()</i> <i>lablib.core.gui.widgets.layout_manager.QLayoutManagedWidget</i> method), 240	(py-	<i>get_esr()</i> ( <i>pylablib.devices.AWG.specific.InstekAFG2225</i> method), 422	
<i>get_element_position()</i> <i>lablib.core.gui.widgets.param_table.ParamTable</i> method), 250	(py-	<i>get_esr()</i> ( <i>pylablib.devices.AWG.specific.RigolDG1000</i> method), 445	
<i>get_element_position()</i> <i>lablib.core.gui.widgets.param_table.StatusTable</i> method), 258	(py-	<i>get_esr()</i> ( <i>pylablib.devices.AWG.specific.RSInstekAFG21000</i> method), 433	
<i>get_EMCCD_gain()</i> <i>lablib.devices.Andor.AndorSDK2.AndorSDK2Camera</i> method), 468	(py-	<i>get_esr()</i> ( <i>pylablib.devices.AWG.specific.TektronixAFG1000</i> method), 439	
<i>get_enabled_channels()</i>	(py-	<i>get_esr()</i> ( <i>pylablib.devices.Cryocon.base.Cryocon1x</i> method), 537	
		<i>get_esr()</i> ( <i>pylablib.devices.Cryomagnetics.base.LM500</i> method), 541	
		<i>get_esr()</i> ( <i>pylablib.devices.Cryomagnetics.base.LM510</i> method), 544	
		<i>get_esr()</i> ( <i>pylablib.devices.Lakeshore.base.Lakeshore218</i> method), 595	
		<i>get_esr()</i> ( <i>pylablib.devices.Lakeshore.base.Lakeshore370</i> method), 600	
		<i>get_esr()</i> ( <i>pylablib.devices.Sirah.Matisse.SirahMatisse</i> method), 741	
		<i>get_esr()</i> ( <i>pylablib.devices.Tektronix.base.DPO2000</i> method), 763	
	(py-	<i>get_esr()</i> ( <i>pylablib.devices.Tektronix.base.ITektronixScope</i>	

method), 753

get\_esr() (pylablib.devices.Tektronix.base.TDS2000 method), 757

get\_esr() (pylablib.devices.Thorlabs.serial.FW method), 814

get\_esr() (pylablib.devices.Thorlabs.serial.FWv1 method), 817

get\_esr() (pylablib.devices.Thorlabs.serial.MDT69xA method), 820

get\_esr() (pylablib.devices.Thorlabs.serial.ThorlabsSerialInterface method), 810

get\_etalon\_lock\_status() (pylablib.devices.M2.solstis.Solstis method), 616

get\_ethernet\_parameters() (pylablib.devices.Newport.picomotor.Picomotor8742 method), 629

get\_exec\_counter() (pylablib.core.thread.controller.QTaskThread method), 309

get\_exec\_counter() (pylablib.core.thread.controller.QThreadController method), 300

get\_executable() (in module pylablib.core.utils.module), 386

get\_export\_clock\_terminal() (pylablib.devices.NI.daq.NIDAQ method), 622

get\_exposure() (pylablib.devices.AlliedVision.Bonito.BonitoIMAQCamera method), 458

get\_exposure() (pylablib.devices.AlliedVision.Bonito.IBonitoCamera method), 451

get\_exposure() (pylablib.devices.Andor.AndorSDK2.AndorSDK2Camera method), 470

get\_exposure() (pylablib.devices.Andor.AndorSDK3.AndorSDK3Camera method), 481

get\_exposure() (pylablib.devices.Basler.pylon.BaslerPylonCamera method), 515

get\_exposure() (pylablib.devices.DCAM.DCAM.DCAMCamera method), 551

get\_exposure() (pylablib.devices.HighFinesse.wlm.WLM method), 558

get\_exposure() (pylablib.devices.interface.camera.IExposureCamera method), 847

get\_exposure() (pylablib.devices.PCO.SC2.PCOS2Camera method), 647

get\_exposure() (pylablib.devices.Photometrics.pvcam.PvcamCamera method), 663

get\_exposure() (pylablib.devices.PhotonFocus.PhotonFocus.IPhotonFocusCamera method), 673

get\_exposure() (pylablib.devices.PhotonFocus.PhotonFocus.PhotonFocusBitFlowCamera method), 696

get\_exposure() (pylablib.devices.PhotonFocus.PhotonFocus.PhotonFocusIMAQCamera method), 680

get\_exposure() (pylablib.devices.PhotonFocus.PhotonFocus.PhotonFocusSiSoCamera method), 689

get\_exposure() (pylablib.devices.PrincetonInstruments.picam.PicamCamera method), 713

get\_exposure() (pylablib.devices.Thorlabs.TLCamera.ThorlabsTLCamera method), 773

get\_exposure() (pylablib.devices.uc480.uc480.UC480Camera method), 870

get\_exposure\_control\_mode() (pylablib.devices.AlliedVision.Bonito.BonitoIMAQCamera method), 458

get\_exposure\_control\_mode() (pylablib.devices.AlliedVision.Bonito.IBonitoCamera method), 451

get\_exposure\_mode() (pylablib.devices.HighFinesse.wlm.WLM method), 558

get\_ext() (pylablib.core.fileio.location.LocationName method), 183

get\_ext\_trigger\_parameters() (pylablib.devices.Andor.AndorSDK2.AndorSDK2Camera method), 469

get\_ext\_trigger\_parameters() (pylablib.devices.DCAM.DCAM.DCAMCamera method), 551

get\_ext\_trigger\_parameters() (pylablib.devices.Thorlabs.TLCamera.ThorlabsTLCamera method), 771

get\_external\_input\_modes() (pylablib.devices.Attocube.anc300.ANC300 method), 504

get\_fan\_mode() (pylablib.devices.Andor.AndorSDK2.AndorSDK2Camera method), 469

get\_fan\_mode() (pylablib.devices.Photometrics.pvcam.PvcamCamera method), 663

get\_fast\_kinetic\_mode\_parameters() (py-

<i>lablib.devices.Andor.AndorSDK2.AndorSDK2Camera</i> <i>method</i> ), 470	791	<i>get_frame()</i> ( <i>pylablib.devices.Thorlabs.TLCamera.ThorlabsTLCamera</i> <i>method</i> ), 771
<i>get_fast_scan_status()</i> ( <i>py-</i> <i>lablib.devices.M2.solstis.Solstis</i> <i>method</i> ), 619		<i>get_frame_delay()</i> ( <i>py-</i> <i>lablib.devices.PCO.SC2.PCOSC2Camera</i> <i>method</i> ), 647
<i>get_fastpiezo_ctl_params()</i> ( <i>py-</i> <i>lablib.devices.Sirah.Matisse.SirahMatisse</i> <i>method</i> ), 739		<i>get_frame_format()</i> ( <i>py-</i> <i>lablib.devices.AlliedVision.Bonito.BonitoIMAQCamera</i> <i>method</i> ), 458
<i>get_fastpiezo_ctl_status()</i> ( <i>py-</i> <i>lablib.devices.Sirah.Matisse.SirahMatisse</i> <i>method</i> ), 739		<i>get_frame_format()</i> ( <i>py-</i> <i>lablib.devices.AlliedVision.Bonito.IBonitoCamera</i> <i>method</i> ), 453
<i>get_fastpiezo_position()</i> ( <i>py-</i> <i>lablib.devices.Sirah.Matisse.SirahMatisse</i> <i>method</i> ), 739		<i>get_frame_format()</i> ( <i>py-</i> <i>lablib.devices.Andor.AndorSDK2.AndorSDK2Camera</i> <i>method</i> ), 473
<i>get_file_creation_time()</i> ( <i>in module py-</i> <i>lablib.core.utils.files</i> ), 361		<i>get_frame_format()</i> ( <i>py-</i> <i>lablib.devices.Andor.AndorSDK3.AndorSDK3Camera</i> <i>method</i> ), 483
<i>get_file_modification_time()</i> ( <i>in module py-</i> <i>lablib.core.utils.files</i> ), 361		<i>get_frame_format()</i> ( <i>py-</i> <i>lablib.devices.Basler.pylon.BaslerPylonCamera</i> <i>method</i> ), 517
<i>get_filesystem_path()</i> ( <i>py-</i> <i>lablib.core.fileio.location.FolderFileSystemDataLocation</i> <i>method</i> ), 187		<i>get_frame_format()</i> ( <i>py-</i> <i>lablib.devices.BitFlow.BitFlow.BitFlowCamera</i> <i>method</i> ), 528
<i>get_filesystem_path()</i> ( <i>py-</i> <i>lablib.core.fileio.location.IFileSystemDataLocation</i> <i>method</i> ), 185		<i>get_frame_format()</i> ( <i>py-</i> <i>lablib.devices.BitFlow.BitFlow.BitFlowFrameGrabber</i> <i>method</i> ), 524
<i>get_filesystem_path()</i> ( <i>py-</i> <i>lablib.core.fileio.location.PrefixedFileSystemDataLocation</i> <i>method</i> ), 187		<i>get_frame_format()</i> ( <i>py-</i> <i>lablib.devices.DCAM.DCAM.DCAMCamera</i> <i>method</i> ), 553
<i>get_filesystem_path()</i> ( <i>py-</i> <i>lablib.core.fileio.location.SingleFileSystemDataLocation</i> <i>method</i> ), 186		<i>get_frame_format()</i> ( <i>py-</i> <i>lablib.devices.IMAQ.IMAQ.IMAQCamera</i> <i>method</i> ), 570
<i>get_fill_status()</i> ( <i>py-</i> <i>lablib.devices.Cryomagnetics.base.LM500</i> <i>method</i> ), 541		<i>get_frame_format()</i> ( <i>py-</i> <i>lablib.devices.IMAQ.IMAQ.IMAQFrameGrabber</i> <i>method</i> ), 565
<i>get_fill_status()</i> ( <i>py-</i> <i>lablib.devices.Cryomagnetics.base.LM510</i> <i>method</i> ), 545		<i>get_frame_format()</i> ( <i>py-</i> <i>lablib.devices.IMAQdx.IMAQdx.EthernetIMAQdxCamera</i> <i>method</i> ), 585
<i>get_filter_settings()</i> ( <i>py-</i> <i>lablib.devices.Lakeshore.base.Lakeshore218</i> <i>method</i> ), 595		<i>get_frame_format()</i> ( <i>py-</i> <i>lablib.devices.IMAQdx.IMAQdx.IMAQdxCamera</i> <i>method</i> ), 580
<i>get_filter_settings()</i> ( <i>py-</i> <i>lablib.devices.Lakeshore.base.Lakeshore370</i> <i>method</i> ), 600		<i>get_frame_format()</i> ( <i>py-</i> <i>lablib.devices.interface.camera.IAttributeCamera</i> <i>method</i> ), 838
<i>get_fine_tuning_status()</i> ( <i>py-</i> <i>lablib.devices.M2.solstis.Solstis</i> <i>method</i> ), 616		<i>get_frame_format()</i> ( <i>py-</i> <i>lablib.devices.interface.camera.IBinROICamera</i> <i>method</i> ), 857
<i>get_fine_wavelength()</i> ( <i>py-</i> <i>lablib.devices.M2.solstis.Solstis</i> <i>method</i> ), 616		<i>get_frame_format()</i> ( <i>py-</i> <i>lablib.devices.interface.camera.ICamera</i> <i>method</i> ), 833
<i>get_first_empty_column()</i> ( <i>in module py-</i> <i>lablib.core.gui.utils</i> ), 264		<i>get_frame_format()</i> ( <i>py-</i> <i>lablib.devices.interface.camera.IExposureCamera</i> <i>method</i> ), 847
<i>get_first_empty_row()</i> ( <i>in module py-</i> <i>lablib.core.gui.utils</i> ), 264		
<i>get_flipper_parameters()</i> ( <i>py-</i> <i>lablib.devices.Thorlabs.kinesis.MFF</i> <i>method</i> ),		



<code>get_frame_format()</code> lablib.devices.interface.camera.IGrabberAttributeCamera method), 843	(py- <code>get_frame_info_fields()</code> lablib.devices.BitFlow.BitFlow.BitFlowCamera method), 528	(py-
<code>get_frame_format()</code> lablib.devices.interface.camera.IROICamera method), 852	(py- <code>get_frame_info_fields()</code> lablib.devices.BitFlow.BitFlow.BitFlowFrameGrabber method), 524	(py-
<code>get_frame_format()</code> lablib.devices.PCO.SC2.PCOSC2Camera method), 649	(py- <code>get_frame_info_fields()</code> lablib.devices.DCAM.DCAM.DCAMCamera method), 553	(py-
<code>get_frame_format()</code> lablib.devices.Photometrics.pvcam.PvcamCamera method), 666	(py- <code>get_frame_info_fields()</code> lablib.devices.IMAQ.IMAQ.IMAQCamera method), 570	(py-
<code>get_frame_format()</code> lablib.devices.PhotonFocus.PhotonFocus.IPhotonFocusCamera method), 674	(py- <code>get_frame_info_fields()</code> lablib.devices.IMAQ.IMAQ.IMAQFrameGrabber method), 565	(py-
<code>get_frame_format()</code> lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusBitFlowCamera method), 696	(py- <code>get_frame_info_fields()</code> lablib.devices.IMAQdx.IMAQdx.EthernetIMAQdxCamera method), 585	(py-
<code>get_frame_format()</code> lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusIMAQCamera method), 680	(py- <code>get_frame_info_fields()</code> lablib.devices.IMAQdx.IMAQdx.IMAQdxCamera method), 581	(py-
<code>get_frame_format()</code> lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusSiSoCamera method), 689	(py- <code>get_frame_info_fields()</code> lablib.devices.interface.camera.IAttributeCamera method), 838	(py-
<code>get_frame_format()</code> lablib.devices.PrincetonInstruments.picam.PicamCamera method), 714	(py- <code>get_frame_info_fields()</code> lablib.devices.interface.camera.IBinROICamera method), 857	(py-
<code>get_frame_format()</code> lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareCamera method), 728	(py- <code>get_frame_info_fields()</code> lablib.devices.interface.camera.ICamera method), 833	(py-
<code>get_frame_format()</code> lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareFrameGrabber method), 724	(py- <code>get_frame_info_fields()</code> lablib.devices.interface.camera.IExposureCamera method), 848	(py-
<code>get_frame_format()</code> lablib.devices.Thorlabs.TLCamera.ThorlabsTLCamera method), 773	(py- <code>get_frame_info_fields()</code> lablib.devices.interface.camera.IGrabberAttributeCamera method), 843	(py-
<code>get_frame_format()</code> lablib.devices.uc480.uc480.UC480Camera method), 870	(py- <code>get_frame_info_fields()</code> lablib.devices.interface.camera.IROICamera method), 853	(py-
<code>get_frame_info_fields()</code> lablib.devices.AlliedVision.Bonito.BonitoIMAQCamera method), 458	(py- <code>get_frame_info_fields()</code> lablib.devices.PCO.SC2.PCOSC2Camera method), 649	(py-
<code>get_frame_info_fields()</code> lablib.devices.AlliedVision.Bonito.IBonitoCamera method), 453	(py- <code>get_frame_info_fields()</code> lablib.devices.Photometrics.pvcam.PvcamCamera method), 666	(py-
<code>get_frame_info_fields()</code> lablib.devices.Andor.AndorSDK2.AndorSDK2Camera method), 473	(py- <code>get_frame_info_fields()</code> lablib.devices.PhotonFocus.PhotonFocus.IPhotonFocusCamera method), 674	(py-
<code>get_frame_info_fields()</code> lablib.devices.Andor.AndorSDK3.AndorSDK3Camera method), 483	(py- <code>get_frame_info_fields()</code> lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusBitFlowCa method), 696	(py-
<code>get_frame_info_fields()</code> lablib.devices.Basler.pylon.BaslerPylonCamera method), 517	(py- <code>get_frame_info_fields()</code> lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusIMAQCam method), 680	(py-

<code>get_frame_info_fields()</code> <i>lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusSiSoCamera</i> <i>method), 689</i>	<i>(py- get_frame_info_format()</i> <i>lablib.devices.interface.camera.IAttributeCamera</i> <i>method), 838</i>	<i>(py-</i>
<code>get_frame_info_fields()</code> <i>lablib.devices.PrincetonInstruments.picam.PicamCamera</i> <i>method), 715</i>	<i>(py- get_frame_info_format()</i> <i>lablib.devices.interface.camera.IBinROICamera</i> <i>method), 857</i>	<i>(py-</i>
<code>get_frame_info_fields()</code> <i>lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareCamera</i> <i>method), 728</i>	<i>(py- get_frame_info_format()</i> <i>lablib.devices.interface.camera.ICamera</i> <i>method), 833</i>	<i>(py-</i>
<code>get_frame_info_fields()</code> <i>lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareFrameGrabber</i> <i>method), 724</i>	<i>(py- get_frame_info_format()</i> <i>lablib.devices.interface.camera.IExposureCamera</i> <i>method), 848</i>	<i>(py-</i>
<code>get_frame_info_fields()</code> <i>lablib.devices.Thorlabs.TLCamera.ThorlabsTLCamera</i> <i>method), 773</i>	<i>(py- get_frame_info_format()</i> <i>lablib.devices.interface.camera.IGrabberAttributeCamera</i> <i>method), 843</i>	<i>(py-</i>
<code>get_frame_info_fields()</code> <i>lablib.devices.uc480.uc480.UC480Camera</i> <i>method), 870</i>	<i>(py- get_frame_info_format()</i> <i>lablib.devices.interface.camera.IROICamera</i> <i>method), 853</i>	<i>(py-</i>
<code>get_frame_info_format()</code> <i>lablib.devices.AlliedVision.Bonito.BonitoIMAQCamera</i> <i>method), 458</i>	<i>(py- get_frame_info_format()</i> <i>lablib.devices.PCO.SC2.PCOSC2Camera</i> <i>method), 650</i>	<i>(py-</i>
<code>get_frame_info_format()</code> <i>lablib.devices.AlliedVision.Bonito.IBonitoCamera</i> <i>method), 453</i>	<i>(py- get_frame_info_format()</i> <i>lablib.devices.Photometrics.pvcam.PvcamCamera</i> <i>method), 666</i>	<i>(py-</i>
<code>get_frame_info_format()</code> <i>lablib.devices.Andor.AndorSDK2.AndorSDK2Camera</i> <i>method), 473</i>	<i>(py- get_frame_info_format()</i> <i>lablib.devices.PhotonFocus.PhotonFocus.IPhotonFocusCamera</i> <i>method), 674</i>	<i>(py-</i>
<code>get_frame_info_format()</code> <i>lablib.devices.Andor.AndorSDK3.AndorSDK3Camera</i> <i>method), 483</i>	<i>(py- get_frame_info_format()</i> <i>lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusBitFlowCamera</i> <i>method), 696</i>	<i>(py-</i>
<code>get_frame_info_format()</code> <i>lablib.devices.Basler.pylon.BaslerPylonCamera</i> <i>method), 517</i>	<i>(py- get_frame_info_format()</i> <i>lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusIMAQCamera</i> <i>method), 680</i>	<i>(py-</i>
<code>get_frame_info_format()</code> <i>lablib.devices.BitFlow.BitFlow.BitFlowCamera</i> <i>method), 528</i>	<i>(py- get_frame_info_format()</i> <i>lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusSiSoCamera</i> <i>method), 689</i>	<i>(py-</i>
<code>get_frame_info_format()</code> <i>lablib.devices.BitFlow.BitFlow.BitFlowFrameGrabber</i> <i>method), 524</i>	<i>(py- get_frame_info_format()</i> <i>lablib.devices.PrincetonInstruments.picam.PicamCamera</i> <i>method), 715</i>	<i>(py-</i>
<code>get_frame_info_format()</code> <i>lablib.devices.DCAM.DCAM.DCAMCamera</i> <i>method), 553</i>	<i>(py- get_frame_info_format()</i> <i>lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareCamera</i> <i>method), 728</i>	<i>(py-</i>
<code>get_frame_info_format()</code> <i>lablib.devices.IMAQ.IMAQ.IMAQCamera</i> <i>method), 570</i>	<i>(py- get_frame_info_format()</i> <i>lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareFrameGrabber</i> <i>method), 724</i>	<i>(py-</i>
<code>get_frame_info_format()</code> <i>lablib.devices.IMAQ.IMAQ.IMAQFrameGrabber</i> <i>method), 565</i>	<i>(py- get_frame_info_format()</i> <i>lablib.devices.Thorlabs.TLCamera.ThorlabsTLCamera</i> <i>method), 773</i>	<i>(py-</i>
<code>get_frame_info_format()</code> <i>lablib.devices.IMAQdx.IMAQdx.EthernetIMAQdxCamera</i> <i>method), 585</i>	<i>(py- get_frame_info_format()</i> <i>lablib.devices.uc480.uc480.UC480Camera</i> <i>method), 870</i>	<i>(py-</i>
<code>get_frame_info_format()</code> <i>lablib.devices.IMAQdx.IMAQdx.IMAQdxCamera</i> <i>method), 581</i>	<i>(py- get_frame_info_period()</i> <i>lablib.devices.AlliedVision.Bonito.BonitoIMAQCamera</i> <i>method), 458</i>	<i>(py-</i>

<code>get_frame_info_period()</code> <i>lablib.devices.AlliedVision.Bonito.IBonitoCamera</i> method), 453	(py- <i>lablib.devices.Photometrics.pvcam.PvcamCamera</i> method), 666	(py- <i>lablib.devices.Photometrics.pvcam.PvcamCamera</i> method), 666
<code>get_frame_info_period()</code> <i>lablib.devices.Andor.AndorSDK2.AndorSDK2Camera</i> method), 473	(py- <i>lablib.devices.PhotonFocus.PhotonFocus.IPhotonFocusCamera</i> method), 675	(py- <i>lablib.devices.PhotonFocus.PhotonFocus.IPhotonFocusCamera</i> method), 675
<code>get_frame_info_period()</code> <i>lablib.devices.Andor.AndorSDK3.AndorSDK3Camera</i> method), 483	(py- <i>lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusBitFlowCa</i> method), 697	(py- <i>lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusBitFlowCa</i> method), 697
<code>get_frame_info_period()</code> <i>lablib.devices.Basler.pylon.BaslerPylonCamera</i> method), 517	(py- <i>lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusIMAQCam</i> method), 681	(py- <i>lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusIMAQCam</i> method), 681
<code>get_frame_info_period()</code> <i>lablib.devices.BitFlow.BitFlow.BitFlowCamera</i> method), 528	(py- <i>lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusSiSoCame</i> method), 689	(py- <i>lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusSiSoCame</i> method), 689
<code>get_frame_info_period()</code> <i>lablib.devices.BitFlow.BitFlow.BitFlowFrameGrabber</i> method), 524	(py- <i>lablib.devices.PrincetonInstruments.picam.PicamCamera</i> method), 715	(py- <i>lablib.devices.PrincetonInstruments.picam.PicamCamera</i> method), 715
<code>get_frame_info_period()</code> <i>lablib.devices.DCAM.DCAM.DCAMCamera</i> method), 554	(py- <i>lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareCamera</i> method), 729	(py- <i>lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareCamera</i> method), 729
<code>get_frame_info_period()</code> <i>lablib.devices.IMAQ.IMAQ.IMAQCamera</i> method), 570	(py- <i>lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareFrameGrabb</i> method), 724	(py- <i>lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareFrameGrabb</i> method), 724
<code>get_frame_info_period()</code> <i>lablib.devices.IMAQ.IMAQ.IMAQFrameGrabber</i> method), 566	(py- <i>lablib.devices.Thorlabs.TLCamera.ThorlabsTLCamera</i> method), 773	(py- <i>lablib.devices.Thorlabs.TLCamera.ThorlabsTLCamera</i> method), 773
<code>get_frame_info_period()</code> <i>lablib.devices.IMAQdx.IMAQdx.EthernetIMAQdxCamera</i> method), 585	(py- <i>lablib.devices.uc480.uc480.UC480Camera</i> method), 870	(py- <i>lablib.devices.uc480.uc480.UC480Camera</i> method), 870
<code>get_frame_info_period()</code> <i>lablib.devices.IMAQdx.IMAQdx.IMAQdxCamera</i> method), 581	(py- <i>lablib.devices.AlliedVision.Bonito.BonitoIMAQCamera</i> method), 459	(py- <i>lablib.devices.AlliedVision.Bonito.BonitoIMAQCamera</i> method), 459
<code>get_frame_info_period()</code> <i>lablib.devices.interface.camera.IAttributeCamera</i> method), 838	(py- <i>lablib.devices.AlliedVision.Bonito.IBonitoCamera</i> method), 452	(py- <i>lablib.devices.AlliedVision.Bonito.IBonitoCamera</i> method), 452
<code>get_frame_info_period()</code> <i>lablib.devices.interface.camera.IBinROICamera</i> method), 857	(py- <i>lablib.devices.Andor.AndorSDK2.AndorSDK2Camera</i> method), 473	(py- <i>lablib.devices.Andor.AndorSDK2.AndorSDK2Camera</i> method), 473
<code>get_frame_info_period()</code> <i>lablib.devices.interface.camera.ICamera</i> method), 833	(py- <i>lablib.devices.Andor.AndorSDK3.AndorSDK3Camera</i> method), 481	(py- <i>lablib.devices.Andor.AndorSDK3.AndorSDK3Camera</i> method), 481
<code>get_frame_info_period()</code> <i>lablib.devices.interface.camera.IExposureCamera</i> method), 848	(py- <i>lablib.devices.Basler.pylon.BaslerPylonCamera</i> method), 515	(py- <i>lablib.devices.Basler.pylon.BaslerPylonCamera</i> method), 515
<code>get_frame_info_period()</code> <i>lablib.devices.interface.camera.IGrabberAttributeCamera</i> method), 843	(py- <i>lablib.devices.DCAM.DCAM.DCAMCamera</i> method), 554	(py- <i>lablib.devices.DCAM.DCAM.DCAMCamera</i> method), 554
<code>get_frame_info_period()</code> <i>lablib.devices.interface.camera.IROICamera</i> method), 853	(py- <i>lablib.devices.interface.camera.IExposureCamera</i> method), 847	(py- <i>lablib.devices.interface.camera.IExposureCamera</i> method), 847
<code>get_frame_info_period()</code> <i>lablib.devices.PCO.SC2.PCOSC2Camera</i> method), 650	(py- <i>lablib.devices.PCO.SC2.PCOSC2Camera</i> method), 647	(py- <i>lablib.devices.PCO.SC2.PCOSC2Camera</i> method), 647



<code>get_frame_period()</code> <i>lablib.devices.Photometrics.pvcam.PvcamCamera</i> method), 666	(py- <code>get_frame_timings()</code> <i>lablib.devices.PhotonFocus.PhotonFocus.IPhotonFocusCamera</i> method), 673	(py-
<code>get_frame_period()</code> <i>lablib.devices.PhotonFocus.PhotonFocus.IPhotonFocusCamera</i> method), 673	(py- <code>get_frame_timings()</code> <i>lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusBitFlowCamera</i> method), 697	(py-
<code>get_frame_period()</code> <i>lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusBitFlowCamera</i> method), 697	(py- <code>get_frame_timings()</code> <i>lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusIMAQCamera</i> method), 681	(py-
<code>get_frame_period()</code> <i>lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusIMAQCamera</i> method), 681	(py- <code>get_frame_timings()</code> <i>lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusSiSoCamera</i> method), 689	(py-
<code>get_frame_period()</code> <i>lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusSiSoCamera</i> method), 689	(py- <code>get_frame_timings()</code> <i>lablib.devices.PrincetonInstruments.picam.PicamCamera</i> method), 713	(py-
<code>get_frame_period()</code> <i>lablib.devices.PrincetonInstruments.picam.PicamCamera</i> method), 713	(py- <code>get_frame_timings()</code> <i>lablib.devices.Thorlabs.TLCamera.ThorlabsTLCamera</i> method), 771	(py-
<code>get_frame_period()</code> <i>lablib.devices.Thorlabs.TLCamera.ThorlabsTLCamera</i> method), 773	(py- <code>get_frame_timings()</code> <i>lablib.devices.uc480.uc480.UC480Camera</i> method), 867	(py-
<code>get_frame_period()</code> <i>lablib.devices.uc480.uc480.UC480Camera</i> method), 870	(py- <code>get_frames_data()</code> <i>lablib.devices.interface.camera.ChunkBufferManager</i> method), 836	(py-
<code>get_frame_readout_time()</code> <i>lablib.devices.DCAM.DCAM.DCAMCamera</i> method), 552	(py- <code>get_frames_status()</code> <i>lablib.devices.AlliedVision.Bonito.BonitoIMAQCamera</i> method), 459	(py-
<code>get_frame_timings()</code> <i>lablib.devices.AlliedVision.Bonito.BonitoIMAQCamera</i> method), 459	(py- <code>get_frames_status()</code> <i>lablib.devices.AlliedVision.Bonito.IBonitoCamera</i> method), 453	(py-
<code>get_frame_timings()</code> <i>lablib.devices.AlliedVision.Bonito.IBonitoCamera</i> method), 452	(py- <code>get_frames_status()</code> <i>lablib.devices.Andor.AndorSDK2.AndorSDK2Camera</i> method), 473	(py-
<code>get_frame_timings()</code> <i>lablib.devices.Andor.AndorSDK2.AndorSDK2Camera</i> method), 471	(py- <code>get_frames_status()</code> <i>lablib.devices.Andor.AndorSDK3.AndorSDK3Camera</i> method), 484	(py-
<code>get_frame_timings()</code> <i>lablib.devices.Andor.AndorSDK3.AndorSDK3Camera</i> method), 481	(py- <code>get_frames_status()</code> <i>lablib.devices.Basler.pylon.BaslerPylonCamera</i> method), 517	(py-
<code>get_frame_timings()</code> <i>lablib.devices.Basler.pylon.BaslerPylonCamera</i> method), 515	(py- <code>get_frames_status()</code> <i>lablib.devices.BitFlow.BitFlow.BitFlowCamera</i> method), 528	(py-
<code>get_frame_timings()</code> <i>lablib.devices.DCAM.DCAM.DCAMCamera</i> method), 552	(py- <code>get_frames_status()</code> <i>lablib.devices.BitFlow.BitFlow.BitFlowFrameGrabber</i> method), 524	(py-
<code>get_frame_timings()</code> <i>lablib.devices.interface.camera.IExposureCamera</i> method), 847	(py- <code>get_frames_status()</code> <i>lablib.devices.DCAM.DCAM.DCAMCamera</i> method), 554	(py-
<code>get_frame_timings()</code> <i>lablib.devices.PCO.SC2.PCOSC2Camera</i> method), 647	(py- <code>get_frames_status()</code> <i>lablib.devices.IMAQ.IMAQ.IMAQCamera</i> method), 570	(py-
<code>get_frame_timings()</code> <i>lablib.devices.Photometrics.pvcam.PvcamCamera</i> method), 664	(py- <code>get_frames_status()</code> <i>lablib.devices.IMAQ.IMAQ.IMAQFrameGrabber</i> method), 566	(py-

<code>get_frames_status()</code> <i>lablib.devices.IMAQdx.IMAQdx.EthernetIMAQdxCamera</i> <i>method</i> ), 586	(py- <code>get_frames_status()</code> <i>lablib.devices.Thorlabs.TLCamera.ThorlabsTLCamera</i> <i>method</i> ), 773	(py-
<code>get_frames_status()</code> <i>lablib.devices.IMAQdx.IMAQdx.IMAQdxCamera</i> <i>method</i> ), 581	<code>get_frames_status()</code> <i>lablib.devices.uc480.uc480.UC480Camera</i> <i>method</i> ), 868	(py-
<code>get_frames_status()</code> <i>lablib.devices.interface.camera.FrameCounter</i> <i>method</i> ), 836	<code>get_framestamp()</code> <i>lablib.devices.AlliedVision.Bonito.BonitoStatusLineChecker</i> <i>method</i> ), 464	(py-
<code>get_frames_status()</code> <i>lablib.devices.interface.camera.IAttributeCamera</i> <i>method</i> ), 838	<code>get_framestamp()</code> <i>lablib.devices.interface.camera.StatusLineChecker</i> <i>method</i> ), 861	(py-
<code>get_frames_status()</code> <i>lablib.devices.interface.camera.IBinROICamera</i> <i>method</i> ), 857	<code>get_framestamp()</code> <i>lablib.devices.PCO.SC2.StatusLineChecker</i> <i>method</i> ), 653	(py-
<code>get_frames_status()</code> <i>lablib.devices.interface.camera.ICamera</i> <i>method</i> ), 832	<code>get_framestamp()</code> <i>lablib.devices.PhotonFocus.PhotonFocus.StatusLineChecker</i> <i>method</i> ), 702	(py-
<code>get_frames_status()</code> <i>lablib.devices.interface.camera.IExposureCamera</i> <i>method</i> ), 848	<code>get_frequencies()</code> <i>lablib.devices.OZOptics.base.EPC04</i> <i>method</i> ), 637	(py-
<code>get_frames_status()</code> <i>lablib.devices.interface.camera.IGrabberAttributeCamera</i> <i>method</i> ), 843	<code>get_frequency()</code> <i>lablib.devices.Attocube.anc300.ANC300</i> <i>method</i> ), 504	(py-
<code>get_frames_status()</code> <i>lablib.devices.interface.camera.IROICamera</i> <i>method</i> ), 853	<code>get_frequency()</code> <i>lablib.devices.Attocube.anc350.ANC350</i> <i>method</i> ), 508	(py-
<code>get_frames_status()</code> <i>lablib.devices.PCO.SC2.PCOSC2Camera</i> <i>method</i> ), 650	<code>get_frequency()</code> <i>lablib.devices.AWG.generic.GenericAWG</i> <i>method</i> ), 404	(py-
<code>get_frames_status()</code> <i>lablib.devices.Photometrics.pvcam.PvcamCamera</i> <i>method</i> ), 666	<code>get_frequency()</code> <i>lablib.devices.AWG.specific.Agilent33220A</i> <i>method</i> ), 416	(py-
<code>get_frames_status()</code> <i>lablib.devices.PhotonFocus.PhotonFocus.IPhotonFocusCamera</i> <i>method</i> ), 675	<code>get_frequency()</code> <i>lablib.devices.AWG.specific.Agilent33500</i> <i>method</i> ), 410	(py-
<code>get_frames_status()</code> <i>lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusBitFluxCamera</i> <i>method</i> ), 697	<code>get_frequency()</code> <i>lablib.devices.AWG.specific.InstekAFG2000</i> <i>method</i> ), 428	(py-
<code>get_frames_status()</code> <i>lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusIMAQdxCamera</i> <i>method</i> ), 681	<code>get_frequency()</code> <i>lablib.devices.AWG.specific.InstekAFG2225</i> <i>method</i> ), 422	(py-
<code>get_frames_status()</code> <i>lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusSiSoCamera</i> <i>method</i> ), 689	<code>get_frequency()</code> <i>lablib.devices.AWG.specific.RigolDG1000</i> <i>method</i> ), 445	(py-
<code>get_frames_status()</code> <i>lablib.devices.PrincetonInstruments.picam.PicamCamera</i> <i>method</i> ), 715	<code>get_frequency()</code> <i>lablib.devices.AWG.specific.RSInstekAFG21000</i> <i>method</i> ), 433	(py-
<code>get_frames_status()</code> <i>lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareCamera</i> <i>method</i> ), 729	<code>get_frequency()</code> <i>lablib.devices.AWG.specific.TektronixAFG1000</i> <i>method</i> ), 439	(py-
<code>get_frames_status()</code> <i>lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareFrameGrabber</i> <i>method</i> ), 724	<code>get_frequency()</code> <i>lablib.devices.HighFinesse.wlm.WLM</i> <i>method</i> ), 558	(py-

<code>get_frequency()</code> <i>lablib.devices.Ophir.base.VegaPowerMeter</i> method), 641	(py-	<code>get_full_info()</code> <i>lablib.devices.Arduino.base.IArduinoDevice</i> method), 501	(py-
<code>get_frequency()</code> <i>lablib.devices.Sirah.tuner.MatisseTuner</i> method), 743	(py-	<code>get_full_info()</code> <i>lablib.devices.Attocube.anc300.ANC300</i> method), 505	(py-
<code>get_frequency()</code> <i>lablib.devices.Thorlabs.elliptec.ElliptecMotor</i> method), 779	(py-	<code>get_full_info()</code> <i>lablib.devices.Attocube.anc350.ANC350</i> method), 509	(py-
<code>get_full_camera_data()</code> <i>lablib.devices.PCO.SC2.PCOS2Camera</i> method), 646	(py-	<code>get_full_info()</code> <i>lablib.devices.AWG.generic.GenericAWG</i> method), 406	(py-
<code>get_full_coarse_tuning_status()</code> <i>lablib.devices.M2.solstis.Solstis</i> method), 616	(py-	<code>get_full_info()</code> <i>lablib.devices.AWG.specific.Agilent33220A</i> method), 416	(py-
<code>get_full_data()</code> <i>lablib.devices.Topica.ibeam.TopicalIBeam</i> method), 824	(py-	<code>get_full_info()</code> <i>lablib.devices.AWG.specific.Agilent33500</i> method), 410	(py-
<code>get_full_fine_tuning_status()</code> <i>lablib.devices.M2.solstis.Solstis</i> method), 615	(py-	<code>get_full_info()</code> <i>lablib.devices.AWG.specific.InstekAFG2000</i> method), 428	(py-
<code>get_full_info()</code> <i>lablib.core.devio.comm_backend.ICommBackendWrapper</i> method), 160	(py-	<code>get_full_info()</code> <i>lablib.devices.AWG.specific.InstekAFG2225</i> method), 422	(py-
<code>get_full_info()</code> <i>lablib.core.devio.interface.IDevice</i> method), 163	(py-	<code>get_full_info()</code> <i>lablib.devices.AWG.specific.RigolDG1000</i> method), 445	(py-
<code>get_full_info()</code> <i>lablib.core.devio.SCPI.SCPIDevice</i> method), 140	(py-	<code>get_full_info()</code> <i>lablib.devices.AWG.specific.RSInstekAFG21000</i> method), 434	(py-
<code>get_full_info()</code> <i>lablib.devices.AlliedVision.Bonito.BonitoIMAQCamera</i> method), 459	(py-	<code>get_full_info()</code> <i>lablib.devices.AWG.specific.TektronixAFG1000</i> method), 440	(py-
<code>get_full_info()</code> <i>lablib.devices.AlliedVision.Bonito.IBonitoCamera</i> method), 453	(py-	<code>get_full_info()</code> <i>lablib.devices.Basler.pylon.BaslerPylonCamera</i> method), 518	(py-
<code>get_full_info()</code> <i>lablib.devices.Andor.AndorSDK2.AndorSDK2Camera</i> method), 474	(py-	<code>get_full_info()</code> <i>lablib.devices.BitFlow.BitFlow.BitFlowCamera</i> method), 528	(py-
<code>get_full_info()</code> <i>lablib.devices.Andor.AndorSDK3.AndorSDK3Camera</i> method), 484	(py-	<code>get_full_info()</code> <i>lablib.devices.BitFlow.BitFlow.BitFlowFrameGrabber</i> method), 524	(py-
<code>get_full_info()</code> <i>lablib.devices.Arcus.performax.GenericPerformaxStage</i> method), 489	(py-	<code>get_full_info()</code> <i>lablib.devices.Conrad.base.RelayBoard</i> method), 534	(py-
<code>get_full_info()</code> <i>lablib.devices.Arcus.performax.Performax2EXStage</i> method), 495	(py-	<code>get_full_info()</code> <i>lablib.devices.Cryocon.base.Cryocon1x</i> method), 537	(py-
<code>get_full_info()</code> <i>lablib.devices.Arcus.performax.Performax4EXStage</i> method), 492	(py-	<code>get_full_info()</code> <i>lablib.devices.Cryomagnetics.base.LM500</i> method), 541	(py-
<code>get_full_info()</code> <i>lablib.devices.Arcus.performax.PerformaxDMXJSAStage</i> method), 499	(py-	<code>get_full_info()</code> <i>lablib.devices.Cryomagnetics.base.LM510</i> method), 545	(py-

<code>get_full_info()</code> <i>lablib.devices.DCAM.DCAM.DCAMCamera</i> <i>method</i> ), 554	(py-	<code>get_full_info()</code> <i>lablib.devices.Leybold.base.GenericITR</i> <i>method</i> ), 608	(py-
<code>get_full_info()</code> <i>lablib.devices.HighFinesse.wlm.WLM</i> <i>method</i> ), 560	(py-	<code>get_full_info()</code> <i>lablib.devices.Leybold.base.ITR90</i> <i>method</i> ), 610	(py-
<code>get_full_info()</code> <i>lablib.devices.IMAQ.IMAQ.IMAQCamera</i> <i>method</i> ), 571	(py-	<code>get_full_info()</code> <i>lablib.devices.LighthousePhotonics.base.SproutG</i> <i>method</i> ), 613	(py-
<code>get_full_info()</code> <i>lablib.devices.IMAQ.IMAQ.IMAQFrameGrabber</i> <i>method</i> ), 566	(py-	<code>get_full_info()</code> <i>lablib.devices.M2.solstis.Solstis</i> <i>method</i> ), 619	(py-
<code>get_full_info()</code> <i>lablib.devices.IMAQdx.IMAQdx.EthernetIMAQdxCamera</i> <i>method</i> ), 586	(py-	<code>get_full_info()</code> <i>lablib.devices.Newport.picomotor.Picomotor8742</i> <i>method</i> ), 630	(py-
<code>get_full_info()</code> <i>lablib.devices.IMAQdx.IMAQdx.IMAQdxCamera</i> <i>method</i> ), 581	(py-	<code>get_full_info()</code> ( <i>pylablib.devices.NI.daq.NIDAQ</i> <i>method</i> ), 626	(py-
<code>get_full_info()</code> <i>lablib.devices.interface.camera.IAttributeCamera</i> <i>method</i> ), 839	(py-	<code>get_full_info()</code> <i>lablib.devices.Ophir.base.OphirDevice</i> <i>method</i> ), 639	(py-
<code>get_full_info()</code> <i>lablib.devices.interface.camera.IBinROICamera</i> <i>method</i> ), 858	(py-	<code>get_full_info()</code> <i>lablib.devices.Ophir.base.VegaPowerMeter</i> <i>method</i> ), 643	(py-
<code>get_full_info()</code> <i>lablib.devices.interface.camera.ICamera</i> <i>method</i> ), 835	(py-	<code>get_full_info()</code> <i>lablib.devices.OZOptics.base.DD100</i> <i>method</i> ), 635	(py-
<code>get_full_info()</code> <i>lablib.devices.interface.camera.IExposureCamera</i> <i>method</i> ), 848	(py-	<code>get_full_info()</code> <i>lablib.devices.OZOptics.base.EPC04</i> <i>method</i> ), 637	(py-
<code>get_full_info()</code> <i>lablib.devices.interface.camera.IGrabberAttributeCamera</i> <i>method</i> ), 843	(py-	<code>get_full_info()</code> <i>lablib.devices.OZOptics.base.OZOpticsDevice</i> <i>method</i> ), 632	(py-
<code>get_full_info()</code> <i>lablib.devices.interface.camera.IROICamera</i> <i>method</i> ), 853	(py-	<code>get_full_info()</code> <i>lablib.devices.OZOptics.base.TF100</i> <i>method</i> ), 634	(py-
<code>get_full_info()</code> <i>lablib.devices.interface.stage.IMultiaxisStage</i> <i>method</i> ), 863	(py-	<code>get_full_info()</code> <i>lablib.devices.PCO.SC2.PCOSC2Camera</i> <i>method</i> ), 650	(py-
<code>get_full_info()</code> <i>lablib.devices.interface.stage.IStage</i> <i>method</i> ), 862	(py-	<code>get_full_info()</code> <i>lablib.devices.Pfeiffer.base.DPG202</i> <i>method</i> ), 658	(py-
<code>get_full_info()</code> <i>lablib.devices.KJL.base.KJL300</i> <i>method</i> ), 591	(py-	<code>get_full_info()</code> <i>lablib.devices.Pfeiffer.base.TPG260</i> <i>method</i> ), 656	(py-
<code>get_full_info()</code> <i>lablib.devices.Lakeshore.base.Lakeshore218</i> <i>method</i> ), 595	(py-	<code>get_full_info()</code> <i>lablib.devices.Photometrics.pvcam.PvcamCamera</i> <i>method</i> ), 666	(py-
<code>get_full_info()</code> <i>lablib.devices.Lakeshore.base.Lakeshore370</i> <i>method</i> ), 601	(py-	<code>get_full_info()</code> <i>lablib.devices.PhotonFocus.PhotonFocus.IPhotonFocusCamera</i> <i>method</i> ), 675	(py-
<code>get_full_info()</code> <i>lablib.devices.LaserQuantum.base.Finesse</i> <i>method</i> ), 605	(py-	<code>get_full_info()</code> <i>lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusBitFlowCa</i> <i>method</i> ), 697	(py-
		<code>get_full_info()</code>	(py-

`lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusIMAQCamera` (method), 681  
`lablib.devices.Thorlabs.kinesis.MFF` (method), 793  
`get_full_info()` (py- `get_full_info()` (py- `lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusSiSoCamera` (method), 689  
`lablib.devices.Thorlabs.serial.FW` (method), 814  
`get_full_info()` (py- `get_full_info()` (py- `lablib.devices.PhysikInstrumente.base.GenericPIController` (method), 704  
`lablib.devices.Thorlabs.serial.FWv1` (method), 817  
`get_full_info()` (py- `get_full_info()` (py- `lablib.devices.PhysikInstrumente.base.PIE516` (method), 706  
`lablib.devices.Thorlabs.serial.MDT69xA` (method), 820  
`get_full_info()` (py- `get_full_info()` (py- `lablib.devices.PrincetonInstruments.picam.PicamCamera` (method), 715  
`lablib.devices.Thorlabs.serial.ThorlabsSerialInterface` (method), 810  
`get_full_info()` (py- `get_full_info()` (py- `lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareCamera` (method), 729  
`lablib.devices.Thorlabs.TLCamera.ThorlabsTLCamera` (method), 774  
`get_full_info()` (py- `get_full_info()` (py- `lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareFrameGrabber` (method), 724  
`lablib.devices.Toptica.ibeam.TopticaIBeam` (method), 825  
`get_full_info()` (py- `get_full_info()` (py- `lablib.devices.Sirah.Matisse.SirahMatisse` (method), 741  
`lablib.devices.Trinamic.base.TMCM1110` (method), 829  
`get_full_info()` (py- `get_full_info()` (py- `lablib.devices.SmarAct.scu3d.SCU3D` (method), 747  
`lablib.devices.uc480.uc480.UC480Camera` (method), 870  
`get_full_info()` (py- `get_full_status()` (py- `lablib.devices.Tektronix.base.DPO2000` (method), 763  
`lablib.core.devio.comm_backend.ICommBackendWrapper` (method), 161  
`get_full_info()` (py- `get_full_status()` (py- `lablib.devices.Tektronix.base.ITektronixScope` (method), 753  
`lablib.core.devio.interface.IDevice` (method), 162  
`get_full_info()` (py- `get_full_status()` (py- `lablib.devices.Tektronix.base.TDS2000` (method), 757  
`lablib.core.devio.SCPi.SCPIDevice` (method), 140  
`get_full_info()` (py- `get_full_status()` (py- `lablib.devices.Thorlabs.elliptec.ElliptecMotor` (method), 780  
`lablib.devices.AlliedVision.Bonito.BonitoIMAQCamera` (method), 459  
`get_full_info()` (py- `get_full_status()` (py- `lablib.devices.Thorlabs.kinesis.BasicKinesisDevice` (method), 783  
`lablib.devices.AlliedVision.Bonito.IBonitoCamera` (method), 453  
`get_full_info()` (py- `get_full_status()` (py- `lablib.devices.Thorlabs.kinesis.KinesisDevice` (method), 789  
`lablib.devices.Andor.AndorSDK2.AndorSDK2Camera` (method), 474  
`get_full_info()` (py- `get_full_status()` (py- `lablib.devices.Thorlabs.kinesis.KinesisMotor` (method), 799  
`lablib.devices.Andor.AndorSDK3.AndorSDK3Camera` (method), 484  
`get_full_info()` (py- `get_full_status()` (py- `lablib.devices.Thorlabs.kinesis.KinesisPiezoMotor` (method), 803  
`lablib.devices.Arcus.performax.GenericPerformaxStage` (method), 489  
`get_full_info()` (py- `get_full_status()` (py- `lablib.devices.Thorlabs.kinesis.KinesisQuadDetector` (method), 808  
`lablib.devices.Arcus.performax.Performax2EXStage` (method), 495  
`get_full_info()` (py- `get_full_status()` (py-



<code>lablib.devices.Arcus.performax.Performax4EXStage</code> <code>method)</code> , 492	<code>lablib.devices.Cryomagnetics.base.LM500</code> <code>method)</code> , 541
<code>get_full_status()</code> (py- <code>lablib.devices.Arcus.performax.PerformaxDMXJSASStage</code> <code>method)</code> , 499	<code>get_full_status()</code> (py- <code>lablib.devices.Cryomagnetics.base.LM510</code> <code>method)</code> , 545
<code>get_full_status()</code> (py- <code>lablib.devices.Arduino.base.IArduinoDevice</code> <code>method)</code> , 501	<code>get_full_status()</code> (py- <code>lablib.devices.DCAM.DCAM.DCAMCamera</code> <code>method)</code> , 554
<code>get_full_status()</code> (py- <code>lablib.devices.Attocube.anc300.ANC300</code> <code>method)</code> , 505	<code>get_full_status()</code> (py- <code>lablib.devices.HighFinesse.wlm.WLM</code> <code>method)</code> , 560
<code>get_full_status()</code> (py- <code>lablib.devices.Attocube.anc350.ANC350</code> <code>method)</code> , 509	<code>get_full_status()</code> (py- <code>lablib.devices.IMAQ.IMAQ.IMAQCamera</code> <code>method)</code> , 571
<code>get_full_status()</code> (py- <code>lablib.devices.AWG.generic.GenericAWG</code> <code>method)</code> , 406	<code>get_full_status()</code> (py- <code>lablib.devices.IMAQ.IMAQ.IMAQFrameGrabber</code> <code>method)</code> , 566
<code>get_full_status()</code> (py- <code>lablib.devices.AWG.specific.Agilent33220A</code> <code>method)</code> , 416	<code>get_full_status()</code> (py- <code>lablib.devices.IMAQdx.IMAQdx.EthernetIMAQdxCamera</code> <code>method)</code> , 586
<code>get_full_status()</code> (py- <code>lablib.devices.AWG.specific.Agilent33500</code> <code>method)</code> , 410	<code>get_full_status()</code> (py- <code>lablib.devices.IMAQdx.IMAQdx.IMAQdxCamera</code> <code>method)</code> , 581
<code>get_full_status()</code> (py- <code>lablib.devices.AWG.specific.InstekAFG2000</code> <code>method)</code> , 428	<code>get_full_status()</code> (py- <code>lablib.devices.interface.camera.IAttributeCamera</code> <code>method)</code> , 839
<code>get_full_status()</code> (py- <code>lablib.devices.AWG.specific.InstekAFG2225</code> <code>method)</code> , 422	<code>get_full_status()</code> (py- <code>lablib.devices.interface.camera.IBinROICamera</code> <code>method)</code> , 858
<code>get_full_status()</code> (py- <code>lablib.devices.AWG.specific.RigolDG1000</code> <code>method)</code> , 445	<code>get_full_status()</code> (py- <code>lablib.devices.interface.camera.ICamera</code> <code>method)</code> , 835
<code>get_full_status()</code> (py- <code>lablib.devices.AWG.specific.RSInstekAFG21000</code> <code>method)</code> , 434	<code>get_full_status()</code> (py- <code>lablib.devices.interface.camera.IExposureCamera</code> <code>method)</code> , 848
<code>get_full_status()</code> (py- <code>lablib.devices.AWG.specific.TektronixAFG1000</code> <code>method)</code> , 440	<code>get_full_status()</code> (py- <code>lablib.devices.interface.camera.IGrabberAttributeCamera</code> <code>method)</code> , 844
<code>get_full_status()</code> (py- <code>lablib.devices.Basler.pylon.BaslerPylonCamera</code> <code>method)</code> , 518	<code>get_full_status()</code> (py- <code>lablib.devices.interface.camera.IROICamera</code> <code>method)</code> , 853
<code>get_full_status()</code> (py- <code>lablib.devices.BitFlow.BitFlow.BitFlowCamera</code> <code>method)</code> , 529	<code>get_full_status()</code> (py- <code>lablib.devices.interface.stage.IMultiaxisStage</code> <code>method)</code> , 863
<code>get_full_status()</code> (py- <code>lablib.devices.BitFlow.BitFlow.BitFlowFrameGrabber</code> <code>method)</code> , 524	<code>get_full_status()</code> (py- <code>lablib.devices.interface.stage.IStage</code> <code>method)</code> , 862
<code>get_full_status()</code> (py- <code>lablib.devices.Conrad.base.RelayBoard</code> <code>method)</code> , 534	<code>get_full_status()</code> (py- <code>lablib.devices.KJL.base.KJL300</code> <code>method)</code> , 591
<code>get_full_status()</code> (py- <code>lablib.devices.Cryocon.base.Cryocon1x</code> <code>method)</code> , 537	<code>get_full_status()</code> (py- <code>lablib.devices.Lakeshore.base.Lakeshore218</code> <code>method)</code> , 596
<code>get_full_status()</code> (py-	<code>get_full_status()</code> (py-

`lablib.devices.Lakeshore.base.Lakeshore370`  
`method)`, 601  
`get_full_status()` (`py-`  
`lablib.devices.LaserQuantum.base.Finesse`  
`method)`, 605  
`get_full_status()` (`py-`  
`lablib.devices.Leybold.base.GenericITR`  
`method)`, 608  
`get_full_status()` (`py-`  
`lablib.devices.Leybold.base.ITR90` `method)`,  
610  
`get_full_status()` (`py-`  
`lablib.devices.LighthousePhotonics.base.SproutG`  
`method)`, 613  
`get_full_status()` (`py-`  
`lablib.devices.M2.solstis.Solstis` `method)`,  
619  
`get_full_status()` (`py-`  
`lablib.devices.Newport.picomotor.Picomotor8742`  
`method)`, 630  
`get_full_status()` (`py-`  
`lablib.devices.NI.daq.NIDAQ` `method)`, 626  
`get_full_status()` (`py-`  
`lablib.devices.Ophir.base.OphirDevice`  
`method)`, 639  
`get_full_status()` (`py-`  
`lablib.devices.Ophir.base.VegaPowerMeter`  
`method)`, 643  
`get_full_status()` (`py-`  
`lablib.devices.OZOptics.base.DD100` `method)`,  
635  
`get_full_status()` (`py-`  
`lablib.devices.OZOptics.base.EPC04` `method)`,  
637  
`get_full_status()` (`py-`  
`lablib.devices.OZOptics.base.OZOpticsDevice`  
`method)`, 632  
`get_full_status()` (`py-`  
`lablib.devices.OZOptics.base.TF100` `method)`,  
634  
`get_full_status()` (`py-`  
`lablib.devices.PCO.SC2.PCOSC2Camera`  
`method)`, 650  
`get_full_status()` (`py-`  
`lablib.devices.Pfeiffer.base.DPG202` `method)`,  
658  
`get_full_status()` (`py-`  
`lablib.devices.Pfeiffer.base.TPG260` `method)`,  
656  
`get_full_status()` (`py-`  
`lablib.devices.Photometrics.pvcam.PvcamCamera`  
`method)`, 666  
`get_full_status()` (`py-`  
`lablib.devices.PhotonFocus.PhotonFocus.IPhotonFocusCamera`  
`method)`, 675  
`get_full_status()` (`py-`  
`lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusBitFlowCamera`  
`method)`, 697  
`get_full_status()` (`py-`  
`lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusIMAQCamera`  
`method)`, 681  
`get_full_status()` (`py-`  
`lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusSiSoCamera`  
`method)`, 689  
`get_full_status()` (`py-`  
`lablib.devices.PhysikInstrumente.base.GenericPIController`  
`method)`, 704  
`get_full_status()` (`py-`  
`lablib.devices.PhysikInstrumente.base.PIE516`  
`method)`, 707  
`get_full_status()` (`py-`  
`lablib.devices.PrincetonInstruments.picam.PicamCamera`  
`method)`, 715  
`get_full_status()` (`py-`  
`lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareCamera`  
`method)`, 729  
`get_full_status()` (`py-`  
`lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareFrameGrabber`  
`method)`, 724  
`get_full_status()` (`py-`  
`lablib.devices.Sirah.Matisse.SirahMatisse`  
`method)`, 741  
`get_full_status()` (`py-`  
`lablib.devices.SmarAct.scu3d.SCU3D` `method)`,  
747  
`get_full_status()` (`py-`  
`lablib.devices.Tektronix.base.DPO2000`  
`method)`, 763  
`get_full_status()` (`py-`  
`lablib.devices.Tektronix.base.ITektronixScope`  
`method)`, 753  
`get_full_status()` (`py-`  
`lablib.devices.Tektronix.base.TDS2000`  
`method)`, 757  
`get_full_status()` (`py-`  
`lablib.devices.Thorlabs.elliptec.ElliptecMotor`  
`method)`, 780  
`get_full_status()` (`py-`  
`lablib.devices.Thorlabs.kinesis.BasicKinesisDevice`  
`method)`, 783  
`get_full_status()` (`py-`  
`lablib.devices.Thorlabs.kinesis.KinesisDevice`  
`method)`, 789  
`get_full_status()` (`py-`  
`lablib.devices.Thorlabs.kinesis.KinesisMotor`  
`method)`, 799  
`get_full_status()` (`py-`  
`lablib.devices.Thorlabs.kinesis.KinesisPiezoMotor`  
`method)`, 799

- method*), 803
- `get_full_status()` (*py-lablib.devices.Thorlabs.kinesis.KinesisQuadDetector* *method*), 808
- `get_full_status()` (*py-lablib.devices.Thorlabs.kinesis.MFF* *method*), 793
- `get_full_status()` (*py-lablib.devices.Thorlabs.serial.FW* *method*), 814
- `get_full_status()` (*py-lablib.devices.Thorlabs.serial.FWv1* *method*), 817
- `get_full_status()` (*py-lablib.devices.Thorlabs.serial.MDT69xA* *method*), 820
- `get_full_status()` (*py-lablib.devices.Thorlabs.serial.ThorlabsSerialInterface* *method*), 810
- `get_full_status()` (*py-lablib.devices.Thorlabs.TLCamera.ThorlabsTLCamera* *method*), 774
- `get_full_status()` (*py-lablib.devices.Toptica.ibeam.TopicalBeam* *method*), 825
- `get_full_status()` (*py-lablib.devices.Trinamic.base.TMCM1110* *method*), 829
- `get_full_status()` (*py-lablib.devices.uc480.uc480.UC480Camera* *method*), 870
- `get_full_web_status()` (*py-lablib.devices.M2.solstis.Solstis* *method*), 615
- `get_function()` (*py-lablib.devices.AWG.generic.GenericAWG* *method*), 403
- `get_function()` (*py-lablib.devices.AWG.specific.Agilent33220A* *method*), 416
- `get_function()` (*py-lablib.devices.AWG.specific.Agilent33500* *method*), 410
- `get_function()` (*py-lablib.devices.AWG.specific.InstekAFG2000* *method*), 428
- `get_function()` (*py-lablib.devices.AWG.specific.InstekAFG2225* *method*), 422
- `get_function()` (*py-lablib.devices.AWG.specific.RigolDG1000* *method*), 446
- `get_function()` (*py-lablib.devices.AWG.specific.RSInstekAFG21000* *method*), 434
- `get_function()` (*py-lablib.devices.AWG.specific.TektronixAFG1000* *method*), 440
- `get_gain_boost()` (*py-lablib.devices.uc480.uc480.UC480Camera* *method*), 868
- `get_gains()` (*py-lablib.devices.uc480.uc480.UC480Camera* *method*), 868
- `get_gate_polarity()` (*py-lablib.devices.AWG.generic.GenericAWG* *method*), 405
- `get_gate_polarity()` (*py-lablib.devices.AWG.specific.Agilent33220A* *method*), 416
- `get_gate_polarity()` (*py-lablib.devices.AWG.specific.Agilent33500* *method*), 410
- `get_gate_polarity()` (*py-lablib.devices.AWG.specific.InstekAFG2000* *method*), 428
- `get_gate_polarity()` (*py-lablib.devices.AWG.specific.InstekAFG2225* *method*), 422
- `get_gate_polarity()` (*py-lablib.devices.AWG.specific.RigolDG1000* *method*), 446
- `get_gate_polarity()` (*py-lablib.devices.AWG.specific.RSInstekAFG21000* *method*), 434
- `get_gate_polarity()` (*py-lablib.devices.AWG.specific.TektronixAFG1000* *method*), 440
- `get_gauge_control_settings()` (*py-lablib.devices.Pfeiffer.base.TPG260* *method*), 655
- `get_gauge_kind()` (*py-lablib.devices.Pfeiffer.base.TPG260* *method*), 655
- `get_gen_move_parameters()` (*py-lablib.devices.Thorlabs.kinesis.KinesisMotor* *method*), 797
- `get_general_input()` (*py-lablib.devices.Trinamic.base.TMCM1110* *method*), 827
- `get_genicam_info_xml()` (*py-lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusSiSoCamera* *method*), 690
- `get_genicam_info_xml()` (*py-lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareCamera* *method*), 729
- `get_genicam_info_xml()` (*py-lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareFrameGrabber* *method*), 722



<code>get_global_parameter()</code> <i>lablib.devices.Trinamic.base.TMCM1110</i> <i>method</i> ), 827	(py- <code>get_grabber_detector_size()</code> <i>lablib.devices.IMAQ.IMAQ.IMAQCamera</i> <i>method</i> ), 571	(py-
<code>get_global_speed()</code> <i>lablib.devices.Arcus.performax.Performax2EXStage</i> <i>method</i> ), 495	(py- <code>get_grabber_detector_size()</code> <i>lablib.devices.IMAQ.IMAQ.IMAQFrameGrabber</i> <i>method</i> ), 562	(py-
<code>get_global_speed()</code> <i>lablib.devices.Arcus.performax.Performax4EXStage</i> <i>method</i> ), 491	(py- <code>get_grabber_detector_size()</code> <i>lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusBitFlowCa</i> <i>method</i> ), 697	(py-
<code>get_grabber_attribute()</code> <i>lablib.devices.interface.camera.IGrabberAttributeCamera</i> <i>method</i> ), 842	(py- <code>get_grabber_detector_size()</code> <i>lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusIMAQCame</i> <i>method</i> ), 681	(py-
<code>get_grabber_attribute()</code> <i>lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusSiSoCamera</i> <i>method</i> ), 690	(py- <code>get_grabber_detector_size()</code> <i>lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusSiSoCame</i> <i>method</i> ), 690	(py-
<code>get_grabber_attribute()</code> <i>lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareCamera</i> <i>method</i> ), 729	(py- <code>get_grabber_detector_size()</code> <i>lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareCamera</i> <i>method</i> ), 729	(py-
<code>get_grabber_attribute()</code> <i>lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareFrameGrabber</i> <i>method</i> ), 724	(py- <code>get_grabber_detector_size()</code> <i>lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareFrameGrabb</i> <i>method</i> ), 722	(py-
<code>get_grabber_attribute_value()</code> <i>lablib.devices.AlliedVision.Bonito.BonitoIMAQCamera</i> <i>method</i> ), 459	(py- <code>get_grabber_roi()</code> <i>lablib.devices.AlliedVision.Bonito.BonitoIMAQCamera</i> <i>method</i> ), 459	(py-
<code>get_grabber_attribute_value()</code> <i>lablib.devices.IMAQ.IMAQ.IMAQCamera</i> <i>method</i> ), 571	(py- <code>get_grabber_roi()</code> <i>lablib.devices.BitFlow.BitFlow.BitFlowCamera</i> <i>method</i> ), 529	(py-
<code>get_grabber_attribute_value()</code> <i>lablib.devices.IMAQ.IMAQ.IMAQFrameGrabber</i> <i>method</i> ), 561	(py- <code>get_grabber_roi()</code> <i>lablib.devices.BitFlow.BitFlow.BitFlowFrameGrabber</i> <i>method</i> ), 522	(py-
<code>get_grabber_attribute_value()</code> <i>lablib.devices.interface.camera.IGrabberAttributeCamera</i> <i>method</i> ), 842	(py- <code>get_grabber_roi()</code> <i>lablib.devices.IMAQ.IMAQ.IMAQCamera</i> <i>method</i> ), 571	(py-
<code>get_grabber_attribute_value()</code> <i>lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusIMAQCamera</i> <i>method</i> ), 681	(py- <code>get_grabber_roi()</code> <i>lablib.devices.IMAQ.IMAQ.IMAQFrameGrabber</i> <i>method</i> ), 562	(py-
<code>get_grabber_attribute_value()</code> <i>lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusSiSoCamera</i> <i>method</i> ), 690	(py- <code>get_grabber_roi()</code> <i>lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusBitFlowCa</i> <i>method</i> ), 697	(py-
<code>get_grabber_attribute_value()</code> <i>lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareCamera</i> <i>method</i> ), 729	(py- <code>get_grabber_roi()</code> <i>lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusIMAQCame</i> <i>method</i> ), 681	(py-
<code>get_grabber_attribute_value()</code> <i>lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareFrameGrabber</i> <i>method</i> ), 724	(py- <code>get_grabber_roi()</code> <i>lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusSiSoCame</i> <i>method</i> ), 690	(py-
<code>get_grabber_detector_size()</code> <i>lablib.devices.AlliedVision.Bonito.BonitoIMAQCamera</i> <i>method</i> ), 459	(py- <code>get_grabber_roi()</code> <i>lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareCamera</i> <i>method</i> ), 729	(py-
<code>get_grabber_detector_size()</code> <i>lablib.devices.BitFlow.BitFlow.BitFlowCamera</i> <i>method</i> ), 529	(py- <code>get_grabber_roi()</code> <i>lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareFrameGrabb</i> <i>method</i> ), 722	(py-
<code>get_grabber_detector_size()</code> <i>lablib.devices.BitFlow.BitFlow.BitFlowFrameGrabber</i> <i>method</i> ), 522	(py- <code>get_grabber_roi_limits()</code> <i>lablib.devices.AlliedVision.Bonito.BonitoIMAQCamera</i> <i>method</i> ), 459	(py-

<code>get_grabber_roi_limits()</code>	(py- lablib.devices.BitFlow.BitFlow.BitFlowCamera method), 529	lablib.core.gui.value_handling.LabelValueHandler method), 271
<code>get_grabber_roi_limits()</code>	(py- lablib.devices.BitFlow.BitFlow.BitFlowFrameGrabber method), 522	<code>get_handler()</code> (py- lablib.core.gui.value_handling.LineEditValueHandler method), 270
<code>get_grabber_roi_limits()</code>	(py- lablib.devices.IMAQ.IMAQ.IMAQCamera method), 571	<code>get_handler()</code> (py- lablib.core.gui.value_handling.ProgressBarValueHandler method), 276
<code>get_grabber_roi_limits()</code>	(py- lablib.devices.IMAQ.IMAQ.IMAQFrameGrabber method), 563	<code>get_handler()</code> (py- lablib.core.gui.value_handling.PropertyValueHandler method), 267
<code>get_grabber_roi_limits()</code>	(py- lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusBitFlowCamera method), 697	<code>get_handler()</code> (py- lablib.core.gui.value_handling.PushButtonValueHandler method), 273
<code>get_grabber_roi_limits()</code>	(py- lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusIMAQCamera method), 682	<code>get_handler()</code> (py- lablib.core.gui.value_handling.StandardValueHandler method), 268
<code>get_grabber_roi_limits()</code>	(py- lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusSiSoCamera method), 690	<code>get_handler()</code> (py- lablib.core.gui.value_handling.ToolButtonValueHandler method), 274
<code>get_grabber_roi_limits()</code>	(py- lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareCamera method), 729	<code>get_handler()</code> (py- lablib.core.gui.value_handling.VirtualValueHandler method), 267
<code>get_grabber_roi_limits()</code>	(py- lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareFrameGrabber method), 722	<code>get_handler()</code> (py- lablib.core.gui.widgets.container.IQContainer method), 202
<code>get_gui_controller()</code> (in module lablib.core.thread.controller), 314	py- lablib.core.thread.threadprop), 320	<code>get_handler()</code> (py- lablib.core.gui.widgets.container.IQWidgetContainer method), 207
<code>get_gui_thread()</code> (in module lablib.core.thread.threadprop), 320	py- lablib.core.gui.value_handling), 281	<code>get_handler()</code> (py- lablib.core.gui.widgets.container.QContainer method), 203
<code>get_gui_values()</code> (in module lablib.core.gui.value_handling), 281	py- lablib.devices.Basler.pylon.BaslerPylonCamera.BaslerPylonCamera method), 516	<code>get_handler()</code> (py- lablib.core.gui.widgets.container.QDialogContainer method), 219
<code>get_handle()</code> (pylablib.devices.Basler.pylon.BaslerPylonCamera.BaslerPylonCamera method), 516	(py- lablib.core.gui.value_handling.CheckboxValueHandler method), 272	<code>get_handler()</code> (py- lablib.core.gui.widgets.container.QFrameContainer method), 215
<code>get_handler()</code>	(py- lablib.core.gui.value_handling.ComboBoxValueHandler method), 275	<code>get_handler()</code> (py- lablib.core.gui.widgets.container.QGroupBoxContainer method), 223
<code>get_handler()</code>	(py- lablib.core.gui.value_handling.GUIValues method), 278	<code>get_handler()</code> (py- lablib.core.gui.widgets.container.QScrollAreaContainer method), 230
<code>get_handler()</code>	(py- lablib.core.gui.value_handling.IBoolValueHandler method), 271	<code>get_handler()</code> (py- lablib.core.gui.widgets.container.QScrollAreaContainer.QContainer method), 227
<code>get_handler()</code>	(py- lablib.core.gui.value_handling.ISingleValueHandler method), 269	<code>get_handler()</code> (py- lablib.core.gui.widgets.container.QTabContainer method), 232
<code>get_handler()</code>	(py- lablib.core.gui.value_handling.IValueHandler method), 266	<code>get_handler()</code> (py- lablib.core.gui.widgets.container.QWidgetContainer method), 211
<code>get_handler()</code>	(py- lablib.core.gui.value_handling.IValueHandler method), 266	<code>get_handler()</code> (py- lablib.core.gui.widgets.container.QWidgetContainer method), 211

`pylablib.core.gui.widgets.param_table.ParamTable`  
`method()`, 250  
`get_handler()` (`py-`  
`pylablib.core.gui.widgets.param_table.StatusTable`  
`method()`, 258  
`get_hardware_id()` (`py-`  
`pylablib.devices.Attocube.anc350.ANC350`  
`method()`, 507  
`get_head_info()` (`py-`  
`pylablib.devices.Ophir.base.VegaPowerMeter`  
`method()`, 641  
`get_help()` (`pylablib.devices.PhysikInstrumente.base.GenericPICController`  
`method()`, 703  
`get_help()` (`pylablib.devices.PhysikInstrumente.base.PIE516`  
`method()`, 707  
`get_high_level()` (`py-`  
`pylablib.devices.Cryomagnetics.base.LM500`  
`method()`, 541  
`get_high_level()` (`py-`  
`pylablib.devices.Cryomagnetics.base.LM510`  
`method()`, 545  
`get_home_offset()` (`py-`  
`pylablib.devices.Thorlabs.elliptec.ElliptecMotor`  
`method()`, 779  
`get_homing_parameters()` (`py-`  
`pylablib.devices.Thorlabs.kinesis.KinesisMotor`  
`method()`, 797  
`get_horizontal_offset()` (`py-`  
`pylablib.devices.Tektronix.base.DPO2000`  
`method()`, 764  
`get_horizontal_offset()` (`py-`  
`pylablib.devices.Tektronix.base.ITektronixScope`  
`method()`, 750  
`get_horizontal_offset()` (`py-`  
`pylablib.devices.Tektronix.base.TDS2000`  
`method()`, 757  
`get_horizontal_span()` (`py-`  
`pylablib.devices.Tektronix.base.DPO2000`  
`method()`, 764  
`get_horizontal_span()` (`py-`  
`pylablib.devices.Tektronix.base.ITektronixScope`  
`method()`, 750  
`get_horizontal_span()` (`py-`  
`pylablib.devices.Tektronix.base.TDS2000`  
`method()`, 757  
`get_hsspeed()` (`py-`  
`pylablib.devices.Andor.AndorSDK2.AndorSDK2Camera`  
`method()`, 468  
`get_hsspeed_frequency()` (`py-`  
`pylablib.devices.Andor.AndorSDK2.AndorSDK2Camera`  
`method()`, 468  
`get_id()` (`pylablib.core.devio.SCPI.SCPIDevice`  
`method()`, 139  
`get_id()` (`pylablib.devices.AWG.generic.GenericAWG`  
`method()`, 406  
`get_id()` (`pylablib.devices.AWG.specific.Agilent33220A`  
`method()`, 416  
`get_id()` (`pylablib.devices.AWG.specific.Agilent33500`  
`method()`, 410  
`get_id()` (`pylablib.devices.AWG.specific.InstekAFG2000`  
`method()`, 428  
`get_id()` (`pylablib.devices.AWG.specific.InstekAFG2225`  
`method()`, 422  
`get_id()` (`pylablib.devices.AWG.specific.RigolDG1000`  
`method()`, 446  
`get_id()` (`pylablib.devices.AWG.specific.RSInstekAFG21000`  
`method()`, 434  
`get_id()` (`pylablib.devices.AWG.specific.TektronixAFG1000`  
`method()`, 440  
`get_id()` (`pylablib.devices.Cryocon.base.Cryocon1x`  
`method()`, 537  
`get_id()` (`pylablib.devices.Cryomagnetics.base.LM500`  
`method()`, 542  
`get_id()` (`pylablib.devices.Cryomagnetics.base.LM510`  
`method()`, 545  
`get_id()` (`pylablib.devices.Lakeshore.base.Lakeshore218`  
`method()`, 596  
`get_id()` (`pylablib.devices.Lakeshore.base.Lakeshore370`  
`method()`, 601  
`get_id()` (`pylablib.devices.Newport.picomotor.Picomotor8742`  
`method()`, 628  
`get_id()` (`pylablib.devices.PhysikInstrumente.base.GenericPICController`  
`method()`, 703  
`get_id()` (`pylablib.devices.PhysikInstrumente.base.PIE516`  
`method()`, 707  
`get_id()` (`pylablib.devices.Sirah.Matisse.SirahMatisse`  
`method()`, 741  
`get_id()` (`pylablib.devices.Tektronix.base.DPO2000`  
`method()`, 764  
`get_id()` (`pylablib.devices.Tektronix.base.ITektronixScope`  
`method()`, 753  
`get_id()` (`pylablib.devices.Tektronix.base.TDS2000`  
`method()`, 757  
`get_id()` (`pylablib.devices.Thorlabs.serial.FW`  
`method()`, 814  
`get_id()` (`pylablib.devices.Thorlabs.serial.FWv1`  
`method()`, 817  
`get_id()` (`pylablib.devices.Thorlabs.serial.MDT69xA`  
`method()`, 821  
`get_id()` (`pylablib.devices.Thorlabs.serial.ThorlabsSerialInterface`  
`method()`, 810  
`get_imag_part_ft()` (in module `py-`  
`pylablib.core.dataproc.fourier`), 120  
`get_image_indexing()` (`py-`  
`pylablib.devices.AlliedVision.Bonito.BonitoIMACCamera`  
`method()`, 459  
`get_image_indexing()` (`py-`  
`pylablib.devices.AlliedVision.Bonito.IBonitoCamera`

<i>method</i> ), 453	<i>method</i> ), 667
<code>get_image_indexing()</code> <i>lablib.devices.Andor.AndorSDK2.AndorSDK2Camera</i> <i>method</i> ), 474	<code>get_image_indexing()</code> <i>lablib.devices.PhotonFocus.PhotonFocus.IPhotonFocusCamera</i> <i>method</i> ), 675
<code>get_image_indexing()</code> <i>lablib.devices.Andor.AndorSDK3.AndorSDK3Camera</i> <i>method</i> ), 484	<code>get_image_indexing()</code> <i>lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusBitFlowCamera</i> <i>method</i> ), 698
<code>get_image_indexing()</code> <i>lablib.devices.Basler.pylon.BaslerPylonCamera</i> <i>method</i> ), 518	<code>get_image_indexing()</code> <i>lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusIMAQCamera</i> <i>method</i> ), 682
<code>get_image_indexing()</code> <i>lablib.devices.BitFlow.BitFlow.BitFlowCamera</i> <i>method</i> ), 529	<code>get_image_indexing()</code> <i>lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusSiSoCamera</i> <i>method</i> ), 690
<code>get_image_indexing()</code> <i>lablib.devices.BitFlow.BitFlow.BitFlowFrameGrabber</i> <i>method</i> ), 524	<code>get_image_indexing()</code> <i>lablib.devices.PrincetonInstruments.picam.PicamCamera</i> <i>method</i> ), 715
<code>get_image_indexing()</code> <i>lablib.devices.DCAM.DCAM.DCAMCamera</i> <i>method</i> ), 554	<code>get_image_indexing()</code> <i>lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareCamera</i> <i>method</i> ), 730
<code>get_image_indexing()</code> <i>lablib.devices.IMAQ.IMAQ.IMAQCamera</i> <i>method</i> ), 571	<code>get_image_indexing()</code> <i>lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareFrameGrabber</i> <i>method</i> ), 725
<code>get_image_indexing()</code> <i>lablib.devices.IMAQ.IMAQ.IMAQFrameGrabber</i> <i>method</i> ), 566	<code>get_image_indexing()</code> <i>lablib.devices.Thorlabs.TLCamera.ThorlabsTLCamera</i> <i>method</i> ), 774
<code>get_image_indexing()</code> <i>lablib.devices.IMAQdx.IMAQdx.EthernetIMAQdxCamera</i> <i>method</i> ), 586	<code>get_image_indexing()</code> <i>lablib.devices.uc480.uc480.UC480Camera</i> <i>method</i> ), 870
<code>get_image_indexing()</code> <i>lablib.devices.IMAQdx.IMAQdx.IMAQdxCamera</i> <i>method</i> ), 581	<code>get_image_mode_parameters()</code> <i>lablib.devices.Andor.AndorSDK2.AndorSDK2Camera</i> <i>method</i> ), 472
<code>get_image_indexing()</code> <i>lablib.devices.interface.camera.IAttributeCamera</i> <i>method</i> ), 839	<code>get_imported_modules()</code> (in module <i>pylablib.core.utils.module</i> ), 386
<code>get_image_indexing()</code> <i>lablib.devices.interface.camera.IBinROICamera</i> <i>method</i> ), 858	<code>get_index()</code> ( <i>pylablib.core.dataproc.table_wrap.Array1DWrapper</i> <i>method</i> ), 127
<code>get_image_indexing()</code> <i>lablib.devices.interface.camera.ICamera</i> <i>method</i> ), 832	<code>get_index()</code> ( <i>pylablib.core.dataproc.table_wrap.Array2DWrapper</i> <i>method</i> ), 131
<code>get_image_indexing()</code> <i>lablib.devices.interface.camera.IExposureCamera</i> <i>method</i> ), 848	<code>get_index()</code> ( <i>pylablib.core.dataproc.table_wrap.DataFrame2DWrapper</i> <i>method</i> ), 133
<code>get_image_indexing()</code> <i>lablib.devices.interface.camera.IGrabberAttributeCamera</i> <i>method</i> ), 844	<code>get_index()</code> ( <i>pylablib.core.dataproc.table_wrap.I1DWrapper</i> <i>method</i> ), 126
<code>get_image_indexing()</code> <i>lablib.devices.interface.camera.IROICamera</i> <i>method</i> ), 853	<code>get_index()</code> ( <i>pylablib.core.dataproc.table_wrap.I2DWrapper</i> <i>method</i> ), 129
<code>get_image_indexing()</code> <i>lablib.devices.PCO.SC2.PCOSC2Camera</i> <i>method</i> ), 650	<code>get_index()</code> ( <i>pylablib.core.dataproc.table_wrap.Series1DWrapper</i> <i>method</i> ), 128
<code>get_image_indexing()</code> <i>lablib.devices.Photometrics.pvcam.PvcamCamera</i> <i>method</i> ), 202	<code>get_index_values()</code> (py- <i>lablib.core.gui.widgets.combo_box.ComboBox</i> <i>method</i> ), 198
	<code>get_indicator()</code> (py- <i>lablib.core.gui.value_handling.GUIValues</i> <i>method</i> ), 280
	<code>get_indicator()</code> (py- <i>lablib.core.gui.widgets.container.IQContainer</i> <i>method</i> ), 202





615 *method*), 428

`get_layout_container()` (in module `pylablib.core.gui.utils`), 263 *method*), 422

`get_layout_shape()` (`pylablib.core.gui.widgets.container.IQWidgetContainer` *method*), 446

*method*), 207 `get_load()` (`pylablib.devices.AWG.specific.RigolDG1000` *method*), 434

`get_layout_shape()` (`pylablib.core.gui.widgets.container.QDialogContainer` *method*), 440

*method*), 219 `get_load()` (`pylablib.devices.AWG.specific.TektronixAFG1000` *method*), 386

`get_layout_shape()` (`pylablib.core.gui.widgets.container.QFrameContainer` *method*), 215 `get_loaded_package_modules()` (in module `pylablib.core.utils.module`), 388

`get_layout_shape()` (`pylablib.core.gui.widgets.container.QGroupBoxContainer` *method*), 223 `get_local_addr()` (in module `pylablib.core.utils.net`), 388

*method*), 223 `get_local_hostname()` (in module `pylablib.core.utils.net`), 388

`get_layout_shape()` (`pylablib.core.gui.widgets.container.QScrollAreaContainer` *method*), 227 `get_local_name()` (`pylablib.core.utils.net.ClientSocket` *method*), 390

*method*), 227 `get_location()` (in module `pylablib.core.fileio.location`), 188

`get_layout_shape()` (`pylablib.core.gui.widgets.container.QWidgetContainer` *method*), 211 `get_low_level()` (`pylablib.devices.Cryomagnetics.base.LM500` *method*), 541

`get_layout_shape()` (`pylablib.core.gui.widgets.layout_manager.IQLayoutManagedWidget` *method*), 239 `get_low_level()` (`pylablib.devices.Cryomagnetics.base.LM510` *method*), 545

`get_layout_shape()` (`pylablib.core.gui.widgets.layout_manager.QLayoutManagedWidget` *method*), 241 `get_mandatory_args()` (`pylablib.core.dataproc.callable.FunctionCallable` *method*), 106

`get_layout_shape()` (`pylablib.core.gui.widgets.param_table.ParamTable` *method*), 250 `get_mandatory_args()` (`pylablib.core.dataproc.callable.ICallable` *method*), 104

`get_layout_shape()` (`pylablib.core.gui.widgets.param_table.StatusTable` *method*), 258 `get_mandatory_args()` (`pylablib.core.dataproc.callable.JoinedCallable` *method*), 105

`get_level()` (`pylablib.devices.Cryomagnetics.base.LM500` *method*), 540 `get_mandatory_args()` (`pylablib.core.dataproc.callable.MethodCallable` *method*), 108

`get_level()` (`pylablib.devices.Cryomagnetics.base.LM510` *method*), 545 `get_mandatory_args()` (`pylablib.core.dataproc.callable.MultiplexedCallable` *method*), 105

`get_library_name()` (in module `pylablib.core.utils.module`), 386 `get_manual_output()` (`pylablib.devices.Thorlabs.kinesis.KinesisQuadDetector` *method*), 806

`get_library_path()` (in module `pylablib.core.utils.module`), 386 `get_matching_paths()` (`pylablib.core.utils.dictionary.Dictionary` *method*), 333

`get_limit_switch_parameters()` (`pylablib.devices.Thorlabs.kinesis.KinesisMotor` *method*), 797 `get_matching_paths()` (`pylablib.core.utils.dictionary.DictionaryPointer` *method*), 338

`get_limit_switches_parameters()` (`pylablib.devices.Trinamic.base.TMCM1110` *method*), 828 `get_matching_paths()` (`pylablib.core.utils.dictionary.FilterTree` *method*), 355

`get_load()` (`pylablib.devices.AWG.generic.GenericAWG` *method*), 403 `get_matching_paths()` (`pylablib.devices.AWG.specific.Agilent33220A` *method*), 416

`get_load()` (`pylablib.devices.AWG.specific.Agilent33500` *method*), 410 `get_matching_paths()` (`pylablib.devices.AWG.specific.InstekAFG2000` *method*), 400

- lablib.core.utils.dictionary.PrefixTree* method), 347
- get\_matching\_subtree()* (*pylablib.core.utils.dictionary.Dictionary* method), 334
- get\_matching\_subtree()* (*pylablib.core.utils.dictionary.DictionaryPointer* method), 338
- get\_matching\_subtree()* (*pylablib.core.utils.dictionary.FilterTree* method), 355
- get\_matching\_subtree()* (*pylablib.core.utils.dictionary.PrefixTree* method), 347
- get\_max\_attenuation()* (*pylablib.devices.OZOptics.base.DD100* method), 635
- get\_max\_gains()* (*pylablib.devices.uc480.uc480.UC480Camera* method), 868
- get\_max\_prefix()* (*pylablib.core.utils.dictionary.Dictionary* method), 327
- get\_max\_prefix()* (*pylablib.core.utils.dictionary.DictionaryPointer* method), 338
- get\_max\_prefix()* (*pylablib.core.utils.dictionary.FilterTree* method), 355
- get\_max\_prefix()* (*pylablib.core.utils.dictionary.PrefixTree* method), 347
- get\_max\_vsspeed()* (*pylablib.devices.Andor.AndorSDK2.AndorSDK2Camera* method), 468
- get\_measurement\_filter()* (*pylablib.devices.Pfeiffer.base.TPG260* method), 655
- get\_measurement\_interval()* (*pylablib.devices.HighFinesse.wlm.WLM* method), 559
- get\_metadata\_mode()* (*pylablib.devices.PCO.SC2.PCOSC2Camera* method), 649
- get\_method\_kind()* (in module *pylablib.core.gui.value\_handling*), 265
- get\_microstep\_resolution()* (*pylablib.devices.Trinamic.base.TMCM1110* method), 828
- get\_min\_attenuation()* (*pylablib.devices.OZOptics.base.DD100* method), 635
- get\_min\_shutter\_times()* (*pylablib.devices.Andor.AndorSDK2.AndorSDK2Camera* method), 469
- get\_missed\_frames\_status()* (*pylablib.devices.Andor.AndorSDK3.AndorSDK3Camera* method), 482
- get\_mode()* (*pylablib.devices.Attocube.anc300.ANC300* method), 503
- get\_mode()* (*pylablib.devices.Cryomagnetics.base.LM500* method), 540
- get\_mode()* (*pylablib.devices.Cryomagnetics.base.LM510* method), 545
- get\_mode()* (*pylablib.devices.OZOptics.base.EPC04* method), 637
- get\_mode\_parameters()* (*pylablib.devices.BitFlow.BitFlow.CameraFileEditor* method), 532
- get\_motor\_info()* (*pylablib.devices.Thorlabs.elliptec.ElliptecMotor* method), 779
- get\_motor\_type()* (*pylablib.devices.Newport.picomotor.Picomotor8742* method), 629
- get\_multi\_track\_mode\_parameters()* (*pylablib.devices.Andor.AndorSDK2.AndorSDK2Camera* method), 471
- get\_names()* (*pylablib.core.dataproc.table\_wrap.Array2DWrapper.Column* method), 130
- get\_names()* (*pylablib.core.dataproc.table\_wrap.DataFrame2DWrapper* method), 132
- get\_new\_frames\_range()* (*pylablib.devices.interface.camera.FrameCounter* method), 836
- get\_new\_images\_range()* (*pylablib.devices.AlliedVision.Bonito.BonitoIMAQCamera* method), 460
- get\_new\_images\_range()* (*pylablib.devices.AlliedVision.Bonito.IBonitoCamera* method), 454
- get\_new\_images\_range()* (*pylablib.devices.Andor.AndorSDK2.AndorSDK2Camera* method), 474
- get\_new\_images\_range()* (*pylablib.devices.Andor.AndorSDK3.AndorSDK3Camera* method), 484
- get\_new\_images\_range()* (*pylablib.devices.Basler.pylon.BaslerPylonCamera* method), 518
- get\_new\_images\_range()* (*pylablib.devices.BitFlow.BitFlow.BitFlowCamera* method), 529
- get\_new\_images\_range()* (*pylablib.devices.BitFlow.BitFlow.BitFlowFrameGrabber* method), 525
- get\_new\_images\_range()* (*pylablib.devices.DCAM.DCAM.DCAMCamera* method), 469

<i>method</i> ), 554		<i>method</i> ), 730	
get_new_images_range() lablib.devices.IMAQ.IMAQ.IMAQCamera <i>method</i> ), 571	(py-	get_new_images_range() lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareFrameGrabber <i>method</i> ), 725	(py-
get_new_images_range() lablib.devices.IMAQ.IMAQ.IMAQFrameGrabber <i>method</i> ), 566	(py-	get_new_images_range() lablib.devices.Thorlabs.TLCamera.ThorlabsTLCamera <i>method</i> ), 774	(py-
get_new_images_range() lablib.devices.IMAQdx.IMAQdx.EthernetIMAQdxCamera <i>method</i> ), 586	(py-	get_new_images_range() lablib.devices.uc480.uc480.UC480Camera <i>method</i> ), 870	(py-
get_new_images_range() lablib.devices.IMAQdx.IMAQdx.IMAQdxCamera <i>method</i> ), 581	(py-	get_noise_filter_mode() lablib.devices.PCO.SC2.PCOSC2Camera <i>method</i> ), 648	(py-
get_new_images_range() lablib.devices.interface.camera.IAttributeCamera <i>method</i> ), 839	(py-	get_number_of_channels() lablib.devices.Cryocon.base.Cryocon1x <i>method</i> ), 536	(py-
get_new_images_range() lablib.devices.interface.camera.IBinROICamera <i>method</i> ), 858	(py-	get_number_of_channels() lablib.devices.Thorlabs.kinesis.BasicKinesisDevice <i>method</i> ), 783	(py-
get_new_images_range() lablib.devices.interface.camera.ICamera <i>method</i> ), 833	(py-	get_number_of_channels() lablib.devices.Thorlabs.kinesis.KinesisDevice <i>method</i> ), 789	(py-
get_new_images_range() lablib.devices.interface.camera.IExposureCamera <i>method</i> ), 848	(py-	get_number_of_channels() lablib.devices.Thorlabs.kinesis.KinesisMotor <i>method</i> ), 799	(py-
get_new_images_range() lablib.devices.interface.camera.IGrabberAttributeCamera <i>method</i> ), 844	(py-	get_number_of_channels() lablib.devices.Thorlabs.kinesis.KinesisPiezoMotor <i>method</i> ), 804	(py-
get_new_images_range() lablib.devices.interface.camera.IROICamera <i>method</i> ), 853	(py-	get_number_of_channels() lablib.devices.Thorlabs.kinesis.KinesisQuadDetector <i>method</i> ), 808	(py-
get_new_images_range() lablib.devices.PCO.SC2.PCOSC2Camera <i>method</i> ), 650	(py-	get_number_of_channels() lablib.devices.Thorlabs.kinesis.MFF <i>method</i> ), 793	(py-
get_new_images_range() lablib.devices.Photometrics.pvcam.PvcamCamera <i>method</i> ), 667	(py-	get_oamp() lablib.devices.Andor.AndorSDK2.AndorSDK2Camera <i>method</i> ), 468	(py-
get_new_images_range() lablib.devices.PhotonFocus.PhotonFocus.IPhotonFocusCamera <i>method</i> ), 675	(py-	get_oamp_desc() lablib.devices.Andor.AndorSDK2.AndorSDK2Camera <i>method</i> ), 468	(py-
get_new_images_range() lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusBitFlowCamera <i>method</i> ), 698	(py-	get_offset() lablib.devices.Attocube.anc300.ANC300 <i>method</i> ), 504	(py-
get_new_images_range() lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusIMAQCamera <i>method</i> ), 682	(py-	get_offset() lablib.devices.Attocube.anc350.ANC350 <i>method</i> ), 508	(py-
get_new_images_range() lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusSiCCamera <i>method</i> ), 690	(py-	get_offset() lablib.devices.AWG.generic.GenericAWG <i>method</i> ), 404	(py-
get_new_images_range() lablib.devices.PrincetonInstruments.picam.PicamCamera <i>method</i> ), 715	(py-	get_offset() lablib.devices.AWG.specific.Agilent33220A <i>method</i> ), 416	(py-
get_new_images_range() lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareCamera <i>method</i> ), 730	(py-	get_offset() lablib.devices.AWG.specific.Agilent33500 <i>method</i> ), 410	(py-
		get_offset() lablib.devices.AWG.specific.InstekAFG2000 <i>method</i> ), 428	(py-
		get_offset() lablib.devices.AWG.specific.InstekAFG2225 <i>method</i> ), 420	(py-
		get_offset() lablib.devices.AWG.specific.RigolDG1000 <i>method</i> ), 420	(py-



method), 446

get\_offset() (pylablib.devices.AWG.specific.RSInstekAFG21000method), 612

method), 432

get\_offset() (pylablib.devices.AWG.specific.TektronixAFG1000method), 440

method), 825

get\_operation\_mode() (py- lablib.devices.Thorlabs.kinesis.KinesisQuadDetector method), 806

get\_output\_power() (py- lablib.devices.LighthousePhotonics.base.SproutG method), 612

get\_output\_range() (py- lablib.devices.AWG.generic.GenericAWG method), 404

get\_output\_range() (py- lablib.devices.AWG.specific.Agilent33220A method), 416

get\_output\_range() (py- lablib.devices.AWG.specific.Agilent33500 method), 410

get\_output\_range() (py- lablib.devices.AWG.specific.InstekAFG2000 method), 428

get\_output\_range() (py- lablib.devices.AWG.specific.InstekAFG2225 method), 422

get\_output\_range() (py- lablib.devices.AWG.specific.RigolDG1000 method), 446

get\_output\_range() (py- lablib.devices.AWG.specific.RSInstekAFG21000 method), 434

get\_output\_range() (py- lablib.devices.AWG.specific.TektronixAFG1000 method), 440

get\_output\_setpoint() (py- lablib.devices.LaserQuantum.base.Finesse method), 605

get\_output\_setpoint() (py- lablib.devices.LighthousePhotonics.base.SproutG method), 612

get\_output\_status() (py- lablib.devices.LaserQuantum.base.Finesse method), 604

get\_output\_trigger\_slope() (py- lablib.devices.AWG.generic.GenericAWG method), 405

get\_output\_trigger\_slope() (py- lablib.devices.AWG.specific.Agilent33220A method), 416

get\_output\_trigger\_slope() (py- lablib.devices.AWG.specific.Agilent33500 method), 411

get\_output\_trigger\_slope() (py- lablib.devices.AWG.specific.InstekAFG2000 method), 428

get\_output\_trigger\_slope() (py- lablib.devices.AWG.specific.InstekAFG2225 method), 422

get\_output\_trigger\_slope() (py- lablib.devices.AWG.specific.RigolDG1000 method), 446

get\_output\_trigger\_slope() (py- lablib.devices.AWG.specific.RSInstekAFG21000 method), 434

get\_output\_trigger\_slope() (py- lablib.devices.AWG.specific.TektronixAFG1000 method), 440

get\_output\_power() (py- lablib.devices.LaserQuantum.base.Finesse method), 605

get\_output\_power() (py- lablib.devices.LighthousePhotonics.base.SproutG method), 612

`pylablib.devices.AWG.specific.RigolDG1000`  
`method`), 446

`get_output_trigger_slope()` (`py-`  
`pylablib.devices.AWG.specific.RSInstekAFG21000`  
`method`), 434

`get_output_trigger_slope()` (`py-`  
`pylablib.devices.AWG.specific.TektronixAFG1000`  
`method`), 440

`get_package_version()` (`in module py-`  
`pylablib.core.utils.module`), 385

`get_path()` (`pylablib.core.fileio.location.LocationName`  
`method`), 183

`get_path()` (`pylablib.core.utils.dictionary.Dictionary`  
`method`), 332

`get_path()` (`pylablib.core.utils.dictionary.DictionaryPointer`  
`method`), 335

`get_path()` (`pylablib.core.utils.dictionary.FilterTree`  
`method`), 355

`get_path()` (`pylablib.core.utils.dictionary.PrefixTree`  
`method`), 347

`get_pcount()` (`pylablib.devices.Thorlabs.serial.FW`  
`method`), 813

`get_pcount()` (`pylablib.devices.Thorlabs.serial.FWv1`  
`method`), 816

`get_peakdet_kernel()` (`in module py-`  
`pylablib.core.dataproc.feature`), 109

`get_peer_args()` (`py-`  
`pylablib.core.utils.ipc.IIPChannel` `method`),  
383

`get_peer_args()` (`py-`  
`pylablib.core.utils.ipc.PipeIPChannel` `method`),  
383

`get_peer_args()` (`py-`  
`pylablib.core.utils.ipc.SharedMemIPChannel`  
`method`), 384

`get_peer_args()` (`py-`  
`pylablib.core.utils.ipc.SharedMemIPCTable`  
`method`), 385

`get_peer_name()` (`py-`  
`pylablib.core.utils.net.ClientSocket` `method`),  
390

`get_phase()` (`pylablib.devices.AWG.generic.GenericAWG`  
`method`), 404

`get_phase()` (`pylablib.devices.AWG.specific.Agilent33220A`  
`method`), 417

`get_phase()` (`pylablib.devices.AWG.specific.Agilent33500`  
`method`), 411

`get_phase()` (`pylablib.devices.AWG.specific.InstekAFG2000`  
`method`), 428

`get_phase()` (`pylablib.devices.AWG.specific.InstekAFG2225`  
`method`), 423

`get_phase()` (`pylablib.devices.AWG.specific.RigolDG1000`  
`method`), 446

`get_phase()` (`pylablib.devices.AWG.specific.RSInstekAFG21000`  
`method`), 434

`get_phase()` (`pylablib.devices.AWG.specific.TektronixAFG1000`  
`method`), 440

`get_pid_parameters()` (`py-`  
`pylablib.devices.Thorlabs.kinesis.KinesisQuadDetector`  
`method`), 806

`get_piezoet_ctl_status()` (`py-`  
`pylablib.devices.Sirah.Matisse.SirahMatisse`  
`method`), 738

`get_piezoet_drive_params()` (`py-`  
`pylablib.devices.Sirah.Matisse.SirahMatisse`  
`method`), 738

`get_piezoet_feedback_params()` (`py-`  
`pylablib.devices.Sirah.Matisse.SirahMatisse`  
`method`), 738

`get_piezoet_feedforward_params()` (`py-`  
`pylablib.devices.Sirah.Matisse.SirahMatisse`  
`method`), 738

`get_piezoet_position()` (`py-`  
`pylablib.devices.Sirah.Matisse.SirahMatisse`  
`method`), 738

`get_pixel_correction_parameters()` (`py-`  
`pylablib.devices.Thorlabs.TLCamera.ThorlabsTLCamera`  
`method`), 771

`get_pixel_distance()` (`py-`  
`pylablib.devices.Photometrics.pvcam.PvcamCamera`  
`method`), 663

`get_pixel_rate()` (`py-`  
`pylablib.devices.PCO.SC2.PCOS2Camera`  
`method`), 647

`get_pixel_rate()` (`py-`  
`pylablib.devices.uc480.uc480.UC480Camera`  
`method`), 867

`get_pixel_rates_range()` (`py-`  
`pylablib.devices.uc480.uc480.UC480Camera`  
`method`), 867

`get_pixel_size()` (`py-`  
`pylablib.devices.Andor.AndorSDK2.AndorSDK2Camera`  
`method`), 467

`get_pixel_size()` (`py-`  
`pylablib.devices.Photometrics.pvcam.PvcamCamera`  
`method`), 663

`get_pixel_size()` (`py-`  
`pylablib.devices.PrincetonInstruments.picam.PicamCamera`  
`method`), 713

`get_points_number()` (`py-`  
`pylablib.devices.Tektronix.base.DPO2000`  
`method`), 764

`get_points_number()` (`py-`  
`pylablib.devices.Tektronix.base.ITektronixScope`  
`method`), 751

`get_points_number()` (`py-`  
`pylablib.devices.Tektronix.base.TDS2000`  
`method`), 757

<code>get_port_index()</code> (in module <code>pylablib.devices.PhotonFocus.PhotonFocus</code> ), 670	<code>lablib.core.fileio.savefile.IBinaryOutputFileFormat</code> method), 193
<code>get_position()</code> (py- <code>lablib.devices.Arcus.performax.Performax2EXStage</code> method), 495	<code>get_preamble()</code> (py- <code>lablib.core.fileio.savefile.TableBinaryOutputFileFormat</code> method), 193
<code>get_position()</code> (py- <code>lablib.devices.Arcus.performax.Performax4EXStage</code> method), 490	<code>get_preamp()</code> (pylablib.devices.Andor.AndorSDK2.AndorSDK2Camera method), 468
<code>get_position()</code> (py- <code>lablib.devices.Arcus.performax.PerformaxDMXJS4Stage</code> method), 497	<code>get_preamp_gain()</code> (py- <code>lablib.devices.Andor.AndorSDK2.AndorSDK2Camera</code> method), 468
<code>get_position()</code> (py- <code>lablib.devices.Attocube.anc350.ANC350</code> method), 508	<code>get_precision()</code> (py- <code>lablib.devices.Attocube.anc350.ANC350</code> method), 508
<code>get_position()</code> (py- <code>lablib.devices.Attocube.anc350.ANC350</code> method), 508	<code>get_precision_mode()</code> (py- <code>lablib.devices.HighFinesse.wlm.WLM</code> method), 559
<code>get_position()</code> (py- <code>lablib.devices.Newport.picomotor.Picomotor8742</code> method), 629	<code>get_pressure()</code> (pylablib.devices.KJL.base.KJL300 method), 590
<code>get_position()</code> (py- <code>lablib.devices.PhysikInstrumente.base.PIE516</code> method), 706	<code>get_pressure()</code> (py- <code>lablib.devices.Leybold.base.GenericITR</code> method), 608
<code>get_position()</code> (py- <code>lablib.devices.Thorlabs.elliptec.ElliptecMotor</code> method), 780	<code>get_pressure()</code> (py- <code>lablib.devices.Leybold.base.ITR90</code> method), 610
<code>get_position()</code> (py- <code>lablib.devices.Thorlabs.kinesis.KinesisMotor</code> method), 796	<code>get_pressure()</code> (py- <code>lablib.devices.Pfeiffer.base.DPG202</code> method), 658
<code>get_position()</code> (py- <code>lablib.devices.Thorlabs.kinesis.KinesisPiezoMotor</code> method), 801	<code>get_pressure()</code> (py- <code>lablib.devices.Pfeiffer.base.TPG260</code> method), 655
<code>get_position()</code> (py- <code>lablib.devices.Thorlabs.serial.FW</code> method), 813	<code>get_prev_len()</code> (in module <code>py- lablib.core.dataproc.fourier</code> ), 117
<code>get_position()</code> (py- <code>lablib.devices.Thorlabs.serial.FWv1</code> method), 816	<code>get_probe_attenuation()</code> (py- <code>lablib.devices.Tektronix.base.DPO2000</code> method), 764
<code>get_position()</code> (py- <code>lablib.devices.Trinamic.base.TMCM1110</code> method), 827	<code>get_probe_attenuation()</code> (py- <code>lablib.devices.Tektronix.base.ITektronixScope</code> method), 751
<code>get_position_lower_limit()</code> (py- <code>lablib.devices.PhysikInstrumente.base.PIE516</code> method), 706	<code>get_probe_attenuation()</code> (py- <code>lablib.devices.Tektronix.base.TDS2000</code> method), 757
<code>get_position_upper_limit()</code> (py- <code>lablib.devices.PhysikInstrumente.base.PIE516</code> method), 706	<code>get_program_files_folder()</code> (in module <code>py- lablib.devices.utils.load_lib</code> ), 873
<code>get_power()</code> (pylablib.devices.Ophir.base.VegaPowerMeter method), 641	<code>get_progress()</code> (py- <code>lablib.core.thread.callsync.QCallResultSynchronizer</code> method), 281
<code>get_preamble()</code> (py- <code>lablib.core.fileio.dict_entry.ExternalNumpyDictionaryEntry</code> method), 175	<code>get_progress()</code> (py- <code>lablib.core.thread.callsync.QDirectResultSynchronizer</code> method), 282
<code>get_preamble()</code> (py- <code>lablib.core.fileio.dict_entry.IExternalFileDictionaryEntry</code> method), 175	<code>get_props()</code> (in module <code>pylablib.core.utils.general</code> ), 374
<code>get_preamble()</code> (py- <code>lablib.devices.HighFinesse.wlm.WLM</code> method), 559	<code>get_pulse_mode()</code> (py- <code>lablib.devices.HighFinesse.wlm.WLM</code> method), 559

`get_pulse_width()` (pylablib.devices.AWG.generic.GenericAWG method), 404  
`get_pulse_width()` (pylablib.devices.AWG.specific.Agilent33220A method), 417  
`get_pulse_width()` (pylablib.devices.AWG.specific.Agilent33500 method), 411  
`get_pulse_width()` (pylablib.devices.AWG.specific.InstekAFG2000 method), 428  
`get_pulse_width()` (pylablib.devices.AWG.specific.InstekAFG2225 method), 423  
`get_pulse_width()` (pylablib.devices.AWG.specific.RigolDG1000 method), 446  
`get_pulse_width()` (pylablib.devices.AWG.specific.RSInstekAFG21000 method), 434  
`get_pulse_width()` (pylablib.devices.AWG.specific.TektronixAFG1000 method), 438  
`get_python_folder()` (in module pylablib.core.utils.module), 386  
`get_ramp_symmetry()` (pylablib.devices.AWG.generic.GenericAWG method), 404  
`get_ramp_symmetry()` (pylablib.devices.AWG.specific.Agilent33220A method), 417  
`get_ramp_symmetry()` (pylablib.devices.AWG.specific.Agilent33500 method), 411  
`get_ramp_symmetry()` (pylablib.devices.AWG.specific.InstekAFG2000 method), 429  
`get_ramp_symmetry()` (pylablib.devices.AWG.specific.InstekAFG2225 method), 423  
`get_ramp_symmetry()` (pylablib.devices.AWG.specific.RigolDG1000 method), 446  
`get_ramp_symmetry()` (pylablib.devices.AWG.specific.RSInstekAFG21000 method), 434  
`get_ramp_symmetry()` (pylablib.devices.AWG.specific.TektronixAFG1000 method), 440  
`get_random_track_mode_parameters()` (pylablib.devices.Andor.AndorSDK2.AndorSDK2Camera method), 471  
`get_range()` (pylablib.devices.Ophir.base.VegaPowerMeter method), 478  
`get_range()` (pylablib.devices.Ophir.base.VegaPowerMeter method), 642  
`get_range_idx()` (pylablib.devices.Ophir.base.VegaPowerMeter method), 642  
`get_range_indices()` (in module pylablib.core.dataproc.utils), 137  
`get_range_info()` (pylablib.devices.Ophir.base.VegaPowerMeter method), 642  
`get_read_mode()` (pylablib.devices.Andor.AndorSDK2.AndorSDK2Camera method), 471  
`get_read_mode()` (pylablib.devices.HighFinesse.wlm.WLM method), 558  
`get_readings()` (pylablib.devices.Thorlabs.kinesis.KinesisQuadDetector method), 806  
`get_readout_mode()` (pylablib.devices.Photometrics.pvcam.PvcamCamera method), 663  
`get_readout_speed()` (pylablib.devices.DCAM.DCAM.DCAMCamera method), 551  
`get_readout_time()` (pylablib.devices.Andor.AndorSDK2.AndorSDK2Camera method), 471  
`get_readout_time()` (pylablib.devices.Photometrics.pvcam.PvcamCamera method), 664  
`get_real_part_ft()` (in module pylablib.core.dataproc.fourier), 120  
`get_refcell_position()` (pylablib.devices.Sirah.Matisse.SirahMatisse method), 739  
`get_refcell_waveform()` (pylablib.devices.Sirah.Matisse.SirahMatisse method), 737  
`get_refcell_waveform_params()` (pylablib.devices.Sirah.Matisse.SirahMatisse method), 739  
`get_reference_cavity_lock_status()` (pylablib.devices.M2.solstis.Solstis method), 617  
`get_region()` (in module pylablib.core.dataproc.image), 121  
`get_region_sum()` (in module pylablib.core.dataproc.image), 121  
`get_relative_position()` (in module pylablib.core.gui.utils), 264  
`get_relative_rectangle()` (in module pylablib.core.gui.utils), 264

`get_relay()` (`pylablib.devices.Conrad.base.RelayBoard` method), 534  
`get_relay_setpoints()` (`pylablib.devices.KJL.base.KJL300` method), 590  
`get_reload_order()` (in module `pylablib.core.utils.module`), 386  
`get_remote_hostname()` (in module `pylablib.core.utils.net`), 388  
`get_resistance()` (`pylablib.devices.Lakeshore.base.Lakeshore370` method), 599  
`get_roi()` (`pylablib.devices.AlliedVision.Bonito.BonitoIMAQCamera` method), 460  
`get_roi()` (`pylablib.devices.AlliedVision.Bonito.IBonitoCamera` method), 451  
`get_roi()` (`pylablib.devices.Andor.AndorSDK2.AndorSDK2Camera` method), 472  
`get_roi()` (`pylablib.devices.Andor.AndorSDK3.AndorSDK3Camera` method), 482  
`get_roi()` (`pylablib.devices.Basler.pylon.BaslerPylonCamera` method), 515  
`get_roi()` (`pylablib.devices.BitFlow.BitFlow.BitFlowCamera` method), 529  
`get_roi()` (`pylablib.devices.BitFlow.BitFlow.BitFlowFrameGrabber` method), 522  
`get_roi()` (`pylablib.devices.DCAM.DCAM.DCAMCamera` method), 552  
`get_roi()` (`pylablib.devices.IMAQ.IMAQ.IMAQCamera` method), 571  
`get_roi()` (`pylablib.devices.IMAQ.IMAQ.IMAQFrameGrabber` method), 562  
`get_roi()` (`pylablib.devices.IMAQdx.IMAQdx.EthernetIMAQdxCamera` method), 586  
`get_roi()` (`pylablib.devices.IMAQdx.IMAQdx.IMAQdxCamera` method), 579  
`get_roi()` (`pylablib.devices.interface.camera.IBinROICamera` method), 856  
`get_roi()` (`pylablib.devices.interface.camera.IROICamera` method), 851  
`get_roi()` (`pylablib.devices.PCO.SC2.PCOSC2Camera` method), 648  
`get_roi()` (`pylablib.devices.Photometrics.pvcam.PvcamCamera` method), 664  
`get_roi()` (`pylablib.devices.PhotonFocus.PhotonFocus.IBinROICamera` method), 673  
`get_roi()` (`pylablib.devices.PhotonFocus.PhotonFocus.PhotonFocusCamera` method), 698  
`get_roi()` (`pylablib.devices.PhotonFocus.PhotonFocus.PhotonFocusIMAQCamera` method), 682  
`get_roi()` (`pylablib.devices.PhotonFocus.PhotonFocus.PhotonFocusSiSocCamera` method), 690  
`get_roi()` (`pylablib.devices.PrincetonInstruments.picam.PicamCamera` method), 713  
`get_roi()` (`pylablib.devices.SiliconSoftware.fgrab.SiliconSoftwareCamera` method), 730  
`get_roi()` (`pylablib.devices.SiliconSoftware.fgrab.SiliconSoftwareFrameGrabber` method), 722  
`get_roi()` (`pylablib.devices.Thorlabs.TLCAmera.ThorlabsTLCAmera` method), 772  
`get_roi()` (`pylablib.devices.uc480.uc480.UC480Camera` method), 869  
`get_roi_limits()` (`pylablib.devices.AlliedVision.Bonito.BonitoIMAQCamera` method), 460  
`get_roi_limits()` (`pylablib.devices.AlliedVision.Bonito.IBonitoCamera` method), 451  
`get_roi_limits()` (`pylablib.devices.Andor.AndorSDK2.AndorSDK2Camera` method), 472  
`get_roi_limits()` (`pylablib.devices.Andor.AndorSDK3.AndorSDK3Camera` method), 483  
`get_roi_limits()` (`pylablib.devices.Basler.pylon.BaslerPylonCamera` method), 515  
`get_roi_limits()` (`pylablib.devices.BitFlow.BitFlow.BitFlowCamera` method), 529  
`get_roi_limits()` (`pylablib.devices.BitFlow.BitFlow.BitFlowFrameGrabber` method), 522  
`get_roi_limits()` (`pylablib.devices.DCAM.DCAM.DCAMCamera` method), 552  
`get_roi_limits()` (`pylablib.devices.IMAQ.IMAQ.IMAQCamera` method), 571  
`get_roi_limits()` (`pylablib.devices.IMAQ.IMAQ.IMAQFrameGrabber` method), 562  
`get_roi_limits()` (`pylablib.devices.IMAQdx.IMAQdx.EthernetIMAQdxCamera` method), 586  
`get_roi_limits()` (`pylablib.devices.IMAQdx.IMAQdx.IMAQdxCamera` method), 579  
`get_roi_limits()` (`pylablib.devices.interface.camera.IBinROICamera` method), 856  
`get_roi_limits()` (`pylablib.devices.interface.camera.IROICamera` method), 852  
`get_roi_limits()` (`pylablib.devices.PCO.SC2.PCOSC2Camera` method), 648  
`get_roi_limits()` (`pylablib.devices.Photometrics.pvcam.PvcamCamera` method), 664  
`get_roi_limits()` (`pylablib.devices.PhotonFocus.PhotonFocus.IBinROICamera` method), 673  
`get_roi_limits()` (`pylablib.devices.PhotonFocus.PhotonFocus.PhotonFocusCamera` method), 698  
`get_roi_limits()` (`pylablib.devices.PhotonFocus.PhotonFocus.PhotonFocusIMAQCamera` method), 682  
`get_roi_limits()` (`pylablib.devices.PhotonFocus.PhotonFocus.PhotonFocusSiSocCamera` method), 690  
`get_roi_limits()` (`pylablib.devices.PrincetonInstruments.picam.PicamCamera` method), 713



`lablib.devices.Photometrics.pvcam.PvcamCamera` `method`), 750  
`method`), 664 `get_selected_channel()` (py-  
`get_roi_limits()` (py- `lablib.devices.Tektronix.base.TDS2000`  
`lablib.devices.PhotonFocus.PhotonFocus.IPhotonFocusCamera` `method`), 757  
`method`), 673 `get_sensor_curve_index()` (py-  
`get_roi_limits()` (py- `lablib.devices.Lakeshore.base.Lakeshore218`  
`lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusBitFlowCamera` `method`), 594  
`method`), 698 `get_sensor_kind()` (py-  
`get_roi_limits()` (py- `lablib.devices.Cryocon.base.Cryocon1x`  
`lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusIMAQCamera` `method`), 536  
`method`), 682 `get_sensor_mode()` (py-  
`get_roi_limits()` (py- `lablib.devices.Thorlabs.serial.FW` `method`),  
`lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusSiSoCamera`  
`method`), 690 `get_sensor_power()` (py-  
`get_roi_limits()` (py- `lablib.devices.Lakeshore.base.Lakeshore370`  
`lablib.devices.PrincetonInstruments.picam.PicamCamera` `method`), 599  
`method`), 713 `get_sensor_reading()` (py-  
`get_roi_limits()` (py- `lablib.devices.Cryocon.base.Cryocon1x`  
`lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareCamera` `method`), 536  
`method`), 730 `get_sensor_reading()` (py-  
`get_roi_limits()` (py- `lablib.devices.Lakeshore.base.Lakeshore218`  
`lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareFrameGrabber` `method`), 594  
`method`), 722 `get_sensor_type()` (py-  
`get_roi_limits()` (py- `lablib.devices.Lakeshore.base.Lakeshore218`  
`lablib.devices.Thorlabs.TLCamera.ThorlabsTLCamera` `method`), 593  
`method`), 772 `get_sensor_voltage()` (py-  
`get_roi_limits()` (py- `lablib.devices.Attocube.anc350.ANC350`  
`lablib.devices.uc480.uc480.UC480Camera` `method`), 508  
`method`), 869 `get_serial_parameter()` (py-  
`get_roi_parameters()` (in module py- `lablib.devices.AlliedVision.Bonito.BonitoIMAQCamera`  
`lablib.devices.Photometrics.pvcam`), 669 `method`), 460  
`get_scale()` (`pylablib.devices.Thorlabs.elliptec.ElliptecMotor` `serial_parameter()` (py-  
`method`), 779 `lablib.devices.AlliedVision.Bonito.IBonitoCamera`  
`get_scale()` (`pylablib.devices.Thorlabs.kinesis.KinesisMotor` `method`), 451  
`method`), 795 `get_serial_params()` (py-  
`get_scale_units()` (py- `lablib.devices.AlliedVision.Bonito.BonitoIMAQCamera`  
`lablib.devices.Thorlabs.kinesis.KinesisMotor` `method`), 460  
`method`), 795 `get_serial_params()` (py-  
`get_scan_params()` (py- `lablib.devices.IMAQ.IMAQ.IMAQCamera`  
`lablib.devices.Sirah.Matisse.SirahMatisse` `method`), 572  
`method`), 740 `get_serial_params()` (py-  
`get_scan_position()` (py- `lablib.devices.IMAQ.IMAQ.IMAQFrameGrabber`  
`lablib.devices.Sirah.Matisse.SirahMatisse` `method`), 564  
`method`), 740 `get_serial_params()` (py-  
`get_scan_status()` (py- `lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusIMAQCamera`  
`lablib.devices.Sirah.Matisse.SirahMatisse` `method`), 682  
`method`), 740 `get_settings()` (py-  
`get_screenshot()` (in module py- `lablib.core.devio.comm_backend.ICommBackendWrapper`  
`lablib.core.gui.utils`), 264 `method`), 161  
`get_selected_channel()` (py- `get_settings()` (py-  
`lablib.devices.Tektronix.base.DPO2000` `lablib.core.devio.interface.IDevice` `method`),  
`method`), 764 162  
`get_selected_channel()` (py- `get_settings()` (py-  
`lablib.devices.Tektronix.base.ITektronixScope` `lablib.core.devio.SCPIDevice` `method`),

140 *method*), 434

*get\_settings()* (py- *get\_settings()* (py-  
*lablib.devices.AlliedVision.Bonito.BonitoIMAQCamera* *lablib.devices.AWG.specific.TektronixAFG1000*  
*method*), 460 *method*), 440

*get\_settings()* (py- *get\_settings()* (py-  
*lablib.devices.AlliedVision.Bonito.IBonitoCamera* *lablib.devices.Basler.pylon.BaslerPylonCamera*  
*method*), 454 *method*), 518

*get\_settings()* (py- *get\_settings()* (py-  
*lablib.devices.Andor.AndorSDK2.AndorSDK2Camera* *lablib.devices.BitFlow.BitFlow.BitFlowCamera*  
*method*), 474 *method*), 529

*get\_settings()* (py- *get\_settings()* (py-  
*lablib.devices.Andor.AndorSDK3.AndorSDK3Camera* *lablib.devices.BitFlow.BitFlow.BitFlowFrameGrabber*  
*method*), 484 *method*), 525

*get\_settings()* (py- *get\_settings()* (py-  
*lablib.devices.Arcus.performax.GenericPerformaxStage* *lablib.devices.Conrad.base.RelayBoard*  
*method*), 489 *method*), 535

*get\_settings()* (py- *get\_settings()* (py-  
*lablib.devices.Arcus.performax.Performax2EXStage* *lablib.devices.Cryocon.base.Cryocon1x*  
*method*), 495 *method*), 537

*get\_settings()* (py- *get\_settings()* (py-  
*lablib.devices.Arcus.performax.Performax4EXStage* *lablib.devices.Cryomagnetics.base.LM500*  
*method*), 492 *method*), 542

*get\_settings()* (py- *get\_settings()* (py-  
*lablib.devices.Arcus.performax.PerformaxDMXJSASStage* *lablib.devices.Cryomagnetics.base.LM510*  
*method*), 499 *method*), 545

*get\_settings()* (py- *get\_settings()* (py-  
*lablib.devices.Arduino.base.IArduinoDevice* *lablib.devices.DCAM.DCAM.DCAMCamera*  
*method*), 502 *method*), 554

*get\_settings()* (py- *get\_settings()* (py-  
*lablib.devices.Attocube.anc300.ANC300* *lablib.devices.HighFinesse.wlm.WLM* *method*),  
*method*), 505 560

*get\_settings()* (py- *get\_settings()* (py-  
*lablib.devices.Attocube.anc350.ANC350* *lablib.devices.IMAQ.IMAQ.IMAQCamera*  
*method*), 509 *method*), 572

*get\_settings()* (py- *get\_settings()* (py-  
*lablib.devices.AWG.generic.GenericAWG* *lablib.devices.IMAQ.IMAQ.IMAQFrameGrabber*  
*method*), 406 *method*), 566

*get\_settings()* (py- *get\_settings()* (py-  
*lablib.devices.AWG.specific.Agilent33220A* *lablib.devices.IMAQdx.IMAQdx.EthernetIMAQdxCamera*  
*method*), 417 *method*), 586

*get\_settings()* (py- *get\_settings()* (py-  
*lablib.devices.AWG.specific.Agilent33500* *lablib.devices.IMAQdx.IMAQdx.IMAQdxCamera*  
*method*), 411 *method*), 582

*get\_settings()* (py- *get\_settings()* (py-  
*lablib.devices.AWG.specific.InstekAFG2000* *lablib.devices.interface.camera.IAttributeCamera*  
*method*), 429 *method*), 839

*get\_settings()* (py- *get\_settings()* (py-  
*lablib.devices.AWG.specific.InstekAFG2225* *lablib.devices.interface.camera.IBinROICamera*  
*method*), 423 *method*), 858

*get\_settings()* (py- *get\_settings()* (py-  
*lablib.devices.AWG.specific.RigolDG1000* *lablib.devices.interface.camera.ICamera*  
*method*), 446 *method*), 835

*get\_settings()* (py- *get\_settings()* (py-  
*lablib.devices.AWG.specific.RSInstekAFG21000* *lablib.devices.interface.camera.IExposureCamera*  
*method*), 446 *method*), 835

*method*), 849

`get_settings()` (*py-lablib.devices.interface.camera.IGrabberAttributeCamera* *method*), 844

`get_settings()` (*py-lablib.devices.interface.camera.IROICamera* *method*), 853

`get_settings()` (*py-lablib.devices.interface.stage.IMultiaxisStage* *method*), 863

`get_settings()` (*py-lablib.devices.interface.stage.IStage* *method*), 862

`get_settings()` (*pylablib.devices.KJL.base.KJL300* *method*), 591

`get_settings()` (*py-lablib.devices.Lakeshore.base.Lakeshore218* *method*), 596

`get_settings()` (*py-lablib.devices.Lakeshore.base.Lakeshore370* *method*), 601

`get_settings()` (*py-lablib.devices.LaserQuantum.base.Finesse* *method*), 606

`get_settings()` (*py-lablib.devices.Leybold.base.GenericITR* *method*), 608

`get_settings()` (*py-lablib.devices.Leybold.base.ITR90* *method*), 610

`get_settings()` (*py-lablib.devices.LighthousePhotonics.base.SproutG* *method*), 613

`get_settings()` (*pylablib.devices.M2.solstis.Solstis* *method*), 620

`get_settings()` (*py-lablib.devices.Newport.picomotor.Picomotor8742* *method*), 631

`get_settings()` (*pylablib.devices.NI.daq.NIDAQ* *method*), 626

`get_settings()` (*py-lablib.devices.Ophir.base.OphirDevice* *method*), 639

`get_settings()` (*py-lablib.devices.Ophir.base.VegaPowerMeter* *method*), 643

`get_settings()` (*py-lablib.devices.OZOptics.base.DD100* *method*), 636

`get_settings()` (*py-lablib.devices.OZOptics.base.EPC04* *method*), 638

`get_settings()` (*py-lablib.devices.OZOptics.base.OZOpticsDevice* *method*), 632

`get_settings()` (*py-lablib.devices.OZOptics.base.TF100* *method*), 634

`get_settings()` (*py-lablib.devices.PCO.SC2.PCOSC2Camera* *method*), 650

`get_settings()` (*py-lablib.devices.Pfeiffer.base.DPG202* *method*), 658

`get_settings()` (*py-lablib.devices.Pfeiffer.base.TPG260* *method*), 656

`get_settings()` (*py-lablib.devices.Photometrics.pvcam.PvcamCamera* *method*), 667

`get_settings()` (*py-lablib.devices.PhotonFocus.PhotonFocus.IPhotonFocusCamera* *method*), 675

`get_settings()` (*py-lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusBitFlowCa* *method*), 698

`get_settings()` (*py-lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusIMAQCam* *method*), 682

`get_settings()` (*py-lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusSiSoCame* *method*), 690

`get_settings()` (*py-lablib.devices.PhysikInstrumente.base.GenericPIController* *method*), 704

`get_settings()` (*py-lablib.devices.PhysikInstrumente.base.PIE516* *method*), 707

`get_settings()` (*py-lablib.devices.PrincetonInstruments.picam.PicamCamera* *method*), 716

`get_settings()` (*py-lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareCamera* *method*), 730

`get_settings()` (*py-lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareFrameGrabb* *method*), 725

`get_settings()` (*py-lablib.devices.Sirah.Matisse.SirahMatisse* *method*), 741

`get_settings()` (*py-lablib.devices.SmarAct.scu3d.SCU3D* *method*), 747

`get_settings()` (*py-lablib.devices.Tektronix.base.DPO2000* *method*), 764

`get_settings()` (*py-lablib.devices.Tektronix.base.ITektronixScope* *method*), 764



`method)`, 753  
`get_settings()` (`py-lablib.devices.Tektronix.base.TDS2000` `method)`, 757  
`get_settings()` (`py-lablib.devices.Thorlabs.elliptec.ElliptecMotor` `method)`, 780  
`get_settings()` (`py-lablib.devices.Thorlabs.kinesis.BasicKinesisDevice` `method)`, 784  
`get_settings()` (`py-lablib.devices.Thorlabs.kinesis.KinesisDevice` `method)`, 789  
`get_settings()` (`py-lablib.devices.Thorlabs.kinesis.KinesisMotor` `method)`, 799  
`get_settings()` (`py-lablib.devices.Thorlabs.kinesis.KinesisPiezoMotor` `method)`, 804  
`get_settings()` (`py-lablib.devices.Thorlabs.kinesis.KinesisQuadDetector` `method)`, 808  
`get_settings()` (`py-lablib.devices.Thorlabs.kinesis.MFF` `method)`, 793  
`get_settings()` (`py-lablib.devices.Thorlabs.serial.FW` `method)`, 814  
`get_settings()` (`py-lablib.devices.Thorlabs.serial.FWv1` `method)`, 817  
`get_settings()` (`py-lablib.devices.Thorlabs.serial.MDT69xA` `method)`, 821  
`get_settings()` (`py-lablib.devices.Thorlabs.serial.ThorlabsSerialInterface` `method)`, 810  
`get_settings()` (`py-lablib.devices.Thorlabs.TLCamera.ThorlabsTLCamera` `method)`, 774  
`get_settings()` (`py-lablib.devices.Toptica.ibeam.TopicalBeam` `method)`, 825  
`get_settings()` (`py-lablib.devices.Trinamic.base.TMCM1110` `method)`, 829  
`get_settings()` (`py-lablib.devices.uc480.uc480.UC480Camera` `method)`, 871  
`get_shape()` (in `module py-lablib.core.utils.array_utils`), 321  
`get_shutter()` (`py-lablib.devices.Andor.AndorSDK2.AndorSDK2Camera` `method)`, 469  
`get_shutter()` (`py-lablib.devices.Andor.AndorSDK3.AndorSDK3Camera` `method)`, 480  
`get_shutter_parameters()` (`py-lablib.devices.Andor.AndorSDK2.AndorSDK2Camera` `method)`, 469  
`get_shutter_status()` (`py-lablib.devices.LaserQuantum.base.Finesse` `method)`, 605  
`get_shutter_status()` (`py-lablib.devices.LighthousePhotonics.base.SproutG` `method)`, 612  
`get_single_track_mode_parameters()` (`py-lablib.devices.Andor.AndorSDK2.AndorSDK2Camera` `method)`, 471  
`get_single_value()` (`py-lablib.core.gui.value_handling.CheckboxValueHandler` `method)`, 272  
`get_single_value()` (`py-lablib.core.gui.value_handling.ComboBoxValueHandler` `method)`, 274  
`get_single_value()` (`py-lablib.core.gui.value_handling.IBoolValueHandler` `method)`, 272  
`get_single_value()` (`py-lablib.core.gui.value_handling.ISingleValueHandler` `method)`, 269  
`get_single_value()` (`py-lablib.core.gui.value_handling.LabelValueHandler` `method)`, 270  
`get_single_value()` (`py-lablib.core.gui.value_handling.LineEditValueHandler` `method)`, 269  
`get_single_value()` (`py-lablib.core.gui.value_handling.ProgressBarValueHandler` `method)`, 275  
`get_single_value()` (`py-lablib.core.gui.value_handling.PushButtonValueHandler` `method)`, 273  
`get_single_value()` (`py-lablib.core.gui.value_handling.ToolButtonValueHandler` `method)`, 274  
`get_slowpiezo_ctl_params()` (`py-lablib.devices.Sirah.Matisse.SirahMatisse` `method)`, 739  
`get_slowpiezo_ctl_status()` (`py-lablib.devices.Sirah.Matisse.SirahMatisse` `method)`, 739  
`get_slowpiezo_position()` (`py-lablib.devices.Sirah.Matisse.SirahMatisse` `method)`, 739  
`get_software_version()` (`py-lablib.devices.Pfeiffer.base.DPG202` `method)`, 658

`get_speed_mode()` (pylablib.devices.Thorlabs.serial.FW method), 813  
`get_stage()` (pylablib.devices.Thorlabs.kinesis.KinesisMotor method), 795  
`get_state()` (pylablib.devices.Thorlabs.kinesis.MFF method), 791  
`get_stats()` (in module pylablib.core.thread.profile), 317  
`get_status()` (pylablib.devices.Andor.AndorSDK2.AndorSDK2Camera method), 467  
`get_status()` (pylablib.devices.Andor.AndorSDK3.AndorSDK3Camera method), 482  
`get_status()` (pylablib.devices.Arcus.performax.Performax2EXStage method), 495  
`get_status()` (pylablib.devices.Arcus.performax.Performax4EXStage method), 491  
`get_status()` (pylablib.devices.Arcus.performax.PerformaxDMXJSAStage method), 498  
`get_status()` (pylablib.devices.Attocube.anc350.ANC350 method), 507  
`get_status()` (pylablib.devices.Basler.pylon.BaslerPylonCamera method), 516  
`get_status()` (pylablib.devices.BitFlow.BitFlow.BitFlowCamera method), 527  
`get_status()` (pylablib.devices.BitFlow.BitFlow.BitFlowFrameGrabber method), 523  
`get_status()` (pylablib.devices.DCAM.DCAM.DCAMCamera method), 552  
`get_status()` (pylablib.devices.PhotonFocus.PhotonFocus.PhotonFocus method), 695  
`get_status()` (pylablib.devices.SmarAct.scu3d.SCU3D method), 746  
`get_status()` (pylablib.devices.Thorlabs.elliptec.ElliptecMotor method), 779  
`get_status()` (pylablib.devices.Thorlabs.kinesis.KinesisMotor method), 795  
`get_status()` (pylablib.devices.Thorlabs.kinesis.KinesisPiezoMotor method), 801  
`get_status()` (pylablib.devices.Thorlabs.kinesis.MFF method), 791  
`get_status()` (pylablib.devices.Thorlabs.TLCamera.ThorlabsTLCamera method), 771  
`get_status_line()` (in module pylablib.devices.PCO.SC2), 653  
`get_status_line_mode()` (pylablib.devices.PCO.SC2.PCOSC2Camera method), 649  
`get_status_line_position()` (in module pylablib.devices.PhotonFocus.PhotonFocus), 701  
`get_status_line_roi()` (in module pylablib.devices.interface.camera), 862  
`get_status_lines()` (in module pylablib.devices.AlliedVision.Bonito), 464  
`get_status_lines()` (in module pylablib.devices.PCO.SC2), 653  
`get_status_lines()` (in module pylablib.devices.PhotonFocus.PhotonFocus), 701  
`get_status_n()` (pylablib.devices.Arcus.performax.Performax2EXStage method), 495  
`get_status_n()` (pylablib.devices.Arcus.performax.Performax4EXStage method), 491  
`get_status_n()` (pylablib.devices.Arcus.performax.PerformaxDMXJSAStage method), 498  
`get_status_n()` (pylablib.devices.Attocube.anc350.ANC350 method), 507  
`get_status_n()` (pylablib.devices.Thorlabs.kinesis.KinesisMotor method), 795  
`get_status_n()` (pylablib.devices.Thorlabs.kinesis.KinesisPiezoMotor method), 801  
`get_status_n()` (pylablib.devices.Thorlabs.kinesis.MFF method), 795  
`get_string_filter()` (in module pylablib.core.utils.string), 398  
`get_sublayout()` (pylablib.core.gui.widgets.container.QDialogContainer method), 207  
`get_sublayout()` (pylablib.core.gui.widgets.container.QDialogContainer method), 219  
`get_sublayout()` (pylablib.core.gui.widgets.container.QFrameContainer method), 215  
`get_sublayout()` (pylablib.core.gui.widgets.container.QGroupBoxContainer method), 223  
`get_sublayout()` (pylablib.core.gui.widgets.container.QScrollAreaContainer method), 227  
`get_sublayout()` (pylablib.core.gui.widgets.container.QWidgetContainer method), 211  
`get_sublayout()` (pylablib.core.gui.widgets.layout\_manager.IQLayoutManagedWidget method), 239  
`get_sublayout()` (pylablib.core.gui.widgets.layout\_manager.QLayoutManagedWidget method), 241  
`get_sublayout()` (pylablib.core.gui.widgets.param\_table.ParamTable method), 464

<i>method</i> ), 250		655	
get_sublayout() lablib.core.gui.widgets.param_table.StatusTable <i>method</i> ), 259	(py-	get_switcher_mode() lablib.devices.HighFinesse.wlm.WLM <i>method</i> ), 558	(py-
get_sublayout_kind() lablib.core.gui.widgets.container.IQWidgetContainer <i>method</i> ), 207	(py-	get_system_info() lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusSiSoCamera <i>method</i> ), 690	(py-
get_sublayout_kind() lablib.core.gui.widgets.container.QDialogContainer <i>method</i> ), 219	(py-	get_system_info() lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareCamera <i>method</i> ), 730	(py-
get_sublayout_kind() lablib.core.gui.widgets.container.QFrameContainer <i>method</i> ), 215	(py-	get_system_info() lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareFrameGrabber <i>method</i> ), 722	(py-
get_sublayout_kind() lablib.core.gui.widgets.container.QGroupBoxContainer <i>method</i> ), 223	(py-	get_system_status() lablib.devices.M2.solstis.Solstis <i>method</i> ), 615	(py-
get_sublayout_kind() lablib.core.gui.widgets.container.QScrollAreaContainer.QComboBoxWidget <i>method</i> ), 227	(py-	get_table_line() lablib.devices.OxfordOxfordMicro.savefile.CSVTableOutputFileFormat <i>method</i> ), 191	(py-
get_sublayout_kind() lablib.core.gui.widgets.container.QWidgetContainer <i>method</i> ), 211	(py-	get_target_position() lablib.devices.Attocube.anc350.ANC350 <i>method</i> ), 508	(py-
get_sublayout_kind() lablib.core.gui.widgets.layout_manager.IQLayoutManagerWidget <i>method</i> ), 239	(py-	get_target_position() lablib.devices.PhysikInstrumente.base.PIE516 <i>method</i> ), 706	(py-
get_sublayout_kind() lablib.core.gui.widgets.layout_manager.QLayoutManagerWidget <i>method</i> ), 241	(py-	get_temperature() lablib.devices.Andor.AndorSDK2.AndorSDK2Camera <i>method</i> ), 467	(py-
get_sublayout_kind() lablib.core.gui.widgets.param_table.ParamTable <i>method</i> ), 250	(py-	get_temperature() lablib.devices.Andor.AndorSDK3.AndorSDK3Camera <i>method</i> ), 481	(py-
get_sublayout_kind() lablib.core.gui.widgets.param_table.StatusTable <i>method</i> ), 259	(py-	get_temperature() lablib.devices.Cryocon.base.Cryocon1x <i>method</i> ), 536	(py-
get_subsampling() lablib.devices.uc480.uc480.UC480Camera <i>method</i> ), 869	(py-	get_temperature() lablib.devices.Lakeshore.base.Lakeshore218 <i>method</i> ), 594	(py-
get_supported_baudrates() lablib.devices.Ophir.base.VegaPowerMeter <i>method</i> ), 642	(py-	get_temperature() lablib.devices.Lakeshore.base.Lakeshore370 <i>method</i> ), 599	(py-
get_supported_binning_modes() lablib.devices.Photometrics.pvcam.PvcamCamera <i>method</i> ), 665	(py-	get_temperature() lablib.devices.PCO.SC2.PCOSC2Camera <i>method</i> ), 646	(py-
get_supported_binning_modes() lablib.devices.uc480.uc480.UC480Camera <i>method</i> ), 869	(py-	get_temperature() lablib.devices.Photometrics.pvcam.PvcamCamera <i>method</i> ), 663	(py-
get_supported_subsampling_modes() lablib.devices.uc480.uc480.UC480Camera <i>method</i> ), 868	(py-	get_temperature_range() lablib.devices.Andor.AndorSDK2.AndorSDK2Camera <i>method</i> ), 467	(py-
get_switch_settings() lablib.devices.Pfeiffer.base.TPG260 <i>method</i> ), 655	(py-	get_temperature_setpoint() lablib.devices.Andor.AndorSDK2.AndorSDK2Camera <i>method</i> ), 467	(py-
get_switch_status() lablib.devices.Pfeiffer.base.TPG260 <i>method</i> ),	(py-	get_temperature_setpoint() lablib.devices.Andor.AndorSDK3.AndorSDK3Camera	(py-

method), 481

get\_temperature\_setpoint() (py-lablib.devices.Photometrics.pvcam.PvcamCamera method), 663

get\_temperature\_status() (py-lablib.devices.Andor.AndorSDK2.AndorSDK2Camera method), 467

get\_temperatures() (py-lablib.devices.LaserQuantum.base.Finesse method), 604

get\_temperatures() (py-lablib.devices.Toptica.ibeam.TopticalBeam method), 825

get\_terascan\_status() (py-lablib.devices.M2.solstis.Solstis method), 618

get\_thinet\_ctl\_params() (py-lablib.devices.Sirah.Matisse.SirahMatisse method), 738

get\_thinet\_ctl\_status() (py-lablib.devices.Sirah.Matisse.SirahMatisse method), 738

get\_thinet\_error\_signal() (py-lablib.devices.Sirah.Matisse.SirahMatisse method), 738

get\_thinet\_power() (py-lablib.devices.Sirah.Matisse.SirahMatisse method), 736

get\_timeout() (py-lablib.core.devio.comm\_backend.FT232DeviceBackend method), 151

get\_timeout() (py-lablib.core.devio.comm\_backend.IDeviceCommBackend method), 144

get\_timeout() (py-lablib.core.devio.comm\_backend.NetworkDeviceBackend method), 153

get\_timeout() (py-lablib.core.devio.comm\_backend.PyUSBDeviceBackend method), 156

get\_timeout() (py-lablib.core.devio.comm\_backend.RecordedDeviceBackend method), 159

get\_timeout() (py-lablib.core.devio.comm\_backend.SerialDeviceBackend method), 148

get\_timeout() (py-lablib.core.devio.comm\_backend.VisaDeviceBackend method), 146

get\_timeout() (pylablib.core.utils.net.ClientSocket method), 389

get\_timestamp\_clock\_frequency() (py-lablib.devices.Thorlabs.TLCamera.ThorlabsTLCamera method), 771

get\_top\_parent() (in module py-lablib.core.gui.utils), 263

get\_transfer\_info() (py-lablib.devices.DCAM.DCAM.DCAMCamera method), 553

get\_trigger\_input() (py-lablib.devices.Attocube.anc300.ANC300 method), 504

get\_trigger\_interleave() (py-lablib.devices.PhotonFocus.PhotonFocus.IPhotonFocusCamera method), 673

get\_trigger\_interleave() (py-lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusBitFlowCamera method), 698

get\_trigger\_interleave() (py-lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusIMACamera method), 682

get\_trigger\_interleave() (py-lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusSiSoCamera method), 691

get\_trigger\_level() (py-lablib.devices.Tektronix.base.DPO2000 method), 764

get\_trigger\_level() (py-lablib.devices.Tektronix.base.ITektronixScope method), 750

get\_trigger\_level() (py-lablib.devices.Tektronix.base.TDS2000 method), 757

get\_trigger\_level\_limits() (py-lablib.devices.Andor.AndorSDK2.AndorSDK2Camera method), 469

get\_trigger\_mode() (py-lablib.devices.Andor.AndorSDK2.AndorSDK2Camera method), 469

get\_trigger\_mode() (py-lablib.devices.Andor.AndorSDK3.AndorSDK3Camera method), 480

get\_trigger\_mode() (py-lablib.devices.DCAM.DCAM.DCAMCamera method), 551

get\_trigger\_mode() (py-lablib.devices.PCO.SC2.PCOSC2Camera method), 646

get\_trigger\_mode() (py-lablib.devices.Photometrics.pvcam.PvcamCamera method), 664

get\_trigger\_mode() (py-lablib.devices.Tektronix.base.DPO2000 method), 764

get\_trigger\_mode() (py-lablib.devices.Tektronix.base.ITektronixScope method), 750

get\_trigger\_mode() (py-

<i>lablib.devices.Tektronix.base.TDS2000</i> <i>method</i> ), 757	<i>lablib.devices.AWG.specific.RSInstekAFG21000</i> <i>method</i> ), 435
<i>get_trigger_mode()</i> ( <i>py-</i> <i>lablib.devices.Thorlabs.serial.FW</i> method), 813	<i>get_trigger_source()</i> ( <i>py-</i> <i>lablib.devices.AWG.specific.TektronixAFG1000</i> method), 441
<i>get_trigger_mode()</i> ( <i>py-</i> <i>lablib.devices.Thorlabs.serial.FWv1</i> method), 816	<i>get_trigger_state()</i> ( <i>py-</i> <i>lablib.devices.Tektronix.base.DPO2000</i> method), 764
<i>get_trigger_mode()</i> ( <i>py-</i> <i>lablib.devices.Thorlabs.TLCamera.ThorlabsTLCamera</i> method), 771	<i>get_trigger_state()</i> ( <i>py-</i> <i>lablib.devices.Tektronix.base.ITektronixScope</i> method), 750
<i>get_trigger_slope()</i> ( <i>py-</i> <i>lablib.devices.AWG.generic.GenericAWG</i> method), 405	<i>get_trigger_state()</i> ( <i>py-</i> <i>lablib.devices.Tektronix.base.TDS2000</i> method), 758
<i>get_trigger_slope()</i> ( <i>py-</i> <i>lablib.devices.AWG.specific.Agilent33220A</i> method), 417	<i>get_type()</i> ( <i>pylablib.core.dataproc.table_wrap.Array1DWrapper</i> method), 127
<i>get_trigger_slope()</i> ( <i>py-</i> <i>lablib.devices.AWG.specific.Agilent33500</i> method), 411	<i>get_type()</i> ( <i>pylablib.core.dataproc.table_wrap.Array2DWrapper</i> method), 131
<i>get_trigger_slope()</i> ( <i>py-</i> <i>lablib.devices.AWG.specific.InstekAFG2000</i> method), 429	<i>get_type()</i> ( <i>pylablib.core.dataproc.table_wrap.DataFrame2DWrapper</i> method), 133
<i>get_trigger_slope()</i> ( <i>py-</i> <i>lablib.devices.AWG.specific.InstekAFG2225</i> method), 423	<i>get_type()</i> ( <i>pylablib.core.dataproc.table_wrap.I1DWrapper</i> method), 126
<i>get_trigger_slope()</i> ( <i>py-</i> <i>lablib.devices.AWG.specific.RigolDG1000</i> method), 446	<i>get_type()</i> ( <i>pylablib.core.dataproc.table_wrap.I2DWrapper</i> method), 129
<i>get_trigger_slope()</i> ( <i>py-</i> <i>lablib.devices.AWG.specific.RSInstekAFG21000</i> method), 435	<i>get_type()</i> ( <i>pylablib.core.dataproc.table_wrap.IGenWrapper</i> method), 125
<i>get_trigger_slope()</i> ( <i>py-</i> <i>lablib.devices.AWG.specific.TektronixAFG1000</i> method), 441	<i>get_type()</i> ( <i>pylablib.core.dataproc.table_wrap.Series1DWrapper</i> method), 128
<i>get_trigger_source()</i> ( <i>py-</i> <i>lablib.devices.AWG.generic.GenericAWG</i> method), 405	<i>get_type()</i> ( <i>pylablib.devices.Cryomagnetics.base.LM500</i> method), 540
<i>get_trigger_source()</i> ( <i>py-</i> <i>lablib.devices.AWG.specific.Agilent33220A</i> method), 417	<i>get_type()</i> ( <i>pylablib.devices.Cryomagnetics.base.LM510</i> method), 545
<i>get_trigger_source()</i> ( <i>py-</i> <i>lablib.devices.AWG.specific.Agilent33500</i> method), 411	<i>get_units()</i> ( <i>pylablib.devices.Leybold.base.GenericITR</i> method), 608
<i>get_trigger_source()</i> ( <i>py-</i> <i>lablib.devices.AWG.specific.InstekAFG2000</i> method), 429	<i>get_units()</i> ( <i>pylablib.devices.Leybold.base.ITR90</i> method), 610
<i>get_trigger_source()</i> ( <i>py-</i> <i>lablib.devices.AWG.specific.InstekAFG2225</i> method), 423	<i>get_units()</i> ( <i>pylablib.devices.Ophir.base.VegaPowerMeter</i> method), 641
<i>get_trigger_source()</i> ( <i>py-</i> <i>lablib.devices.AWG.specific.RigolDG1000</i> method), 447	<i>get_units()</i> ( <i>pylablib.devices.Pfeiffer.base.TPG260</i> method), 654
<i>get_trigger_source()</i> ( <i>py-</i>	<i>get_update()</i> ( <i>pylablib.devices.Leybold.base.GenericITR</i> method), 607
	<i>get_update()</i> ( <i>pylablib.devices.Leybold.base.ITR90</i> method), 610
	<i>get_usb_device_info()</i> ( <i>in module py-</i> <i>lablib.devices.Arcus.performax</i> ), 488
	<i>get_usb_devices_number()</i> ( <i>in module py-</i> <i>lablib.devices.Attocube.anc350</i> ), 506
	<i>get_usb_devices_number()</i> ( <i>in module py-</i> <i>lablib.devices.Newport.picomotor</i> ), 627
	<i>get_value()</i> ( <i>in module py-</i> <i>lablib.core.utils.ctypes_wrap</i> ), 321
	<i>get_value()</i> ( <i>pylablib.core.gui.value_handling.CheckboxValueHandler</i> method), 273



`get_value()` (pylablib.core.gui.value\_handling.ComboBoxValueHandler (pylablib.core.gui.widgets.container.QWidgetContainer method), 275  
`get_value()` (pylablib.core.gui.value\_handling.GUIValues method), 280  
`get_value()` (pylablib.core.gui.value\_handling.IBoolValueHandler (pylablib.core.gui.widgets.edit.TextEdit method), 272  
`get_value()` (pylablib.core.gui.value\_handling.IIndicatorHandler (pylablib.core.gui.widgets.label.EnumLabel method), 276  
`get_value()` (pylablib.core.gui.value\_handling.ISingleValueHandler (pylablib.core.gui.widgets.label.NumLabel method), 269  
`get_value()` (pylablib.core.gui.value\_handling.IValueHandler (pylablib.core.gui.widgets.label.TextLabel method), 265  
`get_value()` (pylablib.core.gui.value\_handling.LabelIndicatorHandler (pylablib.core.gui.widgets.param\_table.ParamTable method), 277  
`get_value()` (pylablib.core.gui.value\_handling.LabelValueHandler (pylablib.core.gui.widgets.param\_table.StatusTable method), 271  
`get_value()` (pylablib.core.gui.value\_handling.LineEditValueHandler (pylablib.core.thread.callsync.QCallResultSynchronizer method), 270  
`get_value()` (pylablib.core.gui.value\_handling.ProgressBarValueHandler (pylablib.core.thread.callsync.QDirectResultSynchronizer method), 276  
`get_value()` (pylablib.core.gui.value\_handling.PropertyValueHandler (pylablib.core.thread.synchronizing.QThreadNotifier method), 267  
`get_value()` (pylablib.core.gui.value\_handling.PushButtonValueHandler (pylablib.devices.Andor.AndorSDK3.AndorSDK3Attribute method), 273  
`get_value()` (pylablib.core.gui.value\_handling.StandardIndicatorHandler (pylablib.devices.Attocube.anc350.ANC350 method), 277  
`get_value()` (pylablib.core.gui.value\_handling.StandardValueHandler (pylablib.devices.Basler.pylon.BaslerPylonAttribute method), 268  
`get_value()` (pylablib.core.gui.value\_handling.ToolButtonValueHandler (pylablib.devices.DCAM.DCAM.DCAMAttribute method), 274  
`get_value()` (pylablib.core.gui.value\_handling.VirtualValueHandler (pylablib.devices.IMAQdx.IMAQdx.IMAQdxAttribute method), 266  
`get_value()` (pylablib.core.gui.widgets.button.ToggleButton method), 198  
`get_value()` (pylablib.core.gui.widgets.combo\_box.ComboBox method), 199  
`get_value()` (pylablib.core.gui.widgets.container.IQContainer method), 202  
`get_value()` (pylablib.core.gui.widgets.container.IQWidgetContainer method), 207  
`get_value()` (pylablib.core.gui.widgets.container.QContainer method), 204  
`get_value()` (pylablib.core.gui.widgets.container.QDialogContainer method), 219  
`get_value()` (pylablib.core.gui.widgets.container.QFrameContainer method), 215  
`get_value()` (pylablib.core.gui.widgets.container.QGroupBoxContainer method), 223  
`get_value()` (pylablib.core.gui.widgets.container.QScrollAreaContainer method), 230  
`get_value()` (pylablib.core.gui.widgets.container.QScrollAreaContainer method), 227  
`get_value()` (pylablib.core.gui.widgets.container.QTabContainer method), 232  
`get_value()` (pylablib.core.gui.value\_handling.CheckboxValueHandler (pylablib.core.gui.value\_handling.ComboBoxValueHandler method), 275  
`get_value()` (pylablib.core.gui.value\_handling.GUIValues method), 280  
`get_value()` (pylablib.core.gui.value\_handling.IBoolValueHandler (pylablib.core.gui.widgets.edit.TextEdit method), 272  
`get_value()` (pylablib.core.gui.value\_handling.IIndicatorHandler (pylablib.core.gui.widgets.label.EnumLabel method), 276  
`get_value()` (pylablib.core.gui.value\_handling.ISingleValueHandler (pylablib.core.gui.widgets.label.NumLabel method), 269  
`get_value()` (pylablib.core.gui.value\_handling.IValueHandler (pylablib.core.gui.widgets.label.TextLabel method), 265  
`get_value()` (pylablib.core.gui.value\_handling.LabelIndicatorHandler (pylablib.core.gui.widgets.param\_table.ParamTable method), 277  
`get_value()` (pylablib.core.gui.value\_handling.LabelValueHandler (pylablib.core.gui.widgets.param\_table.StatusTable method), 271  
`get_value()` (pylablib.core.gui.value\_handling.LineEditValueHandler (pylablib.core.thread.callsync.QCallResultSynchronizer method), 270  
`get_value()` (pylablib.core.gui.value\_handling.ProgressBarValueHandler (pylablib.core.thread.callsync.QDirectResultSynchronizer method), 276  
`get_value()` (pylablib.core.gui.value\_handling.PropertyValueHandler (pylablib.core.thread.synchronizing.QThreadNotifier method), 267  
`get_value()` (pylablib.core.gui.value\_handling.PushButtonValueHandler (pylablib.devices.Andor.AndorSDK3.AndorSDK3Attribute method), 273  
`get_value()` (pylablib.core.gui.value\_handling.StandardIndicatorHandler (pylablib.devices.Attocube.anc350.ANC350 method), 277  
`get_value()` (pylablib.core.gui.value\_handling.StandardValueHandler (pylablib.devices.Basler.pylon.BaslerPylonAttribute method), 268  
`get_value()` (pylablib.core.gui.value\_handling.ToolButtonValueHandler (pylablib.devices.DCAM.DCAM.DCAMAttribute method), 274  
`get_value()` (pylablib.core.gui.value\_handling.VirtualValueHandler (pylablib.devices.IMAQdx.IMAQdx.IMAQdxAttribute method), 266  
`get_value()` (pylablib.devices.Pfeiffer.base.DPG202 method), 657  
`get_value()` (pylablib.devices.Photometrics.pvcam.PvcamAttribute method), 661  
`get_value()` (pylablib.devices.PhotonFocus.PhotonFocus.PFCamAttribute method), 671  
`get_value()` (pylablib.devices.PrincetonInstruments.picam.PicamAttribute method), 711  
`get_value()` (pylablib.devices.SiliconSoftware.fgrab.FGrabAttribute method), 720  
`get_value()` (pylablib.core.gui.value\_handling.IQContainer method), 202  
`get_value()` (pylablib.core.gui.value\_handling.IQWidgetContainer method), 207  
`get_value()` (pylablib.core.gui.value\_handling.QContainer method), 204  
`get_value()` (pylablib.core.gui.value\_handling.QDialogContainer method), 219  
`get_value()` (pylablib.core.gui.value\_handling.QFrameContainer method), 215  
`get_value()` (pylablib.core.gui.value\_handling.QGroupBoxContainer method), 223  
`get_value()` (pylablib.core.gui.value\_handling.QScrollAreaContainer method), 230  
`get_value()` (pylablib.core.gui.value\_handling.QScrollAreaContainer method), 227  
`get_value()` (pylablib.core.gui.value\_handling.QTabContainer method), 232  
`get_value()` (pylablib.core.gui.value\_handling.CheckboxValueHandler (pylablib.core.gui.value\_handling.ComboBoxValueHandler method), 275  
`get_value()` (pylablib.core.gui.value\_handling.GUIValues method), 280  
`get_value()` (pylablib.core.gui.value\_handling.IBoolValueHandler (pylablib.core.gui.widgets.edit.TextEdit method), 272  
`get_value()` (pylablib.core.gui.value\_handling.IIndicatorHandler (pylablib.core.gui.widgets.label.EnumLabel method), 276  
`get_value()` (pylablib.core.gui.value\_handling.ISingleValueHandler (pylablib.core.gui.widgets.label.NumLabel method), 269  
`get_value()` (pylablib.core.gui.value\_handling.IValueHandler (pylablib.core.gui.widgets.label.TextLabel method), 265  
`get_value()` (pylablib.core.gui.value\_handling.LabelIndicatorHandler (pylablib.core.gui.widgets.param\_table.ParamTable method), 277  
`get_value()` (pylablib.core.gui.value\_handling.LabelValueHandler (pylablib.core.gui.widgets.param\_table.StatusTable method), 271  
`get_value()` (pylablib.core.gui.value\_handling.LineEditValueHandler (pylablib.core.thread.callsync.QCallResultSynchronizer method), 270  
`get_value()` (pylablib.core.gui.value\_handling.ProgressBarValueHandler (pylablib.core.thread.callsync.QDirectResultSynchronizer method), 276  
`get_value()` (pylablib.core.gui.value\_handling.PropertyValueHandler (pylablib.core.thread.synchronizing.QThreadNotifier method), 267  
`get_value()` (pylablib.core.gui.value\_handling.PushButtonValueHandler (pylablib.devices.Andor.AndorSDK3.AndorSDK3Attribute method), 273  
`get_value()` (pylablib.core.gui.value\_handling.StandardIndicatorHandler (pylablib.devices.Attocube.anc350.ANC350 method), 277  
`get_value()` (pylablib.core.gui.value\_handling.StandardValueHandler (pylablib.devices.Basler.pylon.BaslerPylonAttribute method), 268  
`get_value()` (pylablib.core.gui.value\_handling.ToolButtonValueHandler (pylablib.devices.DCAM.DCAM.DCAMAttribute method), 274  
`get_value()` (pylablib.core.gui.value\_handling.VirtualValueHandler (pylablib.devices.IMAQdx.IMAQdx.IMAQdxAttribute method), 266  
`get_value()` (pylablib.devices.Pfeiffer.base.DPG202 method), 657  
`get_value()` (pylablib.devices.Photometrics.pvcam.PvcamAttribute method), 661  
`get_value()` (pylablib.devices.PhotonFocus.PhotonFocus.PFCamAttribute method), 671  
`get_value()` (pylablib.devices.PrincetonInstruments.picam.PicamAttribute method), 711  
`get_value()` (pylablib.devices.SiliconSoftware.fgrab.FGrabAttribute method), 720

<code>get_value_changed_signal()</code> <i>lablib.core.gui.value_handling.ISingleValueHandler</i> method), 269	(py- <code>get_value_changed_signal()</code> <i>lablib.core.gui.widgets.container.QTabContainer</i> method), 232	(py-
<code>get_value_changed_signal()</code> <i>lablib.core.gui.value_handling.IValueHandler</i> method), 266	(py- <code>get_value_changed_signal()</code> <i>lablib.core.gui.widgets.container.QWidgetContainer</i> method), 211	(py-
<code>get_value_changed_signal()</code> <i>lablib.core.gui.value_handling.LabelValueHandler</i> method), 271	(py- <code>get_value_changed_signal()</code> <i>lablib.core.gui.widgets.param_table.ParamTable</i> method), 250	(py-
<code>get_value_changed_signal()</code> <i>lablib.core.gui.value_handling.LineEditValueHandler</i> method), 270	(py- <code>get_value_changed_signal()</code> <i>lablib.core.gui.widgets.param_table.StatusTable</i> method), 259	(py-
<code>get_value_changed_signal()</code> <i>lablib.core.gui.value_handling.ProgressBarValueHandler</i> method), 276	(py- <code>get_value_sync()</code> <i>lablib.core.thread.callsync.QCallResultSynchronizer</i> method), 281	(py-
<code>get_value_changed_signal()</code> <i>lablib.core.gui.value_handling.PropertyValueHandler</i> method), 268	(py- <code>get_value_sync()</code> <i>lablib.core.thread.callsync.QDirectResultSynchronizer</i> method), 282	(py-
<code>get_value_changed_signal()</code> <i>lablib.core.gui.value_handling.PushButtonValueHandler</i> method), 273	(py- <code>get_value_sync()</code> <i>lablib.core.thread.synchronizing.QThreadNotifier</i> method), 318	(py-
<code>get_value_changed_signal()</code> <i>lablib.core.gui.value_handling.StandardValueHandler</i> method), 268	(py- <code>get_variable()</code> <i>lablib.core.thread.controller.QTaskThread</i> method), 309	(py-
<code>get_value_changed_signal()</code> <i>lablib.core.gui.value_handling.ToolButtonValueHandler</i> method), 274	(py- <code>get_variable()</code> <i>lablib.core.thread.controller.QThreadController</i> method), 298	(py-
<code>get_value_changed_signal()</code> <i>lablib.core.gui.value_handling.VirtualValueHandler</i> method), 267	(py- <code>get_variable()</code> <i>lablib.core.utils.ipc.SharedMemIPCTable</i> method), 385	(py-
<code>get_value_changed_signal()</code> <i>lablib.core.gui.widgets.container.IQContainer</i> method), 202	(py- <code>get_velocity()</code> <i>lablib.devices.PhysikInstrumente.base.PIE516</i> method), 705	(py-
<code>get_value_changed_signal()</code> <i>lablib.core.gui.widgets.container.IQWidgetContainer</i> method), 207	(py- <code>get_velocity()</code> <i>lablib.devices.Thorlabs.elliptec.ElliptecMotor</i> method), 779	(py-
<code>get_value_changed_signal()</code> <i>lablib.core.gui.widgets.container.QContainer</i> method), 204	(py- <code>get_velocity_factor()</code> <i>lablib.devices.Trinamic.base.TMCM1110</i> method), 828	(py-
<code>get_value_changed_signal()</code> <i>lablib.core.gui.widgets.container.QDialogContainer</i> method), 219	(py- <code>get_velocity_parameters()</code> <i>lablib.devices.Newport.picomotor.Picomotor8742</i> method), 630	(py-
<code>get_value_changed_signal()</code> <i>lablib.core.gui.widgets.container.QFrameContainer</i> method), 215	(py- <code>get_velocity_parameters()</code> <i>lablib.devices.Thorlabs.kinesis.KinesisMotor</i> method), 796	(py-
<code>get_value_changed_signal()</code> <i>lablib.core.gui.widgets.container.QGroupBoxContainer</i> method), 223	(py- <code>get_velocity_parameters()</code> <i>lablib.devices.Trinamic.base.TMCM1110</i> method), 828	(py-
<code>get_value_changed_signal()</code> <i>lablib.core.gui.widgets.container.QScrollAreaContainer</i> method), 230	(py- <code>get_vertical_position()</code> <i>lablib.devices.Tektronix.base.DPO2000</i> method), 764	(py-
<code>get_value_changed_signal()</code> <i>lablib.core.gui.widgets.container.QScrollAreaContainer.QComboBoxWidget</i> method), 227	(py- <code>get_vertical_position()</code> <i>lablib.devices.Tektronix.base.ITektronixScope</i> method), 750	(py-

<code>get_vertical_position()</code> <i>lablib.devices.Tektronix.base.TDS2000 method</i> ), 758	(py- <code>get_vsspeed()</code> <i>lablib.devices.Andor.AndorSDK2.AndorSDK2Camera method</i> ), 468	(py-
<code>get_vertical_span()</code> <i>lablib.devices.Tektronix.base.DPO2000 method</i> ), 764	(py- <code>get_vsspeed_period()</code> <i>lablib.devices.Andor.AndorSDK2.AndorSDK2Camera method</i> ), 468	(py-
<code>get_vertical_span()</code> <i>lablib.devices.Tektronix.base.ITektronixScope method</i> ), 750	(py- <code>get_warning_status()</code> <i>lablib.devices.LighthousePhotonics.base.SproutG method</i> ), 612	(py-
<code>get_vertical_span()</code> <i>lablib.devices.Tektronix.base.TDS2000 method</i> ), 758	(py- <code>get_waveform()</code> <i>lablib.devices.OZOptics.base.EPC04 method</i> ), 637	(py-
<code>get_voltage()</code> <i>lablib.devices.Attocube.anc300.ANC300 method</i> ), 503	(py- <code>get_wavelength()</code> <i>lablib.devices.HighFinesse.wlm.WLM method</i> ), 558	(py-
<code>get_voltage()</code> <i>lablib.devices.Attocube.anc350.ANC350 method</i> ), 508	(py- <code>get_wavelength()</code> <i>lablib.devices.Ophir.base.VegaPowerMeter method</i> ), 642	(py-
<code>get_voltage()</code> <i>lablib.devices.PhysikInstrumente.base.PIE516 method</i> ), 705	(py- <code>get_wavelength()</code> <i>lablib.devices.OZOptics.base.TF100 method</i> ), 633	(py-
<code>get_voltage()</code> <i>lablib.devices.Thorlabs.serial.MDT69xA method</i> ), 820	(py- <code>get_wavelength_correction()</code> <i>lablib.devices.OZOptics.base.TF100 method</i> ), 633	(py-
<code>get_voltage_input_parameters()</code> <i>lablib.devices.NI.daq.NIDAQ method</i> ), 623	(py- <code>get_wavelength_info()</code> <i>lablib.devices.Ophir.base.VegaPowerMeter method</i> ), 641	(py-
<code>get_voltage_lower_limit()</code> <i>lablib.devices.PhysikInstrumente.base.PIE516 method</i> ), 705	<code>get_wfmpre()</code> ( <i>pylablib.devices.Tektronix.base.DPO2000 method</i> ), 765	
<code>get_voltage_output_buffer_fill()</code> <i>lablib.devices.NI.daq.NIDAQ method</i> ), 625	<code>get_wfmpre()</code> ( <i>pylablib.devices.Tektronix.base.ITektronixScope method</i> ), 752	
<code>get_voltage_output_channels()</code> <i>lablib.devices.NI.daq.NIDAQ method</i> ), 625	<code>get_wfmpre()</code> ( <i>pylablib.devices.Tektronix.base.TDS2000 method</i> ), 758	
<code>get_voltage_output_clock_parameters()</code> <i>(pylablib.devices.NI.daq.NIDAQ method)</i> ), 627	<code>get_widget()</code> ( <i>pylablib.core.gui.value_handling.GUIValues method</i> ), 278	
<code>get_voltage_output_parameters()</code> <i>lablib.devices.NI.daq.NIDAQ method</i> ), 625	<code>get_widget()</code> ( <i>pylablib.core.gui.widgets.container.IQContainer method</i> ), 202	
<code>get_voltage_outputs()</code> <i>lablib.devices.NI.daq.NIDAQ method</i> ), 625	<code>get_widget()</code> ( <i>pylablib.core.gui.widgets.container.IQWidgetContainer method</i> ), 207	
<code>get_voltage_pattern()</code> <i>lablib.devices.Attocube.anc300.ANC300 method</i> ), 504	<code>get_widget()</code> ( <i>pylablib.core.gui.widgets.container.QContainer method</i> ), 204	
<code>get_voltage_range()</code> <i>lablib.devices.Thorlabs.serial.MDT69xA method</i> ), 820	<code>get_widget()</code> ( <i>pylablib.core.gui.widgets.container.QDialogContainer method</i> ), 219	
<code>get_voltage_setpoint()</code> <i>lablib.devices.PhysikInstrumente.base.PIE516 method</i> ), 705	<code>get_widget()</code> ( <i>pylablib.core.gui.widgets.container.QGroupBoxContainer method</i> ), 223	
<code>get_voltage_upper_limit()</code> <i>lablib.devices.PhysikInstrumente.base.PIE516 method</i> ), 705	<code>get_widget()</code> ( <i>pylablib.core.gui.widgets.container.QScrollAreaContainer method</i> ), 230	
<code>get_voltages()</code> <i>lablib.devices.OZOptics.base.EPC04 method</i> ), 636	<code>get_widget()</code> ( <i>pylablib.core.gui.widgets.container.QScrollAreaContainer method</i> ), 227	
	<code>get_widget()</code> ( <i>pylablib.core.gui.widgets.container.QTabContainer method</i> ), 233	
	<code>get_widget()</code> ( <i>pylablib.core.gui.widgets.container.QWidgetContainer</i>	



method), 211  
 get\_widget() (pylablib.core.gui.widgets.param\_table.ParamTablemethod), 582  
 method), 248  
 get\_widget() (pylablib.core.gui.widgets.param\_table.StatusTablemethod), 839  
 method), 259  
 get\_widget\_location() (in module py-  
 lablib.core.gui.utils), 263  
 get\_window\_ft\_func() (in module py-  
 lablib.core.dataproc.specfunc), 125  
 get\_window\_func() (in module py-  
 lablib.core.dataproc.specfunc), 125  
 get\_work\_hours() (py-  
 lablib.devices.LaserQuantum.base.Finesse  
 method), 604  
 get\_work\_hours() (py-  
 lablib.devices.LighthousePhotonics.base.SproutG  
 method), 612  
 get\_work\_hours() (py-  
 lablib.devices.Toptica.ibeam.TopicalBeam  
 method), 824  
 get\_x\_column() (in module py-  
 lablib.core.dataproc.utils), 135  
 get\_y\_column() (in module py-  
 lablib.core.dataproc.utils), 135  
 getargsfrom() (in module py-  
 lablib.core.utils.functions), 371  
 getattr\_call() (in module py-  
 lablib.core.utils.functions), 372  
 getattr\_multivar() (in module py-  
 lablib.core.utils.general), 374  
 getdefault() (in module py-  
 lablib.core.utils.funcargparse), 369  
 grab() (pylablib.devices.AlliedVision.Bonito.BonitoIMAQCamera  
 method), 460  
 grab() (pylablib.devices.AlliedVision.Bonito.IBonitoCamera  
 method), 454  
 grab() (pylablib.devices.Andor.AndorSDK2.AndorSDK2Camera  
 method), 474  
 grab() (pylablib.devices.Andor.AndorSDK3.AndorSDK3Camera  
 method), 484  
 grab() (pylablib.devices.Basler.pylon.BaslerPylonCamera  
 method), 518  
 grab() (pylablib.devices.BitFlow.BitFlow.BitFlowCamera  
 method), 529  
 grab() (pylablib.devices.BitFlow.BitFlow.BitFlowFrameGrabber  
 method), 525  
 grab() (pylablib.devices.DCAM.DCAM.DCAMCamera  
 method), 554  
 grab() (pylablib.devices.IMAQ.IMAQ.IMAQCamera  
 method), 572  
 grab() (pylablib.devices.IMAQ.IMAQ.IMAQFrameGrabber  
 method), 566  
 grab() (pylablib.devices.IMAQdx.IMAQdx.EthernetIMAQdxCamera  
 method), 586  
 grab() (pylablib.devices.IMAQdx.IMAQdx.IMAQdxCamera  
 method), 582  
 grab() (pylablib.devices.interface.camera.IAttributeCamera  
 method), 839  
 grab() (pylablib.devices.interface.camera.IBinROICamera  
 method), 858  
 grab() (pylablib.devices.interface.camera.ICamera  
 method), 834  
 grab() (pylablib.devices.interface.camera.IExposureCamera  
 method), 849  
 grab() (pylablib.devices.interface.camera.IGrabberAttributeCamera  
 method), 844  
 grab() (pylablib.devices.interface.camera.IROICamera  
 method), 854  
 grab() (pylablib.devices.PCO.SC2.PCOSC2Camera  
 method), 650  
 grab() (pylablib.devices.Photometrics.pvcam.PvcamCamera  
 method), 667  
 grab() (pylablib.devices.PhotonFocus.PhotonFocus.IPhotonFocusCamera  
 method), 675  
 grab() (pylablib.devices.PhotonFocus.PhotonFocus.PhotonFocusBitFlow  
 method), 698  
 grab() (pylablib.devices.PhotonFocus.PhotonFocus.PhotonFocusIMAQC  
 method), 682  
 grab() (pylablib.devices.PhotonFocus.PhotonFocus.PhotonFocusSiSoCa  
 method), 691  
 grab() (pylablib.devices.PrincetonInstruments.picam.PicamCamera  
 method), 716  
 grab() (pylablib.devices.SiliconSoftware.fgrab.SiliconSoftwareCamera  
 method), 730  
 grab() (pylablib.devices.SiliconSoftware.fgrab.SiliconSoftwareFrameGra  
 method), 725  
 grab() (pylablib.devices.Thorlabs.TLCamera.ThorlabsTLCamera  
 method), 774  
 grab() (pylablib.devices.uc480.uc480.UC480Camera  
 method), 871  
 grab\_continuous() (py-  
 lablib.devices.Tektronix.base.DPO2000  
 method), 765  
 grab\_continuous() (py-  
 lablib.devices.Tektronix.base.ITektronixScope  
 method), 749  
 grab\_continuous() (py-  
 lablib.devices.Tektronix.base.TDS2000  
 method), 758  
 grab\_single() (py-  
 lablib.devices.Tektronix.base.DPO2000  
 method), 765  
 grab\_single() (py-  
 lablib.devices.Tektronix.base.ITektronixScope  
 method), 749  
 grab\_single() (py-  
 lablib.devices.Tektronix.base.TDS2000  
 method), 758

grabber\_info (pylablib.devices.AlliedVision.Bonito.TDeviceInfo attribute), 450

grabber\_info (pylablib.devices.PhotonFocus.PhotonFocus.TDeviceInfo attribute), 671

GrabberClass (pylablib.devices.AlliedVision.Bonito.BonitoMAQCamera attribute), 456

GrabberClass (pylablib.devices.AlliedVision.Bonito.IBonitoCamera attribute), 451

GrabberClass (pylablib.devices.PhotonFocus.PhotonFocus.IPhotonFocusCamera attribute), 672

GrabberClass (pylablib.devices.PhotonFocus.PhotonFocus.PhotonFocusPylibCameras attribute), 694

GrabberClass (pylablib.devices.PhotonFocus.PhotonFocus.PhotonFocusPylibCameras attribute), 678

GrabberClass (pylablib.devices.PhotonFocus.PhotonFocus.PhotonFocusPylibCameras attribute), 687

gui\_thread\_method() (in module pylablib.core.thread.controller), 293

gui\_values\_path (pylablib.core.gui.widgets.container.TChild attribute), 200

GUIValues (class in pylablib.core.gui.value\_handling), 278

GUIValues.IndicatorsSet (class in pylablib.core.gui.value\_handling), 279

## H

hamming\_w() (in module pylablib.core.dataproc.specfunc), 125

hamming\_w\_ft() (in module pylablib.core.dataproc.specfunc), 125

hann\_w() (in module pylablib.core.dataproc.specfunc), 125

hann\_w\_ft() (in module pylablib.core.dataproc.specfunc), 125

has\_arg() (pylablib.core.dataproc.callable.FunctionCallable method), 106

has\_arg() (pylablib.core.dataproc.callable.ICallable method), 104

has\_arg() (pylablib.core.dataproc.callable.JoinedCallable method), 105

has\_arg() (pylablib.core.dataproc.callable.MethodCallable method), 107

has\_arg() (pylablib.core.dataproc.callable.MultiplexedCallable method), 105

has\_calls() (pylablib.core.thread.callsync.QQueueLengthEnforcer method), 288

has\_calls() (pylablib.core.thread.callsync.QQueueScheduler method), 286

has\_calls() (pylablib.core.thread.callsync.QQueueSizeLimitScheduler method), 290

has\_entry() (pylablib.core.utils.dictionary.Dictionary method), 327

high\_pass\_filter() (in module pylablib.core.dataproc.filters), 111

high\_thresh (pylablib.devices.Pfeiffer.base.TTPG260SwitchSettings attribute), 654

high\_value (pylablib.devices.Lakeshore.base.TLakeshore218AnalogSettings attribute), 592

high\_value (pylablib.devices.Lakeshore.base.TLakeshore370AnalogSettings attribute), 598

home() (pylablib.devices.Arcus.performax.Performax2EXStage method), 495

home() (pylablib.devices.Arcus.performax.Performax4EXStage method), 491

home() (pylablib.devices.Arcus.performax.PerformaxDMXJSASStage method), 498

home() (pylablib.devices.OZOptics.base.DD100 method), 635

home() (pylablib.devices.OZOptics.base.TF100 method), 633

home() (pylablib.devices.Thorlabs.elliptec.ElliptecMotor method), 779

home() (pylablib.devices.Thorlabs.kinesis.KinesisMotor method), 795

home\_direction (pylablib.devices.Thorlabs.kinesis.THomeParams attribute), 785

hour (pylablib.devices.uc480.uc480.TTimestamp attribute), 865

hw\_swapped (pylablib.devices.PrincetonInstruments.picam.TROIConstraints attribute), 709

hw\_swapped (pylablib.devices.Thorlabs.kinesis.TLimitSwitchParams attribute), 785

hw\_swapped (pylablib.devices.Thorlabs.kinesis.TLimitSwitchParams attribute), 785

hw\_swapped (pylablib.devices.Thorlabs.kinesis.TLimitSwitchParams attribute), 785

hw\_type (pylablib.devices.Thorlabs.kinesis.TDeviceInfo attribute), 785

- attribute), 781
- hw\_ver (pylablib.devices.Thorlabs.elliptec.TDeviceInfo attribute), 777
- hw\_ver (pylablib.devices.Thorlabs.kinesis.TDeviceInfo attribute), 781
- ## I
- I (pylablib.devices.Sirah.Matisse.TFastpiezoCtlParameters attribute), 735
- I (pylablib.devices.Sirah.Matisse.TSlowpiezoCtlParameters attribute), 735
- I (pylablib.devices.Sirah.Matisse.TThinCtlParameters attribute), 734
- i (pylablib.devices.Thorlabs.kinesis.TQuadDetectorPIDParams attribute), 805
- i () (pylablib.core.dataproc.transform.Indexed2DTransform method), 135
- i () (pylablib.core.dataproc.transform.LinearTransform method), 134
- i () (pylablib.core.devio.interface.CombinedParameterClass method), 167
- i () (pylablib.core.devio.interface.EnumParameterClass method), 166
- i () (pylablib.core.devio.interface.FunctionParameterClass method), 167
- i () (pylablib.core.devio.interface.ICheckingParameterClass method), 164
- i () (pylablib.core.devio.interface.IEnumParameterClass method), 165
- i () (pylablib.core.devio.interface.IParameterClass method), 163
- i () (pylablib.core.devio.interface.RangeParameterClass method), 164
- I1DWrapper (class in pylablib.core.dataproc.table\_wrap), 126
- I1DWrapper.Accessor (class in pylablib.core.dataproc.table\_wrap), 126
- I2DWrapper (class in pylablib.core.dataproc.table\_wrap), 129
- IArduinoDevice (class in pylablib.devices.Arduino.base), 500
- IArduinoDevice.NoParameterCaller (class in pylablib.devices.Arduino.base), 501
- IAttributeCamera (class in pylablib.devices.interface.camera), 837
- IAttributeCamera.NoParameterCaller (class in pylablib.devices.interface.camera), 837
- IBinaryOutputFileFormat (class in pylablib.core.fileio.savefile), 193
- IBinROICamera (class in pylablib.devices.interface.camera), 856
- IBinROICamera.NoParameterCaller (class in pylablib.devices.interface.camera), 856
- IBonitoCamera (class in pylablib.devices.AlliedVision.Bonito), 450
- IBonitoCamera.NoParameterCaller (class in pylablib.devices.AlliedVision.Bonito), 452
- IBoolValueHandler (class in pylablib.core.gui.value\_handling), 271
- ICallable (class in pylablib.core.dataproc.callable), 104
- ICallable.NamesBoundCall (class in pylablib.core.dataproc.callable), 104
- ICamera (class in pylablib.devices.interface.camera), 831
- ICamera.NoParameterCaller (class in pylablib.devices.interface.camera), 834
- ICheckingParameterClass (class in pylablib.core.devio.interface), 163
- ICommBackendWrapper (class in pylablib.core.devio.comm\_backend), 160
- ICommBackendWrapper.NoParameterCaller (class in pylablib.core.devio.comm\_backend), 160
- id (pylablib.core.utils.ipc.TPipeMsg attribute), 383
- id (pylablib.devices.Newport.picomotor.TDeviceInfo attribute), 628
- id (pylablib.devices.Ophir.base.TDeviceInfo attribute), 640
- IDataLocation (class in pylablib.core.fileio.location), 184
- IDevice (class in pylablib.core.devio.interface), 162
- IDevice.NoParameterCaller (class in pylablib.core.devio.interface), 162
- IDeviceCommBackend (class in pylablib.core.devio.comm\_backend), 142
- IDictionaryEntry (class in pylablib.core.fileio.dict\_entry), 170
- idreg (pylablib.devices.BitFlow.BitFlow.TDeviceInfo attribute), 521
- idx (pylablib.devices.BitFlow.BitFlow.TDeviceInfo attribute), 521
- IEnumParameterClass (class in pylablib.core.devio.interface), 164
- IExposureCamera (class in pylablib.devices.interface.camera), 847
- IExposureCamera.NoParameterCaller (class in pylablib.devices.interface.camera), 847
- IExternalFileDictionaryEntry (class in pylablib.core.fileio.dict\_entry), 174
- IExternalTableDictionaryEntry (class in pylablib.core.fileio.dict\_entry), 172
- IFileSystemDataLocation (class in pylablib.core.fileio.location), 185
- IGenWrapper (class in pylablib.core.dataproc.table\_wrap), 125
- IGrabberAttributeCamera (class in py-

- [lablib.devices.interface.camera](#)), 841
- [IGrabberAttributeCamera.NoParameterCaller](#) (class in [pylablib.devices.interface.camera](#)), 842
- [IIndex](#) (class in [pylablib.core.utils.indexing](#)), 381
- [IIndicatorHandler](#) (class in [pylablib.core.gui.value\\_handling](#)), 276
- [IInputFileFormat](#) (class in [pylablib.core.fileio.loadfile](#)), 176
- [IIPCChannel](#) (class in [pylablib.core.utils.ipc](#)), 383
- [iir\\_apply\\_complex\(\)](#) (in module [pylablib.core.dataproc.iir\\_transform](#)), 121
- [ilabels](#) ([pylablib.devices.Andor.AndorSDK3.AndorSDK3Attribute](#) attribute), 478
- [ilabels](#) ([pylablib.devices.Basler.pylon.BaslerPylonAttribute](#) attribute), 513
- [ilabels](#) ([pylablib.devices.DCAM.DCAM.DCAMAttribute](#) attribute), 549
- [ilabels](#) ([pylablib.devices.IMAQdx.IMAQdx.IMAQdxAttribute](#) attribute), 577
- [ilabels](#) ([pylablib.devices.Photometrics.pvcam.PvcamAttribute](#) attribute), 660
- [ilabels](#) ([pylablib.devices.PhotonFocus.PhotonFocus.PFCAAttribute](#) attribute), 671
- [ilabels](#) ([pylablib.devices.PrincetonInstruments.picam.PicamAttribute](#) attribute), 711
- [ilabels](#) ([pylablib.devices.SiliconSoftware.fgrab.FGrabAttribute](#) attribute), 720
- [IMAQCamera](#) (class in [pylablib.devices.IMAQ.IMAQ](#)), 568
- [IMAQCamera.NoParameterCaller](#) (class in [pylablib.devices.IMAQ.IMAQ](#)), 568
- [IMAQdxAttribute](#) (class in [pylablib.devices.IMAQdx.IMAQdx](#)), 576
- [IMAQdxCamera](#) (class in [pylablib.devices.IMAQdx.IMAQdx](#)), 578
- [IMAQdxCamera.NoParameterCaller](#) (class in [pylablib.devices.IMAQdx.IMAQdx](#)), 580
- [IMAQFrameGrabber](#) (class in [pylablib.devices.IMAQ.IMAQ](#)), 561
- [IMAQFrameGrabber.NoParameterCaller](#) (class in [pylablib.devices.IMAQ.IMAQ](#)), 565
- [implemented](#) ([pylablib.devices.Andor.AndorSDK3.AndorSDK3Attribute](#) attribute), 477
- [implemented](#) ([pylablib.devices.Basler.pylon.BaslerPylonAttribute](#) attribute), 512
- [IMultiaxisStage](#) (class in [pylablib.devices.interface.stage](#)), 863
- [IMultiaxisStage.NoParameterCaller](#) (class in [pylablib.devices.interface.stage](#)), 863
- [in\\_use](#) ([pylablib.devices.uc480.uc480.TCameraInfo](#) attribute), 864
- [inc](#) ([pylablib.devices.Basler.pylon.BaslerPylonAttribute](#) attribute), 513
- [inc](#) ([pylablib.devices.IMAQdx.IMAQdx.IMAQdxAttribute](#) attribute), 577
- [inc](#) ([pylablib.devices.Photometrics.pvcam.PvcamAttribute](#) attribute), 660
- [inc](#) ([pylablib.devices.PrincetonInstruments.picam.PicamAttribute](#) attribute), 710
- [inc](#) ([pylablib.devices.SiliconSoftware.fgrab.FGrabAttribute](#) attribute), 720
- [inc\(\)](#) ([pylablib.devices.interface.camera.FrameNotifier](#) method), 836
- [ind](#) ([pylablib.core.gui.value\\_handling.GUIValues.IndicatorsSet](#) attribute), 279
- [index](#) ([pylablib.devices.Attocube.anc350.ANC350.Reply](#) attribute), 507
- [index](#) ([pylablib.devices.Attocube.anc350.ANC350.Telegram](#) attribute), 506
- [index\(\)](#) ([pylablib.core.dataproc.feature.Baseline](#) method), 108
- [index\(\)](#) ([pylablib.core.dataproc.feature.Peak](#) method), 109
- [index\(\)](#) ([pylablib.core.devio.interface.TRawParameterValue](#) method), 168
- [index\(\)](#) ([pylablib.core.gui.utils.TWidgetLocation](#) method), 263
- [index\(\)](#) ([pylablib.core.gui.value\\_handling.GUIValues.IndicatorsSet](#) method), 279
- [index\(\)](#) ([pylablib.core.gui.widgets.container.TChild](#) method), 200
- [index\(\)](#) ([pylablib.core.gui.widgets.container.TTimer](#) method), 199
- [index\(\)](#) ([pylablib.core.gui.widgets.container.TTimerEvent](#) method), 199
- [index\(\)](#) ([pylablib.core.gui.widgets.param\\_table.ParamTable.ParamRow](#) method), 243
- [index\(\)](#) ([pylablib.core.gui.widgets.param\\_table.StatusTable.ParamRow](#) method), 252
- [index\(\)](#) ([pylablib.core.thread.callsync.QScheduledCall.Callback](#) method), 283
- [index\(\)](#) ([pylablib.core.thread.callsync.TDefaultCallInfo](#) method), 284
- [index\(\)](#) ([pylablib.core.thread.controller.QTaskThread.TBatchJob](#) method), 302
- [index\(\)](#) ([pylablib.core.thread.controller.QTaskThread.TCommand](#) method), 302
- [index\(\)](#) ([pylablib.core.thread.multicast\\_pool.TMulticast](#) method), 315
- [index\(\)](#) ([pylablib.core.utils.dictionary.DictionaryDiff](#) method), 334
- [index\(\)](#) ([pylablib.core.utils.dictionary.DictionaryIntersection](#) method), 335
- [index\(\)](#) ([pylablib.core.utils.files.FolderList](#) method), 364
- [index\(\)](#) ([pylablib.core.utils.ipc.TPipeMsg](#) method), 383



- `index()` (`pylablib.core.utils.ipc.TShmemVarDesc` method), 384
- `index()` (`pylablib.core.utils.observer_pool.ObserverPool.Observer` method), 393
- `index()` (`pylablib.core.utils.string.TConversionClass` method), 398
- `index()` (`pylablib.devices.AlliedVision.Bonito.TDeviceInfo` method), 450
- `index()` (`pylablib.devices.AlliedVision.Bonito.TStatusLine` method), 464
- `index()` (`pylablib.devices.Andor.AndorSDK2.TAcqProgress` method), 466
- `index()` (`pylablib.devices.Andor.AndorSDK2.TCycleTimings` method), 466
- `index()` (`pylablib.devices.Andor.AndorSDK2.TDeviceInfo` method), 465
- `index()` (`pylablib.devices.Andor.AndorSDK3.TDeviceInfo` method), 478
- `index()` (`pylablib.devices.Andor.AndorSDK3.TFrameInfo` method), 479
- `index()` (`pylablib.devices.Andor.AndorSDK3.TMissedFramesStatus` method), 479
- `index()` (`pylablib.devices.Attocube.anc300.TDeviceInfo` method), 502
- `index()` (`pylablib.devices.Basler.pylon.TCameraInfo` method), 511
- `index()` (`pylablib.devices.Basler.pylon.TDeviceInfo` method), 513
- `index()` (`pylablib.devices.BitFlow.BitFlow.TDeviceInfo` method), 521
- `index()` (`pylablib.devices.Conrad.base.RelayBoard.TMessage` method), 534
- `index()` (`pylablib.devices.DCAM.DCAM.TDeviceInfo` method), 550
- `index()` (`pylablib.devices.DCAM.DCAM.TFrameInfo` method), 550
- `index()` (`pylablib.devices.HighFinesse.wlm.TDeviceInfo` method), 557
- `index()` (`pylablib.devices.IMAQ.IMAQ.TDeviceInfo` method), 561
- `index()` (`pylablib.devices.IMAQdx.IMAQdx.TCameraInfo` method), 576
- `index()` (`pylablib.devices.IMAQdx.IMAQdx.TDeviceInfo` method), 578
- `index()` (`pylablib.devices.interface.camera.TAcqTimings` method), 846
- `index()` (`pylablib.devices.interface.camera.TAxisROILimit` method), 851
- `index()` (`pylablib.devices.interface.camera.TFrameInfo` method), 831
- `index()` (`pylablib.devices.interface.camera.TFramePosition` method), 831
- `index()` (`pylablib.devices.interface.camera.TFrameSize` method), 830
- `index()` (`pylablib.devices.interface.camera.TFramesStatus` method), 830
- `index()` (`pylablib.devices.interface.camera.TStatusLineDescription` method), 861
- `index()` (`pylablib.devices.KJL.base.TKJL300DeviceInfo` method), 590
- `index()` (`pylablib.devices.Lakeshore.base.TLakeshore218AnalogSettings` method), 592
- `index()` (`pylablib.devices.Lakeshore.base.TLakeshore218CurveHeader` method), 593
- `index()` (`pylablib.devices.Lakeshore.base.TLakeshore218FilterSettings` method), 593
- `index()` (`pylablib.devices.Lakeshore.base.TLakeshore370AnalogSettings` method), 598
- `index()` (`pylablib.devices.Lakeshore.base.TLakeshore370FilterSettings` method), 599
- `index()` (`pylablib.devices.Lakeshore.base.TLakeshore370RangeSettings` method), 598
- `index()` (`pylablib.devices.LaserQuantum.base.TDeviceInfo` method), 603
- `index()` (`pylablib.devices.LaserQuantum.base.TTemperatures` method), 604
- `index()` (`pylablib.devices.LaserQuantum.base.TWorkHours` method), 604
- `index()` (`pylablib.devices.Leybold.base.TDeviceInfo` method), 607
- `index()` (`pylablib.devices.Leybold.base.TITR90Status` method), 609
- `index()` (`pylablib.devices.Leybold.base.TUpdateValue` method), 607
- `index()` (`pylablib.devices.LighthousePhotonics.base.TDeviceInfo` method), 611
- `index()` (`pylablib.devices.LighthousePhotonics.base.TWorkHours` method), 611
- `index()` (`pylablib.devices.Newport.picomotor.TDeviceInfo` method), 628
- `index()` (`pylablib.devices.NI.daq.TDeviceInfo` method), 620
- `index()` (`pylablib.devices.NI.daq.TVoltageOutputClockParameters` method), 621
- `index()` (`pylablib.devices.Ophir.base.TDeviceInfo` method), 640
- `index()` (`pylablib.devices.Ophir.base.THeadInfo` method), 640
- `index()` (`pylablib.devices.Ophir.base.TRangeInfo` method), 641
- `index()` (`pylablib.devices.Ophir.base.TWavelengthInfo` method), 640
- `index()` (`pylablib.devices.PCO.SC2.TCameraStatus` method), 644
- `index()` (`pylablib.devices.PCO.SC2.TDeviceInfo` method), 644
- `index()` (`pylablib.devices.PCO.SC2.TFrameInfo` method), 645

[index\(\) \(pylablib.devices.PCO.SC2.TInternalBufferStatus method\), 645](#)  
[index\(\) \(pylablib.devices.PCO.SC2.TStatusLine method\), 653](#)  
[index\(\) \(pylablib.devices.Pfeiffer.base.TTPG260GaugeCameraSettings method\), 654](#)  
[index\(\) \(pylablib.devices.Pfeiffer.base.TTPG260SwitchSettings method\), 654](#)  
[index\(\) \(pylablib.devices.Photometrics.pvcam.TDeviceInfo method\), 661](#)  
[index\(\) \(pylablib.devices.Photometrics.pvcam.TFrameInfo method\), 661](#)  
[index\(\) \(pylablib.devices.Photometrics.pvcam.TReadoutInfo method\), 662](#)  
[index\(\) \(pylablib.devices.PhotonFocus.PhotonFocus.TCameraInfo method\), 670](#)  
[index\(\) \(pylablib.devices.PhotonFocus.PhotonFocus.TDeviceInfo method\), 671](#)  
[index\(\) \(pylablib.devices.PrincetonInstruments.picam.TCameraInfo method\), 708](#)  
[index\(\) \(pylablib.devices.PrincetonInstruments.picam.TDeviceInfo method\), 711](#)  
[index\(\) \(pylablib.devices.PrincetonInstruments.picam.TFrameInfo method\), 712](#)  
[index\(\) \(pylablib.devices.PrincetonInstruments.picam.TROIConstraint method\), 709](#)  
[index\(\) \(pylablib.devices.SiliconSoftware.fgrab.TAppletInfo method\), 718](#)  
[index\(\) \(pylablib.devices.SiliconSoftware.fgrab.TBoardInfo method\), 718](#)  
[index\(\) \(pylablib.devices.SiliconSoftware.fgrab.TDeviceInfo method\), 720](#)  
[index\(\) \(pylablib.devices.SiliconSoftware.fgrab.TFrameInfo method\), 721](#)  
[index\(\) \(pylablib.devices.SiliconSoftware.fgrab.TFullAppletInfo method\), 719](#)  
[index\(\) \(pylablib.devices.Sirah.Matisse.TFastpiezoCtlParameters method\), 735](#)  
[index\(\) \(pylablib.devices.Sirah.Matisse.TPiezoetDriveParameters method\), 734](#)  
[index\(\) \(pylablib.devices.Sirah.Matisse.TPiezoetFeedbackParameters method\), 734](#)  
[index\(\) \(pylablib.devices.Sirah.Matisse.TPiezoetFeedforwardParameters method\), 734](#)  
[index\(\) \(pylablib.devices.Sirah.Matisse.TRefcellWaveformParameters method\), 735](#)  
[index\(\) \(pylablib.devices.Sirah.Matisse.TScanMode method\), 736](#)  
[index\(\) \(pylablib.devices.Sirah.Matisse.TScanParameters method\), 736](#)  
[index\(\) \(pylablib.devices.Sirah.Matisse.TSlowpiezoCtlParameters method\), 735](#)  
[index\(\) \(pylablib.devices.Sirah.Matisse.TThinCtlParameters method\), 734](#)  
[index\(\) \(pylablib.devices.SmarAct.scu3d.TDeviceInfo method\), 745](#)  
[index\(\) \(pylablib.devices.Tektronix.base.TTriggerParameters method\), 748](#)  
[index\(\) \(pylablib.devices.Thorlabs.elliptec.ElliptecMotor.CommData method\), 778](#)  
[index\(\) \(pylablib.devices.Thorlabs.elliptec.TDeviceInfo method\), 777](#)  
[index\(\) \(pylablib.devices.Thorlabs.elliptec.TMotorInfo method\), 777](#)  
[index\(\) \(pylablib.devices.Thorlabs.kinesis.BasicKinesisDevice.CommData method\), 782](#)  
[index\(\) \(pylablib.devices.Thorlabs.kinesis.BasicKinesisDevice.CommShort method\), 782](#)  
[index\(\) \(pylablib.devices.Thorlabs.kinesis.KinesisDevice.CommData method\), 788](#)  
[index\(\) \(pylablib.devices.Thorlabs.kinesis.KinesisDevice.CommShort method\), 788](#)  
[index\(\) \(pylablib.devices.Thorlabs.kinesis.KinesisMotor.CommData method\), 798](#)  
[index\(\) \(pylablib.devices.Thorlabs.kinesis.KinesisMotor.CommShort method\), 798](#)  
[index\(\) \(pylablib.devices.Thorlabs.kinesis.KinesisPiezoMotor.CommData method\), 802](#)  
[index\(\) \(pylablib.devices.Thorlabs.kinesis.KinesisPiezoMotor.CommShort method\), 802](#)  
[index\(\) \(pylablib.devices.Thorlabs.kinesis.KinesisQuadDetector.CommData method\), 807](#)  
[index\(\) \(pylablib.devices.Thorlabs.kinesis.KinesisQuadDetector.CommShort method\), 807](#)  
[index\(\) \(pylablib.devices.Thorlabs.kinesis.MFF.CommData method\), 792](#)  
[index\(\) \(pylablib.devices.Thorlabs.kinesis.MFF.CommShort method\), 792](#)  
[index\(\) \(pylablib.devices.Thorlabs.kinesis.TDeviceInfo method\), 781](#)  
[index\(\) \(pylablib.devices.Thorlabs.kinesis.TFlipperParameters method\), 790](#)  
[index\(\) \(pylablib.devices.Thorlabs.kinesis.TGenMoveParams method\), 785](#)  
[index\(\) \(pylablib.devices.Thorlabs.kinesis.THomeParams method\), 785](#)  
[index\(\) \(pylablib.devices.Thorlabs.kinesis.TJogParams method\), 784](#)  
[index\(\) \(pylablib.devices.Thorlabs.kinesis.TKCubeTrigIOParams method\), 786](#)  
[index\(\) \(pylablib.devices.Thorlabs.kinesis.TKCubeTrigPosParams method\), 786](#)  
[index\(\) \(pylablib.devices.Thorlabs.kinesis.TLimitSwitchParams method\), 785](#)  
[index\(\) \(pylablib.devices.Thorlabs.kinesis.TPZMotorDriveParams method\), 786](#)  
[index\(\) \(pylablib.devices.Thorlabs.kinesis.TPZMotorJogParams method\), 787](#)

[index\(\)](#) (`pylablib.devices.Thorlabs.kinesis.TQuadDetectorOutputParameter` attribute), 392  
[index\(\)](#) (`pylablib.devices.Thorlabs.kinesis.TQuadDetectorPIDParameter` attribute), 392  
[index\(\)](#) (`pylablib.devices.Thorlabs.kinesis.TQuadDetectorReadingAttribute` attribute), 719  
[index\(\)](#) (`pylablib.devices.Thorlabs.kinesis.TQuadDetectorSetpoint` attribute), 805  
[index\(\)](#) (`pylablib.devices.Thorlabs.kinesis.TVelocityParameter` attribute), 784  
[index\(\)](#) (`pylablib.devices.Thorlabs.TLCamera.TDeviceInfo` attribute), 769  
[index\(\)](#) (`pylablib.devices.Thorlabs.TLCamera.TFrameInfo` attribute), 770  
[index\(\)](#) (`pylablib.devices.Toptica.ibeam.TDeviceInfo` attribute), 823  
[index\(\)](#) (`pylablib.devices.Toptica.ibeam.TTemperatures` attribute), 824  
[index\(\)](#) (`pylablib.devices.Toptica.ibeam.TWorkHours` attribute), 823  
[index\(\)](#) (`pylablib.devices.Trinamic.base.TMCM1110.ReplyData` attribute), 827  
[index\(\)](#) (`pylablib.devices.uc480.uc480.TAcquiredFramesStatus` attribute), 865  
[index\(\)](#) (`pylablib.devices.uc480.uc480.TCameraInfo` attribute), 864  
[index\(\)](#) (`pylablib.devices.uc480.uc480.TDeviceInfo` attribute), 865  
[index\(\)](#) (`pylablib.devices.uc480.uc480.TFrameInfo` attribute), 866  
[index\(\)](#) (`pylablib.devices.uc480.uc480.TTimestamp` attribute), 865  
[index\(\)](#) (`pylablib.devices.utils.load_lib.TLibraryCloseResult` attribute), 874  
[index\(\)](#) (`pylablib.devices.utils.load_lib.TLibraryOpenResult` attribute), 874  
[index\\_to\\_value\(\)](#) (`pylablib.core.gui.widgets.combo_box.ComboBox` attribute), 198  
[Indexed2DTransform](#) (class in `pylablib.core.dataproc.transform`), 134  
[indicator](#) (`pylablib.core.gui.widgets.param_table.ParamTable.ParamRow` attribute), 243  
[indicator](#) (`pylablib.core.gui.widgets.param_table.StatusTable.ParamRow` attribute), 252  
[indicator\\_handler](#) (`pylablib.core.gui.widgets.param_table.ParamTable.ParamRow` attribute), 243  
[indicator\\_handler](#) (`pylablib.core.gui.widgets.param_table.StatusTable.ParamRow` attribute), 252  
[individual](#) (`pylablib.core.utils.dictionary.DictionaryIntersection` attribute), 335  
[infinite\\_list](#) (class in `pylablib.core.utils.numerical`), 392  
[infinite\\_list.counter](#) (class in `pylablib.core.utils.numerical`), 392  
[info](#) (`pylablib.devices.SiliconSoftware.fgrab.TFullAppletInfo` attribute), 719  
[init\\_amp\\_mode\(\)](#) (`pylablib.devices.Andor.AndorSDK2.AndorSDK2Camera` attribute), 469  
[init\\_result](#) (`pylablib.devices.utils.load_lib.TLibraryOpenResult` attribute), 874  
[initial\\_guess\(\)](#) (`pylablib.core.dataproc.fitting.Fitter` attribute), 117  
[InlineTable](#) (class in `pylablib.core.fileio.loadfile_utils`), 182  
[InlineTableDictionaryEntry](#) (class in `pylablib.core.fileio.dict_entry`), 171  
[insert\(\)](#) (`pylablib.core.dataproc.table_wrap.Array1DWrapper` attribute), 127  
[insert\(\)](#) (`pylablib.core.dataproc.table_wrap.Array2DWrapper.ColumnA` attribute), 130  
[insert\(\)](#) (`pylablib.core.dataproc.table_wrap.Array2DWrapper.RowA` attribute), 130  
[insert\(\)](#) (`pylablib.core.dataproc.table_wrap.DataFrame2DWrapper.Col` attribute), 132  
[insert\(\)](#) (`pylablib.core.dataproc.table_wrap.DataFrame2DWrapper.Row` attribute), 132  
[insert\\_column\(\)](#) (`pylablib.core.gui.widgets.container.IQWidgetContainer` attribute), 207  
[insert\\_column\(\)](#) (`pylablib.core.gui.widgets.container.QDialogContainer` attribute), 219  
[insert\\_column\(\)](#) (`pylablib.core.gui.widgets.container.QFrameContainer` attribute), 215  
[insert\\_column\(\)](#) (`pylablib.core.gui.widgets.container.QGroupBoxContainer` attribute), 223  
[insert\\_column\(\)](#) (`pylablib.core.gui.widgets.container.QScrollAreaContainer.QContai` attribute), 227  
[insert\\_column\(\)](#) (`pylablib.core.gui.widgets.container.QWidgetContainer` attribute), 211  
[insert\\_column\(\)](#) (`pylablib.core.gui.widgets.layout_manager.IQLayoutManagedWidget` attribute), 240  
[insert\\_column\(\)](#) (`pylablib.core.gui.widgets.layout_manager.QLayoutManagedWidget` attribute), 241  
[insert\\_column\(\)](#) (`pylablib.core.gui.widgets.param_table.ParamTable` attribute), 250

`insert_column()` (pylablib.core.gui.widgets.param\_table.StatusTable method), 259  
`insert_layout_column()` (in module pylablib.core.gui.utils), 264  
`insert_layout_row()` (in module pylablib.core.gui.utils), 264  
`insert_option()` (pylablib.core.gui.widgets.combo\_box.ComboBox method), 199  
`insert_row()` (pylablib.core.gui.widgets.container.IQWidgetContainer method), 207  
`insert_row()` (pylablib.core.gui.widgets.container.QDialogContainer method), 219  
`insert_row()` (pylablib.core.gui.widgets.container.QFrameContainer method), 215  
`insert_row()` (pylablib.core.gui.widgets.container.QGroupBoxContainer method), 223  
`insert_row()` (pylablib.core.gui.widgets.container.QScrollAreaContainer method), 227  
`insert_row()` (pylablib.core.gui.widgets.container.QWidgetContainer method), 211  
`insert_row()` (pylablib.core.gui.widgets.layout\_manager.IQLayoutManagerWidget method), 240  
`insert_row()` (pylablib.core.gui.widgets.layout\_manager.QLayoutManagerWidget method), 241  
`insert_row()` (pylablib.core.gui.widgets.param\_table.ParamTable method), 250  
`insert_row()` (pylablib.core.gui.widgets.param\_table.StatusTable method), 259  
`insert_status_line()` (in module pylablib.devices.interface.camera), 861  
`install_if_older()` (in module pylablib.core.utils.module), 386  
`InstekAFG2000` (class in pylablib.devices.AWG.specific), 426  
`InstekAFG2000.NoParameterCaller` (class in pylablib.devices.AWG.specific), 426  
`InstekAFG2225` (class in pylablib.devices.AWG.specific), 420  
`InstekAFG2225.NoParameterCaller` (class in pylablib.devices.AWG.specific), 421  
`int2bits()` (in module pylablib.core.utils.strpack), 401  
`int2bytes()` (in module pylablib.core.utils.strpack), 401  
`integer_distance()` (in module pylablib.core.utils.numerical), 391  
`IntegerFormatter` (class in pylablib.core.gui.formatter), 261  
`integrate()` (in module pylablib.core.dataproc.filters), 111  
`interface` (pylablib.devices.IMAQ.IMAQ.TDeviceInfo attribute), 561  
`interface` (pylablib.devices.PCO.SC2.TDeviceInfo attribute), 644  
`interface` (pylablib.devices.PrincetonInstruments.picam.TCameraInfo attribute), 709  
`interface` (pylablib.devices.PrincetonInstruments.picam.TDeviceInfo attribute), 711  
`interpolate1D()` (in module pylablib.core.dataproc.interpolate), 122  
`interpolate1D_func()` (in module pylablib.core.dataproc.interpolate), 122  
`interpolate2D()` (in module pylablib.core.dataproc.interpolate), 122  
`interpolate_trace()` (in module pylablib.core.dataproc.interpolate), 123  
`interpolateND()` (in module pylablib.core.dataproc.interpolate), 123  
`InterruptException`, 320  
`InterruptExceptionStop`, 320  
`IOBase` (pylablib.core.dataproc.image.ROI class method), 121  
`IOBaseContainer` (pylablib.core.dataproc.utils.Range method), 136  
`IOBaseWidget` (pylablib.core.dataproc.fourier), 119  
`IOLayoutManagerWidget` (in module pylablib.core.utils.general), 375  
`Indexed2DTransform` (pylablib.core.dataproc.transform.Indexed2DTransform method), 135  
`LinearTransform` (pylablib.core.dataproc.transform.LinearTransform method), 134  
`io1_oper_mode` (pylablib.devices.Thorlabs.kinesis.TFlipperParameters attribute), 790  
`io1_pulse_width` (pylablib.devices.Thorlabs.kinesis.TFlipperParameters attribute), 790  
`io1_sig_mode` (pylablib.devices.Thorlabs.kinesis.TFlipperParameters attribute), 790  
`io2_oper_mode` (pylablib.devices.Thorlabs.kinesis.TFlipperParameters attribute), 790  
`io2_pulse_width` (pylablib.devices.Thorlabs.kinesis.TFlipperParameters attribute), 790  
`io2_sig_mode` (pylablib.devices.Thorlabs.kinesis.TFlipperParameters attribute), 790  
`io_status` (pylablib.devices.uc480.uc480.TFrameInfo attribute), 866  
`IOObjectCall` (class in pylablib.core.utils.functions), 372  
`IOObjectProperty` (class in pylablib.core.utils.functions), 372  
`IOOutputFileFormat` (class in pylablib.core.fileio.savefile), 190



<code>IParameterClass</code>	(class in <code>py-lablib.core.devio.interface</code> ), 163	<code>lablib.devices.IMAQdx.IMAQdx.IMAQdxCamera</code>	method), 582
<code>IPhotonFocusCamera</code>	(class in <code>py-lablib.devices.PhotonFocus.PhotonFocus</code> ), 672	<code>is_acquisition_setup()</code>	( <code>py-lablib.devices.interface.camera.IAttributeCamera</code> method), 839
<code>IPhotonFocusCamera.NoParameterCaller</code>	(class in <code>py-lablib.devices.PhotonFocus.PhotonFocus</code> ), 674	<code>is_acquisition_setup()</code>	( <code>py-lablib.devices.interface.camera.IBinROICamera</code> method), 858
<code>IQContainer</code>	(class in <code>py-lablib.core.gui.widgets.container</code> ), 200	<code>is_acquisition_setup()</code>	( <code>py-lablib.devices.interface.camera.ICamera</code> method), 831
<code>IQLayoutManagedWidget</code>	(class in <code>py-lablib.core.gui.widgets.layout_manager</code> ), 238	<code>is_acquisition_setup()</code>	( <code>py-lablib.devices.interface.camera.IExposureCamera</code> method), 849
<code>IQWidgetContainer</code>	(class in <code>py-lablib.core.gui.widgets.container</code> ), 205	<code>is_acquisition_setup()</code>	( <code>py-lablib.devices.interface.camera.IGrabberAttributeCamera</code> method), 844
<code>IROICamera</code>	(class in <code>py-lablib.devices.interface.camera</code> ), 851	<code>is_acquisition_setup()</code>	( <code>py-lablib.devices.interface.camera.IROICamera</code> method), 854
<code>IROICamera.NoParameterCaller</code>	(class in <code>py-lablib.devices.interface.camera</code> ), 852	<code>is_acquisition_setup()</code>	( <code>py-lablib.devices.interface.camera.IROICamera</code> method), 854
<code>is_acquisition_setup()</code>	( <code>py-lablib.devices.AlliedVision.Bonito.BonitoIMAQCamera</code> method), 460	<code>is_acquisition_setup()</code>	( <code>py-lablib.devices.PCO.SC2.PCOSC2Camera</code> method), 651
<code>is_acquisition_setup()</code>	( <code>py-lablib.devices.AlliedVision.Bonito.IBonitoCamera</code> method), 454	<code>is_acquisition_setup()</code>	( <code>py-lablib.devices.Photometrics.pvcam.PvcamCamera</code> method), 667
<code>is_acquisition_setup()</code>	( <code>py-lablib.devices.Andor.AndorSDK2.AndorSDK2Camera</code> method), 474	<code>is_acquisition_setup()</code>	( <code>py-lablib.devices.PhotonFocus.PhotonFocus.IPhotonFocusCamera</code> method), 676
<code>is_acquisition_setup()</code>	( <code>py-lablib.devices.Andor.AndorSDK3.AndorSDK3Camera</code> method), 484	<code>is_acquisition_setup()</code>	( <code>py-lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusBitFlowCa</code> method), 698
<code>is_acquisition_setup()</code>	( <code>py-lablib.devices.Basler.pylon.BaslerPylonCamera</code> method), 518	<code>is_acquisition_setup()</code>	( <code>py-lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusIMAQCam</code> method), 682
<code>is_acquisition_setup()</code>	( <code>py-lablib.devices.BitFlow.BitFlow.BitFlowCamera</code> method), 530	<code>is_acquisition_setup()</code>	( <code>py-lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusSiSoCamer</code> method), 691
<code>is_acquisition_setup()</code>	( <code>py-lablib.devices.BitFlow.BitFlow.BitFlowFrameGrabber</code> method), 525	<code>is_acquisition_setup()</code>	( <code>py-lablib.devices.PrincetonInstruments.picam.PicamCamera</code> method), 716
<code>is_acquisition_setup()</code>	( <code>py-lablib.devices.DCAM.DCAM.DCAMCamera</code> method), 555	<code>is_acquisition_setup()</code>	( <code>py-lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareCamera</code> method), 730
<code>is_acquisition_setup()</code>	( <code>py-lablib.devices.IMAQ.IMAQ.IMAQCamera</code> method), 572	<code>is_acquisition_setup()</code>	( <code>py-lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareFrameGrabb</code> method), 725
<code>is_acquisition_setup()</code>	( <code>py-lablib.devices.IMAQ.IMAQ.IMAQFrameGrabber</code> method), 567	<code>is_acquisition_setup()</code>	( <code>py-lablib.devices.Thorlabs.TLCamera.ThorlabsTLCamera</code> method), 774
<code>is_acquisition_setup()</code>	( <code>py-lablib.devices.IMAQdx.IMAQdx.EthernetIMAQdxCamera</code> method), 587	<code>is_acquisition_setup()</code>	( <code>py-lablib.devices.uc480.uc480.UC480Camera</code> method), 871
<code>is_acquisition_setup()</code>	( <code>py-</code>	<code>is_ascending()</code>	(in module <code>py-</code>

<code>lablib.core.dataproc.utils</code> ), 135	<code>lablib.devices.AWG.generic.GenericAWG</code>
<code>is_ascii()</code> ( <code>pylablib.core.devio.data_format.DataFormat</code>	<code>method</code> ), 404
<code>method</code> ), 162	<code>is_burst_enabled()</code> ( <code>py-</code>
<code>is_batch_job_running()</code> ( <code>py-</code>	<code>lablib.devices.AWG.specific.Agilent33220A</code>
<code>lablib.core.thread.controller.QTaskThread</code>	<code>method</code> ), 417
<code>method</code> ), 304	<code>is_burst_enabled()</code> ( <code>py-</code>
<code>is_bool_array()</code> ( <code>in module py-</code>	<code>lablib.devices.AWG.specific.Agilent33500</code>
<code>lablib.core.utils.indexing</code> ), 381	<code>method</code> ), 411
<code>is_branch_path()</code> ( <code>py-</code>	<code>is_burst_enabled()</code> ( <code>py-</code>
<code>lablib.core.utils.dictionary.Dictionary</code> <code>method</code> ),	<code>lablib.devices.AWG.specific.InstekAFG2000</code>
327	<code>method</code> ), 429
<code>is_branch_path()</code> ( <code>py-</code>	<code>is_burst_enabled()</code> ( <code>py-</code>
<code>lablib.core.utils.dictionary.DictionaryPointer</code>	<code>lablib.devices.AWG.specific.InstekAFG2225</code>
<code>method</code> ), 339	<code>method</code> ), 423
<code>is_branch_path()</code> ( <code>py-</code>	<code>is_burst_enabled()</code> ( <code>py-</code>
<code>lablib.core.utils.dictionary.FilterTree</code> <code>method</code> ),	<code>lablib.devices.AWG.specific.RigolDG1000</code>
356	<code>method</code> ), 447
<code>is_branch_path()</code> ( <code>py-</code>	<code>is_burst_enabled()</code> ( <code>py-</code>
<code>lablib.core.utils.dictionary.PrefixTree</code> <code>method</code> ),	<code>lablib.devices.AWG.specific.RSInstekAFG21000</code>
348	<code>method</code> ), 435
<code>is_branch_valid()</code> ( <code>py-</code>	<code>is_burst_enabled()</code> ( <code>py-</code>
<code>lablib.core.fileio.dict_entry.DictEntryParser</code>	<code>lablib.devices.AWG.specific.TektronixAFG1000</code>
<code>method</code> ), 169	<code>method</code> ), 441
<code>is_branch_valid()</code> ( <code>py-</code>	<code>is_call_done()</code> ( <code>py-</code>
<code>lablib.core.fileio.dict_entry.ExpandedContainerDictionaryEntry</code>	<code>pylablib.core.thread.callsync.QCallResultSynchronizer</code>
<code>class method</code> ), 176	<code>method</code> ), 281
<code>is_branch_valid()</code> ( <code>py-</code>	<code>is_call_done()</code> ( <code>py-</code>
<code>lablib.core.fileio.dict_entry.ExternalBinTableDictionaryEntry</code>	<code>pylablib.core.thread.callsync.QDirectResultSynchronizer</code>
<code>class method</code> ), 174	<code>method</code> ), 282
<code>is_branch_valid()</code> ( <code>py-</code>	<code>is_CFR_enabled()</code> ( <code>py-</code>
<code>lablib.core.fileio.dict_entry.ExternalNumpyDictionaryEntry</code>	<code>lablib.devices.PhotonFocus.PhotonFocus.IPhotonFocusCamera</code>
<code>class method</code> ), 176	<code>method</code> ), 673
<code>is_branch_valid()</code> ( <code>py-</code>	<code>is_CFR_enabled()</code> ( <code>py-</code>
<code>lablib.core.fileio.dict_entry.ExternalTextTableDictionaryEntry</code>	<code>pylablib.devices.PhotonFocus.PhotonFocus.PhotonFocusBitFlowCa</code>
<code>class method</code> ), 173	<code>method</code> ), 698
<code>is_branch_valid()</code> ( <code>py-</code>	<code>is_CFR_enabled()</code> ( <code>py-</code>
<code>lablib.core.fileio.dict_entry.IDictionaryEntry</code>	<code>lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusIMACam</code>
<code>class method</code> ), 170	<code>method</code> ), 682
<code>is_branch_valid()</code> ( <code>py-</code>	<code>is_CFR_enabled()</code> ( <code>py-</code>
<code>lablib.core.fileio.dict_entry.IDictionaryEntry</code>	<code>lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusSiSoCamer</code>
<code>method</code> ), 170	<code>method</code> ), 691
<code>is_branch_valid()</code> ( <code>py-</code>	<code>is_channel_enabled()</code> ( <code>py-</code>
<code>lablib.core.fileio.dict_entry.IExternalFileDictionaryEntry</code>	<code>lablib.devices.Tektronix.base.DPO2000</code>
<code>class method</code> ), 175	<code>method</code> ), 765
<code>is_branch_valid()</code> ( <code>py-</code>	<code>is_channel_enabled()</code> ( <code>py-</code>
<code>lablib.core.fileio.dict_entry.IExternalTableDictionaryEntry</code>	<code>lablib.devices.Tektronix.base.ITektronixScope</code>
<code>class method</code> ), 172	<code>method</code> ), 750
<code>is_branch_valid()</code> ( <code>py-</code>	<code>is_channel_enabled()</code> ( <code>py-</code>
<code>lablib.core.fileio.dict_entry.InlineTableDictionaryEntry</code>	<code>lablib.devices.Tektronix.base.TDS2000</code>
<code>class method</code> ), 172	<code>method</code> ), 758
<code>is_branch_valid()</code> ( <code>py-</code>	<code>is_channel_enabled()</code> ( <code>py-</code>
<code>lablib.core.fileio.dict_entry.ITableDictionaryEntry</code>	<code>lablib.devices.Toptica.ibeam.TopticaIBeam</code>
<code>class method</code> ), 171	<code>method</code> ), 824
<code>is_burst_enabled()</code> ( <code>py-</code>	<code>is_command</code> ( <code>pylablib.devices.Andor.AndorSDK3.AndorSDK3Attribute</code>

<a href="#">attribute</a> ), 478	<a href="#">PFCantAttribute</a> class method), 172
<a href="#">is_command</a> ( <a href="#">pylablib.devices.PhotonFocus.PhotonFocus.PFCantAttribute</a> attribute), 671	<a href="#">pylablib.core.fileio.dict_entry.ITableDictionaryEntry</a> class method), 171
<a href="#">is_connected</a> () ( <a href="#">pylablib.core.utils.net.ClientSocket</a> method), 390	<a href="#">is_descending</a> () ( <a href="#">in module py-lablib.core.dataproc.utils</a> ), 135
<a href="#">is_connected</a> () ( <a href="#">py-lablib.devices.Attocube.anc350.ANC350</a> method), 507	<a href="#">is_dict_entry_branch</a> () ( <a href="#">in module py-lablib.core.fileio.dict_entry</a> ), 169
<a href="#">is_continuous</a> () ( <a href="#">py-lablib.devices.Tektronix.base.DPO2000</a> method), 765	<a href="#">is_dictionary</a> () ( <a href="#">in module py-lablib.core.utils.dictionary</a> ), 326
<a href="#">is_continuous</a> () ( <a href="#">py-lablib.devices.Tektronix.base.ITektronixScope</a> method), 749	<a href="#">is_dictionary</a> () ( <a href="#">py-lablib.core.utils.dictionary.Dictionary</a> static method), 326
<a href="#">is_continuous</a> () ( <a href="#">py-lablib.devices.Tektronix.base.TDS2000</a> method), 758	<a href="#">is_dictionary</a> () ( <a href="#">py-lablib.core.utils.dictionary.DictionaryPointer</a> static method), 339
<a href="#">is_convertible</a> () ( <a href="#">in module py-lablib.core.utils.string</a> ), 399	<a href="#">is_dictionary</a> () ( <a href="#">py-lablib.core.utils.dictionary.FilterTree</a> static method), 356
<a href="#">is_cooler_on</a> () ( <a href="#">py-lablib.devices.Andor.AndorSDK2.AndorSDK2Camera</a> method), 467	<a href="#">is_dictionary</a> () ( <a href="#">py-lablib.core.utils.dictionary.PrefixTree</a> static method), 348
<a href="#">is_cooler_on</a> () ( <a href="#">py-lablib.devices.Andor.AndorSDK3.AndorSDK3Camera</a> method), 480	<a href="#">is_diffuser_in</a> () ( <a href="#">py-lablib.devices.Ophir.base.VegaPowerMeter</a> method), 642
<a href="#">is_data_valid</a> () ( <a href="#">py-lablib.core.fileio.dict_entry.DictEntryBuilder</a> method), 169	<a href="#">is_drift_compensation_enabled</a> () ( <a href="#">py-lablib.devices.PhysikInstrumente.base.PIE516</a> method), 705
<a href="#">is_data_valid</a> () ( <a href="#">py-lablib.core.fileio.dict_entry.ExpandedContainerDictionaryEntry</a> class method), 176	<a href="#">is_enabled</a> () ( <a href="#">pylablib.devices.Arcus.performax.Performax2EXStage</a> method), 495
<a href="#">is_data_valid</a> () ( <a href="#">py-lablib.core.fileio.dict_entry.ExternalBinTableDictionaryEntry</a> class method), 174	<a href="#">is_enabled</a> () ( <a href="#">pylablib.devices.Arcus.performax.Performax4EXStage</a> method), 490
<a href="#">is_data_valid</a> () ( <a href="#">py-lablib.core.fileio.dict_entry.ExternalNumpyDictionaryEntry</a> class method), 176	<a href="#">is_enabled</a> () ( <a href="#">pylablib.devices.Arcus.performax.PerformaxDMXJSStage</a> method), 497
<a href="#">is_data_valid</a> () ( <a href="#">py-lablib.core.fileio.dict_entry.ExternalTextTableDictionaryEntry</a> class method), 173	<a href="#">is_enabled</a> () ( <a href="#">pylablib.devices.Attocube.anc300.ANC300</a> method), 503
<a href="#">is_data_valid</a> () ( <a href="#">py-lablib.core.fileio.dict_entry.IDictionaryEntry</a> class method), 170	<a href="#">is_enabled</a> () ( <a href="#">pylablib.devices.Attocube.anc350.ANC350</a> method), 507
<a href="#">is_data_valid</a> () ( <a href="#">py-lablib.core.fileio.dict_entry.IDictionaryEntry</a> method), 170	<a href="#">is_enabled</a> () ( <a href="#">pylablib.devices.Lakeshore.base.Lakeshore218</a> method), 593
<a href="#">is_data_valid</a> () ( <a href="#">py-lablib.core.fileio.dict_entry.IExternalFileDictionaryEntry</a> class method), 175	<a href="#">is_enabled</a> () ( <a href="#">pylablib.devices.LaserQuantum.base.Finesse</a> method), 605
<a href="#">is_data_valid</a> () ( <a href="#">py-lablib.core.fileio.dict_entry.IExternalTableDictionaryEntry</a> class method), 172	<a href="#">is_enabled</a> () ( <a href="#">pylablib.devices.LighthousePhotonics.base.SproutG</a> method), 612
<a href="#">is_data_valid</a> () ( <a href="#">py-lablib.core.fileio.dict_entry.InlineTableDictionaryEntry</a> method), 642	<a href="#">is_enabled</a> () ( <a href="#">pylablib.devices.Pfeiffer.base.TPG260</a> method), 655
	<a href="#">is_enabled</a> () ( <a href="#">pylablib.devices.Topptica.ibeam.ToppticaIBeam</a> method), 824
	<a href="#">is_fastpiezo_locked</a> () ( <a href="#">py-lablib.devices.Sirah.Matisse.SirahMatisse</a> method), 739
	<a href="#">is_filter_in</a> () ( <a href="#">py-lablib.devices.Ophir.base.VegaPowerMeter</a> method), 642

[is\\_frame\\_transfer\\_enabled\(\)](#) (pylablib.devices.Andor.AndorSDK2.AndorSDK2Camera method), 471  
[is\\_free\(\)](#) (pylablib.core.fileio.location.FolderFileSystemDataLocation method), 188  
[is\\_free\(\)](#) (pylablib.core.fileio.location.IDataLocation method), 184  
[is\\_free\(\)](#) (pylablib.core.fileio.location.IFileSystemDataLocation method), 185  
[is\\_free\(\)](#) (pylablib.core.fileio.location.OpenedFileLocation method), 185  
[is\\_free\(\)](#) (pylablib.core.fileio.location.PrefixedFileSystemDataLocation method), 187  
[is\\_free\(\)](#) (pylablib.core.fileio.location.SingleFileSystemDataLocation method), 186  
[is\\_grabbing\(\)](#) (pylablib.devices.Tektronix.base.DPO2000 method), 765  
[is\\_grabbing\(\)](#) (pylablib.devices.Tektronix.base.ITektronixScope method), 749  
[is\\_grabbing\(\)](#) (pylablib.devices.Tektronix.base.TDS2000 method), 758  
[is\\_gui\\_running\(\)](#) (in module pylablib.core.thread.threadprop), 320  
[is\\_gui\\_thread\(\)](#) (in module pylablib.core.thread.threadprop), 320  
[is\\_handled\\_widget\(\)](#) (in module pylablib.core.gui.value\_handling), 276  
[is\\_homed\(\)](#) (pylablib.devices.Thorlabs.kinesis.KinesisMotor method), 796  
[is\\_homing\(\)](#) (pylablib.devices.Thorlabs.kinesis.KinesisMotor method), 796  
[is\\_in\\_controlled\(\)](#) (pylablib.core.thread.controller.QTaskThread method), 309  
[is\\_in\\_controlled\(\)](#) (pylablib.core.thread.controller.QThreadController method), 300  
[is\\_integer\(\)](#) (in module pylablib.core.gui.formatter), 261  
[is\\_layout\\_column\\_empty\(\)](#) (in module pylablib.core.gui.utils), 264  
[is\\_layout\\_row\\_empty\(\)](#) (in module pylablib.core.gui.utils), 264  
[is\\_leaf\\_path\(\)](#) (pylablib.core.utils.dictionary.Dictionary method), 327  
[is\\_leaf\\_path\(\)](#) (pylablib.core.utils.dictionary.DictionaryPointer method), 339  
[is\\_leaf\\_path\(\)](#) (pylablib.core.utils.dictionary.FilterTree method), 378  
[is\\_linear\(\)](#) (in module pylablib.core.dataproc.utils), 135  
[is\\_looping\(\)](#) (pylablib.devices.Basler.pylon.BaslerPylonCamera.Schema method), 107  
[is\\_mandatory\\_arg\(\)](#) (pylablib.core.dataproc.callable.FunctionCallable method), 107  
[is\\_mandatory\\_arg\(\)](#) (pylablib.core.dataproc.callable.ICallable method), 104  
[is\\_mandatory\\_arg\(\)](#) (pylablib.core.dataproc.callable.JoinedCallable method), 106  
[is\\_mandatory\\_arg\(\)](#) (pylablib.core.dataproc.callable.MethodCallable method), 108  
[is\\_mandatory\\_arg\(\)](#) (pylablib.core.dataproc.callable.MultiplexedCallable method), 105  
[is\\_measurement\\_running\(\)](#) (pylablib.devices.HighFinesse.wlm.WLM method), 558  
[is\\_metadata\\_enabled\(\)](#) (pylablib.devices.Andor.AndorSDK3.AndorSDK3Camera method), 481  
[is\\_metadata\\_enabled\(\)](#) (pylablib.devices.Photometrics.pvcam.PvcamCamera method), 663  
[is\\_metadata\\_enabled\(\)](#) (pylablib.devices.PrincetonInstruments.picam.PicamCamera method), 713  
[is\\_moving\(\)](#) (pylablib.devices.Arcus.performax.Performax2EXStage method), 495  
[is\\_moving\(\)](#) (pylablib.devices.Arcus.performax.Performax4EXStage method), 491  
[is\\_moving\(\)](#) (pylablib.devices.Arcus.performax.PerformaxDMXJSStage method), 498  
[is\\_moving\(\)](#) (pylablib.devices.Attocube.anc300.ANC300 method), 505  
[is\\_moving\(\)](#) (pylablib.devices.Attocube.anc350.ANC350 method), 507  
[is\\_moving\(\)](#) (pylablib.devices.Newport.picomotor.Picomotor8742 method), 630  
[is\\_moving\(\)](#) (pylablib.devices.SmarAct.scu3d.SC3D method), 747  
[is\\_moving\(\)](#) (pylablib.devices.Thorlabs.kinesis.KinesisMotor method), 796  
[is\\_moving\(\)](#) (pylablib.devices.Thorlabs.kinesis.KinesisPiezoMotor method), 801  
[is\\_moving\(\)](#) (pylablib.devices.Trinamic.base.TMCM1110 method), 801

*method*), 829  
 is\_online\_enabled() (pylablib.devices.PhysikInstrumente.base.GenericPIController *method*), 429  
*method*), 703  
 is\_online\_enabled() (pylablib.devices.PhysikInstrumente.base.PIE516 *method*), 707  
 is\_opened() (pylablib.core.devio.comm\_backend.FT232RLDeviceBackend *method*), 150  
 is\_opened() (pylablib.core.devio.comm\_backend.ICommBackendWrapper *method*), 161  
 is\_opened() (pylablib.core.devio.comm\_backend.IDeviceCommBackend *method*), 143  
 is\_opened() (pylablib.core.devio.comm\_backend.NetworkDeviceBackend *method*), 153  
 is\_opened() (pylablib.core.devio.comm\_backend.PyUSBDeviceBackend *method*), 156  
 is\_opened() (pylablib.core.devio.comm\_backend.RecordedDeviceBackend *method*), 158  
 is\_opened() (pylablib.core.devio.comm\_backend.SerialDeviceBackend *method*), 148  
 is\_opened() (pylablib.core.devio.comm\_backend.VisaDeviceBackend *method*), 145  
 is\_opened() (pylablib.core.devio.interface.IDevice *method*), 162  
 is\_opened() (pylablib.core.devio.SCPI.SCPIDevice *method*), 141  
 is\_opened() (pylablib.devices.AlliedVision.Bonito.BonitoIMAQCamera *method*), 460  
 is\_opened() (pylablib.devices.AlliedVision.Bonito.IBonitoCamera *method*), 454  
 is\_opened() (pylablib.devices.Andor.AndorSDK2.AndorSDK2Camera *method*), 467  
 is\_opened() (pylablib.devices.Andor.AndorSDK3.AndorSDK3Camera *method*), 479  
 is\_opened() (pylablib.devices.Arcus.performax.GenericPerformaxStage *method*), 488  
 is\_opened() (pylablib.devices.Arcus.performax.Performax2EXStage *method*), 495  
 is\_opened() (pylablib.devices.Arcus.performax.Performax4EXStage *method*), 493  
 is\_opened() (pylablib.devices.Arcus.performax.PerformaxDMXISASStage *method*), 499  
 is\_opened() (pylablib.devices.Arduino.base.IArduinoDevice *method*), 502  
 is\_opened() (pylablib.devices.Attocube.anc300.ANC300 *method*), 506  
 is\_opened() (pylablib.devices.Attocube.anc350.ANC350 *method*), 509  
 is\_opened() (pylablib.devices.AWG.generic.GenericAWGs *method*), 406  
 is\_opened() (pylablib.devices.AWG.specific.Agilent33220A *method*), 417  
 is\_opened() (pylablib.devices.AWG.specific.Agilent33500s *method*), 411  
 is\_opened() (pylablib.devices.AWG.specific.InstekAFG2000 *method*), 429  
 is\_opened() (pylablib.devices.AWG.specific.InstekAFG2225 *method*), 423  
 is\_opened() (pylablib.devices.AWG.specific.RigolDG1000 *method*), 447  
 is\_opened() (pylablib.devices.AWG.specific.RSInstekAFG21000 *method*), 435  
 is\_opened() (pylablib.devices.AWG.specific.TektronixAFG1000 *method*), 441  
 is\_opened() (pylablib.devices.Basler.pylon.BaslerPylonCamera *method*), 514  
 is\_opened() (pylablib.devices.BitFlow.BitFlow.BitFlowCamera *method*), 530  
 is\_opened() (pylablib.devices.BitFlow.BitFlow.BitFlowFrameGrabber *method*), 521  
 is\_opened() (pylablib.devices.Conrad.base.RelayBoard *method*), 535  
 is\_opened() (pylablib.devices.Cryocon.base.CryoconIx *method*), 537  
 is\_opened() (pylablib.devices.Cryomagnetics.base.LM500 *method*), 542  
 is\_opened() (pylablib.devices.Cryomagnetics.base.LM510 *method*), 545  
 is\_opened() (pylablib.devices.DCAM.DCAM.DCAMCamera *method*), 550  
 is\_opened() (pylablib.devices.HighFinesse.wlm.WLM *method*), 557  
 is\_opened() (pylablib.devices.IMAQ.IMAQ.IMAQCamera *method*), 572  
 is\_opened() (pylablib.devices.IMAQ.IMAQ.IMAQFrameGrabber *method*), 561  
 is\_opened() (pylablib.devices.IMAQdx.IMAQdx.EthernetIMAQdxCamera *method*), 587  
 is\_opened() (pylablib.devices.IMAQdx.IMAQdx.IMAQdxCamera *method*), 578  
 is\_opened() (pylablib.devices.interface.camera.IAttributeCamera *method*), 839  
 is\_opened() (pylablib.devices.interface.camera.IBinROICamera *method*), 858  
 is\_opened() (pylablib.devices.interface.camera.ICamera *method*), 835  
 is\_opened() (pylablib.devices.interface.camera.IExposureCamera *method*), 849  
 is\_opened() (pylablib.devices.interface.camera.IGrabberAttributeCamera *method*), 844  
 is\_opened() (pylablib.devices.interface.camera.IROICamera *method*), 854  
 is\_opened() (pylablib.devices.interface.stage.IMultiaxisStage *method*), 864  
 is\_opened() (pylablib.devices.interface.stage.IStage *method*), 862  
 is\_opened() (pylablib.devices.KJL.base.KJL300 *method*), 862





<code>is_output_enabled()</code> (py-lablib.devices.AWG.specific.InstekAFG2000 method), 429	<code>is_running()</code> (pylablib.devices.BitFlow.BitFlow.BitFlowFrameGrabber method), 523
<code>is_output_enabled()</code> (py-lablib.devices.AWG.specific.InstekAFG2225 method), 423	<code>is_running()</code> (pylablib.devices.NI.daq.NIDAQ method), 623
<code>is_output_enabled()</code> (py-lablib.devices.AWG.specific.RigolDG1000 method), 447	<code>is_running()</code> (pylablib.devices.PhotonFocus.PhotonFocus.PhotonFocus method), 695
<code>is_output_enabled()</code> (py-lablib.devices.AWG.specific.RSInstekAFG21000 method), 435	<code>is_sequence()</code> (in module py-lablib.core.utils.funcargparse), 369
<code>is_output_enabled()</code> (py-lablib.devices.AWG.specific.TektronixAFG1000 method), 441	<code>is_servo_enabled()</code> (py-lablib.devices.PhysikInstrumente.base.PIE516 method), 705
<code>is_path_valid()</code> (in module py-lablib.core.utils.files), 362	<code>is_shutter_opened()</code> (py-lablib.devices.LaserQuantum.base.Finesse method), 605
<code>is_peer_closed()</code> (py-lablib.core.utils.ipc.SharedMemIPCTable method), 385	<code>is_slice()</code> (in module pylablib.core.utils.indexing), 381
<code>is_peer_connected()</code> (py-lablib.core.utils.ipc.SharedMemIPCTable method), 385	<code>is_status_line_enabled()</code> (py-lablib.devices.AlliedVision.Bonito.BonitoIMAQCamera method), 460
<code>is_pixel_correction_enabled()</code> (py-lablib.devices.PCO.SC2.PCOS2Camera method), 648	<code>is_status_line_enabled()</code> (py-lablib.devices.AlliedVision.Bonito.IBonitoCamera method), 452
<code>is_range()</code> (in module pylablib.core.utils.indexing), 381	<code>is_status_line_enabled()</code> (py-lablib.devices.PhotonFocus.PhotonFocus.IPhotonFocusCamera method), 673
<code>is_running()</code> (pylablib.core.gui.widgets.container.IQContainer method), 201	<code>is_status_line_enabled()</code> (py-lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusBitFlowCamera method), 698
<code>is_running()</code> (pylablib.core.gui.widgets.container.IQWidgetContainer method), 207	<code>is_status_line_enabled()</code> (py-lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusIMAQCamera method), 683
<code>is_running()</code> (pylablib.core.gui.widgets.container.QContainer method), 204	<code>is_status_line_enabled()</code> (py-lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusSiSoCamera method), 691
<code>is_running()</code> (pylablib.core.gui.widgets.container.QDialogContainer method), 219	<code>is_stopping()</code> (py-lablib.core.gui.widgets.container.IQContainer method), 202
<code>is_running()</code> (pylablib.core.gui.widgets.container.QFrameContainer method), 215	<code>is_stopping()</code> (py-lablib.core.gui.widgets.container.IQWidgetContainer method), 208
<code>is_running()</code> (pylablib.core.gui.widgets.container.QGroupBoxContainer method), 223	<code>is_stopping()</code> (py-lablib.core.gui.widgets.container.QContainer method), 204
<code>is_running()</code> (pylablib.core.gui.widgets.container.QScrollAreaContainer method), 230	<code>is_stopping()</code> (py-lablib.core.gui.widgets.container.QContainedWidget method), 204
<code>is_running()</code> (pylablib.core.gui.widgets.container.QScrollAreaContainer method), 227	<code>is_stopping()</code> (py-lablib.core.gui.widgets.container.QDialogContainer method), 219
<code>is_running()</code> (pylablib.core.gui.widgets.container.QTabContainer method), 233	<code>is_stopping()</code> (py-lablib.core.gui.widgets.container.QWidgetContainer method), 211
<code>is_running()</code> (pylablib.core.gui.widgets.container.QWidgetContainer method), 211	<code>is_stopping()</code> (py-lablib.core.gui.widgets.container.QFrameContainer method), 215
<code>is_running()</code> (pylablib.core.gui.widgets.param_table.ParamTable method), 250	<code>is_stopping()</code> (py-lablib.core.gui.widgets.container.QGroupBoxContainer method), 223
<code>is_running()</code> (pylablib.core.gui.widgets.param_table.StatusTable method), 259	<code>is_stopping()</code> (py-lablib.core.gui.widgets.container.QScrollAreaContainer method), 230
<code>is_running()</code> (pylablib.devices.BitFlow.BitFlow.BitFlowCameraBufferManager method), 527	

<i>method</i> ), 230	<i>method</i> ), 208
is_stopping() (py-lablib.core.gui.widgets.container.QScrollAreaContainer.QComboBoxWidget <i>method</i> ), 227	is_timer_running() (py-lablib.core.gui.widgets.container.QContainer <i>method</i> ), 204
is_stopping() (py-lablib.core.gui.widgets.container.QTabContainer <i>method</i> ), 233	is_timer_running() (py-lablib.core.gui.widgets.container.QDialogContainer <i>method</i> ), 219
is_stopping() (py-lablib.core.gui.widgets.container.QWidgetContainer <i>method</i> ), 211	is_timer_running() (py-lablib.core.gui.widgets.container.QFrameContainer <i>method</i> ), 215
is_stopping() (py-lablib.core.gui.widgets.param_table.ParamTable <i>method</i> ), 250	is_timer_running() (py-lablib.core.gui.widgets.container.QGroupBoxContainer <i>method</i> ), 223
is_stopping() (py-lablib.core.gui.widgets.param_table.StatusTable <i>method</i> ), 259	is_timer_running() (py-lablib.core.gui.widgets.container.QScrollAreaContainer <i>method</i> ), 230
is_switcher_channel_enabled() (py-lablib.devices.HighFinesse.wlm.WLM <i>method</i> ), 559	is_timer_running() (py-lablib.core.gui.widgets.container.QScrollAreaContainer.QContainer <i>method</i> ), 227
is_switcher_channel_shown() (py-lablib.devices.HighFinesse.wlm.WLM <i>method</i> ), 559	is_timer_running() (py-lablib.core.gui.widgets.container.QTabContainer <i>method</i> ), 233
is_sync_output_enabled() (py-lablib.devices.AWG.generic.GenericAWG <i>method</i> ), 403	is_timer_running() (py-lablib.core.gui.widgets.container.QWidgetContainer <i>method</i> ), 211
is_sync_output_enabled() (py-lablib.devices.AWG.specific.Agilent33220A <i>method</i> ), 417	is_timer_running() (py-lablib.core.gui.widgets.param_table.ParamTable <i>method</i> ), 250
is_sync_output_enabled() (py-lablib.devices.AWG.specific.Agilent33500 <i>method</i> ), 411	is_timer_running() (py-lablib.core.gui.widgets.param_table.StatusTable <i>method</i> ), 259
is_sync_output_enabled() (py-lablib.devices.AWG.specific.InstekAFG2000 <i>method</i> ), 429	is_trigger_output_enabled() (py-lablib.devices.AWG.generic.GenericAWG <i>method</i> ), 405
is_sync_output_enabled() (py-lablib.devices.AWG.specific.InstekAFG2225 <i>method</i> ), 423	is_trigger_output_enabled() (py-lablib.devices.AWG.specific.Agilent33220A <i>method</i> ), 417
is_sync_output_enabled() (py-lablib.devices.AWG.specific.RigolDG1000 <i>method</i> ), 447	is_trigger_output_enabled() (py-lablib.devices.AWG.specific.Agilent33500 <i>method</i> ), 411
is_sync_output_enabled() (py-lablib.devices.AWG.specific.RSInstekAFG21000 <i>method</i> ), 435	is_trigger_output_enabled() (py-lablib.devices.AWG.specific.InstekAFG2000 <i>method</i> ), 429
is_sync_output_enabled() (py-lablib.devices.AWG.specific.TektronixAFG1000 <i>method</i> ), 441	is_trigger_output_enabled() (py-lablib.devices.AWG.specific.InstekAFG2225 <i>method</i> ), 423
is_target_reached() (py-lablib.devices.Attocube.anc350.ANC350 <i>method</i> ), 508	is_trigger_output_enabled() (py-lablib.devices.AWG.specific.RigolDG1000 <i>method</i> ), 447
is_timer_running() (py-lablib.core.gui.widgets.container.IQContainer <i>method</i> ), 200	is_trigger_output_enabled() (py-lablib.devices.AWG.specific.RSInstekAFG21000 <i>method</i> ), 435
is_timer_running() (py-lablib.core.gui.widgets.container.IQWidgetContainer <i>method</i> ), 208	is_trigger_output_enabled() (py-lablib.devices.AWG.specific.TektronixAFG1000 <i>method</i> ), 441



- method*), 441
- `is_unprintable_character()` (in module `pylablib.core.fileio.loadfile_utils`), 182
- `is_velocity_control_enabled()` (pylablib.devices.PhysikInstrumente.base.PIE516 *method*), 705
- `is_wait_done()` (pylablib.devices.interface.camera.FrameCounter *method*), 835
- `is_wavemeter_connected()` (pylablib.devices.M2.solstis.Solstis *method*), 615
- `is_wavemeter_lock_on()` (pylablib.devices.M2.solstis.Solstis *method*), 615
- `ISingleValueHandler` (class in `pylablib.core.gui.value_handling`), 269
- `ISkipableNotifier` (class in `pylablib.core.thread.notifier`), 316
- `ispan()` (pylablib.core.dataproc.image.ROI *method*), 121
- `IStage` (class in `pylablib.devices.interface.stage`), 862
- `IStage.NoParameterCaller` (class in `pylablib.devices.interface.stage`), 862
- `ITableDictionaryEntry` (class in `pylablib.core.fileio.dict_entry`), 171
- `ITektronixScope` (class in `pylablib.devices.Tektronix.base`), 748
- `ITektronixScope.NoParameterCaller` (class in `pylablib.devices.Tektronix.base`), 752
- `ItemAccessor` (class in `pylablib.core.utils.dictionary`), 360
- `items()` (pylablib.core.utils.dictionary.Dictionary *method*), 328
- `items()` (pylablib.core.utils.dictionary.DictionaryPointer *method*), 339
- `items()` (pylablib.core.utils.dictionary.FilterTree *method*), 356
- `items()` (pylablib.core.utils.dictionary.PrefixTree *method*), 348
- `iter_sublayout_items()` (pylablib.core.gui.widgets.container.IQWidgetContainer *method*), 208
- `iter_sublayout_items()` (pylablib.core.gui.widgets.container.QDialogContainer *method*), 219
- `iter_sublayout_items()` (pylablib.core.gui.widgets.container.QFrameContainer *method*), 215
- `iter_sublayout_items()` (pylablib.core.gui.widgets.container.QGroupBoxContainer *method*), 223
- `iter_sublayout_items()` (pylablib.core.gui.widgets.container.QScrollAreaContainer *method*), 227
- `iter_sublayout_items()` (pylablib.core.gui.widgets.container.QWidgetContainer *method*), 211
- `iter_sublayout_items()` (pylablib.core.gui.widgets.layout\_manager.IQLayoutManagedWidget *method*), 239
- `iter_sublayout_items()` (pylablib.core.gui.widgets.layout\_manager.QLayoutManagedWidget *method*), 241
- `iter_sublayout_items()` (pylablib.core.gui.widgets.param\_table.ParamTable *method*), 250
- `iter_sublayout_items()` (pylablib.core.gui.widgets.param\_table.StatusTable *method*), 259
- `iteritems()` (pylablib.core.utils.dictionary.Dictionary *method*), 328
- `iteritems()` (pylablib.core.utils.dictionary.DictionaryPointer *method*), 339
- `iteritems()` (pylablib.core.utils.dictionary.FilterTree *method*), 356
- `iteritems()` (pylablib.core.utils.dictionary.PrefixTree *method*), 348
- `iterkeys()` (pylablib.core.utils.dictionary.Dictionary *method*), 330
- `iterkeys()` (pylablib.core.utils.dictionary.DictionaryPointer *method*), 339
- `iterkeys()` (pylablib.core.utils.dictionary.FilterTree *method*), 356
- `iterkeys()` (pylablib.core.utils.dictionary.PrefixTree *method*), 348
- `internodes()` (pylablib.core.utils.dictionary.Dictionary *method*), 330
- `internodes()` (pylablib.core.utils.dictionary.DictionaryPointer *method*), 340
- `internodes()` (pylablib.core.utils.dictionary.FilterTree *method*), 357
- `internodes()` (pylablib.core.utils.dictionary.PrefixTree *method*), 349
- `intervalvalues()` (pylablib.core.utils.dictionary.Dictionary *method*), 329
- `intervalvalues()` (pylablib.core.utils.dictionary.DictionaryPointer *method*), 340
- `intervalvalues()` (pylablib.core.utils.dictionary.FilterTree *method*), 357
- `intervalvalues()` (pylablib.core.utils.dictionary.PrefixTree *method*), 349
- `ITextInputFileFormat` (class in `pylablib.core.fileio.loadfile`), 176
- `ITextOutputFileFormat` (class in `pylablib.core.fileio.savefile`), 190
- `ITR90` (class in `pylablib.devices.Leybold.base`), 609
- `ITR90.NoParameterCaller` (class in `pylablib.devices.Leybold.base`), 609

`lablib.devices.Leybold.base`), 609  
`IValueHandler` (class in `py-lablib.core.gui.value_handling`), 265  
`ivalues` (`pylablib.devices.Andor.AndorSDK3.AndorSDK3Attribute` attribute), 477  
`ivalues` (`pylablib.devices.Basler.pylon.BaslerPylonAttribute` attribute), 513  
`ivalues` (`pylablib.devices.DCAM.DCAM.DCAMAttribute` attribute), 549  
`ivalues` (`pylablib.devices.IMAQdx.IMAQdx.IMAQdxAttribute` attribute), 577  
`ivalues` (`pylablib.devices.Photometrics.pvcam.PvcamAttribute` attribute), 660  
`ivalues` (`pylablib.devices.PhotonFocus.PhotonFocus.PFCamAttribute` attribute), 671  
`ivalues` (`pylablib.devices.PrincetonInstruments.picam.PicamAttribute` attribute), 711  
`ivalues` (`pylablib.devices.SiliconSoftware.fgrab.FGrabAttribute` attribute), 720  
  
**J**  
`job` (`pylablib.core.thread.controller.QTaskThread.TBatchJob` attribute), 302  
`jog` () (`pylablib.devices.Arcus.performax.Performax2EXStage` method), 495  
`jog` () (`pylablib.devices.Arcus.performax.Performax4EXStage` method), 491  
`jog` () (`pylablib.devices.Arcus.performax.PerformaxDMXJSASStage` method), 498  
`jog` () (`pylablib.devices.Attocube.anc300.ANC300` method), 505  
`jog` () (`pylablib.devices.Attocube.anc350.ANC350` method), 509  
`jog` () (`pylablib.devices.Newport.picomotor.Picomotor8742` method), 629  
`jog` () (`pylablib.devices.Thorlabs.kinesis.KinesisMotor` method), 796  
`jog` () (`pylablib.devices.Thorlabs.kinesis.KinesisPiezoMotor` method), 801  
`jog` () (`pylablib.devices.Trinamic.base.TMCM1110` method), 827  
`JoinedCallable` (class in `py-lablib.core.dataproc.callable`), 105  
`JoinedCallable.NamesBoundCall` (class in `py-lablib.core.dataproc.callable`), 106  
`jspan` () (`pylablib.core.dataproc.image.ROI` method), 121  
  
**K**  
`kernel` (`pylablib.core.dataproc.feature.Peak` attribute), 109  
`keyPressEvent` () (`py-lablib.core.gui.widgets.edit.NumEdit` method), 234  
`keyPressEvent` () (`py-lablib.core.gui.widgets.edit.TextEdit` method), 234  
`keys` () (`pylablib.core.utils.dictionary.DictionaryPointer` method), 340  
`keys` () (`pylablib.core.utils.dictionary.FilterTree` method), 357  
`keys` () (`pylablib.core.utils.dictionary.PrefixTree` method), 349  
`kind` (`pylablib.devices.Andor.AndorSDK3.AndorSDK3Attribute` attribute), 477  
`kind` (`pylablib.devices.Basler.pylon.BaslerPylonAttribute` attribute), 512  
`kind` (`pylablib.devices.DCAM.DCAM.DCAMAttribute` attribute), 549  
`kind` (`pylablib.devices.IMAQdx.IMAQdx.IMAQdxAttribute` attribute), 577  
`kind` (`pylablib.devices.interface.camera.TStatusLineDescription` attribute), 861  
`kind` (`pylablib.devices.Photometrics.pvcam.PvcamAttribute` attribute), 660  
`kind` (`pylablib.devices.PhotonFocus.PhotonFocus.PFCamAttribute` attribute), 670  
`kind` (`pylablib.devices.PrincetonInstruments.picam.PicamAttribute` attribute), 709  
`kind` (`pylablib.devices.SiliconSoftware.fgrab.FGrabAttribute` attribute), 719  
`KinesisDevice` (class in `py-lablib.devices.Thorlabs.kinesis`), 787  
`KinesisDevice.CommData` (class in `py-lablib.devices.Thorlabs.kinesis`), 787  
`KinesisDevice.CommShort` (class in `py-lablib.devices.Thorlabs.kinesis`), 788  
`KinesisDevice.NoParameterCaller` (class in `py-lablib.devices.Thorlabs.kinesis`), 788  
`KinesisMotor` (class in `py-lablib.devices.Thorlabs.kinesis`), 794  
`KinesisMotor.CommData` (class in `py-lablib.devices.Thorlabs.kinesis`), 798  
`KinesisMotor.CommShort` (class in `py-lablib.devices.Thorlabs.kinesis`), 798  
`KinesisMotor.NoParameterCaller` (class in `py-lablib.devices.Thorlabs.kinesis`), 798  
`KinesisPiezoMotor` (class in `py-lablib.devices.Thorlabs.kinesis`), 801  
`KinesisPiezoMotor.CommData` (class in `py-lablib.devices.Thorlabs.kinesis`), 802  
`KinesisPiezoMotor.CommShort` (class in `py-lablib.devices.Thorlabs.kinesis`), 802  
`KinesisPiezoMotor.NoParameterCaller`

- (class in *pylablib.devices.Thorlabs.kinesis*), 803
- KinesisQuadDetector (class in *pylablib.devices.Thorlabs.kinesis*), 806
- KinesisQuadDetector.CommData (class in *pylablib.devices.Thorlabs.kinesis*), 807
- KinesisQuadDetector.CommShort (class in *pylablib.devices.Thorlabs.kinesis*), 807
- KinesisQuadDetector.NoParameterCaller (class in *pylablib.devices.Thorlabs.kinesis*), 807
- kinetic\_cycle\_time (pylablib.devices.Andor.AndorSDK2.TCycleTimings attribute), 466
- KJL300 (class in *pylablib.devices.KJL.base*), 590
- KJL300.NoParameterCaller (class in *pylablib.devices.KJL.base*), 591
- KJLBackendError, 589
- KJLError, 589
- ## L
- label (pylablib.core.gui.widgets.param\_table.ParamTable.ParamRow attribute), 243
- label (pylablib.core.gui.widgets.param\_table.StatusTable.ParamRow attribute), 252
- label (pylablib.core.utils.string.TConversionClass attribute), 398
- LabelIndicatorHandler (class in *pylablib.core.gui.value\_handling*), 277
- labels (pylablib.devices.Andor.AndorSDK3.AndorSDK3Attribute attribute), 477
- labels (pylablib.devices.Basler.pylon.BaslerPylonAttribute attribute), 513
- labels (pylablib.devices.DCAM.DCAM.DCAMAttribute attribute), 549
- labels (pylablib.devices.IMAQdx.IMAQdx.IMAQdxAttribute attribute), 577
- labels (pylablib.devices.Photometrics.pvcam.PvcamAttribute attribute), 660
- labels (pylablib.devices.PhotonFocus.PhotonFocus.PFCamAttribute attribute), 671
- labels (pylablib.devices.PrincetonInstruments.picam.PicamAttribute attribute), 711
- labels (pylablib.devices.SiliconSoftware.fgrab.FGrabAttribute attribute), 720
- LabelValueHandler (class in *pylablib.core.gui.value\_handling*), 270
- Lakeshore218 (class in *pylablib.devices.Lakeshore.base*), 593
- Lakeshore218.NoParameterCaller (class in *pylablib.devices.Lakeshore.base*), 595
- Lakeshore370 (class in *pylablib.devices.Lakeshore.base*), 599
- Lakeshore370.NoParameterCaller (class in *pylablib.devices.Lakeshore.base*), 600
- LakeshoreBackendError, 592
- LakeshoreError, 592
- laser (pylablib.devices.LighthousePhotonics.base.TWorkHours attribute), 612
- laser\_enabled (pylablib.devices.LaserQuantum.base.TWorkHours attribute), 604
- laser\_on (pylablib.devices.Toptica.ibeam.TWorkHours attribute), 824
- laser\_threshold (pylablib.devices.LaserQuantum.base.TWorkHours attribute), 604
- LaserQuantumBackendError, 603
- LaserQuantumError, 603
- latching\_trigger() (in module *pylablib.core.dataproc.feature*), 110
- layout (pylablib.core.gui.utils.TWidgetLocation attribute), 263
- left (pylablib.devices.interface.camera.TFramePosition attribute), 831
- level (pylablib.devices.Tektronix.base.TTriggerParameters attribute), 748
- LeyboldBackendError, 606
- LeyboldError, 606
- LibraryController (class in *pylablib.devices.Andor.AndorSDK2*), 465
- LibraryController (class in *pylablib.devices.Andor.AndorSDK3*), 476
- LibraryController (class in *pylablib.devices.Basler.pylon*), 510
- LibraryController (class in *pylablib.devices.DCAM.DCAM*), 548
- LibraryController (class in *pylablib.devices.Photometrics.pvcam*), 659
- LibraryController (class in *pylablib.devices.PhotonFocus.PhotonFocus*), 669
- LibraryController (class in *pylablib.devices.PrincetonInstruments.picam*), 708
- LibraryController (class in *pylablib.devices.SmarAct.scu3d*), 745
- LibraryController (class in *pylablib.devices.Thorlabs.TLCamera*), 769
- LibraryController (class in *pylablib.devices.utils.load\_lib*), 874
- LighthousePhotonicsBackendError, 611
- LighthousePhotonicsError, 611
- limit (pylablib.devices.Lakeshore.base.TLakeshore218CurveHeader attribute), 593
- limit() (pylablib.core.dataproc.image.ROI method), 121

limit_errors_enabled()	(py-	lablib.devices.Arcus.performax.Performax2EXStage	static method), 800
limit_errors_enabled()	(py-	lablib.devices.Arcus.performax.Performax4EXStage	static method), 804
limit_switch	(pylablib.devices.Thorlabs.kinesis.THomeParams	attribute), 785	lablib.devices.Thorlabs.kinesis.KinesisQuadDetector
limit_to_range()	(in module	py-	lablib.devices.Thorlabs.kinesis.MFF
LimitError, 262			static method), 793
LinearTransform	(class in	py-	list_dir() (in module pylablib.core.utils.files), 364
LineEditValueHandler	(class in	py-	list_dir_recursive() (in module py-
list_applets()	(in module	py-	lablib.core.utils.files), 365
list_backend_resources()	(in module	py-	list_kinesis_devices() (in module py-
list_boards()	(in module	py-	lablib.devices.Thorlabs.kinesis), 781
list_cameras()	(in module	py-	list_opened_files() (py-
list_cameras()	(in module	py-	lablib.devices.Basler.pylon), 511
list_cameras()	(in module	py-	lablib.core.fileio.location.FolderFileSystemDataLocation
list_cameras()	(in module	py-	method), 188
list_cameras()	(in module	py-	list_opened_files() (py-
list_cameras()	(in module	py-	lablib.devices.BitFlow.BitFlow), 521
list_cameras()	(in module	py-	lablib.core.fileio.location.IDataLocation
list_cameras()	(in module	py-	method), 185
list_cameras()	(in module	py-	list_opened_files() (py-
list_cameras()	(in module	py-	lablib.devices.IMAQ.IMAQ), 561
list_cameras()	(in module	py-	lablib.core.fileio.location.OpenedFileLocation
list_cameras()	(in module	py-	method), 185
list_cameras()	(in module	py-	list_opened_files() (py-
list_cameras()	(in module	py-	lablib.devices.IMAQdx.IMAQdx), 576
list_cameras()	(in module	py-	lablib.core.fileio.location.PrefixedFileSystemDataLocation
list_cameras()	(in module	py-	method), 187
list_cameras()	(in module	py-	list_opened_files() (py-
list_cameras()	(in module	py-	lablib.devices.PCO.SC2), 644
list_cameras()	(in module	py-	lablib.core.fileio.location.SingleFileSystemDataLocation
list_cameras()	(in module	py-	method), 186
list_cameras()	(in module	py-	list_resources() (py-
list_cameras()	(in module	py-	lablib.devices.Photometrics.pvcam), 659
list_cameras()	(in module	py-	lablib.core.devio.comm_backend.FT232DeviceBackend
list_cameras()	(in module	py-	static method), 151
list_cameras()	(in module	py-	list_resources() (py-
list_cameras()	(in module	py-	lablib.devices.PrincetonInstruments.picam), 709
list_cameras()	(in module	py-	lablib.core.devio.comm_backend.IDeviceCommBackend
list_cameras()	(in module	py-	static method), 144
list_cameras()	(in module	py-	list_resources() (py-
list_cameras()	(in module	py-	lablib.devices.Thorlabs.TLCamera), 769
list_cameras()	(in module	py-	lablib.core.devio.comm_backend.NetworkDeviceBackend
list_cameras()	(in module	py-	static method), 154
list_devices()	(in module	py-	list_resources() (py-
list_devices()	(in module	py-	lablib.devices.NI.daq), 621
list_devices()	(in module	py-	lablib.core.devio.comm_backend.PyUSBDeviceBackend
list_devices()	(in module	py-	static method), 156
list_devices()	(in module	py-	list_resources() (py-
list_devices()	(in module	py-	lablib.devices.SmarAct.scu3d), 745
list_devices()	(in module	py-	lablib.core.devio.comm_backend.RecordedDeviceBackend
list_devices()	(in module	py-	static method), 159
list_devices()	(in module	py-	list_resources() (py-
list_devices()	(in module	py-	lablib.devices.Thorlabs.kinesis.BasicKinesisDevice
list_devices()	(in module	py-	static method), 149
list_devices()	(in module	py-	lablib.devices.Thorlabs.kinesis.KinesisDevice
list_devices()	(in module	py-	static method), 789
list_devices()	(in module	py-	list_resources() (py-
list_devices()	(in module	py-	lablib.core.devio.comm_backend.VisaDeviceBackend



- static method*), 145
- `list_usb_performax_devices()` (in module `pylablib.devices.Arcus.performax`), 488
- `listen()` (in module `pylablib.core.utils.net`), 391
- `ListIndex` (class in `pylablib.core.utils.indexing`), 382
- `ListIndexNoSlice` (class in `pylablib.core.utils.indexing`), 382
- `LM500` (class in `pylablib.devices.Cryomagnetics.base`), 540
- `LM500.NoParameterCaller` (class in `pylablib.devices.Cryomagnetics.base`), 541
- `LM510` (class in `pylablib.devices.Cryomagnetics.base`), 544
- `LM510.NoParameterCaller` (class in `pylablib.devices.Cryomagnetics.base`), 544
- `load()` (in module `pylablib.core.utils.strdump`), 396
- `load()` (`pylablib.core.utils.strdump.StrDumper` method), 396
- `load()` (`pylablib.devices.BitFlow.BitFlow.CameraFileEditor` method), 532
- `load_bin()` (in module `pylablib.core.fileio.loadfile`), 180
- `load_bin_desc()` (in module `pylablib.core.fileio.loadfile`), 180
- `load_csv()` (in module `pylablib.core.fileio.loadfile`), 179
- `load_csv_desc()` (in module `pylablib.core.fileio.loadfile`), 179
- `load_dict()` (in module `pylablib.core.fileio.loadfile`), 180
- `load_file()` (`pylablib.core.fileio.dict_entry.ExternalNumPyDictEntry` class method), 175
- `load_file()` (`pylablib.core.fileio.dict_entry.IExternalFileDictEntry` class method), 175
- `load_generic()` (in module `pylablib.core.fileio.loadfile`), 181
- `load_lib()` (in module `pylablib.devices.utils.load_lib`), 873
- `load_logfile()` (in module `pylablib.core.devio.backend_logger`), 142
- `load_par()` (in module `pylablib`), 875
- `loads()` (in module `pylablib.core.utils.strdump`), 396
- `loads()` (`pylablib.core.utils.strdump.StrDumper` method), 396
- `loc` (`pylablib.core.fileio.location.LocationFile` attribute), 184
- `LocationFile` (class in `pylablib.core.fileio.location`), 184
- `LocationName` (class in `pylablib.core.fileio.location`), 183
- `lock()` (`pylablib.core.devio.comm_backend.FT232DeviceBackend` method), 152
- `lock()` (`pylablib.core.devio.comm_backend.ICommBackendWrapper` method), 161
- `lock()` (`pylablib.core.devio.comm_backend.IDeviceCommBackend` method), 143
- `lock()` (`pylablib.core.devio.comm_backend.NetworkDeviceBackend` method), 154
- `lock()` (`pylablib.core.devio.comm_backend.PyUSBDeviceBackend` method), 157
- `lock()` (`pylablib.core.devio.comm_backend.RecordedDeviceBackend` method), 159
- `lock()` (`pylablib.core.devio.comm_backend.SerialDeviceBackend` method), 149
- `lock()` (`pylablib.core.devio.comm_backend.VisaDeviceBackend` method), 145
- `lock()` (`pylablib.core.devio.SCPI.SCPIDevice` method), 141
- `lock()` (`pylablib.devices.Arduino.base.IArduinoDevice` method), 502
- `lock()` (`pylablib.devices.Attocube.anc300.ANC300` method), 506
- `lock()` (`pylablib.devices.Attocube.anc350.ANC350` method), 509
- `lock()` (`pylablib.devices.AWG.generic.GenericAWG` method), 407
- `lock()` (`pylablib.devices.AWG.specific.Agilent33220A` method), 417
- `lock()` (`pylablib.devices.AWG.specific.Agilent33500` method), 411
- `lock()` (`pylablib.devices.AWG.specific.InstekAFG2000` method), 429
- `lock()` (`pylablib.devices.AWG.specific.InstekAFG2225` method), 423
- `lock()` (`pylablib.devices.AWG.specific.RigolDG1000` method), 447
- `lock()` (`pylablib.devices.AWG.specific.RSInstekAFG21000` method), 435
- `lock()` (`pylablib.devices.AWG.specific.TektronixAFG1000` method), 441
- `lock()` (`pylablib.devices.Conrad.base.RelayBoard` method), 535
- `lock()` (`pylablib.devices.Cryocon.base.Cryocon1x` method), 537
- `lock()` (`pylablib.devices.Cryomagnetics.base.LM500` method), 542
- `lock()` (`pylablib.devices.Cryomagnetics.base.LM510` method), 545
- `lock()` (`pylablib.devices.KJL.base.KJL300` method), 591
- `lock()` (`pylablib.devices.Lakeshore.base.Lakeshore218` method), 596
- `lock()` (`pylablib.devices.Lakeshore.base.Lakeshore370` method), 601
- `lock()` (`pylablib.devices.LaserQuantum.base.Finesse` method), 606
- `lock()` (`pylablib.devices.Leybold.base.GenericITR` method), 608

<code>lock()</code> ( <code>pylablib.devices.Leybold.base.ITR90</code> method), 610	<code>lock()</code> ( <code>pylablib.devices.Thorlabs.serial.ThorlabsSerialInterface</code> method), 811
<code>lock()</code> ( <code>pylablib.devices.LighthousePhotonics.base.SproutLock</code> method), 613	<code>lock()</code> ( <code>pylablib.devices.Toptica.ibeam.TopticaIBeam</code> method), 826
<code>lock()</code> ( <code>pylablib.devices.Newport.picomotor.Picomotor8742</code> method), 631	<code>lock()</code> ( <code>pylablib.devices.Trinamic.base.TMCM1110</code> method), 829
<code>lock()</code> ( <code>pylablib.devices.Ophir.base.OphirDevice</code> method), 639	<code>lock_etalon()</code> ( <code>pylablib.devices.M2.solstis.Solstis</code> method), 616
<code>lock()</code> ( <code>pylablib.devices.Ophir.base.VegaPowerMeter</code> method), 643	<code>lock_reference_cavity()</code> ( <code>pylablib.devices.M2.solstis.Solstis</code> method), 617
<code>lock()</code> ( <code>pylablib.devices.OZOptics.base.DD100</code> method), 636	<code>lock_wavemeter()</code> ( <code>pylablib.devices.M2.solstis.Solstis</code> method), 615
<code>lock()</code> ( <code>pylablib.devices.OZOptics.base.EPC04</code> method), 638	<code>locking()</code> ( <code>pylablib.core.devio.comm_backend.FT232DeviceBackend</code> method), 152
<code>lock()</code> ( <code>pylablib.devices.OZOptics.base.OZOpticsDevice</code> method), 633	<code>locking()</code> ( <code>pylablib.core.devio.comm_backend.ICommBackendWrapper</code> method), 161
<code>lock()</code> ( <code>pylablib.devices.OZOptics.base.TF100</code> method), 634	<code>locking()</code> ( <code>pylablib.core.devio.comm_backend.IDeviceCommBackend</code> method), 143
<code>lock()</code> ( <code>pylablib.devices.Pfeiffer.base.DPG202</code> method), 658	<code>locking()</code> ( <code>pylablib.core.devio.comm_backend.NetworkDeviceBackend</code> method), 154
<code>lock()</code> ( <code>pylablib.devices.Pfeiffer.base.TPG260</code> method), 657	<code>locking()</code> ( <code>pylablib.core.devio.comm_backend.PyUSBDeviceBackend</code> method), 157
<code>lock()</code> ( <code>pylablib.devices.PhysikInstrumente.base.GenericPIController</code> method), 704	<code>locking()</code> ( <code>pylablib.core.devio.comm_backend.RecordedDeviceBackend</code> method), 159
<code>lock()</code> ( <code>pylablib.devices.PhysikInstrumente.base.PIE516</code> method), 707	<code>locking()</code> ( <code>pylablib.core.devio.comm_backend.SerialDeviceBackend</code> method), 149
<code>lock()</code> ( <code>pylablib.devices.Sirah.Matisse.SirahMatisse</code> method), 741	<code>locking()</code> ( <code>pylablib.core.devio.comm_backend.VisaDeviceBackend</code> method), 145
<code>lock()</code> ( <code>pylablib.devices.Tektronix.base.DPO2000</code> method), 765	<code>locking()</code> ( <code>pylablib.core.devio.SCPI.SCPIDevice</code> method), 141
<code>lock()</code> ( <code>pylablib.devices.Tektronix.base.ITektronixScope</code> method), 753	<code>locking()</code> ( <code>pylablib.devices.Arduino.base.IArduinoDevice</code> method), 502
<code>lock()</code> ( <code>pylablib.devices.Tektronix.base.TDS2000</code> method), 758	<code>locking()</code> ( <code>pylablib.devices.Attocube.anc300.ANC300</code> method), 506
<code>lock()</code> ( <code>pylablib.devices.Thorlabs.elliptec.ElliptecMotor</code> method), 780	<code>locking()</code> ( <code>pylablib.devices.Attocube.anc350.ANC350</code> method), 510
<code>lock()</code> ( <code>pylablib.devices.Thorlabs.kinesis.BasicKinesisDevice</code> method), 784	<code>locking()</code> ( <code>pylablib.devices.AWG.generic.GenericAWG</code> method), 407
<code>lock()</code> ( <code>pylablib.devices.Thorlabs.kinesis.KinesisDevice</code> method), 789	<code>locking()</code> ( <code>pylablib.devices.AWG.specific.Agilent33220A</code> method), 417
<code>lock()</code> ( <code>pylablib.devices.Thorlabs.kinesis.KinesisMotor</code> method), 800	<code>locking()</code> ( <code>pylablib.devices.AWG.specific.Agilent33500</code> method), 411
<code>lock()</code> ( <code>pylablib.devices.Thorlabs.kinesis.KinesisPiezoMotor</code> method), 804	<code>locking()</code> ( <code>pylablib.devices.AWG.specific.InstekAFG2000</code> method), 429
<code>lock()</code> ( <code>pylablib.devices.Thorlabs.kinesis.KinesisQuadDetector</code> method), 809	<code>locking()</code> ( <code>pylablib.devices.AWG.specific.InstekAFG2225</code> method), 423
<code>lock()</code> ( <code>pylablib.devices.Thorlabs.kinesis.MFF</code> method), 793	<code>locking()</code> ( <code>pylablib.devices.AWG.specific.RigolDG1000</code> method), 447
<code>lock()</code> ( <code>pylablib.devices.Thorlabs.serial.FW</code> method), 814	<code>locking()</code> ( <code>pylablib.devices.AWG.specific.RSInstekAFG21000</code> method), 435
<code>lock()</code> ( <code>pylablib.devices.Thorlabs.serial.FWv1</code> method), 818	<code>locking()</code> ( <code>pylablib.devices.AWG.specific.TektronixAFG1000</code> method), 441
<code>lock()</code> ( <code>pylablib.devices.Thorlabs.serial.MDT69xA</code> method), 821	

locking() (pylablib.devices.Conrad.base.RelayBoard method), 535

locking() (pylablib.devices.Cryocon.base.Cryocon1x method), 537

locking() (pylablib.devices.Cryomagnetics.base.LM500 method), 542

locking() (pylablib.devices.Cryomagnetics.base.LM510 method), 545

locking() (pylablib.devices.KJL.base.KJL300 method), 591

locking() (pylablib.devices.Lakeshore.base.Lakeshore218 method), 596

locking() (pylablib.devices.Lakeshore.base.Lakeshore370 method), 601

locking() (pylablib.devices.LaserQuantum.base.Finesse method), 606

locking() (pylablib.devices.Leybold.base.GenericITR method), 608

locking() (pylablib.devices.Leybold.base.ITR90 method), 610

locking() (pylablib.devices.LighthousePhotonics.base.Sprout method), 613

locking() (pylablib.devices.Newport.picomotor.Picomotor8742 method), 631

locking() (pylablib.devices.Ophir.base.OphirDevice method), 639

locking() (pylablib.devices.Ophir.base.VegaPowerMeter method), 643

locking() (pylablib.devices.OZOptics.base.DD100 method), 636

locking() (pylablib.devices.OZOptics.base.EPC04 method), 638

locking() (pylablib.devices.OZOptics.base.OZOpticsDevice method), 633

locking() (pylablib.devices.OZOptics.base.TF100 method), 634

locking() (pylablib.devices.Pfeiffer.base.DPG202 method), 658

locking() (pylablib.devices.Pfeiffer.base.TPG260 method), 657

locking() (pylablib.devices.PhysikInstrumente.base.GenericPICO method), 704

locking() (pylablib.devices.PhysikInstrumente.base.PIE116 method), 707

locking() (pylablib.devices.Sirah.Matisse.SirahMatisse method), 741

locking() (pylablib.devices.Tektronix.base.DPO2000 method), 765

locking() (pylablib.devices.Tektronix.base.ITektronixScope method), 753

locking() (pylablib.devices.Tektronix.base.TDS2000 method), 758

locking() (pylablib.devices.Thorlabs.elliptec.ElliptecMotor method), 781

locking() (pylablib.devices.Thorlabs.kinesis.BasicKinesisDevice method), 784

locking() (pylablib.devices.Thorlabs.kinesis.KinesisDevice method), 789

locking() (pylablib.devices.Thorlabs.kinesis.KinesisMotor method), 800

locking() (pylablib.devices.Thorlabs.kinesis.KinesisPiezoMotor method), 804

locking() (pylablib.devices.Thorlabs.kinesis.KinesisQuadDetector method), 809

locking() (pylablib.devices.Thorlabs.kinesis.MFF method), 793

locking() (pylablib.devices.Thorlabs.serial.FW method), 814

locking() (pylablib.devices.Thorlabs.serial.FWv1 method), 818

locking() (pylablib.devices.Thorlabs.serial.MDT69xA method), 821

locking() (pylablib.devices.Thorlabs.serial.ThorlabsSerialInterface method), 811

locking() (pylablib.devices.Toptica.ibeam.TopticaIBeam method), 826

locking() (pylablib.devices.Trinamic.base.TMCM1110 method), 829

lockpoint (pylablib.devices.Sirah.Matisse.TFastpiezoCilParameters attribute), 735

log() (pylablib.core.devio.backend\_logger.BackendLogger method), 142

logerror() (in module pylablib.core.devio.comm\_backend), 142

loop (pylablib.core.gui.widgets.container.TTimerEvent attribute), 200

loop (pylablib.devices.Thorlabs.elliptec.TMotorInfo attribute), 777

lorentzian\_k() (in module pylablib.core.dataproc.specfunc), 124

low\_pass\_filter() (in module pylablib.core.dataproc.filters), 111

low\_thresh (pylablib.devices.Pfeiffer.base.TTPG260SwitchSettings attribute), 654

low\_thresh (pylablib.devices.Lakeshore.base.TLakeshore218AnalogSettings attribute), 592

low\_value (pylablib.devices.Lakeshore.base.TLakeshore370AnalogSettings attribute), 598

lower\_limit (pylablib.devices.Sirah.Matisse.TRefcellWaveformParameters attribute), 735

lower\_limit (pylablib.devices.Sirah.Matisse.TScanParameters attribute), 736

M (Index)

m (pylablib.core.thread.controller.QTaskThread attribute), 301

M2CommunicationError, 614

M2Error, 614

[make\\_comment\\_line\(\)](#) (pylablib.core.fileio.savefile.CSVTableOutputFileFormat method), 192  
[make\\_comment\\_line\(\)](#) (pylablib.core.fileio.savefile.DictionaryOutputFileFormat method), 193  
[make\\_comment\\_line\(\)](#) (pylablib.core.fileio.savefile.ITextOutputFileFormat method), 191  
[make\\_flat\\_nametuple\(\)](#) (in module pylablib.core.utils.general), 375  
[make\\_prop\\_line\(\)](#) (pylablib.core.fileio.savefile.CSVTableOutputFileFormat method), 192  
[make\\_prop\\_line\(\)](#) (pylablib.core.fileio.savefile.DictionaryOutputFileFormat method), 193  
[make\\_prop\\_line\(\)](#) (pylablib.core.fileio.savefile.ITextOutputFileFormat method), 191  
[make\\_savetime\\_line\(\)](#) (pylablib.core.fileio.savefile.CSVTableOutputFileFormat method), 192  
[make\\_savetime\\_line\(\)](#) (pylablib.core.fileio.savefile.DictionaryOutputFileFormat method), 193  
[make\\_savetime\\_line\(\)](#) (pylablib.core.fileio.savefile.ITextOutputFileFormat method), 191  
[make\\_sequence\(\)](#) (in module pylablib.core.utils.funcargparse), 369  
[man\\_value](#) (pylablib.devices.Lakeshore.base.TLakeshore218AnalogSensor attribute), 592  
[man\\_value](#) (pylablib.devices.Lakeshore.base.TLakeshore370AnalogSensor attribute), 598  
[mandatory\\_args\\_num\(\)](#) (pylablib.core.utils.functions.FunctionSignature method), 370  
[manufacturer](#) (pylablib.devices.PhotonFocus.PhotonFocus.TC2000 attribute), 670  
[manufacturer](#) (pylablib.devices.uc480.uc480.TDeviceInfo attribute), 865  
[map\\_container\(\)](#) (in module pylablib.core.utils.general), 374  
[map\\_dict\\_keys\(\)](#) (in module pylablib.core.utils.general), 375  
[map\\_dict\\_values\(\)](#) (in module pylablib.core.utils.general), 375  
[map\\_self\(\)](#) (pylablib.core.utils.dictionary.Dictionary method), 332  
[map\\_self\(\)](#) (pylablib.core.utils.dictionary.DictionaryPoint method), 340  
[map\\_self\(\)](#) (pylablib.core.utils.dictionary.FilterTree method), 357  
[map\\_self\(\)](#) (pylablib.core.utils.dictionary.PrefixTree method), 349  
[mark\\_unscheduled\(\)](#) (pylablib.core.thread.controller.QTaskThread.Job method), 303  
[match\(\)](#) (pylablib.core.utils.dictionary.FilterTree method), 352  
[MatisseTuner](#) (class in pylablib.devices.Sirah.tuner), 743  
[max](#) (pylablib.devices.Andor.AndorSDK3.AndorSDK3Attribute attribute), 477  
[max](#) (pylablib.devices.Basler.pylon.BaslerPylonAttribute attribute), 512  
[max](#) (pylablib.devices.DCAM.DCAM.DCAMAttribute attribute), 549  
[max](#) (pylablib.devices.IMAQdx.IMAQdx.IMAQdxAttribute attribute), 577  
[max](#) (pylablib.devices.interface.camera.TAxisROILimit attribute), 851  
[max](#) (pylablib.devices.Photometrics.pvcam.PvcamAttribute attribute), 660  
[max](#) (pylablib.devices.PhotonFocus.PhotonFocus.PFCamAttribute attribute), 671  
[max](#) (pylablib.devices.PrincetonInstruments.picam.PicamAttribute attribute), 710  
[max](#) (pylablib.devices.SiliconSoftware.fgrab.FGrabAttribute attribute), 720  
[max\\_args\\_num\(\)](#) (pylablib.core.utils.functions.FunctionSignature method), 370  
[max\\_velocity](#) (pylablib.devices.Thorlabs.kinesis.TJogParams attribute), 784  
[max\\_velocity](#) (pylablib.devices.Thorlabs.kinesis.TVelocityParams attribute), 784  
[max\\_voltage](#) (pylablib.devices.Thorlabs.kinesis.TPZMotorDriveParams attribute), 786  
[maxbin](#) (pylablib.devices.interface.camera.TAxisROILimit attribute), 851  
[MDT69xA.NoParameterCaller](#) (class in pylablib.devices.Thorlabs.serial), 820  
[measure\\_capacitance\(\)](#) (pylablib.devices.Attocube.anc300.ANC300 method), 503  
[measure\\_level\(\)](#) (pylablib.devices.Cryomagnetics.base.LM500 method), 540  
[measure\\_level\(\)](#) (pylablib.devices.Cryomagnetics.base.LM510 method), 546  
[median\\_filter\(\)](#) (in module pylablib.core.dataproc.filters), 111  
[merge\(\)](#) (in module pylablib.core.dataproc.utils), 136  
[merge\(\)](#) (pylablib.core.utils.dictionary.Dictionary method), 357



method), 331

merge() (pylablib.core.utils.dictionary.DictionaryPointer method), 341

merge() (pylablib.core.utils.dictionary.FilterTree method), 358

merge() (pylablib.core.utils.dictionary.PrefixTree method), 350

merge() (pylablib.core.utils.functions.FunctionSignature static method), 371

merge\_dicts() (in module pylablib.core.utils.general), 375

messageID (pylablib.devices.Thorlabs.kinesis.BasicKinesisDevice.CommData attribute), 782

messageID (pylablib.devices.Thorlabs.kinesis.BasicKinesisDevice.CommShort attribute), 782

messageID (pylablib.devices.Thorlabs.kinesis.KinesisDevice.CommData attribute), 788

messageID (pylablib.devices.Thorlabs.kinesis.KinesisDevice.CommShort attribute), 788

messageID (pylablib.devices.Thorlabs.kinesis.KinesisMotor.CommData attribute), 798

messageID (pylablib.devices.Thorlabs.kinesis.KinesisMotor.CommShort attribute), 798

messageID (pylablib.devices.Thorlabs.kinesis.KinesisPiezoMotor.CommData attribute), 802

messageID (pylablib.devices.Thorlabs.kinesis.KinesisPiezoMotor.CommShort attribute), 802

messageID (pylablib.devices.Thorlabs.kinesis.KinesisQuadData (pylablib.devices.Thorlabs.kinesis.QuadData attribute), 807

messageID (pylablib.devices.Thorlabs.kinesis.KinesisQuadData (pylablib.devices.Thorlabs.kinesis.QuadData attribute), 807

messageID (pylablib.devices.Thorlabs.kinesis.MFF.CommData attribute), 792

messageID (pylablib.devices.Thorlabs.kinesis.MFF.CommShort attribute), 792

MethodCallable (class in pylablib.core.dataproc.callable), 107

MethodCallable.NamesBoundCall (class in pylablib.core.dataproc.callable), 107

MethodObjectCall (class in pylablib.core.utils.functions), 372

MethodObjectProperty (class in pylablib.core.utils.functions), 373

MFF (class in pylablib.devices.Thorlabs.kinesis), 790

MFF.CommData (class in pylablib.devices.Thorlabs.kinesis), 791

MFF.CommShort (class in pylablib.devices.Thorlabs.kinesis), 792

MFF.NoParameterCaller (class in pylablib.devices.Thorlabs.kinesis), 792

millisecond (pylablib.devices.uc480.uc480.TTimestamp attribute), 866

min (pylablib.devices.Andor.AndorSDK3.AndorSDK3Attribute attribute), 477

min (pylablib.devices.Basler.pylon.BaslerPylonAttribute attribute), 512

min (pylablib.devices.DCAM.DCAM.DCAMAttribute attribute), 549

min (pylablib.devices.IMAQdx.IMAQdx.IMAQdxAttribute attribute), 577

min (pylablib.devices.interface.camera.TAxisROILimit attribute), 851

min (pylablib.devices.Photometrics.pvcam.PvcamAttribute attribute), 660

min (pylablib.devices.PhotonFocus.PhotonFocus.PFCamAttribute attribute), 671

min (pylablib.devices.PrincetonInstruments.picam.PicamAttribute attribute), 710

min (pylablib.devices.SiliconSoftware.fgrab.FGrabAttribute attribute), 719

min\_run\_time (pylablib.core.thread.controller.QTaskThread.TBatchJob attribute), 302

min\_velocity (pylablib.devices.Thorlabs.kinesis.TJogParams attribute), 784

min\_velocity (pylablib.devices.Thorlabs.kinesis.TVelocityParams attribute), 784

minute (pylablib.devices.uc480.uc480.TTimestamp attribute), 866

MissingGUIHandlerError, 278

motorCover (pylablib.devices.Thorlabs.kinesis.TDeviceInfo attribute), 781

quadData (pylablib.devices.Lakeshore.base.TLakeshore218AnalogSettings attribute), 592

quadData (pylablib.devices.Lakeshore.base.TLakeshore370AnalogSettings attribute), 598

quadData (pylablib.devices.Ophir.base.TWavelengthInfo attribute), 640

short (pylablib.devices.Sirah.Matisse.TRefcellWaveformParameters attribute), 735

short (pylablib.devices.Sirah.Matisse.TScanParameters attribute), 736

mode (pylablib.devices.Thorlabs.kinesis.TJogParams attribute), 784

mode (pylablib.devices.Thorlabs.kinesis.TPZMotorJogParams attribute), 787

model (pylablib.devices.Basler.pylon.TCameraInfo attribute), 511

model (pylablib.devices.Basler.pylon.TDeviceInfo attribute), 514

model (pylablib.devices.BitFlow.BitFlow.TDeviceInfo attribute), 521

model (pylablib.devices.DCAM.DCAM.TDeviceInfo attribute), 550

model (pylablib.devices.HighFinesse.wlm.TDeviceInfo attribute), 557

model (pylablib.devices.IMAQdx.IMAQdx.TCameraInfo attribute), 576

model (pylablib.devices.IMAQdx.IMAQdx.TDeviceInfo attribute), 576

attribute), 578  
 model (pylablib.devices.NI.daq.TDeviceInfo attribute), 620  
 model (pylablib.devices.PCO.SC2.TDeviceInfo attribute), 644  
 model (pylablib.devices.PhotonFocus.PhotonFocus.TDeviceInfo attribute), 672  
 model (pylablib.devices.PrincetonInstruments.picam.TCameraInfo attribute), 709  
 model (pylablib.devices.PrincetonInstruments.picam.TDeviceInfo attribute), 711  
 model (pylablib.devices.Thorlabs.TLCamera.TDeviceInfo attribute), 770  
 model (pylablib.devices.uc480.uc480.TCameraInfo attribute), 864  
 model (pylablib.devices.uc480.uc480.TDeviceInfo attribute), 865  
 model\_no (pylablib.devices.Thorlabs.elliptec.TDeviceInfo attribute), 777  
 model\_no (pylablib.devices.Thorlabs.kinesis.TDeviceInfo attribute), 781  
 module (pylablib.devices.Trinamic.base.TMCM1110.ReplyData attribute), 827  
 month (pylablib.devices.uc480.uc480.TTimestamp attribute), 866  
 motor (pylablib.devices.Thorlabs.elliptec.TMotorInfo attribute), 777  
 mousePressEvent () (pylablib.core.gui.widgets.label.EnumLabel method), 236  
 mousePressEvent () (pylablib.core.gui.widgets.label.NumLabel method), 237  
 mousePressEvent () (pylablib.core.gui.widgets.label.TextLabel method), 236  
 move\_by () (pylablib.devices.Arcus.performax.Performax2EXStage method), 496  
 move\_by () (pylablib.devices.Arcus.performax.Performax4EXStage method), 490  
 move\_by () (pylablib.devices.Arcus.performax.PerformaxDMXJSASStage method), 497  
 move\_by () (pylablib.devices.Attocube.anc300.ANC300 method), 505  
 move\_by () (pylablib.devices.Attocube.anc350.ANC350 method), 508  
 move\_by () (pylablib.devices.Newport.picomotor.Picomotor8742 method), 629  
 move\_by () (pylablib.devices.PhysikInstrumente.base.PIE516 method), 706  
 move\_by () (pylablib.devices.SmarAct.scu3d.SCU3D method), 746  
 move\_by () (pylablib.devices.Thorlabs.elliptec.ElliptecMotor method), 780  
 move\_by () (pylablib.devices.Thorlabs.kinesis.KinesisMotor method), 796  
 move\_by () (pylablib.devices.Thorlabs.kinesis.KinesisPiezoMotor method), 801  
 move\_by () (pylablib.devices.Trinamic.base.TMCM1110 method), 827  
 move\_by\_steps () (pylablib.devices.Attocube.anc350.ANC350 method), 508  
 move\_dir () (in module pylablib.core.utils.files), 365  
 move\_file () (in module pylablib.core.utils.files), 363  
 move\_macrostep () (pylablib.devices.SmarAct.scu3d.SCU3D method), 746  
 move\_to () (pylablib.core.utils.dictionary.DictionaryPointer method), 335  
 move\_to () (pylablib.devices.Arcus.performax.Performax2EXStage method), 496  
 move\_to () (pylablib.devices.Arcus.performax.Performax4EXStage method), 490  
 move\_to () (pylablib.devices.Arcus.performax.PerformaxDMXJSASStage method), 497  
 move\_to () (pylablib.devices.Attocube.anc350.ANC350 method), 508  
 move\_to () (pylablib.devices.Newport.picomotor.Picomotor8742 method), 629  
 move\_to () (pylablib.devices.PhysikInstrumente.base.PIE516 method), 706  
 move\_to () (pylablib.devices.Thorlabs.elliptec.ElliptecMotor method), 780  
 move\_to () (pylablib.devices.Thorlabs.kinesis.KinesisMotor method), 796  
 move\_to () (pylablib.devices.Thorlabs.kinesis.KinesisPiezoMotor method), 801  
 move\_to () (pylablib.devices.Trinamic.base.TMCM1110 method), 827  
 move\_to\_state () (pylablib.devices.Thorlabs.kinesis.MFF method), 791  
 move\_up () (pylablib.core.utils.dictionary.DictionaryPointer method), 335  
 multi\_scale\_peakdet () (in module pylablib.core.dataproc.feature), 109  
 MulticastPool (class in pylablib.core.thread.multicast\_pool), 315  
 MultiplexedCallable (class in pylablib.core.dataproc.callable), 104  
 MultiplexedCallable.NamesBoundCall (class in pylablib.core.dataproc.callable), 105  
 multiplied () (pylablib.core.dataproc.transform.Indexed2DTransform method), 135  
 multiplied () (pylablib.core.dataproc.transform.LinearTransform method), 134  
 muxaddr () (in module py-

[lablib.devices.Newport.picomotor](#)), 627  
[muxaddr\(\)](#) (in module [pylablib.devices.Thorlabs.elliptec](#)), 776  
[muxaxis\(\)](#) (in module [pylablib.devices.Attocube.anc300](#)), 502  
[muxaxis\(\)](#) (in module [pylablib.devices.interface.stage](#)), 863  
[muxcall\(\)](#) (in module [pylablib.core.utils.general](#)), 380  
[muxchan\(\)](#) (in module [pylablib.devices.Toptica.ibeam](#)), 823  
[muxchannel\(\)](#) (in module [pylablib.devices.HighFinesse.wlm](#)), 557  
[muxchannel\(\)](#) (in module [pylablib.devices.Tektronix.base](#)), 748  
[muxchannel\(\)](#) (in module [pylablib.devices.Thorlabs.kinesis](#)), 787

## N

[name](#) ([pylablib.core.fileio.location.LocationFileAttribute](#)), 184  
[name](#) ([pylablib.core.gui.widgets.container.TChildAttribute](#)), 200  
[name](#) ([pylablib.core.gui.widgets.container.TTimerAttribute](#)), 199  
[name](#) ([pylablib.core.utils.files.TempFileAttribute](#)), 363  
[name](#) ([pylablib.devices.Andor.AndorSDK3.AndorSDK3AttributeAttribute](#)), 477  
[name](#) ([pylablib.devices.Basler.pylon.BaslerPylonAttributeAttribute](#)), 512  
[name](#) ([pylablib.devices.Basler.pylon.TCameraInfoAttribute](#)), 511  
[name](#) ([pylablib.devices.Basler.pylon.TDeviceInfoAttribute](#)), 514  
[name](#) ([pylablib.devices.DCAM.DCAM.DCAMAttributeAttribute](#)), 549  
[name](#) ([pylablib.devices.IMAQdx.IMAQdx.IMAQdxAttributeAttribute](#)), 577  
[name](#) ([pylablib.devices.IMAQdx.IMAQdx.TCameraInfoAttribute](#)), 576  
[name](#) ([pylablib.devices.Lakeshore.base.TLakeshore218CurveHeaderAttribute](#)), 593  
[name](#) ([pylablib.devices.NI.daq.TDeviceInfoAttribute](#)), 620  
[name](#) ([pylablib.devices.Ophir.base.TDeviceInfoAttribute](#)), 640  
[name](#) ([pylablib.devices.Ophir.base.THeadInfoAttribute](#)), 640  
[name](#) ([pylablib.devices.Photometrics.pvcam.PvcamAttributeAttribute](#)), 660  
[name](#) ([pylablib.devices.PhotonFocus.PhotonFocus.PFCameraAttributeAttribute](#)), 670  
[name](#) ([pylablib.devices.PrincetonInstruments.picam.PicamAttributeAttribute](#)), 709  
[name](#) ([pylablib.devices.PrincetonInstruments.picam.TCameraInfoAttribute](#)), 709  
[name](#) ([pylablib.devices.PrincetonInstruments.picam.TDeviceInfoAttribute](#)), 711  
[name](#) ([pylablib.devices.SiliconSoftware.fgrab.FGrabAttributeAttribute](#)), 719  
[name](#) ([pylablib.devices.SiliconSoftware.fgrab.TAppletInfoAttribute](#)), 718  
[name](#) ([pylablib.devices.SiliconSoftware.fgrab.TBoardInfoAttribute](#)), 718  
[name](#) ([pylablib.devices.SiliconSoftware.fgrab.TFullAppletInfoAttribute](#)), 719  
[name](#) ([pylablib.devices.Thorlabs.TLCamera.TDeviceInfoAttribute](#)), 770  
[NamedUIDGenerator](#) (class in [pylablib.core.utils.general](#)), 378  
[nchannels](#) ([pylablib.devices.Thorlabs.kinesis.TDeviceInfoAttribute](#)), 781  
[ncycles](#) ([pylablib.devices.Thorlabs.kinesis.TKCubeTrigPosParamsAttribute](#)), 786  
[ndim\(\)](#) ([pylablib.core.dataproc.table\\_wrap.Array1DWrapperMethod](#)), 127  
[ndim\(\)](#) ([pylablib.core.dataproc.table\\_wrap.Array2DWrapperMethod](#)), 131  
[ndim\(\)](#) ([pylablib.core.dataproc.table\\_wrap.DataFrame2DWrapperMethod](#)), 133  
[ndim\(\)](#) ([pylablib.core.dataproc.table\\_wrap.I1DWrapperMethod](#)), 126  
[ndim\(\)](#) ([pylablib.core.dataproc.table\\_wrap.I2DWrapperMethod](#)), 129  
[ndim\(\)](#) ([pylablib.core.dataproc.table\\_wrap.IGenWrapperMethod](#)), 125  
[ndim\(\)](#) ([pylablib.core.dataproc.table\\_wrap.Series1DWrapperMethod](#)), 128  
[NetworkDeviceBackend](#) (class in [pylablib.core.devio.comm\\_backend](#)), 152  
[new\\_backend\(\)](#) (in module [pylablib.core.devio.comm\\_backend](#)), 159  
[new\\_frame\(\)](#) ([pylablib.devices.Thorlabs.TLCamera.ThorlabsTLCameraMethod](#)), 771  
[new\\_messages\\_number\(\)](#) ([pylablib.core.thread.controller.QTaskThreadMethod](#)), 309  
[new\\_messages\\_number\(\)](#) ([pylablib.core.thread.controller.QThreadControllerMethod](#)), 294  
[new\\_overflow\(\)](#) ([pylablib.devices.Andor.AndorSDK3.AndorSDK3Camera.BufferManagerMethod](#)), 482  
[NewportBackendError](#), 627  
[NewportError](#), 627  
[next\(\)](#) ([pylablib.core.utils.general.AccessIteratorMethod](#)), 380  
[next\(\)](#) ([pylablib.core.utils.numerical.infinite\\_list.counter](#)

`method`), 392  
NIDAQ (class in `pylablib.devices.NI.daq`), 621  
NIDAQ.NoParameterCaller (class in `pylablib.devices.NI.daq`), 626  
NIDAQmxError, 620  
NIError, 620  
`no_stopping()` (`pylablib.core.thread.controller.QTaskThread` method), 309  
`no_stopping()` (`pylablib.core.thread.controller.QThreadController` method), 295  
NoControllerThreadError, 319  
`nodes()` (`pylablib.core.utils.dictionary.Dictionary` method), 330  
`nodes()` (`pylablib.core.utils.dictionary.DictionaryPointer` method), 341  
`nodes()` (`pylablib.core.utils.dictionary.FilterTree` method), 358  
`nodes()` (`pylablib.core.utils.dictionary.PrefixTree` method), 350  
NoMessageThreadError, 320  
NoParameterError, 265  
`normalize_channel_name()` (`pylablib.devices.Tektronix.base.DPO2000` method), 765  
`normalize_channel_name()` (`pylablib.devices.Tektronix.base.ITektronixScope` method), 749  
`normalize_channel_name()` (`pylablib.devices.Tektronix.base.TDS2000` method), 758  
`normalize_fourier_transform()` (in module `pylablib.core.dataproc.fourier`), 118  
`normalize_path()` (in module `pylablib.core.utils.dictionary`), 325  
`normalize_path()` (in module `pylablib.core.utils.files`), 362  
`normalize_path_entry()` (in module `pylablib.core.utils.dictionary`), 325  
`notes` (`pylablib.devices.Thorlabs.kinesis.TDeviceInfo` attribute), 781  
`notify()` (`pylablib.core.thread.callsync.QCallResultSynchronizer` method), 282  
`notify()` (`pylablib.core.thread.callsync.QDirectResultSynchronizer` method), 283  
`notify()` (`pylablib.core.thread.callsync.QDummyResultSynchronizer` method), 282  
`notify()` (`pylablib.core.thread.notifier.ISkippableNotifier` method), 317  
`notify()` (`pylablib.core.thread.synchronizing.QMultiThreadNotifier` method), 319  
`notify()` (`pylablib.core.thread.synchronizing.QThreadNotifier` method), 318  
`notify()` (`pylablib.core.utils.observer_pool.ObserverPool` method), 393  
`notify()` (`pylablib.devices.PCO.SC2.PCOS2Camera.ScheduleLooper` method), 646  
`notify_exec_point()` (`pylablib.core.thread.controller.QTaskThread` method), 309  
`notify_exec_point()` (`pylablib.core.thread.controller.QThreadController` method), 299  
`notifying_state()` (`pylablib.core.thread.callsync.QCallResultSynchronizer` method), 282  
`notifying_state()` (`pylablib.core.thread.callsync.QDirectResultSynchronizer` method), 283  
`notifying_state()` (`pylablib.core.thread.notifier.ISkippableNotifier` method), 317  
`notifying_state()` (`pylablib.core.thread.synchronizing.QThreadNotifier` method), 318  
`nrois` (`pylablib.devices.PrincetonInstruments.picam.TROIConstraints` attribute), 709  
`num_bk` (`pylablib.devices.Thorlabs.kinesis.TKCubeTrigPosParams` attribute), 786  
`num_fw` (`pylablib.devices.Thorlabs.kinesis.TKCubeTrigPosParams` attribute), 786  
NumberLimit (class in `pylablib.core.gui.limiter`), 262  
NumEdit (class in `pylablib.core.gui.widgets.edit`), 234  
NumLabel (class in `pylablib.core.gui.widgets.label`), 237  
NumpyIndex (class in `pylablib.core.utils.indexing`), 382

## O

`obj_prop()` (in module `pylablib.core.utils.functions`), 374  
ObserverPool (class in `pylablib.core.utils.observer_pool`), 392  
ObserverPool.Observer (class in `pylablib.core.utils.observer_pool`), 392  
`obtain()` (in module `pylablib.core.utils.rpyc_utils`), 395  
`obtain()` (`pylablib.core.utils.rpyc_utils.DeviceService` method), 395  
`obtain()` (`pylablib.core.utils.rpyc_utils.SocketTunnelService` method), 394  
`off_thresh` (`pylablib.devices.Pfeiffer.base.TTPG260GaugeControlSettings` attribute), 654  
`offset` (`pylablib.core.utils.ipc.TShmemVarDesc` attribute), 384  
`offset` (`pylablib.devices.Thorlabs.TLCamera.TFrameInfo` attribute), 770



offset\_distance (pylablib.devices.Thorlabs.kinesis.THomeParams attribute), 785  
 on\_connect () (pylablib.core.utils.rpyc\_utils.DeviceService method), 394  
 on\_connect () (pylablib.core.utils.rpyc\_utils.SocketTunnelService method), 394  
 on\_disconnect () (pylablib.core.utils.rpyc\_utils.DeviceService method), 394  
 on\_disconnect () (pylablib.core.utils.rpyc\_utils.SocketTunnelService method), 394  
 on\_finish () (pylablib.core.thread.controller.QTaskThread method), 304  
 on\_finish () (pylablib.core.thread.controller.QThreadController method), 295  
 on\_overflow () (pylablib.devices.Andor.AndorSDK3.AndorSDK3Camera.BufferManager method), 482  
 on\_start () (pylablib.core.thread.controller.QTaskThread method), 304  
 on\_start () (pylablib.core.thread.controller.QThreadController method), 295  
 on\_thresh (pylablib.devices.Pfeiffer.base.TTPG260GaugeControl attribute), 654  
 opcode (pylablib.devices.Attocube.anc350.ANC350.Telegram attribute), 506  
 open () (pylablib.core.devio.comm\_backend.FT232DeviceBackend method), 150  
 open () (pylablib.core.devio.comm\_backend.ICommBackendWrapper method), 161  
 open () (pylablib.core.devio.comm\_backend.IDeviceCommBackend method), 143  
 open () (pylablib.core.devio.comm\_backend.NetworkDeviceBackend method), 153  
 open () (pylablib.core.devio.comm\_backend.PyUSBDeviceBackend method), 155  
 open () (pylablib.core.devio.comm\_backend.RecordedDeviceBackend method), 158  
 open () (pylablib.core.devio.comm\_backend.SerialDeviceBackend method), 148  
 open () (pylablib.core.devio.comm\_backend.VisaDeviceBackend method), 145  
 open () (pylablib.core.devio.interface.IDevice method), 162  
 open () (pylablib.core.devio.SCPI.SCPIDevice method), 141  
 open () (pylablib.core.fileio.location.FolderFileSystemDataLocation method), 188  
 open () (pylablib.core.fileio.location.IDataLocation method), 184  
 open () (pylablib.core.fileio.location.IFileSystemDataLocation method), 185  
 open () (pylablib.core.fileio.location.LocationFile method), 184  
 open () (pylablib.core.fileio.location.OpenedFileLocation method), 185  
 open () (pylablib.core.fileio.location.PrefixedFileSystemDataLocation method), 187  
 open () (pylablib.core.fileio.location.SingleFileSystemDataLocation method), 186  
 open () (pylablib.devices.AlliedVision.Bonito.BonitoIMAQCamera method), 460  
 open () (pylablib.devices.AlliedVision.Bonito.IBonitoCamera method), 451  
 open () (pylablib.devices.Andor.AndorSDK2.AndorSDK2Camera method), 467  
 open () (pylablib.devices.Andor.AndorSDK2.LibraryController method), 465  
 open () (pylablib.devices.Andor.AndorSDK3.AndorSDK3Camera method), 479  
 open () (pylablib.devices.Andor.AndorSDK3.LibraryController method), 476  
 open () (pylablib.devices.Arcus.performax.GenericPerformaxStage method), 488  
 open () (pylablib.devices.Arcus.performax.Performax2EXStage method), 496  
 open () (pylablib.devices.Arcus.performax.Performax4EXStage method), 493  
 open () (pylablib.devices.Arcus.performax.PerformaxDMXJSStage method), 499  
 open () (pylablib.devices.Arduino.base.IArduinoDevice method), 502  
 open () (pylablib.devices.Attocube.anc300.ANC300 method), 503  
 open () (pylablib.devices.Attocube.anc350.ANC350 method), 510  
 open () (pylablib.devices.AWG.generic.GenericAWG method), 407  
 open () (pylablib.devices.AWG.specific.Agilent33220A method), 417  
 open () (pylablib.devices.AWG.specific.Agilent33500 method), 411  
 open () (pylablib.devices.AWG.specific.InstekAFG2000 method), 429  
 open () (pylablib.devices.AWG.specific.InstekAFG2225 method), 423  
 open () (pylablib.devices.AWG.specific.RigolDG1000 method), 447  
 open () (pylablib.devices.AWG.specific.RSInstekAFG21000 method), 435  
 open () (pylablib.devices.AWG.specific.TektronixAFG1000 method), 441  
 open () (pylablib.devices.Basler.pylon.BaslerPylonCamera method), 514  
 open () (pylablib.devices.Basler.pylon.LibraryController method), 511

`open()` (`pylablib.devices.BitFlow.BitFlow.BitFlowCamera``open()` (`pylablib.devices.LighthousePhotonics.base.SproutG`  
`method`), 530 `method`), 613  
`open()` (`pylablib.devices.BitFlow.BitFlow.BitFlowFrameGrabber``open()` (`pylablib.devices.M2.solstis.Solstis`  
`method`), 521 `method`), 614  
`open()` (`pylablib.devices.Conrad.base.RelayBoard``open()` (`pylablib.devices.Newport.picomotor.Picomotor8742`  
`method`), 533 `method`), 631  
`open()` (`pylablib.devices.Cryocon.base.Cryocon1x``open()` (`pylablib.devices.NI.daq.NIDAQ`  
`method`), 537 `method`), 621  
`open()` (`pylablib.devices.Cryomagnetics.base.LM500``open()` (`pylablib.devices.Ophir.base.OphirDevice`  
`method`), 542 `method`), 640  
`open()` (`pylablib.devices.Cryomagnetics.base.LM510``open()` (`pylablib.devices.Ophir.base.VegaPowerMeter`  
`method`), 546 `method`), 643  
`open()` (`pylablib.devices.DCAM.DCAM.DCAMCamera``open()` (`pylablib.devices.OZOptics.base.DD100`  
`method`), 550 `method`), 636  
`open()` (`pylablib.devices.DCAM.DCAM.LibraryController``open()` (`pylablib.devices.OZOptics.base.EPC04`  
`method`), 548 `method`), 638  
`open()` (`pylablib.devices.HighFinesse.wlm.WLM``open()` (`pylablib.devices.OZOptics.base.OZOpticsDevice`  
`method`), 557 `method`), 633  
`open()` (`pylablib.devices.IMAQ.IMAQ.IMAQCamera``open()` (`pylablib.devices.OZOptics.base.TF100`  
`method`), 572 `method`), 634  
`open()` (`pylablib.devices.IMAQ.IMAQ.IMAQFrameGrabber``open()` (`pylablib.devices.PCO.SC2.PCOSC2Camera`  
`method`), 561 `method`), 645  
`open()` (`pylablib.devices.IMAQdx.IMAQdx.EthernetIMAQdxCamera``open()` (`pylablib.devices.Pfeiffer.base.DPG202`  
`method`), 587 `method`), 659  
`open()` (`pylablib.devices.IMAQdx.IMAQdx.IMAQdxCamera``open()` (`pylablib.devices.Pfeiffer.base.TPG260`  
`method`), 578 `method`), 657  
`open()` (`pylablib.devices.interface.camera.IAttributeCamera``open()` (`pylablib.devices.Photometrics.pvcam.LibraryController`  
`method`), 839 `method`), 659  
`open()` (`pylablib.devices.interface.camera.IBinROICamera``open()` (`pylablib.devices.Photometrics.pvcam.PvcamCamera`  
`method`), 858 `method`), 662  
`open()` (`pylablib.devices.interface.camera.ICamera``open()` (`pylablib.devices.PhotonFocus.PhotonFocus.IPhotonFocusCamera`  
`method`), 835 `method`), 672  
`open()` (`pylablib.devices.interface.camera.IExposureCamera``open()` (`pylablib.devices.PhotonFocus.PhotonFocus.LibraryController`  
`method`), 849 `method`), 669  
`open()` (`pylablib.devices.interface.camera.IGrabberAttributeCamera``open()` (`pylablib.devices.PhotonFocus.PhotonFocus.PhotonFocusBitFlow`  
`method`), 844 `method`), 694  
`open()` (`pylablib.devices.interface.camera.IROICamera``open()` (`pylablib.devices.PhotonFocus.PhotonFocus.PhotonFocusIMAQC`  
`method`), 854 `method`), 678  
`open()` (`pylablib.devices.interface.stage.IMultiaxisStage``open()` (`pylablib.devices.PhotonFocus.PhotonFocus.PhotonFocusSiSoCan`  
`method`), 864 `method`), 687  
`open()` (`pylablib.devices.interface.stage.IStage``open()` (`pylablib.devices.PhysikInstrumente.base.GenericPIController`  
`method`), 862 `method`), 703  
`open()` (`pylablib.devices.KJL.base.KJL300`  
`method`), 591 `open()` (`pylablib.devices.PhysikInstrumente.base.PIE516`  
`method`), 707  
`open()` (`pylablib.devices.Lakeshore.base.Lakeshore218``open()` (`pylablib.devices.PrincetonInstruments.picam.LibraryController`  
`method`), 596 `method`), 708  
`open()` (`pylablib.devices.Lakeshore.base.Lakeshore370``open()` (`pylablib.devices.PrincetonInstruments.picam.PicamCamera`  
`method`), 601 `method`), 712  
`open()` (`pylablib.devices.LaserQuantum.base.Finesse``open()` (`pylablib.devices.SiliconSoftware.fgrab.SiliconSoftwareCamera`  
`method`), 606 `method`), 730  
`open()` (`pylablib.devices.Leybold.base.GenericITR``open()` (`pylablib.devices.SiliconSoftware.fgrab.SiliconSoftwareFrameGra`  
`method`), 608 `method`), 721  
`open()` (`pylablib.devices.Leybold.base.ITR90`  
`method`), 610 `open()` (`pylablib.devices.Sirah.Matisse.SirahMatisse`  
`method`), 741  
`open()` (`pylablib.devices.SmarAct.scu3d.LibraryController`  
`method`), 610 `open()` (`pylablib.devices.SmarAct.scu3d.LibraryController`  
`method`), 610

- method*), 745
  - `open()` (*pylablib.devices.SmarAct.scu3d.SCU3D* *method*), 746
  - `open()` (*pylablib.devices.Tektronix.base.DPO2000* *method*), 765
  - `open()` (*pylablib.devices.Tektronix.base.ITektronixScope* *method*), 753
  - `open()` (*pylablib.devices.Tektronix.base.TDS2000* *method*), 758
  - `open()` (*pylablib.devices.Thorlabs.elliptec.ElliptecMotor* *method*), 781
  - `open()` (*pylablib.devices.Thorlabs.kinesis.BasicKinesisDevice* *method*), 784
  - `open()` (*pylablib.devices.Thorlabs.kinesis.KinesisDevice* *method*), 789
  - `open()` (*pylablib.devices.Thorlabs.kinesis.KinesisMotor* *method*), 800
  - `open()` (*pylablib.devices.Thorlabs.kinesis.KinesisPiezoMotor* *method*), 804
  - `open()` (*pylablib.devices.Thorlabs.kinesis.KinesisQuadDetector* *method*), 809
  - `open()` (*pylablib.devices.Thorlabs.kinesis.MFF* *method*), 793
  - `open()` (*pylablib.devices.Thorlabs.serial.FW* *method*), 814
  - `open()` (*pylablib.devices.Thorlabs.serial.FWv1* *method*), 818
  - `open()` (*pylablib.devices.Thorlabs.serial.MDT69xA* *method*), 821
  - `open()` (*pylablib.devices.Thorlabs.serial.ThorlabsSerialInterface* *method*), 809
  - `open()` (*pylablib.devices.Thorlabs.TLCamera.LibraryController* *method*), 769
  - `open()` (*pylablib.devices.Thorlabs.TLCamera.ThorlabsTLCamera* *method*), 770
  - `open()` (*pylablib.devices.Toptica.ibeam.TopticalBeam* *method*), 824
  - `open()` (*pylablib.devices.Trinamic.base.TMCM1110* *method*), 826
  - `open()` (*pylablib.devices.uc480.uc480.UC480Camera* *method*), 867
  - `open()` (*pylablib.devices.utils.load\_lib.LibraryController* *method*), 874
  - `open_loop_out` (*pylablib.devices.Thorlabs.kinesis.TQuadDetectorOutputParams* *attribute*), 806
  - `open_result` (*pylablib.devices.utils.load\_lib.TLibraryOpenResult* *attribute*), 874
  - `opened` (*pylablib.core.fileio.location.LocationFile* *attribute*), 184
  - `OpenedFileLocation` (*class in pylablib.core.fileio.location*), 185
  - `OphirBackendError`, 638
  - `OphirDevice` (*class in pylablib.devices.Ophir.base*), 639
  - `OphirDevice.NoParameterCaller` (*class in pylablib.devices.Ophir.base*), 639
  - `OphirError`, 638
  - `opid` (*pylablib.devices.utils.load\_lib.TLibraryOpenResult* *attribute*), 874
  - `order_to_pos()` (*in module pylablib.core.gui.formatter*), 261
  - `overflows` (*pylablib.devices.Andor.AndorSDK3.TMissedFramesStatus* *attribute*), 479
  - `overruns` (*pylablib.devices.PCO.SC2.TInternalBufferStatus* *attribute*), 645
  - `oversamp` (*pylablib.devices.Sirah.Matisse.TPiezoetDriveParameters* *attribute*), 734
  - `oversamp` (*pylablib.devices.Sirah.Matisse.TRefcellWaveformParameters* *attribute*), 735
  - `OZOpticsBackendError`, 631
  - `OZOpticsDevice` (*class in pylablib.devices.OZOptics.base*), 632
  - `OZOpticsDevice.NoParameterCaller` (*class in pylablib.devices.OZOptics.base*), 632
  - `OZOpticsError`, 631
- ## P
- `P` (*pylablib.devices.Sirah.Matisse.TPiezoetFeedbackParameters* *attribute*), 734
  - `P` (*pylablib.devices.Sirah.Matisse.TSlowpiezoCtlParameters* *attribute*), 735
  - `P` (*pylablib.devices.Sirah.Matisse.TThinectlParameters* *attribute*), 734
  - `p` (*pylablib.devices.Thorlabs.kinesis.TQuadDetectorPIDParams* *attribute*), 805
  - `pack_int()` (*in module pylablib.core.utils.strpack*), 401
  - `pack_uint()` (*in module pylablib.core.utils.strpack*), 401
  - `pad_borders()` (*pylablib.core.gui.widgets.param\_table.ParamTable* *method*), 243
  - `pad_borders()` (*pylablib.core.gui.widgets.param\_table.StatusTable* *method*), 259
  - `pad_trace()` (*in module pylablib.core.dataproc.utils*), 137
  - `params` (*pylablib.devices.Leybold.base.TDeviceInfo* *attribute*), 607
  - `param1` (*pylablib.devices.Thorlabs.kinesis.BasicKinesisDevice.CommShort* *attribute*), 782
  - `param1` (*pylablib.devices.Thorlabs.kinesis.KinesisDevice.CommShort* *attribute*), 788
  - `param1` (*pylablib.devices.Thorlabs.kinesis.KinesisMotor.CommShort* *attribute*), 798
  - `param1` (*pylablib.devices.Thorlabs.kinesis.KinesisPiezoMotor.CommShort* *attribute*), 803

param1 (pylablib.devices.Thorlabs.kinesis.KinesisQuadDetector.CommShortdata (py-  
 attribute), 807  
 param1 (pylablib.devices.Thorlabs.kinesis.MFF.CommShort static method), 542  
 attribute), 792  
 param2 (pylablib.devices.Thorlabs.kinesis.BasicKinesisDevice.CommShortdata (py-  
 attribute), 782  
 param2 (pylablib.devices.Thorlabs.kinesis.KinesisDevice.CommShortdata (py-  
 attribute), 788  
 param2 (pylablib.devices.Thorlabs.kinesis.KinesisMotor.CommShortstatic method), 596  
 attribute), 798  
 param2 (pylablib.devices.Thorlabs.kinesis.KinesisPiezoMotor.CommShortdata (py-  
 attribute), 803  
 param2 (pylablib.devices.Thorlabs.kinesis.KinesisQuadDetector.CommShortdata (py-  
 attribute), 807  
 param2 (pylablib.devices.Thorlabs.kinesis.MFF.CommShort static method), 741  
 attribute), 792  
 parameter\_range\_error() (in module py- lablib.devices.Tektronix.base.DPO2000 static  
 lablib.core.utils.funcargparse), 368 method), 765  
 parameter\_value\_error() (in module py- parse\_array\_data() (py-  
 lablib.core.utils.funcargparse), 368 lablib.devices.Tektronix.base.ITektronixScope  
 ParamTable (class in py- static method), 753  
 lablib.core.gui.widgets.param\_table), 242 parse\_array\_data() (py-  
 ParamTable.ParamRow (class in py- lablib.devices.Tektronix.base.TDS2000 static  
 lablib.core.gui.widgets.param\_table), 243 method), 758  
 parse\_array\_data() (py- parse\_array\_data() (py-  
 lablib.core.devio.SCPIDevice static lablib.devices.Thorlabs.serial.FW static  
 method), 141 method), 814  
 parse\_array\_data() (py- parse\_array\_data() (py-  
 lablib.devices.AWG.generic.GenericAWG static method), 407 lablib.devices.Thorlabs.serial.FWv1 static  
 static method), 407 method), 818  
 parse\_array\_data() (py- parse\_array\_data() (py-  
 lablib.devices.AWG.specific.Agilent33220A static method), 417 lablib.devices.Thorlabs.serial.MDT69xA  
 static method), 417 static method), 821  
 parse\_array\_data() (py- parse\_array\_data() (py-  
 lablib.devices.AWG.specific.Agilent33500 static method), 411 lablib.devices.Thorlabs.serial.ThorlabsSerialInterface  
 static method), 411 static method), 811  
 parse\_array\_data() (py- parse\_dict\_line() (in module py-  
 lablib.devices.AWG.specific.InstekAFG2000 lablib.core.fileio.loadfile\_utils), 182  
 static method), 429 parse\_float() (in module py-  
 parse\_array\_data() (py- lablib.core.gui.formatter), 260  
 lablib.devices.AWG.specific.InstekAFG2225 static method), 423 parse\_metainfo\_v1() (in module py-  
 static method), 423 lablib.devices.Photometrics.pvcam), 669  
 parse\_array\_data() (py- parse\_metainfo\_v3() (in module py-  
 lablib.devices.AWG.specific.RigolDG1000 static method), 447 lablib.devices.Photometrics.pvcam), 669  
 static method), 447 parse\_stored\_table\_data() (in module py-  
 parse\_array\_data() (py- lablib.core.fileio.dict\_entry), 170  
 lablib.devices.AWG.specific.RSInstekAFG21000 static method), 435 part (pylablib.devices.Photometrics.pvcam.TDeviceInfo  
 static method), 435 attribute), 661  
 parse\_array\_data() (py- partition\_list() (in module py-  
 lablib.devices.AWG.specific.TektronixAFG1000 static method), 441 lablib.core.utils.general), 376  
 static method), 441 pass\_result (pylablib.core.thread.callsync.QScheduledCall.Callback  
 parse\_array\_data() (py- attribute), 283  
 lablib.devices.Cryocon.base.Cryocon1x static passed() (pylablib.core.utils.general.Countdown  
 method), 537 method), 379



passed() (pylablib.core.utils.general.Timer method), 379  
 path (pylablib.devices.SiliconSoftware.fgrab.TFullAppletInfo attribute), 719  
 paths() (pylablib.core.utils.dictionary.Dictionary method), 330  
 paths() (pylablib.core.utils.dictionary.DictionaryPointer method), 341  
 paths() (pylablib.core.utils.dictionary.FilterTree method), 358  
 paths() (pylablib.core.utils.dictionary.PrefixTree method), 350  
 paths\_equal() (in module pylablib.core.utils.files), 362  
 pause() (pylablib.core.thread.controller.QTaskThread.Job method), 303  
 pausing\_acquisition() (py-lablib.devices.AlliedVision.Bonito.BonitoIMAQCamera method), 460  
 pausing\_acquisition() (py-lablib.devices.AlliedVision.Bonito.IBonitoCamera method), 454  
 pausing\_acquisition() (py-lablib.devices.Andor.AndorSDK2.AndorSDK2Camera method), 474  
 pausing\_acquisition() (py-lablib.devices.Andor.AndorSDK3.AndorSDK3Camera method), 485  
 pausing\_acquisition() (py-lablib.devices.Basler.pylon.BaslerPylonCamera method), 518  
 pausing\_acquisition() (py-lablib.devices.BitFlow.BitFlow.BitFlowCamera method), 530  
 pausing\_acquisition() (py-lablib.devices.BitFlow.BitFlow.BitFlowFrameGrabber method), 525  
 pausing\_acquisition() (py-lablib.devices.DCAM.DCAM.DCAMCamera method), 555  
 pausing\_acquisition() (py-lablib.devices.IMAQ.IMAQ.IMAQCamera method), 572  
 pausing\_acquisition() (py-lablib.devices.IMAQ.IMAQ.IMAQFrameGrabber method), 567  
 pausing\_acquisition() (py-lablib.devices.IMAQdx.IMAQdx.EthernetIMAQdxCamera method), 587  
 pausing\_acquisition() (py-lablib.devices.IMAQdx.IMAQdx.IMAQdxCamera method), 582  
 pausing\_acquisition() (py-lablib.devices.interface.camera.IAttributeCamera method), 840  
 pausing\_acquisition() (py-lablib.devices.interface.camera.IBinROICamera method), 858  
 pausing\_acquisition() (py-lablib.devices.interface.camera.ICamera method), 832  
 pausing\_acquisition() (py-lablib.devices.interface.camera.IExposureCamera method), 849  
 pausing\_acquisition() (py-lablib.devices.interface.camera.IGrabberAttributeCamera method), 844  
 pausing\_acquisition() (py-lablib.devices.interface.camera.IROICamera method), 854  
 pausing\_acquisition() (py-lablib.devices.PCO.SC2.PCOSC2Camera method), 651  
 pausing\_acquisition() (py-lablib.devices.Photometrics.pvcam.PvcamCamera method), 667  
 pausing\_acquisition() (py-lablib.devices.PhotonFocus.PhotonFocus.IPhotonFocusCamera method), 676  
 pausing\_acquisition() (py-lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusBitFlowCamera method), 698  
 pausing\_acquisition() (py-lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusIMAQCamera method), 683  
 pausing\_acquisition() (py-lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusSiSoCamera method), 691  
 pausing\_acquisition() (py-lablib.devices.PrincetonInstruments.picam.PicamCamera method), 716  
 pausing\_acquisition() (py-lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareCamera method), 731  
 pausing\_acquisition() (py-lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareFrameGrabber method), 725  
 pausing\_acquisition() (py-lablib.devices.Thorlabs.TLCamera.ThorlabsTLCamera method), 774  
 pausing\_acquisition() (py-lablib.devices.uc480.uc480.UC480Camera method), 871  
 PCOSC2Camera (class in pylablib.devices.PCO.SC2), 645  
 PCOSC2Camera.BufferManager (class in pylablib.devices.PCO.SC2), 646  
 PCOSC2Camera.NoParameterCaller (class in py-

- pylablib.devices.PCO.SC2*), 649
- PCOSC2Camera.ScheduleLooper* (class in *pylablib.devices.PCO.SC2*), 646
- Peak* (class in *pylablib.core.dataproc.feature*), 108
- peaks\_sum\_func()* (in module *pylablib.core.dataproc.feature*), 109
- Performax2EXStage* (class in *pylablib.devices.Arcus.performax*), 493
- Performax2EXStage.NoParameterCaller* (class in *pylablib.devices.Arcus.performax*), 493
- Performax4EXStage* (class in *pylablib.devices.Arcus.performax*), 489
- Performax4EXStage.NoParameterCaller* (class in *pylablib.devices.Arcus.performax*), 492
- PerformaxDMXJSASStage* (class in *pylablib.devices.Arcus.performax*), 497
- PerformaxDMXJSASStage.NoParameterCaller* (class in *pylablib.devices.Arcus.performax*), 498
- period* (*pylablib.core.gui.widgets.container.TTimer* attribute), 199
- PFCamAttribute* (class in *pylablib.devices.PhotonFocus.PhotonFocus*), 670
- PfeifferBackendError*, 653
- PfeifferError*, 653
- phase* (*pylablib.devices.Sirah.Matisse.TPiezoetFeedbackParameters* attribute), 734
- phase* (*pylablib.devices.Sirah.Matisse.TPiezoetFeedforwardParameters* attribute), 735
- PhotonFocusBitFlowCamera* (class in *pylablib.devices.PhotonFocus.PhotonFocus*), 694
- PhotonFocusBitFlowCamera.BufferManager* (class in *pylablib.devices.PhotonFocus.PhotonFocus*), 695
- PhotonFocusBitFlowCamera.NoParameterCaller* (class in *pylablib.devices.PhotonFocus.PhotonFocus*), 695
- PhotonFocusIMAQCamera* (class in *pylablib.devices.PhotonFocus.PhotonFocus*), 678
- PhotonFocusIMAQCamera.NoParameterCaller* (class in *pylablib.devices.PhotonFocus.PhotonFocus*), 678
- PhotonFocusSiSoCamera* (class in *pylablib.devices.PhotonFocus.PhotonFocus*), 686
- PhotonFocusSiSoCamera.NoParameterCaller* (class in *pylablib.devices.PhotonFocus.PhotonFocus*), 687
- PhysikInstrumenteBackendError*, 702
- PhysikInstrumenteError*, 702
- PicamAttribute* (class in *pylablib.devices.PrincetonInstruments.picam*), 709
- PicamCamera* (class in *pylablib.devices.PrincetonInstruments.picam*), 712
- PicamCamera.NoParameterCaller* (class in *pylablib.devices.PrincetonInstruments.picam*), 714
- Picomotor8742* (class in *pylablib.devices.Newport.picomotor*), 628
- Picomotor8742.NoParameterCaller* (class in *pylablib.devices.Newport.picomotor*), 630
- PIE516* (class in *pylablib.devices.PhyisikInstrumente.base*), 705
- PIE516.NoParameterCaller* (class in *pylablib.devices.PhyisikInstrumente.base*), 706
- pip\_install()* (in module *pylablib.core.utils.module*), 386
- PipeIPCChannel* (class in *pylablib.core.utils.ipc*), 383
- pixelclock* (*pylablib.devices.Thorlabs.TLCamera.TFrameInfo* attribute), 770
- pixeltype* (*pylablib.devices.Andor.AndorSDK3.TFrameInfo* attribute), 479
- pixeltype* (*pylablib.devices.DCAM.DCAM.TFrameInfo* attribute), 550
- pixeltype* (*pylablib.devices.Thorlabs.TLCamera.TFrameInfo* attribute), 770
- place\_widget\_at\_location()* (in module *pylablib.core.gui.utils*), 264
- platform* (*pylablib.devices.SiliconSoftware.fgrab.TFullAppletInfo* attribute), 719
- points* (*pylablib.devices.Lakeshore.base.TLakeshore218FilterSettings* attribute), 593
- poke()* (*pylablib.core.thread.controller.QTaskThread* method), 309
- poke()* (*pylablib.core.thread.controller.QThreadController* method), 299
- polynomial()* (in module *pylablib.core.utils.numerical*), 392
- pop()* (*pylablib.core.utils.dictionary.Dictionary* method), 328
- pop()* (*pylablib.core.utils.dictionary.DictionaryPointer* method), 341
- pop()* (*pylablib.core.utils.dictionary.FilterTree* method), 358
- pop()* (*pylablib.core.utils.dictionary.PrefixTree* method), 350



`props` (*pylablib.devices.Basler.pylon.TDeviceInfo* attribute), 514

`pstep` (*pylablib.devices.interface.camera.TAxisROILimit* attribute), 851

`psu` (*pylablib.devices.LaserQuantum.base.TTemperatures* attribute), 604

`psu` (*pylablib.devices.LaserQuantum.base.TWorkHours* attribute), 604

`pulse` (*pylablib.devices.Thorlabs.elliptec.TDeviceInfo* attribute), 777

`PushButtonValueHandler` (class in *pylablib.core.gui.value\_handling*), 273

`pval()` (in module *pylablib.core.devio.interface*), 168

`PvcamAttribute` (class in *pylablib.devices.Photometrics.pvcam*), 659

`PvcamCamera` (class in *pylablib.devices.Photometrics.pvcam*), 662

`PvcamCamera.NoParameterCaller` (class in *pylablib.devices.Photometrics.pvcam*), 665

*pylablib* (module), 875

*pylablib.core* (module), 402

*pylablib.core.dataproc* (module), 138

*pylablib.core.dataproc.callable* (module), 104

*pylablib.core.dataproc.feature* (module), 108

*pylablib.core.dataproc.filters* (module), 110

*pylablib.core.dataproc.fitting* (module), 114

*pylablib.core.dataproc.fourier* (module), 117

*pylablib.core.dataproc.iir\_transform* (module), 120

*pylablib.core.dataproc.image* (module), 121

*pylablib.core.dataproc.interpolate* (module), 122

*pylablib.core.dataproc.specfunc* (module), 124

*pylablib.core.dataproc.table\_wrap* (module), 125

*pylablib.core.dataproc.transform* (module), 134

*pylablib.core.dataproc.utils* (module), 135

*pylablib.core.devio* (module), 168

*pylablib.core.devio.backend\_logger* (module), 141

*pylablib.core.devio.base* (module), 142

*pylablib.core.devio.comm\_backend* (module), 142

*pylablib.core.devio.data\_format* (module), 161

*pylablib.core.devio.interface* (module), 162

*pylablib.core.devio.SCPi* (module), 138

*pylablib.core.fileio* (module), 197

*pylablib.core.fileio.datafile* (module), 168

*pylablib.core.fileio.dict\_entry* (module), 169

*pylablib.core.fileio.loadfile* (module), 176

*pylablib.core.fileio.loadfile\_utils* (module), 182

*pylablib.core.fileio.location* (module), 183

*pylablib.core.fileio.parse\_csv* (module), 188

*pylablib.core.fileio.savefile* (module), 190

*pylablib.core.fileio.table\_stream* (module), 196

*pylablib.core.gui* (module), 281

*pylablib.core.gui.formatter* (module), 260

*pylablib.core.gui.limiter* (module), 262

*pylablib.core.gui.utils* (module), 263

*pylablib.core.gui.value\_handling* (module), 265

*pylablib.core.gui.widgets* (module), 260

*pylablib.core.gui.widgets.button* (module), 197

*pylablib.core.gui.widgets.combo\_box* (module), 198

*pylablib.core.gui.widgets.container* (module), 199

*pylablib.core.gui.widgets.edit* (module), 233

*pylablib.core.gui.widgets.label* (module), 236

*pylablib.core.gui.widgets.layout\_manager* (module), 238

*pylablib.core.gui.widgets.param\_table* (module), 242

*pylablib.core.thread* (module), 321

*pylablib.core.thread.callsync* (module), 281

*pylablib.core.thread.controller* (module), 292

*pylablib.core.thread.multicast\_pool* (module), 315

*pylablib.core.thread.notifier* (module), 316

*pylablib.core.thread.profile* (module), 317

*pylablib.core.thread.synchronizing* (module), 318

*pylablib.core.thread.threadprop* (module), 319

*pylablib.core.thread.utils* (module), 321

[pylablib.core.utils \(module\)](#), 402  
[pylablib.core.utils.array\\_utils \(module\)](#), 321  
[pylablib.core.utils.ctypes\\_wrap \(module\)](#), 321  
[pylablib.core.utils.dictionary \(module\)](#), 325  
[pylablib.core.utils.files \(module\)](#), 361  
[pylablib.core.utils.funcargparse \(module\)](#), 368  
[pylablib.core.utils.functions \(module\)](#), 369  
[pylablib.core.utils.general \(module\)](#), 374  
[pylablib.core.utils.indexing \(module\)](#), 381  
[pylablib.core.utils.ipc \(module\)](#), 383  
[pylablib.core.utils.library\\_parameters \(module\)](#), 385  
[pylablib.core.utils.module \(module\)](#), 385  
[pylablib.core.utils.nbtools \(module\)](#), 387  
[pylablib.core.utils.net \(module\)](#), 387  
[pylablib.core.utils.numerical \(module\)](#), 391  
[pylablib.core.utils.observer\\_pool \(module\)](#), 392  
[pylablib.core.utils.py3 \(module\)](#), 393  
[pylablib.core.utils.rpyc\\_utils \(module\)](#), 394  
[pylablib.core.utils.strdump \(module\)](#), 395  
[pylablib.core.utils.string \(module\)](#), 396  
[pylablib.core.utils.strpack \(module\)](#), 401  
[pylablib.core.utils.units \(module\)](#), 401  
[pylablib.devices \(module\)](#), 875  
[pylablib.devices.AlliedVision \(module\)](#), 465  
[pylablib.devices.AlliedVision.Bonito \(module\)](#), 450  
[pylablib.devices.Andor \(module\)](#), 487  
[pylablib.devices.Andor.AndorSDK2 \(module\)](#), 465  
[pylablib.devices.Andor.AndorSDK3 \(module\)](#), 476  
[pylablib.devices.Andor.atcore\\_features \(module\)](#), 486  
[pylablib.devices.Andor.base \(module\)](#), 486  
[pylablib.devices.Arcus \(module\)](#), 500  
[pylablib.devices.Arcus.base \(module\)](#), 487  
[pylablib.devices.Arcus.performax \(module\)](#), 488  
[pylablib.devices.Arduino \(module\)](#), 502  
[pylablib.devices.Arduino.base \(module\)](#), 500  
[pylablib.devices.Attocube \(module\)](#), 510  
[pylablib.devices.Attocube.anc300 \(module\)](#), 502  
[pylablib.devices.Attocube.anc350 \(module\)](#), 506  
[pylablib.devices.Attocube.base \(module\)](#), 510  
[pylablib.devices.AWG \(module\)](#), 450  
[pylablib.devices.AWG.generic \(module\)](#), 402  
[pylablib.devices.AWG.specific \(module\)](#), 408  
[pylablib.devices.Basler \(module\)](#), 520  
[pylablib.devices.Basler.pylon \(module\)](#), 510  
[pylablib.devices.BitFlow \(module\)](#), 533  
[pylablib.devices.BitFlow.BitFlow \(module\)](#), 520  
[pylablib.devices.Conrad \(module\)](#), 535  
[pylablib.devices.Conrad.base \(module\)](#), 533  
[pylablib.devices.Cryocon \(module\)](#), 539  
[pylablib.devices.Cryocon.base \(module\)](#), 535  
[pylablib.devices.Cryomagnetics \(module\)](#), 548  
[pylablib.devices.Cryomagnetics.base \(module\)](#), 539  
[pylablib.devices.DCAM \(module\)](#), 556  
[pylablib.devices.DCAM.DCAM \(module\)](#), 548  
[pylablib.devices.HighFinesse \(module\)](#), 561  
[pylablib.devices.HighFinesse.wlm \(module\)](#), 557  
[pylablib.devices.IMAQ \(module\)](#), 576  
[pylablib.devices.IMAQ.IMAQ \(module\)](#), 561  
[pylablib.devices.IMAQ.niimaq\\_attrtypes \(module\)](#), 576  
[pylablib.devices.IMAQdx \(module\)](#), 589  
[pylablib.devices.IMAQdx.IMAQdx \(module\)](#), 576  
[pylablib.devices.interface \(module\)](#), 864  
[pylablib.devices.interface.camera \(module\)](#), 830  
[pylablib.devices.interface.stage \(module\)](#), 862  
[pylablib.devices.KJL \(module\)](#), 592  
[pylablib.devices.KJL.base \(module\)](#), 589  
[pylablib.devices.Lakeshore \(module\)](#), 603  
[pylablib.devices.Lakeshore.base \(module\)](#), 592  
[pylablib.devices.LaserQuantum \(module\)](#), 606  
[pylablib.devices.LaserQuantum.base \(module\)](#), 603  
[pylablib.devices.Leybold \(module\)](#), 611  
[pylablib.devices.Leybold.base \(module\)](#), 606  
[pylablib.devices.LighthousePhotonics \(module\)](#), 613



[pylablib.devices.LighthousePhotonics.base \(module\)](#), 611  
[pylablib.devices.M2 \(module\)](#), 620  
[pylablib.devices.M2.solstis \(module\)](#), 614  
[pylablib.devices.Newport \(module\)](#), 631  
[pylablib.devices.Newport.base \(module\)](#), 627  
[pylablib.devices.Newport.picomotor \(module\)](#), 627  
[pylablib.devices.NI \(module\)](#), 627  
[pylablib.devices.NI.daq \(module\)](#), 620  
[pylablib.devices.Ophir \(module\)](#), 644  
[pylablib.devices.Ophir.base \(module\)](#), 638  
[pylablib.devices.OZOptics \(module\)](#), 638  
[pylablib.devices.OZOptics.base \(module\)](#), 631  
[pylablib.devices.PCO \(module\)](#), 653  
[pylablib.devices.PCO.SC2 \(module\)](#), 644  
[pylablib.devices.Pfeiffer \(module\)](#), 659  
[pylablib.devices.Pfeiffer.base \(module\)](#), 653  
[pylablib.devices.Photometrics \(module\)](#), 669  
[pylablib.devices.Photometrics.pvcam \(module\)](#), 659  
[pylablib.devices.PhotonFocus \(module\)](#), 702  
[pylablib.devices.PhotonFocus.PhotonFocus \(module\)](#), 669  
[pylablib.devices.PhysikInstrumente \(module\)](#), 708  
[pylablib.devices.PhysikInstrumente.base \(module\)](#), 702  
[pylablib.devices.PrincetonInstruments \(module\)](#), 718  
[pylablib.devices.PrincetonInstruments.picam \(module\)](#), 708  
[pylablib.devices.SiliconSoftware \(module\)](#), 733  
[pylablib.devices.SiliconSoftware.fgrab \(module\)](#), 718  
[pylablib.devices.Sirah \(module\)](#), 745  
[pylablib.devices.Sirah.Matisse \(module\)](#), 733  
[pylablib.devices.Sirah.tuner \(module\)](#), 743  
[pylablib.devices.SmarAct \(module\)](#), 748  
[pylablib.devices.SmarAct.scu3d \(module\)](#), 745  
[pylablib.devices.Tektronix \(module\)](#), 769  
[pylablib.devices.Tektronix.base \(module\)](#), 748  
[pylablib.devices.Thorlabs \(module\)](#), 823  
[pylablib.devices.Thorlabs.base \(module\)](#), 776  
[pylablib.devices.Thorlabs.elliptec \(module\)](#), 776  
[pylablib.devices.Thorlabs.kinesis \(module\)](#), 781  
[pylablib.devices.Thorlabs.misc \(module\)](#), 809  
[pylablib.devices.Thorlabs.serial \(module\)](#), 809  
[pylablib.devices.Thorlabs.TLCamera \(module\)](#), 769  
[pylablib.devices.Toptica \(module\)](#), 826  
[pylablib.devices.Toptica.base \(module\)](#), 823  
[pylablib.devices.Toptica.ibeam \(module\)](#), 823  
[pylablib.devices.Trinamic \(module\)](#), 830  
[pylablib.devices.Trinamic.base \(module\)](#), 826  
[pylablib.devices.uc480 \(module\)](#), 873  
[pylablib.devices.uc480.uc480 \(module\)](#), 864  
[pylablib.devices.utils \(module\)](#), 875  
[pylablib.devices.utils.load\\_lib \(module\)](#), 873  
[pylablib.widgets \(module\)](#), 875  
[PyUSBDeviceBackend \(class in pylablib.core.devio.comm\\_backend\)](#), 155

## Q

[QCallResultSynchronizer \(class in pylablib.core.thread.callsync\)](#), 281  
[QContainer \(class in pylablib.core.gui.widgets.container\)](#), 202  
[QDialogContainer \(class in pylablib.core.gui.widgets.container\)](#), 217  
[QDirectCallScheduler \(class in pylablib.core.thread.callsync\)](#), 285  
[QDirectResultSynchronizer \(class in pylablib.core.thread.callsync\)](#), 282  
[QDummyResultSynchronizer \(class in pylablib.core.thread.callsync\)](#), 282  
[QFrameContainer \(class in pylablib.core.gui.widgets.container\)](#), 213  
[QGroupBoxContainer \(class in pylablib.core.gui.widgets.container\)](#), 221  
[QLayoutManagedWidget \(class in pylablib.core.gui.widgets.layout\\_manager\)](#), 240  
[QLockNotifier \(class in pylablib.core.thread.synchronizing\)](#), 319  
[QMulticastThreadCallScheduler \(class in pylablib.core.thread.callsync\)](#), 291  
[QMultiQueueScheduler \(class in pylablib.core.thread.callsync\)](#), 290  
[QMultiThreadNotifier \(class in pylablib.core.thread.synchronizing\)](#), 318

[QQueueLengthLimitScheduler \(class in py-lablib.core.thread.callsync\)](#), 287  
[QQueueScheduler \(class in py-lablib.core.thread.callsync\)](#), 285  
[QQueueSizeLimitScheduler \(class in py-lablib.core.thread.callsync\)](#), 288  
[QScheduledCall \(class in py-lablib.core.thread.callsync\)](#), 283  
[QScheduledCall.Callback \(class in py-lablib.core.thread.callsync\)](#), 283  
[QScheduler \(class in pylablib.core.thread.callsync\)](#), 284  
[QScrollAreaContainer \(class in py-lablib.core.gui.widgets.container\)](#), 224  
[QScrollAreaContainer.QContainedWidget \(class in pylablib.core.gui.widgets.container\)](#), 225  
[QTabContainer \(class in py-lablib.core.gui.widgets.container\)](#), 231  
[QTaskThread \(class in py-lablib.core.thread.controller\)](#), 301  
[QTaskThread.CommandAccess \(class in py-lablib.core.thread.controller\)](#), 307  
[QTaskThread.Job \(class in py-lablib.core.thread.controller\)](#), 302  
[QTaskThread.TBatchJob \(class in py-lablib.core.thread.controller\)](#), 302  
[QTaskThread.TCommand \(class in py-lablib.core.thread.controller\)](#), 302  
[QThreadCallScheduler \(class in py-lablib.core.thread.callsync\)](#), 290  
[QThreadController \(class in py-lablib.core.thread.controller\)](#), 293  
[QThreadControllerThread \(class in py-lablib.core.thread.controller\)](#), 292  
[QThreadNotifier \(class in py-lablib.core.thread.synchronizing\)](#), 318  
[query \(\) \(pylablib.devices.Arcus.performax.GenericPerformanceStage method\)](#), 488  
[query \(\) \(pylablib.devices.Arcus.performax.Performax2EXStage method\)](#), 496  
[query \(\) \(pylablib.devices.Arcus.performax.Performax4EXStage method\)](#), 493  
[query \(\) \(pylablib.devices.Arcus.performax.PerformaxDMXISAStage method\)](#), 499  
[query \(\) \(pylablib.devices.Arduino.base.IArduinoDevice method\)](#), 501  
[query \(\) \(pylablib.devices.Attocube.anc300.ANC300 method\)](#), 503  
[query \(\) \(pylablib.devices.Conrad.base.RelayBoard method\)](#), 534  
[query \(\) \(pylablib.devices.KJL.base.KJL300 method\)](#), 590  
[query \(\) \(pylablib.devices.LaserQuantum.base.Finesse method\)](#), 604  
[query \(\) \(pylablib.devices.LighthousePhotonics.base.SproutG method\)](#), 612  
[query \(\) \(pylablib.devices.M2.solstis.Solstis method\)](#), 615  
[query \(\) \(pylablib.devices.Newport.picomotor.Picomotor8742 method\)](#), 628  
[query \(\) \(pylablib.devices.Ophir.base.OphirDevice method\)](#), 639  
[query \(\) \(pylablib.devices.Ophir.base.VegaPowerMeter method\)](#), 643  
[query \(\) \(pylablib.devices.OZOptics.base.DD100 method\)](#), 636  
[query \(\) \(pylablib.devices.OZOptics.base.EPC04 method\)](#), 636  
[query \(\) \(pylablib.devices.OZOptics.base.OZOpticsDevice method\)](#), 632  
[query \(\) \(pylablib.devices.OZOptics.base.TF100 method\)](#), 634  
[query \(\) \(pylablib.devices.Pfeiffer.base.DPG202 method\)](#), 657  
[query \(\) \(pylablib.devices.Pfeiffer.base.TPG260 method\)](#), 654  
[query \(\) \(pylablib.devices.PhysikInstrumente.base.GenericPIController method\)](#), 703  
[query \(\) \(pylablib.devices.PhysikInstrumente.base.PIE516 method\)](#), 707  
[query \(\) \(pylablib.devices.Thorlabs.elliptec.ElliptecMotor method\)](#), 778  
[query \(\) \(pylablib.devices.Thorlabs.kinesis.BasicKinesisDevice method\)](#), 783  
[query \(\) \(pylablib.devices.Thorlabs.kinesis.KinesisDevice method\)](#), 789  
[query \(\) \(pylablib.devices.Thorlabs.kinesis.KinesisMotor method\)](#), 800  
[query \(\) \(pylablib.devices.Thorlabs.kinesis.KinesisPiezoMotor method\)](#), 804  
[query \(\) \(pylablib.devices.Thorlabs.kinesis.KinesisQuadDetector method\)](#), 809  
[query \(\) \(pylablib.devices.Thorlabs.kinesis.MFF method\)](#), 793  
[query \(\) \(pylablib.devices.Toptica.ibeam.TopticaIBeam method\)](#), 824  
[query \(\) \(pylablib.devices.Trinamic.base.TMCM1110 method\)](#), 827  
[query\\_axis \(\) \(pylablib.devices.PhysikInstrumente.base.GenericPIController method\)](#), 703  
[query\\_axis \(\) \(pylablib.devices.PhysikInstrumente.base.PIE516 method\)](#), 707  
[query\\_camera\\_name \(\) \(in module py-lablib.devices.PhotonFocus.PhotonFocus\)](#), 670  
[queue \(\) \(pylablib.devices.Basler.pylon.BaslerPylonCamera.BufferManager method\)](#), 516

`queue_buffer()` (py-lablib.devices.Andor.AndorSDK3.AndorSDK3Camera.BufferManager method), 481  
`quit_sync()` (pylablib.core.thread.controller.QThreadControllerThread method), 292  
`QWidgetContainer` (class in py-lablib.core.gui.widgets.container), 209  
**R**  
`ramp_down` (pylablib.devices.Thorlabs.elliptec.TMotorInfo attribute), 777  
`ramp_up` (pylablib.devices.Thorlabs.elliptec.TMotorInfo attribute), 777  
`Range` (class in pylablib.core.dataproc.utils), 136  
`RangeParameterClass` (class in py-lablib.core.devio.interface), 164  
`ranges` (pylablib.devices.Ophir.base.TRangeInfo attribute), 641  
`rate` (pylablib.devices.NI.daq.TVoltageOutputClockParameters attribute), 621  
`rate` (pylablib.devices.Sirah.Matisse.TPiezoetDriveParameters attribute), 734  
`read()` (pylablib.core.devio.comm\_backend.FT232DeviceBackend method), 151  
`read()` (pylablib.core.devio.comm\_backend.IDeviceCommBackend method), 144  
`read()` (pylablib.core.devio.comm\_backend.NetworkDeviceBackend method), 153  
`read()` (pylablib.core.devio.comm\_backend.PyUSBDeviceBackend method), 156  
`read()` (pylablib.core.devio.comm\_backend.RecordedDeviceBackend method), 158  
`read()` (pylablib.core.devio.comm\_backend.SerialDeviceBackend method), 148  
`read()` (pylablib.core.devio.comm\_backend.VisaDeviceBackend method), 146  
`read()` (pylablib.core.devio.SCPISCPIDevice method), 140  
`read()` (pylablib.core.fileio.loadfile.BinaryTableInputFileFormatter method), 178  
`read()` (pylablib.core.fileio.loadfile.CSVTableInputFileFormat method), 177  
`read()` (pylablib.core.fileio.loadfile.DictionaryInputFileFormat method), 178  
`read()` (pylablib.core.fileio.loadfile.IInputFileFormat method), 176  
`read()` (pylablib.core.fileio.loadfile.ITextInputFileFormat method), 177  
`read()` (pylablib.devices.Andor.AndorSDK3.AndorSDK3Camera.BufferManager method), 481  
`read()` (pylablib.devices.AWG.generic.GenericAWG method), 407  
`read()` (pylablib.devices.AWG.specific.Agilent33220A method), 418  
`read()` (pylablib.devices.AWG.specific.Agilent33500 method), 412  
`read()` (pylablib.devices.AWG.specific.Agilent33500 method), 412  
`read()` (pylablib.devices.AWG.specific.Agilent33220A method), 418  
`read()` (pylablib.devices.AWG.specific.Agilent33500 method), 412  
`read()` (pylablib.devices.AWG.specific.InstekAFG2000 method), 429  
`read()` (pylablib.devices.AWG.specific.InstekAFG2225 method), 424  
`read()` (pylablib.devices.AWG.specific.RigolDG1000 method), 447  
`read()` (pylablib.devices.AWG.specific.RSInstekAFG21000 method), 435  
`read()` (pylablib.devices.AWG.specific.TektronixAFG1000 method), 441  
`read()` (pylablib.devices.Cryocon.base.Cryocon1x method), 538  
`read()` (pylablib.devices.Cryomagnetics.base.LM500 method), 542  
`read()` (pylablib.devices.Cryomagnetics.base.LM510 method), 546  
`read()` (pylablib.devices.Lakeshore.base.Lakeshore218 method), 596  
`read()` (pylablib.devices.Lakeshore.base.Lakeshore370 method), 601  
`read()` (pylablib.devices.NI.daq.NIDAQ method), 624  
`read()` (pylablib.devices.Sirah.Matisse.SirahMatisse method), 741  
`read()` (pylablib.devices.Tektronix.base.DPO2000 method), 765  
`read()` (pylablib.devices.Tektronix.base.ITektronixScope method), 754  
`read()` (pylablib.devices.Tektronix.base.TDS2000 method), 759  
`read()` (pylablib.devices.Thorlabs.serial.FW method), 814  
`read()` (pylablib.devices.Thorlabs.serial.FWv1 method), 818  
`read()` (pylablib.devices.Thorlabs.serial.MDT69xA method), 821  
`read()` (pylablib.devices.Thorlabs.serial.ThorlabsSerialInterface method), 811  
`read_binary_array_data()` (py-lablib.core.devio.SCPISCPIDevice method), 141  
`read_binary_array_data()` (py-lablib.devices.AWG.generic.GenericAWG method), 407  
`read_binary_array_data()` (py-lablib.devices.AWG.specific.Agilent33220A method), 418  
`read_binary_array_data()` (py-lablib.devices.AWG.specific.Agilent33500 method), 412  
`read_binary_array_data()` (py-lablib.devices.AWG.specific.InstekAFG2000 method), 430



<code>read_binary_array_data()</code> (py-lablib.devices.AWG.specific.InstekAFG2225 method), 424	<code>lablib.core.fileio.loadfile_utils</code> ), 182
<code>read_binary_array_data()</code> (py-lablib.devices.AWG.specific.RigolDG1000 method), 447	<code>read_directly</code> (py-lablib.devices.PrincetonInstruments.picam.PicamAttribute attribute), 710
<code>read_binary_array_data()</code> (py-lablib.devices.AWG.specific.RSInstekAFG21000 method), 435	<code>read_in_aux_port()</code> (py-lablib.devices.Andor.AndorSDK2.AndorSDK2Camera method), 469
<code>read_binary_array_data()</code> (py-lablib.devices.AWG.specific.TektronixAFG1000 method), 441	<code>read_multichar_term()</code> (py-lablib.core.devio.comm_backend.FT232DeviceBackend method), 151
<code>read_binary_array_data()</code> (py-lablib.devices.Cryocon.base.Cryocon1x method), 538	<code>read_multichar_term()</code> (py-lablib.core.devio.comm_backend.NetworkDeviceBackend method), 153
<code>read_binary_array_data()</code> (py-lablib.devices.Cryomagnetics.base.LM500 method), 542	<code>read_multichar_term()</code> (py-lablib.core.devio.comm_backend.PyUSBDeviceBackend method), 156
<code>read_binary_array_data()</code> (py-lablib.devices.Cryomagnetics.base.LM510 method), 546	<code>read_multichar_term()</code> (py-lablib.core.devio.comm_backend.SerialDeviceBackend method), 148
<code>read_binary_array_data()</code> (py-lablib.devices.Lakeshore.base.Lakeshore218 method), 596	<code>read_multiple_images()</code> (py-lablib.devices.AlliedVision.Bonito.BonitoIMAQCamera method), 461
<code>read_binary_array_data()</code> (py-lablib.devices.Lakeshore.base.Lakeshore370 method), 601	<code>read_multiple_images()</code> (py-lablib.devices.AlliedVision.Bonito.IBonitoCamera method), 454
<code>read_binary_array_data()</code> (py-lablib.devices.Sirah.Matisse.SirahMatisse method), 742	<code>read_multiple_images()</code> (py-lablib.devices.Andor.AndorSDK2.AndorSDK2Camera method), 475
<code>read_binary_array_data()</code> (py-lablib.devices.Tektronix.base.DPO2000 method), 766	<code>read_multiple_images()</code> (py-lablib.devices.Andor.AndorSDK3.AndorSDK3Camera method), 486
<code>read_binary_array_data()</code> (py-lablib.devices.Tektronix.base.ITektronixScope method), 754	<code>read_multiple_images()</code> (py-lablib.devices.Basler.pylon.BaslerPylonCamera method), 519
<code>read_binary_array_data()</code> (py-lablib.devices.Tektronix.base.TDS2000 method), 759	<code>read_multiple_images()</code> (py-lablib.devices.BitFlow.BitFlow.BitFlowCamera method), 530
<code>read_binary_array_data()</code> (py-lablib.devices.Thorlabs.serial.FW method), 815	<code>read_multiple_images()</code> (py-lablib.devices.BitFlow.BitFlow.BitFlowFrameGrabber method), 525
<code>read_binary_array_data()</code> (py-lablib.devices.Thorlabs.serial.FWv1 method), 818	<code>read_multiple_images()</code> (py-lablib.devices.DCAM.DCAM.DCAMCamera method), 556
<code>read_binary_array_data()</code> (py-lablib.devices.Thorlabs.serial.MDT69xA method), 821	<code>read_multiple_images()</code> (py-lablib.devices.IMAQ.IMAQ.IMAQCamera method), 572
<code>read_binary_array_data()</code> (py-lablib.devices.Thorlabs.serial.ThorlabsSerialInterface method), 811	<code>read_multiple_images()</code> (py-lablib.devices.IMAQ.IMAQ.IMAQFrameGrabber method), 567
<code>read_columns()</code> (in module py-lablib.core.fileio.parse_csv), 189	<code>read_multiple_images()</code> (py-lablib.devices.IMAQdx.IMAQdx.EthernetIMAQdxCamera method), 587
<code>read_dict_and_comments()</code> (in module py-	<code>read_multiple_images()</code> (py-lablib.devices.IMAQdx.IMAQdx.IMAQdxCamera

<i>method</i> ), 582	<i>method</i> ), 766
read_multiple_images() (py-lablib.devices.interface.camera.IAttributeCamera <i>method</i> ), 840	read_multiple_sweeps() (py-lablib.devices.Tektronix.base.ITektronixScope <i>method</i> ), 752
read_multiple_images() (py-lablib.devices.interface.camera.IBinROICamera <i>method</i> ), 859	read_multiple_sweeps() (py-lablib.devices.Tektronix.base.TDS2000 <i>method</i> ), 759
read_multiple_images() (py-lablib.devices.interface.camera.ICamera <i>method</i> ), 834	read_newest_image() (py-lablib.devices.AlliedVision.Bonito.BonitoIMAQCamera <i>method</i> ), 461
read_multiple_images() (py-lablib.devices.interface.camera.IExposureCamera <i>method</i> ), 849	read_newest_image() (py-lablib.devices.AlliedVision.Bonito.IBonitoCamera <i>method</i> ), 455
read_multiple_images() (py-lablib.devices.interface.camera.IGrabberAttributeCamera <i>method</i> ), 845	read_newest_image() (py-lablib.devices.Andor.AndorSDK2.AndorSDK2Camera <i>method</i> ), 475
read_multiple_images() (py-lablib.devices.interface.camera.IROICamera <i>method</i> ), 854	read_newest_image() (py-lablib.devices.Andor.AndorSDK3.AndorSDK3Camera <i>method</i> ), 485
read_multiple_images() (py-lablib.devices.PCO.SC2.PCOSC2Camera <i>method</i> ), 651	read_newest_image() (py-lablib.devices.Basler.pylon.BaslerPylonCamera <i>method</i> ), 519
read_multiple_images() (py-lablib.devices.Photometrics.pvcam.PvcamCamera <i>method</i> ), 665	read_newest_image() (py-lablib.devices.BitFlow.BitFlow.BitFlowCamera <i>method</i> ), 530
read_multiple_images() (py-lablib.devices.PhotonFocus.PhotonFocus.IPhotonFocusCamera <i>method</i> ), 676	read_newest_image() (py-lablib.devices.BitFlow.BitFlow.BitFlowFrameGrabber <i>method</i> ), 525
read_multiple_images() (py-lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusBitFlowCamera <i>method</i> ), 699	read_newest_image() (py-lablib.devices.DCAM.DCAM.DCAMCamera <i>method</i> ), 555
read_multiple_images() (py-lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusIMAQCamera <i>method</i> ), 683	read_newest_image() (py-lablib.devices.IMAQ.IMAQ.IMAQCamera <i>method</i> ), 573
read_multiple_images() (py-lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusSiSoCamera <i>method</i> ), 691	read_newest_image() (py-lablib.devices.IMAQ.IMAQ.IMAQFrameGrabber <i>method</i> ), 567
read_multiple_images() (py-lablib.devices.PrincetonInstruments.picam.PicamCamera <i>method</i> ), 717	read_newest_image() (py-lablib.devices.IMAQdx.IMAQdx.EthernetIMAQdxCamera <i>method</i> ), 587
read_multiple_images() (py-lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareCamera <i>method</i> ), 731	read_newest_image() (py-lablib.devices.IMAQdx.IMAQdx.IMAQdxCamera <i>method</i> ), 582
read_multiple_images() (py-lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareFrameGrabber <i>method</i> ), 725	read_newest_image() (py-lablib.devices.interface.camera.IAttributeCamera <i>method</i> ), 840
read_multiple_images() (py-lablib.devices.Thorlabs.TLCamera.ThorlabsTLCamera <i>method</i> ), 772	read_newest_image() (py-lablib.devices.interface.camera.IBinROICamera <i>method</i> ), 859
read_multiple_images() (py-lablib.devices.uc480.uc480.UC480Camera <i>method</i> ), 872	read_newest_image() (py-lablib.devices.interface.camera.ICamera <i>method</i> ), 834
read_multiple_sweeps() (py-lablib.devices.Tektronix.base.DPO2000 <i>method</i> ), 766	read_newest_image() (py-lablib.devices.interface.camera.IExposureCamera <i>method</i> ), 849

<i>method</i> ), 850	<i>method</i> ), 519
read_newest_image() lablib.devices.interface.camera.IGrabberAttributeCamera <i>method</i> ), 845	(py- read_oldest_image() lablib.devices.BitFlow.BitFlow.BitFlowCamera <i>method</i> ), 530
read_newest_image() lablib.devices.interface.camera.IROICamera <i>method</i> ), 854	(py- read_oldest_image() lablib.devices.BitFlow.BitFlow.BitFlowFrameGrabber <i>method</i> ), 526
read_newest_image() lablib.devices.PCO.SC2.PCOSC2Camera <i>method</i> ), 651	(py- read_oldest_image() lablib.devices.DCAM.DCAM.DCAMCamera <i>method</i> ), 555
read_newest_image() lablib.devices.Photometrics.pvcam.PvcamCamera <i>method</i> ), 667	(py- read_oldest_image() lablib.devices.IMAQ.IMAQ.IMAQCamera <i>method</i> ), 573
read_newest_image() lablib.devices.PhotonFocus.PhotonFocus.IPhotonFocusCamera <i>method</i> ), 676	(py- read_oldest_image() lablib.devices.IMAQ.IMAQ.IMAQFrameGrabber <i>method</i> ), 567
read_newest_image() lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusBitFlowCamera <i>method</i> ), 699	(py- read_oldest_image() lablib.devices.IMAQdx.IMAQdx.EthernetIMAQdxCamera <i>method</i> ), 587
read_newest_image() lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusIMAQCamera <i>method</i> ), 683	(py- read_oldest_image() lablib.devices.IMAQdx.IMAQdx.IMAQdxCamera <i>method</i> ), 583
read_newest_image() lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusSiSoCamera <i>method</i> ), 691	(py- read_oldest_image() lablib.devices.interface.camera.IAttributeCamera <i>method</i> ), 840
read_newest_image() lablib.devices.PrincetonInstruments.picam.PicamCamera <i>method</i> ), 716	(py- read_oldest_image() lablib.devices.interface.camera.IBinROICamera <i>method</i> ), 859
read_newest_image() lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareCamera <i>method</i> ), 731	(py- read_oldest_image() lablib.devices.interface.camera.ICamera <i>method</i> ), 834
read_newest_image() lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareFrameGrabber <i>method</i> ), 726	(py- read_oldest_image() lablib.devices.interface.camera.IExposureCamera <i>method</i> ), 850
read_newest_image() lablib.devices.Thorlabs.TLCamera.ThorlabsTLCamera <i>method</i> ), 775	(py- read_oldest_image() lablib.devices.interface.camera.IGrabberAttributeCamera <i>method</i> ), 845
read_newest_image() lablib.devices.uc480.uc480.UC480Camera <i>method</i> ), 871	(py- read_oldest_image() lablib.devices.interface.camera.IROICamera <i>method</i> ), 855
read_oldest_image() lablib.devices.AlliedVision.Bonito.BonitoIMAQCamera <i>method</i> ), 461	(py- read_oldest_image() lablib.devices.PCO.SC2.PCOSC2Camera <i>method</i> ), 651
read_oldest_image() lablib.devices.AlliedVision.Bonito.IBonitoCamera <i>method</i> ), 455	(py- read_oldest_image() lablib.devices.Photometrics.pvcam.PvcamCamera <i>method</i> ), 668
read_oldest_image() lablib.devices.Andor.AndorSDK2.AndorSDK2Camera <i>method</i> ), 475	(py- read_oldest_image() lablib.devices.PhotonFocus.PhotonFocus.IPhotonFocusCamera <i>method</i> ), 676
read_oldest_image() lablib.devices.Andor.AndorSDK3.AndorSDK3Camera <i>method</i> ), 485	(py- read_oldest_image() lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusBitFlowCa <i>method</i> ), 699
read_oldest_image() lablib.devices.Basler.pylon.BaslerPylonCamera	(py- read_oldest_image() lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusIMAQCam

method), 683

read\_oldest\_image() (py- readable (pylablib.devices.IMAQdx.IMAQdx.IMAQdxAttribute  
lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusSISOCamera (pylablib.devices.Photometrics.pvcam.PvcamAttribute  
method), 692 attribute), 660

read\_oldest\_image() (py- readable (pylablib.devices.PhotonFocus.PhotonFocus.PFCamAttribute  
lablib.devices.PrincetonInstruments.picam.PicamCamera attribute), 670

method), 716 ReadChangeLock (class in pylablib.core.thread.utils),  
321

read\_oldest\_image() (py- read\_change\_lock() (pylablib.core.thread.utils.ReadChangeLock  
lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareFrameGrabber method), 321

method), 731 read\_line() (pylablib.core.devio.comm\_backend.FT232DeviceBackend  
method), 151

read\_oldest\_image() (py- read\_line() (pylablib.core.devio.comm\_backend.IDeviceCommBackend  
lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareFrameGrabber method), 144

method), 726 read\_line() (pylablib.core.devio.comm\_backend.IDeviceCommBackend  
method), 144

read\_oldest\_image() (py- read\_line() (pylablib.core.devio.comm\_backend.NetworkDeviceBackend  
lablib.devices.Thorlabs.TLCamera.ThorlabsTLCamera method), 153

method), 775 read\_line() (pylablib.core.devio.comm\_backend.PyUSBDeviceBackend  
method), 156

read\_oldest\_image() (py- read\_line() (pylablib.core.devio.comm\_backend.RecordedDeviceBackend  
lablib.devices.uc480.uc480.UC480Camera method), 158

method), 871 read\_line() (pylablib.core.devio.comm\_backend.SerialDeviceBackend  
method), 148

read\_raw\_data() (py- read\_line() (pylablib.core.devio.comm\_backend.VisaDeviceBackend  
lablib.devices.Tektronix.base.DPO2000 method), 146

method), 766 read\_lines() (pylablib.core.devio.comm\_backend.FT232DeviceBackend  
method), 152

read\_raw\_data() (py- read\_lines() (pylablib.core.devio.comm\_backend.IDeviceCommBackend  
lablib.devices.Tektronix.base.ITektronixScope method), 144

method), 752 read\_lines() (pylablib.core.devio.comm\_backend.IDeviceCommBackend  
method), 144

read\_raw\_data() (py- read\_lines() (pylablib.core.devio.comm\_backend.IDeviceCommBackend  
lablib.devices.Tektronix.base.TDS2000 method), 154

method), 759 read\_lines() (pylablib.core.devio.comm\_backend.IDeviceCommBackend  
method), 154

read\_sweep() (pylablib.devices.Tektronix.base.DPO2000 method), 154

method), 766 read\_lines() (pylablib.core.devio.comm\_backend.IDeviceCommBackend  
method), 154

read\_sweep() (pylablib.devices.Tektronix.base.ITektronixScope method), 157

method), 752 read\_lines() (pylablib.core.devio.comm\_backend.IDeviceCommBackend  
method), 157

read\_sweep() (pylablib.devices.Tektronix.base.TDS2000 method), 159

method), 759 read\_lines() (pylablib.core.devio.comm\_backend.IDeviceCommBackend  
method), 159

read\_table() (in module py- read\_lines() (pylablib.core.devio.comm\_backend.IDeviceCommBackend  
lablib.core.fileio.parse\_csv), 190 method), 149

read\_trigger() (py- read\_lines() (pylablib.core.devio.comm\_backend.IDeviceCommBackend  
lablib.devices.AlliedVision.Bonito.BonitoIMAQCamera method), 146

method), 461 reason (pylablib.devices.Attocube.anc350.ANC350.Reply  
attribute), 507

read\_trigger() (py- reboot() (pylablib.devices.PCO.SC2.PCOS2Camera  
lablib.devices.IMAQ.IMAQ.IMAQCamera method), 645

method), 573 reboot() (pylablib.devices.PCO.SC2.PCOS2Camera  
method), 645

read\_trigger() (py- reboot() (pylablib.devices.PCO.SC2.PCOS2Camera  
lablib.devices.IMAQ.IMAQ.IMAQFrameGrabber method), 824

method), 563 reconnect() (pylablib.core.devio.SCPI.SCPIDevice  
method), 138

read\_trigger() (py- reconnect() (pylablib.core.devio.SCPI.SCPIDevice  
lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusSISOCamera (pylablib.devices.AWG.generic.GenericAWG  
method), 683 method), 407

readable (pylablib.devices.Andor.AndorSDK3.AndorSDK3Attribute method), 418

attribute), 477 reconnect() (pylablib.devices.AWG.generic.GenericAWG  
method), 412

readable (pylablib.devices.Basler.pylon.BaslerPylonAttribute method), 412

attribute), 512 reconnect() (pylablib.devices.AWG.generic.GenericAWG  
method), 412

readable (pylablib.devices.DCAM.DCAM.DCAMAttribute reconnect() (pylablib.devices.AWG.generic.GenericAWG  
attribute), 549 method), 430

[reconnect\(\)](#) ([pylablib.devices.AWG.specific.InstekAFG2225](#) [recv\\_all\(\)](#) ([pylablib.core.utils.net.ClientSocket](#) [method](#)), 424 [recv\\_all\(\)](#) ([pylablib.core.utils.net.ClientSocket](#) [method](#)), 390  
[reconnect\(\)](#) ([pylablib.devices.AWG.specific.RigolDG1000](#) [recv\\_comm\(\)](#) ([pylablib.devices.Thorlabs.elliptec.ElliptecMotor](#) [method](#)), 447 [recv\\_comm\(\)](#) ([pylablib.devices.Thorlabs.elliptec.ElliptecMotor](#) [method](#)), 778  
[reconnect\(\)](#) ([pylablib.devices.AWG.specific.RSInstekAFG221000](#) [recv\\_comm\(\)](#) ([pylablib.devices.Thorlabs.kinesis.BasicKinesisDevice](#) [method](#)), 436 [recv\\_comm\(\)](#) ([pylablib.devices.Thorlabs.kinesis.BasicKinesisDevice](#) [method](#)), 782  
[reconnect\(\)](#) ([pylablib.devices.AWG.specific.TektronixAFG1000](#) [recv\\_comm\(\)](#) ([pylablib.devices.Thorlabs.kinesis.KinesisDevice](#) [method](#)), 442 [recv\\_comm\(\)](#) ([pylablib.devices.Thorlabs.kinesis.KinesisDevice](#) [method](#)), 789  
[reconnect\(\)](#) ([pylablib.devices.Cryocon.base.Cryocon1x](#) [recv\\_comm\(\)](#) ([pylablib.devices.Thorlabs.kinesis.KinesisMotor](#) [method](#)), 538 [recv\\_comm\(\)](#) ([pylablib.devices.Thorlabs.kinesis.KinesisMotor](#) [method](#)), 800  
[reconnect\(\)](#) ([pylablib.devices.Cryomagnetics.base.LM500](#) [recv\\_comm\(\)](#) ([pylablib.devices.Thorlabs.kinesis.KinesisPiezoMotor](#) [method](#)), 542 [recv\\_comm\(\)](#) ([pylablib.devices.Thorlabs.kinesis.KinesisPiezoMotor](#) [method](#)), 804  
[reconnect\(\)](#) ([pylablib.devices.Cryomagnetics.base.LM510](#) [recv\\_comm\(\)](#) ([pylablib.devices.Thorlabs.kinesis.KinesisQuadDetector](#) [method](#)), 546 [recv\\_comm\(\)](#) ([pylablib.devices.Thorlabs.kinesis.KinesisQuadDetector](#) [method](#)), 809  
[reconnect\(\)](#) ([pylablib.devices.Lakeshore.base.Lakeshore218v](#) [recv\\_comm\(\)](#) ([pylablib.devices.Thorlabs.kinesis.MFF](#) [method](#)), 597 [recv\\_comm\(\)](#) ([pylablib.devices.Thorlabs.kinesis.MFF](#) [method](#)), 793  
[reconnect\(\)](#) ([pylablib.devices.Lakeshore.base.Lakeshore370v](#) [recv\\_decilen\(\)](#) ([pylablib.core.utils.net.ClientSocket](#) [method](#)), 602 [recv\\_decilen\(\)](#) ([pylablib.core.utils.net.ClientSocket](#) [method](#)), 390  
[reconnect\(\)](#) ([pylablib.devices.Sirah.Matisse.SirahMatisse](#) [recv\\_delimiter\(\)](#) ([pylablib.core.utils.net.ClientSocket](#) [method](#)), 742 [recv\\_delimiter\(\)](#) ([pylablib.core.utils.net.ClientSocket](#) [method](#)), 390  
[reconnect\(\)](#) ([pylablib.devices.Tektronix.base.DPO2000](#) [recv\\_fixedlen\(\)](#) ([pylablib.core.utils.net.ClientSocket](#) [method](#)), 766 [recv\\_fixedlen\(\)](#) ([pylablib.core.utils.net.ClientSocket](#) [method](#)), 390  
[reconnect\(\)](#) ([pylablib.devices.Tektronix.base.ITektronixScope](#) [recv\\_JSON\(\)](#) ([in module pylablib.core.utils.net](#)), 390 [recv\\_numpy\(\)](#) ([pylablib.core.utils.ipc.IPCChannel](#) [method](#)), 383  
[reconnect\(\)](#) ([pylablib.devices.Thorlabs.serial.FW](#) [recv\\_numpy\(\)](#) ([pylablib.core.utils.ipc.PipeIPCChannel](#) [method](#)), 815 [recv\\_numpy\(\)](#) ([pylablib.core.utils.ipc.PipeIPCChannel](#) [method](#)), 384  
[reconnect\(\)](#) ([pylablib.devices.Thorlabs.serial.FWv1](#) [recv\\_numpy\(\)](#) ([pylablib.core.utils.ipc.SharedMemIPCChannel](#) [method](#)), 818 [recv\\_numpy\(\)](#) ([pylablib.core.utils.ipc.SharedMemIPCChannel](#) [method](#)), 384  
[reconnect\(\)](#) ([pylablib.devices.Thorlabs.serial.MDT69xA](#) [refresh\\_acquisition\(\)](#) ([pylablib.devices.IMAQdx.IMAQdx.EthernetIMAQdxCamera](#) [method](#)), 821 [refresh\\_acquisition\(\)](#) ([pylablib.devices.IMAQdx.IMAQdx.EthernetIMAQdxCamera](#) [method](#)), 588  
[reconnect\(\)](#) ([pylablib.devices.Thorlabs.serial.ThorlabsSerialInterface](#) [refresh\\_acquisition\(\)](#) ([pylablib.devices.IMAQdx.IMAQdx.IMAQdxCamera](#) [method](#)), 811 [refresh\\_acquisition\(\)](#) ([pylablib.devices.IMAQdx.IMAQdx.IMAQdxCamera](#) [method](#)), 580  
[RecordedDeviceBackend](#) ([class in pylablib.core.devio.comm\\_backend](#)), 157 [register\(\)](#) ([pylablib.devices.Basler.pylon.BaslerPylonCamera.BufferManager](#) [method](#)), 515  
[rectangle\\_k\(\)](#) ([in module pylablib.core.dataproc.specfunc](#)), 124 [regular\\_grid\\_from\\_scatter\(\)](#) ([in module pylablib.core.dataproc.interpolate](#)), 123  
[rectangle\\_w\(\)](#) ([in module pylablib.core.dataproc.specfunc](#)), 124 [relative\\_path\(\)](#) ([in module pylablib.core.utils.files](#)), 362  
[rectangle\\_w\\_ft\(\)](#) ([in module pylablib.core.dataproc.specfunc](#)), 125 [RelayBoard](#) ([class in pylablib.devices.Conrad.base](#)), 533  
[recursive\\_map\(\)](#) ([in module pylablib.core.utils.general](#)), 375 [RelayBoard.NoParameterCaller](#) ([class in pylablib.devices.Conrad.base](#)), 534  
[recv\(\)](#) ([pylablib.core.utils.ipc.IPCChannel](#) [method](#)), 383 [RelayBoard.TMessage](#) ([class in pylablib.devices.Conrad.base](#)), 533  
[recv\(\)](#) ([pylablib.core.utils.ipc.PipeIPCChannel](#) [method](#)), 384 [release\(\)](#) ([pylablib.core.thread.synchronizing.QLockNotifier](#) [method](#)), 319  
[recv\(\)](#) ([pylablib.core.utils.ipc.SharedMemIPCChannel](#) [method](#)), 384 [relevant](#) ([pylablib.devices.PrincetonInstruments.picam.PicamAttribute](#) [attribute](#)), 710  
[recv\(\)](#) ([pylablib.core.utils.net.ClientSocket](#) [method](#)), 390 [reload\\_all\(\)](#) ([in module pylablib](#)), 875  
[recv\\_ack\(\)](#) ([pylablib.core.utils.net.ClientSocket](#) [method](#)), 390



[reload\\_package\\_modules\(\) \(in module py-lablib.core.utils.module\), 386](#)  
[rem\(\) \(pylablib.core.utils.functions.AttrObjectProperty method\), 373](#)  
[rem\(\) \(pylablib.core.utils.functions.IObjectProperty method\), 373](#)  
[rem\(\) \(pylablib.core.utils.functions.MethodObjectProperty method\), 373](#)  
[remap\\_axes\(\) \(pylablib.devices.Arcus.performax.GenericPerformax2EXStage method\), 489](#)  
[remap\\_axes\(\) \(pylablib.devices.Arcus.performax.Performax2EXStage method\), 496](#)  
[remap\\_axes\(\) \(pylablib.devices.Arcus.performax.Performax4EXStage method\), 493](#)  
[remap\\_axes\(\) \(pylablib.devices.Arcus.performax.PerformaxDMX16ASStage method\), 499](#)  
[remap\\_axes\(\) \(pylablib.devices.Attocube.anc300.ANC300 method\), 506](#)  
[remap\\_axes\(\) \(pylablib.devices.Attocube.anc350.ANC350 method\), 510](#)  
[remap\\_axes\(\) \(pylablib.devices.interface.stage.IMultiaxisStage method\), 863](#)  
[remap\\_axes\(\) \(pylablib.devices.Newport.picomotor.Picomotor8742 method\), 631](#)  
[remap\\_axes\(\) \(pylablib.devices.PhysikInstrumente.base.GenericPicoController method\), 704](#)  
[remap\\_axes\(\) \(pylablib.devices.PhysikInstrumente.base.PIE516 method\), 707](#)  
[remap\\_axes\(\) \(pylablib.devices.SmarAct.scu3d.SCUC3D method\), 747](#)  
[remap\\_axes\(\) \(pylablib.devices.Thorlabs.kinesis.KinesisDevice method\), 790](#)  
[remap\\_axes\(\) \(pylablib.devices.Thorlabs.kinesis.KinesisMotor method\), 800](#)  
[remap\\_axes\(\) \(pylablib.devices.Thorlabs.kinesis.KinesisPicoMotor method\), 804](#)  
[remap\\_axes\(\) \(pylablib.devices.Thorlabs.kinesis.MFF method\), 794](#)  
[remote\\_call\(\) \(in module py-lablib.core.thread.controller\), 293](#)  
[remove\\_batch\\_job\(\) \(py-lablib.core.thread.controller.QTaskThread method\), 304](#)  
[remove\\_child\(\) \(py-lablib.core.gui.widgets.container.IQContainer method\), 201](#)  
[remove\\_child\(\) \(py-lablib.core.gui.widgets.container.IQWidgetContainer method\), 205](#)  
[remove\\_child\(\) \(py-lablib.core.gui.widgets.container.QContainer method\), 204](#)  
[remove\\_child\(\) \(py-lablib.core.gui.widgets.container.QDialogContainer method\), 227](#)  
[remove\\_child\(\) \(py-lablib.core.gui.widgets.container.QFrameContainer method\), 215](#)  
[remove\\_child\(\) \(py-lablib.core.gui.widgets.container.QGroupBoxContainer method\), 223](#)  
[remove\\_child\(\) \(py-lablib.core.gui.widgets.container.QScrollAreaContainer method\), 230](#)  
[remove\\_child\(\) \(py-lablib.core.gui.widgets.container.QScrollAreaContainer.QContainer method\), 227](#)  
[remove\\_child\(\) \(py-lablib.core.gui.widgets.container.QTabContainer method\), 233](#)  
[remove\\_child\(\) \(py-lablib.core.gui.widgets.container.QWidgetContainer method\), 212](#)  
[remove\\_child\(\) \(py-lablib.core.gui.widgets.param\\_table.ParamTable method\), 248](#)  
[remove\\_child\(\) \(py-lablib.core.gui.widgets.param\\_table.StatusTable method\), 259](#)  
[remove\\_dir\(\) \(in module pylablib.core.utils.files\), 363](#)  
[remove\\_dir\\_if\\_empty\(\) \(in module py-lablib.core.utils.files\), 363](#)  
[remove\\_exception\\_hook\(\) \(in module py-lablib.core.thread.controller\), 292](#)  
[remove\\_handler\(\) \(py-lablib.core.gui.value\\_handling.GUIValues method\), 278](#)  
[remove\\_indicator\\_handler\(\) \(py-lablib.core.gui.value\\_handling.GUIValues method\), 279](#)  
[remove\\_job\(\) \(pylablib.core.thread.controller.QTaskThread method\), 303](#)  
[remove\\_layout\\_element\(\) \(py-lablib.core.gui.widgets.container.IQWidgetContainer method\), 208](#)  
[remove\\_layout\\_element\(\) \(py-lablib.core.gui.widgets.container.QDialogContainer method\), 219](#)  
[remove\\_layout\\_element\(\) \(py-lablib.core.gui.widgets.container.QFrameContainer method\), 216](#)  
[remove\\_layout\\_element\(\) \(py-lablib.core.gui.widgets.container.QGroupBoxContainer method\), 223](#)  
[remove\\_layout\\_element\(\) \(py-lablib.core.gui.widgets.container.QScrollAreaContainer.QContainer method\), 227](#)

---

<code>remove_layout_element()</code>	(py-lablib.core.gui.widgets.container.QWidgetContainer method), 212	<code>lablib.core.gui.value_handling.CheckboxValueHandler</code> method), 272
<code>remove_layout_element()</code>	(py-lablib.core.gui.widgets.layout_manager.IQLayoutManagedWidget method), 239	<code>repr_single_value()</code> (py-lablib.core.gui.value_handling.ComboBoxValueHandler method), 275
<code>remove_layout_element()</code>	(py-lablib.core.gui.widgets.layout_manager.QLayoutManagedWidget method), 241	<code>repr_single_value()</code> (py-lablib.core.gui.value_handling.IBoolValueHandler method), 271
<code>remove_layout_element()</code>	(py-lablib.core.gui.widgets.param_table.ParamTable method), 250	<code>repr_single_value()</code> (py-lablib.core.gui.value_handling.ISingleValueHandler method), 269
<code>remove_layout_element()</code>	(py-lablib.core.gui.widgets.param_table.StatusTable method), 259	<code>repr_single_value()</code> (py-lablib.core.gui.value_handling.LabelValueHandler method), 271
<code>remove_longest_term()</code>	(in module py-lablib.core.devio.comm_backend), 144	<code>repr_single_value()</code> (py-lablib.core.gui.value_handling.LineEditValueHandler method), 270
<code>remove_observer()</code>	(py-lablib.core.utils.observer_pool.ObserverPool method), 393	<code>repr_single_value()</code> (py-lablib.core.gui.value_handling.ProgressBarValueHandler method), 276
<code>remove_path()</code>	(py-lablib.core.utils.general.StreamFileLogger method), 380	<code>repr_single_value()</code> (py-lablib.core.gui.value_handling.PushButtonValueHandler method), 273
<code>remove_shortcut()</code>	(py-lablib.core.utils.dictionary.PrefixShortcutTree method), 360	<code>repr_single_value()</code> (py-lablib.core.gui.value_handling.ToolButtonValueHandler method), 274
<code>remove_status_line()</code>	(in module py-lablib.devices.interface.camera), 861	<code>repr_value()</code> (pylablib.core.gui.value_handling.CheckboxValueHandler method), 273
<code>remove_status_line()</code>	(in module py-lablib.devices.PhotonFocus.PhotonFocus), 701	<code>repr_value()</code> (pylablib.core.gui.value_handling.ComboBoxValueHandler method), 275
<code>remove_stop_notifier()</code>	(py-lablib.core.thread.controller.QTaskThread method), 310	<code>repr_value()</code> (pylablib.core.gui.value_handling.GUIValues method), 280
<code>remove_stop_notifier()</code>	(py-lablib.core.thread.controller.QThreadController method), 300	<code>repr_value()</code> (pylablib.core.gui.value_handling.IBoolValueHandler method), 272
<code>remove_tab()</code>	(pylablib.core.gui.widgets.container.QTabContainer method), 231	<code>repr_value()</code> (pylablib.core.gui.value_handling.ISingleValueHandler method), 269
<code>remove_widget()</code>	(py-lablib.core.gui.widgets.param_table.ParamTable method), 244	<code>repr_value()</code> (pylablib.core.gui.value_handling.IValueHandler method), 265
<code>remove_widget()</code>	(py-lablib.core.gui.widgets.param_table.StatusTable method), 259	<code>repr_value()</code> (pylablib.core.gui.value_handling.LabelIndicatorHandler method), 277
<code>removed</code>	(pylablib.core.utils.dictionary.DictionaryDiff attribute), 334	<code>repr_value()</code> (pylablib.core.gui.value_handling.LabelValueHandler method), 271
<code>reopen()</code>	(pylablib.devices.Arduino.base.IArduinoDevice method), 501	<code>repr_value()</code> (pylablib.core.gui.value_handling.LineEditValueHandler method), 270
<code>rep</code>	(pylablib.core.utils.string.TConversionClass attribute), 398	<code>repr_value()</code> (pylablib.core.gui.value_handling.ProgressBarValueHandler method), 276
<code>repr</code>	(pylablib.devices.Basler.pylon.BaslerPylonAttribute attribute), 513	<code>repr_value()</code> (pylablib.core.gui.value_handling.PropertyValueHandler method), 268
<code>repr_single_value()</code>	(py-lablib.core.gui.value_handling.LabelValueHandler method), 271	<code>repr_value()</code> (pylablib.core.gui.value_handling.PushButtonValueHandler method), 273
		<code>repr_value()</code> (pylablib.core.gui.value_handling.StandardValueHandler method), 268
		<code>repr_value()</code> (pylablib.core.gui.value_handling.ToolButtonValueHandler method), 274

repr\_value() (pylablib.core.gui.value\_handling.VirtualValueHandler attribute), 614  
 repr\_value() (pylablib.core.gui.widgets.button.ToggleButton attribute), 621  
 repr\_value() (pylablib.core.gui.widgets.combo\_box.ComboBox attribute), 736  
 repr\_value() (pylablib.core.gui.widgets.edit.NumEdit attribute), 762  
 repr\_value() (pylablib.core.gui.widgets.label.EnumLabel attribute), 749  
 repr\_value() (pylablib.core.gui.widgets.label.NumLabel attribute), 756  
 request\_stop() (pylablib.core.thread.controller.QTaskThread method), 310  
 request\_stop() (pylablib.core.thread.controller.QThreadController method), 299  
 requires\_symmetric\_roi() (pylablib.devices.PCO.SC2.PCOSC2Camera method), 648  
 reraise() (in module pylablib.core.devio.comm\_backend), 142  
 reraise() (pylablib.core.utils.general.RetryOnException.ExceptionCatcher method), 377  
 ReraiseError (pylablib.core.devio.SCPI.SCPIDevice attribute), 138  
 ReraiseError (pylablib.devices.AWG.generic.GenericAWG attribute), 403  
 ReraiseError (pylablib.devices.AWG.specific.Agilent33220A attribute), 415  
 ReraiseError (pylablib.devices.AWG.specific.Agilent33500 attribute), 409  
 ReraiseError (pylablib.devices.AWG.specific.InstekAFG2000 attribute), 426  
 ReraiseError (pylablib.devices.AWG.specific.InstekAFG2225 attribute), 421  
 ReraiseError (pylablib.devices.AWG.specific.RigolDG1000 attribute), 444  
 ReraiseError (pylablib.devices.AWG.specific.RSInstekAFG21000 method), 461  
 ReraiseError (pylablib.devices.AWG.specific.TektronixAFG1000 attribute), 438  
 ReraiseError (pylablib.devices.Cryocon.base.Cryocon1x attribute), 536  
 ReraiseError (pylablib.devices.Cryomagnetics.base.LM500 attribute), 540  
 ReraiseError (pylablib.devices.Cryomagnetics.base.LM510 attribute), 544  
 ReraiseError (pylablib.devices.Lakeshore.base.Lakeshore218 attribute), 593  
 ReraiseError (pylablib.devices.Lakeshore.base.Lakeshore370 attribute), 599  
 ReraiseError (pylablib.devices.M2.solstis.Solstis attribute), 614  
 ReraiseError (pylablib.devices.NI.daq.NIDAQ attribute), 621  
 ReraiseError (pylablib.devices.Sirah.Matisse.SirahMatisse attribute), 736  
 ReraiseError (pylablib.devices.Tektronix.base.DPO2000 attribute), 762  
 ReraiseError (pylablib.devices.Tektronix.base.ITektronixScope attribute), 749  
 ReraiseError (pylablib.devices.Tektronix.base.TDS2000 attribute), 756  
 ReraiseError (pylablib.devices.Thorlabs.serial.FW attribute), 813  
 ReraiseError (pylablib.devices.Thorlabs.serial.FWv1 attribute), 817  
 ReraiseError (pylablib.devices.Thorlabs.serial.MDT69xA attribute), 820  
 ReraiseError (pylablib.devices.Thorlabs.serial.ThorlabsSerialInterface attribute), 809  
 res\_range (pylablib.devices.Lakeshore.base.TLakeshore370RangeSetting attribute), 598  
 rescale() (pylablib.core.dataproc.utils.Range method), 136  
 reset() (in module pylablib.core.thread.profile), 317  
 reset() (pylablib.core.dataproc.filters.RunningDebounceFilter method), 114  
 reset() (pylablib.core.dataproc.filters.RunningDecimationFilter method), 114  
 reset() (pylablib.core.devio.SCPI.SCPIDevice method), 139  
 reset() (pylablib.core.utils.general.Countdown method), 378  
 reset() (pylablib.core.utils.general.Timer method), 379  
 reset() (pylablib.core.utils.general.UIDGenerator method), 378  
 reset() (pylablib.devices.AlliedVision.Bonito.BonitoIMAQCamera method), 461  
 reset() (pylablib.devices.Andor.AndorSDK3.AndorSDK3Camera.Buffer method), 481  
 reset() (pylablib.devices.AWG.generic.GenericAWG method), 407  
 reset() (pylablib.devices.AWG.specific.Agilent33220A method), 418  
 reset() (pylablib.devices.AWG.specific.Agilent33500 method), 412  
 reset() (pylablib.devices.AWG.specific.InstekAFG2000 method), 430  
 reset() (pylablib.devices.AWG.specific.InstekAFG2225 method), 424  
 reset() (pylablib.devices.AWG.specific.RigolDG1000 method), 448



`reset()` (`pylablib.devices.AWG.specific.RSInstekAFG21000` method), 818  
`method`), 436 `reset()` (`pylablib.devices.Thorlabs.serial.MDT69xA`  
`reset()` (`pylablib.devices.AWG.specific.TektronixAFG1000` method), 821  
`method`), 442 `reset()` (`pylablib.devices.Thorlabs.serial.ThorlabsSerialInterface`  
`reset()` (`pylablib.devices.BitFlow.BitFlow.BitFlowCamera.BufferManager` method), 811  
`method`), 527 `reset()` (`pylablib.devices.Thorlabs.TLCAmera.ThorlabsTLCAmera.RingL`  
`reset()` (`pylablib.devices.BitFlow.BitFlow.BitFlowFrameGrabber.BufferManager` method), 770  
`method`), 522 `reset_api()` (in module `pylablib.devices.PCO.SC2`),  
`reset()` (`pylablib.devices.Cryocon.base.Cryocon1x` method), 538 644  
`method`), 538 `reset_board()` (`py-`  
`reset()` (`pylablib.devices.Cryomagnetics.base.LM500` lablib.devices.Arduino.base.IArduinoDevice  
`method`), 543 method), 501  
`reset()` (`pylablib.devices.Cryomagnetics.base.LM510` `reset_error()` (`py-`  
`method`), 546 lablib.devices.Pfeiffer.base.TPG260 method),  
`reset()` (`pylablib.devices.IMAQ.IMAQ.IMAQCamera` 656  
`method`), 573 `reset_overflows_counter()` (`py-`  
`reset()` (`pylablib.devices.IMAQ.IMAQ.IMAQFrameGrabber` lablib.devices.Andor.AndorSDK3.AndorSDK3Camera  
`method`), 561 method), 482  
`reset()` (`pylablib.devices.IMAQdx.IMAQdx.EthernetIMAQdxCamera` `reset()` (`py-`  
`method`), 588 lablib.core.gui.widgets.container.QScrollAreaContainer.QContain  
`reset()` (`pylablib.devices.IMAQdx.IMAQdx.IMAQdxCamera` method), 225  
`method`), 578 `restart()` (in module `pylablib.core.utils.general`), 381  
`reset()` (`pylablib.devices.interface.camera.FrameCounter` `restart()` (`pylablib.devices.OZOptics.base.DD100`  
`method`), 835 method), 636  
`reset()` (`pylablib.devices.interface.camera.FrameNotifier` `restart()` (`pylablib.devices.OZOptics.base.OZOpticsDevice`  
`method`), 836 method), 632  
`reset()` (`pylablib.devices.KJL.base.KJL300` method), 590 `restart()` (`pylablib.devices.OZOptics.base.TF100`  
`method`), 590 method), 634  
`reset()` (`pylablib.devices.Lakeshore.base.Lakeshore218` `restart_app()` (in module `py-`  
`method`), 597 lablib.core.thread.controller), 315  
`reset()` (`pylablib.devices.Lakeshore.base.Lakeshore370` `restart_batch_job()` (`py-`  
`method`), 602 lablib.core.thread.controller.QTaskThread  
`reset()` (`pylablib.devices.Newport.picomotor.Picomotor8742` method), 304  
`method`), 628 `restart_lib()` (in module `py-`  
`reset()` (`pylablib.devices.NI.daq.NIDAQ` method), 622 lablib.devices.Andor.AndorSDK2), 465  
`reset()` (`pylablib.devices.Ophir.base.VegaPowerMeter` `restart_lib()` (in module `py-`  
`method`), 641 lablib.devices.Andor.AndorSDK3), 477  
`reset()` (`pylablib.devices.PCO.SC2.PCOS2Camera.ScheduLib` `restart_lib()` (in module `py-`  
`method`), 646 lablib.devices.Basler.pylon), 511  
`reset()` (`pylablib.devices.PhotonFocus.PhotonFocus.PhotonFocusBitFlowCamera.BufferManager` module `py-`  
`method`), 695 lablib.devices.DCAM.DCAM), 548  
`reset()` (`pylablib.devices.PhotonFocus.PhotonFocus.PhotonFocusIMAQCamera` `reset()` (`py-`  
`method`), 683 lablib.devices.Newport.picomotor.Picomotor8742  
`reset()` (`pylablib.devices.Sirah.Matisse.SirahMatisse` method), 628  
`method`), 742 `retrieve()` (`pylablib.devices.Basler.pylon.BaslerPylonCamera.BufferM`  
`reset()` (`pylablib.devices.Tektronix.base.DPO2000` method), 516  
`method`), 766 `retry_clean_dir()` (in module `py-`  
`reset()` (`pylablib.devices.Tektronix.base.ITektronixScope` lablib.core.utils.files), 367  
`method`), 754 `retry_copy()` (in module `pylablib.core.utils.files`),  
`reset()` (`pylablib.devices.Tektronix.base.TDS2000` 366  
`method`), 760 `retry_copy_dir()` (in module `py-`  
`reset()` (`pylablib.devices.Thorlabs.serial.FW` method), 815 lablib.core.utils.files), 366  
`reset()` (`pylablib.devices.Thorlabs.serial.FWv1` `retry_ensure_dir()` (in module `py-`  
`method`), 815 lablib.core.utils.files), 366

[retry\\_move\(\)](#) (in module [pylablib.core.utils.files](#)), [366](#)  
[retry\\_move\\_dir\(\)](#) (in module [pylablib.core.utils.files](#)), [366](#)  
[retry\\_remove\(\)](#) (in module [pylablib.core.utils.files](#)), [366](#)  
[retry\\_remove\\_dir\(\)](#) (in module [pylablib.core.utils.files](#)), [367](#)  
[retry\\_remove\\_dir\\_if\\_empty\(\)](#) (in module [pylablib.core.utils.files](#)), [367](#)  
[retry\\_wait\(\)](#) (in module [pylablib.core.utils.general](#)), [377](#)  
[RetryOnException](#) (class in [pylablib.core.utils.general](#)), [376](#)  
[RetryOnException.ExceptionCatcher](#) (class in [pylablib.core.utils.general](#)), [377](#)  
[revision\\_number](#) ([pylablib.devices.HighFinesse.wlm.TDeviceInfo](#) attribute), [557](#)  
[RigoldG1000](#) (class in [pylablib.devices.AWG.specific](#)), [444](#)  
[RigoldG1000.NoParameterCaller](#) (class in [pylablib.devices.AWG.specific](#)), [444](#)  
[rise\\_speed](#) ([pylablib.devices.Sirah.Matisse.TScanParameters](#) attribute), [736](#)  
[rng](#) ([pylablib.devices.Ophir.base.TWavelengthInfo](#) attribute), [641](#)  
[ROI](#) (class in [pylablib.core.dataproc.image](#)), [121](#)  
[roi](#) ([pylablib.devices.interface.camera.TStatusLineDescription](#) attribute), [861](#)  
[rom\\_version](#) ([pylablib.devices.Ophir.base.TDeviceInfo](#) attribute), [640](#)  
[rotated2d\(\)](#) ([pylablib.core.dataproc.transform.Indexed2DTransform](#) method), [134](#)  
[rotated2d\(\)](#) ([pylablib.core.dataproc.transform.LinearTransform](#) method), [134](#)  
[round\\_significant\(\)](#) (in module [pylablib.core.utils.numerical](#)), [391](#)  
[route](#) ([pylablib.devices.Thorlabs.kinesis.TQuadDetectorOutputParameters](#) attribute), [806](#)  
[RSInstekAFG21000](#) (class in [pylablib.devices.AWG.specific](#)), [432](#)  
[RSInstekAFG21000.NoParameterCaller](#) (class in [pylablib.devices.AWG.specific](#)), [432](#)  
[run\(\)](#) ([pylablib.core.thread.controller.QTaskThread](#) method), [304](#)  
[run\(\)](#) ([pylablib.core.thread.controller.QThreadController](#) method), [295](#)  
[run\(\)](#) ([pylablib.core.thread.controller.QThreadControllerThread](#) method), [292](#)  
[run\\_as\\_batch\\_job\(\)](#) ([pylablib.core.thread.controller.QTaskThread](#) method), [304](#)  
[run\\_device\\_service\(\)](#) (in module [pylablib.core.utils.rpyc\\_utils](#)), [395](#)  
[running\(\)](#) ([pylablib.core.thread.controller.QTaskThread](#) method), [310](#)  
[running\(\)](#) ([pylablib.core.thread.controller.QThreadController](#) method), [299](#)  
[running\(\)](#) ([pylablib.core.utils.general.Countdown](#) method), [378](#)  
[RunningDebounceFilter](#) (class in [pylablib.core.dataproc.filters](#)), [114](#)  
[RunningDecimationFilter](#) (class in [pylablib.core.dataproc.filters](#)), [113](#)

## S

[same](#) ([pylablib.core.utils.dictionary.DictionaryDiff](#) attribute), [334](#)  
[samps\\_per\\_chan](#) ([pylablib.devices.NI.daq.TVoltageOutputClockParameters](#) attribute), [621](#)  
[save\(\)](#) ([pylablib.devices.BitFlow.BitFlow.CameraFileEditor](#) method), [532](#)  
[save\\_bin\(\)](#) (in module [pylablib.core.fileio.savefile](#)), [195](#)  
[save\\_bin\\_desc\(\)](#) (in module [pylablib.core.fileio.savefile](#)), [195](#)  
[save\\_csv\(\)](#) (in module [pylablib.core.fileio.savefile](#)), [194](#)  
[save\\_csv\\_desc\(\)](#) (in module [pylablib.core.fileio.savefile](#)), [194](#)  
[save\\_dict\(\)](#) (in module [pylablib.core.fileio.savefile](#)), [195](#)  
[save\\_file\(\)](#) ([pylablib.core.fileio.dict\\_entry.ExternalNumpyDictionaryEntry](#) method), [175](#)  
[save\\_file\(\)](#) ([pylablib.core.fileio.dict\\_entry.IExternalFileDictionaryEntry](#) method), [175](#)  
[save\\_generic\(\)](#) (in module [pylablib.core.fileio.savefile](#)), [196](#)  
[save\\_parameters\(\)](#) ([pylablib.devices.Newport.picomotor.Picomotor8742](#) method), [628](#)  
[save\\_preset\(\)](#) ([pylablib.devices.OZOptics.base.EPC04](#) method), [637](#)  
[scan\\_both\\_motors\(\)](#) ([pylablib.devices.Sirah.tuner.MatisseTuner](#) method), [744](#)  
[scan\\_both\\_motors\\_quick\(\)](#) ([pylablib.devices.Sirah.tuner.MatisseTuner](#) method), [744](#)  
[scan\\_centered\(\)](#) ([pylablib.devices.Sirah.tuner.MatisseTuner](#) method), [743](#)  
[scan\\_devices\(\)](#) ([pylablib.devices.Newport.picomotor.Picomotor8742](#) method), [629](#)

`scan_quick()` (`pylablib.devices.Sirah.tuner.MatisseTuner` method), 546  
`scan_quick_centered()` (`pylablib.devices.Sirah.tuner.MatisseTuner` method), 743  
`scan_steps()` (`pylablib.devices.Sirah.tuner.MatisseTuner` method), 743  
`schedule()` (`pylablib.core.thread.callsync.QDirectCallScheduler` method), 285  
`schedule()` (`pylablib.core.thread.callsync.QMulticastThreadCallScheduler` method), 291  
`schedule()` (`pylablib.core.thread.callsync.QMultiQueueScheduler` method), 290  
`schedule()` (`pylablib.core.thread.callsync.QQueueLengthLimitScheduler` method), 288  
`schedule()` (`pylablib.core.thread.callsync.QQueueScheduler` method), 286  
`schedule()` (`pylablib.core.thread.callsync.QQueueSizeLimitScheduler` method), 290  
`schedule()` (`pylablib.core.thread.callsync.QScheduler` method), 285  
`schedule()` (`pylablib.core.thread.callsync.QThreadCallScheduler` method), 291  
`schedule()` (`pylablib.core.thread.controller.QTaskThread.Job` method), 302  
`schedule_multiple_queues()` (in module `pylablib.core.thread.callsync`), 290  
`scheduled` (`pylablib.devices.PCO.SC2.TInternalBufferStatus` attribute), 645  
`scheduled_max` (`pylablib.devices.PCO.SC2.TInternalBufferStatus` attribute), 645  
`scheduler` (`pylablib.core.thread.controller.QTaskThread.TCommand` attribute), 302  
`SCPIDevice` (class in `pylablib.core.devio.SCPIDevice`), 138  
`SCPIDevice.NoParameterCaller` (class in `pylablib.core.devio.SCPIDevice`), 140  
`SCU3D` (class in `pylablib.devices.SmarAct.scu3d`), 746  
`SCU3D.NoParameterCaller` (class in `pylablib.devices.SmarAct.scu3d`), 747  
`search_frequency()` (`pylablib.devices.Thorlabs.elliptec.ElliptecMotor` method), 779  
`second` (`pylablib.devices.uc480.uc480.TTimestamp` attribute), 866  
`section()` (`pylablib.core.devio.backend_logger.BackendLogger` method), 142  
`section()` (`pylablib.core.devio.comm_backend.RecordedDeviceBackend` method), 158  
`select_channel()` (`pylablib.devices.Cryomagnetics.base.LM500` method), 540  
`select_channel()` (`pylablib.devices.Cryomagnetics.base.LM510` method), 546  
`select_channel()` (`pylablib.devices.Lakeshore.base.Lakeshore370` method), 599  
`select_channel()` (`pylablib.devices.Tektronix.base.DPO2000` method), 767  
`select_channel()` (`pylablib.devices.Tektronix.base.ITektronixScope` method), 751  
`select_channel()` (`pylablib.devices.Tektronix.base.TDS2000` method), 760  
`select_channel()` (`pylablib.devices.AWG.generic.GenericAWG` method), 403  
`select_current_channel()` (`pylablib.devices.AWG.specific.Agilent33220A` method), 418  
`select_current_channel()` (`pylablib.devices.AWG.specific.Agilent33500` method), 412  
`select_current_channel()` (`pylablib.devices.AWG.specific.InstekAFG2000` method), 430  
`select_current_channel()` (`pylablib.devices.AWG.specific.InstekAFG2225` method), 424  
`select_current_channel()` (`pylablib.devices.AWG.specific.RigolDG1000` method), 448  
`select_current_channel()` (`pylablib.devices.AWG.specific.RSInstekAFG21000` method), 436  
`select_current_channel()` (`pylablib.devices.AWG.specific.TektronixAFG1000` method), 442  
`send()` (`pylablib.core.thread.multicast_pool.MulticastPool` method), 316  
`send()` (`pylablib.core.utils.ipc.IIPCCChannel` method), 383  
`send()` (`pylablib.core.utils.ipc.PipeIPCCChannel` method), 383  
`send()` (`pylablib.core.utils.ipc.SharedMemIPCCChannel` method), 384  
`send()` (`pylablib.core.utils.net.ClientSocket` method), 390  
`send()` (`pylablib.core.utils.net.ClientSocket` method), 390  
`send_comm()` (`pylablib.devices.Thorlabs.elliptec.ElliptecMotor` method), 778  
`send_comm()` (`pylablib.devices.Thorlabs.kinesis.BasicKinesisDevice` method), 782  
`send_comm()` (`pylablib.devices.Thorlabs.kinesis.KinesisDevice` method), 782

*method*), 790  
 send\_comm() (pylablib.devices.Thorlabs.kinesis.KinesisMotor *method*), 800  
 send\_comm() (pylablib.devices.Thorlabs.kinesis.KinesisPiezoMotor *method*), 804  
 send\_comm() (pylablib.devices.Thorlabs.kinesis.KinesisQuadDetector *method*), 809  
 send\_comm() (pylablib.devices.Thorlabs.kinesis.MFF *method*), 794  
 send\_comm\_data() (pylablib.devices.Thorlabs.kinesis.BasicKinesisDevice *method*), 782  
 send\_comm\_data() (pylablib.devices.Thorlabs.kinesis.KinesisDevice *method*), 790  
 send\_comm\_data() (pylablib.devices.Thorlabs.kinesis.KinesisMotor *method*), 800  
 send\_comm\_data() (pylablib.devices.Thorlabs.kinesis.KinesisPiezoMotor *method*), 804  
 send\_comm\_data() (pylablib.devices.Thorlabs.kinesis.KinesisQuadDetector *method*), 809  
 send\_comm\_data() (pylablib.devices.Thorlabs.kinesis.MFF *method*), 794  
 send\_command() (pylablib.devices.Leybold.base.GenericITR *method*), 607  
 send\_command() (pylablib.devices.Leybold.base.ITR90 *method*), 610  
 send\_decllen() (pylablib.core.utils.net.ClientSocket *method*), 390  
 send\_delimiter() (pylablib.core.utils.net.ClientSocket *method*), 390  
 send\_fixedlen() (pylablib.core.utils.net.ClientSocket *method*), 390  
 send\_interrupt() (pylablib.core.thread.controller.QTaskThread *method*), 310  
 send\_interrupt() (pylablib.core.thread.controller.QThreadController *method*), 298  
 send\_message() (pylablib.core.thread.controller.QTaskThread *method*), 310  
 send\_message() (pylablib.core.thread.controller.QThreadController *method*), 298  
 send\_multicast() (pylablib.core.thread.controller.QTaskThread *method*), 310  
 send\_multicast() (pylablib.core.thread.controller.QThreadController *method*), 297  
 send\_multicast\_sync() (pylablib.core.thread.controller.QTaskThread *method*), 311  
 send\_multicast\_sync() (pylablib.core.thread.controller.QThreadController *method*), 297  
 send\_numpy() (pylablib.core.utils.ipc.IIPCChannel *method*), 383  
 send\_numpy() (pylablib.core.utils.ipc.PipeIPCChannel *method*), 384  
 send\_numpy() (pylablib.core.utils.ipc.SharedMemIPCChannel *method*), 384  
 send\_software\_trigger() (pylablib.devices.AlliedVision.Bonito.BonitoIMAQCamera *method*), 461  
 send\_software\_trigger() (pylablib.devices.Andor.AndorSDK2.AndorSDK2Camera *method*), 470  
 send\_software\_trigger() (pylablib.devices.DCAM.DCAM.DCAMCamera *method*), 551  
 send\_software\_trigger() (pylablib.devices.IMAQ.IMAQ.IMAQCamera *method*), 573  
 send\_software\_trigger() (pylablib.devices.IMAQ.IMAQ.IMAQFrameGrabber *method*), 563  
 send\_software\_trigger() (pylablib.devices.PCO.SC2.PCOS2Camera *method*), 646  
 send\_software\_trigger() (pylablib.devices.Photometrics.pvcam.PvcamCamera *method*), 664  
 send\_software\_trigger() (pylablib.devices.PhotonFocus.PhotonFocus.PhotonFocusIMAQCamera *method*), 684  
 send\_software\_trigger() (pylablib.devices.Thorlabs.TLCamera.ThorlabsTLCamera *method*), 771  
 send\_sync() (pylablib.core.thread.controller.QTaskThread *method*), 311  
 send\_sync() (pylablib.core.thread.controller.QThreadController *method*), 298  
 sens\_id (pylablib.devices.uc480.uc480.TCameraInfo *attribute*), 864  
 sensor (pylablib.devices.Leybold.base.TDeviceInfo *attribute*), 607  
 sensor (pylablib.devices.PCO.SC2.TDeviceInfo *attribute*), 644



`serial` (`pylablib.devices.Attocube.anc300.TDeviceInfo` attribute), 502  
`serial` (`pylablib.devices.Basler.pylon.TCameraInfo` attribute), 511  
`serial` (`pylablib.devices.Basler.pylon.TDeviceInfo` attribute), 514  
`serial` (`pylablib.devices.Lakeshore.base.TLakeshore218CaseHeader` attribute), 593  
`serial` (`pylablib.devices.LaserQuantum.base.TDeviceInfo` attribute), 604  
`serial` (`pylablib.devices.LighthousePhotonics.base.TDeviceInfo` attribute), 611  
`serial` (`pylablib.devices.Ophir.base.TDeviceInfo` attribute), 640  
`serial` (`pylablib.devices.Ophir.base.THeadInfo` attribute), 640  
`serial` (`pylablib.devices.Photometrics.pvcam.TDeviceInfo` attribute), 661  
`serial` (`pylablib.devices.Toptica.ibeam.TDeviceInfo` attribute), 823  
`serial_flush()` (`pylablib.devices.AlliedVision.Bonito.BonitoIMAQCamera` method), 461  
`serial_flush()` (`pylablib.devices.IMAQ.IMAQ.IMAQCamera` method), 573  
`serial_flush()` (`pylablib.devices.IMAQ.IMAQ.IMAQFrameGrabber` method), 564  
`serial_flush()` (`pylablib.devices.PhotonFocus.PhotonFocus.PhotonFocusIMAQCamera` method), 684  
`serial_no` (`pylablib.devices.Thorlabs.elliptec.TDeviceInfo` attribute), 777  
`serial_no` (`pylablib.devices.Thorlabs.kinesis.TDeviceInfo` attribute), 781  
`serial_number` (`pylablib.devices.AlliedVision.Bonito.TDeviceInfo` attribute), 450  
`serial_number` (`pylablib.devices.Andor.AndorSDK2.TDeviceInfo` attribute), 465  
`serial_number` (`pylablib.devices.Andor.AndorSDK3.TDeviceInfo` attribute), 478  
`serial_number` (`pylablib.devices.DCAM.DCAM.TDeviceInfo` attribute), 550  
`serial_number` (`pylablib.devices.HighFinesse.wlm.TDeviceInfo` attribute), 557  
`serial_number` (`pylablib.devices.IMAQ.IMAQ.TDeviceInfo` attribute), 561  
`serial_number` (`pylablib.devices.IMAQdx.IMAQdx.TCameraInfo` attribute), 576  
`serial_number` (`pylablib.devices.IMAQdx.IMAQdx.TDeviceInfo` attribute), 578  
`serial_number` (`pylablib.devices.NI.daq.TDeviceInfo` attribute), 620  
`serial_number` (`pylablib.devices.PCO.SC2.TDeviceInfo` attribute), 644  
`serial_number` (`pylablib.devices.PhotonFocus.PhotonFocus.TDeviceInfo` attribute), 672  
`serial_number` (`pylablib.devices.PrincetonInstruments.picam.TCameraInfo` attribute), 709  
`serial_number` (`pylablib.devices.PrincetonInstruments.picam.TDeviceInfo` attribute), 711  
`serial_number` (`pylablib.devices.Thorlabs.TLCamera.TDeviceInfo` attribute), 770  
`serial_number` (`pylablib.devices.uc480.uc480.TCameraInfo` attribute), 864  
`serial_number` (`pylablib.devices.uc480.uc480.TDeviceInfo` attribute), 865  
`serial_query()` (`pylablib.devices.AlliedVision.Bonito.BonitoIMAQCamera` method), 461  
`serial_query()` (`pylablib.devices.AlliedVision.Bonito.IBonitoCamera` method), 451  
`serial_read()` (`pylablib.devices.AlliedVision.Bonito.BonitoIMAQCamera` method), 461  
`serial_read()` (`pylablib.devices.IMAQ.IMAQ.IMAQCamera` method), 573  
`serial_read()` (`pylablib.devices.IMAQ.IMAQ.IMAQFrameGrabber` method), 564  
`serial_read()` (`pylablib.devices.PhotonFocus.PhotonFocus.PhotonFocusIMAQCamera` method), 684  
`serial_readline()` (`pylablib.devices.AlliedVision.Bonito.BonitoIMAQCamera` method), 462  
`serial_readline()` (`pylablib.devices.IMAQ.IMAQ.IMAQCamera` method), 573  
`serial_readline()`

`lablib.devices.IMAQ.IMAQ.IMAQFrameGrabber` `method`), 837  
`method`), 564 `set_all_attribute_values()` (py-  
`serial_readline()` (py- `lablib.devices.Photometrics.pvcam.PvcamCamera`  
`lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusIMAQCamera` `method`), 663  
`method`), 684 `set_all_attribute_values()` (py-  
`serial_write()` (py- `lablib.devices.PhotonFocus.PhotonFocus.IPhotonFocusCamera`  
`lablib.devices.AlliedVision.Bonito.BonitoIMAQCamera` `method`), 672  
`method`), 462 `set_all_attribute_values()` (py-  
`serial_write()` (py- `lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusBitFlowCamera`  
`lablib.devices.IMAQ.IMAQ.IMAQCamera` `method`), 699  
`method`), 574 `set_all_attribute_values()` (py-  
`serial_write()` (py- `lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusIMAQCamera`  
`lablib.devices.IMAQ.IMAQ.IMAQFrameGrabber` `method`), 684  
`method`), 564 `set_all_attribute_values()` (py-  
`serial_write()` (py- `lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusSiSoCamera`  
`lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusIMAQCamera` `method`), 692  
`method`), 684 `set_all_attribute_values()` (py-  
`SerialDeviceBackend` (class in py- `lablib.devices.PrincetonInstruments.picam.PicamCamera`  
`lablib.core.devio.comm_backend`), 147 `method`), 712  
`Series1DWrapper` (class in py- `set_all_frequencies()` (py-  
`lablib.core.dataproc.table_wrap`), 127 `lablib.devices.OZOptics.base.EPC04` `method`),  
`Series1DWrapper.Accessor` (class in py- 637  
`lablib.core.dataproc.table_wrap`), 128 `set_all_grabber_attribute_values()` (py-  
`set()` (pylablib.core.utils.functions.AttrObjectProperty `lablib.devices.interface.camera.IGrabberAttributeCamera`  
`method`), 373 `method`), 842  
`set()` (pylablib.core.utils.functions.IObjectProperty `set_all_grabber_attribute_values()` (py-  
`method`), 373 `lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusSiSoCamera`  
`set()` (pylablib.core.utils.functions.MethodObjectProperty `method`), 692  
`method`), 373 `set_all_grabber_attribute_values()` (py-  
`set_acquisition_mode()` (py- `lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareCamera`  
`lablib.devices.Andor.AndorSDK2.AndorSDK2Camera` `method`), 731  
`method`), 470 `set_all_grabber_attribute_values()` (py-  
`set_active_channel()` (py- `lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareFrameGrabber`  
`lablib.devices.HighFinesse.wlm.WLM` `method`), `method`), 721  
559 `set_all_indicators()` (py-  
`set_addr()` (pylablib.devices.Newport.picomotor.Picomotor8742 `lablib.core.gui.value_handling.GUIValues`  
`method`), 629 `method`), 280  
`set_all_attribute_values()` (py- `set_all_indicators()` (py-  
`lablib.devices.Andor.AndorSDK3.AndorSDK3Camera` `lablib.core.gui.widgets.container.IQContainer`  
`method`), 480 `method`), 202  
`set_all_attribute_values()` (py- `set_all_indicators()` (py-  
`lablib.devices.Basler.pylon.BaslerPylonCamera` `lablib.core.gui.widgets.container.IQWidgetContainer`  
`method`), 514 `method`), 208  
`set_all_attribute_values()` (py- `set_all_indicators()` (py-  
`lablib.devices.DCAM.DCAM.DCAMCamera` `lablib.core.gui.widgets.container.QContainer`  
`method`), 551 `method`), 204  
`set_all_attribute_values()` (py- `set_all_indicators()` (py-  
`lablib.devices.IMAQdx.IMAQdx.EthernetIMAQdxCamera` `lablib.core.gui.widgets.container.QDialogContainer`  
`method`), 588 `method`), 219  
`set_all_attribute_values()` (py- `set_all_indicators()` (py-  
`lablib.devices.IMAQdx.IMAQdx.IMAQdxCamera` `lablib.core.gui.widgets.container.QFrameContainer`  
`method`), 579 `method`), 216  
`set_all_attribute_values()` (py- `set_all_indicators()` (py-  
`lablib.devices.interface.camera.IAttributeCamera` `lablib.core.gui.widgets.container.QGroupBoxContainer`

<i>method</i> ), 224	<i>method</i> ), 212
set_all_indicators() (py- lablib.core.gui.widgets.container.QScrollAreaContainer <i>method</i> ), 230	set_all_values() (py- lablib.core.gui.widgets.param_table.ParamTable <i>method</i> ), 248
set_all_indicators() (py- lablib.core.gui.widgets.container.QScrollAreaContainer.QComboBoxWidget <i>method</i> ), 227	set_all_values() (py- lablib.core.gui.widgets.param_table.StatusTable <i>method</i> ), 259
set_all_indicators() (py- lablib.core.gui.widgets.container.QTabContainer <i>method</i> ), 233	set_all_voltages() (py- lablib.devices.OZOptics.base.EPC04 <i>method</i> ), 636
set_all_indicators() (py- lablib.core.gui.widgets.container.QWidgetContainer <i>method</i> ), 212	set_amp_mode() (py- lablib.devices.Andor.AndorSDK2.AndorSDK2Camera <i>method</i> ), 468
set_all_indicators() (py- lablib.core.gui.widgets.param_table.ParamTable <i>method</i> ), 248	set_amplitude() (py- lablib.devices.AWG.generic.GenericAWG <i>method</i> ), 403
set_all_indicators() (py- lablib.core.gui.widgets.param_table.StatusTable <i>method</i> ), 259	set_amplitude() (py- lablib.devices.AWG.specific.Agilent33220A <i>method</i> ), 418
set_all_relays() (py- lablib.devices.Conrad.base.RelayBoard <i>method</i> ), 534	set_amplitude() (py- lablib.devices.AWG.specific.Agilent33500 <i>method</i> ), 412
set_all_values() (py- lablib.core.gui.value_handling.GUIValues <i>method</i> ), 280	set_amplitude() (py- lablib.devices.AWG.specific.InstekAFG2000 <i>method</i> ), 430
set_all_values() (py- lablib.core.gui.widgets.container.IQContainer <i>method</i> ), 202	set_amplitude() (py- lablib.devices.AWG.specific.InstekAFG2225 <i>method</i> ), 421
set_all_values() (py- lablib.core.gui.widgets.container.IQWidgetContainer <i>method</i> ), 208	set_amplitude() (py- lablib.devices.AWG.specific.RigolDG1000 <i>method</i> ), 448
set_all_values() (py- lablib.core.gui.widgets.container.QContainer <i>method</i> ), 204	set_amplitude() (py- lablib.devices.AWG.specific.RSInstekAFG21000 <i>method</i> ), 436
set_all_values() (py- lablib.core.gui.widgets.container.QDialogContainer <i>method</i> ), 220	set_amplitude() (py- lablib.devices.AWG.specific.TektronixAFG1000 <i>method</i> ), 442
set_all_values() (py- lablib.core.gui.widgets.container.QFrameContainer <i>method</i> ), 216	set_analog_output_value() (py- lablib.devices.Lakeshore.base.Lakeshore218 <i>method</i> ), 595
set_all_values() (py- lablib.core.gui.widgets.container.QGroupBoxContainer <i>method</i> ), 224	set_analog_output_value() (py- lablib.devices.Lakeshore.base.Lakeshore370 <i>method</i> ), 600
set_all_values() (py- lablib.core.gui.widgets.container.QScrollAreaContainer <i>method</i> ), 230	set_attenuation() (py- lablib.devices.OZOptics.base.DD100 <i>method</i> ), 635
set_all_values() (py- lablib.core.gui.widgets.container.QScrollAreaContainer.QComboBoxWidget <i>method</i> ), 228	set_attribute_value() (py- lablib.devices.Andor.AndorSDK3.AndorSDK3Camera <i>method</i> ), 480
set_all_values() (py- lablib.core.gui.widgets.container.QTabContainer <i>method</i> ), 233	set_attribute_value() (py- lablib.devices.Basler.pylon.BaslerPylonCamera <i>method</i> ), 514
set_all_values() (py- lablib.core.gui.widgets.container.QWidgetContainer <i>method</i> ), 212	set_attribute_value() (py- lablib.devices.DCAM.DCAM.DCAMCamera <i>method</i> ), 514

method), 551

set\_attribute\_value() (py- lablib.devices.IMAQdx.IMAQdx.EthernetIMAQdxCamera method), 588

set\_attribute\_value() (py- lablib.devices.IMAQdx.IMAQdx.IMAQdxCamera method), 579

set\_attribute\_value() (py- lablib.devices.interface.camera.IAttributeCamera method), 837

set\_attribute\_value() (py- lablib.devices.Photometrics.pvcam.PvcamCamera method), 662

set\_attribute\_value() (py- lablib.devices.PhotonFocus.PhotonFocus.IPhotonFocusCamera method), 672

set\_attribute\_value() (py- lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusBitFlowCamera method), 699

set\_attribute\_value() (py- lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusIMAQCamera method), 684

set\_attribute\_value() (py- lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusSiSoCamera method), 692

set\_attribute\_value() (py- lablib.devices.PrincetonInstruments.picam.PicamCamera method), 712

set\_axis() (pylablib.devices.PhysikInstrumente.base.GenericPIController method), 703

set\_axis() (pylablib.devices.PhysikInstrumente.base.PIE516 method), 707

set\_axis\_correction() (py- lablib.devices.Attocube.anc300.ANC300 method), 505

set\_axis\_dir() (py- lablib.devices.SmarAct.scu3d.SCUD method), 746

set\_axis\_parameter() (py- lablib.devices.PhysikInstrumente.base.GenericPIController method), 704

set\_axis\_parameter() (py- lablib.devices.PhysikInstrumente.base.PIE516 method), 708

set\_axis\_parameter() (py- lablib.devices.Trinamic.base.TMCM1110 method), 827

set\_axis\_speed() (py- lablib.devices.Arcus.performax.Performax2EXStage method), 496

set\_axis\_speed() (py- lablib.devices.Arcus.performax.Performax4EXStage method), 491

set\_axis\_speed() (py- lablib.devices.Arcus.performax.PerformaxDMXJSAStage method), 498

set\_baudrate() (py- lablib.devices.Arcus.performax.Performax2EXStage method), 496

set\_baudrate() (py- lablib.devices.Arcus.performax.Performax4EXStage method), 490

set\_baudrate() (py- lablib.devices.Ophir.base.VegaPowerMeter method), 642

set\_binning() (py- lablib.devices.uc480.uc480.UC480Camera method), 869

set\_camera\_alignment() (py- lablib.devices.PCO.SC2.PCOSC2Camera method), 649

set\_black\_level\_offset() (py- lablib.devices.AlliedVision.Bonito.BonitoIMAQCamera method), 462

set\_black\_level\_offset() (py- lablib.devices.AlliedVision.Bonito.IBonitoCamera method), 452

set\_black\_level\_offset() (py- lablib.devices.PhotonFocus.PhotonFocus.IPhotonFocusCamera method), 674

set\_black\_level\_offset() (py- lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusBitFlowCamera method), 699

set\_black\_level\_offset() (py- lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusIMAQCamera method), 684

set\_black\_level\_offset() (py- lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusSiSoCamera method), 692

set\_black\_level\_offset() (py- lablib.devices.AWG.generic.GenericAWG method), 405

set\_black\_level\_offset() (py- lablib.devices.AWG.specific.Agilent33220A method), 418

set\_black\_level\_offset() (py- lablib.devices.AWG.specific.Agilent33500 method), 412

set\_black\_level\_offset() (py- lablib.devices.AWG.specific.InstekAFG2000 method), 430

set\_black\_level\_offset() (py- lablib.devices.AWG.specific.InstekAFG2225 method), 424

set\_black\_level\_offset() (py- lablib.devices.AWG.specific.RigolDG1000 method), 448

set\_burst\_mode() (py- lablib.devices.AWG.specific.RigolDG1000 method), 448



<i>lablib.devices.AWG.specific.RSInstekAFG21000</i> <i>method</i> ), 436	<i>lablib.core.gui.widgets.container.QDialogContainer</i> <i>method</i> ), 220
<i>set_burst_mode()</i> (py- <i>lablib.devices.AWG.specific.TektronixAFG1000</i> <i>method</i> ), 442	<i>set_column_stretch()</i> (py- <i>lablib.core.gui.widgets.container.QFrameContainer</i> <i>method</i> ), 216
<i>set_burst_ncycles()</i> (py- <i>lablib.devices.AWG.generic.GenericAWG</i> <i>method</i> ), 405	<i>set_column_stretch()</i> (py- <i>lablib.core.gui.widgets.container.QGroupBoxContainer</i> <i>method</i> ), 224
<i>set_burst_ncycles()</i> (py- <i>lablib.devices.AWG.specific.Agilent33220A</i> <i>method</i> ), 418	<i>set_column_stretch()</i> (py- <i>lablib.core.gui.widgets.container.QScrollAreaContainer.QContain</i> <i>method</i> ), 228
<i>set_burst_ncycles()</i> (py- <i>lablib.devices.AWG.specific.Agilent33500</i> <i>method</i> ), 412	<i>set_column_stretch()</i> (py- <i>lablib.core.gui.widgets.container.QWidgetContainer</i> <i>method</i> ), 212
<i>set_burst_ncycles()</i> (py- <i>lablib.devices.AWG.specific.InstekAFG2000</i> <i>method</i> ), 430	<i>set_column_stretch()</i> (py- <i>lablib.core.gui.widgets.layout_manager.IQLayoutManagedWidget</i> <i>method</i> ), 239
<i>set_burst_ncycles()</i> (py- <i>lablib.devices.AWG.specific.InstekAFG2225</i> <i>method</i> ), 424	<i>set_column_stretch()</i> (py- <i>lablib.core.gui.widgets.layout_manager.QLayoutManagedWidget</i> <i>method</i> ), 241
<i>set_burst_ncycles()</i> (py- <i>lablib.devices.AWG.specific.RigolDG1000</i> <i>method</i> ), 448	<i>set_column_stretch()</i> (py- <i>lablib.core.gui.widgets.param_table.ParamTable</i> <i>method</i> ), 250
<i>set_burst_ncycles()</i> (py- <i>lablib.devices.AWG.specific.RSInstekAFG21000</i> <i>method</i> ), 436	<i>set_column_stretch()</i> (py- <i>lablib.core.gui.widgets.param_table.StatusTable</i> <i>method</i> ), 259
<i>set_burst_ncycles()</i> (py- <i>lablib.devices.AWG.specific.TektronixAFG1000</i> <i>method</i> ), 442	<i>set_container()</i> (py- <i>lablib.core.dataproc.table_wrap.Array2DWrapper</i> <i>method</i> ), 129
<i>set_by_name()</i> (py- <i>lablib.core.gui.widgets.container.QTabContainer</i> <i>method</i> ), 231	<i>set_control_mode()</i> (py- <i>lablib.devices.Cryomagnetics.base.LM510</i> <i>method</i> ), 544
<i>set_calibration_factor()</i> (py- <i>lablib.devices.Pfeiffer.base.TPG260</i> <i>method</i> ), 655	<i>set_cooler()</i> ( <i>pylablib.devices.Andor.AndorSDK2.AndorSDK2Camera</i> <i>method</i> ), 467
<i>set_camera_id()</i> (py- <i>lablib.devices.uc480.uc480.UC480Camera</i> <i>method</i> ), 867	<i>set_cooler()</i> ( <i>pylablib.devices.Andor.AndorSDK3.AndorSDK3Camera</i> <i>method</i> ), 480
<i>set_channel_power()</i> (py- <i>lablib.devices.Toptica.ibeam.TopticaIBeam</i> <i>method</i> ), 825	<i>set_coupling()</i> (py- <i>lablib.devices.Tektronix.base.DPO2000</i> <i>method</i> ), 767
<i>set_clear_cycles()</i> (py- <i>lablib.devices.Photometrics.pvcam.PvcamCamera</i> <i>set</i> <i>method</i> ), 664	<i>set_coupling()</i> (py- <i>lablib.devices.Tektronix.base.ITektronixScope</i> <i>method</i> ), 751
<i>set_clear_mode()</i> (py- <i>lablib.devices.Photometrics.pvcam.PvcamCamera</i> <i>set</i> <i>method</i> ), 664	<i>set_coupling()</i> (py- <i>lablib.devices.Tektronix.base.TDS2000</i> <i>method</i> ), 760
<i>set_color_mode()</i> (py- <i>lablib.devices.uc480.uc480.UC480Camera</i> <i>method</i> ), 867	<i>set_cursor_order()</i> (py- <i>lablib.core.gui.widgets.edit.NumEdit</i> <i>method</i> ), 235
<i>set_column_stretch()</i> (py- <i>lablib.core.gui.widgets.container.IQWidgetContainer</i> <i>method</i> ), 208	<i>set_curve()</i> ( <i>pylablib.devices.Lakeshore.base.Lakeshore218</i> <i>method</i> ), 594
<i>set_column_stretch()</i> (py-	<i>set_curve_header()</i> (py- <i>lablib.devices.Lakeshore.base.Lakeshore218</i> <i>method</i> ), 594
	<i>set_custom_steps()</i> (py-

<i>lablib.core.gui.widgets.edit.NumEdit</i> method), 235	<i>lablib.core.devio.comm_backend.ICommBackendWrapper</i> method), 161
<i>set_data_format()</i> ( <i>py-lablib.devices.Tektronix.base.DPO2000</i> method), 767	<i>set_device_variable()</i> ( <i>py-lablib.core.devio.interface.IDevice</i> method), 163
<i>set_data_format()</i> ( <i>py-lablib.devices.Tektronix.base.ITektronixScope</i> method), 751	<i>set_device_variable()</i> ( <i>py-lablib.core.devio.SCPIDevice</i> method), 141
<i>set_data_format()</i> ( <i>py-lablib.devices.Tektronix.base.TDS2000</i> method), 760	<i>set_device_variable()</i> ( <i>py-lablib.devices.AlliedVision.Bonito.BonitoIMAQCamera</i> method), 462
<i>set_data_pts_range()</i> ( <i>py-lablib.devices.Tektronix.base.DPO2000</i> method), 767	<i>set_device_variable()</i> ( <i>py-lablib.devices.AlliedVision.Bonito.IBonitoCamera</i> method), 455
<i>set_data_pts_range()</i> ( <i>py-lablib.devices.Tektronix.base.ITektronixScope</i> method), 751	<i>set_device_variable()</i> ( <i>py-lablib.devices.Andor.AndorSDK2.AndorSDK2Camera</i> method), 475
<i>set_data_pts_range()</i> ( <i>py-lablib.devices.Tektronix.base.TDS2000</i> method), 760	<i>set_device_variable()</i> ( <i>py-lablib.devices.Andor.AndorSDK3.AndorSDK3Camera</i> method), 485
<i>set_default_addr()</i> ( <i>py-lablib.devices.Thorlabs.elliptec.ElliptecMotor</i> method), 778	<i>set_device_variable()</i> ( <i>py-lablib.devices.Arcus.performax.GenericPerformaxStage</i> method), 489
<i>set_default_channel()</i> ( <i>py-lablib.devices.HighFinesse.wlm.WLM</i> method), 558	<i>set_device_variable()</i> ( <i>py-lablib.devices.Arcus.performax.Performax2EXStage</i> method), 496
<i>set_default_channel()</i> ( <i>py-lablib.devices.Thorlabs.kinesis.KinesisDevice</i> method), 787	<i>set_device_variable()</i> ( <i>py-lablib.devices.Arcus.performax.Performax4EXStage</i> method), 493
<i>set_default_channel()</i> ( <i>py-lablib.devices.Thorlabs.kinesis.KinesisMotor</i> method), 800	<i>set_device_variable()</i> ( <i>py-lablib.devices.Arcus.performax.PerformaxDMXJSASage</i> method), 500
<i>set_default_channel()</i> ( <i>py-lablib.devices.Thorlabs.kinesis.KinesisPiezoMotor</i> method), 805	<i>set_device_variable()</i> ( <i>py-lablib.devices.Arduino.base.IArduinoDevice</i> method), 502
<i>set_default_channel()</i> ( <i>py-lablib.devices.Thorlabs.kinesis.MFF</i> method), 794	<i>set_device_variable()</i> ( <i>py-lablib.devices.Attocube.anc300.ANC300</i> method), 506
<i>set_defect_correct_mode()</i> ( <i>py-lablib.devices.DCAM.DCAM.DCAMCamera</i> method), 552	<i>set_device_variable()</i> ( <i>py-lablib.devices.Attocube.anc350.ANC350</i> method), 510
<i>set_device_number()</i> ( <i>py-lablib.devices.Arcus.performax.GenericPerformaxStage</i> method), 488	<i>set_device_variable()</i> ( <i>py-lablib.devices.AWG.generic.GenericAWG</i> method), 407
<i>set_device_number()</i> ( <i>py-lablib.devices.Arcus.performax.Performax2EXStage</i> method), 496	<i>set_device_variable()</i> ( <i>py-lablib.devices.AWG.specific.Agilent33220A</i> method), 418
<i>set_device_number()</i> ( <i>py-lablib.devices.Arcus.performax.Performax4EXStage</i> method), 493	<i>set_device_variable()</i> ( <i>py-lablib.devices.AWG.specific.Agilent33500</i> method), 412
<i>set_device_number()</i> ( <i>py-lablib.devices.Arcus.performax.PerformaxDMXJSASage</i> method), 499	<i>set_device_variable()</i> ( <i>py-lablib.devices.AWG.specific.InstekAFG2000</i> method), 430
<i>set_device_variable()</i> ( <i>py-</i>	<i>set_device_variable()</i> ( <i>py-</i>

<code>lablib.devices.AWG.specific.InstekAFG2225</code> <code>method)</code> , 424	<code>lablib.devices.interface.camera.IBinROICamera</code> <code>method)</code> , 859
<code>set_device_variable()</code> (py- <code>lablib.devices.AWG.specific.RigolDG1000</code> <code>method)</code> , 448	<code>set_device_variable()</code> (py- <code>lablib.devices.interface.camera.ICamera</code> <code>method)</code> , 835
<code>set_device_variable()</code> (py- <code>lablib.devices.AWG.specific.RSInstekAFG21000</code> <code>method)</code> , 436	<code>set_device_variable()</code> (py- <code>lablib.devices.interface.camera.IExposureCamera</code> <code>method)</code> , 850
<code>set_device_variable()</code> (py- <code>lablib.devices.AWG.specific.TektronixAFG1000</code> <code>method)</code> , 442	<code>set_device_variable()</code> (py- <code>lablib.devices.interface.camera.IGrabberAttributeCamera</code> <code>method)</code> , 845
<code>set_device_variable()</code> (py- <code>lablib.devices.Basler.pylon.BaslerPylonCamera</code> <code>method)</code> , 519	<code>set_device_variable()</code> (py- <code>lablib.devices.interface.camera.IROICamera</code> <code>method)</code> , 855
<code>set_device_variable()</code> (py- <code>lablib.devices.BitFlow.BitFlow.BitFlowCamera</code> <code>method)</code> , 531	<code>set_device_variable()</code> (py- <code>lablib.devices.interface.stage.IMultiaxisStage</code> <code>method)</code> , 864
<code>set_device_variable()</code> (py- <code>lablib.devices.BitFlow.BitFlow.BitFlowFrameGrabber</code> <code>method)</code> , 526	<code>set_device_variable()</code> (py- <code>lablib.devices.interface.stage.IStage</code> <code>method)</code> , 863
<code>set_device_variable()</code> (py- <code>lablib.devices.Conrad.base.RelayBoard</code> <code>method)</code> , 535	<code>set_device_variable()</code> (py- <code>lablib.devices.KJL.base.KJL300</code> <code>method)</code> , 591
<code>set_device_variable()</code> (py- <code>lablib.devices.Cryocon.base.Cryocon1x</code> <code>method)</code> , 538	<code>set_device_variable()</code> (py- <code>lablib.devices.Lakeshore.base.Lakeshore218</code> <code>method)</code> , 597
<code>set_device_variable()</code> (py- <code>lablib.devices.Cryomagnetics.base.LM500</code> <code>method)</code> , 543	<code>set_device_variable()</code> (py- <code>lablib.devices.Lakeshore.base.Lakeshore370</code> <code>method)</code> , 602
<code>set_device_variable()</code> (py- <code>lablib.devices.Cryomagnetics.base.LM510</code> <code>method)</code> , 546	<code>set_device_variable()</code> (py- <code>lablib.devices.LaserQuantum.base.Finesse</code> <code>method)</code> , 606
<code>set_device_variable()</code> (py- <code>lablib.devices.DCAM.DCAM.DCAMCamera</code> <code>method)</code> , 555	<code>set_device_variable()</code> (py- <code>lablib.devices.Leybold.base.GenericITR</code> <code>method)</code> , 609
<code>set_device_variable()</code> (py- <code>lablib.devices.HighFinesse.wlm.WLM</code> <code>method)</code> , 560	<code>set_device_variable()</code> (py- <code>lablib.devices.Leybold.base.ITR90</code> <code>method)</code> , 610
<code>set_device_variable()</code> (py- <code>lablib.devices.IMAQ.IMAQ.IMAQCamera</code> <code>method)</code> , 574	<code>set_device_variable()</code> (py- <code>lablib.devices.LighthousePhotonics.base.SproutG</code> <code>method)</code> , 613
<code>set_device_variable()</code> (py- <code>lablib.devices.IMAQ.IMAQ.IMAQFrameGrabber</code> <code>method)</code> , 567	<code>set_device_variable()</code> (py- <code>lablib.devices.M2.solstis.Solstis</code> <code>method)</code> , 620
<code>set_device_variable()</code> (py- <code>lablib.devices.IMAQdx.IMAQdx.EthernetIMAQdxCamera</code> <code>method)</code> , 588	<code>set_device_variable()</code> (py- <code>lablib.devices.Newport.picomotor.Picomotor8742</code> <code>method)</code> , 631
<code>set_device_variable()</code> (py- <code>lablib.devices.IMAQdx.IMAQdx.IMAQdxCamera</code> <code>method)</code> , 583	<code>set_device_variable()</code> (py- <code>lablib.devices.NI.daq.NIDAQ</code> <code>method)</code> , 626
<code>set_device_variable()</code> (py- <code>lablib.devices.interface.camera.IAttributeCamera</code> <code>method)</code> , 840	<code>set_device_variable()</code> (py- <code>lablib.devices.Ophir.base.OphirDevice</code> <code>method)</code> , 640
<code>set_device_variable()</code> (py-	<code>set_device_variable()</code> (py- <code>lablib.devices.Ophir.base.VegaPowerMeter</code>

<i>method</i> ), 643	<i>method</i> ), 742
<code>set_device_variable()</code> (py-lablib.devices.OZOptics.base.DD100 <i>method</i> ), 636	<code>set_device_variable()</code> (py-lablib.devices.SmarAct.scu3d.SC3D <i>method</i> ), 747
<code>set_device_variable()</code> (py-lablib.devices.OZOptics.base.EPC04 <i>method</i> ), 638	<code>set_device_variable()</code> (py-lablib.devices.Tektronix.base.DPO2000 <i>method</i> ), 767
<code>set_device_variable()</code> (py-lablib.devices.OZOptics.base.OZOpticsDevice <i>method</i> ), 633	<code>set_device_variable()</code> (py-lablib.devices.Tektronix.base.ITektronixScope <i>method</i> ), 754
<code>set_device_variable()</code> (py-lablib.devices.OZOptics.base.TF100 <i>method</i> ), 634	<code>set_device_variable()</code> (py-lablib.devices.Tektronix.base.TDS2000 <i>method</i> ), 760
<code>set_device_variable()</code> (py-lablib.devices.PCO.SC2.PCOSC2Camera <i>method</i> ), 651	<code>set_device_variable()</code> (py-lablib.devices.Thorlabs.elliptec.ElliptecMotor <i>method</i> ), 781
<code>set_device_variable()</code> (py-lablib.devices.Pfeiffer.base.DPG202 <i>method</i> ), 659	<code>set_device_variable()</code> (py-lablib.devices.Thorlabs.kinesis.BasicKinesisDevice <i>method</i> ), 784
<code>set_device_variable()</code> (py-lablib.devices.Pfeiffer.base.TPG260 <i>method</i> ), 657	<code>set_device_variable()</code> (py-lablib.devices.Thorlabs.kinesis.KinesisDevice <i>method</i> ), 790
<code>set_device_variable()</code> (py-lablib.devices.Photometrics.pvcam.PvcamCamera <i>method</i> ), 668	<code>set_device_variable()</code> (py-lablib.devices.Thorlabs.kinesis.KinesisMotor <i>method</i> ), 800
<code>set_device_variable()</code> (py-lablib.devices.PhotonFocus.PhotonFocus.IPhotonFocusCamera <i>method</i> ), 677	<code>set_device_variable()</code> (py-lablib.devices.Thorlabs.kinesis.KinesisPiezoMotor <i>method</i> ), 805
<code>set_device_variable()</code> (py-lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusBitFluxCamera <i>method</i> ), 699	<code>set_device_variable()</code> (py-lablib.devices.Thorlabs.kinesis.KinesisQuadDetector <i>method</i> ), 809
<code>set_device_variable()</code> (py-lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusIMACamera <i>method</i> ), 684	<code>set_device_variable()</code> (py-lablib.devices.Thorlabs.kinesis.MFF <i>method</i> ), 794
<code>set_device_variable()</code> (py-lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusSiSoCamera <i>method</i> ), 692	<code>set_device_variable()</code> (py-lablib.devices.Thorlabs.serial.FW <i>method</i> ), 815
<code>set_device_variable()</code> (py-lablib.devices.PhysikInstrumente.base.GenericPIController <i>method</i> ), 704	<code>set_device_variable()</code> (py-lablib.devices.Thorlabs.serial.FWv1 <i>method</i> ), 818
<code>set_device_variable()</code> (py-lablib.devices.PhysikInstrumente.base.PIE516 <i>method</i> ), 708	<code>set_device_variable()</code> (py-lablib.devices.Thorlabs.serial.MDT69xA <i>method</i> ), 822
<code>set_device_variable()</code> (py-lablib.devices.PrincetonInstruments.picam.PicamCamera <i>method</i> ), 716	<code>set_device_variable()</code> (py-lablib.devices.Thorlabs.serial.ThorlabsSerialInterface <i>method</i> ), 811
<code>set_device_variable()</code> (py-lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareCamera <i>method</i> ), 731	<code>set_device_variable()</code> (py-lablib.devices.Thorlabs.TLCamera.ThorlabsTLCamera <i>method</i> ), 775
<code>set_device_variable()</code> (py-lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareFrameGrabber <i>method</i> ), 726	<code>set_device_variable()</code> (py-lablib.devices.Toptica.ibeam.TopicalIBeam <i>method</i> ), 826
<code>set_device_variable()</code> (py-lablib.devices.Sirah.Matisse.SirahMatisse <i>method</i> ), 742	<code>set_device_variable()</code> (py-lablib.devices.Trinamic.base.TMCM1110 <i>method</i> ), 844

<i>method</i> ), 829		<i>set_duty_cycle()</i>	(py-
<i>set_device_variable()</i>	(py-	<i>lablib.devices.AWG.specific.Agilent33500</i>	
<i>lablib.devices.uc480.uc480.UC480Camera</i>		<i>method</i> ), 412	
<i>method</i> ), 871		<i>set_duty_cycle()</i>	(py-
<i>set_diffuser()</i>	(py-	<i>lablib.devices.AWG.specific.InstekAFG2000</i>	
<i>lablib.devices.Ophir.base.VegaPowerMeter</i>		<i>method</i> ), 430	
<i>method</i> ), 642		<i>set_duty_cycle()</i>	(py-
<i>set_digital_gain()</i>	(py-	<i>lablib.devices.AWG.specific.InstekAFG2225</i>	
<i>lablib.devices.AlliedVision.Bonito.BonitoIMAQCamera</i>		<i>method</i> ), 424	
<i>method</i> ), 462		<i>set_duty_cycle()</i>	(py-
<i>set_digital_gain()</i>	(py-	<i>lablib.devices.AWG.specific.RigolDG1000</i>	
<i>lablib.devices.AlliedVision.Bonito.IBonitoCamera</i>		<i>method</i> ), 448	
<i>method</i> ), 452		<i>set_duty_cycle()</i>	(py-
<i>set_digital_output()</i>	(py-	<i>lablib.devices.AWG.specific.RSInstekAFG21000</i>	
<i>lablib.devices.Arcus.performax.Performax2EXStage</i>		<i>method</i> ), 436	
<i>method</i> ), 496		<i>set_duty_cycle()</i>	(py-
<i>set_digital_output()</i>	(py-	<i>lablib.devices.AWG.specific.TektronixAFG1000</i>	
<i>lablib.devices.Arcus.performax.Performax4EXStage</i>		<i>method</i> ), 442	
<i>method</i> ), 492		<i>set_edge_trigger_coupling()</i>	(py-
<i>set_digital_output()</i>	(py-	<i>lablib.devices.Tektronix.base.DPO2000</i>	
<i>lablib.devices.Arcus.performax.PerformaxDMXJSASStage</i>		<i>method</i> ), 767	
<i>method</i> ), 500		<i>set_edge_trigger_coupling()</i>	(py-
<i>set_digital_output_register()</i>	(py-	<i>lablib.devices.Tektronix.base.ITektronixScope</i>	
<i>lablib.devices.Arcus.performax.Performax2EXStage</i>		<i>method</i> ), 749	
<i>method</i> ), 496		<i>set_edge_trigger_coupling()</i>	(py-
<i>set_digital_output_register()</i>	(py-	<i>lablib.devices.Tektronix.base.TDS2000</i>	
<i>lablib.devices.Arcus.performax.Performax4EXStage</i>		<i>method</i> ), 760	
<i>method</i> ), 492		<i>set_edge_trigger_slope()</i>	(py-
<i>set_digital_output_register()</i>	(py-	<i>lablib.devices.Tektronix.base.DPO2000</i>	
<i>lablib.devices.Arcus.performax.PerformaxDMXJSASStage</i>		<i>method</i> ), 767	
<i>method</i> ), 500		<i>set_edge_trigger_slope()</i>	(py-
<i>set_digital_outputs()</i>	(py-	<i>lablib.devices.Tektronix.base.ITektronixScope</i>	
<i>lablib.devices.NI.daq.NIDAQ</i> <i>method</i> ), 624		<i>method</i> ), 749	
<i>set_diode_power_lowlevel()</i>	(py-	<i>set_edge_trigger_slope()</i>	(py-
<i>lablib.devices.Sirah.Matisse.SirahMatisse</i>		<i>lablib.devices.Tektronix.base.TDS2000</i>	
<i>method</i> ), 736		<i>method</i> ), 760	
<i>set_direct_index_action()</i>	(py-	<i>set_edge_trigger_source()</i>	(py-
<i>lablib.core.gui.widgets.combo_box.ComboBox</i>		<i>lablib.devices.Tektronix.base.DPO2000</i>	
<i>method</i> ), 198		<i>method</i> ), 767	
<i>set_display_channel()</i>	(py-	<i>set_edge_trigger_source()</i>	(py-
<i>lablib.devices.Pfeiffer.base.TPG260</i> <i>method</i> ),		<i>lablib.devices.Tektronix.base.ITektronixScope</i>	
655		<i>method</i> ), 749	
<i>set_display_resolution()</i>	(py-	<i>set_edge_trigger_source()</i>	(py-
<i>lablib.devices.Pfeiffer.base.TPG260</i> <i>method</i> ),		<i>lablib.devices.Tektronix.base.TDS2000</i>	
655		<i>method</i> ), 760	
<i>set_display_units()</i>	(py-	<i>set_EMCCD_gain()</i>	(py-
<i>lablib.devices.Cryocon.base.Cryocon1x</i>		<i>lablib.devices.Andor.AndorSDK2.AndorSDK2Camera</i>	
<i>method</i> ), 536		<i>method</i> ), 468	
<i>set_duty_cycle()</i>	(py-	<i>set_enabled()</i>	(py-
<i>lablib.devices.AWG.generic.GenericAWG</i>		<i>lablib.core.gui.widgets.param_table.ParamTable</i>	
<i>method</i> ), 404		<i>method</i> ), 248	
<i>set_duty_cycle()</i>	(py-	<i>set_enabled()</i>	(py-
<i>lablib.devices.AWG.specific.Agilent33220A</i>		<i>lablib.core.gui.widgets.param_table.StatusTable</i>	
<i>method</i> ), 418		<i>method</i> ), 259	





<i>lablib.core.gui.widgets.edit.NumEdit</i> method), 235	<i>lablib.devices.interface.camera.IExposureCamera</i> method), 850
<i>set_formatter()</i> ( <i>py-lablib.core.gui.widgets.label.NumLabel</i> method), 237	<i>set_frame_format()</i> ( <i>py-lablib.devices.interface.camera.IGrabberAttributeCamera</i> method), 845
<i>set_frame_delay()</i> ( <i>py-lablib.devices.PCO.SC2.PCOSC2Camera</i> method), 647	<i>set_frame_format()</i> ( <i>py-lablib.devices.interface.camera.IROICamera</i> method), 855
<i>set_frame_format()</i> ( <i>py-lablib.devices.AlliedVision.Bonito.BonitoIMAQCamera</i> method), 462	<i>set_frame_format()</i> ( <i>py-lablib.devices.PCO.SC2.PCOSC2Camera</i> method), 652
<i>set_frame_format()</i> ( <i>py-lablib.devices.AlliedVision.Bonito.IBonitoCamera</i> method), 455	<i>set_frame_format()</i> ( <i>py-lablib.devices.Photometrics.pvcam.PvcamCamera</i> method), 668
<i>set_frame_format()</i> ( <i>py-lablib.devices.Andor.AndorSDK2.AndorSDK2Camera</i> method), 475	<i>set_frame_format()</i> ( <i>py-lablib.devices.PhotonFocus.PhotonFocus.IPhotonFocusCamera</i> method), 677
<i>set_frame_format()</i> ( <i>py-lablib.devices.Andor.AndorSDK3.AndorSDK3Camera</i> method), 485	<i>set_frame_format()</i> ( <i>py-lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusBitFlowCa</i> method), 699
<i>set_frame_format()</i> ( <i>py-lablib.devices.Basler.pylon.BaslerPylonCamera</i> method), 519	<i>set_frame_format()</i> ( <i>py-lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusIMAQCam</i> method), 685
<i>set_frame_format()</i> ( <i>py-lablib.devices.BitFlow.BitFlow.BitFlowCamera</i> method), 531	<i>set_frame_format()</i> ( <i>py-lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusSiSoCamer</i> method), 692
<i>set_frame_format()</i> ( <i>py-lablib.devices.BitFlow.BitFlow.BitFlowFrameGrabber</i> method), 526	<i>set_frame_format()</i> ( <i>py-lablib.devices.PrincetonInstruments.picam.PicamCamera</i> method), 716
<i>set_frame_format()</i> ( <i>py-lablib.devices.DCAM.DCAM.DCAMCamera</i> method), 555	<i>set_frame_format()</i> ( <i>py-lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareCamera</i> method), 731
<i>set_frame_format()</i> ( <i>py-lablib.devices.IMAQ.IMAQ.IMAQCamera</i> method), 574	<i>set_frame_format()</i> ( <i>py-lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareFrameGrabb</i> method), 726
<i>set_frame_format()</i> ( <i>py-lablib.devices.IMAQ.IMAQ.IMAQFrameGrabber</i> method), 567	<i>set_frame_format()</i> ( <i>py-lablib.devices.Thorlabs.TLCamera.ThorlabsTLCamera</i> method), 775
<i>set_frame_format()</i> ( <i>py-lablib.devices.IMAQdx.IMAQdx.EthernetIMAQdxCamera</i> method), 588	<i>set_frame_format()</i> ( <i>py-lablib.devices.uc480.uc480.UC480Camera</i> method), 872
<i>set_frame_format()</i> ( <i>py-lablib.devices.IMAQdx.IMAQdx.IMAQdxCamera</i> method), 583	<i>set_frame_info_format()</i> ( <i>py-lablib.devices.AlliedVision.Bonito.BonitoIMAQCamera</i> method), 463
<i>set_frame_format()</i> ( <i>py-lablib.devices.interface.camera.IAttributeCamera</i> method), 840	<i>set_frame_info_format()</i> ( <i>py-lablib.devices.AlliedVision.Bonito.IBonitoCamera</i> method), 455
<i>set_frame_format()</i> ( <i>py-lablib.devices.interface.camera.IBinROICamera</i> method), 859	<i>set_frame_info_format()</i> ( <i>py-lablib.devices.Andor.AndorSDK2.AndorSDK2Camera</i> method), 475
<i>set_frame_format()</i> ( <i>py-lablib.devices.interface.camera.ICamera</i> method), 833	<i>set_frame_info_format()</i> ( <i>py-lablib.devices.Andor.AndorSDK3.AndorSDK3Camera</i> method), 485
<i>set_frame_format()</i> ( <i>py-</i>	<i>set_frame_info_format()</i> ( <i>py-</i>

<i>lablib.devices.Basler.pylon.BaslerPylonCamera</i> method), 519	<i>lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusIMAQCamera</i> method), 685
<i>set_frame_info_format()</i> (py- <i>lablib.devices.BitFlow.BitFlow.BitFlowCamera</i> method), 531	<i>set_frame_info_format()</i> (py- <i>lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusSiSoCamera</i> method), 692
<i>set_frame_info_format()</i> (py- <i>lablib.devices.BitFlow.BitFlow.BitFlowFrameGrabber</i> method), 526	<i>set_frame_info_format()</i> (py- <i>lablib.devices.PrincetonInstruments.picam.PicamCamera</i> method), 717
<i>set_frame_info_format()</i> (py- <i>lablib.devices.DCAM.DCAM.DCAMCamera</i> method), 555	<i>set_frame_info_format()</i> (py- <i>lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareCamera</i> method), 732
<i>set_frame_info_format()</i> (py- <i>lablib.devices.IMAQ.IMAQ.IMAQCamera</i> method), 574	<i>set_frame_info_format()</i> (py- <i>lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareFrameGrabber</i> method), 726
<i>set_frame_info_format()</i> (py- <i>lablib.devices.IMAQ.IMAQ.IMAQFrameGrabber</i> method), 568	<i>set_frame_info_format()</i> (py- <i>lablib.devices.Thorlabs.TLCamera.ThorlabsTLCamera</i> method), 775
<i>set_frame_info_format()</i> (py- <i>lablib.devices.IMAQdx.IMAQdx.EthernetIMAQdxCamera</i> method), 588	<i>set_frame_info_format()</i> (py- <i>lablib.devices.uc480.uc480.UC480Camera</i> method), 872
<i>set_frame_info_format()</i> (py- <i>lablib.devices.IMAQdx.IMAQdx.IMAQdxCamera</i> method), 583	<i>set_frame_info_period()</i> (py- <i>lablib.devices.AlliedVision.Bonito.BonitoIMAQCamera</i> method), 463
<i>set_frame_info_format()</i> (py- <i>lablib.devices.interface.camera.IAttributeCamera</i> method), 840	<i>set_frame_info_period()</i> (py- <i>lablib.devices.AlliedVision.Bonito.IBonitoCamera</i> method), 455
<i>set_frame_info_format()</i> (py- <i>lablib.devices.interface.camera.IBinROICamera</i> method), 859	<i>set_frame_info_period()</i> (py- <i>lablib.devices.Andor.AndorSDK2.AndorSDK2Camera</i> method), 476
<i>set_frame_info_format()</i> (py- <i>lablib.devices.interface.camera.ICamera</i> method), 833	<i>set_frame_info_period()</i> (py- <i>lablib.devices.Andor.AndorSDK3.AndorSDK3Camera</i> method), 485
<i>set_frame_info_format()</i> (py- <i>lablib.devices.interface.camera.IExposureCamera</i> method), 850	<i>set_frame_info_period()</i> (py- <i>lablib.devices.Basler.pylon.BaslerPylonCamera</i> method), 520
<i>set_frame_info_format()</i> (py- <i>lablib.devices.interface.camera.IGrabberAttributeCamera</i> method), 845	<i>set_frame_info_period()</i> (py- <i>lablib.devices.BitFlow.BitFlow.BitFlowCamera</i> method), 531
<i>set_frame_info_format()</i> (py- <i>lablib.devices.interface.camera.IROICamera</i> method), 855	<i>set_frame_info_period()</i> (py- <i>lablib.devices.BitFlow.BitFlow.BitFlowFrameGrabber</i> method), 526
<i>set_frame_info_format()</i> (py- <i>lablib.devices.PCO.SC2.PCOSC2Camera</i> method), 652	<i>set_frame_info_period()</i> (py- <i>lablib.devices.DCAM.DCAM.DCAMCamera</i> method), 556
<i>set_frame_info_format()</i> (py- <i>lablib.devices.Photometrics.pvcam.PvcamCamera</i> method), 668	<i>set_frame_info_period()</i> (py- <i>lablib.devices.IMAQ.IMAQ.IMAQCamera</i> method), 574
<i>set_frame_info_format()</i> (py- <i>lablib.devices.PhotonFocus.PhotonFocus.IPhotonFocusCamera</i> method), 677	<i>set_frame_info_period()</i> (py- <i>lablib.devices.IMAQ.IMAQ.IMAQFrameGrabber</i> method), 568
<i>set_frame_info_format()</i> (py- <i>lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusBitFlowFrameGrabber</i> method), 700	<i>set_frame_info_period()</i> (py- <i>lablib.devices.IMAQdx.IMAQdx.EthernetIMAQdxCamera</i> method), 588
<i>set_frame_info_format()</i> (py-	<i>set_frame_info_period()</i> (py-



<i>lablib.devices.IMAQdx.IMAQdx.IMAQdxCamera</i> <i>method</i> ), 583		<i>lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusSiSoCamera</i> <i>method</i> ), 693
<i>set_frame_info_period()</i> <i>lablib.devices.interface.camera.IAttributeCamera</i> <i>method</i> ), 841	(py-	<i>set_frame_merge()</i> <i>lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareCamera</i> <i>method</i> ), 732
<i>set_frame_info_period()</i> <i>lablib.devices.interface.camera.IBinROICamera</i> <i>method</i> ), 860	(py-	<i>set_frame_merge()</i> <i>lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareFrameGrabber</i> <i>method</i> ), 722
<i>set_frame_info_period()</i> <i>lablib.devices.interface.camera.ICamera</i> <i>method</i> ), 833	(py-	<i>set_frame_period()</i> <i>lablib.devices.AlliedVision.Bonito.BonitoIMAQCamera</i> <i>method</i> ), 463
<i>set_frame_info_period()</i> <i>lablib.devices.interface.camera.IExposureCamera</i> <i>method</i> ), 850	(py-	<i>set_frame_period()</i> <i>lablib.devices.AlliedVision.Bonito.IBonitoCamera</i> <i>method</i> ), 452
<i>set_frame_info_period()</i> <i>lablib.devices.interface.camera.IGrabberAttributeCamera</i> <i>method</i> ), 846	(py-	<i>set_frame_period()</i> <i>lablib.devices.Andor.AndorSDK2.AndorSDK2Camera</i> <i>method</i> ), 470
<i>set_frame_info_period()</i> <i>lablib.devices.interface.camera.IROICamera</i> <i>method</i> ), 855	(py-	<i>set_frame_period()</i> <i>lablib.devices.Andor.AndorSDK3.AndorSDK3Camera</i> <i>method</i> ), 481
<i>set_frame_info_period()</i> <i>lablib.devices.PCO.SC2.PCOSC2Camera</i> <i>method</i> ), 652	(py-	<i>set_frame_period()</i> <i>lablib.devices.Basler.pylon.BaslerPylonCamera</i> <i>method</i> ), 515
<i>set_frame_info_period()</i> <i>lablib.devices.Photometrics.pvcam.PvcamCamera</i> <i>method</i> ), 668	(py-	<i>set_frame_period()</i> <i>lablib.devices.PCO.SC2.PCOSC2Camera</i> <i>method</i> ), 647
<i>set_frame_info_period()</i> <i>lablib.devices.PhotonFocus.PhotonFocus.IPhotonFocusCamera</i> <i>method</i> ), 677	(py-	<i>set_frame_period()</i> <i>lablib.devices.PhotonFocus.PhotonFocus.IPhotonFocusCamera</i> <i>method</i> ), 673
<i>set_frame_info_period()</i> <i>lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusBitFlowCamera</i> <i>method</i> ), 700	(py-	<i>set_frame_period()</i> <i>lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusBitFlowCamera</i> <i>method</i> ), 700
<i>set_frame_info_period()</i> <i>lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusIMAQCamera</i> <i>method</i> ), 685	(py-	<i>set_frame_period()</i> <i>lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusIMAQCamera</i> <i>method</i> ), 685
<i>set_frame_info_period()</i> <i>lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusSiSoCamera</i> <i>method</i> ), 692	(py-	<i>set_frame_period()</i> <i>lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusSiSoCamera</i> <i>method</i> ), 693
<i>set_frame_info_period()</i> <i>lablib.devices.PrincetonInstruments.picam.PicamCamera</i> <i>method</i> ), 717	(py-	<i>set_frame_period()</i> <i>lablib.devices.uc480.uc480.UC480Camera</i> <i>method</i> ), 867
<i>set_frame_info_period()</i> <i>lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareCamera</i> <i>method</i> ), 732	(py-	<i>set_frameskip_behavior()</i> <i>lablib.devices.uc480.uc480.UC480Camera</i> <i>method</i> ), 868
<i>set_frame_info_period()</i> <i>lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareFrameGrabber</i> <i>method</i> ), 726	(py-	<i>set_frequency()</i> <i>lablib.devices.Attocube.anc300.ANC300</i> <i>method</i> ), 504
<i>set_frame_info_period()</i> <i>lablib.devices.Thorlabs.TLCamera.ThorlabsTLCamera</i> <i>method</i> ), 775	(py-	<i>set_frequency()</i> <i>lablib.devices.Attocube.anc350.ANC350</i> <i>method</i> ), 508
<i>set_frame_info_period()</i> <i>lablib.devices.uc480.uc480.UC480Camera</i> <i>method</i> ), 872	(py-	<i>set_frequency()</i> <i>lablib.devices.AWG.generic.GenericAWG</i> <i>method</i> ), 404
<i>set_frame_merge()</i>	(py-	<i>set_frequency()</i>

<i>lablib.devices.AWG.specific.Agilent33220A</i> <i>method</i> ), 418	<i>lablib.devices.AWG.specific.TektronixAFG1000</i> <i>method</i> ), 442
<i>set_frequency()</i> ( <i>py-</i> <i>lablib.devices.AWG.specific.Agilent33500</i> <i>method</i> ), 412	<i>set_gain_boost()</i> ( <i>py-</i> <i>lablib.devices.uc480.uc480.UC480Camera</i> <i>method</i> ), 868
<i>set_frequency()</i> ( <i>py-</i> <i>lablib.devices.AWG.specific.InstekAFG2000</i> <i>method</i> ), 430	<i>set_gains()</i> ( <i>pylablib.devices.uc480.uc480.UC480Camera</i> <i>method</i> ), 868
<i>set_frequency()</i> ( <i>py-</i> <i>lablib.devices.AWG.specific.InstekAFG2225</i> <i>method</i> ), 424	<i>set_gate_polarity()</i> ( <i>py-</i> <i>lablib.devices.AWG.generic.GenericAWG</i> <i>method</i> ), 405
<i>set_frequency()</i> ( <i>py-</i> <i>lablib.devices.AWG.specific.RigolDG1000</i> <i>method</i> ), 448	<i>set_gate_polarity()</i> ( <i>py-</i> <i>lablib.devices.AWG.specific.Agilent33220A</i> <i>method</i> ), 419
<i>set_frequency()</i> ( <i>py-</i> <i>lablib.devices.AWG.specific.RSInstekAFG21000</i> <i>method</i> ), 436	<i>set_gate_polarity()</i> ( <i>py-</i> <i>lablib.devices.AWG.specific.Agilent33500</i> <i>method</i> ), 413
<i>set_frequency()</i> ( <i>py-</i> <i>lablib.devices.AWG.specific.TektronixAFG1000</i> <i>method</i> ), 442	<i>set_gate_polarity()</i> ( <i>py-</i> <i>lablib.devices.AWG.specific.InstekAFG2000</i> <i>method</i> ), 430
<i>set_frequency()</i> ( <i>py-</i> <i>lablib.devices.OZOptics.base.EPC04</i> <i>method</i> ), 637	<i>set_gate_polarity()</i> ( <i>py-</i> <i>lablib.devices.AWG.specific.InstekAFG2225</i> <i>method</i> ), 424
<i>set_frequency()</i> ( <i>py-</i> <i>lablib.devices.Thorlabs.elliptec.ElliptecMotor</i> <i>method</i> ), 779	<i>set_gate_polarity()</i> ( <i>py-</i> <i>lablib.devices.AWG.specific.RigolDG1000</i> <i>method</i> ), 448
<i>set_func_variable()</i> ( <i>py-</i> <i>lablib.core.thread.controller.QTaskThread</i> <i>method</i> ), 311	<i>set_gate_polarity()</i> ( <i>py-</i> <i>lablib.devices.AWG.specific.RSInstekAFG21000</i> <i>method</i> ), 436
<i>set_func_variable()</i> ( <i>py-</i> <i>lablib.core.thread.controller.QThreadController</i> <i>method</i> ), 298	<i>set_gate_polarity()</i> ( <i>py-</i> <i>lablib.devices.AWG.specific.TektronixAFG1000</i> <i>method</i> ), 442
<i>set_function()</i> ( <i>py-</i> <i>lablib.devices.AWG.generic.GenericAWG</i> <i>method</i> ), 403	<i>set_general_output()</i> ( <i>py-</i> <i>lablib.devices.Trinamic.base.TMCM1110</i> <i>method</i> ), 827
<i>set_function()</i> ( <i>py-</i> <i>lablib.devices.AWG.specific.Agilent33220A</i> <i>method</i> ), 418	<i>set_global_parameter()</i> ( <i>py-</i> <i>lablib.devices.Trinamic.base.TMCM1110</i> <i>method</i> ), 827
<i>set_function()</i> ( <i>py-</i> <i>lablib.devices.AWG.specific.Agilent33500</i> <i>method</i> ), 412	<i>set_global_speed()</i> ( <i>py-</i> <i>lablib.devices.Arcus.performax.Performax2EXStage</i> <i>method</i> ), 496
<i>set_function()</i> ( <i>py-</i> <i>lablib.devices.AWG.specific.InstekAFG2000</i> <i>method</i> ), 430	<i>set_global_speed()</i> ( <i>py-</i> <i>lablib.devices.Arcus.performax.Performax4EXStage</i> <i>method</i> ), 491
<i>set_function()</i> ( <i>py-</i> <i>lablib.devices.AWG.specific.InstekAFG2225</i> <i>method</i> ), 424	<i>set_grabber_attribute_value()</i> ( <i>py-</i> <i>lablib.devices.AlliedVision.Bonito.BonitoIMAQCamera</i> <i>method</i> ), 463
<i>set_function()</i> ( <i>py-</i> <i>lablib.devices.AWG.specific.RigolDG1000</i> <i>method</i> ), 448	<i>set_grabber_attribute_value()</i> ( <i>py-</i> <i>lablib.devices.IMAQ.IMAQ.IMAQCamera</i> <i>method</i> ), 574
<i>set_function()</i> ( <i>py-</i> <i>lablib.devices.AWG.specific.RSInstekAFG21000</i> <i>method</i> ), 436	<i>set_grabber_attribute_value()</i> ( <i>py-</i> <i>lablib.devices.IMAQ.IMAQ.IMAQFrameGrabber</i> <i>method</i> ), 562
<i>set_function()</i> ( <i>py-</i>	<i>set_grabber_attribute_value()</i> ( <i>py-</i> <i>lablib.devices.interface.camera.IGrabberAttributeCamera</i>

<a href="#">method</a> ), 842	<a href="#">method</a> ), 779
<a href="#">set_grabber_attribute_value()</a> ( <a href="#">py-lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusIMAQCamera</a> ), 685	<a href="#">set_horizontal_offset()</a> ( <a href="#">py-lablib.devices.Tektronix.base.DPO2000</a> ), 767
<a href="#">set_grabber_attribute_value()</a> ( <a href="#">py-lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusSiSoCamera</a> ), 693	<a href="#">set_horizontal_offset()</a> ( <a href="#">py-lablib.devices.Tektronix.base.ITektronixScope</a> ), 750
<a href="#">set_grabber_attribute_value()</a> ( <a href="#">py-lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareCamera</a> ), 732	<a href="#">set_horizontal_offset()</a> ( <a href="#">py-lablib.devices.Tektronix.base.TDS2000</a> ), 760
<a href="#">set_grabber_attribute_value()</a> ( <a href="#">py-lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareFrameGrabber</a> ), 726	<a href="#">set_horizontal_span()</a> ( <a href="#">py-lablib.devices.Tektronix.base.DPO2000</a> ), 767
<a href="#">set_grabber_roi()</a> ( <a href="#">py-lablib.devices.AlliedVision.Bonito.BonitoIMAQCamera</a> ), 463	<a href="#">set_horizontal_span()</a> ( <a href="#">py-lablib.devices.Tektronix.base.ITektronixScope</a> ), 750
<a href="#">set_grabber_roi()</a> ( <a href="#">py-lablib.devices.BitFlow.BitFlow.BitFlowCamera</a> ), 531	<a href="#">set_horizontal_span()</a> ( <a href="#">py-lablib.devices.Tektronix.base.TDS2000</a> ), 760
<a href="#">set_grabber_roi()</a> ( <a href="#">py-lablib.devices.BitFlow.BitFlow.BitFlowFrameGrabber</a> ), 522	<a href="#">set_image_indexing()</a> ( <a href="#">py-lablib.devices.AlliedVision.Bonito.BonitoIMAQCamera</a> ), 463
<a href="#">set_grabber_roi()</a> ( <a href="#">py-lablib.devices.IMAQ.IMAQ.IMAQCamera</a> ), 574	<a href="#">set_image_indexing()</a> ( <a href="#">py-lablib.devices.AlliedVision.Bonito.IBonitoCamera</a> ), 455
<a href="#">set_grabber_roi()</a> ( <a href="#">py-lablib.devices.IMAQ.IMAQ.IMAQFrameGrabber</a> ), 562	<a href="#">set_image_indexing()</a> ( <a href="#">py-lablib.devices.Andor.AndorSDK2.AndorSDK2Camera</a> ), 476
<a href="#">set_grabber_roi()</a> ( <a href="#">py-lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusBitFlowCamera</a> ), 700	<a href="#">set_image_indexing()</a> ( <a href="#">py-lablib.devices.Andor.AndorSDK3.AndorSDK3Camera</a> ), 486
<a href="#">set_grabber_roi()</a> ( <a href="#">py-lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusIMAQCamera</a> ), 685	<a href="#">set_image_indexing()</a> ( <a href="#">py-lablib.devices.Basler.pylon.BaslerPylonCamera</a> ), 520
<a href="#">set_grabber_roi()</a> ( <a href="#">py-lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusSiSoCamera</a> ), 693	<a href="#">set_image_indexing()</a> ( <a href="#">py-lablib.devices.BitFlow.BitFlow.BitFlowCamera</a> ), 531
<a href="#">set_grabber_roi()</a> ( <a href="#">py-lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareCamera</a> ), 732	<a href="#">set_image_indexing()</a> ( <a href="#">py-lablib.devices.BitFlow.BitFlow.BitFlowFrameGrabber</a> ), 526
<a href="#">set_grabber_roi()</a> ( <a href="#">py-lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareFrameGrabber</a> ), 722	<a href="#">set_image_indexing()</a> ( <a href="#">py-lablib.devices.DCAM.DCAM.DCAMCamera</a> ), 556
<a href="#">set_hardware_id()</a> ( <a href="#">py-lablib.devices.Attocube.anc350.ANC350</a> ), 507	<a href="#">set_image_indexing()</a> ( <a href="#">py-lablib.devices.IMAQ.IMAQ.IMAQCamera</a> ), 575
<a href="#">set_high_level()</a> ( <a href="#">py-lablib.devices.Cryomagnetics.base.LM500</a> ), 541	<a href="#">set_image_indexing()</a> ( <a href="#">py-lablib.devices.IMAQ.IMAQ.IMAQFrameGrabber</a> ), 568
<a href="#">set_high_level()</a> ( <a href="#">py-lablib.devices.Cryomagnetics.base.LM510</a> ), 546	<a href="#">set_image_indexing()</a> ( <a href="#">py-lablib.devices.IMAQdx.IMAQdx.EthernetIMAQdxCamera</a> ), 588
<a href="#">set_home_offset()</a> ( <a href="#">py-lablib.devices.Thorlabs.elliptec.ElliptecMotor</a> ), 507	<a href="#">set_image_indexing()</a> ( <a href="#">py-lablib.devices.IMAQdx.IMAQdx.IMAQdxCamera</a> ), 588

<i>method</i> ), 583	<i>method</i> ), 198
set_image_indexing() lablib.devices.interface.camera.IAttributeCamera <i>method</i> ), 841	(py- set_indicator() lablib.core.gui.value_handling.GUIValues <i>method</i> ), 280
set_image_indexing() lablib.devices.interface.camera.IBinROICamera <i>method</i> ), 860	(py- set_indicator() lablib.core.gui.widgets.container.IQContainer <i>method</i> ), 202
set_image_indexing() lablib.devices.interface.camera.ICamera <i>method</i> ), 832	(py- set_indicator() lablib.core.gui.widgets.container.IQWidgetContainer <i>method</i> ), 208
set_image_indexing() lablib.devices.interface.camera.IExposureCamera <i>method</i> ), 850	(py- set_indicator() lablib.core.gui.widgets.container.QContainer <i>method</i> ), 204
set_image_indexing() lablib.devices.interface.camera.IGrabberAttributeCamera <i>method</i> ), 846	(py- set_indicator() lablib.core.gui.widgets.container.QDialogContainer <i>method</i> ), 220
set_image_indexing() lablib.devices.interface.camera.IROICamera <i>method</i> ), 855	(py- set_indicator() lablib.core.gui.widgets.container.QFrameContainer <i>method</i> ), 216
set_image_indexing() lablib.devices.PCO.SC2.PCOSC2Camera <i>method</i> ), 652	(py- set_indicator() lablib.core.gui.widgets.container.QGroupBoxContainer <i>method</i> ), 224
set_image_indexing() lablib.devices.Photometrics.pvcam.PvcamCamera <i>method</i> ), 668	(py- set_indicator() lablib.core.gui.widgets.container.QScrollAreaContainer <i>method</i> ), 230
set_image_indexing() lablib.devices.PhotonFocus.PhotonFocus.IPhotonFocusCamera <i>method</i> ), 677	(py- set_indicator() lablib.core.gui.widgets.container.QScrollAreaContainer.QContainer <i>method</i> ), 228
set_image_indexing() lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusBitFlowCamera <i>method</i> ), 700	(py- set_indicator() lablib.core.gui.widgets.container.QTabContainer <i>method</i> ), 233
set_image_indexing() lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusIMACamera <i>method</i> ), 685	(py- set_indicator() lablib.core.gui.widgets.container.QWidgetContainer <i>method</i> ), 212
set_image_indexing() lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusSiSoCamera <i>method</i> ), 693	(py- set_indicator() lablib.core.gui.widgets.param_table.ParamTable <i>method</i> ), 248
set_image_indexing() lablib.devices.PrincetonInstruments.picam.PicamCamera <i>method</i> ), 717	(py- set_indicator() lablib.core.gui.widgets.param_table.StatusTable <i>method</i> ), 260
set_image_indexing() lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareCamera <i>method</i> ), 732	(py- set_interval() lablib.devices.Cryomagnetics.base.LM500 <i>method</i> ), 540
set_image_indexing() lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareFrameGrabber <i>method</i> ), 727	(py- set_interval() lablib.devices.Cryomagnetics.base.LM510 <i>method</i> ), 546
set_image_indexing() lablib.devices.Thorlabs.TLCamera.ThorlabsTLCamera <i>method</i> ), 775	(py- set_limiter() lablib.core.gui.widgets.edit.NumEdit <i>method</i> ), 234
set_image_indexing() lablib.devices.uc480.uc480.UC480Camera <i>method</i> ), 872	(py- set_limiter() lablib.core.gui.widgets.label.NumLabel <i>method</i> ), 237
set_index_values() lablib.core.gui.widgets.combo_box.ComboBox	(py- set_load() (pylablib.devices.AWG.generic.GenericAWG <i>method</i> ), 403

<code>set_load()</code> ( <code>pylablib.devices.AWG.specific.Agilent33220A</code> method), 419	<code>lablib.devices.PCO.SC2.PCOSC2Camera</code> method), 648
<code>set_load()</code> ( <code>pylablib.devices.AWG.specific.Agilent33500</code> method), 413	<code>set_offset()</code> ( <code>pylablib.devices.Attocube.anc300.ANC300</code> method), 504
<code>set_load()</code> ( <code>pylablib.devices.AWG.specific.InstekAFG2000</code> method), 430	<code>set_offset()</code> ( <code>pylablib.devices.Attocube.anc350.ANC350</code> method), 508
<code>set_load()</code> ( <code>pylablib.devices.AWG.specific.InstekAFG2225</code> method), 425	<code>set_offset()</code> ( <code>pylablib.devices.AWG.generic.GenericAWG</code> method), 404
<code>set_load()</code> ( <code>pylablib.devices.AWG.specific.RigolDG1000</code> method), 448	<code>set_offset()</code> ( <code>pylablib.devices.AWG.specific.Agilent33220A</code> method), 419
<code>set_load()</code> ( <code>pylablib.devices.AWG.specific.RSInstekAFG21000</code> method), 436	<code>set_offset()</code> ( <code>pylablib.devices.AWG.specific.Agilent33500</code> method), 413
<code>set_load()</code> ( <code>pylablib.devices.AWG.specific.TektronixAFG1000</code> method), 442	<code>set_offset()</code> ( <code>pylablib.devices.AWG.specific.InstekAFG2000</code> method), 430
<code>set_low_level()</code> ( <code>pylablib.devices.Cryomagnetics.base.LM500</code> method), 541	<code>set_offset()</code> ( <code>pylablib.devices.AWG.specific.InstekAFG2225</code> method), 420
<code>set_low_level()</code> ( <code>pylablib.devices.Cryomagnetics.base.LM510</code> method), 546	<code>set_offset()</code> ( <code>pylablib.devices.AWG.specific.RigolDG1000</code> method), 448
<code>set_manual_output()</code> ( <code>pylablib.devices.Thorlabs.kinesis.KinesisQuadDetector</code> method), 806	<code>set_offset()</code> ( <code>pylablib.devices.AWG.specific.RSInstekAFG21000</code> method), 436
<code>set_measurement_filter()</code> ( <code>pylablib.devices.Pfeiffer.base.TPG260</code> method), 655	<code>set_offset()</code> ( <code>pylablib.devices.AWG.specific.TektronixAFG1000</code> method), 442
<code>set_measurement_interval()</code> ( <code>pylablib.devices.HighFinesse.wlm.WLM</code> method), 559	<code>set_operation_mode()</code> ( <code>pylablib.devices.Thorlabs.kinesis.KinesisQuadDetector</code> method), 806
<code>set_metadata_mode()</code> ( <code>pylablib.devices.PCO.SC2.PCOSC2Camera</code> method), 649	<code>set_options()</code> ( <code>pylablib.core.gui.widgets.combo_box.ComboBox</code> method), 199
<code>set_microstep_resolution()</code> ( <code>pylablib.devices.Trinamic.base.TMCM1110</code> method), 828	<code>set_options()</code> ( <code>pylablib.core.gui.widgets.label.EnumLabel</code> method), 236
<code>set_mode()</code> ( <code>pylablib.devices.Attocube.anc300.ANC300</code> method), 503	<code>set_out_aux_port()</code> ( <code>pylablib.devices.Andor.AndorSDK2.AndorSDK2Camera</code> method), 469
<code>set_mode()</code> ( <code>pylablib.devices.Cryomagnetics.base.LM500</code> method), 540	<code>set_out_of_range()</code> ( <code>pylablib.core.gui.widgets.combo_box.ComboBox</code> method), 198
<code>set_mode()</code> ( <code>pylablib.devices.Cryomagnetics.base.LM510</code> method), 546	<code>set_out_of_range()</code> ( <code>pylablib.core.gui.widgets.label.EnumLabel</code> method), 236
<code>set_mode()</code> ( <code>pylablib.devices.OZOptics.base.EPC04</code> method), 637	<code>set_output_mode()</code> ( <code>pylablib.devices.LighthousePhotonics.base.SproutG</code> method), 612
<code>set_mode_parameters()</code> ( <code>pylablib.devices.BitFlow.BitFlow.CameraFileEditor</code> method), 532	<code>set_output_parameters()</code> ( <code>pylablib.devices.Thorlabs.kinesis.KinesisQuadDetector</code> method), 807
<code>set_motor_type()</code> ( <code>pylablib.devices.Newport.picomotor.Picomotor8742</code> method), 629	<code>set_output_polarity()</code> ( <code>pylablib.devices.AWG.generic.GenericAWG</code> method), 403
<code>set_names()</code> ( <code>pylablib.core.dataproc.table_wrap.Array2DWrapper.ColumnAccessor</code> method), 130	<code>set_output_polarity()</code> ( <code>pylablib.devices.AWG.specific.Agilent33220A</code> method), 419
<code>set_names()</code> ( <code>pylablib.core.dataproc.table_wrap.DataFrame2DWrapper.ColumnAccessor</code> method), 132	<code>set_output_polarity()</code> ( <code>pylablib.devices.AWG.specific.Agilent33500</code> method), 413
<code>set_noise_filter_mode()</code> ( <code>pylablib.devices.AWG.specific.Agilent33220A</code> method), 419	



<i>method</i> ), 413		<i>method</i> ), 413	
<code>set_output_polarity()</code> ( <i>pylablib.devices.AWG.specific.InstekAFG2000</i> <i>method</i> ), 431		<code>set_output_trigger_slope()</code> ( <i>pylablib.devices.AWG.specific.InstekAFG2000</i> <i>method</i> ), 431	
<code>set_output_polarity()</code> ( <i>pylablib.devices.AWG.specific.InstekAFG2225</i> <i>method</i> ), 425		<code>set_output_trigger_slope()</code> ( <i>pylablib.devices.AWG.specific.InstekAFG2225</i> <i>method</i> ), 425	
<code>set_output_polarity()</code> ( <i>pylablib.devices.AWG.specific.RigolDG1000</i> <i>method</i> ), 448		<code>set_output_trigger_slope()</code> ( <i>pylablib.devices.AWG.specific.RigolDG1000</i> <i>method</i> ), 448	
<code>set_output_polarity()</code> ( <i>pylablib.devices.AWG.specific.RSInstekAFG21000</i> <i>method</i> ), 436		<code>set_output_trigger_slope()</code> ( <i>pylablib.devices.AWG.specific.RSInstekAFG21000</i> <i>method</i> ), 437	
<code>set_output_polarity()</code> ( <i>pylablib.devices.AWG.specific.TektronixAFG1000</i> <i>method</i> ), 442		<code>set_output_trigger_slope()</code> ( <i>pylablib.devices.AWG.specific.TektronixAFG1000</i> <i>method</i> ), 443	
<code>set_output_power()</code> ( <i>pylablib.devices.LaserQuantum.base.Finesse</i> <i>method</i> ), 605		<code>set_overflow_behavior()</code> ( <i>pylablib.devices.Andor.AndorSDK3.AndorSDK3Camera</i> <i>method</i> ), 482	
<code>set_output_power()</code> ( <i>pylablib.devices.LighthousePhotonics.base.SproutG</i> <i>method</i> ), 612		<code>set_pcount()</code> ( <i>pylablib.devices.Thorlabs.serial.FW</i> <i>method</i> ), 813	
<code>set_output_range()</code> ( <i>pylablib.devices.AWG.generic.GenericAWG</i> <i>method</i> ), 404		<code>set_phase()</code> ( <i>pylablib.devices.AWG.generic.GenericAWG</i> <i>method</i> ), 404	
<code>set_output_range()</code> ( <i>pylablib.devices.AWG.specific.Agilent33220A</i> <i>method</i> ), 419		<code>set_phase()</code> ( <i>pylablib.devices.AWG.specific.Agilent33220A</i> <i>method</i> ), 419	
<code>set_output_range()</code> ( <i>pylablib.devices.AWG.specific.Agilent33500</i> <i>method</i> ), 413		<code>set_phase()</code> ( <i>pylablib.devices.AWG.specific.Agilent33500</i> <i>method</i> ), 413	
<code>set_output_range()</code> ( <i>pylablib.devices.AWG.specific.Agilent33500</i> <i>method</i> ), 413		<code>set_phase()</code> ( <i>pylablib.devices.AWG.specific.InstekAFG2000</i> <i>method</i> ), 431	
<code>set_output_range()</code> ( <i>pylablib.devices.AWG.specific.InstekAFG2000</i> <i>method</i> ), 431		<code>set_phase()</code> ( <i>pylablib.devices.AWG.specific.InstekAFG2225</i> <i>method</i> ), 425	
<code>set_output_range()</code> ( <i>pylablib.devices.AWG.specific.InstekAFG2225</i> <i>method</i> ), 425		<code>set_phase()</code> ( <i>pylablib.devices.AWG.specific.RigolDG1000</i> <i>method</i> ), 448	
<code>set_output_range()</code> ( <i>pylablib.devices.AWG.specific.RigolDG1000</i> <i>method</i> ), 448		<code>set_phase()</code> ( <i>pylablib.devices.AWG.specific.RSInstekAFG21000</i> <i>method</i> ), 437	
<code>set_output_range()</code> ( <i>pylablib.devices.AWG.specific.RSInstekAFG21000</i> <i>method</i> ), 436		<code>set_phase()</code> ( <i>pylablib.devices.AWG.specific.TektronixAFG1000</i> <i>method</i> ), 443	
<code>set_output_range()</code> ( <i>pylablib.devices.AWG.specific.TektronixAFG1000</i> <i>method</i> ), 442		<code>set_pid_parameters()</code> ( <i>pylablib.devices.Thorlabs.kinesis.KinesisQuadDetector</i> <i>method</i> ), 806	
<code>set_output_trigger_slope()</code> ( <i>pylablib.devices.AWG.generic.GenericAWG</i> <i>method</i> ), 407		<code>set_piezoet_ctl_status()</code> ( <i>pylablib.devices.Sirah.Matisse.SirahMatisse</i> <i>method</i> ), 738	
<code>set_output_trigger_slope()</code> ( <i>pylablib.devices.AWG.specific.Agilent33220A</i> <i>method</i> ), 419		<code>set_piezoet_drive_params()</code> ( <i>pylablib.devices.Sirah.Matisse.SirahMatisse</i> <i>method</i> ), 738	
<code>set_output_trigger_slope()</code> ( <i>pylablib.devices.AWG.specific.Agilent33500</i> <i>method</i> ), 419		<code>set_piezoet_feedback_params()</code> ( <i>pylablib.devices.Sirah.Matisse.SirahMatisse</i> <i>method</i> ), 738	
		<code>set_piezoet_feedforward_params()</code> ( <i>pylablib.devices.Sirah.Matisse.SirahMatisse</i> <i>method</i> ), 738	
		<code>set_piezoet_position()</code> ( <i>pylablib.devices.Sirah.Matisse.SirahMatisse</i> <i>method</i> ), 738	

`method`), 738  
`set_pixel_rate()` (py-  
`lablib.devices.PCO.SC2.PCOSC2Camera`  
`method`), 647  
`set_pixel_rate()` (py-  
`lablib.devices.uc480.uc480.UC480Camera`  
`method`), 867  
`set_points_number()` (py-  
`lablib.devices.Tektronix.base.DPO2000`  
`method`), 767  
`set_points_number()` (py-  
`lablib.devices.Tektronix.base.ITektronixScope`  
`method`), 751  
`set_points_number()` (py-  
`lablib.devices.Tektronix.base.TDS2000`  
`method`), 760  
`set_position()` (py-  
`lablib.devices.Thorlabs.serial.FW` `method`),  
813  
`set_position()` (py-  
`lablib.devices.Thorlabs.serial.FWv1` `method`),  
816  
`set_position_lower_limit()` (py-  
`lablib.devices.PhysikInstrumente.base.PIE516`  
`method`), 706  
`set_position_reference()` (py-  
`lablib.devices.Arcus.performax.Performax2EXStage`  
`method`), 497  
`set_position_reference()` (py-  
`lablib.devices.Arcus.performax.Performax4EXStage`  
`method`), 490  
`set_position_reference()` (py-  
`lablib.devices.Arcus.performax.PerformaxDMXJSASStage`  
`method`), 497  
`set_position_reference()` (py-  
`lablib.devices.Newport.picomotor.Picomotor8742`  
`method`), 629  
`set_position_reference()` (py-  
`lablib.devices.Thorlabs.kinesis.KinesisMotor`  
`method`), 796  
`set_position_reference()` (py-  
`lablib.devices.Thorlabs.kinesis.KinesisPiezoMotor`  
`method`), 801  
`set_position_reference()` (py-  
`lablib.devices.Trinamic.base.TMCM1110`  
`method`), 827  
`set_position_upper_limit()` (py-  
`lablib.devices.PhysikInstrumente.base.PIE516`  
`method`), 706  
`set_precision()` (py-  
`lablib.devices.Attocube.anc350.ANC350`  
`method`), 508  
`set_precision_mode()` (py-  
`lablib.devices.HighFinesse.wlm.WLM` `method`),  
559  
`set_probe_attenuation()` (py-  
`lablib.devices.Tektronix.base.DPO2000`  
`method`), 767  
`set_probe_attenuation()` (py-  
`lablib.devices.Tektronix.base.ITektronixScope`  
`method`), 751  
`set_probe_attenuation()` (py-  
`lablib.devices.Tektronix.base.TDS2000`  
`method`), 760  
`set_props()` (in module `pylablib.core.utils.general`),  
374  
`set_pulse_mode()` (py-  
`lablib.devices.HighFinesse.wlm.WLM` `method`),  
559  
`set_pulse_width()` (py-  
`lablib.devices.AWG.generic.GenericAWG`  
`method`), 404  
`set_pulse_width()` (py-  
`lablib.devices.AWG.specific.Agilent33220A`  
`method`), 419  
`set_pulse_width()` (py-  
`lablib.devices.AWG.specific.Agilent33500`  
`method`), 413  
`set_pulse_width()` (py-  
`lablib.devices.AWG.specific.InstekAFG2000`  
`method`), 431  
`set_pulse_width()` (py-  
`lablib.devices.AWG.specific.InstekAFG2225`  
`method`), 425  
`set_pulse_width()` (py-  
`lablib.devices.AWG.specific.RigolDG1000`  
`method`), 449  
`set_pulse_width()` (py-  
`lablib.devices.AWG.specific.RSInstekAFG21000`  
`method`), 437  
`set_pulse_width()` (py-  
`lablib.devices.AWG.specific.TektronixAFG1000`  
`method`), 438  
`set_ramp_symmetry()` (py-  
`lablib.devices.AWG.generic.GenericAWG`  
`method`), 404  
`set_ramp_symmetry()` (py-  
`lablib.devices.AWG.specific.Agilent33220A`  
`method`), 419  
`set_ramp_symmetry()` (py-  
`lablib.devices.AWG.specific.Agilent33500`  
`method`), 413  
`set_ramp_symmetry()` (py-  
`lablib.devices.AWG.specific.InstekAFG2000`  
`method`), 431  
`set_ramp_symmetry()` (py-  
`lablib.devices.AWG.specific.InstekAFG2225`  
`method`), 425

<code>set_ramp_symmetry()</code> <i>pylablib.devices.AWG.specific.RigolDG1000</i> <i>method</i> ), 449	<code>set_roi()</code> ( <i>pylablib.devices.IMAQdx.IMAQdx.EthernetIMAQdxCamera</i> <i>method</i> ), 588
<code>set_ramp_symmetry()</code> <i>pylablib.devices.AWG.specific.RSInstekAFG21000</i> <i>method</i> ), 437	<code>set_roi()</code> ( <i>pylablib.devices.IMAQdx.IMAQdx.IMAQdxCamera</i> <i>method</i> ), 579
<code>set_ramp_symmetry()</code> <i>pylablib.devices.AWG.specific.TektronixAFG1000</i> <i>method</i> ), 443	<code>set_roi()</code> ( <i>pylablib.devices.interface.camera.IBinROICamera</i> <i>method</i> ), 856
<code>set_range_idx()</code> <i>pylablib.devices.Ophir.base.VegaPowerMeter</i> <i>method</i> ), 642	<code>set_roi()</code> ( <i>pylablib.devices.interface.camera.IROICamera</i> <i>method</i> ), 852
<code>set_read_mode()</code> <i>pylablib.devices.Andor.AndorSDK2.AndorSDK2Camera</i> <i>method</i> ), 471	<code>set_roi()</code> ( <i>pylablib.devices.PCO.SC2.PCOSC2Camera</i> <i>method</i> ), 648
<code>set_read_mode()</code> <i>pylablib.devices.HighFinesse.wlm.WLM</i> <i>method</i> ), 558	<code>set_roi()</code> ( <i>pylablib.devices.Photometrics.pvcam.PvcamCamera</i> <i>method</i> ), 664
<code>set_readout_mode()</code> <i>pylablib.devices.Photometrics.pvcam.PvcamCamera</i> <i>method</i> ), 663	<code>set_roi()</code> ( <i>pylablib.devices.PhotonFocus.PhotonFocus.IPhotonFocusCa</i> <i>method</i> ), 673
<code>set_readout_speed()</code> <i>pylablib.devices.DCAM.DCAM.DCAMCamera</i> <i>method</i> ), 551	<code>set_roi()</code> ( <i>pylablib.devices.PhotonFocus.PhotonFocus.PhotonFocusBit</i> <i>method</i> ), 700
<code>set_refcell_position()</code> <i>pylablib.devices.Sirah.Matisse.SirahMatisse</i> <i>method</i> ), 739	<code>set_roi()</code> ( <i>pylablib.devices.PhotonFocus.PhotonFocus.PhotonFocusIMA</i> <i>method</i> ), 685
<code>set_refcell_waveform_params()</code> <i>pylablib.devices.Sirah.Matisse.SirahMatisse</i> <i>method</i> ), 739	<code>set_roi()</code> ( <i>pylablib.devices.PhotonFocus.PhotonFocus.PhotonFocusSiS</i> <i>method</i> ), 693
<code>set_relay()</code> ( <i>pylablib.devices.Conrad.base.RelayBoard</i> <i>method</i> ), 534	<code>set_roi()</code> ( <i>pylablib.devices.PrincetonInstruments.picam.PicamCamera</i> <i>method</i> ), 713
<code>set_relay_setpoints()</code> <i>pylablib.devices.KJL.base.KJL300</i> <i>method</i> ), 590	<code>set_roi()</code> ( <i>pylablib.devices.SiliconSoftware.fgrab.SiliconSoftwareCame</i> <i>method</i> ), 732
<code>set_roi()</code> ( <i>pylablib.devices.AlliedVision.Bonito.BonitoIMAQCamera</i> <i>method</i> ), 463	<code>set_roi()</code> ( <i>pylablib.devices.SiliconSoftware.fgrab.SiliconSoftwareFrame</i> <i>method</i> ), 722
<code>set_roi()</code> ( <i>pylablib.devices.AlliedVision.Bonito.IBonitoCamera</i> <i>method</i> ), 451	<code>set_roi()</code> ( <i>pylablib.devices.Thorlabs.TLCCamera.ThorlabsTLCCamera</i> <i>method</i> ), 772
<code>set_roi()</code> ( <i>pylablib.devices.Andor.AndorSDK2.AndorSDK2Camera</i> <i>method</i> ), 472	<code>set_roi()</code> ( <i>pylablib.devices.uc480.uc480.UC480Camera</i> <i>method</i> ), 869
<code>set_roi()</code> ( <i>pylablib.devices.Andor.AndorSDK3.AndorSDK3Camera</i> <i>method</i> ), 482	<code>set_row_stretch()</code> ( <i>py-</i> <i>pylablib.core.gui.widgets.container.IQWidgetContainer</i> <i>method</i> ), 208
<code>set_roi()</code> ( <i>pylablib.devices.Basler.pylon.BaslerPylonCamera</i> <i>method</i> ), 515	<code>set_row_stretch()</code> ( <i>py-</i> <i>pylablib.core.gui.widgets.container.QDialogContainer</i> <i>method</i> ), 220
<code>set_roi()</code> ( <i>pylablib.devices.BitFlow.BitFlow.BitFlowCamera</i> <i>method</i> ), 531	<code>set_row_stretch()</code> ( <i>py-</i> <i>pylablib.core.gui.widgets.container.QFrameContainer</i> <i>method</i> ), 216
<code>set_roi()</code> ( <i>pylablib.devices.BitFlow.BitFlow.BitFlowFrameGrabber</i> <i>method</i> ), 522	<code>set_row_stretch()</code> ( <i>py-</i> <i>pylablib.core.gui.widgets.container.QGroupBoxContainer</i> <i>method</i> ), 224
<code>set_roi()</code> ( <i>pylablib.devices.DCAM.DCAM.DCAMCamera</i> <i>method</i> ), 552	<code>set_row_stretch()</code> ( <i>py-</i> <i>pylablib.core.gui.widgets.container.QScrollAreaContainer.QContai</i> <i>method</i> ), 228
<code>set_roi()</code> ( <i>pylablib.devices.IMAQ.IMAQ.IMAQCamera</i> <i>method</i> ), 575	<code>set_row_stretch()</code> ( <i>py-</i> <i>pylablib.core.gui.widgets.container.QWidgetContainer</i> <i>method</i> ), 212
<code>set_roi()</code> ( <i>pylablib.devices.IMAQ.IMAQ.IMAQFrameGrabber</i> <i>method</i> ), 575	<code>set_row_stretch()</code> ( <i>py-</i> <i>pylablib.core.gui.widgets.layout_manager.IQLayoutManagedWidge</i> <i>method</i> ), 239
	<code>set_row_stretch()</code> ( <i>py-</i> <i>pylablib.core.gui.widgets.layout_manager.QLayoutManagedWidge</i> <i>method</i> ), 239



<i>method</i> ), 241		<i>set_single_value()</i> (py-	
<i>set_row_stretch()</i> (py-		<i>lablib.core.gui.value_handling.ISingleValueHandler</i>	
<i>lablib.core.gui.widgets.param_table.ParamTable</i>		<i>method</i> ), 269	
<i>method</i> ), 250		<i>set_single_value()</i> (py-	
<i>set_row_stretch()</i> (py-		<i>lablib.core.gui.value_handling.LabelValueHandler</i>	
<i>lablib.core.gui.widgets.param_table.StatusTable</i>		<i>method</i> ), 270	
<i>method</i> ), 260		<i>set_single_value()</i> (py-	
<i>set_scale()</i> ( <i>pylablib.devices.Thorlabs.elliptec.ElliptecMotor</i>		<i>lablib.core.gui.value_handling.LineEditValueHandler</i>	
<i>method</i> ), 779		<i>method</i> ), 270	
<i>set_scan_params()</i> (py-		<i>set_single_value()</i> (py-	
<i>lablib.devices.Sirah.Matisse.SirahMatisse</i>		<i>lablib.core.gui.value_handling.ProgressBarValueHandler</i>	
<i>method</i> ), 740		<i>method</i> ), 275	
<i>set_scan_position()</i> (py-		<i>set_single_value()</i> (py-	
<i>lablib.devices.Sirah.Matisse.SirahMatisse</i>		<i>lablib.core.gui.value_handling.PushButtonValueHandler</i>	
<i>method</i> ), 740		<i>method</i> ), 273	
<i>set_scan_status()</i> (py-		<i>set_single_value()</i> (py-	
<i>lablib.devices.Sirah.Matisse.SirahMatisse</i>		<i>lablib.core.gui.value_handling.ToolButtonValueHandler</i>	
<i>method</i> ), 740		<i>method</i> ), 274	
<i>set_sensor_curve_index()</i> (py-		<i>set_slowpiezo_ctl_params()</i> (py-	
<i>lablib.devices.Lakeshore.base.Lakeshore218</i>		<i>lablib.devices.Sirah.Matisse.SirahMatisse</i>	
<i>method</i> ), 594		<i>method</i> ), 739	
<i>set_sensor_kind()</i> (py-		<i>set_slowpiezo_ctl_status()</i> (py-	
<i>lablib.devices.Cryocon.base.Cryocon1x</i>		<i>lablib.devices.Sirah.Matisse.SirahMatisse</i>	
<i>method</i> ), 536		<i>method</i> ), 739	
<i>set_sensor_mode()</i> (py-		<i>set_slowpiezo_position()</i> (py-	
<i>lablib.devices.Thorlabs.serial.FW</i>	<i>method</i> ),	<i>lablib.devices.Sirah.Matisse.SirahMatisse</i>	
813		<i>method</i> ), 739	
<i>set_sensor_type()</i> (py-		<i>set_span()</i> ( <i>pylablib.devices.KJL.base.KJL300</i>	
<i>lablib.devices.Lakeshore.base.Lakeshore218</i>		<i>method</i> ), 591	
<i>method</i> ), 594		<i>set_speed_mode()</i> (py-	
<i>set_sensor_voltage()</i> (py-		<i>lablib.devices.Thorlabs.serial.FW</i>	<i>method</i> ),
<i>lablib.devices.Attocube.anc350.ANC350</i>		813	
<i>method</i> ), 508		<i>set_status_line_mode()</i> (py-	
<i>set_serial_parameter()</i> (py-		<i>lablib.devices.PCO.SC2.PCOS2Camera</i>	
<i>lablib.devices.AlliedVision.Bonito.BonitoIMACamera</i>		<i>method</i> ), 648	
<i>method</i> ), 463		<i>set_subsampling()</i> (py-	
<i>set_serial_parameter()</i> (py-		<i>lablib.devices.uc480.uc480.UC480Camera</i>	
<i>lablib.devices.AlliedVision.Bonito.IBonitoCamera</i>		<i>method</i> ), 869	
<i>method</i> ), 451		<i>set_supported_channels()</i> (py-	
<i>set_shutter()</i> (py-		<i>lablib.devices.Thorlabs.kinesis.KinesisMotor</i>	
<i>lablib.devices.Andor.AndorSDK3.AndorSDK3Camera</i>		<i>method</i> ), 795	
<i>method</i> ), 480		<i>set_switcher_mode()</i> (py-	
<i>set_shutter()</i> (py-		<i>lablib.devices.HighFinesse.wlm.WLM</i>	<i>method</i> ),
<i>lablib.devices.LaserQuantum.base.Finesse</i>		559	
<i>method</i> ), 605		<i>set_temperature()</i> (py-	
<i>set_single_value()</i> (py-		<i>lablib.devices.Andor.AndorSDK2.AndorSDK2Camera</i>	
<i>lablib.core.gui.value_handling.CheckboxValueHandler</i>		<i>method</i> ), 467	
<i>method</i> ), 272		<i>set_temperature()</i> (py-	
<i>set_single_value()</i> (py-		<i>lablib.devices.Andor.AndorSDK3.AndorSDK3Camera</i>	
<i>lablib.core.gui.value_handling.ComboBoxValueHandler</i>		<i>method</i> ), 481	
<i>method</i> ), 275		<i>set_temperature()</i> (py-	
<i>set_single_value()</i> (py-		<i>lablib.devices.Photometrics.pvcam.PvcamCamera</i>	
<i>lablib.core.gui.value_handling.IBoolValueHandler</i>		<i>method</i> ), 663	
<i>method</i> ), 272		<i>set_thinet_ctl_params()</i> (py-	

<code>lablib.devices.Sirah.Matisse.SirahMatisse</code> <code>method</code> ), 738	<code>set_trigger_level()</code> <code>lablib.devices.Tektronix.base.TDS2000</code> <code>method</code> ), 761	(py-
<code>set_thinnet_ctl_status()</code> <code>lablib.devices.Sirah.Matisse.SirahMatisse</code> <code>method</code> ), 738	<code>set_trigger_mode()</code> <code>lablib.devices.Andor.AndorSDK2.AndorSDK2Camera</code> <code>method</code> ), 469	(py-
<code>set_timeout()</code> <code>lablib.core.devio.comm_backend.FT232DeviceBackend</code> <code>method</code> ), 150	<code>set_trigger_mode()</code> <code>lablib.devices.Andor.AndorSDK3.AndorSDK3Camera</code> <code>method</code> ), 480	(py-
<code>set_timeout()</code> <code>lablib.core.devio.comm_backend.IDeviceCommBackend</code> <code>method</code> ), 143	<code>set_trigger_mode()</code> <code>lablib.devices.DCAM.DCAM.DCAMCamera</code> <code>method</code> ), 551	(py-
<code>set_timeout()</code> <code>lablib.core.devio.comm_backend.NetworkDeviceBackend</code> <code>method</code> ), 153	<code>set_trigger_mode()</code> <code>lablib.devices.PCO.SC2.PCOS2Camera</code> <code>method</code> ), 646	(py-
<code>set_timeout()</code> <code>lablib.core.devio.comm_backend.PyUSBDeviceBackend</code> <code>method</code> ), 156	<code>set_trigger_mode()</code> <code>lablib.devices.Photometrics.pvcam.PvcamCamera</code> <code>method</code> ), 664	(py-
<code>set_timeout()</code> <code>lablib.core.devio.comm_backend.RecordedDeviceBackend</code> <code>method</code> ), 159	<code>set_trigger_mode()</code> <code>lablib.devices.Tektronix.base.DPO2000</code> <code>method</code> ), 767	(py-
<code>set_timeout()</code> <code>lablib.core.devio.comm_backend.SerialDeviceBackend</code> <code>method</code> ), 148	<code>set_trigger_mode()</code> <code>lablib.devices.Tektronix.base.ITektronixScope</code> <code>method</code> ), 750	(py-
<code>set_timeout()</code> <code>lablib.core.devio.comm_backend.VisaDeviceBackend</code> <code>method</code> ), 146	<code>set_trigger_mode()</code> <code>lablib.devices.Tektronix.base.TDS2000</code> <code>method</code> ), 761	(py-
<code>set_timeout()</code> <code>lablib.core.utils.general.Countdown</code> <code>method</code> ), 379	<code>set_trigger_mode()</code> <code>lablib.devices.Thorlabs.serial.FW</code> <code>method</code> ), 813	(py-
<code>set_timeout()</code> ( <code>pylablib.core.utils.net.ClientSocket</code> <code>method</code> ), 389	<code>set_trigger_mode()</code> <code>lablib.devices.Thorlabs.serial.FWv1</code> <code>method</code> ), 816	(py-
<code>set_timeout()</code> ( <code>pylablib.devices.M2.solstis.Solstis</code> <code>method</code> ), 614	<code>set_trigger_mode()</code> <code>lablib.devices.Thorlabs.TLCamera.ThorlabsTLCamera</code> <code>method</code> ), 771	(py-
<code>set_trigger_input()</code> <code>lablib.devices.Attocube.anc300.ANC300</code> <code>method</code> ), 504	<code>set_trigger_slope()</code> <code>lablib.devices.AWG.generic.GenericAWG</code> <code>method</code> ), 405	(py-
<code>set_trigger_interleave()</code> <code>lablib.devices.PhotonFocus.PhotonFocus.IPhotonFocusCamera</code> <code>method</code> ), 673	<code>set_trigger_slope()</code> <code>lablib.devices.AWG.specific.Agilent33220A</code> <code>method</code> ), 419	(py-
<code>set_trigger_interleave()</code> <code>lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusBitFluxGenerator</code> <code>method</code> ), 700	<code>set_trigger_slope()</code> <code>lablib.devices.AWG.specific.Agilent33500</code> <code>method</code> ), 413	(py-
<code>set_trigger_interleave()</code> <code>lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusIMACalibrator</code> <code>method</code> ), 685	<code>set_trigger_slope()</code> <code>lablib.devices.AWG.specific.InstekAFG2000</code> <code>method</code> ), 431	(py-
<code>set_trigger_interleave()</code> <code>lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusSiSoCalibrator</code> <code>method</code> ), 693	<code>set_trigger_slope()</code> <code>lablib.devices.AWG.specific.InstekAFG2225</code> <code>method</code> ), 425	(py-
<code>set_trigger_level()</code> <code>lablib.devices.Tektronix.base.DPO2000</code> <code>method</code> ), 767	<code>set_trigger_slope()</code> <code>lablib.devices.AWG.specific.RigolDG1000</code> <code>method</code> ), 449	(py-
<code>set_trigger_level()</code> <code>lablib.devices.Tektronix.base.ITektronixScope</code> <code>method</code> ), 750		

---

`set_trigger_slope()` (pylablib.devices.AWG.specific.RSInstekAFG21000 method), 437  
`set_trigger_slope()` (pylablib.devices.AWG.specific.TektronixAFG1000 method), 443  
`set_trigger_source()` (pylablib.devices.AWG.generic.GenericAWG method), 405  
`set_trigger_source()` (pylablib.devices.AWG.specific.Agilent33220A method), 419  
`set_trigger_source()` (pylablib.devices.AWG.specific.Agilent33500 method), 413  
`set_trigger_source()` (pylablib.devices.AWG.specific.InstekAFG2000 method), 431  
`set_trigger_source()` (pylablib.devices.AWG.specific.InstekAFG2225 method), 425  
`set_trigger_source()` (pylablib.devices.AWG.specific.RigolDG1000 method), 449  
`set_trigger_source()` (pylablib.devices.AWG.specific.RSInstekAFG21000 method), 437  
`set_trigger_source()` (pylablib.devices.AWG.specific.TektronixAFG1000 method), 443  
`set_units()` (pylablib.devices.Leybold.base.ITR90 method), 609  
`set_units()` (pylablib.devices.Pfeiffer.base.TPG260 method), 654  
`set_value()` (pylablib.core.gui.value\_handling.CheckboxValueHandler method), 273  
`set_value()` (pylablib.core.gui.value\_handling.ComboBoxValueHandler method), 275  
`set_value()` (pylablib.core.gui.value\_handling.GUIValueHandler method), 280  
`set_value()` (pylablib.core.gui.value\_handling.IBoolValueHandler method), 272  
`set_value()` (pylablib.core.gui.value\_handling.IIndicatorValueHandler method), 276  
`set_value()` (pylablib.core.gui.value\_handling.ISingleValueHandler method), 269  
`set_value()` (pylablib.core.gui.value\_handling.IValueHandler method), 265  
`set_value()` (pylablib.core.gui.value\_handling.LabelIndicatorValueHandler method), 277  
`set_value()` (pylablib.core.gui.value\_handling.LabelValueHandler method), 271  
`set_value()` (pylablib.core.gui.value\_handling.LineEditValueHandler method), 270  
`set_value()` (pylablib.core.gui.value\_handling.ProgressBarValueHandler method), 276  
`set_value()` (pylablib.core.gui.value\_handling.PropertyValueHandler method), 267  
`set_value()` (pylablib.core.gui.value\_handling.PushButtonValueHandler method), 273  
`set_value()` (pylablib.core.gui.value\_handling.StandardIndicatorHandler method), 277  
`set_value()` (pylablib.core.gui.value\_handling.StandardValueHandler method), 268  
`set_value()` (pylablib.core.gui.value\_handling.ToolButtonValueHandler method), 274  
`set_value()` (pylablib.core.gui.value\_handling.VirtualValueHandler method), 266  
`set_value()` (pylablib.core.gui.widgets.button.ToggleButton method), 198  
`set_value()` (pylablib.core.gui.widgets.combo\_box.ComboBox method), 199  
`set_value()` (pylablib.core.gui.widgets.container.IQContainer method), 202  
`set_value()` (pylablib.core.gui.widgets.container.IQWidgetContainer method), 208  
`set_value()` (pylablib.core.gui.widgets.container.QContainer method), 204  
`set_value()` (pylablib.core.gui.widgets.container.QDialogContainer method), 220  
`set_value()` (pylablib.core.gui.widgets.container.QFrameContainer method), 216  
`set_value()` (pylablib.core.gui.widgets.container.QGroupBoxContainer method), 224  
`set_value()` (pylablib.core.gui.widgets.container.QScrollAreaContainer method), 230  
`set_value()` (pylablib.core.gui.widgets.container.QScrollAreaContainer method), 228  
`set_value()` (pylablib.core.gui.widgets.container.QTabContainer method), 233  
`set_value()` (pylablib.core.gui.widgets.container.QWidgetContainer method), 212  
`set_value()` (pylablib.core.gui.widgets.edit.NumEdit method), 236  
`set_value()` (pylablib.core.gui.widgets.edit.TextEdit method), 234  
`set_value()` (pylablib.core.gui.widgets.label.EnumLabel method), 237  
`set_value()` (pylablib.core.gui.widgets.label.NumLabel method), 238  
`set_value()` (pylablib.core.gui.widgets.label.TextLabel method), 236  
`set_value()` (pylablib.core.gui.widgets.param\_table.ParamTable method), 248  
`set_value()` (pylablib.core.gui.widgets.param\_table.StatusTable method), 260  
`set_value()` (pylablib.devices.Andor.AndorSDK3.AndorSDK3Attribute method), 478

`set_value()` (`pylablib.devices.Attocube.anc350.ANC350` `method`), 248  
`method`), 507 `set_visible()` (`py-`  
`set_value()` (`pylablib.devices.Basler.pylon.BaslerPylonAttribute` `lablib.core.gui.widgets.param_table.StatusTable`  
`method`), 513 `method`), 260  
`set_value()` (`pylablib.devices.DCAM.DCAM.DCAMAttribute` `set_voltage()` (`py-`  
`method`), 550 `lablib.devices.Attocube.anc300.ANC300`  
`set_value()` (`pylablib.devices.IMAQdx.IMAQdx.IMAQdxAttribute` `method`), 503  
`method`), 578 `set_voltage()` (`py-`  
`set_value()` (`pylablib.devices.Photometrics.pvcam.PvcamAttribute` `lablib.devices.Attocube.anc350.ANC350`  
`method`), 661 `method`), 508  
`set_value()` (`pylablib.devices.PhotonFocus.PhotonFocus.PFCamAttribute` (`py-`  
`method`), 671 `lablib.devices.OZOptics.base.EPC04` `method`),  
`set_value()` (`pylablib.devices.PrincetonInstruments.picam.PicamAttribute`  
`method`), 711 `set_voltage()` (`py-`  
`set_value()` (`pylablib.devices.SiliconSoftware.fgrab.FGrabAttribute` `lablib.devices.PhysikInstrumente.base.PIE516`  
`method`), 720 `method`), 705  
`set_value_labels()` (`py-` `set_voltage()` (`py-`  
`lablib.core.gui.widgets.button.ToggleButton` `lablib.devices.Thorlabs.serial.MDT69xA`  
`method`), 198 `method`), 820  
`set_variable()` (`py-` `set_voltage_lower_limit()` (`py-`  
`lablib.core.thread.controller.QTaskThread` `lablib.devices.PhysikInstrumente.base.PIE516`  
`method`), 311 `method`), 705  
`set_variable()` (`py-` `set_voltage_outputs()` (`py-`  
`lablib.core.thread.controller.QThreadController` `lablib.devices.NI.daq.NIDAQ` `method`), 625  
`method`), 298 `set_voltage_pattern()` (`py-`  
`set_variable()` (`py-` `lablib.devices.Attocube.anc300.ANC300`  
`lablib.core.utils.ipc.SharedMemIPCTable` `method`), 504  
`method`), 385 `set_voltage_upper_limit()` (`py-`  
`set_velocity()` (`py-` `lablib.devices.PhysikInstrumente.base.PIE516`  
`lablib.devices.PhysikInstrumente.base.PIE516` `method`), 705  
`method`), 705 `set_vsspeed()` (`py-`  
`set_velocity()` (`py-` `lablib.devices.Andor.AndorSDK2.AndorSDK2Camera`  
`lablib.devices.Thorlabs.elliptec.ElliptecMotor` `method`), 468  
`method`), 779 `set_wait_callback()` (`py-`  
`set_vertical_position()` (`py-` `lablib.core.utils.net.ClientSocket` `method`),  
`lablib.devices.Tektronix.base.DPO2000` 389  
`method`), 767 `set_waveform()` (`py-`  
`set_vertical_position()` (`py-` `lablib.devices.OZOptics.base.EPC04` `method`),  
`lablib.devices.Tektronix.base.ITektronixScope` 637  
`method`), 750 `set_wavelength()` (`py-`  
`set_vertical_position()` (`py-` `lablib.devices.Ophir.base.VegaPowerMeter`  
`lablib.devices.Tektronix.base.TDS2000` `method`), 642  
`method`), 761 `set_wavelength()` (`py-`  
`set_vertical_span()` (`py-` `lablib.devices.OZOptics.base.TF100` `method`),  
`lablib.devices.Tektronix.base.DPO2000` 633  
`method`), 767 `set_wavelength_correction()` (`py-`  
`set_vertical_span()` (`py-` `lablib.devices.OZOptics.base.TF100` `method`),  
`lablib.devices.Tektronix.base.ITektronixScope` 633  
`method`), 750 `set_xarg_name()` (`py-`  
`set_vertical_span()` (`py-` `lablib.core.dataproc.fitting.Fitter` `method`),  
`lablib.devices.Tektronix.base.TDS2000` 115  
`method`), 761 `set_zero()` (`pylablib.devices.KJL.base.KJL300`  
`set_visible()` (`py-` `method`), 591  
`lablib.core.gui.widgets.param_table.ParamTable` `setattr_call()` (`in` `module` `py-`

[lablib.core.utils.functions\), 372](#)  
[setbp \(\) \(in module pylablib\), 875](#)  
[setbp \(\) \(in module pylablib.core.utils.general\), 380](#)  
[setdefault \(\) \(pylablib.core.utils.dictionary.Dictionary method\), 328](#)  
[setdefault \(\) \(pylablib.core.utils.dictionary.DictionaryPointer method\), 342](#)  
[setdefault \(\) \(pylablib.core.utils.dictionary.FilterTree method\), 358](#)  
[setdefault \(\) \(pylablib.core.utils.dictionary.ItemAccess method\), 361](#)  
[setdefault \(\) \(pylablib.core.utils.dictionary.PrefixTree method\), 350](#)  
[setpoint \(pylablib.devices.Sirah.Matisse.TFastpiezoCtlParameters attribute\), 735](#)  
[setpoint \(pylablib.devices.Sirah.Matisse.TSlowpiezoCtlParameters attribute\), 735](#)  
[setpoint \(pylablib.devices.Sirah.Matisse.TThinnetCtlParameters attribute\), 734](#)  
[settle\\_time \(pylablib.devices.Lakeshore.base.TLakeshore370FilterStub device attribute\), 599](#)  
[setup \(\) \(pylablib.core.gui.widgets.container.IQContainers method\), 200](#)  
[setup \(\) \(pylablib.core.gui.widgets.container.IQWidgetContainer method\), 205](#)  
[setup \(\) \(pylablib.core.gui.widgets.container.QContainer method\), 204](#)  
[setup \(\) \(pylablib.core.gui.widgets.container.QDialogContainer method\), 220](#)  
[setup \(\) \(pylablib.core.gui.widgets.container.QFrameContainer method\), 216](#)  
[setup \(\) \(pylablib.core.gui.widgets.container.QGroupBoxContainer method\), 221](#)  
[setup \(\) \(pylablib.core.gui.widgets.container.QScrollAreaContainer method\), 229](#)  
[setup \(\) \(pylablib.core.gui.widgets.container.QScrollAreaContainer.QCloudWidget method\), 228](#)  
[setup \(\) \(pylablib.core.gui.widgets.container.QTabContainer method\), 233](#)  
[setup \(\) \(pylablib.core.gui.widgets.container.QWidgetContainer method\), 212](#)  
[setup \(\) \(pylablib.core.gui.widgets.layout\\_manager.IQLayoutManager method\), 238](#)  
[setup \(\) \(pylablib.core.gui.widgets.layout\\_manager.QLayoutManager method\), 241](#)  
[setup \(\) \(pylablib.core.gui.widgets.param\\_table.ParamTable method\), 242](#)  
[setup \(\) \(pylablib.core.gui.widgets.param\\_table.StatusTable method\), 251](#)  
[setup \(\) \(pylablib.devices.Thorlabs.TLCamera.ThorlabsTLCamera method\), 771](#)  
[setup\\_accum\\_mode \(\) \(pylablib.devices.Andor.AndorSDK2.AndorSDK2Camera method\), 470](#)  
[setup\\_acquisition \(\) \(pylablib.devices.AlliedVision.Bonito.BonitoIMAQCamera method\), 463](#)  
[setup\\_acquisition \(\) \(pylablib.devices.AlliedVision.Bonito.IBonitoCamera method\), 451](#)  
[setup\\_acquisition \(\) \(pylablib.devices.Andor.AndorSDK2.AndorSDK2Camera method\), 472](#)  
[setup\\_acquisition \(\) \(pylablib.devices.Andor.AndorSDK3.AndorSDK3Camera method\), 482](#)  
[setup\\_acquisition \(\) \(pylablib.devices.Basler.pylon.BaslerPylonCamera method\), 516](#)  
[setup\\_acquisition \(\) \(pylablib.devices.BitFlow.BitFlow.BitFlowCamera method\), 531](#)  
[setup\\_acquisition \(\) \(pylablib.devices.BitFlow.BitFlow.BitFlowFrameGrabber method\), 523](#)  
[setup\\_acquisition \(\) \(pylablib.devices.DCAM.DCAM.DCAMCamera method\), 552](#)  
[setup\\_acquisition \(\) \(pylablib.devices.IMAQ.IMAQ.IMAQCamera method\), 575](#)  
[setup\\_acquisition \(\) \(pylablib.devices.IMAQ.IMAQ.IMAQFrameGrabber method\), 564](#)  
[setup\\_acquisition \(\) \(pylablib.devices.IMAQdx.IMAQdx.EthernetIMAQdxCamera method\), 589](#)  
[setup\\_acquisition \(\) \(pylablib.devices.IMAQdx.IMAQdx.IMAQdxCamera method\), 574](#)  
[setup\\_acquisition \(\) \(pylablib.devices.interface.camera.IAttributeCamera method\), 841](#)  
[setup\\_acquisition \(\) \(pylablib.devices.interface.camera.IBinROICamera method\), 860](#)  
[setup\\_acquisition \(\) \(pylablib.devices.interface.camera.ICamera method\), 831](#)  
[setup\\_acquisition \(\) \(pylablib.devices.interface.camera.IExposureCamera method\), 850](#)  
[setup\\_acquisition \(\) \(pylablib.devices.interface.camera.IGrabberAttributeCamera method\), 846](#)  
[setup\\_acquisition \(\) \(pylablib.devices.interface.camera.IROICamera method\), 855](#)



setup_acquisition() <i>lablib.devices.PCO.SC2.PCOSC2Camera</i> <i>method</i> ), 647	(py-	setup_clock() <i>method</i> ), 622	(pylablib.devices.NI.daq.NIDAQ
setup_acquisition() <i>lablib.devices.Photometrics.pvcam.PvcamCamera</i> <i>method</i> ), 665	(py-	setup_cont_mode() <i>method</i> ), 470	(py-
setup_acquisition() <i>lablib.devices.PhotonFocus.PhotonFocus.IPhotonFocusCamera</i> <i>method</i> ), 677	(py-	setup_cooldown() <i>method</i> ), 152	(py-
setup_acquisition() <i>lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusBitFlowCamera</i> <i>method</i> ), 695	(py-	setup_cooldown() <i>method</i> ), 143	(py-
setup_acquisition() <i>lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusIMACamera</i> <i>method</i> ), 686	(py-	setup_cooldown() <i>method</i> ), 154	(py-
setup_acquisition() <i>lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusSiSoCamera</i> <i>method</i> ), 693	(py-	setup_cooldown() <i>method</i> ), 157	(py-
setup_acquisition() <i>lablib.devices.PrincetonInstruments.picam.PicamCamera</i> <i>method</i> ), 714	(py-	setup_cooldown() <i>method</i> ), 159	(py-
setup_acquisition() <i>lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareCamera</i> <i>method</i> ), 732	(py-	setup_cooldown() <i>method</i> ), 149	(py-
setup_acquisition() <i>lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareFrameGrabber</i> <i>method</i> ), 723	(py-	setup_cooldown() <i>method</i> ), 146	(py-
setup_acquisition() <i>lablib.devices.Thorlabs.TLCamera.ThorlabsTLCamera</i> <i>method</i> ), 771	(py-	setup_current() <i>method</i> ), 828	(py-
setup_acquisition() <i>lablib.devices.uc480.uc480.UC480Camera</i> <i>method</i> ), 868	(py-	setup_drive() <i>method</i> ), 802	(py-
setup_analog_output() <i>lablib.devices.Lakeshore.base.Lakeshore218</i> <i>method</i> ), 594	(py-	setup_edge_trigger() <i>method</i> ), 767	(py-
setup_analog_output() <i>lablib.devices.Lakeshore.base.Lakeshore370</i> <i>method</i> ), 600	(py-	setup_edge_trigger() <i>method</i> ), 750	(py-
setup_autocalibration() <i>lablib.devices.HighFinesse.wlm.WLM</i> <i>method</i> ), 560	(py-	setup_edge_trigger() <i>method</i> ), 761	(py-
setup_camlink_pixel_format() <i>lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusSiSoCamera</i> <i>method</i> ), 693	(py-	setup_ethernet() <i>method</i> ), 629	(py-
setup_camlink_pixel_format() <i>lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareCamera</i> <i>method</i> ), 732	(py-	setup_ext_trigger() <i>method</i> ), 469	(py-
setup_camlink_pixel_format() <i>lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareFrameGrabber</i> <i>method</i> ), 722	(py-	setup_ext_trigger() <i>method</i> ), 551	(py-
setup_channel_range() <i>lablib.devices.Lakeshore.base.Lakeshore370</i> <i>method</i> ), 599	(py-	setup_ext_trigger() <i>method</i> ), 771	(py-
		setup_fast_kinetic_mode() <i>method</i> ), 771	(py-

<code>lablib.devices.Andor.AndorSDK2.AndorSDK2Camera</code> <code>method</code> ), 470	<code>lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusSiSoCamera</code> <code>method</code> ), 693
<code>setup_filter()</code> ( <code>py-</code> <code>lablib.devices.Lakeshore.base.Lakeshore218</code> <code>method</code> ), 595	<code>setup_multi_track_mode()</code> ( <code>py-</code> <code>lablib.devices.Andor.AndorSDK2.AndorSDK2Camera</code> <code>method</code> ), 471
<code>setup_filter()</code> ( <code>py-</code> <code>lablib.devices.Lakeshore.base.Lakeshore370</code> <code>method</code> ), 600	<code>setup_name()</code> ( <code>pylablib.core.gui.widgets.container.IQContainer</code> <code>method</code> ), 200
<code>setup_flipper()</code> ( <code>py-</code> <code>lablib.devices.Thorlabs.kinesis.MFF</code> <code>method</code> ), 791	<code>setup_name()</code> ( <code>pylablib.core.gui.widgets.container.IQWidgetContainer</code> <code>method</code> ), 208
<code>setup_func()</code> ( <code>in</code> <code>module</code> <code>py-</code> <code>lablib.core.utils.ctypes_wrap</code> ), 321	<code>setup_name()</code> ( <code>pylablib.core.gui.widgets.container.QContainer</code> <code>method</code> ), 204
<code>setup_gauge_control()</code> ( <code>py-</code> <code>lablib.devices.Pfeiffer.base.TPG260</code> <code>method</code> ), 656	<code>setup_name()</code> ( <code>pylablib.core.gui.widgets.container.QDialogContainer</code> <code>method</code> ), 220
<code>setup_gen_move()</code> ( <code>py-</code> <code>lablib.devices.Thorlabs.kinesis.KinesisMotor</code> <code>method</code> ), 797	<code>setup_name()</code> ( <code>pylablib.core.gui.widgets.container.QFrameContainer</code> <code>method</code> ), 216
<code>setup_homing()</code> ( <code>py-</code> <code>lablib.devices.Thorlabs.kinesis.KinesisMotor</code> <code>method</code> ), 797	<code>setup_name()</code> ( <code>pylablib.core.gui.widgets.container.QGroupBoxContainer</code> <code>method</code> ), 224
<code>setup_image_mode()</code> ( <code>py-</code> <code>lablib.devices.Andor.AndorSDK2.AndorSDK2Camera</code> <code>method</code> ), 472	<code>setup_name()</code> ( <code>pylablib.core.gui.widgets.container.QScrollAreaContainer</code> <code>method</code> ), 230
<code>setup_jog()</code> ( <code>pylablib.devices.Thorlabs.kinesis.KinesisMotor</code> <code>method</code> ), 797	<code>setup_name()</code> ( <code>pylablib.core.gui.widgets.container.QScrollAreaContainer</code> <code>method</code> ), 228
<code>setup_jog()</code> ( <code>pylablib.devices.Thorlabs.kinesis.KinesisPiezoMotor</code> <code>method</code> ), 802	<code>setup_name()</code> ( <code>pylablib.core.gui.widgets.container.QTabContainer</code> <code>method</code> ), 233
<code>setup_kcube_trigio()</code> ( <code>py-</code> <code>lablib.devices.Thorlabs.kinesis.KinesisMotor</code> <code>method</code> ), 797	<code>setup_name()</code> ( <code>pylablib.core.gui.widgets.container.QWidgetContainer</code> <code>method</code> ), 212
<code>setup_kcube_trigpos()</code> ( <code>py-</code> <code>lablib.devices.Thorlabs.kinesis.KinesisMotor</code> <code>method</code> ), 798	<code>setup_name()</code> ( <code>pylablib.core.gui.widgets.param_table.ParamTable</code> <code>method</code> ), 251
<code>setup_kinetic_mode()</code> ( <code>py-</code> <code>lablib.devices.Andor.AndorSDK2.AndorSDK2Camera</code> <code>method</code> ), 470	<code>setup_name()</code> ( <code>pylablib.core.gui.widgets.param_table.StatusTable</code> <code>method</code> ), 260
<code>setup_limit_switch()</code> ( <code>py-</code> <code>lablib.devices.Thorlabs.kinesis.KinesisMotor</code> <code>method</code> ), 797	<code>setup_pixel_correction()</code> ( <code>py-</code> <code>lablib.devices.Thorlabs.TLCamera.ThorlabsTLCamera</code> <code>method</code> ), 771
<code>setup_limit_switches()</code> ( <code>py-</code> <code>lablib.devices.Trinamic.base.TMCM1110</code> <code>method</code> ), 828	<code>setup_random_track_mode()</code> ( <code>py-</code> <code>lablib.devices.Andor.AndorSDK2.AndorSDK2Camera</code> <code>method</code> ), 471
<code>setup_max_baudrate()</code> ( <code>py-</code> <code>lablib.devices.PhotonFocus.PhotonFocus.IPhotonFocusCamera</code> <code>method</code> ), 672	<code>setup_serial_params()</code> ( <code>py-</code> <code>lablib.devices.AlliedVision.Bonito.BonitoIMAQCamera</code> <code>method</code> ), 463
<code>setup_max_baudrate()</code> ( <code>py-</code> <code>lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusBitFluxCamera</code> <code>method</code> ), 700	<code>setup_serial_params()</code> ( <code>py-</code> <code>lablib.devices.IMAQ.IMAQ.IMAQCamera</code> <code>method</code> ), 575
<code>setup_max_baudrate()</code> ( <code>py-</code> <code>lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusIMAQCamera</code> <code>method</code> ), 686	<code>setup_serial_params()</code> ( <code>py-</code> <code>lablib.devices.IMAQ.IMAQ.IMAQFrameGrabber</code> <code>method</code> ), 564
<code>setup_max_baudrate()</code> ( <code>py-</code> <code>lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusIMAQCamera</code> <code>method</code> ), 686	<code>setup_serial_params()</code> ( <code>py-</code> <code>lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusIMAQCamera</code> <code>method</code> ), 686
<code>setup_max_baudrate()</code> ( <code>py-</code> <code>lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusBitFluxCamera</code> <code>method</code> ), 700	<code>setup_shutter()</code> ( <code>py-</code> <code>lablib.devices.Andor.AndorSDK2.AndorSDK2Camera</code> <code>method</code> ), 469
<code>setup_max_baudrate()</code> ( <code>py-</code> <code>lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusIMAQCamera</code> <code>method</code> ), 686	<code>setup_single_track_mode()</code> ( <code>py-</code> <code>lablib.devices.Andor.AndorSDK2.AndorSDK2Camera</code> <code>method</code> ), 471
<code>setup_max_baudrate()</code> ( <code>py-</code>	<code>setup_switch()</code> ( <code>py-</code>

[lablib.devices.Pfeiffer.base.TPG260](#) (*method*), [shutdown\(\)](#) (*pylablib.devices.Photometrics.pvcam.LibraryController* *method*), [659](#)  
[655](#)  
[setup\\_task\(\)](#) (*pylablib.core.thread.controller.QTaskThread* *method*), [304](#)  
[shutdown\(\)](#) (*pylablib.devices.PhotonFocus.PhotonFocus.LibraryController* *method*), [669](#)  
[setup\\_terascan\(\)](#) (*pylablib.devices.M2.solstis.Solstis* *method*), [708](#)  
[617](#)  
[shutdown\(\)](#) (*pylablib.devices.SmarAct.scu3d.LibraryController* *method*), [745](#)  
[setup\\_velocity\(\)](#) (*pylablib.devices.Newport.picomotor.Picomotor8742* *method*), [769](#)  
[630](#)  
[shutdown\(\)](#) (*pylablib.devices.Thorlabs.TLCamera.LibraryController* *method*), [769](#)  
[setup\\_velocity\(\)](#) (*pylablib.devices.Thorlabs.kinesis.KinesisMotor* *method*), [796](#)  
[signature\(\)](#) (*pylablib.core.utils.functions.FunctionSignature* *method*), [370](#)  
[setup\\_velocity\(\)](#) (*pylablib.devices.Trinamic.base.TMCM1110* *method*), [828](#)  
[SilenceException](#) (*class in pylablib.core.utils.general*), [377](#)  
[setup\\_voltage\\_output\\_clock\(\)](#) (*pylablib.devices.NI.daq.NIDAQ* *method*), [626](#)  
[SiliconSoftwareCamera](#) (*class in pylablib.devices.SiliconSoftware.fgrab*), [727](#)  
[sfglob\(\)](#) (*in module pylablib.core.utils.string*), [398](#)  
[SiliconSoftwareCamera.NoParameterCaller](#) (*class in pylablib.devices.SiliconSoftware.fgrab*), [727](#)  
[sfregex\(\)](#) (*in module pylablib.core.utils.string*), [398](#)  
[shape\(\)](#) (*pylablib.core.dataproc.table\_wrap.Array1DWrapper* *method*), [127](#)  
[SiliconSoftwareFrameGrabber](#) (*class in pylablib.devices.SiliconSoftware.fgrab*), [721](#)  
[shape\(\)](#) (*pylablib.core.dataproc.table\_wrap.Array2DWrapper* *method*), [131](#)  
[SiliconSoftwareFrameGrabber.NoParameterCaller](#) (*class in pylablib.devices.SiliconSoftware.fgrab*), [723](#)  
[shape\(\)](#) (*pylablib.core.dataproc.table\_wrap.DataFrame2DWrapper* *method*), [133](#)  
[single\\_op\(\)](#) (*pylablib.core.devio.comm\_backend.FT232DeviceBackend* *method*), [150](#)  
[shape\(\)](#) (*pylablib.core.dataproc.table\_wrap.I1DWrapper* *method*), [126](#)  
[single\\_op\(\)](#) (*pylablib.core.devio.comm\_backend.SerialDeviceBackend* *method*), [148](#)  
[shape\(\)](#) (*pylablib.core.dataproc.table\_wrap.I2DWrapper* *method*), [129](#)  
[SingleFileSystemDataLocation](#) (*class in pylablib.core.fileio.location*), [186](#)  
[shape\(\)](#) (*pylablib.core.dataproc.table\_wrap.IGenWrapper* *method*), [125](#)  
[SirahMatisse](#) (*class in pylablib.devices.Sirah.Matisse*), [736](#)  
[shape\(\)](#) (*pylablib.core.dataproc.table\_wrap.Series1DWrapper* *method*), [128](#)  
[SirahMatisse.NoParameterCaller](#) (*class in pylablib.devices.Sirah.Matisse*), [740](#)  
[SharedMemIPCCChannel](#) (*class in pylablib.core.utils.ipc*), [384](#)  
[size](#) (*pylablib.core.utils.ipc.TShmemVarDesc* *attribute*), [384](#)  
[SharedMemIPCTable](#) (*class in pylablib.core.utils.ipc*), [384](#)  
[shifted\(\)](#) (*pylablib.core.dataproc.transform.Indexed2DTransform* *method*), [135](#)  
[shifted\(\)](#) (*pylablib.devices.Andor.AndorSDK3.TFrameInfo* *attribute*), [479](#)  
[shifted\(\)](#) (*pylablib.core.dataproc.transform.LinearTransform* *method*), [134](#)  
[size\(\)](#) (*pylablib.devices.uc480.uc480.TFrameInfo* *attribute*), [866](#)  
[show\\_value\(\)](#) (*pylablib.core.gui.widgets.edit.NumEdit* *method*), [236](#)  
[size\(\)](#) (*pylablib.core.dataproc.image.ROI* *method*), [121](#)  
[show\\_value\(\)](#) (*pylablib.core.gui.widgets.edit.TextEdit* *method*), [234](#)  
[method\), \[328\]\(#\)  
\[shutdown\\(\\)\]\(#\) \(\*pylablib.devices.Andor.AndorSDK2.LibraryController\* \*method\*\), \[465\]\(#\)  
\[DictionaryPointer\]\(#\) \(\*pylablib.core.utils.dictionary.DictionaryPointer\* \*method\*\), \[342\]\(#\)  
\[shutdown\\(\\)\]\(#\) \(\*pylablib.devices.Andor.AndorSDK3.LibraryController\* \*method\*\), \[476\]\(#\)  
\[method\\), \\[358\\]\\(#\\)  
\\[shutdown\\\(\\\)\\]\\(#\\) \\(\\*pylablib.devices.Basler.pylon.LibraryController\\* \\*method\\*\\), \\[511\\]\\(#\\)  
\\[FilterTree\\]\\(#\\) \\(\\*pylablib.core.utils.dictionary.PrefixTree\\* \\*method\\*\\), \\[350\\]\\(#\\)  
\\[shutdown\\\(\\\)\\]\\(#\\) \\(\\*pylablib.devices.DCAM.DCAM.LibraryController\\* \\*method\\*\\), \\[548\\]\\(#\\)  
\\[call\\\\_sync\\]\\(#\\) \\(\\*pylablib.core.thread.callsync.QScheduledCall\\* \\*method\\*\\), \\[284\\]\\(#\\)\]\(#\)](#)



[skipped \(pylablib.devices.Andor.AndorSDK3.TMissedFramesStatusmethod\), 822](#)  
[attribute\), 479](#)      [sleep \(\) \(pylablib.devices.Thorlabs.serial.ThorlabsSerialInterface](#)  
[skipped \(pylablib.devices.interface.camera.TFramesStatus](#)      [method\), 811](#)  
[attribute\), 830](#)      [sliding\\_average \(\) \(in module py-](#)  
[skipped \(\) \(pylablib.core.thread.callsync.QCallResultSynchronizerlablib.core.dataproc.filters\), 111](#)  
[method\), 281](#)      [sliding\\_filter \(\) \(in module py-](#)  
[skipped \(\) \(pylablib.core.thread.callsync.QDirectResultSynchronizerlablib.core.dataproc.filters\), 111](#)  
[method\), 282](#)      [slope \(pylablib.devices.Tektronix.base.TTriggerParameters](#)  
[attribute\), 748](#)  
[SkippedCallError, 320](#)      [slow\\_piezo\\_tune\\_to \(\) \(py-](#)  
[sleep \(\) \(pylablib.core.devio.SCPI.SCPIDevice](#)      [lablib.devices.Sirah.tuner.MatisseTuner](#)  
[method\), 139](#)      [method\), 744](#)  
[sleep \(\) \(pylablib.core.thread.controller.QTaskThread](#)  
[method\), 311](#)      [snap \(\) \(pylablib.devices.AlliedVision.Bonito.BonitoIMAQCamera](#)  
[sleep \(\) \(pylablib.core.thread.controller.QThreadController](#)      [method\), 464](#)  
[method\), 295](#)      [snap \(\) \(pylablib.devices.AlliedVision.Bonito.IBonitoCamera](#)  
[sleep \(\) \(pylablib.devices.AWG.generic.GenericAWG](#)      [method\), 456](#)  
[method\), 407](#)      [snap \(\) \(pylablib.devices.Andor.AndorSDK2.AndorSDK2Camera](#)  
[sleep \(\) \(pylablib.devices.AWG.specific.Agilent33220A](#)      [method\), 476](#)  
[method\), 419](#)      [snap \(\) \(pylablib.devices.Andor.AndorSDK3.AndorSDK3Camera](#)  
[sleep \(\) \(pylablib.devices.AWG.specific.Agilent33500](#)      [method\), 486](#)  
[method\), 413](#)      [snap \(\) \(pylablib.devices.Basler.pylon.BaslerPylonCamera](#)  
[sleep \(\) \(pylablib.devices.AWG.specific.InstekAFG2000](#)      [method\), 520](#)  
[method\), 431](#)      [snap \(\) \(pylablib.devices.BitFlow.BitFlow.BitFlowCamera](#)  
[sleep \(\) \(pylablib.devices.AWG.specific.InstekAFG2225](#)      [method\), 532](#)  
[method\), 425](#)      [snap \(\) \(pylablib.devices.BitFlow.BitFlow.BitFlowFrameGrabber](#)  
[sleep \(\) \(pylablib.devices.AWG.specific.RigolDG1000](#)      [method\), 526](#)  
[method\), 449](#)      [snap \(\) \(pylablib.devices.DCAM.DCAM.DCAMCamera](#)  
[sleep \(\) \(pylablib.devices.AWG.specific.RSInstekAFG21000](#)      [method\), 556](#)  
[method\), 437](#)      [snap \(\) \(pylablib.devices.IMAQ.IMAQ.IMAQCamera](#)  
[sleep \(\) \(pylablib.devices.AWG.specific.TektronixAFG1000](#)      [method\), 575](#)  
[method\), 443](#)      [snap \(\) \(pylablib.devices.IMAQ.IMAQ.IMAQFrameGrabber](#)  
[sleep \(\) \(pylablib.devices.Cryocon.base.Cryocon1x](#)      [method\), 568](#)  
[method\), 538](#)      [snap \(\) \(pylablib.devices.IMAQdx.IMAQdx.EthernetIMAQdxCamera](#)  
[sleep \(\) \(pylablib.devices.Cryomagnetics.base.LM500](#)      [method\), 589](#)  
[method\), 543](#)      [snap \(\) \(pylablib.devices.IMAQdx.IMAQdx.IMAQdxCamera](#)  
[sleep \(\) \(pylablib.devices.Cryomagnetics.base.LM510](#)      [method\), 583](#)  
[method\), 547](#)      [snap \(\) \(pylablib.devices.interface.camera.IAttributeCamera](#)  
[sleep \(\) \(pylablib.devices.Lakeshore.base.Lakeshore218](#)      [method\), 841](#)  
[method\), 597](#)      [snap \(\) \(pylablib.devices.interface.camera.IBinROICamera](#)  
[sleep \(\) \(pylablib.devices.Lakeshore.base.Lakeshore370](#)      [method\), 860](#)  
[method\), 602](#)      [snap \(\) \(pylablib.devices.interface.camera.ICamera](#)  
[sleep \(\) \(pylablib.devices.Sirah.Matisse.SirahMatisse](#)      [method\), 834](#)  
[method\), 742](#)      [snap \(\) \(pylablib.devices.interface.camera.IExposureCamera](#)  
[sleep \(\) \(pylablib.devices.Tektronix.base.DPO2000](#)      [method\), 851](#)  
[method\), 768](#)      [snap \(\) \(pylablib.devices.interface.camera.IGrabberAttributeCamera](#)  
[sleep \(\) \(pylablib.devices.Tektronix.base.ITektronixScope](#)      [method\), 846](#)  
[method\), 754](#)      [snap \(\) \(pylablib.devices.interface.camera.IROICamera](#)  
[sleep \(\) \(pylablib.devices.Tektronix.base.TDS2000](#)      [method\), 855](#)  
[method\), 761](#)      [snap \(\) \(pylablib.devices.PCO.SC2.PCOS2Camera](#)  
[sleep \(\) \(pylablib.devices.Thorlabs.serial.FW method\),](#)      [method\), 652](#)  
[815](#)      [snap \(\) \(pylablib.devices.Photometrics.pvcam.PvcamCamera](#)  
[sleep \(\) \(pylablib.devices.Thorlabs.serial.FWv1](#)      [method\), 668](#)  
[method\), 818](#)      [snap \(\) \(pylablib.devices.PhotonFocus.PhotonFocus.IPhotonFocusCamera](#)  
[sleep \(\) \(pylablib.devices.Thorlabs.serial.MDT69xA](#)      [method\), 677](#)

`snap()` (`pylablib.devices.PhotonFocus.PhotonFocus.PhotonFocusBitFlowCamera` attribute), 708  
`method`), 700  
`source` (`pylablib.devices.Thorlabs.kinesis.KinesisPiezoMotor.CommData` attribute), 802  
`snap()` (`pylablib.devices.PhotonFocus.PhotonFocus.PhotonFocusImageCamera` attribute), 686  
`method`), 686  
`source` (`pylablib.devices.Thorlabs.kinesis.KinesisPiezoMotor.CommShortData` attribute), 803  
`snap()` (`pylablib.devices.PhotonFocus.PhotonFocus.PhotonFocusSiCImageCamera` attribute), 693  
`method`), 693  
`source` (`pylablib.devices.Thorlabs.kinesis.KinesisQuadDetector.CommData` attribute), 807  
`snap()` (`pylablib.devices.PrincetonInstruments.picam.PicamCamera` attribute), 717  
`method`), 717  
`source` (`pylablib.devices.Thorlabs.kinesis.KinesisQuadDetector.CommShortData` attribute), 807  
`snap()` (`pylablib.devices.SiliconSoftware.fgrab.SiliconSoftwareCamera` attribute), 733  
`method`), 733  
`source` (`pylablib.devices.Thorlabs.kinesis.MFF.CommData` attribute), 792  
`snap()` (`pylablib.devices.SiliconSoftware.fgrab.SiliconSoftwareFrameGrabber` attribute), 727  
`method`), 727  
`source` (`pylablib.devices.Thorlabs.kinesis.MFF.CommShortData` attribute), 792  
`snap()` (`pylablib.devices.Thorlabs.TLCamera.ThorlabsTLCamera` attribute), 776  
`method`), 776  
`speed_freq` (`pylablib.devices.Photometrics.pvcam.TReadoutInfo` attribute), 662  
`snap()` (`pylablib.devices.uc480.uc480.UC480Camera` attribute), 872  
`method`), 872  
`speed_idx` (`pylablib.devices.Photometrics.pvcam.TReadoutInfo` attribute), 662  
`sock` (`pylablib.core.utils.net.ClientSocket` attribute), 389  
`SocketError`, 387  
`SocketTimeout`, 388  
`SocketTunnelService` (class in `pylablib.core.utils.rpyc_utils`), 394  
`software_version` (`pylablib.devices.Andor.AndorSDK3.TDeviceInfo` attribute), 479  
`software_version` (`pylablib.devices.LaserQuantum.base.TDeviceInfo` attribute), 604  
`software_version` (`pylablib.devices.SiliconSoftware.fgrab.TDeviceInfo` attribute), 720  
`Solstis` (class in `pylablib.devices.M2.solstis`), 614  
`Solstis.NoParameterCaller` (class in `pylablib.devices.M2.solstis`), 619  
`sort_by()` (in module `pylablib.core.dataproc.utils`), 136  
`sort_set_by_list()` (in module `pylablib.core.utils.general`), 376  
`source` (`pylablib.devices.Lakeshore.base.TLakeshore218AnalogSettings` attribute), 593  
`source` (`pylablib.devices.Lakeshore.base.TLakeshore370AnalogSettings` attribute), 599  
`source` (`pylablib.devices.Tektronix.base.TTriggerParameters` attribute), 748  
`source` (`pylablib.devices.Thorlabs.kinesis.BasicKinesisDevice.CommData` attribute), 782  
`source` (`pylablib.devices.Thorlabs.kinesis.BasicKinesisDevice.CommShortData` attribute), 782  
`source` (`pylablib.devices.Thorlabs.kinesis.KinesisDevice.CommData` attribute), 788  
`source` (`pylablib.devices.Thorlabs.kinesis.KinesisDevice.CommShortData` attribute), 788  
`source` (`pylablib.devices.Thorlabs.kinesis.KinesisMotor.CommData` attribute), 798  
`source` (`pylablib.devices.Thorlabs.kinesis.KinesisMotor.CommShortData` attribute), 798  
`split_in_groups()` (in module `pylablib.core.utils.general`), 376  
`split_into_bins()` (in module `pylablib.core.dataproc.filters`), 113  
`split_path()` (in module `pylablib.core.utils.dictionary`), 325  
`split_units()` (in module `pylablib.core.utils.units`), 401  
`SproutG` (class in `pylablib.devices.LighthousePhotonics.base`), 612  
`SproutG.NoParameterCaller` (class in `pylablib.devices.LighthousePhotonics.base`), 612  
`src` (`pylablib.core.thread.multicast_pool.TMulticast` attribute), 315  
`sstep` (`pylablib.devices.interface.camera.TAxisROILimit` attribute), 851  
`StandardIndicatorHandler` (class in `pylablib.core.gui.value_handling`), 276  
`StandardValueHandler` (class in `pylablib.core.gui.value_handling`), 268  
`start` (`pylablib.core.dataproc.utils.Range` attribute), 136  
`start` (`pylablib.core.gui.widgets.container.TTimerEvent` attribute), 200  
`start` (in module `pylablib.core.thread.profile`), 317  
`start()` (`pylablib.core.devio.backend_logger.BackendLogger` attribute), 141  
`start()` (`pylablib.core.devio.comm_backend.RecordedDeviceBackend` attribute), 158  
`start()` (`pylablib.core.gui.widgets.container.IQContainer` attribute), 201  
`start()` (`pylablib.core.gui.widgets.container.IQWidgetContainer` attribute), 201  
`start()` (`pylablib.core.gui.widgets.container.QContainer` attribute), 204

`start()` (`pylablib.core.gui.widgets.container.QDialogContainer` method), 220  
`start()` (`pylablib.core.gui.widgets.container.QFrameContainer` method), 216  
`start()` (`pylablib.core.gui.widgets.container.QGroupBoxContainer` method), 224  
`start()` (`pylablib.core.gui.widgets.container.QScrollAreaContainer` method), 230  
`start()` (`pylablib.core.gui.widgets.container.QScrollAreaContainer.QContainedWidget` method), 228  
`start()` (`pylablib.core.gui.widgets.container.QTabContainer` method), 233  
`start()` (`pylablib.core.gui.widgets.container.QWidgetContainer` method), 212  
`start()` (`pylablib.core.gui.widgets.param_table.ParamTable` method), 251  
`start()` (`pylablib.core.gui.widgets.param_table.StatusTable` method), 260  
`start()` (`pylablib.core.thread.controller.QTaskThread` method), 312  
`start()` (`pylablib.core.thread.controller.QThreadController` method), 299  
`start()` (`pylablib.devices.NI.daq.NIDAQ` method), 623  
`start_acquisition()` (`pylablib.devices.AlliedVision.Bonito.BonitoIMAQCamera` method), 464  
`start_acquisition()` (`pylablib.devices.AlliedVision.Bonito.IBonitoCamera` method), 456  
`start_acquisition()` (`pylablib.devices.Andor.AndorSDK2.AndorSDK2Camera` method), 472  
`start_acquisition()` (`pylablib.devices.Andor.AndorSDK3.AndorSDK3Camera` method), 482  
`start_acquisition()` (`pylablib.devices.Basler.pylon.BaslerPylonCamera` method), 516  
`start_acquisition()` (`pylablib.devices.BitFlow.BitFlow.BitFlowCamera` method), 532  
`start_acquisition()` (`pylablib.devices.BitFlow.BitFlow.BitFlowFrameGrabber` method), 523  
`start_acquisition()` (`pylablib.devices.DCAM.DCAM.DCAMCamera` method), 552  
`start_acquisition()` (`pylablib.devices.IMAQ.IMAQ.IMAQCamera` method), 575  
`start_acquisition()` (`pylablib.devices.IMAQ.IMAQ.IMAQFrameGrabber` method), 565  
`start_acquisition()` (`pylablib.devices.IMAQdx.IMAQdx.EthernetIMAQdxCamera` method), 589  
`start_acquisition()` (`pylablib.devices.IMAQdx.IMAQdx.IMAQdxCamera` method), 580  
`start_acquisition()` (`pylablib.devices.interface.camera.IAttributeCamera` method), 841  
`start_acquisition()` (`pylablib.devices.interface.camera.IBinROICamera` method), 860  
`start_acquisition()` (`pylablib.devices.interface.camera.ICamera` method), 831  
`start_acquisition()` (`pylablib.devices.interface.camera.IExposureCamera` method), 851  
`start_acquisition()` (`pylablib.devices.interface.camera.IGrabberAttributeCamera` method), 846  
`start_acquisition()` (`pylablib.devices.interface.camera.IROICamera` method), 856  
`start_acquisition()` (`pylablib.devices.PCO.SC2.PCOSC2Camera` method), 647  
`start_acquisition()` (`pylablib.devices.Photometrics.pvcam.PvcamCamera` method), 665  
`start_acquisition()` (`pylablib.devices.PhotonFocus.PhotonFocus.IPhotonFocusCamera` method), 677  
`start_acquisition()` (`pylablib.devices.PhotonFocus.PhotonFocus.PhotonFocusBitFlowCamera` method), 700  
`start_acquisition()` (`pylablib.devices.PhotonFocus.PhotonFocus.PhotonFocusIMAQCamera` method), 686  
`start_acquisition()` (`pylablib.devices.PhotonFocus.PhotonFocus.PhotonFocusSiSoCamera` method), 693  
`start_acquisition()` (`pylablib.devices.PrincetonInstruments.picam.PicamCamera` method), 714  
`start_acquisition()` (`pylablib.devices.SiliconSoftware.fgrab.SiliconSoftwareCamera` method), 733  
`start_acquisition()` (`pylablib.devices.SiliconSoftware.fgrab.SiliconSoftwareFrameGrabber` method), 723  
`start_acquisition()` (`pylablib.devices.Thorlabs.TLCamera.ThorlabsTLCamera` method), 771  
`start_acquisition()` (`pylablib.devices.Thorlabs.TLCamera.ThorlabsTLCamera` method), 771

<code>lablib.devices.uc480.uc480.UC480Camera</code>	<code>method</code> ), 220
<code>method</code> ), 868	<code>start_timer()</code> (py-
<code>start_batch_job()</code> (py-	<code>lablib.core.gui.widgets.container.QFrameContainer</code>
<code>lablib.core.thread.controller.QTaskThread</code>	<code>method</code> ), 216
<code>method</code> ), 304	<code>start_timer()</code> (py-
<code>start_bk (pylablib.devices.Thorlabs.kinesis.TKCubeTrigPosParam</code>	<code>lablib.core.gui.widgets.container.QGroupBoxContainer</code>
<code>attribute</code> ), 786	<code>method</code> ), 224
<code>start_degas()</code> (py-	<code>start_timer()</code> (py-
<code>lablib.devices.Leybold.base.ITR90</code>	<code>method</code> ),
609	<code>lablib.core.gui.widgets.container.QScrollAreaContainer</code>
<code>start_fast_scan()</code> (py-	<code>method</code> ), 230
<code>lablib.devices.M2.solstis.Solstis</code>	<code>method</code> ),
618	<code>lablib.core.gui.widgets.container.QScrollAreaContainer.QContai</code>
<code>start_fill()</code> (pylablib.devices.Cryomagnetics.base.LM500	<code>method</code> ), 228
<code>method</code> ), 541	<code>start_timer()</code> (py-
<code>start_fill()</code> (pylablib.devices.Cryomagnetics.base.LM510	<code>method</code> ), 233
<code>method</code> ), 547	<code>start_timer()</code> (py-
<code>start_fw (pylablib.devices.Thorlabs.kinesis.TKCubeTrigPosParam</code>	<code>lablib.core.gui.widgets.container.QWidgetContainer</code>
<code>attribute</code> ), 786	<code>method</code> ), 212
<code>start_link()</code> (pylablib.devices.M2.solstis.Solstis	<code>start_timer()</code> (py-
<code>method</code> ), 615	<code>lablib.core.gui.widgets.param_table.ParamTable</code>
<code>start_loop()</code> (pylablib.devices.Andor.AndorSDK3.AndorSDK3Camera	<code>method</code> ), 212
<code>method</code> ), 481	<code>start_timer()</code> (py-
<code>start_loop()</code> (pylablib.devices.Basler.pylon.BaslerPylonCamera	<code>lablib.core.gui.widgets.param_table.StatusTable</code>
<code>method</code> ), 516	<code>method</code> ), 260
<code>start_loop()</code> (pylablib.devices.BitFlow.BitFlow.BitFlowCamera	<code>method</code> ), 312
<code>method</code> ), 527	<code>lablib.core.thread.controller.QTaskThread</code>
<code>start_loop()</code> (pylablib.devices.BitFlow.BitFlow.BitFlowCamera	<code>attribute</code> ), 312
<code>method</code> ), 523	<code>FrameContainer (pylablib.devices</code>
<code>method</code> ), 523	<code>FrameContainer (pylablib.devices</code>
<code>start_loop()</code> (pylablib.devices.PhotonFocus.PhotonFocus	<code>attribute</code> ), 293
<code>method</code> ), 695	<code>PhotonFocus.PhotonFocus</code>
<code>start_measurement()</code> (py-	<code>PhotonFocus.PhotonFocus</code>
<code>lablib.devices.Cryomagnetics.base.LM500</code>	<code>PhotonFocus.PhotonFocus</code>
<code>method</code> ), 540	<code>PhotonFocus.PhotonFocus</code>
<code>start_measurement()</code> (py-	<code>status (pylablib.devices.PCO.SC2.TCameraStatus</code>
<code>lablib.devices.Cryomagnetics.base.LM510</code>	<code>attribute</code> ), 644
<code>method</code> ), 547	<code>status (pylablib.devices.Trinamic.base.TMCM1110.ReplyData</code>
<code>start_measurement()</code> (py-	<code>attribute</code> ), 827
<code>lablib.devices.HighFinesse.wlm.WLM</code>	<code>status (pylablib.devices.uc480.uc480.TCameraInfo</code>
557	<code>attribute</code> ), 864
<code>start_terascan()</code> (py-	<code>status_bits (pylablib.devices.Attocube.anc350.ANC350</code>
<code>lablib.devices.M2.solstis.Solstis</code>	<code>attribute</code> ), 507
617	<code>status_bits (pylablib.devices.Thorlabs.kinesis.KinesisDevice</code>
<code>start_timer()</code> (py-	<code>attribute</code> ), 787
<code>lablib.core.gui.widgets.container.IQContainer</code>	<code>status_bits (pylablib.devices.Thorlabs.kinesis.KinesisMotor</code>
<code>method</code> ), 200	<code>attribute</code> ), 800
<code>start_timer()</code> (py-	<code>status_bits (pylablib.devices.Thorlabs.kinesis.KinesisPiezoMotor</code>
<code>lablib.core.gui.widgets.container.IQWidgetContainer</code>	<code>attribute</code> ), 805
<code>method</code> ), 208	<code>status_bits (pylablib.devices.Thorlabs.kinesis.MFF</code>
<code>start_timer()</code> (py-	<code>attribute</code> ), 794
<code>lablib.core.gui.widgets.container.QContainer</code>	<code>StatusLineChecker (class in py-</code>
<code>method</code> ), 204	<code>lablib.devices.interface.camera</code> ), 861
<code>start_timer()</code> (py-	<code>StatusLineChecker (class in py-</code>
<code>lablib.core.gui.widgets.container.QDialogContainer</code>	<code>lablib.devices.PCO.SC2</code> ), 653
<code>method</code> ), 702	<code>StatusLineChecker (class in py-</code>
	<code>lablib.devices.PhotonFocus.PhotonFocus</code> ),



StatusTable (class in py-	stop () (pylablib.core.thread.controller.QTaskThread
lablib.core.gui.widgets.param_table), 251	method), 312
StatusTable.ParamRow (class in py-	stop () (pylablib.core.thread.controller.QThreadController
lablib.core.gui.widgets.param_table), 252	method), 299
step (pylablib.devices.DCAM.DCAM.DCAMAttribute	stop () (pylablib.core.utils.general.Countdown
attribute), 549	method), 378
step_bk (pylablib.devices.Thorlabs.kinesis.TKCubeTrigPosParam	stop () (pylablib.devices.Arcus.performax.Performax2EXStage
attribute), 786	method), 497
step_fw (pylablib.devices.Thorlabs.kinesis.TKCubeTrigPosParam	stop () (pylablib.devices.Arcus.performax.Performax4EXStage
attribute), 786	method), 491
step_size (pylablib.devices.Thorlabs.kinesis.TJogParams	stop () (pylablib.devices.Arcus.performax.PerformaxDMX.JSASage
attribute), 785	method), 498
step_size_bk (pylablib.devices.Thorlabs.kinesis.TPZMotorJogParam	stop () (pylablib.devices.Attocube.anc300.ANC300
attribute), 787	method), 505
step_size_fw (pylablib.devices.Thorlabs.kinesis.TPZMotorJogParam	stop () (pylablib.devices.Attocube.anc350.ANC350
attribute), 787	method), 509
step_voltage () (py-	stop () (pylablib.devices.Newport.picomotor.Picomotor8742
lablib.devices.OZOptics.base.EPC04 method),	method), 630
637	stop () (pylablib.devices.NI.daq.NIDAQ method), 623
stitched_scan () (py-	stop () (pylablib.devices.PhysikInstrumente.base.PIE516
lablib.devices.Sirah.tuner.MatisseTuner	method), 706
method), 744	stop () (pylablib.devices.SmarAct.scu3d.SCU3D
stop (pylablib.core.dataproc.utils.Range attribute), 136	method), 747
stop (pylablib.core.gui.widgets.container.TTimerEvent	stop () (pylablib.devices.Thorlabs.kinesis.KinesisMotor
attribute), 200	method), 796
stop () (in module pylablib.core.thread.profile), 317	stop () (pylablib.devices.Thorlabs.kinesis.KinesisPiezoMotor
stop () (pylablib.core.devio.backend_logger.BackendLogger	method), 802
method), 142	stop () (pylablib.devices.Trinamic.base.TMCM1110
stop () (pylablib.core.devio.comm_backend.RecordedDeviceBackend	method), 828
method), 158	stop_acquisition () (py-
stop () (pylablib.core.gui.widgets.container.IQContainer	lablib.devices.AlliedVision.Bonito.BonitoIMAQCamera
method), 201	method), 464
stop () (pylablib.core.gui.widgets.container.IQWidgetContainer	stop_acquisition () (py-
method), 208	lablib.devices.AlliedVision.Bonito.IBonitoCamera
stop () (pylablib.core.gui.widgets.container.QContainer	method), 456
method), 204	stop_acquisition () (py-
stop () (pylablib.core.gui.widgets.container.QDialogContainer	lablib.devices.Andor.AndorSDK2.AndorSDK2Camera
method), 220	method), 472
stop () (pylablib.core.gui.widgets.container.QFrameContainer	stop_acquisition () (py-
method), 216	lablib.devices.Andor.AndorSDK3.AndorSDK3Camera
stop () (pylablib.core.gui.widgets.container.QGroupBoxContainer	method), 482
method), 224	stop_acquisition () (py-
stop () (pylablib.core.gui.widgets.container.QScrollAreaContainer	lablib.devices.Basler.pylon.BaslerPylonCamera
method), 231	method), 516
stop () (pylablib.core.gui.widgets.container.QScrollAreaContainer.QContainedWidget	stop_acquisition () (py-
method), 228	lablib.devices.BitFlow.BitFlow.BitFlowCamera
stop () (pylablib.core.gui.widgets.container.QTabContainer	method), 532
method), 233	stop_acquisition () (py-
stop () (pylablib.core.gui.widgets.container.QWidgetContainer	lablib.devices.BitFlow.BitFlow.BitFlowFrameGrabber
method), 212	method), 523
stop () (pylablib.core.gui.widgets.param_table.ParamTable	stop_acquisition () (py-
method), 251	lablib.devices.DCAM.DCAM.DCAMCamera
stop () (pylablib.core.gui.widgets.param_table.StatusTable	method), 552
method), 260	stop_acquisition () (py-

`lablib.devices.IMAQ.IMAQ.IMAQCamera`  
`method)`, 575

`stop_acquisition()` (py-  
`lablib.devices.IMAQ.IMAQ.IMAQFrameGrabber`  
`method)`, 565

`stop_acquisition()` (py-  
`lablib.devices.IMAQdx.IMAQdx.EthernetIMAQdxCamera`  
`method)`, 589

`stop_acquisition()` (py-  
`lablib.devices.IMAQdx.IMAQdx.IMAQdxCamera`  
`method)`, 580

`stop_acquisition()` (py-  
`lablib.devices.interface.camera.IAttributeCamera`  
`method)`, 841

`stop_acquisition()` (py-  
`lablib.devices.interface.camera.IBinROICamera`  
`method)`, 860

`stop_acquisition()` (py-  
`lablib.devices.interface.camera.ICamera`  
`method)`, 832

`stop_acquisition()` (py-  
`lablib.devices.interface.camera.IExposureCamera`  
`method)`, 851

`stop_acquisition()` (py-  
`lablib.devices.interface.camera.IGrabberAttributeCamera`  
`method)`, 846

`stop_acquisition()` (py-  
`lablib.devices.interface.camera.IROICamera`  
`method)`, 856

`stop_acquisition()` (py-  
`lablib.devices.PCO.SC2.PCOSC2Camera`  
`method)`, 647

`stop_acquisition()` (py-  
`lablib.devices.Photometrics.pvcam.PvcamCamera`  
`method)`, 665

`stop_acquisition()` (py-  
`lablib.devices.PhotonFocus.PhotonFocus.IPhotonFocusCamera`  
`method)`, 677

`stop_acquisition()` (py-  
`lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusBitFlowCamera`  
`method)`, 701

`stop_acquisition()` (py-  
`lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusIMAQCamera`  
`method)`, 686

`stop_acquisition()` (py-  
`lablib.devices.PhotonFocus.PhotonFocus.PhotonFocusSiSoCamera`  
`method)`, 694

`stop_acquisition()` (py-  
`lablib.devices.PrincetonInstruments.picam.PicamCamera`  
`method)`, 714

`stop_acquisition()` (py-  
`lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareCamera`  
`method)`, 733

`stop_acquisition()` (py-  
`lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareFrameGrabber`  
`method)`, 723

`stop_acquisition()` (py-  
`lablib.devices.Thorlabs.TLCamera.ThorlabsTLCamera`  
`method)`, 772

`stop_acquisition()` (py-  
`lablib.devices.uc480.uc480.UC480Camera`  
`method)`, 868

`stop_all_controllers()` (in module py-  
`lablib.core.thread.controller)`, 315

`stop_all_operation()` (py-  
`lablib.devices.M2.solstis.Solstis`  
`method)`, 619

`stop_app()` (in module py-  
`lablib.core.thread.controller)`, 315

`stop_batch_job()` (py-  
`lablib.core.thread.controller.QTaskThread`  
`method)`, 304

`stop_coarse_tuning()` (py-  
`lablib.devices.M2.solstis.Solstis`  
`method)`, 616

`stop_controller()` (in module py-  
`lablib.core.thread.controller)`, 314

`stop_degas()` (pylablib.devices.Leybold.base.ITR90  
`method)`, 609

`stop_fast_scan()` (py-  
`lablib.devices.M2.solstis.Solstis`  
`method)`, 619

`stop_fine_tuning()` (py-  
`lablib.devices.M2.solstis.Solstis`  
`method)`, 616

`stop_grabbing()` (py-  
`lablib.devices.Tektronix.base.DPO2000`  
`method)`, 768

`stop_grabbing()` (py-  
`lablib.devices.Tektronix.base.ITektronixScope`  
`method)`, 749

`stop_grabbing()` (py-  
`lablib.devices.Tektronix.base.TDS2000`  
`method)`, 761

`stop_loop()` (pylablib.devices.Andor.AndorSDK3.AndorSDK3Camera.L  
`method)`, 481

`stop_loop()` (pylablib.devices.Basler.pylon.BaslerPylonCamera.Schedu  
`method)`, 516

`stop_loop()` (pylablib.devices.BitFlow.BitFlow.BitFlowCamera.BufferM  
`method)`, 527

`stop_loop()` (pylablib.devices.BitFlow.BitFlow.BitFlowFrameGrabber.L  
`method)`, 523

`stop_loop()` (pylablib.devices.PhotonFocus.PhotonFocus.PhotonFocusL  
`method)`, 695

`stop_lower()` (pylablib.devices.Sirah.Matisse.TSscanMode  
`attribute)`, 736

`stop_measurement()` (py-  
`lablib.devices.HighFinesse.wlm.WLM`  
`method)`,

557  
 stop\_mode (pylablib.devices.Thorlabs.kinesis.TJogParams attribute), 785  
 stop\_scan\_web () (pylablib.devices.M2.solstis.Solstis method), 619  
 stop\_terascan () (pylablib.devices.M2.solstis.Solstis method), 618  
 stop\_timer () (pylablib.core.gui.widgets.container.IQContainer method), 200  
 stop\_timer () (pylablib.core.gui.widgets.container.IQWidgetContainer method), 208  
 stop\_timer () (pylablib.core.gui.widgets.container.QContainer method), 204  
 stop\_timer () (pylablib.core.gui.widgets.container.QDialogContainer method), 220  
 stop\_timer () (pylablib.core.gui.widgets.container.QFrameContainer method), 216  
 stop\_timer () (pylablib.core.gui.widgets.container.QGroupBoxContainer method), 224  
 stop\_timer () (pylablib.core.gui.widgets.container.QScrollAreaContainer method), 231  
 stop\_timer () (pylablib.core.gui.widgets.container.QScrollAreaContainer method), 228  
 stop\_timer () (pylablib.core.gui.widgets.container.QTabContainer method), 233  
 stop\_timer () (pylablib.core.gui.widgets.container.QWidgetContainer method), 212  
 stop\_timer () (pylablib.core.gui.widgets.param\_table.ParamTable method), 251  
 stop\_timer () (pylablib.core.gui.widgets.param\_table.StatusTable method), 260  
 stop\_upper (pylablib.devices.Sirah.Matisse.TScanMode attribute), 736  
 store\_axis\_parameter () (pylablib.devices.Trinamic.base.TMCM1110 method), 827  
 store\_defaults () (pylablib.devices.Arcus.performax.GenericPerformaxStage method), 488  
 store\_defaults () (pylablib.devices.Arcus.performax.Performax2EXStage method), 497  
 store\_defaults () (pylablib.devices.Arcus.performax.Performax4EXStage method), 493  
 store\_defaults () (pylablib.devices.Arcus.performax.PerformaxDMXJSStage method), 500  
 store\_parameters () (pylablib.devices.Thorlabs.elliptec.ElliptecMotor method), 779  
 store\_settings () (pylablib.devices.Thorlabs.serial.FW method), 813  
 str\_to\_float () (in module pylablib.core.gui.formatter), 261  
 StrDumper (class in pylablib.core.utils.strdump), 395  
 StreamFileLogger (class in pylablib.core.utils.general), 379  
 strerror (pylablib.core.thread.threadprop.TimeoutThreadError attribute), 320  
 StreamReader (pylablib.core.utils.net.SocketError attribute), 388  
 StreamWriter (pylablib.core.utils.net.SocketTimeout attribute), 388  
 strindex (pylablib.devices.Andor.AndorSDK3.TFrameInfo attribute), 479  
 strcasecmp () (in module pylablib.core.utils.string), 396  
 str\_idx () (in module pylablib.core.utils.indexing), 381  
 str\_idx () (class in pylablib.core.utils.string), 397  
 strprep () (in module pylablib.core.utils.ctypes\_wrap), 324  
 subcolumn () (pylablib.core.dataproc.table\_wrap.Array1DWrapper method), 126  
 subcolumn () (pylablib.core.dataproc.table\_wrap.I1DWrapper method), 126  
 subcolumn () (pylablib.core.dataproc.table\_wrap.Series1DWrapper method), 128  
 subscribe\_commsync () (pylablib.core.thread.controller.QTaskThread method), 306  
 subscribe\_direct () (pylablib.core.thread.controller.QTaskThread method), 312  
 subscribe\_direct () (pylablib.core.thread.controller.QThreadController method), 296  
 subscribe\_direct () (pylablib.core.thread.multicast\_pool.MulticastPool method), 315  
 subscribe\_sync () (pylablib.core.thread.controller.QTaskThread method), 312  
 subscribe\_sync () (pylablib.core.thread.controller.QThreadController method), 295  
 subtable () (pylablib.core.dataproc.table\_wrap.Array2DWrapper method), 131  
 subtable () (pylablib.core.dataproc.table\_wrap.DataFrame2DWrapper method), 133  
 subtable () (pylablib.core.dataproc.table\_wrap.I2DWrapper method), 129  
 subtract\_baseline () (in module pylablib.core.dataproc.feature), 108

[success\\_wait\(\)](#) ([pylablib.core.thread.callsync.QCallResultSynchronizer](#) method), 282
 [sync\\_variable\(\)](#) ([pylablib.core.thread.controller.QTaskThread](#) method), 313

[success\\_wait\(\)](#) ([pylablib.core.thread.callsync.QDirectResultSynchronizer](#) method), 283
 [sync\\_variable\(\)](#) ([pylablib.core.thread.controller.QThreadController](#) method), 299

[success\\_wait\(\)](#) ([pylablib.core.thread.notifier.ISkippableNotifier](#) method), 317
 [sync\\_with\\_ai\(\)](#) ([pylablib.devices.NI.daq.TVoltageOutputClockParameter](#) attribute), 621

[success\\_wait\(\)](#) ([pylablib.core.thread.synchronizing.QThreadNotifier](#) method), 318
 [system\(\)](#) ([pylablib.devices.Photometrics.pvcam.TDeviceInfo](#) attribute), 661

[sum\(\)](#) ([pylablib.devices.Thorlabs.kinesis.TQuadDetectorReadings](#) attribute), 805
 [system\\_info\(\)](#) ([pylablib.devices.SiliconSoftware.fgrab.TDeviceInfo](#) attribute), 720

[sw\\_kind\(\)](#) ([pylablib.devices.Thorlabs.kinesis.TLimitSwitchParams](#) attribute), 785
 [table\\_entry\\_builder\(\)](#) ([in module py-lablib.core.fileio.dict\\_entry](#)), 174

[sw\\_position\\_ccw](#) ([pylablib.devices.Thorlabs.kinesis.TLimitSwitchParams](#) attribute), 785
 [TableBinaryOutputFileFormat](#) ([class in py-lablib.core.fileio.savefile](#)), 193

[sw\\_position\\_cw](#) ([pylablib.devices.Thorlabs.kinesis.TLimitSwitchParams](#) attribute), 785
 [TableStreamFile](#) ([class in py-lablib.core.fileio.table\\_stream](#)), 196

[swver](#) ([pylablib.devices.KJL.base.TKJL300DeviceInfo](#) attribute), 590
 [TAcqProgress](#) ([class in py-lablib.devices.Andor.AndorSDK2](#)), 466

[swver](#) ([pylablib.devices.Leybold.base.TDeviceInfo](#) attribute), 607
 [TAcqTimings](#) ([class in py-lablib.devices.interface.camera](#)), 846

[sync\\_controller\(\)](#) ([in module py-lablib.core.thread.controller](#)), 314
 [TAcquiredFramesStatus](#) ([class in py-lablib.devices.uc480.uc480](#)), 865

[sync\\_exec\\_point\(\)](#) ([pylablib.core.thread.controller.QTaskThread](#) method), 313
 [tag](#) ([pylablib.core.thread.multicast\\_pool.TMulticast](#) attribute), 315

[sync\\_exec\\_point\(\)](#) ([pylablib.core.thread.controller.QThreadController](#) method), 300
 [tags](#) ([pylablib.devices.SiliconSoftware.fgrab.TFullAppletInfo](#) attribute), 719

[sync\\_phase\(\)](#) ([pylablib.devices.AWG.generic.GenericAWG](#) method), 404
 [TAppletInfo](#) ([class in py-lablib.devices.SiliconSoftware.fgrab](#)), 718

[sync\\_phase\(\)](#) ([pylablib.devices.AWG.specific.Agilent33220A](#) method), 419
 [TAxisROIILimit](#) ([class in py-lablib.devices.interface.camera](#)), 851

[sync\\_phase\(\)](#) ([pylablib.devices.AWG.specific.Agilent33500](#) method), 413
 [TBoardInfo](#) ([class in py-lablib.devices.SiliconSoftware.fgrab](#)), 718

[sync\\_phase\(\)](#) ([pylablib.devices.AWG.specific.InstekAFG2000](#) method), 431
 [TCameraInfo](#) ([class in py-lablib.devices.Basler.pylon](#)), 511

[sync\\_phase\(\)](#) ([pylablib.devices.AWG.specific.InstekAFG2225](#) method), 425
 [TCameraInfo](#) ([class in py-lablib.devices.IMAQdx.IMAQdx](#)), 576

[sync\\_phase\(\)](#) ([pylablib.devices.AWG.specific.RigolDG1000](#) method), 444
 [TCameraInfo](#) ([class in py-lablib.devices.PhotonFocus.PhotonFocus](#)), 670

[sync\\_phase\(\)](#) ([pylablib.devices.AWG.specific.RSInstekAFG21000](#) method), 437
 [TCameraInfo](#) ([class in py-lablib.devices.uc480.uc480](#)), 864

[sync\\_phase\(\)](#) ([pylablib.devices.AWG.specific.TektronixAFG1000](#) method), 443
 [TCameraStatus](#) ([class in py-lablib.devices.PCO.SC2](#)), 644

[sync\\_stop\(\)](#) ([pylablib.core.thread.controller.QTaskThread](#) method), 313
 [TChild](#) ([class in py-lablib.core.gui.widgets.container](#)), 200

[sync\\_stop\(\)](#) ([pylablib.core.thread.controller.QThreadController](#) method), 299
 [TConversionClass](#) ([class in py-lablib.core.utils.string](#)), 398



TCycleTimings (class in <i>pylablib.devices.Andor.AndorSDK2</i> ), 465	py- TDeviceInfo (class in <i>pylablib.devices.Thorlabs.TLCamera</i> ), 769
TDefaultCallInfo (class in <i>pylablib.core.thread.callsync</i> ), 284	py- TDeviceInfo (class in <i>pylablib.devices.Toptica.ibeam</i> ), 823
TDeviceInfo (class in <i>pylablib.devices.AlliedVision.Bonito</i> ), 450	py- TDeviceInfo (class in <i>pylablib.devices.uc480.uc480</i> ), 865
TDeviceInfo (class in <i>pylablib.devices.Andor.AndorSDK2</i> ), 465	py- TDS2000 (class in <i>pylablib.devices.Tektronix.base</i> ), 755
TDeviceInfo (class in <i>pylablib.devices.Andor.AndorSDK3</i> ), 478	py- TDS2000.NoParameterCaller (class in <i>pylablib.devices.Tektronix.base</i> ), 755
TDeviceInfo (class in <i>pylablib.devices.Attocube.anc300</i> ), 502	py- TektronixAFG1000 (class in <i>pylablib.devices.AWG.specific</i> ), 438
TDeviceInfo (class in <i>pylablib.devices.Basler.pylon</i> ), 513	TektronixAFG1000.NoParameterCaller (class in <i>pylablib.devices.AWG.specific</i> ), 438
TDeviceInfo (class in <i>pylablib.devices.BitFlow.BitFlow</i> ), 521	TektronixBackendError, 748
TDeviceInfo (class in <i>pylablib.devices.DCAM.DCAM</i> ), 550	TektronixError, 748
TDeviceInfo (class in <i>pylablib.devices.HighFinesse.wlm</i> ), 557	temp_library_parameters() (in module <i>pylablib.core.utils.library_parameters</i> ), 385
TDeviceInfo (class in <i>pylablib.devices.IMAQ.IMAQ</i> ), 561	temp_open() ( <i>pylablib.devices.Andor.AndorSDK2.LibraryController</i> method), 465
TDeviceInfo (class in <i>pylablib.devices.IMAQdx.IMAQdx</i> ), 578	temp_open() ( <i>pylablib.devices.Andor.AndorSDK3.LibraryController</i> method), 476
TDeviceInfo (class in <i>pylablib.devices.LaserQuantum.base</i> ), 603	temp_open() ( <i>pylablib.devices.Basler.pylon.LibraryController</i> method), 511
TDeviceInfo (class in <i>pylablib.devices.Leybold.base</i> ), 606	temp_open() ( <i>pylablib.devices.DCAM.DCAM.LibraryController</i> method), 548
TDeviceInfo (class in <i>pylablib.devices.LighthousePhotonics.base</i> ), 611	temp_open() ( <i>pylablib.devices.Photometrics.pvcam.LibraryController</i> method), 659
TDeviceInfo (class in <i>pylablib.devices.Newport.picomotor</i> ), 627	temp_open() ( <i>pylablib.devices.PhotonFocus.PhotonFocus.LibraryController</i> method), 669
TDeviceInfo (class in <i>pylablib.devices.NI.daq</i> ), 620	temp_open() ( <i>pylablib.devices.PrincetonInstruments.picam.LibraryController</i> method), 708
TDeviceInfo (class in <i>pylablib.devices.Ophir.base</i> ), 640	temp_open() ( <i>pylablib.devices.SmarAct.scu3d.LibraryController</i> method), 745
TDeviceInfo (class in <i>pylablib.devices.PCO.SC2</i> ), 644	temp_open() ( <i>pylablib.devices.Thorlabs.TLCamera.LibraryController</i> method), 769
TDeviceInfo (class in <i>pylablib.devices.Photometrics.pvcam</i> ), 661	temp_open() ( <i>pylablib.devices.utils.load_lib.LibraryController</i> method), 875
TDeviceInfo (class in <i>pylablib.devices.PhotonFocus.PhotonFocus</i> ), 671	TempFile (class in <i>pylablib.core.utils.files</i> ), 362
TDeviceInfo (class in <i>pylablib.devices.PrincetonInstruments.picam</i> ), 711	test_columns_line() (in module <i>pylablib.core.fileio.loadfile_utils</i> ), 182
TDeviceInfo (class in <i>pylablib.devices.SiliconSoftware.fgrab</i> ), 720	test_row_type() (in module <i>pylablib.core.fileio.loadfile_utils</i> ), 182
TDeviceInfo (class in <i>pylablib.devices.SmarAct.scu3d</i> ), 745	test_savetime_comment() (in module <i>pylablib.core.fileio.loadfile_utils</i> ), 182
TDeviceInfo (class in <i>pylablib.devices.Thorlabs.elliptec</i> ), 776	TextEdit (class in <i>pylablib.core.gui.widgets.edit</i> ), 233
TDeviceInfo (class in <i>pylablib.devices.Thorlabs.kinesis</i> ), 781	TextLabel (class in <i>pylablib.core.gui.widgets.label</i> ), 236
	TF100 (class in <i>pylablib.devices.OZOptics.base</i> ), 633
	TF100.NoParameterCaller (class in <i>pylablib.devices.OZOptics.base</i> ), 633
	TFastpiezoCtlParameters (class in <i>pylablib.devices.Sirah.Matisse</i> ), 735
	TFlipperParameters (class in <i>pylablib.devices.Thorlabs.kinesis</i> ), 790

TFrameInfo	(class in py-lablib.devices.Andor.AndorSDK3), 479	lablib.devices.Sirah.Matisse.SirahMatisse method), 738
TFrameInfo	(class in pylablib.devices.DCAM.DCAM), 550	thinet_wait_move() (py-lablib.devices.Sirah.Matisse.SirahMatisse method), 737
TFrameInfo	(class in py-lablib.devices.interface.camera), 831	THomeParams (class in py-lablib.devices.Thorlabs.kinesis), 785
TFrameInfo	(class in pylablib.devices.PCO.SC2), 645	ThorlabsBackendError, 776
TFrameInfo	(class in py-lablib.devices.Photometrics.pvcam), 661	ThorlabsError, 776
TFrameInfo	(class in py-lablib.devices.PrincetonInstruments.picam), 711	ThorlabsSerialInterface (class in py-lablib.devices.Thorlabs.serial), 809
TFrameInfo	(class in py-lablib.devices.SiliconSoftware.fgrab), 720	ThorlabsSerialInterface.NoParameterCaller (class in pylablib.devices.Thorlabs.serial), 810
TFrameInfo	(class in py-lablib.devices.Thorlabs.TLCamera), 770	ThorlabsTimeoutError, 776
TFrameInfo	(class in pylablib.devices.uc480.uc480), 866	ThorlabsTLCamera (class in py-lablib.devices.Thorlabs.TLCamera), 770
TFramePosition	(class in py-lablib.devices.interface.camera), 830	ThorlabsTLCamera.NoParameterCaller (class in pylablib.devices.Thorlabs.TLCamera), 773
TFrameSize	(class in py-lablib.devices.interface.camera), 830	ThorlabsTLCamera.RingBuffer (class in py-lablib.devices.Thorlabs.TLCamera), 770
TFramesStatus	(class in py-lablib.devices.interface.camera), 830	ThreadError, 319
TFullAppletInfo	(class in py-lablib.devices.SiliconSoftware.fgrab), 718	time_left() (pylablib.core.thread.controller.QTaskThread.Job method), 303
TGenMoveParams	(class in py-lablib.devices.Thorlabs.kinesis), 785	time_left() (pylablib.core.utils.general.Countdown method), 378
THeadInfo	(class in pylablib.devices.Ophir.base), 640	time_left() (pylablib.core.utils.general.Timer method), 379
thinet_clear_errors()	(py-lablib.devices.Sirah.Matisse.SirahMatisse method), 737	time_passed() (py-lablib.core.utils.general.Countdown method), 379
thinet_get_position()	(py-lablib.devices.Sirah.Matisse.SirahMatisse method), 737	TimeoutError (pylablib.devices.AlliedVision.Bonito.BonitoIMAQCamera attribute), 456
thinet_get_range()	(py-lablib.devices.Sirah.Matisse.SirahMatisse method), 737	TimeoutError (pylablib.devices.AlliedVision.Bonito.IBonitoCamera attribute), 452
thinet_get_status()	(py-lablib.devices.Sirah.Matisse.SirahMatisse method), 737	TimeoutError (pylablib.devices.Andor.AndorSDK2.AndorSDK2Camera attribute), 467
thinet_get_status_n()	(py-lablib.devices.Sirah.Matisse.SirahMatisse method), 737	TimeoutError (pylablib.devices.Andor.AndorSDK3.AndorSDK3Camera attribute), 479
thinet_home()	(py-lablib.devices.Sirah.Matisse.SirahMatisse method), 738	TimeoutError (pylablib.devices.Basler.pylon.BaslerPylonCamera attribute), 514
thinet_is_moving()	(py-lablib.devices.Sirah.Matisse.SirahMatisse method), 737	TimeoutError (pylablib.devices.BitFlow.BitFlow.BitFlowCamera attribute), 527
thinet_move_to()	(py-lablib.devices.Sirah.Matisse.SirahMatisse method), 737	TimeoutError (pylablib.devices.BitFlow.BitFlow.BitFlowFrameGrabber attribute), 521
thinet_stop()	(py-lablib.devices.Sirah.Matisse.SirahMatisse method), 737	TimeoutError (pylablib.devices.DCAM.DCAM.DCAMCamera attribute), 550
		TimeoutError (pylablib.devices.IMAQ.IMAQ.IMAQCamera attribute), 569
		TimeoutError (pylablib.devices.IMAQ.IMAQ.IMAQFrameGrabber attribute), 561
		TimeoutError (pylablib.devices.IMAQdx.IMAQdx.EthernetIMAQdxCam attribute), 584
		TimeoutError (pylablib.devices.IMAQdx.IMAQdx.IMAQdxCamera attribute), 578

TimeoutError (pylablib.devices.interface.camera.IAttributeCameraattribute), 213  
     attribute), 838      TimerUIDGenerator (py-  
 TimeoutError (pylablib.devices.interface.camera.IBinROICameralablib.core.gui.widgets.container.QGroupBoxContainer  
     attribute), 856      attribute), 221  
 TimeoutError (pylablib.devices.interface.camera.ICameraattribute), 831      TimerUIDGenerator (py-  
     attribute), 831      lablib.core.gui.widgets.container.QScrollAreaContainer  
 TimeoutError (pylablib.devices.interface.camera.IExposureCameraattribute), 229  
     attribute), 847      TimerUIDGenerator (py-  
 TimeoutError (pylablib.devices.interface.camera.IGrabberAttributeCinecamlablib.core.gui.widgets.container.QScrollAreaContainer.QContain  
     attribute), 842      attribute), 225  
 TimeoutError (pylablib.devices.interface.camera.IROICameraattribute), 852      TimerUIDGenerator (py-  
     attribute), 852      lablib.core.gui.widgets.container.QTabContainer  
 TimeoutError (pylablib.devices.PCO.SC2.PCOSC2Cameraattribute), 231  
     attribute), 645      TimerUIDGenerator (py-  
 TimeoutError (pylablib.devices.Photometrics.pvcam.PvcamCameraattribute), 209  
     attribute), 662      attribute), 209  
 TimeoutError (pylablib.devices.PhotonFocus.PhotonFocus.IPhotonFocusCameraattribute), 674      TimerUIDGenerator (py-  
     attribute), 674      lablib.core.gui.widgets.param\_table.ParamTable  
 TimeoutError (pylablib.devices.PhotonFocus.PhotonFocus.PhotonFocusBinFLWCameraattribute), 695      TimerUIDGenerator (py-  
     attribute), 695      attribute), 252  
 TimeoutError (pylablib.devices.PhotonFocus.PhotonFocus.PhotonFocusIMAQCameraattribute), 678      TimerUIDGenerator (py-  
     attribute), 678      attribute), 252  
 TimeoutError (pylablib.devices.PhotonFocus.PhotonFocus.PhotonFocusSiliconDevicesattribute), 687      TimerUIDGenerator (py-  
     attribute), 687      attribute), 721  
 TimeoutError (pylablib.devices.PrincetonInstruments.picam.PicamCameraattribute), 712      TimerUIDGenerator (py-  
     attribute), 712      attribute), 866  
 TimeoutError (pylablib.devices.SiliconSoftware.fgrab.SiliconSoftwareCameraattribute), 727      TimerUIDGenerator (py-  
     attribute), 727      attribute), 727  
 TimeoutError (pylablib.devices.SiliconSoftware.fgrab.SiliconSoftwareCameraGrabberattribute), 721      TimerUIDGenerator (py-  
     attribute), 721      attribute), 721  
 TimeoutError (pylablib.devices.Thorlabs.TLCamera.ThorlabsTLCameraattribute), 770      TimerUIDGenerator (py-  
     attribute), 770      attribute), 866  
 TimeoutError (pylablib.devices.uc480.uc480.UC480Cameraattribute), 867      TimerUIDGenerator (py-  
     attribute), 867      attribute), 712  
 TimeoutThreadError, 319      timestamp\_end\_ns (py-  
 Timer (class in pylablib.core.utils.general), 379      timestamp\_end\_ns (py-  
 timer (pylablib.core.gui.widgets.container.TTimerattribute), 199      timestamp\_end\_ns (py-  
     attribute), 199      attribute), 661  
 timer (pylablib.core.gui.widgets.container.TTimerEventattribute), 200      timestamp\_long (py-  
     attribute), 200      attribute), 721  
 TimerUIDGenerator (py-  
     lablib.core.gui.widgets.container.IQContainerattribute), 200      timestamp\_start (py-  
     attribute), 200      attribute), 712  
 TimerUIDGenerator (py-  
     lablib.core.gui.widgets.container.IQWidgetContainerattribute), 206      timestamp\_start\_ns (py-  
     attribute), 206      attribute), 661  
 TimerUIDGenerator (py-  
     lablib.core.gui.widgets.container.QContainerattribute), 202      timestamp\_us (pylablib.devices.DCAM.DCAM.TFrameInfo  
     attribute), 202      attribute), 550  
 TimerUIDGenerator (py-  
     lablib.core.gui.widgets.container.QDialogContainerattribute), 217      timing() (in module pylablib.core.utils.general), 380  
 TimerUIDGenerator (py-  
     lablib.core.gui.widgets.container.QFrameContainerattribute), 217      TInternalBufferStatus (class in py-  
     attribute), 217      lablib.devices.PCO.SC2), 644  
 TimerUIDGenerator (py-  
     lablib.core.gui.widgets.container.QFrameContainerattribute), 217      TITR90Status (class in py-  
     attribute), 217      lablib.devices.Leybold.base), 609

TJogParams (class in py-lablib.devices.Thorlabs.kinesis), 784  
 TKCubeTrigIOParams (class in py-lablib.devices.Thorlabs.kinesis), 785  
 TKCubeTrigPosParams (class in py-lablib.devices.Thorlabs.kinesis), 786  
 TKJL300DeviceInfo (class in py-lablib.devices.KJL.base), 590  
 TLakeshore218AnalogSettings (class in py-lablib.devices.Lakeshore.base), 592  
 TLakeshore218CurveHeader (class in py-lablib.devices.Lakeshore.base), 593  
 TLakeshore218FilterSettings (class in py-lablib.devices.Lakeshore.base), 593  
 TLakeshore370AnalogSettings (class in py-lablib.devices.Lakeshore.base), 598  
 TLakeshore370FilterSettings (class in py-lablib.devices.Lakeshore.base), 599  
 TLakeshore370RangeSettings (class in py-lablib.devices.Lakeshore.base), 598  
 TLibraryCloseResult (class in py-lablib.devices.utils.load\_lib), 874  
 TLibraryOpenResult (class in py-lablib.devices.utils.load\_lib), 874  
 TLimitSwitchParams (class in py-lablib.devices.Thorlabs.kinesis), 785  
 TMCM1110 (class in pylablib.devices.Trinamic.base), 826  
 TMCM1110.NoParameterCaller (class in py-lablib.devices.Trinamic.base), 829  
 TMCM1110.ReplyData (class in py-lablib.devices.Trinamic.base), 827  
 TMissedFramesStatus (class in py-lablib.devices.Andor.AndorSDK3), 479  
 TMotorInfo (class in py-lablib.devices.Thorlabs.elliptec), 777  
 TMulticast (class in py-lablib.core.thread.multicast\_pool), 315  
 to\_alias() (pylablib.core.devio.interface.EnumParameterClass method), 167  
 to\_alias() (pylablib.core.devio.interface.FunctionParameterClass method), 167  
 to\_alias() (pylablib.core.devio.interface.ICheckingParameterClass method), 164  
 to\_alias() (pylablib.core.devio.interface.IEnumParameterClass method), 165  
 to\_alias() (pylablib.core.devio.interface.RangeParameterClass method), 164  
 to\_alias() (pylablib.core.devio.interface.RangeParameterClass method), 164  
 to\_callable() (in module py-lablib.core.dataproc.callable), 108  
 to\_desc() (pylablib.core.devio.data\_format.DataFormat method), 162  
 to\_dict() (in module pylablib.core.utils.general), 375  
 to\_dict() (pylablib.core.fileio.dict\_entry.ExpandedContainerDictionaryEntry method), 176  
 to\_dict() (pylablib.core.fileio.dict\_entry.ExternalBinTableDictionaryEntry method), 174  
 to\_dict() (pylablib.core.fileio.dict\_entry.ExternalNumpyDictionaryEntry method), 176  
 to\_dict() (pylablib.core.fileio.dict\_entry.ExternalTextTableDictionaryEntry method), 173  
 to\_dict() (pylablib.core.fileio.dict\_entry.IDictionaryEntry method), 170  
 to\_dict() (pylablib.core.fileio.dict\_entry.IExternalFileDictionaryEntry method), 174  
 to\_dict() (pylablib.core.fileio.dict\_entry.IExternalTableDictionaryEntry method), 172  
 to\_dict() (pylablib.core.fileio.dict\_entry.InlineTableDictionaryEntry method), 172  
 to\_dict() (pylablib.core.fileio.dict\_entry.ITableDictionaryEntry method), 171  
 to\_double\_index() (in module py-lablib.core.utils.indexing), 382  
 to\_Pa() (pylablib.devices.Pfeiffer.base.TPG260 method), 654  
 to\_pairs\_list() (in module py-lablib.core.utils.general), 375  
 to\_path() (pylablib.core.fileio.location.LocationName method), 183  
 to\_predicate() (in module py-lablib.core.utils.general), 374  
 to\_range() (in module pylablib.core.utils.indexing), 381  
 to\_range() (in module pylablib.core.utils.string), 400  
 to\_string() (in module pylablib.core.utils.string), 399  
 to\_string() (pylablib.core.fileio.location.LocationName method), 183  
 to\_struct() (pylablib.core.utils.ctypes\_wrap.CStructWrapper method), 324  
 to\_value() (pylablib.core.devio.interface.EnumParameterClass method), 167  
 to\_value() (pylablib.core.devio.interface.FunctionParameterClass method), 167  
 to\_value() (pylablib.core.devio.interface.ICheckingParameterClass method), 164  
 to\_value() (pylablib.core.devio.interface.IEnumParameterClass method), 165  
 to\_value() (pylablib.core.devio.interface.RangeParameterClass method), 164  
 to\_value() (pylablib.core.devio.interface.RangeParameterClass method), 164  
 TToggleButton (class in py-lablib.core.gui.widgets.button), 197  
 ToolButtonValueHandler (class in py-lablib.core.gui.value\_handling), 274  
 tooltip (pylablib.devices.Basler.pylon.BaslerPylonAttribute attribute), 512  
 tooltip (pylablib.devices.IMAQdx.IMAQdx.IMAQdxAttribute attribute), 577



top (*pylablib.devices.interface.camera.TFramePosition* attribute), 831  
 toploopSlot() (in module *pylablib.core.thread.controller*), 292  
 topological\_order() (in module *pylablib.core.utils.general*), 376  
 TopticaBackendError, 823  
 TopticaError, 823  
 TopticaIBeam (class in *pylablib.devices.Toptica.ibeam*), 824  
 TopticaIBeam.NoParameterCaller (class in *pylablib.devices.Toptica.ibeam*), 825  
 touch() (in module *pylablib.core.utils.files*), 361  
 TPG260 (class in *pylablib.devices.Pfeiffer.base*), 654  
 TPG260.NoParameterCaller (class in *pylablib.devices.Pfeiffer.base*), 656  
 TPiezoetDriveParameters (class in *pylablib.devices.Sirah.Matisse*), 734  
 TPiezoetFeedbackParameters (class in *pylablib.devices.Sirah.Matisse*), 734  
 TPiezoetFeedforwardParameters (class in *pylablib.devices.Sirah.Matisse*), 734  
 TPipeMsg (class in *pylablib.core.utils.ipc*), 383  
 TPZMotorDriveParams (class in *pylablib.devices.Thorlabs.kinesis*), 786  
 TPZMotorJogParams (class in *pylablib.devices.Thorlabs.kinesis*), 787  
 TQuadDetectorOutputParams (class in *pylablib.devices.Thorlabs.kinesis*), 806  
 TQuadDetectorPIDParams (class in *pylablib.devices.Thorlabs.kinesis*), 805  
 TQuadDetectorReadings (class in *pylablib.devices.Thorlabs.kinesis*), 805  
 TQuadDetectorSetpoint (class in *pylablib.devices.Thorlabs.kinesis*), 805  
 TRangeInfo (class in *pylablib.devices.Ophir.base*), 641  
 transfer() (in module *pylablib.core.utils.rpyc\_utils*), 394  
 transfer() (*pylablib.core.utils.rpyc\_utils.DeviceService* method), 395  
 transfer() (*pylablib.core.utils.rpyc\_utils.SocketTunnelService* method), 394  
 transfer\_missed (pylablib.devices.uc480.uc480.TAcquiredFramesStatus attribute), 865  
 transit\_time (*pylablib.devices.Thorlabs.kinesis.TFlipperParameters* attribute), 790  
 translate\_string\_filter() (in module *pylablib.core.utils.string*), 397  
 travel (*pylablib.devices.Thorlabs.elliptec.TDeviceInfo* attribute), 777  
 TRawParameterValue (class in *pylablib.core.devio.interface*), 168  
 TReadoutInfo (class in *pylablib.devices.Photometrics.pvcam*), 661  
 TRefcellWaveformParameters (class in *pylablib.devices.Sirah.Matisse*), 735  
 trig1\_mode (*pylablib.devices.Thorlabs.kinesis.TKCubeTrigIOParams* attribute), 786  
 trig1\_pol (*pylablib.devices.Thorlabs.kinesis.TKCubeTrigIOParams* attribute), 786  
 trig2\_mode (*pylablib.devices.Thorlabs.kinesis.TKCubeTrigIOParams* attribute), 786  
 trig2\_pol (*pylablib.devices.Thorlabs.kinesis.TKCubeTrigIOParams* attribute), 786  
 trigger() (*pylablib.core.utils.general.Countdown* method), 378  
 trim\_frames() (in module *pylablib.devices.interface.camera*), 835  
 trim\_frames\_range() (pylablib.devices.interface.camera.FrameCounter method), 836  
 TrinamicBackendError, 826  
 TrinamicError, 826  
 TROIConstraints (class in *pylablib.devices.PrincetonInstruments.picam*), 709  
 truncate\_roi\_axis() (in module *pylablib.devices.interface.camera*), 851  
 truncate\_trace() (in module *pylablib.core.dataproc.fourier*), 117  
 truncate\_value() (pylablib.devices.Andor.AndorSDK3.AndorSDK3Attribute method), 478  
 truncate\_value() (pylablib.devices.Basler.pylon.BaslerPylonAttribute method), 513  
 truncate\_value() (pylablib.devices.IMAQdx.IMAQdx.IMAQdxAttribute method), 577  
 truncate\_value() (pylablib.devices.Photometrics.pvcam.PvcamAttribute method), 660  
 truncate\_value() (pylablib.devices.PhotonFocus.PhotonFocus.PFCamAttribute method), 671  
 truncate\_value() (pylablib.devices.PrincetonInstruments.picam.PicamAttribute method), 711  
 truncate\_value() (pylablib.devices.SiliconSoftware.fgrab.FGGrabAttribute method), 720  
 TScanMode (class in *pylablib.devices.Sirah.Matisse*), 735  
 TScanParameters (class in *pylablib.devices.Sirah.Matisse*), 736  
 TShmemVarDesc (class in *pylablib.core.utils.ipc*), 384

TSlowpiezoCtlParameters (class in py-  
     lablib.devices.Sirah.Matisse), 735  
 TStatusLine (class in py-  
     lablib.devices.AlliedVision.Bonito), 464  
 TStatusLine (class in pylablib.devices.PCO.SC2),  
     652  
 TStatusLineDescription (class in py-  
     lablib.devices.interface.camera), 860  
 TTemperatures (class in py-  
     lablib.devices.LaserQuantum.base), 604  
 TTemperatures (class in py-  
     lablib.devices.Toptica.ibeam), 824  
 TThinnetCtlParameters (class in py-  
     lablib.devices.Sirah.Matisse), 734  
 TTimer (class in pylablib.core.gui.widgets.container),  
     199  
 TTimerEvent (class in py-  
     lablib.core.gui.widgets.container), 199  
 TTimestamp (class in pylablib.devices.uc480.uc480),  
     865  
 TTPG260GaugeControlSettings (class in py-  
     lablib.devices.Pfeiffer.base), 654  
 TTPG260SwitchSettings (class in py-  
     lablib.devices.Pfeiffer.base), 653  
 TTriggerParameters (class in py-  
     lablib.devices.Tektronix.base), 748  
 tune\_etalon() (pylablib.devices.M2.solstis.Solstis  
     method), 616  
 tune\_laser\_resonator() (py-  
     lablib.devices.M2.solstis.Solstis method),  
     616  
 tune\_reference\_cavity() (py-  
     lablib.devices.M2.solstis.Solstis method),  
     617  
 tune\_to() (pylablib.devices.Sirah.tuner.MatisseTuner  
     method), 744  
 tunnel\_recv() (py-  
     lablib.core.utils.rpyc\_utils.DeviceService  
     method), 395  
 tunnel\_recv() (py-  
     lablib.core.utils.rpyc\_utils.SocketTunnelService  
     method), 394  
 tunnel\_send() (py-  
     lablib.core.utils.rpyc\_utils.DeviceService  
     method), 395  
 tunnel\_send() (py-  
     lablib.core.utils.rpyc\_utils.SocketTunnelService  
     method), 394  
 tup() (pylablib.core.dataproc.image.ROI method), 121  
 tup() (pylablib.core.dataproc.utils.Range method), 136  
 tup() (pylablib.core.utils.ctypes\_wrap.CStructWrapper  
     method), 325  
 tup() (pylablib.core.utils.indexing.Index method), 382  
 tup() (pylablib.core.utils.indexing.ListIndex method),  
     382  
 tup() (pylablib.core.utils.indexing.ListIndexNoSlice  
     method), 382  
 tup() (pylablib.core.utils.indexing.NumpyIndex  
     method), 382  
 tup\_struct() (pylablib.core.utils.ctypes\_wrap.CStructWrapper  
     class method), 325  
 TUpdateValue (class in py-  
     lablib.devices.Leybold.base), 607  
 TVelocityParams (class in py-  
     lablib.devices.Thorlabs.kinesis), 784  
 TVoltageOutputClockParameters (class in py-  
     lablib.devices.NI.daq), 621  
 TWavelengthInfo (class in py-  
     lablib.devices.Ophir.base), 640  
 TWidgetLocation (class in pylablib.core.gui.utils),  
     263  
 TWorkHours (class in py-  
     lablib.devices.LaserQuantum.base), 604  
 TWorkHours (class in py-  
     lablib.devices.LighthousePhotonics.base),  
     611  
 TWorkHours (class in pylablib.devices.Toptica.ibeam),  
     823  
 type (pylablib.devices.IMAQdx.IMAQdx.TCameraInfo  
     attribute), 576  
 type (pylablib.devices.Ophir.base.THeadInfo attribute),  
     640  
 type (pylablib.devices.PhotonFocus.PhotonFocus.TCameraInfo  
     attribute), 670

## U

UC480Camera (class in pylablib.devices.uc480.uc480),  
     866  
 UC480Camera.NoParameterCaller (class in py-  
     lablib.devices.uc480.uc480), 869  
 uid (pylablib.devices.SiliconSoftware.fgrab.TFullAppletInfo  
     attribute), 719  
 UIDGenerator (class in pylablib.core.utils.general),  
     377  
 unescape\_string() (in module py-  
     lablib.core.utils.string), 399  
 unittest\_result (py-  
     lablib.devices.utils.load\_lib.TLibraryCloseResult  
     attribute), 874  
 unique\_slices() (in module py-  
     lablib.core.dataproc.utils), 136  
 unit (pylablib.devices.DCAM.DCAM.DCAMAttribute  
     attribute), 549  
 units (pylablib.devices.Basler.pylon.BaslerPylonAttribute  
     attribute), 513  
 units (pylablib.devices.IMAQdx.IMAQdx.IMAQdxAttribute  
     attribute), 577  
 unity() (in module pylablib.core.utils.numerical), 392

[unload\\_all\(\) \(in module pylablib\), 875](#)  
[unload\\_package\\_modules\(\) \(in module pylablib.core.utils.module\), 386](#)  
[unlock\(\) \(pylablib.core.devio.comm\\_backend.FT232RLDeviceBackend\), 152](#)  
[unlock\(\) \(pylablib.core.devio.comm\\_backend.ICommBackendWrapper\), 161](#)  
[unlock\(\) \(pylablib.core.devio.comm\\_backend.IDeviceCommBackend\), 143](#)  
[unlock\(\) \(pylablib.core.devio.comm\\_backend.NetworkDeviceBackend\), 155](#)  
[unlock\(\) \(pylablib.core.devio.comm\\_backend.PyUSBDeviceBackend\), 157](#)  
[unlock\(\) \(pylablib.core.devio.comm\\_backend.RecordedDeviceBackend\), 159](#)  
[unlock\(\) \(pylablib.core.devio.comm\\_backend.SerialDeviceBackend\), 149](#)  
[unlock\(\) \(pylablib.core.devio.comm\\_backend.VisaDeviceBackend\), 145](#)  
[unlock\(\) \(pylablib.core.devio.SCPI.SCPIDevice\), 141](#)  
[unlock\(\) \(pylablib.devices.Arduino.base.IArduinoDevice\), 502](#)  
[unlock\(\) \(pylablib.devices.Attocube.anc300.ANC300\), 506](#)  
[unlock\(\) \(pylablib.devices.Attocube.anc350.ANC350\), 510](#)  
[unlock\(\) \(pylablib.devices.AWG.generic.GenericAWG\), 407](#)  
[unlock\(\) \(pylablib.devices.AWG.specific.Agilent33220A\), 419](#)  
[unlock\(\) \(pylablib.devices.AWG.specific.Agilent33500\), 413](#)  
[unlock\(\) \(pylablib.devices.AWG.specific.InstekAFG2000\), 431](#)  
[unlock\(\) \(pylablib.devices.AWG.specific.InstekAFG2225\), 425](#)  
[unlock\(\) \(pylablib.devices.AWG.specific.RigolDG1000\), 449](#)  
[unlock\(\) \(pylablib.devices.AWG.specific.RSInstekAFG21000\), 437](#)  
[unlock\(\) \(pylablib.devices.AWG.specific.TektronixAFG1000\), 443](#)  
[unlock\(\) \(pylablib.devices.Conrad.base.RelayBoard\), 535](#)  
[unlock\(\) \(pylablib.devices.Cryocon.base.Cryocon1x\), 538](#)  
[unlock\(\) \(pylablib.devices.Cryomagnetics.base.LM500\), 543](#)  
[unlock\(\) \(pylablib.devices.Cryomagnetics.base.LM510\), 547](#)  
[unlock\(\) \(pylablib.devices.KJL.base.KJL300 method\), 592](#)  
[unlock\(\) \(pylablib.devices.Lakeshore.base.Lakeshore218\), 597](#)  
[unlock\(\) \(pylablib.devices.Lakeshore.base.Lakeshore370 method\), 602](#)  
[unlock\(\) \(pylablib.devices.LaserQuantum.base.Finesse\), 606](#)  
[unlock\(\) \(pylablib.devices.Leybold.base.GenericITR\), 609](#)  
[unlock\(\) \(pylablib.devices.Leybold.base.ITR90\), 611](#)  
[unlock\(\) \(pylablib.devices.LighthousePhotonics.base.SproutG\), 613](#)  
[unlock\(\) \(pylablib.devices.Newport.picomotor.Picomotor8742\), 631](#)  
[unlock\(\) \(pylablib.devices.Ophir.base.OphirDevice\), 640](#)  
[unlock\(\) \(pylablib.devices.Ophir.base.VegaPowerMeter\), 643](#)  
[unlock\(\) \(pylablib.devices.OZOptics.base.DD100\), 636](#)  
[unlock\(\) \(pylablib.devices.OZOptics.base.EPC04\), 638](#)  
[unlock\(\) \(pylablib.devices.OZOptics.base.OZOpticsDevice\), 633](#)  
[unlock\(\) \(pylablib.devices.OZOptics.base.TF100\), 635](#)  
[unlock\(\) \(pylablib.devices.Pfeiffer.base.DPG202\), 659](#)  
[unlock\(\) \(pylablib.devices.Pfeiffer.base.TPG260\), 657](#)  
[unlock\(\) \(pylablib.devices.PhysikInstrumente.base.GenericPIController\), 705](#)  
[unlock\(\) \(pylablib.devices.PhysikInstrumente.base.PIE516\), 708](#)  
[unlock\(\) \(pylablib.devices.Sirah.Matisse.SirahMatisse\), 742](#)  
[unlock\(\) \(pylablib.devices.Tektronix.base.DPO2000\), 768](#)  
[unlock\(\) \(pylablib.devices.Tektronix.base.ITektronixScope\), 754](#)  
[unlock\(\) \(pylablib.devices.Tektronix.base.TDS2000\), 761](#)  
[unlock\(\) \(pylablib.devices.Thorlabs.elliptec.ElliptecMotor\), 781](#)  
[unlock\(\) \(pylablib.devices.Thorlabs.kinesis.BasicKinesisDevice\), 784](#)  
[unlock\(\) \(pylablib.devices.Thorlabs.kinesis.KinesisDevice\), 790](#)  
[unlock\(\) \(pylablib.devices.Thorlabs.kinesis.KinesisMotor\), 800](#)  
[unlock\(\) \(pylablib.devices.Thorlabs.kinesis.KinesisPiezoMotor\), 805](#)  
[unlock\(\) \(pylablib.devices.Thorlabs.kinesis.KinesisQuadDetector\), 809](#)  
[unlock\(\) \(pylablib.devices.Thorlabs.kinesis.MFF\), 809](#)

method), 794

unlock() (pylablib.devices.Thorlabs.serial.FW method), 815

unlock() (pylablib.devices.Thorlabs.serial.FWv1 method), 818

unlock() (pylablib.devices.Thorlabs.serial.MDT69xA method), 822

unlock() (pylablib.devices.Thorlabs.serial.ThorlabsSerialInterface method), 811

unlock() (pylablib.devices.Toptica.ibeam.TopticaIBeam method), 826

unlock() (pylablib.devices.Trinamic.base.TMCM1110 method), 830

unlock\_all() (pylablib.devices.Sirah.tuner.MatisseTuner method), 744

unlock\_etalon() (pylablib.devices.M2.solstis.Solstis method), 616

unlock\_reference\_cavity() (pylablib.devices.M2.solstis.Solstis method), 617

unpack\_int() (in module pylablib.core.utils.strpack), 401

unpack\_numpy\_u12bit() (in module pylablib.core.utils.strpack), 401

unpack\_uint() (in module pylablib.core.utils.strpack), 401

unread() (pylablib.devices.interface.camera.TFramesStatus attribute), 830

unschedule() (pylablib.core.thread.callsync.QQueueLengthLimitScheduler method), 288

unschedule() (pylablib.core.thread.callsync.QQueueScheduler method), 286

unschedule() (pylablib.core.thread.callsync.QQueueSizeLimitScheduler method), 290

unschedule() (pylablib.core.thread.controller.QTaskThread.Job method), 303

unsubscribe() (pylablib.core.thread.controller.QTaskThread method), 313

unsubscribe() (pylablib.core.thread.controller.QThreadController method), 297

unsubscribe() (pylablib.core.thread.multicast\_pool.MulticastPool method), 316

unwrap\_mod\_data() (in module pylablib.core.dataproc.utils), 137

unzip\_file() (in module pylablib.core.utils.files), 368

unzip\_folder() (in module pylablib.core.utils.files), 368

update() (pylablib.core.utils.dictionary.Dictionary method), 331

update() (pylablib.core.utils.dictionary.DictionaryPointer method), 342

update() (pylablib.core.utils.dictionary.FilterTree method), 358

update() (pylablib.core.utils.dictionary.PrefixTree method), 350

update\_acquired\_frames() (pylablib.devices.interface.camera.FrameCounter method), 835

update\_attribute\_value() (pylablib.devices.PhotonFocus.PhotonFocus.IPhotonFocusCamera method), 672

update\_attribute\_value() (pylablib.devices.PhotonFocus.PhotonFocus.PhotonFocusBitFlowCamera method), 701

update\_attribute\_value() (pylablib.devices.PhotonFocus.PhotonFocus.PhotonFocusIMACamera method), 686

update\_attribute\_value() (pylablib.devices.PhotonFocus.PhotonFocus.PhotonFocusSiSoCamera method), 694

update\_available\_axes() (pylablib.devices.Attocube.anc300.ANC300 method), 503

update\_fit\_parameters() (pylablib.core.dataproc.fitting.Fitter method), 115

update\_fixed\_parameters() (pylablib.core.dataproc.fitting.Fitter method), 115

update\_full\_data() (pylablib.devices.PCO.SC2.PCOSC2Camera method), 646

update\_indicators() (pylablib.core.gui.value\_handling.GUIValues method), 280

update\_indicators() (pylablib.core.gui.widgets.container.IQContainer method), 202

update\_indicators() (pylablib.core.gui.widgets.container.IQWidgetContainer method), 208

update\_indicators() (pylablib.core.gui.widgets.container.QContainer method), 204

update\_indicators() (pylablib.core.gui.widgets.container.QDialogContainer method), 220

update\_indicators() (pylablib.core.gui.widgets.container.QFrameContainer method), 216

update\_indicators() (pylablib.core.gui.widgets.container.QGroupBoxContainer method), 224



update_indicators()	(py- lablib.core.gui.widgets.container.QScrollAreaContainer method), 231	update_value() (py- lablib.core.gui.value_handling.GUIValues method), 281
update_indicators()	(py- lablib.core.gui.widgets.container.QScrollAreaContainer.QComboBoxWidget method), 228	update_value() (py- lablib.core.gui.widgets.container.IQContainer method), 202
update_indicators()	(py- lablib.core.gui.widgets.container.QTabContainer method), 233	update_value() (py- lablib.core.gui.widgets.container.IQWidgetContainer method), 208
update_indicators()	(py- lablib.core.gui.widgets.container.QWidgetContainer method), 212	update_value() (py- lablib.core.gui.widgets.container.QContainer method), 204
update_indicators()	(py- lablib.core.gui.widgets.param_table.ParamTable method), 248	update_value() (py- lablib.core.gui.widgets.container.QDialogContainer method), 220
update_indicators()	(py- lablib.core.gui.widgets.param_table.StatusTable method), 260	update_value() (py- lablib.core.gui.widgets.container.QFrameContainer method), 216
update_limits()	(py- lablib.devices.Andor.AndorSDK3.AndorSDK3Attribute method), 478	update_value() (py- lablib.core.gui.widgets.container.QGroupBoxContainer method), 224
update_limits()	(py- lablib.devices.Basler.pylon.BaslerPylonAttribute method), 513	update_value() (py- lablib.core.gui.widgets.container.QScrollAreaContainer method), 231
update_limits()	(py- lablib.devices.DCAM.DCAM.DCAMAttribute method), 549	update_value() (py- lablib.core.gui.widgets.container.QScrollAreaContainer.QContainer method), 228
update_limits()	(py- lablib.devices.IMAQdx.IMAQdx.IMAQdxAttribute method), 577	update_value() (py- lablib.core.gui.widgets.container.QTabContainer method), 233
update_limits()	(py- lablib.devices.Photometrics.pvcam.PvcamAttribute method), 660	update_value() (py- lablib.core.gui.widgets.container.QWidgetContainer method), 212
update_limits()	(py- lablib.devices.PhotonFocus.PhotonFocus.PFCamAttribute method), 671	update_value() (py- lablib.core.gui.widgets.param_table.ParamTable method), 248
update_limits()	(py- lablib.devices.PrincetonInstruments.picam.PicamAttribute method), 711	update_value() (py- lablib.core.gui.widgets.param_table.StatusTable method), 260
update_limits()	(py- lablib.devices.SiliconSoftware.fgrab.FGrabAttribute method), 720	updated() (pylablib.core.utils.dictionary.Dictionary method), 331
update_properties()	(py- lablib.devices.Andor.AndorSDK3.AndorSDK3Attribute method), 478	updated() (pylablib.core.utils.dictionary.DictionaryPointer method), 342
update_reports()	(py- lablib.devices.M2.solstis.Solstis method), 615	updated() (pylablib.core.utils.dictionary.FilterTree method), 359
update_status()	(py- lablib.core.thread.controller.QTaskThread method), 305	updated() (pylablib.core.utils.dictionary.PrefixShortcutTree method), 360
update_status_line()	(py- lablib.core.gui.widgets.param_table.StatusTable method), 252	updated() (pylablib.core.utils.dictionary.PrefixTree method), 351
		upper_limit (pylablib.devices.Sirah.Matisse.TRefcellWaveformParameter attribute), 735
		upper_limit (pylablib.devices.Sirah.Matisse.TScanParameters attribute), 736
		usb_version (pylablib.devices.uc480.uc480.TDeviceInfo attribute), 736

attribute), 865  
 use\_parameters() (in module py-  
 lablib.core.devio.interface), 168  
 use\_xarg() (pylablib.core.dataproc.fitting.Fitter  
 method), 115  
 user\_name (pylablib.devices.Basler.pylon.TCameraInfo  
 attribute), 511  
 user\_name (pylablib.devices.Basler.pylon.TDeviceInfo  
 attribute), 514  
 using\_channel() (py-  
 lablib.devices.Thorlabs.kinesis.KinesisDevice  
 method), 787  
 using\_channel() (py-  
 lablib.devices.Thorlabs.kinesis.KinesisMotor  
 method), 800  
 using\_channel() (py-  
 lablib.devices.Thorlabs.kinesis.KinesisPiezoMotor  
 method), 805  
 using\_channel() (py-  
 lablib.devices.Thorlabs.kinesis.MFF method),  
 794  
 using\_default\_addr() (py-  
 lablib.devices.Thorlabs.elliptec.ElliptecMotor  
 method), 778  
 using\_device() (py-  
 lablib.core.devio.interface.CombinedParameterClass  
 method), 168  
 using\_device() (py-  
 lablib.core.devio.interface.EnumParameterClass  
 method), 167  
 using\_device() (py-  
 lablib.core.devio.interface.FunctionParameterClass  
 method), 167  
 using\_device() (py-  
 lablib.core.devio.interface.ICheckingParameterClass  
 method), 164  
 using\_device() (py-  
 lablib.core.devio.interface.IEnumParameterClass  
 method), 166  
 using\_device() (py-  
 lablib.core.devio.interface.IParameterClass  
 method), 163  
 using\_device() (py-  
 lablib.core.devio.interface.RangeParameterClass  
 method), 164  
 using\_layout() (py-  
 lablib.core.gui.widgets.container.IQWidgetContainer  
 method), 209  
 using\_layout() (py-  
 lablib.core.gui.widgets.container.QDialogContainer  
 method), 220  
 using\_layout() (py-  
 lablib.core.gui.widgets.container.QFrameContainer  
 method), 216  
 using\_layout() (py-  
 lablib.core.gui.widgets.container.QGroupBoxContainer  
 method), 224  
 using\_layout() (py-  
 lablib.core.gui.widgets.container.QScrollAreaContainer.QContain  
 method), 228  
 using\_layout() (py-  
 lablib.core.gui.widgets.container.QWidgetContainer  
 method), 213  
 using\_layout() (py-  
 lablib.core.gui.widgets.layout\_manager.IQLayoutManagedWidget  
 method), 238  
 using\_layout() (py-  
 lablib.core.gui.widgets.layout\_manager.QLayoutManagedWidget  
 method), 241  
 using\_layout() (py-  
 lablib.core.gui.widgets.param\_table.ParamTable  
 method), 251  
 using\_layout() (py-  
 lablib.core.gui.widgets.param\_table.StatusTable  
 method), 260  
 using\_method() (in module py-  
 lablib.core.utils.general), 374  
 using\_new\_sublayout() (py-  
 lablib.core.gui.widgets.container.IQWidgetContainer  
 method), 209  
 using\_new\_sublayout() (py-  
 lablib.core.gui.widgets.container.QDialogContainer  
 method), 220  
 using\_new\_sublayout() (py-  
 lablib.core.gui.widgets.container.QFrameContainer  
 method), 217  
 using\_new\_sublayout() (py-  
 lablib.core.gui.widgets.container.QGroupBoxContainer  
 method), 224  
 using\_new\_sublayout() (py-  
 lablib.core.gui.widgets.container.QScrollAreaContainer.QContain  
 method), 228  
 using\_new\_sublayout() (py-  
 lablib.core.gui.widgets.container.QWidgetContainer  
 method), 213  
 using\_new\_sublayout() (py-  
 lablib.core.gui.widgets.layout\_manager.IQLayoutManagedWidget  
 method), 239  
 using\_new\_sublayout() (py-  
 lablib.core.gui.widgets.layout\_manager.QLayoutManagedWidget  
 method), 241  
 using\_new\_sublayout() (py-  
 lablib.core.gui.widgets.param\_table.ParamTable  
 method), 243  
 using\_new\_sublayout() (py-  
 lablib.core.gui.widgets.param\_table.StatusTable  
 method), 260  
 using\_timeout() (py-

<code>lablib.core.devio.comm_backend.FT232DeviceBackend</code> <code>method</code> ), 152	<code>lablib.devices.Cryomagnetics.base.LM500</code> <code>method</code> ), 543
<code>using_timeout()</code> <code>lablib.core.devio.comm_backend.IDeviceCommBackend</code> <code>method</code> ), 144	<code>(py- using_write_buffer()</code> <code>lablib.devices.Cryomagnetics.base.LM510</code> <code>method</code> ), 547
<code>using_timeout()</code> <code>lablib.core.devio.comm_backend.NetworkDeviceBackend</code> <code>method</code> ), 155	<code>(py- using_write_buffer()</code> <code>lablib.devices.Lakeshore.base.Lakeshore218</code> <code>method</code> ), 597
<code>using_timeout()</code> <code>lablib.core.devio.comm_backend.PyUSBDeviceBackend</code> <code>method</code> ), 157	<code>(py- using_write_buffer()</code> <code>lablib.devices.Lakeshore.base.Lakeshore370</code> <code>method</code> ), 602
<code>using_timeout()</code> <code>lablib.core.devio.comm_backend.RecordedDeviceBackend</code> <code>method</code> ), 159	<code>(py- using_write_buffer()</code> <code>lablib.devices.Sirah.Matisse.SirahMatisse</code> <code>method</code> ), 742
<code>using_timeout()</code> <code>lablib.core.devio.comm_backend.SerialDeviceBackend</code> <code>method</code> ), 149	<code>(py- using_write_buffer()</code> <code>lablib.devices.Tektronix.base.DPO2000</code> <code>method</code> ), 768
<code>using_timeout()</code> <code>lablib.core.devio.comm_backend.VisaDeviceBackend</code> <code>method</code> ), 147	<code>(py- using_write_buffer()</code> <code>lablib.devices.Tektronix.base.ITektronixScope</code> <code>method</code> ), 754
<code>using_timeout()</code> <code>lablib.core.utils.net.ClientSocket</code> 389	<code>(py- using_write_buffer()</code> <code>lablib.devices.Tektronix.base.TDS2000</code> <code>method</code> ), 761
<code>using_write_buffer()</code> <code>lablib.core.devio.SCPI.SCPIDevice</code> 139	<code>(py- using_write_buffer()</code> <code>lablib.devices.Thorlabs.serial.FW</code> 815
<code>using_write_buffer()</code> <code>lablib.devices.AWG.generic.GenericAWG</code> <code>method</code> ), 407	<code>(py- using_write_buffer()</code> <code>lablib.devices.Thorlabs.serial.FWv1</code> 818
<code>using_write_buffer()</code> <code>lablib.devices.AWG.specific.Agilent33220A</code> <code>method</code> ), 419	<code>(py- using_write_buffer()</code> <code>lablib.devices.Thorlabs.serial.MDT69xA</code> <code>method</code> ), 822
<code>using_write_buffer()</code> <code>lablib.devices.AWG.specific.Agilent33500</code> <code>method</code> ), 413	<code>(py- using_write_buffer()</code> <code>lablib.devices.Thorlabs.serial.ThorlabsSerialInterface</code> <code>method</code> ), 811
<code>using_write_buffer()</code> <code>lablib.devices.AWG.specific.InstekAFG2000</code> <code>method</code> ), 431	<b>V</b>
<code>using_write_buffer()</code> <code>lablib.devices.AWG.specific.InstekAFG2225</code> <code>method</code> ), 425	<code>value (pylablib.core.devio.interface.TRawParameterValue</code> <code>attribute)</code> , 168
<code>using_write_buffer()</code> <code>lablib.devices.AWG.specific.RigolDG1000</code> <code>method</code> ), 449	<code>value (pylablib.core.thread.multicast_pool.TMulticast</code> <code>attribute)</code> , 315
<code>using_write_buffer()</code> <code>lablib.devices.AWG.specific.RSInstekAFG21000</code> <code>method</code> ), 437	<code>value (pylablib.devices.Leybold.base.TUpdateValue</code> <code>attribute)</code> , 607
<code>using_write_buffer()</code> <code>lablib.devices.AWG.specific.TektronixAFG1000</code> <code>method</code> ), 443	<code>value (pylablib.devices.Trinamic.base.TMCM11110.ReplyData</code> <code>attribute)</code> , 827
<code>using_write_buffer()</code> <code>lablib.devices.Cryocon.base.Cryocon1x</code> <code>method</code> ), 538	<code>value_access (pylablib.devices.PrincetonInstruments.picam.PicamAttr</code> <code>attribute)</code> , 710
<code>using_write_buffer()</code>	<code>value_changed</code> <code>lablib.core.gui.widgets.button.ToggleButton</code> <code>attribute)</code> , 198
	<code>value_changed</code> <code>lablib.core.gui.widgets.combo_box.ComboBox</code> <code>attribute)</code> , 199

value_changed (pylablib.core.gui.widgets.edit.NumEdit attribute), 235	VegaPowerMeter (class in py-lablib.devices.Ophir.base), 641
value_changed (pylablib.core.gui.widgets.edit.TextEdit attribute), 234	VegaPowerMeter.NoParameterCaller (class in py-lablib.devices.Ophir.base), 642
value_changed (pylablib.core.gui.widgets.label.EnumLabel attribute), 237	velocity (pylablib.devices.Thorlabs.kinesis.THomeParams attribute), 785
value_changed (pylablib.core.gui.widgets.label.NumLabel attribute), 238	velocity (pylablib.devices.Thorlabs.kinesis.TPZMotorDriveParams attribute), 786
value_changed (pylablib.core.gui.widgets.label.TextLabel attribute), 236	velocity (pylablib.devices.Thorlabs.kinesis.TPZMotorJogParams attribute), 787
value_entered (pylablib.core.gui.widgets.edit.NumEdit attribute), 235	vendor (pylablib.devices.Basler.pylon.TCameraInfo attribute), 511
value_entered (pylablib.core.gui.widgets.edit.TextEdit attribute), 234	vendor (pylablib.devices.Basler.pylon.TDeviceInfo attribute), 514
value_handler (pylablib.core.gui.widgets.param_table.ParamTable attribute), 243	vendor (pylablib.devices.DCAM.DCAM.TDeviceInfo attribute), 550
value_handler (pylablib.core.gui.widgets.param_table.StatusTable attribute), 252	vendor (pylablib.devices.IMAQdx.IMAQdx.TCameraInfo attribute), 576
value_to_index() (pylablib.core.gui.widgets.combo_box.ComboBox method), 198	vendor (pylablib.devices.IMAQdx.IMAQdx.TDeviceInfo attribute), 578
values (pylablib.devices.Andor.AndorSDK3.AndorSDK3Attribute attribute), 477	vendor (pylablib.devices.Photometrics.pvcam.TDeviceInfo attribute), 661
values (pylablib.devices.Basler.pylon.BaslerPylonAttribute attribute), 513	version (pylablib.devices.AlliedVision.Bonito.TDeviceInfo attribute), 450
values (pylablib.devices.DCAM.DCAM.DCAMAttribute attribute), 549	version (pylablib.devices.Attocube.anc300.TDeviceInfo attribute), 502
values (pylablib.devices.IMAQdx.IMAQdx.IMAQdxAttribute attribute), 577	version (pylablib.devices.IMAQdx.IMAQdx.TCameraInfo attribute), 576
values (pylablib.devices.Photometrics.pvcam.PvcamAttribute attribute), 660	version (pylablib.devices.LighthousePhotonics.base.TDeviceInfo attribute), 611
values (pylablib.devices.PhotonFocus.PhotonFocus.PFCamAttribute attribute), 671	version (pylablib.devices.PhotonFocus.PhotonFocus.TCameraInfo attribute), 670
values (pylablib.devices.PrincetonInstruments.picam.PicamAttribute attribute), 711	version (pylablib.devices.SiliconSoftware.fgrab.TFullAppletInfo attribute), 719
values (pylablib.devices.SiliconSoftware.fgrab.FGGrabAttribute attribute), 720	version (pylablib.devices.Toptica.ibeam.TDeviceInfo attribute), 823
values() (pylablib.core.utils.dictionary.Dictionary method), 329	viewitems() (pylablib.core.utils.dictionary.Dictionary method), 328
values() (pylablib.core.utils.dictionary.DictionaryPointer method), 342	viewitems() (pylablib.core.utils.dictionary.DictionaryPointer method), 342
values() (pylablib.core.utils.dictionary.FilterTree method), 359	viewitems() (pylablib.core.utils.dictionary.FilterTree method), 359
values() (pylablib.core.utils.dictionary.PrefixTree method), 351	viewitems() (pylablib.core.utils.dictionary.PrefixTree method), 351
values() (pylablib.core.utils.dictionary.PrefixTree method), 351	viewkeys() (pylablib.core.utils.dictionary.Dictionary method), 329
	viewkeys() (pylablib.core.utils.dictionary.DictionaryPointer method), 343
	viewkeys() (pylablib.core.utils.dictionary.FilterTree method), 359
	viewkeys() (pylablib.core.utils.dictionary.PrefixTree method), 351
	viewvalues() (pylablib.core.utils.dictionary.Dictionary method), 329

[viewvalues\(\) \(pylablib.core.utils.dictionary.DictionaryPrinter method\), 343](#)  
[viewvalues\(\) \(pylablib.core.utils.dictionary.FilterTree method\), 360](#)  
[viewvalues\(\) \(pylablib.core.utils.dictionary.PrefixTree method\), 352](#)  
[virtual\\_gui\\_values\(\) \(in module pylablib.core.gui.value\\_handling\), 281](#)  
[VirtualIndicatorHandler \(in module pylablib.core.gui.value\\_handling\), 276](#)  
[VirtualValueHandler \(class in pylablib.core.gui.value\\_handling\), 266](#)  
[VisaDeviceBackend \(class in pylablib.core.devio.comm\\_backend\), 144](#)  
[visibility \(pylablib.devices.Basler.pylon.BaslerPylonAttribute attribute\), 512](#)  
[visibility \(pylablib.devices.IMAQdx.IMAQdx.IMAQdxAttribute attribute\), 577](#)

**W**

[wait\(\) \(pylablib.core.devio.SCPIDevice method\), 139](#)  
[wait\(\) \(pylablib.core.thread.callsync.QCallResultSynchronizer method\), 282](#)  
[wait\(\) \(pylablib.core.thread.callsync.QDirectResultSynchronizer method\), 282](#)  
[wait\(\) \(pylablib.core.thread.notifier.ISkippableNotifier method\), 317](#)  
[wait\(\) \(pylablib.core.thread.synchronizing.QMultiThreadNotifier method\), 318](#)  
[wait\(\) \(pylablib.core.thread.synchronizing.QThreadNotifier method\), 318](#)  
[wait\(\) \(pylablib.devices.AWG.generic.GenericAWG method\), 408](#)  
[wait\(\) \(pylablib.devices.AWG.specific.Agilent33220A method\), 419](#)  
[wait\(\) \(pylablib.devices.AWG.specific.Agilent33500 method\), 414](#)  
[wait\(\) \(pylablib.devices.AWG.specific.InstekAFG2000 method\), 431](#)  
[wait\(\) \(pylablib.devices.AWG.specific.InstekAFG2225 method\), 425](#)  
[wait\(\) \(pylablib.devices.AWG.specific.RigolDG1000 method\), 449](#)  
[wait\(\) \(pylablib.devices.AWG.specific.RSInstekAFG21000 method\), 437](#)  
[wait\(\) \(pylablib.devices.AWG.specific.TektronixAFG1000 method\), 443](#)  
[wait\(\) \(pylablib.devices.Cryocon.base.Cryocon1x method\), 538](#)  
[wait\(\) \(pylablib.devices.Cryomagnetics.base.LM500 method\), 543](#)  
[wait\(\) \(pylablib.devices.Cryomagnetics.base.LM510 method\), 547](#)

[wait\(\) \(pylablib.devices.interface.camera.FrameNotifier method\), 836](#)  
[wait\(\) \(pylablib.devices.Lakeshore.base.Lakeshore218 method\), 597](#)  
[wait\(\) \(pylablib.devices.Lakeshore.base.Lakeshore370 method\), 602](#)  
[wait\(\) \(pylablib.devices.Sirah.Matisse.SirahMatisse method\), 742](#)  
[wait\(\) \(pylablib.devices.Tektronix.base.DPO2000 method\), 768](#)  
[wait\(\) \(pylablib.devices.Tektronix.base.ITektronixScope method\), 754](#)  
[wait\(\) \(pylablib.devices.Tektronix.base.TDS2000 method\), 761](#)  
[wait\(\) \(pylablib.devices.Thorlabs.serial.FW method\), 815](#)  
[wait\(\) \(pylablib.devices.Thorlabs.serial.FWv1 method\), 819](#)  
[wait\(\) \(pylablib.devices.Thorlabs.serial.MDT69xA method\), 822](#)  
[wait\(\) \(pylablib.devices.Thorlabs.serial.ThorlabsSerialInterface method\), 811](#)  
[wait\\_dev\(\) \(pylablib.core.devio.SCPIDevice method\), 139](#)  
[wait\\_dev\(\) \(pylablib.devices.AWG.generic.GenericAWG method\), 408](#)  
[wait\\_dev\(\) \(pylablib.devices.AWG.specific.Agilent33220A method\), 420](#)  
[wait\\_dev\(\) \(pylablib.devices.AWG.specific.Agilent33500 method\), 414](#)  
[wait\\_dev\(\) \(pylablib.devices.AWG.specific.InstekAFG2000 method\), 431](#)  
[wait\\_dev\(\) \(pylablib.devices.AWG.specific.InstekAFG2225 method\), 425](#)  
[wait\\_dev\(\) \(pylablib.devices.AWG.specific.RigolDG1000 method\), 449](#)  
[wait\\_dev\(\) \(pylablib.devices.AWG.specific.RSInstekAFG21000 method\), 437](#)  
[wait\\_dev\(\) \(pylablib.devices.AWG.specific.TektronixAFG1000 method\), 443](#)  
[wait\\_dev\(\) \(pylablib.devices.Cryocon.base.Cryocon1x method\), 538](#)  
[wait\\_dev\(\) \(pylablib.devices.Cryomagnetics.base.LM500 method\), 543](#)  
[wait\\_dev\(\) \(pylablib.devices.Cryomagnetics.base.LM510 method\), 547](#)  
[wait\\_dev\(\) \(pylablib.devices.Lakeshore.base.Lakeshore218 method\), 597](#)  
[wait\\_dev\(\) \(pylablib.devices.Lakeshore.base.Lakeshore370 method\), 602](#)  
[wait\\_dev\(\) \(pylablib.devices.Sirah.Matisse.SirahMatisse method\), 742](#)  
[wait\\_dev\(\) \(pylablib.devices.Tektronix.base.DPO2000 method\), 768](#)



<code>wait_dev()</code> ( <code>pylablib.devices.Tektronix.base.ITektronixScope</code> method), 754	<code>lablib.devices.IMAQ.IMAQ.IMAQFrameGrabber</code> method), 568
<code>wait_dev()</code> ( <code>pylablib.devices.Tektronix.base.TDS2000</code> method), 761	<code>wait_for_frame()</code> ( <code>py-</code> <code>lablib.devices.IMAQdx.IMAQdx.EthernetIMAQdxCamera</code> method), 589
<code>wait_dev()</code> ( <code>pylablib.devices.Thorlabs.serial.FW</code> method), 815	<code>wait_for_frame()</code> ( <code>py-</code> <code>lablib.devices.IMAQdx.IMAQdx.IMAQdxCamera</code> method), 583
<code>wait_dev()</code> ( <code>pylablib.devices.Thorlabs.serial.FWv1</code> method), 819	<code>wait_for_frame()</code> ( <code>py-</code> <code>lablib.devices.interface.camera.IAttributeCamera</code> method), 841
<code>wait_dev()</code> ( <code>pylablib.devices.Thorlabs.serial.MDT69xA</code> method), 822	<code>wait_for_frame()</code> ( <code>py-</code> <code>lablib.devices.interface.camera.IBinROICamera</code> method), 860
<code>wait_dev()</code> ( <code>pylablib.devices.Thorlabs.serial.ThorlabsSerialInterface</code> method), 812	<code>wait_for_frame()</code> ( <code>py-</code> <code>lablib.devices.interface.camera.ICamera</code> method), 832
<code>wait_done()</code> ( <code>pylablib.devices.interface.camera.FrameCounter</code> method), 836	<code>wait_for_frame()</code> ( <code>py-</code> <code>lablib.devices.interface.camera.IExposureCamera</code> method), 851
<code>wait_for_any_message()</code> ( <code>py-</code> <code>lablib.core.thread.controller.QTaskThread</code> method), 314	<code>wait_for_frame()</code> ( <code>py-</code> <code>lablib.devices.interface.camera.IGrabberAttributeCamera</code> method), 846
<code>wait_for_any_message()</code> ( <code>py-</code> <code>lablib.core.thread.controller.QThreadController</code> method), 294	<code>wait_for_frame()</code> ( <code>py-</code> <code>lablib.devices.AlliedVision.Bonito.BonitoIMAQCamera</code> method), 464
<code>wait_for_fine_tuning()</code> ( <code>py-</code> <code>lablib.devices.M2.solstis.Solstis</code> method), 615	<code>wait_for_frame()</code> ( <code>py-</code> <code>lablib.devices.AlliedVision.Bonito.IBonitoCamera</code> method), 456
<code>wait_for_frame()</code> ( <code>py-</code> <code>lablib.devices.Andor.AndorSDK2.AndorSDK2Camera</code> method), 476	<code>wait_for_frame()</code> ( <code>py-</code> <code>lablib.devices.Andor.AndorSDK2.AndorSDK2Camera</code> method), 668
<code>wait_for_frame()</code> ( <code>py-</code> <code>lablib.devices.Andor.AndorSDK3.AndorSDK3Camera</code> method), 486	<code>wait_for_frame()</code> ( <code>py-</code> <code>lablib.devices.Andor.AndorSDK3.AndorSDK3Camera</code> method), 677
<code>wait_for_frame()</code> ( <code>py-</code> <code>lablib.devices.Andor.AndorSDK3.AndorSDK3Camera</code> method), 482	<code>wait_for_frame()</code> ( <code>py-</code> <code>lablib.devices.Andor.AndorSDK3.AndorSDK3Camera</code> method), 701
<code>wait_for_frame()</code> ( <code>py-</code> <code>lablib.devices.Basler.pylon.BaslerPylonCamera</code> method), 520	<code>wait_for_frame()</code> ( <code>py-</code> <code>lablib.devices.Andor.AndorSDK3.AndorSDK3Camera</code> method), 701
<code>wait_for_frame()</code> ( <code>py-</code> <code>lablib.devices.BitFlow.BitFlow.BitFlowCamera</code> method), 532	<code>wait_for_frame()</code> ( <code>py-</code> <code>lablib.devices.Basler.pylon.BaslerPylonCamera</code> method), 686
<code>wait_for_frame()</code> ( <code>py-</code> <code>lablib.devices.BitFlow.BitFlow.BitFlowFrameGrabber</code> method), 526	<code>wait_for_frame()</code> ( <code>py-</code> <code>lablib.devices.BitFlow.BitFlow.BitFlowCamera</code> method), 694
<code>wait_for_frame()</code> ( <code>py-</code> <code>lablib.devices.DCAM.DCAM.DCAMCamera</code> method), 556	<code>wait_for_frame()</code> ( <code>py-</code> <code>lablib.devices.BitFlow.BitFlow.BitFlowFrameGrabber</code> method), 717
<code>wait_for_frame()</code> ( <code>py-</code> <code>lablib.devices.IMAQ.IMAQ.IMAQCamera</code> method), 575	<code>wait_for_frame()</code> ( <code>py-</code> <code>lablib.devices.DCAM.DCAM.DCAMCamera</code> method), 733
<code>wait_for_frame()</code> ( <code>py-</code> <code>lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareCamera</code> method), 727	<code>wait_for_frame()</code> ( <code>py-</code> <code>lablib.devices.DCAM.DCAM.DCAMCamera</code> method), 733
	<code>wait_for_frame()</code> ( <code>py-</code> <code>lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareFrameGrabber</code> method), 727
	<code>wait_for_frame()</code> ( <code>py-</code> <code>lablib.devices.SiliconSoftware.fgrab.SiliconSoftwareFrameGrabber</code> method), 727

`lablib.devices.Thorlabs.TLCamera.ThorlabsTLCamera.wait_for_stop()` (py-  
 method), 776  
`wait_for_frame()` (py-  
`lablib.devices.Thorlabs.TLCamera.ThorlabsTLCamera.wait_for_frame()` (py-  
 method), 771  
`wait_for_frame()` (py-  
`lablib.devices.uc480.uc480.UC480Camera.wait_for_frame()` (py-  
 method), 872  
`wait_for_grabbing()` (py-  
`lablib.devices.Tektronix.base.DPO2000.wait_for_grabbing()` (py-  
 method), 768  
`wait_for_grabbing()` (py-  
`lablib.devices.Tektronix.base.ITektronixScope.wait_for_grabbing()` (py-  
 method), 749  
`wait_for_grabbing()` (py-  
`lablib.devices.Tektronix.base.TDS2000.wait_for_grabbing()` (py-  
 method), 761  
`wait_for_home()` (py-  
`lablib.devices.Thorlabs.kinesis.KinesisMotor.wait_for_home()` (py-  
 method), 796  
`wait_for_keypress()` (in module py-  
`lablib.core.utils.general`), 381  
`wait_for_measurement()` (py-  
`lablib.devices.Cryomagnetics.base.LM500.wait_for_measurement()` (py-  
 method), 540  
`wait_for_measurement()` (py-  
`lablib.devices.Cryomagnetics.base.LM510.wait_for_measurement()` (py-  
 method), 547  
`wait_for_message()` (py-  
`lablib.core.thread.controller.QTaskThread.wait_for_message()` (py-  
 method), 314  
`wait_for_message()` (py-  
`lablib.core.thread.controller.QThreadController.wait_for_message()` (py-  
 method), 294  
`wait_for_report()` (py-  
`lablib.devices.M2.solstis.Solstis.wait_for_report()` (py-  
 method), 615  
`wait_for_sample()` (py-  
`lablib.devices.NI.daq.NIDAQ.wait_for_sample()` (py-  
 method), 624  
`wait_for_scan()` (py-  
`lablib.devices.Newport.picomotor.Picomotor8742.wait_for_scan()` (py-  
 method), 629  
`wait_for_status()` (py-  
`lablib.devices.SmarAct.scu3d.SC3D.wait_for_status()` (py-  
 method), 746  
`wait_for_status()` (py-  
`lablib.devices.Thorlabs.kinesis.KinesisMotor.wait_for_status()` (py-  
 method), 795  
`wait_for_status()` (py-  
`lablib.devices.Thorlabs.kinesis.KinesisPiezoMotor.wait_for_status()` (py-  
 method), 801  
`wait_for_status()` (py-  
`lablib.devices.Thorlabs.kinesis.MFF.wait_for_status()` (py-  
 method), 791  
`wait_for_stop()` (py-  
`lablib.devices.Thorlabs.kinesis.KinesisMotor.wait_for_stop()` (py-  
 method), 796  
`wait_for_sync()` (py-  
`lablib.core.thread.controller.QTaskThread.wait_for_sync()` (py-  
 method), 314  
`wait_for_sync()` (py-  
`lablib.core.thread.controller.QThreadController.wait_for_sync()` (py-  
 method), 294  
`wait_for_terascan_update()` (py-  
`lablib.devices.M2.solstis.Solstis.wait_for_terascan_update()` (py-  
 method), 618  
`wait_move()` (pylablib.devices.Arcus.performax.Performax2EXStage  
 method), 497  
`wait_move()` (pylablib.devices.Arcus.performax.Performax4EXStage  
 method), 491  
`wait_move()` (pylablib.devices.Arcus.performax.PerformaxDMXJSStage  
 method), 498  
`wait_move()` (pylablib.devices.Attocube.anc300.ANC300  
 method), 505  
`wait_move()` (pylablib.devices.Attocube.anc350.ANC350  
 method), 508  
`wait_move()` (pylablib.devices.Newport.picomotor.Picomotor8742  
 method), 630  
`wait_move()` (pylablib.devices.SmarAct.scu3d.SC3D  
 method), 747  
`wait_move()` (pylablib.devices.Thorlabs.kinesis.KinesisMotor  
 method), 796  
`wait_move()` (pylablib.devices.Thorlabs.kinesis.KinesisPiezoMotor  
 method), 802  
`wait_move()` (pylablib.devices.Trinamic.base.TMCM1110  
 method), 829  
`wait_scan()` (pylablib.devices.Sirah.Matisse.SirahMatisse  
 method), 740  
`wait_start()` (pylablib.devices.interface.camera.FrameCounter  
 method), 835  
`wait_sync()` (pylablib.core.devio.SCPI.SCPIDevice  
 method), 139  
`wait_sync()` (pylablib.devices.AWG.generic.GenericAWG  
 method), 408  
`wait_sync()` (pylablib.devices.AWG.specific.Agilent33220A  
 method), 420  
`wait_sync()` (pylablib.devices.AWG.specific.Agilent33500  
 method), 414  
`wait_sync()` (pylablib.devices.AWG.specific.InstekAFG2000  
 method), 432  
`wait_sync()` (pylablib.devices.AWG.specific.InstekAFG2225  
 method), 426  
`wait_sync()` (pylablib.devices.AWG.specific.RigolDG1000  
 method), 449  
`wait_sync()` (pylablib.devices.AWG.specific.RSInstekAFG21000  
 method), 437  
`wait_sync()` (pylablib.devices.AWG.specific.TektronixAFG1000  
 method), 443

`wait_sync()` (`pylablib.devices.Cryocon.base.Cryocon1x` attribute), 644  
`method`), 538  
`wait_sync()` (`pylablib.devices.Cryomagnetics.base.LM500` `wheelEvent()` (`pylablib.core.gui.widgets.combo_box.ComboBox`  
`method`), 543 `method`), 198  
`wait_sync()` (`pylablib.devices.Cryomagnetics.base.LM510` `widget` (`pylablib.core.gui.widgets.container.TChild` at-  
`method`), 547 `tribute`), 200  
`wait_sync()` (`pylablib.devices.Lakeshore.base.Lakeshore218` `widget` (`pylablib.core.gui.widgets.param_table.ParamTable.ParamRow`  
`method`), 597 `attribute`), 243  
`wait_sync()` (`pylablib.devices.Lakeshore.base.Lakeshore370` `widget` (`pylablib.core.gui.widgets.param_table.StatusTable.ParamRow`  
`method`), 602 `attribute`), 252  
`wait_sync()` (`pylablib.devices.Sirah.Matisse.SirahMatisse` `width` (`pylablib.core.dataproc.feature.Baseline` at-  
`method`), 742 `tribute`), 108  
`wait_sync()` (`pylablib.devices.Tektronix.base.DPO2000` `width` (`pylablib.core.dataproc.feature.Peak` attribute),  
`method`), 768 109  
`wait_sync()` (`pylablib.devices.Tektronix.base.ITektronixScope` `width` (`pylablib.devices.interface.camera.TFrameSize`  
`method`), 755 `attribute`), 830  
`wait_sync()` (`pylablib.devices.Tektronix.base.TDS2000` `width` (`pylablib.devices.Thorlabs.kinesis.TKCubeTrigPosParams`  
`method`), 761 `attribute`), 786  
`wait_sync()` (`pylablib.devices.Thorlabs.serial.FW` `window` (`pylablib.devices.Lakeshore.base.TLakeshore218FilterSettings`  
`method`), 815 `attribute`), 593  
`wait_sync()` (`pylablib.devices.Thorlabs.serial.FWv1` `window` (`pylablib.devices.Lakeshore.base.TLakeshore370FilterSettings`  
`method`), 819 `attribute`), 599  
`wait_sync()` (`pylablib.devices.Thorlabs.serial.MDT69xA` `with_traceback()` (py-  
`method`), 822 `lablib.core.devio.base.DeviceError` `method`),  
142  
`wait_sync()` (`pylablib.devices.Thorlabs.serial.ThorlabsSerialInterface` `with_traceback()` (py-  
`method`), 812 `lablib.core.devio.comm_backend.DeviceBackendError`  
`method`), 142  
`wait_until()` (`pylablib.core.thread.controller.QTaskThread` `with_traceback()` (py-  
`method`), 314  
`wait_until()` (`pylablib.core.thread.controller.QThreadController` `lablib.core.devio.comm_backend.DeviceFT232Error`  
`method`), 294 `method`), 150  
`wait_until()` (`pylablib.core.thread.synchronizing.QMultiThreadNotifier` `with_traceback()` (py-  
`method`), 319 `lablib.core.devio.comm_backend.DeviceNetworkError`  
`method`), 152  
`waiting()` (`pylablib.core.thread.callsync.QCallResultSynchronizer` `method`), 152  
`method`), 282 `with_traceback()` (py-  
`waiting()` (`pylablib.core.thread.callsync.QDirectResultSynchronizer` `lablib.core.devio.comm_backend.DeviceRecordedError`  
`method`), 283 `method`), 157  
`waiting()` (`pylablib.core.thread.notifier.ISkipableNotifier` `with_traceback()` (py-  
`method`), 317 `lablib.core.devio.comm_backend.DeviceSerialError`  
`method`), 147  
`waiting()` (`pylablib.core.thread.synchronizing.QThreadNotifier` `with_traceback()` (py-  
`method`), 318  
`waiting_state()` (py- `lablib.core.devio.comm_backend.DeviceUSBError`  
`lablib.core.thread.callsync.QCallResultSynchronizer` `method`), 155  
`method`), 282 `with_traceback()` (py-  
`waiting_state()` (py- `lablib.core.devio.comm_backend.DeviceVisaError`  
`lablib.core.thread.callsync.QDirectResultSynchronizer` `method`), 144  
`method`), 283 `with_traceback()` (py-  
`waiting_state()` (py- `lablib.core.gui.limiter.LimitError` `method`),  
`lablib.core.thread.notifier.ISkipableNotifier` 262  
`method`), 317 `with_traceback()` (py-  
`waiting_state()` (py- `lablib.core.gui.value_handling.MissingGUIHandlerError`  
`lablib.core.thread.synchronizing.QThreadNotifier` `method`), 278  
`method`), 318 `with_traceback()` (py-  
`walk_dir()` (in module `pylablib.core.utils.files`), 364  
`warnings` (`pylablib.devices.PCO.SC2.TCameraStatus` `method`), 265  
`method`), 265



<code>with_traceback()</code> <code>lablib.core.thread.threadprop.DuplicateControllerThreadError</code> <code>method), 319</code>	(py- <code>with_traceback()</code> <code>lablib.devices.Arduino.base.ArduinoError</code> <code>method), 500</code>	(py- <code>with_traceback()</code> <code>lablib.devices.Arduino.base.ArduinoError</code> <code>method), 500</code>
<code>with_traceback()</code> <code>lablib.core.thread.threadprop.InterruptException</code> <code>method), 320</code>	(py- <code>with_traceback()</code> <code>lablib.devices.Attocube.base.AttocubeBackendError</code> <code>method), 510</code>	(py- <code>with_traceback()</code> <code>lablib.devices.Attocube.base.AttocubeBackendError</code> <code>method), 510</code>
<code>with_traceback()</code> <code>lablib.core.thread.threadprop.InterruptExceptionStop</code> <code>method), 320</code>	(py- <code>with_traceback()</code> <code>lablib.devices.Attocube.base.AttocubeError</code> <code>method), 510</code>	(py- <code>with_traceback()</code> <code>lablib.devices.Attocube.base.AttocubeError</code> <code>method), 510</code>
<code>with_traceback()</code> <code>lablib.core.thread.threadprop.NoControllerThreadError</code> <code>method), 319</code>	(py- <code>with_traceback()</code> <code>lablib.devices.AWG.generic.GenericAWGBackendError</code> <code>method), 402</code>	(py- <code>with_traceback()</code> <code>lablib.devices.AWG.generic.GenericAWGBackendError</code> <code>method), 402</code>
<code>with_traceback()</code> <code>lablib.core.thread.threadprop.NoMessageThreadError</code> <code>method), 320</code>	(py- <code>with_traceback()</code> <code>lablib.devices.AWG.generic.GenericAWGError</code> <code>method), 402</code>	(py- <code>with_traceback()</code> <code>lablib.devices.AWG.generic.GenericAWGError</code> <code>method), 402</code>
<code>with_traceback()</code> <code>lablib.core.thread.threadprop.SkippedCallError</code> <code>method), 320</code>	(py- <code>with_traceback()</code> <code>lablib.devices.BitFlow.BitFlow.BitFlowError</code> <code>method), 520</code>	(py- <code>with_traceback()</code> <code>lablib.devices.BitFlow.BitFlow.BitFlowError</code> <code>method), 520</code>
<code>with_traceback()</code> <code>lablib.core.thread.threadprop.ThreadError</code> <code>method), 319</code>	(py- <code>with_traceback()</code> <code>lablib.devices.BitFlow.BitFlow.BitFlowTimeoutError</code> <code>method), 521</code>	(py- <code>with_traceback()</code> <code>lablib.devices.BitFlow.BitFlow.BitFlowTimeoutError</code> <code>method), 521</code>
<code>with_traceback()</code> <code>lablib.core.thread.threadprop.TimeoutThreadError</code> <code>method), 320</code>	(py- <code>with_traceback()</code> <code>lablib.devices.Conrad.base.ConradBackendError</code> <code>method), 533</code>	(py- <code>with_traceback()</code> <code>lablib.devices.Conrad.base.ConradBackendError</code> <code>method), 533</code>
<code>with_traceback()</code> <code>lablib.core.utils.net.SocketError</code> <code>method), 388</code>	(py- <code>with_traceback()</code> <code>lablib.devices.Conrad.base.ConradError</code> <code>method), 533</code>	(py- <code>with_traceback()</code> <code>lablib.devices.Conrad.base.ConradError</code> <code>method), 533</code>
<code>with_traceback()</code> <code>lablib.core.utils.net.SocketTimeout</code> <code>method), 388</code>	(py- <code>with_traceback()</code> <code>lablib.devices.Cryocon.base.CryoconBackendError</code> <code>method), 535</code>	(py- <code>with_traceback()</code> <code>lablib.devices.Cryocon.base.CryoconBackendError</code> <code>method), 535</code>
<code>with_traceback()</code> <code>lablib.devices.AlliedVision.Bonito.BonitoError</code> <code>method), 450</code>	(py- <code>with_traceback()</code> <code>lablib.devices.Cryocon.base.CryoconError</code> <code>method), 535</code>	(py- <code>with_traceback()</code> <code>lablib.devices.Cryocon.base.CryoconError</code> <code>method), 535</code>
<code>with_traceback()</code> <code>lablib.devices.Andor.base.AndorError</code> <code>method), 486</code>	(py- <code>with_traceback()</code> <code>lablib.devices.Cryomagnetics.base.CryomagneticsBackendError</code> <code>method), 540</code>	(py- <code>with_traceback()</code> <code>lablib.devices.Cryomagnetics.base.CryomagneticsBackendError</code> <code>method), 540</code>
<code>with_traceback()</code> <code>lablib.devices.Andor.base.AndorFrameTransferError</code> <code>method), 487</code>	(py- <code>with_traceback()</code> <code>lablib.devices.Cryomagnetics.base.CryomagneticsError</code> <code>method), 539</code>	(py- <code>with_traceback()</code> <code>lablib.devices.Cryomagnetics.base.CryomagneticsError</code> <code>method), 539</code>
<code>with_traceback()</code> <code>lablib.devices.Andor.base.AndorNotSupportedError</code> <code>method), 487</code>	(py- <code>with_traceback()</code> <code>lablib.devices.interface.camera.DefaultFrameTransferError</code> <code>method), 830</code>	(py- <code>with_traceback()</code> <code>lablib.devices.interface.camera.DefaultFrameTransferError</code> <code>method), 830</code>
<code>with_traceback()</code> <code>lablib.devices.Andor.base.AndorTimeoutError</code> <code>method), 487</code>	(py- <code>with_traceback()</code> <code>lablib.devices.KJL.base.KJLBackendError</code> <code>method), 590</code>	(py- <code>with_traceback()</code> <code>lablib.devices.KJL.base.KJLBackendError</code> <code>method), 590</code>
<code>with_traceback()</code> <code>lablib.devices.Arcus.base.ArcusBackendError</code> <code>method), 487</code>	(py- <code>with_traceback()</code> <code>lablib.devices.KJL.base.KJLError</code> <code>method), 589</code>	(py- <code>with_traceback()</code> <code>lablib.devices.KJL.base.KJLError</code> <code>method), 589</code>
<code>with_traceback()</code> <code>lablib.devices.Arcus.base.ArcusError</code> <code>method), 487</code>	(py- <code>with_traceback()</code> <code>lablib.devices.Lakeshore.base.LakeshoreBackendError</code> <code>method), 592</code>	(py- <code>with_traceback()</code> <code>lablib.devices.Lakeshore.base.LakeshoreBackendError</code> <code>method), 592</code>
<code>with_traceback()</code> <code>lablib.devices.Arduino.base.ArduinoBackendError</code> <code>method), 500</code>	(py- <code>with_traceback()</code> <code>lablib.devices.Lakeshore.base.LakeshoreError</code> <code>method), 592</code>	(py- <code>with_traceback()</code> <code>lablib.devices.Lakeshore.base.LakeshoreError</code> <code>method), 592</code>



method), 370  
 writable (pylablib.devices.Andor.AndorSDK3.AndorSDK3Attribute  
 attribute), 477  
 writable (pylablib.devices.Basler.pylon.BaslerPylonAttribute  
 attribute), 512  
 writable (pylablib.devices.DCAM.DCAM.DCAMAttribute  
 attribute), 549  
 writable (pylablib.devices.IMAQdx.IMAQdx.IMAQdxAttribute  
 attribute), 577  
 writable (pylablib.devices.Photometrics.pvcam.PvcamAttribute  
 attribute), 660  
 writable (pylablib.devices.PhotonFocus.PhotonFocus.PhotonFocusAttribute  
 attribute), 670  
 writable (pylablib.devices.PrincetonInstruments.picam.PicamAttribute  
 attribute), 710  
 write() (pylablib.core.devio.comm\_backend.FT232DeviceBackend  
 method), 151  
 write() (pylablib.core.devio.comm\_backend.IDeviceCommunicationBackend  
 method), 144  
 write() (pylablib.core.devio.comm\_backend.NetworkDeviceBackend  
 method), 154  
 write() (pylablib.core.devio.comm\_backend.PyUSBDeviceBackend  
 method), 156  
 write() (pylablib.core.devio.comm\_backend.RecordedDeviceBackend  
 method), 158  
 write() (pylablib.core.devio.comm\_backend.SerialDeviceBackend  
 method), 148  
 write() (pylablib.core.devio.comm\_backend.VisaDeviceBackend  
 method), 146  
 write() (pylablib.core.devio.SCPI.SCPIDevice write()  
 method), 139  
 write() (pylablib.core.fileio.savefile.CSVTableOutputFileFormat  
 method), 192  
 write() (pylablib.core.fileio.savefile.DictionaryOutputFileFormat  
 method), 193  
 write() (pylablib.core.fileio.savefile.IBinaryOutputFileFormat  
 method), 193  
 write() (pylablib.core.fileio.savefile.IOutputFileFormat  
 method), 190  
 write() (pylablib.core.fileio.savefile.ITextOutputFileFormat  
 method), 191  
 write() (pylablib.core.fileio.savefile.TableBinaryOutputFileFormat  
 method), 194  
 write() (pylablib.core.utils.general.StreamFileLogger  
 method), 380  
 write() (pylablib.devices.AWG.generic.GenericAWG  
 method), 408  
 write() (pylablib.devices.AWG.specific.Agilent33220A  
 method), 420  
 write() (pylablib.devices.AWG.specific.Agilent33500  
 method), 414  
 write() (pylablib.devices.AWG.specific.InstekAFG2000  
 method), 432  
 write() (pylablib.devices.AWG.specific.InstekAFG2225  
 method), 426  
 write() (pylablib.devices.AWG.specific.RigolDG1000  
 method), 449  
 write() (pylablib.devices.AWG.specific.RSInstekAFG21000  
 method), 437  
 write() (pylablib.devices.AWG.specific.TektronixAFG1000  
 method), 443  
 write() (pylablib.devices.Cryocon.base.Cryocon1x  
 method), 539  
 write() (pylablib.devices.Cryomagnetics.base.LM500  
 method), 543  
 write() (pylablib.devices.Cryomagnetics.base.LM510  
 method), 547  
 write() (pylablib.devices.Lakeshore.base.Lakeshore218  
 method), 597  
 write() (pylablib.devices.Lakeshore.base.Lakeshore370  
 method), 602  
 write() (pylablib.devices.Sirah.Matisse.SirahMatisse  
 method), 742  
 write() (pylablib.devices.Tektronix.base.DPO2000  
 method), 768  
 write() (pylablib.devices.Tektronix.base.ITektronixScope  
 method), 755  
 write() (pylablib.devices.Tektronix.base.TDS2000  
 method), 761  
 write() (pylablib.devices.Thorlabs.serial.FW method),  
 815  
 write() (pylablib.devices.Thorlabs.serial.FWv1  
 method), 819  
 write() (pylablib.devices.Thorlabs.serial.MDT69xA  
 method), 822  
 write() (pylablib.devices.Thorlabs.serial.ThorlabsSerialInterface  
 method), 812  
 write\_comments() (py-  
 lablib.core.fileio.savefile.CSVTableOutputFileFormat  
 method), 192  
 write\_comments() (py-  
 lablib.core.fileio.savefile.DictionaryOutputFileFormat  
 method), 193  
 write\_comments() (py-  
 lablib.core.fileio.savefile.ITextOutputFileFormat  
 method), 191  
 write\_data() (pylablib.core.fileio.savefile.CSVTableOutputFileFormat  
 method), 192  
 write\_data() (pylablib.core.fileio.savefile.DictionaryOutputFileFormat  
 method), 193  
 write\_data() (pylablib.core.fileio.savefile.IBinaryOutputFileFormat  
 method), 193  
 write\_data() (pylablib.core.fileio.savefile.IOutputFileFormat  
 method), 190  
 write\_data() (pylablib.core.fileio.savefile.ITextOutputFileFormat  
 method), 191  
 write\_data() (pylablib.core.fileio.savefile.TableBinaryOutputFileFormat  
 method), 194

[write\\_file\(\)](#) (pylablib.core.fileio.savefile.CSVTableOutputFileFormat (pylablib.core.fileio.savefile.CSVTableOutputFileFormat attribute), 192  
[write\\_file\(\)](#) (pylablib.core.fileio.savefile.DictionaryOutputFileFormat (pylablib.core.fileio.savefile.DictionaryOutputFileFormat attribute), 193  
[write\\_file\(\)](#) (pylablib.core.fileio.savefile.IBinaryOutputFileFormat (pylablib.core.fileio.savefile.IBinaryOutputFileFormat attribute), 193  
[write\\_file\(\)](#) (pylablib.core.fileio.savefile.IOutputFileFormat (pylablib.core.fileio.savefile.IOutputFileFormat attribute), 190  
[write\\_file\(\)](#) (pylablib.core.fileio.savefile.ITextOutputFileFormat (pylablib.core.fileio.savefile.ITextOutputFileFormat attribute), 191  
[write\\_file\(\)](#) (pylablib.core.fileio.savefile.TableBinaryOutputFileFormat (pylablib.core.fileio.savefile.TableBinaryOutputFileFormat attribute), 194  
[write\\_header\(\)](#) (pylablib.core.utils.general.StreamFileLogger (pylablib.core.utils.general.StreamFileLogger attribute), 380  
[write\\_line\(\)](#) (pylablib.core.fileio.savefile.CSVTableOutputFileFormat (pylablib.core.fileio.savefile.CSVTableOutputFileFormat static method), 192  
[write\\_line\(\)](#) (pylablib.core.fileio.savefile.DictionaryOutputFileFormat (pylablib.core.fileio.savefile.DictionaryOutputFileFormat static method), 193  
[write\\_line\(\)](#) (pylablib.core.fileio.savefile.ITextOutputFileFormat (pylablib.core.fileio.savefile.ITextOutputFileFormat static method), 191  
[write\\_multiple\\_rows\(\)](#) (pylablib.core.fileio.table\_stream.TableStreamFile (pylablib.core.fileio.table\_stream.TableStreamFile attribute), 197  
[write\\_props\(\)](#) (pylablib.core.fileio.savefile.CSVTableOutputFileFormat (pylablib.core.fileio.savefile.CSVTableOutputFileFormat attribute), 192  
[write\\_props\(\)](#) (pylablib.core.fileio.savefile.DictionaryOutputFileFormat (pylablib.core.fileio.savefile.DictionaryOutputFileFormat attribute), 193  
[write\\_props\(\)](#) (pylablib.core.fileio.savefile.ITextOutputFileFormat (pylablib.core.fileio.savefile.ITextOutputFileFormat attribute), 191  
[write\\_row\(\)](#) (pylablib.core.fileio.table\_stream.TableStreamFile (pylablib.core.fileio.table\_stream.TableStreamFile attribute), 197  
[write\\_savetime\(\)](#) (pylablib.core.fileio.savefile.CSVTableOutputFileFormat (pylablib.core.fileio.savefile.CSVTableOutputFileFormat attribute), 192  
[write\\_savetime\(\)](#) (pylablib.core.fileio.savefile.DictionaryOutputFileFormat (pylablib.core.fileio.savefile.DictionaryOutputFileFormat attribute), 193  
[write\\_savetime\(\)](#) (pylablib.core.fileio.savefile.ITextOutputFileFormat (pylablib.core.fileio.savefile.ITextOutputFileFormat attribute), 191  
[write\\_text\\_lines\(\)](#) (pylablib.core.fileio.table\_stream.TableStreamFile (pylablib.core.fileio.table\_stream.TableStreamFile attribute), 197  
[wrng](#) (pylablib.devices.PrincetonInstruments.picam.TROIConstraints (pylablib.devices.PrincetonInstruments.picam.TROIConstraints attribute), 709

## X

[xbins](#) (pylablib.devices.PrincetonInstruments.picam.TROIConstraints (pylablib.devices.PrincetonInstruments.picam.TROIConstraints attribute), 709